



# U.S. Department of Health & Human Services

---

**Food and Drug Administration**

## SAVE REQUEST

**USER:** (kml)  
**FOLDER:** K032115 - 437 pages  
**COMPANY:** POLYGANICS BV (POLYGANICS)  
**PRODUCT:** CUFF, NERVE (JXI)  
**SUMMARY:** Product: NEUROLAC NERVE GUIDE MODELS NG01-15/03, NG01-020/03, NG01 025/03, NG01  
**DATE REQUESTED:** Jun 7, 2016  
**DATE PRINTED:** Jun 7, 2016  
**Note:** Printed



OCT 10 2003

Page 1 of 2

Neurolac® Nerve Guide  
Polyganics BV

Traditional 510(k) Premarket Notification



K032115

510(k)

Summary of Safety and Effectiveness

**Submitter:** Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands  
[www.polyganics.com](http://www.polyganics.com)

**Contact Person:** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Tel : +31 50 588 6588  
Fax : +31 50 588 6599  
Mobile : +31 653 211 303  
E-mail : [hak@polyganics.com](mailto:hak@polyganics.com)

**Date Prepared:** May 20, 2003

**General Provisions:** Trade Name: Neurolac® Nerve guide  
Common Name: Nerve guide  
Classification Name: Nerve Cuff, 21 CFR 882.5275  
Device Classification: Class II

**Predicate Devices:**

- Neurotube™ Neuroregen L.L.C. K983007
- NeuroGen™ Integra Life Sciences Corp. K011168

**Performance Standards** For the Nerve Cuff performance, the FDA, under section 514 of the Food, Drug and Cosmetic Act, has not established standards.

**Indications for Use** The Neurolac® nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.



---

**Device Description**

Neurolac® is designed to be a flexible and transparent resorbable poly(DL-lactide-co-ε-caprolactone) tube to provide a protective environment for peripheral nerve regeneration after injury and to create a conduit to guide axonal growth across a nerve gap.

Neurolac® nerve guides are provided sterile in Tyvek pouch packages and retainer in a variety of sizes.

---

**Performance Data:**

The safety and effectiveness of the Neurolac nerve guides have been demonstrated via data collected from design verification tests and analyses. The design verification testing consisted of the following:

- In vitro suture retention testing
- In vitro degradation testing
- In vivo nerve function recovery

---

**Summary of Substantial Equivalence**

The design, fundamental technology and intended use (safety and efficacy) featured with the Neurolac® Nerve Guide are substantially equivalent to those featured with the competitor devices Neurotube™ (ref. 510(k) 983007; Neuroregen L.L.C.) and the NeuroGen™ Nerve Guide (ref. 510(k) 011168; Integra Life Sciences Corporation).

Biocompatibility, mechanical and physical property testing, in vitro degradation testing, and performance testing in an animal model provide reasonable scientific evidence that Neurolac® nerve guide is substantially equivalent to the predicate devices. Evaluation of the Polyganics Neurolac® Nerve guide based on biocompatibility testing, animal tests, results from literature and the comparison of the Neurolac® nerve guide with its predicate devices, shows that the Neurolac® nerve guide is safe for implantation.



OCT 1 0 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Jan Bart Hak, Ph.D.  
Manager, Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands

Re: K032115  
Trade/Device Name: Neurolac® Nerve Guide  
Regulation Number: 21 CFR 882.5275  
Regulation Name: Nerve cuff  
Regulatory Class: II  
Product Code: JXI  
Dated: July 3, 2003  
Received: July 17, 2003

Dear Dr. Hak:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 - Jan Bart Hak, Ph.D.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

*Miriam C. Provost*  
for

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and Radiological Health

Enclosure





K032115/A1

Polyganics BV  
L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands

Telephone (31) 50 588 65 88  
Telefax (+31) 50 588 65 99  
E-mail [mail@polyganics.com](mailto:mail@polyganics.com)  
Internet [www.polyganics.com](http://www.polyganics.com)

Mrs. M. Shulman  
Food and Drug Administration  
Center for Devices and Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850  
USA

Your ref. : 510k# K032115  
Our ref. : Neurolac Nerve guide  
Subject : Correspondence address

2003 JUL -4 P 2:30

Groningen, July 30, 2003

Dear Mrs. Shulman,


In relation to the premarket notification of our product Neurolac nerve guide with the 510(k) number K032115, we have noted that you use our US agent's address for correspondence.

I would greatly appreciate if it would be possible to send any future correspondence directly to me so that we can address any issue timely and efficiently. Our address is (as indicated on the submission cover sheet):

Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands

Tel: +31 50 588 6588  
Fax: +31 50 588 6599  
E-mail: [hak@polyganics.com](mailto:hak@polyganics.com)

Sincerely yours,

  
Jan Bart Hak  
Manager Clinical and Regulatory Affairs

All our orders and deliveries are subject to our General Terms and Conditions, registered at the Chamber of Commerce in Groningen under number 02067151. On request a copy will be sent free of charge.

Bank account  
ABN AMRO Groningen 51.35.29.195  
Commercial Reg. No. 02067151  
VAT No. NL80.85.30.525.B01

SKS-8 17

Memorandum

Date: 10/22/03

From: DMC (HFZ-401)

Subject: Premarket Notification Number(s): K032115/A2

To: Division Director: NE / DGRND

The attached information has been received by the 510(k) DMC on the above referenced 510(k) submission(s). Since a final decision has been rendered, this record is officially closed.

Please review the attached document and return it to the DMC, with one of the statements checked below.

Information does not change the status of the 510(k); no other action required by the DMC; please add to image file. (Prepare K-25) THIS DOES NOT APPLY TO TRANSFER OF OWNERSHIP. PLEASE BRING ANY TRANSFER OF OWNERSHIP TO POS.

Additional information requires a new 510(k); however, the information submitted is incomplete; (Notify company to submit a new 510(k); [Prepare the K30 Letter on the LAN]

No response necessary (e.g., hard copy of fax for the truthful and accuracy statement, 510(k) statement, change of address, phone number, or fax number).

**CLIA CATEGORIZATION refers to laboratory test system devices reviewed by the Division of Clinical Laboratory Devices (HFZ-440)**

Information requires a CLIA CATEGORIZATION; the complexity may remain the same as the original 510(k) or may change as a result of the additional information (Prepare a CAT letter)

Additional information requires a CLIA CATEGORIZATION; however, the information submitted is incomplete; (call or fax firm)

No response necessary

This information should be returned to the DMC within 10 working days from the date of this memorandum.

Reviewed by: Daino B. Barbera

Date: 10/23/03

DA/B

Draft #2 : 9/8/99  
Draft #3: 1/3/00  
Draft #4: 3/7/03

OCT 23

nmj



Polyganics BV  
L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands

Telephone (31) 50 588 65 88  
Telefax (+31) 50 588 65 99  
E-mail mail@polyganics.com  
Internet www.polyganics.com

K032115/A<sup>2</sup>

David B. Berkowitz, Ph.D., V.M.D.  
Center for Devices and Radiological Health  
Food and Drug Administration  
9200 Corporate Blvd. HFZ-410  
Rockville, MD 20850  
USA

Your ref. : K032115  
Our ref. : 03-061  
Subject : Nerve Cuff, Neurolac


Groningen, Oct 8, 2003

Dear Mr. Berkowitz,

Please find enclosed the hardcopies of the files, which I've sent to you by electronic mail on Oct 8, 2003. This package includes the "510k file" and the instructions for use.

Of the 510k file, page 6, 8, 15, 17, 18 and 32 (this page is part of the 510k summary) are adjusted according to your recommendation.

Yours sincerely,

  
Jan-Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs

FDI/AMRO  
2003 OCT 22 P 3:18

SK-33

All our orders and deliveries are subject to our General Terms and Conditions, registered at the Chamber of Commerce in Groningen under number 02067151. On request a copy will be sent free of charge.

Bank account  
ABN AMRO Groningen 51.35.29.95  
Commercial Reg. No. 02067151  
VAT No. NL80.85.30.525.B01



**510(k)  
PREMARKET NOTIFICATION  
TABLE OF CONTENTS**

<u>Section</u>	<u>Description</u>	<u>Page</u>
FDA ADMINISTRATIVE FORMS .....		3
	Submission Cover Sheet .....	4
	Indications for Use Form .....	8
APPLICANT STATEMENTS .....		9
	Truthful And Accurate Statement .....	10
	Substantial Equivalence Terminology Statement .....	11
Section 1: GENERAL INFORMATION .....		12
	Applicant .....	12
	Contact Person .....	12
	Device Name .....	12
	Device Classification .....	12
	Manufacturing Facility .....	12
	Sterilization Facility (EtO) .....	13
	Performance Standards / Special Controls .....	13
	Purpose of Premarket Notification .....	13
	Predicate Devices .....	13
	Software Validation and Certification .....	13
	Kit Certification .....	13
	Class III Certification .....	14
	510(k) Summary .....	14
Section 2: DEVICE DESCRIPTION .....		15
	Drawings .....	15
	Intended Use .....	15
	Device Description .....	15
Section 3: PROPOSED LABELING .....		17
	Subject Device Labeling .....	17
	Intended Use of the Subject Device .....	17
	Promotional Materials .....	17
Section 4: COMPARATIVE INFORMATION .....		18
	Background .....	18
	Intended Use .....	18
	Subject Device .....	18
	Predicate Device .....	18
	Device Characteristics .....	18
Section 5: PERFORMANCE VERIFICATION .....		20
	Performance Verification .....	20
	Suture retention testing .....	21
	In vitro degradation testing .....	22
	Experimental details .....	22
	Nerve function recovery: sciatic nerve model .....	24
	Neurolac nerve guide versus autologous nerve graft .....	25



---

Functional nerve recovery after bridging a 15 mm nerve gap with a biodegradable nerve guide .....	26
Section 6: BIOCOMPATIBILITY .....	28
Background .....	28
Testing .....	28
Section 7: STERILIZATION AND PYROGENICITY INFORMATION .....	29
Introduction .....	29
Sterilization Method .....	29
Validation Method .....	29
Sterility Assurance Level .....	29
Pyrogen Test Method .....	29
Section 8: PACKAGING AND SHELF-LIFE INFORMATION .....	30
Introduction .....	30
Packaging Description .....	30
Product Shelf – Life .....	30
APPENDIX A: 510(K) SUMMARY OF SAFETY AND EFFECTIVNESS .....	31
APPENDIX B: DEVICE DRAWING .....	34
APPENDIX C: DEVICE LABELING .....	35
APPENDIX D: PREDICATE DEVICE NEUROTUBE™ .....	39
APPENDIX E: PREDICATE DEVICE NEUROGEN™ .....	40
APPENDIX F: PERFORMANCE TEST REPORTS .....	41
APPENDIX G: BIOCOMPATIBILITY DATA .....	47



---

## FDA ADMINISTRATIVE FORMS

Premarket Submission Cover Sheet

Indications for Use Form



<b>Submission Cover Sheet</b>				
Date of Submission: <b>July 3, 2003</b>			FDA Document Number:	
<b>Section A Type of Submission</b>				
<b>PMA</b> <input type="radio"/> Original Submission <input type="radio"/> Modules Submission <input type="radio"/> Amendment <input type="radio"/> Report <input type="radio"/> Report Amendment	<b>PMA Supplement</b> <input type="radio"/> Regular <input type="radio"/> Special <input type="radio"/> Panel Track <input type="radio"/> 30-day Supplement <input type="radio"/> 30-day Notice <input type="radio"/> 135-day Supplement <input type="radio"/> Real-time Review <input type="radio"/> Amendment to PMA Supplement	<b>PDP</b> <input type="radio"/> Presubmission summary <input type="radio"/> Original PDP <input type="radio"/> Notice of intent to start clinical trials <input type="radio"/> Intention to submit Notice of Completion <input type="radio"/> Notice of Completion <input type="radio"/> Amendment to PDP <input type="radio"/> Report	<b>510(k)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> <u>Traditional</u> <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated <input type="checkbox"/> Additional Information <input type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated	<b>Meeting</b> <input type="checkbox"/> Pre-IDE meeting <input type="checkbox"/> Pre-PMA meeting <input type="checkbox"/> Pre-PDP meeting <input type="checkbox"/> 180-day meeting <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="radio"/> Original Submission <input type="radio"/> Amendment <input type="radio"/> Supplement	<b>Humanitarian Device Exemption</b> <input type="radio"/> Original submission <input type="radio"/> Amendment <input type="radio"/> Supplement <input type="radio"/> Report	<b>Class II Exemption</b> <input type="radio"/> Original submission <input type="radio"/> Additional information	<b>Evaluation of Automatic Class III Designation</b> <input type="radio"/> Original submission <input type="radio"/> Additional information	<b>Other Submission</b> Describe submission:
<b>Section B Applicant or Sponsor</b>				
Company / Institution name: <b>Polyganics BV</b>		Establishment registration number: <b>Not registered.</b>		
Division name (if applicable): <b>Not Applicable</b>		Phone number (include area code): <b>+31 50 588 6588</b>		
Street address: <b>L.J. Zielstraweg 1</b>		FAX number (include area code): <b>+31 50 588 6599</b>		
City: <b>Groningen</b>	State/Province: <b>Groningen</b>	Country: <b>The Netherlands</b>	Zip/Postal Code: <b>9713-GX</b>	
Contact name: <b>Jan-Bart Hak</b>				
Contact title: <b>Manager Clinical and Regulatory Affairs</b>		Contact e-mail address: <b>hak@polyganics.com</b>		
<b>Section C Submission Correspondent (if different from above)</b>				
Company / Institution name: not applicable		Establishment registration number: not applicable		
Division name (if applicable): not applicable		Phone number (include area code): not applicable		
Street address: not applicable		FAX number (include area code): not applicable		
City: not applicable	State/Province: not applicable	Country: not applicable	Zip/Postal Code: not applicable	
Contact name: not applicable				
Contact title: not applicable		Contact e-mail address: not applicable		




Section D1 Reason for Submission – PMA, PDP, or HDE		
<ul style="list-style-type: none"> <li><input type="checkbox"/> New device</li> <li><input type="checkbox"/> Withdrawal</li> <li><input type="checkbox"/> Additional or expanded indications</li> <li><input type="checkbox"/> Licensing agreement</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in design, component, or specifications:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Software</li> <li><input type="checkbox"/> Color Additive</li> <li><input type="checkbox"/> Material</li> <li><input type="checkbox"/> Specifications</li> <li><input type="checkbox"/> Other (specify below):</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Location Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Manufacturer</li> <li><input type="checkbox"/> Sterilizer</li> <li><input type="checkbox"/> Packager</li> <li><input type="checkbox"/> Distributor</li> </ul> </li> </ul>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Process Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Manufacturing</li> <li><input type="checkbox"/> Sterilization</li> <li><input type="checkbox"/> Packaging</li> </ul> </li> <li><input type="checkbox"/> Other (specify below):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Labeling Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Indications</li> <li><input type="checkbox"/> Instructions</li> <li><input type="checkbox"/> Performance characteristics</li> <li><input type="checkbox"/> Shelf Life</li> <li><input type="checkbox"/> Trade Name</li> </ul> </li> <li><input type="checkbox"/> Other (specify below):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Report submissions:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Annual or periodic</li> <li><input type="checkbox"/> Post-approval study</li> <li><input type="checkbox"/> Adverse reaction</li> <li><input type="checkbox"/> Device defect</li> <li><input type="checkbox"/> Amendment</li> </ul> </li> </ul>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Response to FDA correspondence:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Request for applicant hold</li> <li><input type="checkbox"/> Request for removal of applicant hold</li> <li><input type="checkbox"/> Request for extension</li> <li><input type="checkbox"/> Request to remove or add manufacturing site</li> </ul> </li> <li><input type="checkbox"/> Other reason (specify):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in ownership</li> <li><input type="checkbox"/> Change in correspondent</li> </ul>	
Section D2 Reason for Submission – IDE		
<ul style="list-style-type: none"> <li><input type="checkbox"/> New device</li> <li><input type="checkbox"/> Addition of institution</li> <li><input type="checkbox"/> Expansion/extension of study</li> <li><input type="checkbox"/> IRB certification</li> <li><input type="checkbox"/> Request hearing</li> <li><input type="checkbox"/> Request waiver</li> <li><input type="checkbox"/> Termination of Study</li> <li><input type="checkbox"/> Withdrawal of application</li> <li><input type="checkbox"/> Unanticipated adverse effect</li> <li><input type="checkbox"/> Notification of emergency use</li> <li><input type="checkbox"/> Compassionate use request</li> <li><input type="checkbox"/> Treatment IDE</li> <li><input type="checkbox"/> Continuing availability request</li> <li><input type="checkbox"/> Other reason (specify):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Correspondent</li> <li><input type="checkbox"/> Design</li> <li><input type="checkbox"/> Informed Consent</li> <li><input type="checkbox"/> Manufacturer</li> <li><input type="checkbox"/> Manufacturing process</li> <li><input type="checkbox"/> Protocol – feasibility</li> <li><input type="checkbox"/> Protocol – other</li> <li><input type="checkbox"/> Sponsor</li> </ul> </li> <li><input type="checkbox"/> Report Submission:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Current investigator</li> <li><input type="checkbox"/> Annual progress</li> <li><input type="checkbox"/> Site waiver limit reached</li> <li><input type="checkbox"/> Final</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Response to FDA letter concerning:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Conditional approval</li> <li><input type="checkbox"/> Deemed approved</li> <li><input type="checkbox"/> Deficient final report</li> <li><input type="checkbox"/> Deficient progress report</li> <li><input type="checkbox"/> Deficient investigator report</li> <li><input type="checkbox"/> Disapproval</li> <li><input type="checkbox"/> Request extension of time to respond to FDA</li> <li><input type="checkbox"/> Request meeting</li> </ul> </li> </ul>
Section D3 Reason for Submission – 510(k)		
<ul style="list-style-type: none"> <li><input checked="" type="checkbox"/> New device</li> <li><input type="checkbox"/> Additional or expanded indications</li> <li><input type="checkbox"/> Other reason (specify)</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in technology</li> <li><input type="checkbox"/> Change in design</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in materials</li> <li><input type="checkbox"/> Change in manufacturing process</li> </ul>



<b>Section E Additional Information on 510(k) Submissions</b>					
Product codes of devices to which substantial equivalence is claimed:				Summary of, or statement concerning, safety and effectiveness data:	
1 <b>JXI</b>	2	3	4	<input type="checkbox"/> 510(k) summary attached <input type="checkbox"/> 510(k) statement	
Information on devices to which substantial equivalence is claimed:					
510(k) Number		Trade or proprietary or model name		Manufacturer	
1. #K983007		Neurotube™		Neuroregen L.L.C.	
2. #K011168		NeuroGen™ Nerve Guide		Integra Life Sciences Corp.	
<b>Section F Product Information – Applicable to All Applications</b>					
Common or usual or classification name: Nerve Cuff					
Trade or proprietary or model name				Model number	
Neurolac® Nerve Guide				NG01-015/03	
Neurolac® Nerve Guide				NG01-020/03	
Neurolac® Nerve Guide				NG01-025/03	
Neurolac® Nerve Guide				NG01-030/03	
FDA document numbers of all prior related submissions (regardless of outcome):					
1 N.A.	2	3	4	5	6
Data included in submission: <input type="checkbox"/> Laboratory testing <input type="checkbox"/> Animal trials <input type="checkbox"/> Human trials					
<b>Section G Product Classification – Applicable to All Applications</b>					
Product code: <b>JXI</b>		C.F.R. section: <b>21 CFR 882.5275</b>		Device class: <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
				Device class: <input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
				Device class: <input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
Classification panel: <b>Neurology (Neurological Therapeutic Devices)</b>					
<b>Indications (from labeling):</b> The Neurolac nerve guide is indicated for the reconstruction of peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.					



		FDA Document Number:	
<b>Section H Manufacturing / Packaging / Sterilization Sites</b>			
<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer	
		<input type="checkbox"/> Contract sterilizer <input type="checkbox"/> Repackager / relabeler	
Company/institution name: <b>Polyganics BV</b>		Establishment registration number: <b>Not registered</b>	
Division name (if applicable):		Phone number (include area code): <b>+31 50 588 6588</b>	
Street address: <b>L.J. Zielstraweg 1</b>		Fax number (include area code): <b>+31 50 588 6599</b>	
City: <b>Groningen</b>	State/Province: <b>Groningen</b>	Country: <b>The Netherlands</b>	Zip/Postal Code: <b>9713-GX</b>
Contact Name: <b>(see Section B)</b>			
Contact Title: <b>(see Section B)</b>		Contact e-mail address: <b>(see Section B)</b>	
<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer	
		<input checked="" type="checkbox"/> Contract sterilizer (EtO) <input type="checkbox"/> Repackager / relabeler	



Company/institution name:		Establishment registration number:	
Division name (if applicable):		Phone number (include area code): ( )	
Street address:		Fax number (include area code): ( )	
City:	State/Province:	Country:	Zip/Postal Code:
Contact Name:			
Contact Title:		Contact e-mail address:	



---

**Indications for Use Form**

510(k) Number: \_\_\_\_\_

Device Name: **Neurolac® Nerve Guide**

**Indications for Use:**

The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use \_\_\_\_\_ OR Over-The-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109) (Optional Format 1-2-96)

\_\_\_\_\_  
(Division Sign-Off)

510(k) Number \_\_\_\_\_



## APPLICANT STATEMENTS

Truth and Accuracy Certification

Substantial Equivalence Terminology Statement



---

## PREMARKET NOTIFICATION

### Truthful And Accurate Statement

[As Required by 21 CFR 807.87(k)]

I certify that, in my capacity as a Regulatory Affairs representative of Polyganics BV, I believe to the best of my knowledge, that all data and information submitted in the pre-market notification are truthful and accurate and that no material fact has been omitted.

---

J.B. Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

---

Date

---

Premarket Notification [510(k)] Number

## **Substantial Equivalence Terminology Statement**

### **USE OF THE TERM "SUBSTANTIALLY EQUIVALENT"**

The use of the term "substantially equivalent" as used herein is intended to be a determination of substantial equivalence under the Federal Food, Drug, and Cosmetic Act, as Amended, and relates to the fact that the product can be marketed without premarket approval or reclassification. Such a determination is not intended to have any bearing whatsoever on the resolution of patent infringement suits or other patent matters.



---

## Section 1: GENERAL INFORMATION

**Applicant** Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

FDA Establishment Registration Number: Not Registered yet

---

**Contact Person** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

Tel: +31 50 588 6588  
Fax: +31 50 588 6599  
E-mail: [hak@polyganics.com](mailto:hak@polyganics.com)

---

Device Name	Trade Name	Common Name	Classification Name
	Neurolac Nerve Guide	Nerve Guide	Nerve Cuff 21 CFR 882.5275

---

**Device Classification**

**Classification Name:** "Nerve Cuff" (Ref. Codes of Federal regulations, title 21 – Food and Drugs, Part 882 – Neurological devices, subpart F – Neurological Therapeutic devices, Sec 882.5275)

**Class:** Class II

**Description:** (a) Identification. A nerve cuff is a tubular silicone rubber sheath used to encase a nerve for aid in repairing the nerve (e.g., to prevent ingrowth of scar tissue) and for capping the end of the nerve to prevent the formation of neuroma (tumors).

---

**Manufacturing Facility** Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

FDA Establishment Registration Number: not registered

---



**Sterilization Facility (EtO)**

(b) (4)

**Performance Standards / Special Controls**

No performance standards are indicated for this product.

**Purpose of Premarket Notification**

The reason for this premarket notification is to inform the Food and Drug Administration (FDA) of our intent to market the Polyganics Neurolac Nerve Guide.

**Predicate Devices**

The following table provides information on the predicate devices.

Trade Name	Manufacturer	510k #	Concurrence Date	Substantial Equivalence
Neurotube™	Neuroregen L.L.C.	K983007 Traditional APPENDIX D	03/22/1999	The tube provides an optimal environment for longitudinal nerve axon growth of the peripheral nerve. For single use only in patients with a peripheral nerve injury where the nerve gap is more than or equal to 8 mm, but less than or equal to 3 cm.
NeuroGen™ Nerve Guide	Integra Life Sciences Corp.	K011168 Traditional APPENDIX E	06/22/2001	NeuroGen™ Nerve guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**Software Validation and Certification**

This section is not applicable, since the Neurolac Nerve Guide does not utilize software in the performance of its intended use.

**Kit Certification**

This section is not applicable.



---

**Class III Certification**

This section is not applicable. The Neurolac Nerve Guide is a Class II device

---

**510(k) Summary**

A summary of safety and effectiveness for the Neurolac Nerve Guide is provided in APPENDIX A.

---

## Section 2: DEVICE DESCRIPTION

**Drawings** Please refer to APPENDIX B of this 510(k) premarket notification for an engineering drawing of the subject device.


---

**Intended Use** The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.  
There are no known contraindications.

For FDA administrative purposes, the intended use of the Polyganics Neurolac Nerve Guide is also documented in a separate form that can be found at the beginning of the 510(k) notification in the section entitled "FDA Administrative Forms"

---

**Device Description** The Neurolac nerve guide is composed of the bioresorbable copolyester poly(DL-lactide- $\epsilon$ -caprolactone). The Neurolac nerve guide provides guidance and protection to regenerating axons.



**Figure 1.** Example of the Neurolac nerve guide

The Neurolac nerve guide elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the tube by fibrous tissue. Degradation of the Neurolac nerve guide occurs through hydrolysis leading to gradual reduction of molecular weight. The Neurolac nerve guide retains its initial mechanical properties up to 8 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and  $\omega$ -hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. Animal studies demonstrated that a Neurolac nerve guide is resorbed within 16 months.

The Neurolac nerve guide inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The Neurolac nerve guide is indicated for single-use.

The Polyganics Nerve Guide further consists of:



- 
- Packaging
  - Labeling
  - IFU

---

### Section 3: PROPOSED LABELING

**Subject Device Labeling**      The following proposed labeling for the subject device Neurolac Nerve Guide is provided in APPENDIX C:

- Outer label (Carton)
- Inner label (Pouch)
- Pre-printed carton text and graphics
- Instructions for Use

---

**Intended Use of the Subject Device**      The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

---

**Promotional Materials**      At present, no promotional materials are available for the subject device.



## Section 4: COMPARATIVE INFORMATION

**Background** The Polyganics Neurolac nerve guide is a biodegradable tube for the repair of transected peripheral nerves.

**Intended Use** The Indications for Use for the subject and predicate devices are described in the table below.

Subject Device	Indication for Use
Neurolac	The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

Predicate Device	Indication for Use
Neurotube™ APPENDIX D	The Neurotube is intended for single use in patients with an injury to a peripheral nerve, in which the nerve gap is more than or equal to 8 mm but less than or equal to 3 cm. The nerve gap may be created primarily at the time of injury or created secondarily at the time of exploration of failed primary repair
NeuroGen™ APPENDIX E	NeuroGen Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**Device Characteristics** The technological characteristics (i.e., design, dimensions material etc.) of the subject Neurolac Nerve Guide and the predicate devices are presented in the table below. The table provides a comparison, demonstrating that the Neurolac is substantially equivalent to the currently marketed predicate device.



**Device Characteristics of the Subject and Predicate devices**

Characteristics General	Subject Device	Predicate Devices	
	Neurolac®	Neurotube™	NeuroGen™
<b>510(k) Reference</b>	<b>This 510(k)</b>	<b>K983007</b>	<b>K011168</b>
<b>Sterile</b>	Sterile device	Sterile device	Sterile device
<b>Single Use</b>	Single-use	Single-use	Single-use
<b>Contents packaging</b>	Nerve guide and instructions for use	Nerve guide and instructions for use	Nerve guide and instructions for use
<b>Length</b>	(b)(4) Product specs		
<b>Inner Diameter (In.)</b>			
<b>Material</b>			
<b>Biodegradable</b>	Yes	Yes	Yes
<b>Animal derived</b>	No	Yes	No
<b>Indication</b>	Peripheral Nerve discontinuity	Peripheral Nerve discontinuity	Peripheral Nerve discontinuity
<b>Transparent</b>	Yes	No	No



## Section 5: PERFORMANCE VERIFICATION

### Performance Verification


The performance testing on the Neurolac Nerve Guide was conducted to verify that the meets performance characteristics. The following performance tests were conducted:

	Subject/Test	Report
1	(b) (4)	
2		
3		
4		
5		

Testing was performed on finished sterile products. The test results were all favorable for the Neurolac Nerve guide. The complete performance test reports have been attached as APPENDIX F.




(b)(4) Testing






(b)(4) Testing






(b)(4) Testing






(b)(4) Testing






(b)(4) Testing





(b)(4) Testing





(b)(4) Testing






---

## Section 6: BIOCOMPATIBILITY

(b)(4) Testing

**Background**

**Testing**



---

## Section 7: STERILIZATION AND PYROGENICITY INFORMATION

Introduction (b) (4)

**Sterilization  
Method**

**Validation  
Method**

**Sterility As-  
surance Level**

**Pyrogen Test  
Method**

---

## Section 8: PACKAGING AND SHELF-LIFE INFORMATION

**Introduction** The Neurolac nerve guides are packaged in a tray (mechanical protection) and subsequently sealed in a Tyvek/PET-PE pouch. The packaging does comply with standards ISO11607.

---

**Packaging Description** The Neurolac nerve guides are packaged in a tray serving as a mechanical protection.

The tray together with a Neurolac nerve guide are packaged in a Tyvek/PET-PE pouch. The pouch carries a pouch label on the PET-PE layer of the pouch.

The pouch with product in the tray is packaged together with the instructions for use in a carton box. The carton box carries the carton box label.

---

**Product Shelf – Life** The Neurolac nerve guide has a shelf-life (Use By date) of 12 months.

## **APPENDIX A: 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS**

---

**510(k)**  
**Summary of Safety and Effectiveness**

---

**Submitter:** Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands  
[www.polyganics.com](http://www.polyganics.com)

**Contact Person:** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Tel : +31 50 588 6588  
Fax : +31 50 588 6599  
Mobile : +31 653 211 303  
E-mail : [hak@polyganics.com](mailto:hak@polyganics.com)

**Date Prepared:** May 20, 2003

---

**General Provisions:** Trade Name: Neurolac® Nerve guide  
Common Name: Nerve guide  
Classification Name: Nerve Cuff, 21 CFR 882.5275  
Device Classification: Class II

---

**Predicate Devices:**

• Neurotube™	Neuroregen L.L.C.	K983007
• NeuroGen™	Integra Life Sciences Corp.	K011168

---

**Performance Standards** For the Nerve Cuff performance, the FDA, under section 514 of the Food, Drug and Cosmetic Act, has not established standards.

---

**Indications for Use** The Neurolac® nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

---



---

**Device Description**

Neurolac® is designed to be a flexible and transparent resorbable poly(DL-lactide-co-ε-caprolactone) tube to provide a protective environment for peripheral nerve regeneration after injury and to create a conduit to guide axonal growth across a nerve gap.

Neurolac® nerve guides are provided sterile in Tyvek pouch packages and retainer in a variety of sizes.

---

**Performance Data:**

The safety and effectiveness of the Neurolac nerve guides have been demonstrated via data collected from design verification tests and analyses. The design verification testing consisted of the following:

- In vitro suture retention testing
- In vitro degradation testing
- In vivo nerve function recovery

---

**Summary of Substantial Equivalence**


The design, fundamental technology and intended use (safety and efficacy) featured with the Neurolac® Nerve Guide are substantially equivalent to those featured with the competitor devices Neurotube™ (ref. 510(k) 983007; Neuroregen L.L.C.) and the NeuroGen™ Nerve Guide (ref. 510(k) 011168; Integra Life Sciences Corporation).

Biocompatibility, mechanical and physical property testing, in vitro degradation testing, and performance testing in an animal model provide reasonable scientific evidence that Neurolac® nerve guide is substantially equivalent to the predicate devices. Evaluation of the Polyganics Neurolac® Nerve guide based on biocompatibility testing, animal tests, results from literature and the comparison of the Neurolac® nerve guide with its predicate devices, shows that the Neurolac® nerve guide is safe for implantation.

---

## APPENDIX B: DEVICE DRAWING

(b) (4)



## **APPENDIX C: DEVICE LABELING**

**Subject device: Neurolac Nerve Guide**

**DRAFT**

## DEVICE LABELING


### Carton, DRAFT

<b>REF</b> NG01-02B03 <b>LOT</b> NSA2003 04 0701 2004-04	<b>NEUROLAC®</b>	<b>STERILE</b> Sterile	 <b>POLYGANICS</b> Polymeric innovations in tissue recovery
		<b>STERILE EO</b> Sterilized with ethylen oxide gas	
	3cm	<b>LOT</b> Lot No.	
	I.D. 2.0 mm	 Use By	
	<b>Bioresorbable peripheral nerve guide</b>	 For one use only	
<b>Bioresorbierbaren peripheren nervenleitschiene</b>	 Attention, see instructions for Use		
<b>Bioresorbierbare perifere zenuwgeleider</b>	<b>REF</b> Catalog No.		
<b>Guide de nerf périphérique bioresorbable</b>	<b>Nonpyrogenic.</b> Do not use open or damaged packages. Store in a cool, dry place at or below 4 °C. Caution: Federal (USA) law restricts this device to sale by or on the order of a physician		
<b>Guía del nervio periférico bioabsorbible</b>	 <b>POLYGANICS</b> Polymeric innovations in tissue recovery		
	L.J. Zielstraweg 1 9713 GX Groningen The Netherlands		

## DEVICE LABELING

### Inner label pouch, DRAFT

**NEUROLAC®**



3cm I.D. 2.0 mm

Bioresorbable peripheral nerve guide  
Bioresorbierbaren peripheren nervenleitschiene  
Bioresorbeerbare perifere zenuwgeleider  
Guide de nerf périphérique biorésorbable  
Guía del nervio periférico bioabsorbible

**STERILE**  
**STERILE EO**

REF **NG01-020/03**

**LOT** **NGA2003 04 2501**

 **2004-04**

Nonpyrogenic. Do not use open or damaged packages. Store in a cool, dry place at or below 4 °C. **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician



L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands

NGA L101-001/03

## **DEVICE LABELING**

**Instructions for Use, DRAFT**

**PG002a 2003-01-09 MM024 IFU EN version 19**



## **APPENDIX D: PREDICATE DEVICE NEUROTUBE™**

**Neuroregen L.L.C.**

**510(k): K983007**



---

**APPENDIX E: PREDICATE DEVICE NEUROGEN™**

**Integra Life Sciences Corporation**

**510(k): K011168**

## **APPENDIX F: PERFORMANCE TEST REPORTS**

**Suture retention testing**

**In vitro degradation testing**

**Nerve function recovery: sciatic nerve model**

**Neurolac nerve guide versus autologous nerve graft**

**Functional nerve recovery after bridging a 15 mm nerve gap with a Neurolac nerve guide**



---

## Suture retention testing



---

## In vitro degradation testing



---

## Nerve function recovery: sciatic nerve model



---

## Neurolac nerve guide versus autologous nerve graft



---

**Functional nerve recovery after bridging a 15 mm nerve gap with a Neurolac nerve guide**

## **APPENDIX G: BIOCOMPATIBILITY DATA**

**Cytotoxicity**  
**Irritation**  
**Sensitization**  
**Hemocompatibility**  
**Acute systemic toxicity**  
**Pyrogenicity**  
**Mutagenicity /genotoxicity**  
**Sub chronic toxicity**  
**Carcinogenicity**  
**Chronic toxicity**  
**Reproductive toxicity**  
**Implantation**



---

## Cytotoxicity



---

**Irritation**



**Sensitization**



---

## Hemocompatibility



---

---

## Acute systemic toxicity



---

## Pyrogenicity



---

---

## Mutagenicity /genotoxicity



---

## Sub chronic toxicity



---

**Carcinogenicity**

**Chronic toxicity**

**Reproductive toxicity**



---

## Implantation



Instructions for Use, English

**STERILE.** Sterilized with ethylene oxide gas. For single use only. Do not autoclave.  
**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

**Neuroloc® peripheral nerve guide**

**Description**

The Neuroloc nerve guide is composed of the bioresorbable copolyester poly(DL-lactide-ε-caprolactone). The Neuroloc nerve guide provides guidance and protection to regenerating axons. The Neuroloc nerve guide elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the tube by fibrous tissue. Degradation of the Neuroloc nerve guide occurs through hydrolysis leading to gradual reduction of molecular weight. The Neuroloc nerve guide retains its initial mechanical properties up to 8 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and ε-hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. Animal studies demonstrated that a Neuroloc nerve guide is resorbed within 16 months.

The Neuroloc nerve guide inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The Neuroloc nerve guide is indicated for single-use.

**Indications**

The Neuroloc nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

**Contraindications**

There are no known contraindications.

**Warnings**

- The Neuroloc nerve guide is for single use only. Do not resterilize or re-use. Structural integrity and/or function may be impaired through cleaning, resterilization, or re-use and may cause adverse patient reactions. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from re-use of (or any part of) the Neuroloc nerve guide.
- Sterile unless package has been opened or damaged.
- Discard open unused nerve guides;
- The Neuroloc nerve guide should only be used by physicians who are trained in nerve defect repair techniques. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel. The physician should consult recent literature on current medical practice on peripheral nerve repair.
- Nerve regeneration may be suboptimal in elderly, malnourished or debilitated patients or in patients suffering from cancer, anaemia, obesity, diabetes, infection or other conditions which may delay wound healing, infected wounds, or moderate tissue inflammatory response characteristic of foreign body response.
- Precautions**
- Use prior to "Use by date";
- Store in dark, dry place at or below 4°C (39°F);
- Do not expose the nerve guide to organic solvents (e.g. chloroform, acetone);
- Do not use absorbable sutures for fixation of the nerve stumps into the nerve guide;



Refer to accompanying instructions for use.

Use by

Catalog number

For single use only

Lot number

Sterile product

Sterilized by ethylene oxide

Polyganics BV Telephone +31 50 588 6588  
 L.J. Zielsraweg 1 Fax +31 50 588 6599  
 9713 GX Groningen  
 The Netherlands

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**  
 There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Polyganics product(s) described in this publication. Under no circumstances shall Polyganics BV be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Polyganics to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Polyganics BV printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Polyganics BV will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Neuroloc® is a registered trademark of Polyganics



**NOTE:** It is recommended that the nerve ends are pulled into the tube for at least 3 mm for optimal nerve regeneration.

- Fill the tube with heparinized saline, using a solution containing 1000 units of heparin per 100 ml of normal saline (Fig 1.5).
- Subsequently, use the same procedure, to pull the distal nerve stump into the nerve guide.
- A minimum space of 5 mm should be left between the nerve ends in the nerve guide.
- Fill any remaining space with heparinized saline (Fig 1.6) by injecting along the nerve into the lumen of the tube or by penetrating the tube (not the nerve).




Figure 1. Schematic representation of suture technique for suturing the nerve ends into the nerve guide.

**CAUTION:** Ensure that no blood enters the nerve guide lumen since this may hinder nerve recovery.

**CAUTION:** The nerve guide should be implanted and sutured with all joints in an extended position as to assure that no tension occurs on the proximal or distal nerve end when joints are being mobilized.

Close the wound and splint to prevent kinking for the first three postoperative weeks. Long-term compression of the nerve guide should be avoided. Patients may be administered oral antibiotics for the first post-operative week.

Dispose contaminated implantation and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.

60



OCT 10 2003

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Jan Bart Hak, Ph.D.  
Manager, Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands

Re: K032115  
Trade/Device Name: Neurolac® Nerve Guide  
Regulation Number: 21 CFR 882.5275  
Regulation Name: Nerve cuff  
Regulatory Class: II  
Product Code: JXI  
Dated: July 3, 2003  
Received: July 17, 2003

Dear Dr. Hak:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4659. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/edrh/dsma/dsmamain.html>

Sincerely yours.

*Miriam C. Provost*  
For Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



**Indications for Use Form**

510(k) Number:     K032115    

Device Name:     Neurolac® Nerve Guide    

Indications for Use:

The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use  OR Over-The-Counter Use   
(Per 21 CFR 801.109) (Optional Format 1-2-96)

\_\_\_\_\_  
(Division Sign-Off)

510(k) Number \_\_\_\_\_

    Miriam C. Provost      
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

July 21, 2003

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

POLYGANICS BV  
C/O VIKTOR J. NICKOLSON  
33 RIVER ST. UNIT 813  
HOBOKEN, NJ 07030  
ATTN: VIKTOR J. NICKOLSON

510(k) Number: K032115  
Received: 17-JUL-2003  
Product: NEUROLAC NERVE GUIDE  
MODELS NG01-15/03,  
NG01-020/03, NG01  
025/03, NG01-030/03

The Food and Drug Administration (FDA), Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. **YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.**

The Act, as amended by the Medical Device User Fee and Modernization Act of 2002 (MDUFMA)(Public Law 107-250), authorizes FDA to collect user fees for premarket notification submissions. (For more information on MDUFMA, you may refer to our website at <http://www.fda.gov/oc/mdufma>).

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (DMC)(HFZ-401) at the above letterhead address. Correspondence sent to any address other than the one above will not be considered as part of your official premarket notification submission. Also, please note the new Blue Book Memorandum regarding Fax and E-mail Policy entitled, "Fax and E-Mail Communication with Industry about Premarket Files Under Review". Please refer to this guidance for information on current fax and e-mail practices at [www.fda.gov/cdrh/ode/a02-01.html](http://www.fda.gov/cdrh/ode/a02-01.html).

You should be familiar with the manual entitled, "Premarket Notification 510(k) Regulatory Requirements for Medical Devices" available from DSMICA. If you have other procedural or policy questions, or want information on how to check on the status of your submission, please contact DSMICA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsmamain.html> or me at (301)594-1190.

Sincerely yours,

Marjorie Shulman  
Supervisory Consumer Safety Officer  
Office of Device Evaluation  
Center for Devices and Radiological Health

Food and Drug Administration  
Center for Devices and  
Radiological Health  
Office of Device Evaluation  
Document Mail Center (HFZ-401)  
9200 Corporate Blvd.  
Rockville, Maryland 20850

July 09, 2003

POLYGANICS BV  
C/O VIKTOR J. NICKOLSON  
33 RIVER ST. UNIT 813  
HOBOKEN, NJ 07030  
ATTN: VIKTOR J. NICKOLSON

510(k) Number: K032115  
Received: 09-JUL-2003  
Product: NEUROLAC NERVE GUIDE  
User Fee ID Number: 83121-15/03,  
NG01-020/03,  
NG01-025/03,

The Food and Drug Administration (FDA) Center for Devices and Radiological Health (CDRH), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

The Act, as amended by the Medical Device User Fee and Modernization Act of 2002 (MDUFMA) (Public Law 107-250), specifies that a submission shall be considered incomplete and shall not be accepted for filing until fees have been paid (Section 738(f)). Our records indicate that you have not submitted the user fee payment information and therefore your 510(k) cannot be filed and has been placed on hold. The payment information we need in order to begin the review of your 510(k) includes, the user fees cover sheet with the payment ID faxed to the Office of Financial Management at (301) 827-9213 and a check mailed to:

By Regular Mail

By Private Courier (e.g., Fed Ex, UPS, etc.)

-----  
Food and Drug Administration  
P.O. Box 956733  
St. Louis, MO 63195-6733.

-----  
U.S. Bank  
956733  
1005 Convention Plaza  
St. Louis, MO 63101  
(314) 418-4983

The check should be made out to the Food and Drug Administration referencing the payment identification number, and a copy of the User Fee Cover sheet should be included with the check. A copy of the Medical Device User Fee Cover Sheet should also be faxed to CDRH at (301) 594-2977 referencing the 510(k) number if you have not already sent it in with your 510(k) submission. After the FDA has been notified of the receipt of your user fee payment, your 510(k) will be filed and the review will begin. If payment has not been received within 30 days, your 510(k) will be deleted from the system. Additional information on user fees and how to submit your user fee payment may be found at <http://www.fda.gov/oc/mdufma>.

Please note that since your 510(k) has not been reviewed, additional information may be required during the review process and the file may be placed on hold once again. If you are unsure as to whether or not you need to file an application with FDA or what type of application to file, you should first telephone the Division of Small Manufacturers, International and Consumer Assistance (DSMICA), for guidance at (301)443-6597 or its toll-free number (800)638-2041, or contact them at their Internet address <http://www.fda.gov/cdrh/dsmamain.html>, or you may submit a 513(g) request to the Document Mail Center at the address above. If you have any questions concerning the contents of this letter, you may contact me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman  
Consumer Safety Officer  
Office of Device Evaluation  
Center for Devices and  
Radiological Health



Polyganics BV  
L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands

Telephone (31) 50 588 65 88  
Telefax (+31) 50 588 65 99  
E-mail [mail@polyganics.com](mailto:mail@polyganics.com)  
Internet [www.polyganics.com](http://www.polyganics.com)

United States Agent Notification  
Information processing and Office Automation Branch  
Office of Compliance  
Center for Devices and Radiological Health  
9200 Corporate Blvd, HFZ-308  
Rockville, MD 20850-4015  
USA

RECEIVED  
2003 JUL -9 A 11: 09  
FDA/CDRH/ODE/PMO

Our ref. : 03-029/PG002a Neurolac  
Subject : U.S. Agent Notification

Groningen, July 3, 2003

Dear Sir, Madam,

I am providing the following information to comply with the U.S. agent requirement for foreign establishments.

Establishment information

Name: Polyganics BV.  
Registration #: Not yet registered, 510k process ongoing.  
Street address: L.J. Zielstraweg 1  
City: Groningen  
Zip code: 9713 GX  
Country: The Netherlands  
Official Correspondent: Jan Bart Hak, Ph.D.  
Phone: +31 50 588 6588  
Fax: +31 50 588 6599  
E-mail: [hak@polyganics.com](mailto:hak@polyganics.com)  
Web: [www.polyganics.com](http://www.polyganics.com)

All our orders and deliveries are subject to our General Terms and Conditions, registered at the Chamber of Commerce in Groningen under number 02067151. On request a copy will be sent free of charge.


Bank account  
ABN AMRO Groningen 51.35.29.195  
Commercial Reg. No. 02067151  
VAT No. NL80.85.30.525.B01

39  
SIFA  
TF-30  
11

United States Agent Information

Name: Viktor J. Nickolson  
Street address: 33 River St. Unit 813  
City: Hoboken  
State: NJ  
Zip code: 07030  
Country: United States of America  
Phone: +1 201 656 4123  
E-mail: [victorj.nickolson@mac.com](mailto:victorj.nickolson@mac.com)

Yours sincerely,

A handwritten signature in black ink, appearing to read 'Jan Bart Hak', written over a horizontal line.

Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory affairs

**510(k)  
PREMARKET NOTIFICATION  
TABLE OF CONTENTS**

<b><u>Section</u></b>	<b><u>Description</u></b>	<b><u>Page</u></b>
FDA ADMINISTRATIVE FORMS .....		3
	Submission Cover Sheet .....	4
	Indications for Use Form .....	8
APPLICANT STATEMENTS .....		9
	Truthful And Accurate Statement .....	10
	Substantial Equivalence Terminology Statement.....	11
Section 1: GENERAL INFORMATION .....		12
	Applicant .....	12
	Contact Person .....	12
	Device Name.....	12
	Device Classification .....	12
	Manufacturing Facility .....	12
	Sterilization Facility (EtO).....	13
	Performance Standards / Special Controls .....	13
	Purpose of Premarket Notification .....	13
	Predicate Devices .....	13
	Software Validation and Certification .....	13
	Kit Certification .....	13
	Class III Certification .....	14
	510(k) Summary.....	14
Section 2: DEVICE DESCRIPTION.....		15
	Drawings .....	15
	Intended Use.....	15
	Device Description .....	15
Section 3: PROPOSED LABELING.....		17
	Subject Device Labeling.....	17
	Intended Use of the Subject Device.....	17
	Promotional Materials .....	17
Section 4: COMPARATIVE INFORMATION .....		18
	Background .....	18
	Intended Use.....	18
	Subject Device .....	18
	Predicate Device .....	18
	Device Characteristics .....	18
Section 5: PERFORMANCE VERIFICATION .....		20
	Performance Verification .....	20
	Suture retention testing .....	21
	In vitro degradation testing.....	22
	Experimental details.....	22
	Nerve function recovery: sciatic nerve model .....	24
	Neurolac nerve guide versus autologous nerve graft.....	25

---

Functional nerve recovery after bridging a 15 mm nerve gap with a biodegradable nerve guide .....	26
Section 6: BIOCOMPATIBILITY .....	28
Background .....	28
Testing .....	28
Section 7: STERILIZATION AND PYROGENICITY INFORMATION .....	29
Introduction .....	29
Sterilization Method.....	29
Validation Method .....	29
Sterility Assurance Level.....	29
Pyrogen Test Method.....	29
Section 8: PACKAGING AND SHELF-LIFE INFORMATION.....	30
Introduction .....	30
Packaging Description .....	30
Product Shelf – Life.....	30
APPENDIX A: 510(K) SUMMARY OF SAFETY AND EFFECTIVINESS.....	31
APPENDIX B: DEVICE DRAWING .....	34
APPENDIX C: DEVICE LABELING.....	35
APPENDIX D: PREDICATE DEVICE NEUROTUBE™ .....	39
APPENDIX E: PREDICATE DEVICE NEUROGEN™ .....	40
APPENDIX F: PERFORMANCE TEST REPORTS .....	41
APPENDIX G: BIOCOMPATIBILITY DATA .....	47

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION  <b>MEDICAL DEVICE USER FEE COVER SHEET</b>	PAYMENT IDENTIFICATION NUMBER: (b)(4)  Write the Payment Identification Number on your check.
---	---

**See Instructions Before Completing This Cover Sheet**

A completed cover sheet must accompany each original premarket application or supplement listed in Box 3 of this cover sheet. Other premarket application types do not require the use of this cover sheet; see list in the instructions. Payment instructions and fee rates can be found at the following website: <http://www.fda.gov/oc/mdufma>. The following three actions must be taken to properly submit your premarket application and fee payment:

- FAX a copy of this completed cover sheet to the Food and Drug Administration at (301) 827-9213 before payment is sent.
- Include a copy of this completed cover sheet with the check made payable to the Food and Drug Administration and mail them to the Food and Drug Administration, P.O. Box 956733, St. Louis, MO 63195-6733. (Note: in no case should payment be submitted with the premarket application.) Also remember that the Payment Identification Number must be written on the check.  
If you prefer to send a check by a courier, the courier may deliver the check and cover sheet to: US Bank, 956733, 1005 Convention Plaza, St. Louis, MO 63101. (Note: This address is for courier delivery only. Contact the US Bank at 314-418-4821 if you have any questions concerning courier delivery.)
- Include a copy of this completed cover sheet in volume one of the premarket application when submitting to the Food and Drug Administration at either the CBER or CDRH Document Mail Center.

1. COMPANY NAME AND ADDRESS (Include name, street address, city, state, country, and post office code)  POLYGANICS BV L.J. ZIELSTRAWEG 1 GRONINGEN, 9713 GX NL  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	2. CONTACT NAME JAN-BART HAK  2.1 E-MAIL ADDRESS hak@polyganics.com  2.2 TELEPHONE NUMBER (Include area code) +31 50 588 6588  2.3 FACSIMILE (FAX) NUMBER (Include area code) +31 50 588 6599
---	---

3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <http://www.fda.gov/oc/mdufma> )

Select an application type: <input checked="" type="checkbox"/> Premarket notification (510(k)); except for third party reviews <input type="checkbox"/> Biologic License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR)	3.1 Select one of the types below: <input checked="" type="checkbox"/> Original Application  Supplement Types: <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)
--	--

4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status.)

YES, I meet the small business criteria and have submitted the required qualifying documents to FDA  NO, I am not a small business

4.1 If Yes, please enter your Small Business Decision Number:

5. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.

<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates, parents, and partner firms	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population
<input type="checkbox"/> This biologic application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially

6. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA).)

YES  NO

7. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION

(b)(4)



<b>Submission Cover Sheet</b>				
Date of Submission: <b>July 3, 2003</b>			FDA Document Number:	
<b>Section A</b>		<b>Type of Submission</b>		
<b>PMA</b>	<b>PMA Supplement</b>	<b>PDP</b>	<b>510(k)</b>	<b>Meeting</b>
<input type="checkbox"/> Original Submission <input type="checkbox"/> Modules Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<input type="checkbox"/> Regular <input type="checkbox"/> Special <input type="checkbox"/> Panel Track <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA Supplement	<input type="checkbox"/> Presubmission summary <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of intent to start clinical trials <input type="checkbox"/> Intention to submit Notice of Completion <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP <input type="checkbox"/> Report	<input checked="" type="checkbox"/> Original Submission <input checked="" type="checkbox"/> <u>Traditional</u> <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated <input type="checkbox"/> Additional Information <input type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated	<input type="checkbox"/> Pre-IDE meeting <input type="checkbox"/> Pre-PMA meeting <input type="checkbox"/> Pre-PDP meeting <input type="checkbox"/> 180-day meeting <input type="checkbox"/> Other (specify):
<b>IDE</b>	<b>Humanitarian Device Exemption</b>	<b>Class II Exemption</b>	<b>Evaluation of Automatic Class III Designation</b>	<b>Other Submission</b>
<input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<input type="checkbox"/> Original submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report	<input type="checkbox"/> Original submission <input type="checkbox"/> Additional information	<input type="checkbox"/> Original submission <input type="checkbox"/> Additional information	Describe submission:
<b>Section B</b>		<b>Applicant or Sponsor</b>		
Company / Institution name: <b>Polyganics BV</b>		Establishment registration number: <b>Not registered.</b>		
Division name (if applicable): <b>Not Applicable</b>		Phone number (include area code): <b>+31 50 588 6588</b>		
Street address: <b>L.J. Zielstraweg 1</b>		FAX number (include area code): <b>+31 50 588 6599</b>		
City: <b>Groningen</b>	State/Province: <b>Groningen</b>	Country: <b>The Netherlands</b>	Zip/Postal Code: <b>9713-GX</b>	
Contact name: <b>Jan-Bart Hak</b>				
Contact title: <b>Manager Clinical and Regulatory Affairs</b>		Contact e-mail address: <b>hak@polyganics.com</b>		
<b>Section C</b>		<b>Submission Correspondent (if different from above)</b>		
Company / Institution name: not applicable		Establishment registration number: not applicable		
Division name (if applicable): not applicable		Phone number (include area code): not applicable		
Street address: not applicable		FAX number (include area code): not applicable		
City: not applicable	State/Province: not applicable	Country: not applicable	Zip/Postal Code: not applicable	
Contact name: not applicable				
Contact title: not applicable		Contact e-mail address: not applicable		



<b>Section D1 Reason for Submission – PMA, PDP, or HDE</b>		
<ul style="list-style-type: none"> <li><input type="checkbox"/> New device</li> <li><input type="checkbox"/> Withdrawal</li> <li><input type="checkbox"/> Additional or expanded indications</li> <li><input type="checkbox"/> Licensing agreement</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in design, component, or specifications:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Software</li> <li><input type="checkbox"/> Color Additive</li> <li><input type="checkbox"/> Material</li> <li><input type="checkbox"/> Specifications</li> <li><input type="checkbox"/> Other (specify below):</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Location Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Manufacturer</li> <li><input type="checkbox"/> Sterilizer</li> <li><input type="checkbox"/> Packager</li> <li><input type="checkbox"/> Distributor</li> </ul> </li> </ul>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Process Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Manufacturing</li> <li><input type="checkbox"/> Sterilization</li> <li><input type="checkbox"/> Packaging</li> </ul> </li> <li><input type="checkbox"/> Other (specify below):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Labeling Change:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Indications</li> <li><input type="checkbox"/> Instructions</li> <li><input type="checkbox"/> Performance characteristics</li> <li><input type="checkbox"/> Shelf Life</li> <li><input type="checkbox"/> Trade Name</li> </ul> </li> <li><input type="checkbox"/> Other (specify below):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Report submissions:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Annual or periodic</li> <li><input type="checkbox"/> Post-approval study</li> <li><input type="checkbox"/> Adverse reaction</li> <li><input type="checkbox"/> Device defect</li> <li><input type="checkbox"/> Amendment</li> </ul> </li> </ul>
<ul style="list-style-type: none"> <li><input type="checkbox"/> Response to FDA correspondence:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Request for applicant hold</li> <li><input type="checkbox"/> Request for removal of applicant hold</li> <li><input type="checkbox"/> Request for extension</li> <li><input type="checkbox"/> Request to remove or add manufacturing site</li> </ul> </li> <li><input type="checkbox"/> Other reason (specify):</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in ownership</li> <li><input type="checkbox"/> Change in correspondent</li> </ul>	
<b>Section D2 Reason for Submission – IDE</b>		
<ul style="list-style-type: none"> <li><input type="checkbox"/> New device</li> <li><input type="checkbox"/> Addition of institution</li> <li><input type="checkbox"/> Expansion/extension of study</li> <li><input type="checkbox"/> IRB certification</li> <li><input type="checkbox"/> Request hearing</li> <li><input type="checkbox"/> Request waiver</li> <li><input type="checkbox"/> Termination of Study</li> <li><input type="checkbox"/> Withdrawal of application</li> <li><input type="checkbox"/> Unanticipated adverse effect</li> <li><input type="checkbox"/> Notification of emergency use</li> <li><input type="checkbox"/> Compassionate use request</li> <li><input type="checkbox"/> Treatment IDE</li> <li><input type="checkbox"/> Continuing availability request</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Correspondent</li> <li><input type="checkbox"/> Design</li> <li><input type="checkbox"/> Informed Consent</li> <li><input type="checkbox"/> Manufacturer</li> <li><input type="checkbox"/> Manufacturing process</li> <li><input type="checkbox"/> Protocol – feasibility</li> <li><input type="checkbox"/> Protocol – other</li> <li><input type="checkbox"/> Sponsor</li> </ul> </li> <li><input type="checkbox"/> Report Submission:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Current investigator</li> <li><input type="checkbox"/> Annual progress</li> <li><input type="checkbox"/> Site waiver limit reached</li> <li><input type="checkbox"/> Final</li> </ul> </li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Response to FDA letter concerning:                             <ul style="list-style-type: none"> <li><input type="checkbox"/> Conditional approval</li> <li><input type="checkbox"/> Deemed approved</li> <li><input type="checkbox"/> Deficient final report</li> <li><input type="checkbox"/> Deficient progress report</li> <li><input type="checkbox"/> Deficient investigator report</li> <li><input type="checkbox"/> Disapproval</li> <li><input type="checkbox"/> Request extension of time to respond to FDA</li> <li><input type="checkbox"/> Request meeting</li> </ul> </li> </ul>
<b>Section D3 Reason for Submission – 510(k)</b>		
<ul style="list-style-type: none"> <li><input type="checkbox"/> New device</li> <li><input type="checkbox"/> Additional or expanded indications</li> <li><input type="checkbox"/> Other reason (specify)</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in technology</li> <li><input type="checkbox"/> Change in design</li> </ul>	<ul style="list-style-type: none"> <li><input type="checkbox"/> Change in materials</li> <li><input type="checkbox"/> Change in manufacturing process</li> </ul>



<b>Section E Additional Information on 510(k) Submissions</b>					
Product codes of devices to which substantial equivalence is claimed:				Summary of, or statement concerning, safety and effectiveness data:	
1 <b>JXI</b>	2	3	4	<input checked="" type="checkbox"/> 510(k) summary attached <input type="checkbox"/> 510(k) statement	
Information on devices to which substantial equivalence is claimed:					
510(k) Number		Trade or proprietary or model name		Manufacturer	
1. #K983007		Neurotube™		Neuroregen L.L.C.	
2. #K011168		NeuroGen™ Nerve Guide		Integra Life Sciences Corp.	
<b>Section F Product Information – Applicable to All Applications</b>					
Common or usual or classification name: Nerve Cuff					
Trade or proprietary or model name				Model number	
Neurolac® Nerve Guide				NG01-015/03	
Neurolac® Nerve Guide				NG01-020/03	
Neurolac® Nerve Guide				NG01-025/03	
Neurolac® Nerve Guide				NG01-030/03	
FDA document numbers of all prior related submissions (regardless of outcome):					
1 N.A.	2	3	4	5	6
Data included in submission: <input type="checkbox"/> Laboratory testing <input type="checkbox"/> Animal trials <input type="checkbox"/> Human trials					
<b>Section G Product Classification – Applicable to All Applications</b>					
Product code: <b>JXI</b>		C.F.R. section: <b>21 CFR 882.5275</b>		Device class: <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
				Device class: <input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
				Device class: <input type="checkbox"/> Class I <input type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified	
Classification panel: <b>Neurology (Neurological Therapeutic Devices)</b>					
<b>Indications (from labeling):</b> The Neurolac nerve guide is indicated for the reconstruction of peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve. The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.					


## **FDA ADMINISTRATIVE FORMS**

Premarket Submission Cover Sheet

Indications for Use Form



			FDA Document Number:		
<b>Section H Manufacturing / Packaging / Sterilization Sites</b>					
<input checked="" type="radio"/> Original <input type="radio"/> Add <input type="radio"/> Delete		<input checked="" type="radio"/> Manufacturer <input type="radio"/> Contract manufacturer		<input type="radio"/> Contract sterilizer <input type="radio"/> Repackager / relabeler	
Company/institution name: <b>Polyganics BV</b>			Establishment registration number: <b>Not registered</b>		
Division name (if applicable):			Phone number (include area code): <b>+31 50 588 6588</b>		
Street address: <b>L.J. Zielstraweg 1</b>			Fax number (include area code): <b>+31 50 588 6599</b>		
City: <b>Groningen</b>		State/Province: <b>Groningen</b>		Country: <b>The Netherlands</b>	
				Zip/Postal Code: <b>9713-GX</b>	
Contact Name: <b>(see Section B)</b>					
Contact Title: <b>(see Section B)</b>			Contact e-mail address: <b>(see Section B)</b>		
<input checked="" type="radio"/> Original <input type="radio"/> Add <input type="radio"/> Delete		<input type="radio"/> Manufacturer <input type="radio"/> Contract manufacturer		<input checked="" type="radio"/> Contract sterilizer (EtO) <input type="radio"/> Repackager / relabeler	



<input type="radio"/> Original <input type="radio"/> Add <input type="radio"/> Delete		<input type="radio"/> Manufacturer <input type="radio"/> Contract manufacturer		<input type="radio"/> Contract sterilizer <input type="radio"/> Repackager / relabeler	
Company/institution name:			Establishment registration number:		
Division name (if applicable):			Phone number (include area code): ( )		
Street address:			Fax number (include area code): ( )		
City:		State/Province:		Country:	
				Zip/Postal Code:	
Contact Name:					
Contact Title:			Contact e-mail address:		



---

**Indications for Use Form**

510(k) Number: \_\_\_\_\_

Device Name: **Neurolac® Nerve Guide**

**Indications for Use:**

The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use \_\_\_\_\_  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_

(Optional Format 1-2-96)

\_\_\_\_\_  
(Division Sign-Off)

510(k) Number \_\_\_\_\_

## **APPLICANT STATEMENTS**

Truth and Accuracy Certification

Substantial Equivalence Terminology Statement


---

## PREMARKET NOTIFICATION

### Truthful And Accurate Statement


[As Required by 21 CFR 807.87(k)]

I certify that, in my capacity as a Regulatory Affairs representative of Polyganics BV, I believe to the best of my knowledge, that all data and information submitted in the pre-market notification are truthful and accurate and that no material fact has been omitted.



---

J.B. Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands



---

Date

---

Premarket Notification [510(k)] Number

---

## **Substantial Equivalence Terminology Statement**

### **USE OF THE TERM “SUBSTANTIALLY EQUIVALENT**

The use of the term “substantially equivalent” as used herein is intended to be a determination of substantial equivalence under the Federal Food, Drug, and Cosmetic Act, as Amended, and relates to the fact that the product can be marketed without premarket approval or reclassification. Such a determination is not intended to have any bearing whatsoever on the resolution of patent infringement suits or other patent matters.

---

## Section 1: GENERAL INFORMATION

**Applicant** Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

FDA Establishment Registration Number: Not Registered yet

---

**Contact Person** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

Tel: +31 50 588 6588  
Fax: +31 50 588 6599  
E-mail: [hak@polyganics.com](mailto:hak@polyganics.com)

---

Device Name	Trade Name	Common Name	Classification Name
	Neurolac Nerve Guide	Nerve Guide	Nerve Cuff 21 CFR 882.5275

---

**Device Classification**

**Classification Name:** Classification name is "Nerve Cuff" (Ref. Codes of Federal regulations, title 21 – Food and Drugs, Part 882 – Neurological devices, subpart F – Neurological Therapeutic devices, Sec 882.5275)

**Class:** Class II

**Description:** (a) Identification. A nerve cuff is a tubular silicone rubber sheath used to encase a nerve for aid in repairing the nerve (e.g., to prevent ingrowth of scar tissue) and for capping the end of the nerve to prevent the formation of neuroma (tumors).

---

**Manufacturing Facility** Polyganics BV  
L.J. Zielstraweg 1  
9713-GX Groningen  
The Netherlands

FDA Establishment Registration Number: not registered

---



**Sterilization Facility (EtO)**

(b)(4)

**Performance Standards / Special Controls**

No performance standards are indicated for this product.

**Purpose of Premarket Notification**

The reason for this premarket notification is to inform the Food and Drug Administration (FDA) of our intent to market the Polyganics Neurolac Nerve Guide.

**Predicate Devices**

The following table provides information on the predicate devices.

Trade Name	Manufacturer	510k #	Concurrence Date	Substantial Equivalence
Neurotube™	Neuroregen L.L.C.	K983007 Traditional APPENDIX D	03/22/1999	The tube provides an optimal environment for longitudinal nerve axon growth of the peripheral nerve. For single use only in patients with a peripheral nerve injury where the nerve gap is more than or equal to 8 mm, but less than or equal to 3 cm.
NeuroGen™ Nerve Guide	Integra Life Sciences Corp.	K011168 Traditional APPENDIX E	06/22/2001	NeuroGen™ Nerve guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**Software Validation and Certification**

This section is not applicable, since the Neurolac Nerve Guide does not utilize software in the performance of its intended use.

**Kit Certification**

This section is not applicable.

---

**Class III Certification**

This section is not applicable. The Neurolac Nerve Guide is a Class II device

---

**510(k) Summary**

A summary of safety and effectiveness for the Neurolac Nerve Guide is provided in APPENDIX A.


## Section 2: DEVICE DESCRIPTION

**Drawings** Please refer to APPENDIX B of this 510(k) premarket notification for an engineering drawing of the subject device.

**Intended Use** The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.  
The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

For FDA administrative purposes, the intended use of the Polyganics Neurolac Nerve Guide is also documented in a separate form that can be found at the beginning of the 510(k) notification in the section entitled "FDA Administrative Forms"

**Device Description** The Neurolac nerve guide is composed of the bioresorbable copolyester poly(DL-lactide- $\epsilon$ -caprolactone). The Neurolac nerve guide provides guidance and protection to regenerating axons, and prevent ingrowth of fibrous tissue into the nerve gap during nerve regeneration from the proximal to the distal nerve stump of the transected nerve.



**Figure 1.** Example of the Neurolac nerve guide

The Neurolac nerve guide elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the tube by fibrous tissue. Degradation of the Neurolac nerve guide occurs through hydrolysis leading to gradual reduction of molecular weight. The Neurolac nerve guide retains its initial mechanical properties up to 8 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and  $\omega$ -hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. Animal studies demonstrated that a Neurolac nerve guide is resorbed within 16 months.

The Neurolac nerve guide inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The Neurolac nerve guide is indi-

---

cated for single-use.

The Polyganics Nerve Guide further consists of:

- Packaging
- Labeling
- IFU

---

## Section 3: PROPOSED LABELING

**Subject Device Labeling** The following proposed labeling for the subject device Neurolac Nerve Guide is provided in APPENDIX C:

- Outer label (Carton)
  - Inner label (Pouch)
  - Pre-printed carton text and graphics
  - Instructions for Use
- 

**Intended Use of the Subject Device** The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

---

**Promotional Materials** At present, no promotional materials are available for the subject device.

## Section 4: COMPARATIVE INFORMATION

**Background** The Polyganics Neurolac nerve guide is a biodegradable tube for the repair of transected peripheral nerves.

**Intended Use** The Indications for Use for the subject and predicate devices are described in the table below.

Subject Device	Indication for Use
Neurolac	The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve. The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

Predicate Device	Indication for Use
Neurotube™ APPENDIX D	The Neurotube is intended for single use in patients with an injury to a peripheral nerve, in which the nerve gap is more than or equal to 8 mm but less than or equal to 3 cm. The nerve gap may be created primarily at the time of injury or created secondarily at the time of exploration of failed primary repair
NeuroGen™ APPENDIX E	NeuroGen Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**Device Characteristics** The technological characteristics (i.e., design, dimensions material etc.) of the subject Neurolac Nerve Guide and the predicate devices are presented in the table below. The table provides a comparison, demonstrating that the Neurolac is substantially equivalent to the currently marketed predicate device.

**Device Characteristics of the Subject and Predicate devices**

Characteristics General	Subject Device	Predicate Devices	
	Neurolac®	Neurotube™	NeuroGen™
<b>510(k) Reference</b>	<b>This 510(k)</b>	<b>K983007</b>	<b>K011168</b>
<b>Sterile</b>	Sterile device	Sterile device	Sterile device
<b>Single Use</b>	Single-use	Single-use	Single-use
<b>Contents packaging</b>	Nerve guide and instructions for use	Nerve guide and instructions for use	Nerve guide and instructions for use
<b>Length</b>	(b)(4)		
<b>Inner Diameter (In.)</b>	(b)(4)		
<b>Material</b>	(b)(4)		
<b>Biodegradable</b>	Yes	Yes	Yes
<b>Animal derived</b>	No	Yes	No
<b>Indication</b>	Peripheral Nerve discontinuity	Peripheral Nerve discontinuity	Peripheral Nerve discontinuity
<b>Transparent</b>	Yes	No	No


## Section 5: PERFORMANCE VERIFICATION

**Performance Verification** The performance testing on the Neurolac Nerve Guide was conducted to verify that the meets performance characteristics. The following performance tests were conducted:


	Subject/Test	Report
1	(b)(4)	
2		
3		
4		
5		

Testing was performed on finished sterile products. The test results were all favorable for the Neurolac Nerve guide. The complete performance test reports have been attached as APPENDIX F.


(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4) Testing




(b)(4) Testing



<sup>1</sup> SFI = 0: normal; SFI = -100: total impairment.

(b)(4) Testing




---

## Section 6: BIOCOMPATIBILITY

Background

(b) (4)



Testing

---

## Section 7: STERILIZATION AND PYROGENICITY INFORMATION

Introduction

(b) (4)

Sterilization  
Method

Validation  
Method

Sterility As-  
surance Level

Pyrogen Test  
Method

---

## Section 8: PACKAGING AND SHELF-LIFE INFORMATION

**Introduction** The Neurolac nerve guides are packaged in a tray (mechanical protection) and subsequently sealed in a Tyvek/PET-PE pouch. The packaging does comply with standards ISO11607.

---

**Packaging Description** The Neurolac nerve guides are packaged in a tray serving as a mechanical protection.

The tray together with a Neurolac nerve guide are packaged in a Tyvek/PET-PE pouch. The pouch carries a pouch label on the PET-PE layer of the pouch.

The pouch with product in the tray is packaged together with the instructions for use in a carton box. The carton box carries the carton box label.

---

**Product Shelf – Life** The Neurolac nerve guide has a shelf-life (Use By date) of 12 months.

---

## **APPENDIX A: 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS**

---

**510(k)**  
**Summary of Safety and Effectiveness**

---

**Submitter:** Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands  
[www.polyganics.com](http://www.polyganics.com)

**Contact Person:** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Tel : +31 50 588 6588  
Fax : +31 50 588 6599  
Mobile : +31 653 211 303  
E-mail : [hak@polyganics.com](mailto:hak@polyganics.com)

**Date Prepared:** May 20, 2003

---

**General Provisions:** Trade Name: Neurolac® Nerve guide  
Common Name: Nerve guide  
Classification Name: Nerve Cuff, 21 CFR 882.5275  
Device Classification: Class II

---

**Predicate Devices:**

• Neurotube™	Neuroregen L.L.C.	K983007
• NeuroGen™	Integra Life Sciences Corp.	K011168

---

**Performance Standards** For the Nerve Cuff performance, the FDA, under section 514 of the Food, Drug and Cosmetic Act, has not established standards.

---

**Indications for Use** The Neurolac® nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.  
The use of a Neurolac nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

---

---

**Device Description**

Neurolac® is designed to be a flexible and transparent resorbable poly(DL-lactide-co-ε-caprolactone) tube to provide a protective environment for peripheral nerve regeneration after injury and to create a conduit to guide axonal growth across a nerve gap.

Neurolac® nerve guides are provided sterile in Tyvek pouch packages and retainer in a variety of sizes.

---

**Performance Data:**

The safety and effectiveness of the Neurolac nerve guides have been demonstrated via data collected from design verification tests and analyses. The design verification testing consisted of the following:

- In vitro suture retention testing
- In vitro degradation testing
- In vivo nerve function recovery

---

**Summary of Substantial Equivalence**


The design, fundamental technology and intended use (safety and efficacy) featured with the Neurolac® Nerve Guide are substantially equivalent to those featured with the competitor devices Neurotube™ (ref. 510(k) 983007; Neuroregen L.L.C.) and the NeuroGen™ Nerve Guide (ref. 510(k) 011168; Integra Life Sciences Corporation).

Biocompatibility, mechanical and physical property testing, in vitro degradation testing, and performance testing in an animal model provide reasonable scientific evidence that Neurolac® nerve guide is substantially equivalent to the predicate devices. Evaluation of the Polyganics Neurolac® Nerve guide based on biocompatibility testing, animal tests, results from literature and the comparison of the Neurolac® nerve guide with its predicate devices, shows that the Neurolac® nerve guide is safe for implantation.

---

## APPENDIX B: DEVICE DRAWING

(b) (4)




## **APPENDIX C: DEVICE LABELING**

**Subject device: Neurolac Nerve Guide**

**DRAFT**

## DEVICE LABELING

### Carton, DRAFT



-  Sterile
-  Sterilized with ethylen oxide gas
-  Lot No.
-  Use By
-  For one use only
-  Attention, see Instructions for Use
-  Catalog No.

Nonpyrogenic. Do not use open or damaged packages. Store in a cool, dry place at or below 4 °C. **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician.




L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands



## DEVICE LABELING

### Inner label pouch, DRAFT



Bioresorbable peripheral nerve guide  
Bioresorbierbaren peripheren nervenleitschiene  
Bioresorbeerbare perifere zenuwgeleider  
Guide de nerf périphérique biorésorbable  
Guía del nervio periférico bioabsorbible

STERILE  
STERILE EO  
REF NG01-020/03  
LOT NGA2003 04 2501  
2004-04

Nonpyrogenic. Do not use open or damaged packages. Store in a cool, dry place at or below 4 °C. **Caution:** Federal (USA) law restricts this device to sale by or on the order of a physician



POLYGANICS

L.J. Zielstraweg 1  
9713 GX Groningen  
The Netherlands

## **DEVICE LABELING**

**Instructions for Use, DRAFT**

**PG002a 2003-01-09 MM024 IFU EN version 19**



**Instructions for Use, English**

**STERILE. Sterilized with ethylene oxide gas. For single use only. Do not autoclave.**  
**CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.**

**Neuroloc® peripheral nerve guide**

The Neuroloc® nerve guide is composed of the bioresorbable copolyester poly(DL-lactide-ε-caprolactone). The Neuroloc nerve guide provides guidance and protection to regenerating axons, and prevents ingrowth of fibrous tissue into the nerve gap during nerve regeneration from the proximal to the distal nerve stump of a transected nerve.

The Neuroloc nerve guide elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the tube by fibrous tissue. Degradation of the Neuroloc nerve guide occurs through hydrolysis leading to gradual reduction of molecular weight. The Neuroloc nerve guide attains its initial mechanical properties up to 8 weeks, whereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and ω-hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. Animal studies demonstrated that a Neuroloc nerve guide is absorbed within 16 months.

The Neuroloc nerve guide inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The Neuroloc nerve guide is indicated for single-use.

**indications**

The Neuroloc nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

**Contraindications**

The use of a Neuroloc nerve guide is contraindicated in patients with known hypersensitivity or allergies to its components.

**Warnings**

The Neuroloc nerve guide is for single use only. Do not resterilize or re-use. Structural integrity and/or function may be impaired through cleaning, resterilization, or re-use and may cause adverse patient reactions. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from re-use of (or any part of) the Neuroloc nerve guide.

Sterile unless package has been opened or damaged. Discard open unused nerve guides.

The Neuroloc nerve guide should only be used by physicians who are trained in nerve defect repair techniques. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel. The physician should consult recent literature on current medical practice on peripheral nerve repair.

Persons with allergic reactions to the biodegradable material or breakdown products may suffer an allergic response to this implant.

Nerve regeneration may be suboptimal in elderly, malnourished or debilitated patients or in patients suffering from cancer, anaemia, obesity, diabetes, infection or other conditions which may delay wound healing, infected wounds, or moderate tissue inflammatory response characteristic of foreign body response.

**Precautions**

- Use prior to "Use by date";
- Store in dark, dry place at or below 4°C (39°F);
- Do not expose the nerve guide to organic solvents (e.g. chloroform, acetone);
- Do not use absorbable sutures for fixation of the nerve stumps into the nerve guide;
- Avoid crushing, crimping, kinking or other damage due to application of surgical instruments such as forceps, needles, and scissors or during handling of the device;
- Avoid tension on the nerve ends;
- Prevent compression and/or kinking of the Neuroloc after the procedure. The use of a protective splint is recommended.

**Adverse effects**

- Adverse events associated with the use of a Neuroloc® nerve guide may include but are not limited to:
  - Failure to provide adequate nerve regeneration at sites where too much tension or compression occurs;
  - Failure to provide adequate/complete nerve regeneration
  - Transitory/local irritation
  - Infection
  - Allergy
  - Delayed wound healing

**Opening of the package**

The pouch is opened in such a way that the tray remains sterile. The tray can be opened by sliding the lid. By clamping the nerve guide at one of its ends between a pair of tweezers, it can be taken from the tray. The lid contains a ruler that may be used as a reference to estimate the gap length or nerve stump diameter.

**Surgical Procedure**

1. Surgically expose the injured nerve.
  2. Resect the injured segment distally and proximally until a nerve stump is identified with no residual intrafascicular scarring.
- NOTE: Do not crush the nerve stumps as this can cause extrusion of intra-fascicular components.**
3. Measure the length of the defect with all joints in an extended position.
  4. If the gap length is between 0 and 20 mm, the injured nerve can be reconstructed with a Neuroloc nerve guide.
  5. Select the Neuroloc nerve guide with the proper internal diameter.
- NOTE: It is essential that the internal nerve guide diameter is slightly larger than the diameter of the transected nerve to guarantee optimal nerve regeneration.**
6. Cut the selected nerve guide with a pair of scissors or a knife so that the nerve guide is 1 cm longer than the nerve gap.

Under some circumstances immobilization of the nerve ends, as to avoid tension on the nerve ends, may be employed at the discretion of the surgeon. To secure adequate fixation of the nerve ends in the nerve guide, the accepted surgical technique of flat, square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon, is required.

**Suturing Technique**

1. Place the Neuroloc nerve guide in warm saline (37°C) for approximately 1 minute before implantation. This will make the tube more flexible and ease needle passage during suturing.
2. Suture the Neuroloc nerve guide by passing the suture (8-0 suture) first through the tube from the outside to the inside and then transversely and superficially through the epi-

neurium and back through the tube from the inside to the outside, after which a tie is made (Fig. 1.1-1.3).

When positioning optimization of the nerve ends in the nerve guide is required, it is recommended to place a second suture in the same nerve end (Fig. 1.4).

Pull the proximal nerve stump into the nerve guide.  
**NOTE: It is recommended that the nerve ends are pulled into the tube for at least 3 mm for optimal nerve regeneration.**

5. Fill the tube with heparinized saline, using a solution containing 1000 units of heparin per 100 ml of normal saline (Fig 1.5).
6. Subsequently, use the same procedure, to pull the distal nerve stump into the nerve guide.
7. A minimum space of 5 mm should be left between the nerve ends in the nerve guide.
8. Fill any remaining space with heparinized saline (Fig 1.6) by injecting along the nerve into the lumen of the tube or by penetrating the tube (not the nerve).




Figure 1. Schematic representation of suture technique for suturing the nerve ends into the nerve guide.

**CAUTION: Ensure that no blood enters the nerve guide lumen since this may hinder nerve recovery.**  
**CAUTION: The nerve guide should be implanted and sutured with all joints in an extended position as to assure that no tension occurs on the proximal or distal nerve end when joints are being mobilized.**

9. Close the wound and splint to prevent kinking for the first three postoperative weeks. Long-term compression of the nerve guide should be avoided. Patients may be administered oral antibiotics for the first post-operative week.

Dispose contaminated implantation and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.



= Refer to accompanying instructions for use.



= Use by

= Catalog number

= For one use only

= Lot number

= Sterile product

= Sterilized by ethylene oxide



**Polyganics BV** Telephone +31 50 588 6588  
 L.J. Zielstraweg 1 Fax +31 50 588 6599  
 9713 GX Groningen  
 The Netherlands

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**  
 There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Polyganics product(s) described in this publication. Under no circumstances shall Polyganics BV be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Polyganics to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Polyganics BV printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Polyganics BV will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Neuroloc® is a registered trademark of Polyganics

73

---

## **APPENDIX D: PREDICATE DEVICE NEUROTUBE™**

**Neuroregen L.L.C.**

**510(k): K983007**

3/22/99

K983007

PREMARKET NOTIFICATION [510(K)] SUMMARY  
(as required by 21CFR 807.92(a))

Submitted by: John E. Barham, Managing Member, Neuroregen L.L.C.  
43 N. Bond Street, Bel Air, MD 21014  
Phone 410 838-8090 Fax 410 838-8092

Date: August 28, 1998

Trade Name: Neurotube<sup>™</sup>

Common Name: Nerve Conduit

Classification Name: Nerve Cuff (per 21 CFR Section 882.5275)

Equivalent Device: Silicone Nerve Cuffs

Description: The Neurotube is a woven, flexible, polyglycolic acid tube which has been heat treated to achieve a configuration corrugated externally for wall strength. The tube is 2.3 mm in diameter and 4 cm in length.

Intended Use: The tube provides an optimal environment for longitudinal nerve axon growth of the peripheral nerve. For single use only in patients with a peripheral nerve injury where the nerve gap is more than or equal to 8 mm, but less than or equal to 3 cm.

Technological Characteristics Compared to Predicate Device: The Neurotube is fabricated from a bioresorbable material in contrast to nonresorbable silicone and thus precludes the need for a second surgery. Both the silicone and bioresorbable products are tubular in design to facilitate nerve regeneration.

Clinical Data: Clinical trials for the Neurotube were conducted in the United States over a period of three and one half years to support a determination of substantial equivalence.. A total of 98 subjects were enrolled at five clinical trial sites. One hundred two nerve reconstructions were evaluated. There were 56 in the control group using classic end-to-end nerve graft repairs and 46 received the Neurotube. Subjects were given sensory evaluations at 3, 6, 9, and 12 month intervals. Static and moving sensory discrimination tests were performed. Results were equivalent between the two groups. The only adverse effects reported were delayed healing of a skin closure and skin separation with partial extrusion of the Neurotube. The study indicated that a single stage, biodegradable, polyglycolic acid conduit can be used as an alternative to a nerve graft or a biodurable nerve tube.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAR 22 1999

Mr. John E. Barham  
Managing Member  
Neuroregen, L.L.C.  
43 North Bond Street  
Bel Air, Maryland 21014

Re: K983007  
Trade Name: Neurotube™  
Regulatory Class: II  
Product Code: JXI  
Dated: December 18, 1998  
Received: December 22, 1998

Dear Mr. Barham:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 – Mr. John E. Barham

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4595. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Celia M. Witten, Ph.D., M.D.  
Director  
Division of General and  
Restorative Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K983007

510 (k) Number K983007

Device Name: Neurotube™

Indications for Use:


The Neurotube is intended for single use in patients with an injury to a peripheral nerve, in which the nerve gap is more than or equal to 8 mm but less than or equal to 3 cm. The nerve gap may be created primarily at the time of injury or created secondarily at the time of exploration of failed primary repair.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

  
(Division Sign-Off)  
Division of General Restorative Devices

510(k) Number

K983007

JUN 22 2001

NeuroGen™ Nerve Guide

510(K) SUMMARY

**Submitter's name and address:**

Integra LifeSciences Corporation  
105 Morgan Lane  
Plainsboro, NJ 08536 USA

**Contact person and telephone number:**

Judith E. O'Grady  
Senior Vice President, Regulatory Affairs, Quality Assurance and Clinical Research  
(609) 275-0500

**Date Summary was prepared:**

April 16, 2001

**Name of the device:**

Proprietary Name: NeuroGen™  
Common Name: Nerve Guide  
Classification Name: Nerve Cuff, Product Code 84JXI

**Substantial Equivalence:**

NeuroGen™ Nerve Guide is substantially equivalent in function and intended use to the following products which have been cleared to market under Premarket Notifications 510(k): Salumedica™ Nerve Cuff, NeuroTube® and Fastube™ Nerve Cuff.

**Device Description:**

NeuroGen™ Nerve Guide is an implant designed for repair of peripheral nerve discontinuities. NeuroGen™ Nerve Guide provides a protective environment for peripheral nerve repair after injury. NeuroGen™ Nerve Guide is designed to be an interface between the nerve and surrounding tissue and to create a conduit for axonal growth across a nerve gap.

NeuroGen™ Nerve Guide is flexible to accommodate movement of joint and associated tendons while retaining its shape and is resistant to occlusive forces from surrounding tissue. When hydrated, NeuroGen™ is an easy to handle, soft, pliable, nonfriable, porous collagen tube. NeuroGen™ Nerve Guides are provided sterile in double blister packages in a variety of sizes.

**Intended Use:**

NeuroGen™ Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**Safety**

Biocompatibility studies have demonstrated NeuroGen™ Nerve Guides to be: noncytotoxic, nonpyrogenic, nonirritating, and nonsensitizing. The following studies were conducted:

- a) Cytotoxicity
- b) Irritation / Intracutaneous Reactivity
- c) Sensitization
- d) Acute Systemic Toxicity
- e) Subchronic Toxicity
- f) Chronic Toxicity
- g) Genotoxicity
- h) Implantation
- i) Hemolysis

**Performance Characteristics:**

The effectiveness of NeuroGen™ Nerve Guide to repair peripheral nerve discontinuities was studied in a long-term (3.5 years) primate study and in rodent animal models. The studies demonstrated that NeuroGen™ Nerve Guides are substantially equivalent to nerve graft, direct suture and silicone tubes.

The mechanical and physical characteristics of the NeuroGen™ Nerve Guides were evaluated in a series of tests. These tests were conducted to ensure that the NeuroGen™ Nerve Guides possess the mechanical properties (suture retention and mechanical compression) as well as physical properties (porosity and permeability) that determine their suitability for use in the human body. Testing has demonstrated that the nerve guides are able to hold a suture, resist repeated compression from surrounding tissues, have a porous outer surface and tube wall, and allow the passage of molecules of specific size through the tube wall.

**Technological Characteristics Compared to Predicate Devices:**

NeuroGen™ Nerve Guide is a tubular device which is equivalent to the predicate devices, Salumedia™ Nerve Cuff, NeuroTube® and Fastube™ Nerve Cuff in its design for repair of peripheral nerve discontinuities. Like the predicate devices, NeuroGen™ is provided sterile, for single use only. The NeuroGen™ Nerve Guide is manufactured from a bioresorbable material, as is one of the predicate devices, NeuroTube®. NeuroGen™ Nerve Guide meets ISO 10993 requirements for Biocompatibility testing.

**Conclusion**

NeuroGen™ Nerve Guide is indicated for the repair of nerve discontinuities where gap closure can be achieved by flexion of the extremity. NeuroGen™ Nerve Guide is flexible to accommodate movement of joint and associated tendons while retaining its shape and is resistant to occlusive forces from surrounding tissue.

Biocompatibility studies have demonstrated NeuroGen™ Nerve Guide to be non-cytotoxic, non-sensitizing, non-toxic and non-mutagenic. Extensive, long-term evaluations in primates demonstrates NeuroGen™ Nerve Guide to be biocompatible and provides an environment for axonal growth.

Valid scientific evidence through substantial testing of descriptive characteristics, Biocompatibility, mechanical and physical property testing and extensive performance testing in a primate model, provide reasonable assurance that NeuroGen™ Nerve Guide is safe and effective under the proposed conditions of use, and substantially equivalent to its predicate devices.



JUN 22 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Judith E. O'Grady, RN, MSN  
Senior Vice President, Regulatory Affairs,  
Quality Assurance and Clinical Affairs  
Integra Life Sciences Corporation  
105 Morgan Lane  
Plainsboro, New Jersey 08536

Re: K011168  
Trade/Device Name: NeuroGen Nerve Guide  
Regulation Number: 882.5275  
Regulatory Class: II  
Product Code: JXI  
Dated: April 16, 2001  
Received: April 17, 2001

Dear Ms. O'Grady:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Judith E. O'Grady, RN, MSN

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

for Mark N. Milkman

Celia M. Witten, Ph.D., M.D.

Director

Division of General, Restorative  
and Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

510(k) Number:

K011168

Device Name: NeuroGen™ Nerve Guide

**Indications for Use**

NeuroGen™ Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)**

(Concurrence of CDRH, Office of Device Evaluation (ODE))

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-the-Counter Use \_\_\_\_\_

*for Mark A. Milburn*  
\_\_\_\_\_  
(Division Sign-Off)

Division of General, Restorative  
and Neurological Devices

510(k) Number \_\_\_\_\_

K011168

## **APPENDIX E: PREDICATE DEVICE NEUROGEN™**

**Integra Life Sciences Corporation**

**510(k): K011168**

JUN 22 2001

## NeuroGen™ Nerve Guide

## 510(K) SUMMARY

**Submitter's name and address:**

Integra LifeSciences Corporation  
105 Morgan Lane  
Plainsboro, NJ 08536 USA

**Contact person and telephone number:**

Judith E. O'Grady  
Senior Vice President, Regulatory Affairs, Quality Assurance and Clinical Research  
(609) 275-0500

**Date Summary was prepared:**

April 16, 2001

**Name of the device:**

Proprietary Name: NeuroGen™  
Common Name: Nerve Guide  
Classification Name: Nerve Cuff, Product Code 84JXI

**Substantial Equivalence:**

NeuroGen™ Nerve Guide is substantially equivalent in function and intended use to the following products which have been cleared to market under Premarket Notifications 510(k): Salumedica™ Nerve Cuff, NeuroTube® and Fastube™ Nerve Cuff.

**Device Description:**

NeuroGen™ Nerve Guide is an implant designed for repair of peripheral nerve discontinuities. NeuroGen™ Nerve Guide provides a protective environment for peripheral nerve repair after injury. NeuroGen™ Nerve Guide is designed to be an interface between the nerve and surrounding tissue and to create a conduit for axonal growth across a nerve gap.

NeuroGen™ Nerve Guide is flexible to accommodate movement of joint and associated tendons while retaining its shape and is resistant to occlusive forces from surrounding tissue. When hydrated, NeuroGen™ is an easy to handle, soft, pliable, nonfriable, porous collagen tube. NeuroGen™ Nerve Guides are provided sterile in double blister packages in a variety of sizes.

**Intended Use:**

NeuroGen™ Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

D0001

## Safety

Biocompatibility studies have demonstrated NeuroGen™ Nerve Guides to be: noncytotoxic, nonpyrogenic, nonirritating, and nonsensitizing. The following studies were conducted:

- a) Cytotoxicity
- b) Irritation / Intracutaneous Reactivity
- c) Sensitization
- d) Acute Systemic Toxicity
- e) Subchronic Toxicity
- f) Chronic Toxicity
- g) Genotoxicity
- h) Implantation
- i) Hemolysis

## Performance Characteristics:

The effectiveness of NeuroGen™ Nerve Guide to repair peripheral nerve discontinuities was studied in a long-term (3.5 years) primate study and in rodent animal models. The studies demonstrated that NeuroGen™ Nerve Guides are substantially equivalent to nerve graft, direct suture and silicone tubes.

The mechanical and physical characteristics of the NeuroGen™ Nerve Guides were evaluated in a series of tests. These tests were conducted to ensure that the NeuroGen™ Nerve Guides possess the mechanical properties (suture retention and mechanical compression) as well as physical properties (porosity and permeability) that determine their suitability for use in the human body. Testing has demonstrated that the nerve guides are able to hold a suture, resist repeated compression from surrounding tissues, have a porous outer surface and tube wall, and allow the passage of molecules of specific size through the tube wall.

## Technological Characteristics Compared to Predicate Devices:

NeuroGen™ Nerve Guide is a tubular device which is equivalent to the predicate devices, Salumedia™ Nerve Cuff, NeuroTube® and Fastube™ Nerve Cuff in its design for repair of peripheral nerve discontinuities. Like the predicate devices, NeuroGen™ is provided sterile, for single use only. The NeuroGen™ Nerve Guide is manufactured from a bioresorbable material, as is one of the predicate devices, NeuroTube®. NeuroGen™ Nerve Guide meets ISO 10993 requirements for Biocompatibility testing.

## Conclusion

NeuroGen™ Nerve Guide is indicated for the repair of nerve discontinuities where gap closure can be achieved by flexion of the extremity. NeuroGen™ Nerve Guide is flexible to accommodate movement of joint and associated tendons while retaining its shape and is resistant to occlusive forces from surrounding tissue.

Biocompatibility studies have demonstrated NeuroGen™ Nerve Guide to be non-cytotoxic, non-sensitizing, non-toxic and non-mutagenic. Extensive, long-term evaluations in primates demonstrates NeuroGen™ Nerve Guide to be biocompatible and provides an environment for axonal growth.

Valid scientific evidence through substantial testing of descriptive characteristics, Biocompatibility, mechanical and physical property testing and extensive performance testing in a primate model, provide reasonable assurance that NeuroGen™ Nerve Guide is safe and effective under the proposed conditions of use, and substantially equivalent to its predicate devices.

D0002



JUN 22 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Judith E. O'Grady, RN, MSN  
Senior Vice President, Regulatory Affairs,  
Quality Assurance and Clinical Affairs  
Integra Life Sciences Corporation  
105 Morgan Lane  
Plainsboro, New Jersey 08536

Re: K011168  
Trade/Device Name: NeuroGen Nerve Guide  
Regulation Number: 882.5275  
Regulatory Class: II  
Product Code: JXI  
Dated: April 16, 2001  
Received: April 17, 2001

Dear Ms. O'Grady:


We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

for Celia M. Witten, Ph.D., M.D.

Director  
Division of General, Restorative  
and Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number: K011168

Device Name: NeuroGen™ Nerve Guide

**Indications for Use**

NeuroGen™ Nerve Guide is indicated for repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity.

**PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)**

(Concurrence of CDRH, Office of Device Evaluation (ODE))  
Prescription Use X OR Over-the-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109)

*for Mark A. Milburn*  
\_\_\_\_\_  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices  
510(k) Number K011168

## **APPENDIX F: PERFORMANCE TEST REPORTS**

**Suture retention testing**

**In vitro degradation testing**

**Nerve function recovery: sciatic nerve model**

**Neurolac nerve guide versus autologous nerve graft**


**Functional nerve recovery after bridging a 15 mm nerve gap with a Neurolac nerve guide**



---

## Suture retention testing

(b)(4) Testing



**CONTENTS**

1 INTRODUCTION..... 3

2 INCENTIVE..... 3

3 OBJECTIVE..... 3

4 EXPERIMENTAL DETAILS ..... 3

    4.1 PREPARATION OF SUTURED NERVE GUIDES ..... 3

    4.2 SCREENING EXPERIMENTS ..... 4


5 RESULTS ..... 4

(b)(4) Testing ..... 4


6 CONCLUSION..... 5

REFERENCES ..... 5


(b)(4) Testing



(b)(4) Testing




(b)(4) Testing




## **In vitro degradation testing**

(b)(4) Testing




(b)(4) Testing




**Contents**

	<b>Summary .....</b>	<b>2</b>
<b>1</b>	<b>Purpose.....</b>	<b>4</b>
<b>2</b>	<b>References.....</b>	<b>4</b>
<b>3</b>	<b>Compliance.....</b>	<b>4</b>
<b>4</b>	<b>Acceptance criteria.....</b>	<b>4</b>
<b>5</b>	<b>Material and methods .....</b>	<b>4</b>
5.1	Test samples .....	4
5.2	Testing schedule .....	5
5.3	Methods.....	5
(b)(4)Testing	.....	5
(b)(4)Testing	.....	5
(b)(4)Testing	.....	6
<b>6</b>	<b>Results and discussion.....</b>	<b>7</b>
6.1	Visual observations .....	7
6.2	Dimensional changes .....	7
6.3	Water uptake and weight loss .....	8
6.3.1	Intrinsic viscosity .....	8
6.3.2	Nuclear Magnetic Resonance (NMR).....	9
6.3.3	Differential Scanning Calorimetry (DSC).....	11
6.4	Tensile testing: mechanical properties.....	12
<b>7</b>	<b>Conclusion.....</b>	<b>13</b>
<b>8</b>	<b>References.....</b>	<b>14</b>
<b>9</b>	<b>APPENDIX 1.....</b>	<b>15</b>


(b)(4) Testing




(b)(4) Testing




(b)(4)Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4) Testing




(b)(4) Testing



## **Nerve function recovery: sciatic nerve model**

(b)(4) Testing




(b)(4) Testing




**Contents**

<b>1</b>	<b>Purpose</b> .....	<b>4</b>
<b>2</b>	<b>References</b> .....	<b>4</b>
<b>3</b>	<b>Compliance</b> .....	<b>4</b>
<b>4</b>	<b>Implant specimens</b> .....	<b>4</b>
4.1	Test article.....	4
4.2	Control articles.....	5
<b>5</b>	<b>Animals and implantation</b> .....	<b>5</b>
5.1	Justification of the test system.....	5
5.2	Justification of the test periods.....	5
5.3	Animals used in the test.....	5
5.4	Animal care and maintenance.....	5
5.5	Surgical implantation.....	5
5.6	Functional assessment.....	6
<b>6</b>	<b>Evaluation</b> .....	<b>7</b>
<b>7</b>	<b>Discussion and conclusions</b> .....	<b>8</b>
	<b>References</b> .....	<b>9</b>


(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4) Testing




(b)(4) Testing




by  
in  
in  
re  
a  
g  
er  
e.  
e  
of  
d  
-  
e  
-  
y  
s  
e  
e

(b)(4) Testing



(b)(4) Testing







---

## Neurolac nerve guide versus autologous nerve graft

(b)(4) Testing



(b)(4) Testing




**Contents**


<b>1</b>	<b>Purpose</b> .....	<b>4</b>
<b>2</b>	<b>References</b> .....	<b>4</b>
<b>3</b>	<b>Compliance</b> .....	<b>4</b>
<b>4</b>	<b>Implant specimens</b> .....	<b>4</b>
4.1	Test article.....	4
4.2	Control groups.....	5
<b>5</b>	<b>Animals and implantation</b> .....	<b>5</b>
5.1	Justification of test system.....	5
5.2	Justification of the test periods.....	5
5.3	Animals used in the test.....	5
5.4	Animal care and maintenance.....	5
5.5	Surgical implantation.....	5
5.6	Functional assessment.....	6
5.6.1	Test schedule.....	6
5.6.2	Electromyography.....	6
5.6.3	Video Analysis of Walking.....	6
<b>6</b>	<b>Evaluation</b> .....	<b>6</b>
6.1.1	EMG.....	6
6.1.2	Evaluation of Gait.....	8
<b>7</b>	<b>Discussion and conclusions</b> .....	<b>9</b>
	<b>References</b> .....	<b>9</b>

(b)(4) Testing


(b)(4)Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4)Testing




(b)(4)Testing




---

## Functional nerve recovery after bridging a 15 mm nerve gap with a Neurolac nerve guide

(b)(4) Testing



(b)(4)Testing




**Contents**


<b>1</b>	<b>Purpose</b> .....	<b>4</b>
<b>2</b>	<b>References</b> .....	<b>4</b>
<b>3</b>	<b>Compliance</b> .....	<b>4</b>
<b>4</b>	<b>Implant specimens</b> .....	<b>4</b>
4.1	Test article.....	4
4.2	Control articles.....	5
<b>5</b>	<b>Animals and implantation</b> .....	<b>5</b>
5.1	Justification of test system.....	5
5.2	Justification of the test periods.....	5
5.3	Animals used in the test.....	5
5.4	Animal care and maintenance.....	5
5.5	Surgical implantation.....	5
5.6	Functional assessment.....	6
5.6.1	Test schedule.....	6
5.6.2	Walking-track analysis.....	6
5.6.3	Withdrawal test.....	6
5.6.4	Conduction velocity (CV).....	7
5.6.5	Electromyography (EMG).....	7
<b>6</b>	<b>Evaluation</b> .....	<b>8</b>
6.1.1	Sciatic function index (SF1).....	8
6.1.2	Stance factor.....	8
6.1.3	Withdrawal test.....	8
6.1.4	Conduction velocity (CV).....	8
6.1.5	Electromyography (EMG).....	8
<b>7</b>	<b>Discussion and conclusions</b> .....	<b>9</b>
	<b>References</b> .....	<b>11</b>




(b)(4)Testing

A large black rectangular redaction box covering the majority of the page content, starting below the header and extending to the bottom of the page.


(b)(4) Testing




(b)(4)Testing




(b)(4) Testing




(b)(4)Testing




(b)(4) Testing



(b)(4) Testing



(b)(4) Testing



---

## **APPENDIX G: BIOCOMPATIBILITY DATA**

**Cytotoxicity**  
**Irritation**  
**Sensitization**  
**Hemocompatibility**  
**Acute systemic toxicity**  
**Pyrogenicity**  
**Mutagenicity /genotoxicity**  
**Sub chronic toxicity**  
**Carcinogenicity**  
**Chronic toxicity**  
**Reproductive toxicity**  
**Implantation**

## Cytotoxicity

























---

## Irritation





























## Sensitization









































































---

## Hemocompatibility







































## Acute systemic toxicity















































## Pyrogenicity

























## **Mutagenicity /genotoxicity**

































































































## **Sub chronic toxicity**




























**Carcinogenicity**

**Chronic toxicity**

**Reproductive toxicity**

(b)(4) Testing



Contents

Study summary ..... 2

1 Purpose ..... 3

2 References ..... 3

3 Compliance ..... 3

4 Device and manufacturing process characteristics ..... 3

    4.1 Test article specifications ..... 3

    4.2 Manufacturing process ..... 3

        4.2.1 Materials ..... 3

        4.2.2 Purification and Polymerization ..... 4

        4.2.3 Dissolving and filtration ..... 4

        4.2.4 Dipcoating ..... 4

        4.2.5 Extraction and drying ..... 4

        4.2.6 Packaging ..... 4

        4.2.7 Sterilization ..... 4

        (b) (4) ..... 4

        ..... 4

        ..... 4

        ..... 4

        ..... 4

        ..... 5

        ..... 5

        ..... 5

        ..... 5

        ..... 5

        ..... 6

        ..... 7

        ..... 7

        ..... 8

        ..... 8

        ..... 8


        ..... 8

        ..... 8


5 Conclusion ..... 9

Study summary


(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4)Testing



(b)(4) Testing



(b)(4) Testing



(b)(4)Testing

## Implantation

(b)(4) Testing

(b)(4) Testing

**Contents**

**Study summary** ..... 2

**1 Purpose** ..... 3

**2 References** ..... 3

**3 Compliance** ..... 3

**4 Implant specimens** ..... 3

    4.1 Test article ..... 3

    4.2 Control articles ..... 4

(b)(4) Testing ..... 4

..... 4

..... 4

..... 4

..... 4

..... 4

..... 5

**6 Evaluation** ..... 5

    6.1 Post-operative assessment of animal behaviour ..... 5

    6.2 Macroscopic appearance of the test specimen ..... 5

    6.3 Macroscopic and histological assessment of surrounding tissue ..... 5


**7 Final evaluation** ..... 6

**8 Literature** ..... 6


(b)(4) Testing

(b)(4) Testing

(b)(4) Testing



(b)(4) Testing




349

(b)(4) Testing



(b)(4) Testing




**Contents**


	<b>Study summary</b> .....	<b>2</b>
<b>1</b>	<b>Purpose</b> .....	<b>3</b>
<b>2</b>	<b>References</b> .....	<b>3</b>
<b>3</b>	<b>Compliance</b> .....	<b>3</b>
<b>4</b>	<b>Implant specimens</b> .....	<b>3</b>
4.1	Test article .....	3
4.2	Control articles .....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
(b) (4)	.....	4
<b>6</b>	<b>Evaluation</b> .....	<b>5</b>
(b) (4)	.....	5
(b) (4)	.....	5
(b) (4)	.....	7
<b>7</b>	<b>Final evaluation</b> .....	<b>10</b>
<b>8</b>	<b>Conclusion</b> .....	<b>11</b>
	<b>References</b> .....	<b>11</b>

**Study summary**


(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4)Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4) Testing




(b)(4) Testing







(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing




(b)(4) Testing



(b)(4) Testing



(b)(4) Testing



From: Reviewer(s) - Name(s) David B. Berkowitz MA  
Subject: 510(k) Number K032115  
To: The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept).
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.
- Other (e.g., exempt by regulation, not a device, duplicate, etc.)

Is this device subject to Section 522 Postmarket Surveillance?  YES  NO  
 Is this device subject to the Tracking Regulation?  YES  NO  
 Was clinical data necessary to support the review of this 510(k)?  YES  NO  
 Is this a prescription device?  YES  NO  
 Was this 510(k) reviewed by a Third Party?  YES  NO  
 Special 510(k)?  YES  NO  
 Abbreviated 510(k)? Please fill out form on H Drive 510k/boilers  YES  NO

Truthful and Accurate Statement  Requested  Enclosed  
 A 510(k) summary OR  A 510(k) statement  
 The required certification and summary for class III devices  
 The indication for use form

Combination Product Category (Please see algorithm on H drive 510k/Boilers) N

Animal Tissue Source  YES  NO Material of Biological Origin  YES  NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):  
 No Confidentiality  Confidentiality for 90 days  Continued Confidentiality exceeding 90

Predicate Product Code with class: \_\_\_\_\_ Additional Product Code(s) with panel (optional): \_\_\_\_\_

JXI Nerve Cuff class II  
 Review: Steph Ruelo PR JB 10/9/03  
 (Branch Chief) (Branch Code) (Date)  
 Final Review: Miriam C. Provost 10/10/03  
 (Division Director) (Date)



**5 1 0 (K) M E M O R A N D U M**

**TO:** K032115

**FROM:** David B. Berkowitz, Veterinarian  
ODE/DGRND/Plastic and Reconstructive Surgery Devices Branch

**DATE:** October 7, 2003

**SUBJ:** Neurolac® Nerve Guide  
Polyganics BV  
Jan Bart Hak, Manager of clinical and Regulatory Affairs (See Aug. 4, letter)  
31 50 588 6588

**Recommendation:** S.E.

Procode: JXI

Class: II

Regulation Number: 882.5275

Regulation Name: Nerve Cuff

**REVIEW:**

Device Description: This is a co-polyester of poly(DL-lactide-ε-caprolactone) tube used to guide regenerating axons and to exclude ingrowth of fibrous tissue during regeneration.

**1. Comparison of the Intended Use/Indications of the Subject Device and Predicate(s)  
Subject Device**

For the reconstruction of a peripheral nerve discontinuity up to 20mm in patients who have sustained a complete division of a nerve.

**Predicate devices**

K983007 Neurotube Neuroregen LLC For single use in patients with a peripheral nerve injury where the nerve gap is more than or equal to 8 mm, but less than 3 cm. The tube is made of polyglycolic acid.

K011168 NeuroGen™ Nerve Guide – For the repair of peripheral nerve discontinuities where gap closure can be achieved by flexion of the extremity. The tube is made of bovine collagen.

**Discussion of whether the intended use/indications are the same**

The K011168 closure is limited to cases in which closure can be achieved by flexion. The current device is not restricted to flexion. I do not know the significance of this restriction. The current device can be used for up to a 20 mm gap, well within the range of the predicates.

**2. Comparison of the Technological Characteristics (Design, Materials, Sizes, Shapes, etc.) of the Subject Device and Predicate(s)**

**Subject Device**

The chemical composition is provided. This is a copolyester of poly(DL-lactide-ε-caprolactone) tube. The tube retains strength for 8 weeks, and then rapidly degrades. (b) (4)

[Redacted]

[Redacted] (b) (4)

**Predicate Devices**

K983007 Neurotube™ Nerve Cuff – Made from polyglycolic acid. 2.3 mm diameter and 4 cm long. Full absorption in 6 months. The corrugated configuration said to prevent collapse of tube from soft tissue pressure.

K0111168 – Made from bovine collagen. 2 – 4 cm, with i.d. 2, 4, 5, 6, and 7 mm.

**Discussion of whether the subject device has a significant change in technological characteristics.**


The materials are different for each device. Other predicates include silicone and collagen tubes. Equivalence rests on the safety and effectiveness similarities.

**3. Comparative Data (in vitro, animal and/or clinical)**


**Safety Data - Subject Device**

(b) (4)

(b)(4) Performance Data




(b)(4)Performance Data



**Safety Data - Predicate Devices**

(b)(4)Performance Data




**Effectiveness Data – Subject Device**


(b)(4)Performance Data



(b)(4) Performance Data




(b)(4)Performance Data




**Effectiveness Data - Predicate Devices**

(b)(4)Performance Data



**Discussion of whether the data demonstrate that the subject device is as safe and effective as the predicate(s)**

(b)(4)Performance Data




**4. Does the product contain drugs or biologicals?**

- a. **If yes, what drug(s)/biologic(s): No**  
**Combination Product Code: N**

**5. Sterilization**


(b) (4)



**6. Is the Labeling Adequate?**

(OTC and/or Prescription) Rx  
Package Insert (page ) 17 and appendix C

(b) (4)



**7. Claims**

(b) (4)



**8. Has sponsor provided all administrative requirements?**

- Truthful and Accurate Statement 10
- 510(k) Summary or Statement 32
- Indication for Use Page 8

**9. Analysis of the Equivalence of the Subject and Predicates**

(b) (4)





THE 510(K) DOCUMENTATION FORMS ARE AVAILABLE ON THE LAN UNDER 510(K) BOILERPLATES TITLED "DOCUMENTATION" AND MUST BE FILLED OUT WITH EVERY FINAL DECISION (SE, NSE, NOT A DEVICE, ETC.).

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K032115

Reviewer: David Berkowitz  
 Division/Branch: DGRND/PRSR  
 Device Name: Neurolac Nerve Guide  
 Product To Which Compared (510(K) Number If Known): \_\_\_\_\_

	YES	NO	
1. Is Product A Device	<input checked="" type="checkbox"/>		If NO = Stop
2. Is Device Subject To 510(k)?	<input checked="" type="checkbox"/>		If NO = Stop
3. Same Indication Statement?	<input checked="" type="checkbox"/>		If YES = Go To 5
4. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NE
5. Same Technological Characteristics?		<input checked="" type="checkbox"/>	If YES = Go To 7
6. Could The New Characteristics Affect Safety Or Effectiveness?	<input checked="" type="checkbox"/>		If YES = Go To 8
7. Descriptive Characteristics Precise Enough?			If NO = Go To 10 If YES = Stop SE
8. New Types Of Safety Or Effectiveness Questions?		<input checked="" type="checkbox"/>	If YES = Stop NE
9. Accepted Scientific Methods Exist?	<input checked="" type="checkbox"/>		If NO = Stop NE
10. Performance Data Available?	<input checked="" type="checkbox"/>		If NO = Request Data
11. Data Demonstrate Equivalence?	<input checked="" type="checkbox"/>		Final Decision:

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

1. Intended Use: *Reconstruction of a peripheral nerve discontinuity up to 20 mm in patients with a complete nerve lesion.*
2. Device Description: Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Is the device life-supporting or life sustaining? Is the device implanted (short-term or long-term)? Does the device design use software? Is the device sterile? Is the device for single use? Is the device over-the-counter or prescription use? Does the device contain drug or biological product as a component? Is this device a kit? Provide a summary about the device's design, materials, physical properties and toxicology profile if important.

EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

1. Explain why not a device:
2. Explain why not subject to 510(k):
3. How does the new indication differ from the predicate device's indication:
4. Explain why there is or is not a new effect or safety or effectiveness issue: *The material is new for this type of device.*
5. Describe the new technological characteristics:  
*The device is made from poly(DL-lactid-ε-caprolactone).*
6. Explain how new characteristics could or could not affect safety or effectiveness: *May interfere with neural anastomosis.*
7. Explain how descriptive characteristics are not precise enough:
8. Explain new types of safety or effectiveness questions raised or why questions are not new:
9. Explain why existing scientific methods can not be used:
10. Explain what performance data is needed:  
*Device assist the rejoining of nerves.*
11. Explain how the performance data demonstrates that the device is or is not substantially equivalent:  
*Performs as well as predicates.*

ATTACH ADDITIONAL SUPPORTING INFORMATION

# Internal Administrative Form

	YES	I
1. Did the firm request expedited review?		X
2. Did we grant expedited review?		
3. Have you verified that the Document is labeled Class III for GMP purposes?		
4. If, not, has POS been notified?	X	
5. Is the product a device?	X	
6. Is the device exempt from 510(k) by regulation or policy?	X	
7. Is the device subject to review by CDRH?		X
8. Are you aware that this device has been the subject of a previous NSE decision?		
9. If yes, does this new 510(k) address the NSE issue(s), (e.g., performance data)?		
10. Are you aware of the submitter being the subject of an integrity investigation?		X
11. If, yes, consult the ODE Integrity Officer.		
12. Has the ODE Integrity Officer given permission to proceed with the review? (Blue Book Memo #191-2 and Federal Register 90N0332, September 10, 1991.		


**Berkowitz, David**

---

**From:** Hak [hak@polyganics.com]  
**Sent:** Wednesday, October 08, 2003 4:22 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application

Dear Mr. Berkowitz,

(b) (4)




Kind regards,

Jan-Bart Hak

----- Original Message -----

**From:** Berkowitz, David  
**To:** 'Hak'  
**Sent:** Tuesday, October 07, 2003 2:44 PM  
**Subject:** RE: 510(k) application

(b) (4)




-----Original Message-----


**From:** Hak [mailto:hak@polyganics.com]  
**Sent:** Friday, October 03, 2003 9:55 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application

Dear Mr. Berkowitz,

(b) (4)



(b) (4)



Kind regards,


Jan-Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
tel:+31 50 588 6586  
gsm: +31 653 211 303

The above information is intended only for the person or entity to whom it is addressed and may contain confidential and/or privileged information. Any review, retransmission, dissemination of, or taking action in reliance upon this information by others than the intended recipient is prohibited. If you are not the intended recipient, please return this e-mail to the sender and delete it from any computer system.

----- Original Message -----

**From:** Berkowitz, David  
**To:** 'Hak'  
**Sent:** Thursday, October 02, 2003 1:10 PM  
**Subject:** RE: 510(k) application

(b) (4)



-----Original Message-----

**From:** Hak [mailto:hak@polyganics.com]  
**Sent:** Thursday, October 02, 2003 2:46 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application

Dear Mr. Berkowitz,


I think e-mail is a perfect to address any issue quickly. I looking forward to your questions.

Kind regards,

Jan-Bart Hak, Ph.D.

Manager Clinical and Regulatory Affairs  
Polyganics BV  
tel: +31 50 588 6586  
Mobile: +31 653 211 303

(b) (4)



----- Original Message -----

**From:** Berkowitz, David  
**To:** 'hak@polyganics.com'  
**Sent:** Wednesday, October 01, 2003 7:22 PM  
**Subject:** 510(k) application

(b) (4)




**Berkowitz, David**

---

**From:** Hak [hak@polygenics.com]  
**Sent:** Wednesday, October 08, 2003 7:48 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application

Dear Mr. Berkowitz,

(b) (4)




Kind regards,

Jan-Bart Hak

----- Original Message -----

**From:** Berkowitz, David  
**To:** 'Hak'  
**Sent:** Tuesday, October 07, 2003 2:44 PM  
**Subject:** RE: 510(k) application

(b) (4)




-----Original Message-----

**From:** Hak [mailto:hak@polygenics.com]  
**Sent:** Friday, October 03, 2003 9:55 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application


Dear Mr. Berkowitz,

(b) (4)



24

(b) (4)



Kind regards,


Jan-Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
tel:+31 50 588 6586  
gsm: +31 653 211 303

The above information is intended only for the person or entity to whom it is addressed and may contain confidential and/or privileged information. Any review, retransmission, dissemination of, or taking action in reliance upon this information by others than the intended recipient is prohibited. If you are not the intended recipient, please return this e-mail to the sender and delete it from any computer system.

----- Original Message -----

**From:** Berkowitz, David  
**To:** 'Hak'  
**Sent:** Thursday, October 02, 2003 1:10 PM  
**Subject:** RE: 510(k) application

(b) (4)



-----Original Message-----

**From:** Hak [mailto:hak@polyganics.com]  
**Sent:** Thursday, October 02, 2003 2:46 AM  
**To:** Berkowitz, David  
**Subject:** Re: 510(k) application

Dear Mr. Berkowitz,

I think e-mail is a perfect to address any issue quickly. I looking forward to your questions.

Kind regards,

Jan-Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Polyganics BV  
tel: +31 50 588 6586  
Mobile: +31 653 211 303

The above information is intended only for the person or entity to whom it is addressed and may contain confidential and/or privileged information. Any review, retransmission, dissemination of, or taking action in reliance upon this information by others than the intended recipient is prohibited. If you are not the intended recipient, please return this e-mail to the sender and delete it from any computer system.

----- Original Message -----

**From:** Berkowitz, David  
**To:** 'hak@polyganics.com'  
**Sent:** Wednesday, October 01, 2003 7:22 PM  
**Subject:** 510(k) application

(b) (4)





---

**Indications for Use Form**

510(k) Number: \_\_\_\_\_

Device Name: **Neurolac® Nerve Guide**

**Indications for Use:**

The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use \_\_\_\_\_ OR Over-The-Counter Use \_\_\_\_\_  
(Per 21 CFR 801.109) (Optional Format 1-2-96)

\_\_\_\_\_  
(Division Sign-Off)

510(k) Number \_\_\_\_\_



---

**510(k)**  
**Summary of Safety and Effectiveness**

---

**Submitter:** Polyganics BV  
L.J. Zielstraweg 1  
9713 GX, Groningen  
The Netherlands  
[www.polyganics.com](http://www.polyganics.com)

**Contact Person:** Jan Bart Hak, Ph.D.  
Manager Clinical and Regulatory Affairs  
Tel : +31 50 588 6588  
Fax : +31 50 588 6599  
Mobile : +31 653 211 303  
E-mail : [hak@polyganics.com](mailto:hak@polyganics.com)

**Date Prepared:** May 20, 2003

---

**General Provisions:** Trade Name: Neurolac® Nerve guide  
Common Name: Nerve guide  
Classification Name: Nerve Cuff, 21 CFR 882.5275  
Device Classification: Class II

---

**Predicate Devices:**

- Neurotube™ Neuroregen L.L.C. K983007
- NeuroGen™ Integra Life Sciences Corp. K011168

---

**Performance Standards** For the Nerve Cuff performance, the FDA, under section 514 of the Food, Drug and Cosmetic Act, has not established standards.

---

**Indications for Use** The Neurolac® nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

---



---

**Device Description**

Neurolac® is designed to be a flexible and transparent resorbable poly(DL-lactide-co-ε-caprolactone) tube to provide a protective environment for peripheral nerve regeneration after injury and to create a conduit to guide axonal growth across a nerve gap.

Neurolac® nerve guides are provided sterile in Tyvek pouch packages and retainer in a variety of sizes.

---

**Performance Data:**

The safety and effectiveness of the Neurolac nerve guides have been demonstrated via data collected from design verification tests and analyses. The design verification testing consisted of the following:

- In vitro suture retention testing
- In vitro degradation testing
- In vivo nerve function recovery

---

**Summary of Substantial Equivalence**

The design, fundamental technology and intended use (safety and efficacy) featured with the Neurolac® Nerve Guide are substantially equivalent to those featured with the competitor devices Neurotube™ (ref. 510(k) 983007; Neuroregen L.L.C.) and the NeuroGen™ Nerve Guide (ref. 510(k) 011168; Integra Life Sciences Corporation).

Biocompatibility, mechanical and physical property testing, in vitro degradation testing, and performance testing in an animal model provide reasonable scientific evidence that Neurolac® nerve guide is substantially equivalent to the predicate devices. Evaluation of the Polyganics Neurolac® Nerve guide based on biocompatibility testing, animal tests, results from literature and the comparison of the Neurolac® nerve guide with its predicate devices, shows that the Neurolac® nerve guide is safe for implantation.



---

### Section 3: PROPOSED LABELING

**Subject Device Labeling**      The following proposed labeling for the subject device Neurolac Nerve Guide is provided in APPENDIX C:

- Outer label (Carton)
- Inner label (Pouch)
- Pre-printed carton text and graphics
- Instructions for Use

---

**Intended Use of the Subject Device**      The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

---

**Promotional Materials**      At present, no promotional materials are available for the subject device.

---

## Polyganics' Neurolac® bioresorbable nerve guide

Instructions for Use, English

**STERILE.** Sterilized with ethylene oxide gas. For single use only. Do not autoclave.  
**CAUTION:** Federal (USA) law restricts this device to sale by or on the order of a physician.

### Neurolac® peripheral nerve guide

#### Description

The Neurolac nerve guide is composed of the bioresorbable copolyester poly(DL-lactide-ε-caprolactone). The Neurolac nerve guide provides guidance and protection to regenerating axons. The Neurolac nerve guide elicits a minimal acute inflammatory reaction of the surrounding tissue, which is followed by gradual encapsulation of the tube by fibrous tissue. Degradation of the Neurolac nerve guide occurs through hydrolysis leading to gradual reduction of molecular weight. The Neurolac nerve guide retains its initial mechanical properties up to 8 weeks, thereafter rapid loss of mechanical strength and gradual mass loss occur. The final degradation products, lactic acid and ω-hydroxy hexanoic acid, are resorbed, metabolized and excreted by the body. Animal studies demonstrated that a Neurolac nerve guide is resorbed within 16 months.

The Neurolac nerve guide inner diameter is indicated on the label, and is packed in a tray placed in a Tyvek pouch. The Neurolac nerve guide is indicated for single-use.

#### Indications

The Neurolac nerve guide is indicated for the reconstruction of a peripheral nerve discontinuity up to 20 mm in patients who have sustained a complete division of a nerve.

#### Contraindications

There are no known contraindications.

#### Warnings

- The Neurolac nerve guide is for single use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired through cleaning, resterilization, or re-use and may cause adverse patient reactions. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from re-use of (or any part of) the Neurolac nerve guide.
  - Sterile unless package has been opened or damaged.
  - Discard open unused nerve guides.
  - The Neurolac nerve guide should only be used by physicians who are trained in nerve defect repair techniques. Accordingly, Polyganics BV will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel. The physician should consult recent literature on current medical practice on peripheral nerve repair.
  - Nerve regeneration may be suboptimal in elderly, malnourished or debilitated patients or in patients suffering from cancer, anaemia, obesity, diabetes, infection or other conditions which may delay wound healing, injected wounds, or moderate to severe inflammatory response characteristic of foreign body response.
- Precautions**
- Use prior to "Use by date".
  - Store in dark, dry place at or below 4°C (39°F).
  - Do not expose the nerve guide to organic solvents (e.g. chloroform, acetone).
  - Do not use absorbable sutures for fixation of the nerve stumps into the nerve guide.



- Avoid crushing, crimping, kinking or other damage due to application of surgical instruments such as forceps, needle, holders and scissors or during handling of the device.
- Avoid tension on the nerve ends.
- Prevent compression and/or kinking of the Neurolac after the procedure. The use of a protective splint is recommended.

#### Adverse effects

- Adverse events associated with the use of a Neurolac nerve guide may include but are not limited to:
- Failure to provide adequate nerve regeneration at sites where too much tension or compression occurs;
  - Failure to provide adequate/complete nerve regeneration;
  - Transitory local irritation;
  - Infection;
  - Allergy;
  - Delayed wound healing.

#### Opening of the package

The pouch is opened in such a way that the tray remains sterile. The tray can be opened by sliding the lid. By clasping the nerve guide at one of its ends between a pair of tweezers, it can be taken from the tray. The lid contains a ruler that may be used as a reference to estimate the gap length or nerve stump diameter.

#### Surgical Procedure

- Surgically expose the injured nerve.
  - Resect the injured segment distally and proximally until a nerve stump is identified with no residual intrafascicular scarring.
- NOTE: Do not crush the nerve stumps as this can cause extrusion of intra-fascicular components.**
- Measure the length of the defect with all joints in an extended position.
  - If the gap length is between 0 and 20 mm, the injured nerve can be reconstructed with a Neurolac nerve guide.
  - Select the Neurolac nerve guide with the proper internal diameter.

**NOTE: It is essential that the internal nerve guide diameter is slightly larger than the diameter of the transected nerve to guarantee optimal nerve regeneration.**

- Cut the selected nerve guide with a pair of scissors or a knife so that the nerve guide is 1 cm longer than the nerve gap.

Under some circumstances immobilization of the nerve ends, as to avoid tension on the nerve ends, may be employed at the discretion of the surgeon. To secure adequate fixation of the nerve ends in the nerve guide, the accepted surgical technique of flat square ties with additional throws as warranted by surgical circumstance and the experience of the surgeon, is required.

#### Suturing Technique

- Place the Neurolac nerve guide in warm saline (37°C) for approximately 1 minute before implantation. This will make the tube more flexible and ease needle passage during suturing.
- Suture the Neurolac nerve guide by passing the suture (8-0 suture) first through the tube from the outside to the inside and then transversally and superficially through the epineurium and back through the tube from the inside to the outside, after which a tie is made (Fig. 1.1-1.3).
- When positioning optimization of the nerve ends in the nerve guide is required, it is recommended to place a second suture in the same nerve end (Fig. 1.4).
- Pull the proximal nerve stump into the nerve guide.

## Polyganics' Neurolac® bioresorbable nerve guide

- NOTE: It is recommended that the nerve ends are pulled into the tube for at least 3 mm for optimal nerve regeneration.**
- Fill the tube with heparinized saline, using a solution containing 1000 units of heparin per 100 ml of normal saline (Fig 1.5).
  - Subsequently, use the same procedure, to pull the distal nerve stump into the nerve guide.
  - A minimum space of 5 mm should be left between the nerve ends in the nerve guide.
  - Fill any remaining space with heparinized saline (Fig 1.6) by injecting along the nerve into the lumen of the tube or by penetrating the tube (not the nerve).




Figure 1. Schematic representation of suture technique for suturing the nerve ends into the nerve guide.

**CAUTION:** Ensure that no blood enters the nerve guide lumen since this may hinder nerve recovery.  
**CAUTION:** The nerve guide should be implanted and sutured with all joints in an extended position as to assure that no tension occurs on the proximal or distal nerve end when joints are being mobilized.

- Close the wound and splint to prevent kinking for the first three postoperative weeks. Long-term compression of the nerve guide should be avoided. Patients may be administered oral antibiotics for the first post-operative week.

Dispose contaminated implantation and packaging materials utilizing standard hospital procedures and universal precautions for bio-hazardous waste.



- Refer to accompanying instructions for Use.
- Use by
- Catalog number
- For single use only
- Lot number
- Sterile product
- Sterilized by ethylene oxide

Polyganics BV  
 C.J. Ziebellweg 1  
 9713 GK Groningen  
 The Netherlands  
 Telephone +31 50 568 6588  
 Fax +31 50 568 6599

**DISCLAIMER OF WARRANTY AND LIMITATION OF REMEDY**  
 There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Polyganics product(s) described in this publication. Under no circumstances shall Polyganics BV be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Polyganics to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Polyganics BV printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Polyganics BV will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

Neurolac® is a registered trademark of Polyganics