



## Dental Morelli Ltda.

### SECTION 6

#### 510(k) SUMMARY

<b>Proprietary Name</b>	Edgewise Ceramic Brackets; Roth Ceramic Brackets
<b>Date Prepared</b>	April 15, 2013
<b>Submitter</b>	DENTAL MORELLI LTDA Alameda Jundiaí, 230 – Jardim Saira - Sorocaba CEP: 18085-090 Brazil Telephone: 55 (15) 3238-8200
<b>Official Contact</b>	Tara Conrad TechLink International Consulting 18851 NE 29 <sup>th</sup> Avenue Suite 720 Aventura, FL 33180 TEL- (305) 377-0077
<b>Common Name</b>	Orthodontic Ceramic Brackets
<b>Trade Name</b>	Edgewise Ceramic Brackets; Roth Ceramic Brackets
<b>Classification</b>	Class II
<b>Product Code</b>	NJM
<b>Classification Panel</b>	Dental
<b>Regulation Numbers</b>	21 CFR 872.5470
<b>Substantial Equivalence</b>	K102803 Clarity Advanced Ceramic Brackets

AUG 27 2013

#### Description of Proposed Device

Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The bracket is uncoated. The Morelli bracket is not self-bonding and requires a primer and adhesive to bond the bracket to the teeth.



## Dental Morelli Ltda.

### Indications for Use

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

### Device Comparison Table

	Edgewise Ceramic Brackets by Dental Morelli	Roth Ceramic Brackets by Dental Morelli	Clarity Advance Ceramic Brackets by 3M Unitek Corporation K102803
Indications for use	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Clarity Advanced Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.
Target Population	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction
Anatomical Site	Teeth	Teeth	Teeth
Location of use	Use only by professional orthodontists	Use only by professional orthodontists	Use only by professional orthodontists
Materials	Aluminum oxide	Aluminum oxide	Aluminum oxide
Biocompatibility	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Compatibility with the environment and other devices	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Sterility	Non-sterile	Non-sterile	Non-sterile
Maxillary In-out (mm)	0.94	0.6-1.2	0.53-.089
Maxillary Torque	0	-7 to +8	-7 to +17
Maxillary Angulation	0	0 to +12	0 to +8
Mandibular In-out (mm)	0.94	0.6 – 1.2	0.51-1.14
Mandibular Torque	0	-22 to 0	-17 to 0
Slot	0.022"	0.022"	0.022"



Dental Morelli Ltda.

---

### **Substantial Equivalence**

Both the non-clinical data and the biocompatibility evaluation indicate that Edgewise and Roth Ceramic Brackets are safe and effective for their intended use in orthodontic treatment and perform as well as the predicate. The subject and predicate devices have the same intended use, indications for use, compositions, device design and performance.

### **Device Material and Design**

The body of the subject and predicate devices are composed of ceramic. The Edgewise and Roth Ceramic Bracket is not coated and does not have a liner. The Morelli Brackets have rounded edges and corners. The Edgewise and Roth Ceramic Bracket is not built to facilitate debonding.

### **Conclusion**

This premarket notification is being submitted to request clearance for the Edgewise and Roth Ceramic Brackets. The analysis on the Edgewise Ceramic Brackets demonstrates substantial equivalence to the 3M Unitek Corporation predicate.



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

August 27, 2013

Dental Morelli Limited  
C/O Ms. Tara Conrad  
Regulatory Affairs Manager  
Techlink International Consultants  
18851 NE 29<sup>th</sup> Avenue Suite 720  
AVENTURA FL 33180

Re: K131197  
Trade/Device Name: Edgewise Ceramic Brackets; Roth Ceramic Brackets  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: NJM  
Dated: May 15, 2013  
Received: June 4, 2013

Dear Ms Conrad:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**Mary S. Runner -S**

Kwame Ulmer M.S.  
Acting Division Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K131197

**Indications for Use Statement**

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

Prescription Use  (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use \_\_\_\_\_ (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEED

Sheena A. Green-5  
2013.08.27 14:29:58 -0400

for M. Susan Runner, DDS, MA

(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K131197

Concurrence of CDRH, Office of Device Evaluation  
(ODE) Page 1 of 1



Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center – WO66-G609  
Silver Spring, MD 20993-0002

August 27, 2013

Dental Morelli Limited  
C/O Ms. Tara Conrad  
Regulatory Affairs Manager  
Techlink International Consultants  
18851 NE 29<sup>th</sup> Avenue Suite 720  
AVENTURA FL 33180

Re: K131197  
Trade/Device Name: Edgewise Ceramic Brackets; Roth Ceramic Brackets  
Regulation Number: 21 CFR 872.5470  
Regulation Name: Orthodontic Plastic Bracket  
Regulatory Class: II  
Product Code: NJM  
Dated: May 15, 2013  
Received: June 4, 2013

Dear Ms Conrad:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you, however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

Mary S. Runner -S

Kwame Ulmer M.S.  
Acting Division Director  
Division of Anesthesiology, General Hospital,  
Respiratory, Infection Control and  
Dental Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Concurrence & Template History Page**

[THIS PAGE IS INCLUDED IN IMAGE COPY ONLY]

Full Submission Number: K131197

For Office of Compliance Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=318](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=318)

For Office of Surveillance and Biometrics Contact Information:

[http://insideportlets.fda.gov:9010/portal/page?\\_pageid=197,415881&\\_dad=portal&\\_schema=PORTAL&org=423](http://insideportlets.fda.gov:9010/portal/page?_pageid=197,415881&_dad=portal&_schema=PORTAL&org=423)

<b>Digital Signature Concurrence Table</b>	
Reviewer Sign-Off	Myra Browne
Branch Chief Sign-Off	Sheena A. Green for M. Susan Runner, DDS, MA
Division Sign-Off	 <p>Mary S. Runner -S 2013.08.27 14:52:59 -04'00'</p>

Template Name: K1(A) – SE after 1996

Template History:

Date of Update	By	Description of Update
7/27/09	Brandi Stuart	Added Updates to Boiler Table
8/7/09	Brandi Stuart	Updated HFZ Table
1/11/10	Diane Garcia	Liability/Warranty sentence added at bottom of 1 <sup>st</sup> page
10/4/11	M. McCabe Janicki	Removed IFU sheet and placed in Forms
9/25/12	Edwena Jones	Added digital signature format
12/12/12	M. McCabe Janicki	Added an extra line between letter signature block and the word "Enclosure". Also, added a missing digit in 4-digit extension on letterhead zip code: "002" should be "0002".
4/2/2013	M. McCabe Janicki	Edited sentence that starts "If you desire specific advice for your device on our labeling regulation (21 CFR Part 801)..." Replaced broken Compliance link with general link to DSMICA.
4/12/2013	Margaret McCabe Janicki	Fixed a typo: Paragraph 1, final sentence, "We remind you, however; that device labeling must be truthful..." Replaced incorrect semicolon with a comma.

K131197

**Indications for Use Statement**

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

Prescription Use  (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use  (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEED

Sheena A. Green - 5  
2013.08.27 14:29:58 - 0f00

for M. Susan Runner, DDS, MA

(Division Sign-Off)

Division of Anesthesiology, General Hospital  
Infection Control, Dental Devices

510(k) Number: K131197

Concurrence of CDRH, Office of Device Evaluation  
(ODE) Page 1 of 1



Dental Morelli Ltda.

K/31/97  
FDA CDRH DMC

APR 26 2013

Received

---

**CONFIDENTIAL**

**COVER LETTER**

April 15, 2013

Document Mail Center W066-G609  
Center for Devices and Radiological Health Food and Drug  
Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002 - USA

RE: 510(k) Traditional Pre Market Notification Request

**COMPANY NAME AND ADDRESS**

DENTAL MORELLI LTDA  
Alameda Jundiá, 230 -  
Jardim Saira - Sorocaba  
CEP: 18085-090  
Brazil

Telephone: 55 (15) 3238-8200

**CONSULTANT NAME AND ADDRESS**

TechLink International Consultants  
18851 NE 29<sup>th</sup> Avenue  
Suite 720  
Aventura, Florida 33180

Telephone: (305) 377-0077

Primary Contact: Tara Conrad  
Secondary Contact: Lilian Llull

**Attention: Document Control Clerk**

**Device Trade Name:**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

According to Section 510 (k) of the Federal Food, Drug and Cosmetic Act, as amended (ACT), Dental Morelli Ltda proposes to introduce the Orthocontic Ceramic

---



## Dental Morelli Ltda.

Brackets into interstate commerce for commercial distribution and hereby requests 510 (k) clearance by the Food and Drug Administration, as required by law.

The following information on these products (according to 21 CFR 807.87) is submitted for your consideration.

**Common Name:** Orthodontic Ceramic Brackets  
**Trade Name:** Edgewise Ceramic Brackets; Roth Ceramic Brackets  
**Classification:** Class II  
**Product Code:** NJM  
**Classification Panel:** Dental  
**Regulation Numbers:** 21 CFR 872.5470  
**Substantial Equivalence:** K102803 Clarity Advanced Ceramic Brackets

The eCopy is an exact duplicate of the paper copy.

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain components derived from a tissue or other biologic source?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device provided sterile?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device intended for single use?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device a reprocessed single use device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, does this device type require reprocessed validation data?	N/A	N/A
Does the device contain a drug?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain a biologic?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device use software?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the submission include clinical information?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device implanted?	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## Dental Morelli Ltda.

---

The predicate device and proposed device have the same intended use, same indications, and are designed with the same technological characteristics. A complete list of indications and a comparison table are included in Section 13. Indications for Use Statement can be found in Section 5. Labeling specifications are detailed in section 14. Considering our intent to market these devices, all information submitted is confidential; with the exception of the 510(k) summary.

Dental Morelli has not submitted prior applications for the same device. However, there are other manufacturers who have submitted 510(k) applications for similar devices. These similar devices are approved by the FDA. In particular the Clarity Advanced Ceramic Brackets that has been used as a predicate.

We are confident that this information will be sufficient for you to reach a favorable decision. However, please do not hesitate to contact me if you have any additional questions, concerns, or if you feel that I can be of any further assistance.

Best regards,

A handwritten signature in blue ink, appearing to read 'Tara Conrad'.

Tara Conrad, Biomedical Engineer  
Regulatory Affairs Manager

(305) 377-0077 – ph  
taraconrad@techlinkusa.net

---



**Confidential**

## **Table of Contents**

Section	Page	
1	Table of Contents	1
2	Form 3514 CDRH Premarket Review Submission Cover Sheet	2
3	Form 3601 Medical Device User Fee Cover Sheet	7
4	Cover Letter	8
5	Indications for Use	11
6	510(K) Summary	12
7	Truthful and Accurate Statement	15
8	Class III Summary	16
9	Financial Disclosure	17
10	Declaration of Conformity	18
11	Executive Summary	23
12	Device Description	27
13	Substantial Equivalence Discussion	38
14	Proposed Labeling	47
15	Sterilization	55
16	Biocompatibility	56
17	Software Validation	57
18	Electromagnetic Compatibility	58
19	Bench Test	59
20	Animal Tests	82
21	Clinical Tests	83

**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**

Date of Submission 4/18/2013	User Fee Payment ID Number <b>(b)(4)Trade Secret</b>	FDA Submission Document Number (if known)
---------------------------------	---	---

**SECTION A TYPE OF SUBMISSION**

<b>PMA</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<b>PMA &amp; HDE Supplement</b> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<b>PDP</b> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<b>510(k)</b> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<b>Meeting</b> <input type="checkbox"/> Pre-510(K) Meeting <input type="checkbox"/> Pre-IDE Meeting <input type="checkbox"/> Pre-PMA Meeting <input type="checkbox"/> Pre-PDP Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Other (specify):
<b>IDE</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<b>Humanitarian Device Exemption (HDE)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<b>Class II Exemption Petition</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Evaluation of Automatic Class III Designation (De Novo)</b> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<b>Other Submission</b> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Dental Morelli	Establishment Registration Number (if known) Pending		
Division Name (if applicable)	Phone Number (including area code) 55 (15) 3238-8200		
Street Address Alameda Jundiai, 230/250	FAX Number (including area code)		
City Sorocaba	State / Province Sao Paulo	ZIP/Postal Code 18085-090	Country Brazil

**(b)(4)Trade Secret Process**

**SECTION C APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name TechLink International Consulting	Phone Number (including area code) 305-377-0077		
Division Name (if applicable)	FAX Number (including area code)		
Street Address 18851 NE 29th Ave 720	FAX Number (including area code)		
City Aventura	State / Province Florida	ZIP Code 33180	Country USA
Contact Name Tara Conrad			
Contact Title Regulatory Affaris Manager	Contact E-mail Address taraconrad@techlinkusa.net		

**SECTION D1**

**REASON FOR APPLICATION - PMA, PDP, OR HDE**

- New Device
- Withdrawal
- Additional or Expanded Indications
- Request for Extension
- Post-approval Study Protocol
- Request for Applicant Hold
- Request for Removal of Applicant Hold
- Request to Remove or Add Manufacturing Site

- Change in design, component, or specification:
  - Software / Hardware
  - Color Additive
  - Material
  - Specifications
  - Other (*specify below*)

- Location change:
  - Manufacturer
  - Sterilizer
  - Packager

- Process change:
  - Manufacturing     Packaging
  - Sterilization
  - Other (*specify below*)

- Labeling change:
  - Indications
  - Instructions
  - Performance Characteristics
  - Shelf Life
  - Trade Name
  - Other (*specify below*)

- Report Submission:
  - Annual or Periodic
  - Post-approval Study
  - Adverse Reaction
  - Device Defect
  - Amendment

- Response to FDA correspondence:

- Change in Ownership
- Change in Correspondent
- Change of Applicant Address

- Other Reason (*specify*):

**SECTION D2**

**REASON FOR APPLICATION - IDE**

- New Device
- New Indication
- Addition of Institution
- Expansion / Extension of Study
- IRB Certification
- Termination of Study
- Withdrawal of Application
- Unanticipated Adverse Effect
- Notification of Emergency Use
- Compassionate Use Request
- Treatment IDE
- Continued Access

- Change in:
  - Correspondent / Applicant
  - Design / Device
  - Informed Consent
  - Manufacturer
  - Manufacturing Process
  - Protocol - Feasibility
  - Protocol - Other
  - Sponsor

- Response to FDA Letter Concerning:
  - Conditional Approval
  - Deemed Approved
  - Deficient Final Report
  - Deficient Progress Report
  - Deficient Investigator Report
  - Disapproval
  - Request Extension of Time to Respond to FDA
  - Request Meeting
  - Request Hearing

- Report submission:
  - Current Investigator
  - Annual Progress Report
  - Site Waiver Report
  - Final

- Other Reason (*specify*):

**SECTION D3**

**REASON FOR SUBMISSION - 510(k)**

- New Device

- Additional or Expanded Indications

- Change in Technology

- Other Reason (*specify*):

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed

1	NJM	2		3		4	
5		6		7		8	

Summary of, or statement concerning, safety and effectiveness information

- 510 (k) summary attached  
 510 (k) statement

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1	K102803	Clarity Advanced Ceramic Brackets	3M Unitek Corporation
2			
3			
4			
5			
6			

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name

Orthodontic Ceramic Brackets

	Trade or Proprietary or Model Name for This Device	Model Number
1	Edgewise Ceramic Brackets	1
2	Roth Ceramic Brackets	2
3		3
4		4
5		5

FDA document numbers of all prior related submissions (regardless of outcome)

1	2	3	4	5	6
7	8	9	10	11	12

Data Included in Submission

- Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code NJM	C.F.R. Section (if applicable) 872.5470	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel Dental		

Indications (from labeling)

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

**Note:** Submission of the information entered in Section H does not affect the need to submit device establishment registration.

FDA Document Number (if known)

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler		
Company / Institution Name			Establishment Registration Number		
Division Name (if applicable)			Phone Number (including area code)		
Street Address			FAX Number (including area code)		
City		State / Province	ZIP Code	Country	
Contact Name		Contact Title		Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler		
Company / Institution Name			Establishment Registration Number		
Division Name (if applicable)			Phone Number (including area code)		
Street Address			FAX Number (including area code)		
City		State / Province	ZIP Code	Country	
Contact Name		Contact Title		Contact E-mail Address	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Contract Manufacturer <input type="checkbox"/> Repackager / Relabeler		
Company / Institution Name			Establishment Registration Number		
Division Name (if applicable)			Phone Number (including area code)		
Street Address			FAX Number (including area code)		
City		State / Province	ZIP Code	Country	
Contact Name		Contact Title		Contact E-mail Address	

## SECTION I

## UTILIZATION OF STANDARDS

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	N/A	FDA	Format for Traditional and Abbreviated 510 (k)'s	N/A	8/12/2005
2	ISO 10993-1	International Standards Organization	Biological evaluation of medical Part 1: Evaluation and testing		01/01/2009
3	ISO 14971	International Standards Organization	Medical Devices- application of risk management to medical		01/01/2007
4	ISO 15223-1	International Standards Organization	Medical devices- Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General Requirements		01/01/2012
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

**Public reporting burden for this collection of information** is estimated to average 0.5 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
1350 Piccard Drive, Room 400  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

Form Approved OMB No. 0910-511 Expiration Date February 28, 2013. See Instructions for OMB Statement

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>		PAYMENT IDENTIFICATION NUMBER: (b)(4)Trade Write the Payment Identification number on your check.
A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <a href="http://www.fda.gov/oc/mdufma/cover sheet.html">http://www.fda.gov/oc/mdufma/cover sheet.html</a>		
1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  DENTAL MORELLI LTDA. ALAMEDA JUNDIA,230/250 JD.SAIRA SOROCABA / SP 18085-090 BR  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN)	2. CONTACT NAME Roger Morelli  2.1 E-MAIL ADDRESS assuntosregulatorios@morelli.com.br  2.2 TELEPHONE NUMBER (include Area code) 55-15-32388200  2.3 FACSIMILE (FAX) NUMBER (Include Area code)	
3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/oc/mdufma">http://www.fda.gov/oc/mdufma</a> ) <u>Select an application type:</u> <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> Annual Fee for Periodic Reporting (APR) <input type="checkbox"/> 30-Day Notice		
3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER <u>3.2 Select one of the types below</u> <input checked="" type="checkbox"/> Original Application <u>Supplement Types:</u> <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)		
4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input checked="" type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input type="checkbox"/> NO, I am not a small business 4.1 If Yes, please enter your Small Business Decision Number: SBD135206		
5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/cdrh/mdufma">http://www.fda.gov/cdrh/mdufma</a> for additional information)		
6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION. <input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only <input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially		
7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO		
PAPERWORK REDUCTION ACT STATEMENT Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to the address below.  Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive, 4th Floor Rockville, MD 20850 [Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]		
8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION (b)		25-Feb-2013

Form FDA 3601 (01/2007)

["Close Window"](#) [Print Cover sheet](#)



Dental Morelli Ltda.

---

**CONFIDENTIAL**

**COVER LETTER**

April 15, 2013

Document Mail Center W066-G609  
Center for Devices and Radiological Health Food and Drug  
Administration  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002 – USA

RE: 510(k) Traditional Pre Market Notification Request

**COMPANY NAME AND ADDRESS**

DENTAL MORELLI LTDA  
Alameda Jundiaí, 230 –  
Jardim Saira - Sorocaba  
CEP: 18085-090  
Brazil

Telephone: 55 (15) 3238-8200

**CONSULTANT NAME AND ADDRESS**

TechLink International Consultants  
18851 NE 29<sup>th</sup> Avenue  
Suite 720  
Aventura, Florida 33180

Telephone: (305) 377-0077

Primary Contact: Tara Conrad  
Secondary Contact: Lilian Llull

**Attention: Document Control Clerk**

**Device Trade Name:**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

According to Section 510 (k) of the Federal Food, Drug and Cosmetic Act, as amended (ACT), Dental Morelli Ltda proposes to introduce the Orthodontic Ceramic



## Dental Morelli Ltda.

Brackets into interstate commerce for commercial distribution and hereby requests 510 (k) clearance by the Food and Drug Administration, as required by law.

The following information on these products (according to 21 CFR 807.87) is submitted for your consideration.

**Common Name:** Orthodontic Ceramic Brackets  
**Trade Name:** Edgewise Ceramic Brackets; Roth Ceramic Brackets  
**Classification:** Class II  
**Product Code:** NJM  
**Classification Panel:** Dental  
**Regulation Numbers:** 21 CFR 872.5470  
**Substantial Equivalence:** K102803 Clarity Advanced Ceramic Brackets

The eCopy is an exact duplicate of the paper copy.

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain components derived from a tissue or other biologic source?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device provided sterile?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device intended for single use?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device a reprocessed single use device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, does this device type require reprocessed validation data?	N/A	N/A
Does the device contain a drug?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain a biologic?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device use software?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the submission include clinical information?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device implanted?	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## Dental Morelli Ltda.

---

The predicate device and proposed device have the same intended use, same indications, and are designed with the same technological characteristics. A complete list of indications and a comparison table are included in Section 13. Indications for Use Statement can be found in Section 5. Labeling specifications are detailed in section 14. Considering our intent to market these devices, all information submitted is confidential; with the exception of the 510(k) summary.

Dental Morelli has not submitted prior applications for the same device. However, there are other manufacturers who have submitted 510(k) applications for similar devices. These similar devices are approved by the FDA. In particular the Clarity Advanced Ceramic Brackets that has been used as a predicate.

We are confident that this information will be sufficient for you to reach a favorable decision. However, please do not hesitate to contact me if you have any additional questions, concerns, or if you feel that I can be of any further assistance.

Best regards,

Tara Conrad, Biomedical Engineer  
Regulatory Affairs Manager

(305) 377-0077 - ph  
taraconrad@techlinkusa.net

## Indications for Use Statement

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

Perscription Use  (Part 21 CFR 801 Subpart D) AND/OR

Over-The-Counter Use  (21 CFR 801 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEED

---

Concurrence of CDRH, Office of Device Evaluation  
(ODE) Page 1 of 1



## Dental Morelli Ltda.

---

### SECTION 6

#### 510(k) SUMMARY

<b>Proprietary Name</b>	Edgewise Ceramic Brackets; Roth Ceramic Brackets
<b>Date Prepared</b>	April 15, 2013
<b>Submitter</b>	DENTAL MORELLI LTDA Alameda Jundiaí, 230 – Jardim Saira - Sorocaba CEP: 18085-090 Brazil Telephone: 55 (15) 3238-8200
<b>Official Contact</b>	Tara Conrad TechLink International Consulting 18851 NE 29 <sup>th</sup> Avenue Suite 720 Aventura, FL 33180 TEL- (305) 377-0077
<b>Common Name</b>	Orthodontic Ceramic Brackets
<b>Trade Name</b>	Edgewise Ceramic Brackets; Roth Ceramic Brackets
<b>Classification</b>	Class II
<b>Product Code</b>	NJM
<b>Classification Panel</b>	Dental
<b>Regulation Numbers</b>	21 CFR 872.5470
<b>Substantial Equivalence</b>	K102803 Clarity Advanced Ceramic Brackets

#### Description of Proposed Device

Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The bracket is uncoated. The Morelli bracket is not self-bonding and requires a primer and adhesive to bond the bracket to the teeth.



## Dental Morelli Ltda.

### Indications for Use

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

### Device Comparison Table

	Edgewise Ceramic Brackets by Dental Morelli	Roth Ceramic Brackets by Dental Morelli	Clarity Advance Ceramic Brackets by 3M Unitek Corporation K102803
Indications for use	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Clarity Advanced Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.
Target Population	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction
Anatomical Site	Teeth	Teeth	Teeth
Location of use	Use only by professional orthodontists	Use only by professional orthodontists	Use only by professional orthodontists
Materials	Aluminum oxide	Aluminum oxide	Aluminum oxide
Biocompatibility	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Compatibility with the environment and other devices	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Sterility	Non-sterile	Non-sterile	Non-sterile
Maxillary In-out (mm)	0.94	0.6-1.2	0.53-.089
Maxillary Torque	0	-7 to +8	-7 to +17
Maxillary Angulation	0	0 to +12	0 to +8
Mandibular In-out (mm)	0.94	0.6 - 1.2	0.51-1.14
Mandibular Torque	0	-22 to 0	-17 to 0
Slot	0.022"	0.022"	0.022"



---

### **Substantial Equivalence**

Both the non-clinical data and the biocompatibility evaluation indicate that Edgewise and Roth Ceramic Brackets are safe and effective for their intended use in orthodontic treatment and perform as well as the predicate. The subject and predicate devices have the same intended use, indications for use, compositions, device design and performance.

### **Device Material and Design**

The body of the subject and predicate devices are composed of ceramic. The Edgewise and Roth Ceramic Bracket is not coated and does not have a liner. The Morelli Brackets have rounded edges and corners. The Edgewise and Roth Ceramic Bracket is not built to facilitate debonding.

### **Conclusion**

This premarket notification is being submitted to request clearance for the Edgewise and Roth Ceramic Brackets. The analysis on the Edgewise Ceramic Brackets demonstrates substantial equivalence to the 3M Unitek Corporation predicate.



Dental Morelli Ltda.

---

**CONFIDENTIAL**

**SECTION 7**

**Premarket Notification Truthful and Accurate Statement**

[As Required by 21 CFR 807.87(k)]

In my capacity as Director of Dental Morelli Ltda., I certify that to the best of my knowledge all data and information submitted in the premarket notification for the Edgewise and Roth Ceramic Brackets is truthful and accurate and that no material fact has been omitted.

*Roger Morelli*

---

Name: Roger Morelli

Position: Director

Date: 3/20/13

(Premarket Notification 510(k) number pending)



**Dental Morelli Ltda.**

---

**CONFIDENTIAL**

**SECTION 8**

**CLASS III SUMMARY AND CERTIFICATION**

Edgewise and Roth Ceramic Brackets are class II devices.

This section does not apply.



**Dental Morelli Ltda.**

---

**CONFIDENTIAL**

**SECTION 9**

**FINANCIAL DISCLOSURE**

Edgewise and Roth Ceramic Brackets are class II devices.

This submission is a traditional 510(k). A Clinical Trial was not conducted. The requirement for financial certification or disclosure as described in 21 CFR 807.87(i) does not apply to this submission.

**Confidential**

**Section 10**

**Declaration of Conformity**

Edgewise Ceramic Brackets

Roth Ceramic Brackets

Compliance with Performance Standards:

- ISO 10993-1:2009 - Biological evaluation of medical Part 1: Evaluation and testing
- ISO 14971:2007 - Medical devices - Application of risk management to medical
- ISO 15223-1:2012 - Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied Part 1: General requirements

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 15223-1: Medical devices- Symbols to be used with medical device labels, labelling, and information to be supplied -Part 1: Ge

**Please answer the following questions**

Yes    No

Is this standard recognized by FDA <sup>2</sup>? .....    

FDA Recognition number <sup>3</sup> ..... #5-73

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....    

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....       
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....    

Does this standard include acceptance criteria? .....       
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....       
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard?.....       
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....    

Were deviations or adaptations made beyond what is specified in the FDA SIS?.....       
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....       
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard?.....       
If yes, was the guidance document followed in preparation of this 510k? .....    

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 14971 Medical Devices -Application of risk management to medical devices

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....      

FDA Recognition number <sup>3</sup> ..... # 5-40

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....      

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....         
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....      

Does this standard include acceptance criteria? .....         
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....         
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard?.....         
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....      

Were deviations or adaptations made beyond what is specified in the FDA SIS?.....         
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....         
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard?.....         
If yes, was the guidance document followed in preparation of this 510k? .....      

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 10993-1:2009 Biological evaluation of medical devices --Part 1: Evaluation and testing within a risk management process

**Please answer the following questions**

Yes    No

Is this standard recognized by FDA <sup>2</sup>? .....    

FDA Recognition number <sup>3</sup> ..... #2-156

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....    

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....       
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....    

Does this standard include acceptance criteria? .....       
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....       
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard?.....       
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....    

Were deviations or adaptations made beyond what is specified in the FDA SIS?.....       
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....       
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard?.....       
If yes, was the guidance document followed in preparation of this 510k? .....    

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

BS EN ISO 27020: 2010: Dentistry. Brackets and tubes or use in orthodontics

**Please answer the following questions**

Yes      No

Is this standard recognized by FDA <sup>2</sup>? .....  Yes       No

FDA Recognition number <sup>3</sup> ..... # \_\_\_\_\_

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....  Yes       No

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....  Yes       No  
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....  Yes       No

Does this standard include acceptance criteria? .....  Yes       No  
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....  Yes       No  
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard?.....  Yes       No  
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....  Yes       No

Were deviations or adaptations made beyond what is specified in the FDA SIS?.....  Yes       No  
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....  Yes       No  
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard?.....  Yes       No  
If yes, was the guidance document followed in preparation of this 510k? .....  Yes       No

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)



**Dental Morelli Ltda.**

---

**CONFIDENTIAL**

**SECTION 11**

**EXECUTIVE SUMMARY**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

**Description of Proposed Device**

Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The brackets are uncoated. The Morelli brackets are not self-bonding and require a primer and adhesive to bond the bracket to the teeth. The Edgewise and Roth Ceramic Bracket are composed of aluminum oxide. These devices are designed for orthodontic use only. They are single use devices and are not to be reused.

Prescriptions	Model	Torque (°)		Angle (°)		In-out (mm)		With hook	No hook	Slot
		S	I	S	I	S	I			
Roth		+5	-	+12	-	0,7	-	-	1	.022"
		+8	-	+9	-	1,2	-	-	2	
		0	0	0	0	1,2	1,2	-	1,2	
		-2	-11	+9	+7	0,6	0,6	3	-	
		-2	-11	+13	+7	0,6	0,6	3	-	
		-7	-17	0	0	0,6	0,6	4	-	
		-7	-22	0	0	0,6	0,6	5	-	
Edgewise		0		0		0,94		-	1	.022"

Raw Material

(b)(4) Trade Secret Process

**Indications for Use**

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.

**Substantial Equivalence**

Both the non-clinical data and the biocompatibility evaluation indicate that Edgewise and Roth Ceramic Brackets are safe and effective for they intended use in orthodontic treatment and perform as well or better than the predicated. The subject and predicate devices have the same intended use, indications for use, compositions, device design and performance. Further details can be observed in the table below:

	Edgewise Ceramic Brackets by Dental Morelli	Roth Ceramic Brackets by Dental Morelli	Clarity Advance Ceramic Brackets by 3M Unitek Corporation K102803
Indications for use	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Clarity Advanced Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.
Target Population	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction
Anatomical Site	Teeth	Teeth	Teeth
Location of use	Use only by professional orthodontists	Use only by professional orthodontists	Use only by professional orthodontists
Materials	Aluminum oxide	Aluminum oxide	Aluminum oxide
Biocompatibility	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Compatibility with the environment and other devices	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets
Sterility	Non-sterile	Non-sterile	Non-sterile
Maxillary In-out (mm)	0.94	0.6-1.2	0.53-.089
Maxillary Torque	0	-7 to +8	-7 to +17
Maxillary Angulation	0	0 to +12	0 to +8
Mandibular In-out (mm)	0.94	0.6 – 1.2	0.51-1.14
Mandibular Torque	0	-22 to 0	-17 to 0
Slot	0.022"	0.022"	0.022"



**Dental Morelli Ltda.**

---

**Summary of Performance Testing**

The Edgewise and Roth Ceramic Brackets materials were evaluated according to ISO 10993-1. The results of this test confirmed that the Morelli Ceramic Brackets met the biomaterial compatibility requirements. A detailed test report can be seen in Section 20 of this submission.



Dental Morelli Ltda.


---

**CONFIDENTIAL**

**SECTION 12**

**Device Description**

Edgewise Ceramic Bracket  
Roth Ceramic Bracket





## Dental Morelli Ltda.

---

### Description of Proposed Device

Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The bracket is uncoated. The Morelli bracket is not self-bonding and requires a primer and adhesive to bond the bracket to the teeth. (b)

(4)Trade  
Secret  
Process

The Morelli bracket is manufactured using an injection-molding process. This process produces a stronger bracket than machining the bracket. (b)(4)Trade Secret Process


rounded edges on the brackets instead of square edges seen on machined brackets. The Edgewise Ceramic Bracket is composed of aluminum oxide.

(b)(4)Trade Secret Process



Dental Morelli Ltda.

(b)(4)Trade Secret Process



Physical and Chemical Properties of the Material:

(b)(4)Trade Secret Process



**Material**

(b)(4)Trade Secret Process





---

**Packaging**

The requirements of brackets and tubes are differentiated by a colored stripe on the packaging:

Technique	Color stripe
Roth	
Edgewise	

A more detailed explanation can be found in section 15-Proposed Labeling.




Dental Morelli Ltda.

---

**Cleaning Process**

(b)(4)Trade Secret Process





Dental Morelli Ltda.

# MSDS

	<b>Dental Morelli Ltda.</b>	SDS	
	Safety Data Sheet	1 / 2	Rev. / Date /

### Identification of the Product

1. Identification of the substance/preparation and of the company/undertaking;

1.1- Identification of the substance/preparation;

#### Aluminium Oxide Ceramic Products.

1.2- Identification of the of the company / undertaking;

#### Dental Morelli Ltda.

Alameda Jundiá, 230 / 250

Jardim Saira

ZIP 18085-090 – Sorocaba – São Paulo – Brazil.

Technical responsible: Eng. Roger. Roger Morelli

Emergency telephone number of the company and / or official advisory body in accordance with Article 12 of Directive 88/379/EEC.

In Brazil:

+55 (15) 0800-141255

+55 (15) 3238-8200

In Europe:

EUROPEAN REPRESENTATIVE

Nuno Flores

Al. Bonifácio Lázaro Lozano, 3 - Piso 0 - C

2780-125 Oeiras - Portugal

[info.morelli@euroconexao.com](mailto:info.morelli@euroconexao.com)

Tel. +351 214439292 / Fax +351 214439294

In USA:

Mr. Yesid Arias Urrea

[ariasint@bellsouth.net](mailto:ariasint@bellsouth.net)

Phone: 1-954 2362788

2. Composition/information on ingredients;

Aluminum oxide -  $Al_2O_3$ .

CAS no. 1344-28-1.

3. Hazards identification;

The information of this document if they refer to the product in your supply condition. This product is not considered dangerous, however it is possible the accidental inhalation or ingestion of parts or fragments during the use.

Eventual grinding or cut operations can generate dangerous powders.



# Dental Morelli Ltda.

SDS

## Safety Data Sheet

2 / 2

Rev. / Date  
/

4. First-aid measures;

In case of inhalation seek a doctor.

In the case of powder inhalation, to move the patient for place with fresh air and to seek a doctor.

5. Fire-fighting measures;

Not inflammable material.

6. Accidental release measures;

To clean the place.

**To observe the directives and effective laws.**

7. Handling and storage;

7.1 Handling;

Handling foreseen by capable and qualified professionals.

7.2 Storage;

To maintain in the original packing until the moment of use.

8. Exposure controls/personal protection;

Use good work practices; gloves and masks.

9. Physical and chemical properties;

Solid material.

Translucent color.

Insoluble in water.

Odorless.

10. Stability and reactivity;

Stable material in conditions of supply.

11. Toxicological information;

Not hazardous material.

12. Ecological information;

Any ecological effect is not known.

**To observe the directives and effective laws.**

13. Disposal considerations;

Use good biosafety practices.

To observe the directives and effective laws.

14. Transport information;

Transports in a safety way for do not damage the package.

15. Regulatory information;

Attempt to good work practices.

Material intended for use by qualified professionals.

16. Other information.

The document's information is based on the knowledge obtained until the present date and it is intended to a description of the product regarding its safety measures.

This description does not ensure the product's properties or characteristics.

The Dental Morelli Ltd., has rights reserved on doing technical alterations in technical specifications, with no prior advertising.



Dental Morelli Ltda.

---

# Engineering Drawings







Dental Morelli Ltda.

---

**CONFIDENTIAL**

**SECTION 13**

**SUBSTANTIAL EQUIVALENCE DISCUSSION**

Edgewise Ceramic Brackets

Roth Ceramic Brackets

**CLAIMING SUBSTANTIAL EQUIVALENCE WITH:**

**Device Name:** Clarity Advance Ceramic Brackets

**Manufacturer:** 3M Unitek Corporation

**K Number:** K102803

	Edgewise Ceramic Brackets by Dental Morelli	Roth Ceramic Brackets by Dental Morelli	Clarity Advance Ceramic Brackets by 3M Unitek Corporation K102803	Differences
Indications for use	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	Clarity Advanced Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.	
Target Population	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction	Patients in need of teeth alignment correction	
Anatomical Site	Teeth	Teeth	Teeth	
Location of use	Use only by professional orthodontists	Use only by professional orthodontists	Use only by professional orthodontists	
Materials	Aluminum oxide	Aluminum oxide	Aluminum oxide	
Biocompatibility	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	
Compatibility with the environment and other devices	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	Aluminum oxide is medical grade and is accepted for ceramic brackets	
Sterility	Non-sterile	Non-sterile	Non-sterile	
Maxillary In-out (mm)	0.94	0.6-1.2	0.53-.089	The differences do not propose and safety or effectiveness concerns. Differences are expected as each manufacturer designs their devices to meet as many patients needs
Maxillary Torque	0	-7 to +8	-7 to +17	
Maxillary Angulation	0	0 to +12	0 to +8	
Mandibular In-out (mm)	0.94	0.6 - 1.2	0.51-1.14	
Mandibular Torque	0	-22 to 0	-17 to 0	



				as possible. The differences are between the accepted values allowed for these types of devices.
Slot	0.022"	0.022"	0.022"	

The difference between the prescriptions is summarized in angulations, torque, in-out of the brackets that per action promoted by orthodontic arches results in correction of malocclusion. The orthodontist chooses the prescription that gives best results.

Dental Morelli chose manufactured ceramic brackets in the prescriptions already developed: Straight-Wire Prescription (Edgewise) and Roth Prescription.

Ceramic Brackets follow constructive characteristics of torque, angle, in-out described in BS EN ISO 27020:2010 Dentistry – Brackets and tubes for use in orthodontics.

Therefore, the differences between the prescriptions subject and predicate devices are torque, angle, in-out and length.

The application of the Brackets prescriptions is the same, the orthodontist decide what will be the best prescription to be applied to the treatment of each patient.

The maxillary in-out, maxillary torque, maxillary angulation, mandibular in-out and mandibular torque all have slight variances from the predicate. The small variances do not affect the safety or effectiveness of the subject devices. The subject devices are well within the standard of acceptance for the roth and edgewise ceramic brackets.



Dental Morelli Ltda.

---

**Predicate Device Information**

K102803 Clarity Advanced Ceramic Brackets

# Clarity™ ADVANCED Ceramic Brackets

## A Technical Perspective

by Nicole Wagner, Bill Wyllie, and Glenys Thorstenson



Nicole Wagner is a Senior Technical Service Engineer at 3M Unitek. She received her B.S. in Chemistry from the State University of New

York at Stony Brook. Her M.S. and Ph.D. are in Mechanical Engineering from the University of Minnesota, where her research focused on synthesis, characterization, and reaction modeling of hard, wear-resistant ceramic materials. She has been at 3M since 2007, joining 3M Unitek in 2010.



Bill Wyllie is a Product Development Specialist at 3M Unitek. He is active in the development of new aesthetic materials

for brackets, and has worked on a variety of product development teams such as SmartClip™ and Clarity™ SL Self-Ligating Brackets and Forsus™ Class II Correctors. He received his B.S. in Materials and Metallurgical Engineering from the University of Michigan and his M.S. and Ph.D. in Materials Engineering from Rensselaer Polytechnic Institute. He has been at 3M Unitek since 1997.



Glenys Thorstenson received her B.S. in Materials Science and Engineering from Michigan State University and her

Ph.D. in Biomedical Engineering from the University of North Carolina at Chapel Hill. Her dissertation focused on the resistance to sliding of novel orthodontic bracket systems. She has been at 3M Unitek since 2003.


### Introduction

As more patients seek a more aesthetic orthodontic treatment, orthodontists still demand functionality in their orthodontic appliances. Through various discussions with orthodontists and assistants, the 3M Unitek product development team assessed that ceramic brackets need to maintain the characteristics of aesthetics, small physical size, strength, predictable debonding, and a design that is comfortable to patients. New Clarity™ ADVANCED Ceramic Brackets incorporate these features into a revolutionary design to give orthodontists the aesthetics and efficiency they require.




### Bracket Material and Design

Advances in materials, manufacturing technologies, and bracket design have enabled new levels of performance in the Clarity brand of aesthetic brackets. Clarity ADVANCED Ceramic Brackets are made of polycrystalline alumina, which consists of small crystals, called “grains” (Figure 1A-B). As the size of these grains decreases, the strength of the ceramic material increases (Figure 2). Clarity ADVANCED Brackets are made of the same material



**Figure 1A-B** Average alumina grain size of (A) 15 µm (Clarity™ Metal-Reinforced Ceramic Bracket) and (B) 0.9 µm (Clarity™ ADVANCED Ceramic Bracket).




**Figure 2** As average grain size decreases, strength of material increases.



as that used in Clarity™ SL Self-Ligating Brackets, which is a finer-grained ceramic than the leading polycrystalline ceramic brackets. In addition, the finer grain size of the ceramic material in Clarity™ ADVANCED Ceramic Brackets improves its inherent material strength as compared to the material used in Clarity brackets. Therefore, as seen in tie-wing crush strength testing, while overall smaller in size, the strength of the Clarity ADVANCED brackets is comparable to Clarity brackets. Also, since the material is the same as that used in Clarity SL brackets, the material is proven to resist staining to various staining agents throughout the course of treatment. In addition, the translucent material of the Clarity ADVANCED brackets blends with the color of various tooth shades.

Clarity ADVANCED brackets are fabricated by an injection-molding process. This method permits the creation of smooth, rounded corners designed to reduce binding and notching at the bracket slot corners. Binding is an element of friction that contributes to the resistance to sliding when the archwire is in contact with the corners of the bracket slot. It is impacted by the materials and geometries of the archwires and brackets, and does not depend on the force applied by the ligature<sup>1</sup>.

Another factor that contributes to friction is notching, which is the resistance to sliding when the bracket permanently deforms the archwire. Most often, notching is due to the ligature force and occurs on the lingual side of the archwire. However, notching can also occur on the occlusal or gingival sides<sup>2</sup>. Images of bracket slot corners of Clarity ADVANCED brackets and other ceramic brackets that are currently on the market are shown in Figure 3A-D. The bracket slot corners of the Clarity ADVANCED brackets appear to be more rounded and smooth compared to the other ceramic brackets.



**Figure 3A-D** Bracket slot corners of (A) 3M Unitek Clarity™ ADVANCED, (B) American Radiance™, (C) GAC Mystique™, and (D) Ormco Inspire ICE™ Brackets.

To allow for an increased inter-bracket distance, Clarity ADVANCED brackets are designed to have small mesio-distal dimensions. For reduced occlusal interference, Clarity ADVANCED brackets have small occlusal-gingival dimensions. In addition, Clarity ADVANCED brackets feature a generous under-tie-wing area to allow for both single- and double-ligation. Of the doctors surveyed during a ligation study (3M Unitek), 91% found that the Clarity ADVANCED brackets easily accommodate double-ligation.

## Bonding Base and Predictable Debonding


The bonding base of Clarity ADVANCED brackets has been designed with a tooth-specific anatomy, similar to Clarity and Clarity SL brackets, to contour to the shape of patients' teeth and provide a better fit to each tooth. Also, similar to Clarity and Clarity SL brackets, Clarity ADVANCED brackets have a micro-crystalline surface on the base to create a mechanical bond with the adhesive (Figure 4).



**Figure 4** Clarity™ ADVANCED Ceramic Bracket bonding base with stress-concentrator.

To maintain the predictable debonding that orthodontists are accustomed to with the Clarity brand of brackets, the new Clarity ADVANCED ceramic brackets also feature the proprietary stress-concentrator vertically along their bracket base (Figure 4). When debonding a bracket, the adhesive first breaks at the edge of the bracket, initiating a crack that continues through the adhesive layer along the bracket base, after which the stress-concentrator collapses the bracket vertically in half<sup>3</sup>. With a mesial-distal 'rocking' motion, first towards the half of the bracket where the adhesive first broke, then towards the other half, the Clarity ADVANCED bracket can be removed entirely from the tooth. The recommended debonding tool is the same as that used for Clarity SL brackets, namely the Unitek™ Self-Ligating Bracket Debonding Instrument. To remove a bracket, this instrument is inserted in the labial side of the bracket with the instrument blade along the vertical

center slot and its ledges seated on the tie-wings. Using the mesial-distal squeeze debonding technique, Clarity™ ADVANCED Ceramic Brackets can be debonded on or off the archwire (Figure 5A-B). When debonding on the archwire, the ligature supports the collapsed bracket halves. Care should be taken to grasp and hold the collapsed bracket when debonding off the archwire.



**Figure 5A-B** Debonding Clarity™ ADVANCED Ceramic Brackets using the Unitek™ Self-Ligating Bracket Debonding Instrument either (A) on or (B) off the archwire.

## Patient Comfort

Clarity ADVANCED brackets are designed to provide enhanced patient comfort. By using an injection-molding process, smooth, rounded corners are created. The dome-shaped design and rounded bi-directional ball hooks are intended to further improve patient comfort.

The low profile design of Clarity ADVANCED brackets aims to provide patients with enhanced comfort. In addition, the low profile of lower anterior Clarity ADVANCED brackets reduces occlusal interference, giving orthodontists more flexibility to use ceramic brackets on a patient's lower arch. Clarity ADVANCED brackets have an in/out dimension that is compatible with that of Victory Series™ Low Profile Brackets.

## Conclusions

Clarity ADVANCED brackets are a new generation of ceramic brackets with both aesthetics and functionality. The aesthetics of the bracket are enabled by its translucent fine-grained alumina material and low-profile design. The smooth, rounded features of the bracket can both reduce binding and notching during treatment and assist with increasing patient comfort. With these features and the predictable debonding that remains a key feature of the Clarity brand, the Clarity ADVANCED bracket system provides an excellent aesthetic solution for both patients and orthodontists.

## References

1. Thorstenson GA, Kusy RP, "Effect of archwire size and material on the resistance to sliding of self-ligating brackets with second-order angulation in the dry state" *Am J Orthod Dentofacial Orthop* 2002; 122: 295-305.
2. Articulo LC, Kusy K, Saunders CR, Kusy RP, "Influence of ceramic and stainless steel brackets on the notching of archwires during clinical treatment" *Eur J Orthod* 2000; 22: 409-425.
3. Hansen J. "Ceramic Orthodontic Bracket with Debonding Channel" US Patent Number 5439379, 8 Aug 1995.



## Clarity™ ADVANCED Ceramic Brackets Parts List

**CLARITY™** | ADVANCED  
advanced ceramic brackets

Maxillary			IN/OUT		0.018 in. Slot		0.022 in. Slot	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Central	+17	+4	0.026	0.66	006-201	006-202	006-301	006-302
Lateral	+10	+8	0.035	0.89	006-205	006-206	006-305	006-306
Cuspid HK	0	+8	0.021	0.53	006-209	006-210	006-309	006-310
Universal Bicuspid	-7	0	0.029	0.74	006-211		006-311	
Mandibular			IN/OUT		0.018 in. Slot		0.022 in. Slot	
Tooth	Torque	Angulation	in.	mm	L	R	L	R
Anterior	-6	0	0.045	1.14	006-250		006-350	
Cuspid HK	0	+3	0.020	0.51	006-253	006-254	006-353	006-354
1st Bicuspid	-12	+2	0.030	0.76	006-255	006-256	006-355	006-356
2nd Bicuspid	-17	+2	0.034	0.86	006-259	006-260	006-359	006-360
Adding Prefix "3" is APC™ II Adhesive System								
Adding Prefix "5" is APC™ PLUS Adhesive System								

Single Patient Kits	0.018 in. Slot	0.022 in. Slot
U/L 5x5 Clarity™ ADVANCED Brackets Cuspid HK	006-100	006-110
U/L 3x3 Clarity ADVANCED Brackets Cuspid HK	006-103	006-113
U 5x5 Clarity ADVANCED Brackets Cuspid HK	006-105	006-115
U 3x3 Clarity ADVANCED Brackets Cuspid HK	006-108	006-118
Adding Prefix "3" is APC II Adhesive System		
Adding Prefix "5" is APC PLUS Adhesive System		

**APC™ II**  
Available in APC™ II Adhesive  
Coated Appliance System

**APC™ PLUS**  
Available in APC™ PLUS Adhesive  
Coated Appliance System

**MBT™**  
Versatile+ Appliance System

Page 45 of 84

Traditional 510(k) Dental Morelli Ceramic Brackets

**3M** Unitek

# Clarity™ ADVANCED Ceramic Brackets Practice Marketing Tools

**CLARITY™ | ADVANCED**  
advanced ceramic brackets

Description	Part Number
Bracket Typodont	600-236
20x Demonstration Model	600-237
Patient Brochure - Large	021-118
Patient Brochure - Small	021-119
Lenticular - Comparison With Metal Braces	012-255
Desktop Display (11x8.5 in.) - Group	014-536
Desktop Display (8.5x11 in.) - Blonde Woman	014-537
Desktop Display (8.5x11 in.) - Brunette Woman	014-538
Desktop Display (8.5x11 in.) - Dark Hair Woman	014-539
Desktop Display (8.5x11 in.) - Male	014-540
Office Poster (24x18 in.) - Group	014-541
Office Poster (18x24 in.) - Blonde Woman	014-542
Office Poster (18x24 in.) - Brunette Woman	014-543
Office Poster (18x24 in.) - Dark Hair Woman	014-544
Office Poster (18x24 in.) - Male	014-546

For product and model images to use in your practice marketing, contact 3M Unitek at 800-423-4588 and ask for extension 4303.

## Desktop Displays/Office Posters



**3M**

**3M Unitek**  
**Orthodontic Products**  
2724 South Peck Road  
Monrovia, CA 91016 USA  
www.3MUnitek.com

In U.S. and Puerto Rico: 1-800-423-4588 • 626-574-4000  
In Canada: 1-800-443-1661  
Technical Helpline: 1-800-265-1943 • 626-574-4577  
CE Hotline: 1-800-852-1990 x4649 • 626-574-4649  
Outside these areas, contact your local representative.

APC, Clarity, and MBT are trademarks of 3M.  
Please recycle. Printed in USA.  
© 2012 3M. All rights reserved.  
012-261-1 1204

Page 46 of 84

Traditional 510(k) Dental Morelli Ceramic Brackets



Dental Morelli Ltda.

---

**CONFIDENTIAL**  
**SECTION 14**  
**PROPOSED LABELING**

Edgewise Ceramic Brackets

Roth Ceramic Brackets

The requirements of brackets and tubes are differentiated by a colored stripe on the packaging:

Technique	Color stripe
Roth	
Edgewise	

Presentation forms of health product (including packaging, quantity and dimensions).


Product	Packaging Type	Quantity	Dimension
Kit	box	20 units	105mm (length) x 62mm (width) x 10mm (height)
Loose	blister	05 units	78mm (length) x 43mm (width) x 8mm (height)

NOTE – The Kit is composed by specific brackets for each tooth; loose are similar parts used for replacement.

The requirements of brackets and tubes are differentiated by a colored stripe on the packaging:

Technique	Color stripe
Roth	
Edgewise	


The packaging of the kits for the positioning of the parts:



1. Device Description- Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The bracket is uncoated. The Morelli bracket is not self-bonding and requires a primer and adhesive to bond the bracket to the teeth.

2. Rx Only

3. Indications for use: Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.



<b>ICONS USED BY DENTAL MORELLI ON EXPANSION SCREWS COMPOSITE</b>			
Legal company name and Organization address	Product registration number at the Ministry of Health / ANVISA	Code bar for product identification	
CE Mark for medical devices class II and IIa* 	Data referring to the European Representative	Code bar for manufacturing lot identification	
Lot number for traceability* 	European Representative 	Product description and instructions for use in three languages	
For exclusive use by the qualified professional** 	Attention! Read the instructions for use * 	Product expiration date (year/month)* 	06-2011
Do not use the product if the packaging has been torn* 	Product code related to the catalogue* REF	Non-reusable product* 	

<p>Non-sterile product*</p> 	<p>Recycle</p> 	<p>Manufactured by</p> 
---	--	--

**The symbols employed vary according to the need under consideration of the technical characteristics of each product.**

\* Symbols employed in compliance with Harmonized Standards named herein (DIN EN 980 and ABNT NBR ISO 15223).


\*\* Symbols adapted to complete the given information with the purpose of minimizing eventual risks inherent to the product use.

# Edgewise Ceramic



*Aesthetic line*

The orthodontic "Ceramic" brackets by Morelli incorporate beauty, functionality and strength into a sophisticated and sleek device.




The brackets are manufactured from a polycrystalline ceramic. The brackets are translucent, chemically stable, biocompatible, do not wear out or lose their color.

The design has been specially created to provide great comfort to the patient, emphasizing the rounded forms and low profile.

The base is produced of ridges in order to anatomically adhere to the cervico occlusal. This design allows for greater adhesion to the enamel without the undesirable removal effects.

The Ceramic brackets offer a high degree of sliding, due the slot with rounding in the ends and the exclusive polishing process for ceramics. Thus, the free movement is ensured even under a important misalignment of teeth.




# Edgewise Ceramic

Aesthetic line



Strength  
Beauty  
Discretion

Length	3,2 mm	3,2 mm	3,2 mm	2,9 mm	3,6 mm	3,6 mm	2,9 mm	3,2 mm	3,2 mm	3,2 mm
Torq	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/out	1,20 mm									
Replacements 022" x .030"	10.18.013	10.18.013	10.18.013	10.18.012	10.18.011	10.18.011	10.18.012	10.18.013	10.18.013	10.18.013



Length	3,2 mm	3,2 mm	3,2 mm	2,9 mm	2,9 mm	2,9 mm	2,9 mm	3,2 mm	3,2 mm	3,2 mm
Torq	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
Angle	0°	0°	0°	0°	0°	0°	0°	0°	0°	0°
In/out	1,20 mm									
Replacements 022" x .030"	10.18.013	10.18.013	10.18.013	10.18.012	10.18.012	10.18.012	10.18.012	10.18.013	10.18.013	10.18.013

The Morelli Ceramic Brackets are available through Morelli or over 150 distributors/re-salers throughout Brazil.



SISTEMA DE QUALIDADE  
CERTIFICADO POR NORMAS  
INTERNACIONAIS E  
REGULAMENTOS DO  
MINISTÉRIO DA SAÚDE

DENTAL MORELLI LTDA  
Alameda Jundiá, 230/250  
Sorocaba - SP - Brasil - CEP: 18085-090  
Serviço de atendimento ao cliente  
(55) 15 - 3238-8200 / 0800 12-1455



# Edgewise Ceramic

Aesthetic line



## Rx Only

## D evice Description

Edgewise and Roth Ceramic Brackets are intended to be bonded to teeth, upon which an orthodontic wire is used to move the teeth to new positions. The Morelli bracket consists of a translucent alumina body. The bracket is uncoated. The Morelli bracket is not self-bonding and requires a primer and adhesive to bond the bracket to the teeth.

## I ndications for use

Edgewise and Roth Ceramic Brackets are intended for use in orthodontic treatments. The brackets are affixed to teeth so that pressure can be exerted on the teeth.



Dental Morelli Ltda.

---

**CONFIDENTIAL**

**SECTION 15**

**STERILIZATION AND SHELF LIFE**

(b)(4)Trade Secret Process





Dental Morelli Ltda.


---

**CONFIDENTIAL**

**SECTION 16**

**BIOCOMPATIBILITY**

(b)(4)Trade Secret Process





Dental Morelli Ltda.


---

**CONFIDENTIAL**

**SECTION 17**

**SOFTWARE VALIDATION**

(b)(4)Trade Secret Process





**Confidential**

**Section 18**

**Electromagnetic Safety**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

This section does not apply.



Dental Morelli Ltda.

---

**CONFIDENTIAL**

**SECTION 19**

**BENCH TESTING**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets















































**Dental Morelli Ltda.**

---

**CONFIDENTIAL**

**SECTION 20**

**ANIAML TESTING**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

The Edgewise and Roth Ceramic Brackets did not involve animal testing. This section does not apply.



**CONFIDENTIAL**

**SECTION 21**

**CLINICAL TRIALS**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

This is a traditional 510(k) premarket clearance. Certification of Compliance under 42 U.S.C. does not apply.

The Edgewise and Roth Ceramic Brackets did not involve clinical trials. This section does not apply.



DEPARTMENT OF HEALTH AND HUMAN SERVICES  
Food and Drug Administration

**Certification of Compliance, under 42 U.S.C. § 282(j)(5)(B), with Requirements of ClinicalTrials.gov Data Bank (42 U.S.C. § 282(j))**

(For submission with an application/submission, including amendments, supplements, and resubmissions, under §§ 505, 515, 520(m), or 510(k) of the Federal Food, Drug, and Cosmetic Act or § 351 of the Public Health Service Act.)

**SPONSOR / APPLICANT / SUBMITTER INFORMATION**

1. NAME OF SPONSOR/APPLICANT/SUBMITTER Dental Morelli LTDA	2. DATE OF THE APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES Feb 6, 2013
3. ADDRESS (Number, Street, State, and ZIP Code) Alameda Jundiai, 230- Jardim Saira-Sorocaba CEP: 18085-090 Brazil	4. TELEPHONE AND FAX NUMBERS (Include Area Code) (Tel.) 55 15 3238-8200 (Fax)

**PRODUCT INFORMATION**

5. **FOR DRUGS/BIOLOGICS:** Include Any/All Available Established, Proprietary and/or Chemical/Biochemical/Blood/Cellular/Gene Therapy Product Name(s)  
**FOR DEVICES:** Include Any/All Common or Usual Name(s), Classification, Trade or Proprietary or Model Name(s) and/or Model Number(s)  
(Attach extra pages as necessary)

Edgewise Ceramic Brackets

Roth Ceramic Brackets

**APPLICATION / SUBMISSION INFORMATION**

6. TYPE OF APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES  
 IND     NDA     ANDA     BLA     PMA     HDE     510(k)     PDP     Other

7. INCLUDE IND/NDA/ANDA/BLA/PMA/HDE/510(k)/PDP/OTHER NUMBER (If number previously assigned)  
N/A

8. SERIAL NUMBER ASSIGNED TO APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES  
N/A

**CERTIFICATION STATEMENT / INFORMATION**

9. CHECK ONLY ONE OF THE FOLLOWING BOXES (See instructions for additional information and explanation)

A. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, do not apply because the application/submission which this certification accompanies does not reference any clinical trial.

B. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, do not apply to any clinical trial referenced in the application/submission which this certification accompanies.

C. I certify that the requirements of 42 U.S.C. § 282(j), Section 402(j) of the Public Health Service Act, enacted by 121 Stat. 823, Public Law 110-85, apply to one or more of the clinical trials referenced in the application/submission which this certification accompanies and that those requirements have been met.

10. IF YOU CHECKED BOX C, IN NUMBER 9, PROVIDE THE NATIONAL CLINICAL TRIAL (NCT) NUMBER(S) FOR ANY "APPLICABLE CLINICAL TRIAL(S)," UNDER 42 U.S.C. § 282(j)(1)(A)(i), SECTION 402(j)(1)(A)(i) OF THE PUBLIC HEALTH SERVICE ACT, REFERENCED IN THE APPLICATION/SUBMISSION WHICH THIS CERTIFICATION ACCOMPANIES (Attach extra pages as necessary)

NCT Number(s):

The undersigned declares, to the best of her/his knowledge, that this is an accurate, true, and complete submission of information. I understand that the failure to submit the certification required by 42 U.S.C. § 282(j)(5)(B), section 402(j)(5)(B) of the Public Health Service Act, and the knowing submission of a false certification under such section are prohibited acts under 21 U.S.C. § 331, section 301 of the Federal Food, Drug, and Cosmetic Act.  
**Warning:** A willfully and knowingly false statement is a criminal offense, U.S. Code, title 18, section 1001.

11. SIGNATURE OF SPONSOR/APPLICANT/SUBMITTER OR AN AUTHORIZED REPRESENTATIVE (Sign) 	12. NAME AND TITLE OF THE PERSON WHO SIGNED IN NO. 11 (Name) <u>Lilian Clull</u> (Title) <u>Senior Partner</u>
13. ADDRESS (Number, Street, State, and ZIP Code) (of person identified in Nos. 11 and 12) <u>1885 NE 29th Ave Suite 720</u> <u>Aventura FL 33180</u>	14. TELEPHONE AND FAX NUMBERS (Include Area Code) (Tel.) <u>305-377-0077</u> (Fax)
15. DATE OF CERTIFICATION <u>3/28/13</u>	

**Barlow, Lenny \***

---

**From:** Barlow, Lenny \*  
**Sent:** Thursday, August 29, 2013 3:45 PM  
**To:** 'taraconrad@techlinkusa.net'  
**Cc:** DCCLetters  
**Subject:** k131197 Correspondence  
**Attachments:** k131197.pdf

<b>Tracking:</b>	<b>Recipient</b>	<b>Delivery</b>
	'taraconrad@techlinkusa.net' DCCLetters	Delivered: 8/29/2013 3:45 PM



# COVER SHEET MEMORANDUM

Food and Drug Administration  
Office of Device Evaluation &  
Office of In Vitro Diagnostics and  
Radiological Health

**From:** Reviewer Name Myra E. Browne  
**Subject:** 510(k) Number K131197  
**To:** The Record

**Please list CTS decision code:** SE - Substantially Equivalent

- Refused to Accept (Note: this is considered the first review cycle. See [screening checklist](#).)
- Hold (Additional Information or Telephone Hold)
- Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.)

Please complete the following for a final clearance decision (i.e. SE, SE with Limitations, etc.)	YES	NO
Indications for Use Page ( <i>Attach IFU</i> )	X	
510(k) Summary or 510(k) Statement ( <i>Attach Summary or Statement</i> )	X	
Truthful and Accurate Statement ( <i>Must be present for a Final Decision</i> )	X	
Is the device Class III?		X
Does firm reference standards? (If yes, please attach <a href="#">Form 3654</a> .)		X
Is this a combination product?		X
Is this a reprocessed single use device? (See <a href="#">Guidance for Industry and FDA Staff - MDUFMA - Validation Data in 510(k)s for Reprocessed Single-Use Medical Devices</a> .)		X
Is this device intended for pediatric use only?		X
Is this a prescription device? (If both prescription & OTC, check both boxes.)	X	
Is clinical data necessary to support the review of this 510(k)?		X
For United States based clinical studies only, did the application include a completed Form FDA 3674, Certification with Requirements of <a href="#">ClinicalTrials.gov</a> Data Bank? (If study was conducted in the United States and Form FDA 3674 was not included or was incomplete, then applicant must be contacted to obtain completed form.)		X
Does this device include an Animal Tissue Source?		X
All Pediatric Patients age <= 21		X
Neonate/Newborn (Birth to 28 days)		X
Infant (29 days to < 2 years)		X
Child (2 years to <12 years)	X	
Adolescent (12 years to <18 years)	X	
Transitional Adolescent A (18 years to <21 years); Special considerations are being given to this group, different from adults age >= 21 (different device design or testing, different protocol procedures, etc.)		X
Transitional Adolescent B (18 years to <21 years); No special considerations compared to adults >= 21 years)		X

Nanotechnology		×
Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance)		×

<b>Regulation Number:</b>	872.5470
<b>Class:</b>	II
<b>Product Code:</b>	NJM
<b>Additional Product Codes:</b>	

**Digital Signature Concurrence Table**  
(Not all signatures may be required)

Branch Chief Sign-Off	<p align="center">Sheena A. Green -S 2013.08.27 14:27:54 -04'00'</p>
Division Sign-Off	<p align="center">Mary S. Runner -S 2013.08.27 14:49:28 -04'00'</p> <p><i>Susan Runner DDS, MA</i></p>

K131197/9001



Dental Morelli Ltda.

---

**CONFIDENTIAL**

**COVER LETTER**

April 15, 2013

Food and Drug Administration,  
Center for Devices and Radiological Health  
Document Mail Center W066-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002 – USA

FDA CDRH DMC

JUN 04 2013

Received

RE: 510(k) Traditional Pre Market Notification Request

**COMPANY NAME AND ADDRESS**

DENTAL MORELLI LTDA  
Alameda Jundiaí, 230 –  
Jardim Saira - Sorocaba  
CEP: 18085-090  
Brazil

Telephone: 55 (15) 3238-8200

**CONSULTANT NAME AND ADDRESS**

TechLink International Consulting  
18851 NE 29<sup>th</sup> Avenue  
Suite 720  
Aventura, Florida 33180

Telephone: (305) 377-0077

Primary Contact: Tara Conrad  
Secondary Contact: Lilian Llull

**Attention: Document Control Clerk**

**Device Trade Name:**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

According to Section 510 (k) of the Federal Food, Drug and Cosmetic Act, as amended (ACT), Dental Morelli Ltda proposes to introduce the Orthodontic Ceramic

---



## Dental Morelli Ltda.

Brackets into interstate commerce for commercial distribution and hereby requests 510 (k) clearance by the Food and Drug Administration, as required by law.

The following information on these products (according to 21 CFR 807.87) is submitted for your consideration.

**Common Name:** Orthodontic Ceramic Brackets  
**Trade Name:** Edgewise Ceramic Brackets; Roth Ceramic Brackets  
**Classification:** Class II  
**Product Code:** NJM  
**Classification Panel:** Dental  
**Regulation Numbers:** 21 CFR 872.5470  
**Substantial Equivalence:** K102803 Clarity Advanced Ceramic Brackets

The eCopy is an exact duplicate of the paper copy.

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain components derived from a tissue or other biologic source?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device provided sterile?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device intended for single use?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device a reprocessed single use device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, does this device type require reprocessed validation data?	N/A	N/A
Does the device contain a drug?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain a biologic?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device use software?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the submission include clinical information?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device implanted?	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## Dental Morelli Ltda.

---

The predicate device and proposed device have the same intended use, same indications, and are designed with the same technological characteristics. A complete list of indications and a comparison table are included in Section 13. Indications for Use Statement can be found in Section 5. Labeling specifications are detailed in section 14. Considering our intent to market these devices, all information submitted is confidential; with the exception of the 510(k) summary.

Dental Morelli has not submitted prior applications for the same device. However, there are other manufacturers who have submitted 510(k) applications for similar devices. These similar devices are approved by the FDA. In particular the Clarity Advanced Ceramic Brackets that has been used as a predicate.

We are confident that this information will be sufficient for you to reach a favorable decision. However, please do not hesitate to contact me if you have any additional questions, concerns, or if you feel that I can be of any further assistance.

Best regards,

A handwritten signature in blue ink, appearing to read 'Tara Conrad', is written over a light blue circular scribble.

Tara Conrad, Biomedical Engineer  
Regulatory Affairs Manager

(305) 377-0077 – ph  
taraconrad@techlinkusa.net

---







Dental Morelli Ltda.

---

**CONFIDENTIAL**

**COVER LETTER**

April 15, 2013

Food and Drug Administration  
Center for Devices and Radiological Health  
Document Mail Center W066-G609  
10903 New Hampshire Avenue  
Silver Spring, MD 20993-0002 – USA

RE: 510(k) Traditional Pre Market Notification Request

**COMPANY NAME AND ADDRESS**

DENTAL MORELLI LTDA  
Alameda Jundiaí, 230 –  
Jardim Saira - Sorocaba  
CEP: 18085-090  
Brazil

Telephone: 55 (15) 3238-8200

**CONSULTANT NAME AND ADDRESS**

TechLink International Consulting  
18851 NE 29<sup>th</sup> Avenue  
Suite 720  
Aventura, Florida 33180

Telephone: (305) 377-0077

Primary Contact: Tara Conrad  
Secondary Contact: Lilian Llull

**Attention: Document Control Clerk**

**Device Trade Name:**

Edgewise Ceramic Brackets  
Roth Ceramic Brackets

According to Section 510 (k) of the Federal Food, Drug and Cosmetic Act, as amended (ACT), Dental Morelli Ltda proposes to introduce the Orthodontic Ceramic



## Dental Morelli Ltda.

Brackets into interstate commerce for commercial distribution and hereby requests 510 (k) clearance by the Food and Drug Administration, as required by law.

The following information on these products (according to 21 CFR 807.87) is submitted for your consideration.

**Common Name:** Orthodontic Ceramic Brackets  
**Trade Name:** Edgewise Ceramic Brackets; Roth Ceramic Brackets  
**Classification:** Class II  
**Product Code:** NJM  
**Classification Panel:** Dental  
**Regulation Numbers:** 21 CFR 872.5470  
**Substantial Equivalence:** K102803 Clarity Advanced Ceramic Brackets

The eCopy is an exact duplicate of the paper copy.

Question	YES	NO
Is the device intended for prescription use (21 CFR 801 Subpart D)?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain components derived from a tissue or other biologic source?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device provided sterile?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device intended for single use?	<input checked="" type="checkbox"/>	<input type="checkbox"/>
Is the device a reprocessed single use device?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
If yes, does this device type require reprocessed validation data?	N/A	N/A
Does the device contain a drug?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device contain a biologic?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the device use software?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Does the submission include clinical information?	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Is the device implanted?	<input type="checkbox"/>	<input checked="" type="checkbox"/>



## Dental Morelli Ltda.

---

The predicate device and proposed device have the same intended use, same indications, and are designed with the same technological characteristics. A complete list of indications and a comparison table are included in Section 13. Indications for Use Statement can be found in Section 5. Labeling specifications are detailed in section 14. Considering our intent to market these devices, all information submitted is confidential; with the exception of the 510(k) summary.

Dental Morelli has not submitted prior applications for the same device. However, there are other manufacturers who have submitted 510(k) applications for similar devices. These similar devices are approved by the FDA. In particular the Clarity Advanced Ceramic Brackets that has been used as a predicate.

We are confident that this information will be sufficient for you to reach a favorable decision. However, please do not hesitate to contact me if you have any additional questions, concerns, or if you feel that I can be of any further assistance.

Best regards,

A handwritten signature in blue ink, appearing to read 'Tara Conrad', is written over a light blue circular scribble.

Tara Conrad, Biomedical Engineer  
Regulatory Affairs Manager

(305) 377-0077 – ph  
taraconrad@techlinkusa.net