



# U.S. Department of Health & Human Services

---

**Food and Drug Administration**

## SAVE REQUEST

**USER:** (kml)  
**FOLDER:** K140300 - 269 pages  
**COMPANY:** ROTATION MEDICAL, INC. (ROTAMEDI)  
**PRODUCT:** MESH, SURGICAL, COLLAGEN, ORTHOPAEDICS, REINFORCEMENT OF TENDON (OWY)  
**SUMMARY:** Product: COLLAGEN TENDON SHEET-DDI (CTS-DDI)

**DATE REQUESTED:** Nov 4, 2015

**DATE PRINTED:** Nov 4, 2015

**Note:** Printed





## 510(k) Summary

### Applicant Information

**Applicant Name:** Rotation Medical, Inc.  
**Applicant Address:** 15350 25<sup>th</sup> Avenue North, Suite 100  
 Plymouth, MN 55447  
**Telephone:** 763-746-7502  
**Fax:** 763-746-7501  
**Contact Person:** Jeff Sims  
 Vice President, Clinical Programs and Regulatory Affairs  
**Date Prepared:** March 19, 2014

### Name of Device

**Device Common Name:** Tendon Protector  
**Device Trade Name:** Collagen Tendon Sheet-DDI  
**Device Classification Name:** Mesh, Surgical, Collagen, Orthopaedics,  
 Reinforcement of Tendon  
 878.3300  
 Class II  
 OWY  
**Secondary Classification Code** ORQ

### Legally Marketed Devices to Which Substantial Equivalence is Claimed

**Predicate Device(s):** Collagen Tendon Sheet-D, K122048  
 Rotation Medical, Inc.

### Description of the Device

Collagen Tendon Sheet-DDI is a resorbable type I collagen matrix that provides a layer of collagen over injured tendons. Collagen Tendon Sheet-DDI is designed to provide a layer between the tendon and the surrounding tissue during healing. When hydrated, Collagen Tendon Sheet-DDI is an easy-to-use, soft, pliable, nonfriable, porous collagen sheet. Collagen Tendon Sheet-DDI is provided sterile, non-pyrogenic, for single use only, in a variety of sizes, and is packaged preloaded in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration.



### **Intended Use**

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### **Summary/Comparison of Technical Characteristics**

Collagen Tendon Sheet-DDI is exactly the same product as its predicate, Collagen Tendon Sheet-D (K122048) in regard to all aspects of the collagen scaffold. Only the packaging has changed between the current device and the predicate device. Collagen Tendon Sheet-DDI is packaged in a cartridge to facilitate arthroscopic delivery and positioning of the collagen matrix.

Past safety, biocompatibility and mechanical characterization tests of the predicate product are directly applicable to the current product. Additional tests to confirm the performance characteristics of the packaging change were performed. To confirm performance of the device in the new packaging; suture pull-out testing, hydrothermal transition temperature, and endotoxin testing were evaluated. To confirm the performance of the cartridge assembly itself biocompatibility, simulated use and mechanical integrity tests were evaluated.

In summary, the Collagen Tendon Sheet-DDI device passed all applicable FDA Blue Book Memorandum G95-1 and ISO 10993-1 testing for the biological evaluation of medical devices. Testing was conducted in accordance with FDA's Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh.

### **Conclusion**

The purpose of this 510(k) application was to notify the Food and Drug Administration of proposed modifications to the packaging of the previously cleared Collagen Tendon Sheet-D (K122048). Collagen Tendon Sheet-DDI and Collagen Tendon Sheet-D are exactly the same device excepting their respective packaging. The proposed packaging modifications for the new device raise no new questions regarding safety and effectiveness. Therefore, the proposed Collagen Tendon Sheet-DDI is substantially equivalent to the predicate Collagen Tendon Sheet-D.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

March 26, 2014

Rotation Medical Incorporated  
Mr. Jeff Sims  
Vice President, Clinical Programs and Regulatory Affairs  
15350 25<sup>th</sup> Avenue North, Suite 100  
Plymouth, Minnesota 55447

Re: K140300  
Trade/Device Name: Collagen Tendon Sheet-DDI  
Regulation Number: 21 CFR 878.3300  
Regulation Name: Surgical Mesh  
Regulatory Class: Class II  
Product Code: OWY, ORQ  
Dated: February 3, 2014  
Received: February 6, 2014

Dear Mr. Sims:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**David Krause -S**

for Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K140300

Device Name  
Collagen Tendon Sheet-DDI

*Indications for Use (Describe)*

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

Peter L. Hudson -S

2014.03.25 14:50:46 -04'00'

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
10903 New Hampshire Avenue  
Document Control Center - WO66-G609  
Silver Spring, MD 20993-0002

March 26, 2014

Rotation Medical Incorporated  
Mr. Jeff Sims  
Vice President, Clinical Programs and Regulatory Affairs  
15350 25<sup>th</sup> Avenue North, Suite 100  
Plymouth, Minnesota 55447

Re: K140300  
Trade/Device Name: Collagen Tendon Sheet-DDI  
Regulation Number: 21 CFR 878.3300  
Regulation Name: Surgical Mesh  
Regulatory Class: Class II  
Product Code: OWY, ORQ  
Dated: February 3, 2014  
Received: February 6, 2014

Dear Mr. Sims:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration. Please note: CDRH does not evaluate information related to contract liability warranties. We remind you; however, that device labeling must be truthful and not misleading.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); medical device reporting (reporting of medical

device-related adverse events) (21 CFR 803); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). For questions regarding the reporting of adverse events under the MDR regulation (21 CFR Part 803), please go to <http://www.fda.gov/MedicalDevices/Safety/ReportaProblem/default.htm> for the CDRH's Office of Surveillance and Biometrics/Division of Postmarket Surveillance.

You may obtain other general information on your responsibilities under the Act from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 796-7100 or at its Internet address <http://www.fda.gov/MedicalDevices/ResourcesforYou/Industry/default.htm>.

Sincerely yours,

**David Krause -S**

for Binita S. Ashar, M.D., M.B.A., F.A.C.S.  
Acting Director  
Division of Surgical Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

**Indications for Use**

510(k) Number (if known)  
K140300

Device Name  
Collagen Tendon Sheet-DDI

Indications for Use (Describe)

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

Type of Use (Select one or both, as applicable)

Prescription Use (Part 21 CFR 801 Subpart D)

Over-The-Counter Use (21 CFR 801 Subpart C)

**PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON A SEPARATE PAGE IF NEEDED.**

**FOR FDA USE ONLY**

Concurrence of Center for Devices and Radiological Health (CDRH) (Signature)

**Peter L. Hudson -S**

**2014.03.25 14:50:46 -04'00'**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 79 hours per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov

*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB number."*

**Collins, Virginia \***

---

**From:** Collins, Virginia \*  
**Sent:** Thursday, March 27, 2014 4:25 PM  
**To:** 'jsims@rotationmedical.com'  
**Cc:** DCCLetters  
**Subject:** K140300 SE Letter  
**Attachments:** K140300.pdf

**Tracking:**

**Recipient**

'jsims@rotationmedical.com'

DCCLetters

**Delivery**

Delivered: 3/27/2014 4:25 PM



K140300

February 3, 2014

FDA CDRH DMC

VIA FEDERAL EXPRESS

Document Mail Center (WO66-G609)  
Center for Devices and Radiological Health  
Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

FEB 06 2014

Received

Re: 510(k) Premarket Notification (Traditional)  
Collagen Tendon Sheet-DDI (b) (4)

Note: This submission includes an electronic copy of the 510(k) submission per *eCopy Program for Medical Device Submissions*. The eCopy is an exact duplicate of the paper copy.

Dear Sir or Madam:

Pursuant to 21 CFR Part 807, Subpart E, Premarket Notification Procedures, Section 807.81, Rotation Medical, Inc. is submitting this 510(k) Premarket Notification for its Collagen Tendon Sheet-DDI. Rotation Medical has determined that Collagen Tendon Sheet-DDI is substantially equivalent to current legally marketed tendon protector devices and intends to manufacture and market the device.

(b) (4)

<b>Trade Name or Proprietary Name:</b>	Collagen Tendon Sheet-DDI
<b>Common or Usual Name:</b>	Tendon Protector
<b>Device Classification Name:</b>	Mesh, Surgical
<b>Regulation Number:</b>	878.3300
<b>Product Code:</b>	FTM
<b>Device Class:</b>	Class II
<b>Name and Address of Manufacturer:</b>	Rotation Medical, Inc. 15350 25 <sup>th</sup> Avenue North, Suite 100 Plymouth, MN 55447

---

15350 25<sup>th</sup> Avenue No • Suite 100 • Plymouth MN 55447 • 763.746.7500



**Establishment Registration No.:** 3009351468

**Name, Address, and Telephone Number of Contact Person:** Jeff Sims  
 Vice President, Clinical Programs and Regulatory Affairs  
 Rotation Medical, Inc.  
 15350 25<sup>th</sup> Avenue North, Suite 100  
 Plymouth, MN 55447  
 Tel: 763.746.7502  
 Fax: 763.746.7501  
 E-mail: jsims@rotationmedical.com

Below is the recommended "Design and Use of the Device" questions per the *Guidance for Industry and FDA Staff: Format for Traditional and Abbreviated 510(k)s*.

Question	Yes	No
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X	
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?		X
Is the device provided sterile?	X	
Is the device intended for single use?	X	
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		n/a
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?	X	

The purpose of this application is to gain clearance for Collagen Tendon Sheet-DDI, an additional model of the predicate device, Rotation Medical's Collagen Tendon Sheet-D (K122048). Collagen Tendon Sheet-DDI is exactly the same device as its predicate except it is

(b) (4)

The indication for use for Collagen Tendon Sheet-DDI is the same as the predicate Collagen Tendon Sheet-D:

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.



Arthroscopic tendon repair procedures constitute a growing percentage of all such procedures.

(b) (4)

If you have any questions, please do not hesitate to contact me. Thank you all in advance for your work toward the clearance of this device.

Sincerely,

A handwritten signature in blue ink, appearing to read "Jeff Sims", with a stylized flourish at the end.

Jeff Sims  
Vice President, Clinical Programs and Regulatory Affairs

Enclosures (submitted including an electronic copy)

**Premarket Notification  
510(k) Application**

**For**

**Collagen Tendon Sheet-DDI**

(b) (4)

**February 3, 2014**

**Rotation Medical, Inc.  
15350 25<sup>th</sup> Avenue North, Suite 100  
Plymouth, Minnesota 55447**

**Premarket Notification, 510(k) Application  
for  
Collagen Tendon Sheet-DDI**

(b) (4)

**Table of Contents**

1.0	Medical Device User Fee Cover Sheet (Form FDA 3601).....	Page 4
2.0	CDRH Premarket Review Submission Cover Sheet .....	6
3.0	510(k) Cover Letter.....	11
4.0	Indications for Use Statement.....	14
5.0	510(k) Summary .....	15
6.0	Truthful and Accuracy Statement.....	17
7.0	Class III Summary and Certification .....	18
8.0	Financial Certification or Disclosure Statement.....	19
9.0	Declarations of Conformity and Summary Reports .....	20
10.0	Executive Summary.....	21
11.0	Device Description.....	23
12.0	Substantial Equivalence Discussion .....	30
13.0	Proposed Labeling.....	39
14.0	Sterilization and Shelf Life .....	40
15.0	Biocompatibility.....	43
16.0	Software .....	44
17.0	Electromagnetic Compatibility and Electrical Safety.....	45
18.0	Performance Testing – Bench .....	46
19.0	Performance Testing – Animal.....	48
20.0	Performance Testing – Clinical .....	49

## APPENDICES

Appendix	A	(b) (4)
Appendix	B	Biocompatibility Test Reports
Appendix	C	Proposed Label and Instructions for Use
Appendix	D	Shelf Life Study Report
Appendix	E	Design Verification Records
Appendix	F	Standards Data Reports (Form FDA 3654)
Appendix	G	Simulated Use and Cartridge Design Integrity Test Reports

DEPARTMENT OF HEALTH AND HUMAN SERVICES FOOD AND DRUG ADMINISTRATION <b>MEDICAL DEVICE USER FEE COVER SHEET</b>	PAYMENT IDENTIFICATION NUMBER: (b) (4) Write the Payment Identification number on your check.
---	---

A completed cover sheet must accompany each original application or supplement subject to fees. If payment is sent by U.S. mail or courier, please include a copy of this completed form with payment. Payment and mailing instructions can be found at: <http://www.fda.gov/oc/mdufma/coversheet.html>

1. COMPANY NAME AND ADDRESS (include name, street address, city state, country, and post office code)  ROTATION MEDICAL, INC. FKA DENALI MEDICAL, INC. 15350 25th AVENUE N., SUITE 100 PLYMOUTH MN 55447 US  1.1 EMPLOYER IDENTIFICATION NUMBER (EIN) (b) (4)	2. CONTACT NAME Jeff Sims 2.1 E-MAIL ADDRESS jsims@rotationmedical.com 2.2 TELEPHONE NUMBER (include Area code) 763-746-7502 2.3 FACSIMILE (FAX) NUMBER (Include Area code) 763-746-7501
--	---

3. TYPE OF PREMARKET APPLICATION (Select one of the following in each column; if you are unsure, please refer to the application descriptions at the following web site: <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/ucm345263.htm</a> Select an application type: <input checked="" type="checkbox"/> Premarket notification(510(k)); except for third party <input type="checkbox"/> 513(g) Request for Information <input type="checkbox"/> Biologics License Application (BLA) <input type="checkbox"/> Premarket Approval Application (PMA) <input type="checkbox"/> Modular PMA <input type="checkbox"/> Product Development Protocol (PDP) <input type="checkbox"/> Premarket Report (PMR) <input type="checkbox"/> 30-Day Notice	3.1 Select a center <input checked="" type="checkbox"/> CDRH <input type="checkbox"/> CBER 3.2 Select one of the types below <input checked="" type="checkbox"/> Original Application Supplement Types: <input type="checkbox"/> Efficacy (BLA) <input type="checkbox"/> Panel Track (PMA, PMR, PDP) <input type="checkbox"/> Real-Time (PMA, PMR, PDP) <input type="checkbox"/> 180-day (PMA, PMR, PDP)
--	---

4. ARE YOU A SMALL BUSINESS? (See the instructions for more information on determining this status) <input type="checkbox"/> YES, I meet the small business criteria and have submitted the required qualifying documents to FDA <input checked="" type="checkbox"/> NO, I am not a small business	4.1 If Yes, please enter your Small Business Decision Number:
--	---

5. FDA WILL NOT ACCEPT YOUR SUBMISSION IF YOUR COMPANY HAS NOT PAID AN ESTABLISHMENT REGISTRATION FEE THAT IS DUE TO FDA. HAS YOUR COMPANY PAID ALL ESTABLISHMENT REGISTRATION FEES THAT ARE DUE TO FDA? <input checked="" type="checkbox"/> YES (All of our establishments have registered and paid the fee, or this is our first device, and we will register and pay the fee within 30 days of FDA's approval/clearance of this device.) <input type="checkbox"/> NO (If "NO," FDA will not accept your submission until you have paid all fees due to FDA. This submission will not be processed; see <a href="http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/RegistrationandListing/ucm053165.htm">http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/HowtoMarketYourDevice/RegistrationandListing/ucm053165.htm</a> for additional information)	
---	--

6. IS THIS PREMARKET APPLICATION COVERED BY ANY OF THE FOLLOWING USER FEE EXCEPTIONS? IF SO, CHECK THE APPLICABLE EXCEPTION.	
<input type="checkbox"/> This application is the first PMA submitted by a qualified small business, including any affiliates <input type="checkbox"/> This biologics application is submitted under section 351 of the Public Health Service Act for a product licensed for further manufacturing use only	<input type="checkbox"/> The sole purpose of the application is to support conditions of use for a pediatric population <input type="checkbox"/> The application is submitted by a state or federal government entity for a device that is not to be distributed commercially

7. IS THIS A SUPPLEMENT TO A PREMARKET APPLICATION FOR WHICH FEES WERE WAIVED DUE TO SOLE USE IN A PEDIATRIC POPULATION THAT NOW PROPOSES CONDITION OF USE FOR ANY ADULT POPULATION? (If so, the application is subject to the fee that applies for an original premarket approval application (PMA)). <input type="checkbox"/> YES <input checked="" type="checkbox"/> NO
---

**PAPERWORK REDUCTION ACT STATEMENT**  
 Public reporting burden for this collection of information is estimated to average 18 minutes per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information,

including suggestions for reducing this burden, to the address below.

Department of Health and Human Services, Food and Drug Administration, Office of Chief Information Officer, 1350 Piccard Drive,  
4th Floor Rockville, MD 20850

[Please do NOT return this form to the above address, except as it pertains to comments on the burden estimate.]

8. USER FEE PAYMENT AMOUNT SUBMITTED FOR THIS PREMARKET APPLICATION

(b) (4)

27-Jan-2014

Form FDA 3601 (01-2007)

["Close Window"](#) [Print Cover sheet](#)

**CDRH PREMARKET REVIEW SUBMISSION COVER SHEET**

Form Approval  
 OMB No. 0910-0120  
 Expiration Date: December 31, 2013  
 See PRA Statement on page 5.

Date of Submission  
 February 3, 2014

User Fee Payment ID Number

(b) (4)

FDA Submission Document Number (if known)

**SECTION A**

**TYPE OF SUBMISSION**

<p><b>PMA</b></p> <input type="checkbox"/> Original Submission <input type="checkbox"/> Premarket Report <input type="checkbox"/> Modular Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment <input type="checkbox"/> Licensing Agreement	<p><b>PMA &amp; HDE Supplement</b></p> <input type="checkbox"/> Regular (180 day) <input type="checkbox"/> Special <input type="checkbox"/> Panel Track (PMA Only) <input type="checkbox"/> 30-day Supplement <input type="checkbox"/> 30-day Notice <input type="checkbox"/> 135-day Supplement <input type="checkbox"/> Real-time Review <input type="checkbox"/> Amendment to PMA & HDE Supplement <input type="checkbox"/> Other	<p><b>PDP</b></p> <input type="checkbox"/> Original PDP <input type="checkbox"/> Notice of Completion <input type="checkbox"/> Amendment to PDP	<p><b>510(k)</b></p> <input checked="" type="checkbox"/> Original Submission: <input checked="" type="checkbox"/> Traditional <input type="checkbox"/> Special <input type="checkbox"/> Abbreviated (Complete section I, Page 5) <input type="checkbox"/> Additional Information <input type="checkbox"/> Third Party	<p><b>Request for Feedback</b></p> <input type="checkbox"/> Pre-Submission <input type="checkbox"/> Informational Meeting <input type="checkbox"/> Submission Issue Meeting <input type="checkbox"/> Day 100 Meeting <input type="checkbox"/> Agreement Meeting <input type="checkbox"/> Determination Meeting <input type="checkbox"/> Study Risk Determination <input type="checkbox"/> Other (specify):
<p><b>IDE</b></p> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement	<p><b>Humanitarian Device Exemption (HDE)</b></p> <input type="checkbox"/> Original Submission <input type="checkbox"/> Amendment <input type="checkbox"/> Supplement <input type="checkbox"/> Report <input type="checkbox"/> Report Amendment	<p><b>Class II Exemption Petition</b></p> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<p><b>Evaluation of Automatic Class III Designation (De Novo)</b></p> <input type="checkbox"/> Original Submission <input type="checkbox"/> Additional Information	<p><b>Other Submission</b></p> <input type="checkbox"/> 513(g) <input type="checkbox"/> Other (describe submission):

Have you used or cited Standards in your submission?  Yes  No (If Yes, please complete Section I, Page 5)

**SECTION B**

**SUBMITTER, APPLICANT OR SPONSOR**

Company / Institution Name Rotation Medical, Inc.	Establishment Registration Number (if known) 3009351468		
Division Name (if applicable)	Phone Number (including area code) 763.746.7502		
Street Address 15350 25th Avenue N. Suite 100	FAX Number (including area code) 763.746.7501		
City Plymouth	State / Province MN	ZIP/Postal Code 55447	Country USA
Contact Name Jeff Sims			
Contact Title Vice President, Clinical Programs and Regulatory Affairs		Contact E-mail Address jsims@rotationmedical.com	

**SECTION C**

**APPLICATION CORRESPONDENT (e.g., consultant, if different from above)**

Company / Institution Name			
Division Name (if applicable)	Phone Number (including area code)		
Street Address	FAX Number (including area code)		
City	State / Province	ZIP Code	Country
Contact Name			
Contact Title		Contact E-mail Address	

**SECTION D1**

**REASON FOR APPLICATION - PMA, PDP, OR HDE**

- New Device
- Withdrawal
- Additional or Expanded Indications
- Request for Extension
- Post-approval Study Protocol
- Request for Applicant Hold
- Request for Removal of Applicant Hold
- Request to Remove or Add Manufacturing Site

- Change in design, component, or specification:
  - Software / Hardware
  - Color Additive
  - Material
  - Specifications
  - Other (*specify below*)

- Location change:
  - Manufacturer
  - Sterilizer
  - Packager

- Process change:
  - Manufacturing     Packaging
  - Sterilization
  - Other (*specify below*)

- Labeling change:
  - Indications
  - Instructions
  - Performance Characteristics
  - Shelf Life
  - Trade Name
  - Other (*specify below*)

- Report Submission:
  - Annual or Periodic
  - Post-approval Study
  - Adverse Reaction
  - Device Defect
  - Amendment

- Response to FDA correspondence:

- Change in Ownership
- Change in Correspondent
- Change of Applicant Address

- Other Reason (*specify*):

**SECTION D2**

**REASON FOR APPLICATION - IDE**

- New Device
- New Indication
- Addition of Institution
- Expansion / Extension of Study
- IRB Certification
- Termination of Study
- Withdrawal of Application
- Unanticipated Adverse Effect
- Notification of Emergency Use
- Compassionate Use Request
- Treatment IDE
- Continued Access

- Change in:
  - Correspondent / Applicant
  - Design / Device
  - Informed Consent
  - Manufacturer
  - Manufacturing Process
  - Protocol - Feasibility
  - Protocol - Other
  - Sponsor

- Response to FDA Letter Concerning:
  - Conditional Approval
  - Deemed Approved
  - Deficient Final Report
  - Deficient Progress Report
  - Deficient Investigator Report
  - Disapproval
  - Request Extension of Time to Respond to FDA
  - Request Meeting
  - Request Hearing

- Report submission:
  - Current Investigator
  - Annual Progress Report
  - Site Waiver Report
  - Final

- Other Reason (*specify*):

**SECTION D3**

**REASON FOR SUBMISSION - 510(k)**

- New Device

- Additional or Expanded Indications

- Change in Technology

- Other Reason (*specify*):

**SECTION E ADDITIONAL INFORMATION ON 510(K) SUBMISSIONS**

Product codes of devices to which substantial equivalence is claimed

1	FTM	2		3		4	
5		6		7		8	

Summary of, or statement concerning, safety and effectiveness information

- 510 (k) summary attached  
 510 (k) statement

Information on devices to which substantial equivalence is claimed (if known)

	510(k) Number	Trade or Proprietary or Model Name	Manufacturer
1	K122048	Collagen Tendon Sheet-D	Rotation Medical, Plymouth MN
2			
3			
4			
5			
6			

**SECTION F PRODUCT INFORMATION - APPLICATION TO ALL APPLICATIONS**

Common or usual name or classification name

Surgical Mesh

	Trade or Proprietary or Model Name for This Device	Model Number
1	Collagen Tendon Sheet-DDI (CTS-DDI)	1 TBD
2		2
3		3
4		4
5		5

FDA document numbers of all prior related submissions (regardless of outcome)

1	K131632	2		3		4		5		6	
7		8		9		10		11		12	

Data Included in Submission

- Laboratory Testing       Animal Trials       Human Trials

**SECTION G PRODUCT CLASSIFICATION - APPLICATION TO ALL APPLICATIONS**

Product Code FTM	C.F.R. Section (if applicable) 878.3300	Device Class <input type="checkbox"/> Class I <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Class III <input type="checkbox"/> Unclassified
Classification Panel General and Plastic Surgery		

Indications (from labeling)

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

**Note:** Submission of the information entered in Section H does not affect the need to submit device establishment registration.

FDA Document Number (if known)

**SECTION H MANUFACTURING / PACKAGING / STERILIZATION SITES RELATING TO A SUBMISSION**

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number <b>(b)(4)</b>	<input type="checkbox"/> Manufacturer <input checked="" type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name <b>(b)(4)</b>		Establishment Registration Number <b>(b)(4)</b>	
Division Name (if applicable)		Phone Number (including area code) <b>(b)(4)</b>	
Street Address <b>(b)(4)</b>		FAX Number (including area code) <b>(b)(4)</b>	
City <b>(b)(4)</b>		State / Province <b>(b)</b>	ZIP Code <b>(b)</b>
		Country <b>(b)(4)</b>	
Contact Name <b>(b)(4)</b>	Contact Title <b>(b)(4)</b>	Contact E-mail Address <b>(b)(4)</b>	

<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number <b>(b)(4)</b>	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input checked="" type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name <b>(b)(4)</b>		Establishment Registration Number <b>(b)(4)</b>	
Division Name (if applicable)		Phone Number (including area code) <b>(b)(4)</b>	
Street Address <b>(b)(4)</b>		FAX Number (including area code) <b>(b)(4)</b>	
City <b>(b)(4)</b> Trade		State / Province <b>(b)</b>	ZIP Code <b>(b)(4)</b>
		Country <b>(b)(4)</b>	
Contact Name <b>(b)(4)</b>	Contact Title <b>(b)(4)</b>	Contact E-mail Address <b>(b)(4)</b>	

<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete	Facility Establishment Identifier (FEI) Number	<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract Manufacturer	<input type="checkbox"/> Contract Sterilizer <input type="checkbox"/> Repackager / Relabeler
Company / Institution Name		Establishment Registration Number	
Division Name (if applicable)		Phone Number (including area code)	
Street Address		FAX Number (including area code)	
City		State / Province	ZIP Code
		Country	
Contact Name	Contact Title	Contact E-mail Address	

## SECTION I

## UTILIZATION OF STANDARDS

**Note:** Complete this section if your application or submission cites standards or includes a "Declaration of Conformity to a Recognized Standard" statement.

	Standards No.	Standards Organization	Standards Title	Version	Date
1	11135-1	ISO	Sterilization of health care products - Ethylene oxide	1st version	05/01/2007
2	10993	ISO	Biological evaluation of medical devices	3rd edition	08/01/2003
3					
4					
5					
6					
7					

**Please include any additional standards to be cited on a separate page.**

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF ADDRESS BELOW.\***

The burden time for this collection of information is estimated to average 0.5 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
 Food and Drug Administration  
 Office of Chief Information Officer  
 Paperwork Reduction Act (PRA) Staff  
 1350 Piccard Drive, Room 400  
 Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*



February 3, 2014

VIA FEDERAL EXPRESS

Document Mail Center (WO66-G609)  
Center for Devices and Radiological Health  
Food and Drug Administration  
10903 New Hampshire Avenue  
Silver Spring, Maryland 20993-0002

Re: 510(k) Premarket Notification (Traditional)  
Collagen Tendon Sheet-DDI (b) (4)

Note: This submission includes an electronic copy of the 510(k) submission per *eCopy Program for Medical Device Submissions*. The eCopy is an exact duplicate of the paper copy.

Dear Sir or Madam:

Pursuant to 21 CFR Part 807, Subpart E, Premarket Notification Procedures, Section 807.81, Rotation Medical, Inc. is submitting this 510(k) Premarket Notification for its Collagen Tendon Sheet-DDI. Rotation Medical has determined that Collagen Tendon Sheet-DDI is substantially equivalent to current legally marketed tendon protector devices and intends to manufacture and market the device.

(b) (4)

<b>Trade Name or Proprietary Name:</b>	Collagen Tendon Sheet-DDI
<b>Common or Usual Name:</b>	Tendon Protector
<b>Device Classification Name:</b>	Mesh, Surgical
<b>Regulation Number:</b>	878.3300
<b>Product Code:</b>	FTM
<b>Device Class:</b>	Class II
<b>Name and Address of Manufacturer:</b>	Rotation Medical, Inc. 15350 25 <sup>th</sup> Avenue North, Suite 100 Plymouth, MN 55447



**Establishment Registration No.:** 3009351468

**Name, Address, and Telephone Number of Contact Person:** Jeff Sims  
 Vice President, Clinical Programs and Regulatory Affairs  
 Rotation Medical, Inc.  
 15350 25<sup>th</sup> Avenue North, Suite 100  
 Plymouth, MN 55447  
 Tel: 763.746.7502  
 Fax: 763.746.7501  
 E-mail: jsims@rotationmedical.com

Below is the recommended “Design and Use of the Device” questions per the *Guidance for Industry and FDA Staff: Format for Traditional and Abbreviated 510(k)s*.

Question	Yes	No
Is the device intended for prescription use (21 CFR 801 Subpart D)?	X	
Is the device intended for over-the-counter use (21 CFR 807 Subpart C)?		X
Is the device provided sterile?	X	
Is the device intended for single use?	X	
Is the device a reprocessed single use device?		X
If yes, does this device type require reprocessed validation data?		n/a
Does the device contain a drug?		X
Does the device contain a biologic?		X
Does the device use software?		X
Does the submission include clinical information?		X
Is the device implanted?	X	

The purpose of this application is to gain clearance for Collagen Tendon Sheet-DDI, an additional model of the predicate device, Rotation Medical’s Collagen Tendon Sheet-D (K122048). Collagen Tendon Sheet-DDI is exactly the same device as its predicate except it is

(b) (4)

The indication for use for Collagen Tendon Sheet-DDI is the same as the predicate Collagen Tendon Sheet-D:

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.



Arthroscopic tendon repair procedures constitute a growing percentage of all such procedures.

(b) (4)

If you have any questions, please do not hesitate to contact me. Thank you all in advance for your work toward the clearance of this device.

Sincerely,

A handwritten signature in blue ink, appearing to read "Jeff Sims", with a stylized flourish at the end.

Jeff Sims  
Vice President, Clinical Programs and Regulatory Affairs

Enclosures (submitted including an electronic copy)

## Indications for Use

510(k) Number (if known): \_\_\_\_\_

Device Name: Collagen Tendon Sheet-DDI

### Indications for Use:

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

Prescription Use X  
(Part 21 CFR 801 Subpart D)

AND/OR

Over-The-Counter Use \_\_\_\_\_  
(21 CFR 807 Subpart C)

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)



## 510(k) Summary

### Applicant Information

<b>Applicant Name:</b>	Rotation Medical, Inc.
<b>Applicant Address:</b>	15350 25 <sup>th</sup> Avenue North, Suite 100 Plymouth, MN 55447
<b>Telephone:</b>	763-746-7502
<b>Fax:</b>	763-746-7501
<b>Contact Person:</b>	Jeff Sims Vice President, Clinical Programs and Regulatory Affairs
<b>Date Prepared:</b>	February 3, 2014

### Name of Device

<b>Device Common Name:</b>	Tendon Protector
<b>Device Trade Name:</b>	Collagen Tendon Sheet-DDI
<b>Device Classification Name:</b>	Mesh, Surgical 878.3300 Class II FTM

### Legally Marketed Devices to Which Substantial Equivalence is Claimed

<b>Predicate Device(s):</b>	Collagen Tendon Sheet-D, K122048 Rotation Medical, Inc.
-----------------------------	--

### Description of the Device

Collagen Tendon Sheet-DDI is a resorbable type I collagen matrix that provides a layer of collagen over injured tendons. Collagen Tendon Sheet-DDI is designed to provide a layer between the tendon and the surrounding tissue during healing. When hydrated, Collagen Tendon Sheet-DDI is an easy-to-use, soft, pliable, nonfriable, porous collagen sheet. Collagen Tendon Sheet-DDI is provided sterile, non-pyrogenic, for single use only, in a variety of sizes, and is packaged preloaded in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration.



### **Intended Use**

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### **Summary/Comparison of Technical Characteristics**

Collagen Tendon Sheet-DDI is exactly the same product as its predicate, Collagen Tendon Sheet-D (K122048) in regard to all aspects of the collagen scaffold. Only the packaging has changed between the current device and the predicate device. Collagen Tendon Sheet-DDI is packaged in a cartridge to facilitate arthroscopic delivery and positioning of the collagen matrix.

No new device performance or characterization data is included in this submission. Past safety, biocompatibility and mechanical characterization tests of the predicate product are directly applicable to the current product. Additional tests to confirm the performance characteristics of the packaging change were performed.


In summary, the Collagen Tendon Sheet-DDI device passed all applicable FDA Blue Book Memorandum G95-1 and ISO 10993-1 testing for the biological evaluation of medical devices. Testing was conducted in accordance with FDA's Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh.


### **Conclusion**

The purpose of this 510(k) application was to notify the Food and Drug Administration of proposed modifications to the packaging of the previously cleared Collagen Tendon Sheet-D (K122048). Collagen Tendon Sheet-DDI and Collagen Tendon Sheet-D are exactly the same device excepting their respective packaging. The proposed packaging modifications for the new device raise no new questions regarding safety and effectiveness. Therefore, the proposed Collagen Tendon Sheet-DDI is substantially equivalent to the predicate Collagen Tendon Sheet-D.

## 6. TRUTHFUL AND ACCURACY STATEMENT

I certify that, in my capacity as Vice President, Clinical Programs and Regulatory Affairs for Rotation Medical, Inc., I believe to the best of my knowledge, that all data and information submitted in the Premarket Notification are truthful and accurate and that no material fact has been omitted.

  
\_\_\_\_\_  
Jeff Sims,  
Vice President, Clinical Programs and Regulatory Affairs  
Rotation Medical, Inc.

  
\_\_\_\_\_  
Date

## **7. Class III Summary and Certification**

Class III Summary and Certification is not applicable to this device.

## **8. Financial Certification or Disclosure Statement**

Since clinical studies were not performed in support of this application, no financial certifications or disclosure statements are applicable.

## **9. Declarations of Conformity and Summary Reports**

### **9.1 Sterilization**

Sterilization is conducted in accordance with ISO 11135: *Sterilization of Health Care Products – Ethylene Oxide, 2007*

### **9.2 Biocompatibility**

The biocompatibility of the finished product was tested according to ISO 10993-1: *Biological Evaluation of Medical Devices – Part 1 Evaluation and Testing, 2003*

### **9.3 FDA Guidance Documents**

FDA Guidance Document entitled “Guidance for the Preparation of a Premarket Notification Application for a Surgical Mesh,” issued on March 2, 1999, was used in the development and testing of the product.

Standards Data Reports (Form FDA 3654) for the above referenced sterilization and biocompatibility standards are included in Appendix F.


## 10. EXECUTIVE SUMMARY

### 10.1 Brief Description of Device

Collagen Tendon Sheet-DDI is exactly the same product as its predicate; each are a bioabsorbable implant device that provides a layer of collagen over injured tendons. Collagen Tendon Sheet-DDI is designed to provide a layer of collagen between a flat tendon and the surrounding tissue during tendon healing. After hydration, Collagen Tendon Sheet-DDI is an easy-to-handle, soft, pliable, nonfriable, porous collagen sheet. Collagen Tendon Sheet-DDI is provided sterile, non-pyrogenic, for single use only, and is packaged in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration.

### 10.2 Design Philosophy

(b) (4)




In addition to the direct predicate, Rotation Medical's Collagen Tendon Sheet (K112423), and

(b) (4)



protective environment for tissue healing. Other tendon protector products successfully used for the protection of tendon injuries include (b) (4)



(b) (4)



Arthroscopic tendon repair procedures constitute a growing percentage of all such procedures.

Collagen Tendon Sheet-DDI (b) (4)



### 10.3 Device Comparison

Rotation Medical's Collagen Tendon Sheet-DDI is substantially equivalent to Collagen Tendon Sheet-D (K122048). Collagen Tendon Sheet-DDI is exactly the same device as its predicate, Collagen Tendon Sheet-D (K122048), in regard to all aspects of the collagen scaffold. Only the packaging has changed between the current device and the predicate device. Collagen Tendon Sheet-DDI is packaged preloaded in a cartridge to facilitate arthroscopic delivery and positioning of the collagen matrix.

#### **10.4 Summary of Performance Testing**

The results of the safety and biocompatibility testing that were previously submitted in support of the predicate device, Collage Tendon Sheet-D, are directly applicable to Collagen Tendon Sheet-DDI because the two devices are identical. In like manner, the results of the material and mechanical characterization of the predicate device are also directly applicable to the current device.

The materials used in the delivery instrument cartridge head assembly have passed additional biocompatibility and simulated-use and integrity safety tests (Section 15 and 20 respectively). The proposed packaging modifications to the predicate device raise no new questions regarding safety and effectiveness. Therefore, the proposed Collagen Tendon Sheet-DDI is substantially equivalent to the predicate Collagen Tendon Sheet-D.

## 11. DEVICE DESCRIPTION

### 11.1 Description


Collagen Tendon Sheet-DDI, like its predicate, is a resorbable type I collagen matrix that provides a layer of collagen over injured tendons. Collagen Tendon Sheet-DDI is designed to provide a layer of collagen between a flat tendon and the surrounding tissue during tendon healing. After hydration, Collagen Tendon Sheet-DDI is an easy-to-use, soft, pliable, nonfriable, porous collagen sheet. Collagen Tendon Sheet-DDI is provided sterile, non-pyrogenic, for single use only, in a variety of sizes, and is packaged in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration.

### 11.2 Intended Use

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### 11.3 Contract Development and Manufacture

(b) (4)




### 11.4 Collagen Materials

All of the materials in Collagen Tendon Sheet-DDI are identical to the predicate Collagen Tendon Sheet-D. As often stated within this current submission, Collagen Tendon Sheet-DDI and the predicate Collagen Tendon Sheet-D are the exact same device excepting their respective packaging. (b) (4)




### 11.5 Performance Specifications/Design Requirements

The design specifications of Collagen Tendon Sheet-DDI are the same as the predicate Collagen Tendon Sheet-D. Only the packaging has changed between the two products. Therefore the associated shelf life of the device has also changed; (b) (4)



The key design parameters and final product specifications of the Collagen Tendon Sheet-DDI are summarized in Table 11-1, (b) (4)



. The design has been verified and validated through the design control process.

**Table 11-1 Summary of Product Specifications**

Parameter	Product Specification
<b>Performance</b>	Management and protection of tendon injuries in which there has been no substantial loss of tendon tissue
<b>Chemical and Physical Properties</b>	
(b) (4)	[REDACTED]
[REDACTED]	[REDACTED] ECH: ≤ 9 mg
<b>Biological Properties</b>	
Biocompatibility	Biocompatible (Pass FDA G95-1 and ISO 10993)
Pyrogenicity <sup>†</sup>	Non-pyrogenic (≤ 0.5 EU/ml)
<b>Stability</b>	
(b) (4)	[REDACTED]
<b>Packaging</b>	
(b) (4)	(b) (4)


(b) (4)

## **11.6 Manufacturing Process**

### **11.6.1 Manufacturing Flowchart**

The manufacturing flowchart is shown on the following page. There has been no change in the flowchart as a result of the packaging change between the processing of Collagen Tendon Sheet-DDI and the predicate Collagen Tendon Sheet-D. The packaging change between the two products is described in Section 11.6.2 below.

**Manufacturing Flowchart for Collagen Tendon Sheet-DDI**



### 11.6.2 Description of the Manufacturing Process

The only change between the processing of Collagen Tendon Sheet-DDI and the predicate Collagen Tendon Sheet-D is the (b) (4)


All other aspects of the predicate device manufacturing are applicable to the current device manufacturing and will not be repeated here.

However, since the packaging of the device has changed, this change will be elaborated upon.

(b) (4)

**Table 11-2 Description of the manufacturing process** (b) (4)

(b) (4)




(b) (4)



Note, the actual loading of the collagen scaffold in the cartridge assembly is done in a clean room environment with gloved hands.

(b) (4)



## **12. SUBSTANTIAL EQUIVALENCE DISCUSSION**

### **12.1 Predicate Devices**

The data presented in this 510(k) demonstrate that Collagen Tendon Sheet-DDI is substantially equivalent to the following predicate devices:

Collagen Tendon Sheet-D (K122048)  
Rotation Medical, Inc. Plymouth, MN

### **12.2 Substantial Equivalence Comparison Table**

The Substantial Equivalence Comparison Table is presented on the following page (Table 12-1). As the table shows, the technological characteristics of the Collagen Tendon Sheet-DDI are identical and therefore are substantially equivalent to the predicate device referenced above. Indeed they are the same device, only the packaging has changed between the two products as reflected in the last line of Table 12-1.


Shelf life assessments of the Collagen Tendon Sheet-DDI device in its new packaging are sufficient to demonstrate that the new device is not impacted by the change in packaging. These assessments and the protocol for future assessments are submitted in Appendix D and Section 14 respectively. Collagen Tendon Sheet-DDI meets its design specifications and is equivalent to its predicate device. Further the ongoing lot acceptance tests are performed on product after packaging and sterilization, these tests too will continually verify the packaging change has no effect on the performance characteristics of the device.

**Table 12-1 Substantial Equivalence Comparison Chart with Comparative Data**

	<b>Collagen Tendon Sheet-DDI (This submission)</b>	<b>Collagen Tendon Sheet-D K122048</b>
<b>Indications for use</b>	Management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.	Management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.
<b>Material</b>	(b)(4) Trade Secret Process-Product Specs	
<b>Source</b>		
<b>Form</b>		
<b>Color</b>		
<b>Dimensions</b>		
<b>Thickness</b>		
<b>Density</b>		
<b>Mechanical Strength</b>		
<b>St</b>		
<b>Biocompatibility</b>		
<b>(b)(4) Trade Secret</b>	(b)(4) Trade Secret P	
<b>(b)(4) Trade Secret</b>		
<b>Sterility</b>	Sterile, SAL 10 <sup>-6</sup> ETO sterilization	Sterile, SAL 10 <sup>-6</sup> ETO sterilization
<b>Pyrogenicity</b>	Non-pyrogenic (≤ 0.5 EU/ml)	Non-pyrogenic (≤ 0.5 EU/ml)
<b>Single use / Reuse</b>	(b)(4) Trade Secret Process-Product Specs	
<b>Packaging</b>		

### 12.3 Substantial Equivalence Decision Flowchart

The following Substantial Equivalence Decision-making Flowchart was used to assess Collagen Tendon Sheet-DDI against the predicate Collagen Tendon Sheet-D.



## 12.4 Description of Substantial Equivalence Flowchart Decisions

### **Does the new device have the same indication statement?**

Yes. Collagen Tendon Sheet-DDI and Collagen Tendon Sheet-D have the same indication statement, which is: the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### **New device has same intended use and may be substantially equivalent.**

Yes. Collagen Tendon Sheet-DDI and Collagen Tendon Sheet-D have the same intended use, which is to provide a protective environment for tissue repair and healing.

### **Does the new device have the same technological characteristics, e.g., design, materials, etc.?**

Yes. There is no difference in the technological, material or processing characteristics of the implanted devices. Indeed they are the same device. Only the packaging has changed between the two products. The only difference is that Collagen Tendon Sheet-DDI is (b) (4)

### **Could the new characteristics affect safety and effectiveness?**

Yes. Although the packaging change was not likely to affect either the safety or performance of the device, a complete set of shelf life performance tests was undertaken on three lots of Collagen Tendon Sheet-DDI to demonstrate such was indeed the case. See Sections 12.5 below for more information on the packaging change and Section 14 for more information on the associated shelf life testing protocol.

### **Substantial equivalence**

The results of the data presented demonstrate that Collagen Tendon Sheet-DDI is substantially equivalent to its predicate Collagen Tendon Sheet-D.

## 12.5 Detailed Description of Substantial Equivalence

The overall product comparisons were presented in Table 12-1 above. In this section, a detailed comparison of the key material characterization test results between the subject device and the predicate device will be discussed:

- Collagen Tendon Sheet-DDI – Subject Device
- Collagen Tendon Sheet-D (K122048) – Predicate Device

### 12.5.1 Intended Use

The intended use of the Collagen Tendon Sheet-DDI and its predicate is identical, which is: the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### 12.5.2 Form, Appearance, and Dimensions

Collagen Tendon Sheet-DDI and its predicate device, Collagen Tendon Sheet-D, are exactly the same device prior to Collagen Tendon Sheet-DDI being loaded into the delivery instrument cartridge for packaging. The manufacturing process for preloading the device in the cartridge is described in Section 11.6.2. (b) (4)






Figure 12-1 Collagen Tendon Sheet-DDI

Figure 12-2 Collagen Tendon Sheet-D





The size and dimensions of the current product and its predicate are identical.

	Collagen Tendon Sheet-DDI	Collagen Tendon Sheet-D
Form	(b)(4) Trade Secret Process-Product Specs	
Color		
Dimensions		
Thickness		

The Collagen Tendon Sheet-DDI is provided in the same sizes and thickness as the predicate Collagen Tendon Sheet-D. Dimensions and thickness are verified on all devices in manufacturing prior to final packaging.

### 12.5.3 Material Composition

The Collagen Tendon Sheet-DDI and its predicate Collagen Tendon Sheet-D are both manufactured from the identical material, (b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-Product Specs

	Collagen Tendon Sheet-DDI	Collagen Tendon Sheet-D
	(b)(4) Trade Secret Process-Product Specs	(b)(4) Trade Secret Process-Product Specs

### 12.5.4 Physical Properties

(b)(4) Trade Secret

Process-Product Specs

(b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-

(b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-Product Specs

### 12.5.5 Physico-chemical Properties

(b)(4) Trade Secret Process-Product

(b)(4) Trade Secret Process-Product Specs

	Collagen Tendon Sheet-DDI	Collagen Tendon Sheet-D
(b)(4) Trade Secret Process-Product Specs	(b)(4) Trade Secret Process-	(b)(4) Trade Secret

(b)(4) Trade Secret Process-Product Specs

### 12.5.6 Sterility

The subject device and its predicates are all terminally sterilized using Ethylene oxide at a sterility assurance level of  $10^{-6}$ .

	Collagen Tendon Sheet-DDI	Collagen Tendon Sheet-D
<b>Sterilization Method</b>	Ethylene oxide	Ethylene oxide

### 12.5.7 (b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-Product Specs  
(b)(4) Trade Secret Process-Product Specs

(b)(4) Trade Secret Process-Product Specs

the Collagen Tendon Sheet-D, and therefore Collagen Tendon Sheet-DDI, [REDACTED]

[REDACTED]

	Collagen Tendon Sheet-D	Collagen Tendon Sheet
(b)(4) Trade Secret Process-Product Specs	[REDACTED]	[REDACTED]

(b)(4) Trade Secret Process-Product Specs [REDACTED]

(b)(4) Trade Secret Process-Product Specs [REDACTED]

## 12.6 Conclusions of Comparative Analysis

Based on the comparative analysis, Collagen Tendon Sheet-DDI has been shown to be substantially equivalent to its predicate, Collagen Tendon Sheet-D; indeed they are exactly the same device, only the packaging has changed. Additional support is provided throughout this submission.

### 13. PROPOSED LABELING

The product label and instructions for use of the Collagen Tendon Sheet-DDI are provided in Appendix C. Because of the similarities in intended use, materials, and product characteristics, the label and instructions for use of the predicate device Collagen Tendon Sheet-D were used as a guide in developing the insert for Collagen Tendon Sheet-DDI.

(b)(4) Trade Secret Process-Product Specs



## 14. STERILIZATION AND SHELF LIFE

### 14.1 Sterilization

The method of sterilization is ethylene oxide sterilization for both Collagen Tendon Sheet-DDI and its predicate, Collagen Tendon Sheet-D. The sterilization method has been validated in accordance with ISO 11135-1, Sterilization of health care products - Ethylene oxide. The sterility assurance level of the device is  $10^{-6}$ .

Pursuant to ISO requirements, the ethylene oxide and ethylene chlorohydrins residual must meet the following acceptance criteria:

Ethylene oxide:  $\leq 4$  mg  
Ethylene chlorohydrin  $\leq 9$  mg

The final packaging used to maintain sterility is a dual Tyvek sterile barrier tray-in-tray design, typical of sterile medical device packaging.

### 14.2 Shelf Life

Product stability assessments examine the stability of the product over time as packaged and stored under accelerated or real time aging conditions. (b)(4) Trade Secret Process-Product Specs. The results of the tests at each time point must meet overall product specifications.

In order to assess shelf life in real time, Rotation Medical (b)(4) Trade Secret Process-Product Specs. The protocol when (b)(4) Trade Secret Process-Product Specs. Table 14-1 summarizes the on-going shelf life test protocol for the device.

(b)(4) Trade Secret Process-Product Specs.

Table 14-1 Self Life Testing Protocol

(b)(4) Trade Secret Process-Product Specs	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED] [REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED] [REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED] [REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]

(b)(4) Trade Secret Process-Product Specs

Table 14-2 Package Integrity Protocol Testing

(b)(4) Trade Secret Process-Product Specs	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]
[REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]

(b)(4) Trade Secret Process-Product Specs

[REDACTED] The test results to date are provided in Appendix D.

## 15. BIOCOMPATIBILITY

The extensive biocompatibility testing done in support of the predicate product, Collagen Tendon Sheet-D, is directly applicable to Collagen Tendon Sheet-DDI and was not repeated for the current application.

The following table summarizes the additional biocompatibility tests and assessments made for the cartridge head assembly. The associated final reports are provided in Appendix B. The additional biocompatibility tests were performed (b)(4) Trade Secret Process-Product Specs to established protocols in conformance with ISO 10993-1, 2003 for parts that contact the tissue and bone less than 24 hours. (b)(4) Trade Secret Process-Product Specs

Table 15-1. (b)(4) Trade Secret Process-Product Specs

Test	Test Method/Model or Assessment	Results
(b)(4) Trade Secret Process-Product Specs		

### Conclusion

Collagen Tendon Sheet-DDI is exactly the same device as its predicate Collagen Tendon Sheet-D, excepting their respective packaging, and like its predicate, the current device is biocompatible and safe for human implantation. The addition of the acute cartridge components have also been shown to be appropriately biocompatible.

## **16. SOFTWARE**

The device does not contain software; therefore this section is not applicable to this submission.

## **17. ELECTROMAGNETIC COMPATIBILITY AND ELECTRICAL SAFETY**

The device does not include an electronic component; therefore this section is not applicable to this submission.



Parameter	Product Specification / Design Input	Design Output
(b)(4) Trade Secret		
Residual EO and ECH	EO: ≤ 4 mg ECH: ≤ 9 mg	(b)(4) Trade Secret
<b>Biological Properties</b>		
Biocompatibility	Biocompatible (Pass FDA G95-1 and ISO 10993)	(b)(4) Trade Secret Process-Product Specs
Pyrogenicity	Non-pyrogenic (≤ 0.5 EU/ml)	(b)(4) Trade Secret
<b>Stability</b>		
Shelf Life	(b)(4) Trade Secret Process-Product Specs	
<b>Packaging</b>		
Double peel packages	(b)(4) Trade Secret Process-Product Specs	

## 18.2 Design Verification Test Methods and Results

Collagen Tendon Sheet-DDI met all of the product specifications and design requirements set out for the product.

The design verification data sheets that support the summary results above for Collagen Tendon Sheet-DDI are included in Appendix E. (b)(4) Trade Secret Process-Product Specs

[Redacted]

## 19. PERFORMANCE TESTING - ANIMAL

There was no additional animal performance testing done in support of Collagen Tendon Sheet-DDI. The animal testing completed in support of the root predicate product, Collagen Tendon Sheet, is applicable to Collagen Tendon Sheet-DDI since the product specifications and design requirements of both products are the same, (b)(4) Trade Secret Process-Product Specs and the new packaging for Collagen Tendon Sheet-DDI.

## 20. PERFORMANCE TESTING - CLINICAL

### 20.1 Clinical Experience with Collagen Tendon Sheet

Rotation Medical Inc. has not yet made its Collagen Tendon Sheet products commercially available, (b)(4) Trade Secret Process-Product Specs



### 20.2 Simulated-Use and Clinical Performance Testing

Rotation Medical has completed a series of simulated use and integrity tests demonstrating the (b)(4) Trade Secret Process-Product Specs Collagen Tendon Sheet-DDI. (b)(4) Trade Secret Process-Product Specs



In the first test (Report Number 2324 following Protocol Number 2304 included in Appendix G) three areas were assessed, namely: (b)(4) Trade Secret Process-Product Specs



(b)(4) Trade Secret Process-Product Specs



(b)(4) Trade Secret Process-Product Specs



(b)(4) Trade Secret Process-Product Specs



(b)(4) Trade Secret Process- Testing






Figure 20-1

(b)(4) Trade Secret Process- Testing



(b)(4) Trade Secret Process- Testing



(b)(4) Trade Secret Process- Testing






Table 20-2

(b)(4) Trade Secret Process- Testing




### 20.3 Conclusions from the Simulated-Use and Clinical Performance Testing

Rotation Medical (b)(4) Trade Secret Process- Testing



Additionally, the functional integrity and safety of the cartridge (b)(4) Trade Secret Process- Testing



Appendix A – (b)(4) Trade Secret Process- Testing


(b)(4) Trade Secret Process- Testing

(b)(4) Trade Secret Process- Testing


Figure A1. Rotation Medical Arthroscopic  
Delivery Instrumentation for Collagen Tendon Sheet-DDI

(b)(4) Trade Secret Process- Testing


(b)(4) Trade Secret Process- Testing




(b)(4) Trade Secret Process- Testing



(b)(4) Trade Secret Process- Testing



(b)(4) Trade Secret Process- Testing



(b)(4) Trade Secret Process- Testing

Instrument	Nominal Length	Nominal Width
(b)(4) Trade Secret Process- Testing		

**Table 1 Dimensions of Acute Surgical Instrumentation**

(b)(4) Trade Secret Process- Testing





























































































































































# Collagen Tendon Sheet-DDI

REF XXXXXX

Quantity: 1 Membrane

Size: XXXXXXXXX

LOT XXXXXXXXXXX

 YYYY-MM  
EXPIRATION DATE

 DO NOT RE-USE

**STERILEEO** Non-pyrogenic

**R<sub>x</sub>** only

 ATTENTION.  
SEE INSTRUCTIONS FOR USE

(b)(4) Trade Secret  
Process- Product  
Rotation Medical, Inc., 15350 25<sup>th</sup> Avenue North, Suite 100,  
Plymouth, MN 55447 USA

(b)(4) Trade Secret Process- Product Specs

CTS-DDI label draft

## Collagen Tendon Sheet-DDI (CTS-DDI)

### Description

Collagen Tendon Sheet-DDI (CTS-DDI) is a bioabsorbable implant device that provides a layer of collagen over injured tendons. CTS-DDI is designed to provide a layer of collagen between a flat tendon and the surrounding tissue. After hydration, CTS-DDI is an easy-to-handle, soft, pliable, nonfriable, porous collagen sheet. CTS-DDI is provided sterile, non-pyrogenic, for single use only, in a variety of sizes, and is packaged pre-loaded in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration.

### Indications for Use

Collagen Tendon Sheet-DDI is indicated for the management and protection of tendon injuries in which there has been no substantial loss of tendon tissue.

### Contraindications

Collagen Tendon Sheet-DDI is not designed, sold, or intended for use except as described in the indications for use and is contraindicated in the following situations:

- Collagen Tendon Sheet-DDI is not indicated to replace damaged tendon or to reinforce the strength of any tendon repair.
- Collagen Tendon Sheet-DDI is not indicated for patients with a known history of hypersensitivity to bovine-derived materials.

### Instructions for Use

1. Follow standard procedures for treatment of the injured tendon.
2. Determine the tendon width in millimeters (mm) using a suitable measuring instrument.
3. Select a CTS-DDI size that is slightly smaller than the width of the tendon.
4. Pre-hydrate CTS-DDI, in the cartridge, in sterile saline for at least 2 minutes.
5. After hydration, use the Rotation Medical Delivery Instrument and accessories to assist in positioning the scaffold over the tendon with one end overlapping the tendon insertion. Reference Delivery Instrument IFU.
6. Secure the CTS-DDI to the tendon and bone. Ensure that CTS-DDI is in good contact with the tendon.

7. Thoroughly irrigate the surgical site and close the incision in the standard fashion.
8. Application of the Collagen Tendon Sheet-DDI does not modify the postoperative treatment. The surgeon must determine motion and strength requirements according to standard practice.

### Warnings

- Do not reuse or re-sterilize.
- Do not use if the product package is damaged or opened.

### Precautions

- Collagen Tendon Sheet-DDI should not be applied until bleeding and infection are controlled.

### Storage

Store at room temperature. Avoid excessive heat or humidity.

### How Supplied

Collagen Tendon Sheet-DDI is supplied sterile for single-use, and is packaged pre-loaded in a cartridge, for use with the Rotation Medical Delivery Instrument, in a dual sterile seal tray-in-tray configuration. Contents of the package are guaranteed sterile and non-pyrogenic unless the package is opened or damaged. The CTS-DDI product and packaging do not contain natural rubber latex.

### Caution

Federal (USA) law restricts this device to sale by or on the order of a physician.

### Symbols Used on Labeling



See Instructions for Use



Expiration Date



Do not reuse after opening



Lot Number



Method of sterilization – ethylene oxide

(b)(4) Trade Secret Process-

Rotation Medical, Inc.

15350 25<sup>th</sup> Ave. N., Plymouth, MN 55447 USA

(b)

(4)

TT

## Rotation Medical Delivery Instrument

### Description

The Rotation Medical Delivery Instrument consists of the delivery instrument and a “clamshell” which is used to facilitate attachment of the cartridge containing the Rotation Medical Scaffold onto the end of the delivery instrument. The Rotation Medical Delivery Instrument is provided sterile, non-pyrogenic, for single-use only in a dual sterile seal tray-in-tray configuration.

### Indications for Use

The Scaffold Delivery Instrument is indicated for arthroscopic delivery of the Rotation Medical Scaffold into the subacromial space.

### Contraindications

The Scaffold Delivery Instrument is not designed, sold, or intended for use except as described in the indications for use and is contraindicated in the following situation:

- The Scaffold Delivery Instrument is not indicated for use with scaffolds manufactured by any company other than Rotation Medical.

### Instructions for Use

1. Place the cartridge containing the Rotation Medical Scaffold into the clamshell with the slot in the introducer tube engaging with the nub on the blue side of the clamshell and snap the clamshell closed.
2. Submerge the clamshell containing the cartridge in sterile saline for at least 2 minutes to hydrate the scaffold.
3. After hydration, snap the clamshell onto the end of the delivery instrument with the blue side of the clamshell facing down.
4. Pinch the wings of the clamshell and pull the clamshell, including the introducer tube, off of the delivery instrument.
5. Insert the guidewire into the small hole in the shaft of the delivery instrument and advance the cartridge down the guidewire until the cartridge has passed through the arthroscopic portal into the subacromial space. The cartridge is fully advanced to the desired location when the red button on the back of the delivery instrument is pushed out.
6. Press the trigger release button on the side of the delivery instrument and squeeze the trigger to

retract the sheath, which allows the scaffold to flatten out over the tendon.

7. While holding the handle of the delivery instrument, align the scaffold with the direction of the tendon and maintain that position and alignment while the scaffold is attached to the tendon.
8. After the scaffold has been securely attached to the tendon, pull back on the delivery instrument to remove it from the arthroscopic portal. The Nitinol fingers will easily bend to allow removal.

### Warnings

- Do not reuse or re-sterilize.
- Do not use if the product package is damaged or opened.

### Precautions

- The cartridge should not be delivered into the subacromial space unless it is securely snapped onto the end of the delivery instrument.

### Storage

Store at room temperature.

### How Supplied

The Scaffold Delivery Instrument is supplied sterile for single-use in a dual sterile seal tray-in-tray configuration. Contents of the package are guaranteed sterile and non-pyrogenic unless the package is opened or damaged. The Scaffold Delivery Instrument and packaging do not contain natural rubber latex.

### Caution

Federal (USA) law restricts this device to sale by or on the order of a physician.

### Symbols Used on Labeling



See Instructions for Use



Expiration Date



Do not reuse after opening



Lot Number



Method of sterilization – ethylene oxide

### Manufactured by:

Rotation Medical, Inc.  
15350 25<sup>th</sup> Ave. N., Plymouth, MN 55447 USA



























































































































Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 11135-1, Sterilization of health care products - Ethylene oxide - Part 1, 2007

**Please answer the following questions**

Yes    No

Is this standard recognized by FDA <sup>2</sup>? .....    

FDA Recognition number <sup>3</sup> ..... # 14-331

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....    

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....       
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....    

Does this standard include acceptance criteria? .....       
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....       
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....       
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....    

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....       
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....       
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....       
If yes, was the guidance document followed in preparation of this 510(k)? .....    

Title of guidance: \_\_\_\_\_

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], [www.fda.gov/cdrh/stdsprog.html](http://www.fda.gov/cdrh/stdsprog.html)

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and address of the test laboratory or

certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at [www.fda.gov/cdrh/guidance.html](http://www.fda.gov/cdrh/guidance.html)

**EXTENT OF STANDARD CONFORMANCE  
SUMMARY REPORT TABLE**

STANDARD TITLE  
ISO 11135-1:2007 Sterilization of health care products - Ethylene oxide - Part 1

**CONFORMANCE WITH STANDARD SECTIONS\***

SECTION NUMBER Annex B	SECTION TITLE Conservative determination of lethal rate of the sterilization process, Overkill	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
---------------------------	---	--

TYPE OF DEVIATION OR OPTION SELECTED \*  
Annex B is the option selected.

DESCRIPTION  
PCDs and biological indicators were used in fractional cycle analyses; sterilization for the device was then subjected to a full cycle at least double the half cycle.

JUSTIFICATION  
The methods of sterilization assurance are recognized by the standard and implemented by the contract sterilizer, (b)(4) Trade Secret, (b)(4) using established procedures in their facilities, (b)(4) Trade Secret, certified as a Contract Sterilization Service using ISO 11135-1.

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED \*

DESCRIPTION

JUSTIFICATION

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED \*

DESCRIPTION

JUSTIFICATION

\* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.

\* Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.

**Paperwork Reduction Act Statement**

Public reporting burden for this collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
1350 Piccard Drive, Room 400  
Rockville, MD 20850

*An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number.*

Department of Health and Human Services  
Food and Drug Administration  
**STANDARDS DATA REPORT FOR 510(k)s**  
(To be filled in by applicant)

This report and the Summary Report Table are to be completed by the applicant when submitting a 510(k) that references a national or international standard. A separate report is required for each standard referenced in the 510(k).

TYPE OF 510(K) SUBMISSION

Traditional       Special       Abbreviated

STANDARD TITLE <sup>1</sup>

ISO 10993-1:2003, Biologic Evaluation of medical devices - Part 1

**Please answer the following questions**

Yes    No

Is this standard recognized by FDA <sup>2</sup>? .....    

FDA Recognition number <sup>3</sup> ..... #2-117

Was a third party laboratory responsible for testing conformity of the device to this standard identified in the 510(k)? .....    

Is a summary report <sup>4</sup> describing the extent of conformance of the standard used included in the 510(k)? .....       
If no, complete a summary report table.

Does the test data for this device demonstrate conformity to the requirements of this standard as it pertains to this device? .....    

Does this standard include acceptance criteria? .....       
If no, include the results of testing in the 510(k).

Does this standard include more than one option or selection of tests? .....       
If yes, report options selected in the summary report table.

Were there any deviations or adaptations made in the use of the standard? .....       
If yes, were deviations in accordance with the FDA supplemental information sheet (SIS) <sup>5</sup>? .....    

Were deviations or adaptations made beyond what is specified in the FDA SIS? .....       
If yes, report these deviations or adaptations in the summary report table.

Were there any exclusions from the standard? .....       
If yes, report these exclusions in the summary report table.

Is there an FDA guidance <sup>6</sup> that is associated with this standard? .....       
If yes, was the guidance document followed in preparation of this 510(k)? .....    

Title of guidance: Use of International Standard ISO-10993, 'Biological Evaluation of Medical Devices Part 1: Evaluation...'

<sup>1</sup> The formatting convention for the title is: [SDO] [numeric identifier] [title of standard] [date of publication]

<sup>2</sup> Authority [21 U.S.C. 360d], <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/Standards/default.htm>

<sup>3</sup> <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>4</sup> The summary report should include: any adaptations used to adapt to the device under review (for example, alternative test methods); choices made when options or a selection of methods are described; deviations from the standard; requirements not applicable to the device; and the name and

address of the test laboratory or certification body involved in conformance assessment to this standard. The summary report includes information on all standards utilized during the development of the device.

<sup>5</sup> The supplemental information sheet (SIS) is additional information which is necessary before FDA recognizes the standard. Found at <http://www.accessdata.fda.gov/scripts/cdrh/cfdocs/cfStandards/search.cfm>

<sup>6</sup> The online search for CDRH Guidance Documents can be found at <http://www.fda.gov/MedicalDevices/DeviceRegulationandGuidance/GuidanceDocuments/default.htm>

**EXTENT OF STANDARD CONFORMANCE  
SUMMARY REPORT TABLE**

STANDARD TITLE  
ISO 10993-1:2003, Biologic Evaluation of medical devices - Part 1

**CONFORMANCE WITH STANDARD SECTIONS\***

SECTION NUMBER 10993-5:2009	SECTION TITLE (b) (4)	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
--------------------------------	--------------------------	--

TYPE OF DEVIATION OR OPTION SELECTED \*  
Tests performed by (b) (4) in accordance with the standard and established professional business practices.

DESCRIPTION

(b) (4)

JUSTIFICATION

(b)(4) Trade Secret Process- Product Specs

SECTION NUMBER 10993-10:2010	SECTION TITLE (b) (4)	CONFORMANCE? <input checked="" type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
---------------------------------	--------------------------	--

TYPE OF DEVIATION OR OPTION SELECTED \*  
Tests performed (b) (4) in accordance with the standard and established professional business practices.

DESCRIPTION

(b)(4) Trade Secret Process- Product Specs

JUSTIFICATION

(b)(4) Trade Secret Process- Product Specs

SECTION NUMBER	SECTION TITLE	CONFORMANCE? <input type="checkbox"/> Yes <input type="checkbox"/> No <input type="checkbox"/> N/A
----------------	---------------	---

TYPE OF DEVIATION OR OPTION SELECTED \*

DESCRIPTION

JUSTIFICATION

\* For completeness list all sections of the standard and indicate whether conformance is met. If a section is not applicable (N/A) an explanation is needed under "justification." Some standards include options, so similar to deviations, the option chosen needs to be described and adequately justified as appropriate for the subject device. Explanation of all deviations or description of options selected when following a standard is required under "type of deviation or option selected," "description" and "justification" on the report. More than one page may be necessary.

♦ Types of deviations can include an exclusion of a section in the standard, a deviation brought out by the FDA supplemental information sheet (SIS), a deviation to adapt the standard to the device, or any adaptation of a section.

This section applies only to requirements of the Paperwork Reduction Act of 1995.

**\*DO NOT SEND YOUR COMPLETED FORM TO THE PRA STAFF EMAIL ADDRESS BELOW.\***


The burden time for this collection of information is estimated to average 1 hour per response, including the time to review instructions, search existing data sources, gather and maintain the data needed and complete and review the collection of information. Send comments regarding this burden estimate or any other aspect of this information collection, including suggestions for reducing this burden, to:

Department of Health and Human Services  
Food and Drug Administration  
Office of Chief Information Officer  
Paperwork Reduction Act (PRA) Staff  
PRAStaff@fda.hhs.gov


*"An agency may not conduct or sponsor, and a person is not required to respond to, a collection of information unless it displays a currently valid OMB control number."*

<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2304 Rev A</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2304 Rev A</b>
---	------------------------------	-----------------------------------


(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




ms

(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num</b> 2324 <b>Rev A</b>
---	------------------------------	--

(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num</b> 2324 <b>Rev A</b>
---	------------------------------	--

(b)(4) Trade Secret Process- Design Specs Testing




<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num</b> 2324 <b>Rev A</b>
---	------------------------------	--


(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




Design Verification (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2324 Rev A
Report: Performance of delivery system, cartridge, and scaffold		

Design Verification (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2304 Rev A
--------------------------	------------------------------	--------------------------


(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




Design Verification (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2324 Rev A
--------------------------------	------------------------------	--------------------------


(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




Design Verification (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2324 Rev A
--------------------------------	------------------------------	--------------------------

(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2337 Rev A</b>
---	------------------------------	-----------------------------------


(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing



l  
f

e


(b)(4) Trade Secret Process- Design Specs Testing



(b)(4) Trade Secret Process- Design Specs Testing




(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




DESIGN VERIFICATION (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2309 Rev C
--------------------------------	------------------------------	--------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing




<b>DESIGN VERIFICATION (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num 2309 Rev C</b>
---	------------------------------	-----------------------------------

(b)(4) Trade Secret Process- Design Specs Testing



Design Verification (DV)	<b>Rotation Medical, Inc</b>	Doc Num 2364 Rev A
--------------------------------	------------------------------	--------------------------


(b)(4) Trade Secret Process- Design Specs Testing



e


<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num</b> 2364 <b>Rev A</b>
---	------------------------------	--

(b)(4) Trade Secret Process- Design Specs Testing



<b>Design Verification (DV)</b>	<b>Rotation Medical, Inc</b>	<b>Doc Num</b> 2364 <b>Rev A</b>
---	------------------------------	--

(b)(4) Trade Secret Process- Design Specs Testing





# COVER SHEET MEMORANDUM

Food and Drug Administration  
Office of Device Evaluation &  
Office of In Vitro Diagnostics and  
Radiological Health

**NOTE: This form is REQUIRED for holds and for final decisions.**

Reviewer Name Maegen Colehour, MS

510(k) Number K140300

Please list CTS decision code: SE - Substantially Equivalent

Hold (Additional Information or Telephone Hold)      Hold Date

Final Decision (SE, SE with Limitations, NSE (select code below), Withdrawn, etc.)

Incomplete Response - Convert Supplement to Amendment (attach email sent to firm)

Add to File

(review staff should follow the instructions and complete the memo/routing sheet at:  
[http://eroom.fda.gov/eRoom/CDRH3/CDRHPreMarketNotification510kProgram/0\\_3bbaZ](http://eroom.fda.gov/eRoom/CDRH3/CDRHPreMarketNotification510kProgram/0_3bbaZ). DCC should refer to that documentation for the close-out code and mail any provided letter.)

The remainder of this form must be filled out for close-outs only

**Class:** II  
**Regulation Number:** 878.3300  
**Product Code:** OWY  
**Additional Product Codes:** ORQ

Please complete the following for a final clearance decision (i.e, SE, SE with Limitations, etc.)	YES	NO
Indications for Use Page (Attach IFU)	X	
510(k) Summary or 510(k) Statement (Attach Summary)	X	
Truthful and Accurate Statement (Must be present for a Final Decision)	X	
Is the device Class III?		X
Is this a combination product?		X
Is this device intended for pediatric use only?		X
Is this a prescription device? (If both prescription & OTC, check both boxes.)	X	
Is clinical data necessary to support the review of this 510(k)?		X
For United States based clinical studies only, did the application include a completed Form FDA 3674, Certification with Requirements of ClinicalTrials.gov Data Bank? (If study was conducted in the United States and Form FDA 3674 was not included or was incomplete, then applicant must be contacted to obtain completed form.)		
Does this device include an Animal Tissue Source?		X
All Pediatric Patients age <= 21		X

DCC  
3.27

Neonate/Newborn (Birth to 28 days)		×
Infant (29 days to < 2 years)		×
Child (2 years to <12 years)		×
Adolescent (12 years to <18 years)		×
Transitional Adolescent A (18 years to <21 years); Special considerations are being given to this group, different from adults age >= 21 (different device design or testing, different protocol procedures, etc.)		×
Transitional Adolescent B (18 years to <21 years); No special considerations compared to adults >= 21 years)		×
Is this device subject to the Tracking Regulation? (Medical Device Tracking Guidance)		×

**Digital Signature Concurrence Table**

(Not all signatures may be required)

Branch Chief Sign-Off	Peter L. Hudson -S 2014.03.25 14:45:18 -04'00'
Division Sign-Off	David Krause -S 2014.03.26 16:17:32 -04'00'