USER: GRAY, ILKA K (ixg)

FOLDER: K984638 - 83 pages (FOI:08003844)

COMPANY: I-FLOW CORP. (IFLOW)

PRODUCT: SET, ADMINISTRATION, INTRAVASCULAR
(FPA)

SUMMARY: Product: PARAGON BOLUS ACCESSORY
SET

DATE REQUESTED: Fri Jan 15 24:00:00 2010

DATE PRINTED: Tue Nov 16 14:09:44 2010

Note: Releasable Version

Table of Contents

510KSUM - 7 pages	1
FOLDER - ADMINISTRATION INTRAVASCULAR SET - 74 pages	8

2/9/99

I-FLOW
CORPORATION

20202 Windrow Rd.
Lake Forest, CA 92630
(800) 448-3569 (714) 206-2700
Fax (714) 206-2600

K984638

SUMMARY OF SAFETY AND EFFECTIVENESS

December 30, 1998

Trade Name: Paragon Bolus Accessory Set

Common Name: Bolus Accessory

Classification Name: Set, Administration, Intravascular

All questions and/or comments concerning this document should be made to:

Robert J. Bard, Esq., R.A.C.
Vice President of Regulatory and Legal Affairs

I-Flow Corporation
20202 Windrow Drive
Lake Forest, CA 92630

Telephone: 949.206.2700
Fax: 949.206.2600

1.0 GENERAL INFORMATION

1.1 Purpose of Submission

- 1.1.1 This submission is intended to notify the Federal Food and Drug Administration that I-Flow Corporation intends to market an accessory device for the Paragon Infusion System (K923875), originally identified as the SideKick 50 Plus and SideKick 100 Plus. This new accessory will be known as the Paragon Bolus Accessory Set, hereafter identified as the Bolus Accessory. No change will be made to the existing Paragon pump or administration sets.
- 1.1.2 Trade Name: Paragon Bolus Accessory Set
- 1.1.3 Common Name: Bolus Accessory Set
- 1.1.4 Classification Name: Set, Administration, Intravascular
- 1.1.5 Product Code: 80 FPA
- 1.1.6 Device Classification: Class II, 880.5440
- 1.1.7 Classification Panel: General Hospital and Personal Use Device

1.2 Statement of Equivalence

- 1.2.1 The Bolus Accessory is substantially equivalent to Patient Control Module (K884505) marketed by Baxter Healthcare Corporation and the I-Flow Bolus Dispenser (K935811).

2.0 PHYSICAL SPECIFICATIONS AND DESCRIPTIONS

2.1 Description of the Bolus Accessory

Wherever the Paragon pump and administration set are mentioned, they may be replaced by any constant 6 psi pressure system.

- 2.1.1 The Bolus Accessory may connect to any Paragon administration set to deliver fixed boluses of medication upon demand by the patient or healthcare provider.
- 2.1.2 The Bolus Accessory consists of plastic housing, medication reservoir, bolus button activator and wrist bands.
- 2.1.3 The bolus button allows patient controlled administration of medication as needed.

2.2 Product Configuration

- 2.2.1 The Bolus Accessory is available in 0.5 ml bolus volume.

2.3 Components and Materials

- 2.3.1 All fluid path components of the Bolus Accessory are in conformance with ISO 10993 Part 1.

2.4 Power Requirements

- 2.4.1 The Bolus Accessory is a mechanical device that requires no external power.

3.0 OPERATIONAL SPECIFICATIONS AND DESCRIPTIONS

3.1 Standard Operating Conditions:

Bolus Volume:	0.5 ml
Refill Time:	variable
Priming/Residual Volume:	<=4 ml
Operating Temperature:	90 ± 2°F
Calibration Solution:	0.9% NaCl
Operating Pressure:	6.0 psi pressure source
Head Height:	0"
Accuracy:	bolus volume: ±10% at 95% confidence interval at the identified lockout times.

3.2 Performance Data: Testing occurred at standard operating conditions. All models performed within the specified accuracy when tested at nominal conditions.

3.3 Safety / Alarm Functions

3.3.1 This device contains no alarms or indicators.

3.3.2 The non-linear refill adds additional patient safety if the bolus button is activated prior to the lockout time.

4.0 BIOLOGICAL SPECIFICATIONS

4.1 Biological testing is in conformance with ISO 10993 Part 1 for all fluid path components.

5.0 CHEMICAL AND DRUG SPECIFICATIONS

5.1 Compatibility

5.1.1 There are no specific drugs referenced in the labeling for the Bolus Accessory.

5.1.2 The Bolus Accessory is intended for general purpose drugs and pain medication.

6.0 INTENDED USE

6.1 The Bolus Accessory is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider.

6.2 The routes of administration intravenous, epidural, intramuscular and subcutaneous.

6.3 The Bolus Accessory is not intended for continuous delivery.

6.4 The Bolus Accessory is single patient use only.

6.5 The Bolus Accessory is suitable for use as an ambulatory device and is intended for use in the hospital, home environment or alternative care sites.

6.6 No testing has been conducted to determine the efficacy of Bolus Accessory for the delivery of blood, blood products, lipids or fat emulsions. The Bolus Accessory is not intended for the delivery of blood, blood products, lipids or fat emulsions.

7.0 STANDARDS

7.1 There are currently no standards established for mechanical PCA infusion devices.

8.0 PACKAGING

8.1 Packaging is suitable for radiation or ETO sterilization.

9.0 STERILIZATION

9.1 The method of sterilization is gamma radiation (cobalt 60).

10.0 COMPARISON TO LEGALLY MARKETED DEVICES

10.1 The Bolus Accessory has the same intended use as the predicate Baxter Pain Control Module and the I-Flow Bolus Dispenser. The Bolus Accessory has similar bolus volumes and lockout times as its predicate devices.

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB - 9 1999

Robert J. Bard, Esq., R.A.C.
Vice President Regulatory and Legal Affairs
I-Flow Corporation
20202 Window Drive
Lake Forest, California 92630

Re: K984638
Trade Name: Paragon Bolus Accessory Set
Regulatory Class: II
Product Code: FPA
Dated: December 30, 1998
Received: December 31, 1998

Dear Mr. Bard:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Pre-market Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your pre-market notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Bard

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K984638

Device Name: Paragon Bolus Accessory Set

Indications for Use:

1. The Paragon Bolus Accessory Set is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider. The routes of administration are intravenous, epidural, intramuscular and subcutaneous.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUED ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use
(Per 21 CFR 801.109)

Patricia Ciccardi Over-The-Counter Use
(Division Sign-Off)
Division of Dental, Infection Control,
and General Hospital Devices

(Optional Format 1-2-96)

510(k) Number K984638

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

FEB - 9 1999

Robert J. Bard, Esq., R.A.C.
Vice President Regulatory and Legal Affairs
I-Flow Corporation
20202 Window Drive
Lake Forest, California 92630

Re: K984638
Trade Name: Paragon Bolus Accessory Set
Regulatory Class: II
Product Code: FPA
Dated: December 30, 1998
Received: December 31, 1998

Dear Mr. Bard:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Mr. Bard

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

510(k) Number (if known): K984638

Device Name: Paragon Bolus Accessory Set

Indications for Use:

1. The Paragon Bolus Accessory Set is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider. The routes of administration are intravenous, epidural, intramuscular and subcutaneous.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUED ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use
(Per 21 CFR 801.109)

Palmer Curran Over-The-Counter Use **3**
 (Division Sign-Off)
 Division of Dental, Infection Control,
 and General Hospital Devices
 510(k) Number K984638
 (Optional Format 1-2-96)

DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service -
Food And Drug Administration

Memorandum

TO: Reviewer(s) - Name(s) HUNG TRINIT

FROM: 510(k) Number K 984638

The Record - It is my recommendation that the subject 510(k) Notification:

- Refused to accept.
- Requires additional information (other than refuse to accept).
- Accepted for review _____.
- Is substantially equivalent to marketed devices.
- NOT substantially equivalent to marketed devices.

De Novo Classification Candidate? YES NO

Other (e.g., exempt by regulation, not a device, duplicate, etc.)

- Is this device subject to Postmarket Surveillance? YES NO
- Is this device subject to the Tracking Regulation? YES NO
- Was clinical data necessary to support the review of this 510(k)? YES NO
- Is it a prescription device? YES NO
- Was this 510(k) reviewed by a Third Party? YES NO
- Special 510(k)? YES NO
- Abbreviated 510(k)? YES NO

This 510(k) contains:

Truthful and Accurate Statement Requested Enclosed
(required for originals received 3-14-95 and after)

A 510(k) summary OR A 510(k) statement

The required certification and summary for class III devices

The indication for use form (required for originals received 1-1-96 and after)

Material of Biological Origin YES NO

The submitter requests under 21 CFR 807.95 (doesn't apply for SEs):

Confidentiality Confidentiality for 90 days Continued Confidentiality exceeding 90 days

Predicate Product Code with class: _____ Additional Product Code(s) with panel (optional): _____

80 FPA CLASS II

Palutia Cuervo

(Branch Code)

(Date)

6/12/99 2/19/99

Screening Checklist For all Premarket Notification 510(k) Submissions

Device Name: <i>Paragon Bolus Accessory Set</i>						K984638	
Submitter (Company): <i>I-Flow Corp</i>							
Items which should be included <i>(circle missing & needed information)</i>	SPECIAL		ABBREVIATED		TRADITIONAL		✓ IF ITEM IS NEEDED AND IS MISSING
	YES	NO	YES	NO	YES	NO	
1. Cover Letter clearly identifies Submission as:							
a) "Special 510(k): Device Modification"							
b) "Abbreviated 510(k)"							
c) Traditional 510(k)							
	GO TO # 2,4		GO TO # 3,4,5		GO TO # 4,5		
2. "SPECIALS" - ONLY FOR MODIFICATIONS TO MANUFACTURER'S OWN CLASS II, III OR RESERVED CLASS I DEVICE							
a) Name & 510(k) number of legally marketed (unmodified) predicate device							
b) STATEMENT - INTENDED USE AND INDICATIONS FOR USE OF MODIFIED DEVICE AS DESCRIBED IN ITS LABELING HAVE NOT CHANGED*							
c) STATEMENT - FUNDAMENTAL SCIENTIFIC TECHNOLOGY OF THE MODIFIED DEVICE HAS NOT CHANGED*							
d) Design Control Activities Summary							
i) Identification of Risk Analysis method(s) used to assess the impact of the modification on the device and its components, and the results of the analysis							
ii) Based on the Risk Analysis, an identification of the verification and/or validation activities required, including methods or tests used and acceptance criteria to be applied							
iii) A declaration of conformity with design controls. The declaration of conformity should include:							
1) A statement signed by the individual responsible, that, as required by the risk analysis, all verification and validation activities were performed by the designated individual(s) and the results demonstrated that the predetermined acceptance criteria were met							
2) A statement signed by the individual responsible, that manufacturing facility is in conformance with design control procedure requirements as specified in 21 CFR 820.30 and the records are available for review.							

5

	SPECIALS		ABBREVIATED		TRADITIONAL		✓ IF ITEM IS NEEDED AND IS MISSING
	YES	NO	YES	NO	YES	NO	
3. ABBREVIATED 510(K): SPECIAL CONTROLS/CONFORMANCE TO RECOGNIZED STANDARDS							
a) For a submission, which relies on a guidance document and/or special control(s), a summary report that describes how the guidance and/or special control(s) was used to address the risks associated with the particular device type							
b) If a manufacturer elects to use an alternate approach to address a particular risk, sufficient detail should be provided to justify that approach.							
c) For a submission, which relies on a recognized standard, a declaration of conformity to the standard. The declaration should include the following:							
i) An identification of the applicable recognized consensus standards that were met							
ii) A specification, for each consensus standard, that all requirements were met, except for inapplicable requirements or deviations noted below							
iii) An identification, for each consensus standard, of any way(s) in which the standard may have been adapted for application to the device under review, e.g., an identification of an alternative series of tests that were performed							
iv) An identification, for each consensus standard, of any requirements that were not applicable to the device							
v) A specification of any deviations from each applicable standard that were applied							
vi) A specification of the differences that may exist, if any, between the tested device and the device to be marketed and a justification of the test results in these areas of difference							
vii) Name/address of test laboratory/certification body involved in determining the conformance of the device with applicable consensus standards and a reference to any accreditations for those organizations							
d) Data/information to address issues not covered by guidance documents, special controls, and/or recognized standards							

Ce

REVISED:3/14/95

THE 510(K) DOCUMENTATION FORMS ARE AVAILABLE ON THE LAN UNDER 510(K)
BOILERPLATES TITLED "DOCUMENTATION" AND MUST BE FILLED OUT WITH
EVERY FINAL DECISION (SE, NSE, NOT A DEVICE, ETC.).

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K _____

Reviewer: _____

Division/Branch: _____

Device Name: _____

Product To Which Compared (510(K) Number If Known): _____

YES NO

	YES	NO	
1. Is Product A Device			If NO = Stop
2. Is Device Subject To 510(k)?			If NO = Stop
3. Same Indication Statement?			If YES = Go To 5
4. Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NE
5. Same Technological Characteristics?			If YES = Go To 7
6. Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 8
7. Descriptive Characteristics Precise Enough?			If NO = Go To 10 If YES = Stop SE
8. New Types Of Safety Or Effectiveness Questions?			If YES = Stop NE
9. Accepted Scientific Methods Exist?			If NO = Stop NE
10. Performance Data Available?			If NO = Request Data
11. Data Demonstrate Equivalence?			Final Decision:

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

8

1. Intended Use:
2. Device Description: Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Is the device life-supporting or life sustaining? Is the device implanted (short-term or long-term)? Does the device design use software? Is the device sterile? Is the device for single use? Is the device for home use or prescription use? Does the device contain drug or biological product as a component? Is this device a kit? Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

EXPLANATIONS TO "YES" AND "NO" ANSWERS TO QUESTIONS ON PAGE 1 AS NEEDED

1. Explain why not a device:
2. Explain why not subject to 510(k):
3. How does the new indication differ from the predicate device's indication:
4. Explain why there is or is not a new effect or safety or effectiveness issue:
5. Describe the new technological characteristics:
6. Explain how new characteristics could or could not affect safety or effectiveness:
7. Explain how descriptive characteristics are not precise enough:
8. Explain new types of safety or effectiveness questions raised or why the questions are not new:
9. Explain why existing scientific methods can not be used:
10. Explain what performance data is needed:
11. Explain how the performance data demonstrates that the device is or is not substantially equivalent:

ATTACH ADDITIONAL SUPPORTING INFORMATION

**MEMO TO THE RECORD
510(K) REVIEW**

DATE: 2/1/99
FROM: Hung Trinh *HT 2/3/99*
DOCUMENT: K984638
COMPANY NAME: I-Flow Corp
DEVICE NAME: Paragon Bollus Accessory set

OFFICE: HFZ-480
DIVISION: DDIG/GHDB

Contact point: Robert Bard, VC Regulatory and Legal Affairs
949-206-2700
949-206-2600 (fx)

NARRATIVE DEVICE DESCRIPTION

1. **SUMMARY DESCRIPTION OF THE SUBMISSION UNDER REVIEW:**
The firm intends to extend the administration set product line for the Paragon Infusion System (K923875), originally identified as the SideKick 50 Plus and SideKick 100 Plus, to include a bolus accessory set.
2. **INTENDED USE:**
The device is intended to allow patient controlled bolus delivery.
3. **DEVICE DESCRIPTION:**
 - A. **Life-supporting or life-sustaining:** no
 - B. **Implant (short-term or long-term):** no
 - C. **Is the device sterile?** yes
If yes, is sterility information provided? yes
Gamma radiation (cobalt 60) or (b)(4), (b)(5) Gamma radiation validation methodology is by ANSI/AAMI ST32-1191/EN 552; dosage from 25 to 35 kGy, (b)(4), (b)(5) sterilization complies with ANSI/AAMI/ISO (b)(4), (b)(5)
 - D. **Is the device for single use?** Single patient use
 - E. **Is the device for prescription use?** yes
If yes, is prescription labeling included? yes
 - F. **Is the device for home use or portable?** yes
 - G. **Does the device contain drug or biological product as a component?** no
 - H. **Is this device a kit?** no
 - I. **Software-driven:** no
 - J. **Electrically Operated:** no
 - K. **Applicable standards to which conformance has been demonstrated (e.g., IEC, ANSI, ASTM, etc.):** no
 - L. **Device(s) to which equivalence is claimed, manufacturer, and 510(k) number or preamendment status:**
Baxter Basal/Bolus Infusor (K884505)
I-Flow Bolus Dispenser (K935811)
 - M. **Submission provides comparative specifications** yes
comparative in vitro data no
performance data yes

10

animal testing no
clinical testing no
biocompatibility testing no (same material as predicate)

N. Provide a statement of how the device is either similar to and/or different from other marketed devices, plus data (if necessary) to support the statement. Provide a summary about the devices design, materials, physical properties and toxicology profile if important.

The device consists of a plastic housing, medication reservoir, bolus button activator, and wrist band. The device can be connected to any Paragon administration set to deliver fixed boluses of medication upon demand by the patient or healthcare provider. This device **does not** deliver continuous medication. The refill time of the medication reservoir is determined by the flow control orifice of the Paragon administration set, and therefore, is inversely proportional to the flow rate (see page 16 of 19).

The Baxter Patient Control module is very similar to the subject device in that it is composed of a plastic housing, medication reservoir, bolus button activator, and wrist bands. It also requires a pressure source to fill the medication reservoir.

Comparative specs:

Comparison Element	Paragon Bolus Accessory Set	Baxter Patient Control Module	I-Flow Bolus Dispenser
Bolus Volume	0.5 ml; $\pm 10\%$ @ 95% confidence interval	0.5 ml	0-1.0 ml
Bolus Lockout time	3.6, 9, 18, 35 and 70 min; $\pm 15\%$ @ 95% confidence interval	6, 15, 60 min	15, 30,, 60, 120 min
Pump Volume	100 ml	65 ml	30
Pressure source	Mechanical spring energy of the pump	Strain energy of elastomeric membranes	Mechanical Spring
Fluid Reservoir	PVC drug bag	Elastomeric membrane	PVC drug bag or polypropyplene syringe

A more detailed comparison is available on page 22 of 23.

O. Does the submission include a summary of safety and effectiveness information upon which an equivalence determination is based? yes

P. RECOMMENDATION:

I believe that this device is equivalent to: 80 FPA

Classification should be based on:

11

880.5440 Intravascular (IV) administration set

Class: II

If the device is substantially equivalent to a class III device, does the submission include: (1) certification that a reasonable search of all information known, or otherwise available, about the generic type of device has been performed and (2) a summary description of the types of safety and effectiveness problems associated with the type of device and a citation to the literature, or other sources of information, upon which they have based the description? n/a

12

"SUBSTANTIAL EQUIVALENCE" (SE) DECISION MAKING DOCUMENTATION

K984638

Reviewer: Hung Trinh

Division/Branch: DDIGD/GHDB

Device Name: Paragon Bolus Administration set

Product To Which Compared (510(K) Number If Known):

YES NO

1.	Is Product A Device	X		If NO = Stop
2.	Is Device Subject To 510(k)?	X		If NO = Stop
3.	Same Indication Statement?	X		If YES = Go To 5
4.	Do Differences Alter The Effect Or Raise New Issues of Safety Or Effectiveness?			If YES = Stop NE
5.	Same Technological Characteristics?	X		If YES = Go To 7
6.	Could The New Characteristics Affect Safety Or Effectiveness?			If YES = Go To 8
7.	Descriptive Characteristics Precise Enough?		X	If NO = Go To 10 If YES = Stop SE
8.	New Types Of Safety Or Effectiveness Questions?			If YES = Stop NE
9.	Accepted Scientific Methods Exist?			If NO = Stop NE
10.	Performance Data Available?	X		If NO = Request Data
11.	Data Demonstrate Equivalence?	X		Final Decision: SE

Note: In addition to completing the form on the LAN, "yes" responses to questions 4, 6, 8, and 11, and every "no" response requires an explanation.

B

20202 Windrow Dr
Lake Forest, CA 92630
Tel: 949.206.2670
Fax: 949.206.2603

Fax

To: Hung Trinh **From:** Robert J Bard, Esq.

Fax: 301.480.3002 **Pages:** 6

Phone: 301.594.1287 x 130 **Date:** 02/02/99

Re: K984063, K984146, K984638 **CC:**

Urgent **For Review** **Please Comment** **Please Reply** **Please Recycle**

● Comments: The documents accompanying this facsimile transmission contains information which may be legally privileged and confidential. The information is intended only for the use of the recipient named below. If you have received this facsimile in error, please immediately notify us by telephone to arrange for return of the original documents to us. Any disclosure, copying, distribution or the taking of any action in reliance on the contents of this faxed information is strictly prohibited.

4

I-FLOW
CORPORATION

20202 Windrow Drive
Lake Forest, CA 92630
(800) 448-3569 (949) 206-2700
Fax (949) 206-2600

February 2, 1999

VIA FACSIMILE

Hung Trinh
Food and Drug Administration
Center for Devices and Radiological Health
Office Device Evaluation
9200 Corporate Blvd.
Rockville, Maryland 20850

Re: K984063
K984146
K984638

Dear Mr. Trinh:

Pursuant to our conversation of February 2, 1999, I have made the following changes to the above referenced premarket notifications.

K984063 The **Indications for Use** description has been changed to a single statement. Lines 3, 4 and 5 have been deleted and Line 2 has been incorporated into Line 1.

No other changes have been made to premarket notification K984063.

K984146 A supplemental submission was made on January 13, 1999 to include synovial infusions as an additional indication for use. The basis of this supplemental submission was the finding of a predicate device (the McKinley Outbound K982256) that included (b)(4), (b)(5) as an indication for use.

I-Flow has agreed to remove the January 13 revision from our submission based on your statement that the McKinley will not be allowed to have an indication for use that includes (b)(4), (b)(5). Based on your statement, I-Flow believes the Office of Device Evaluation will require McKinley to modify their Indications for Use Statement and that the posting on the FDA 510(k) webpage will be changed.

15


The **Indications for Use** description has been changed to a single statement including both Lines 1 and 2.

K984638 The **Indications for Use** description has been changed to a single statement including both Lines 1 and 2.

No other changes have been made to premarket notification K984638.

If you have any additional issues specific to the above identified premarket notifications, either to the information provided or in general, I can be reached at 949.206.2670 or 800.206.2700.

Sincerely,

Robert J. Bard, Esq., R.A.C.
Vice President, Regulatory and Legal Affairs

16

510(k) Number (if known): K984638

Device Name: Paragon Bolus Accessory Set

Indications for Use:

1. The Paragon Bolus Accessory Set is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider. The routes of administration are intravenous, epidural, intramuscular and subcutaneous.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUED ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over-The-Counter Use

17

(Optional Format 1-2-96)

Food and Drug Administration
Center for Devices and
Radiological Health
Office of Device Evaluation
Document Mail Center (HFZ-401)
9200 Corporate Blvd.
Rockville, Maryland 20850

December 31, 1998

I-FLOW CORP.
20202 WINDROW DR.
LAKE FOREST, CA 92630
ATTN: ROBERT J. BARD

510(k) Number: K984638
Received: 31-DEC-1998
Product: PARAGON BOLUS
ACCESSORY SET

The Center for Devices and Radiological Health (CDRH), Office of Device Evaluation (ODE), has received the Premarket Notification you submitted in accordance with Section 510(k) of the Federal Food, Drug, and Cosmetic Act (Act) for the above referenced product. We have assigned your submission a unique 510(k) number that is cited above. Please refer prominently to this 510(k) number in any future correspondence that relates to this submission. We will notify you when the processing of your premarket notification has been completed or if any additional information is required. YOU MAY NOT PLACE THIS DEVICE INTO COMMERCIAL DISTRIBUTION UNTIL YOU RECEIVE A LETTER FROM FDA ALLOWING YOU TO DO SO.

On January 1, 1996, FDA began requiring that all 510(k) submitters provide on a separate page and clearly marked "Indication For Use" the indication for use of their device. If you have not included this information on a separate page in your submission, please complete the attached and amend your 510(k) as soon as possible. Also if you have not included your 510(k) Summary or 510(k) Statement, or your Truthful and Accurate Statement, please do so as soon as possible. There may be other regulations or requirements affecting your device such as Postmarket Surveillance (Section 522(a)(1) of the Act) and the Device Tracking regulation (21 CFR Part 821). Please contact the Division of Small Manufacturers Assistance (DSMA) at the telephone or web site below for more information.

Please remember that all correspondence concerning your submission MUST be sent to the Document Mail Center (HFZ-401) at the above letterhead address. Correspondence sent to any address other than the Document Mail Center will not be considered as part of your official premarket notification submission. Because of equipment and personnel limitations, we cannot accept telefaxed material as part of your official premarket notification submission, unless specifically requested of you by an FDA official. Any telefaxed material must be followed by a hard copy to the Document Mail Center (HFZ-401).

You should be familiar with the manual entitled, "Premarket Notification 510(k) Regulatory Requirements for Medical Devices" available from DSMA. If you have other procedural or policy questions, or want information on how to check on the status of your submission (after 90 days from the receipt date), please contact DSMA at (301) 443-6597 or its toll-free number (800) 638-2041, or at their Internet address <http://www.fda.gov/cdrh/dsmamain.html> or me at (301) 594-1190.

Sincerely yours,

Marjorie Shulman
Consumer Safety Officer
Premarket Notification Staff
Office of Device Evaluation
Center for Devices and Radiological Health

18

15984638

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH Premarket Submission Cover Sheet

Date of Submission: 12/30/98 FDA Document Number: ω

Section A	Type of Submission		
<input checked="" type="checkbox"/> 510(k)	<input type="checkbox"/> IDE	<input type="checkbox"/> PMA	<input type="checkbox"/> PMA Supplement - Regular
<input type="checkbox"/> 510(k) Add'l information	<input type="checkbox"/> IDE Amendment	<input type="checkbox"/> PMA Amendment	<input type="checkbox"/> PMA Supplement - Special
	<input type="checkbox"/> IDE Supplement	<input type="checkbox"/> PMA Report	<input type="checkbox"/> PMA Supplement - 30 day
	<input type="checkbox"/> IDE Report		<input type="checkbox"/> PMA Supplement - Panel Track

Section B1 Reason for Submission — 510(k)s Only

New device Additional or expanded indications Change in technology, design, materials, or manufacturing process

Other reason (specify):

Section B2 Reason for Submission — PMAs Only

<input type="checkbox"/> New device	<input type="checkbox"/> Change in design, component, or specification:	<input type="checkbox"/> Location change:
<input type="checkbox"/> Withdrawal	<input type="checkbox"/> Software	<input type="checkbox"/> Manufacturer
<input type="checkbox"/> Additional or expanded indications	<input type="checkbox"/> Color Additive	<input type="checkbox"/> Sterilizer
<input type="checkbox"/> Licensing agreement	<input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Packager
<input type="checkbox"/> Labeling change:	<input type="checkbox"/> Process change:	<input type="checkbox"/> Report submission:
<input type="checkbox"/> Indications	<input type="checkbox"/> Manufacturer	<input type="checkbox"/> Annual or periodic
<input type="checkbox"/> Instructions	<input type="checkbox"/> Sterilizer	<input type="checkbox"/> Post-approval study
<input type="checkbox"/> Performance Characteristics	<input type="checkbox"/> Packager	<input type="checkbox"/> Adverse reaction
<input type="checkbox"/> Shelf life		<input type="checkbox"/> Device defect
<input type="checkbox"/> Trade name	<input type="checkbox"/> Response to FDA correspondence (specify below)	<input type="checkbox"/> Amendment
<input type="checkbox"/> Other (specify below)	<input type="checkbox"/> Request for applicant hold	
<input type="checkbox"/> Change in ownership	<input type="checkbox"/> Request for removal of applicant hold	
<input type="checkbox"/> Change in correspondent	<input type="checkbox"/> Request for extension	
<input type="checkbox"/> Other reason (specify):	<input type="checkbox"/> Request to remove or add manufacturing site	

Section B3 Reason for Submission — IDEs Only

<input type="checkbox"/> New device	<input type="checkbox"/> Change in:	<input type="checkbox"/> Response to FDA letter concerning:
<input type="checkbox"/> Addition of institution	<input type="checkbox"/> Correspondent	<input type="checkbox"/> Conditional approval
<input type="checkbox"/> Expansion / extension of study	<input type="checkbox"/> Design	<input type="checkbox"/> Deemed approved
<input type="checkbox"/> IRB certification	<input type="checkbox"/> Informed consent	<input type="checkbox"/> Deficient final report
<input type="checkbox"/> Request hearing	<input type="checkbox"/> Manufacturer	<input type="checkbox"/> Deficient progress report
<input type="checkbox"/> Request waiver	<input type="checkbox"/> Manufacturing	<input type="checkbox"/> Deficient investigator report
<input type="checkbox"/> Termination of study	<input type="checkbox"/> Protocol - feasibility	<input type="checkbox"/> Disapproval
<input type="checkbox"/> Withdrawal of application	<input type="checkbox"/> Protocol- other	<input type="checkbox"/> Request extension of time to respond to FDA
<input type="checkbox"/> Unanticipated adverse effect	<input type="checkbox"/> Sponsor	<input type="checkbox"/> Request meeting
<input type="checkbox"/> Emergency use:	<input type="checkbox"/> Report submission:	<input type="checkbox"/> IOL submissions only:
<input type="checkbox"/> Notification of emergency use	<input type="checkbox"/> Current investigator	<input type="checkbox"/> Change in IOL style
<input type="checkbox"/> Additional information	<input type="checkbox"/> Annual progress	<input type="checkbox"/> Request for protocol waiver
<input type="checkbox"/> Other reason (specify):	<input type="checkbox"/> Site waiver limit reached	
	<input type="checkbox"/> Final	

skg

19
HO
II

FDA Document Number:

Section C

Product Classification

Product code: 80 FPA	C.F.R. Section: 80.5440	Device class: <input type="checkbox"/> Class I <input type="checkbox"/> Class III <input checked="" type="checkbox"/> Class II <input type="checkbox"/> Unclassified
Classification panel: General Hospital and Personal Use Device		

Section D

Information on 510(k) Submissions

Product codes of devices to which substantial equivalence is claimed:				Summary of, or statement concerning, safety and effectiveness data: <input checked="" type="checkbox"/> 510(k) summary attached <input type="checkbox"/> 510(k) statement
1 80 FPA	2 80 FRN	3	4	
5	6	7	8	

Information on devices to which substantial equivalence is claimed:

510(k) Number	Trade or proprietary or model name	Manufacturer
1 K884505	1 Patient Control Module	1 Baxter Healthcare Corp.
2 K935811	2 Bolus Dispenser	2 I-Flow Corp.
3	3	3
4	4	4
5	5	5
6	8	8

Section E

Product Information — Applicable to All Applications

Common or usual name or classification name:
 Set, Administration, Intravascular

Trade or proprietary or model name	Model number
1 Paragon Bolus Accessory Set	1
2	2
3	3
4	4
5	5
6	6

FDA document numbers of all prior related submissions (regardless of outcome):

1	2	3	4	5	6
7	8	9	10	11	12

Data included in submission: Laboratory testing Animal trials Human trials


Indications (from labeling):

The Paragon Bolus Accessory Set is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider. The routes of administration are intravenous, epidural, intramuscular and subcutaneous.

20


				FDA Document Number:			
Section F Manufacturing / Packaging / Sterilization Sites							
<input checked="" type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		FDA establishment registration number: 2026095		<input checked="" type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer		<input type="checkbox"/> Contract sterilizer <input type="checkbox"/> Repackager / relabeler	
Company / Institution name: I-Flow Corporation							
Division name (if applicable):				Phone number (include area code): (949) 206-2700 ext. 2670			
Street address: 20202 Windrow Drive				FAX number (include area code): (949) 206-2603			
City: Lake Forest		State / Province: CA		Country: U.S.A.		ZIP / Postal Code: 92630	
Contact name: Robert J. Bard, Esq., R.A.C.							
Contact title: Vice President of Regulatory and Legal Affairs							
<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		FDA establishment registration number:		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer		<input type="checkbox"/> Contract sterilizer <input type="checkbox"/> Repackager / relabeler	
Company / Institution name:							
Division name (if applicable):				Phone number (include area code): ()			
Street address:				FAX number (include area code): ()			
City:		State / Province:		Country:		ZIP / Postal Code:	
Contact name:							
Contact title:							
<input type="checkbox"/> Original <input type="checkbox"/> Add <input type="checkbox"/> Delete		FDA establishment registration number:		<input type="checkbox"/> Manufacturer <input type="checkbox"/> Contract manufacturer		<input type="checkbox"/> Contract sterilizer <input type="checkbox"/> Repackager / relabeler	
Company / Institution name:							
Division name (if applicable):				Phone number (include area code): ()			
Street address:				FAX number (include area code): ()			
City:		State / Province:		Country:		ZIP / Postal Code:	
Contact name:							
Contact title:							

21

				FDA Document Number:	
Section G Applicant or Sponsor					
Company / Institution name: I-Flow Corporation				FDA establishment registration number: 2026095	
Division name (if applicable):				Phone number (include area code): (949) 206-2700 ext. 2670	
Street address: 20202 Windrow Drive				FAX number (include area code): (949) 206-2603	
City: Lake Forest	State / Province: CA	Country: U.S.A.	ZIP / Postal Code: 92630		
Signature: 					
Name: Robert J. Bard, Esq., R.A.C.					
Title: Vice President of Regulatory and Legal Affairs					
Section H Submission correspondent (if different from above)					
Company / Institution name:					
Division name (if applicable):				Phone number (include area code): ()	
Street address:				FAX number (include area code): ()	
City:	State / Province:	Country:	ZIP / Postal Code:		
Contact name:					
Contact title:					

Your voluntary completion of this Premarket Submission Cover Sheet will not affect any FDA decision concerning your submission, but will help FDA's Center for Devices and Radiological Health process your submission more efficiently. The information you provide should apply *only* to a single accompanying submission. Please do not send cover sheets for any previous submissions. See the instructions for additional information on completing the cover sheet. If you have a question concerning completion of the cover sheet, please contact the Division of Small Manufacturers Assistance at (800) 638-2041 or (301) 443-6597.

Handwritten mark

20202 Windrow Drive
 Lake Forest, CA 92630
 (800) 448-3569 (949) 206-2700
 Fax (949) 206-2600

Premarket Notification – 510(k)

Via Federal Express
 December 30, 1998

Food and Drug Administration
 Center for Devices and Radiological Health
 Office of Device Evaluation
 Document Mail Center HFZ – 401
 9200 Corporate Blvd.
 Rockville, Maryland 20850

SKJ
 RECEIVED
 31 DEC 30 10 50
 FDA/CDRH/ODE/DWG

Reviewing Staff:

In accordance with §510(k) of the Federal Food, Drug, and Cosmetic Act and in conformance with Title 21 CFR §807.81, I-Flow Corporation is submitting this premarket notification for the *Paragon Bolus Accessory Set* prior to the introduction into interstate commerce for commercial distribution.

I-Flow intends to market an accessory device for the Paragon Infusion System (K923875), originally identified as the SideKick 50 Plus and SideKick 100 Plus. This new accessory will be known as the *Paragon Bolus Accessory Set*, hereafter identified as the *Bolus Accessory*. The *Bolus Accessory* is substantially equivalent to the Patient Control Module (K884505) marketed by Baxter Healthcare Corporation and the I-Flow Bolus Dispenser (K935811).

All questions and/or comments concerning this document should be made to:
 Robert J. Bard, Esq., R.A.C.
 Vice President Regulatory and Legal Affairs

I-Flow Corporation
 20202 Windrow Drive
 Lake Forest, CA 92630
 Telephone: 949.206.2700
 Fax: 949.206.2600

Sincerely,

Robert J. Bard, Esq., R.A.C.
 Vice President Regulatory and Legal Affairs

23

I-FLOW
CORPORATION

20202 Windrow Rd.
Lake Forest, CA 92630
(800) 448-3569 (714) 206-2700
Fax (714) 206-2600

**PREMARKET NOTIFICATION
TRUTHFUL AND ACCURATE STATEMENT
(As required by 21 CFR 807.87(j))**

I certify that, in my capacity as the Vice President of Regulatory and Legal Affairs of I-Flow Corporation, I believe to the best of my knowledge, that all data and information submitted in the premarket notification for the Paragon Bolus Accessory Set are truthful and accurate and that no material fact has been omitted.

Signature

Robert J Bard, Vice President of Regulatory and Legal Affairs

Name	Title
I-Flow Corporation	12/30/98
Company	Dated

Premarket Notification (510(k) Number)

24

I-FLOW CORPORATION

20202 Windrow Rd.
Lake Forest, CA 92630
(800) 448-3569 (714) 206-2700
Fax (714) 206-2600

510(k) Number (if known): _____

Device Name: Paragon Bolus Accessory Set

Indications for Use:

1. The Paragon Bolus Accessory Set is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider.
2. The routes of administration are intravenous, epidural, intramuscular and subcutaneous.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUED ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

25

Prescription Use _____
(Per 21 CFR 801.109)

OR

Over-The-Counter Use _____

(Optional Format 1-2-96)

TABLE OF CONTENTS

1.0 GENERAL INFORMATION Page 1

2.0 PHYSICAL SPECIFICATIONS AND DESCRIPTION Page 2

3.0 OPERATIONS SPECIFICATIONS AND DESCRIPTION..... Page 4

4.0 BIOLOGICAL SPECIFICATIONS..... Page 6

Chart #1 Page 7

Chart #2 Page 8

Chart #3 Page 9

Chart #4 Page 10

Chart #5 Page 11

Chart #6a Page 12

Chart #6b Page 13

5.0 CHEMICAL AND DRUG SPECIFICATIONS..... Page 14

6.0 INTENDED USE Page 14

7.0 LABELS AND LABELING Page 14

8.0 STANDARDS..... Page 16

9.0 PACKAGING Page 16

10.0 STERILIZATION INFORMATION..... Page 16

11.0 COMPARISON TO LEGALLY MARKETED DEVICES..... Page 16

Appendix A – Paragon Bolus Accessory Set Drawings

Appendix B – Paragon Bolus Accessory Set Labeling

Appendix C – Predicate Labeling

- **Baxter Patient Control Module**

Appendix D – Reference Labeling

- **I-Flow Paragon Infusion System**

Appendix E – Summary of Safety and Effectiveness

26

1.0 GENERAL INFORMATION

1.1 Purpose of Submission

- 1.1.1 This submission is intended to notify the Federal Food and Drug Administration that I-Flow Corporation intends to market an accessory device for the Paragon Infusion System (K923875), originally identified as the SideKick 50 Plus and SideKick 100 Plus. This new accessory will be known as the Paragon Bolus Accessory Set, hereafter identified as the Bolus Accessory. No change will be made to the existing Paragon pump or administration sets.
- 1.1.2 Trade Name: Paragon Bolus Accessory Set
- 1.1.3 Common Name: Bolus Accessory Set
- 1.1.4 Classification Name: Set, Administration, Intravascular
- 1.1.5 Product Code: 80 FPA
- 1.1.6 Device Classification: Class II, 880.5440
- 1.1.7 Classification Panel: General Hospital and Personal Use Device

1.2 Statement of Equivalence

- 1.2.1 The Bolus Accessory is substantially equivalent to Patient Control Module (K884505) marketed by Baxter Healthcare Corporation and the I-Flow Bolus Dispenser (K935811).

1.3 Overview of the Existing Paragon Infusion System

- 1.3.1 The Paragon Infusion System consists of two main components, the Paragon pump and administration set. This premarket notification proposes a new accessory for use with the Paragon Infusion System, i.e. the Bolus Accessory. No change will be made to the existing Paragon pump or administration sets.
- 1.3.2 The Paragon Pump
- 1.3.2.1 The Paragon pump consists of two cylindrical shells. The top half of the pump has (b)(4) (b)(4) (b)(4) (b)(4) of the bottom half of the pump.
- 1.3.2.2 The top incorporates a (b)(4) (b)(4) which applies a load to the pliable drug bag. The load is applied to the drug bag by way of (b)(4) (b)(4) mechanism. The (b)(4) (b)(4) creates a near constant pressure in the drug bag.
- 1.3.2.3 The bottom half of the pump is slotted to allow for positioning of the administration set.
- 1.3.2.4 When the top and bottom halves of the pump are fully (b)(4) (b)(4) together, the (b)(4) (b)(4) contacts the drug bag and acts as the pressurizing element.
- 27
D

1.3.3 The Paragon Administration Set

1.3.3.1 The Paragon administration sets consist of a PVC drug bag attached to the administration line.

1.3.3.2 Each administration set has either flow control tubing or orifice to regulate the flow rate to the patient.

1.3.3.2.1 The delivery time characteristic is derived from the flow rate of the device which is in turn approximated by (b)(4), (b)(5)

(b)(4), (b)(5)
(b)(4), (b)(5)

1.3.3.2.2

(b)(4), (b)(5)

(b)(4), (b)(5)

equation provides an approximation of the actual delivery time.

2.0 PHYSICAL SPECIFICATIONS AND DESCRIPTIONS

2.1 Description of the Bolus Accessory

Wherever the Paragon pump and administration set are mentioned, they may be replaced by any constant 6 psi pressure system. See Appendix A for drawings.

2.1.1 The Bolus Accessory may connect to any Paragon administration set to deliver fixed boluses of medication upon demand by the patient or healthcare provider.

2.1.2 The Bolus Accessory consists of plastic housing, medication reservoir, bolus button activator and wrist bands.

2.1.3 The bolus button allows patient controlled administration of medication as needed.

(b)(4), (b)(5)

(b)(4), (b)(5)

Top View

Side View

2.2 Bolus Refill Cycle (i.e. Lockout Time)

2.2.1 The Paragon pump is the pressure source pushing fluid into the inlet port of the Bolus Accessory.


2.2.2 The flow control orifice of the Paragon administration set establishes a lockout period during which the patient cannot receive another full bolus of medication.

2.2.2.1 The time required to fill the bolus cavity is called the “lockout time” or “refill cycle”.

2.2.2.2 Lockout times vary according to the flow rate of the Paragon administration set. See table below.

Flow Rate of Paragon Set (ml/hr)	Lockout Time (min.)
0.5	70.0
1.0	35.0
2.0	18.0
4.0	9.0
10	3.6

The lockout time is approximated as follows:

2.2.3 As fluid flows into the bolus cavity, the diaphragm is continuously displaced until the cavity is full.

2.2.3.1 As a safety feature, the bolus cavity is filled non-linearly, slower in the beginning and more rapidly towards the end of the refill cycle. See chart #6 for the pressure profile.

2.2.3.2 The bolus reservoir is 0.5 ml volume.

2.3 Bolus Activation

2.3.1 On the opposite side of the diaphragm is a bolus button. The button moves with the diaphragm as it is displaced.

2.3.2 At any time while the bolus cavity is filling, the button can be activated by the user. Only when the bolus cavity is full will the patient receive the entire bolus.

2.3.3 When the button is activated, the fluid in the bolus cavity is pushed out and forced through a check valve. At this point, the refill cycle begins.

29

2.4 Product Configuration

Bolus Accessory Model

2.4.1 B000000: 0.5 ml bolus volume.

2.5 Components and Materials

No change will be made to the Paragon pump or administration set. The Paragon pump and administration set components remain the same as K923875.

The Bolus Accessory is a disposable device intended for single patient use.

2.5.1 Non-fluid path components (bolus component only)

2.5.1.1 Bolus Button: (b)(4), (b)(5) (b)(4), (b)(5)

2.5.1.2 Outer Housing: (b)(4), (b)(5) (b)(4), (b)(5)

2.5.1.3 Luer Cap: (b)(4), (b)(5) (b)(4), (b)(5)
(b)(4), (b)(5)

2.5.2 Fluid path components (bolus component only)

2.5.2.1 Inner Housing: (b)(4), (b)(5) (b)(4), (b)(5)

2.5.2.2 Diaphragm: (b)(4), (b)(5) (b)(4), (b)(5)

2.5.2.3 Check Valve: (b)(4), (b)(5) (b)(4), (b)(5)

2.5.2.4 Male Luer Adapter: (b)(4), (b)(5) (b)(4), (b)(5) or equivalent.

2.5.2.5 Female Luer Adapter: (b)(4), (b)(5) (b)(4), (b)(5) or equivalent.

2.5.2.6 Tubing (Make-up): (b)(4), (b)(5) (b)(4), (b)(5)

2.5.2.7 Solvent Bonding: (b)(4), (b)(5) (b)(4), (b)(5)
(b)(4), (b)(4), (b)(5) or equivalent.

2.6 Power Requirements

2.6.1 The Bolus Accessory is a mechanical device that requires no external power to operate.

3.0 OPERATIONAL SPECIFICATIONS AND DESCRIPTIONS

3.1 Standard Operating Conditions:

Bolus Volume: 0.5 ml

Refill Time: variable

Priming/Residual Volume: <=4 ml

Operating Temperature: 90 ± 2°F

Calibration Solution: 0.9% NaCl

Operating Pressure: 6.0 psi pressure source

Head Height: 0"

Accuracy: bolus volume: ±10% at 95% confidence interval at the identified lockout times.

30

3.2 Performance Data: (b)(4), (b)(5) (b)(4), (b)(5)

(b)(4), (b)(5)

(b)(4), (b)(5)

31

(b)(4), (b)(5)

3.2.5

(b)(4), (b)(5)

3.2.6

3.3 **Safety/Alarm Functions**

- 3.3.1 The Bolus Accessory will not be recommended for any application that exceeds the minimum internal pressure of the system.
- 3.3.2 If for any reason the patient needs to stop his or her infusions, each administration set is supplied with a pinch clamp to stop the infusion.
- 3.3.3 This device contains no alarms or indicators for flow other than visual.
- 3.3.4 This device contains no alarms or indicators to detect air in line or an occlusion; however, each set may include an integrated air-eliminating filter.
- 3.3.5 The non-linear refill adds additional patient safety if the bolus button is activated prior to the lockout time.

4.0 **BIOLOGICAL SPECIFICATIONS**


- 4.1 Biological testing is in conformance with ISO 10993 Part 1 for all fluid path components.
- 4.2 The Bolus Accessory is categorized as follows:
 - 4.2.1 Device Category: External Communicating Device.
 - 4.2.2 Body Contact: Blood Path, Indirect
 - 4.2.3 Contact Duration: Prolonged

32

Chart #1
 Bolus Accessory
Bolus Volume Accuracy
 various lockout times

Bolus #	Volume (ml)				
	3.6 min. lockout	9 min. lockout	18 min. lockout	35 min. lockout	70 min. lockout
1	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
2	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
3	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
4	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
5	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
6	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
7	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)
8	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)	(b)(4), (b)(5)


Min
 Max
 N
 Average
 Std. Dev.

3/23


Chart #2
Bolus Accessory
Bolus Volume Accuracy
 3.6 and 70 min. lockout

Bolus #	3.6 min. lockout	70 min. lockout
1	(b)(4), (b)(5)	(b)(4), (b)(5)
2	(b)(4), (b)(5)	(b)(4), (b)(5)
3	(b)(4), (b)(5)	(b)(4), (b)(5)
4	(b)(4), (b)(5)	(b)(4), (b)(5)
5	(b)(4), (b)(5)	(b)(4), (b)(5)
6	(b)(4), (b)(5)	(b)(4), (b)(5)
7	(b)(4), (b)(5)	(b)(4), (b)(5)
8	(b)(4), (b)(5)	(b)(4), (b)(5)
9	(b)(4), (b)(5)	(b)(4), (b)(5)
10	(b)(4), (b)(5)	(b)(4), (b)(5)
11	(b)(4), (b)(5)	(b)(4), (b)(5)
12	(b)(4), (b)(5)	(b)(4), (b)(5)
13	(b)(4), (b)(5)	(b)(4), (b)(5)
14	(b)(4), (b)(5)	(b)(4), (b)(5)
15	(b)(4), (b)(5)	(b)(4), (b)(5)
16	(b)(4), (b)(5)	(b)(4), (b)(5)
17	(b)(4), (b)(5)	(b)(4), (b)(5)
Min		
Max		
N		
Average		
Std. Dev.		

34

Chart #3
Bolus Accessory
Bolus Volume Accuracy
Life Test


(b)(4)

(b)(4)

35

Chart #4
Bolus Accessory
Bolus Volume Accuracy
Life Test


(b)(4)

(b)(4)

36

Chart #5
Bolus Accessory
Bolus Refill Time Accuracy
Life Test


(b)(4)

(b)(4)

37

Chart #6a
Paragon Bolus Accessory Set
Pressure Profile
During Refill Cycle
3.6 min. lockout

304705


(b)(4)

(b)(4)

58

Chart #6b
Baxter Patient Control Module
Pressure Profile
During Refill Cycle

(b)(4)

(b)(4)

BT

5.0 CHEMICAL AND DRUG SPECIFICATIONS

5.1 Compatibility

- 5.1.1 There are no specific drugs referenced in the labeling for the Bolus Accessory.
- 5.1.2 The Bolus Accessory is intended for general purpose drugs and pain medication.

5.2 Drug Stability

- 5.2.1 There are no specific drugs referenced in the labeling for the Bolus Accessory.

6.0 INTENDED USE

- 6.1 The Bolus Accessory is intended to deliver fixed boluses of medication upon demand by the patient or healthcare provider.
- 6.2 The routes of administration intravenous, epidural, intramuscular and subcutaneous.
- 6.3 The Bolus Accessory is not intended for continuous delivery.
- 6.4 The Bolus Accessory is single patient use only.
- 6.5 The Bolus Accessory is suitable for use as an ambulatory device and is intended for use in the hospital, home environment or alternative care sites.
- 6.6 No testing has been conducted to determine the efficacy of Bolus Accessory for the delivery of blood, blood products, lipids or fat emulsions. The Bolus Accessory is not intended for the delivery of blood, blood products, lipids or fat emulsions.

7.0 LABELS AND LABELING

- 7.1 I-Flow Corporation believes the proposed labels and labeling, where appropriate, meets the requirements of 21 CFR Part 801 as it relates to a determination of intended use and adequate directions for use.
- 7.2 The Bolus Accessory Directions for Use labeling:
 - 7.2.1 Provides comprehensive directions for preparation and use for the Bolus Accessory.
 - 7.2.2 Describes the routes of administration as it relates to intended use.
 - 7.2.3 Contains warning information.
 - 7.2.4 Contains the prescription statement required under 801.109 (b)(1).
 - 7.2.5 Includes the specifications of the Bolus Accessory. The specifications include the priming volume, residual volume, accuracy and operating conditions.
- 7.3 Identification labels and labeling
 - 7.3.1 I-Flow has developed product identification labeling for the Bolus Accessory. Refer to Appendix B for examples.
- 7.4 Packaging labels
 - 7.4.1 Contains the prescription statement required under 801.109 (b)(1).

40

7.5 Appendix C contains predicate labeling for the Patient Control Module marketed by Baxter Healthcare Corporation.

8.0 STANDARDS

8.1 There are currently no standards established for mechanical PCA infusion devices.

9.0 PACKAGING

9.1 The Bolus Accessory is packaged in a sealed Tyvek pouch, 12 pouches per case.

9.2 Packaging is suitable for radiation or ETO sterilization.

9.3 Package aging tests have been conducted on the Tyvek pouch. The results of (b) (4) testing has determined that the Tyvek pouches used to package the disposable Bolus Accessory maintains sterility in excess of three years.

10.0 STERILIZATION

10.1 The method of sterilization is gamma radiation (cobalt 60).

10.2 Sterilization validation methodology is by ANSI/AAMI ST32-1991 / EN552 Method 1 for gamma radiation.

10.2.1 (b) (4)

10.3 The sterile product under review here will have a sterilization assurance level (SAL) of 10^{-6} . Under AAMI Method 1 for Gamma sterilized product, no sterility test is required.

10.4 The Bolus Accessory is labeled pyrogen free and is tested for pyrogens using either the USP Rabbit Pyrogen Test or LAL test methods.

10.4.1 I-Flow products have been validated for LAL testing.

10.4.2 Either method may be used.

11.0 COMPARISON TO LEGALLY MARKETED DEVICES

See Table 1 that follows this section for more specific information.

11.1 Intended Use

11.1.1 The Bolus Accessory, the Baxter Patient Control Module and the I-Flow Bolus Dispenser have the same intended use:

11.1.1.1 To deliver fixed boluses of medication upon demand by the patient or healthcare provider.

11.2 Descriptions of the Paragon Bolus Accessory

Wherever the Paragon pump and administration set are mentioned, they may be replaced by any constant 6 psi pressure system.

11.2.1 The Bolus Accessory may connect to any Paragon administration set to deliver fixed boluses of medication upon demand by the patient or healthcare provider. The Bolus Accessory does not deliver continuous medication.

11.2.2 The Bolus Accessory consists of plastic housing, medication reservoir, bolus button activator and wrist bands.

41

11.2.3 Bolus Refill Cycle (i.e. Lockout Time)

11.2.3.1 The Paragon pump is the pressure source pushing fluid into the inlet port of the Bolus Accessory.

11.2.3.2 The flow control orifice of the Paragon administration set establishes a lockout period during which the patient cannot receive another full bolus of medication.

11.2.3.3 The time required to fill the bolus cavity is called the "lockout time" or "refill cycle".

11.2.3.4 Lockout times vary according to the flow rate of the Paragon administration set. See table below.

Flow Rate of Paragon Set (ml/hr)	Lockout Time (min.)
0.5	70.0
1.0	35.0
2.0	18.0
4.0	9.0
10	3.6

11.2.3.5 As fluid flows into the bolus cavity, the diaphragm is continuously displaced until the cavity is full.

11.2.3.5.1 The bolus cavity is filled non-linearly, slower in the beginning and more rapidly towards the end of the refill cycle.

11.2.3.5.2 The bolus reservoir is 0.5 ml volume.

11.2.4 Bolus Activation

11.2.4.1 On the opposite side of the diaphragm is a bolus button. The button moves with the diaphragm as it is displaced.

11.2.4.2 At any time while the bolus cavity is filling, the button can be activated by the user. Only when the bolus cavity is full will the patient receive the entire bolus.

11.2.4.3 When the button is activated, the fluid in the bolus cavity is pushed out and forced through a check valve. At this point, the refill cycle begins.

11.3 Comparison to the Baxter Patient Control Module.

11.3.1 The Baxter Patient Control Module is very similar to the I-Flow Bolus Accessory.

11.3.2 The Patient Control Module may connect to any Baxter Infusor to deliver fixed boluses of medication upon demand by the patient or healthcare provider.

11.3.3 The Patient Control Module consists of plastic housing, medication reservoir, bolus button activator and wrist bands.

11.3.4 Bolus Refill Cycle (i.e. Lockout Time)

42

- 11.3.4.1 The Baxter Infusor is the pressure source pushing fluid into the inlet port of the Patient Control Module.
- 11.3.4.2 The time required to fill the bolus cavity is called the “lockout time” or “refill cycle”.
- 11.3.4.3 Both the Patient Control Module and the Bolus Accessory determine the lockout time by the flow control orifice of the device they attach to. During the lockout time, the patient cannot receive another full bolus of medication.
- 11.3.4.4 Lockout times vary according to the flow rate of the Baxter Infusor. See table below.

Flow Rate of Baxter Infusor (ml/hr)	Lockout Time (min.)
0.5	60
2.0	15
5.0	6

- 11.3.4.4.1 The bolus cavity is filled linearly as opposed to non-linearly for the Bolus Accessory. See Charts #6a and #6b.
- 11.3.4.4.2 The bolus reservoir is the same as the Bolus Accessory (i.e. 0.5 ml volume).

11.3.5 Bolus Activation

- 11.3.5.1 For both the Baxter device and the Bolus Accessory, at any time while the bolus cavity is filling, the button can be activated by the user. Only when the bolus cavity is full will the patient receive the entire bolus.
- 11.3.5.2 When the button is activated, the fluid in the bolus cavity is pushed out and forced through a check valve. At this point, the refill cycle begins.

11.4 Comparison to the I-Flow Bolus Dispenser

- 11.4.1 The I-Flow Bolus Dispenser is a stand alone device that does not need to connect to an infusion pump such as the Paragon pump for the Bolus Accessory and the Baxter Infusor for the Patient Control Module.
- 11.4.2 The Bolus Dispenser consists of a reusable plastic dispenser, disposable medication reservoir (plastic syringe or mini-bag) and disposable administration set.

11.5 Specifications

- 11.5.1 The Bolus Accessory, the Baxter Patient Control Module and the I-Flow Bolus Dispenser have similar bolus volumes and lockout times. See Table 1.

11.6 Materials

- 11.6.1 All fluid path materials of the Bolus Accessory are in conformance with ISO 10993 Part 1.

- 11.7 Based upon the data presented in this section 11.0 and Table 1, I-Flow Corporation has determined that the Bolus Accessory is substantially equivalent to the named predicate devices. 43

Table 1
Comparison to Legally Marketed Devices

Comparison Element	Paragon Bolus Accessory Set (subject device)	SE ¹ Baxter Patient Control Module (K884505)	SE ¹ I-Flow Bolus Dispenser (K935811)
Intended Use	To deliver fixed boluses of medication upon demand by the patient or healthcare provider.	To deliver fixed boluses of medication upon demand by the patient or healthcare provider.	To deliver fixed boluses of medication upon demand by the patient or healthcare provider.
Route of Administration	Intravenous, epidural, intramuscular and subcutaneous	Intravenous, epidural, intra-arterial and subcutaneous	Intravenous, epidural, intramuscular and subcutaneous
Contraindications	Not intended for delivery of blood, blood products, lipids or fat emulsions.		Not intended for delivery of blood or blood products.
Reuse Capability	Disposable, Single Patient Only	Disposable, Single Patient Only	Disposable, Single Patient Only
Description			
Bolus Volume	0.5 ml	0.5 ml	0.25, 0.5 and 1.0 ml
Bolus Lockout Times	3.6, 9, 18, 35 and 70 min.	variable (6 to 60 min.)	15, 30, 60 and 120 min.
Pump Type	Available for use with any constant 6 psi pressure system such as the Paragon pump and administration set.	Available for use with Baxter Infusor pump	vacuum
Pump Volume	100 ml	65 ml	30 ml
Power Requirements	None	None	None
Pressure Source	Mechanical spring energy of the Paragon Pump	Strain energy of elastomeric membranes	vacuum
Fluid Reservoir	PVC drug bag	Elastomeric membranes	PVC drug bag or polypropylene syringe
Safety / Alarm Functions	The Bolus Accessory attaches to fixed flow rate tubing which prevents fluid runaway conditions. The non-linear refill adds additional patient safety if the bolus button is activated prior to the lockout time.	The Patient Control Module attaches to fixed flow rate tubing which prevents fluid runaway conditions.	

AF

¹SE = Substantially Equivalent

Comparison Element	Paragon Bolus Accessory Set (subject device)	SE ¹ Baxter Basal/Bolus Infusor (K984505)	SE ¹ I-Flow Bolus Dispenser (K935811)
Non-fluid Path Components	(b)(4)	(b)(4)	(b)(4)
Outer Bolus Housing	(b)(4)	(b)(4)	(b)(4)
Bolus Button	(b)(4)	(b)(4)	(b)(4)
Luer Cap	(b)(4)	(b)(4)	(b)(4)
Pinch Clamp	(b)(4)	(b)(4)	(b)(4)
Fluid Path Components	(b)(4)	(b)(4)	(b)(4)
Inner Housing	(b)(4)	(b)(4)	(b)(4)
Diaphragm (bolus reservoir)	(b)(4)	(b)(4)	(b)(4)
Check Valve	(b)(4)	(b)(4)	(b)(4)
Luer Adapters	(b)(4)	(b)(4)	(b)(4)
Tubing (make-up)	(b)(4)	(b)(4)	(b)(4)
Packaging (sterile pouch)	(b)(4)	(b)(4)	(b)(4)
Sterilization	Gamma 80 FPA	ETO 80 FPA	Gamma or ETO 80 FPA
Product Code			

45

¹SE = Substantially Equivalent

Appendix A
Paragon Bolus Accessory Set
Drawings

44

Appendix B
Paragon Bolus Accessory Set
Labeling

50

PARAGON BOLUS ACCESSORY SET

Directions for Use

NOMENCLATURE

1. Bolus Accessory Set ①
2. Luer Lock ②
3. PARAGON Administration Set ③
4. Bolus Button ④

INTENDED USE

The Paragon Bolus Accessory Set, in combination with a Paragon Administration Set, is intended to allow patient controlled bolus delivery. The bolus component of the administration set enables fixed boluses to be delivered upon demand by the patient or healthcare provider. The routes of administration include intravenous, epidural, intramuscular and subcutaneous.

DO NOT USE IF PACKAGE HAS BEEN OPENED OR IS DAMAGED OR IF EITHER PROTECTOR CAP IS NOT IN PLACE. THE PARAGON BOLUS ACCESSORY SET IS STERILE AND NON-PYROGENIC.


DO NOT RESTERILIZE. SINGLE PATIENT USE.

CAUTION

1. Not for blood, blood products, lipids or fat emulsions delivery. It is recommended that the Bolus Accessory Set be changed in accordance with established guidelines.
2. Epidural administration of analgesics is limited to use with indwelling catheters specifically indicated for either short-term or long-term analgesic drug delivery. Do not use Y adapter with epidural delivery.
3. If the device is to be used for epidural analgesic drug administration, it should be labeled to differentiate from other routes of administration. When using this device for epidural drug administration, make certain only drugs recommended for this route of administration are used.
4. The Bolus Accessory Set does not provide basal flow rate. It should be used in conjunction with another infusion line providing a continuous delivery in order to keep the IV line patent.
5. Refer to the Directions for Use for the Paragon Administration Set for additional instructions.

THE BOLUS ACCESSORY SET

The Paragon Bolus Accessory Set is designed to allow the patient or healthcare provider to administer a 0.5 ml bolus on demand. The refill time for the bolus depends on the Paragon Administration Set to which the Bolus Accessory Set is attached. See the Refill Time Table for specific refill times.

PRIMING THE ADMINISTRATION SET

1. Prime the administration set first.
2. Attach the Bolus Accessory Set to the administration set.
3. Using appropriate aseptic technique, remove the cap from the Luer lock at the end of the set. Open the clamp on the administration set tubing. The medication will flow into the bolus chamber on the Bolus Accessory Set.
4. When the bolus button fills (button rises to top of housing), press the Bolus button to prime. Repeat until air is purged down stream from the bolus housing, which may take several pushes.
5. Pinch the clamp closed and replace the cap.

STARTING THE INFUSION

1. Attach the Bolus Accessory Set to the infusion site. Secure the bolus housing against the skin, using either the strap provided or tape.
2. The Bolus button may be attached to the patient's wrist using the strap provided.
3. Opening the clamp on the administration set. A bolus can be delivered immediately.

BOLUS ACTIVATION

1. To receive a bolus of medication, press firmly down on the bolus button until it stops. ④
2. At any time during the infusion, the bolus button can be pushed to deliver a bolus of medication.
3. The next full bolus will be available after the refill time noted in the Refill Time Table.
4. Pressing the bolus button prior to the end of the refill time will result in a partial dose.

IMPORTANT

1. Only PARAGON administration sets distributed by I-Flow Corporation are authorized for use with this product. I-Flow Corporation accepts no responsibility for performance, or the liability for damages, caused by the misuse of this product when used with unauthorized administration sets.
2. This product uses DEHP plasticized PVC. Certain solutions may be incompatible with the PVC material used in the Bolus Accessory Set. Consult the drug package insert and other available sources of information for a more thorough understanding of possible incompatibility problems.

51

The PARAGON Bolus Accessory Set Specifications

Bolus Volume: 0.5 ml

Refill Time: 3.6 min. - 70 min., depending on the attached Paragon Administration Set.

Accuracy:

Bolus Volume: $\pm 10\%$ at 95% confidence interval at the identified lockout time.

Priming volume: Allow 1 ml for loss during priming.

Residual volume: Approximately 4 ml (including the Paragon Administration Set)

<u>Paragon Administration Set</u>	<u>Refill Time</u>
0.5 ml/hr	70 min.
1.0 ml/hr	35 min.
2.0 ml/hr	18 min.
4.0 ml/hr	9 min.
10.0 ml/hr	3.6 min.

NOTES

- Actual refill times may vary from the specified range due to:
 - viscosity and/or drug concentration.
 - temperatures above or below the operating conditions.
 - the positioning of the *PARAGON* Infuser above or below the infusion site.

52

CAUTION: Federal (U.S.A.) law restricts this device to sale by or on the order of a healthcare professional.

For Customer Service
Call: 1.800.448.3569
949.206.2700

European Representative:
MPS Medical Product Service GmgH
Borngasse 20, 35619 Braunfels, Germany

A PRODUCT OF

I-FLOW CORPORATION
LAKE FOREST, CA 92630
U.S.A.

1302268A
12/98

A PRODUCT OF / EINE PRODUKT VON / UN PRODUIT DE / UN PRODUCTO DE CONTENTS / INHALT /
CONTENU / CONTENIDO: 1

I-FLOW CORPORATION, LAKE FOREST, CA U.S.A.

REF B000000

Paragon Bolus Accessory Set

0.5 ml bolus volume

STERILE

LOT

SEE DIRECTIONS FOR USE.
CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A HEALTHCARE PROFESSIONAL.

Manufactured by / Hersteller von /
Fabrique par / Fabricado por:
I-Flow Corporation
Lake Forest, CA 92630 U.S.A.

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REF B000000

PART NO. 500XXXX

I-FLOW CORPORATION, LAKE FOREST, CA U.S.A.

Paragon Bolus Accessory Set

0.5 ml bolus volume

LOT

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A HEALTHCARE PROFESSIONAL. SEE DIRECTIONS FOR USE.

Manufactured by / Hersteller von /
Fabrique par / Fabricado por:
I-Flow Corporation
Lake Forest, CA 92630 U.S.A.

European Representative / Europäische Vertretung /
Représentant pour l'Europe / Representante Europeo:
MPS Medical Product Service GmbH
Borngasse 20, 35619 Braunfels, Germany

130XXXXA

54

Appendix C
Predicate Labeling

55

Infusor Patient Control Module

Description

The Patient Control Module is an accessory for Baxter's Infusors. When the medication demand button of the module is depressed, a small quantity of drug solution is delivered, nominally 0.5 mL. There is no infusion between depressions of the medication demand button.

Depression of the medication demand button empties the 0.5 mL reservoir of the Module. Upon release of the button, the reservoir again begins to fill with medication. The time to fill the Module's reservoir is determined by the Infusor selected.

For example:

Infusor	Flow Rate	Approximate Module Fill Time
2C1070/2C1071	2 mL/hr	15 minutes
2C1073	5 mL/hr	6 minutes
2C1080	0.5 mL/hr	60 minutes
2C1954	0.5 mL/hr	60 minutes
2C1955	0.5 mL/hr	15 minutes

Note: All times reflect use of the Infusor as directed in the individual direction sheet.

Indications and Usage

The Patient Control Module is intended for use with patients in both the hospital and home environment. It is convenient for use by ambulatory or nonambulatory patients. **This accessory can only be used with Baxter's Infusors.**

Warnings

- To ensure proper performance the device should be used as directed.
- This accessory **must** be used with Baxter's Infusors to ensure safe operation.
- Do not use if white shipping insert is not affixed to the medication demand button when received.
- Remove white shipping insert prior to connection to patient. **Failure to remove shipping insert will cause continuous infusion.**

Not intended for use for continuous infusion.

Directions For Filling

Use aseptic technique throughout entire procedure

A. Standard Method

- Select the Infusor with the desired flow rate. Module filling time determines the maximum number of 0.5 mL doses available per hour.
- Prepare the Infusor per its enclosed direction sheet. The drug concentration should be set so that the unit delivers the maximum dose/hr prescribed by the physician at either the 6, 15, or 60 minute dose interval. The patient can then use the Module to titrate the drug within the prescribed limits.

- When the Infusor is ready for use, remove the tip protector and winged Luer cap from the Module. Do not remove white shipping insert at this time. Place the male Luer of Infusor into the female Luer of the Module. **Note: Avoid touch contamination of the Luer fittings.**
- The Module will be primed by the flow of the Infusor. Priming time will depend upon the flow rate of the Infusor selected.
- Observe the reservoir filling through the clear back plate. **When the liquid starts to flow from the male Luer, remove the white shipping insert. Failure to remove shipping insert will cause continuous infusion.** When the reservoir is filled, depress and release the medication demand button to complete the priming operation and expel any air in the reservoir.
- Replace the winged Luer cap onto the male Luer of the Module.
- Package the entire unit in the Infusor Dispenser Bag for transport to the patient.
- Store in a clean area, protected from sunlight. Refer to drug manufacturer's package insert for storage requirements.
- Optional Method
 - Proceed per standard method above, sections A.1 and A.2.
 - When the Infusor is ready for use, remove the tip protector and winged Luer cap from the Module. **Remove the white shipping insert from the Patient Control Module. Failure to remove shipping insert will cause continuous infusion.** Using a syringe containing a minimum of 2 mL of solution or diluent, attach the syringe to the female Luer of the Module. While filling the Module reservoir, cycle the Module medication demand button as needed to eliminate air in the reservoir or tubing. Disconnect the syringe from the female Luer of the Module.
 - Connect the male Luer of Infusor to the female Luer of the Module.
 - Proceed per standard method above, sections A.6 through A.8. **Do not store above 38°C (100°F) or in direct sunlight.**

Baxter Healthcare Corporation

Deerfield, IL 60015 USA

Printed in USA


©Copyright 1985, 1988, 1992, Baxter Healthcare Corporation. All rights reserved.

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Warning: Failure to remove shipping insert will cause continuous infusion.

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Infusor Patient Control Module

Sterile, nonpyrogenic fluid path.

Do not use if tip protectors are not in place.

The Patient Control Module is intended for use in both the hospital and home environment by ambulatory or nonambulatory patients.

For use with Baxter's infusors only.

Not intended for use for continuous infusion.

Failure to remove shipping insert will cause continuous infusion.

See accompanying directions for use.

Caution: Federal (USA) law restricts this device to sale by or on order of a physician.

Do not store above 38°C (100°F) or in direct sunlight.

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
L016311R/
04-11-93/
11-90/20-06-90

56

PORTABLE INFUSION SYSTEMS WITH ELASTOMERIC TECHNOLOGY

Model	Number	Flow Rate	Lockout Time	Bolus Dose	Capacity	Units Per Case
Singleday Infusor	2C1071	2 mL/hr	15 min*	0.5 mL*	65 mL	6
12-hour Infusor	2C1073	5 mL/hr	6 min*	0.5 mL*	65 mL	6
Multiday Infusor	2C1080	0.5 mL/hr	-	-	65 mL	6
Basal-Bolus 15	2C1955	0.5 mL/hr	15 min	0.5 mL	65 mL	6
Patient Control Module	2C1079	-	-	0.5 mL	0.5 mL	12

*When used with Patient Control Module.

1. Carlson RW, Sikic BI. Continuous infusion or bolus injection in cancer chemotherapy. *Ann Intern Med.* 1983;99:823-833.

Baxter

Baxter Healthcare Corporation, I.V. Systems Division, 1425 Lake Cook Road, Deerfield, IL 60015

59
/

U.S. Food and Drug Administration

DETAILED INFORMATION

Device Classification Name: PUMP, INFUSION

Regulation Number: 880.5725

510(k) Number: K884505

Device Name: BASAL/BOLUS INFUSORS

Applicant: BAXTER HEALTHCARE CORP.

ROUTE 120 & WILSON ROAD

ROUND LAKE, IL 60073

Contact: PATRICIA S BARSANTI

Product Code: FRN

Date Received: 10/25/88

Decision Date: 07/14/89

Decision: Substantially Equivalent

Classification Advisory Committee: General Hospital

Review Advisory Committee: General Hospital

Statement/Summary/Purged Indicator: Purged, no summary or statement

[Return For Another Search](#)

[CDRH Home Page](#)

[FDA HOME PAGE](#)

[COMMENTS](#)

(Database Updated April 6, 1998)

Appendix D
Reference Labeling
Paragon Infusion System

Jc
62

PARAGON USER'S GUIDE

Directions for Use

NOMENCLATURE

- 1. PARAGON Infuser
- 2. Fluid Level Indicator
- 3. 1.2 micron air-eliminating filter
- 4. PARAGON Administration Set
- 5. Flow Rate Label
- 6. Luer Lock
- 7. End of Infusion Indicator

CAUTION

- Do not use the administration set if the sterile pouch is opened or damaged. If either protective cap is missing or not in place, the sterility of the administration sets is no longer guaranteed.
- Not for blood or blood products delivery. It is recommended that the administration set be changed every 24-48 hours or in accordance with CDC guidelines or institutional policies.
- Do not restaink administration set. Administration sets are intended for single patient use only. The fluid pathway is sterile and nonpyrogenic.

INTRODUCTION

The PARAGON is a drug delivery system consisting of a reusable mechanical infuser and specially designed administration sets. The PARAGON provides precise delivery of medications requiring slow and continuous infusions, such as chemotherapeutics and analgesics. The PARAGON also infuses medications which require faster delivery, such as antibiotics.

THE PARAGON ADMINISTRATION SET

Administration sets are made of PVC. Each set is approximately 127 cm long. A 1.2 micron air-eliminating filter is built into all administration sets. The flow rate at which the drug flows to the patient is controlled by a flow restrictor built into the end of the set. Flow rates for each set are printed on the air-eliminating filter label.

FILLING THE PARAGON IV BAG - USE ASEPTIC TECHNIQUE

- Remove the IV bag with attached administration set from its package.
 - Move the flow clamp next to the filling valve and close the clamp.
 - Fill a sterile syringe with the solution to be dispensed into the IV bag.
 - Connect the tip of the syringe to the filling valve and inject the solution into the IV bag. Refill the syringe and repeat if necessary.
- NOTE:** The PARAGON Infuser is designed to hold a total of 100 ml of fluid. The maximum fill volume is 110 ml. If the amount of fluid exceeds 110 ml, it may be difficult to engage the threads on the top and bottom of the PARAGON infuser.
- Remove air from the IV bag by aspirating with a syringe attached to the filling valve. Squeezing the sides of the IV bag when pulling back on the syringe will aid in removing the air.
 - Be certain to replace the cap on the filling valve.
 - Do not place labels on the IV bag. Labels may be wrapped around the set.

LOADING THE IV BAG INTO THE PARAGON INFUSER

- Twist open the top and bottom halves of the PARAGON Infuser.
- Before placing the IV bag into the PARAGON Infuser, slide the thin portion of the administration set through the slot found on the bottom of the infuser.
- Center the bag in the bottom and press all around the edge of the bag to fully seat the bag in the bottom. Make sure there are no wrinkles in the bag.
- Pull gently on the thick portion of the tubing so that it is fully extended and seated at the bottom of the slot.
- Twist the top and bottom halves of the PARAGON Infuser together until they meet.

PRIMING THE ADMINISTRATION SET

- Using aseptic technique, remove the cap from the luer lock at the end of the set. Open the clamp on the IV tubing. The medication will flow toward the end of the luer lock.
- Confirm that fluid is flowing by observing the formation of a drop at the end of the luer lock. It may take 10 minutes for a drop to form when priming the 0.5 ml/hr set.
- Finish the clamp closed and replace the cap.

STARTING THE INFUSION

- Attach the administration set to the IV site. Secure the connection against the skin.
- Start the infusion by opening the clamp on the administration set. The infusion will begin immediately.

THE FLUID LEVEL INDICATOR

- The window with the markings on the side of the infuser is used to estimate how far the infusion has progressed.
- When the PARAGON IV bag is filled to its capacity of 100-110 ml, the top of the pressure plate will be aligned with the top round marker.
- As the infusion progresses, the plate will move to the bottom marker indicating the bag is nearly empty.

THE END OF THE INFUSION

The infusion is complete when at least three (out of the six) small blue dots appear through the bottom of the PARAGON Infuser.

THE CARRYING CASE

The carrying case can be worn on a belt, over the shoulder, or around the waist.

- Place the PARAGON infuser in the carrying case so that the bottom of the infuser can be seen through the clear plastic window.
- Lift the Velcro strap and slide the administration set down so that the set exits the carrying case at the side window opening. Close the strap. (Positioning the infuser in this way allows for the viewing of the Fluid Level Indicator.)
- The front flap of the carrying case lifts up to reveal a clear plastic window, allowing for the viewing of the End of Infusion Indicator.
- If necessary, a small lock can be placed through the larger of the two holes on the zipper, and then through the cloth loop on the side of the carrying case. (This may discourage tampering with the infuser during an infusion.)

CARE OF THE PARAGON

The PARAGON Infuser is durable and is intended to be used for repeated drug deliveries. After each patient use, the exposed surfaces, except the threads, may be wiped clean using isopropyl alcohol or a 10% bleach solution.

NOTE: Do not submerge the PARAGON Infuser in a bleach solution. After cleaning, if the PARAGON is difficult to twist together, place a small drop of lubricating ointment (such as K-Y™ Jelly) on a small section of the threads on the bottom of the infuser. Twist the top of the infuser onto the bottom to spread out the ointment.

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Handwritten marks and scribbles at the bottom right of the page.

IMPORTANT

1. Only administration sets distributed by I-Flow Corporation are authorized for use with this product. I-Flow Corporation accepts no responsibility for performance, or the liability for damages, caused by the misuse of this product when used with unauthorized administration sets.
2. This product uses DEHP plasticized PVC. Certain solutions may be incompatible with the PVC material used in the IV administration set. Consult the drug package insert and other available sources of information for a more thorough understanding of possible incompatibility problems.
3. This device contains natural rubber latex. Individuals with known natural rubber latex sensitivities should not use this product.

The PARAGON Infuser Specifications

Size: 5.8 cm high; 10.2 cm in diameter
 Weight: 260 Gms
 Flow Rates: 0.5, 1, 2, 4, and 10 ml/hr
 Delivery accuracy: Accuracy is at $\pm 10\%$ at 95% confidence interval.
 Priming volume: 1.5 ml
 Residual volume: 5 ml or less

NOTES

1. The infusion rates for each administration set are indicated on the administration set label.
2. Actual infusion rates may vary from the specified range due to:
 - viscosity and/or drug concentration.
 - temperatures above or below the operating conditions.
 - the positioning of the PARAGON Infuser above or below the IV site.
3. The Paragon Drug Delivery System has been calibrated using Normal Saline (NS) as the diluent and skin contact temperature (32°C, 90°F) as the operating environment. When using NS and skin temperature the Paragon System will flow at the specified nominal rate. The use of other diluents or operating temperatures other than the above will affect the nominal flow rate. For example, if 5% Dextrose (D5W) is used as the final diluent, the Paragon System will flow at 10% below the nominal rate due to higher solution viscosity.

CAUTION: FEDERAL LAW (U.S.A.) RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A HEALTHCARE PROFESSIONAL.

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Manufactured by
 I-Flow Corp
 Lake Forest, CA 92630
 Printed in the U.S.A.

U.S. and Foreign Patents Pending

European Representative
 MPS Medical Product Service GmbH
 Borgasse 20, 35619 Braunes, Germany

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