

1 Synergo system consists of two separate components,  
2 the Synergo hyperthermia device and the Synergo kit.  
3 The Synergo hyperthermia device is a computer-  
4 controlled device that monitors and controls  
5 treatment and contains the user interface. The  
6 Synergo kit contains single use components required  
7 to perform each Synergo treatment, namely the  
8 disposable catheter tubing set and two vials of the  
9 mitomycin C, both of which are packaged sterile. The  
10 catheter tubing set consists of a specialized urinary  
11 catheter used to deliver treatment and tubing to  
12 circulate the drug solution in and out of the bladder  
13 during treatment.

14           The vials of mitomycin C are provided in  
15 the clinical dosage specified in the Synergo system  
16 labeling. Mitomycin C is provided to Medical  
17 Enterprises in final package form from the drug  
18 manufacturer, Bedford Laboratories. Mitomycin C is  
19 legally marketed in the United States for a different  
20 indication and route of administration and is sold by  
21 Bedford Laboratories under the trade name mitomycin  
22 for Injection, USP.

23           By copackaging the drug with the device  
24 disposable components, the firm is meeting FDA's  
25 requirements for mutually confirming labeling.

1 Specifically, the labeling insert package with the  
2 Synergo kit is intended to provide the user with  
3 adequate directions for use of both the device and  
4 drug components.

5           The sponsor proposes that the Synergo  
6 system be indicated to deliver heat transurethrally  
7 by means of radio frequency energy to the urinary  
8 bladder walls for the treatment of superficial  
9 transitional cell carcinoma of the bladder or STCCB,  
10 concomitant with intravesical installing of Mitomycin  
11 for Injection, USP. The combination of Synergo and  
12 mitomycin C is intended for the prophylactic  
13 treatment of recurrence in patients following  
14 endoscopic removal of Ta and T1 and grades 1 through  
15 3 STCCB. Synergo and mitomycin C treatment is  
16 clinically indicated for STCCB patients of  
17 intermediate and high-risk.

18           FDA has conducted a very comprehensive  
19 review of this PMA. Since this is a device/drug  
20 combination product, our review has spanned to two  
21 centers, the Center for Devices and Radiological  
22 Health and the Center for Drug Evaluation and  
23 Research. Additionally, the Office of Regional  
24 Operations was involved to perform bioresearch  
25 monitoring inspections which are FDA's clinical data

1 audits.

2 I'd like to take a few moments to  
3 acknowledge all the individuals that continued to the  
4 review of this PMA. In addition to myself as Team  
5 Leader, this slide lists the clinical reviewers as  
6 well as Janine Morris, who until recently was Chief  
7 of the Urology and Lithotripsy Devices Branch.

8 The reviewers listed on this slide reviewed  
9 the statistical information, pharmacokinetic data and  
10 post-approval study plan.

11 These individuals reviewed the preclinical  
12 testing information and patient labeling.

13 And these individuals, they reviewed the  
14 sterilization information and manufacturing process  
15 or conducted or oversaw the bioresearch monitoring  
16 inspections.

17 The principle of operation of this  
18 combination product has already been presented in  
19 detail by the sponsor. Briefly, Synergo treatment  
20 involves hyperthermia of the internal urinary bladder  
21 wall with circulation of an intravesical dose of the  
22 chemotherapy drug, mitomycin C. The hyperthermia  
23 device heats the bladder wall using radio frequency  
24 energy, specifically in the microwave range,  
25 delivered from a catheter-based antenna.

1           In addition to the antenna, the specialized  
2 transurethral catheter contains five thermocouples,  
3 three in the patient's bladder and two along the  
4 catheter shaft to provide temperature feedback to the  
5 system, as well as internal lumens for the continuous  
6 circulation of the drug solution between the bladder  
7 and an external heat exchanger. The treating  
8 physician adjusts both the radio frequency power  
9 output and drug circulation flow rate to maintain the  
10 temperature of the bladder wall in the range of 42  
11 plus or minus 2 degrees Celsius.

12           Treatment is initiated following  
13 transurethral resection of all tumors. A complete  
14 course of treatment consists of eight weekly  
15 inductive sessions followed by four monthly  
16 maintenance sessions. Each treatment session lasts  
17 60 minutes. AS previously stated, the physician  
18 controls the hyperthermia system to maintain a  
19 temperature of 42 plus or minute 2 degrees Celsius at  
20 the bladder wall. Prior to the session, the  
21 patient's bladder is drained, local anesthesia is  
22 given along the urethra, and the treatment catheter  
23 is inserted into the bladder. Mitomycin C is  
24 administered in 2 consecutive instillations of 20  
25 milligrams in 50 mL distilled water, each for a dwell

1 time of 30 minutes. At the end of the second 30-  
2 minute period, hyperthermia is terminated, the  
3 patient's bladder is emptied, and the catheter is  
4 withdrawn.

5 This slide lists the preclinical tests that  
6 were performed to evaluate the safety and functional  
7 performance of the Synergo device. The sterilization  
8 and manufacturing information was reviewed from both  
9 a chemistry, manufacturing and controls perspective  
10 in the Center for Drug Evaluation and Research, and a  
11 good manufacturing perspective in the Center for  
12 Devices and Radiological Health.

13 FDA has no major concerns regarding the  
14 preclinical data provided. As noted in the FDA  
15 Executive Summary, minor issues remain regarding the  
16 claimed shelf life of the Synergo kit.

17 There are several sources of clinical data  
18 provided for review in the PMA, all of which were  
19 collected at sites in Europe and Israel. The table  
20 presented here lists the data sources that were  
21 collected in prospective clinical studies as well as  
22 the supportive clinical data. I would like to bring  
23 your attention to Study 101.1, highlighted in yellow  
24 in this table. As you heard, this is a randomized  
25 clinical trial intended as a pivotal study dataset in

1 support of the safety and effectiveness of the  
2 Synergo system. In this study, eligible subjects  
3 were randomized to either Synergo treatment  
4 consisting of hyperthermia plus mitomycin C or  
5 mitomycin C alone. The main effectiveness analysis  
6 for this study is comparison of the two-year rate of  
7 recurrence of bladder cancer between study arms.  
8 Safety is evaluated by analyzing the incidence and  
9 severity of adverse events.

10           The remaining sources of clinical data were  
11 submitted to provide supportive data regarding the  
12 safety and/or effectiveness of Synergo treatment.  
13 For various reasons, these studies are not intended  
14 to provide pivotal evidence of the safety and  
15 effectiveness of this combination product.  
16 Specifically, the human pharmacokinetic study was  
17 conducted to evaluate the systemic levels of  
18 mitomycin C during Synergo treatment. Study 102.1  
19 reports the results of an unplanned interim analysis  
20 of an ongoing study comparing Synergo treatment to  
21 Bacillus Calmette-Guerin or BCG. And, Study 101.4 is  
22 an exploratory study performed to assess the use of  
23 the Synergo system in a different patient than  
24 specified in the PMA's proposed indication.

25           The other clinical data sources provided in

1 the PMA, and listed on the law row are comparisons to  
2 historical mitomycin C and BCG data reported in the  
3 literature, commercial data collected on the Synergo  
4 system and reported as the European Prophylactic  
5 Study and data from Bladder Salvage and European  
6 Ablation Patients, collected at European sites to  
7 evaluate different indications than posed in the PMA.

8           FDA's review comments regarding this  
9 clinical information will be presented later by Drs.  
10 Kane and Li.

11           To understand how the pivotal study data  
12 and patients were managed, it is important that I  
13 briefly discuss the chronology and progress of study  
14 events.

15           Study 101.1 did not originate as a planned  
16 study to support a marketing application in the U.S.  
17 The protocol for the study was originated in 1993.  
18 In 1997, the study sponsorship was transferred from  
19 one of the study sites to Medical Enterprises or MEL.  
20 That same year, MEL's investors contracted with SMS,  
21 a British firm, to create case report forms and  
22 conduct monitoring visits. During the monitoring  
23 visits conducted in 1997, the case report forms were  
24 retrospectively completed by transcribing existing  
25 data from the patients' medical records. After the

1 1997 monitoring visits, the case report forms were  
2 completed prospectively.

3           Following study completion and preliminary  
4 discussions with FDA, MEL submitted this PMA in  
5 August of 2001. FDA conducted its in-depth review of  
6 the PMA in the summer and fall of 2001, and  
7 additional in-depth review cycles occurred in 2002  
8 and 2004 following the firm's submission of  
9 additional information. Due to the unconventional  
10 way in which the case report forms were developed,  
11 FDA determined during the 2004 review cycle that  
12 additional review of the PMA should not proceed until  
13 the clinical data were audited and bioresearch  
14 monitoring inspections.

15           Based on the firm's availability, FDA  
16 conducted these inspections of all sites in late  
17 2005. These inspections verified that the data  
18 contained in the case report forms and PMA data  
19 listings are accurate reflections of the  
20 investigator's source documents. However,  
21 limitations to the use of retrospectively created and  
22 transcribed case report forms include the potential  
23 for relevant, pre-1997 study data not being collected  
24 in the patient records and the potential for recall  
25 bias.

1           This slide outlines the remainder of FDA's  
2 presentation. I would now like to introduce  
3 Dr. Herrera, who will give a clinical overview of  
4 bladder cancer and current treatment options.

5           DR. HERRERA: Good morning. I'm Hector  
6 Herrera, the urologist for the Urology Devices  
7 Branch. I'm going to talk about bladder cancer. All  
8 the Panel members are experts in this subject. So  
9 this presentation is given for the benefit of the  
10 public not familiar with bladder cancer.

11           More than 60,000 cases annually are  
12 diagnosed in the United States. The main causes is  
13 tobacco use, industrial carcinogens and in general,  
14 the aging of the population.

15           The most common pathological subtype is the  
16 transitional cell carcinoma or TCC. It's observed in  
17 over 90 percent of the tumors. Squamous cell  
18 carcinoma is 5 percent and adenocarcinomas  
19 approximately 1 percent.

20           The prognosis factors for the bladder  
21 cancer are two, the stage and grade. Stage is how  
22 far the disease has spread, and the grade is the  
23 appearance of these cells to the microscopic  
24 examination.

25           The tumor stages of the non-muscle invasive

1 bladder cancer is Ta, T1 and Tis or CIS. The Ta is  
2 invading the urothelium, the T1 invading the  
3 urothelium and the lamina propria and Tis is confined  
4 to the mucosa but this tumor is highly malignant and  
5 aggressive and as I say before is called also CIS.

6 We have the representation of the bladder  
7 wall in this slide. Here we have the urine in  
8 contact with the mucosa, then urothelium followed by  
9 the lamina propria, the muscle cells, the fat and the  
10 peritoneum. We have here what is described as TIS,  
11 is a very superficial tumor confined only to the  
12 mucosa but is very aggressive. The Ta invades the  
13 urothelium and the T1 invades the mucosa, urothelium  
14 and doesn't go into the three layers of muscles that  
15 involve the bladder.

16 The tumor grade as I said before is the  
17 appearance of the cells to the microscope and is  
18 divided in well differentiated, moderately  
19 differentiated and poorly differentiated.

20 We can see an example of what is a well  
21 differentiated, we have here in the bottom, the  
22 normal cells of the bladder and on top is the  
23 increased number of cells but as you're noticing  
24 here, the architecture is almost the same.

25 In the moderately differentiated, we notice

1 that we have a lot of the architecture but we can see  
2 perfectly well each cell.

3 In the poorly differentiated, we cannot  
4 distinguish anymore the cells. That's why it's  
5 called poorly differentiated.

6 The treatment options, the main one is the  
7 transurethral resection of the bladder tumor. Then  
8 there is the intravesical chemotherapy or  
9 immunotherapy or intravesical laser ablation and  
10 photodynamic therapy.

11 Basically taken to the operating room and  
12 the bladder tumor is resected. This is called TURBT  
13 or transurethral resection. It's essential for the  
14 diagnostic of the tumor but you need to be deep  
15 enough to get a sampling of the first layer of  
16 muscles.

17 The most commonly employed intravesical  
18 agents in the USA is the Bacillus Calmette-Guerin or  
19 BCG and the mitomycin C.

20 The mechanism of action, BCG is  
21 inflammatory host response with release of special  
22 proteins, that they are called cytokines. The  
23 mitomycin C is antineoplastic, inhibits DNA  
24 synthesis. Both are used as adjuvant or maintenance  
25 fashion.

1           The mitomycin C usually is given between 20  
2 to 60 milligrams instilled in the bladder and  
3 dissolve in distilled water most commonly.

4           And the BCG is the first line of treatment  
5 for the Tis. It's effective as prophylaxis for  
6 recurrences but the optimal dosing has not yet been  
7 established.

8           There are side effects for both of them.  
9 The mitomycin C, the rate is 22 to 24 percent with  
10 multiple does with or without maintenance. The main  
11 side effects is urgency, frequency, urethral  
12 infection and bladder contracture.

13           The BCG, 38 percent with induction and 57  
14 percent with induction plus maintenance. Side  
15 effects, fever, chills, malaise, altered liver  
16 functions and sepsis.

17           The recurrence and progression rates after  
18 resection, the recurrences are decreased by 31  
19 percent with BCG maintenance and by 18 percent with  
20 mitomycin C. The progression rate estimate in all  
21 patient risk groups is 8 percent for BCG and 4  
22 percent for mitomycin C.

23           The predicting factor, they are divided in  
24 three, the low risk, the intermediate risk and high-  
25 risk, according to the European Urology Association

1 and is based on the number of tumors, the size of  
2 them, the prior recurrence, the T category, tumor  
3 grade and the presence of Tis. Thank you very much.

4 DR. KANE: Good morning. I'm Robert Kane.  
5 I'm a Medical Oncologist, and I'm a Clinical Reviewer  
6 in the Center for Drug Evaluation, CDER. I will be  
7 presenting the clinical review of the Synergo  
8 application this morning.

9 I will discuss superficial bladder cancer  
10 briefly to relate Dr. Herrera's discussion to the  
11 Synergo application. I will then discuss regulatory  
12 considerations for marketing. I will discuss certain  
13 of the Synergo studies, 101.1 in particular, and  
14 102.1, and conclude with some regulatory concerns.

15 For non-muscle invasive disease,  
16 superficial bladder cancer, as you have heard,  
17 primary therapy is cystoscopic complete removal.  
18 Pathologic staging and grading is extremely important  
19 to established prognosis and additional treatment  
20 options, which may include adjuvant drug therapy  
21 intravesically into the bladder and you have heard  
22 that BCG is FDA approved for this indication.  
23 Mitomycin C, Interferon and other agents have been  
24 employed.

25 You have also heard the population is

1 heterogeneous and thus outcomes vary. Recurrence is  
2 common. Progression is less common and not  
3 apparently altered by our therapies to date.  
4 Prognostic factors affect this, and these have been  
5 discussed as well. And thus, standard of care has  
6 been complete resection, observation, and repeat of  
7 the procedure, and recurrence after TUR alone has  
8 been estimated anywhere from 10 to 70 percent by 2  
9 years based on the characteristics of the patients  
10 examined. And, adjuvant treatment is indicated for  
11 the higher risk patients.

12 I will abbreviate some of the comments from  
13 these slides because you already know them and have  
14 heard them enough times I think.

15 With intravesical treatment with mitomycin  
16 C, either as immediate following cystoscopy and TUR  
17 or maintenance, there appears to be an overall  
18 suggestion of a reduction of about 18 percent in the  
19 probability of recurrence over the 2 year time  
20 interval. Again, there's considerable variation in  
21 this information.

22 While mitomycin C is not FDA approved for  
23 the treatment or prophylaxis of superficial bladder  
24 cancer, it has wide usage for this condition. And  
25 you've also heard already about the availability of

1 BCG under these conditions, and I think we would all  
2 agree that much additional improvement in the  
3 management of this condition would be very desirable.

4           The study endpoints that have usually been  
5 examined for this condition include recurrence and  
6 recurrence by the two-year time interval has commonly  
7 been used. Progression, indicating advancement of  
8 the T stage, as mentioned, is generally unusual  
9 unless there's high grade disease present.

10           Recurrence-free interval or time to  
11 recurrence, time to event endpoints may be utilized.  
12 Long-term outcomes, however, have not been  
13 systematically studied for this condition.

14           Turning now for a moment to consider  
15 requirements for marketing approval, we look for  
16 reasonable assurance based on valid scientific  
17 evidence demonstrating effectiveness with acceptable  
18 safety. And the source of this evidence should be  
19 adequate and controlled investigations. In addition,  
20 other types of studies can provide evidence as well.

21           Now, this slide indicates what we might  
22 classify as ideal characteristics for a single study  
23 to support effectiveness. It should be multicenter.  
24 It should be adequately powered. There should be  
25 consistent efficacy across multiple endpoints and key

1 study subsets. The statistical results should be  
2 very persuasive, and the study should be clinically  
3 compelling for the benefit to the risk assessment.

4 Turning to the Synergo application, this  
5 represents a device and drug combination product as  
6 has already been described. And the proposed  
7 indication is for use for the prophylactic treatment  
8 of recurrence following endoscopic removal, and  
9 patients should be considered to be at intermediate  
10 or high-risk for this to be applicable.

11 I will be commenting on the following four  
12 studies submitted by the sponsor.

13 With respect to Study 101.4, titled the  
14 tumor ablation study, patients enrolled in this study  
15 were those who were not fully resectable. It's a  
16 single-arm evaluation of the Synergo device plus  
17 mitomycin C.

18 Study 102.1 is a randomized comparison of  
19 Synergo plus mitomycin C versus BCG alone without the  
20 use of hyperthermia. Patients had to have definitive  
21 TURBT, and this study began in 2002. It consists of  
22 six weekly inductions and then six months of  
23 maintenance once a month. Four different types of  
24 BCG are being used in this study. The planned  
25 enrollment according to the applicant is 300, and

1 based on a data analysis performed in April 2007, at  
2 this time, we've been informed about 90 evaluable  
3 patients accrued over 5 years at 10 different sites.  
4 And the Synergo plus mitomycin C arm is of interest  
5 for safety.

6           The applicant has identified Study 101.1 as  
7 their pivotal study, a comparative study of  
8 intravesical mitomycin C instillation or mitomycin C  
9 and local hyperthermia for the prophylaxis of a  
10 recurrence of superficial transitional cell bladder  
11 carcinoma. As you have heard, the study was  
12 conducted at three sites. This began as a research  
13 collaboration among the sites. Enrollment between  
14 1994 and 1999, and during this time interval, 83  
15 patients were accrued and randomized. You've also  
16 heard case report forms were initiated in 1997.

17           And according to the plan, patients were to  
18 be randomized one to one between Synergo hyperthermia  
19 at 42 degrees plus mitomycin C which was preheated  
20 versus mitomycin C alone which was administered at  
21 room temperature, and 2 successive instillations for  
22 each patient were used. Patients received each over  
23 30 minutes projected dwell time, and thus a total  
24 exposure time was to be 60 minutes for each  
25 treatment. All patients were to receive eight weekly

1 cycles and then monthly times four. Cystoscopy was  
2 to be performed after the initial eight weeks and  
3 then quarterly.

4           And patients on the Synergo arm, since the  
5 catheter was retained in the bladder, had the bladder  
6 drained via the catheter after the one hour  
7 hyperthermia treatments, and in the protocol, the  
8 treatment could occur between a time interval of 40  
9 to 75 minutes. So it was considered reasonable.  
10 Patients on mitomycin C alone were instructed to void  
11 spontaneously after the total of one hour, and  
12 treatment was to being 20 to 40 days after the  
13 definitive resection, with quarterly follow-up to  
14 include cystoscopy, cytology, biopsy of any  
15 suspicious areas until two years or recurrence. If  
16 recurrence occurred, patients then went off study.

17           In the 1993 protocol, several endpoints  
18 were stated without development of the hierarchy or  
19 specific identification of what was to be the primary  
20 endpoint. And these endpoints included disease-free  
21 interval, the tumor recurrence rate, including  
22 evaluation of the number of recurrent tumors at two  
23 years post-treatment, progression of stage,  
24 progression of grade, occurrence of CIS, occurrence  
25 of distant metastases. In a submission by MEL in

1 2001, recurrence-free survival was indicated as the  
2 primary endpoint. In the MEL submission of January  
3 of this year, in preparation for this meeting, the  
4 primary endpoint was identified as time to  
5 recurrence.

6 Now, the inclusion criteria have also been  
7 discussed and I will not repeat these for you.

8 Looking at the demographics of the enrolled  
9 population, as noted, a total of 83 patients  
10 randomized, 41 and 42 respectively. There was a  
11 slight increase in the proportion of the patients of  
12 age 66 or greater in the control arm. Overall, for  
13 the first episode, 37 percent of patients presented  
14 in that fashion and 40 percent were stage Ta and the  
15 majority of patients were grade 2 or grade 3 in the  
16 enrollment.

17 Now, this slide indicates an overview of  
18 study results, and through the next several slides I  
19 will try to reconstruct as the sponsor has also  
20 presented the various analysis populations that have  
21 been considered. Of the total of 83 initially, 8  
22 patients in total discontinued before the first time  
23 cystoscopy follow-up that was to take place. A total  
24 of seven of these were on the Synergo plus mitomycin  
25 C arm, one on the control arm. Subtracting these 8

1 patients, the applicant has identified a population  
2 considered valid for analysis or protocol consisting  
3 of 75 patients, 35 and 40. Just to note, in summary,  
4 54 patients completed all therapy, that is the 8  
5 induction and 4 maintenance treatments, 65 percent in  
6 total, somewhat evenly distributed.

7           Now, I want to spend a moment talking about  
8 the applicant's description of these early  
9 discontinuations, a total of eight patients, as we  
10 mentioned seven on the Synergo arm, one on the  
11 control arm. The first six, the applicant has  
12 identified as early withdraws or dropouts and the  
13 final two are described as protocol deviations.  
14 These were also patients in whom consent was  
15 withdrawn for the reasons indicated, and after the  
16 patients were withdrawn, they received mitomycin C  
17 alone. So keep this total of eight patients in mind.  
18 I have combined them for the subsequent slides.

19           Now, going back to the per-protocol or  
20 valid for analysis population as identified by the  
21 sponsor, the total of 83, and removing these 8 early  
22 discontinuations, the applicant provided us their  
23 information according to the Kaplan-Meier curves of  
24 the estimate of the two-year recurrence rates, and  
25 that is shown on this slide, 17.1, 61.6, with a log-

1 rank P value test of significance. And again, this  
2 analysis group includes the total of the eight  
3 patients who discontinued before first follow-up.

4 Now, you've heard also a description of  
5 additional analysis populations. This table attempts  
6 to summarize and present to you two additional  
7 analysis populations as described by the applicant.  
8 These analyses adjust both for randomization errors  
9 and add back that two protocol deviation patients but  
10 not the six early dropouts. So if we look at the  
11 total of 83 and remove the 6 early dropouts,  
12 discontinuations, you have 77 patients. And you get  
13 slightly different results for the analysis based --  
14 because of this additional issue of the problem with  
15 randomization.

16 As you've heard, 5 pairs or 10 patients, 12  
17 percent, had their randomization assignments  
18 switched, and apparently these 5 pairs symmetrically  
19 received the opposite treatment, and I think we don't  
20 completely understand just the mechanics of how this  
21 occurred, at least we haven't. The sponsor  
22 identified an evaluable group which they've called  
23 randomized intended. I think they merely mean  
24 analyzed on the basis of treatment that should have  
25 been given, and that gives you these two differences

1 in results for the two-year recurrence rate. And  
2 then the other group is randomized as treated, that  
3 is the analysis population is for patients evaluated  
4 based on the actual treatment as given in which the  
5 10 patients were switched.

6 So putting these three populations  
7 together, this table indicates the applicant's  
8 estimated two-year recurrence rate for the two  
9 populations that adjust for the randomization errors  
10 and exclude the six early dropouts and the protocol  
11 population, the lower row here, shows the initial  
12 group I had defined for you minus eight patients.  
13 And again, the log-rank procedure was used for the  
14 analysis test, statistical test.

15 Now, this table adds one additional  
16 category. This is an attempt to get to an ITT  
17 population. All patients as they were assigned, that  
18 is as the treatment was intended to be, and this also  
19 uses a form of worst case censoring. The mechanism  
20 of censoring was not prespecified in the protocol.  
21 So with a total of 83 patients, and again using the  
22 applicant's calculation for their estimated 2 year  
23 recurrence rate, one has a 2 year recurrence rate of  
24 38 percent on the Synergo arm versus 51 percent for  
25 the control arm of mitomycin C alone, and application

1 of the log-rank test does not yield a value less than  
2 .05.

3 In this worst case censoring analysis, as  
4 the sponsor did mention, the six dropouts were added  
5 back, and the assumption was made that one control  
6 patient was recurrence free and the five Synergo  
7 patients who had dropped out had disease recurrence,  
8 and that's why this is termed a worst case censoring  
9 application.

10 Now, this represents from the January 2008  
11 applicant information. The analysis of a per-  
12 protocol group. This is not intention to treat.  
13 This is the 75 patient subgroups of the entire study.  
14 The lower curve is mitomycin C alone, the blue line  
15 here, and the number of patients on each arm are  
16 indicated here, and this is to capture the two-year  
17 experience.

18 This curve represents the data analysis if  
19 one used the intention to treat population combined  
20 with that description I gave you of worst case  
21 censoring. So all patients are added back and  
22 censoring convention applied as described here. The  
23 blue line represents again the mitomycin C control  
24 group. There's an early decline in the Synergo group  
25 due to the dropouts, the curves cross and then they

1 are roughly parallel but the mitomycin C line is  
2 lower.

3 Now, turning to secondary endpoints, the  
4 applicant has provided information. No stage  
5 progression was identified. No CIS subsequently  
6 developed, and a late follow-up of those patients  
7 available occurred, and no overall survival  
8 difference was identified between the two groups of  
9 patients. The applicant has reported to us a total  
10 of 14 deaths. No bladder cancer deaths. Five are  
11 described as due to tumors in other organs. We do  
12 not have pathology reports to review for that. There  
13 were three patients in this January summary to the  
14 FDA described to have metastases on the mitomycin C  
15 arm based on scans, but we have no pathology to  
16 examine here.

17 Now, this is the overall survival analysis  
18 submitted earlier. This is not the intention to  
19 treat population, but this is the -- and this is a  
20 shorter time interval than what you had seen earlier  
21 from the applicant. And one would not ordinarily  
22 expect to see a difference in overall survival due to  
23 the nature of this condition.

24 Turning to adverse events, there were more  
25 reported in total on the hyperthermia arm. We noted

1 that pain occurred in 40 percent of patients versus 0  
2 percent on the control arm. There was more dysuria  
3 reported, and this bladder wall thermal reaction,  
4 which we'll talk a little bit more in a moment.

5 Pain on the Synergo therapy was reported in  
6 17 of the 42 patients who received the hyperthermia.  
7 It was described as pain or intolerance to treatment  
8 in 10 of 42, and bladder spasms attributed to the  
9 catheter placement in another 7. And of this total  
10 17 patients, 7 were rated as mild, 7 moderate and 3  
11 were rated as severe. One Synergo patient, as well,  
12 developed a false urinary passage attributed to the  
13 larger size of the catheter that's necessary for this  
14 therapy.

15 The posterior wall thermal tissue reaction  
16 as mentioned related to the radio frequency antenna  
17 is visible by cystoscopy. The investigators  
18 described the condition as consisting of tissue  
19 necrosis in 23 of the 42 patients, 55 percent, and it  
20 was graded as mild, moderate or severe, in these  
21 numbers. A posterior wall ulcer was reported in 1 of  
22 the 42. Now, as described, none of these were  
23 determined to be transmural and no interventions were  
24 described as necessary to manage the condition.

25 Looking at another of the applicant's

1 studies for just a moment, the EPP, European  
2 Prophylactic Patient study, is an ongoing study of  
3 the commercial use in practice of this device. It is  
4 a single-arm study, and we are informed of 186  
5 patients enrolled, 122 evaluated for efficacy. On  
6 this single-arm study, the estimated two-year  
7 recurrence rate is identified by the applicant as 32  
8 percent. While we are reluctant to do cross-study  
9 comparisons, we did observe that in Study 102.1, the  
10 BCG alone control arm has a two-year recurrence rate  
11 of 32 percent.

12           Now, this is a study, commercial use in  
13 practice. With regard to the studies discussed here,  
14 Study 102.1 addresses a different question than 101.1  
15 but the safety does appear to be similar.  
16 Approximately one-third of the patients are accrued  
17 at this time, and we do await the full study report.  
18 Study 101.1 is a prospective study of 83 patients, 42  
19 of whom received the Synergo hyperthermia treatment  
20 for the proposed approval indication. The general  
21 treatment plan, the endpoints to your follow-up are  
22 reasonable choices for this disease. Adverse events  
23 are greater on the Synergo arm including all AEs,  
24 pain, dysuria, and this posterior bladder wall tissue  
25 reaction.

1           Now, in examining the study, the protocol,  
2 the information submitted and the amendments between  
3 2001 and now with additional data, we do want to  
4 indicate what we term clinical limitations regarding  
5 the sponsor's pivotal study. At the time of  
6 randomization, on the randomization form, the  
7 patient's baseline characteristics were disclosed.  
8 The information was available. A misrandomization  
9 occurred. At one site, the BIMO auditor examined  
10 eight patients and four were misrandomized, and  
11 you've heard that a total of 5 pairs of patients have  
12 been identified in which the treatment assignments  
13 were switched, 10 of 83 or 12 percent of the  
14 population.

15           The original study protocol lacked a  
16 prespecified definition for the primary endpoint, or  
17 for an interim efficacy analysis performed after 64  
18 patients had been enrolled, or the decision to stop  
19 at 83. Also a plan for analysis of the missing data  
20 and censoring conventions were lacking.

21           With respect to data collection in the  
22 study, case report forms to document study events  
23 were not introduced until three to four years after  
24 the study started. Two-thirds of the patients had  
25 been enrolled by that point. And while the source

1 records as audited by the BIMO examiner appeared to  
2 be transcribed accurately, there's no assurance that  
3 additional vital study information was not lost by  
4 the failure to document prospectively for all  
5 patients.

6 Pathology information regarding baseline  
7 and follow-up biopsy results was not prospectively  
8 and consistently captured in a structured format to  
9 assure accuracy and completeness of the information.  
10 For example, in the pathology reports available,  
11 baseline description regarding the presence of or  
12 involvement of the muscular layer of the bladder was  
13 inconsistently reported.

14 Also with respect to pathology, there was  
15 no central pathology review performed. And in the  
16 applicant's submission of January 2008, amendment 16,  
17 it states specifically no pathology report is  
18 available for 4 of 29 of the total reported  
19 recurrences.

20 Now, blinding has been commented on  
21 earlier, and patients and investigators knew their  
22 treatment plan. We agree, blinding of the treatment  
23 of itself was not practical. We do point out that no  
24 independent assessment or adjudication of the  
25 endpoint of tumor recurrence occurred in this study.

1           The dwell time in the bladder for the  
2 intravesical mitomycin C was not consistently  
3 recorded. Part of this issue had to do with the late  
4 availability of the case report forms. But treatment  
5 differences may have resulted from that variable, and  
6 we're not fully able to assess that.

7           For the Synergo patients, mitomycin C was  
8 removed by the catheter, whereas for the control  
9 patients, they were advised to void per schedule.

10           We also noted that the follow-up cystoscopy  
11 exams were not performed consistently at study sites.  
12 They may have been performed at local sites for  
13 patient convenience. We're not entirely clear of  
14 what information was captured or forwarded or how  
15 those study sites may have been instructed to perform  
16 the examinations.

17           Concomitant medications were not  
18 prospectively recorded and, therefore, concurrent use  
19 of other medications which may inform further  
20 regarding efficacy or safety is uncertain.

21           The potential for bias was not minimized in  
22 the applicant's pivotal study. A substantial  
23 asymmetry in early dropouts on the Synergo arm, seven  
24 versus 1, suggests additional adverse features of  
25 this therapy may be present. The strength of the

1 efficacy conclusions vary with the choice of the  
2 analysis population. Study 102.1 shares some of the  
3 similar concerns as we have identified in 101.1.

4           Now, of interest, the Study 101.1 was  
5 reported in the Journal of Clinical Oncology in 2003  
6 by site investigators, and they described it as a  
7 multicentric study comparing intravesical  
8 chemotherapy alone with local hyperthermia for  
9 prophylaxis of recurrence at the three sites  
10 indicated and consisting of the 83 patients  
11 randomized. And the authors' conclusions in this  
12 published analysis was that, "These results were  
13 preliminary and need to be confirmed by larger,  
14 prospective, multicenter studies. The combined  
15 treatment was more cumbersome and requires a larger  
16 catheter. However, the reduction in the proportion  
17 of recurrences at 24 months in favor of thermo-  
18 chemotherapy encourages further clinical  
19 investigations."

20           Thank you.

21           DR. LI: Good morning, ladies and  
22 gentlemen, distinguished Panel members and guests.  
23 My name is Xuefeng Li, a Statistician in the Center  
24 for Devices and Radiological Health. I'm going to  
25 give you a brief overview of this PMA from the

1 statistical perspective.

2           This is the outline for my presentation.  
3 First, I will give you a brief introduction of the  
4 pivotal study design and the effectiveness results.  
5 Then I'll talk about the limitations of the pivotal  
6 study. After that, I will talk about other relevant  
7 data and studies for the Synergo device and their  
8 limitations. Finally, I will conclude with a  
9 summary.

10           The pivotal design Study 101.1 is a  
11 randomized, multicenter, unblinded study. Since  
12 previous speakers have talked about this design in  
13 great details, so I will not repeat them here. Note  
14 that the sponsor originally planned to enroll 150  
15 patients. An intimate look was taken when 64  
16 patients were enrolled and the study was stopped  
17 early when 83 patients have had had 2-year data.

18           This table gives a patient accountability  
19 and also lists the analysis populations. Since you  
20 have also heard it before, I will not repeat it here.

21           So this table summarizes effectiveness  
22 results for the primary effective endpoints based on  
23 different patient populations. In the all study  
24 patient population, missing patients were treated as  
25 failure in the Synergo group but success in the

1 control group. Except for this worst-case scenario,  
2 the Synergo group appears to have statistical lower  
3 two-year probability of recurrence than the control  
4 group. Note that those P-values are not adjusted for  
5 the interim look or the early stopping.

6 This table summarizes the results for all  
7 secondary endpoints.

8 Now, I would like to discuss some  
9 limitations associated with the pivotal trial.  
10 Actually, Dr. Kane has mentioned most of them as  
11 clinical limitations. Since they also have  
12 statistical impact, I will repeat them here briefly.  
13 First, this trial is an unblinded study. An  
14 independent, blinded review committee to adjudicate  
15 tumor related endpoints was not implemented in this  
16 study. Second, the sponsor listed six endpoints but  
17 did not clearly specify which ones are primarily and  
18 which ones are secondary in the protocol. Third, the  
19 primary hypothesis was not clearly stated, and the  
20 statement about the analysis method of the primary  
21 endpoint was not included in the protocol.

22 However, also primary endpoint and the  
23 corresponding hypothesis could be inferred from the  
24 weight of statement about the sample size, first, how  
25 missing data would be handled was not pre-specified

1 in the protocol. Finally, the interim analysis was  
2 severely lacking in detail. The adverse -- function,  
3 the -- of earlier stop for success and the sample  
4 size recalculation method were not stated in the  
5 protocol.

6           When considering these limitations  
7 together, it is very difficult to evaluate the  
8 potential bias introduced by them.

9           As you have seen, the sponsor also provided  
10 the long-term data and other relevant data and  
11 analysis for Synergo treatment.

12           So next I will briefly go through these  
13 analyses and discuss some limitations associated with  
14 them.

15           Here are the results of the long-term data  
16 of Study 101.1. Regarding the recurrence rate, the  
17 Synergo group has lower rates than control at both 5  
18 years and 10 years. In terms of the overall  
19 survival, there are nine deaths in the control group  
20 and five deaths in the Synergo group until the year  
21 2006. None of these deaths were due to the device.

22           It should be noted that the analyses of the  
23 long-term data were not pre-specified. The sponsor  
24 used the as-treated population and the missing data  
25 was not.

1           Here the descriptive data from Study 102.1,  
2 the sponsor expected to enroll 150 patients per arm.  
3 At the interim look, there are 51 Synergo patients  
4 and 53 control patients from Study 102.1, which are  
5 about 1/3 of the total patients expected to be  
6 enrolled. The 2-year recurrence rate for the Synergo  
7 group is 16.9 percent and for control group is 31.7  
8 percent. No progression was observed at the interim  
9 look.

10           There are some limitations with the data  
11 from the Study 102.1. First, the data is from an  
12 unplanned interim look and the study is still  
13 ongoing. Second, the Study 102.1 and the Study 101.1  
14 have different treatment regimen as previously  
15 discussed by Dr. Kane. Third, it is difficult to  
16 fully verify the patient comparability.

17           The sponsor also provided a combined  
18 analysis that pooled the data from Studies 101.1 and  
19 102.1 together. The 2-year recurrence rate of the  
20 combined Synergo group is 17.1 percent based on the  
21 per protocol group.

22           Again, this combined analysis was not a  
23 planned analysis, and the patient comparability  
24 cannot be fully verified.

25           The sponsor also conducted the literature

1 review of control treatments in Studies 101.1 and  
2 102.1. Based on meta-analysis, the sponsor reported  
3 a 2-year recurrence rate of 41.5 percent for the MMC  
4 treatment and 35.5 percent for the BCG treatment.

5           Again, this literature review is an  
6 unplanned analysis, and it is very difficult to  
7 verify the poolability across studies. In addition,  
8 the literature review is -- subject to publication  
9 bias. Here the publication bias means that positive  
10 results have a better chance of being published.

11           To summarize, the pivotal study appears to  
12 show that the Synergo treatment had a smaller two-  
13 year recurrence rate than that of the control.  
14 However, there are some limitations associated with  
15 this study as I have mentioned before. When these  
16 limitations are considered collectively, it  
17 potentially impairs the ability to interpret the  
18 statistical results and increase the uncertainty of  
19 the state of the conclusions. Other analyses  
20 provided by the sponsor appear to show supportive  
21 evidence of the safety and effectiveness of the  
22 Synergo device, but they all have their own  
23 limitations as mentioned before.

24           This ends my presentation. Thank you very  
25 much. Now, I turn the podium to Dr. Wei.

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 DR. WEI: Good morning, distinguished Panel  
2 members and members of audience. I'm Dr. Shaokui  
3 Wei. I'm an Epidemiologist in Division of Post-  
4 Market Surveillance, Office of Surveillance and  
5 Biometrics.

6 Today I will talk about post-approval study  
7 as have been proposed for the Synergo hyperthermia  
8 system submitted by Medical Enterprise. This  
9 presentation is based on the original protocol  
10 submitted to the FDA on January 4, 2008, that differ  
11 slightly from what the sponsor presented earlier.

12 Before we talk about the post-approval  
13 study, we need to clarify that the discussion of a  
14 post-approval study prior to a formal recommendation  
15 on the approvability of the PMA should not be  
16 interpreted to mean FDA is suggesting the Panel find  
17 the device approval. The plan to conduct the post-  
18 approval study does not decrease the threshold of  
19 evidence required to find the device approval. The  
20 pre-market data submitted to the Agency and discussed  
21 today must stand on its own in demonstrating a  
22 reasonable assurance of the safety and the  
23 effectiveness in order for the device to be found  
24 approvable.

25 This is the outline of my presentation. I

1 will start by describing the general principles and  
2 objectives of the post-approval studies. Here are  
3 the two general principles on the originals for the  
4 post-approval studies. The first one is to evaluate  
5 device performance and the potential device-related  
6 problems in the broader population over an extended  
7 period of time after pre-market establishment of  
8 reasonable device safety and effectiveness. Second,  
9 post-approval studies should not be used to evaluate  
10 unresolved issues from the pre-market phase that are  
11 important to the initial establishment of the  
12 reasonable assurance of the device safety and  
13 effectiveness.

14           These are the general objectives of the  
15 post-approval studies. The post-approval studies are  
16 conducted to gather post-market information on the  
17 longer-term performance. As we all know, pre-market  
18 clinical data are collected from patients that are  
19 highly selected and treated by the community-based  
20 physician. In contrast, the device be permitted to  
21 be on the market, patient that receive the device are  
22 less respected and the physicians who treated those  
23 patients are not limited to a physician that are  
24 participating in the pre-market study. Therefore, a  
25 post-approval study are also conducted to obtain the

1 data on (1) community performance; (2) effectiveness  
2 of training programs; (3) performance of the device  
3 in subgroups of the population. Additionally, some  
4 of the rare adverse events that were not observed in  
5 the pre-market may be present in the post-market  
6 phase as observation period attained and the patient  
7 population -- Therefore, post-approval study can be  
8 designed to monitor such adverse events.

9 Another reason for conducting post-approval  
10 studies to address its use and concern that the Panel  
11 members may -- based on their observation and  
12 experience.

13 Now, I will present the overview of the  
14 sponsor's proposal followed by assessment of the  
15 proposed study.

16 The Panel's proposal is to conduct a  
17 prospective, single-arm study to evaluate the safety  
18 of the device compared with the pre-market data  
19 101.1. The main study endpoint will be the frequency  
20 of the anticipated adverse events, and the secondary  
21 endpoint will be any unanticipated adverse events.  
22 And, they plan to conduct null and alternative  
23 hypothesis test for each adverse event. They propose  
24 to include 211 patients in the intermediate and high-  
25 risk STCCB, including transurethral resection of

1 their tumor under the 12-month follow-up period with  
2 assessment every 3 months. They claim that the 211  
3 subjects will allow for 80 percent of the power.

4 The next two slides I will present more  
5 detail about the primary endpoint and the study  
6 hypothesis.

7 As I mentioned earlier, the primary  
8 endpoint is frequency of anticipated adverse event.  
9 The adverse events are listed here, and they are  
10 posterior wall thermal reaction, urethral stenosis,  
11 hematuria, false route, hypotonic bladder, reduced  
12 bladder capacity, urinary tract infection and  
13 necrosis in other areas of the bladder.

14 And for each of adverse event risk, a null  
15 or alternative hypothesis will be conducted. The  
16 null and alternative hypotheses are formulated as  
17 follows: the null hypothesis is  $P$  greater than or  
18 equal to the  $\pi$  plus delta; alternative hypothesis  $P$   
19 less than  $\pi$  plus delta. Where the  $\pi$  represent the  
20 current rate of adverse event as established in the  
21 pivotal Study 101.1, and the delta is the largest  
22 acceptable difference between the study rate and  
23 the -- rate. The sponsor proposes the delta to be 10  
24 percent for proportions over 50 percent and 5 percent  
25 when adverse event occurrence rate is small, for

1 example, a rare event. For each test -- will be  
2 claimed if the upper limit of 90 percent -- if less  
3 than the pi plus the delta.

4 Now, I would like to move onto the  
5 assessment of the sponsor's proposal. The first  
6 concern we have related to the basic study design,  
7 comparison group, the sponsor proposes to use a pre-  
8 market study group and a target control group. FDA  
9 believes this proposal is limited because, as I  
10 mentioned earlier, pre-market study of performance is  
11 highly selected patient population and are conducted  
12 in highly specialized center, and not represent or be  
13 comparable with post-market experience. And all the  
14 pre-market data came from outside U.S. We believe  
15 the most appropriate post-market comparator may be  
16 the existing of standard of care.

17 Next, I will present more pre-market data  
18 to make one more argument.

19 This table shows adverse events the sponsor  
20 plans to evaluate in the post-market study. Please  
21 note the pre-market safety profile of the device  
22 different by the study. Therefore, it is difficult  
23 to determine which pre-market study is the best  
24 comparator. Additionally, there are some frequent  
25 adverse events in the pre-market study that the

1 sponsor did not include in the original PSA proposal,  
2 such as pain which was observed in 41 percent in the  
3 Synergo patients in Study 101. There was no pain in  
4 the control Study 102. Fifty-three percent of the  
5 Synergo patient suffer with pain versus two percent  
6 in the control.

7           Starting -- with the study design, the  
8 study endpoint, the sponsor only considered the  
9 evaluation of the safety. There is no plan to assess  
10 post-market effectiveness. The long-term  
11 effectiveness in large and diverse patient population  
12 is still not clear.

13           Additionally, please note that although  
14 both 101 and 102 studies show Synergo has a reduced  
15 two-year recurrence rate, the deficiencies in these  
16 studies, addressed by Dr. Kane and Dr. Li in the  
17 previous FDA presentation, make the recurrence data  
18 difficult to interpret. There is also the question  
19 about the long-term survival. The pivotal Study 101  
20 show that there was no significant difference in the  
21 overall survival rate between the two treatment  
22 groups. There was a total of 14 deaths, 9 in the  
23 control group, 5 in the Synergo, reported since the  
24 start of the clinical study in 1999 until later in  
25 2006. None of the deaths were known to be due to

1 STCCB.

2           The third limitation of study design is the  
3 length of follow-up. Current proposal only considers  
4 one year of follow-up. Since STCCB has a long  
5 waiting period, slow progression, most cancer  
6 recurrence may not occur in the one-year follow-up  
7 and some adverse event, such as urethral stricture  
8 develop slowly and appear after one year of -- So  
9 one year of follow-up may not be sufficient to assess  
10 long-term performance of this device.

11           Our first concern relates to the  
12 statistical analysis plan. They plan to conduct a  
13 hypothesis test for each anticipated adverse events.  
14 There's a 10 percent delta for the adverse event with  
15 50 percent frequency and a delta of 5 percent for the  
16 rare adverse event. The definitions of a common  
17 versus a rare event are not clear. The rationale for  
18 the selection of dealt is not presented.

19           Specifically, here are chosen delta for  
20 each adverse event. Those at the present are less  
21 than 10 percent of frequency in the pre-market Study  
22 101. If you compare with the post-market rate there  
23 was a 5 percent delta. The posterior wall reaction  
24 if you compare is a 10 percent delta.

25           Now, I will present the issues that we

1 would like the Panel members to consider this  
2 afternoon when addressing the Panel questions.

3           First, we would like to hear your comments  
4 on the comparison group. As mentioned earlier, the  
5 sponsor's proposal to use a pre-market study group as  
6 a historic control group. Please discuss what is the  
7 most appropriate comparison group for the post-  
8 approval study? We believe the pre-market group is  
9 not appropriate as a comparison group for the post-  
10 approval study and the -- of the more appropriate  
11 comparison group may be the standard of care.

12           Second, we would like to hear your comment  
13 on the study endpoint. The proposed study only  
14 includes the evaluation of the device safety but it  
15 does not include a plan to assess the post-market  
16 effectiveness. Should the effectiveness be studied  
17 post-market? If so, what endpoint should be studied?  
18 Should it be the recurrence rate, patient survival or  
19 both?

20           And thirdly, the current proposal consider  
21 only one-year follow-up. We would like you to  
22 discuss whether one year of follow is sufficient to  
23 assess the long-term performance of this device.  
24 Please consider that STCCB has a long waiting period  
25 and a slow progression, and the some of the adverse

1 events develop slowly.

2           Finally, please discuss and make  
3 recommendations on the definition for rare versus  
4 common adverse events and the rationale for the  
5 chosen deltas.

6           This concludes my presentation. The floor  
7 is now open to the questions for the FDA. Thank you.

8           DR. TALAMINI: I'd like to thank the FDA  
9 speakers for their presentations and open up the  
10 Panel to questions specifically for the FDA. We have  
11 15 minutes scheduled prior to lunch for that, and  
12 we're just about on time. So I'd ask Panel members,  
13 if you have questions specifically to the FDA  
14 regarding their presentations. Dr. Connor, perhaps I  
15 would invite you specifically since there's so much  
16 statistical --

17           DR. CONNOR: (Off mic.)

18           DR. TALAMINI: In the meantime,  
19 Dr. Donatucci?

20           DR. DONATUCCI: Yeah. While you're  
21 reviewing your notes, I just have a question for  
22 Dr. Kane. You presented a slide I believe and I've  
23 been looking for it and haven't found it, about  
24 progression, differences between BCG and mitomycin C.  
25 If I recall, there was an eight percent and four

1 percent --

2 DR. KANE: That was Dr. Herrera.

3 DR. DONATUCCI: Who was that? I'm sorry.

4 DR. KANE: Dr. Herrera.

5 DR. DONATUCCI: Dr. Herrera, okay. Maybe  
6 that's why I can't find it. Thank you.

7 DR. CONNOR: So I guess my first question,  
8 maybe this is the question of the day, so it's good  
9 to get FDA's opinion. I'm not used to having to have  
10 the microphone this close. One of the questions  
11 that's continually put out are are there biases in  
12 this study, and we recognize, you know, that the  
13 study design and implementation was not ideal. You  
14 know, we don't know the consequences, the biases of  
15 this study. Do those in FDA who have thought long  
16 and hard about the consequences of these biases have  
17 an opinion, and let me know if this is a question I  
18 can ask, if those biases are enough to overcome what  
19 looks like to be an effective treatment? I mean, is  
20 the concern that biases are so huge that it doesn't  
21 look like this is as effective as at least in the  
22 efficacy data looks like the treatment is?

23 MS. BROGDON: I'd like to respond to that.

24 DR. CONNOR: Okay.

25 MS. BROGDON: This is one of the issues

1 that we really wanted the Panel's input on. I think  
2 it wouldn't be appropriate for us to give an opinion.

3 DR. CONNOR: Okay.

4 DR. TALAMINI: Perhaps I could ask a very  
5 practical question, and I couldn't quite tell, and  
6 this is to FDA specifically. It sounds as if the  
7 sponsors have expressed a change in their plans for a  
8 post-market study. Was the FDA aware of that? And,  
9 has the FDA taken that into consideration in this  
10 current presentation?

11 DR. BAXLEY: The post-approval study that  
12 we presented was sent to us as part of our PMA review  
13 in 2007, and that was the only one submitted  
14 officially to the PMA. I think recently they've had  
15 some change in what they would like to propose, and  
16 that was not formally reviewed by us. That came in  
17 the Panel package that was sent to you all. So we've  
18 looked at it but we haven't formally reviewed it.

19 DR. TALAMINI: Thank you.

20 DR. BAXLEY: We weren't able to present  
21 anything else.

22 DR. TALAMINI: Other questions for the FDA?  
23 Dr. Connor. Okay. Let's go around this way.  
24 Dr. Redman.

25 DR. REDMAN: I have several but I want a

1 clarification that we can ask them questions later on  
2 after --

3 DR. TALAMINI: We indeed can ask questions  
4 later on.

5 DR. REDMAN: One questions involves slide  
6 60, including all the patients with a signed and  
7 tended to treat. Has anybody done the numbers with  
8 just dropping out the three patients that never  
9 received any treatment which probably reflects the  
10 consent process and not the treatment effect in the  
11 worst-case scenario, and if those three patients are  
12 dropped, is the worst-case scenario those that  
13 included patients who received any treatment though  
14 may have stopped treatment prematurely, does that  
15 treatment effect still hold? Should I go through  
16 them or --

17 DR. TALAMINI: Why don't you go through  
18 them and then we'll allow the FDA to figure which  
19 ones they can answer quickly and which might require  
20 more analysis for later today.

21 DR. REDMAN: Slide 70 looked at possible  
22 data being supplied to the randomization. I'm  
23 assuming this was not a stratified randomization.  
24 Pre-study stratified by grade or Ta, T1, in that  
25 regard. Okay. However, even if that didn't occur,

1 my memory is that there was no major difference in  
2 the major prognostic factors. It still held up. So  
3 even though it wasn't stratified for that, it's  
4 unfortunate, but there was still no difference in  
5 that?

6 MS. BROGDON: Excuse me. I think it would  
7 be good if we could get the sponsor's responses on  
8 the record rather than having nodding of heads or  
9 shaking of heads.

10 DR. TALAMINI: So would you prefer that the  
11 sponsors -- that we broaden this set of questioning  
12 now to beyond just the FDA or have that occur later  
13 on?

14 MS. BROGDON: I think it's up to you.

15 DR. TALAMINI: Perhaps we could have an  
16 acknowledgment from the sponsor in response to those  
17 questions.

18 DR. O'DONNELL: Can you repeat the  
19 question?

20 DR. REDMAN: No pre-study stratification.

21 DR. O'DONNELL: Correct. No pre-study  
22 stratification.

23 DR. REDMAN: But your data showed that the  
24 major prognostic factors, there were no differences  
25 between the groups, if my memory's right.

1 DR. O'DONNELL: Right. They were evenly  
2 balanced, and there was no difference in the  
3 prognosis facts between the two groups.

4 DR. REDMAN: This is nexus to the FDA. On  
5 slide 72, just clarify for me, you've seen the CRFs  
6 and you're comfortable with the CRFs but there was a  
7 comment made that even though there was retrospective  
8 data collected on those CRFs, that there was a fear  
9 that some additional vital study information could  
10 not be lost in that process. What was the potential  
11 vital study information that you were concerned  
12 about?

13 DR. KANE: I think if we knew it, we would  
14 not be concerned about it. Our observation was  
15 simply that we were unsure what may not have been  
16 captured, written into a hospital record, that may  
17 have occurred and I guess corollary to that is the  
18 fact that some of the follow-up cystoscopies were not  
19 performed at the study sites, and we're looking at a  
20 great deal of reconstructed information.

21 DR. REDMAN: I accept that, but the second  
22 point you made about the cystoscopies not being done  
23 at the study site, were those cystoscopies during the  
24 first two years of follow-up that were protocol  
25 specific that were done --

1 DR. KANE: Yes.

2 DR. REDMAN: -- outside --

3 DR. KANE: Yes.

4 DR. REDMAN: -- and do you know the number  
5 of those for each arm that were performed not at the  
6 study site?

7 DR. KANE: No, I don't. No, I don't. Just  
8 the comment that that --

9 DR. REDMAN: Just the comment that it  
10 occurred.

11 DR. KANE: That is occurred.

12 DR. REDMAN: Okay. And the original biopsy  
13 material, and I guess this would go to the sponsor,  
14 too, my concern is that even though some of the ones  
15 that did not have muscle on the biopsy material, were  
16 there also blind biopsies done to confirm that there  
17 was no CIS and was that equal between groups and  
18 again that I don't expect you to be able to answer at  
19 this point in time, but it would be interesting to  
20 know.

21 DR. TALAMINI: I think the sponsor has a  
22 response at least to one of those specifics regarding  
23 the studies. So we'll let them respond.

24 MS. STEIN: I just wanted to comment the  
25 issue that Dr. Kane mentioned regarding the

1 cystoscopies were not performed at the site. I'm not  
2 sure exactly where he got that information from  
3 because all cystoscopies were performed only at the  
4 sites by the PIs, by the principal investigators and  
5 recorded both in the hospital source documentation  
6 and on the CRFs. Okay. I just wanted to clarify  
7 that.

8 DR. TALAMINI: Thank you. Dr. Redman,  
9 additional? Dr. Dahm?

10 DR. DAHM: One of the gentlemen that gave  
11 the statistics presentation, I would ask to comment  
12 on what I found as a footnote to slide 86 where it  
13 talks about not adjusted for interim look or early  
14 stopping. I guess you're suggesting that the fact  
15 that the study was stopped early for benefit may be a  
16 potential source of bias. I think that relates to  
17 recent studies in the methodology literature, that  
18 that's just something that we should look for, too.  
19 Is this something you can comment on?

20 DR. LI: We know that if we do an interim  
21 analysis, that -- rate would be inflated. So it  
22 would not be under 5 percent, might be like 10  
23 percent, something like that. So we have to pre-  
24 specify -- function to control the overall -- rate  
25 once we have interim analysis that we could study for

1 early claim of success. So the sponsor did pre-  
2 specify a method for doing that. So we're not sure  
3 whether -- rate is controlled or not. Actually, I  
4 have a backup slide.

5 DR. DAHM: Do I understand correctly that  
6 this would exaggerate the observed treatment effect,  
7 that that would be the direction of the bias if you  
8 stopped the trial early for benefit which is common  
9 by the way, which is very common in these kind of  
10 studies? But is that the concern?

11 DR. LI: Yes. So here these slides lists  
12 two, why we used the method for control of -- rate  
13 when the trial has an interim look. The first one is  
14 O'Brian-Fleming -- rate. If we take a look at we  
15 have 64 patients out of the 150 patients are randomly  
16 expected, so the P-value should be under .0013 in  
17 order to claim study success at this interim look.  
18 The second method is called Pocock boundary. So, for  
19 example, if we take an interim look, and we have 64  
20 patients. So the P-value should be under .0277 in  
21 order to claim early success.

22 So the sponsor P-value, as-treated  
23 populations, is well below this O'Brian-Fleming  
24 boundary. So it's still statistically significant  
25 but for as intended population, the P-values are

1 above this O'Brian-Fleming boundary. So when --

2 DR. DAHM: Can you repeat your last  
3 sentence again?

4 DR. LI: If we use the first method, the  
5 sponsor's analysis based on as intended populations  
6 will not be statistically significant. That means  
7 the sponsor P-value, if we check the slides, my  
8 slides 85 or 86, actually 86, we see that -- excuse  
9 me. Look at the P-values. These two P-values are  
10 above .0013, the threshold from the O'Brian-Fleming  
11 boundary. So the first two analysis populations  
12 would give statistical insignificant results. For  
13 the last two, as-treated populations, we still have  
14 statistical significant.

15 DR. CONNOR: Can I speak to something  
16 Dr. Dahm asked?

17 DR. TALAMINI: Yeah, from the Panel's point  
18 of view, we're still supposed to be asking questions.

19 DR. CONNOR: I think it's a brief  
20 clarifying thing maybe. Your question was whether  
21 the interim analysis would bias the treatment effect  
22 estimate. So the answer is the treatment effect we  
23 observe is not biased by looking early. The  
24 probability of them making a type one error, the  
25 probability of them, you know, saying the effect,

1 that is inflated, but the size of the effect we  
2 estimate is not inflated by this. And that's about  
3 all we can talk about.

4 DR. TALAMINI: So I want to make sure that  
5 we keep to asking the FDA questions regarding their  
6 presentations as much as we can while still being  
7 fair to the sponsors to allow them to respond. So  
8 other -- does that get to your question?

9 DR. DAHM: Yes, it does. Thank you.

10 DR. DONATUCCI: Can I? I'm confused now.  
11 I wasn't confused going in. Now, I'm really  
12 confused. You know, I haven't had statistics since  
13 college but my understanding of a P-value was that  
14 either you prove the known hypotheses or you don't.  
15 So it's either by chance or not by -- it's just a  
16 percentage of probability that this is a chance  
17 observance rather than a treatment effect. I've  
18 never seen before anybody say specifically -- this  
19 boundary thing is completely new to me. So you're  
20 saying that, and I used to dismiss anything less than  
21 .01 as meaningless, the difference of being  
22 meaningless, but you're now telling me that you need  
23 to see, in order for that treatment effect to be real  
24 and not just by chance, in an interim analysis you  
25 need to see like Pocock's hypothesis, three

1 statistically -- three significant figures? You  
2 know, so less than .01 is meaningless to me now.  
3 It -- I'm lost. Help me here.

4 DR. LI: Yes. My idea is that you  
5 frequently check the data. By chance you may observe  
6 success earlier but this will inflate the -- rate.  
7 That means actually there's no treatment effect but  
8 you claim there is a treatment effect if you still  
9 use .05 to control the -- rate. So we have to --

10 DR. DONATUCCI: Lower the P-value.

11 DR. LI: -- lower the P-value, lower the  
12 bar.

13 DR. DONATUCCI: That's pretty low. That's  
14 the first time I heard of that. Thank you.

15 DR. TALAMINI: So, Dr. Connor, perhaps  
16 you'd be willing to -- as I understand it, as a  
17 layman, this is the statistical bar that you need to  
18 say that you can stop the study early or evaluate  
19 midway through. Is that correct?

20 DR. CONNOR: Right. You can imagine a  
21 situation where every time a patient comes in, we  
22 look at the data again, and so we would be looking at  
23 the data 150 times in a 150 patient trial. If ever  
24 that P-value got below .05, the sponsor goes, we win.  
25 Obviously the type one error rate is much higher than

1 .05 in that. So I think that example may make sense.  
2 So these values are a way to quote spend, your 5  
3 percent error. The sponsor says we have, you know,  
4 or the agreement between FDA and the sponsor says we  
5 have 5 percent error rate. We're willing to make an  
6 error rate of, you know, 1 percent early on and then  
7 have a 4 percent error rate later. So it's a way to  
8 spread out that 5 percent chance of saying something  
9 is effective when it actually isn't.

10 DR. TALAMINI: So on this particular  
11 question, does the sponsor have a comment or  
12 response? Yes. Yes, sir.

13 DR. MARCOVICH: Just a couple of minor  
14 comments. First of all, at the time of this, this  
15 was called by the Data Safety Monitoring Board, and  
16 they weren't aware of the randomization or it was  
17 discovered much later. The second thing is that the  
18 overall recurrence rates in the 101 study is  
19 essentially identical to the 102 study. So this does  
20 not appear to be just a fluke. Those curves are, you  
21 know, right on, and the mitomycin recurrence rates  
22 are right on with the expectation of the literature.  
23 So there doesn't appear to be anything unusual in the  
24 results in the 101 study.

25 DR. TALAMINI: Thank you. Are there

1 questions about the FDA presentation from the Panel?

2 Yes, sir.

3 DR. BHUTANI: I have questions for the FDA.  
4 I refer to slide 13 which is on the chronology. In  
5 2001, the PMA was submitted to the FDA. What  
6 feedback was provided to the sponsor regarding the  
7 study design or problems, if any, or about  
8 statistical concerns and then I want to know what the  
9 in-depth review cycles are from 2001 to 2004. I  
10 presume they are with the FDA reviewing the data, and  
11 then lastly in 2005, when FDA conducted bioresearch  
12 monitoring inspections, and found the CRFs,  
13 retrospectively collected CRFs to be acceptable, was  
14 any feedback provided to the sponsor about potential  
15 problems of vital information being missed that could  
16 become an issue if the device is brought for a PMA to  
17 the FDA in the future? Thank you.

18 DR. BAXLEY: Let me start with the in-depth  
19 review cycles. Those were times when the ball was in  
20 FDA's court for our review, and the first review  
21 cycle was with review of the original PMA. The next  
22 two were with review of additional information we  
23 received, that we requested, and they were mainly  
24 with how the study was conducted and what the results  
25 are and additional information on case report forms

1 and pathology reports, detailed information. The  
2 bioresearch monitoring inspection was just to audit  
3 the data and to see if it's factual and it was. But  
4 the concern raised here is a review concern just on  
5 how the study's conducted, not an audit issue.

6 DR. TALAMINI: Yes, sir. Dr. Layton.

7 DR. LAYTON: I have a couple of questions  
8 relative to some of the earlier slides. One of them  
9 was the preclinical review and during the preclinical  
10 review you said that you had no concerns with the  
11 device other than the shelf life. Has that been  
12 taken care of? And the preclinical review that you  
13 did do relative to that, is there any comments  
14 relative to the failure modes or the adverse events?

15 DR. BAXLEY: Our first several review  
16 cycles were also heavily involved with the  
17 preclinical review, and additional information was  
18 provided in those two cycles that adequately  
19 addressed, you know, the concerns we had. So we're  
20 looking at that in depth during the early phase. And  
21 there was hazard analysis information in there that  
22 looked at failure modes and reliability and that's  
23 since been resolved.

24 DR. LAYTON: All right. The other question  
25 that I have was relative to, and I think this is what

1 I heard, during the clinical data sources slide, you  
2 made a comment, I thought, of unplanned studies, and  
3 I didn't know what you meant by that.

4 DR. BAXLEY: I think I used unplanned  
5 interim analysis of studies.

6 DR. LAYTON: Okay. All right. Thank you.

7 DR. BAXLEY: But not unplanned studies.

8 DR. LAYTON: All right. Thank you.

9 DR. TALAMINI: Seeing no other questions, I  
10 think there's at least one or two questions still on  
11 the table. Dr. Redman had the question about the  
12 worst-case scenario re-analysis that I hope the FDA  
13 group can address or can answer for us later on. I  
14 think that's the only one that's still on the table.  
15 Is that correct?

16 Okay. So given that, we will plan to break  
17 for lunch. I remind the Panel members that there  
18 should be no discussion of the PMA during the break  
19 amongst themselves, with the sponsor or with the  
20 public. We will reconvene again in this room at 1:00  
21 p.m. Please take any personal belongings you may  
22 want with you at this time. The room will be secured  
23 by FDA staff during the lunch break. You will not be  
24 allowed back into the room until we reconvene. Thank  
25 you. See you at 1:00.

1                   (Whereupon, at 12:00 p.m., a luncheon  
2 recess was taken.)  
3  
4  
5  
6  
7  
8  
9  
10  
11  
12  
13  
14  
15  
16  
17  
18  
19  
20  
21  
22  
23  
24  
25

A F T E R N O O N   S E S S I O N

(1:05 p.m.)

1  
2  
3 DR. TALAMINI: I'd like to call the Panel  
4 back to order, it being it looks like about 1:05.  
5 And I'd like to take up, if we could with the  
6 questions still on the table just prior to lunch that  
7 Dr. Redman asked the FDA to clarify, and I think we  
8 do have a response. Is that correct?

9 DR. BAXLEY: Yeah. The information that  
10 was asked about in that question, FDA didn't have  
11 presented to us. In the meantime, the firm has  
12 provided that information and a copy has been handed  
13 out to the Panel of a Kaplan-Meier curve. If you  
14 have any questions, I would refer to the company.

15 DR. TALAMINI: Okay. With that being said,  
16 I would ask the company to come forward to address  
17 any detailed issues raised in the morning session  
18 that the company has been asked to address after the  
19 lunch break and perhaps also this new data analysis.  
20 So sponsors?

21 MS. DEUTSCH: Okay. I will just quickly go  
22 through the re-analysis that we provided.

23 DR. TALAMINI: Can everybody hear her?

24 UNIDENTIFIED SPEAKER: No.

25 MS. DEUTSCH: Sorry. What we were

1 requested to do was provide an analysis of the worst-  
2 case scenario excluding three patients from the  
3 analysis and this is what we have in the Kaplan-Meier  
4 curve that you have in front of you. The total  
5 number of patients instead of 83 is 80. Estimated  
6 from the Kaplan-Meier curves, you can see that  
7 there's a 25.6 recurrence rate in the Synergo group  
8 and nearly 60 percent in the mitomycin group, the --  
9 test. We still have a significant at P equals 0.037.

10 DR. TALAMINI: Thank you. Other responses  
11 from the sponsors?

12 DR. GROSSMAN: I think we've addressed most  
13 of the major concerns that we were asked this  
14 morning. A couple of minor points. The number of  
15 adverse events recorded prior to '97 and after '97  
16 were very similar, suggesting that the doubts, the  
17 review of the case before us, retrospective case  
18 report forms was not a problem. The overall  
19 presentation of the data of the 101 and 102 studies  
20 again show consistency in the data, the overall  
21 experience with mitomycin and BCG are very similar.  
22 This is in general a very high-risk group. Half the  
23 patients were high-risk overall. So this puts the  
24 treated population at very high risk for recurrence,  
25 and I think the data stands for themselves. We, of

1 course, would be happy to address any other concerns  
2 that you have.

3 DR. TALAMINI: Thank you very much.

4 So now the Panel enters into a phase of  
5 public discussion, both questions and discussion  
6 which eventually will funnel down to answering the  
7 specific sets of questions that the Panel has before  
8 them for consideration but before getting to those  
9 specific questions, we'd like to entertain any other  
10 more general discussion or questions that would open  
11 up such discussion from the panel. So I'll open it  
12 up to any Panel member. Yes, ma'am.

13 MS. STOKES: I have a question of FDA. The  
14 introduction of the lit review of controls, what  
15 significance does that play? The slide number is 96.

16 DR. TALAMINI: Perhaps you could repeat the  
17 question for us.

18 MS. STOKES: Slide number 96 refers to lit  
19 review of controls, and my question is what is the  
20 significance of introducing lit review, and is it a  
21 necessary support of evidence?

22 DR. LI: So can you repeat the question  
23 again? I'm sorry.

24 MS. STOKES: If you look at slide 96,  
25 you've introduced the lit review of the controls, and

1 my question is, what is the significant of that  
2 information?

3 DR. LI: Well, the sponsor wanted to show  
4 that data they observed in Study 101.1 and 102.1,  
5 only for control arms, they're consistent with those  
6 studies in the literature. That's their purpose, to  
7 show this literature review. So I'm a statistician.  
8 I don't know how to evaluate this particular  
9 information. So I just wanted to point out that  
10 there are some limitations with the literature review  
11 from a statistical perspective.

12 MS. STOKES: Okay. Thank you.

13 DR. LI: You're welcome.

14 DR. TALAMINI: Other questions from the  
15 Panel? Dr. Connor.

16 DR. CONNOR: So this is a question  
17 regarding the choice of control group and precisely  
18 what was done in the control group. So I notice that  
19 the MMC in the first study, 101.1, the MMC was  
20 preheated in the treatment group but not administered  
21 at room temperature in the control group. Presumably  
22 increasing the temperature of the bladder contributes  
23 to the efficacy benefit we see. I wondered if  
24 there's any information whatsoever on what the  
25 treatment effect would be just from preheating the

1 MMC but not using your device. Or why for instance  
2 or why perhaps that in the control group, the drug  
3 wasn't preheated similar to like it was in the  
4 treatment group?

5 DR. TALAMINI: So this is a question for  
6 the sponsor?

7 DR. CONNOR: Yes.

8 DR. O'DONNELL: So the traditional way we  
9 give intravesical chemotherapy is we either make it  
10 up in our clinic or we have our pharmacy group  
11 prepare it and send it down to us. It comes at room  
12 temperature and we give it. So that's basically the  
13 standard of how we do it. You can imagine that the  
14 bladder, as any body cavity, is a huge heat sink. If  
15 you put something in it, it will quickly equilibrate  
16 to the body temperature, whereas it may take a few  
17 minutes to make that equilibration. I'm not aware  
18 that anyone has been able to show, in fact, there  
19 were earlier attempts to preheat liquids and put them  
20 in the bladder and nothing ever came of that from a  
21 clinical standpoint, and it probably is due to the  
22 fact that simply, you know, small volumes in a big  
23 cavity, boom, it's going to equilibrate rather  
24 quickly and not make any difference.

25 DR. CONNOR: Can I follow that along then

1 and ask why then you did preheat MMC for patients in  
2 the treatment group?

3 DR. GROSSMAN: That was done to prevent a  
4 significant altered sensation because the bladder was  
5 being heated by the microwave, and we didn't want to  
6 put in cold mitomycin to increase the rate of bladder  
7 spasms. So we were trying to get the mitomycin up to  
8 temperature to the bladder and then, in fact, during  
9 the treatment, the mitomycin is extracted and cooled  
10 and circulated so it doesn't overheat but it was just  
11 done initially to equalize things and minimize  
12 thermal shock.

13 DR. CONNOR: Okay. Thank you.

14 DR. TALAMINI: So I have two questions  
15 which will sound naïve, but I think they're probably  
16 important for the public record. The first is what  
17 is the evidence that this is not simply giving  
18 systemic doses of mitomycin C and that that is the  
19 effect that we're seeing here, and the microwave and  
20 so forth is simply an epi phenomenon. That's  
21 question number 1.

22 Question number 2, if you could speak to  
23 the discrepancy between the survival data and the  
24 recurrence data, and what you might predict, and it's  
25 not fair to predict, but what you might predict could

1 be the survival effect of this treatment in the long  
2 run.

3 DR. O'DONNELL: Okay. I'll speak first to  
4 the idea of a systemic administration of  
5 chemotherapy. One of the things we've learned from  
6 at least the clinical experience is that chemotherapy  
7 when it is given appropriately for advanced bladder  
8 cancer, which includes disease still in existence in  
9 the bladder, is reasonably good at dealing with  
10 deeply invasive disease but actually is relatively  
11 poor at dealing with surface disease such as  
12 carcinoma in situ and papillary disease. The same  
13 can be said with radiation, and it may very well be  
14 the case that the blood supply to the areas closest  
15 to the mucosa that basically incorporates with the  
16 urine is a source of concentration gradient and so  
17 you lose that effect. So we don't really -- even  
18 when the case is when it's given for another reason,  
19 we don't really see an effect on superficial disease.

20 The second point is that there is actually  
21 another device-related treatment for enhancing  
22 mitomycin C administration. Actually, there have  
23 been two parallel kinds of studies. One was an  
24 attempt, in the U.S., to increase the dosage and give  
25 it better concentration gradient and with that, you

1 saw a corresponding increase in serum levels and  
2 improvement in mitomycin effect. So there seems to  
3 be a relationship between the ability to give a  
4 better concentration gradient, dissolution of the  
5 tissue and tissue levels were measured in that one  
6 and, in fact, to give it a clinical effect.

7           The same thing was seen from an  
8 electromotive therapy, use -- pheresis to try drive  
9 things through. Again, the same epi phenomenon of --  
10 well, I shouldn't say epi phenomenon. The same  
11 phenomenon of an increased serum level translating  
12 into an increased tissue level translating into a  
13 beneficial clinical effect. So these in aggregate  
14 would suggest that what you're getting is a tissue  
15 penetration effect, a better tissue penetration with  
16 a hyperthermia and the medication, rather not with  
17 one or the other by itself. The levels of mitomycin  
18 that are going in are nowhere near even those  
19 absorbed, the level that you would consider for  
20 systemic does given for other reasons like indication  
21 of pancreatic cancer or something like that.

22           And the second point was about survival.  
23 One thing to understand about the recurrence and  
24 survival rates of bladder cancer is that -- I showed  
25 you that recurrence curve which is an exponential

1 decay. The progression rate on the other hand is  
2 more of a linear phenomenon, and it's an order of  
3 magnitude different. And so in order to see any  
4 differences in survival on progression over time, you  
5 need large studies, highly risk event prone patients  
6 and it's extremely rare to see any data for  
7 progression and survival data in the literature. The  
8 only results that have really come out have been in  
9 large med analysis with 2500 patients, that have been  
10 able to see a small difference for BCG progression.

11 DR. TALAMINI: So then in follow-up, to ask  
12 it a slightly different way, there are certainly  
13 diseases where preventing recurrences does not help  
14 long-term survival and is therefore just torturing  
15 the patient in the interim. So what is the evidence  
16 in bladder cancer to the contrary?

17 DR. O'DONNELL: You have to ask yourself is  
18 there a benefit in preventing the patient from coming  
19 in multiple times to have a procedure under  
20 anesthesia, enduring the psychological harm of  
21 knowing they have cancer or cancer has recurred, and  
22 then you have the suggestion of maybe there is a  
23 benefit but we don't -- you saw even a difference in  
24 this group but the numbers were too small to save any  
25 statistical matters on it. So that will probably

1 remain an unanswered question but the other benefits  
2 of preventing the necessary treatments, the  
3 anesthesia, enduring the psychological aspects of  
4 having cancer, having to go in for recurrences and  
5 living under that umbrella, still remain.

6 DR. TALAMINI: Dr. Donatucci.

7 DR. DONATUCCI: To follow-up on what you  
8 just said because for me, I mean discussion really  
9 has to do with survival, and my knowledge of bladder  
10 cancer is limited, despite the fact I'm a urologist,  
11 because that's not what I do, but my memory is that  
12 patients with invasive disease generally tend to show  
13 up with it when they present. Patients with  
14 superficial disease usually remain superficial, and  
15 the hard part is trying to determine ahead of time  
16 which of the very few people are going to go from  
17 superficial to invasive disease. Now, we do know  
18 that patients with CIS and T1 disease high grade are  
19 the group that are likely to progress, and just to  
20 follow your thought one step further and take it back  
21 to the public representative who spoke earlier,  
22 that's the group that may actually have prophylactic  
23 cystectomy and that is -- we talk about adverse  
24 events and what the adverse event comparative issues  
25 are, but one adverse event we don't talk about is

1 taking a patient who's given a prophylactic  
2 cystectomy because of their high-risk of progression  
3 to muscle invasive disease, and there's not that many  
4 of them, but that is clearly a group of patients who  
5 would benefit.

6           But I do agree with what you were saying  
7 earlier that for most patients with superficial  
8 disease, the fact that you reduce recurrences does  
9 not translate into any appreciable diminution of risk  
10 of progression. And so we are talking about  
11 essentially what Dr. O'Donnell says, reducing the  
12 number of returns for cystoscopies and TURs, et  
13 cetera, et cetera, and it's a -- we're not talking  
14 survival. We're talking about just in psychologic  
15 and local symptomatology and then you have to say  
16 exactly whether the cure is worse than the risk of  
17 leaving them uncured.

18           DR. TALAMINI: Does the sponsor have a  
19 response?

20           DR. GROSSMAN: Yeah. I'd like to address  
21 the issue because it's really very, very critical.  
22 BCG, of course, is approved for recurrence. In the  
23 large Southwest Oncology Group trial, BCG maintenance  
24 versus no maintenance, there was a large significant  
25 difference in recurrence rate. There was a smaller

1 but still significant difference as with maintenance  
2 compared to no maintenance, in the rate of, I forget  
3 exactly the term, but the requirement to switch  
4 therapy to something else which was reasonably  
5 aggressive, largely cystectomy. So there were fewer  
6 rates of significant treatment failures with the  
7 maintenance arm, and there was even a trend. I think  
8 the P-value was .08 for survival which is really very  
9 impressive because the study was a modest size. All  
10 patients came in with non-muscle invasive disease and  
11 many of these deaths were actually from other causes,  
12 not all from bladder cancer.

13           So the odds are stacked against you for  
14 showing improvement in survival. So is preventing a  
15 recurrence good? Yes. Bladder cancer is the most  
16 expensive cancer to treat overall, and most of the  
17 costs are actually driven by recurrent disease and  
18 transurethral resections and repeated  
19 hospitalizations. This is an older population, and  
20 every time you go to the operating room and get an  
21 anesthetic, it's not good for you. There are small  
22 but there are real risks associated with  
23 transurethral resections. So preventing recurrence  
24 is a good thing. Prevent recurrence is a sure way of  
25 preventing progression, and if you had a large enough

1 study I think you would almost certainly show that  
2 aggressive intravesical treatment compared to  
3 transurethral resection alone is going to decrease  
4 recurrence, decrease progression, decrease survival.  
5 It hasn't been done just because the numbers are very  
6 small. They're stacked against you, and nobody has  
7 had large enough studies and especially with long  
8 enough follow-up to reliably answer that question.

9 DR. DONATUCCI: May I ask Dr. Kane to  
10 respond because slide 39 you presented data that said  
11 with superficial transitional cell carcinoma of the  
12 bladder, it's a heterogeneous population, outcomes  
13 vary, recurrence is common, progression is less  
14 common and not altered by therapy. So that's in  
15 direct contradiction to what you just said about  
16 preventing recurrence, prevent progression.

17 DR. O'DONNELL: Would you like me to answer  
18 that? I'm sorry, or address that.

19 DR. DONATUCCI: Somebody.

20 DR. O'DONNELL: So progression survival are  
21 the hardest terms to come to grip with because of the  
22 small event rate over time. As I mentioned, you need  
23 large studies, high-risk patients, for long periods  
24 of time.

25 Just to give you an example of how the

1 field is evolving, when the first guidelines came out  
2 in 1999 from the AUA and the 2002 guidelines from the  
3 EAU, progression was not even on the table because  
4 there essentially was not data. The data really came  
5 out just about two or three years ago with large meta  
6 analyses, mostly from the European ERTC Group which  
7 looked at -- there were two other meta analyses that  
8 backed it up, that looked at all accumulated data  
9 with BCG and discovered, for instance, that in  
10 studies that had BCG plus at least a year of  
11 maintenance, there was a real difference. It was a  
12 relative difference in progression, reduction by 27  
13 percent.

14 DR. DONATUCCI: All right. But what's the  
15 absolute number? Twenty-seven percent.

16 DR. O'DONNELL: It went from 13 percent to  
17 9 percent. It was small amounts.

18 DR. DONATUCCI: Right, but I mean that  
19 could be two patients.

20 DR. O'DONNELL: It depends on your sample  
21 size obviously.

22 DR. DONATUCCI: Right.

23 DR. O'DONNELL: So the point is that you  
24 can see a difference with aggressive therapy. To ask  
25 any study to demonstrate that requires you to make an

1 enormous investment in healthcare resources to ask  
2 that question.

3 DR. DONATUCCI: So what you're saying, if I  
4 understand it, is that frankly while this may occur,  
5 if it takes that many patients, it's not that big of  
6 problem.

7 DR. O'DONNELL: It's a --

8 DR. DONATUCCI: Not an individual  
9 perspective but from a healthcare perspective.

10 DR. O'DONNELL: Relatively speaking, it's a  
11 small percentage of patients.

12 DR. DONATUCCI: Relatively speaking, it's a  
13 small number of patients that progress from  
14 superficial to invasive disease.

15 DR. O'DONNELL: About 20 percent of those  
16 patients that present superficial disease that enter  
17 into the muscle invasive pool of patients, yes.

18 DR. TALAMINI: So having raised the  
19 question, let me rein the question in a little bit.  
20 In this device, we're looking at recurrence and not  
21 overall survival. So it was a little bit of an  
22 unfair question to ask but for me I thought it was  
23 important to put within the framework of the overall  
24 disease and how people do.

25 DR. O'DONNELL: Certainly.

1 DR. TALAMINI: Yes, sir.

2 DR. WITJES: Maybe I could just add some of  
3 the figures because I was co-author on that -- and  
4 the -- analysis was 4500 patients and we had a  
5 difference in progression of 12 versus 9 percent,  
6 11.8 percent. So indeed it's a limited number of  
7 patients that have progression and the relative  
8 reduction is around 25 percent.

9 DR. TALAMINI: Other questions for the  
10 Panel? Dr. Dahm?

11 DR. DAHM: I have two questions. One  
12 relates to understanding the device better. There's  
13 a rendering -- artistic rendering of the device,  
14 slide 24 of the sponsor's presentation that suggests  
15 that the catheter is kind of in a central position in  
16 the bladder but I would suspect, unless there's a  
17 specific mechanism to put it there, it actually in  
18 many patients will be lying on the posterior wall of  
19 the bladder and that would explain some of the side  
20 effects. Is that a correct understanding? So to a  
21 certain extent this picture may be misleading just  
22 with that regard.

23 DR. WITJES: Well, the misleading part is  
24 that the tip of the catheter is not as large as it's  
25 suggested here.

1 DR. DAHM: Yeah.

2 DR. WITJES: You put the catheter in and  
3 you empty the bladder. You pull the catheter out of  
4 the bladder so that you sort of position the catheter  
5 at the point you'd like --

6 DR. DAHM: Do you put traction on the  
7 catheter?

8 DR. WITJES: We don't put traction on, but  
9 it's quite a heavy catheter. So it's positioned at  
10 the bladder neck, and then obviously what you'll say  
11 is the tip of the catheter is close to the posterior  
12 wall. Now, you have some additional -- that's a  
13 mechanical irritation of the mucosa and of course,  
14 the radio frequency antenna is closest to that part.  
15 That explains the posterior wall --

16 DR. DAHM: Okay.

17 DR. WITJES: -- thermal reaction,  
18 absolutely.

19 DR. DAHM: Okay. The second question,  
20 we're asked to look at 101 as the pivotal study, and  
21 I took from Dr. O'Donnell's presentation as well as  
22 from Dr. Grossman's presentation that most people  
23 that we're targeting this device for with  
24 intermediate and high grade bladder cancer would  
25 be -- in this country would be getting BCG rather

1 than mitomycin C. And I notice that in your second  
2 study, that is your control. Would BCG not be a  
3 better control, a more appropriate control to look at  
4 the effectiveness of this, of this device?

5 DR. WITJES: Well, if you refer to Study  
6 101, yes, we realize that this study was thought of  
7 in the beginning of the nineties, when BCG was  
8 actually just at the table. The initial BCG studies  
9 are from that period.

10 DR. DAHM: Okay.

11 DR. WITJES: And I think the choice in  
12 those days when also mitomycin C was still advised  
13 for even progression in the EAU, AUA guidelines in  
14 1999. So I think as to the choice in those has been  
15 Synergo versus mitomycin C, and obviously when we  
16 spoke about the new trial in the beginning of this  
17 century, which was for BCG, so that's I think the  
18 development of the knowledge of how you treat bladder  
19 cancer.

20 DR. TALAMINI: Dr. Lippert.

21 DR. LIPPERT: I also wanted to ask about  
22 the posterior part because in that same diagram I see  
23 little prongs coming out interiorally, laterally and  
24 posteriorally. Is it -- so it's not really true that  
25 it's really lying on the posterior bladder wall that

1 you have more of a tissue effect on the bladder wall?  
2 I mean does it make sense that the whole bladder is  
3 heated equally?

4 DR. WITJES: Well, we have five  
5 thermocouples. Two are in the urethra and three,  
6 they're very floppy things. So they just press  
7 against the bladder wall. They don't damage that at  
8 all --

9 DR. LIPPERT: But they tend to go  
10 posteriorally?

11 DR. WITJES: The thermal antennas? No,  
12 they are from side, left side and right side. So  
13 they are let's say more or less in different spaces  
14 in the bladder. That's not only on the backside.  
15 Obviously we don't measure the temperature in the  
16 whole bladder. You have sort of an average of the  
17 five thermocouples you have, and indeed there are  
18 maybe at the spot, which is closest to the radio  
19 frequency wire, that you might have a little bit  
20 higher temperature.

21 DR. LIPPERT: So the posterior radio  
22 frequency wire is posterior?

23 DR. WITJES: Yeah, because the patients are  
24 attached to the device and they're laying down  
25 which --

1 DR. LIPPERT: So because they're supined,  
2 the wire, it falls down to the posterior bladder  
3 wall?

4 DR. WITJES: More or less, yeah, but that's  
5 not the intended. It's just so light that it is not  
6 due to gravity going -- it's not that all the  
7 thermocouples are in the back of the bladder. They  
8 are spread through the bladder.

9 DR. LIPPERT: Right.

10 DR. O'DONNELL: If I may, just a point of  
11 clarification. The microwave energy doesn't come  
12 from the thermocouples at all.

13 DR. LIPPERT: Right.

14 DR. O'DONNELL: They are temperature probes  
15 only. So the antenna is in the middle part of the  
16 catheter itself. So those aren't -- that is fixed by  
17 the catheter position which may by gravity or by  
18 position tend to go towards the posterior wall.

19 DR. LIPPERT: That's what I -- okay. And  
20 then the question he asked about the -- well, the  
21 literature review, you looked at BCG for comparison  
22 and you looked at mitomycin C, but you didn't look at  
23 any comparisons with BCG with Interferon.

24 DR. O'DONNELL: I've done a lot of work  
25 with BCG with Interferon and to date there has only

1 been one BCG+Interferon versus BCG randomized trial.  
2 We actually just presented it at this year's AUA, and  
3 we don't have the details of that. We haven't even  
4 sorted it out yet. I would just tell you that the  
5 result of that presentation was that there was no  
6 different BCG and BCG with Interferon for first time  
7 patients. So at the moment, we don't see that BCG  
8 with Interferon is a different or superior control  
9 versus BCG for another new therapy. That's all we  
10 can say at this point, even though for recurrences,  
11 BCG failures, we at least have the phase two data  
12 that we see reasonably good effects of  
13 BCG+Interferon.

14 DR. LIPPERT: But you're saying that would  
15 be a different patient group than this that we're  
16 looking at, intermediate and high-risk, because you  
17 said those are first time.

18 DR. O'DONNELL: Okay. So not really a  
19 different patient group in terms of the intermediate  
20 high-risk groups by the EAU criteria but those  
21 patients that we've at least focused most of the data  
22 on for BCG with Interferon, are patients that have  
23 failed BCG at least once, whereas the application  
24 here for the Synergo group is even for up front  
25 patients with intermediate and high-risk group as

1 even a first time treatment. And it does represent,  
2 you know, a new therapy option for the high-risk  
3 patients which up to now have only been BCG or BCG  
4 containing regimens.

5 DR. TALAMINI: Dr. Bhutani.

6 DR. BHUTANI: Yes, I have a few questions.  
7 One is on is there any data -- we're talking about  
8 survival. Is there any data before, for example, BCG  
9 was -- became involved, about the progression of  
10 transitional superficial bladder cancer in terms of  
11 mortality and the progression rate without treatment.  
12 Of course, at this time, it may not be ethical to do  
13 a study like this when there is treatment but when  
14 there were limited treatments out there, any older  
15 studies on this?

16 DR. O'DONNELL: The best study I'm aware of  
17 was for a relative uniform group of patients that had  
18 carcinoma in situ, and this is, as you all may know,  
19 one of the clear indications for BCG use in the  
20 United States. In fact, it's acknowledged by  
21 everyone essentially in the world, this is the number  
22 one treatment.

23 The studies from the Mayo Clinic was before  
24 the BCG era and they recorded a progression rate of  
25 about 7 percent per year regardless of when they got

1 their cystectomies. In general, the death rate is 50  
2 percent of the progression rate because half the  
3 patients that have muscle invasive disease die of it.  
4 That gives you just a general idea of what a rate  
5 would be.

6 Now, there have been recently an attempt to  
7 define risk tables of progression and recurrence  
8 based on the EORTCs, a risk calculator which is now  
9 available online, and that does give you an idea of  
10 what the expected progression rate is. It's  
11 generally about 5 to 10 times lower than the  
12 recurrence rate, and can give you a ballpark of what  
13 relative differences in the order of magnitudes are  
14 between recurrence and progression.

15 DR. BHUTANI: Thanks.

16 DR. GROSSMAN: The other point is that T1  
17 disease is associated with the 30 percent risk of  
18 progression to muscle invasive disease. So that's  
19 not high enough to say that everybody that walks in  
20 the door with T1 should have an immediate cystectomy,  
21 but it's enough to say they deserve an aggressive  
22 course of therapy, and if they don't respond promptly  
23 and still have T1 disease, then you should have a  
24 very short fuse for proceeding with a cystectomy, not  
25 continuing with therapy indefinitely.

1 DR. DONATUCCI: Let me just ask you further  
2 about that. So you give an aggressive course of  
3 therapy. They don't have recurrence disease, but is  
4 there evidence that actually, other than preventing  
5 the recurrence when you go back and biopsy, does the  
6 percent progress change? Does the percentage change?  
7 Is it still 30 percent? In other words, you just  
8 said 30 percent of T1 disease progress.

9 DR. GROSSMAN: Yes.

10 DR. DONATUCCI: And therefore they deserve  
11 aggressive therapy.

12 DR. GROSSMAN: Yes.

13 DR. DONATUCCI: And you reduce their  
14 percentage of recurrence but does that actually -- is  
15 there data now to say and then subsequently you have  
16 also reduced their chance of progression from 30  
17 percent to a lower number?

18 DR. GROSSMAN: There's data with T1  
19 patients in this series. There's data with T1  
20 patients in other series which show that BCG is  
21 better than mitomycin for decreasing the recurrence  
22 rate of T1 disease. Similarly, with this trial which  
23 included high-risk patients, many of which had T1,  
24 again treatment with Synergo was better than  
25 mitomycin in preventing disease recurrence but having

1 large numbers of T1s in a single series showing that  
2 treatment is better than no treatment, I don't think  
3 that has been looked at. So it's basically one  
4 treatment versus something else.

5 DR. DONATUCCI: I don't want to leave you  
6 with the impression that I'm not in favor of treating  
7 these patients. I'm just trying to point out that  
8 the data, the evidence just isn't there. And I'm not  
9 saying you shouldn't treatment but the evidence isn't  
10 there.

11 DR. GROSSMAN: Right, and I wish the data  
12 was stronger but that's why there's been several meta  
13 analyses done which show that BCG does prevent  
14 progression, and the problem is you need large  
15 numbers of patients and long follow-up and even big  
16 large randomized trials frequently don't address all  
17 the questions that we like to rigorously address, and  
18 I'm as skeptical about meta analysis as anybody else  
19 but frequently that's what we have to resort to, to  
20 try and answer some of these difficult questions  
21 which haven't been rigorously addressed by the  
22 existing literature.

23 DR. TALAMINI: Yes, sir.

24 DR. BHUTANI: My follow-up question is  
25 again regarding, we've talked about survival data or

1 lack thereof in superficial bladder cancer. What I'd  
2 like to understand is short of survival, 30 percent  
3 of people will likely progress to muscle invasive  
4 bladder cancer. Then they may need cystectomy or  
5 other therapies that are more aggressive. I'd like  
6 to understand the morbidity to the patient when they  
7 have a muscle invasive cancer, their quality of life,  
8 sexual, urinary function and so on and what happens?

9 DR. TALAMINI: Just to clarify. Do you  
10 mean between these two therapies, cystectomy versus  
11 what's being --

12 DR. BHUTANI: In the sense that we are  
13 talking about whether we should treat superficial  
14 bladder cancer at all because we, you know, BCG  
15 they're saying prevents progression but it doesn't  
16 affect survival or we don't know if it affects  
17 survival but if a patient say, didn't have BCG and  
18 his tumor or her -- his tumor would be more likely to  
19 progress, then what would be the effect without any  
20 treatment on a patient with superficial bladder  
21 cancer and morbidity and quality of life?

22 DR. O'DONNELL: Okay. So getting back to  
23 the survival issue, again we're encumbered by the  
24 fact that it's a small event rate. I will tell you  
25 that this year's AUA and at the European Association

1 National Meeting, there actually was for the first  
2 time data presented on a large randomized trial of  
3 Epirubicin and TURBT versus BCG and TURBT and there  
4 was a survival advantage. Now, I know you talk  
5 about, well, how many percentage of patients does  
6 that affect, and my answer there really would be that  
7 it's not so much the absolute number of patients,  
8 it's the relative benefit you may get from it. So I  
9 would still say that if you look at chemotherapy  
10 trials and so forth that, you know, a relative  
11 benefit of 25 percent is considered clinically  
12 significant, and that's the level at least we're  
13 working on with progression.

14 I'd also make a point that, you know, of  
15 older folks and superficial bladder cancer, let's  
16 take a more common cancer like prostate cancer and we  
17 have been at it for, you know, a decade with hundreds  
18 of thousands of patients getting radical  
19 prostectomies. We are just now beginning to see in  
20 our national statistics a small decrease in the death  
21 rate of prostate cancer. It's too early to say  
22 whether radical prostectomy is responsible for that,  
23 but that's the kind of level of evidence you're  
24 asking for, to get these kinds of answers. And so,  
25 you know, the real bottom line answer is that

1 recurrence is a major problem and when patients do  
2 progress to bladder cancer, they have their bladders  
3 removed, most of them are impotent, if they're men,  
4 over 50 percent. You have to wear a bag. At least  
5 two-thirds of those that receive a cystectomy in the  
6 United States are wearing an appliance and a bag.  
7 Granted, they adapt to that. There is a 30 percent  
8 complication rate in the early 3 months. There's  
9 another 30 percent later on. There's a 2 to 3  
10 percent mortality in high volume centers. It's  
11 double that in the community center. All of these  
12 are very major issues for patients to undergo. So  
13 even if you could prevent a 25 percent progression  
14 relative, I would say that that -- and that  
15 translates into only half the benefit for survival,  
16 that's still a benefit for the whole group. Again,  
17 relative is still important for the disease  
18 population that you're focusing on.

19 DR. TALAMINI: So let me just take the  
20 privilege of the Chair and again sort of try and draw  
21 back from the discussion of long term survival  
22 because that's not the question we're being asked  
23 today, and we do have a lot of things that we need to  
24 cover about what we're being asked today. So with  
25 that in mind, Dr. Marcovich.

1 DR. MARCOVICH: This is more of a  
2 methodology question, another one with the practice  
3 patterns were at the time this started, but would any  
4 of these patients have gotten a single mitomycin at  
5 the time of resection and, if so, how would that have  
6 affected the outcomes?

7 DR. WITJES: The answer very simply is no  
8 because at the time this study was started, the 101,  
9 the 101 -- was not an item at all and still is not in  
10 the United States.

11 DR. TALAMINI: Dr. Redman.

12 DR. GROSSMAN: Could I just supplement  
13 that? The AUA has looked at that issue with the  
14 current guidelines. What's not so clear in the  
15 published literature but is very clear when you talk  
16 to the investigators who have done the studies is  
17 that a single dose of mitomycin is optimal for low  
18 risk patients. It is not nearly as effective in  
19 preventing recurrence in intermediate and high-risk  
20 patients. So this would not be the target population  
21 where a suitable dose of intravesical chemotherapy is  
22 likely to afford a long-term response.

23 DR. TALAMINI: So again, let me take the  
24 privilege of the Chair because we're just about out  
25 of time for generalized discussion, and there are two

1 issues still on the table that I think we need to get  
2 settled if we can. One is the issue of cystoscopy  
3 and the study and it sounded before lunch like there  
4 was a difference of opinion or dispute, and I wonder  
5 whether between the FDA and the sponsor, that still  
6 exists or if that's a settled question. So that's  
7 one thing about whether the cystoscopies were, in  
8 fact, done by the same investigators in the study.  
9 That's number one.

10           And then, number two, is this issue of the  
11 post-market study. It sounds as if the sponsors now  
12 have proposed a different plan and since that's one  
13 of the specific questions that we're going to  
14 address, I wonder if the Panel feels that we need to  
15 know a little bit more about that. I do. So perhaps  
16 the sponsors could address those two issues.

17           MS. STERN: I'll address the first issue,  
18 but I'm basically just going to repeat what I said  
19 before that all the cystoscopies were done by the  
20 principal investigators at each of the  
21 investigational sites. I don't know where the FDA  
22 got other information but maybe they can provide  
23 that.

24           DR. TALAMINI: Response from the FDA?

25           DR. KANE: I did not bring all of my

1 supporting documentation with me. So I'm not able to  
2 give you a numerical answer to that question. It --  
3 I looked over lunchtime through the published report  
4 and through the protocol and at least one of these  
5 sources indicated that patients had the option to  
6 have procedures performed locally in follow-up, but  
7 more than that, I can't give you a specific answer at  
8 this time.

9 I did find an additional reference to  
10 Dr. Redman's question about what kind of vital  
11 information might be missing and in the sponsor's  
12 submission of January of this year, for example, and  
13 I think this is page 13 where the comment is made  
14 that pathology reports were not available on four  
15 patients of the recurrence, the comment is made that  
16 these reports were not transcribed to the case report  
17 form and therefore were not available for later  
18 follow-up, and that would be one example of  
19 information that was not able to be carried forward.

20 DR. TALAMINI: Okay. Thank you.

21 DR. KANE: I think also that you may have  
22 perceived somewhat different numerical consequences  
23 of this issue of progression, and I think you have to  
24 sort out patients by their baseline risk group. It  
25 may be 30 percent of T1s, but it's not anywhere near

1 that for Tas. Many more patients typically present  
2 with Tas. So we're focusing on a subgroup of the  
3 overall group. The study, as information is  
4 available, did not demonstrate a difference in  
5 progression. What type of sample size or population  
6 might be necessary to be able to examine that  
7 question, I think that's a whole different question  
8 to open the door for, that we don't have available  
9 for you.

10 DR. TALAMINI: And the post-market study.  
11 Sorry, did you have a response to that?

12 DR. WITJES: Well, we've been thinking  
13 about the different locations, and we thought the  
14 same, that sometimes patients are within one hospital  
15 looked at, at the OR, and then the other time at the  
16 outpatients which is the different location but the  
17 same hospital. And with regard to the missing  
18 pathology, there were two finally patients with  
19 missing pathology, one in the Synergo arm and one in  
20 the mitomycin C arm.

21 DR. TALAMINI: And that post-market study  
22 changes or differences from what was in the original  
23 packet?

24 DR. O'DONNELL: So as you probably gathered  
25 from the description of the older study, it was

1 relatively convoluted with non-inferiority for eight  
2 different variables with highly prevalent adverse  
3 events and low prevalence adverse events and was  
4 attempting to compare it with historical control  
5 which was the proper control. We realize that that  
6 was going to be untenable, and we agree that it's  
7 really an unworkable situation and provides data that  
8 you have to wonder what it's really telling you. And  
9 so we envisioned doing a much more simplified study  
10 where we get a reasonable number of patients and  
11 provide the kinds of rigorous adverse event reporting  
12 so that it can go into the labeling, examine,  
13 whatever, the actual usage among the subgroups.  
14 Where it would be used in the U.S. would again  
15 provide very relevant patient population data so that  
16 would, in fact, make it clear where this treatment is  
17 indeed working. And I want to reiterate the sponsor  
18 is very open to working with the FDA and with the  
19 Panel to work on the details of this, to provide the  
20 proper data acquisition so that you get the data that  
21 you think is critical for this kind of study, and  
22 that we remain very cooperative and very open to this  
23 kind of suggestion. We just presented what we think  
24 is a reasonable and very simplified, very pertinent  
25 type of post-marketing strategy, post-marketing

1 study, and we take that as basically the basis from  
2 which to go further on.

3 DR. TALAMINI: Okay. Dr. Redman.

4 DR. REDMAN: Yeah, this question I didn't  
5 think was relevant, but now I'm thinking it is  
6 relevant because I thought I had known that that's  
7 what I'm going to ask Dr. Grossman and Dr. O'Donnell  
8 and maybe the AUA has this information. I read the  
9 AUA Guidelines. This is mitomycin C versus BCG as a  
10 standard of care in this country for T1, Ta. You  
11 know, my sense was from my group of urologists and  
12 Southwest Oncology Group that it's not 100 percent  
13 BCG use, because if it is, this becomes relevant.  
14 That's why I'm asking.

15 DR. GROSSMAN: No, it is not. I still  
16 treat patient with mitomycin. For the really high-  
17 risk patients, the AUA Guidelines recommend BCG and I  
18 agree, BCG is appropriate currently for T1 and for  
19 carcinoma in situ. For other patients, mitomycin is  
20 a reasonable choice, and I do use mitomycin C in some  
21 of my patients.

22 DR. REDMAN: Do we get data from the AUA or  
23 Medicare use or anything about what is -- I mean is  
24 it known?

25 DR. O'DONNELL: In the United States, BCG

1 is used overall two to one versus all forms of  
2 chemotherapy, and mitomycin is the most common form  
3 of chemotherapy used in the United States. Just to  
4 clarify what Dr. Grossman said is that the AUA and  
5 EUA Guidelines specify the use of mitomycin as a  
6 course of therapy as was used in the control arm for  
7 patients with intermediate risk disease or BCG. So  
8 that's where the option is, and the reason for that  
9 is if you fail mitomycin, you often will respond to  
10 BCG, whereas the converse is not necessarily the  
11 case. But for the high-risk group, which is the  
12 highest risk of progression, only BCG or cystectomy  
13 is currently in the Guidelines from both national  
14 groups.

15 DR. TALAMINI: We are running dangerously  
16 close to being short on time. Dr. Lippert, another  
17 discussion question?

18 DR. LIPPERT: A quick question. In your  
19 pivotal study comparing heated versus unheated  
20 mitomycin treatments, was there a difference in serum  
21 white count?

22 DR. KOREN: We did not have any case of  
23 mito suppression in 101 and as mentioned before, the  
24 levels are far too low to cause anything which was  
25 described in this.

1 DR. LIPPERT: You mentioned the levels, but  
2 you never mentioned the white count.

3 DR. KOREN: No mito suppression, no.

4 DR. LIPPERT: Thank you.

5 DR. TALAMINI: So, Dr. Connor, I saw you  
6 had your hand up with a question, but I would also  
7 ask you to serve the Panel as we move towards  
8 beginning to address these specific questions being  
9 the statistical expert on the Panel, if you could  
10 discuss your opinion of the data that you've seen  
11 thus far with respect to the statistical issues.

12 DR. CONNOR: So, and this gets to the broad  
13 question of efficacy versus adverse event rate versus  
14 biases that are potentially here. Okay. So I think  
15 one thing that was brought up by FDA that I think if  
16 I were sitting in your chair, I would have said, so  
17 I'm going to say something minor, defend is too  
18 strong a word, but it was brought up that seven  
19 patients randomized to the treatment group didn't  
20 receive therapy and only one patient randomized, the  
21 control didn't receive the therapy but because of  
22 this mix up then, in fact, those numbers are really  
23 five versus three which isn't quite so bad. I'm  
24 sorry, didn't receive therapy but withdrew from the  
25 study. So seven versus one withdrew from the study.

1 If we look at how this switching occurred, it's five  
2 treatment versus three control withdrew from the  
3 study and if we actually look at patients who  
4 received at least one treatment course, only two  
5 patients randomized to your treatment versus three on  
6 the control group, elected to remove themselves from  
7 the study. So I think that speaks positively at  
8 least that there's not a difference in patients who  
9 didn't like being treated by this. So I think that  
10 didn't come up clearly, but something I was able to  
11 count from this.

12 I think some of the biases, and this was  
13 the question I asked, that I was told was our  
14 question to answer, not FDA's question, you know,  
15 there are definitely biases here. This isn't the  
16 ideal trial, but I've never seen an ideal trial.  
17 This may be more or further from the ideal trial than  
18 we would like which is why we're a Panel in the first  
19 place. But the treatment effect is pretty big, and I  
20 think that definitely helps overcome some of these  
21 biases. So I think I would punt the question you  
22 asked to me to the clinicians on the Panel and ask in  
23 particular if you did the post-approval study, and I  
24 think that the way that you're proposing to do it is  
25 better than the first way because I think the fairest

1 thing you can do is identify the risks and precisely  
2 quantify the risks which your study does much better  
3 than the first study in a way that clinicians are  
4 going to understand. So the clinicians and their  
5 patients can discuss and say, here are the risks in a  
6 clear way and here are the benefits.

7           So I think my key question that I have to  
8 ask of the other clinicians up here is what is the  
9 risk or what are the cons of approving a treatment,  
10 and I ask this without alluding to what my beliefs  
11 are, but is the risk or the cons of approving a  
12 treatment that looks like it works with quantified  
13 risks, some of which we have seen, and it sounds like  
14 we all have better quantified risks, so that doctors  
15 and patients can discuss a treatment that is  
16 available to them and decide if that treatment is  
17 right to them, and I say that because it looks like  
18 the efficacy is pretty good and understanding that  
19 it's smaller than we think, but to me personally  
20 probably not. The biases aren't enough to overcome  
21 what looks to be efficacious in particular and  
22 improves quality of life, if not survival, which I  
23 understand is important. So what would the risk be  
24 of letting this get to market and letting patients  
25 and doctors make their own decisions? So I think

1 that to me is the bigger question right now than any  
2 huge statistical concern I have that would be a  
3 problem.

4 DR. TALAMINI: And I would remind the Panel  
5 that we are instructed to only consider this based  
6 upon the data that we have and not data that we would  
7 get from a post-market study.

8 So with that discussion, I'd like to move  
9 towards addressing the specific four questions that  
10 the FDA has put before us. Each of you has a copy of  
11 them in written form in your packet, and if we could  
12 please show the first question up on the slide.

13 So I'm going to read from what's in front  
14 of me because I can't look back and talk in the  
15 microphone. I trust that they're the same.

16 Question Number 1, Effectiveness. Under 21  
17 C.F.R. 860.7(e)(1), effectiveness is defined as  
18 reasonable assurance, based upon valid scientific  
19 evidence, that, in a significant portion of the  
20 population, the use of the device for its intended  
21 uses and conditions of use, when accompanied by  
22 adequate directions for use and warnings against  
23 unsafe use, will provide clinically significant  
24 results.

25 The PMA for the Synergo system presents

Free State Reporting, Inc.  
1378 Cape Saint Claire Road  
Annapolis, MD 21409  
(410) 974-0947

1 clinical data from a single pivotal trial, Study 101,  
2 and several additional supporting clinical data  
3 sources, the review of which present the following  
4 challenges in assessing the effectiveness of this  
5 combination product.

6           Significant limitations exist in the design  
7 and conduct of Study 101.1, which when considered  
8 collectively, potentially impair the ability to  
9 interpret the results and increase the uncertainty of  
10 the state conclusions including multiple sources of  
11 bias; the absence of structured methods for obtaining  
12 pathology samples and recording pathology  
13 information; potential variation in mitomycin C  
14 exposure between the study groups; and reliance on a  
15 small, limited study population to perform the  
16 risk/benefit analysis and generalize the study  
17 results to the general U.S. population.

18           The supporting clinical data sources were  
19 not designed to provide stand-alone evidence of the  
20 safety and effectiveness of the Synergo system for  
21 the proposed indication.

22           So the question then, considering the  
23 design and conduct of Study 101.1 and the supporting  
24 clinical data sources, please discuss whether the  
25 clinical data in the PMA provide reasonable assurance