

FDA Advisory Panel

Drug-Eluting Stent Safety Profile

Stent Thrombosis



Campbell Rogers, M.D.

Chief Technology Officer



Agenda



CYPHER® Stent Benefits

Campbell Rogers, MD

Chief Technology Officer
Cordis Corporation

Safety Profile: Death and MI

Dennis Donohoe, MD

Worldwide Vice President
Clinical Research & Regulatory Affairs
Cordis Corporation

Stent thrombosis


Laura Mauri, MD, MSc

Chief Scientific Officer
Harvard Clinical Research Institute



Conclusions

Campbell Rogers, MD

Freedom From TLR Through 4-Years



Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS




Catheterization and Cardiovascular Interventions 64:163-168 (2005)

Coronary Artery Perforation by Cutting Balloon Resulting in Dissecting Subepicardial Hematoma and Avulsion of the Vasculature

Vu H. Quan,¹ MD, James R. Stone,² MD, PhD, Gregory S. Couper,³ MD, and Campbell Rogers,^{1*} MD

Cordis
a Johnson & Johnson company




Hematoma

Hemorrhagic Infarct

Restenosis and Late Mortality


Consecutive Patients



**2,272 patients undergoing BMS
implantation between 1992-1996**

Schühlen H et al. Am Heart J 2004;147:3174.


Diabetics



**603 diabetic patients undergoing
PTCA between 1987-1995**


Van Belle E et al. Circulation 2001;103:1218.

Clinical Presentation of BMS Restenosis



**10% of BMS
In-Stent Restenosis Cases
Presented as an MI**

Nayak AK., et al., *Circ* 2006;70:1026-9.



Clinical Consequences of Acute Coronary Syndrome in BMS Restenosis

Acute Coronary Syndrome May Occur With In-Stent Restenosis and Is Associated With Adverse Outcomes (The PRESTO Trial)

Abid R. Assali, MD^a, Ali Moustapha, MD^b, Stefano Sdringola, MD^{d,*}, Ali E. Denktas, MD^c, James T. Willerson, MD^d, David R. Holmes, Jr., MD^c, and Richard W. Smalling, MD, PhD^d


(Am J Cardiol 2006;98:729-733)

Table 3
Clinical events


| | Stable Angina (n = 617) | ACS (n = 824) | p Value |
|---|----------------------------|------------------|--------------|
| In-hospital events | | | |
| Any event | 4 (1%) | 3 (0.3%) | 0.4 |
| Myocardial infarction | 2 (0.3%) | 2 (0.2%) | 0.66 |
| Target vessel revascularization | 2 (0.3%) | 1 (0.1%) | 1.00 |
| Death | 0 (0%) | 0 (0%) | N/A |
| Follow-up events | | | |
| Any event | 135 (22%) | 285 (35%) | <0.001 |
| Myocardial infarction | 13 (2%) | 16 (2%) | 0.83 |
| Target vessel revascularization | 130 (21%) | 271 (33%) | <0.001 |
| Death | 3 (0.5%) | 18 (2%) | 0.008 |
| Restenosis defined by $\geq 50\%$ narrowing* | 42 (42%) | 78 (56%) | 0.043 |
| Restenosis defined by $\geq 50\%$ loss of gain* | 50 (51%) | 90 (64%) | 0.033 |

*Restenosis measurements are from a subset of the angiographic sub-study patients (n = 99, stable angina group; n = 140, ACS group).

Any Event Before Nine Months by Unstable Angina Status




History of Stent thrombosis



1. Schatz et al. Circulation.1991;83:148; 2. Fischman et al. N Engl J Med. 1994;331:1496; 3. Colombo et al. Circulation.1995;91:1676; 4. Schömig et al. Circulation.1994;90:2716; 5. Leon et al. N Engl J Med. 1998;339:1665;

Clinical DES Data

- Data accurately analyzed and reported
 - Need for independence and transparency



Independent

All analyses shown done by HCRI

Clinical DES Data

- Data accurately analyzed and reported
 - Need for independence and transparency
- A common set of definitions
 - ARC definitions

Clinical DES Data

- Data accurately analyzed and reported
 - Need for independence and transparency
- A common set of definitions
 - ARC definitions
- Patient based clinical end-points
 - Vigilance for safety signals

Clinical DES Data



- Data accurately analyzed and reported
 - Need for independence and transparency
- A common set of definitions
 - ARC definitions
- Patient based clinical end-points
 - Vigilance for safety signals
- Opinions of experts with deep domain expertise

Ralph D'Agostino, Ph.D.

Harvard/BU

Director Biostatistics

Elazer Edelman, M.D., Ph.D

Harvard/MIT

Vascular biology and DES

Daniel Simon, M.D.

Case Western

Antiplatelet therapies

Frederic Resnic, M.D., MSc


Harvard

Registry and Outcomes

Evidence on CYPHER[®] Stent




1. Wide range of levels of evidence will be seen by panel



Evidence on CYPHER[®] Stent


1. Wide range of levels of evidence will be seen by panel
2. Need to consider each device as a discrete entity
 - Different drugs
 - Different polymers
 - Different rate and duration of drug delivery



Evidence on CYPHER[®] Stent



- 1. Wide range of levels of evidence will be seen by panel**
- 2. Need to consider each device as a discrete entity**
 - Different drugs
 - Different polymers
 - Different rate and duration of drug delivery
- 3. Strong evidence of Safety and Efficacy of the CYPHER[®] Stent**
 - Wide variety of clinical settings
 - Data from more 45,000 patients enrolled in clinical trials



Overview



- Patient level pooled analysis of RCTs in patients treated with the CYPHER[®] Stent Sirolimus-eluting Stent vs. BMS through 4-year follow-up demonstrates no significant differences in Death or MI
- Patients treated with CYPHER[®] Stent and BMS have a similar overall risk of stent thrombosis over 4-years
- Although early, late and very late events occur in both arms, there are more events before year one for BMS and more events after year one for CYPHER[®] Stent
- Commitment to continued research and education to improve patient outcomes

Review Of Safety Data

Dennis Donohoe, MD
Vice President, Clinical Research & Regulatory Affairs



Patient Level Pooled Analysis of 4 RCTs




| | RAVEL* | SIRIUS* | C-SIRIUS | E-SIRIUS |
|--------------------------------|--|---|---|----------------------------------|
| Study Type | Prospective, Multi-Center, Blinded, Randomized | | | |
| # of Patients | 238 (120 CYPHER®, 118 BMS) | 1,058 (533 CYPHER®, 525 BMS) | 100 (50 CYPHER®, 50 BMS) | 352 (175 CYPHER®, 177 BMS) |
| Lesion Type | Single <i>de novo</i> lesion in native coronary artery | | | |
| RVD | ≥2.5 to ≤3.5 mm | | ≥2.5 to ≤3.0 mm | |
| Lesion Length | Lesion had to be covered with a single 18 mm stent | 15 to 30 mm in length coverable with 2 stents | 15 to 32 mm in length coverable with 2 stents | |
| Aspirin | Indefinitely | | | |
| Clopidogrel or Ticlopidine | 2 months | 3 months | 2 months | 2 months |
| Compliance to 4-year follow-up | CYPHER® – 94.2% BMS – 94.1% | CYPHER® – 96.8% BMS – 97.0% | CYPHER® – 98.0% BMS – 98.0% | CYPHER® – 97.1% BMS – 98.3% |

Summary of Safety Events

- Analysis of 4-year follow-up data across the 4 RCTs was conducted
- No significant differences noted in death and non-fatal MI combined, death or MI rates
- Multiple subgroup analyses were conducted and a significant finding was noted only in the diabetic subgroup

Cumulative Incidence of Myocardial Infarction*: 0 – 1,440 Days (4-Years)




| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 847 | 832 | 807 | 779 | 741 |
| Bx Velocity | 870 | 837 | 824 | 806 | 782 | 740 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

* Non -Q wave CK levels greater than 2 times normal with elevated CKMB


Cumulative Incidence of Death: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 870 | 863 | 842 | 817 | 776 |
| Bx Velocity | 870 | 863 | 857 | 843 | 824 | 781 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS


Cumulative Incidence of Cardiac Death: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 870 | 863 | 842 | 817 | 776 |
| Bx Velocity | 870 | 863 | 857 | 843 | 824 | 781 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Cumulative Incidence of Non-Cardiac Death: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 870 | 863 | 842 | 817 | 776 |
| Bx Velocity | 870 | 863 | 857 | 843 | 824 | 781 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

4-Year Mortality in CYPHER® Stent Trials



| Mortality | CYPHER® Stent | Control | p-Value |
|----------------------|------------------|-----------------|--------------|
| Cardiac | | | |
| RAVEL (n=238) | 3 (2.7%) | 5 (4.6%) | 0.493 |
| SIRIUS (n=1058) | 17 (3.4%) | 12 (2.4%) | 0.452 |
| C-SIRIUS (n=100) | 1 (2.0%) | 1 (2.0%) | 1.00 |
| E-SIRIUS (n=352) | 8 (4.7%) | 5 (2.9%) | 0.409 |
| Non-Cardiac | | | |
| RAVEL (n=238) | 10 (8.9%) | 2 (1.9%) | 0.034 |
| SIRIUS (n=1058) | 15 (3.0%) | 13 (2.6%) | 0.849 |
| C-SIRIUS (n=100) | 1 (2.0%) | 2 (4.1%) | 1.00 |
| E-SIRIUS (n=352) | 2 (1.2%) | 5 (2.9%) | 0.449 |


Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Exploratory Post-Hoc Subgroup Analyses



- Multiple subgroups within these 4 RCTs were evaluated for mortality and MI rates
 - Multiple stents, overlapping stents, long lesions (>20mm), small vessels (2.5mm stent), diabetics and non-diabetics
- There were no pre-specified subgroup analyses or sub-randomizations in these four trials
- All subgroups demonstrated similar mortality and MI rates except for the diabetic subgroup

Summary of Contemporary Diabetic Mortality Data



*Lee T et al., Am J Cardiol, 2006; 98:718-721

\$DIABETES: Sabaté M., et al., ESC 2006; Oral Presentation.
 SCORPIUS: Baumgart D., et al., TCT 2006; Oral Presentation.
 Letter from Don Baim, M.D.

DECODE: Chan C., et al., AHA 2005; Oral Presentation.
 4 RCTs CYPHER: Internal Data, Cordis Corporation.

Diabetic Mortality

- The difference at 4-years in post hoc subgroup analysis appears to be an anomalous finding
 - This difference was driven by the SIRIUS study ($P=0.037$) and at 5-years the difference no longer significant ($P=0.21$)
 - BMS-treatment group had an unusually good outcome compared to published data
 - A Cox regression analysis of these data indicates that presence of diabetes **lowered** the risk of death
- Three randomized prospective trials in diabetics did not support this observation

Conclusions: Safety

- Patient level pooled analysis of randomized controlled multi-center studies in patients treated with the CYPHER® Sirolimus-eluting Stent vs. BMS through 4-year follow-up demonstrates no significant differences in:
 - Death (cardiac and non-cardiac)
 - MI
 - Death and non-fatal MI

Pooled Analysis Of Stent Thrombosis Using The Academic Research Consortium (ARC) Definition

Laura Mauri, MD, MSc

Chief Scientific Officer

Harvard Clinical Research Institute, Boston, MA

Stent Thrombosis




- Academic Research Consortium Definition
- Results of Blinded Adjudication
 - Overall results
 - Relationship of target lesion revascularization to subsequent stent thrombosis
 - Relationship of stent thrombosis to clinical endpoints (death, myocardial infarction)
 - Relationship to antiplatelet therapy
- Conclusions

Is There An Optimum Definition For Stent Thrombosis



- **Protocol Definition**
 - Most restrictive
 - Excluded intervening TLR, potential bias
- **Definite ARC**
 - More reliable to distinguish mechanism, but may miss some stent thrombosis events
- **Possible ARC**
 - Driven by unexplained death
 - Best captured by clinical endpoints
- **Definite + Probable ARC**
 - Balances sensitivity and specificity


Cumulative Incidence of Stent Thrombosis (Protocol Definition): 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 866 | 858 | 836 | 811 | 770 |
| Bx Velocity | 870 | 860 | 853 | 839 | 820 | 777 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS


Cumulative Incidence of Definite ARC Stent Thrombosis: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 867 | 859 | 837 | 811 | 769 |
| Bx Velocity | 870 | 861 | 853 | 838 | 818 | 775 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Cumulative Incidence of Definite or Probable ARC Stent Thrombosis: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 866 | 858 | 835 | 809 | 768 |
| Bx Velocity | 870 | 856 | 848 | 834 | 813 | 772 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS


Thrombosis Incidence Analysis: ARC Definite or Probable



| ARC Definite or Probable Stent Thrombosis | SES (N=878 Patients) | BMS (N=870 Patients) |
|--|-------------------------|-------------------------|
| Acute Thrombosis (0-1) | 0.0% (0/878) | 0.0% (0/870) |
| Sub Acute Thrombosis (2-30) | 0.4% (4/877) | 0.3% (3/870) |
| Late Thrombosis (31-360) | 0.1% (1/874) | 1.0% (8/865) |
| Very Late Thrombosis (361-1440) | 0.9% (8/848) | 0.5% (4/843) |
| Any Thrombosis (0-1440) | 1.5% (13/848) | 1.8% (15/843) |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS





Cumulative Incidence of Any ARC Stent Thrombosis: 0 – 1,440 Days (4-Years)



| # Entered | 0 D | 180 D | 360 D (1-yr) | 720 D (2-yr) | 1080 D (3-yr) | 1440 D (4-yr) |
|-------------|-----|-------|--------------|--------------|---------------|---------------|
| Sirolimus | 878 | 866 | 858 | 835 | 809 | 768 |
| Bx Velocity | 870 | 856 | 848 | 834 | 813 | 772 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Cumulative Incidence of Stent Thrombosis to Latest Follow-up (4-5 Years, 4 Trials)



— SES - - - BMS

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Prior TLR and Stent Thrombosis: 4-Year Data



| | CYPHER® Stent ST N=30 | BMS ST N=28 |
|--|--------------------------|----------------|
| Prior TLR | 0 | 10 |
| Time from TLR to ST, days, median (min,max) | N/A | 269 (15, 1141) |
| ST Classification | | |
| Definite/Probable (N=13 vs. 15) | 0 | 6 |
| Possible (N=17 vs. 13) | 0 | 4 |
| TLR Procedure | 0 | 10 |
| BMS only | 0 | 0 |
| Any SES | 0 | 1 |
| PTCA only | 0 | 2 |
| Brachytherapy/PTCA | 0 | 5 |
| Brachytherapy/BMS Stent | 0 | 2 |

Pooled Data from RAVEL, SIRIUS, E-SIRIUS, and C-SIRIUS

Relationship of Stent Thrombosis (Definite and Probable) to Clinical Endpoints



| | SES (N=13) | BMS (N=15) |
|-----------------------|---------------|---------------|
| Death | 4 | 5 |
| Myocardial Infarction | 13 | 13 |
| Fatal MI | 4 | 4 |
| Q Wave MI | 8 | 5 |
| Non-Q Wave MI | 5 | 8 |

Similar mortality observed for SES and BMS thrombosis

Antiplatelet or Warfarin Therapy at Time of Event for Patients with Stent Thrombosis (0-1,440 days)




| | SES | BMS |
|---|----------------|----------------|
| Patients with <u>Any ARC</u> Stent Thrombosis | 30 | 28 |
| Yes | | |
| ASA, Clopidogrel/Ticlopidine | 27% (8) | 43% (12) |
| ASA only | 37% (11) | 18% (5) |
| Clopidogrel only | 1 | 0 |
| Coumadin | 1 | 1 |
| No* | 6% (2) | 11% (3) |
| <i>*Stopped at 7, 17, 58, & 403 days prior to ST, with 1 unknown stop date)</i> | | |
| Unknown** | 23% (7) | 25% (7) |
| <i>** Last known dose between 82 – 1325 days prior to ST – 5 unknown stop dates</i> | | |

Stent Thrombosis



- Patients treated with CYPHER® Stent and BMS have a similar overall risk of stent thrombosis over 4-years, and to last available follow-up beyond 4-years
- Although early, late and very late events occur in both arms, there are more events before year one for BMS and more events after year one for CYPHER® Stent
- Evaluation of varying frequencies over time is limited by small numbers of events (proportional hazards assumption was not rejected)
- Patients with BMS were more likely to have ST if they had TLR
- Although rare, patients with TLR following CYPHER® Stent did not have an increased frequency of stent thrombosis
- Clinical outcomes following ST were similar for CYPHER® Stent and BMS

Time Course of Events in SES Does Not Support a Single Simple Hypothesis



Final Conclusion

- The CYPHER[®] Stent has demonstrated impressive, sustained benefit in reducing the need for repeat revascularizations
- No difference in overall risk of stent thrombosis
 - No significant difference in death, and death or MI
 - Temporal distribution of stent thrombosis may vary between CYPHER[®] Stents and BMS
- Cordis will continue to work with the FDA to:
 - Provide physician and patient education
 - Generate the appropriate data to understand better how to reduce the risk of stent thrombosis

Recommendations and Commitments

- Endorse ACC/AHA/SCAI PCI Guidelines regarding dual antiplatelet therapy (up to 12 months for suitable patients)
- Educate on the need for dual antiplatelet therapy
 - Primary cardiologists, gastroenterologists and dentists
- Exploring patient programs to enhance compliance
 - Financial assistance, education, awareness
- Define optimal stent procedural technique

Recommendations and Commitments



- **Extend follow-up of 3 SIRIUS trials to 8-years**
- **Coordinate the extended follow-up of 10 RCTs (n=4,500 patients) to 5-years**
- **Conduct appropriately powered CYPHER® Stent PMS study with randomization to 2 durations of dual APT pending recommendations from panel**
- **Continue to examine the mechanisms of thrombosis and to optimize device design**