

1 flexion-extension radiographs, this patient had a
2 considerable range of motion, so it's not interfering
3 with motion. Next.

4 This is an HO Class of II, a little more,
5 as you can see on the far left of the A/P, but based
6 on flexion-extension this patient still had a very
7 high degree of range of motion at the operative level.
8 And again, this is 24 months post-operatively. Next.

9 And this just demonstrates the progression
10 of HO in a single patient. So if you take a patient
11 immediately post-op and follow them through two years
12 post-op, this is what the HO, how it presents itself.
13 And you can see it appearing at six months or six
14 weeks post-operatively, and then you get some
15 densification of the HO and by 24 months, it's very
16 evidence on the A/P film. Next.

17 And this is just a CT scan demonstrating
18 the location of the HO not within the disc space, but
19 it tends to be in the peri-annular region adjacent to
20 the disc and not in the psoas. This happens to be
21 four years post-op. Next.

22 And this is the Visual Analog Scale, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 data of interest here. HO and non-HO cases pre-
2 operatively, very similar based on VAS scores by 24
3 month post-operatively using the Wilcoxon Rank Sum.
4 There was no statistical differences between the two
5 groups in VAS. Next. And ODI, similar findings.
6 Pre-op, nearly identical and at 24 months, no
7 differences between the two groups. Next.

8 And this is the flexion-extension range of
9 motion. Interestingly, if you were to group these out
10 into HO and non-HO cases based on the pre-operative
11 plain films, it's a little over 6 degrees of motion
12 for both treatments. But then interestingly, at the
13 24 month post-operative period, the HO cases had more,
14 a higher range of motion at the operative level, not
15 statistical, but pretty close compared to the non-HO
16 cases. So it doesn't appear to be -- actually, the
17 range of motion is higher with the incidence of HO.
18 Next. And that is pretty much what I have for that.

19 CHAIRPERSON YASZEMSKI: Thank you, Mr.
20 Cunningham. Thank you, Dr. McAfee. Dr. Kim?

21 DR. KIM: Just going back to the question
22 of why there is a proportion of people that don't get

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 better. A proportion of them probably is coming from
2 the facet joint. If that's the case, if you looked at
3 that, and should reconsider looking at the facet joint
4 more closely prior to having somebody undergo this
5 procedure. Can one of the sponsors comment on that?

6 DR. MCAFEE: I'll give it a shot after
7 much deliberation. Paul McAfee. With any interbody
8 fusion device, you tend to unload the facets. You're
9 increasing the disc space height. The digitized
10 results of all the series both at 045 and L5-S1 did
11 show statistically that the increase in disc space
12 height was better for the SB Charite group versus the
13 BAK.

14 In addition, from the immediate six week
15 visit film to two years, it turns out that the
16 maintenance of the height was also better in the SB
17 Charite group. In other words, there was slightly
18 more subsidence with the BAKs. So we thought the main
19 thing, the main purpose of the Charite, was to unload
20 the facets and I can start with slide 247 if you would
21 like. Looking at the baboons, I know it's only a six
22 month follow-up, but the facet joints were normal at

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 sacrifice. Secondly, you know --

2 DR. KIM: Sorry to interrupt.

3 DR. MCAFEE: Yes.

4 DR. KIM: I guess what I was trying to get
5 at is if you had the opportunity to sub-stratify those
6 patients that persisted with back pain and compared
7 them to the Charite group that did not have back pain
8 and you just looked at the facet joint, would there be
9 a difference that you know of?

10 DR. MCAFEE: We didn't look at that, but
11 in a way we're selecting out patients that have
12 problems with the facet joints pre-operatively.
13 Remember in our workup, if the patient has some
14 element of mechanical back pain, we would tend to get
15 posterior facet joint blocks and if that facet joint
16 is a pain indicator, then they are selected out of the
17 study.

18 CHAIRPERSON YASZEMSKI: Thanks, Dr.
19 McAfee. Dr. Kim, does that answer your question?

20 DR. KIM: It does.

21 CHAIRPERSON YASZEMSKI: Thank you. Dr.
22 Naidu?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. NAIDU: I agree with the rest of the
2 group. Pain is a subjective measure. I think the
3 only reason that the Charite group probably had more
4 pain is probably the patient population itself was a
5 more active group. From the data presented, the
6 Charite group had a significantly lower body mass
7 index and, in general, more active. So I mean, it is
8 a subjective measure, but those are the only two
9 reasons I can think of attributing it to.

10 CHAIRPERSON YASZEMSKI: All right. Thank
11 you, Dr. Naidu. Dr. Witten, the Panel's discussion on
12 Aim 3 is that pain being a subjective measure, there
13 really are no concerns that the Charite group had this
14 percentage of people who still had continued pain, and
15 that this reflects the general treatment of low back
16 pain be it by nonsurgical or surgical methods other
17 than disc replacement.

18 Have we adequately answered this question
19 to FDA's satisfaction?

20 DR. WITTEN: Yes, thank you.

21 CHAIRPERSON YASZEMSKI: Thank you. Let's
22 move on to Question 4. Within the Charite group, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 mean range of motion and flexion-extension at the
2 treated level at three, six, 12 and 24 months was 4.9,
3 6.0, 7.0 and 7.4 degrees, respectively. Lateral
4 bending and axial rotation range of motion were not
5 reported in this investigation. Please, comment on
6 the sponsor's claim that the Charite permits "near
7 physiological segment movement with up to 15 degrees
8 bending and flexion-extension and a similar degree of
9 lateral bending and axial rotation to the natural
10 disc." Dr. Kim?

11 DR. KIM: I think the sponsor used very
12 good preclinical data to show that the disc does
13 achieve near physiologic motion at the time of
14 implantation, but the results of the clinical study
15 clearly show that the range of motion changes and is
16 variable being as low as 0 and as high as 22 degrees.

17 There is a table that the FDA put together
18 to try to make a correlation between outcome and that
19 range of motion, and there really isn't a
20 statistically significant correlation, although there
21 is a trend toward better results if you have 5 to 9
22 degrees range of motion.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Given that, it's hard to decide whether or
2 not this is significant. In addition, this range of
3 motion question is going to require long-term follow-
4 up, because one of the advantages is potentially
5 decreasing adjacent segment disease and we may not see
6 that for five to 10 years. So I don't think that this
7 study has the ability to claim that it maintains
8 physiologic motion and that that motion is the key to
9 success.

10 Going onto the second question as to
11 whether or not there is an equal amount of lateral
12 bending and axial rotation, they show that clearly in
13 the preclinical studies, but not in the clinical
14 study, because they only looked at flexion-extension.
15 The design is symmetric though and if something is
16 moving 7 degrees in the flexion-extension plane, I
17 have no problem assuming that in the lateral bending
18 plane, in the actual rotation, we'll get similar
19 degrees of motion.

20 So I'm not too worried about the comment
21 on the lateral bending and axial rotation, but I think
22 the difficulty lies in the significance of the range

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in motion with clinical outcome.

2 CHAIRPERSON YASZEMSKI: Thank you, Dr.
3 Kim. Dr. Naidu?

4 DR. NAIDU: The normal range of motion
5 cited in the literature in the PMA provided at the L5-
6 S1 is about 9 degrees of flexion, 5 degrees of
7 extension with a fairly large standard deviation of
8 plus or minus 5 degrees, and I certainly think that
9 the sponsor's claim that 4.9 degrees is physiologic in
10 flexion-extension plane is valid, although this
11 validity isn't confined to the flexion-extension plane
12 only. Obviously, rotation is questionable, at best,
13 at this point. That's it.

14 CHAIRPERSON YASZEMSKI: Thank you, Dr.
15 Naidu. Dr. Kirkpatrick?

16 DR. KIRKPATRICK: No additional comment.

17 CHAIRPERSON YASZEMSKI: Thanks, Dr.
18 Kirkpatrick. Dr. Blumenstein?

19 DR. BLUMENSTEIN: Could I ask George Chu
20 a question here?

21 CHAIRPERSON YASZEMSKI: Yes.

22 DR. BLUMENSTEIN: In the, I think you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 called it, addendum, I didn't understand the two
2 tables that were presented, the difference between the
3 two tables.

4 DR. CHU: Which two, the range of motion,
5 the histogram?

6 DR. BLUMENSTEIN: Yes.

7 DR. CHU: Or the 2 by --

8 DR. BLUMENSTEIN: Yes, the 2 by 8 tables.

9 I just want to make sure I understand. I mean, I see
10 that one of them is repeated in the question that
11 we're addressing at this point.

12 DR. CHU: The histogram is based on the
13 available data for the randomized Charite patients.
14 It's about 175.

15 DR. BLUMENSTEIN: Yes.

16 DR. CHU: So from the histogram, it looks
17 like the range of motion at 24 months post-op is
18 equally distributed about 10 percent among the
19 different range. And for the table in the Panel
20 draft, it's a 2 by 8 table, it basically just tries to
21 see the general association between the range of
22 motion and the outcome at 24 months. So the general

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 association test, the P-value is not significant.

2 DR. BLUMENSTEIN: Yes. But there was a
3 second table in your addendum.

4 DR. CHU: Which table?

5 DR. BLUMENSTEIN: I received an addendum.

6 DR. CHU: Okay. I don't have that one.

7 Can you show me that one? Yes, the second table is,
8 basically, the assessed output for the statistical
9 test.

10 DR. BLUMENSTEIN: May I ask, do we have
11 that table in the FDA? That was part of your
12 presentation, wasn't it, Dr. Chu?

13 DR. CHU: No.

14 DR. BLUMENSTEIN: Do we have that table on
15 a slide, that everybody could see it?

16 DR. CHU: Yes, actually the main point
17 here, just looking at the general association test
18 between this success outcome and the range of motion.

19 DR. BLUMENSTEIN: I mean, the titles of
20 the two tables seemed the same to me and not the
21 second smaller table underneath the big one, but the
22 second page of tables.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. CHU: The second table is just
2 collapsed the last four column of the bigger table
3 into just one column.

4 DR. BLUMENSTEIN: Oh, okay. I'm sorry,
5 sorry to cause that. Then I have no comment.

6 CHAIRPERSON YASZEMSKI: Okay. Thank you.
7 Dr. Besser?

8 DR. BESSER: While the ranges reported
9 were slightly less than have been reported in the
10 literature for physiological changes, they are well
11 within the range of normal and I guess only time will
12 tell whether, in fact, implanting a device that gives
13 this extra range of motion prevents adjacent segments
14 from needing fusing and future surgery, etcetera.

15 CHAIRPERSON YASZEMSKI: Okay.

16 DR. BESSER: No other comments.

17 CHAIRPERSON YASZEMSKI: Thanks, Dr.
18 Besser. Ms. Maher?

19 MS. MAHER: No comment.

20 CHAIRPERSON YASZEMSKI: Ms. Luckner?

21 MS. LUCKNER: No comment.

22 CHAIRPERSON YASZEMSKI: Thank you. Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Diaz?

2 DR. DIAZ: No comment.

3 CHAIRPERSON YASZEMSKI: Okay. Dr. Mabrey?

4 DR. MABREY: I have no comment.

5 CHAIRPERSON YASZEMSKI: Dr. Finnegan?

6 DR. FINNEGAN: No comment.

7 CHAIRPERSON YASZEMSKI: Thank you. Dr.
8 Witten, we have talked about the ranges of motion. We
9 generally feel that they are within the range of
10 normal, that the flexion-extension numbers are in the
11 physiologic range and less so for the other modes of
12 motion. The range of motion link to clinical
13 improvement shows a trend, but has not been met. In
14 general, the Panel doesn't have any concerns on this
15 issue.

16 And have we answered and discussed it
17 adequately?

18 DR. WITTEN: Yes.

19 CHAIRPERSON YASZEMSKI: Thank you. Let's
20 move on to number 5. Do clinical data provide
21 reasonable assurance of safety? Dr. Finnegan?

22 DR. FINNEGAN: Well, I would have to say

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that if this device was going to do its main purpose
2 over a short-term, that is a two to three year period,
3 and then would basically be physiologically non-
4 functioning, then the data does suggest that this is
5 probably safe. Unfortunately, this device is designed
6 for a much longer period of time and I do not think
7 that there is data present at the present time to say
8 that it is, in fact -- there is reasonable assurance
9 that it is safe for the lifetime that it is predicted
10 to be necessary for.

11 CHAIRPERSON YASZEMSKI: Okay. Thank you.
12 I'm going to go around to Dr. Kim, but before I leave,
13 based on that, Dr. Finnegan, I'm going to back to you
14 with Question 7 and ask what you think. Dr. Kim?

15 DR. KIM: I would agree with that. I
16 think the sponsors have done an excellent job in
17 providing an honest assessment of their device, and it
18 is absolutely clear that in the two year period that
19 this device is safe, but, once again, I agree with Dr.
20 Finnegan. This is a complex device. It's the first
21 of its kind and designed to last for a long time, and
22 we can't get at that question until we wait.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 CHAIRPERSON YASZEMSKI: Thanks, Dr. Kim.
2 Dr. Naidu?

3 DR. NAIDU: I totally concur with Dr.
4 Finnegan and Dr. Kim.

5 CHAIRPERSON YASZEMSKI: Thank you, Dr.
6 Naidu. Dr. Kirkpatrick?

7 DR. KIRKPATRICK: I concur and nothing
8 further to add.

9 CHAIRPERSON YASZEMSKI: Thank you, Dr.
10 Kirkpatrick. Dr. Blumenstein?

11 DR. BLUMENSTEIN: I concur.

12 CHAIRPERSON YASZEMSKI: Thank you. Dr.
13 Besser?

14 DR. BESSER: I concur.

15 CHAIRPERSON YASZEMSKI: Thank you. Ms.
16 Maher?

17 MS. MAHER: I would just urge the Panel to
18 remember that we also have to look at least burdensome
19 as we're figuring out how to evaluate the safety and
20 effectiveness of this device.

21 CHAIRPERSON YASZEMSKI: Thank you, Ms.
22 Maher. Ms. Luckner?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MS. LUCKNER: I concur with the earlier
2 statements.

3 CHAIRPERSON YASZEMSKI: Thank you. As we
4 come around to Dr. Diaz, Mr. Melkerson, if you have a
5 comment, may I go to Dr. Diaz and then --

6 MR. MELKERSON: Actually, just a question.

7 CHAIRPERSON YASZEMSKI: Okay. Go ahead.

8 MR. MELKERSON: This is to the Chair
9 himself. Being that Question 5 and 6 are related to
10 safety and effectiveness, do you want to hear the
11 other public speakers before you answer this question?

12 CHAIRPERSON YASZEMSKI: I think that we
13 can probably hear them afterwards and then incorporate
14 their thoughts when we get to voting if that would be
15 okay. Dr. Witten, is that acceptable to you? Thank
16 you. Thank you, Mr. Melkerson. Dr. Diaz?

17 DR. DIAZ: Nothing additional.

18 CHAIRPERSON YASZEMSKI: Thank you. Dr.
19 Mabrey?

20 DR. MABREY: Well, as has been pointed out
21 by other Panel Members, I mean, this is a complex
22 device. It's brand new and it's going to be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 eventually implanted by a lot more than the original
2 surgeons. So I think we're looking at two levels of
3 safety. One is is it safe to implant, and I think
4 over the first two years you have demonstrated that
5 with your trained surgeons that it is safe to implant.
6 And then the second question is is the device itself
7 safe over a long period of time, and I don't think two
8 years is long enough.

9 So I do have one question, and I would
10 address this to Dr. Blumenthal or Dr. McAfee. At your
11 training centers, the only analogous situation I can
12 come up with is the experience with another company's
13 foray into minimally invasive hip surgery and
14 restricting access to that to those who have been
15 trained at the company's facility.

16 One comment that I have heard from the
17 trainers there is that there is a training of the
18 trainers that goes on, and I guess my question is do
19 the clinicians feel that with more experience, that
20 your initial training of those people who will be
21 implanting the device, does that become easier and
22 have you learned to avoid some of the major problems

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 you encountered during the initial phase of this
2 study?

3 DR. MCAFEE: I'll start, because one of
4 the main concerns is neurologic problems, so there has
5 been a great advance in instrumentation and, honestly,
6 the key with any anterior interbody device is to keep
7 it in the midline, so that newer instruments over the
8 last two and half years happen to be called the
9 centerline instruments, but they keep the implant in
10 the midline and they prevent the surgeon from going
11 into the lateral recess and causing a neurologic
12 problem.

13 Secondly, if you can put up -- start with
14 maybe slide 493. We did an analysis. We wanted to
15 see what the effect of the training was. The
16 training, you know, we had a perfect opportunity to do
17 that, because we had a cohort of five cases from each
18 group that were training cases, and we could compare
19 how well the training patients did with the rest.

20 Well, the idea is whether there is a
21 surgeon volume effect. Is the data good enough to
22 show a surgeon volume effect, and this is from

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Birkmeyer's New England Journal of Medicine and
2 "JAMA," a lead article in December. And the key was
3 there are 10 different operations that all show a
4 surgeon volume effect, coronary artery bypass, grafts,
5 aortic valve replacements, surgeons that did more and
6 had a higher volume who had lower complications, but
7 not a single spine procedure was in this group. Next
8 slide, please.

9 So we looked at actually four areas.
10 First, we looked at the 71 training cases versus the
11 randomized cases. Then the next analysis was we
12 looked at the four highest enrolling sites that all
13 did more than 40 procedures versus the 11 remaining
14 sites that didn't do as many cases. Now, the key is
15 that all four groups fulfilled the FDA's success
16 criterion of greater than 25 percent improvement in
17 the Oswestry, as well as no neurologic progression, no
18 return to the OR and no major complications. Next
19 slide.

20 But there is a definite volume effect, and
21 if you look at the training cases, which are on the
22 left, the surgery time was larger in the person's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 first five cases versus the enrolled. The length of
2 stay was slightly longer for those patients, and the
3 overall number of complications was higher as well.
4 Next slide.

5 You know, we looked at all the parameters,
6 but I'm just showing the ones that are significant.
7 And then the high enrolling sites versus the low
8 enrolling sites, the surgery time was much less for
9 those sites that did more than 40 cases. The length
10 of stay was less and the device failure incidence was
11 lower. So in summary, surgeon volume really did have
12 an effect. And then the last slide.

13 I think the key is to learn from the
14 European experience, really almost memorize the IDE
15 prospective randomized trial data. The key to
16 avoiding complications are to identify a vascular
17 access surgeon, go to company-sponsored courses to
18 learn the specific instrumentation.

19 But what is even more important is what we
20 haven't talked about and that is the model of the
21 Scoliosis Research Society, which is something like
22 the Spine Arthroplasty Society, what we'll hear from

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the public comment period, is they are going to the
2 forefront of continuous reassessment of results and
3 having a surgeon group take some responsibility for
4 setting the bench mark.

5 CHAIRPERSON YASZEMSKI: Thank you, Dr.
6 McAfee. Dr. Witten, we have discussed the issue of
7 safety and the consensus of the Panel is that over the
8 study period in the short-term, this device is safe,
9 and questions remain, of course, over the long period,
10 because the data is not gathered yet and it's a new
11 device without precedent.

12 Have we discussed this to the satisfaction
13 of FDA?

14 DR. WITTEN: Yes, thanks.

15 CHAIRPERSON YASZEMSKI: Thanks, Dr.
16 Witten. We're going to move on to number 6. Do
17 clinical data provide reasonable assurance of
18 effectiveness? Dr. Diaz, can you lead off with this
19 one?

20 DR. DIAZ: Yes. The question that we're
21 being asked and I will read a little bit of the
22 definition the FDA wants us to adhere to for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 effectiveness. "A device is considered effective if,
2 when used in a significant portion of a targeted
3 population for the intended use and under the
4 conditions of use, it provides significant results for
5 that population."

6 To answer that question, I would look at
7 the effect, clinical effect, on various aspects of the
8 individual's life. When I operate on a patient, the
9 question the patient asks me very frequently is will
10 I get better? Will my pain get better? Will I be
11 able to get back to work? And probably more
12 important, will I be able to get back to play? Work
13 in a back pain patient group is not always what they
14 want to do, but play certainly is what they want to
15 get back to do.

16 And as we have heard, the back problem,
17 spine progressive degeneration is a dynamic problem.
18 Coming from the Rust Belt in Detroit, when I talk to
19 patients about spine surgery, I tell them that it is
20 like dealing with rust. All patients in Detroit can
21 understand rust. If you get it on the fender and you
22 clean it, you patch it, you fix it, it will show up on

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the door and a few years later, it will show up on the
2 fender. So they can relate to the idea that perhaps
3 this is not the only time that we'll see Dr. Diaz to
4 take care of their problem.

5 So in regard to the effectiveness, does
6 the device provide pain relief? Yes, it does.
7 Perhaps not as well as it did for the BAK, but
8 significantly equal. Does it restore function? We
9 believe it does, based on the anatomical and on the
10 mechanical presentations given. Does it allow
11 patients to get back to their usual activities? The
12 answer is yes. And to my personal liking, I was very
13 pleased to see that patients could get back to very
14 active function very quickly. Because within the week
15 after surgery, they could be doing a lot of things
16 that I keep my patients from doing when I fuse them.

17 If I fuse a patient, I really keep them
18 sedentary for a long time. I don't like that. I like
19 to be able to get people up and moving very quickly
20 and I think this device provides for that opportunity.
21 Do they return to work? Yes, perhaps they do. Maybe
22 not as much as I would like them to see. But I don't

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 think we will ever be able to get back patients to get
2 back to work as much as we would like them to do that.

3 And does it prevent adjacent level
4 disease? I don't think that this single device will
5 be the answer for preventing adjacent level disease,
6 but I think it delays it, which I think is a very
7 important achievement. So in my mind, I believe for
8 the intended use that the device was proposed in the
9 population as targeted with the possible applications
10 as provided, it does fulfill the requirements of
11 effectiveness under the FDA guidelines.

12 CHAIRPERSON YASZEMSKI: Thank you, Dr.
13 Diaz. Dr. Mabrey?

14 DR. MABREY: I fully concur with Dr. Diaz'
15 comments.

16 CHAIRPERSON YASZEMSKI: Thank you. Dr.
17 Finnegan?

18 DR. FINNEGAN: I agree.

19 CHAIRPERSON YASZEMSKI: Thank you. Dr.
20 Kim?

21 DR. KIM: I agree as well.

22 CHAIRPERSON YASZEMSKI: Thank you. Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Naidu?

2 DR. NAIDU: Same here.

3 CHAIRPERSON YASZEMSKI: Thank you. Dr.
4 Kirkpatrick?

5 DR. KIRKPATRICK: Well, we're not going to
6 get a dissertation today. I think there are some
7 concerns about effectiveness, but I think by the FDA's
8 definition, I would agree with Dr. Diaz.

9 CHAIRPERSON YASZEMSKI: Thank you. Dr.
10 Blumenstein?

11 DR. BLUMENSTEIN: I believe the device has
12 been shown to be not inferior to the standard control
13 device that was used in the trial.

14 CHAIRPERSON YASZEMSKI: Thank you. Dr.
15 Besser?

16 DR. BESSER: I agree with almost
17 everything Dr. Diaz said, other than I don't think we
18 have any evidence to support that it will, in fact,
19 delay adjacent segment disease. We hope, we'll see,
20 nothing now.

21 CHAIRPERSON YASZEMSKI: Thank you. Ms.
22 Maher?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MS. MAHER: Well, I would agree with Dr.
2 Diaz. I would also remind the Panel that this product
3 has been on the market since 1987 in Europe and there
4 have been 7,000 cases or implantations. So when you
5 talk about only having a two year follow-up, as we did
6 on the previous question, I want you to remember that
7 there actually is a much longer history outside the
8 U.S. There is two years of good data from within the
9 U.S.

10 CHAIRPERSON YASZEMSKI: Thank you. Ms.
11 Luckner?

12 MS. LUCKNER: I concur with Dr. Diaz.

13 CHAIRPERSON YASZEMSKI: Thanks very much.
14 Dr. Witten, the Panel feels that the device as
15 presented is effective. Have we adequately discussed
16 this?

17 DR. WITTEN: Yes.

18 CHAIRPERSON YASZEMSKI: Thank you. We're
19 going to move on to number 7. Number 7, if you
20 recommend approvability for this PMA, do you recommend
21 a post-approval study? If so, please, discuss what
22 types of end points would be useful for an updated

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 label and recommend the duration of such a study. Dr.
2 Kirkpatrick, would you lead this one, please?

3 DR. KIRKPATRICK: I think he asked me
4 because I already provided a list of suggestions. I
5 would like to see mobility testing data for the
6 complete reference, rather than just a two paragraph
7 summary. And I did get a little bit more of it in the
8 presentation today, but I would like to be able to
9 review the data, and I think the FDA would like to be
10 able to review that as well.

11 I think a little added study in the
12 biomechanics lab of demonstrating that facet stresses
13 or strains or some other element of a facet function
14 is either unchanged or minimally changed after the
15 insertion of the disc. I did not get that out of my
16 read of the PMA and again, as I have repeatedly said,
17 if you have that data and I missed it, please, tell me
18 what page to find it on.

19 The third one is I still think that the
20 wear data to 50 million cycles would be more
21 appropriate. I did have concerns whether the curve on
22 the wear data may have been accelerating over time.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Your curve on, I think it was, the weight of the
2 specimen versus cycles was drawn as a curve going
3 downward, which means it is accelerating as you get to
4 the end of that 10 million cycles. I would like to
5 know what it does in the next four decades of the 50
6 million cycle testing.

7 I think we need an acceptable rationale
8 for not testing the response to submicron particles
9 more extensively. Number 5 that is on the list, we
10 can exclude, because I clarified that in my
11 presentation and you clarified it in your presentation
12 that the osteointegration studies are not relevant to
13 the device you are presenting, so you can eliminate
14 number five. However, I would like to know what the
15 rationale is behind the long-term fixation of the
16 device. Is it just the pegs or do you expect there to
17 be some bone implant interface adherence?

18 6, I would like a clerical -- clarifying
19 of the neurologic rating scale that you used, so I can
20 understand how these statistics were applied to a
21 qualitative physical exam. 7, I think that was
22 handled by the FDA as far as stratifying the range of

1 motion, and it appears that outcome is not
2 significantly different, although there were trends.
3 So 7, I'll leave it up to the Panel whether we think
4 we need to go further with that.

5 8, indication groups, especially the ones
6 that did have known facet changes at the implantation.
7 I would like to know if stratifying those out would
8 make a difference in either your BAK group or your
9 charity group. I would suggest again based upon the
10 literature as well as my presentation that the concept
11 is if we're preserving motion, we need to demonstrate
12 that. And if people lose motion, I would like to know
13 if that resulted in a difference in their other
14 measures of effectiveness, such as the VAS, the ODI
15 and that sort of thing.

16 So number 9 would say include those 0 to
17 5 degrees as failures and see if that correlates to
18 clinical failure. And also, if you called them
19 failures, what would happen to your statistics on the
20 study success. 11 you can eliminate based upon the
21 discussion of HO and the presentation you did in
22 answer to one of our questions. I'm sorry, that's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 number 10. Number 11 gets back to the facet issue.
2 Can you do axial imaging at 24 months and look at the
3 comparison between your pre-op and your axial imaging
4 at 24 months and tell us whether there are facet
5 changes.

6 Number 12, adjacent segment degeneration,
7 I think, should be looked at. We obviously have the
8 x-rays stored on computer data and that should be
9 something that could be doable. And then 13 is
10 perhaps the most difficult and that is, I think, the
11 follow-up should be extended to five years to get to
12 some semblance of a number of where we would see
13 adjacent segment, so that we can back up the rationale
14 that we are preserving the adjacent segment from the
15 standpoint of the philosophy of the disc replacement.
16 Thank you.

17 CHAIRPERSON YASZEMSKI: Thanks, Dr.
18 Kirkpatrick. May I ask before we move to Dr.
19 Blumenstein, may I ask that among the suggestions you
20 made there are some that seem to be answerable by
21 relooking at the existing preclinical and clinical
22 data and perhaps one or more that may require further

1 study data after approval. Would you care to comment
2 on? It would seem to me that number 13 might, of
3 course, require more study data clinically, and that
4 maybe number 3 would require in-vitro data, but that
5 the others perhaps could be answered by looking at the
6 existing data. Would that be accurate?

7 DR. KIRKPATRICK: I agree that number 3
8 and 13 definitely would require additional work. I
9 think the remaining things, as I recall, are either a
10 discussion of their existing data or an expansion on
11 analysis of their existing data.

12 CHAIRPERSON YASZEMSKI: Okay. Thank you.
13 And then when we have the sponsor summary, I'll ask
14 the sponsors to comment on these. Dr. Blumenstein?

15 DR. WITTEN: Can I just mention one thing?

16 CHAIRPERSON YASZEMSKI: Dr. Witten?

17 DR. WITTEN: Yes, I'll just mention one
18 thing which is that when you get to the vote, you'll
19 have to clarify for us for each of these
20 recommendations whether or not these are things that
21 you would expect to see pre-approval or post-approval,
22 because if it is new data, for example, then that's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 not a condition of approval recommendation. It's a
2 recommendation to put the PMA in approvable form. So
3 I just want a clarification of when you think we need
4 this information.

5 CHAIRPERSON YASZEMSKI: Thank you. Dr.
6 Blumenstein?

7 DR. BLUMENSTEIN: My main concern is the
8 long-term follow-up and I think that has been
9 addressed.

10 CHAIRPERSON YASZEMSKI: Thank you. Dr.
11 Besser?

12 DR. BESSER: I concur.

13 CHAIRPERSON YASZEMSKI: Thank you. Ms.
14 Maher?

15 MS. MAHER: I'm going to sound a little
16 bit like a broken record like I always do and again
17 remind everybody that we do have data back to 1987.
18 We do have a significant patient population outside
19 the United States, and so maybe a post-approval study
20 following the other patients in the study now for
21 longer would be appropriate, but some of the rest of
22 the data may not be necessary.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 CHAIRPERSON YASZEMSKI: Thank you, Ms.
2 Maher. Ms. Luckner?

3 MS. LUCKNER: Nothing else to add.

4 CHAIRPERSON YASZEMSKI: Thank you. Dr.
5 Diaz?

6 DR. DIAZ: I would like to agree with Ms.
7 Maher, because I think there is significant clinical
8 data available in the world literature that indicates
9 the longevity and the effectiveness of this device in
10 the treatment of discogenic disease. The available
11 literature does not answer all the questions that Dr.
12 Kirkpatrick mentioned, but those could be answered on
13 an ongoing type analysis, rather than trying to
14 redesign a new study. I believe that there is
15 sufficient information already available to answer
16 many of these things. And going back to square one,
17 I don't think is necessary.

18 CHAIRPERSON YASZEMSKI: Thank you, Dr.
19 Diaz. Dr. Mabrey?

20 DR. MABREY: Yes, I concur with the plan
21 to go ahead with looking at those individuals that are
22 currently under study and also perhaps extend some of

1 these investigations to those patients who are
2 available in Europe. And particularly, I'm
3 interesting in looking at the possibility of
4 osteolysis at four and five years out. I think you
5 have the potential to continue to look at radiographic
6 data on the original 200 patients here and it's no
7 additional great burden and it would be nice to see
8 data from European studies indicating that there is no
9 osteolysis.

10 CHAIRPERSON YASZEMSKI: Thank you.

11 MS. MAHER: Can I ask for a clarification
12 for that? Were you talking about a post-approval type
13 of look at it for that data?

14 DR. MABREY: Okay. This is my first Panel
15 meeting and so we've been talking about post-approval.

16 CHAIRPERSON YASZEMSKI: And may I
17 interrupt?

18 DR. MABREY: And PMAs and PDPs.

19 CHAIRPERSON YASZEMSKI: May I interrupt,
20 Dr. Mabrey? What I'll suggest is all of these things
21 that we recommend to FDA, we will need to recommend
22 whether they are things that need to be done before

1 the approval, and thus defer the approval, vote that
2 this is non-approvable or whether we would agree that
3 this is an approvable application, and in addition as
4 terms of the approval, we would like them to do
5 further work and follow the patients.

6 DR. MABREY: Okay.

7 CHAIRPERSON YASZEMSKI: And we will just
8 have to make that distinction when we come to voting.

9 DR. MABREY: Okay.

10 CHAIRPERSON YASZEMSKI: And thank you, Ms.
11 Maher, for bringing that clarification up.

12 DR. MABREY: Okay.

13 CHAIRPERSON YASZEMSKI: Dr. Finnegan?

14 DR. FINNEGAN: Now, you've changed my
15 train of thought. I was getting all my thoughts
16 together. One of the problems with looking over the
17 data here is that unfortunately of the 205 or so
18 patients that had the Charite implants, a number of
19 them have actually not reached two years yet. And I
20 think that we have already all pretty much agreed that
21 two years is probably not a safe length of time to
22 follow these patients. So I do think that there needs

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to be a long-term follow-up. I definitely agree with
2 Dr. Kirkpatrick on that.

3 And I know that we are looking at the end
4 plate with the spikes. If there is a change to a
5 coated end plate, then that's obviously going to
6 drastically change the biomechanics on the
7 polyethylene and understanding that the company feels
8 the polyethylene is cross-linked to some degree, it is
9 obviously a random cross-linking. So I think all of
10 those are things that need to be considered.

11 And as well, I do think that the company
12 needs to be at least -- or the sponsor needs to be at
13 least familiar with neurological response to chronic
14 inflammation and be comfortable that that is not going
15 to be a long-term problem. And I certainly agree with
16 Dr. Kirkpatrick on adjacent segments. I think that
17 some kind of study needs to be done on that.

18 CHAIRPERSON YASZEMSKI: Thank you, Dr.
19 Finnegan. Dr. Kim?

20 DR. KIM: I think if we only looked at the
21 U.S. clinical trial data, I would be nervous with just
22 two year results. But I agree with Ms. Maher that we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 have the -- we are fortunate that we have a pretty
2 extensive European experience of 7,000 patients. I
3 think with the excellent U.S. clinical data of two
4 years combined with the European data, it's a very
5 promising device. And based on that, I feel like this
6 device is safe and effective. But there are a lot of
7 questions that remain and I think looking at the
8 existing patients over the long-term, maybe over five
9 or even 10 years, is a reasonable condition.

10 CHAIRPERSON YASZEMSKI: Thank you, Dr.
11 Kim. Dr. Naidu?

12 DR. NAIDU: You know, I would like to
13 listen to the second open public hearings prior to
14 commenting on this.

15 CHAIRPERSON YASZEMSKI: Thank you, Dr.
16 Naidu. Dr. Blumenstein?

17 DR. BLUMENSTEIN: I just get so enamored
18 of the case series data that are likely to come out of
19 Europe, other places like that. There's nothing more
20 valuable than the data that have been invested into
21 this clinical trial in a comparative way and the
22 potential for that data to come out with an unbiased

1 comparison of the results as opposed to the
2 uncontrolled convenient samples that are often
3 published in the literature.

4 CHAIRPERSON YASZEMSKI: Thanks, Dr.
5 Blumenstein. Dr. Witten, as you've heard, we've had
6 a more extensive discussion on this question. And I
7 would like to review it and then to say that when we
8 get around as a Panel to making a recommendation to
9 vote on, we'll consider which of these suggestions
10 might be conditions of approval and which we would
11 want to be done after the approval vote.

12 We've talked about adding mobility testing
13 data and in-vitro study of the facets, wear data to 50
14 million cycles, test response to submicron particles,
15 to consider using data from the existing European
16 studies and we have heard pros and cons about that
17 from several Members of the Panel, osteolysis at four
18 to five years. Many of the Panel Members thought that
19 going out about five years for several of these end
20 points would be appropriate. Dr. Finnegan mentioned
21 the effects of inflammation on the neurologic tissues
22 by chronic inflammation and Dr. Kirkpatrick about

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 adjacent segment.

2 Have we discussed this to your
3 satisfaction?

4 DR. WITTEN: Yes, thank you.

5 CHAIRPERSON YASZEMSKI: Thanks, Dr.
6 Witten. That's going to conclude the discussion on
7 the specific questions that the FDA has posed of us.
8 We are going to move now to the second open public
9 hearing. There are three people who wish to present
10 to the Panel, at this point. These are Dr.
11 Hochschuler, Dr. Van Ooij and Ms. Adams. Dr.
12 Hochschuler will be first with a time of five minutes.

13 DR. HOCHSCHULER: I am Steve Hochschuler.
14 I am a spine surgeon. I am Chairman of the Texas Back
15 Institute, and today I am representing the Spine
16 Arthroplasty Society in my presentation. First, I
17 want to thank you for allowing me to come to the
18 podium today. Secondly, despite having been a spine
19 surgeon for about 28 years now, this is the first time
20 I have been at an FDA Panel meeting. And I have to
21 say as a citizen, I'm very impressed.

22 I have my own bias as to whether this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 should be approved or not approved, but I must say
2 it's almost like a TV show. I'm not sure how you are
3 going to vote an I'm really intrigued. I think it's
4 a wonderful process and I don't want to be
5 supercilious, but I would like to compliment you and
6 thank you.

7 Having said that, I feel as a surgeon that
8 it is our primary responsibility to care for the
9 patient. And with that, I have been charged along
10 with the rest of the Spine Arthroplasty Society to put
11 together a position statement on some of the items you
12 discussed earlier in terms of safety, how do we
13 protect our patients, how do we get better outcomes.
14 With this in mind, despite the fact I usually don't
15 like to read directly, I would like to read this
16 statement to you, since I've got limited time, and
17 then go from there.

18 "Spinal Arthroplasty Society Educational
19 Objectives." "The board of directors of the Spinal
20 Arthroplasty Society has decided to take a unique step
21 in establishing education and training goals for spine
22 surgeons interested in new arthroplasty technologies.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 The ultimate goals of this effort are to improve
2 clinical outcomes and reduce technical complications
3 in patients undergoing surgical treatment utilizing
4 these new technologies by providing a strong
5 educational core of knowledge for surgeons.

6 Traditionally, rigorous patient selection
7 criteria have been required for inclusion in FDA
8 trials. Additionally, investigators are specifically
9 selected by the companies who design the studies based
10 on their reputations and experience. However, when
11 devices are approved for marketing to surgeons in the
12 community, there has been no formal standardization
13 for training these physicians in your use.

14 Training historically has run the gamut
15 from a product introduction by a company
16 representative sitting across the surgeon's desk to a
17 brief course with a lecture in the morning followed by
18 a crowded hands-on training using saw bones models to
19 a comprehensive training program incorporating surgeon
20 education for diagnostic workup, patient selection
21 criteria, management of complications and ample time
22 in a cadaver lab developing familiarity with the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 instrumentation and surgical exposure.

2 Ideally, comprehensive formal training
3 should be followed by proctorship at the training
4 surgeon's hospital for its first case or cases by a
5 teaching surgeon with a high level of expertise. This
6 would serve to close the loop of the surgical
7 proctoring process. Obviously, this level of training
8 is expensive and time consuming, but it offers
9 significant long-term advantages for patients,
10 surgeons, industry and hospitals.

11 For patients, technical complications may
12 be reduced and outcomes improved. For surgeons, their
13 patients' clinical results may be more gratifying and
14 litigation avoided. It is important for industry so
15 that their devices can produce the best results
16 possible. A product may be unjustly criticized for
17 high complications and poor outcomes if surgeons have
18 poor technical skills or employ too broad patient
19 selection criteria.

20 Hospitals also have a vested interest in
21 the training of surgeons. The hospital's mission like
22 that of the surgeon is to ensure the maximum benefit

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to the patient. While the technology of spinal
2 surgery is steadily advancing, clinical safety and
3 outcomes cannot be expected to improve unless the
4 appropriate patient selection and optimal surgical
5 techniques are taught. SAS is prepared to take a pro-
6 active role in addressing surgeon education. A
7 program of organized processes for training surgeons
8 on new devices will incorporate didactic lectures,
9 hands-on training and proctorships.

10 The role of the Society will be to develop
11 guidelines for content of educational programs,
12 identify training centers with adequate facilities and
13 staffing for consistent quality training and organize
14 access to specialists with experience with the
15 specific devices to provide proctorships. Due to
16 liability issues, certification can verify that the
17 surgeon has completed training, but not that he or she
18 has adequate skills. A document will be issued only
19 to verify course attendance and subsequent
20 proctorship.

21 The fact that training is provided through
22 a Society and performed in an organized, standardized

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 format across the country and hopefully the world
2 should enhance the overall quality of care for our
3 patients. All parties concerned recognize the
4 importance of having surgeons properly trained when
5 introducing new technologies. With the rapid
6 developments in spinal implants, SAS has an
7 unprecedented and unique opportunity to play an
8 important role in improving patient care, optimizing
9 the application of new technologies and furthering the
10 development of new implants by increasing the safety
11 of new product introduction and adoption of these
12 standardized training programs. Thank you.

13 CHAIRPERSON YASZEMSKI: Thank you very
14 much, Dr. Hochschuler. Dr. Van Ooij? Dr. Van Ooij is
15 scheduled for 10 minutes.

16 DR. VAN OOIJ: Thank you, Mr. Chairman.
17 I am very honored to be here and to speak to you. I
18 am an orthopedic surgeon, spine surgeon in Maastricht,
19 the Netherlands for 24 years and I'm a member of the
20 SRS, the Scoliosis Research Society and the European
21 Spine Society. I will talk about the other side of
22 the Charite disc prosthesis, so I talk about

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 complications that I see in a cohort of about 500
2 patients that were operated in a neighboring hospital,
3 and another 500 has been operated by countrymen of me
4 in Munich in Germany. So this is a cohort of about 50
5 patients. They are already a little bit more, but 49
6 were evaluated.

7 So if we could have the next slide? Where
8 do I have to press? All right. Oh, yes. So in eight
9 years, I saw 49 patients, 28 women, 28 men, and with
10 a young age, of course, because that is in the
11 indication and there were operations performed as
12 early as in 1989. 20 in the period of the first five
13 years. Then 24 in the second five years. And seven
14 in the last four or five years. The next one. So
15 most patients were operated in one level, of course,
16 let me see, some in two levels. Two levels in 10
17 patients and three levels in two patients. There were
18 a lot of previous operations done, but most had no
19 operations before this.

20 Next one, please. So there were early
21 complications, subluxation of a prosthesis and removal
22 after a few days, some hematoma. In men, there is a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 risk of retrograde ejaculation and erectile
2 dysfunction. If you ask the male people, they
3 sometimes have a dysfunction without retrograde
4 ejaculation. There is one patient that had a urethra
5 lesion with a large urinoma.

6 Next slide. This is the patient with a
7 stint in place and here a large urinoma from a urethra
8 lesion. This was punctured several times and there
9 was pseudomonas involved and probably she has a low
10 grade infect now in one of these prostheses, so this
11 is really a problem. She is in a bed. All these
12 patients have really terrible leg and back pain. They
13 have VAS scores about 8 to 9 mostly.

14 Next patient -- sorry, next slide, please.
15 The leg complications are migration. These mainly are
16 prosthesis uncoated, so there were anterior migration,
17 posterior migration, even the main cause of complaints
18 after a year or a lot of them remainder complaints is
19 disc degeneration at the other levels. In 13 patients
20 it was -- this was not obvious before the operation on
21 plain x-rays and discography, but there were the other
22 ones had more or less degenerated disc, but without

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 pain on discography. Facet joint degeneration is a
2 big problem, I think. In the late situation, we saw
3 it 16 times.

4 Next, please. This is the patient with
5 anterior migration. You can see that in 10 years this
6 is 1991 and this is 2001, this was sliding anteriorly
7 and pressed on the big vessels and we had to remove it
8 and luckily we were successful in doing it without
9 lesions of the vessels, but this has been reported and
10 undoubtedly many times if you hear the conferences and
11 there was a fusion done, and this is the only one of
12 the re-operations that the patient is satisfied.

13 Next, please. There is a big issue, I
14 think, in facet joints arthrosis because it's probably
15 the biomechanical behavior of the prosthesis. I'm
16 very worried about axial rotation that is increased in
17 the prosthesis compared to the normal disc, so you get
18 a big load, I think, on the facet joints. Probably
19 also when it is more anterior located, you must put
20 the prosthesis really posterior to get some kind of
21 motion and these are the facet joints that are really
22 very hypertrophic and are triadic.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Next, please. Subsidence is a big issue
2 in this series. In 17 of these patients the
3 prosthesis was obviously too small. There was some
4 subluxation of the core. One big issue that was not
5 spoken about today is breakage of the metal wire. If
6 you look good at the x-rays, you can see the breakage
7 and the flattening of the polyethylene core and
8 probably also some wear debris. Hyperlordosis should
9 be an issue if you distract the segments, you get
10 easily in hyperlordosis and asymmetrical loading of
11 your facet joints. And I think that the patient that
12 Mr. McAfee demonstrated had, in my mind, really aware
13 of the prosthesis.

14 Next case. So this is the patient, a
15 patient with a subsidence that is seen many times and
16 it can go all the way posteriorly or anteriorly or
17 sideways and this was fused, but the patient keeps on
18 complaining. Probably, I think, that the posterior
19 stabilization and fusion is not the answer, because
20 most people keep complaining because of micro motion
21 in the prosthesis, despite the posterior fusion.

22 Next case. This is a patient that really

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 bothers me with degenerative scoliosis developing
2 after seven years. Some patients have multiple
3 degeneration above the discs in the Charite and they
4 get degenerative scoliosis above or including the
5 prosthesis. I think mainly from axial rotation
6 problems or because of the forces that go through the
7 spine that are blocked by the prosthesis. You can say
8 that these are stones in the shoe. If you have a
9 stone in the shoe, you get pain in all your leg and I
10 doubt that this will mimic the natural motion so
11 intimately that you prevent really motion
12 degeneration.

13 Next. This is a case with a broken ring
14 and a flattening of the core at the posterior side,
15 but this is hard to see on this slide, but you should
16 observe that and look very carefully at it. Next
17 slide. And this is the patient with the wear, which
18 has already been shown by Mr. McAfee, of one of my
19 patients with holes in the bone and scoliosis and the
20 flattening of the core. That is indicative of wear,
21 I think.

22 Next, please. We did 21 or 21 additional

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 operations were done. We did 11 of them. 10 were
2 already performed before. Most had posterior fusions,
3 but most patients without much benefit, so I really
4 would stress that it is not a good solution for the
5 problem. Probably this prosthesis will be a pain
6 source afterwards.

7 Next, please. Many patients in this
8 series, I think, had back surgery, back placement and
9 back sizing, but also in the boot placement and boot
10 sizing that were problems and it seems from this very
11 experienced surgeon that it's not really -- that it is
12 really difficult to do good surgery like this disc
13 prosthesis placement. And I think that it is not
14 behaving as a normal disc. The center of rotation has
15 been talked about. If you really put it posteriorly,
16 it could be well, but then you have the risk of going
17 over the edge and getting a rim fracture. Two
18 patients in my series have a posterior placement.
19 Nobody talks about shock absorption, but Lueck has
20 shown that there is no shock absorption and that the
21 normal forces that go down the disc are not going like
22 normally when you have a disc prosthesis in.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Next, please. And the rotation, I already
2 talked about. So a lot of problems will be seen, I
3 think. Also, in the United States if you wait long
4 enough, two years, in my mind, is far too short to see
5 those problems. Wear will be a big issue in the
6 future. I'm convinced of that, because the forces on
7 the lower spine are very, very high. Revision is
8 dangerous and sometimes impossible and I go to series
9 from surgeons and they say that they couldn't reach
10 the prosthesis, because the vessels were too adherent
11 and the claim of preventing adjacent disc degeneration
12 is not substantiated.

13 Next, please. So that was the end. I
14 want to report an investigation that was presented
15 yesterday in Porto, where I was yesterday, and it was
16 from the Charite group. They sell from East Berlin of
17 Berlin now where it was originated. They did a 17
18 year follow-up of 53 patients out of a group of 71
19 patients. And 60 percent of the segments were fused,
20 really fused, didn't move anything and didn't move at
21 all and most had already bone in it and were really
22 fused. But they were the better one and the patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that still fused still had some motion in it were the
2 bed one, and I think you should look also to that
3 series. That is the most, the longest experience to
4 date, 17 years, and the conclusion of Mr. Brooks here
5 was that there was no indication for a disc prosthesis
6 in the disc disease. Thank you to showing it to you.

7 CHAIRPERSON YASZEMSKI: Thanks, Dr. Van
8 Ooij. May I ask before you leave, we have asked all
9 the speakers to state for the transcriptionist for our
10 record, the Conflict of Interest statement, the three
11 questions and, please, your industry relations, any
12 financial aspects that you might have and the source
13 of funding for your trip here.

14 DR. VAN OOIJ: Yes, thank you. I forgot
15 to name that. I have no personal financial
16 relationships with any industry. Medtronic Company
17 brought me here, provided for the travel. And
18 further, I have no relationship whatsoever.

19 CHAIRPERSON YASZEMSKI: All right. Thank
20 you very much, sir. The third speaker will be Ms. Pam
21 Adams from OSMA. And, Ms. Adams, you are scheduled
22 for five minutes.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MS. ADAMS: Good afternoon. My name is
2 Pamela Adams and I speak here today representing the
3 Orthopedic Surgical Manufacturers Association or OSMA.
4 OSMA, a trade association, with over 30 member
5 companies welcomes this opportunity to provide general
6 comments at today's Panel meeting. OSMA's comments
7 should not be taken as an endorsement of the product
8 being discussed today. We ask instead that our
9 comments be considered during today's Panel
10 deliberations. These comments represent the careful
11 compilation of our member companies' views.

12 I would like, first, to provide a brief
13 introduction and background. OSMA was formed over 45
14 years ago and has worked cooperatively with FDA and
15 the American Academy of Orthopedic Surgeons, the
16 American Society of Testing and Materials and other
17 professional medical societies and standards
18 development bodies. This collaboration has helped to
19 ensure that orthopedic medical products are safe, of
20 uniform high quality and supplied in quantities
21 sufficient to meet national needs.

22 OSMA membership currently includes

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 companies who produce over 85 percent of all
2 orthopedic implants intended for clinical use in the
3 United States. OSMA has a strong and vested interest
4 in ensuring the ongoing availability of safe and
5 effective medical devices. The deliberations of the
6 Panel today and the Panel's recommendations to the FDA
7 will have a direct bearing on the availability of new
8 products.

9 We make these comments to remind the Panel
10 of the regulatory burden that must be met today. We
11 urge the Panel to focus its deliberations on the
12 product safety and effectiveness based on the data
13 provided. As regards reasonable assurance of safety
14 and effectiveness, the FDA is responsible for
15 protecting the American public from drugs, devices,
16 food and cosmetics that are either adulterated or are
17 unsafe or ineffective. However, FDA has another role
18 to foster innovation.

19 The Orthopedic Devices Branch is fortunate
20 to have available a staff of qualified reviewers,
21 including a Board certified orthopedic surgeon to
22 evaluate the types of applications brought before this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Panel. The role of this Panel is also very important
2 to the analysis of the data in the manufacturer's
3 application and to determine the availability of new
4 and innovative products in the U.S. marketplace.

5 Those of you on the Panel have been
6 selected based on your expertise and training. You
7 also bring the view of practicing clinicians who treat
8 patients with commercially available products. OSMA
9 is aware that you have received training from FDA on
10 the law and the regulation and I do not intend to
11 repeat that information today. We do, however, want
12 to emphasize two points that may have a bearing on
13 today's deliberations.

14 Firstly, reasonable assurance of safety
15 and effectiveness and secondly valid scientific
16 evidence. As regards the first point, there is a
17 reasonable assurance that a device is safe when it can
18 be determined that the probable benefits outweigh the
19 probable risks. Some important caveats associated
20 with this over simplified statement include valid
21 scientific evidence and proper labeling and that
22 safety data may be generated in the lab, in animals or

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in humans.

2 There is a reasonable assurance that a
3 device is effective when it provides a clinically
4 significant result. Again, labeling and valid
5 scientific evidence play important roles in this
6 determination. The regulation and the law clearly
7 state that the standard to be met is a reasonable
8 assurance of safety and effectiveness. Reasonable is
9 defined as moderate, fair and inexpensive.

10 As regards the second point, valid
11 scientific evidence, the regulation states that well-
12 controlled investigations shall be the principal means
13 to generate the data used in the effectiveness
14 determination. The following principles are cited in
15 the regulation as being recognized by the scientific
16 community as essentials in a well-controlled
17 investigation, a study protocol, a method of selecting
18 subjects, a method of observation and recording
19 results and comparison of results with a control.

20 In conclusion, the Panel has an important
21 job today. You must listen to the data presented by
22 the sponsor, evaluate the FDA presentations and make

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 a recommendation about the approvability of the
2 sponsor's application. We speak for many applicants
3 when we ask for your careful consideration. Please,
4 keep in mind that the standard is a reasonable
5 assurance balancing the benefits with the risks. The
6 regulatory standard is not proof beyond a shadow of a
7 doubt.

8 When considering making recommendations
9 for further studies, remember that FDA takes these
10 recommendations seriously, often as a consensus of the
11 Panel of a whole and such recommendations may delay
12 the introduction of a useful product or result in
13 burdensome and expensive additional data collection.
14 Therefore, you play an important role in reducing the
15 burden of bringing new products, products that you and
16 your colleagues use in treating patients to the
17 market.

18 Please, be thoughtful in weighing the
19 evidence. Remember that the standard is a reasonable
20 assurance of safety and effectiveness and that there
21 is a legally broad range of valid scientific evidence
22 to support that determination. On behalf of OSMA, I

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 would like to thank the FDA and the Panel for the
2 opportunity to speak today. Our association trusts
3 that its comments are taken in the spirit offered,
4 which is to help the FDA decide whether to make a new
5 product available for use in the U.S. marketplace.
6 Thank you.

7 CHAIRPERSON YASZEMSKI: Thanks very much,
8 Ms. Adams. We're going to break now and then proceed
9 with the summation from both the FDA and the sponsors.
10 It's about 3:51. Let's come back and start at 10
11 minutes past 4, 4:10.

12 (Whereupon, at 3:52 p.m. a recess until
13 4:14 p.m.)

14 CHAIRPERSON YASZEMSKI: Can I ask
15 everybody to take your seats, please? We're going to
16 get started. We're going to ask, at this time, for
17 the FDA and sponsor summations and then we're going to
18 proceed to the voting. And I will first ask FDA. Dr.
19 Witten, would FDA like to add anything, and I would
20 specifically like to ask you to comment on the rules
21 regarding conditions and their effect on the vote?

22 DR. WITTEN: I'm sorry. You're asking me

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to clarify the rules on post-approval commissions and
2 the vote?

3 CHAIRPERSON YASZEMSKI: Yes.

4 DR. WITTEN: Yes. Okay. Thank you for
5 asking me. I will just amplify what I said before the
6 break, which is that if there is a condition that is
7 asking for new data or a new analysis, if the request
8 is for that new data or new analysis to be provided to
9 us after approval to answer some focused, specific
10 question or a series of questions then that, you know,
11 would be what we would consider a condition of
12 approval.

13 If what you are requesting or what the
14 Panel recommends is a condition where you're asking
15 for new data or a new analysis, that you want it
16 provided to us for our review prior to approval, then
17 what that is to us in terms of the vote and the
18 recommendation is a non-approvable recommendation, and
19 what you're providing us with are the recommendations
20 of how to put the application in approvable form.

21 So that's why I had said prior to the
22 break when Dr. Kirkpatrick was going through his list

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that we need to understand whether, for each of these
2 that you may agree on, the Panel is recommending that
3 we have the data in hand to review prior to approval,
4 which would mean you're really making a non-approvable
5 recommendation with a recommendation of how to put it
6 in approvable form versus telling us you would like us
7 to look at a specific, focused question and get some
8 data around those questions after approval, which
9 would be a post-approval condition.

10 And let me just clarify one additional
11 thing, in case the question should come up, which is,
12 you know, under what would we take such an application
13 back to Panel and it would be our option of whether or
14 not we felt that we had, you know, additional issues
15 we wanted to ask the Panel about.

16 CHAIRPERSON YASZEMSKI: Okay. Thank you,
17 Dr. Witten. I would like to ask if the sponsor has
18 any summation comments to make, Mr. Christianson?

19 MR. CHRISTIANSON: Thank you, Dr.
20 Yaszemski. Jack, could I have the first slide,
21 please? We heard some discussion today that several
22 of the Panel Members expressed concerns that we don't

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 know what the long-term safety profile of the Charite
2 Artificial Disc is, and I would just like to remind
3 the Panel that we did conduct a 24 month randomized
4 prospective study per the FDA guidance document on
5 spinal devices, and that has been used for all
6 previous spinal devices and, indeed, all previous
7 orthopedic devices have been approved based on a two
8 year follow-up study.

9 In addition, several people did remind the
10 Panel during the course of the discussion that there
11 are, indeed, long-term follow-up data from Europe.
12 There is a very good case series from Dr. LeMaire and
13 a six year case series from Dr. David, and case series
14 do meet the FDA definition of valid scientific
15 evidence, so the Panel can, indeed, take those series
16 into consideration that are available in the
17 literature.

18 And we also heard some discussion about a
19 post-approval study. Indeed, that's the place that's
20 appropriate and accepted to develop the longer term
21 data in a post-approval study after the device has
22 been approved and the company certainly is amenable to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 conducting a five year follow-up study as Dr.
2 Kirkpatrick recommended in his document. Next slide,
3 please.

4 Sorry, Dr. Kirkpatrick, we didn't know
5 what to title this slide when we put it together, but
6 reviewing your list of recommendations that you passed
7 around, we agree that most of the recommendations on
8 your list are reasonable and we will certainly discuss
9 them with FDA and put the answers together that we can
10 with our existing data.

11 However, I must comment on your
12 recommendation for 50 million cycle testing. The
13 company believes that that is an excessive requirement
14 for testing. For example, for metal on a polyethylene
15 device, that will take at least 15 weeks to conduct,
16 probably longer, and that would potentially delay the
17 approval of the device for a significant period of
18 time.

19 The testing that we have submitted, the 10
20 million cycle testing that is already in our PMA, does
21 represent 80 years of significant bends while listing
22 a 20 kilogram weight, so we do think that we provided

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 adequate long-term mechanical test data. And so if
2 that's an issue that we'll need to negotiate with FDA,
3 I wanted to get that statement on the record. Last
4 slide, please, Jack.

5 And I'll close with the same statements
6 that I made when we closed our Panel presentation.
7 We're presenting a device to you that's got a long
8 clinical history of use in Europe, fully
9 biomechanically characterized, robust, valid,
10 scientific evidence that the device is safe and
11 effective and we, again, ask the Panel to recommend
12 that this device be approved for use in patients in
13 the U.S.

14 CHAIRPERSON YASZEMSKI: Thank you, Mr.
15 Christianson. Ms. Scudiero will now read the three
16 possible Panel recommendation options for pre-market
17 approval applications. Ms. Scudiero?

18 MS. SCUDIERO: These were in the meeting
19 handouts. They are entitled "Panel Recommendation
20 Options for Pre-Market Approval Applications." The
21 medical device amendments to the Federal Food, Drug
22 and Cosmetic Act, as amended by the State Medical

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Devices Act of 1999, allows the Food and Drug
2 Administration to obtain a recommendation from an
3 expert Advisory Panel on designated medical device
4 pre-market approval applications that are filed with
5 the Agency.

6 The PMA must stand on its own merits and
7 your recommendation must be supported by the safety
8 and effectiveness data in the application or by
9 applicable publicly available information. Safety is
10 defined in the Act as "Reasonable assurance based on
11 valid, scientific evidence that the probable benefits
12 of health under the conditions on intended use
13 outweigh any probable risk." Effectiveness is defined
14 as "A reasonable assurance that in a significant
15 portion of the population, the use of the device for
16 its intended uses and conditions of use, when labeled,
17 will provide clinically significant results."

18 Your recommendation options for the vote
19 are as follows: (1) Approval, if there are no
20 conditions attached. (2) Approvable with conditions.
21 The Panel may recommend that the PMA be found
22 approvable subject to specified conditions, such as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 physician or patient education, labeling changes or
2 further analysis of existing data. Prior to voting,
3 all the conditions should be discussed by the Panel.

4 (3) Non-approvable. The Panel may recommend that the
5 PMA is not approvable if the data do not provide a
6 reasonable assurance that the device is safe or if a
7 reasonable assurance has not been given that the
8 device is effective under the conditions of use
9 prescribed, recommended or suggested in the proposed
10 labeling.

11 Following the voting, the Chair will ask
12 each Panel Member to present a brief statement
13 outlining the reasons for their vote, and this became
14 effective June 14, 1999.

15 CHAIRPERSON YASZEMSKI: Thanks, Ms.
16 Scudiero. I would like to make a few comments before
17 we ask for a motion. First, with respect to voting,
18 the eight Panel Members will vote. Our consumer
19 representative, our patient representative, that is,
20 and our industry representative will not vote. I will
21 only vote in the event of a tie.

22 Regarding the motion, the sequence that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 can happen is we have a motion, a second for the
2 motion, discussion and a vote. If that sequence
3 occurs and the vote is for the motion, then we're
4 finished. If the vote is either against the motion or
5 if there is not a second for the motion, then we'll
6 ask for another motion.

7 With that in mind, the lead reviewer for
8 this was Dr. Kirkpatrick and I'm going to ask him to
9 make a motion. Dr. Kirkpatrick?

10 DR. KIRKPATRICK: To borrow from Dr.
11 Hochschuler, I felt a little bit like Simon at the
12 beginning of my discussion and I hope that we can
13 understand each other as far as where we're coming
14 from.

15 A recent editorial in the NAS Journal
16 indicated that I am part of an increasing or a
17 decreasing majority of spine surgeons. The editorial
18 was discussing the fact that there is a number of
19 spine surgeons who will do things on patients that
20 they would never consider for themselves. This
21 reminds me of what the FDA's purpose is and that is,
22 first, protecting the public. As such, some of my

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 comments and my motion will be directed towards that.

2 The second rationale I have for my motion
3 is being a bone setter from Alabama means I have
4 adopted certain habits and customs. One of those
5 customs happens to be watching NASCAR. As many of you
6 know, NASCAR is a race around a track that generally
7 runs between 250 and 500 miles or sometimes, on one
8 occasion a year, 600 miles in length. The design of
9 the tires is specific for the track and the type of
10 racing done and is not expected to exceed the length
11 of time that the gas tank is full. In other words,
12 when you run out of gas, your tires are going to need
13 to be replaced.

14 I think we need to think of the same
15 design rationale as far as a disc replacement. We
16 need to make sure that we are assured of both the
17 safety and effectiveness for the intended length of
18 use. Now, I know that's an onerous thing if we're
19 talking 50 years, and I don't propose that at all.
20 However, I do think that a two year follow-up, in all
21 due respect to Mr. Christianson and his colleagues as
22 far as discussion of precedent, this is an

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 unprecedented device and I don't think two years is
2 adequate.

3 As such, I would recommend or, excuse me,
4 I move that we call this PMA not approvable and that
5 would be my motion.

6 CHAIRPERSON YASZEMSKI: Thank you, Dr.
7 Kirkpatrick. Do we have a second for that motion?
8 Dr. Finnegan?

9 DR. FINNEGAN: Yes.

10 CHAIRPERSON YASZEMSKI: We have a second.
11 Discussion?

12 MS. MAHER: Well, can I lead off the
13 discussion?

14 CHAIRPERSON YASZEMSKI: Ms. Maher?

15 MS. MAHER: As a non-voting member, I can
16 lead this off pretty well. I have to take exception,
17 Dr. Kirkpatrick, to what you're saying, because spinal
18 cages were approved initially with two years follow-up
19 and they also, at the time they came on board, were a
20 first of their kind. And if I actually recall
21 correctly, they even had much less animal data and
22 other data than we have on this product, which, again,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 as Mr. Christianson put up on the slide, we do have
2 data since 1987 showing that it has been used safely
3 and effectively in Europe for many years.

4 So I have some deep concerns that if you
5 tell a company they can't launch something for five
6 years after they have started developing it, we're
7 going to put a stop to new product innovation in the
8 medical device or the orthopedic world. And I'm
9 wondering why you feel that that's more appropriate
10 than having a post-market approval study, a post-
11 approval study where you can follow the devices and
12 look at what's happening.

13 You have got a cohort of patients that
14 already has two years. You can have three more years
15 and you will have the five year data, in which case
16 you'll have the other patients. It will be available
17 for sale. It will be being sold and being used, but
18 I think they have provided adequate evidence that it
19 would be safe and effective, so I have to disagree
20 with you.

21 CHAIRPERSON YASZEMSKI: Thank you, Ms.
22 Maher. Dr. Diaz?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. DIAZ: I also would like to disagree
2 with Dr. Kirkpatrick's statement, because I think we
3 are making a statement that flies against a very large
4 body of evidence. There has been 17 years of use of
5 these devices throughout the world. We are the only
6 country in the industrialized world that does not
7 approve its use yet.

8 To expect to compare a mechanical device
9 like a tire that is running on a NASCAR track to the
10 function of the human body is counter-intuitive. If
11 the good Lord had designed our gas tank to allow our
12 functioning parts to last the exact same time, we
13 would die in perfect physical condition and that is
14 not a reality. We run out of gas at about the same
15 time when all our parts have fallen apart.

16 So I think that the motion is not what I
17 would endorse. I disagree entirely that there is not
18 sufficient evidence to indicate its use. It can be
19 done, I believe, with some continuing monitoring and
20 perhaps longitudinal studies to answer some of the
21 questions, but I believe the experience in France and
22 Germany have already shown the various things that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 will happen.

2 And if we were to say that ultimately, all
3 the device does is delay the occurrence of a natural
4 fusion, as was presented already in a relatively small
5 comment made at the end of the presentation of Dr.
6 Ooij, I believe is his name, from the Netherlands,
7 Ooij, Dr. Ooij. Even if we gain 10 years of extension
8 on the function of a disc, I think we have provided a
9 sufficient opportunity for the individual not to have
10 a fusion that perhaps would occur spontaneously or
11 that may be induced by the introduction of mechanical
12 devices, which we have already approved. So I think
13 that the decision not to allow it is incorrect.

14 CHAIRPERSON YASZEMSKI: Thank you, Dr.
15 Diaz. Dr. Mabrey?

16 DR. MABREY: Well, having trained in North
17 Carolina at the same time that Dr. Kirkpatrick did at
18 Duke, I can certainly share his observations of NASCAR
19 as a NASCAR dad, but I do see an opportunity here to
20 provide additional information after approval. The
21 developers of the device and the clinicians have
22 demonstrated that they have a very good cohort of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 patients. They have gone out of their way to document
2 every complication that has occurred, and I think we
3 have the opportunity to follow that data after
4 approval and look at that data both at four and five
5 years out.

6 John, I agree. I think, you know, maybe
7 50 million cycles isn't an unreasonable number of
8 cycles to go through, but, again, that type of data
9 could be ongoing rather than pre-approval. So I would
10 argue for approval with certain conditions.

11 CHAIRPERSON YASZEMSKI: Thank you, Dr.
12 Mabrey. Dr. Finnegan, you were the seconder.

13 DR. FINNEGAN: Yes, I'm not sure that some
14 of the Panel Members and maybe Ms. Maher will
15 understand that just because we say not approval
16 doesn't mean this is going into the closet. Not
17 approval means that, at the present time, the Panel is
18 not comfortable with all of the data. It does not
19 mean it has to come back to Panel. It just means that
20 certain things have to be done before the FDA makes a
21 decision and that it's quite possible, given the
22 discussion today, that it will not need to come back

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to Panel. But if there are some things that we really
2 feel strongly need to be done before the FDA gives it
3 approval, then by regulation we cannot approve it.

4 Now, one of the biggest concerns is that
5 of those two year follow-up patients, they haven't all
6 reached two year follow-up and, in fact, if I read the
7 numbers correctly, in fact, the latest patient to get
8 this is probably less than 10 months ago. So if you
9 take all of those patients out to two years, you're
10 actually going to have some three and a half or four
11 year data, which will be much more helpful than doing
12 it now when some of the patients haven't reached their
13 two year mark.

14 CHAIRPERSON YASZEMSKI: Thanks, Dr.
15 Finnegan. Dr. Kim?

16 DR. KIM: This is such a difficult topic
17 to vote on because of the complexity of the disease
18 and the fact that this is a brand new product, but I
19 was reading the FDA Modernization Act of 1997 and what
20 that Act basically entails is the spirit of trying to
21 promote innovation, and I think by requiring much
22 longer follow-up, it will deter companies from being

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 able to produce these innovative materials and I think
2 the burden will be too onerous.

3 So I think that the two year clinical
4 data, which is excellent, combined with the long-term
5 follow-up of the European literature, which, as Ms.
6 Maher pointed out, is data that we can use as an FDA
7 Panel to make decisions, are compelling and I would
8 lean toward approval with specific conditions that
9 addresses some of the concerns that we have.

10 CHAIRPERSON YASZEMSKI: Thank you, Dr.
11 Kim. Dr. Naidu?

12 DR. NAIDU: You know, after listening to
13 the presentation from the Netherlands, the physician
14 from the Netherlands as far as device complications,
15 it appears as if device related complications
16 including anterior/posterior migration is less than 1
17 percent. In addition, the sponsor has conducted an
18 excellent study where they have shown a significant
19 improvement in objective outcomes, including the ODI
20 and the VAS.

21 They have also shown that it's non-
22 inferior to BAK, that's the fusion device, and they

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 have also shown that it's at least equivalent and it's
2 not inferior, and I'm not sure as to why we are
3 debating as to approvability of this device. I think
4 it's approvable without any conditions.

5 CHAIRPERSON YASZEMSKI: Thank you, Dr.
6 Naidu. Dr. Blumenstein?

7 DR. BLUMENSTEIN: I can't go along with
8 disapproval. I have to think that there are some
9 conditions that we could put on with an approval with
10 conditions that would satisfy the long-term follow-up
11 requirement.

12 CHAIRPERSON YASZEMSKI: Thank you. Dr.
13 Besser?

14 DR. BESSER: I would also look for
15 approval with conditions. I think we can resolve the
16 issues here post-approval.

17 CHAIRPERSON YASZEMSKI: Thank you, Dr.
18 Besser. Would anybody else like to add commentary?
19 Hearing none, what we're going to do now is vote on
20 this motion, which is for non-approval. I will go
21 around the room and ask everybody to say yes or no for
22 non-approval. If you vote yes, that means you agree

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 and you would like this to be not approved.

2 If this motion passes, then we're finished
3 and our recommendation is non-approval and we'll
4 discuss after that conditions that need to be met to
5 make it approvable. If the motion does not pass, then
6 we will ask Dr. Kirkpatrick if he might entertain a
7 new motion.

8 Let's start, Dr. Diaz, with you.

9 DR. DIAZ: I disagree with the motion.

10 CHAIRPERSON YASZEMSKI: Thank you. Dr.
11 Mabrey?

12 DR. MABREY: Disagree.

13 CHAIRPERSON YASZEMSKI: Dr. Finnegan?

14 DR. FINNEGAN: I agree.

15 CHAIRPERSON YASZEMSKI: Dr. Kim?

16 DR. KIM: I disagree.

17 CHAIRPERSON YASZEMSKI: Dr. Naidu?

18 DR. NAIDU: I disagree.

19 CHAIRPERSON YASZEMSKI: Dr. Kirkpatrick?

20 DR. KIRKPATRICK: I agree.

21 CHAIRPERSON YASZEMSKI: Dr. Blumenstein?

22 DR. BLUMENSTEIN: Disagree.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 CHAIRPERSON YASZEMSKI: Dr. Besser?

2 DR. BESSER: Disagree.

3 CHAIRPERSON YASZEMSKI: The motion does
4 not pass. Our two remaining options are approval or
5 approval with conditions and I will ask, at this time,
6 Dr. Kirkpatrick, would you entertain another motion?

7 DR. KIRKPATRICK: I would be glad to. I
8 would also like to take a moment to recognize the
9 beauty of democracy and the fact that we can agree to
10 disagree, and that we have the freedom to do so at the
11 expense of a number of our countrymen right now.

12 I would suggest, I would like to make the
13 motion that it is approvable with conditions and if
14 that passes, I would like to itemize conditions and
15 take them individually if that's okay with the Chair.

16 CHAIRPERSON YASZEMSKI: Yes, the way that
17 we're going to do it is we're going to go around. If
18 you make a motion for approval with conditions, what
19 we now do is consider the conditions first. I'm
20 sorry, a point of order. Ms. Scudiero just reminded
21 me that I did not ask for a second.

22 Dr. Kirkpatrick has made a motion for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 approval with conditions. Do I have a second? Dr.
2 Besser has seconded the motion. Thank you, Ms.
3 Scudiero.

4 What we'll do prior to voting is have a
5 discussion and ask for conditions. And if someone
6 brings a condition up, we'll discuss that condition
7 and then vote on that condition, and that condition
8 will then either be included or not included. If
9 persons have disagreements with the conditions that
10 get voted in, then they can exercise their
11 disagreement by voting no for the motion when it comes
12 to a vote.

13 So I would like to entertain now if there
14 is a motion for a condition from anybody. Yes, Dr.
15 Kirkpatrick?

16 DR. KIRKPATRICK: In a follow-up to my
17 concern about the length of follow-up, I would suggest
18 that a condition would be that all of the currently
19 enrolled patients, including the -- I can't remember
20 what you termed it, but basically the patients that
21 aren't in the IDE, but the ones that have continued to
22 be done.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 UNIDENTIFIED SPEAKER: The continued
2 access.

3 DR. KIRKPATRICK: The continued access
4 group be followed to the last of the continued access
5 group being a minimum of two years follow-up. That
6 should give us close to five years on most of the IDE
7 patients if I'm remembering correctly on your block of
8 time.

9 CHAIRPERSON YASZEMSKI: Yes. Thank you.

10 DR. KIRKPATRICK: That would be the first
11 of several conditions I would propose.

12 CHAIRPERSON YASZEMSKI: We have a motion
13 for a condition to include all the continued access
14 patients until they have completed two year follow-up.
15 Is there a second for this motion?

16 DR. DIAZ: Second.

17 CHAIRPERSON YASZEMSKI: Dr. Diaz and Dr.
18 Finnegan, we have seconds. But is there discussion on
19 this motion? Dr. Besser?

20 DR. BESSER: I'm questioning as to whether
21 two years is long enough. I'm not sure how many of
22 the people will be out to five years at that two years

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 after the last patient. If, in fact, we're looking
2 for data out to five years, I would like to see the
3 last patient at five years and that would give us even
4 longer data and better data for the rest.

5 CHAIRPERSON YASZEMSKI: Mr. Christianson,
6 could I ask you or a member of your company to comment
7 on this question from Dr. Besser?

8 MR. CHRISTIANSON: Yes, the first patient
9 was enrolled in 2000 and the last continued access
10 patient was enrolled last week. So if we follow that
11 patient through five years, the patient from 2000,
12 someone do the math for me quick, is going to be
13 extensive.

14 UNIDENTIFIED SPEAKER: Nine years out.

15 MR. CHRISTIANSON: So I believe that the
16 entire randomized cohort will be at or beyond five
17 years if we followed the last continued access patient
18 through two years.

19 CHAIRPERSON YASZEMSKI: Okay. Thank you,
20 Mr. Christianson. Dr. Besser, does that answer your
21 question?

22 DR. BESSER: It answers my question, but

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 I'm not sure I'm convinced to shorten that. I would
2 still like -- you know, nine years of data would be
3 great.

4 CHAIRPERSON YASZEMSKI: Okay. Thank you.
5 Others?

6 DR. KIRKPATRICK: May I amend my --

7 CHAIRPERSON YASZEMSKI: No, we have a
8 second already.

9 DR. KIRKPATRICK: Okay.

10 CHAIRPERSON YASZEMSKI: We'll have to vote
11 on it first. Any other discussion on this point?

12 DR. BESSER: I believe you can withdraw a
13 motion.

14 DR. KIRKPATRICK: I don't want to withdraw
15 it.

16 CHAIRPERSON YASZEMSKI: Let's have
17 discussion.

18 DR. KIRKPATRICK: I want to include yours.

19 CHAIRPERSON YASZEMSKI: No, we can make
20 another one. More discussion? Seeing none, Dr. Diaz,
21 you're in the number one position here. I'm going to
22 keep asking you first.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Let's vote, vote on this condition. The
2 condition is to include all the continued access
3 patients until they have completed two years of
4 follow-up. This would be a condition to a motion for
5 approval with conditions. This is the first
6 condition.

7 DR. KIRKPATRICK: That would be a report
8 on all patients once the last of the continued access
9 reaches two years.

10 CHAIRPERSON YASZEMSKI: Yes, that's
11 assumed.

12 DR. KIRKPATRICK: I just wanted to clarify
13 that.

14 CHAIRPERSON YASZEMSKI: That's assumed.
15 Yes, sir, Dr. Diaz?

16 DR. DIAZ: I agree.

17 CHAIRPERSON YASZEMSKI: Dr. Mabrey?

18 DR. MABREY: I agree.

19 CHAIRPERSON YASZEMSKI: Dr. Finnegan?

20 DR. FINNEGAN: I agree.

21 CHAIRPERSON YASZEMSKI: Dr. Kim?

22 DR. KIM: I hate to do this, but wouldn't

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 it be better to follow the IDE patients that are
2 randomized for a total of five years?

3 CHAIRPERSON YASZEMSKI: We can do that as
4 a separate motion. I think we need to finish voting
5 here.

6 DR. KIM: So based on that, I would
7 disagree.

8 CHAIRPERSON YASZEMSKI: Okay. Thank you.
9 Dr. Naidu?

10 DR. NAIDU: I agree.

11 CHAIRPERSON YASZEMSKI: Dr. Kirkpatrick?

12 DR. KIRKPATRICK: Agree.

13 CHAIRPERSON YASZEMSKI: Dr. Blumenstein?

14 DR. BLUMENSTEIN: Disagree.

15 CHAIRPERSON YASZEMSKI: Dr. Besser?

16 DR. BESSER: Disagree.

17 CHAIRPERSON YASZEMSKI: And this motion
18 passes 5 to 3 and so one condition for Dr.
19 Kirkpatrick's motion for approval with conditions is
20 that the currently enrolled patients in the continued
21 access category and all other enrolled patients in the
22 IDE study have follow-up at the time that the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 continued access patients reach two years follow-up.

2 Now, would anybody like to introduce a
3 second condition? Dr. Finnegan?

4 DR. FINNEGAN: It seems to me that this is
5 the ultimate device for device tracking and I would,
6 therefore, like to introduce the condition that this
7 device be tracked.

8 CHAIRPERSON YASZEMSKI: We have a motion
9 to include a condition for device tracking. I would
10 like to ask Dr. Witten to comment on the device
11 tracking condition.

12 DR. WITTEN: Well, I just would like --
13 that term is always really confusion, and so I would
14 like clarification as to what exactly that means,
15 whether it's that we want to be able to track the
16 device to the patients or there are specific data
17 elements we want when it gets implanted. Is this all
18 patients and, for example, maybe I could start with
19 asking what the objective would be and then we could
20 better understand what it is.

21 DR. FINNEGAN: See, the objective would be
22 as this is put in, one of the unfortunate problems at

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 present is when a device is put in and five years
2 later significant complications are known and
3 probably, what am I thinking about, the -- what's the
4 hip that got -- anyway, we have reason to -- yes, the
5 Saltzer Hip, that we have had recent experience that
6 patients are panicking, lawyers are calling everybody
7 to find out if they have got the device in them or not
8 and no one has the answer, which maybe the sponsor
9 thinks would be a good idea.

10 But anyway, what we're looking for is a
11 way that a patient would know what device was in them.
12 The physician would know what patients they had
13 implanted the device in. The sponsor would know that
14 the device was in Patient X, so that when, long-term,
15 something came up, you would know where to go.

16 DR. WITTEN: Okay. So it's to identify
17 the patients and the physicians, and that's for
18 anybody who receives the implant. It's not a data
19 collection mechanism?

20 DR. FINNEGAN: That is correct.

21 DR. WITTEN: Okay.

22 CHAIRPERSON YASZEMSKI: Thank you. We

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 have a motion for device tracking. Is there a second?
2 Dr. Mabrey. Discussion? Dr. Diaz?

3 DR. DIAZ: I would like to just make sure
4 that we call it device ID follow-up rather than
5 tracking, because tracking to me implies something
6 very different. To me that means a responsibility on
7 the corporation to follow every single device that's
8 implanted wherever it happens to end, and I think
9 that's an onerous condition of its approval. I think
10 it is important for the patient and the surgeon to
11 know what device was implanted in whom when and where
12 and leave it at that.

13 CHAIRPERSON YASZEMSKI: Okay. Thanks, Dr.
14 Diaz. May I go out of the order first and I would
15 like commentary from Ms. Maher on this.

16 MS. MAHER: Yes, I would support what Dr.
17 Diaz just said. Device tracking, when you go to the
18 degree as to what that term actually means, is
19 exceedingly burdensome to the industry and in the days
20 of HIPAA is almost going to be impossible to do. It's
21 not one of those things that patients want to be
22 followed and want to be tracked. You know, they move.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 They don't tell you they have moved. Keeping track of
2 where they are is virtually impossible.

3 The device tracking requirements were
4 originally put into place for products where if they
5 were to fail, such as the heart valves, it would be
6 catastrophic to the patient immediately and I don't
7 see that this product necessarily fits that definition
8 of being catastrophic immediately.

9 I like what Dr. Diaz suggested, that we
10 actually train people more and have the labeling
11 require more, that the patient is supposed to know
12 that you have gotten a DePuy Charite Disc. I know
13 many people who have gotten joint replacements and
14 have no idea what joint they had placed in them, which
15 I also find bizarre.

16 But I think that if you go to this Nth
17 degree, you're adding a burden that is almost
18 impossible to meet and I'm not sure I see the benefit
19 of it, especially given what that regulation and law
20 was originally intended for.

21 CHAIRPERSON YASZEMSKI: Thanks. And, Dr.
22 Diaz, may I ask for a clarification on your use of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 term follow-up? Would it be similar to what Ms. Maher
2 has just said? What would you suggest the follow-up
3 be?

4 DR. DIAZ: I think it should be limited
5 only to providing the patient with a name of the
6 device, perhaps an ID number that all of these devices
7 have. The patient would have the name of the surgeon,
8 the place where the surgery was done and the date and
9 leave it at that.

10 CHAIRPERSON YASZEMSKI: Thank you. Now,
11 I haven't yet asked for a second and I would like to--

12 UNIDENTIFIED SPEAKER: You have a second.

13 CHAIRPERSON YASZEMSKI: I have a second?
14 I'm sorry. I did so.

15 DR. DIAZ: This is just a friendly
16 editorial amendment.

17 CHAIRPERSON YASZEMSKI: Yes. May I come
18 first to Dr. Mabrey and then we're going to come to
19 Ms. Luckner?

20 MS. LUCKNER: From the patient's
21 perspective, I think you are asking for patient
22 identification to know what device was implanted and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to know the number. That is totally reasonable from
2 a patient perspective. I do not wish to reveal
3 totally in this room, but I have two knee
4 replacements. I carry in my wallet that I have a knee
5 replacement, so that I have no difficulty with going
6 through airport security systems.

7 Now, I will tell you it does not say on it
8 the manufacturer and I am one of those people that,
9 over here my colleague said, I have no idea what knee
10 replacement I have. Listening to this conversation,
11 you can believe when I return to Toledo, Ohio, I will
12 find out exactly what is in my knees.

13 CHAIRPERSON YASZEMSKI: Thanks, Ms.
14 Luckner. Dr. Mabrey?

15 DR. MABREY: Yes, I agree with Ms. Maher's
16 comments that a tracking type of program would be
17 somewhat onerous and my chief concern is trying to
18 keep in line with all the HIPAA regulations. I think
19 that becomes a quagmire, if I can borrow from another
20 era.

21 However, every one of my total joint
22 patients get a card and they know exactly what implant

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 they have in them and I make sure they have it, but
2 then again, not everybody puts in total joints. I
3 think it's reasonable to provide the patient not only
4 with a card that identifies what implant they have and
5 the date it was implanted, but also a serial number
6 much along the lines of the pacemakers.

7 I believe most pacemaker implants have a
8 serial number associated with them. The manufacturer
9 will be keeping at least a registry of those serial
10 numbers and should there ever be a problem with that
11 group, it seems like it would be a simple matter for
12 the patient then to take the initiative and contact
13 the physician or the company to follow-up on that
14 serial number.

15 It's certainly also helpful. I would love
16 to have every total hip patient and total knee patient
17 in the country carrying around their serial numbers,
18 so that I could call up that company and find out
19 exactly what size implant I'm going to revise, and I
20 don't see that as being too onerous on industry or too
21 onerous on the patients and, certainly, you know,
22 would put everyone's mind at ease.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 CHAIRPERSON YASZEMSKI: Along the lines of
2 what you have just discussed, would it be reasonable
3 to consider having a card come with every prosthesis
4 that has that identification number and then which the
5 surgeon just fills out his or her name, the patient's
6 name, date and hospital?

7 DR. MABREY: Well, I don't put these
8 devices in, but it looks like there's three parts to
9 it and they all come in three separate boxes. But
10 there are peel-off stickers that could go with that.
11 I would only caution you that if they go into your
12 wallet, then after a couple of years the numbers will
13 probably wear off. So I'm not sure how we would
14 handle that.

15 CHAIRPERSON YASZEMSKI: All right.
16 Thanks, Dr. Mabrey. Dr. Blumenstein?

17 DR. BLUMENSTEIN: Well, as a concerned
18 father of teenagers, perhaps we could put the tattoo
19 parlors to work and have them --

20 CHAIRPERSON YASZEMSKI: Thank you, Dr.
21 Blumenstein. Dr. Kirkpatrick?

22 UNIDENTIFIED SPEAKER: They would be too

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 busy to work on your kids.

2 CHAIRPERSON YASZEMSKI: Dr. Kirkpatrick?

3 DR. KIRKPATRICK: May I suggest that
4 implied in this condition would be that the FDA would
5 work with the manufacturer in order to make sure that
6 there is a legal way to do this, because the Joint
7 Registry is already significantly alone working on
8 those problems. FDA is aware of those issues and if
9 it can be done, it can, but if it can't be done, that
10 the Panel would accept that, but we would encourage it
11 to be done. Is that implicit in the motion?

12 DR. FINNEGAN: It is implicit in the
13 motion, but I think, Jay, I don't want to give people
14 sort of a cop-out, because I think Jay is right. I
15 think if you put the serial number on and the patient
16 has the access to the serial number, that doesn't have
17 to have a lot of data and then you can do the same
18 thing that, you know, Mercedes does when their brakes
19 don't work. They say if your car, you know, has this
20 serial number or was between this and this date, then
21 you need to call the company. Here is the 1-800
22 number. I mean, that's a pretty straightforward thing

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 to do.

2 DR. KIRKPATRICK: I agree that it's
3 straightforward, but, believe it or not, that minimal
4 of a data set is significantly complicated in trying
5 to get through other federal agencies as far as
6 whether it's HIPAA compliant.

7 DR. FINNEGAN: But if the patient signs
8 the consent, I don't think it is.

9 DR. KIRKPATRICK: Exactly.

10 DR. FINNEGAN: So if the patient signs, it
11 says that they are quite happy to have the serial
12 number and to have the company know what serial number
13 they have, then I don't think that's --

14 CHAIRPERSON YASZEMSKI: Okay. Thanks.
15 Dr. Mabrey?

16 DR. MABREY: Well, if I could just
17 clarify. I don't even think that we're asking that
18 the company know which patient has it, but just that
19 the company know what implants are out there and have
20 been implanted, and we're placing part of the -- yes,
21 we're putting part of the responsibility on the
22 patient now to look at the serial number and then get

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in contact with industry. I think this then keeps us
2 out of all the problems with HIPAA. Carrying around
3 a serial number that only you know you have and
4 industry having that same serial number, but having no
5 idea who you are, I think that's reasonable, and I
6 think it may be reasonable to at least keep track of
7 which physician put it in.

8 And I know the way industry does the total
9 joints, they certainly know what region it goes into
10 and I know that our distributor keeps track of just
11 about every implant I have put in anyway. I will have
12 to check and see what HIPAA rules we're violating on
13 that when I get back though.

14 CHAIRPERSON YASZEMSKI: All right. Thanks
15 very much. Dr. Finnegan, do you have additional
16 comments? Otherwise, I'll go to Dr. Kim.

17 DR. FINNEGAN: I have no other comments.

18 CHAIRPERSON YASZEMSKI: Okay. Dr. Kim?

19 DR. KIM: I would agree with what Dr.
20 Mabrey said.

21 CHAIRPERSON YASZEMSKI: Thanks. Dr.
22 Naidu?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. NAIDU: Same here.

2 CHAIRPERSON YASZEMSKI: Thanks. Dr.
3 Kirkpatrick, additional comments?

4 DR. KIRKPATRICK: No.

5 CHAIRPERSON YASZEMSKI: Dr. Blumenstein?

6 DR. BLUMENSTEIN: No comments.

7 CHAIRPERSON YASZEMSKI: Dr. Besser?

8 DR. BESSER: No comments.

9 CHAIRPERSON YASZEMSKI: Thanks. Ms.
10 Maher?

11 MS. MAHER: I would just like a little
12 clarification. We're probably talking about lot
13 numbers here, not serial numbers, and I would
14 recommend that we actually, since the FDA now knows
15 from this details conversation that what we really
16 want is for patients to know what device they have had
17 implanted and what lot numbers it was, that we leave
18 it to the FDA and the sponsor to work out the best way
19 to obtain that information.

20 DR. MABREY: I think lot numbers are --

21 CHAIRPERSON YASZEMSKI: Okay. All right.

22 Thank you. Now, before we vote, I want to ask Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Finnegan, because you made the motion and you did
2 start with the words device tracking, but in light of
3 the discussion, would you be okay if the motion did
4 not include those specific words, which impose a
5 certain level of --

6 DR. FINNEGAN: All I want --

7 CHAIRPERSON YASZEMSKI: -- responsibility
8 on the company, but to go along with what we
9 discussed.

10 DR. FINNEGAN: All I want is the patient
11 to be able to know if the implant they have has a
12 problem.

13 CHAIRPERSON YASZEMSKI: Okay. Thank you.
14 I will state the motion then in the form it is after
15 discussion. The motion is that the patients be
16 supplied with the name and lot number of the device,
17 the doctor and hospital and date that it was put in,
18 that the company know only that the device was
19 implanted and that if problems do arise, the company
20 can send out a notice and it would be the patient's
21 responsibility to recognize that they have one of the
22 implants in them that was in the notice.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 We're going to vote on this now. Dr.
2 Diaz?

3 DR. DIAZ: I agree.

4 CHAIRPERSON YASZEMSKI: Dr. Mabrey?

5 DR. MABREY: I agree.

6 CHAIRPERSON YASZEMSKI: Dr. Finnegan?

7 DR. FINNEGAN: I better agree.

8 CHAIRPERSON YASZEMSKI: Dr. Kim?

9 DR. KIM: I agree.

10 CHAIRPERSON YASZEMSKI: Dr. Naidu?

11 DR. NAIDU: I agree.

12 CHAIRPERSON YASZEMSKI: Dr. Kirkpatrick?

13 DR. KIRKPATRICK: Agree.

14 CHAIRPERSON YASZEMSKI: Dr. Blumenstein?

15 DR. BLUMENSTEIN: Agree.

16 CHAIRPERSON YASZEMSKI: Dr. Besser?

17 DR. BESSER: I agree.

18 CHAIRPERSON YASZEMSKI: This motion passes
19 as the second condition of the motion for approval
20 with conditions.

21 We will now move on and ask if there are
22 other conditions that people would like to raise and

1 include in the motion for approval with conditions.

2 Dr. Kirkpatrick?

3 DR. KIRKPATRICK: I need a clarification
4 before I make my motion. When Mr. Christianson said
5 that it would be 15 weeks, there was a lot of mumbling
6 in the background and I assume that means it's longer
7 than that.

8 CHAIRPERSON YASZEMSKI: Mr. Christianson,
9 would you care to comment?

10 DR. KIRKPATRICK: That's on the 50 million
11 cycle tests.

12 MR. CHRISTIANSON: Yes. Thank you for
13 asking that question. When I got back, I was told.
14 I meant to say 15 months. It's not 15 weeks.

15 CHAIRPERSON YASZEMSKI: Thank you. So 50
16 million cycles, 15 months.

17 DR. KIRKPATRICK: In the spirit of that,
18 I would like to suggest a post-approval study that
19 takes the wear data out to 50 million cycles as
20 discussed in my presentation. I would also like to
21 ask if they could do a study looking at the other
22 coupled motion, meaning flexion and extension coupled

1 with lateral bending, and provide a rationale for the
2 length of that testing that is reasonable. I don't
3 know.

4 I don't think that needs to be taken to 50
5 million necessarily, because I think what is going to
6 happen is we need to see what happens after somewhere
7 intermediate range, you know, maybe 5 to 7 million and
8 then change directions and see if that makes more wear
9 debris come off or, if you want to, you can do the
10 coupled motion all the time to make sure such as like
11 a figure 8 motion doesn't make a different wear debris
12 pattern. Is that clear enough?

13 CHAIRPERSON YASZEMSKI: I'll --

14 DR. KIRKPATRICK: First, to summarize it,
15 it's, basically, number one is extending data post-
16 approval for 50 million cycles and studying the effect
17 of coupled motion of flexion-extension with lateral
18 bending, as opposed to axial rotation.

19 CHAIRPERSON YASZEMSKI: Okay. I will ask
20 for a second for this motion.

21 DR. MABREY: Second.

22 CHAIRPERSON YASZEMSKI: Dr. Mabrey has

1 seconded. Discussion? Dr. Diaz?

2 DR. DIAZ: I don't have any comment.

3 CHAIRPERSON YASZEMSKI: No comments? Dr.
4 Mabrey?

5 DR. MABREY: No comments.

6 CHAIRPERSON YASZEMSKI: Dr. Finnegan? Dr.
7 Kim?

8 DR. KIM: It seems so excessive to have to
9 test the device for 15 months, so I would question the
10 need to do that, but I do agree with testing the other
11 motions, but for a reasonable period of time.

12 CHAIRPERSON YASZEMSKI: Okay. What would
13 you consider a reasonable period of time?

14 DR. KIM: The 10 million cycles, which
15 represents 80 years if, in fact, that is correct seems
16 very reasonable to me.

17 CHAIRPERSON YASZEMSKI: Thank you. Thank
18 you, Dr. Kim. Dr. Naidu?

19 DR. NAIDU: I concur with Dr. Kim. It
20 appears as if 50 million cycles will be excessive.

21 CHAIRPERSON YASZEMSKI: Thank you. Dr.
22 Blumenstein?

1 DR. BLUMENSTEIN: Let me see if I can
2 understand this. This is not done inside the body,
3 but is done outside the body? It seems excessive to
4 me.

5 CHAIRPERSON YASZEMSKI: Dr. Besser?

6 DR. BESSER: As the testimony we have
7 heard today, they referred to failed prostheses and
8 the failed procedures, none of them in my memory were
9 because of device problems with particulate matter.
10 There were other reasons why it failed. I also would
11 think that the 50 million cycles is probably excessive
12 and the 10 million cycle data that has already been
13 presented is adequate.

14 I would, however, like to see the multiple
15 modes. I'm wondering whether we can separate this
16 motion into two.

17 CHAIRPERSON YASZEMSKI: What we'll do is
18 vote on this and if it passes, they will both occur
19 and if it doesn't pass, we can entertain another
20 motion for one or the other of them. Ms. Maher?

21 MS. MAHER: Yes, I would just like to take
22 this opportunity again. It seems excessive to me. I