



<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation          Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 1 von 7</b>

1020a079e.doc

*Handwritten notes and signatures:*  
 12.11.2001  
 E  
 14  
 20

## Determination of Octopirox in Shampoos by HPLC

### Summary

Octopirox is determined by reversed phase (RP) HPLC after reduction with Ti(III) chloride. The evaluation is carried out via an external standard.

### Contents

1. Purpose
2. Scope of application
3. Definitions
4. Responsibilities
5. Description
  - 5.1 Task and field of application
  - 5.2 Principle
  - 5.3 Apparatus
  - 5.4 Reagents
  - 5.5 Procedure
  - 5.6 Validation of the method
  - 5.7 Remarks
  - 5.8 Appendix
6. Notes
  - 6.1 Further applicable documents
  - 6.2 Literature references
7. Documentation
8. Updating service
9. Distribution

Erstellt	Geprüft	Freigegeben	gültig ab
DI Hirschen	DI Hirschen	DI Hirschen	02.03.2001

<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 2 von 7</b>

*Working copy, only valid 05.03.01* |

**1. Purpose**

The method is used for the determination of the Octopirox content in shampoos in the range from 0.1 to 5 g/100 g.

**2. Scope of application**

Division Functional Chemicals, RQA

**3. Definitions**

Eluent = mobile phase

**4. Responsibilities**

The determination is carried out by a trained employee. The head of the group shall ensure a procedure in accordance with these operating instructions. Any deviations necessary shall be approved by him and recorded by the employee.



<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 4 von 7</b>

*Working copy, only valid 05.03.01* |

## 5.5 Procedure

### 5.5.1 Sample preparation

- weigh approx. 50 mg, accurately to 0.1 mg, of the shampoo to be analyzed into a 25 ml measuring flask
- add 1 ml mobile phase and 100 µl Ti(III) solution
- leave the sample to react for 5 min, make up to volume with the mobile phase
- filter the sample over a membrane filter
- inject 25 µl of the filtered sample
- inject 25 µl of the external standard solution

### 5.5.2 HPLC conditions

Column : Purospher RP-18, 5 µm, 125 x 4.0 mm  
 Eluent : MeOH/H<sub>2</sub>O = 95/5 (v/v), degassed  
 Flow rate : 0.5 ml/min (pressure: approx. 180 bar)  
 Detector : UV detector  
 Wavelength : 300 nm  
 Injection : 25 µl

### 5.5.3 Result of the determination

$$\text{Octopirox content (g/100 g)} = \frac{A_s \cdot E_{St} \cdot 100}{A_{St} \cdot E_s}$$

In this equation:

$A_s$  = peak area of the Octopirox in the sample  
 $A_{St}$  = peak area of the standard  
 $E_s$  = mass of the sample in mg/ml  
 $E_{St}$  = mass of the standard in mg/ml

<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 5 von 7</b>

*Working copy, only valid 05.03.01*

## 5.6 Validation of the Method

### Precision

The precision is checked by injecting an Octopirox standard solution (0.01 mg/ml) six times.

Area <sub>1</sub>	=	458.85 counts	Area (n = 6)	=	458.25 counts
Area <sub>2</sub>	=	458.47 counts	std	=	0.43 counts
Area <sub>3</sub>	=	458.68 counts	VC	=	0.09%
Area <sub>4</sub>	=	457.82 counts			
Area <sub>5</sub>	=	458.28 counts			
Area <sub>6</sub>	=	459.00 counts			

### Linearity

The linearity is checked by duplicate injection of 5 Octopirox standard solutions in the concentration range from 0.001 to 0.06 mg/ml.

x mg/ml	Area <sub>1</sub> counts	Area <sub>2</sub> counts
0.0584	2550.5	2549.9
0.0292	1284.1	1284.7
0.0103	458.8	458.5
0.0021	92.0	92.8
0.0011	46.2	46.5

Standard deviation of the method: 0.00014 mg/ml

Relative standard deviation of the method: 0.68%

Sensitivity:  $43674 \frac{\text{counts}}{\text{mg/ml}}$

<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation          Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 6 von 7</b>



*Working copy, only valid 05.03.01*

### 5.7 Remarks

The procedure is described for Octopirox contents of approx. 0.5 %. For higher Octopirox contents the sample masses are to be adjusted accordingly. If conventional RP-18 columns are used, a slight increase in the retention time of Octopirox may occur in the chromatography of shampoos. In such cases, it should be ensured that the correct signal is evaluated by standard addition of Octopirox to the appropriate samples.

### 5.8 Appendix

typical chromatogram



<b>Clariant GmbH</b>  Division Functional Chemicals RQA	<b>BMH</b>	<b>EQH-1020-AA-0079</b>
	<b>Standard Operation Procedure</b>	<b>Version 02 (02.03.01)</b>
		<b>Seite 7 von 7</b>

*Working copy, only valid 05.03.01*

**6. Notes**

**6.1 Further applicable documents**

Monitoring of balances (EQH-1110-AA-0018)

**6.2 Literature references**

Research report from Hoechst AG no. 06-Z30-0120-75 of 03.02.1975

**7. Documentation**

The head of the group is responsible for the electronic documentation of the Standard Operation Procedure. The documentation of the paper version is delegated to the archive employee.

**8. Updating service**

The head of the group is responsible for the update of this Standard Operation Procedure.

**9. Distribution**

This Standard Operation Procedure is to via the document management system accessible for all involved employees of RQA.