

LABELING

Caution: Federal Law (USA and Canada) restricts this device to sale by or on the order of a physician.

encore
surgical

Important Information for the Operating Surgeon Ceramic/Ceramic Total Hip System

PRODUCT DESCRIPTION AND IMPLANT MATERIALS

The Ceramic / Ceramic acetabular liners (28mm x 48/64mm and 32mm x 52/64mm diameter sizes) and femoral heads (28mm and 32mm/Neutral, +4, -4) are manufactured from aluminum oxide ceramic (alumina - ISO 6474). The metal acetabular shells used with the ceramic liner are manufactured from Ti6Al4V (ASTM F136) and are porous coated (CP titanium beads, ASTM F67) for cementless press-fit applications. The corresponding metal hip stem components used in the study are from Encore's FOUNDATION[®], Revelation[®] and Linear[®] Hip Systems that are porous-coated, semi-constrained total hip replacements intended for cementless press-fit applications.

INDICATIONS

The Ceramic/Ceramic Total Hip System is indicated for primary total hip arthroplasty in patients with the following conditions:

- 1) inflammatory tissue disorders; and,
- 2) noninflammatory degenerative joint disease including osteoarthritis, post-traumatic arthritis or secondary arthritis, and avascular necrosis

CONTRAINDICATIONS

- 1) infection;
- 2) sepsis;
- 3) osteomyelitis;
- 4) rapid joint destruction or bone absorption apparent on roentgenogram;
- 5) skeletally immature patients and cases where there is a loss of abductor musculature, poor bone stock, poor skin coverage around hip joint which would make the procedure unjustifiable;
- 6) uncooperative patient or a patient with neurologic disorders and incapable of following instructions;
- 7) osteoporosis;
- 8) metabolic disorders which may impair bone formation;
- 9) osteomalacia; and
- 10) distant foci of infections (which may cause a hematogenous spread to the implant site);
- 11) obesity; and,
- 12) foreign body sensitivity

WARNINGS

Improper selection, placement, positioning, and fixation of the implant components may result in unusual stress conditions and a subsequent reduction in service life of the prosthetic implant. The surgeon is to be thoroughly familiar with the implant, instruments, and surgical procedure prior to performing surgery.

PRESS FIT APPLICATIONS

Tight fixation at the time of surgery is critical to the success of the procedure. The femoral component stem must press fit into the femur, which necessitates precise operative technique and the use of specified instruments. Intraoperative fracture of the femur can occur during seating of the prosthesis. Bone stock must be adequate to support the device.

ACETABULAR FIXATION SCREWS

Caution should be exercised when using acetabular shell fixation screws to avoid the perforation of the pelvis and possible rupture of blood vessels.

AMIC ACETABULAR SHELL/LINER

- Fixation screws should be fully seated to assure stable fixation of the shell, and to avoid interference with the ceramic liner which could lead to premature failure/fracture of the component.
- Prior to seating the ceramic liner into the shell component, surgical debris must be cleaned from the interior of the shell. Debris may inhibit the liner from locking into the shell component. Failure to properly seat the liner into the shell can lead to disassociation of the liner from the shell and may cause it to bind in the wrong position, chip, or be damaged.
- The ceramic liner should be placed in the shell by hand. Prior to impacting, take your finger and run along the rim to ensure the liner is even. Gently impact with the impactor provided.
- The ceramic liner should not be implanted if the liner is damaged (e.g. if the liner is dropped on the floor, or if the liner is scratched by an instrument) as this can significantly affect the structural integrity of the components.
- Replace both the ceramic liner and metal acetabular shell if the liner is chipped, cracked, or otherwise damaged during the implant procedure or postoperative timeframe. This is because the acetabular shell taper, once deformed through assembly to its mating ceramic liner, should not be reassembled to another ceramic liner.
- Do not reassemble a ceramic liner and metal acetabular shell once they have been disassembled, due to the deformation incurred by the taper locking mechanism during the initial assembly.
- Do not re-use ceramic liners. Even though the implant may appear undamaged, it may have small defects and internal stress patterns which lead to early failure/fracture of the component.

CERAMIC HEAD

- Do not use ceramic heads with components from other manufacturers because design, material, or tolerance differences may lead to premature device failure. Ceramic heads must only be used with hip stem tapers that correspond with the head tapers. Components of the system have been specifically designed to work together.
- The hip stem taper and head taper should be dry and free of all contamination to ensure proper seating and assembly. Failure to properly seat the head onto the stem can lead to disassociation of the head from the stem or may cause it to be damaged.
- The head component should be placed on the stem taper gently while keeping the head and taper in alignment, and then firmly seated by sharply hitting the head using a soft plastic hammer or impactor.
- Ceramic heads should not be implanted if the head or the cone of the stem is damaged (e.g. if the head is dropped on the floor, if the stem taper is scratched by an instrument, or if the head and stem were assembled and disassembled) as this can significantly affect the structural integrity of the component.
- Do not reassemble a ceramic head and metal stem once they have been disassembled, due to the deformation incurred by the taper locking mechanism during the initial assembly.
- If the ceramic head must be revised for any reason and the hip stem is firmly fixed, the revision should be made with a CoCr head and corresponding polyethylene liner and metal shell.
- Do not re-use ceramic heads. Even though the implant may appear undamaged, it may have small defects and internal stress patterns which lead to early failure of the device.

PRECAUTIONS

The patient must be advised of the limitation of the reconstruction and the need for protection of the implant from full weight bearing until adequate fixation and healing have occurred. Excessive activity and trauma affecting the joint replacement have been implicated with premature failure of the reconstruction by loosening, fracture and/or wear of the prosthetic implants. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult. The patient should be advised to report any related pain, decrease in range of motion, swelling, fever, and unusual incidences.

The patient is to be cautioned to govern activities and protect the replaced joint from unreasonable stresses, and follow the instructions of the physician with respect to follow-up care and treatment.

The patient is to be warned of surgical risks, and made aware of possible adverse effects. The patient is to be warned that the device does not replace normal healthy bone, and that the implant can break or become damaged as a result of strenuous activity or trauma.

Specialized instruments are available and must be used to assure the accurate implantation of prosthetic components.

An implant should never be reused. While it may appear undamaged, imperfections may exist which would reduce the service life of the implant.

While rare, intraoperative fracture or breaking of instruments can occur. Instruments which have experienced extensive use or excessive force are susceptible to fracture. Instruments should be examined for wear prior to surgery.

The safety and effectiveness of the Ceramic/Ceramic Total Hip System has not been established in patients with the following conditions, or undergoing the following procedures

- 1) functional deformity,
- 2) revision procedures where other treatments or devices have failed,
- 3) cemented fixation, and,
- 4) treatment of nonunion, femoral neck and trochanteric fractures of the proximal femur with head involvement

VERSE EFFECTS

Please refer to the Adverse Events Table, below, for the operative site adverse events reported with the use of this device system

Table 1. Time Course Distribution of Operative Site Adverse Events for Ceramic/Ceramic Total Hip System vs. Control Ceramic/Polyethylene Hip System Out to 24 Months Post-Operatively.

Operative Site Adverse Events	Ceramic/Ceramic Original Study Population 250 cases (237 patients) enrolled							Continued Access Population 447 cases (400 patients) enrolled							Ceramic/Polyethylene Control Population 250 cases (242 patients) enrolled						
	Op	6	12	24	Tot	%	Op	6	12	24	Tot	%	Op	6	12	24	Tot	%			
N = cases evaluated	250	220	207	196	-	-	447	343	272	115	-	-	250	216	205	187	-	-			
Revision*	0	1	0	0	1	.4%	0	1	0	1	2	.4%	0	5	0	1	6	2.4%			
Fractured Liner	1	0	0	0	1	.4%	0	0	0	0	0	0%	0	0	0	0	0	0%			
Clicking w/Walking	0	0	0	1	1	.4%	0	0	0	0	0	0%	0	0	0	0	0	0%			
Hip Pain	2	2	1	1	6	2.4%	0	1	0	0	1	.2%	0	0	0	3	3	1.2%			
Dislocation	7	0	0	0	7	2.8%	6	2	0	1	9	2%	10	0	0	0	10	4%			
Fracture femur	8	0	0	0	8	3.2%	6	0	0	0	6	1.3%	2	0	0	0	2	8%			
Trochanteric Bursitis	2	1	0	2	5	2%	0	3	1	1	5	1.1%	2	1	5	2	10	4%			
Wound Infection	2	0	2	1	5	2%	1	1	0	0	2	.4%	2	0	2	0	4	1.6%			
DVT	3	1	0	0	4	1.6%	1	0	0	0	1	.2%	4	0	0	0	4	1.6%			
Hematoma	2	0	1	0	3	1.2%	2	0	0	0	2	.4%	1	0	0	0	1	.4%			
Flexor Tendonitis	0	1	0	0	1	.4%	0	1	0	0	1	.2%	0	0	0	0	0	0%			
Acetabular Cell Tumor	1	0	0	0	1	.4%	0	0	0	0	0	0%	0	0	0	0	0	0%			
Heterotopic Ossification	0	1	0	0	1	.4%	0	1	1	0	2	.4%	0	3	0	0	3	1.2%			
Trochanteric Wire Break	1	0	0	0	1	.4%	4	0	0	0	4	.9%	1	0	0	0	3	1.2%			
Acetabular Loosening	0	0	0	0	0	0%	0	0	0	0	0	0%	1	0	0	1	2	.8%			
Leg Length Inequality	0	0	0	0	0	0%	0	0	0	0	0	0%	2	0	0	0	2	8%			
Fracture Pubic Rami	0	0	0	0	0	0%	0	0	0	0	0	0%	0	0	1	0	1	.4%			
Abductor Weakness	0	0	0	0	0	0%	0	0	0	0	0	0%	1	0	0	0	1	.4%			
Misoriented Shell	0	0	0	0	0	0%	1	0	0	0	1	.2%	0	0	0	0	0	0%			
Hip Infection	0	0	0	0	0	0%	2	0	1	0	3	.67%	0	0	0	0	0	0%			
TOTAL	29	7	4	5	45		23	10	3	3	39		28	9	8	7	52				

Op = intraoperatively; 6 = 6 months; 12 = 12 months; 24 = 24 months postoperative; Tot = total. Table includes all operative site adverse events - device related and 'unrelated'. % = total number of a particular adverse event reported divided by number of cases enrolled. * Revisions, included here, were not included among the operative site adverse events reported and analyzed by the sponsor. However, the adverse events that led to the revisions were reported and analyzed. The sponsor listed and analyzed revisions separately under survivorship.

General adverse effects reported for any total hip replacement surgery include:

- 1) Accelerated wear of the articulating surfaces of acetabular components has been reported following total hip replacement. Such wear may be initiated by particles of cement, metal, or other debris which can cause abrasion of the articulating surfaces. Accelerated wear shortens the useful life of the prosthesis, and leads to early revision surgery to replace the worn prosthetic components
- 2) Metallosis and osteolysis may be implicated from wear debris associated with the use of orthopedic implants
- 3) Peripheral neuropathies have been reported following total joint surgery. Subclinical nerve damage occurs more frequently, possibly the result of surgical trauma.
- 4) Metal sensitivity reactions in patients following joint replacement have been rarely reported. Implantation of foreign material in tissues can result in histological reactions involving macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to, or during the healing process. In some cases, wear debris can initiate the process of histiocytic granuloma formation and consequent osteolysis and loosening of the implant
- 5) Dislocation and subluxation of implant components can result from improper positioning of the components. Muscle and fibrous tissue laxity can also contribute to these conditions
- 6) Implants can loosen or migrate due to trauma or loss of fixation
- 7) Infection can lead to failure of the joint replacement
- 8) While rare, fatigue fracture of the implant can occur as a result of strenuous activity, improper alignment, or duration of service
- 9) Fracture of the femur can occur while press-fitting (seating) the femoral stem component into the prepared femoral canal
- 10) Allergic reactions

Intraoperative and early postoperative complications can include

- 1) femoral or acetabular perforation, or fracture,
- 2) femoral fracture can occur while seating the device,
- 3) damage to blood vessels,
- 4) temporary or permanent nerve damage resulting in pain or numbness of the affected limb,
- 5) undesirable shortening or lengthening of the limb,
- 6) traumatic arthrosis of the knee from intraoperative positioning of the extremity,
- 7) cardiovascular disorders including venous thrombosis, pulmonary embolism, or myocardial infarction,
- 8) gastrointestinal complications,
- 9) genitourinary complications,
- 10) hematoma,
- 11) delayed wound healing,
- 12) infection,
- 13) pain, and,
- 14) death

Late postoperative complications can include

- 1) trochanteric avulsion as a result of excess muscular weakening,
- 2) trochanteric non-union due to inadequate reattachment and/or early weight bearing,
- 3) aggravated problems of the knee or ankle of the affected limb or contralateral extremity by leg length discrepancy, too much femoral medialization, or muscle deficiency,
- 4) femoral fracture by trauma or excessive loading, particularly in the presence of poor bone stock,
- 5) periparticular calcification or ossification, with or without impediment to joint mobility,
- 6) inadequate range of motion due to improper selection or positioning of components, by femoral impingement, and periparticular calcification,
- 7) localized progressive bone resorption (osteolysis)

Study Design

The study was a prospective, multi-center, concurrently controlled clinical trial. A total of 500 cases (479 patients) were enrolled and randomized into 1 of 2 treatment groups. The investigational study group received the ceramic/ceramic hip, and the control group received a standard ceramic/polyethylene hip. Patients were followed until the last patient implanted attained their 2-year evaluation. Study patients consisted of those who had inflammatory tissue disorders (e.g. rheumatoid arthritis, lupus, etc), osteoarthritis, post-traumatic arthritis (secondary arthritis), or avascular necrosis, and less than 70 preoperative Harris Hip Score; and, were candidates for primary total hip replacement. After fulfillment of the original study enrollment, an additional 447 ceramic/ceramic devices were implanted in 400 patients under 'Continued Access'. Safety data for this group is reported in the form of adverse events in Table 1.

Clinical Patient Assessment

Each patient was evaluated at 6, 12, and 24 month post-operative time intervals. At each visit a Harris Hip Score and A/P, M/L radiographs were obtained. Radiographs were evaluated by an independent radiologist. The study endpoints consisted of a comparison of mean HHS scores, complication rates, survival rates, and radiographic failure rates (based on lucencies/migration). Success/failure of the study (i.e., of the investigational ceramic/ceramic device) was based on comparison of the results of the study endpoints for the two treatment groups. Acceptance criteria for success of the ceramic/ceramic hip system required that it would perform within 5 points or 5 percentage points of the control hip system for each of the study endpoints.

Study Population/Demographics

A total of 500 procedures were performed on 479 patients by 17 surgeons at 17 sites in the original clinical trial. A total of 250 investigational ceramic/ceramic devices were implanted into 237 patients in the original clinical trial. Over the same time period, 250 control devices were implanted in 242 patients. Forty-two patients (42) received bilateral implants, 12 received investigational devices in both hips, 9 received control devices in both hips, and 21 patients received 1 investigational and 1 control device. The primary analysis included 5 investigational patients and 8 control patients who were protocol deviations. These patients were analyzed along with the group to which they were originally randomized. The 42 bilateral patients (84 cases) were not included in the efficacy analysis for comparison of mean HHS scores. The bilateral HHS data was reported and analyzed separately, however. The bilateral patients/cases were included in the overall analysis of adverse events, device survival, and radiographic failure rates.

Table 2: Demographics

Category	*Demographics					
	Female		Male		Total	
	Study	Control	Study	Control	Study	Control
Number of cases	112	133	138	117	250	250
Mean Age	55.8	61.9	54.3	59.8	55.0	60.9
Age Range	17-91	19-86	18-87	27-94	17-91	19-94
Mean Preop HHS	42.9	42.7	46.4	46.5	44.8	44.5
HHS Range	8-65	10-68	12-69	12-69	8-69	10-69
Right Hip	53	77	67	58	120	135
Left Hip	59	56	71	59	130	115

*includes bilateral cases, Study = Ceramic/Ceramic, Control = Ceramic/Polyethylene

Device Accountability

The course distribution of device accountability at each follow-up interval out to 2 years is provided in Table 3.

Table 3: Device Accountability

	Ceramic/Ceramic Study Population (n=250)			Ceramic/PE Control Population (n=250)		
	6 mo	1 yr	>2 yrs*	6 mo	1 yr	>2 yrs*
Theoretical Follow-up	250	250	250	250	250	250
ΔDeaths	0	1	6	0	1	3
ΔRevisions	1	0	1	5	0	6
Expected Follow-up	249	248	243	245	244	241
Lost To Follow-Up	0	0	14	0	0	26
Actual Follow-Up %	88.3% 200/249	83.5% 207/248	94.2% (229/243)	87.8% 216/246	84% 205/244	89.2% (215/241)

Theoretical Follow-Up (TFU) = theoretical number of cases available for follow-up; Expected Follow-Up (EFU) = TFU minus deaths and revisions; Actual Follow-Up (AFU) = cases that had clinical data available at specified follow-up interval; Lost to Follow-Up (LTFU) = EFU minus AFU; AFU (%) = AFU / EFU; * includes 2 year or next available evaluation (3 or 4 year); n = number of cases; Δ cumulative from previous follow-up interval; table includes bilateral cases and protocol deviations

Safety and Effectiveness Data

Clinical testing was performed to determine if the ceramic/ceramic device was as safe and effective as the control device. Effectiveness was assessed using the Harris Hip Score. Safety was assessed according to complications, survival of the devices, and radiographic analyses. Predetermined acceptance criteria were established for each endpoint. The two treatment groups: investigational (ceramic liners) and control (polyethylene liners) were compared based on the acceptance criteria that the ceramic/ceramic hip system would perform within 5 points or 5 percentage points for each endpoint (6 percentage points for complications).

The results of the study, presented in Table 4, show that the study success criteria established to demonstrate the safety and effectiveness of the ceramic/ceramic hip system as compared to the control total hip system, were met. There were no increased risks associated with the ceramic insert. The study showed both treatment groups to be equivalent. All scores and percentages reflect the outcome at the 2-year evaluation or if a particular patient missed their 2-year evaluation, the next annual evaluation was used. There were no radiographic failures; therefore, the percentages are not indicated in the following table.

Table 4: Safety and Effectiveness Data

	Study	Control
Mean HHS Scores*	92.36	92.16
% of Patients with Related Complications	6.4%	4.8%
% of Patients with Unrelated Complications	9.2%	14.1%
Survival**	99.6%	97.2%
Radiographic Failures#	0	0

Does not include continued access data, only data from the original study population. * Includes unilateral patients only (i.e. 186 study and 183 control), all other categories include bilateral patients/devices. ** lack of revision or removal, #lack of specified radiolucencies and migration

STERILIZATION

Unless opened or damaged, components are supplied sterile in triple peel pouches or double blister trays. Check all packaging for punctures or other damage. Sterilization of components is by gamma radiation to achieve a Sterility Assurance Level (SAL) of 10⁻⁶. When removing the implant, appropriate aseptic procedures must be observed. The outer barrier and as appropriate, the middle barrier may be opened by non-aseptic personnel. The inner barrier must be retrieved by aseptic personnel and the implant must not come in contact with any objects that might damage the surface.

Warning: Do not autoclave any ceramic components or ceramic heads.

<p><u>in Resterilization</u> - -FOR NON-METALLIC COMPONENTS ONLY! Remove from supplied packaging and wrap in protective wrap</p>	<p>270° F. (132°C.) for 30 minutes Gravity Displacement Autoclave Sterility Assurance Level (SAL) of 10⁻⁶</p>
--	--

WARNING: Protect all machined tapers, coatings and polished surfaces. Standard cleaning procedures cannot be relied upon to remove contamination from plasma coating.

Key to the symbols on the labeling:

	Single use - do not reuse	STERILE R	Sterility symbol R:gamma radiation min.25 kGy
	Expiration Date	non-sterile	Non-sterile symbol :nonsterile
LOT	Lot number	Qty.	Quantity of items in package
			See "Instructions for Use"

For further information regarding the use of the ENCORE® Ceramic / Ceramic Acetabular components, contact your ENCORE® representative or distributor.
 ENCORE MEDICAL, L.P., 9800 Metric Blvd., Austin, TX 78758 USA (Made in the USA)
 FOUNDATION®, LINEAR®, KEYSTONE®, are registered trademarks of Encore Medical Corporation, Austin, TX 78758 USA or its affiliates.

0400-0082


PATIENT LABELING FOR THE CERAMIC/CERAMIC TOTAL HIP SYSTEM

Table of Contents

What is the Ceramic/Ceramic Total Hip System	1
What is the purpose of the device?	1
What are the definitions for the above-mentioned conditions?	2
What are the contraindications of this device?	2
What are the benefits and risks of using this device?	3
What procedure will the surgeon follow?	4
What will I have to do before my surgery?	4
What are the expectations after my surgery?	4
Are there any other devices and procedures that I could choose?.....	5
GLOSSARY OF TERMS	6
Printing Date.....	7

What is the Ceramic/Ceramic Total Hip System?

The Ceramic/Ceramic Total Hip System is an artificial hip joint made up of the following parts: a ceramic (alumina) ball, a ceramic (alumina) liner that fits into a metal shell; and a compatible metal hip stem that is used with these parts.



What is the purpose of the device?

The Ceramic/Ceramic Hip System is indicated for use in total hip joint replacement for the reduction or relief of pain and/or improvement in function. The device may be used in the treatment of the following conditions: inflammatory tissue disorders and non-inflammatory degenerative joint disease, including osteoarthritis, post-traumatic arthritis (secondary arthritis), and avascular necrosis.

What are the definitions for the above-mentioned conditions?

Inflammatory tissue disorders - These disorders (like rheumatoid arthritis) cause the lining covering the joints, referred to as the synovium, to become thickened and inflamed. These conditions allow too much synovial fluid within the joint space causing a chronic inflammation that damages the cartilage. This damage results in cartilage loss, pain and stiffness.

Non-Inflammatory degenerative joint disease - Is a condition that is caused by the breakdown and eventual loss of the cartilage of one or more joints. Cartilage is a protein substance that serves as a "cushion" between the bones of the joints. Non-Inflammatory degenerative joint disease includes the following three conditions listed below:

Osteoarthritis – This is often referred to as degenerative arthritis because it is a “wearing out” condition involving the breakdown of cartilage and bones. When the cartilage breaks down, the bones rub against each other, resulting in pain and stiffness.

Post-traumatic arthritis (secondary arthritis) – This type of arthritis causes inflammation of a joint resulting from an injury. The symptoms of this condition include: breakdown of the bone and rubbery tissue, bleeding in the joint space, increased thickness of the bone, flattening of the joint surface, cartilage separation from the underlying bone, and wearing away (erosion) of the bone.

Avascular Necrosis – This condition results from a loss of blood supply to the hip bones, which can lead to necrosis in portions of these bones. Without a proper supply of blood the bone's structure weakens and may collapse damaging the cartilage.

What are the contraindications of this device?

Reasons for not using this device may include:

1. An infection exists in the hip,
2. There are remote centers of infection that may be spread through the blood supply to the actual implant site,
3. Rapid joint destruction or bone loss is seen on multiple X-rays,
4. The patient's bones are still growing,
5. Patients muscles may be too weak to work properly (e.g., loss of function, fusion, or poor skin coverage around the hip joint),
6. The patient has weak bones from various causes (e.g., metabolic disorders such as Osteoporosis, Paget's disease, and Osteomalacia), which may prevent adequate support for the device.
7. The patient is unable or unwilling to follow proper procedures after joint replacement,
8. The patient weighs too much for implanted device to support the patient's full body weight, or
9. The patient is allergic to the implant materials.

What are the benefits and risks of using this device?

The potential benefits of using the Ceramic/Ceramic Total Hip replacement are the relief of pain and the improved function of your hip joint.

A two year clinical study to evaluate the Ceramic/Ceramic Total Hip System found it to be comparable, in terms of safety and effectiveness, to a marketed ceramic and plastic hip system. The study followed 479 patients who met the indications for use identified above. Two hundred and thirty seven (237) patients were implanted with the Ceramic/Ceramic Hip, and 242 were implanted with the ceramic/plastic hip. Evaluation of safety and effectiveness included comparisons of complications, x-ray results, survivorship (i.e., was the device still in place), function, and pain relief. There were complications reported but they were consistent with a marketed total hip replacement system. All patients had successful x-ray results, and 99.6% of the patients still had their device in place two years after surgery.

There is a risk that this device may not restore the hip to a condition of a normal healthy person. Activities after implanting the artificial hip may be somewhat limited.

With any surgery there are risks and complications that can arise during and after surgery. For example there may be problems with the anesthesia, development of an infection, hematoma, or hemorrhage.

Possible risks and complications associated with any total hip surgery, listed in order of severity, may include:

- death,
- myocardial infarction,
- the development of blood clots that may dislodge and cause decreased blood flow to the lungs (pulmonary embolism),
- deep vein thrombosis,
- allergic reactions,
- nerve injury, (e.g., numbness or pain),
- not being able to urinate due to pain and/or inflammation,
- swelling that continues,
- fracture of the implant,
- fracture of the bone,
- revision of the implant,
- non-union of implant and bone,
- changing position of the implant,

- dislocation of the implant,
- unequal leg lengths,
- muscle weakness
- aggravated problems in the joints of the affected limb, opposite limb, or back due to unequal leg lengths or muscle weakness,
- decreased range of motion,
- local bone resorption,
- heterotopic bone formation,
- bursitis,
- tendonitis,
- cellulitis,
- pain

Use of the ceramic liner includes the risks of the implant changing position, or the liner separating from the metal shell portion of the artificial socket. Also, there may be excessive wearing of the ceramic material, causing joint loosening and impairment of function.

The above list of risks may not be complete, and any or all risks are possible. If complications occur, you may require further hip surgery.

What procedure will the surgeon follow?

Your doctor has determined that you might benefit from a surgical procedure that replaces the hip joint with artificial parts. In this procedure, your surgeon will remove the ball portion (femoral head) of your hip and replace it with a new ceramic ball. The ball part of your new hip is attached to a long metal stem that allows your surgeon to place it firmly in the thigh bone. Your surgeon will replace the socket part (acetabulum) of your hip with a new one that fits with the new ceramic ball. The new socket is made up of the ceramic liner and metal shell.

What will I have to do before my surgery?

You will have to have routine pre-surgical tests and procedures before your replacement hip surgery, including blood tests and X-rays. Your surgeon will perform your hip surgery and take care of you throughout your surgery and healing period after surgery. The surgical procedure for this hip replacement is similar to any other hip replacement.

What are the expectations after my surgery?

You will experience discomfort following the hip surgery and the necessary rehabilitation. There are limits to what you can do after you receive your new hip replacement. Following

the surgery, you may need to use a walker, crutches, or a cane for up to three (3) months, depending upon the speed of your recovery. Your surgeon will tell you when and how much weight you should place on your new hip. After your operation, you will need physical therapy. It will not differ from other hip replacement rehabilitation. Please follow your surgeon's post surgery directions carefully. This will help speed your proper recovery during the rehabilitation period. Your surgeon will schedule the number of postoperative visits based on your ability to heal. At each visit your surgeon will examine you and ask for X-rays of your implant.

After healing has occurred, you may not be able to do the activities you did prior to your surgery (such as playing sports or heavy physical work). The device does not replace normal healthy bone. Your new hip can break or become damaged as a result of strenuous activity, trauma, or even normal use. An artificial hip has a limited expected service life, and may need to be replaced at some time in the future. Unexpected injury to the hip can cause broken bones, implant loosening, or wear of the implants.

You should contact your doctor immediately if you experience any of the following:

- excessive pain (pain that limits and or prevents normal activity),
- decrease in range of motion (post surgical range of motion will be limited; however, range of motion should increase as rehabilitation progresses. If during or after rehabilitation range of motion decreases, your doctor should be contacted.),
- swelling, or
- fever.

Are there any other devices and procedures that I could choose?

You could elect not to have surgery and choose a non-surgical method of treatment of reduced activity and medicine. You could elect to have a total hip replacement with another commercially available implant or implant part. Choices available to you include a metal femoral head and a plastic liner (metal/plastic), a metal femoral head and a metal liner (metal/metal), or a ceramic femoral head and plastic liner (ceramic/plastic). Other alternative operations may be available to you, including fusion of the hip joint. Please talk to your orthopedic doctor about these and other options.

GLOSSARY OF TERMS

Bursitis – inflammation of the bursa, a pad like sac or cavity found in the joint area and acts as a lubricant to the joints.

Cartilage – rubbery type of tissue that pads the joints

Cellulitis – inflammation of the skin around the hip due to infection

Deep vein thrombosis – the development of a blood clot in the vascular system

Dislocation - the ball at the top of the thighbone (the femoral head) is not stable within the ceramic insert

Fusion – joining together

Hematoma – collection of blood under the skin

Hemorrhage – bleeding

Heterotopic bone – formation of bone in areas in which bone is not intended to form
Inflammation – swelling, redness, and pain in tissues

Local bone resorption – the deterioration of bone affected by the surgery

Myocardial infarction – changes that occur to the heart muscle due to the lack of circulating blood (i.e., heart attack)

Necrosis – cell death

Non-union – the bone does not grow into the artificial device possibly leading to loosening

Osteomalacia - Softening of bone, particularly in the sense of bone weakened by demineralization (the loss of mineral) and most notably by the depletion of calcium from bone

Osteoporosis - Thinning of the bones with reduction in bone mass due to depletion of calcium and bone protein

Paget's disease - A chronic bone disorder that typically results in enlarged, deformed bones due to excessive breakdown and formation of bone tissue that can cause bones to weaken and may result in bone pain, arthritis, deformities or fractures

Revision – removal of the man-made devices

Printing Date