

Summary Basis for Regulatory Action

October 21, 2016

From: Ricardo Espinola, Chair of the Review Committee

BLA/ STN#: Anti-Human Globulin (Rabbit/Murine Monoclonal) (Formulated for Automated Testing) STN 125529/0

Applicant Name: Bio-Rad Medical Diagnostics GmbH License No. 1845

Date of Submission: February 27, 2014

MDUFA Goal Date: October 21, 2016

Proprietary Name: IH-Card AHG Anti-IgG,-C3d

Established Name (common or usual name): Anti-Human Globulin (Rabbit/Murine Monoclonal) (Formulated for Automated Testing).

Intended Use: The Anti-Human Globulin reagent is intended to be used as a component to manufacture the IH gel card products. - The final in-vitro products (IH-Cards) that contain Anti-IgG,-C3d are intended to be used on the IH-1000 analyzer for the detection and identification of clinical relevant antibodies, cross matching, and Direct Antiglobulin Testing (DAT), based on the principles of agglutination and gel filtration.

Recommended Action: Approval

Signatory Authorities Action:

Offices Signatory Authority: Jay Epstein, MD, Director, Office of Blood Research and Review

- I concur with the summary review.
- I concur with the summary review and include a separate review to add further analysis.
- I do not concur with the summary review and include a separate review.

Offices Signatory Authority: Mary Malarkey, Director, Office of Compliance and Biologics Quality

- I concur with the summary review.
- I concur with the summary review and include a separate review to add further analysis.
- I do not concur with the summary review and include a separate review.

Material Reviewed/ Consulted - Specific documentation used in developing the SBRA	
	Reviewer Name – Document(s) Date
Clinical Review	Joyce Rockwell November 10, 2015 Ricardo Espinola August 08, 2016
Statistical Review	Zhen Jiang, PhD November 26, 2014 July 14, 2015 October 30, 2015 August 30, 2016
CMC Facility Review	Chad Burger

	December 17, 2014 November 13, 2015
CMC Product Review	Joyce Rockwell June 23, 2014 August 21, 2014 Ricardo Espinola August 08, 2016
Bioresearch Monitoring Review	Bioresearch monitoring inspections were not conducted for this BLA and Efficacy Supplement
Labeling	Dana Jones August 26, 2014
Facility Review	Chad Burger December 17, 2014
Establishment Inspection Report	Chad Burger October 27, 2015
Lot Release Protocols/Testing Plans	Karen Campbell (DBSQC – OCBQ) October 28, 2015 (Laboratory Quality Product Testing Plan) November 17, 2015 (Lot Release Protocol Templates)
Advisory Committee Transcript	Not applicable
Other (list)	Not applicable

1. Introduction

Bio-Rad Medical Diagnostics GmbH (BMD), located in Dreieich, Germany (Establishment Registration Number 9610824) submitted to the FDA, 17 applications to obtain approval for an automated immunohematology test system called the IH-System. The submissions consist of:

- Three Biologics License Applications (BLAs): one Anti-Human Globulin and two Blood Grouping Reagents (BGRs).

- Ten Efficacy Supplements: one Anti-Human Globulin, eight BGRs, and one Reagent Red Blood Cells (RRBCs).
- Four 510(k) premarket notifications for the analyzer, software, control and neutral card.

The following is a list of all submissions associated with the IH-System:

- BMD - BLAs and Efficacy Supplements:
 - Anti-Human Globulin (Rabbit/Murine Monoclonal)(Formulated for Automated Testing), STN 125529/0
 - Anti-Human Globulin (Formulated for Automated Testing), STN 125098/88
 - Blood Grouping Reagent, Anti-B (Murine Monoclonal)(IgG)(Formulated for Automated Testing), STN 125532/0
 - Blood Grouping Reagent, Anti-D (Monoclonal Blend)(Formulated for Automated Testing), STN 125533/0
 - Blood Grouping Reagent, Anti-A (Murine Monoclonal)(Formulated for Automated Testing), STN 125094/113
 - Blood Grouping Reagent, Anti-A,B (Murine Monoclonal)(Formulated for Automated Testing), STN 125096/58
 - Blood Grouping Reagent, Anti-D (Monoclonal)(IgM)(Formulated for Automated Testing), STN 125097/67
 - Blood Grouping Reagent, Anti-E (Monoclonal)(Formulated for Automated Testing), STN 125202/50
 - Blood Grouping Reagent, Anti-e (Monoclonal)(Formulated for Automated Testing), STN 125203/48
 - Blood Grouping Reagent, Anti-K (Monoclonal)(Formulated for Automated Testing), STN 125204/46
 - Blood Grouping Reagent, Anti-c (Monoclonal)(Formulated for Automated Testing), STN 125205/46
 - Blood Grouping Reagent, Anti-C (Monoclonal)(Formulated for Automated Testing), STN 125206/48
 - Reagent Red Blood Cells For Use in Automated Systems, STN 125208/70

- BMD - Companion 510(k) submissions:
 - BK140106 IH-1000 Analyzer System
 - BK140107 IH-COM (data management software)
 - BK140138 IH-Card Neutral
 - BK140139 IH-Card Control

- (b) (4) [Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]

The above submissions were grouped as follows: one group containing Anti-Human Globulin reagents, one group containing ten Blood Grouping Reagents, one group containing eight Reagent Red Blood Cells, and one group containing four 510(k) premarket notifications.

- (b) (4) [Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]

Millipore was approved by the FDA on August 14, 2008 to supplement their license to supply the following FFMUs to Biotest Medical Diagnostic GmbH, U.S. License No. 1798 under a shared manufacturing arrangement: Anti-C (clone MS-24), STN 103858/5046, Anti-c (Clone MS-33), STN 103860/5048, Anti-K (clone MS-56), STN 103864/5045 and Anti-e (clones MS-16/MS-21/MS-63), STN103866/5055/5056/5057. BMD U.S. License No. 1845 subsequently acquired portions of Biotest and acquired ownership of these FFMU products.

The IH-System performs ABO grouping, red blood cell antigen typing, detection and identification of clinically significant red blood cell antibodies, crossmatching, and direct antiglobulin testing, based on the principles of agglutination and gel filtration.

It generates results from individual images that must be verified by visual inspection by a qualified operator.

The IH-System consists of:

- IH-Card: a plastic card, consisting of six microtubes containing the active component, i.e., BGR or AHG, in a buffered (b) (4) gel suspension.
- IH-Anti-D Blend: vialled Anti-D reagent for performing weak D and DVI testing using the IH-AHG Anti-IgG card.
- IH-Cell products: vialled RRBCs (i.e., reverse grouping cells, screening cells, pool cells, and identification panel cells).
- IH-1000 Automated Analyzer System: an automated, blood grouping and antibody test system analyzer for the IH-Cards.
- IH-COM: stand-alone software to be used for data management, and the evaluation and interpretation of assay results. The software is directly linked to the IH-1000 via a bidirectional interface and can also be interfaced with the customer's Laboratory Information System (LIS).
- IH-Card Neutral: a plastic card, consisting of six microtubes. The microtubes contain a neutral gel suspension and do not contain AHG reagent. The IH-Card Neutral is designed for serum grouping with A₂ RRBCs and identification of unexpected antibodies with papainized red cells.
- IH-Card Control: a plastic card, consisting of six microtubes filled with (b) (4) (b) (4) containing buffer, diluent medium, and preservative, and is intended for use as a supplemental control for IH-Cards with monoclonal Blood Grouping Reagent without a control well.
- IH-LISS Rack (Class II Exempt from premarket notification procedures): consists of 10 plastic cards, each with six microtubes, filled with a suspending medium, i.e., modified Low Ionic Strength Solution. The IH-LISS is used for preparing red blood cell suspensions for use with the appropriate IH-Card.

The IH-System is not a first of its kind device. Other manufacturers have been approved and cleared to market manual and automated immunohematology test systems using the column agglutination technique first described by Yves Lapierre in 1985 for the detection of red blood cell agglutination.

2. Background

Meetings with FDA

FDA held a pre-submission (CRMTS # 8105, PTS PS001492) meeting with BMD on October 6, 2011. The discussion items included performance studies design, statistical analysis and data reporting, instrument changes, and submission strategy. The meeting package indicated that the future submissions would include both manual and automated testing methods and instrumentation. Prior to submitting the respective submissions to CBER in February 2014, BMD decided to only submit information and data for automated testing using the IH-1000 Automated Analyzer System.

Marketing History

The design of the IH-Cards is based on the technology transfer from two commercially distributed products manufactured by DiaMed Ltd. and Bio-Rad Laboratories, Inc. The plastic card is used by DiaMed Ltd. (Morat, Switzerland) for the DiaMed ID-Micro Typing System that was introduced to non-US markets in 1988 and is still manufactured and marketed by Bio-Rad in Switzerland to non-US markets. The gel was used by Bio-Rad Laboratories, Inc. in (b) (4) for manufacturing the ScanGel® Cards distributed to non-US markets from the late 1990's to 2014.

In 2005, FDA approved the same Anti-IgG (Rabbit) in vitro substance intended to be used as the active component for the IH-Card for the manufacture of Anti-IgG (Rabbit) (Formulated for Automated Testing) used with BMD's 510(k) cleared TANGO® instrument. However, for the IH-Card AHG Anti-IgG,-C3d, the formulation of the in vitro substance is different from what was approved in 2005.

The Anti-C3d component that is used as the active component for the IH-Card is produced from cell culture [REDACTED] 053A-714.

Device Description and Function

The IH-System is an Immunohematology Test System that consists of an analyzer, software (IH-COM), Anti-Human Globulin (IH-Cards AHG), and supplemental reagents (FDA licensed or cleared) for automated testing. The test principle is based on gel filtration and column agglutination. In gel filtration technique, the gel in the microtube acts as a sieve; after centrifugation of the card, non- agglutinated RBCs settle at the bottom of the microtube while the agglutinated cells are dispersed throughout the gel depending on their size.

The IH-Card is a plastic card composed of six microtubes. Each microtube has an incubation chamber at the top of a long and narrow microtube. Each microtube is filled with a mixture of [REDACTED] gel, buffer, and a specific antibody i.e., anti-IgG, -C3d; polyspecific. The agglutination occurs when the red blood cells sensitized *in vivo* or *in vitro* by human IgG antibodies come in contact with the Anti-IgG, -C3d present in the gel solution. The gel column acts as a sieve that traps agglutinated cells as they pass through the column during the centrifugation of the card. The gel separates agglutinated red blood cells based on size. Non-agglutinated red blood cells form a pellet at the bottom of the microtube. The Anti-IgG,-C3d is contained in two IH-Cards.

Table 1: IH-Cards containing Anti-IgG,-C3d

Card Name	Microtube Contents					
	1	2	3	4	5	6
IH-Card AHG Anti-IgG,-C3d	Anti-IgG,- C3d; polyspecific	Anti-IgG,- C3d; polyspecific	Anti-IgG,- C3d; polyspecific	Anti-IgG,- C3d; polyspecific	Anti-IgG,- C3d; polyspecific	Anti-IgG,- C3d; polyspecific
IH-Card	Anti-A	Anti-B	Anti-A,B	Anti-D	Control	Anti-IgG,-

ABO/Rh(DVI+)				(DVI+)		C3d; polyspecific
--------------	--	--	--	--------	--	----------------------

The information provided in the table above was extracted from the submission.

The IH-1000 is a fully automated high throughput analyzer for gel card technology, which allows performance of different assays. The gel cards, reagents and samples are automatically identified by the barcode reader after being placed on the analyzer. Sample pipetting, reagent pipetting and incubation of reaction, if applicable, are all performed automatically without interaction from the operator. Reactions in the gel microtube are captured by the camera and analyzed by the image evaluation software for grading. The evaluated images are transferred to the IH-COM external data management software for further interpretation and generation of results. Every individual result that is generated from the instrument is reviewed; validated and equivocal results are edited (changed) as needed.

The various assays performed by the IH-System provide test results for blood collection establishments, transfusion services, and hospitals for managing donors and patients.

Chronology

CBER received the Anti-Human Globulin BLA on March 7, 2014. CBER issued a Filing with No Deficiencies Letter on May 1, 2014. CBER subsequently received twenty amendments from BMD in response to fifteen information requests for STN 125529/0. A Complete Response Letter was issued on December 31, 2014. A final amendment dated June 30, 2016 completed BMD’s responses to all outstanding issues associated with the Anti-Human Globulin BLA.

3. Chemistry Manufacturing and Controls (CMC)

All manufacturing is carried out in a controlled environment. The applications were submitted in accordance with the recommendations in FDA’s Guidance for Industry: *“Content and Format of Chemistry, Manufacturing, and Controls Information and Establishment Description Information for a Biological in-Vitro Diagnostic Product”*.

a. Manufacturing Summary

Manufacturing of the IH-Card AHG Anti-IgG,-C3d consists of four main manufacturing stages: production of the Anti-IgG and Anti-C3d components, the (b) (4) [redacted], preparation of the gel (b) (4) [redacted], and filling/packaging of the final product. The manufacturing process for the IH- Cards uses the same methodologies, no matter what the antibody type (IgG or IgM), origin (monoclonal or polyclonal), or specificity (i.e., Anti-IgG or Anti-C3D) of the in vitro substances used. The set specifications for the in-process controls and quality control testing demonstrated product homogeneity, reproducibility, and consistency of the manufacturing process.

Manufacture of the Anti-C3d in vitro substance

The manufacture of the Anti-C3d in vitro substance occurs at Bio-Rad Laboratories in (b) (4) [redacted]. The (b) (4) [redacted] facility has been registered with FDA as a medical device manufacturer since 2009.

The Anti-C3d in vitro substance is derived from the (b) (4) [redacted] cell line 053A-714. (b) (4) [redacted]
[redacted]
[redacted]

(b) (4) [redacted]
[redacted]
[redacted]
[redacted]
[redacted]
[redacted]
[redacted]
[redacted]
[redacted].

Raw materials

The fetal calf serum is cell culture grade, made from raw material collected in the United States, and suitable for the manufacture of biopharmaceutical products (92/118/EEC). Prior to purchase, potential lots of fetal calf serum are evaluated by (b) (4). The desired lot is purchased and accepted based on incoming goods inspection and testing, and information listed on the Certificate of Analysis.

The Anti-C3d in vitro substance manufacturing process has been successfully validated with three lots.

Manufacture of the Anti-IgG in-vitro substance

Polyclonal antibody production consists of (b) (4).

The Anti-IgG component has been used for other BMD Anti-Human Globulin products since 1990. The reactivity of the Anti-IgG is not heavy chain specific, but is directed towards the immunoglobulin light chains of IgG, and thus may also bind to IgM or IgA sensitized red blood cells. There is no reactivity with complement coated red blood cells.

Further processing of the Anti-IgG in vitro substance

The rabbit serum received from (b) (4) is further processed by BMD, and includes the addition of 0.1% sodium azide (preservative) and (b) (4).

(b) (4) of the in vitro substance are tested and found negative for anti-HIV1 and 2, anti-HCV, anti-HBV, HBsAg, and syphilis. The (b) (4) with (b) (4).

(b) (4)

[Redacted]

Production of the gel (b) (4)

The gel (b) (4) is manufactured in a (b) (4)

[Redacted]

Date of Manufacture (DOM) / Expiration Date

The date of manufacture (DOM) is the date when the (b) (4) the cards. The expiration date of the filled IH-Card AHG Anti-IgG,-C3d is 16 months from the DOM.

Filling of the IH-Card

In preparation for filling the IH-Card, the required number of (b) (4)

[Redacted]

The filling is performed using an (b) (4)

[Redacted]

(b) (4)
[Redacted]
[Redacted]
[Redacted]

Labeling and Packaging of IH-Cards

The filling line (b) (4)
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]
[Redacted]

The in-process controls, physical inspections, and testing are performed to ensure all specifications are met throughout the filling, labeling, and packaging operation.

Quality Control Testing

Final serological testing is performed on the filled IH-Cards at least (b) (4)
[Redacted]. This validated hold time allows for proper sedimentation of the gel particles and supernatant, resulting in conditions necessary for performing serological tests. Anti-IgG potency testing uses (b) (4)
[Redacted]
[Redacted]
[Redacted]
[Redacted]

The indirect antiglobulin test and crossmatching are performed using the IH-1000 Analyzer. The following table summarizes the release testing and acceptance criteria for the IH-Card AHG Anti-IgG,-C3d:

Table 2: Release testing of IH-Card AHG Anti-IgG,-C3d

Release Testing: IH-Card AHG Anti-IgG,-C3d			
Positive Specificity	Indirect Antiglobulin Test	(b) (4)	(b) (4)
	Direct Antiglobulin Test	(b) (4)	(b) (4)
	Crossmatch	(b) (4)	(b) (4)
Negative Specificity	Specificity against antigen negative cells	(b) (4)	(b) (4)
	Indirect Antiglobulin Test	(b) (4)	(b) (4)
	Direct Antiglobulin Test	(b) (4)	(b) (4)
	Crossmatch	(b) (4)	(b) (4)

	(b) (4)	
Potency	(b) (4)	(b) (4)
Visual Inspection	(b) (4)	(b) (4)
Bioburden	(b) (4)	(b) (4)

**s means: numerous small clumps.*

The information provided in the table above was extracted from the submission.

Quality control of packaging and labeling is performed after mechanical labeling and packaging occurs. The number of samples taken for inspection is based on a statistical sampling plan. The inspection includes label identity, label position, presence of gel and supernatant, check for air bubbles in the gel, gel level check, sealing area, check for completeness of packaged carton and for the correct version of the Instructions for Use.

Microbiology/Bioburden

The AHG Anti-IgG,-C3d reagent is microbiologically controlled, and as such is not considered labeled as sterile. The manufacturing of the in-vitro products includes the addition of sodium azide as preservative and filtration using a (b) (4) filter. The bioburden test method was qualified in accordance with (b) (4) The proposed sodium azide formulation concentration was shown to have effective anti-microbial properties in accordance with (b) (4)

b. CBER Lot Release

The lot release protocol template was submitted to CBER for review and found to be

acceptable after revisions. A lot release testing plan was developed by CBER and will be used for routine lot release.

c. Facility Description / Inspection

Facility information and data provided in the BLA for the manufacture of the Anti-Human Globulin (AHG) were reviewed by CBER and found to be sufficient and acceptable. The facilities involved in the manufacture of AHG Anti-IgG,-C3d and AHG Anti-IgG are listed in the table below. The activities performed and inspectional histories are noted in the table and are further described in the paragraphs that follow.

Table 3: Manufacturing Facilities for AHG Anti-IgG,-C3d and AHG Anti-IgG

Name/address	FEI number	DUNS number	Inspection/ waiver	Results/ Justification
<i>Final device</i> Manufacturing and Testing Bio-Rad Medical Diagnostics GmbH Industriestr. 1 Dreieich, Hessen, Germany	3002806595	312576506	Surveillance Inspection Pre-License Inspection	Team Biologics March 16 – 24, 2015 VAI CBER October 1 – 10, 2014 VAI
<i>Component</i> Manufacturing Bio-Rad Laboratories, (b) (4) [Redacted]	(b) (4)	(b) (4)	N/A*	N/A
<i>Component</i> Manufacturing (b) (4) [Redacted] [Redacted]	N/A	N/A	N/A*	N/A

--	--	--	--	--

* Due to the nature of this product the *in vitro* substance manufacturer facilities were not required to be inspected.

CBER performed a Pre-License Inspection of the Dreieich, Germany facility from October 1-10, 2014. At the end of the inspection, a Form FDA 483 with seven observations was issued. The firm responded adequately addressing all 483 observations.

Subsequent to the PLI, Team Biologics performed a surveillance inspection of the Bio-Rad Medical Diagnostics GmbH manufacturing facility from March 16-24, 2015. A 483 Form was issued. The corrective actions were found to be acceptable and the inspection was classified as Voluntary Action Indicated (VAI).

d. Environmental Assessment

Bio-Rad Medical Diagnostics GmbH included a request for categorical exclusion from performing an Environmental Assessment under 21CFR Part 25.31(c). The FDA concluded that this request is justified as the manufacturing of this product does not alter significantly the concentration and distribution of naturally occurring substances, and no extraordinary circumstances exist that would require an environmental assessment.

e. Container Closure System

The Anti-Human Globulin is filled into polypropylene plastic cards (IH-Cards) with overall dimensions of 70 x 9 x 53 mm (b) (4). Each gel card has six small columns (micro tubes) integrated into them, which are filled with the Anti-Human Globulin and (b) (4) gel. The opening of the columns is covered with a heat sealing lacquered aluminum foil (b) (4).

Bio-Rad Medical Diagnostics GmbH conducted the container closure integrity testing for the IH-Cards at their Dreieich location. This testing consisted of the (b) (4) IH-Cards at a temperature between 18°C to (b) (4). Then, the cards were tested for serological reactivity, bioburden, and visual inspection for the detection of leaks at 0, 3, 6, 9, 12, 16 (b) (4) of the IH-Cards); all acceptance criteria were met.

4. Analytical Studies

Analytical studies performed included a reproducibility/repeatability study, a lot-to-lot reproducibility study, stability studies (shelf life and on-board), sample aging and anticoagulant studies and an interfering substances study.

Shelf Life Stability

Stability studies were performed on (b) (4) conformance lots of IH-Card AHG Anti-IgG,-C3d. The lots were used to execute a real-time stability study and (b) (4) lot was used in a transport simulation study. Test methods (potency and specificity) used to evaluate the stability of IH-Card AHG Anti-IgG,-C3d are the same serological methods used for release testing.

The stability data at the (b) (4) test interval met specifications; however BMD selected a dating period of 16 months for the IH-Card AHG Anti-IgG,-C3d and IH-Card AHG Anti-IgG when stored at the recommended temperature of 18 – 25 °C.

On-Board stability

BMD provided data that supports the labeling claim for seven day storage of unopened and two hour storage of partially used IH-Cards on the IH-1000 Analyzer. On-board stability testing was performed using one lot of IH-Card AHG Anti-IgG,-C3d with the following parameters:

Table 4: On-Board stability study. Card/Storage Conditions

Test Interval	Card / Storage Conditions
t_{end}	(b) (4) months shelf life
$t_{\text{end}} + 5 \text{ days}$	(b) (4) months shelf life and five days on-board storage conditions
$t_{\text{end}} + (b) (4)$	(b) (4) months shelf life and (b) (4) on-board storage conditions
$t_{\text{end}} + (b) (4) + \text{opened}$	(b) (4) months shelf life and (b) (4) on-board storage

	conditions and (b) (4)	storage of opened cards
--	------------------------	-------------------------

The information provided in this table was extracted from the submission.

The stored cards were tested for specificity (positive and negative samples) and potency. The acceptance criteria were the same as established for quality control release testing. There were no differences between results obtained from unstored cards (t_{end}) and the cards stored at the three storage conditions. All acceptance criteria were met after (b) (4) closed and (b) (4) closed plus (b) (4) opened, when stored at the maximum allowed IH-1000 Analyzer on-board temperature of (b) (4). The testing used cards stored under worst-case conditions, as the shelf life is established at 16 months, and the on-board storage conditions exceeded the labeling claims of up to seven days (closed) and two hours (opened, mimicking re-use of partially used cards.)

Lot-to-Lot Study

A Lot-to-Lot Study was conducted internally at Bio-Rad. Three lots of IH-Card AHG Anti-IgG, -C3d were tested using a precision panel of known red blood cell phenotypes, complement-coated cells, and antibodies. Each sample was tested in duplicate, with two runs per day on five non-consecutive days over a 20 day period, providing 60 data points per sample (i.e., $60 = 3 \text{ lots} \times 2 \text{ duplicates} \times 2 \text{ runs} \times 5 \text{ days}$). The results for all three lots showed 100% agreement.

Reproducibility and Repeatability Study

A reproducibility study was conducted at three sites (two external and one internal) using one lot of IH-Card AHG Anti-IgG, C3d using the same precision panel as used in the lot-to-lot study. Each sample was tested in duplicate, with two runs per day on five non-consecutive test dates over a 20 day period, providing 60 data points per sample (i.e., $60 = 3 \text{ sites} \times 2 \text{ duplicates} \times 2 \text{ runs} \times 5 \text{ days}$). The results for the IH-Card AHG Anti-IgG,-C3d showed 100% agreement.

Limitation Study

Sensitivity of the crossmatch with the IH-Card AHG Anti-IgG,-C3d was evaluated internally by BMD. The evaluation was performed by (b) (4) containing different clinically significant antibodies (i.e., (b) (4) that reacted by the indirect antiglobulin test. (b) (4) of each of the (b) (4) antibody samples were prepared and tested with (b) (4) different ABO-compatible donor samples with heterozygous expression of the corresponding antigen with the IH-Card AHG cards in comparison to a FDA licensed reference, Ortho MTS Anti-IgG Card. The following table indicates the titer endpoint achieved with the investigational cards in comparison with the reference method.

Table 5: Comparison of titer endpoints between the investigational cards and the reference method

Antibody Specificity	IH-Card AHG Anti-IgG,-C3d	Ortho MTS Anti-IgG Card
(b) (4)	16	16
	8	16
	256	256
	128	256
	8	8
	8	8
	512	64
	256	64
	128	64
	64	64

The information provided in this table was extracted from the submission.

The data demonstrates reliable detection of red cell antibodies in the crossmatch assay on the IH-1000 Analyzer with the investigational IH-IH-Card AHG Anti-IgG,-C3d. The investigational IH-Card AHG Anti-IgG,-C3d exceeded titer endpoints in three out of ten tests, equal titer endpoints in five out of ten tests, and one titer endpoint lower in two tests (Anti-D and Anti-K).

Sample aging and anticoagulant studies

Data for supporting the sample age and type, and preservative solutions were obtained from the clinical comparison studies and internal studies conducted by BMD to support the labeling claims regarding specimens used in testing. The data demonstrated that the specimen type and age (as listed below) are suitable for use in testing with the IH-Card AHG Anti-IgG,-C3d. The IH-Card AHG Anti-IgG,-C3d Instructions for Use and the IH-1000 User Manual were revised to reflect the following sample age restrictions:

Table 6: Specimen type and age

Test	
Indirect Antiglobulin Test	<p>Fresh EDTA, citrated, and serum samples are stored at 2 to 8 °C and may be tested up to 10 days post-collection.</p> <p>Frozen serum or plasma that has been stored at -20 °C may be tested:</p> <ul style="list-style-type: none"> • Sodium citrate samples frozen for up to 674 days • EDTA samples frozen for up to 30 days • Serum samples frozen for up to 26 days
Direct Antiglobulin Test	EDTA and cord blood samples may be tested up to five days post-collection and stored at 2 to 8 °C.
Crossmatching	<p>Donor blood collected in CPD, CP2D, and (b) (4) may be tested until the expiration of the unit and stored at 1 to 8 °C</p> <p>Donor blood containing AS-1 or AS-3 may be tested up to 30 days post-collection.</p>

	Recipient samples collected in EDTA may be tested up to 10 days post-collection.
--	--

The information provided in this table was extracted from the submission

Interfering substances

A study for interfering substances (hemolysis, icteric, and lipemic) was conducted internally at BMD. The testing of samples was performed with the IH-Card AHG Anti-IgG,-C3d in comparison to an FDA licensed product. The results obtained with the IH-Card AHG cards were concordant with the results obtained with the reference reagent. The results demonstrated that higher than normal concentrations of triglycerides, bilirubin, or hemoglobin did not have an adverse effect on the performance of the IH-Card AHG Anti-IgG,-C3d.

Samples from the Elderly

More than ^{(b) (4)} EDTA samples from patients older than 80 years of age, were tested with the IH-Card AHG Anti-IgG,C3d using the antibody detection method.

The results were compared with a reference (i.e., FDA licensed Anti-Human Globulin reagent). All results of the antibody detection except for one sample were concordant between the IH-Card AHG Anti-IgG, C3d and the reference. The sample with discordant results was negative with the reference and positive with the IH-Card AHG Anti-IgG,-C3d. The quantity of the sample was insufficient for further resolution testing.

Incubation temperatures used for serological testing

When performing the various assays using the indirect antiglobulin test, the IH-Card AHG Anti-IgG,-C3d is incubated at 37 °C. The IH-1000 Analyzer incubator temperature range of 35 to 38 °C is controlled by the instrument software. A study tested IgG antibodies and complement-coated cells with the IH-Card AHG cards incubated at the minimum and maximum parameters. Data confirmed that the IH-Card AHG cards perform as expected

when the cards are incubated at the IH-1000 Analyzer minimum and maximum incubator temperatures.

5. Performance Studies

a) Clinical Studies

BMD conducted a clinical comparison study to evaluate the IH-Card AHG Anti-IgG,-C3d performance in a clinical setting. The studies were performed at four external testing sites: Puget Sound Blood Center (PSBC) located in Renton, Washington, Miriam Hospital (MH) located in Providence, Rhode Island, Vanderbilt University Medical Center (VUMC) located in Nashville, Tennessee, and LifeSource Testing Laboratory (LSTL) located in Rosemont, Illinois. In addition, internal studies were conducted at the Bio-Rad Research and Development Laboratory in Cressier, Switzerland.

The clinical comparison study evaluated the performance of the IH-1000 Analyzer for antibody detection, antibody identification, direct antiglobulin testing, and antiglobulin crossmatch testing. The study included the use of BMD IH-Cell products (Reagent Red Blood Cells, $0.6 \pm 0.1\%$), quality control material, IH-LISS (diluent), at least two lots of IH-Card AHG Anti-IgG,-C3d, and three different IH-1000 Analyzers. The clinical study consisted of testing the IH-Card AHG Anti-IgG,-C3d in parallel with an FDA-licensed reagent and cleared instruments for concordance. The table shown below shows the comparator test method used by site.

Table 7: Comparator test method by site

Test				
	Immucor Galileo Neo (Capture R)	Ortho ID- MTS Anti- IgG Card (ProVue)	Immucor Galileo Echo (Capture R)	BMD TANGO optimo System

Antibody Identification	Not Performed	Ortho ID-MTS Anti-IgG Card (Manual)	Ortho ID-MTS Anti-IgG Card (Manual)	BMD Anti-IgG (Tube Test)
Antiglobulin crossmatch	Not Performed	Ortho ID-MTS Anti-IgG Card (ProVue)	Ortho ID-MTS Anti-IgG Card (Manual)	Not Performed
Direct Antiglobulin Test	Not Performed	Ortho Anti-IgG,-C3d (Tube Test)	Not Performed	BMD Anti-IgG,-C3d (Tube Test)

Note: internal testing used BMD FDA licensed reagents as the reference.

This table contains information extracted from the submission.

Results from the investigational method and the reference assay were compared for concordance. If the results of the investigational and reference method were discordant, testing was repeated. If the results were still discordant after repeat testing, then a third FDA licensed reagent was used as a reference. If the antibody identification results were discordant after repeat testing, antigen dosage and profile were evaluated as potential cause for discrepancy.

BMD used the following acceptance criteria for the studies:

- Antibody Detection, Antibody Identification, and Direct Antiglobulin Test: The lower bound of the one-sided 95% confidence intervals for the positive percent agreement and the negative percent agreement had to exceed 0.95.
- Antiglobulin Crossmatch: The lower bound of the one-sided 95% confidence intervals for the positive percent agreement and the negative percent agreement had to exceed 0.99.

Some of the IH-Card AHG Anti-IgG, C3d testing results described below in the clinical sites failed to meet the acceptance criteria to support approval. In a submission issue meeting held on January 14, 2016, FDA requested BMD to perform additional testing in an In-House Performance Study. The result verification process would include reviewing and editing (if necessary) of the camera images and assay interpretation. Editing of the results was not limited to indeterminate results. Initially positive or negative reactions may have to be edited after visual review. The edited results of the studies performed at the clinical sites are described below.

Antibody Detection

The detection of unexpected antibodies was conducted using IH-Card AHG Anti-IgG,-C3d and three different IH-Cell products used to screen for red blood cell antibodies. The following table summarizes the number of samples collected and analyzed per site.

Table 8: Number of samples collected and analyzed per site

Site			
PSBC	500	0	0
MH	0	0	800
VUMC	0	902	400
LSTL	0	1,100	0
BMD	137	0	137
Sub-total	637	2,002	1,337
TOTAL	3,976		

This table contains data extracted from the submission.

The performance evaluations for the IH-Card AHG Anti-IgG,-C3d were based on the combined edited results for all three IH-Cell products per IH-Card.

Table 9: Results of IH-Card AHG-Anti-IgG,-C3d for antibody detection

	ANTIBODY DETECTION		
	(Edited Results)		
	Positive Percent Agreement	Negative Percent Agreement	Overall Percent Agreement
	[Lower 95% CB]	[Lower 95% CB]	[Lower 95% CB]
IH-Card AHG Anti-IgG,-C3d	95.42% 146/153 [91.58%]	98.77% 3774/3821 [98.43%]	98.59% 3920/3976 [98.24%]

This table contains information extracted from the submission.

The positive percent agreement did not meet the 95% acceptance criteria when combining results for all screened Reagent Red Blood Cells (RRBCs) on the IH-Card AHG Anti-IgG, -C3d.

A summary of the factors contributing to not meeting the acceptance criteria is listed below:

- Five samples reacted negative with the investigational method and positive with the reference method, and no red cell antibodies were identified.
- In one of the five samples mentioned above, the technologist performing the testing noted that the positive reaction in the Ortho ID-MTS Anti-IgG Card was due to diffuse hemolysis in the sample and not true agglutination so would report as negative. The sample was negative on the IH-1000 on repeat testing and positive on the reference method.

Antibody Identification

Antibody identification testing was performed on all samples with positive or equivocal antibody screens by any method. The identification of unexpected antibodies was conducted

using IH-Card AHG Anti-IgG,-C3d with different IH-Cell products used for red blood cell antibody identification.

Table 10: Summary of the number of samples collected and analyzed per site

Site	IH-Panel 11 (11-cell panel)	IH-Panel 6 Plus (6-cell panel)
PSBC	0	0
MH	48	40
VUMC	71	47
LSTL	3	2
BMD (internal)	84	0
Total	206	89

This table contains data extracted from the submission.

The performance evaluations for the IH-Card AHG Anti-IgG,-C3d were based on the combined initial results for all three IH-Cell products per IH-Card (see Table 11).

Table 11: Results of IH-Card AHG-Anti-IgG,-C3d for antibody identification

	ANTIBODY IDENTIFICATION (Initial Results) (per sample)		
	Positive Percent Agreement [Lower 95% CB]	Negative Percent Agreement [Lower 95% CB]	Overall Percent Agreement [Lower 95% CB]
	IH-Card AHG Anti-IgG,-C3d	97.92% 47/48 [90.49%]	87.45% 216/247 [83.44%]

The table presented above contains information extracted from the submission.

A total of 48 antibodies were identified in the 295 samples tested. Both the investigational and reference methods identified the same antibodies in 40/48 samples. – There were 31 discordant negative results- i.e. the reference method did not identify an antibody but the investigational method either identified an antibody or gave nonspecific positive results. The positive percent agreement and the negative percent agreement did not meet the 95% acceptance criteria.

A summary of the factors contributing to not meeting the acceptance criteria is listed below:

Of the 48 antibody positive samples:

- The number of positive samples (48) is not large enough to achieve statistical significance.
- In 7/48 samples, only the investigational method detected antibodies of the following specificities: anti-E, anti-Fy^a, anti-M, anti-Lu^a, and anti-Jk^b.
- In 1/48 sample – The reference method identified an anti-Le^a, and the investigational method obtained nonspecific positive reactivity.

Of the 31 discordant negative results:

- 7/31 samples – Tested negative with the reference method only, but the investigational method identified the following antibodies. Anti-E, anti-Fy^a, anti-Co^b, anti-M, and anti-Lu^a.
- 15/31 samples – Unidentifiable antibody with the investigational method and negative with the reference method.
- 9/31 samples - Negative with reference method, and the investigational method had one equivocal microtube result with each sample, which precluded an initial interpretation of negative.

Direct Antiglobulin Test

The direct antiglobulin test (DAT) was evaluated using the IH-Card AHG Anti-IgG,-C3d in comparison with a licensed Anti-IgG,-C3d reagent. The comparison testing was performed on a total of 650 samples; 250 at MH and 400 at LSTL.

A total of 65 samples tested positive on the reference method. There was agreement with the investigational method in 58/65 samples.

Table 12: Results for DAT testing with the IH-Card AHG Anti-IgG, -C3d

N = 650 samples

IH-Card AHG Anti-IgG, -C3d		FDA Licensed AHG Anti-IgG, -C3d Reference				Positive % Agreement [one-sided Exact 95% Lower CI]	Negative % Agreement [one-sided Exact 95% Lower CI]	Overall % Agreement [one -sided Exact 95% Lower CI]
		Pos	EQV	Neg	Total			
Investigational Anti-IgG, -C3d 74010 N = 650	Pos	58	0	14	72	89.23% 58/65 [80.72%]	97.61% 571/585 [96.28%]	96.77% 629/650 [95.38%]
	EQV	0	0	0	0			
	Neg	7	0	571	578			
	Total	65	0	585	650			

This table was extracted from the revised statistical summary dated March 4, 2016, which is included in the submission.

The negative percent agreement and the overall agreement exceeded the 95% acceptance criteria with a 95% lower confidence interval (LCI). The positive percent agreement did not meet the acceptance criteria of 95% with a 95% LCI.

A summary of the factors contributing to not meeting the positive percent agreement acceptance criteria is listed below:

- The number of positive samples tested is not large enough to be statistically significant.
- One sample was negative by the reference method.
- One sample was positive by the reference method.
- Four samples were not tested by the reference method.

Antiglobulin Crossmatch

The antiglobulin crossmatch using the IH-Card AHG Anti-IgG,-C3d was performed with 90 samples (34 antibody positive and 56 antibody negative) that were crossmatched with 760 compatible and incompatible donor units. Crossmatches performed per site were 350 at MH and 410 at VUMC. The results are presented in Table 13.

Table 13: Results for the detection of antiglobulin crossmatches with IH-Card AHG Anti-IgG,-C3d

N = 760 crossmatches

IH-Card AHG Anti-IgG, -C3d		FDA Licensed Anti-IgG Reference				Positive % Agreement [one-sided Exact 95% Lower CI]	Negative % Agreement [one-sided Exact 95% Lower CI]	Overall % Agreement [one -sided Exact 95% Lower CI]
		Pos	EQV	Neg	Total			
Investigational Anti-IgG, -C3d 74010 N = 760	Pos	317	0	19	336	95.20% 317/333 [92.79%]	95.55% 408/427 [93.54%]	95.39% 725/760 [93.94%]
	EQV	0	0	0	0			
	Neg	16	0	408	424			
	Total	333	0	427	760			

This table was extracted from the revised statistical summary dated March 4, 2016, which is included in the submission.

The negative and positive percent agreement did not meet the 99% acceptance criteria with a 95% LCI.

Failure to meet the acceptance criteria was due to a high number (18) of false negative results with the reference method when the expected result was positive (incompatible). In all these crossmatches, the investigational method was positive (incompatible) and concordant with the expected results.

Summary of Clinical Studies Results

During the clinical study the test results obtained with the following assays failed to meet the primary endpoint:

- Detection of Unexpected Antibodies with IH-Card AHG Anti-IgG, -C3d

- The positive percent agreement did not meet the 95% acceptance criteria.
- Antibody-identification with IH-Card AHG Anti-IgG, -C3d
 - The positive percent agreement and the negative percent agreement did not meet the 95% acceptance criteria.
- Direct Antiglobulin Test (DAT) with IH-Card AHG Anti-IgG, -C3d
 - The positive percent agreement did not meet the 95% acceptance criteria.
- AHG Crossmatch with IH-Card AHG Anti-IgG, C3d
 - The positive and negative agreement did not meet the 99% acceptance criteria.

Because the assays mentioned above failed to meet the primary endpoint, on January 14, 2016, FDA requested BMD to do additional testing in an In-House Performance Study, to support the approval of the IH-Card AHG Anti-IgG, -C3d for detection and identification of unexpected antibodies, DAT and AHG compatibility testing on the IH-1000. For the in-house performance study BMD used known positive and negative samples and/or contrived samples.

b) In-House Performance Study

The additional studies requested by FDA for the AHG reagent and RRBCs are shown in the table below.

- Green indicates the study endpoint was met in the clinical study and no further study is required
- Red indicates that the study endpoint was not met in the clinical study and that testing of additional samples in the in-house performance study is required to demonstrate that the predetermined study endpoints are met.

Table 14: Additional studies for IH-Card Anti-IgG,-C3d and for RRBCS

Assay	Antibody screening		Antibody Identification		Direct Antiglobulin Testing (DAT)		AHG Crossmatching*	
	PPA	NPA	PPA	NPA	PPA	NPA	PPA	NPA
IH-Card AHG Anti-IgG	Red	Green	Red	Red	Red	Green	Green	Red
IH-Card AHG Anti-IgG, -C3d	Red	Green	Red	Red	Red	Green	Red	Red
Card	IH-Card Anti-IgG	AHG	IH-Card Anti-IgG	AHG -C3d				
	PPA	NPA	PPA	NPA				
IH-Cell Pool	Red	Green	Red	Green				
IH-Cell I-II	Red	Green	Red	Green				
IH-Cell I-II-III	Red	Green	Red	Green				
IH-Panel 11	Red	Red	Red	Red				
IH-Panel 6 Plus	Green	Red	Red	Red				
IH-Panel 11 Papain	Red	Red	NA					

* detection of incompatibility due to IgG antibodies

NPA = Negative Percent Agreement; PPA = Positive Percent Agreement

The information provided in the table above was obtained from the in-house performance study report included in the submission.

Sample Size for Eight Additional In-House Performance Studies

1) Number of samples for RRBCs

- Antibody Screening: 191 antibody positive samples
- Antibody Identification: 425 antibody negative samples and 128 antibody positive samples.

2) Number of samples for AHG Anti-IgG,-C3d crossmatching

- Number of samples by characterization result: 38 antibody positive frozen samples collected in CPD and 9 antibody positive frozen samples collected in sodium citrate. Thirty-nine (39) antibody negative frozen EDTA samples.

- Total number of crossmatches: 299 samples (with 23 antigen positive, ABO compatible donor units collected in sodium citrate) and 301 samples (with 19 ABO compatible donor units collected in sodium citrate)

3) Number of known and/or contrived DAT positive samples

- Forty six (46) contrived samples
- 22 well-characterized samples

Testing Algorithm

After the sample characterization with at least two FDA licensed methods was completed, the samples were tested with the investigational reagents using the same sample aliquot used for the characterization. Before final test interpretation, all results were validated by trained operators.

BMD adopted the same acceptance criteria used for the clinical study.

Results

A) Detection of unexpected antibodies with IH-Card AHG Anti-IgG, -C3d

Sixty-four samples were tested with IH-Cell Pool and IH-Cell I-II. Sixty three (63) of those samples were also tested with IH-Cell I-II-III for the total of 191 antibody screening tests.

The information presented in the tables shown below was obtained from the in-house performance study report included in the submission.

Table 15: Testing results with IH-Cell Pool

Ab screen IH-Card AHG Anti-IgG, - C3d IH-Cell Pool		Expected result	
		Positive	Negative
Investigational method	Positive	64	NA
	EQV	0	
	Negative	0	
Total		64	
% Agreement		100%	
One-sided Exact 95% LCL		95.43%	

Table 16: Testing results with IH-Cell I-II

Ab screen IH-Card AHG Anti-IgG, - C3d IH-Cell I-II		Expected result	
		Positive	Negative
Investigational method	Positive	64	NA
	EQV	0	
	Negative	0	
Total		64	
% Agreement		100%	
One-sided Exact 95% LCL		95.43%	

Table 17: Testing results with IH-Cell I-II-III

Ab screen IH-Card AHG Anti-IgG, - C3d IH-Cell I-II-III		Expected result	
		Positive	Negative
Investigational method	Positive	63	NA
	EQV	0	
	Negative	0	
Total		63	
% Agreement		100%	
One-sided Exact 95% LCL		95.36%	

Conclusion: For each of the three panels, the results for the samples shown in the three tables met the predetermined acceptance criteria.

B) Antibody identification with IH-Card AHG Anti-IgG, -C3d

Table 18: Testing results with IH-Panel 11

AbID IH-Card AHG Anti-IgG, - C3d IH-Panel 11		Expected result	
		Positive	Negative
Investigational method	Positive	64	4 ^a
	EQV	0	0
	Negative	0	208 ^b
Total		64	212
% Agreement		100%	98.11%
One-sided Exact 95% LCL		95.43%	95.73%

Conclusion: The results met the predetermined acceptance criteria

Table 19: Testing results with IH-Panel Plus 6

AbID IH-Card AHG Anti-IgG, - C3d IH-Panel Plus 6		Expected result	
		Positive	Negative
Investigational method	Positive	64	3 ^a
	EQV	0	0
	Negative	0	210 ^{b,c}
Total		64	213
% Agreement		100%	98.59%
One-sided Exact 95% LCL		95.43%	96.40%

Conclusion: The results met the predetermined acceptance criteria

Table 20: Testing results with IH-Panel 11

IH-Card AHG Anti-IgG, -C3d (74010)	Investigational RRBC Results		Negative % Agreement [One-sided Exact 95% Lower CI]
	Pos		
IH-Panel 11 (07100)	Pos	58	96.66% 1765/1826 [95.88 %]
	EQV	3	
	Neg	1765	
	Total	1826	

Conclusion: The results met the predetermined acceptance criteria.

C) Direct Antiglobulin Test (DAT) with IH-Card AHG Anti-IgG, -C3d

Table 21: Testing results

DAT IH-Card AHG Anti-IgG, - C3d		Expected result	
		Positive	Negative
Investigational method	Positive	69	NA
	EQV	0	
	Negative	0	
Total	69		
% Agreement		100%	
One-sided Exact 95% LCL		95.75%	

Conclusion: The results met the predetermined acceptance criteria

D) AHG Crossmatch with IH-Card AHG Anti-IgG, -C3d

Testing included 301 expected compatible and 299 incompatible crossmatches.

Table 22: Testing results

AHG crossmatch IH-Card AHG Anti-IgG, - C3d		Expected result	
		Positive	Negative
Investigational method	Positive	299	0
	EQV	0	0
	Negative	0	301 ^{a,b}
Total		299	301
% Agreement		100%	100%
One-sided Exact 95% LCL		99.00%	99.01%

Conclusion: The results met the predetermined acceptance criteria.

Summary of in-house performance study results

During the in-house performance study, the test results obtained for the following assays met the acceptance criteria:

- Detection of Unexpected Antibodies with IH-Card AHG Anti-IgG, -C3d (positive agreement)

- Antibody-identification with IH-Card AHG Anti-IgG, -C3d (positive and negative agreement).
- Direct Antiglobulin Test (DAT) with IH-Card AHG Anti-IgG, -C3d (positive percent agreement).
- AHG Crossmatch with IH-Card AHG Anti-IgG, C3d (positive and negative agreement).

b) Pediatrics

Data from the clinical studies found that cord blood samples could be properly tested using the IH-Card AHG Anti-IgG,-C3d with the IH-System.

c) Other Special Populations

In the Sample Aging and Anticoagulant analytical studies, cord blood samples and samples from elderly patients were included.

d) Overall Comparability Assessment

The results of the clinical comparison study and product labeling support the conclusion that the IH-Card AHG Anti-IgG,-C3d and the IH-Card AHG Anti-IgG formulated for automated testing on the IH-1000 Analyzer are safe and effective. In addition, the lots manufactured in support of these submissions, demonstrate the reliability of the manufacturing process to consistently produce IH-Cards AHG that meet established specifications, perform as intended, and remain stable throughout its shelf life.

6. Advisory Committee Meeting

The Bio-Rad Medical Diagnostics GmbH IH-System does not include novel technology; therefore, an advisory committee meeting was not held or required.

7. Other Relevant Regulatory Issues

The review committee members from DBCD, DMPQ, DB, DCM, and DBSQC reviewed their specific sections of the BLA and Efficacy Supplement and resolved any issues through information requests and teleconferences with BMD. The Review Team sought the expertise of their respective management, when warranted. No internal or external disagreements were communicated to the regulatory project manager or chairperson. All reviewers recommended approval of the Anti-Human Globulin (Formulated for Automated Testing) BLA and Efficacy Supplement.

No postmarketing commitments are associated with the Anti-Human Globulin (Formulated for Automated Testing) BLA and Efficacy Supplement.

8. Labeling

The labeling for Anti-Human Globulin (Formulated for Automated Testing) complies with Title 21 CFR 610.62, 610.63, 610.64, 610.65, 660.55 and 809.10. In the teleconference held on June 8, 2016, BMD agreed to include in the final labeling all the recommendations provided by FDA. The final labeling for IH-Card AHG Anti-IgG,-C3d and IH-Card AHG Anti-IgG includes the following changes.

- The use of a table format to represent the results from the clinical study and the in-house performance study in the package inserts of AHG under the Specific Performance Characteristics section.
- All package inserts (BGRs, AHGs, and RRBCs) are revised to define visual reading and editing (if applicable) as a required process. All package inserts are revised to include the “edited results” rather than the “initial results”.

Unique Device Identification (UDI) review performed by CBER found the required elements to comply with Title 21 CFR 830. The labeling met the UDI requirements ahead of the September 24, 2016 compliance date for this classification of medical devices.

9. Recommendations and Risk/ Benefit Assessment

a) Recommended Regulatory Action

The review committee members, representing the necessary review disciplines (DBCD, DMPQ, DB, DCM, and DBSQC) recommend approval. These are independent conclusions based on content of the BLA, issues satisfactorily resolved during the review cycle, and concurrence by the respective management. No internal or external disagreements were brought to the attention of the chairperson.

b) Risk/ Benefit Assessment

- The IH-1000 Automated Analyzer and the reagents used by the IH-System, provide potential advantages to support transfusion medicine.
- The clinical benefits using the IH-System include greater patient safety and timely availability of transfusion products to the patient through improved productivity.
- Features that impact patient safety include reduction in errors associated with subjective interpretation due to manual testing, transcription errors, test errors (i.e., using expired reagents or the wrong reagent), and the capability to review of stored test results, if necessary.
- Features that impact timely availability of transfusion products include reduction in hands-on technologist time by automating the process, time required for recording assay reagents, controls, and equipment, as well as turn-around time.

c) Recommendation for Postmarketing Activities

There are no postmarketing commitments associated with the Anti-Human Globulin (Formulated for Automated Testing) BLA or Efficacy Supplement.