

August 11, 2017

W. Ian Lipkin, MD
Director
The Center for Infection and Immunity
Columbia University
722 West 168th St., 17th Floor
New York, NY 10032

Dear Dr. Lipkin:

This letter is in response to your request that the Food and Drug Administration (FDA) issue an Emergency Use Authorization (EUA) for emergency use of The Center for Infection and Immunity, Columbia University's ("Columbia University") CII-ArboViroPlex rRT-PCR assay for the qualitative detection and differentiation of RNA from Zika virus, dengue virus, chikungunya virus, and West Nile virus in serum, and for the qualitative detection of Zika virus RNA in urine (collected alongside a patient-matched serum specimen). The assay is intended for use with specimens collected from individuals meeting Centers for Disease Control and Prevention (CDC) Zika virus clinical criteria (e.g., clinical signs and symptoms associated with Zika virus infection) and/or CDC Zika virus epidemiological criteria (e.g., history of residence in or travel to a geographic region with active Zika transmission at the time of travel, or other epidemiological criteria for which Zika virus testing may be indicated), by laboratories in the United States (U.S.) that are certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. § 263a, to perform high complexity tests, or by similarly qualified non-U.S. laboratories, pursuant to section 564 of the Federal Food, Drug, and Cosmetic Act (the Act) (21 U.S.C. § 360bbb-3).¹ Assay results are for the identification of Zika, dengue, chikungunya, and West Nile viral RNA. Viral RNA is generally detectable in serum during the acute phase of infection and, according to the updated CDC Guidance for U.S. Laboratories Testing for Zika Virus Infection,² up to 14 days in serum and urine (possibly longer in urine), following onset of symptoms, if present. Positive results are indicative of current infection.

On February 26, 2016, pursuant to section 564(b)(1)(C) of the Act (21 U.S.C. § 360bbb-3(b)(1)(C)), the Secretary of Health and Human Services (HHS) determined that there is a significant potential for a public health emergency that has a significant potential to affect national security or the health and security of U.S. citizens living abroad and that involves Zika virus.³ Pursuant to section 564(b)(1) of the Act (21 U.S.C. § 360bbb-3(b)(1)), and on the basis

¹ For ease of reference, this letter will refer to "laboratories in the United States (U.S.) that are certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA), 42 U.S.C. § 263a, to perform high complexity tests, or...similarly qualified non-U.S. laboratories" as "authorized laboratories."

² Available at <http://www.cdc.gov/zika/laboratories/lab-guidance.html> (last updated on July 24, 2017).

³ As amended by the Pandemic and All Hazards Preparedness Reauthorization Act, Pub. L. No. 113-5, under section

of such determination, the Secretary of HHS then declared that circumstances exist justifying the authorization of the emergency use of *in vitro* diagnostic tests for detection of Zika virus and/or diagnosis of Zika virus infection, subject to the terms of any authorization issued under 21 U.S.C. § 360bbb-3(a).⁴

Having concluded that the criteria for issuance of this authorization under section 564(c) of the Act (21 U.S.C. § 360bbb-3(c)) are met, I am authorizing the emergency use of the CII-ArboViroPlex rRT-PCR assay (as described in the Scope of Authorization section of this letter (Section II)) in individuals meeting CDC Zika virus clinical criteria (e.g., clinical signs and symptoms associated with Zika virus infection) and/or CDC Zika virus epidemiological criteria (e.g., history of residence in or travel to a geographic region with active Zika transmission at the time of travel, or other epidemiological criteria for which Zika virus testing may be indicated) (as described in the Scope of Authorization section of this letter (Section II)) for the detection of Zika virus infection by authorized laboratories, subject to the terms of this authorization.

I. Criteria for Issuance of Authorization

I have concluded that the emergency use of the CII-ArboViroPlex rRT-PCR assay for the detection of Zika virus and diagnosis of Zika virus infection in the specified population meets the criteria for issuance of an authorization under section 564(c) of the Act, because I have concluded that:

1. The Zika virus can cause Zika virus infection, a serious or life-threatening disease or condition to humans infected with the virus;
2. Based on the totality of scientific evidence available to FDA, it is reasonable to believe that the CII-ArboViroPlex rRT-PCR assay, when used with the specified instrument(s) and in accordance with the Scope of Authorization, may be effective in detecting Zika virus and diagnosing Zika virus infection, and that the known and potential benefits of the CII-ArboViroPlex rRT-PCR assay for detecting Zika virus and diagnosing Zika virus infection outweigh the known and potential risks of such product; and
3. There is no adequate, approved, and available alternative to the emergency use of the CII-ArboViroPlex rRT-PCR assay for detecting Zika virus and diagnosing Zika virus infection.⁵

II. Scope of Authorization

I have concluded, pursuant to section 564(d)(1) of the Act, that the scope of this authorization is limited to the use of the authorized CII-ArboViroPlex rRT-PCR assay by authorized laboratories for the qualitative detection and differentiation of RNA from Zika virus, dengue virus,

564(b)(1)(C) of the Act, the Secretary may make a determination of a public health emergency, or of a significant potential for a public health emergency.

⁴ HHS. *Determination and Declaration Regarding Emergency Use of in Vitro Diagnostic Tests for Detection of Zika Virus and/or Diagnosis of Zika Virus Infection*. 81 Fed. Reg. 10878 (March 2, 2016).

⁵ No other criteria of issuance have been prescribed by regulation under section 564(c)(4) of the Act.

chikungunya virus, and West Nile virus in serum, and for the qualitative detection of Zika virus RNA in urine (collected alongside a patient-matched serum specimen) in individuals meeting CDC Zika virus clinical criteria (e.g., clinical signs and symptoms associated with Zika virus infection) and/or CDC Zika virus epidemiological criteria (e.g., history of residence in or travel to a geographic region with active Zika transmission at the time of travel, or other epidemiological criteria for which Zika virus testing may be indicated).

The Authorized CII-ArboViroPlex rRT-PCR assay

The CII-ArboViroPlex rRT-PCR assay is a multiplex one-step real-time reverse transcription polymerase chain reaction (rRT-PCR) assay for the qualitative detection and differentiation of RNA from Zika virus, dengue virus, chikungunya virus, and West Nile virus in serum, and other authorized specimen types. The CII-ArboViroPlex rRT-PCR assay can also be used for the qualitative detection of Zika virus RNA in urine when collected alongside a patient-matched serum specimen and other authorized whole blood derived specimen types.

To perform the CII-ArboViroPlex rRT-PCR assay, the RNA is first extracted and purified from the patient specimen. The RNA is then reverse transcribed into cDNA which is amplified using the primer set and detected using the specific probe. The rRT-PCR is performed on the CFX96 Real-Time PCR Detection System (Bio-Rad), or other authorized instruments.

The CII-ArboViroPlex rRT-PCR assay includes the following materials or other authorized materials:

- ZIKV-MIX, DENV-MIX, CHIKV-MIX, WNV-MIX and RP-MIX vials containing primers and probes for the assay targets and internal control
- ZPC, DPC, CP, WPC, HSC, eHSC, NTC vials containing the positive and negative controls used in the assay
- Diluent vial used to reconstitute dried vials

The CII-ArboViroPlex rRT-PCR assay also requires the use of additional materials and ancillary reagents that are not included with the test but are commonly used in clinical laboratories and are described in the authorized CII-ArboViroPlex rRT-PCR assay Instructions for Use.

The CII-ArboViroPlex rRT-PCR assay requires the following control materials, or other authorized control materials; all controls listed below must generate expected results in order for a test to be considered valid, as outlined in the CII-ArboViroPlex rRT-PCR assay Instructions for Use:

- Human Specimen Control: A human cell culture preparation used as an extraction control and positive control for the RNase P primer and probe set that is extracted and tested concurrently with the test specimens.
- Extracted Human Specimen Control (eHSC): Extracted total nucleic acid from a human cell culture preparation known to contain RNase P (eHSC), but negative for viral targets, is used as a control for performance of RNase P primer/probe set and PCR reagent function.
- Positive Controls for viruses: Run with each batch of patient specimens. Monitors for failures of rRT-PCR reagents and reaction conditions.
 - ZIKV Positive Control (ZPC), synthetic *in vitro* transcribed RNA

- DENV Positive Control (DPC), synthetic *in vitro* transcribed RNA
- CHIKV Positive Control (CPC), synthetic *in vitro* transcribed RNA
- WNV Positive Control (WPC), synthetic *in vitro* transcribed RNA
- No Template Control (NTC): Sterile, nuclease-free water—two NTC run with each PCR plate. Monitors for reagent and system contamination.
- RNase P control in clinical samples: All clinical samples and HSC are tested for human RNase P, using the RP primer and probe set, to control for specimen quality and as an indicator that nucleic acid resulted from the extraction process.

The above described CII-ArboViroPlex rRT-PCR assay, when labeled consistently with the labeling authorized by FDA entitled “CII-ArboViroPlex rRT-PCR assay Instructions for Use” (available at <http://www.fda.gov/MedicalDevices/Safety/EmergencySituations/ucm161496.htm>), which may be revised by Columbia University in consultation with, and with concurrence of, the Division of Microbiology Devices (DMD)/Office of In Vitro Diagnostics and Radiological Health (OIR)/Center for Devices and Radiological Health (CDRH), is authorized to be distributed to and used by authorized laboratories under this EUA, despite the fact that it does not meet certain requirements otherwise required by federal law.

The above described CII-ArboViroPlex rRT-PCR assay is authorized to be accompanied by the following information pertaining to the emergency use, which is authorized to be made available to healthcare providers and patients, including pregnant women:

- Fact Sheet for Healthcare Providers: Interpreting CII-ArboViroPlex rRT-PCR Assay Test Results
- Fact Sheet for Patients: Understanding Results from the CII-ArboViroPlex rRT-PCR Assay

As described in Section IV below, Columbia University and its authorized distributors are also authorized to make available additional information relating to the emergency use of the authorized CII-ArboViroPlex rRT-PCR assay that is consistent with, and does not exceed, the terms of this letter of authorization.

I have concluded, pursuant to section 564(d)(2) of the Act, that it is reasonable to believe that the known and potential benefits of the authorized CII-ArboViroPlex rRT-PCR assay in the specified population, when used for detection of Zika virus and to diagnose Zika virus infection and used consistently with the Scope of Authorization of this letter (Section II), outweigh the known and potential risks of such a product.

I have concluded, pursuant to section 564(d)(3) of the Act, based on the totality of scientific evidence available to FDA, that it is reasonable to believe that the authorized CII-ArboViroPlex rRT-PCR assay may be effective in the detection of Zika virus and diagnosis of Zika virus infection, when used consistently with the Scope of Authorization of this letter (Section II), pursuant to section 564(c)(2)(A) of the Act.

FDA has reviewed the scientific information available to FDA, including the information supporting the conclusions described in Section I above, and concludes that the authorized CII-

ArboViroPlex rRT-PCR assay, when used for detection of Zika virus and to diagnose Zika virus infection in the specified population (as described in the Scope of Authorization of this letter (Section II)), meets the criteria set forth in section 564(c) of the Act concerning safety and potential effectiveness.

The emergency use of the authorized CII-ArboViroPlex rRT-PCR assay under this EUA must be consistent with, and may not exceed, the terms of this letter, including the Scope of Authorization (Section II) and the Conditions of Authorization (Section IV). Subject to the terms of this EUA and under the circumstances set forth in the Secretary of HHS's determination described above and the Secretary of HHS's corresponding declaration under section 564(b)(1), the CII-ArboViroPlex rRT-PCR assay described above is authorized to detect Zika virus and diagnose Zika virus infection in individuals meeting CDC Zika virus clinical criteria (e.g., clinical signs and symptoms associated with Zika virus infection) and/or CDC Zika virus epidemiological criteria (e.g., history of residence in or travel to a geographic region with active Zika virus transmissions at the time of travel, or other epidemiological criteria for which Zika virus testing may be indicated).

This EUA will cease to be effective when the HHS declaration that circumstances exist to justify the EUA is terminated under section 564(b)(2) of the Act or when the EUA is revoked under section 564(g) of the Act.

III. Waiver of Certain Requirements

I am waiving the following requirements for the CII-ArboViroPlex rRT-PCR assay during the duration of this EUA:

- Current good manufacturing practice requirements, including the quality system requirements under 21 CFR Part 820 with respect to the design, manufacture, packaging, labeling, storage, and distribution of the CII-ArboViroPlex rRT-PCR assay.
- Labeling requirements for cleared, approved, or investigational devices, including labeling requirements under 21 CFR 809.10 and 21 CFR 809.30, except for the intended use statement (21 CFR 809.10(a)(2), (b)(2)); adequate directions for use (21 U.S.C. 352(f)), (21 CFR 809.10(b)(5), (7), and (8)); any appropriate limitations on the use of the device including information required under 21 CFR 809.10(a)(4); and any available information regarding performance of the device, including requirements under 21 CFR 809.10(b)(12).

IV. Conditions of Authorization

Pursuant to section 564 of the Act, I am establishing the following conditions on this authorization:

Columbia University and Its Authorized Distributor(s)

- A. Columbia University and its authorized distributor(s) will distribute the authorized CII-ArboViroPlex rRT-PCR assay with the authorized labeling only to authorized laboratories. Changes to the authorized labeling may be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- B. Columbia University and its authorized distributor(s) will provide to authorized laboratories the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Healthcare Providers and the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Patients.
- C. Columbia University and its authorized distributor(s) will make available on their websites the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Healthcare Providers and the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Patients.
- D. Columbia University and its authorized distributor(s) will inform authorized laboratories and relevant public health authority(ies) of this EUA, including the terms and conditions herein.
- E. Columbia University and its authorized distributor(s) will ensure that the authorized laboratories using the authorized CII-ArboViroPlex rRT-PCR assay have a process in place for reporting test results to healthcare providers and relevant public health authorities,⁶
- F. Through a process of inventory control, Columbia University and its authorized distributor(s) will maintain records of device usage.
- G. Columbia University and its authorized distributor(s) will collect information on the performance of the test. Columbia University will report to FDA any suspected occurrence of false positive and false negative results and significant deviations from the established performance characteristics of the test of which Columbia University becomes aware.
- H. Columbia University and its authorized distributor(s) are authorized to make available additional information relating to the emergency use of the authorized CII-ArboViroPlex rRT-PCR assay that is consistent with, and does not exceed, the terms of this letter of authorization.

Columbia University

- I. Columbia University will notify FDA of any authorized distributor(s) of the CII-

⁶For questions related to reporting Zika test results to relevant public health authorities, it is recommended that Columbia University, other authorized distributor(s), and authorized laboratories consult with the applicable country, state, or territory health department(s). According to CDC, Zika virus disease is a nationally notifiable condition (see <http://www.cdc.gov/zika/>).

ArboViroPlex rRT-PCR assay, including the name, address, and phone number of any authorized distributor(s).

- J. Columbia University will provide its authorized distributor(s) with a copy of this EUA, and communicate to its authorized distributor(s) any subsequent amendments that might be made to this EUA and its authorized accompanying materials (e.g., Fact Sheets, Instructions for Use).
- K. Columbia University may request changes to the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Healthcare Providers and the authorized CII-ArboViroPlex rRT-PCR assay Fact Sheet for Patients. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- L. Columbia University may request the addition of other instruments for use with the authorized CII-ArboViroPlex rRT-PCR assay. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- M. Columbia University may request the addition of other extraction methods for use with the authorized CII-ArboViroPlex rRT-PCR assay. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- N. Columbia University may request the addition of other specimen types for use with the authorized CII-ArboViroPlex rRT-PCR assay. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- O. Columbia University may request the addition and/or substitution of other control materials for use with the authorized CII-ArboViroPlex rRT-PCR assay. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- P. Columbia University may request the addition and/or substitution of other ancillary reagents and materials for use with the authorized CII-ArboViroPlex rRT-PCR assay. Such requests will be made by Columbia University in consultation with, and require concurrence of, DMD/OIR/CDRH.
- Q. Columbia University will assess traceability⁷ of the CII-ArboViroPlex rRT-PCR assay with FDA-recommended reference material(s). After submission to FDA and DMD/OIR/CDRH's review of and concurrence with the data, Columbia University will update its labeling to reflect the additional testing.
- R. Columbia University will track adverse events and report to FDA under 21 CFR Part 803.

Authorized Laboratories

⁷ *Traceability* refers to tracing analytical sensitivity/reactivity back to a FDA-recommended reference material.

- S. Authorized laboratories will include with reports of the results of the CII-ArboViroPlex rRT-PCR assay the authorized Fact Sheet for Healthcare Providers and the authorized Fact Sheet for Patients. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- T. Authorized laboratories will perform the CII-ArboViroPlex rRT-PCR assay using the NucliSENS easyMAG automated extraction platform (bioMérieux) or with other authorized extraction methods.
- U. Authorized laboratories will perform the CII-ArboViroPlex rRT-PCR assay on the CFX96 Real-Time PCR Detection System (Bio-Rad), or other authorized instruments.
- V. Authorized laboratories will perform the CII-ArboViroPlex rRT-PCR assay for Zika virus, dengue virus, chikungunya virus, and West Nile virus on human serum or other authorized specimen types.
- W. Authorized laboratories will perform the CII-ArboViroPlex rRT-PCR assay for Zika virus on human urine when collected alongside a patient-matched serum specimen and other authorized whole blood derived specimen types.
- X. Authorized laboratories will have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.⁸
- Y. Authorized laboratories will collect information on the performance of the test and report to DMD/OIR/CDRH (via email: CDRH-EUA-Reporting@fda.hhs.gov) and Columbia University any suspected occurrence of false positive or false negative results of which they become aware.
- Z. All laboratory personnel using the test should be appropriately trained in RT-PCR techniques and use appropriate laboratory and personal protective equipment when handling this kit, and use the test in accordance with the authorized labeling.

Columbia University, Its Authorized Distributor(s), and Authorized Laboratories

- AA. Columbia University, its authorized distributor(s), and authorized laboratories, will ensure that any records associated with this EUA are maintained until notified by FDA. Such records will be made available to FDA for inspection upon request.

Conditions Related to Advertising and Promotion

- BB. All advertising and promotional descriptive printed matter relating to the use of the

⁸ For questions related to reporting Zika test results to relevant public health authorities, it is recommended that Columbia University, other authorized distributor(s), and authorized laboratories consult with the applicable country, state, or territory health department(s). According to CDC, Zika virus disease is a nationally notifiable condition (<http://www.cdc.gov/zika/>).

authorized CII-ArboViroPlex rRT-PCR assay shall be consistent with the Fact Sheets and authorized labeling, as well as the terms set forth in this EUA and the applicable requirements set forth in the Act and FDA regulations.

CC. All advertising and promotional descriptive printed matter relating to the use of the authorized CII-ArboViroPlex rRT-PCR assay shall clearly and conspicuously state that:

- This test has not been FDA cleared or approved;
- This test has been authorized by FDA under an EUA for use by authorized laboratories;
- This test has been authorized only for the detection and differentiation of RNA from Zika virus, dengue virus, chikungunya virus, and West Nile virus, not for any other viruses or pathogens; and
- This test is only authorized for the duration of the declaration that circumstances exist justifying the authorization of the emergency use of *in vitro* diagnostic tests for detection of Zika virus and/or diagnosis of Zika virus infection under section 564(b)(1) of the Act, 21 U.S.C. § 360bbb-3(b)(1), unless the authorization is terminated or revoked sooner.

No advertising or promotional descriptive printed matter relating to the use of the authorized CII-ArboViroPlex rRT-PCR assay may represent or suggest that this test is safe or effective for the diagnosis of Zika virus infection.

The emergency use of the authorized CII-ArboViroPlex rRT-PCR assay as described in this letter of authorization must comply with the conditions and all other terms of this authorization.

V. Duration of Authorization

This EUA will be effective until the declaration that circumstances exist justifying the authorization of the emergency use of *in vitro* diagnostic tests for detection of Zika virus and/or diagnosis of Zika virus infection is terminated under section 564(b)(2) of the Act or the EUA is revoked under section 564(g) of the Act.

Sincerely,

Scott Gottlieb, M.D.
Commissioner of Food and Drugs

Enclosures