

ELEVAIR™ Endobronchial Coil System

U.S. Food & Drug Administration
Anesthesiology and Respiratory Therapy Devices Panel
of the Medical Devices Advisory Committee

June 14, 2018

ELEVAIR™ Endobronchial Coil System

Introduction

Julia Anastas, MPH

Vice President, Regulatory Affairs

PneumRx, Inc., a BTG International group company

Unmet Need for Management of Severe Emphysema


- Emphysema is a progressive disease that results in significant debilitation
 - Approximately 3.5 million adults diagnosed in the US^a
 - 1.2 million estimated to have severe emphysema^b
- Patients with emphysema have increasing lack of mobility and limitations in activities of daily living
- The ELEVAIR™ System is intended for patients with significant symptom burden despite optimal medical therapy

^a CDC FastStats: <http://www.cdc.gov/nchs/fastats/default.htm>. Accessed March 8, 2018.


^b Tilert T, et al. *Respir Res.* 2013;14:103; Murphy DE and Panos RJ. *Int J Chron Obstruct Pulmon Dis.* 2013;8:199-208; Mapel DW, et al. *Int J Chron Obstruct Pulmon Dis.* 2011;6:573-581; Hurst JR, et al. *N Engl J Med.* 2010;363(12):1128-1138.

ELEVAIR™ Endobronchial Coil System Overview

- Designed to treat the specific pathophysiology of emphysema
 - Reduce lung hyperinflation, relieve breathlessness, and improve quality of life
- Device consists of Coils and a Delivery System
 - Coils: nitinol shape-memory implants
 - Delivery System: Minimally invasive, bronchoscopic



ELEVAIR™ System Clinical Development Program



Overview of IDE Clinical Program

RENEW

Randomized controlled trial
Treatment vs optimal medical therapy
12-month primary endpoint

RENEW Roll-in

First 2 subjects
at newly trained
clinical sites


Crossover

Eligible RENEW control
subjects receive
Coil treatment


Long-term clinical follow-up

1 2 3 4 5
up to 5 years post-procedure for all Coil-treated subjects

FDA Interactions and PMA Regulatory Timeline



RENEW Protocol Amendment to Lower RV Threshold



Key Conclusions From Clinical Experience

- RENEW met all primary and secondary endpoints based on ITT population
- Greater clinical improvement seen across all endpoints in originally defined population with RV $\geq 225\%$ predicted
- Totality of evidence demonstrates ELEVAIR™ System safety and effectiveness
 - Improves quality of life, exercise capacity, and lung function in patients with severe emphysema (heterogeneous or homogeneous) and severe hyperinflation (RV $\geq 225\%$ predicted)
 - Manageable safety profile through 24 months after implantation
- Data demonstrate a favorable overall benefit-risk profile in the population with severe hyperinflation

Proposed Indications for Use

- The ELEVAIR™ Endobronchial Coil System is indicated for bronchoscopic placement of ELEVAIR Coils in patients with severe emphysema (homogeneous and/or heterogeneous) and severe hyperinflation to improve quality of life, lung function, and exercise capacity

Agenda

**Emphysema
Disease Background**

James Donohue, MD
UNC School of Medicine

RENEW Trial Design

Claire Daugherty, MS
Director, Biostatistics, BTG International, Inc.

**Development and
Effectiveness**

Gerard J. Criner, MD
Founding Chair of Thoracic Medicine and Surgery, Temple University

Safety Profile

David Hahn, MD
Head, PneumRx, Inc.
University of Chicago, Pritzker School of Medicine

Post-Market Plan

Julia Anastas, MPH
Vice President, Regulatory Affairs, PneumRx, Inc.

Patient Preference

A. Brett Hauber, PhD
Vice President, Health Preference Assessment, RTI Health Solutions
University of Washington, School of Pharmacy

Clinical Context

Frank Sciurba, MD
University of Pittsburgh Medical Center

Advisors

Malcolm M. DeCamp, MD

Chair, Cardiothoracic Surgery
School of Medicine and Public Health
University of Wisconsin

Ashley Burns

Director, Clinical Development
PneumRx, Inc.

Scott Berry, PhD

President and Senior Statistical Scientist,
Berry Consultants