


Sponsor Summation

Frank Scieurba, MD


Professor of Medicine and Education
University of Pittsburgh Medical Center

Severe Emphysema Patients With Severe Hyperinflation Have Poor Quality of Life



ELEVAIR™ Endobronchial Coil System Overview

- Designed to treat the specific pathophysiology of emphysema
 - Reduce lung hyperinflation, relieve breathlessness, and improve quality of life
- Device consists of Coils and a Delivery System
 - Coils: nitinol shape-memory implants
 - Delivery System: Minimally invasive, bronchoscopic



Proposed Indications for Use


- The ELEVAIR™ Endobronchial Coil System is indicated for bronchoscopic placement of ELEVAIR Coils in patients with severe emphysema (homogeneous and/or heterogeneous) and severe hyperinflation to improve quality of life, lung function, and exercise capacity

RENEW ITT Baseline Characteristics: Severe and Symptomatic Patient Population

Characteristic ^a	Treatment N=158	Control N=157
Emphysema distribution, % (n)		
Heterogeneous	22.8 (36)	22.9 (36)
Homogeneous	77.2 (122)	77.1 (121)
GOLD status, % (n)		
GOLD 3	24.1 (38)	28.7 (45)
GOLD 4	75.9 (120)	71.3 (112)
FEV ₁ % predicted	25.7 ± 6.3	26.3 ± 6.7
RV % predicted	245.9 ± 39.1	244.5 ± 38.7
6MWT, meters	312.0 ± 79.9	302.7 ± 79.3
mMRC dyspnea scale	2.88 ± 0.74	2.84 ± 0.73
SGRQ, points	60.1 ± 12.8	57.4 ± 14.8
DLCO % predicted	34.1 ± 10.5	34.5 ± 10.7
Number of comorbidities	2.6 ± 2.0	2.3 ± 1.8

^a Mean ±SD, unless otherwise noted.

Effectiveness Outcomes are Consistently Greater in Subjects with RV $\geq 225\%$ Compared to $<225\%$






Effectiveness in both US and Outside US Subjects with RV ≥225%


Endpoint **Between-Group Difference**

▲ OUS, RV ≥225% N=107
● US, RV ≥225% N=128

Poolability test:
All p≥0.35




^a Median difference; ^b Adjusted mean difference.
All endpoints shown are prespecified.



Apparent regional differences in effectiveness were driven by differences in baseline RV by region

Note: Figure was not provided within the PMA; however, underlying information / analysis was included.

Crossover Subject Selection



*N=55 at baseline, N=38 at 12 months

^aAll values reported as non-adjusted medians based on available data at each visit.

Note: Figure was not provided within the PMA; however, underlying information / analysis was included.

Major Complications (MCs) Through 12 Months

RENEW Randomized


Major complication	Subjects, n (%)		P-value ^a
	Treatment N=155	Control N=157	
Total major complication events (95% CI) ^b	54 (34.8) (27.4, 42.9)	30 (19.1) (13.3, 26.1)	0.0021
Death	10 (6.5)	8 (5.1)	NS
Pneumothorax	1 (0.6)	1 (0.6)	NS
Hemoptysis requiring intervention	2 (1.3)	0	NS
COPD exacerbation	18 (11.6)	13 (8.3)	NS
Lower respiratory infections	29 (18.7)	7 (4.5)	<0.0001
Respiratory failure	6 (3.9)	6 (3.8)	NS
Unanticipated bronchoscopy	0	0	N/A

^a Nominal p value By Fisher's exact test unadjusted for multiplicity.

^b By Clopper-Pearson method

Most Frequent Serious Adverse Events Though 12 Months

RENEW Randomized



Pneumonia, terms used: Pneumonia, Pneumonia Bacterial, Pneumonia Staphylococcal, Bronchopulmonary Aspergillosis, Pneumonia Necrotizing, Pneumonia Respiratory Syncytial Virus, Lower Respiratory Tract Infection.

Hemoptysis/Hemorrhage, terms used: Hemoptysis, Post Procedural Hemorrhage, Procedural Hemorrhage, Pulmonary Hemorrhage, Respiratory Tract Hemorrhage.

Note: Figure was not provided within the PMA; however, underlying information / analysis was included.

Post-Market Plan Elements

- Launch based on “Centers of Excellence” model
 - Treatment centers: offer therapy
 - Model treatment centers: offer therapy and host training program
- Comprehensive physician training program
 - Adapted based on post-market plan feedback
- Focused US post-approval study to collect additional safety and effectiveness data
- Ongoing European clinical study program and RENEW to collect long-term safety data
- Ongoing post-market surveillance program

The Study Suggests There is a Group of Patients in the Target Population for Whom the Benefit of a Coil Treatment Likely Outweighs the Risk when Compared to Maximum Medical Therapy

Full Survey Sample



- Benefit > Risks (32%)
- Risks > Benefits (68%)

Survey Sample With RV $\geq 225\%$



- Benefit > Risks (51%)
- Risks > Benefits (49%)

Patient preference for the additional risk of pneumonia requiring hospitalization was extrapolated from 15% to 17.5% (full sample) or 17.3% (RV $\geq 225\%$ predicted).

Benefit - Risk Profile Assessment

Benefits

- Meaningful improvement in
- Lung function
 - Quality of life
 - Exercise capacity



Risks

- Increased risk of
- COPD exacerbation
 - Pneumonia
 - Pneumothorax
 - Bleeding

The ELEVAIR™ System Provides a Meaningful Clinical Tool to Help Severe Emphysema Patients Who Are Seeking Options

