

## **Section 4 – 510(k) Summary**

### **Submitter's Details**

Name Diagnostic Grifols S.A.  
Address Passeig Fluvial, 24, Parets del Valles  
Barcelona, 08150, Spain  
Establishment Registration Number: 3002772505  
Contact Person: Joaquín Alberto Tamparillas  
Technical Director at Diagnostic Industrial Group of Grifols, S.A.  
[joaquin.alberto@grifols.com](mailto:joaquin.alberto@grifols.com)  
Phone: (34) 670-924-632  
Date of Summary: September 5<sup>th</sup>, 2018

### **Name of Device:**

Trade Name: Erytra<sup>®</sup>  
Classification Name: Automated Blood Grouping and Antibody Test System  
Device Class: II  
Product Code: KSZ  
Regulation Number: 21 CFR 864.9175

### **Identification of the Legally Marketed Device (Predicated Device):**

Trade Name: DG Gel 8 Neutral  
Classification Name: Automated Blood Grouping and Antibody Test System  
510(k) Number: BK160100  
Device Class: II  
Product Code: KSZ  
Regulation Number: 21 CFR 864.9175  
Clearance Letter: November 9, 2016

### **Device Description:**

DG Gel 8 Neutral card consists of a plastic support with eight microtubes containing gel in a buffered medium. The microtubes, without antibodies, are used in techniques where antibodies react directly with the red blood cells.

### **Indications for Use:**

The DG Gel 8 Neutral card is for the detection of antibodies to red blood cell antigens of human blood samples and for its use as a control microtube.

For use with the DG Gel System.

For in vitro diagnostic use.

**Comparison to Predicate Device:**

<b>Parameter</b>	<b>Predicate Device DG Gel 8 Neutral (BK160100)</b>	<b>Subject Device DG Gel 8 Neutral</b>
<b>Intended Use</b>	The DG Gel 8 Neutral card is for the detection of antibodies to red blood cell antigens of human blood samples and for its use as a control microtube. For use with the DG Gel System. For in vitro diagnostic use.	The DG Gel 8 Neutral card is for the detection of antibodies to red blood cell antigens of human blood samples and for its use as a control microtube. For use with the DG Gel System. For in vitro diagnostic use.
<b>Regulatory Class</b>	II	II
<b>Product Code</b>	KSZ	KSZ
<b>Regulation number</b>	21 CFR 864.9175	21 CFR 864.9175
<b>Common Name</b>	Automated Blood Grouping and Antibody Test Systems	Automated Blood Grouping and Antibody Test Systems
<b>Principle of the test</b>	Immunoematology gel technique	Immunoematology gel technique
<b>Techniques</b>	Saline techniques: - Determination of the ABO reverse group - ABO compatibility test (immediate spin crossmatch test) - Antibody detection (at cold, room temperature and 37°C) - Use as a Control microtube  Enzymatic techniques: - Screening and/or identification of unexpected antibodies.	Saline techniques: - Determination of the ABO reverse group - ABO compatibility test (immediate spin crossmatch test) - Antibody detection (at cold, room temperature and 37°C) - Use as a Control microtube  Enzymatic techniques: - Screening and/or identification of unexpected antibodies.
<b>Results</b>	Report results as an agglutination grade, absence of agglutination or hemolysis.	Report results as an agglutination grade, absence of agglutination or hemolysis.
<b>Storage conditions</b>	2-25°C	2-25°C
<b>Reagent preparation</b>	Ready to use.	Ready to use.
<b>Number of microtubes</b>	8	8
<b>Instruments for automated method</b>	Erytra or WADiana Compact.	Erytra Eflexis, Erytra or WADiana Compact.

**Performance:**

All required performance tests have been conducted with acceptable results. All performance studies have demonstrated the device is safe and effective when used with Erytra Eflexis.

**Conclusions:**

Diagnostic Grifols S.A. concludes, based on all information submitted and discussed in this submission and in this summary, DG Gel 8 Neutral when used with Erytra Eflexis for the defined indications for use performs as well as or better than the legally marketed predicate device DG Gel 8 Neutral (BK160100) when used with Erytra.