



## 510(k) Summary

This summary of 510(k) safety and effectiveness information is being submitted in accordance to requirements of SMDA 1990 and 21 CFR §807.92.

### Submitter's Details

Name: Bio-Rad Medical Diagnostics GmbH  
Address: Industriestrasse 1, 63303 Dreieich, Germany  
Establishment Registration Number: 9610824

Contact: Dr. Andreas Grosse  
Telephone number: +49 (0) 6103 3130-580  
Fax number: +49 (0) 6103 3130-646  
Date Summary Prepared: April 10, 2019

### Name of Device

Trade or Proprietary Name: TANGO infinity  
Common Name: Automated Blood Grouping and Antibody Test System  
Classification Name: Automated Blood Grouping and Antibody Test System  
510(k) number: BK190312  
Device Class: II  
Product Code: KSZ  
Regulation number: 21 CFR 864.9175

### Identification of the Legally Marketed Device (Predicated Device)

Trade or Proprietary Name: TANGO optimo  
Common Name: Automated Blood Bank System  
Classification Name: Automated Blood Grouping and Antibody Test System  
510(k) number: BK080013  
Device Class: II  
Product Code: KSZ  
Regulation number: 21 CFR 864.9175  
Clearance Letter: July 23, 2008



## DESCRIPTION OF THE DEVICE

TANGO infinity is an automated instrument intended for the *in vitro* serological analysis of human blood specimens utilizing hemagglutination and solid phase methodology. The analyzer is capable of automating the testing process from pipetting sequences of liquids (including diluting samples) through automated and standardized analysis of the end reaction. Process control features on the TANGO infinity include positive sample identification, automated recording of reagents (lot number, expiration date, onboard storage time tracking, sample dilutions, reagent addition and mixing, incubation (room temperature and 37°C), washing, centrifugation, re-suspension, and image analysis. Orders can be generated manually at the TANGO infinity or can be imported by a Laboratory Information System (LIS) interface.

This 510(k) summary applies to the addition of AHG crossmatch testing on TANGO infinity.

## DEVICE COMPARISON TABLE

**Table 1: Substantial Equivalence Comparison**

Parameter	Predicate Device <b>Bio-Rad TANGO optimo</b>	Subject Device <b>Bio-Rad TANGO infinity</b>
Indications for Use Statement	<p>Intended for use in ABO, Rh(D) (including weak D testing) and Rh Phenotype and Kell blood grouping for patient and donor testing, including on receipt, independent confirmation of donor unit labeling. The TANGO optimo Automated Blood Bank Analyzer System uses Anti-Human Globulin Anti-IgG Solidscreen II to detect human red blood cells sensitized with IgG alloantibodies and/or autoantibodies to include tests for antibody screening, antibody identification and the Direct Antiglobulin Test. Also, Blood Grouping Reagent Solidscreen II Anti-D is used to detect weak D and partial D (DVI and DVII) antigens.</p> <p>The TANGO optimo automated system utilizes Blood Grouping Reagents and Anti-Human Globulin Reagents specifically formulated and manufactured by or for Bio-Rad Medical</p>	<p>The TANGO infinity automated blood grouping and antibody test system is an automated instrument intended for the <i>in vitro</i> serological analysis of human blood specimens. In the USA, TANGO infinity is "Rx only" and operates with TANGO infinity Software Version 1.5. It generates results from individual images that must be verified by visual inspection by a qualified person called a validator. In rare cases, where the software is unable to decide, manual editing may be necessary.</p> <p>TANGO infinity utilizes Erytype S agglutination methodology for ABO (forward and reverse blood grouping), Rh(D), Rh Phenotype and Kell blood grouping. TANGO infinity utilizes Solidscreen II solidphase methodology for weak D and partial D testing, antibody screening and identification of red cell alloantibodies, crossmatch, auto control and direct</p>



Parameter	Predicate Device Bio-Rad TANGO optimo	Subject Device Bio-Rad TANGO infinity
	Diagnostics GmbH, and licensed and cleared exclusively for this system.	antiglobulin testing. The TANGO infinity uses reagents and microwell strips approved for the TANGO infinity. Operators of TANGO infinity must be trained by authorized personnel and must be familiar with the content of this TANGO infinity User Manual. The TANGO infinity User Manual is designed to serve as a reference manual for operations and troubleshooting, to assist writing institutional SOPs and to assist training.
Classification	II	Same
510(k) Number	BK080013	BK190312
Product Code	KSZ	Same
Regulation Number	21 CFR 864.9175	Same
Common name	Automated Blood Grouping and Antibody Test System	Same
Tests Performed	<ul style="list-style-type: none"> <li>• Blood Group and Rh(D)</li> <li>• Antigen typing</li> <li>• Antibody Screening</li> <li>• Antibody Identification</li> <li>• Direct Antiglobulin Test</li> <li>• AHG crossmatch</li> </ul>	<p>Same</p> <ul style="list-style-type: none"> <li>• Auto control</li> </ul>
Primary components	<ul style="list-style-type: none"> <li>• Analyzer</li> <li>• Computer</li> <li>• Software</li> <li>• Printer</li> </ul>	Same
Specimen Types	Plasma, Serum and Red Blood Cells	Same
Capability to process STAT samples	Yes	Same
QC procedures implemented	Yes	Same
Barcode Reading	Sample Identification, Reagent Lot No. and Expiration Date	Same
Manual Entry of Sample IDs or Reagent Data	Yes	Same
Sample Loading Random Access	Yes	Same



Parameter	Predicate Device Bio-Rad TANGO optimo	Subject Device Bio-Rad TANGO infinity
Sample loading capacity	144 tubes simultaneously, 12 racks with 12 samples each	120 tubes simultaneously, 10 racks with 12 samples each
Sample Barcode Type	<ul style="list-style-type: none"> <li>• Coda-bar</li> <li>• Code 39</li> <li>• Interleaved 2/5</li> <li>• EAN</li> <li>• (isbt) Code 128 (Subset C)</li> </ul>	Same
Reagent loading capacity	18 positions <ul style="list-style-type: none"> <li>• 3 racks with 4 positions for Reagent Red Blood Cells</li> <li>• 2 Racks with 3 positions for Diluents)</li> </ul>	27 positions <ul style="list-style-type: none"> <li>• 3 racks with 5 positions for Reagent Red Blood Cells</li> <li>• 3 Racks with 4 positions for Diluents</li> </ul>
Reagent Barcode type	Code 128 (Subset C)	Same
Reagent Red Cell Suspension	Yes	Same
Plate loading capacity	10 Plates with 12 strips each	Same
Plate and Strip Transport System	Yes	Same
Pipettor	2 independent pipettor arms	1 pipettor arm
Aspirate Verification	For the whole strip before centrifuge step	After each aspirate step
Strip linear shaking system	Yes	Same
Incubator	Two independent temperature areas 37°C and room temperature	Same
Incubator Temperature Tolerance	Incubator temperature specification at 37 °C is $\pm 2^{\circ}\text{C}$	Same
Incubator capacity	16 positions <ul style="list-style-type: none"> <li>• 5 position for RT</li> <li>• 10 postions for 37 °C</li> <li>• 1 position for balance strip</li> </ul>	Same
Centrifuge	1 Centrifuge; 4 strips capacity	Same
Centrifuge speed	The centrifuge speed is 30 to 900 g depending on the requirements of the test assay	Same
Centrifugation time	The centrifugation time depends on the requirements of the test assay	Same
Orbital shaker	4 strip capacity	Same
Orbital shaking frequency	The orbital shaking frequency is depending on the requirements of the test assay	Same



Parameter	Predicate Device Bio-Rad TANGO optimo	Subject Device Bio-Rad TANGO infinity
Results reading	Yes. CCD camera and Digital Image Analysis.	Same
Test interpretation	Yes. According to predefined rules stated in a definition file.	Same
System solutions and waste containers	<ul style="list-style-type: none"> <li>• System liquid: DI water</li> <li>• Wash buffer</li> <li>• Liquid waste container</li> <li>• Strip waste using conveyer belt</li> </ul>	<ul style="list-style-type: none"> <li>• System liquid: isotonic saline</li> <li>• Wash buffer</li> <li>• Liquid waste container</li> <li>• Strip waste using a bin</li> </ul>
Total speed	12 samples (ABO/ABS 3) per hour, including forward & reverse group	Same
Interfaces	Bidirectional with Laboratory Information System	Same
Useful life	5 years minimum considering a normal function of 230 tests/day and 300 days/year	Same
Computer	TANGO optimo is supplied with a general purpose computer meeting the required specifications can be used	TANGO infinity is supplied with a general purpose computer meeting the required specifications can be used
Operating system	Window XP, Service Pack 3	Windows 7
IH-Com Connectivity	No	Yes

## PERFORMANCE DATA PERCENT AGREEMENT TABLES

Performance Characteristics for crossmatch (Indirect Antiglobulin Test) on the TANGO infinity Analyzer

A multi-center clinical trial, which included testing at two different US clinical sites and an internal site, was conducted to evaluate the performance of Solidscreen II for AHG crossmatch testing on the TANGO infinity Analyzer. The clinical trial included testing of patient and donor samples. The positive and negative percent agreements were calculated for the Anti-Human Globulin Anti-IgG Solidscreen II for AHG crossmatch testing in comparison to the FDA licensed reference methods. Additional internal studies have been performed with well-characterized samples to evaluate the performance of the AHG crossmatch testing on the TANGO infinity Analyzer.

Results of the positive percent agreement and negative percent agreement, with the one-sided Exact 95% Lower Confidence Limit (LCL) are listed in the data table below.



Test	Results from Clinical Trials			
	Negative Agreement		Positive Agreement	
	N	Point Estimate (one-sided Exact 95% LCL)	N	Point Estimate (one-sided Exact 95% LCL)
IAT Crossmatch	539	98.70% (97.57%)	449	99.33% (98.28%)

Test	Results from In-House Study with well-characterized samples			
	Negative Agreement		Positive Agreement	
	N	Point Estimate (one-sided Exact 95% LCL)	N	Point Estimate (one-sided Exact 95% LCL)
IAT Crossmatch	344	100% (99.13%)	320	100% (99.07%)

Agreement between the methods does not imply which method obtained the correct result. The above results do not reflect any discrepancy resolution between the methods.

## CONCLUSION

In support of this premarket notification a multi-center clinical study at three external sites and at Bio-Rad was performed to validate the use of Anti-Human Globulin Anti-IgG Solidscreen II on TANGO infinity for AHG crossmatching. The studies included a number of 1590 crossmatches.

The results obtained in the clinical trial support the conclusion that the testing of Anti-Human Globulin Testing using the Solidscreen II method on the TANGO infinity Analyzer is safe and effective. The yielded results demonstrate that end users, with proper training, could use the TANGO infinity Analyzer to perform crossmatching and that the testing with the specified reagents on the TANGO infinity Analyzer does generate results comparable to established FDA licensed reference reagents and FDA approved predicates.