

## 510(k) Summary

**Date:** September 30, 2020

**Submitter:**

Haemonetics Corporation  
125 Summer Street  
Boston, MA 02110

**Contact:**

Julie A. Ryan  
Sr. Director, Regulatory Affairs  
Phone: 781-690-0575  
Fax: 781-356-3558  
Email: julie.ryan@haemonetics.com

**Device Information:**

*Trade Name:* NexSys PCS® Plasma Collection System with Persona™ Technology  
*Common Name:* Automated Blood Cell Separator  
*Classification Name:* Separator, Automated, Blood Cell, Diagnostic  
*Regulation Number:* 21 CFR 864.9245  
*Review Panel:* Hematology  
*Product Code:* GKT  
*Device Class:* 2

**Device Characteristics Summary:**

The Haemonetics NexSys PCS® Plasma Collection System with Persona™ Technology includes the Persona Plasma Pooling Bottle and ability to integrate with NexLynk DMS® Donor Management System. It is designed for separation of whole blood by centrifugation, collection of plasma, and return of the remaining components to the donor.

Proprietary Persona Technology incorporates the new Persona nomogram which is tailored to each donor's individual characteristics. When Persona is enabled, the target plasma volume is determined using each donor's height, weight, and hematocrit. The plasma collected by the NexSys PCS may be designated for use in therapeutic transfusion or be conserved, used as source plasma, and subsequently fractionated into plasma-derived products.

## Indications for Use:

The NexSys PCS® Plasma Collection System with Persona™ Technology is intended for use as an automated cell separator system and blood component collector in conjunction with single-use sterile disposable sets, with or without saline compensation.

Products that can be collected using the NexSys PCS® system include source plasma and plasma for transfusion.

## Non-Clinical Testing Summary:

The following non-clinical performance testing was submitted in support of a determination of substantial equivalence between the subject and predicate device. A summary of the performance testing is presented below in Table 1.

Test data demonstrates that the device met all performance requirements, and that the subject device is as safe, as effective, and performs as well as or better than the predicate device.

**Table 1: Summary of Performance Studies**

Test Name	Test Report #	Test Intent	Test Result
Software Verification	TR-SOF-100731	To verify that NexSys PCS software functions as intended and meets all design requirements.	Passed

## Clinical Testing Summary:

Clinical testing was submitted in support of a determination of substantial equivalence between the subject and predicate device. A summary of the clinical testing is presented below in Table 2.

Haemonetics conducted a prospective, double-blinded, randomized, controlled, multicenter clinical trial to demonstrate the safety and effectiveness of the NexSys PCS® Plasma Collection System with Persona™ Technology. The data provide evidence that plasmapheresis with the NexSys PCS with Persona Technology is non-inferior to the control with regard to the incidence rates of significant hypotensive adverse events. Secondary analysis with regard to total plasma volume collected per procedure demonstrated that more plasma volume per procedure was collected using the NexSys PCS with Persona Technology.

Clinical data demonstrates that the study met the primary endpoint and that the NexSys PCS Plasma Collection System with Persona Technology is substantially equivalent to the predicate device.

**Table 2: Summary of Clinical Studies**

<b>Test Name</b>	<b>Test Report #</b>	<b>Test Intent</b>	<b>Test Result</b>
IMPACT Clinical Trial	TR-CLN-100467	To demonstrate safety and effectiveness of the NexSys PCS Plasma Collection System with Persona Technology.	The primary endpoint was met.

**Comparison to Predicate:**

The Haemonetics NexSys® PCS Plasma Collection System with Persona™ Technology is substantially equivalent to the Haemonetics NexSys PCS (PCS 300) Plasma Collection System with YES® Technology cleared under BK180185. The NexSys PCS is intended for use in the same operating environment with the same donor/operator population as the predicate device. The indications for use are the same. The manner in which the software protocol operates to process blood and collect plasma is the same. The technological characteristics of the subject device differ from the predicate only in features of the embedded software that enable use of the Persona nomogram and do not impact the clinical functionality of the device. These differences do not render the device non-substantially equivalent because clinical and non-clinical testing has demonstrated that the subject device is as safe and effective as the predicate and the results of verification and validation have not raised different questions of safety and effectiveness than the predicate.

A summary comparison is presented below in Table 3.

**Table 3: Comparison of the NexSys PCS® Plasma Collection System with Persona™ Technology to the Predicate**

	<b>Predicate</b> <b>NexSys PCS® (PCS 300) Plasma Collection System w/ YES® Technology (BK180185)</b>	<b>Subject</b> <b>NexSys PCS® Plasma Collection System with Persona™ Technology</b>
<b>Manufacturer</b>	Haemonetics Corporation	Same
<b>Trade Name</b>	NexSys PCS® (PCS 300) Plasma Collection System w/ YES™ Technology	NexSys PCS® Plasma Collection System w/ Persona™ Technology
<b>Common Name</b>	Automated Blood Cell Separator	Same
<b>Classification Name</b>	Separator, Automated, Blood Cell, Diagnostic	Same
<b>Regulation Number</b>	21 CFR 864.9245	Same
<b>Review Panel</b>	Hematology	Same
<b>Product Code</b>	GKT	Same
<b>Device Class</b>	2	Same
<b>Indications for Use</b>	<p>The PCS 300 Plasma Collection System is intended for use as an automated cell separator system and blood component collector in conjunction with single-use sterile disposable sets, with or without saline compensation.</p> <p>Products that can be collected using the PCS 300 system include source plasma and plasma for transfusion.</p>	<p>The NexSys PCS® Plasma Collection System with Persona™ Technology is intended for use as an automated cell separator system and blood component collector in conjunction with single-use sterile disposable sets, with or without saline compensation.</p> <p>Products that can be collected using the NexSys PCS® system include source plasma and plasma for transfusion.</p>

# HAEMONETICS®

Hardware		
<b>Pumps</b>	Peristaltic pumps, 1 ml per rotation	Same
<b>Effluent Line Sensor</b>	Absorbance optical system (LED beam across transparent tubing) for detection of air/plasma interface and plasma/buffy coat interface	Same
<b>Air Detectors</b>	Ultrasonic	Same
<b>Pressure Sensor (DPM)</b>	Donor Pressure Monitor with interlock to regulate pump speed based on pressure	Same
<b>Wireless Connectivity</b>	Yes	Same
<b>Centrifuge</b>	Nominal speed = 7500 rpm	Same
<b>Bowl Optics</b>	Absorbance optical system for detection of air/plasma interface and plasma/buffy coat interface in the separation bowl	Same
<b>Valves</b>	Pneumatic valves	Same
<b>Plasma Weigher</b>	Fixed front load cell	Same
<b>User Interface</b>	8" color touch screen	Same
<b>Bar Code Reader</b>	Embedded; used for operator, donor, donation, disposable set readings	Embedded or external USB; used for operator, donor, donation, disposable set readings
<b>Donor Display</b>	Digital display on each side of the device, communicates info to donor about the procedure	Same
<b>Anticoagulant (AC) Weigher</b>	Load cell on pole with hook for hanging the AC bag	Same
<b>Status Beacon</b>	Beacon light above touch screen display, indicates status of procedure	Same
Software		
<b>Self-Test</b>	Yes	Same
<b>Plasma Target Selection</b>	Yes, manual and through server if connected	Same, with addition of plasma target programming based on Persona™ nomogram
<b>Modifiable Parameters</b>	Yes, cuff pressure, draw and return speed, max plasma per cycle, saline	Same

# HAEMONETICS®

<b>Express Donor Draw and Return Flow Control</b>	Yes	Same
<b>AC Short Prime</b>	Yes	Same
<b>Disposable Detection</b>	Detection of the installed disposables: DPM, line sensor tubing, and plasma container; detection of line sensor cover and centrifuge cover lock; disposables bar codes can also be scanned	Same
<b>Diagnostics</b>	Manual and automated diagnostics	Same
<b>Notifications</b>	Main and hints screens, individual ID for each notice	Same
<b>Procedure Technical Data</b>	Records data for up to 100 procedures	Records data for up to 10,000 procedures
<b>Phlebotomy Workflow</b>	Yes	Same
<b>User Access Control</b>	Yes	Same
<b>Disposable Sets</b>		
<b>Disposables</b>	Previously-cleared disposable bowls, bottles, and harnesses	Same, with addition of a larger Persona™ plasma pooling bottle

\_\_\_\_\_  
 Julie A. Ryan  
 Sr. Director, Regulatory Affairs  
 Haemonetics Corporation

\_\_\_\_\_  
 Date