510k Number

BK200520

A. Introduction

In accordance with the requirements of 21CFR807.92, the following information provides sufficient detail to understand the basis for a determination of substantial equivalence.

I. Submitter

EKF-diagnostic GmbH Ebendorfer Chaussee 3 Barleben, Germany 39179

Primary Contact: Karen Golomb

Phone: +1 (830) 249-0772 Fax: +1 (830) 249-0851

Email: karengolomb@ekfdiagnostics.com

Date prepared: October 21, 2020

II. Device Name

Proprietary names: DiaSpect Tm

DiaSpect Tm Cuvettes

Common names: As above

Classification: 21 CFR 864.5620 – Automated hemoglobin system

Class II Hematology

Product Code: GKR

III. Predicate Device

The DiaSpect Tm system is substantially equivalent to the currently marketed HemoCue Hb 301 System (BK060048, K061047).

IV. Device Description

The DiaSpect Tm system consists of an analyzer and cuvettes. The DiaSpect Tm analyzer is a spectrophotometric instrument for the measurement of total hemoglobin concentration in unaltered human blood. The DiaSpect Tm Cuvette is injection-molded of poly methyl methacrylate (PMMA) and contains a cavity of $10~\mu L$ volume. The cavity is empty.

V. Indications for Use

The DiaSpect Tm system is intended for the *in vitro* quantitative measurement of total hemoglobin in non-anticoagulated capillary whole blood and venous whole blood drawn in K2-EDTA or lithium heparin tubes in point-of-care settings and in non-anticoagulated capillary whole blood and venous whole blood drawn in K2-EDTA tubes in blood bank settings. The DiaSpect Tm system consists of the DiaSpect Tm analyzer and specifically designed disposable cuvettes. The DiaSpect Tm analyzer is only to be used with DiaSpect Tm Cuvettes. Rx only.

VI. Comparison with Predicate

Predicate device name(s): HemoCue Hb 301 System

Predicate 510(k) number(s): BK060048, K061047

Similarities compared to the predicate device (BK060048,K061047)

Performance	Predicate Device HemoCue Hb 301 System (BK060048, K061047)	Candidate Device DiaSpect Tm system (K172173)
Indications for use	Determine hemoglobin content of whole blood	Same
Analyte	Hemoglobin	Same
Sample preparation	None	Same
Automation	Fully automated assay	Same
Calibration procedure	Factory calibrated	Same
Test methodology	Optical absorbance	Same
Sample volume	10 μL	Same
Cuvette reagent components	None	Same
Result memory	No	Same
Quality Control	Quality control materials are available	Same

Differences compared to the predicate device (BK060048, K061047)

100	Predicate Device	Candidate Device
Performance	HemoCue Hb 301 System	DiaSpect Tm system
	(BK060048, K061047)	(K172173)
Intended Use	The HemoCue Hb 301 System is	The DiaSpect Tm system is intended for
	designed for quantitative point-of-care	the <i>in vitro</i> quantitative measurement of
	whole blood hemoglobin determination	total hemoglobin in non-anticoagulated
	in primary care or blood donation	capillary whole blood and venous whole
	settings using a specially designed	blood drawn in K2-EDTA or lithium
	analyzer, the HemoCue Hb 301	heparin tubes in point-of-care settings and
	Analyzer, and specially designed	in non-anticoagulated capillary whole
	microcuvettes, the HemoCue Hb 301	blood and venous whole blood drawn in
	Microcuvettes. The HemoCue Hb 301	K2-EDTA tubes in blood bank settings.
	System is for In Vitro Diagnostic use	The DiaSpect Tm system consists of the
	only. The HemoCue Hb 301 Analyzer	DiaSpect Tm analyzer and specifically
	is only to be used with HemoCue Hb	designed disposable cuvettes. The
	301 Microcuvettes.	DiaSpect Tm analyzer is only to be used
		with DiaSpect Tm Cuvettes.
Comple type	Conillany outonial on you are whole	Rx only.
Sample type	Capillary, arterial or venous whole blood	Capillary or venous whole blood
Built-In Quality	Auto self-check at start and at regular	Auto self-check confirms each
Control	time intervals	measurement
Expected values		
Adult Males	13.0 - 18.0 g/dL	13.0 - 17.0 g/dL
Adult Females	11.0 - 16.0 g/dL	12.0 – 15.0 g/dL
Children	2 yrs to teenage: 11.0 – 16.0 g/dL	>2 yrs to 21: 11.0 – 15.5 g/dL
Infants	Post-natal: 10.0 – 14.0 g/dL	1 mo. to 2 yrs: 9.4 – 16.5 g/dL
Operating	10 to 40°C (50 to 104°F)	10 to 42°C (50 to 107°F)
temperature		
Cuvette storage	10 – 40°C (50 to 104 °C)	0 – 50°C (32 to 122 °F)
Cuvette	Composed of polystyrene	Composed of poly methyl methacrylate
composition		(PMMA)

Performance	Predicate Device HemoCue Hb 301 System (BK060048, K061047)	Candidate Device DiaSpect Tm system (K172173)
Control Kit	Two concentration levels of liquid	Three concentration levels of liquid
components	controls	controls
Measurement Range	0.0 - 25.6 g/dL	1.2 – 25.5 g/dL
Connectivity	Serial port	Wireless Bluetooth Low EnergyUSB
Measuring Time	3 seconds	1 second
User Interface	Display	• Display
	Beeper	Beeper
	One button	Cuvette holder
Microcuvette	Place on a tray	Slot in
Insertion Technique		
Power Sources	AC adapter	• 3.6 V integrated lithium ion rechargeable
	• 4 x AA batteries	batteries
		USB Adaptor
Dimensions	6.3x5.5x2.8 inch	6x3.5x1.6 inch

VII. Performance Characteristics

1. Analytical Performance

Analytical studies to determine the substantial equivalence of the DiaSpect Tm system in point-of-care settings have been previously cleared under K172173. Precision study summaries from K172173 are included as follows:

a. Precision/Reproducibility:

20-Day Precision

Reproducibility was conducted at three intended use sites over 20 operating days using three DiaSpect Cuvette lots (one lot per site), three DiaSpect Tm Analyzers (one instrument per site), and one lot of DiaSpect Controls at three levels (low, medium and high). Each control set was run in duplicate twice daily for 20 days, by two operators at each site. A total of 160 test results were generated for each control level at each site. SD and %CV for within-run, between-run, between-day, between-operator, and between-site were calculated for each site and all sites combined. Reproducibility results at all test sites were within the defined acceptance criteria.

Table 1 20-Day precision data met the criteria of <7% CV.

Sample	N	Mean	Within- Run SD, %CV	Between- Run SD, %CV	Between- Day SD, %CV	Between- Operator SD, %CV	Between- Site SD, %CV	Total SD, %CV
Low	240	7.99	0.085, 1.06%	0.05, 0.59%	0.04, 0.47%	0, 0%	0.04, 0.45%	0.11, 1.38%
Medium	240	12.58	0.11, 0.8%	0.05, 0.38%	0.03, 0.22%	0, 0%	0.05, 0.47%	0.14, 1.09%
High	240	15.82	0.15, 0.92%	0.06, 0.36%	0.04, 0.27%	0, 0%	0.15, 0.97%	0.22, 1.41%

Single-Day Precision

Single-day precision was performed at a single site using K2-EDTA venous whole blood. Five donors provided five samples (one sample each), some of which were manipulated to increase or decrease hemoglobin levels. The study was performed at a single site using three instruments, three lots of cuvettes, and three operators. Each of three operators ran duplicate tests on all five levels, providing 54 measurements for each level. Single-day precision results were within the defined acceptance criteria.

Table 2 Single-Day precision data met the criteria of <7% CV.

Sample	N	Mean	Within-Run (SD, %CV)	Between-Lot (SD, %CV)	Between- Operator (SD, %CV)	Between- Analyzer (SD, %CV)	Total (SD, %CV)
Level 1	54	4.87	(0.05, 0.97%)	(0.04,0.82%)	(0.02, 0.42%)	(0.11, 2.19%)	(0.13, 2.57%)
Level 2	54	10.19	(0.19, 1.87%)	(0, 0%)	(0.06, .06%)	(0.12, 1.17%)	(0.23, 2.29%)
Level 3	54	13.75	(0.38, 2.77%)	(0.07, 0.48%)	(0.21, 1.5%)	(0,0%)	(0.44, 3.19%)
Level 4	54	17.46	(0.69, 3.93%)	(0,0%)	(0.08, 0.48%)	(0.37, 2.14%)	(0.79, 4.5%)
Level 5	54	22.93	(1.45, 6.33%)	(0,0%)	(0.3, 1.32%)	(0,0%)	(1.48, 6.47%)

b. Linearity

Please refer to submission K172173.

c. Traceability, Stability, Expected values (Controls, Calibrators)

The DiaSpect Control HBT is produced in three concentrations that correspond to three levels of human hemoglobin.

d. Detection Limit:

Please refer to submission K172173.

e. Analytical specificity:

Please refer to submission K172173.

f. Assay cut-off:

Not applicable

g. Stability

Please refer to submission K172173.

h. Calibration

The DiaSpect Tm analyzer is factory calibrated and requires no further calibration.

i. Quality Control

DiaSpect control solutions are available to facilitate compliance with local, state and/or federal regulations or accreditation requirements. The DiaSpect Control HBT is produced in three concentrations that correspond to three known levels of human hemoglobin. Each vial contains 1.9 mL of a solution of a red dye (Rhodamine) in purified water. The reagent does not contain any material of human or animal origin. The DiaSpect Control HBT is 510(k) exempt.

2. Method Comparison Studies

Method comparison studies were conducted with finger stick capillary whole blood and paired K2-EDTA venous whole blood at blood bank settings. A total of 150 capillary and 147 venous blood donors were tested with the DiaSpect Tm and HemoCue Hb 301 at two sites.

DiaSpect Tm vs HemoCue Hb 301 (linear regression)

Site	Matrix	N	Slope (95% CI)	Intercept (95% CI)	r
1	Capillary	75	0.9644 (0.8911 ~ 1.0377)	0.1086 (-1.011 ~ 1.228)	0.951
1	Venous	72	1.1521 (1.0967 ~ 1.2074)	-2.8551 (-3.7028 ~ -2.0074)	0.980
2	Capillary	75	0.9472 (0.8988 ~ 0.9957)	0.3588 (-0.3555 ~ 1.0731)	0.977
2	Venous	75	1.1784 (1.1341 ~ 1.2227)	-2.7653 (-3.4184 ~ -2.1121)	0.987
Total	Capillary	150	0.9541 (0.9134 ~ 0.9948)	0.2619 (-0.3496 ~ 0.8733)	0.967
Total	Venous	147	1.1465 (1.1068 ~ 1.1863)	-2.5299 (-3.4183 ~ -2.1122)	0.978

DiaSpect Tm vs HemoCue Hb 301 (Passing-Bablok regression)

Site #	Matrix	N	Slope (95% CI)	Intercept (95% CI)	r
ALL	Capillary	150	1.000 (0.9512 ~ 1.032)	-0.45 (-1.011 ~ 0.2805)	0.967
ALL	Venous	147	1.140 (1.100 ~ 1.179)	-2.477 (-3.067 ~ -1.895	0.978

The bias of the DiaSpect Tm at hemoglobin acceptance levels for donors was calculated from the regression curves using the Hb 301 as the reference method.

Daman Critaff aug III. 201	DiaSpect Tn	n Capillary	Capillary DiaSpect Tm Ver	
Donor Cutoff vs Hb 301	Hb Level	% Bias	Hb Level	% Bias
Adult Female, 12.5 g/dL	12.2 g/dL	-2.4%	11.8	-5.6%
Adult Male, 13 g/dL	12.7 g/dL	-2.3%	12.4	-4.8%

Using the Hb 301 as the reference method, the agreements for blood donor acceptance at these cutoff levels are as follows:

Capillary All Donors	Hb 301 ≥Cutoff	Hb 301 < Cutoff	
DiaSpect Tm ≥Cutoff	127	2	
DiaSpect Tm < Cutoff	7	14	

Agreement (Score 95% CI)

Accept = 127/134 = 94.8% (89.6~97.4%) Reject = 14/16 = 87.5% (64.0~96.5%)

Venous All Donors	Hb 301 ≥Cutoff	Hb 301 < Cutoff	
DiaSpect Tm ≥Cutoff	120	1	
DiaSpect Tm <cutoff< td=""><td>13</td><td>13</td></cutoff<>	13	13	

Agreement (Score 95% CI)

Accept = 120/133 = 90.2% (84.0~94.2%) Reject = 13/14 = 92.9% (68.5~98.7%)

The negative bias observed for the DiaSpect Tm was resolved by comparing both test methods to a laboratory reference method. In a separate study conducted at one site, 100 capillary and 100 venous samples were tested with the DiaSpect Tm and HemoCue Hb 301 and compared with the same venous samples that were also tested with the Sysmex XP-300. The DiaSpect Tm was found to be substantially equivalent to the Sysmex XP-300 (calibrated to the hemiglobin yanide method, HiCN), while the HemoCue Hb 301 consistently provided results that are higher by approximately 0.6 g/dL. The results are summarized below.

The bias of the DiaSpect Tm at hemoglobin acceptance levels for donors was calculated from the regression curves using the Sysmex XP-300 as the reference method.

Donor Cutoff vs XP-300	DiaSpect Tn	n Capillary	DiaSpect T	Tm Venous
Donor Cuton vs AP-300	Hb Level	% Bias	Hb Level	% Bias
Adult Female, 12.5 g/dL	12.56	0.5%	12.53	0.2%
Adult Male, 13 g/dL	13.05	0.4%	13.07	0.7%

The bias of the Hb 301 at hemoglobin acceptance levels for donors was calculated from the regression curves using the Sysmex XP-300 as the reference method.

Donor Cutoff vs XP-300	Нь 301 С	Capillary	Hb 301	Venous
Donor Cuton vs AF-300	Hb Level	% Bias	Hb Level	% Bias
Adult Female, 12.5 g/dL	13.12	5.0%	13.03	4.2%
Adult Male, 13 g/dL	13.62	4.8%	13.55	4.2%

Using the Sysmex-XP-300 as the reference method, the agreements for blood donor acceptance at these cutoff levels are as follows:

Capillary All Donors	Sysmex XP-300 ≥Cutoff	Sysmex XP-300 < Cutoff
DiaSpect Tm ≥Cutoff	75	3
DiaSpect Tm <cutoff< td=""><td>6</td><td>16</td></cutoff<>	6	16

Agreement (Score 95% CI)

 $Accept = 75/81 = 92.6\% (84.8 \sim 96.6\%)$

Reject = 16/19 = 84.2% (62.4~94.5%)

Venous All Donors	Sysmex XP-300 ≥Cutoff	Sysmex XP-300 < Cutoff
DiaSpect Tm ≥Cutoff	77	0
DiaSpect Tm <cutoff< td=""><td>4</td><td>19</td></cutoff<>	4	19

Agreement (Score 95% CI)

 $Accept = 77/81 = 95.1\% (88.0 \sim 98.1\%)$ Reject = 19/19 = 100.0% (83.2~100.0%)

Matrix Comparison

Please refer to submission K172173.

VII. Conclusion

The DiaSpect Tm data presented and provided is complete and supports the basis for substantial equivalence to the predicate device for use in blood bank settings.