

Individuals using assistive technology may not be able to fully access the information contained in this file. For assistance, please call 800-835-4709 or 240-402-8010, extension 1. CBER Consumer Affairs Branch or send an e-mail to: ocod@fda.hhs.gov and include 508 Accommodation and the title of the document in the subject line of your e-mail.

Anti-Human Globulin Anti-IgG, -C3d Polyspecific (Rabbit/Murine Monoclonal)

MTS™ Anti-IgG, -C3d Card

INSTRUCTIONS FOR USE

REF

MTS084014

Rx ONLY

Intended Use

For *in vitro* diagnostic and laboratory professional use

Qualitative procedure for the detection of IgG or complement bound to red blood cells

For use with the ID-Micro Typing System™

Contains: 6 tests per card

Blood samples from patients and or donors intended for direct antiglobulin testing should be drawn into EDTA to prevent *in vitro* complement binding.

Anti-Human Globulin Anti-IgG, -C3d Polyspecific (Rabbit/Murine Monoclonal) MTS™ Anti-IgG, -C3d Card is used in the investigation of transfusion reactions, autoimmune hemolytic anemia, and hemolytic disease of the newborn in a direct antiglobulin test (DAT).

Anti-Human Globulin Anti-IgG, -C3d Polyspecific (Rabbit/Murine Monoclonal) MTS™ Anti-IgG, -C3d Card is intended for manual, semi-automated or automated column agglutination technology using ORTHO® Workstation, ORTHO Optix™ Reader and the ORTHO VISION® and VISION® Max Analyzers.

Observable Indications

Drying, discoloration, bubbles, crystals, other artifacts, opened or damaged seals may indicate product alteration.

Summary and Explanation of the Test

Anti-Human Globulin was described in 1945 by Coombs, Mourant, and Race.¹ Blood group antibodies of the IgG class, that were previously undetectable, reacted in the direct or indirect antiglobulin test (also known as the Coombs test). Anti-IgG reagents remain important tools for determining the presence or absence of IgG on human red blood cells. The reagent is used in the investigation of hemolytic disease of the newborn (refer to Limitations of Procedure, item 11), transfusion reactions, and autoimmune hemolytic anemia in a direct antiglobulin test (DAT). The DAT detects IgG and/or C3 using either a polyspecific reagent solely or monospecific Anti-IgG and Anti-C3.

Principles of the Procedure

The combination of the antiglobulin reagent incorporated into gel, known as the ID-MTS™ Gel Test² was first described by Dr. Yves Lapierre.³ The MTS™ Anti-IgG, -C3d Card can be used in the direct antiglobulin test.²

Red blood cells that are coated with IgG and/or complement due to *in vivo* sensitization are detected with the direct antiglobulin test.

The MTS™ Anti-IgG, -C3d Card restricts the unbound IgG from moving through the gel during centrifugation. The unbound IgG and complement components do not neutralize the Anti-IgG, -C3d incorporated in the gel.

Red blood cells sensitized with IgG, and/or, C3 react with the corresponding antiglobulin component in the microtube during centrifugation. Strongly positive agglutination reactions produce a line of red blood cells layered at the top of the gel.

Positive reactions will have varying degrees of visible red blood cell agglutinates suspended in the gel. Uncoated (unsensitized) red blood cells are not agglutinated by the Anti-IgG and/or, Anti-C3 and will form a button at the bottom of the microtube.

Reagents

Anti-Human Globulin Anti-IgG, -C3d; Polyspecific (Rabbit/Murine Monoclonal) for the MTS™ Anti-IgG, -C3d Card is prepared from blended pools of sera obtained from rabbits that have been immunized with human IgG. The Anti-C3 antibodies are produced by murine hybridoma cell lines secreting specific complement antibodies which have been grown in mouse

INSTRUCTIONS FOR USE

MTS™ Anti-IgG, -C3d

Storage Requirements

ascites. The sera and ascites are further processed to remove unwanted heterospecific antibodies and are suspended in a buffered gel solution. The reagent meets present potency and specificity requirements of the FDA.

Sodium Azide (0.1% final concentration) is added as a preservative.

Anti-Human Globulin, Anti-IgG,-C3d; Polyspecific (Rabbit/Murine Monoclonal) is suspended in a diluent and buffered gel solution and is contained in the 6 microtubes of the MTS™ Anti-IgG,-C3d Card.

Storage Requirements

Store cards upright at 2–25 °C.

Precautions

- Do not use beyond expiration date.
- Do not freeze or expose cards to excessive heat.
- Use reagents as furnished.

Caution: All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents.

Caution: Sodium azide may be toxic if ingested and may react with lead and copper plumbing to form explosive compounds. If discarded into sink, flush with a large volume of water to prevent azide buildup.

WARNING: Once a gel card is used in testing, it may contain infectious material and should therefore be handled and disposed of as biohazard waste.

- A clear liquid layer should appear on top of the opaque gel in each microtube. Do not use gel cards if the gel matrix is absent or the liquid level in the microtube is at or below the top of the gel matrix. Do not use gel cards that show signs of drying, discoloration, bubbles, crystals, or other artifacts. Do not use cards if foil seals appear damaged or opened.

Note: Refer to the ID-Micro Typing System™ Interpretation Guide⁴ for additional information related to the visual inspection of gel cards before use.

- Do not remove foil seal until ready to use. Foil should be removed immediately before testing or within 1 hour of testing. Once opened, the gel may begin to dry out which could affect test results (refer to Limitations of the Procedure, item 10).
- After removing the foil, visually inspect all gel cards to ensure that residual film does not block the opening of any microtube.

Caution: The pipet tip should not touch the gel card. Erroneous results due to carryover may occur.

- Customers who choose to use commercial antisera in an off-label manner must ensure that the test method is appropriate by validating its intended use.
- Do not use gel cards that have not been shipped in an upright position.

Specimen Collection and Preparation

No special preparation of the patient is required prior to specimen collection. Collect all blood samples using acceptable phlebotomy techniques.

Cell washing is an optional step in specimen preparation. Performance of the device is not compromised by washing cells with saline three times.

Samples for Direct Antiglobulin Test (DAT)

Samples intended for direct antiglobulin testing should be drawn into EDTA to prevent *in vitro* complement binding. Red blood cells should be tested within 24 hours after collection. Some samples such as cord blood, blood stored for extended periods of time, or blood that has been incompletely anticoagulated, may develop fibrin clots or particulates. The fibrin clots or particulates may interfere with the ID-MTS™ Gel Test and cause red blood cell entrapment at the top of the microtube. Testing should be repeated using red blood cells that have been washed to remove the clots or particulates.

Red blood cells that are stored for extended periods of time may become coated *in vitro* with complement and/or globulin proteins. Those samples coated with IgG will then test as DAT positive with this reagent.

Hemolyzed and grossly icteric blood samples may cause difficulty in interpretation, and test results should be used with caution.

INSTRUCTIONS FOR USE

Reagent Preparation

MTS™ Anti-IgG, -C3d

Rouleaux caused by serum or plasma with abnormally high concentrations of protein (such as in patients with multiple myeloma or Waldenstrom's macroglobulinemia or from patients who have received plasma expanders of high molecular weight) may infrequently cause difficulties in the ID-MTS™ Gel Test interpretation.⁴ False positive results or hazy reactions may occur with these samples but are rare. Samples exhibiting rouleaux should be washed several times in saline and retested.⁵ Laboratories are advised to consult their approved procedures.

Reagent Preparation

The MTS™ Anti-IgG,-C3d Card is provided ready to use. Each microtube contains Anti-IgG,-C3d suitable for one test. The gel card is heat-sealed with aluminum foil to preserve the integrity of the reagents. Variations in the liquid and/or gel levels between microtubes may normally be observed. However, do not use cards if the liquid level in the microtube is at or below the top of the gel matrix (refer to Precautions).

Procedure

The procedures identified below are for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures and are advised to consult the appropriate regulatory agencies to determine validation requirements. Refer to ID-Micro Typing System™ Interpretation Guide⁴ and ID-Micro Typing System™ Implementation Guide and Procedures⁶ for additional information.

Materials Provided

Anti-Human Globulin, Anti-IgG,-C3d; Polyspecific (Rabbit/Murine Monoclonal) suspended in a final diluent and buffered gel solution is contained in the 6 microtubes of the MTS™ Anti-IgG,-C3d Card.

Materials Required but Not Provided

For manual gel card processing:

- Quality Control Material known to give the appropriate positive and negative test results for each reagent requiring quality control. Examples include, but are not limited to, Hemo bioscience HBS – C3 Control Cells
- MTS™ Diluent 2
- Pipets: 10 µL and 50 µL
- Pipet Tips
- Test Tubes
- Dispenser pipet capable of delivering 1.0 mL
- Marking Pen
- ORTHO® Workstation
- ORTHO Optix™ Reader

For automated gel card processing with the ORTHO VISION® Analyzer or ORTHO VISION® Max Analyzer:

- Quality Control Material known to give the appropriate positive and negative test results for each reagent requiring quality control
- MTS™ Diluent 2
- ORTHO VISION® Analyzer
- ORTHO VISION® Max Analyzer

Test Procedure

Direct Antiglobulin Test

1. Bring samples and reagents to room temperature (18–25 °C).
2. Visually inspect gel cards before use. Each microtube should have a clear liquid layer on top of opaque gel.

Caution:

Do not use gel cards if the gel matrix is absent or the liquid level in the microtube is at or below the top of the gel matrix. Do not use gel cards that show signs of drying, discoloration, bubbles, crystals, or other artifacts. Do not use cards if foil seals appear damaged or opened.

Note:

Refer to ID-Micro Typing System™ Interpretation Guide⁴ for additional information related to the visual inspection of gel cards before use.

3. Prepare a red blood cell suspension of approximately 0.8% in MTS™ Diluent 2 (e.g., deliver 1.0 mL of MTS™ Diluent 2 into a test tube and pipet 10 µL packed red blood cells into the diluent), mix gently.
4. Label the gel card appropriately.

INSTRUCTIONS FOR USE

MTS™ Anti-IgG, -C3d

Interpretation of Results

- Remove the foil seal from the MTS™ Gel Card or from the individual microtubes to be used for testing. After removing the foil, visually inspect all gel cards to ensure that residual film does not block the opening of any microtube.

Note:

Foil should be removed immediately before testing or within 1 hour of testing. Once opened, the gel may begin to dry out which could affect test results (refer to Limitations of the Procedure, item 10).

- Add 50 µL of red blood cells (cells must be diluted in MTS™ Diluent 2 to approximately 0.8% or be a commercial 0.8% red blood cell in a low ionic strength diluent specifically approved for use in the ID-Micro Typing System™) to each microtube. It is not necessary that the blood come into contact with the gel.

Caution:

The pipet tip should not touch the gel card. Erroneous results due to carryover may occur.

- Centrifuge the prepared cards in the ORTHO® Workstation at the preset conditions installed by the manufacturer.
- After centrifugation, remove the gel card(s) from the centrifuge. Observe, read macroscopically the front and back of each microtube for agglutination and/or hemolysis and record reactions. If either side of the microtube is positive, the reaction is to be considered positive. See Diagram 1.

Interpretation of Results

Refer to ID-Micro Typing System™ Interpretation Guide ⁴ for additional information.

Negative Result : No agglutination and no hemolysis of the red blood cells is a negative test result. A complete sedimentation of all red blood cells is present in the bottom of the microtube. A negative test result indicates the absence of detectable IgG or C3 on the red blood cells.

Positive Result : Agglutination and/or hemolysis of the red blood cells is a positive test result. Red blood cells may remain suspended on the top of the gel or are dispersed throughout the gel in varying degrees. A few red blood cells may form a button in the bottom of the microtube in some positive reactions. A positive test result indicates the presence of IgG and/or C3 on the red blood cells.

Reaction Grading Guide (Use in conjunction with Diagram 1)

0 Negative	Unagglutinated red blood cells form a well-defined button at the bottom of the microtube.
1+ Reaction	Red blood cell agglutinates are observed predominantly in the lower half of the gel microtube. Unagglutinated red blood cells form a button in the bottom of the microtube.
2+ Reaction	Red blood cell agglutinates are dispersed throughout the length of the gel microtube. Few unagglutinated red blood cells may be observed in the bottom of the microtube.
3+ Reaction	The majority of red blood cell agglutinates are trapped in the upper half of the gel microtube.
4+ Reaction	Solid band of red blood cell agglutinates on top of the gel. A few agglutinates may filter into the gel but remain near the predominant band.
Mixed Field	Red blood cell agglutinates at the top of the gel or dispersed throughout the gel microtube accompanied by a button of negative red blood cells in the bottom of the microtube. See Note below.

Note:

Caution must be taken in interpreting a reaction as mixed field. Additional patient history and testing will be necessary for resolution. However, not all mixed cell situations have a sufficient minor population to be detected.

Caution:


Clots, particulates or other artifacts may cause some red blood cells to be entrapped at the top of the gel that may cause an anomalous result in a negative test (refer to Limitations of the Procedure, item 8).

INSTRUCTIONS FOR USE

Stability of Reaction

MTS™ Anti-IgG, -C3d

Diagram 1: Examples of Reaction Grades



Stability of Reaction

For best results, it is recommended that reactions should be read immediately following centrifugation. Interpretations may be affected by the drying out of the gel, hemolysis of the red blood cells, and slanting of the reaction patterns due to storage in a non-upright position.

Quality Control

To confirm the specificity and reactivity of the MTS™ Anti-IgG, -C3d Card, it is recommended that each lot be tested each day of use with known positive and negative samples. Reactivity must be present with the positive sample only. The anti-complement reactivity of this reagent can be assessed by using complement coated cells.¹

Limitations of the Procedure

Refer to ID-Micro Typing System™ Interpretation Guide⁴ for additional information.

1. Strict adherence to the procedures and recommended equipment is essential.
2. Proper centrifuge calibration is particularly important to the performance of the MTS™ Gel Test. The ORTHO® Workstation, ORTHO VISION® Analyzer and ORTHO VISION® Max Analyzer have been exclusively designed to provide the correct time, speed and angle.
3. This card is intended for direct antiglobulin testing only.
4. Not all positive reactions imply the presence of clinically significant antibodies. It is important to distinguish between "nuisance" reactions in which cell bound serum globulins are present, but which have no clinical significance from positive reactions due to clinically significant antibodies.⁷
5. Red blood cells must be suspended in MTS™ Diluent 2 or be a commercial 0.8% red blood cell in a low ionic strength diluent specifically approved for use in the ID-Micro Typing System™.
6. Variations in red blood cell concentration can markedly affect the sensitivity of test results.¹ If red blood cell suspensions are too concentrated, they can give weaker results due to the increase in the antigen/antibody ratio. In addition, red blood cells may fail to completely migrate to the bottom of the microtube and could cause a false positive interpretation. When red blood cells are too low in concentration, they become difficult to visualize, and, in extreme cases, a weak positive can fail to be detected.
7. False positive or false negative test results can occur from bacterial or chemical contamination of test materials, inadequate incubation time or temperature, improper centrifugation, improper storage of materials, or omission of test samples.
8. Anomalous results may be caused by fibrin or other particulate matter in blood samples that could stick to the sides of the microtube.
9. Red blood cells that test as DAT positive should not be used in an indirect antiglobulin test procedure.
10. False-positive results may occur if a card that shows signs of drying is used in testing.
11. Negative direct antiglobulin test results do not necessarily rule out hemolytic disease of the newborn (HDN), especially if ABO incompatibility is suspected.
12. Rouleaux caused by serum or plasma with abnormally high concentrations of protein (such as in patients with multiple myeloma or Waldenstrom's macroglobulinemia or from patients who have received plasma expanders of high molecular weight) may infrequently cause difficulties in the ID-MTS™ Gel Test interpretation.⁴ False positive results or hazy reactions may occur with these samples but are rare. If false positive reactions (e.g., rouleaux, cells coated with immunoglobulins, etc.) occur in the control gel, the blood group cannot be established. Additional testing will be necessary to resolve this false positive reaction. If the control test is positive, the test cells should be washed several times in warm saline and retested.⁵ If the control test again gives a positive reaction, a valid interpretation of the results obtained cannot be made. Laboratories are advised to consult their approved procedures.

INSTRUCTIONS FOR USE

MTS™ Anti-IgG, -C3d

Specific Performance Characteristics

13. Hemolyzed and grossly icteric blood samples may cause difficulty in interpretation, and test results should be used with caution.
14. When using automated instruments, refer to the limitations contained in the operator's manual provided by the device manufacturer.
15. DAT positive test results can occur with red blood cells from patients who are undergoing treatment with monoclonal antibody therapies (CD38-directed cytolytic antibodies). Laboratories are advised to consult their approved procedures when interpreting test results if there is a history of monoclonal antibody treatment.
16. Individual laboratory procedures should be adhered to with respect to blood sample exposure and contact with laboratory cleaning and disinfectant solutions including 10% bleach, 70% isopropyl alcohol (IPA) and quaternary ammonium compounds. Open red blood cell sample exposure and direct interaction with these chemicals may induce a positive reaction in the gel card, promote red cell hemolysis and / or unacceptable physical appearance (i.e., hazy or smeared reactions).
17. The impact of red blood cell concentration due to manual versus automated sample preparation, in conjunction with the defined centrifugation specifications of the ORTHO VISION® / VISION® Max Analyzers versus the ORTHO® Workstation, though within the same tolerance range, can result in increased detectability of very weak positive samples on automated systems versus the manual gel system.

Specific Performance Characteristics

Each lot of MTS™ Anti-IgG, -C3d Card meets FDA requirements. The potency of Anti-IgG and Anti-C3d are verified by tests with red blood cells sensitized with decreasing amounts of Anti-D, Anti-Fy^a, Anti-C3b and Anti-C3d according to methods approved by FDA. Additionally, each lot is tested with a known antibody to ensure Anti-IgG sensitivity of 0.1 IU/mL or greater.

The absence of antibodies to C4 components has been confirmed by methods approved by FDA.⁸

The absence of contaminating heterophile agglutinins has been verified in tests employing group A₁, B, and O red blood cells.

Testing was performed on the ORTHO VISION® / VISION® Max Analyzers, ORTHO Optix™ Reader and ORTHO™ Workstation at four sites (two external and two internal sites: Sites 1-4).

Further investigation into the low positive percent agreement (PPA) values obtained on all platforms identified three major contributing factors specific to use of conventional tube test technology as a comparator for column agglutination technology. The major contributing factors identified were method variability of the tube test, operator variability within the tube test, and subjectivity of reaction grading in the tube test (specifically at site 4, with a greater number of different tube method operators). Interactions between these three factors contributed to a low PPA at patient sites due to increased incidence of positive samples within the patient population compared to donors.

Additional testing of patient samples was performed at a third external site (Test Site 5) and has been excluded from the overall performance analysis for all platforms. This site included samples from patients who were undergoing monoclonal antibody therapies. Identification of specific samples from patients undergoing these therapies was not made available. An investigation was performed on the impact of monoclonal antibody therapies on the gel DAT test. This investigation demonstrated that exposure of DAT negative samples to CD38-directed cytolytic antibodies can lead to unexpected positive reactions in the DAT test as the CD38 antigen is present on red blood cells (Refer to Limitation of Procedures, item 15). It was not able to be concluded whether monoclonal antibody therapies cause equivalent interference in tube technology versus column agglutination technology. It was also not able to be concluded if the high number of positive samples were from patients undergoing this therapy or if the results were due to a combination of the factors related to tube test variability observed at the other test sites in addition to the CD38-directed cytolytic antibody therapy.

Performance Characteristics on ORTHO VISION® / VISION® Max Analyzer

Clinical studies were performed at four sites (two external and two internal sites) that routinely perform immunohematology testing. Random clinical specimens including donors and patients (N=2801) were tested on the ORTHO VISION® / VISION® Max Analyzer.

Random sample results were assessed on a microtube-to-tube test basis using a paired sample comparison between the comparator tube method and the gel card under evaluation. For reaction grades to be concordant between methods either both results had to be negative, or both had to be positive (any reaction grade 1+ through 4+). For the tube test method, a weak reaction was considered positive when comparing to the gel card microtube. The combined results from all sites are summarized in the following table. Percent agreement indicates concordance between the two assays and does not indicate which method gave the correct result. The results below do not reflect testing to resolve initial discrepant results between methods.

INSTRUCTIONS FOR USE

Specific Performance Characteristics

MTS™ Anti-IgG, -C3d

Anti-IgG, -C3d Microtube Concordance

All Sites Combined	Random Clinical Specimens								
	Total			Positive			Negative		
Test	N	Percent Agreement	Lower Bound of 95% CI	N*	Percent Agreement	Lower Bound of 95% CI	N	Percent Agreement	Lower Bound of 95% CI
DAT (IgG C3d) - AHG Poly	2801	96%	95%	46	70%	57%	2755	96%	96%

* The expected frequency of positive DAT results has been reported as 1 in 1,000 to 1 in 14,000 for blood donors and 1 to 15% of hospital patients with the wide variation most likely due to different DAT techniques.⁵ In the clinical studies performed, the frequency of random positive DAT results among the clinical sample set was 1.6%. The lower bound of the 95% confidence interval (CI) for the positive sample set is low due to the small number of random positive samples encountered during the clinical study. The observed variability of the comparator tube method, as discussed at the beginning of the Specific Performance Characteristics section, also had an impact in the calculation of the 95% CI summarized in the discordant investigation discussion below.

The results of positive percent agreement and negative percent agreement for each individual test site is summarized in the following table. The results below do not reflect testing to resolve initial discrepant results between methods.

Anti-IgG, -C3d Microtube Concordance - Individual Test Site Results

Test Site	Sample Type	Total Sample Size	Overall Agree	Overall Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 1 - External	Donor	644	616	96%	1	1	100%	643	615	96%
Site 2 - External	Donor	815	803	99%	5	5	100%	810	798	99%
Site 3 - Internal	Patient	845	813	96%	25	24	96%	820	789	96%
Site 4 - Internal	Patient	497	454	91%	15	2	13%	482	452	94%

There were 115 discrepant results between the MTS™ Anti-IgG, -C3d Card and an FDA licensed Anti-IgG, -C3d Polyspecific tube reagent method across the four sites. Discordant investigation consisted of sample testing with a second, different, FDA licensed tube method followed by repeat testing of the sample using the initial discrepant test method. Nine initial positive tube reagent results and 13 initial negative tube reagent results changed to become concordant when repeated as part of discordant investigation. Four MTS™ Anti-IgG, -C3d Card results changed to become concordant when repeated as part of discordant investigation.

The following table summarizes the additional patient testing results from the third external testing site that have been excluded from the combined Anti-IgG, -C3d Microtube Concordance.

Anti-IgG, -C3d Microtube Concordance - Additional Patient Test Site

Test Site	Sample Type	Total Sample Size	Total Agree	Total Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 5 - External	Patient	673	625	93%	113	92	81%	560	533	95%

Anti-IgG, -C3d Microtube Agreement - Contrived Positive Samples

Due to the small positive sample size, IgG and C3 contrived positive samples (N=205 unique samples) were used and met a 100% point estimate agreement for each ID-MTS test system. The characterization of the contrived samples spanned a range of reactivity from 1+ to 3+.

Contrived Positive Samples VISION			Contrived Positive Samples VISION Max		
Test	N	Percent Agreement	Test	N	Percent Agreement
DAT (IgG C3d) - AHG Poly	205	100%	DAT (IgG C3d) - AHG Poly	205	100%

INSTRUCTIONS FOR USE

MTS™ Anti-IgG, -C3d

Specific Performance Characteristics

Performance Characteristics on ORTHO Optix™ Reader

Clinical studies were performed at four sites (two external and two internal sites) that routinely perform immunohematology testing. Random clinical specimens including donors and patients (N=2793) were tested on the ORTHO Optix™ Reader. Random sample results were assessed on a microtube-to-tube test basis using a paired sample comparison between the comparator tube method and the gel card under evaluation. For reaction grades to be concordant between methods either both results had to be negative, or both had to be positive (any reaction grade 1+ through 4+). For the tube test method, a weak reaction was considered positive when comparing to the gel card microtube. The combined results from all sites are summarized in the following table. Percent agreement indicates concordance between the two assays and does not indicate which method gave the correct result. The results below do not reflect testing to resolve initial discrepant results between methods.

Anti-IgG, -C3d Microtube Concordance

Random Clinical Specimens									
All Sites Combined	Total			Positive			Negative		
Test	N	Percent Agreement	Lower Bound of 95% CI	N*	Percent Agreement	Lower Bound of 95% CI	N	Percent Agreement	Lower Bound of 95% CI
DAT (IgG C3d) - AHG Poly	2793	98%	97%	45	58%	44%	2748	99%	98%

* The expected frequency of positive DAT results has been reported as 1 in 1,000 to 1 in 14,000 for blood donors and 1 to 15% of hospital patients with the wide variation most likely due to different DAT techniques.⁵ In the clinical studies performed, the frequency of random positive DAT results among the clinical sample set was 1.6%. The lower bound of the 95% confidence interval (CI) for the positive sample set is low due to the small number of random positive samples encountered during the clinical study. The observed variability of the comparator tube method, as discussed at the beginning of the Specific Performance Characteristics section, also had an impact in the calculation of the 95% CI summarized in the discordant investigation discussion below.

The results of positive percent agreement and negative percent agreement for each individual test site is summarized in the following table. The results below do not reflect testing to resolve initial discrepant results between methods.

Anti-IgG, -C3d Microtube Concordance - Individual Test Site Results

Test Site	Sample Type	Total Sample Size	Overall Agree	Overall Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 1 - External	Donor	759	752	99%	1	1	100%	758	751	99%
Site 2 - External	Donor	693	690	>99%	2	2	100%	691	688	>99%
Site 3 - Internal	Patient	843	827	98%	27	21	78%	816	806	99%
Site 4 - Internal	Patient	498	467	94%	15	2	13%	483	465	96%

There were 57 discrepant results between the MTS™ Anti-IgG, -C3d Card and an FDA licensed Anti-IgG, -C3d Polyspecific tube reagent method across the four sites. Discordant investigation consisted of sample testing with a second, different, FDA licensed tube method followed by repeat testing of the sample using the initial discrepant test method. Nine initial positive tube reagent results and 11 initial negative tube reagent results changed to become concordant when repeated as part of discordant investigation. No MTS™ Anti-IgG, -C3d Card results changed to become concordant when repeated as part of discordant investigation.

The following table summarizes the additional patient testing results from the third external testing site that have been excluded from the combined Anti-IgG, -C3d Microtube Concordance.

Anti-IgG, -C3d Microtube Concordance - Additional Patient Test Site

Test Site	Sample Type	Total Sample Size	Total Agree	Total Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 5 - External	Patient	671	626	93%	107	79	74%	564	547	97%

INSTRUCTIONS FOR USE

Specific Performance Characteristics

MTS™ Anti-IgG, -C3d

Anti-IgG, -C3d Microtube Agreement - Contrived Positive Samples

Due to the small positive sample size, IgG and C3 contrived positive samples (N=205 unique samples) were used and met a 100% point estimate agreement for each ID-MTS test system. The characterization of the contrived samples spanned a range of reactivity from 1+ to 3+.

Contrived Positive Samples		
Test	N	Percent Agreement
DAT (IgG C3d) - AHG Poly	205	100.0%

Performance Characteristics on ORTHO™ Workstation

Clinical studies were performed at four sites (two external and two internal sites) that routinely perform immunohematology testing. Random clinical specimens including donors and patients (N=2816) were tested on the ORTHO™ Workstation. Random sample results were assessed on a microtube-to-tube test basis using a paired sample comparison between the comparator tube method and the gel card under evaluation. For reaction grades to be concordant between methods either both results had to be negative, or both had to be positive (any reaction grade 1+ through 4+). For the tube test method, a weak reaction was considered positive when comparing to the gel card microtube. The combined results from all sites are summarized in the following table. Percent agreement indicates concordance between the two assays and does not indicate which method gave the correct result. The results below do not reflect testing to resolve initial discrepant results between methods.

Anti-IgG, -C3d Microtube Concordance

Random Clinical Specimens									
All Sites Combined	Total			Positive			Negative		
Test	N	Percent Agreement	Lower Bound of 95% CI	N*	Percent Agreement	Lower Bound of 95% CI	N	Percent Agreement	Lower Bound of 95% CI
DAT (IgG C3d) - AHG Poly	2816	98%	97%	46	57%	43%	2770	99%	98%

* The expected frequency of positive DAT results has been reported as 1 in 1,000 to 1 in 14,000 for blood donors and 1 to 15% of hospital patients with the wide variation most likely due to different DAT techniques.⁵ In the clinical studies performed, the frequency of random positive DAT results among the clinical sample set was 1.6%. The lower bound of the 95% confidence interval (CI) for the positive sample set is low due to the small number of random positive samples encountered during the clinical study. The observed variability of the comparator tube method, as discussed at the beginning of the Specific Performance Characteristics section, also had an impact in the calculation of the 95% CI summarized in the discordant investigation discussion below.

The results of positive percent agreement and negative percent agreement for each individual test site is summarized in the following table. The results below do not reflect testing to resolve initial discrepant results between methods.

Anti-IgG, -C3d Microtube Concordance - Individual Test Site Results

Test Site	Sample Type	Total Sample Size	Overall Agree	Overall Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 1 - External	Donor	760	753	99%	1	1	100%	759	752	99%
Site 2 - External	Donor	695	692	>99%	2	2	100%	693	690	>99%
Site 3 - Internal	Patient	852	835	98%	28	21	75%	824	814	99%
Site 4 - Internal	Patient	509	477	94%	15	2	13%	494	475	96%

There were 59 discrepant results between the MTS™ Anti-IgG, -C3d Card and an FDA licensed Anti-IgG, -C3d Polyspecific tube reagent methods across the four sites. Discordant investigation consisted of sample testing with a second, different, FDA licensed tube method followed by repeat testing of the sample using the initial discrepant test method. Nine initial positive tube reagent results and 12 initial negative tube reagent results changed to become concordant when repeated as part of discordant investigation. Two MTS™ Anti-IgG, -C3d Card results changed to become concordant when repeated as part of discordant investigation.

The following table summarizes the additional patient testing results from the third external testing site that have been excluded from the combined Anti-IgG, -C3d Microtube Concordance.

INSTRUCTIONS FOR USE

MTS™ Anti-IgG, -C3d

References

Anti-IgG, -C3d Microtube Concordance - Additional Patient Test Site

Test Site	Sample Type	Total Sample Size	Total Agree	Total Percent Agreement	Positive Sample Size	Positive Agreement	Positive Percent Agreement	Negative Sample Size	Negative Agreement	Negative Percent Agreement
Site 5 - External	Patient	685	634	93%	115	85	74%	570	549	96%

Anti-IgG, -C3d Microtube Agreement - Contrived Positive Samples

Due to the small positive sample size, IgG and C3 contrived positive samples (N=205 unique samples) were used and met a 100% point estimate agreement for each ID-MTS test system. The characterization of the contrived samples spanned a range of reactivity from 1+ to 3+.

























Contrived Positive Samples		
Test	N	Percent Agreement
DAT (IgG C3d) - AHG Poly	205	100.0%

References

1. Coombs RRA, Mourant AE, Race RR. Detection of weak and "incomplete" Rh agglutinins. A new test. *Lancet*. 1945;ii:15.
2. Malyska H, Weiland D. The gel test. *Laboratory Medicine*. 1994;25:81-85.
3. Lapiere Y, Rigal D, Adam J, et al. The gel test: a new way to detect red cell antigen-antibody reactions. *Transfusion*. 1990;30:109-113.
4. ID-Micro Typing System™ Interpretation Guide (6902201), Ortho Clinical Diagnostics.
5. Cohn CS, Delaney M, Johnson S, et al, eds. Technical Manual. 21st Edition. Bethesda, MD: AABB, 2023.
6. ID-Micro Typing System™ Implementation Guide and Procedures (6902200), Ortho Clinical Diagnostics.
7. Nasongkla M, Hummert J, Chaplin Jr. H. Weak "false positive" direct antiglobulin test C3d. *Transfusion*. 1982;22:273-275.
8. Office of Biologics Research and Review, FDA. Recommended methods for Anti-Human Globulin Evaluation. Docket No. 84S-0182.

Glossary of Symbols

The following symbols may have been used in the labeling of this product.

	Do Not Reuse		Contains Sufficient for "n" Tests		Fragile, Handle with Care.
	Use by or Expiration Date (Day-Month-Year)		<i>In vitro</i> Diagnostic Medical Device		Keep Dry
	Batch Code or Lot Number		Upper Limit of Temperature		This end up
	Serial Number		Lower Limit of Temperature		Do Not Use if Damaged
	Catalog Number or Product Code		Temperature Limitation		Cards
	Caution		Consult instructions for use		Concentration
	Date of Manufacture		Biological Risks		Der Grüne Punkt (the Green Dot). Manufacturer follows certain packaging material waste disposal management regulations.
	Manufacturer		Health Hazards		
	Authorized Representative in the European Community				

INSTRUCTIONS FOR USE

Revision History

MTS™ Anti-IgG, -C3d

Revision History

Date of Revision	Version	Description of Technical Changes*
2025-01-17	7.0	<ul style="list-style-type: none">• Replacement of rabbit polyclonal complement (C3) reagent with murine monoclonal C3 antibody reagents.• Intended Use Section updated for clarity to include:<ul style="list-style-type: none">– "Qualitative procedure for the detection of IgG or complement bound to red blood cells"– Indications for direct antiglobulin testing, patient population, sample type and anticoagulant used.– Description of instruments on which the MTS™ Anti-IgG, -C3d Card test can be used.• Specimen Collection and Preparation: Added information regarding optional cell washing.• Limitations of the Procedure: Added list item 15, 16 and 17.• Specific Performance Characteristics: Updated content and tables.• References:<ul style="list-style-type: none">– Corrected author name from "Mourant EE" to "Mourant AE"– Removed unrelated references– Updated Technical Manual reference• Updated patent statement

* The change bars indicate the position of a technical amendment to the text with respect to the previous version of the document.

MTS™ Anti-IgG, -C3d

INSTRUCTIONS FOR USE

Revision History

Patents: www.quidelortho.com/us/en/patents



Micro Typing Systems, Inc.
an Ortho-Clinical Diagnostics Company
1295 S.W. 29th Avenue
Pompano Beach, Florida 33069

U.S. License No. 1177

MTS, ID-MTS and ID-Micro Typing System are trademarks of
Micro Typing Systems, Inc.
© Micro Typing Systems, Inc., 2008–2025

Ortho Clinical Diagnostics

Micro Typing Systems, Inc.