

July 1, 2021

Linda Staswick
Regulatory Affic Project Manager
Bio-Rad Lak ratory
4000 Alfi An obel Dr
Hercuk CA S 47

Device: B Plex 2200 SARS-CoV-2 IgG

EUA Number. EUA202689

Company: Bio-R Laboratories

Indication: Qual ative and semi-quantitative detection of IgG antibodies to

SA S-CoV a Liuman serum and plasma (dipotassium EDTA, tri otalian EDTA lithium heparin, sodium heparin, and sodium chate). Intende for us as an aid in identifying individuals with an adaptive amune response to SARS-CoV-2, indicating recent or prior infect in. Emergence use of this test is limited to authorized

laboratories

Authorized Laboratories: Laboratories continued and Clinical Laboratory Improvement

Amendments of 1988 (2004) 263a, that meet requirements to perfect mode ate or high complexity tests.

### Dear Ms. Staswick:

This letter is in response to your¹ request that the Food and Drug dministration (FD) issue an Emergency Use Authorization (EUA) for emergency use of your roduct pursuant Section 564 of the Federal Food, Drug, and Cosmetic Act (the Act) (21 U.S. §360bbb).

On February 4, 2020, pursuant to Section 564(b)(1)(C) of the Act, the Secretary of Department of Health and Human Services (HHS) determined that there is a public health emergency that has a significant potential to affect national security or the health and security of United States citizens living abroad, and that involves the virus that causes COVID-19. Pursuant to Section 564 of the Act, and on the basis of such determination, the Secretary of HHS then declared that circumstances exist justifying the authorization of emergency use of in vitro

<sup>&</sup>lt;sup>1</sup> For ease of reference, this letter will use the term "you" and related terms to refer to Bio-Rad Laboratories.

<sup>&</sup>lt;sup>2</sup> For ease of reference, this letter will use the term "your product" to refer to the BioPlex 2200 SARS-CoV-2 IgG for the indication identified above.

diagnostics for detection and/or diagnosis of the virus that causes COVID-19 subject to the terms of any authorization issued under Section 564(a) of the Act.<sup>3</sup>

FDA considered the totality of scientific information available in authorizing the emergency use of your product for the indication above. A summary of the performance information FDA relied upon is included in the "BioPlex 2200 SARS-CoV-2 IgG" Instructions for Use (identified below).

Having concluded that the criteria for issuance of this authorization under Section 564(c) of the Act are met. It is a thorizing the emergency use of your product, described in the Scope of Authorization of this atter (Section II), subject to the terms of this authorization.

### I. Criteria for Sance of Authorization

I have concluded that a temerge we use of your product meets the criteria for issuance of an authorization under section 564(c) of the Act, because I have concluded that:

- 1. The SARS-C V-2 can cause a serious or life-threatening disease or condition, including sever respirator illness, to humans infected by this virus;
- 2. Based on the totality of please evidence available to FDA, it is reasonable to believe that your product may be effective in Lagnos. To recent or prior infection with SARS-CoV-2 by identifying individuals with an adaptive immune response to the virus that causes COVID-19, and that the known and potential isks a your product; and
- 3. There is no adequate, approved, and available a creative to the emergency use of your product. 4

### II. Scope of Authorization

I have concluded, pursuant to Section 564(d)(1) of the Act, the de scop of this authorization is limited to the indication above.

### **Authorized Product Details**

Your product is a bead-based immunoassay intended for the qualitative and emi-qualitative detection of IgG antibodies to SARS-CoV-2 in human serum and plasma (diperation EDTA, tripotassium EDTA, lithium heparin, sodium heparin, and sodium citrate) using the Bio-Rad BioPlex 2200 System (includes BioPlex 2200 System Operation Manual). The product is intended for use as an aid in identifying individuals with an adaptive immune response to SARS-CoV-2, indicating recent or prior infection. At this time, it is unknown for how long antibodies

<sup>&</sup>lt;sup>3</sup> U.S. Department of Health and Human Services, *Determination of a Public Health Emergency and Declaration that Circumstances Exist Justifying Authorizations Pursuant to Section 564(b) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C.* § 360bbb-3. 85 FR 7316 (February 7, 2020).

<sup>&</sup>lt;sup>4</sup> No other criteria of issuance have been prescribed by regulation under Section 564(c)(4) of the Act.

persist following infection and if the presence of antibodies confers protective immunity. Semiquantitative results from your product should not be interpreted as an indication or degree of immunity or protection from reinfection. Testing is limited to laboratories certified under CLIA that meet requirements to perform moderate or high complexity tests.

Your product is an automated assay that employs fluoromagnetic, dyed beads coated with Spike protein S1 subunit (including the Receptor Binding Domain) to identify the presence of IgG antibodies to SARS-CoV-2 virus in a two-step assay format. The BioPlex 2200 System combines an alignot of sample, sample diluent, and bead reagent into a reaction vessel. The mixture is inc at 37°C. After a wash cycle, murine monoclonal anti-human IgG phyco ythrin (PE) is added to the dyed beads, and this mixture is incubated at ess congate is removed in another wash cycle, and the beads are re-suspended in Tables mixture then passes through the detector. The identity of the dyed beads conjugated 37°C. T sheath aid. T is determined b the fluor ace of the dyes, and the amount of antibody captured by the he fly rescence of the attached PE. Raw data are calculated in relative antigen is determined by fluorescence i ensity

Two additional dye beads, a Internal Standard Bead (ISB) and a Serum Verification Bead (SVB), are present in a charaction dixture to verify detector response and the addition of serum or plasma to the reaction lessel, respectively.

The system is calibrated using a set  $_{\rm SIX}$  (6) distinct calibrator vials available separately from you. Five (5) levels of antibody affinity are studed, which are used for establishing a Four-Parameter Logistic calibration curve. The smi-quantitative results are expressed in U/mL. A direct relationship exists between the amount of SARS-  $_{\rm SI}$  oV-2 IgG antibody present in the sample and the RFI calculated by the system. The quantitative ssay results are reported as negative if < 10 U/mL or positive if  $\ge 10$  U.A.

The BioPlex 2200 SARS-CoV-2 IgG includes the following naterials or other authorized materials (as may be requested under Condition M. Flow): lead Set, Conjugate and Sample Diluent.

Your product requires the use of the BioPlex 2200 SARS-CoV a.gG librator Set which is not included with your product kit but is available from you with it "BioPlex 2200 SCRS-CoV-2 IgG Calibrator Set" Instructions for Use. The BioPlex 2200 SAR CoV-2 IgG calibrator Set consists of six vials in a human serum matrix made from defibrinates lasms with added known analyte concentrations of SARS-CoV-2 recombinant antibodies and antibrates derived tom inactivated human disease state plasma.

Your product also requires the use of the BioPlex 2200 SARS-CoV-2 IgG Control Set external control kit which is not included with your product kit but is available from you with the "BioPlex 2200 SARS-CoV-2 IgG Control Set" Instructions for Use, or other authorized control material (as may be requested under Condition M below). Controls are run as outlined in the "BioPlex 2200 SARS-CoV-2 IgG" Instructions for Use and the "BioPlex 2200 SARS-CoV-2 IgG Control Set" Instructions for Use:

- Negative Control: Human serum matrix made from defibrinated/delipidated plasma.
- Positive Control: Human serum matrix made from defibrinated/delipidated plasma,

with known analyte concentrations of SARS-CoV-2 recombinant antibodies and antibodies derived from inactivated human disease state plasma.

Your product also requires the use of additional authorized materials and authorized ancillary reagents that are not included with your product and are described in the authorized labeling.

The labeling entitled "BioPlex 2200 SARS-CoV-2 IgG" Instructions for Use, the "BioPlex 2200 SARS-CoV-2 IgG Calibrator Set" Instructions for Use and the "BioPlex 2200 SARS-CoV-2 IgG Control Set" Instructions for Use (available at <a href="https://www.fda.gov/medical-devices/coror of the disease-2019-covid-19-emergency-use-authorizations-medical-devices/inco-diagn tics-euas">https://www.fda.gov/medical-devices/coror of the disease-2019-covid-19-emergency-use-authorizations-medical-devices/inco-diagn tics-euas</a>), the Product Information Card (PIC) and the following product-specific aform tion potaining to the emergency use is required to be made available as set forth in the condition authorization (Section IV) and are collectively referred to "authorized labeling":

- Fact Short for a althour Providers: Bio-Rad Laboratories BioPlex 2200 SARS-CoV-2
- Fact Sheet or Recipient Bio-Rad Laboratories BioPlex 2200 SARS-CoV-2 IgG

The above described product, we see the companies that the authorized labeling as set forth in the Conditions of Authorization (Section IV), is authorized to be distributed to and used by authorized laboratories under this EUA, despite the fact that it does not meet certain requirements otherwise required by applicable federal law.

I have concluded, pursuant to Section 564(a), the Argenta it is reasonable to believe that the known and potential benefits of your product, where the state of the Scope of Authorization of this letter (Section II), outweigh the known and potential risks of your product.

I have concluded, pursuant to Section 564(d)(3) of the Act based of the total w of scientific evidence available to FDA, that it is reasonable to believe that y of pursuant may be effective in diagnosing recent or prior infection with SARS-CoV-2 by identifying induciduals with an adaptive immune response to the virus that causes COVID-19, when used ansistent with the stope of Authorization of this letter (Section II), pursuant to Section 564(c)(x A) of the Act.

FDA has reviewed the scientific information available to FDA, including the information supporting the conclusions described in Section I above, and concludes that the purpose of the Scope of Authorization of this letter (Section II) meets the criteria set forth in Section 564(c) of the Act concerning safety and potential effectiveness.

The emergency use of your product under this EUA must be consistent with, and may not exceed, the terms of this letter, including the Scope of Authorization (Section II) and the Conditions of Authorization (Section IV). Subject to the terms of this EUA and under the circumstances set forth in the Secretary of HHS's determination under Section 564(b)(1)(C) of the Act described above and the Secretary of HHS's corresponding declaration under Section 564(b)(1) of the Act, your product is authorized for the indication above.

### III. Waiver of Certain Requirements

I am waiving the following requirements for your product during the duration of this EUA:

• Current good manufacturing practice requirements, including the quality system requirements under 21 CFR Part 820 with respect to the design, manufacture, packaging, labeling, storage, and distribution of your product, but excluding Subpart H (Acceptance Activities, 2 CFR 820.80 and 21 CFR 820.86), Subpart I (Nonconforming Product, 21 CFR 820.90), 2 A Supart O (Statistical Techniques, 21 CFR 820.250).

### IV. Contation of Authorization

Pursuant to Se on 5640 of Le Act, I am establishing the following conditions on this authorization:

# Bio-Rad Laborate ies (You, and A thorized Distributor(s)<sup>5</sup>

- A. Your product me, comply with the following labeling requirements under FDA regulations: the intended se state and (21 CFR 809.10(a)(2), (b)(2)); adequate directions for use (21 U GC 52(f)), (21 CFR 809.10(b)(5), (7), and (8)); appropriate limitations on the use of the device auding a formation required under 21 CFR 809.10(a)(4); and any available in fraction reg ding performance of the device, including requirements under 21 CFR 809.10(b) 12).
- B. You and authorized distributor(s) may your roduct available with the authorized labeling to authorized laboratories.
- C. You and authorized distributor(s) must make availage on your scite(s) the Fact Sheet for Healthcare Providers and the Fact Sheet for Recipitats.
- D. You and authorized distributor(s) must include a physical copy of the PIC ord with each shipped BioPlex 2200 SARS-CoV-2 IgG to authorized laboratories and will make the "BioPlex 2200 SARS-CoV-2 IgG" Instructions for Use electronically available with the opportunity to request a copy in paper form, and there such recest, you must promptly provide the requested information without adda ional cost.
- E. You and authorized distributor(s) must inform authorized laboratories and relevant public health authorities of this EUA, including the terms and conditions herein, and any updates made to your product, and authorized labeling.
- F. Through a process of inventory control, you and authorized distributor(s) must maintain records of the authorized laboratories to which they distribute the test and number of tests they distribute.

<sup>&</sup>lt;sup>5</sup> "Authorized Distributor(s)" are identified by you, Bio-Rad Laboratories, in your EUA submission as an entity allowed to distribute your product.

- G. You and authorized distributor(s) must collect information on the performance of your product. You must report to the Division of Microbiology (DMD)/Office of Health Technology 7 (OHT7)-Office of In Vitro Diagnostics and Radiological Health (OIR)/Office of Product Evaluation and Quality (OPEQ)/Center for Devices and Radiological Health (CDRH) (via email:CDRH-EUA-Reporting@fda.hhs.gov) any suspected occurrence of false positive and false negative results and significant deviations from the established performance characteristics of the product of which you become aware.
- H. You and authorized distributor(s) are authorized to make available additional form ion reling to the emergency use of your product that is consistent with, and does not see al., the arms of this letter of authorization.
- I. You are authorised districtor(s) must make available the calibration material (BioPlex 22 to SARS-CoV-2 IgG Calibrator Set) with the "BioPlex 2200 SARS-CoV-2 IgG Calibrator Set instructions for Use and the control material (BioPlex 2200 SARS-CoV-2 IgG Control Set" Instructions for Use of other authorized calibration or control materials (as may be requested under Condition M below at the same time as your product.

## **Bio-Rad Laboratories (You)**

- J. You must notify FDA of any authorized distributor(s) of your product, including the name, address, and phone number of any authorized attributor(s).
- K. You must provide authorized distributor(s) with the second EUA and communicate to authorized distributor(s) any subsequent amount that might be made to this EUA and its authorized accompanying materials (e.g., Fa Sheets).
- L. You must comply with the following requirements Aursus Att SDA regulations: 21 CFR 820 Subpart H (Acceptance Activities, 21 CFR 820.80 and 21 CR 820.86) Subpart I (Nonconforming Product, 21 CFR 820.90), and Subpart C (Statistical Teachinges, 21 CFR 820.250).
- M. You may request changes to this EUA for your product, including a the Scope Authorization (Section II in this letter) or to the authorized labeling, sluding equests to make available additional authorized labeling specific to an authorized distributor. Such additional labeling may use another name for the product but otherwise must be consistent with the authorized labeling, and not exceed the terms of authorization of this letter. Any request for changes to this EUA should be submitted to DMD/OHT7-OIR/OPEQ/CDRH and require appropriate authorization from FDA prior to implementation.
- N. You must evaluate the performance and assess traceability of your product with any

<sup>&</sup>lt;sup>6</sup> Traceability refers to tracing analytical sensitivity/reactivity back to an FDA-recommended reference material.

FDA-recommended reference material(s) or established panel(s) of characterized clinical specimens. After submission to and concurrence with the data by FDA, you must update your labeling to reflect the additional testing. Such labeling updates will be made in consultation with, and require concurrence of, DMD/OHT7-OIR/OPEQ/CDRH.

- O. You must have a process in place to track adverse events, including any occurrence of false results and report to FDA pursuant to 21 CFR Part 803.
- P. You rest have lot release procedures and the lot release procedures, including the study design and statistical power, must ensure that the tests released for distribution have the nical and an extical performance claimed in the authorized labeling.
- Q. If required by FYA, you must submit lot release procedures to FDA, including sampling protects, testing protocols, and acceptance criteria, that you use to release lots of your project for distribution in the U.S. If such lot release procedures are requested by FDA, you must provide it within 48 hours of the request.
- R. If requested by N. A, you cust participate in a National Cancer Institute study on the evaluation of your product. After comission to and concurrence with the data by FDA you will update your product abeling to reflect the additional testing. Such labeling updates will be made in consultation with, an require concurrence of, DMD/OHT7-OIR/OPEQ/CDRH.
- S. If requested by FDA, you must periodically so mit revolots for testing at the National Cancer Institute, or by another gove, a general sesignated by FDA, to confirm continued performance characteristics across lead to the FDA may request records regarding lot release data for tests to be districted already distributed. If such lot release data are requested by FDA, you must provide it within a cours of the request.
- T. You must complete the agreed upon real-time stable ty strayer your poduct. After submission to and concurrence with the data by FDA, you will a date your labeling to reflect the additional testing. Such labeling updates will be made in correct to more its ion with, and require concurrence of, DMD/OHT7-OIR/OPEQ/CDR

#### **Authorized Laboratories**

- U. Authorized laboratories using your product must include with test result reports, all authorized Fact Sheets. Under exigent circumstances, other appropriate methods for disseminating these Fact Sheets may be used, which may include mass media.
- V. Authorized laboratories must use your product as outlined in the authorized labeling. Deviations from the authorized procedures, including the authorized instrument, authorized clinical specimen types, authorized control materials, authorized other

FDA may request, for example, that you perform this study in the event that we receive reports of adverse events concerning your product.

- ancillary reagents and authorized materials required to use your product are not permitted.
- W. Authorized laboratories that receive your product must notify the relevant public health authorities of their intent to run your product prior to initiating testing.
- X. Authorized laboratories using your product must have a process in place for reporting test results to healthcare providers and relevant public health authorities, as appropriate.
- Y. Author 260 boratories must collect information on the performance of your product and report to DMN OHT7-OIR/OPEQ/CDRH (via email: CDRH-EUA
  port v@fd hhs.gov) and to you (TechSupport.USSD@bio-rad.com) any suspected occurre to calse visitive or false negative results and significant deviations from the established performance characteristics of your product of which they become aware.
- Z. All laborates personnel using your product must be appropriately trained in automated immunoaste techniques and se appropriate laboratory and personal protective equipment who handling the kit, and use your product in accordance with the authorized labels. All laboratory personnel using the assay must also be trained in and be familiar with the interpretation of results of the product.

# Bio-Rad Laboratories (You), Authorized Astribut. (s) and Authorized Laboratories

AA. You, authorized distributor(s), an authorized la oratories using your product must ensure that any records associated with this ELL are a intained until otherwise notified by FDA. Such records will be made with the FLA for inspection upon request.

## Conditions Related to Printed Materials, Adverting an Promotion

- BB. All descriptive printed matter, advertising, and proportional caterials relating to the use of your product, shall be consistent with the authors ed laterials as well as the terms set forth in this EUA and meet the requirements set forth it section (2(a), (q)(1), and (r) of the Act, as applicable, and FDA implementing regulation
- CC. No descriptive printed matter, advertising, or promotional material relating to the use of your product may represent or suggest that this test is safe or exective for the detection of SARS-CoV-2.
- DD. All descriptive printed matter, advertising, and promotional materials relating to the use of your product shall clearly and conspicuously state that:
  - This product has not been FDA cleared or approved but has been authorized for emergency use by FDA under an EUA for use by authorized laboratories;
  - This product has been authorized only for detecting the presence of IgG antibodies to SARS-CoV-2, not for any other viruses or pathogens; and

• The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

The emergency use of your product as described in this letter of authorization must comply with the conditions and all other terms of this authorization.

# V. Durati of Auth sization

This EVA will be active at il the declaration that circumstances exist justifying the authorization of the emergency use of in vitro diagnostic tests for detection and/or diagnosis of COVID-19 is rminar a under Section 564(b)(2) of the Act or the EUA is revoked under Section 564(g) of the Act.

sincerely,

RADM inise ... Hinton

Food and tration

Enclosure