

## **FDA Executive Summary**

Prepared for the  
**Spring 2024 review** by the  
FDA's Pediatric Advisory Committee

**Medtronic Activa Neurostimulator for Dystonia Treatment  
H020007**

## Table of Contents

I. INTRODUCTION .....	3
II. ANNUAL DISTRIBUTION NUMBER (ADN) AND US DEVICE DISTRIBUTION DATA .....	3
III. POST-MARKET DATA: MEDICAL DEVICE REPORTS (MDRs).....	4
IV. POST-MARKET LITERATURE REVIEW: SAFETY DATA .....	11
V. REFERENCES OF NOT INCLUDED PAPERS .....	16

## **I. INTRODUCTION**

In accordance with the Pediatric Medical Device Safety and Improvement Act, this review provides a safety update based on the post-market experience with the use of the Medtronic Activa® Dystonia Therapy in pediatric patients since approval in 2003. The purpose of this review is to provide the Pediatric Advisory Committee (PAC) with post-market safety data so the committee can advise the Food and Drug Administration (FDA) on whether they have any new safety concerns and whether they believe that the Humanitarian Device Exemption (HDE) remains appropriately approved for pediatric use.

The Medtronic Activa® Dystonia Therapy system is indicated for unilateral or bilateral stimulation of the internal globus pallidus (GPi) or the subthalamic nucleus (STN) to aid in the management of chronic, intractable (drug refractory) primary dystonia, including generalized and/or segmental dystonia, hemidystonia, and cervical dystonia (torticollis) in patients seven years of age or above. Other Medtronic device models have been approved under the dystonia therapy in pediatric patients' indication for use HDE H020007. For the purposes of this document, Medtronic Activa® Dystonia Therapy describes any device model approved under this HDE (H020007).

This memorandum summarizes the safety data regarding H020007 for the current review period including pre-market clinical data, post-market medical device reporting (MDR) for adverse events, and peer-reviewed literature regarding safety data associated with the device.

At this time, in review of the safety and effectiveness data, FDA believes the HDE remains appropriately approved for pediatric use.

## **II. ANNUAL DISTRIBUTION NUMBER (ADN) AND US DEVICE DISTRIBUTION DATA**

Section 520(m)(6)(A)(ii) of the Federal Food, Drug, and Cosmetic Act (FD&C Act) allows HDEs indicated for pediatric use to be sold for profit as long as the number of devices distributed in any calendar year does not exceed the annual distribution number (ADN). On December 13, 2016, the 21st Century Cures Act (Pub. L. No. 114-255) updated the definition of ADN to be the number of devices "reasonably needed to treat, diagnose, or cure a population of 8,000 individuals in the United States." Based on this definition, FDA calculates the ADN to be 8,000 multiplied by the number of devices reasonably necessary to treat an individual. The Medtronic Activa Dystonia Therapy Kits are composed of only the neurostimulator if used for neurostimulator replacement or include the neurostimulator, extension, lead, and controller for implantation of the entire system. Therefore, the number of kits implanted provides a reasonable representation of the number of individuals treated with the device. No Medtronic Activa Dystonia Kits were sold in the US in the year 2023 (see below). The ADN of 8,000 has not been exceeded in 2023.

<b>Medtronic Dystonia Kit Number</b>	<b>Number of Kits Sold</b>
3307	0
3309	0
3310	0
3317	0
3319	0

3320	0
3330	0
3337	0
3339	0
33TH17	0
33TH19	0
33TH37	0
33TH39	0
33TH40	0
33TH47	0
33TH49	0
33TH57	0
33TH59	0
<b>Total</b>	<b>0</b>

Data timeframe: January 1, 2023 - December 31, 2023

<b>Number of dystonia devices implanted and active implants (in use) in the calendar year 2023</b>	
#devices implanted	508
#active implants	4079
#implants in pediatric patients in the year	76
#active implants in pediatric patients in the year	440

Data timeframe: January 1, 2023 - December 31, 2023

### **III. POST-MARKET DATA: MEDICAL DEVICE REPORTS (MDRs)**

#### **Overview of the MDR Database**

Each year, the FDA receives over 1.4 million MDRs of suspected device-associated deaths, serious injuries, and malfunctions. The database houses MDRs submitted to the FDA by mandatory reporters (manufacturers, importers, and device user facilities) and voluntary reporters such as health care professionals, patients, and consumers. The FDA uses MDRs to monitor device performance, detect potential device-related safety issues, and contribute to benefit-risk assessments of these products. MDR reports can be used effectively to:

- Establish a qualitative snapshot of adverse events for a specific device or device type
- Detect actual or potential device problems used in a “real world” setting, including:
  - Rare, serious, or unexpected adverse events
  - Adverse events that occur during long-term device use
  - Adverse events associated with vulnerable populations
  - Use error

Although MDRs are a valuable source of information, this passive surveillance system has limitations, including: the potential submission of incomplete, inaccurate, untimely, unverified, or biased data. In addition, the incidence or prevalence of an event cannot be determined from this reporting system alone due to potential under-reporting of events and lack of information about frequency of device use. Because of this, MDRs comprise only one of the FDA's several important post-market surveillance data sources.

- MDR data alone cannot be used to establish rates of events, evaluate a change in event rates over time, or compare event rates between devices. The number of reports cannot be interpreted or used in isolation to reach conclusions about the existence, severity, or frequency of problems associated with devices.
- Confirming whether a device caused a specific event can be difficult based solely on information provided in a given report. Establishing a cause-and-effect relationship is especially difficult if circumstances surrounding the event have not been verified or if the device in question has not been directly evaluated.
- MDR data is subject to reporting bias, attributable to potential causes such as reporting practice, increased media attention, and/or other agency regulatory actions.
- MDR data does not represent all known safety information for a reported medical device and should be interpreted in the context of other available information when making device-related or treatment decisions.

### MDRs Associated with the Medtronic Activa Neurostimulator for Dystonia Treatment

The Agency searched the MDR database to identify reports associated with the Medtronic Activa Neurostimulator for Dystonia Treatment entered between September 28, 2022 and September 27, 2023. The reports entered during this timeframe are related to devices implanted between February 24, 2010 and June 28, 2023. The search resulted in the identification of 196 MDRs. For the purpose of this MDR analysis, these 196 MDRs will be referred to as the 2024 PAC data. All of the MDRs were submitted by the manufacturer (N= 196 MDRs). Patient gender information was reported in 154 of the MDRs of which 94 were female and 60 were male patients. The event types by age category are presented in Tables 1a, 1b, and 1c. The number of MDRs in PAC data sets by PAC year are displayed graphically in Chart 1.

**Table 1a. Event types by age category for MDRs included in the 2015, 2016, 2017, and 2018 PAC data sets.**

Event Type	2015 PAC				2016 PAC				2017 PAC				2018 PAC			
	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total
Malfunction	19 (13.9)	91 (66.9)	26 (19.1)	136	22 (15.1)	101 (69.6)	22 (15.1)	145	27 (15.9)	107 (63.3)	35 (20.7)	169	29 (15.5)	136 (72.7)	22 (11.7)	187
Injury	22 (15.2)	84 (58.3)	38 (26.3)	144	34 (18.3)	122 (65.9)	29 (15.6)	185	31 (20.1)	90 (58.4)	33 (21.4)	154	18 (12.1)	102 (68.9)	28 (18.9)	148
Death	1 (50)	1 (50)	0 (0)	2	0 (0)	0 (0)	3 (100)	3	0 (0)	1 (100)	0 (0)	1	6 (75)	2 (25)	0 (0)	8
<b>Total</b>	<b>42 (14.8)</b>	<b>176 (62.4)</b>	<b>64 (22.6)</b>	<b>282</b>	<b>56 (16.8)</b>	<b>223 (66.9)</b>	<b>54 (16.2)</b>	<b>333</b>	<b>58 (17.9)</b>	<b>198 (61.1)</b>	<b>68 (20.9)</b>	<b>324</b>	<b>53 (15.4)</b>	<b>240 (69.9)</b>	<b>50 (14.5)</b>	<b>343</b>


**Table 1b. Event types by age category for MDRs included in the 2019, 2020, 2021, and 2022 PAC data sets.**

Event Type	2019 PAC				2020 PAC				2021 PAC				2022 PAC			
	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total
Malfunction	22 (16.2)	102 (75.5)	11 (8.1)	135	24 (18.6)	98 (75.9)	7 (5.4)	129	9 (12)	50 (66.6)	16 (21.3)	75	8 (8.8)	56 (61.5)	27 (29.7)	91
Injury	19 (21.3)	56 (62.9)	14 (15.7)	89	20 (26.6)	47 (62.6)	8 (10.6)	75	10 (15.1)	37 (56)	19 (28.7)	66	10 (13)	36 (46.8)	31 (40.2)	77
Death	0 (0)	3 (100)	0 (0)	3	0 (0)	0 (0)	0 (0)	0	0 (0)	0 (0)	0 (0)	0	0 (0)	0 (0)	0 (0)	0
<b>Total</b>	<b>41 (18)</b>	<b>161 (70.9)</b>	<b>25 (11)</b>	<b>227</b>	<b>44 (21.5)</b>	<b>145 (71)</b>	<b>15 (7.3)</b>	<b>204</b>	<b>19 (13.4)</b>	<b>87 (61.7)</b>	<b>35 (24.8)</b>	<b>141</b>	<b>18 (10.7)</b>	<b>92 (54.7)</b>	<b>58 (34.5)</b>	<b>168</b>

**Table 1c. Event types by age category for MDRs included in the 2023 and 2024 PAC data sets.**

Event Type	2023 PAC				2024 PAC			
	PEDS (%)	ADULT (%)	UNK (%)	Total	PEDS (%)	ADULT (%)	UNK (%)	Total
Malfunction	20 (19.6)	50 (49.0)	32 (31.4)	102	12 (9.6)	58 (46.7)	54 (43.5)	124
Injury	13 (14.6)	46 (51.7)	30 (33.7)	89	11 (15.2)	32 (44.4)	29 (40.2)	72
Death	0 (0)	0 (0)	0 (0)	0	0 (0)	0 (0)	0 (0)	0
<b>Total</b>	<b>33 (17.3)</b>	<b>96 (50.2)</b>	<b>62 (32.5)</b>	<b>191</b>	<b>23 (11.7)</b>	<b>90 (45.9)</b>	<b>83 (42.3)</b>	<b>196</b>

**Chart 1. The Number of MDRs in Activa PAC data set by year**



Patient age was available in 113 MDRs, which included 23 pediatric reports and 90 adult reports. The patient age was unknown in 83 reports. The number of MDRs that originated in the United States (US)

and outside of the US (OUS) for the 2024 PAC data is presented by age category in Table 2. The majority of MDRs originated from within the US.


**Table 2. The Number of US and OUS MDRs by age category in the 2024 PAC data set**

Reporter Country	Pediatric	Adult	Unknown	Total
US	18	79	44	141
OUS	5	10	38	53
Unknown	0	1	1	2
<b>Total</b>	23	90	83	196

Pediatric MDR Review (N= 23)

The reporting country for the majority of Pediatric MDRs was the United States (N= 18 MDRs) and 5 MDRs were reported from outside the United States. Within the pediatric reports, 11 MDRs were associated with female patients, 10 MDRs were associated with male patients, and 2 MDRs did not report patient gender. Pediatric patient age ranged from 10 years of age to 21 years of age. The average age of the patients in the pediatric reports was 16 years of age. The percentages of pediatric reports within PAC data sets reviewed annually between the 2015 and 2024 datasets ranged from 11% and 22% (see Chart 2).

**Chart 2. Percentage of Pediatric Reports by PAC Dataset Year**



*Time to Event (TTE) for Pediatric MDRs*

In an effort to separate reports for events that occurred zero to 30 days post-implant from those that occurred greater than 30 days post-implant, an analysis of the TTE was conducted on the pediatric MDRs. The TTE was calculated based on the implant date and date of event provided for each report; and was calculable for 19 of the 23 pediatric reports received. Reported problems and event types for pediatric MDRs by TTE are presented in Tables 3 and 4. The range of TTE was from 0 to 2693 days with an average of 819 days and median of 673 days.

There were 5 reports in which the event occurred between zero- and 30-days post-implant procedure and 14 reports in which the event occurred greater than 30 days post-implant procedure (see Table 3 and Table 4).

**Table 3. Reported problems and event types for pediatric MDRs\* in the 2024 PAC data set with TTE ≤ 30 days (n= 5)**

Reported Problem	Injury	Malfunction
Device explanted	2	0
Impedance issue	2	0
Battery charging issue	1	1
Infection	1	0
Discomfort	2	0
Lead break/fracture	1	0
Worsening symptoms	1	0
Electromagnetic Interference	0	0

\* A single MDR may be associated with more than one problem of clinical interest.

**Table 4. Reported problems and event types for pediatric MDRs\* in the 2024 PAC data set with TTE > 30 days (n=14)**

Reported Problem	Injury	Malfunction
Impedance issue	2	2
Battery charging issue	2	6
Device explanted	5	0
Worsening symptoms	0	0
Discomfort	0	0
Infection	2	0
Lead break/fracture	2	0
Electromagnetic Interference	0	0

\* A single MDR may be associated with more than one problem of clinical interest.

All pediatric reports were individually reviewed to identify events that were previously determined to be clinically significant or concerning by CDRH clinicians with input from previous PAC panel members, and to be consistent with prior MDR analyses. The specific adverse events are presented in Table 5 and explained in detail in the appropriate subsections below by the number of unique events. Please note that more than one contributing factor may have been associated with each of the events presented in Table 5.

**Table 5. Clinically concerning pediatric reports\* in the 2024 PAC data set**

Adverse Event	MDR Report Count	Number of Unique events
Battery/Charging issue	10	9
Device explanted	8	7
Device replaced	5	4
Infection	3	3
Lead break/fracture	3	2
Return or worsening of symptoms	1	1
Potential electromagnetic interference	0	0
Cognitive issue	0	0
Stroke	0	0



\* A single MDR may be associated with more than one type of adverse event.

- Battery/Charging Issues (N=10 MDRs, 9 unique events): Reports of battery/charging issues described resolved, unresolved, and unknown outcomes:
  - Resolved (N= 3 unique events)
    - Overdischarge and device position issues resolved with a Physician Recharge Mode action (N= 1)
    - Device communication issue resolved with technical services troubleshooting (N= 1)
    - Implanted stimulator (battery) migration/flipped resolved with surgical intervention (N= 1)
  - Unresolved (N= 3 unique events)
    - Device communication issue/poor coupling. MDR noted that a new recharger was being sent (N= 2)
    - Rapid battery drain and low impedance on 2 electrodes (N= 1)
    - Lost charging unit and device unresponsive (N= 1)
  - Unknown (N= 3 unique events)
    - Overdischarged stimulator (N= 1)
    - Premature battery discharge (less than 1 year) no anomalies found during device evaluation (N= 1)
    - Device position issue; implanted under the muscle (N= 1)
- Device Explant (N= 8 MDRs, 7 unique events) and Device Replacement (N= 5 MDRs, 4 unique events):
  - 3 unique events were associated with explant without replacement described as impedance issue and discomfort (N= 1) and infection (N= 2)
  - 4 unique events note explant and replacement and were associated with
    - Impedance issue and infection (N= 1)
    - Impedance issue and broken lead (N= 1)
    - Battery/charging issue (N=1)
    - Foreign body reaction (N= 1)
- Infection (N= 3 MDRs, 3 unique events):
  - Stimulator/battery implant site infection with battery and leads removed without reported replacement. TTE of 88 days (N= 1)
  - Incision site infection twice, impedance issues, and device replacement. TTE 22 days (N= 1)
  - Infection at extension wire towards base of skull. Entire system removed without replacement due to it being the patient's second infection. TTE 76 days (N=1)
- Return or Worsening of Dystonia Symptoms (N= 1 MDRs, 1 unique events): Reports of worsening of dystonia symptoms described as unresolved:
  - One (1) contact (electrode) broken on implant date and patient discomfort
- Lead break/fracture (N= 3 MDRs, 2 unique events):
  - Unresolved worsening dystonia symptoms and 1 contact broken on implant date. Also described in above bullet (N= 1 unique event)
  - One unique event described impedance issues and kinked lead resolved with replacing all four leads (N= 1 unique event)

## MDR Conclusions


A total of 23 MDRs, reporting 21 unique events, were associated with use of the Dystonia indication of the Medtronic Activa® Dystonia Therapy system in pediatric patients. Device explant/replacement was the most frequently reported pediatric patient problem. The labeling does address the issue and these events are known to occur with use of other neurostimulators. Other reported patient problems are noted in either the device labeling and/or clinical summary.

The most frequently reported device problem was battery/charging issues associated with device were overdischarge and communication issues. These device problems stated in the MDRs are noted in the device labeling or are known device issues with neurostimulator devices in general.

No MDRs associated with pediatric death were reported within the 2024 PAC data.

No new patient or device problems were identified in the 2024 PAC data when compared to PAC data from previous years. The most frequently reported clinically significant or concerning pediatric reports by PAC year are presented in Chart 3. There were no cognitive issues reported in the PAC datasets, and stroke has only been reported in the 2016 dataset thus far.

**Chart 3. Comparison of the number of clinically concerning pediatric reports\* for 2015 – 2024 PAC data sets**



\* A single report may be associated with more than one type of adverse event.

#### IV. POST-MARKET LITERATURE REVIEW: SAFETY DATA

##### Purpose

The objective of this systematic literature review is to provide an update of post-market safety/adverse events (AEs) associated with the use of the Medtronic Activa neurostimulator. This is an update on the systematic assessment of published literature since the 2023 PAC meeting.

Specifically, the systematic review was conducted to address the following question:

- What is the safety of Medtronic Activa neurostimulator device for the treatment of dystonia in the pediatric population?

##### Methods

A literature search was conducted using similar search criteria applied in previous presentations to the PAC:

(medtronic dystonia) OR (medtronic activa deep brain stimulation) OR (medtronic dbs) OR (medtronic activa) OR (activa) OR (soletra) OR (percept) OR (dbs) AND (pediatric) AND (Dystonia).

The search was conducted on November 6, 2023 using two electronic biomedical databases (PubMed and Embase) for the period between November 6, 2021 and November 6, 2022 (dates included). Although the publication period of interest for this review is between 2022 and 2023, the searches were expanded to 2021 to guarantee that no relevant studies were missed. The following inclusion and exclusion criteria were used (Table 6):

**Table 6.**

PICOTS	Inclusion Criteria	Exclusion Criteria
<b>Population</b>	Children aged 7 to <22 with chronic, intractable primary dystonia, including generalized and/or segmental dystonia, hemidystonia, and cervical dystonia (torticollis)	Non-pediatric or combined (pediatric and adult) populations where pediatric and adult populations are not analyzed separately  Not a primary dystonia (secondary or acquired)
<b>Intervention</b>	Medtronic Activa® Dystonia Therapy system  Medtronic Percept™ PC  (both on- and off-label use)	No use of Medtronic device or unknown device
<b>Comparison</b>	<ul style="list-style-type: none"> <li>• Other active treatments or standard of care (e.g., medications,</li> </ul>	No exclusion

PICOTS	Inclusion Criteria	Exclusion Criteria
	occupational/physical therapy, speech therapy, surgery) <ul style="list-style-type: none"> <li>• No comparison group</li> </ul>	
<b>Outcomes</b>	Safety <ol style="list-style-type: none"> <li>1. New safety concerns not listed at the time of HDE approval</li> <li>2. Known/anticipated safety concerns               <ol style="list-style-type: none"> <li>a. Hemiplegia/Hemiparesis</li> <li>b. Worsening of Motor Impairment</li> <li>c. Dysphagia</li> <li>d. Sensory Impairment</li> <li>e. Speech/Language</li> <li>f. Subcutaneous Hemorrhage/Seroma</li> <li>g. Cerebral Spinal Fluid Abnormality</li> <li>h. General*                   <ol style="list-style-type: none"> <li>i. Infection</li> <li>ii. Erosion</li> <li>iii. Lead fractures</li> <li>iv. Hardware Breakage</li> <li>v. Implanted Pulse Generator (IPG) Failure</li> </ol> </li> <li>i. Déjà vu corrected by surgically revised lead placement</li> <li>j. Irritating cough with stimulation ON</li> </ol> </li> <li>3. Other AEs e.g., those similar to AEs recorded with Activa systems approved for Parkinson's disease and Essential Tremor</li> </ol> * Includes adverse events related to the system components	Studies were excluded if they did not report safety outcomes.
<b>Timing</b>	Any	No exclusion
<b>Setting</b>	US and OUS	No exclusion
<b>Study Design</b>	<ul style="list-style-type: none"> <li>• Randomized controlled trials</li> <li>• Cohort studies (prospective/retrospective)</li> <li>• Case-control studies</li> <li>• Cross-sectional studies</li> <li>• Case series and case reports</li> <li>• SLRs, meta-analyses</li> </ul>	Laboratory studies, animal studies, economic and cost-effectiveness analyses

PICOTS	Inclusion Criteria	Exclusion Criteria
		SLRs and meta-analyses for which all included references were published prior to November 6, 2022
<b>Language</b>	Articles published in English	Non-English language articles
<b>Publication Dates</b>	November 7, 2022 to November 6, 2023	Published outside of date range

Abbreviations: AE: Adverse Event; HDE: Humanitarian Device Exemption; IPG: Implanted Pulse Generator; OUS: Outside the US; SLR: Systematic Literature Review; US: United States

## Results

In total, 102 unique records were identified from the database searches and screened at the title/abstract level. After excluding 58 records that were not relevant to the review at the title/abstract level, there were 44 full-text records assessed for eligibility. Of the 44 records retrieved and screened at the full-text level, no studies were relevant to this review update. A list of excluded full texts and their reasons for exclusion is available in Table 7. Figure 1 illustrates the PRISMA diagram of the literature flow.


Three studies of relevant Medtronic devices appeared in the search. These were excluded due to population not of interest (adults with drug-resistant epilepsy)<sup>7</sup> or the outcomes of interest were not reported<sup>8,9</sup>.

**Table 7. Excluded Studies**

Reference	Reason for Exclusion
Cif et al. 2023 <sup>10</sup>	Intervention not of interest
El et al. 2023 <sup>11</sup>	Intervention not of interest
Garofalo et al. 2023 <sup>12</sup>	Intervention not of interest
Grossen et al. 2023 <sup>13</sup>	Intervention not of interest
Hernandez-Martin et al. 2023 <sup>14</sup>	Intervention not of interest
Koy et al. 2023 <sup>15</sup>	Population not of interest
Lumsden et al. 2023A <sup>16</sup>	Outcomes not of interest
Lumsden et al 2023B <sup>17</sup>	Study design not of interest
MacLean et al. 2023 <sup>7</sup>	Outcomes not of interest
McEvoy et al. 2023 <sup>9</sup>	Outcomes not of interest
Novelli et al. 2023 <sup>18</sup>	Study design not of interest
Peltola et a. 2023 <sup>19</sup>	Population not of interest
Singha et al. 2023 <sup>20</sup>	Intervention not of interest
Tunyi et al. 2023 <sup>21</sup>	Intervention not of interest
Vogt et al. 2023 <sup>22</sup>	Intervention not of interest
Wang et al. 2023 <sup>23</sup>	Intervention not of interest

Reference	Reason for Exclusion
Zaman et al. 2023 <sup>8</sup>	Outcomes not of interest
Alkubaisi et al. 2022A <sup>24</sup>	Published outside of the included date range
Alkubaisi et al. 2022B <sup>25</sup>	Published outside of the included date range
Bozkaya et al. 2022 <sup>26</sup>	Study design not of interest
Chaib et al. 2022 <sup>27</sup>	Published outside of the included date range
Dhar et al. 2022 <sup>28</sup>	Published outside of the included date range
Fasano et al. 2022 <sup>29</sup>	Published outside of the included date range
Ferrero-Turrión et al. 2022 <sup>30</sup>	Study design not of interest
Fung et al. 2022 <sup>31</sup>	Published outside of the included date range
Koy et al. 2022 <sup>32</sup>	Published outside of the included date range
Li et al. 2022 <sup>33</sup>	Published outside of the included date range
Malatt et al. 2022 <sup>34</sup>	Published outside of the included date range
Mandarano et al. 2022 <sup>35</sup>	Published outside of the included date range
Munoz et al. 2022 <sup>36</sup>	Published outside of the included date range
Okazaki et al. 2022 <sup>37</sup>	Published outside of the included date range
Salamatova et al. 2022 <sup>38</sup>	Population not of interest
Shalash et al. 2022 <sup>39</sup>	Published outside of the included date range
Srinivasan et al. 2022 <sup>40</sup>	Study design not of interest
Thomas et al. 2022 <sup>41</sup>	Study design not of interest
Villessot et al. 2022 <sup>42</sup>	Published outside of the included date range
D'Hardemare et al. 2021 <sup>43</sup>	Study design not of interest
Gimeno et al. 2021 <sup>44</sup>	Published outside of the included date range
Goswami et al. 2021 <sup>45</sup>	Published outside of the included date range
Rajan et al. 2021 <sup>46</sup>	Published outside of the included date range
Saryyeva et al. 2021 <sup>47</sup>	Study design not of interest
Tai et al. 2021 <sup>48</sup>	Published outside of the included date range
Eggink et al. 2020 <sup>49</sup>	Published outside of the included date range
Skogseid et al. 2018 <sup>50</sup>	Published outside of the included date range

**Figure. 1. Article Retrieval and Selection**



## Literature Review Conclusions

No studies related to safety of the Medtronic Activa® Dystonia Therapy in pediatric patients were identified within the searches of published literature between November 7, 2022, and November 6, 2023. No new conclusions regarding the safety of the Medtronic Activa® Dystonia Therapy in pediatric populations can be drawn at this time based on the available literature.

## SUMMARY

FDA’s Review Team has identified no new safety concerns compared to what was known/anticipated at the time of HDE approval in 2003. Based on the available data, and taking into account the probable benefits and risks, FDA concludes that the HDE remains appropriately approved for pediatric use. FDA will continue routine surveillance including MDR and literature reviews. FDA will provide focused updated safety and use data to the PAC in 2025.

FDA will continue surveillance and will report the following to the PAC in 2025:

- Annual distribution number

- MDR review
- Literature review

## **V. REFERENCES OF NOT INCLUDED PAPERS**

1. Cook DJ, Mulrow CD, Haynes RB. Systematic reviews: synthesis of best evidence for clinical decisions. *Annals of internal medicine*. 1997;126(5):376-380.
2. Cumpston M, Li T, Page MJ, et al. Updated guidance for trusted systematic reviews: a new edition of the Cochrane Handbook for Systematic Reviews of Interventions. *Cochrane Database Syst Rev*. 2019;10:ED000142.
3. Page MJ, McKenzie JE, Bossuyt PM, et al. The PRISMA 2020 statement: an updated guideline for reporting systematic reviews. *Bmj*. Mar 29 2021;372:n71. doi:10.1136/bmj.n71
4. Garritty C, Gartlehner G, Nussbaumer-Streit B, et al. Cochrane Rapid Reviews Methods Group offers evidence-informed guidance to conduct rapid reviews. *Journal of Clinical Epidemiology*. 2021/02/01/ 2021;130:13-22. doi:<https://doi.org/10.1016/j.jclinepi.2020.10.007>
5. Moher D, Liberati A, Tetzlaff J, Altman DG. Preferred reporting items for systematic reviews and meta-analyses: the PRISMA statement. *PLoS Med*. Jul 21 2009;6(7):e1000097. doi:10.1371/journal.pmed.1000097
6. Robinson KA, Chou R, Berkman ND, et al. Twelve recommendations for integrating existing systematic reviews into new reviews: EPC guidance. *J Clin Epidemiol*. Feb 2016;70:38-44. doi:10.1016/j.jclinepi.2015.05.035
7. MacLean JA, Nataraj J, Olaya J, Liker MA, Sanger TD. Deep brain stimulation in an adolescent with hypomyelination with atrophy of the basal ganglia and cerebellum due to a TUBB4A mutation: illustrative case. *J Neurosurg Case Lessons*. Jul 10 2023;6(2)
8. Zaman Z, Straka N, Pinto AL, et al. Deep brain stimulation for medically refractory status dystonicus in UBA5-related disorder. Note. *Movement Disorders*. 2023;38(9):1757-1759.
9. McEvoy SD, Limbrick DD, Raskin JS. Neurosurgical management of non-spastic movement disorders. *Childs Nerv Syst*. Oct 2023;39(10):2887-2898.
10. Cif L, Demailly D, Gehin C, et al. Deep brain stimulation effect in genetic dyskinetic cerebral palsy: The case of ADCY5- related disease. *Mol Genet Metab*. Jan 2023;138(1):106970.
11. El Otmani H, El Moutawakil B, Daghi M, et al. Deep Brain Stimulation for Dystonia: Experience of a Moroccan University Hospital. *Pediatr Neurol*. Nov 2023;148:23-27.
12. Garofalo M, Beudel M, Dijk JM, et al. Elective and Emergency Deep Brain Stimulation in Refractory Pediatric Monogenetic Movement Disorders Presenting with Dystonia: Current Practice Illustrated by Two Cases. *Neuropediatrics*. Feb 2023;54(1):44-52.
13. Grossen A, Shi HH, Schenk M, et al. Deep brain stimulation in a pediatric dystonia patient with cochlear implants and mitochondrial disorder: novel application of a frameless stereotactic system and navigating the anesthesia choice and neurosurgical complexities. Illustrative case. *J Neurosurg Case Lessons*. Aug 7 2023;6(6)
14. Hernandez-Martin E, Vidmark JSL, MacLean JA, Sanger TD. What is the effect of benzodiazepines on deep brain activity? A study in pediatric patients with dystonia. Article. *Frontiers in Neurology*. 2023;14
15. Koy A, Kühn AA, Schiller P, et al. Long-Term Follow-Up of Pediatric Patients with Dyskinetic Cerebral Palsy and Deep Brain Stimulation. Article. *Movement Disorders*. 2023;38(9):1736-1742.
16. Lumsden DE, Cif L, Capuano A, Allen NM. The changing face of reported status dystonicus - A



systematic review. *Parkinsonism Relat Disord.* Jul 2023;112:105438.

doi:10.1016/j.parkreldis.2023.105438

17. Lumsden DE, Tsagkaris S, Perides S, et al. 10-Year Follow-Up of Deep Brain Stimulation in the Management of Childhood Dystonia. Conference Abstract. *Developmental Medicine and Child Neurology.* 2023;65:11.
18. Novelli M, Galosi S, Zorzi G, et al. GNAO1-related movement disorder: An update on phenomenology, clinical course, and response to treatments. *Parkinsonism Relat Disord.* Jun 2023;111:105405.
19. Peltola J, Colon AJ, Pimentel J, et al. Deep Brain Stimulation of the Anterior Nucleus of the Thalamus in Drug-Resistant Epilepsy in the MORE Multicenter Patient Registry. *Neurology.* May 2 2023;100(18):e1852-e1865.
20. Singha S, Dwarakanath S, Yadav R, et al. Deep brain stimulation in pediatric dystonia: calls for therapeutic realism over nihilism. *Childs Nerv Syst.* Oct 25 2023;
21. Tunyi J, Abreu NJ, Tripathi R, et al. Deep Brain Stimulation for the Management of AIFM1-Related Disabling Tremor: A Case Series. Article. *Pediatric Neurology.* 2023;142:47-50.
22. Vogt LM, Yan H, Santyr B, et al. Deep Brain Stimulation for Refractory Status Dystonicus in Children: Multicenter Case Series and Systematic Review. *Ann Neurol.* Sep 15 2023;
23. Wang KL, Li JP, Shan YZ, et al. Centromedian-parafascicular complex deep brain stimulation improves motor symptoms in rapid onset Dystonia-Parkinsonism (DYT12-ATP1A3). Letter. *Brain Stimulation.* 2023;16(5):1310-1312.
24. Alkubaisi A. Neuromodulation With DBS as a Rescue Strategy in Life-Threatening Paediatric Dystonic Storm. Conference Abstract. *Neuromodulation.* 2022;25(6):S1.
25. Alkubaisi A, Sandhu MK, Polyhronopoulos NE, Honey CR. Deep brain stimulation as a rescue for pediatric dystonic storm. Case reports and literature review. Review. *Interdisciplinary Neurosurgery: Advanced Techniques and Case Management.* 2022;30
26. Bozkaya Yilmaz S, Olgaç Dünder N, Basarir G, Gencpinar P, Ersen A. DBS Experience in A Pediatric Neurowilson Case. Conference Abstract. *Movement Disorders.* 2022;37:S119.
27. Chaib H, Schoene-Bake JC, Saryyeva A, Jack T, Hartmann H, Krauss JK. DBS emergency surgery for treatment of dystonic storm associated with rhabdomyolysis and acute colitis in DYT-GNAO1. *Childs Nerv Syst.* Sep 2022;38(9):1821-1824.
28. Dhar D, Holla VV, Kamble N, Yadav R, Srinivas D, Pal PK. Surgical Outcomes in Rare Movement Disorders: A Report of Seventeen Patients from India and Review of Literature. Article. *Tremor and Other Hyperkinetic Movements.* 2022;12(1)
29. Fasano A, Gorodetsky C, Paul D, et al. Local Field Potential-Based Programming: A Proof-of-Concept Pilot Study. *Neuromodulation.* Feb 2022;25(2):271-275.
30. Ferrero-Turrión J, Salazar-Villacorta A, Bescos A, et al. Setting Up a Pediatric Deep Brain Stimulation Unit. Conference Abstract. *Movement Disorders.* 2022;37:S149-S150.
31. Fung EL, Mo CY, Fung ST, et al. Deep brain stimulation in a young child with GNAO1 mutation - Feasible and helpful. *Surg Neurol Int.* 2022;13:285.
32. Koy A, Kühn AA, Huebl J, et al. Quality of Life After Deep Brain Stimulation of Pediatric Patients with Dyskinetic Cerebral Palsy: A Prospective, Single-Arm, Multicenter Study with a Subsequent Randomized Double-Blind Crossover (STIM-CP). *Mov Disord.* Apr 2022;37(4):799-811.
33. Li H, Li D, Yang W, Yan H, Zhao Z, Yang H. Deep Brain Stimulation (DBS) with Subthalamic Nucleus (STN) as Target for Pediatric Patients with PKAN. *World Neurosurg.* Jul 2022;163:e317-e322.
34. Malatt C, Tagliati M. Long-Term Outcomes of Deep Brain Stimulation for Pediatric Dystonia. *Pediatr Neurosurg.* 2022;57(4):225-237.
35. Mandarano R, Danieli A, Petacchi E, et al. Deep Brain Stimulation in childhood-onset dystonia due to brain pathology. A long-term study. *Eur J Paediatr Neurol.* Mar 2022;37:62-67.
36. Munoz Chesta D, Troncoso Schifferli M, Hernandez Chavez M, Ramirez Bruggink K, Aguirre Padilla D. Deep Brain Stimulation in Progressive Generalized Dystonia in Childhood Associated With

- ADAR1 Gene Variant. Article. *Movement Disorders Clinical Practice*. 2022;9(S2):S29-S31.
37. Okazaki Y, Sasaki T, Kawai K, et al. Two Cases of Monozygotic Twins with Early-onset Isolated (DYT1) Dystonia Effectively Treated with Bilateral Globus Pallidus Internus Stimulation. *NMC Case Rep J*. 2022;9:307-312.
38. Salamatova Y, Malaty I, Ghosh S. Pediatric autoimmune Parkinsonism and response to deep brain stimulation. Article. *Child's Nervous System*. 2022;38(1):203-206.
39. Shalash A, Fayed ZY, Hamid E, et al. Outcome of pallidal stimulation of idiopathic generalized dystonia with predominant mobile truncal dystonia: case report. Article. *International Journal of Neuroscience*. 2022;132(5):429-433.
40. Srinivasan HL, Tambirajoo K, Lumsden D, et al. Long-term hardware complications in DBS for dystonia in children: incidence, prevention and management. Conference Abstract. *British Journal of Neurosurgery*. 2022;36(1):131-132.
41. Thomas KF, Krach LE, Feyma T. Functional Improvement in Pediatric Patient with Gabriele-de Vries Syndrome After Deep Brain Stimulator Placement. Conference Abstract. *PM and R*. 2022;14:S179-S180.
42. Villessot M, Demailly D, Chan-Seng E, et al. Unpredicted failure of deep brain stimulation by the impedance measures in a child with severe PANK2- gene related generalized dystonia. Letter. *Brain Stimulation*. 2022;15(4):921-923.
43. D'Hardemare V, Dorison N, Doummar D, et al. Paediatric Deep Brain Stimulation: The Challenge of Secondary Dystonia. Conference Abstract. *Stereotactic and Functional Neurosurgery*. 2021;99(SUPPL 1):71-72.
44. Gimeno H, Polatajko HJ, Cornelius V, Lin JP, Brown RG. Rehabilitation in childhood-onset hyperkinetic movement disorders including dystonia: Treatment change in outcomes across the ICF and feasibility of outcomes for full trial evaluation. Article. *European Journal of Paediatric Neurology*. 2021;33:159-167.
45. Goswami JN, Roy S, Patnaik SK. Pediatric Dystonic Storm. Article. *Neurology: Clinical Practice*. 2021;11(5):E645-E653.
46. Rajan R, Garg K, Saini A, et al. Pallidal deep brain stimulation for KMT2B related dystonia in an Indian patient. Article. *Annals of Indian Academy of Neurology*. 2021;24(4):586-588.
47. Saryyeva A, Schoene-Bake JC, Jack T, Buelmann E, Hartmann H, Krauss JK. Emergency Deep Brain Stimulation in a Girl with Dystonic Storm Due to GNAO1-Associated Hyperkinetic-Dystonic Movement Disorder. Conference Abstract. *Stereotactic and Functional Neurosurgery*. 2021;99(SUPPL 1):56.
48. Tai CH, Lee WT, Tseng SH. Dyt6 dystonia mimicking adolescent idiopathic scoliosis successfully treated by pallidal stimulation. Article. *International Medical Case Reports Journal*. 2021;14:315-321.
49. Eggink H, Toonen RF, van Zijl van Egmond JCME, et al. The effectiveness of deep brain stimulation in dystonia: A patient-centered approach. Article. *Tremor and Other Hyperkinetic Movements*. 2020;10
50. Skogseid IM, Røsby O, Konglund A, et al. Dystonia-deafness syndrome caused by ACTB p.Arg183Trp heterozygosity shows striatal dopaminergic dysfunction and response to pallidal stimulation. Article. *Journal of Neurodevelopmental Disorders*. 2018;10(1)