

Ventura® Interatrial Shunt for Reduction of Heart Failure Hospitalizations in Patients with Heart Failure with Reduced Ejection Fraction (HFrEF)

December 3, 2025

Circulatory Systems Device Panel

V-Wave, a Johnson & Johnson Company



Introduction & Trial Design

William T. Abraham, MD

College of Medicine Distinguished Professor

The Ohio State University

Chief Medical Officer

V-Wave, a Johnson & Johnson Company

Proposed Indication Limited to HFrEF Patients With Focus on Reduced Risk of HF Hospitalization

The Ventura Shunt is indicated for NYHA Class III HF patients who remain symptomatic despite guideline-directed medical therapy, have a LVEF of $\leq 40\%$, and who are judged by a Heart Team to be appropriate for shunt therapy, to reduce the risk of hospitalization for heart failure.

Agenda

Introductory Remarks and RELIEVE-HF Trial Design

Safety and Effectiveness Results

Mechanistic Basis for Differential Effects of Shunt Treatment in HFrEF vs HFpEF

Benefit-Risk Summary and Clinical Perspective

William T. Abraham, MD

College of Medicine Distinguished Professor,
The Ohio State University
Chief Medical Officer, V-Wave, a Johnson & Johnson Company

Gregg Stone, MD

Director of Academic Affairs
Professor of Medicine (Cardiology)
Professor of Population Health Sciences and Policy
Icahn School of Medicine at Mount Sinai, New York

Michael R. Zile, MD

Charles Ezra Daniel Professor of Medicine, Cardiology Division
Medical University of South Carolina

JoAnn Lindenfeld, MD

Professor of Medicine, Division of Cardiovascular Medicine
Vanderbilt University Medical Center


Significant Unmet Need for Therapies that Improve Prognosis in HFrEF Patients

- Significant residual risk of HF hospitalization and mortality despite guideline-directed medical therapy (GDMT)
- Elevated left atrial pressure is the primary cause of these poor clinical outcomes
- Lowering left atrial pressure improves clinical outcomes in HF but is difficult to achieve with medical therapy¹
- No alternative therapies to improve clinical outcomes in NYHA Class III HFrEF patients who are treated with optimal GDMT


V-Wave Ventura Interatrial Shunt

Addresses primary cause of HF hospitalizations -
Elevated left atrial pressure resulting in pulmonary congestion

- Permanently implanted across fossa ovalis
- 5 mm opening
- 14 Fr delivery system
- Performed in cardiac catheterization lab
- Procedure time ~1 hr



Ventura Interatrial Shunt Mechanism of Action



**Shunt flow increases “automatically”
when left atrial pressure rises**

Ventura Interatrial Shunt Granted Breakthrough Device Designation

Criteria	Description (Statute and Guidance)	Ventura Shunt
First Criterion	Device provides for more effective treatment or diagnosis of life-threatening or irreversibly debilitating human disease or condition	<input checked="" type="checkbox"/>
<hr/> Device also meets <u>at least one</u> of the following:		
Second Criterion	a) Represents Breakthrough Technology	<input checked="" type="checkbox"/>
	b) No approved or cleared alternatives exist – OR –	<input checked="" type="checkbox"/>
	b) Offers significant advantages over existing approved or cleared alternatives	<input checked="" type="checkbox"/>
	c) Device availability is in best interest of patients	<input checked="" type="checkbox"/>

Breakthrough Device Designation Granted in 2019

Approval Under Breakthrough Device Designation

“As part of the benefit-risk determination for Breakthrough Devices subject to a PMA, FDA may accept a greater extent of uncertainty of the benefit-risk profile for these devices if appropriate under the circumstances, including that the uncertainty is sufficiently balanced by other factors, such as the probable benefits for patients to have earlier access to the device...and adequate postmarket controls to support premarket approval...”

Well-Executed, Robust, Double-Blind Study

- Multicenter, randomized, double-blind, sham-controlled
 - 114 sites
- Symptomatic HF on optimal (i.e., maximally tolerated) GDMT
 - Assessed by Central Eligibility Committee
- Randomized 508 patients
 - > 95% NYHA Class III
 - Median follow-up 22 months
 - 98.4% follow-up at primary analysis
- Few major protocol deviations
- No confounding interventions

Pre-Study Expectations, Considerations, and Planning

- Anticipated potential treatment benefit for HFrEF and HFpEF
- Included patients across the full spectrum of LVEF
- Acknowledged functional differences between HFrEF and HFpEF phenotypes could affect response
- Randomization stratified by LVEF group
- Prespecified interaction testing between strata

Device and Procedure Safety Clearly Demonstrated

- Primary Safety Endpoint was met
- No device- or procedure-related major adverse cardiovascular or neurologic events (MACNE) through 2 years
- Stroke, MI, and thromboembolic events occurred infrequently and at similar rates in HFrEF Shunt and Control groups

Summary of Effectiveness Findings

- Primary effectiveness endpoint (win ratio) not met
- Interaction observed between LVEF strata for primary endpoint and therefore could not be pooled for analysis of effectiveness
- LVEF stratum separately analyzed → directionally opposite outcomes
- Benefit in HFrEF patients - Harm in HFpEF patients

Totality of Data Support Probable Benefit in HFrEF

	HFrEF Shunt N = 101 pts [Events]	HFrEF Control N = 105 pts [Events]	2yr Nelson-Aalen Hazard Rate Ratio (95% CI)	Reduction (%)	NNT*
Primary Endpoint HF Events [2yr]	76	134	0.49 (0.35, 0.65)	51%	1.0
All-cause death	13	20	0.48 (0.20, 1.06)	52%	6.4
LVAD/HT	1	6	0.15 (0.00, 0.98)	85%	12.7
All HFHs	41	78	0.46 (0.29, 0.68)	54%	1.6
All outpatient WHFs	21	30	0.64 (0.33, 1.17)	36%	7.3

LVAD = Left ventricular assist device

HT = Heart transplant

HFH = Heart failure hospitalization

WHF = Worsening heart failure


NNT = Number needed to treat

* Number of patients needed to treat with the shunt to prevent on average 1 event during 2-year follow-up


Biologic Plausibility: Echo Core Lab Analyses Show Differences in Cardiac Structure and Function

Differences in Left Ventricular Structure and Function Affect Responses to Interatrial Shunt

HFrEF



HFpEF




LEFT VENTRICLE **undergoes REVERSE REMODELING** (\downarrow LVEDV)
in response to Shunt placement

LEFT VENTRICLE **unable to undergo REVERSE REMODELING**
in response to Shunt placement

Differences in Right Ventricular Structure and Function Affect Responses to Interatrial Shunt

HFrEF




Enlarged RV size

Compliant RV

RIGHT VENTRICLE **able to accept** an increase in redistributed blood volume without resulting in changes in right heart size or increased PA pressure

HFpEF



Normal RV size

Non-compliant RV

RIGHT VENTRICLE **NOT able to accept** an increase in redistributed blood volume, resulting in increased right heart size and increased PA pressure

Access to Ventura Shunt Will Benefit HFrEF Patients



HFrEF high rates of HF hospitalization and mortality despite GDMT¹
HF hospitalization and mortality rates increasing since 2012^{1,2}

RELIEVE-HF Control Group →

- 15% annualized mortality rate
- 89% annualized rate of HF Events³
- 52% annualized rate of HF hospitalization

No alternative therapies to reduce risk of HF hospitalization in HFrEF patients

Ventura Interatrial Shunt is needed for HFrEF patients

Conditions of Approval Will Support Safe and Responsible Use and Data Collection*

"As part of the benefit-risk determination for Breakthrough Devices subject to a PMA, FDA may accept a greater extent of uncertainty of the benefit-risk profile for these devices if appropriate under the circumstances, including that the uncertainty is sufficiently balanced by other factors, such as the probable benefits for patients to have earlier access to the device ... and adequate postmarket controls to support premarket approval..."

...FDA intends to weigh the device's impact on patient health, including the probable benefit of earlier access to the device, against the probable risk of harm to patients from the device should subsequent data collection demonstrate that the device is ineffective or unsafe."

- Limited indication
- Requirement for local heart team for patient selection
- Extensive physician training plan
- Gradual and controlled commercial roll-out
- Robust post-approval study
- Registry enrolling all U.S. patients treated with a commercial device


Benefit of Therapy Outweighs the Risk in HFrEF

- ✓ Safety
- ✓ Effectiveness
- ✓ Biological Plausibility
- ✓ Unmet Need / Breakthrough
- ✓ Post-Approval Commitments

RELIEVE-HF Trial Design

Robust, Prespecified Randomization by LVEF

RELIEVE-HF – A Randomized, Double-Blind, Sham-Controlled Trial



Roll-in cohort = unblinded cases used to train sites on the procedure as well as to develop more safety data

Quality Control Assured Throughout the Trial

- Sponsor-independent processes
 - Eligibility committee confirmed patients met eligibility criteria, received maximally-tolerated GDMT, and remained symptomatic
 - Clinical Events Committee (CEC) adjudicated all outcomes
 - Data and Safety Monitoring Board (DSMB) provided trial oversight
 - Echocardiography core laboratory evaluated imaging
 - Data management and biostatistics
- Strict blinding procedures and blinding assessments at 2 timepoints

Key Inclusion Criteria

1. Ischemic or non-ischemic cardiomyopathy with any LVEF and documented HF for at least 6 months
2. NYHA class II, III, or ambulatory IV functional class despite maximally-tolerated class I GDMT and device therapy for HF as assessed by a Central Eligibility Committee
3. HF hospitalization within the prior 12 months and/or elevated (BMI-adjusted) BNP/NT-proBNP – both required for NYHA class II patients
4. 6-minute walk test (6MWT) \geq 100 meters to \leq 450 meters

Key Exclusion Criteria

1. Resting SBP < 90 or > 160 mmHg or intractable heart failure
2. Severe pulmonary hypertension defined as PA systolic pressure > 70 mmHg (or pulmonary vascular resistance > 4.0 WU on right heart catheterization that cannot be reduced by vasodilator therapy)
3. RV dysfunction defined as TAPSE < 12 mm or RVFAC \leq 25% on TTE
4. LVEDD > 8 cm on TTE
5. Atrial septal defect, patent foramen ovale, anomalous pulmonary venous return, corrected congenital heart defect, severe valve lesions

Final Key Exclusion Criteria After Right Heart Catheterization and Transesophageal or Intracardiac Echocardiography Performed

Hemodynamic, heart rhythm, or respiratory instability

- Need for IV vasopressor or inotrope medication
- Malignant arrhythmias
- Acute respiratory distress or hypoxemia
- SBP < 90 or > 160 mmHg
- Cardiac index < 1.5 L/min/m²
- Mean PCWP < 7 mmHg or > 35 mmHg
- Severe pulmonary hypertension*
- RA pressure \geq LA pressure (or PCWP) when left atrial pressure (PCWP) is \geq 7 mmHg

Anatomical anomaly that precludes implanting shunt

- Minimal fossa ovalis thickness > 6 mm or lengths < 10 mm
- Atrial septal defect or patent foramen ovale with > trace shunting
- Atrial septal aneurysm
- Intracardiac thrombus

* PA systolic pressure > 70 mmHg with pulmonary vascular resistance > 4.0 WU on right heart catheterization that cannot be reduced by vasodilator therapy

SBP = systolic blood pressure; PCWP = pulmonary capillary wedge pressure; RA = right atrial; LA = left atrial

Primary Safety Endpoint

- Composite of device- or procedure-related major adverse cardiovascular or neurologic events (MACNE) in shunt arm within 30 days
 - All-cause death, stroke, systemic embolism, need for open cardiac surgery or major endovascular surgical repair
- **Performance goal 11%, analyzed with an exact binomial test**

Primary Effectiveness Endpoint

- Win ratio hierarchical composite of
 1. All-cause death
 2. Heart transplantation or left ventricular assist device (LVAD) implantation (HTLV)
 3. All HF hospitalizations (HFH)
 4. All outpatient worsening HF Events (WHF)
 5. Change in health status through longest blinded follow-up measured by KCCQ-OSS
- **Analyzed by win ratio when last enrolled patient reached 12 months with longest follow-up through 24 months**

Prespecified Secondary Endpoints

1. KCCQ change from baseline to 12 months (ANCOVA)
2. All HFHs adjusted for all-cause mortality through 24 months (joint frailty)
3. Composite death, HTLV, or HFH through 24 months (time-to-first event)
4. Composite death or HFH through 24 months (time-to-first event)
5. Cumulative HFHs through 24 months (Nelson-Aalen estimate)
6. HFH through 24 months (time-to-first event)
7. Primary Effectiveness Endpoint including mortality, LVAD/Transplant, HFH, and WHF treated as an outpatient, but without KCCQ (win ratio)
8. 6MWT change from baseline to 12 months (ANCOVA)

Statistical Considerations for ITT Analysis

- Effectiveness powered for all randomized patients
- Randomization stratified by baseline LVEF (HFrEF vs HFpEF), to assess “poolability” (i.e., consistent treatment effects) using an interaction test
- Study protocol
 - *“The safety and effectiveness of the shunt according to pre-specified LVEF subgroups will be assessed by interaction testing.”¹*
 - *“Primary effectiveness endpoint analysis will be performed on a combined HFrEF and HFpEF population. The homogeneity of the treatment effect will be examined in an analysis of the interaction between treatment effect and the HFrEF/HFpEF subpopulation...”¹*



Safety and Effectiveness Results

Gregg Stone, MD

Director of Academic Affairs


Professor of Medicine (Cardiology)

Professor of Population Health Sciences and Policy

Icahn School of Medicine at Mount Sinai, New York

Results: Patient Disposition and Baseline Characteristics

RELIEVE-HF: Patient Disposition



Selected Baseline Characteristics

	HFrEF Total N = 206		HFpEF Total N = 302		p-value*
	n	%	n	%	
Age, years	70 (62, 76)		75 (68, 80)		< 0.0001
Male	168	81.6%	151	50.0%	< 0.0001
Diabetes mellitus	105	51.0%	144	47.7%	0.47
LVEF (core lab)	30.5% (24.4, 35.3)		55.5% (48.5, 62.3)		< 0.0001
Ischemic cardiomyopathy	129	62.6%	105	34.8%	< 0.0001
NYHA Class III	196	95.1%	294	97.4%	0.19
KCCQ Summary Score	55 (39, 72)		48 (33, 64)		0.002
eGFR < 60 mL/min/1.73m²	150	72.8%	226	74.8%	0.61
NTproBNP, pg/mL	2162 (1069, 3840)		1547 (816, 2700)		0.0003

Nearly all patients were NYHA class III and had multiple comorbidities

Baseline Medications and Electrical Therapies

	HFrEF Total N = 206		HFpEF Total N = 302		p-value*
	n	%	n	%	
Beta-blockers	200	97.1%	247	81.8%	< 0.0001
Renin-angiotensin system inhibitors	188	91.3%	173	57.3%	< 0.0001
ACEi	14	6.8%	56	18.5%	< 0.0001
ARB	15	7.3%	62	20.5%	< 0.0001
ARNi	159	77.2%	55	18.2%	< 0.0001
Mineralocorticoid receptor antagonists	151	73.3%	168	55.6%	< 0.0001
SGLT2 inhibitors	104	50.5%	102	33.8%	0.0002
ICD or CRT-D	184	89.3%	54	17.8%	< 0.0001
CRT-D or CRT-P	92	44.7%	37	12.3%	< 0.0001

Use of GDMT and CRT was high, especially in HFrEF

ACEi = angiotensin converting enzyme inhibitor; ARB = angiotensin II receptor blocker; ARNi = angiotensin receptor-neprilysin inhibitor

* p-value between HFrEF vs HFpEF (not reviewed by FDA)

Results: Implant Procedure Metrics

Procedure Metrics (ITT)

	Shunt N = 250		Control N = 258	
	n	%	n	%
Shunt implant attempt	250	100%	1	0.4%
Shunt implanted successfully	250	100%	1	0.4%
Procedure duration, mins		80 (59, 100)		43 (30, 55)
Fluoroscopy time, mins		14 (10, 20)		4 (2, 7)
Heparin administered, units		9000 (7000, 12000)		-
Activated clotting time, secs		291 (246, 342)		-
Contrast administered, mL		0 (0, 0)		0 (0, 0)
Hospital duration post procedure, days		1 (1, 1)		1 (1, 1)

Shunt implant success was 100%
Procedures were brief and most patients were discharged the next day

Results: Safety

Primary Safety Endpoint at 30 Days was Met (ITT)

Shunt Through 30 Days N = 250		
	Total n	Rate*
Any device/procedure-related MACNE	0	0%
All-cause death	0	0%
Stroke	0	0%
Systemic embolism	0	0%
Need for open cardiac surgery	0	0%
Need for major endovascular surgical repair	0	0%

**30-day Device- or Procedure-Related MACNE = 0% (upper 97.5% CI 1.5%),
< performance goal of 11% (p < 0.0001)**

Analysis of All Shunt-Treated Patients: Primary Safety Endpoint at 30 Days and 2 Years

Shunt-treated patients (250 ITT Shunt + 1 ITT Control + 97 Roll-in)	Shunt Through 30 Days N = 348		Shunt Through 2 Years N = 348	
	Total n	Rate*	Total n	Rate*
Any device/procedure-related MACNE	0	0%	0	0%
All-cause death	0	0%	0	0%
Stroke	0	0%	0	0%
Systemic embolism	0	0%	0	0%
Need for open cardiac surgery	0	0%	0	0%
Need for major endovascular surgical repair	0	0%	0	0%

- 30-day Device- or Procedure-Related MACNE = 0%
- 2-year Device- or Procedure-Related MACNE = 0%

All MACNE Through 2 Years: HFrEF Population

Whether or not device- or procedure-related

	Shunt N = 101		Control N = 105		Hazard Ratio (95% CI)
	Total n	Rate*	Total n	Rate*	
Any MACNE during first 2 years (first event)	16	16.6%	28	32.7%	0.56 (0.30, 1.07)
All-cause death	13	14.3%	20	26.8%	0.63 (0.31, 1.26)
Stroke	3	1.1%	2	2.2%	1.54 (0.26, 9.23)
Systemic embolism	0	0%	0	0%	-
Need for open cardiac surgery	1	1.5%	6	9.0%	0.16 (0.02, 1.32)
Need for major endovascular surgical repair	0	0%	0	0%	-

No differences in all MACNE or its component safety events between Shunt and Control

Additional Safety Endpoints Through 2 Years: HFrEF

Whether or not device- or procedure-related

	Shunt N = 101		Control N = 105		p-value
	Total n	Rate*	Total n	Rate*	
Shunt implant embolization/thrombosis	0	0%	-	-	-
Pericardial effusion/cardiac tamponade	0	0%	0	0%	-
BARC types 3 or 5 bleeding	1	1.0%	3	3.1%	0.33
Cerebrovascular events	4	4.1%	3	3.2%	0.67
CNS infarction (stroke)	3	3.1%	2	2.2%	0.63
CNS hemorrhage (intracerebral or subarachnoid)	0	0%	1	1.2%	0.31
Transient ischemic attack	1	1.0%	1	1.0%	0.98
Myocardial infarction	1	1.1%	3	3.5%	0.32
Systemic embolization	0	0%	0	0%	-
Pulmonary emboli	1	1.7%	0	0%	0.36

No safety issues were identified in Shunt-treated patients with HFrEF

Site-Reported (Blinded) SAEs and Adverse Device Effects: HFrEF Population


	Shunt N = 101	Control N = 105	RR (95% CI)	p-value
All SAE	172	282	0.59 (0.50, 0.69)	< 0.0001
Cardiovascular	102	160	0.62 (0.49, 0.79)	0.0001
Non-cardiovascular	70	122	0.56 (0.42, 0.76)	0.0002
Respiratory	14	22		
Renal urinary	11	19		
Infection	7	17		
Injury	3	11		
Gastrointestinal	13	10		
Neoplasm	2	10		
Other	20	33		
Unexpected adverse device effect	0	0		

Summary: Ventura Shunt has Strong Safety Profile

- Primary 30-day safety endpoint of device- or procedure-related MACNE met: $p < 0.0001$
- No device- or procedure-related MACNE during 2-year follow-up
- Shunt implant success: 100%
- Peri-procedural complications rare and not increased in Shunt group; and no Shunt embolization nor adherent thrombus during 2-year follow-up
- Stroke, MI, or thromboembolic events occurred infrequently and at similar rates in Shunt and blinded Control groups

Results: Effectiveness


Win Ratio Analysis for the Primary Hierarchical Composite Effectiveness Endpoint (ITT)



LVEF Strata versus 14 Exploratory Subgroups

Characteristic	LVEF Strata (HFrEF and HFpEF)	All Other Subgroups* (Age, Sex, BMI, Diabetes, Hypertension, Ischemic/Non, BL NYHA class, BL 6MWT, BL KCCQ, Shunt encapsulation process, US vs Non, Prior COVID-19 Y/N, patent/non-patent)
Stratification variable (at randomization)	YES	NO
Included in “Scientific Rationale for Study Design” section of protocol	YES	NO
Separate sample size calculation based on these groups in protocol	YES	NO
Included in “Topline Analysis” prior to unblinding for primary and secondary endpoints	YES	NO
Included in “Primary Effectiveness Endpoint” section of SAP (section 4.1.2) with an interaction analysis	YES	NO

Primary Effectiveness Outcome by LVEF



Interaction p-value = 0.0146


* Weighted for interim analysis. Confidence intervals are provided to illustrate the variability of the corresponding summary statistic; they are provided for descriptive purposes and should not be used to draw statistical inference

Analyses of the Totality of Effectiveness Data from RELIEVE-HF

Global Statistical Test (GST) Evaluates Totality of Evidence

- Integrates treatment effects from prespecified primary & secondary outcomes
 - HF Events and KCCQ (components of primary endpoint)
 - First and recurrent events
- “Null Hypothesis” = Shunt Treatment does not affect HF outcomes
- “Alternative Hypothesis” = Shunt Treatment affects HF outcomes
- GST accounts only for unique contribution of endpoints & avoids double counting

Totality of Evidence Supports Effectiveness of Ventura Interatrial Shunt in HFrEF



Low likelihood that the positive GST outcome is a chance finding (i.e. Type-I error)
 (See figure 24 of the briefing document)

Limitations of Win Ratio Analysis

Win Ratio

- Only counts one win, loss, or tie per patient pair
- Does not reflect all events patients are experiencing
- Underestimates the total burden of disease

Recurrent Events Analysis

- Cumulative hazard rates for both first and recurrent events over time
- Represents risk of all adverse outcomes
- Appropriately reflects overall burden of disease


Number of HFH and WHF Events Excluded from the Win Ratio Analysis in HFrEF Patients

		Shunt N = 101	Control N = 105
Win Ratio Tier 3 HFH Events	All HFHs	41	78
	HFHs excluded due to death or HT/LVAD	10	40
	% HFHs excluded	24.4%	51.3%
Win Ratio Tier 4 WHF Events	HFHs remaining	31	38
	All WHFs	21	30
	WHDs excluded due to death or HT/LVAD or HFH	5	10
	% WHFs excluded	24.0%	33.3%
	WHDs remaining	16	20

The win ratio analysis concealed 38% of all HF-related events, including nearly twice as many events that occurred in the control group compared with the shunt group


HFrEF Benefit: Reduction of All HF Events

(All-Cause Mortality, Cardiac Transplant or LVAD, Heart Failure Hospitalization, Worsening HF Events)



HFrEF patients \Rightarrow Greatest risk of HF Events & greater unmet clinical need
 HFrEF patients with shunt-treatment \Rightarrow Pronounced reduction in HF Events

Continuous Relationship Between Baseline LVEF and All HF Events



Sensitivity Analysis of Baseline LVEF Ranges > 40% to ≤ 45%

	LVEF > 40% to ≤ 43% N = 32		LVEF > 40% and ≤ 44% N = 38		LVEF > 40% and ≤ 45% N = 46	
	Shunt N = 17	Control N = 15	Shunt N = 20	Control N = 18	Shunt N = 23	Control N = 23
Patient-years of follow-up	26.61	22.54	30.54	26.58	35.73	35.79
Death, HT/LVAD, HFH, WHF, n (%/y)	8 (30.1%)	20 (88.7%)	12 (39.3%)	20 (75.2%)	19 (53.2%)	23 (64.3%)
RR (95% CI)	0.34 (0.14, 0.76)		0.52 (0.25, 1.06)		0.83 (0.44, 1.53)	

HF event rates were reduced by Shunt treatment up to an LVEF of 43%. The effect was attenuated when including higher LVEFs but were not worse even when patients with LVEF up to 45% were included.


Probable HFrEF Benefit: Reduction of All HF Events

	HFrEF Shunt N = 101 pts [Events]	HFrEF Control N = 105 pts [Events]	2yr Nelson-Aalen Hazard Rate Ratio (95% CI)	Reduction (%)	NNT*
Primary Endpoint HF Events [2yr]	76	134	0.49 (0.35, 0.65)	51%	1.0
All-cause death	13	20	0.48 (0.20, 1.06)	52%	6.4
LVAD/HT	1	6	0.15 (0.00, 0.98)	85%	12.7
All HFHs	41	78	0.46 (0.29, 0.68)	54%	1.6
All outpatient WHFs	21	30	0.64 (0.33, 1.17)	36%	7.3

Shunt-treated HFrEF \Rightarrow Consistent relative risk reductions for all primary HF Event components

*Number of patients needed to treat with the shunt to prevent on average 1 event during 2-year follow-up. Confidence intervals are provided to illustrate the variability of the corresponding summary statistic; they are provided for descriptive purposes and should not be used to draw statistical inference.

HFrEF Benefit: Heart Failure Hospitalization



Shunt-treated HFrEF \Rightarrow Marked reduction in HF hospitalization with ongoing benefit at 2 years


HFrEF Benefit: Adjudicated Causes of Death

	Shunt N = 101 pts	Control N = 105 pts
All-cause	13	20
Cardiovascular or unknown	12	14
Cardiovascular	11	12
Unknown	1	2
Non-cardiovascular	1	6

- There were 2 fewer “CV or other” deaths with Shunt treatment, and 5 fewer non-CV deaths with Shunt treatment
- Non-cardiac diseases trigger/worsen heart failure, and mortality rates are increased in HF patients with non-cardiac disorders (especially cancer and infections)*


Numerically fewer deaths in the HFrEF cohort treated with the Shunt

Risk of All-Cause Death or Cardiac Transplant / LVAD Implant in HFrEF (Terminal Events)




Shunt-treated patients had a 48% reduction in death, heart transplant, or LVAD treatment

Nelson-Aalen Cumulative Hazard Rates at 24 Months (HFrEF)



In HFrEF, Shunt treatment was associated with consistent, robust reductions in all event types (Including hospitalizations, outpatient worsening episodes, and terminal outcomes)

Recurrent Event Outcomes by 24 Months in Reduced LVEF $\leq 40\%$ (HFrEF) Stratum



*Non-fatal events, adjusted for competing risk of terminal events

Confidence intervals are provided to illustrate the variability of the corresponding summary statistic; they are provided for descriptive purposes and should not be used to draw statistical inference.


HFrEF Sample Size and Number of Events Were Large Enough to Provide Robust Results

HFrEF	
Total events in 206 pts	
Event Type	
Single Events	
Death	33
Heart transplant/LVAD (HT/LV)	7
Hospitalizations for HF (HFH)	119
All-cause hospitalizations (ACH)	272
Worsening HF out-patient events (WHF)	51
Composite Events	
Terminal events (Death, HT/LV)	40
Primary endpoint HF events (Death, HT/LV, HFH, WHF)	210
All events (Death, HT/LV, ACH, WHF)	363

Symmetric Trimming Analysis* for HF Events and HFHs in HFrEF Patients



- Equal numbers of patients were sequentially removed from both groups, starting with those with the highest number of HF events
- For **all HF Events**, the 95% UCB of < 1 was preserved until 23 pairs of Shunt and Control patients were removed, representing 72% of all events
- For **all HFHs**, the 95% UCB of < 1 was preserved until 13 pairs of Shunt and Control patients were removed, representing 55% of all events

These results show that the treatment effect of the atrial shunt in HFrEF is robust and distributed across the full spectrum of patients experiencing the majority of events—not just restricted to a few outliers



*FDA has reviewed the event data but has not reviewed this analysis, which was performed after Sponsor saw FDA's slide presentation.

Change in KCCQ-OSS During 2-Year Follow-up



All patient groups felt better, with no differences in KCCQ between Shunt and Control, despite marked differences in HF Events

Recent Blinded, Randomized Trials Showing Reduced Death and/or HFH with Minimal Between-Group Change in KCCQ


Study	N	Disease	Intervention	Primary Endpoint		Between-group Change in KCCQ (MCID = 5 pts)
				Clinical Event(s)	HR or RRR	
EMPEROR-Reduced ¹	3,730	HFrEF	Empagliflozin	CV death or HF hospitalization	0.75	+1.7
EMPEROR-Preserved ²	5,988	HFpEF	Empagliflozin	CV death or HF hospitalization	0.79	+1.3
DAPA-HF ³	4,744	HFrEF	Dapagliflozin	Worsening HF or CV death	0.74	+2.8
PARADIGM-HF ⁴	8,399	HFrEF	Sacubitril/ Valsartan	CV death or HF hospitalization	0.80	+1.6
FINEARTS-HF ⁵	6,001	HFpEF	Finerenone	CV death or worsening HF	0.79	+1.6
RELIEVE-HF (Reduced) ⁶	206	HFrEF	Ventura Atrial Shunt	All-cause death, LVAD/HT, all HFHs, all outpatient WHFs	0.49	+0.4
RELIEVE-HF (Preserved) ⁶	302	HFpEF			1.69	-1.7

MCID = minimal clinically important difference


1. Packer M et al. NEJM 2020; 2. Anker SD et al. NEJM 2021; 3. McMurray JJV et al. NEJM 2019; 4. McMurray JJV et al. NEJM 2014; 5. Solomon SD et al. NEJM 2024

6. Stone GW et al. Circulation 2024.

HFrEF: KCCQ Responder Analysis



HFrEF: NYHA Responder Analysis



P-values are provided for descriptive purposes and should not be used to draw statistical inference.

RELIEVE-HF Summary: The Ventura Atrial Shunt has Been Demonstrated to be Safe

- Primary Endpoint was met
 - 30-day MACNE = 0% (N = 250) < performance goal of 11% (p < 0.0001)
 - 2-year MACNE = 0% (N = 348)
- Peri-procedural complications were rare and not increased in the Shunt group; no Shunt embolizations nor adherent thrombus occurred during 2-year follow-up
- Stroke, MI, or thromboembolic events occurred infrequently and at similar rates in the Shunt and blinded Control groups
- SAEs, both CV and non-CV, were less frequent in Shunt-treated patients than in Control group patients with HFrEF

RELIEVE-HF Summary: The Totality of Data Support Ventura Atrial Shunt Effectiveness in HFrEF

	HFrEF Shunt N = 101 pts [Events]	HFrEF Control N = 105 pts [Events]	2yr Nelson-Aalen Hazard Rate Ratio (95% CI)	Reduction (%)	NNT*
Primary Endpoint HF Events [2yr]	76	134	0.49 (0.35, 0.65)	51%	1.0
All-cause death	13	20	0.48 (0.20, 1.06)	52%	6.4
LVAD/HT	1	6	0.15 (0.00, 0.98)	85%	12.7
All HFHs	41	78	0.46 (0.29, 0.68)	54%	1.6
All outpatient WHFs	21	30	0.64 (0.33, 1.17)	36%	7.3

Given absence of safety risk, the data support a positive benefit-risk ratio in HFrEF



Mechanistic Basis for Differential Effects of Shunt Treatment in HFrEF vs HFpEF

Michael R. Zile, MD

Charles Ezra Daniel Professor of Medicine
Cardiology Division, Medical University of
South Carolina


Structural and Functional Remodeling in HF

Hypothesis: Differences in cardiac structure and function between HFrEF and HFpEF provide biologically plausible explanation for differences in responses to interatrial Shunt

Structural Determinants of Mortality and Morbidity

- LV structural remodeling: Δ End diastolic volume
- RV structural remodeling: RV diastolic function
 Δ PA pressure


RELIEVE-HF Echocardiographic Study Consort Diagram: Echo from Baseline vs 12 Months




*Number of patients with paired 12-month and baseline echocardiographic values after 2-step imputation
Zile et al, JACC CV Imaging Published Online Aug 29, 2025

Differences in Left Ventricular Structure and Function in HFrEF vs HFpEF Affect Responses to Interatrial Shunt

HFrEF




HFpEF




LEFT VENTRICLE **undergoes REVERSE REMODELING** (\downarrow LVEDV)
in response to Shunt placement

LEFT VENTRICLE **unable to undergo REVERSE REMODELING**
in response to Shunt placement


Changes in Left Heart Structure in RELIEVE-HF



Relationship Between Change in LV End-Diastolic Volume and Mortality in HFrEF




Relationship Between Change in LV End-Diastolic Volume and Mortality in HFrEF




Differences in RV Structure and Function in HFrEF vs HFpEF Affect Responses to Interatrial Shunt

HFrEF



RIGHT VENTRICLE **able to accept** an increase in redistributed blood volume without resulting in changes in right heart size or increased PA pressure


HFpEF




RIGHT VENTRICLE **NOT able to accept** an increase in redistributed blood volume, resulting in increased right heart size and increased PA pressure

Changes in Right Heart Structure in RELIEVE-HF

RV End Diastolic Area Index




RA Area Index




Change from Baseline to 12 months

Changes in RV Systolic Function in RELIEVE-HF

RV Fractional Area Change




Tricuspid Annular Plane Systolic Excursion



Change from Baseline to 12 months

Changes in PA Systolic Pressure in RELIEVE-HF



Effects of Change from Baseline in PA Pressure on Subsequent All-Cause Mortality

	Analysis Method	Results	p-value
PA Diastolic	Continuous Variable	Hazard Ratio (95% CI) 1.03 (1.01, 1.05)	0.0042
	Categorical Variable	% Change in Mortality Risk	0.0237
	≤ - 2 mmHg	- 14.7%	
	≥ + 2 mmHg	+ 26.7%	
PA Systolic	Continuous Variable	Hazard Ratio (95% CI) 1.02 (1.01, 1.03)	0.0020
	Categorical Variable	% Change in Mortality Risk	0.0198
	≤ - 3 mmHg	- 14.2%	
	≥ + 3 mmHg	+ 23.8%	
PA Mean	Continuous Variable	Hazard Ratio (95% CI) 1.03 (1.01, 1.04)	0.0023
	Categorical Variable	% Change in Mortality Risk	0.0546
	≤ - 2 mmHg	- 16.7%	
	≥ + 2 mmHg	+ 13.9%	

Biologically Plausible Mechanisms for Differences in Responses to Interatrial Shunt

- At baseline, there are critical differences in cardiac structure / function
- After Shunt placement, RV compliance determined the ability to accommodate the LA to RA shunted volume
- Changes in LV remodeling and PA pressures predict subsequent CV mortality and HF morbidity
- RELIEVE-HF: HFrEF patients treated with an interatrial Shunt had improved morbidity / mortality related to structure / function characteristics



Benefit-Risk Summary and Clinical Perspective

JoAnn Lindenfeld, MD

Professor of Medicine, Division of
Cardiovascular Medicine

Vanderbilt University Medical Center

Benefit-Risk Determination Guidance For ANY Device

Factors to Consider When Making Benefit-Risk Determinations in Medical Device Premarket Approval and De Novo Classifications

Guidance for Industry and Food and Drug Administration Staff

Document issued on August 30, 2019.

Document originally issued on March 28, 2012.

This document supersedes "Factors to Consider When Making Benefit-Risk Determinations in Medical Device Premarket Approvals and De Novo Classifications" issued August 24, 2016.

For questions about this document concerning devices regulated by CDRH, contact the Office of Policy at 301-796-5441. For questions about this document concerning CBER-regulated devices, contact the Office of Communication, Outreach and Development (OCOD) by calling 800-835-4709 or 240-402-8010.



U.S. Department of Health and Human Services
Food and Drug Administration
Center for Devices and Radiological Health
Center for Biologics Evaluation and Research

"Do the benefits of the device outweigh the risks for its intended use?"

Benefit

- Is there any evidence of clinical benefit?***
- What is the extent of uncertainty for the benefit?***

Risk

- Are known/probable risks more than minimal?***
- What is the extent of uncertainty for the risks?***

Benefit - Risk

- Do the Benefits outweigh the risks?***
- Do the Benefits outweigh the risks, considering post-market actions?***

Ability to Look Beyond Statistical Significance of Primary Datasets

Assessment of Benefit

1. Is there any evidence of clinical benefit?

Is a clinical benefit demonstrated for the device for this indication (e.g., from any one or more of the primary and/or secondary datasets or from associated real-world evidence)? Benefit may be considered in terms of how a patient feels, functions, survives, or an acceptable surrogate outcome. This information may be collected using validated tools such as quality of life questionnaires. Benefit may also be considered in terms of convenience in managing or diagnosing a disease or condition. Benefit should be considered based on the assessment of the data, whether or not the results are statistically significant. *Select any of the following that demonstrate benefit, and then answer the question in the box below.*

A favorable change in at least 1 clinical assessment that is equal to or greater than seen in the control group

A favorable perfor

A favorab minim

A favorab chang

A favorab patien consid

A favorab clinic

A favorab specif treatm

Acceptabl

Other(s):

None

Q1: Is there any evi



YES → Continue

NO → Move to Q

1. Is there any evidence of clinical benefit?

Is a clinical benefit demonstrated for the device for this indication (e.g., from any one or more of the **primary and/or secondary datasets** or from associated real-world evidence)? Benefit may be considered in terms of how a patient feels, functions, survives, or an acceptable surrogate outcome. This information may be collected using validated tools such as quality of life questionnaires. Benefit may also be considered in terms of convenience in managing or diagnosing a disease or condition. Benefit should be considered based on the assessment of the data, **whether or not the results are statistically significant**. *Select any of the following that demonstrate benefit, and then answer the question in the box below.*

Heart Failure: Is There an Unmet Need?



Joyst Maddox KE et al, Circulation. 2024 Jun 4;
 Van Nuys KE et al., JACC Heart Fail. 2018; 6:401-9;
 Heidenreich PA et al. Circ Heart Fail. 2013; 6:606-19


Modified from Shah KS, et al. J Am Coll Cardiol. 2017;70:2476-86.

Heart Failure is a large & growing clinical burden and
 HFrEF patients have very poor post-admission outcomes

HFrEF Patients Have High Risk of Recurrent HFH and CV Death Despite GDMT


EMPEROR-Reduced

75% NYHA II




RELIEVE-HF HFrEF


96% NYHA III



Half of All 50-Year-Olds on Ideal GDMT Will Die Within 12 Years



Was the Study Well-Designed and Executed?



Well-Executed, Robust, Double-Blind Study

- Multicenter, randomized, double-blind, sham-controlled
 - 114 sites
- Symptomatic HF on optimal (i.e., maximally tolerated) GDMT
 - Assessed by Central Eligibility Committee
- Randomized 508 patients
 - $> 95\%$ NYHA Class III
 - Median follow-up 22 months
 - 98.4% follow-up at primary analysis
- Few major protocol deviations
- No confounding interventions

Can I Safely Treat My Patients With This Device?

Analysis of All Shunt-Treated Patients: Primary Safety Endpoint at 30 Days and 2 Years

40

Shunt-treated patients (250 ITT Shunt + 1 ITT Control + 97 Roll-in)	Shunt Through 30 Days N = 348		Shunt Through 2 Years N = 348	
	Total n	Rate*	Total n	Rate*
Any device/procedure-related MACNE	0	0%	0	0%
All-cause death	0	0%	0	0%
Stroke	0	0%	0	0%
Systemic embolism	0	0%	0	0%
Need for open cardiac surgery	0	0%	0	0%
Need for major endovascular surgical repair	0	0%	0	0%

- 30-day Device- or Procedure-Related MACNE = 0%
- 2-year Device- or Procedure-Related MACNE = 0%

* Event rates are Kaplan-Meier estimates

Additional Safety Endpoints Through 2 Years: HFrEF

Whether or not device- or procedure-related




42

	Shunt N = 101		Control N = 105		p-value
	Total n	Rate*	Total n	Rate*	
Shunt implant embolization/thrombosis	0	0%	-	-	-
Pericardial effusion/cardiac tamponade	0	0%	0	0%	-
BARC types 3 or 5 bleeding	1	1.0%	3	3.1%	0.33
Cerebrovascular events	4	4.1%	3	3.2%	0.67
CNS infarction (stroke)	3	3.1%	2	2.2%	0.63
CNS hemorrhage (intracerebral or subarachnoid)	0	0%	1	1.2%	0.31
Transient ischemic attack	1	1.0%	1	1.0%	0.98
Myocardial infarction	1	1.1%	3	3.5%	0.32
Systemic embolization	0	0%	0	0%	-
Pulmonary emboli	1	1.7%	0	0%	0.36

No safety issues were identified in Shunt-treated patients with HFrEF

* Event rates are Kaplan-Meier estimates


Will My Patients Benefit From This Device?



Is There a Plausible Biological Mechanism to Explain Differences in Treatment Effect?

Differences in Left Ventricular Structure and Function in HFrEF vs HFpEF Affect Responses to Interatrial Shunt

74




LEFT VENTRICLE undergoes REVERSE REMODELING (\downarrow LVEDV) in response to Shunt placement

LEFT VENTRICLE unable to undergo REVERSE REMODELING in response to Shunt placement

Zile et al, JACC CV Imaging Published Online Aug 29, 2025

Differences in RV Structure and Function in HFrEF vs HFpEF Affect Responses to Interatrial Shunt

78




RIGHT VENTRICLE able to accept an increase in redistributed blood volume without resulting in changes in right heart size or increased PA pressure

RIGHT VENTRICLE NOT able to accept an increase in redistributed blood volume, resulting in increased right heart size and increased PA pressure

Zile et al, JACC CV Imaging Published Online Aug 29, 2025

Changes in Left Heart Structure in RELIEVE-HF


75



Zile et al, JACC CV Imaging Published Online Aug 29, 2025

Changes in Right Heart Structure in RELIEVE-HF

79




Zile et al, JACC CV Imaging Published Online Aug 29, 2025

Proposed Indication: Supported by the Data & in the Best Interest of Patients in Need


The Ventura Shunt is indicated for **NYHA Class III** HF patients who remain symptomatic **despite guideline-directed medical therapy**, have a **LVEF of $\leq 40\%$** , and who are judged by a Heart Team to be appropriate for shunt therapy, to **reduce the risk of heart failure hospitalization**.

Totality of Evidence Support Approval


Significant Unmet Need




Well-Executed, Randomized, Sham-Controlled Trial










Convincing Safety Profile



Improvements in All Clinical Measures



Biologically Plausible Mechanism



Favorable Benefit-Risk Profile

Ventura® Interatrial Shunt for Reduction of Heart Failure Hospitalizations in Patients with Heart Failure with Reduced Ejection Fraction (HFrEF)

December 3, 2025


Circulatory Systems Device Panel

V-Wave, a Johnson & Johnson Company


ADDITIONAL SLIDES SHOWN

Quality Assessment of Echo Core Lab - LVEF

Bland-Altman Plot: Inter-
observer variability of LVEF
(Reader 1 vs Reader 2)




Bland-Altman Plot: Intra-
observer variability
(Within Reader 1)



From Table 38: Echocardiographic Data at Baseline and 12 Months in HFrEF (LVEF ≤ 40%)

	Shunt N = 101			Control N = 105			Difference Shunt- Control	ANCOVA Nominal P-value
	Baseline	12 Months	12 Month- Baseline Difference	Baseline	12 Months	12 Month- Baseline Difference		
LV end-diastolic volume index, ml/m ²	98.1 (90.9, 105.3)	94.5 (86.8, 102.1)	-3.7 (-10.2, 2.9)	96.9 (89.8, 103.9)	105.0 (97.2, 112.9)	8.2 (1.3, 15.0)	-11.9 (-21.3, -2.5)	0.01
LV end-systolic volume index, ml/m ²	69.4 (63.1, 75.7)	66.0 (59.4, 72.6)	-3.4 (-9.0, 2.1)	70.0 (63.9, 76.2)	75.6 (68.7, 82.4)	5.5 (-0.4, 11.4)	-8.9 (-17.2, -0.7)	0.03
LV ejection fraction, %	30.0 (28.4, 31.7)	32.2 (30.4, 34.0)	2.2 (0.4, 4.0)	29.2 (27.6, 30.8)	30.5 (28.8, 32.2)	1.3 (-0.4, 3.0)	0.9 (-1.5, 3.3)	0.46
LV global longitudinal strain, %	9.6 (8.9, 10.3)	10.3 (9.6, 11.0)	0.7 (-0.0, 1.5)	9.9 (9.2, 10.6)	9.7 (9.0, 10.4)	-0.2 (-0.9, 0.5)	1.0 (-0.1, 2.0)	0.06
Left atrial volume index, ml/m ²	45.2 (41.1, 49.3)	45.7 (41.5, 49.9)	0.5 (-3.7, 4.8)	40.9 (36.9, 44.9)	47.2 (43.0, 51.3)	6.3 (2.1, 10.5)	-5.8 (-11.8, 0.2)	0.06
E/e	18.2 (16.1, 20.3)	17.6 (15.3, 19.9)	-0.6 (-2.7, 1.6)	16.3 (14.2, 18.4)	19.5 (17.3, 21.6)	3.2 (1.2, 5.2)	-3.7 (-6.7, -0.7)	0.02

LVEF Values at Baseline Reflect Epidemiology




Type I Error Estimate

Three opportunities for success resulting in potential inflation of type 1 error:

1. Overall
2. HFrEF
3. HFpEF


Analysis includes primary and 7 secondaries

Tested 100,000 simulations



Since the null is true for these random scenarios, “Positive” indicates a false positive, or a Type I Error under the Null

Nelson-Aalen Cumulative Hazard Rates at 24 Months (HFrEF)



In HFrEF, Shunt treatment was associated with consistent, robust reductions in all event types
(Including hospitalizations, outpatient worsening episodes, and terminal outcomes)

Table 19: HFrEF Subgroup Medications at Baseline and 12 months

	Baseline		12 months	
	Shunt group (N = 101)	Control group (N = 105)	Shunt group (N = 101)	Control group (N = 105)
Beta-blockers	99 (98.0%)	101 (96.2%)	92 (100%)	91 (96.8%)
Renin-angiotensin system inhibitors	95 (94.1%)	93 (88.6%)	87 (94.6)	84 (89.4%)
-ACEi	7 (6.9%)	7 (6.7%)	7 (7.6%)	7 (7.4%)
-ARB	8 (7.9%)	7 (6.7%)	6 (6.5%)	3 (3.2%)
-ARNi	80 (79.2%)	79 (75.2%)	74 (80.4%)	74 (78.7%)
Mineralocorticoid receptor antagonists	74 (73.3%)	77 (73.3%)	66 (71.7%)	65 (69.1%)
Sodium-glucose cotransporter-2 inhibitors	48 (47.5%)	56 (53.3%)	59 (64.1%)	56 (59.6%)
Vasodilators	8 (7.9%)	13 (12.4%)	10 (10.9%)	12 (12.8%)
- Long-acting nitrates	7 (6.9%)	11 (10.5%)	9 (9.8%)	8 (8.5%)
- Hydralazine	2 (2.0%)	8 (7.6%)	3 (3.3%)	8 (8.5%)
Diuretics	93 (92.1%)	98 (93.3%)	85 (92.4%)	83 (88.3%)
Antiplatelet agents	51 (50.5%)	52 (49.5%)	49 (53.3%)	49 (52.1%)
Chronic oral anticoagulation	63 (62.4%)	54 (51.4%)	59 (64.1%)	56 (59.6%)

Table 13: Post-Procedure Antiplatelet/Anticoagulation Treatment at Implant Procedure Discharge – ITT Population

	Shunt Group (N=250) n (%)	Control Group (N=258) n (%)
Antiplatelet agents, open-label (clinical)	121 (48.4)	132 (51.2)
Antiplatelet agents, study medications ^a	55 (22.0)	63 (24.4)
Chronic oral anticoagulation	158 (63.2)	150 (58.1)

a. Aspirin and clopidogrel (one or both) unless the patient was otherwise taking open-label aspirin and a platelet P2Y12 receptor inhibitor or on anticoagulation due to a clinical indication.

Roll-In TEE (Patency)

	Shunt N = 97	6 Months N = 90	12 Months N = 82	p-value
Studies analyzed				
TEE or TTE	97 (100%)	87 (97%)	75 (91%)	
TEE	86 (89%)	69 (77%)	56 (68%)	
Time to TEE, months	0 [0-0]	6.2 [5.7-6.6]	12.3 [11.9-12.9]	
Results				
Shunt patent	97 (100%)	87 (100%)	72 (100%) ^a	1.000
ΔP, mmHg	4.2 ± 2.9	5.1 ± 3.1	5.1 ± 3.9	0.316
Q, mL/min	1037 ± 385	1124 ± 417	1137 ± 463	0.384
Qp/Qs, mean ± SD			1.22 ± 0.12	

ΔP = mean interatrial pressure gradient; C_d = discharge coefficient; D_{eff} = effective diameter; D_{vc} = vena contracta diameter; Q = flow

a. Three cases with echocardiographic imaging, but no colour Doppler views of the shunt; b. Thrombus seen in left atrial appendage in one patient

ESC Heart Fail. 2024 Oct;11(5):2499-2509. doi: 10.1002/ejhf2.14859. Epub 2024 May 22.

In vivo fluid dynamics of the Ventura interatrial shunt device in patients with heart failure

Michael Pfeiffer 1, John Boehmer 1, John Gorcsan 1, Shunsuke Eguchi 1, Yoshiyuki Orihara 1, Michal Laufer Perl 2, Neal Eigler 3, William T Abraham 4, Julio Nuñez Villota 5, Elizabeth Lee 6, Antoni Bayés-Genís 7, Gil Moravsky 8, Saibal Kar 9, Michael R Zile 10, Richard Holcomb 11, Stefan D Anker 12, Gregg W Stone 13, Josep Rodés-Cabau 14, JoAnn Lindenfeld 15, Jeroen J Bax 16


Patient Sex Distribution

	ITT Total N = 508		HFrEF Total N = 206		HFpEF Total N = 302		p-value*
	n	%	n	%	n	%	
Male	319	62.8%	168	81.6%	151	50.0%	
Female	189	37.2%	38	18.4%	151	50.0%	< 0.0001

Overall Demographics US vs OUS (HFrEF)

	Total N = 206		US N = 97		OUS N = 109	
	n	%	n	%	n	%
Race: White	184	89.3%	76	78.4%	108	99.0%
Race: Black	16	7.8%	15	15.5%	1	0.9%
Ethnicity: Hispanic	25	12.1%	14	14.4%	11	10.1%

LVEF Values at Baseline Reflect Epidemiology



Events in Patients with LVEF ≤ 20 at Baseline

	Patients with LVEF ≤ 20 at Baseline	
	Shunt N = 6	Control N = 8
HFH events	2	17
Death	2	2
HT/LV events	0	2
WHF events	3	2
HF events	7	23

Changes in MR and TR – ITT Population

	Baseline				12 Months			
	Treatment (N = 250)	Control (N = 258)	Difference (95% CI)	p- value	Treatment (N = 250)	Control (N = 258)	Difference (95% CI)	p- value
Mitral Regurgitation								
Moderate or greater	20%	15%	4.9 (-1.7, 11.4)	0.15	11%	10%	0.6 (-5.2, 6.5)	0.83
Moderate	18%	14%	4.8 (-1.5, 11.2)	0.14	10%	10%	-0.3 (-5.9, 5.3)	0.91
Mod – Severe	1%	1%	0.0 (-2.3, 2.5)	1.00	1%	0.5%	1.0 (-1.2, 3.8)	0.36
Severe	0	0	0	N/A	0	0	N/A	N/A
Tricuspid Regurgitation								
Moderate or greater	20%	18%	2.7 (-4.1, 9.6)	0.43	16%	11%	4.8 (-1.7, 11.2)	0.15
Moderate	19%	16%	3.1 (-3.6, 9.7)	0.36	13%	11%	2.8 (-3.3, 9.0)	0.36
Mod – Severe	1%	0.4%	0.8 (-1.1, 3.2)	0.36	2%	0.5%	1.9 (-0.5, 5.1)	0.12
Severe	0	1%	-1.2 (-3.4, 0.4)	0.25	0	0	N/A	N/A

Changes in MR and TR – HFrEF

	Baseline				12 Months			
	Treatment (N = 101)	Control (N = 105)	Difference (95% CI)	p- value	Treatment (N = 101)	Control (N = 105)	Difference (95% CI)	p- value
Mitral Regurgitation								
Moderate or greater	24%	18%	5.7 (-5.4, 16.8)	0.32	14%	10%	4.4 (-5.2, 14.1)	0.37
Moderate	22%	15%	6.5 (-4.0, 17.1)	0.23	10%	10%	0.9 (-8.1, 10.0)	0.84
Mod – Severe	2%	3%	-0.9 (-6.5, 4.6)	1.00	3%	0	3.5 (-1.1, 9.9)	0.25
Severe	0	0	N/A	N/A	0	0	N/A	N/A
Tricuspid Regurgitation								
Moderate or greater	12%	16%	-3.9 (-13.5, 5.6)	0.42	8%	10%	-1.7 (-10.4, 7.0)	0.69
Moderate	12%	14%	-2.0 (-11.4, 7.3)	0.67	6%	10%	-4.1 (-12.2, 4.1)	0.33
Mod – Severe	0	0	N/A	N/A	2%	0	2.3 (-2.3, 8.2)	0.50
Severe	0	2%	-1.9 (-6.8, 2.0)	0.50	0	0	N/A	N/A


Event Rates and Rate Ratio for Baseline Resting Pulmonary Vascular Resistance (HFrEF)

	PVR < 2.0 Wood Units			p-value	PVR ≥ 2.0 Wood Units			p-value
	Shunt N = 44	Control N = 43	RRR (95% CI)		Shunt N = 55	Control N = 60	RRR (95% CI)	
Death, HT/LV, HFH, WHF	22 (30.2%/yr)	42 (68.5%/yr)	0.44 (0.25, 0.76)	0.0021	53 (66.8%/yr)	88 (102.2%/yr)	0.65 (0.46, 0.93)	0.0169


Cardiomyopathy (HFrEF): All Heart Failure Events (All-cause Death, LVAD/HT, All HFHs, All Out-patient Worsening HF Events)



Win Ratio vs GST vs Recurrent Events Analysis



Gail-Simon GST (with KCCQ) p-value: 0.0201



Missing Echo Data

Pentara

Missing Values but has Echo

AVISIT	TRT	HFrEFL	RAA	RVEDA	RVSA
Baseline	CONTROL ARM	0	1	0	0
Baseline	V-WAVE-SHUNT ARM	0	0	0	0
Baseline	CONTROL ARM	1	1	0	0
Baseline	V-WAVE-SHUNT ARM	1	2	0	0
Month 12	CONTROL ARM	0	4	3	3
Month 12	V-WAVE-SHUNT ARM	0	2	2	2
Month 12	CONTROL ARM	1	0	2	2
Month 12	V-WAVE-SHUNT ARM	1	2	4	4

Missing Echo

AVISIT	TRT	HFrEFL	N
Month 12	CONTROL ARM	0	15
Month 12	V-WAVE-SHUNT ARM	0	21
Month 12	CONTROL ARM	1	12
Month 12	V-WAVE-SHUNT ARM	1	14

ACC / AHA Guidelines: HF with Improved EF

“ The ACC/AHA also notes that, while the EF may now be in the "preserved" or "mid-range" category, the clinical trajectory and residual cardiac dysfunction of HFimpEF are distinct from de novo HFpEF, and the evidence base for management aligns with HFrEF rather than HFpEF. This distinction is supported by the observation that improvement in LVEF does not equate to normalization of cardiac structure or function, and ongoing therapy is necessary. ”

Outcomes in Patients with Improved LVEF

HF Events death, transplant, LVAD, all HFH, all WHF	Shunt Hazard Rate (95% CI)	Control Hazard Rate (95% CI)	Hazard Rate Ratio (95% CI)
HFrEF improved LVEF*	0.51 (0.19, 0.95) N = 20	0.79 (0.29, 1.52) N = 19	0.64 (0.17, 2.16)
HFrEF non-improved LVEF	1.04 (0.74, 1.48) N = 81	2.30 (1.56, 3.43) N = 86	0.45 (0.26, 0.77)
HFpEF	1.17 (0.95, 1.51) N = 149	0.69 (0.54, 0.93) N = 153	1.69 (1.14, 2.52)

These findings suggest that HFrEF patients with improved LVEF do not behave like de novo HFpEF patients and continue to benefit from interatrial shunting.


Site Reported SAE Heart Rhythm Events in HFrEF Stratum (MedDra Coded)

Heart Rhythm Event	Shunt (N = 101) N pts (N Events)	Control (N = 105) N pts (N Events)
Atrial Fibrillation / Flutter	6 (6)	5 (5)
Ventricular Tachycardia	8 (11)	9 (12)
Ventricular Fibrillation	0	1 (1)
Arrhythmia - other	1 (1)	1 (1)
Syncope	2 (2)	2 (2)
CRT, new	2	1
Pacemaker, new	1	0
Lead Dislodgement	0	1

RELIEVE-HF HFpEF Patients Higher Risk than REDUCE LAP HF II Patients

Baseline Characteristic	RELIEVE-HF HFpEF	REDUCE LAP HF II ¹
	Total N=302	Total N=578
Age, y	73.8	72.5
Female, %	50.0	61.4
Diabetes, %	47.7	36.3
Ischemic Cardiomyopathy, %	34.8	16.2
Atrial fibrillation, %	67.9	53.1
ARNI, %	18.2	0.0
SGLT2, %	33.8	2.6
NT-proBNP, pg/ml median	1553	444
eGFR, mL/min/1.73 m ²	47.0	56.0
NYHA functional class, %		
II	2.0	21.8
III	97.3	78.2
E/e'	15.8	12.5
LAVI, ml/m ²	40.7	31.6
TAPSE	17.0	20
Cardiac Output, L/min	4.6	5.2
PVR, WU	2.2	1.5

REDUCE LAP-HF HFH + Outpatient WHF Events Substantiate RELIEVE-HF Conclusions in HFpEF



Gustafsson F, et al. JACC Heart Fail., 2024

Adapted from Stone, et al. Circulation 2024

Replicating REDUCE non-responder outcomes gives
credence to the totality of RELIEVE-HF data

Table 11: Additional Safety Endpoints

	Shunt group (N=250)	Control group (N=258)	Relative risk or difference	P value
Secondary safety endpoints:				
MACNE* or BARC types 3 or 5 bleeding at 30 days ¹	2 (0.8%)	-	-	-
BARC types 3 or 5 bleeding at 30 days ¹	2 (0.8%)	1 (0.4%)	2.07 [0.19, 22.85] ²	0.54
MACNE* at 1 year ¹	0 (0.0%)	-	-	-
MACNE* at 2 years ¹	0 (0.0%)	-	-	-
Cerebrovascular events at 2 years, any ¹	11 (5.1%)	6 (2.5%)	1.92 [0.71, 5.18] ²	0.19
CNS infarction (stroke) ^{1,**}	7 (3.3%)	5 (2.1%)	1.46 [0.46, 4.60] ²	0.52
CNS hemorrhage (intracerebral or subarachnoid) ^{1,†}	0 (0.0%)	1 (0.5%)	-	0.33
Transient ischemic attack ¹	4 (1.9%)	1 (0.4%)	4.12 [0.46, 36.91] ²	0.17
Myocardial infarction at 2 years ¹	8 (3.8%)	13 (6.6%)	0.63 [0.26, 1.52] ²	0.30
Systemic embolization events at 2 years ¹	0 (0.0%)	0 (0.0%)	-	-
Pulmonary embolization events at 2 years ¹	2 (1.0%)	0 (0.0%)	-	0.16
Shunt implant embolization at 2 years ¹	0 (0.0%)	-	-	-

* MACNE was device-related or procedure-related.


** The 7 strokes in patients who were treated with the Shunt were classified by the CEC as being due to cerebrovascular disease (n=3), embolic due to atrial fibrillation (n=2) and undetermined (n=2). The 5 strokes in Control group patients who were treated with a placebo-procedure were classified by the CEC as being due to cerebrovascular disease (n=1), embolic due to atrial fibrillation (n=2), subarachnoid hemorrhage (n=1) and undetermined (n=1). Only one stroke occurred within 30 days of randomization, that being in the Control group.

† Does not include 1 additional patient in the placebo group with an ischemic stroke and hemorrhagic transformation.


1. Event rates were number of events (Kaplan-Meier time-to-first event estimates. Not done for MACNE as there were no events.

2. Hazard ratio [95% confidence interval].


Subgroup Analysis (HFrEF): All Heart Failure Events (All-cause Death, LVAD/HT, All HFHs, All Out-patient Worsening HF Events)



Subgroup Analysis (HFpEF): All Heart Failure Events (All-cause Death, LVAD/HT, All HFHs, All Out-patient Worsening HF Events)



Shunt-Treated Roll-in Cohort Observed vs Predicted 12-month Mortality (HFrEF)



HFrEF Roll-in vs Randomized Shunt Patient Baseline Characteristics – Key Differences

Characteristic	Roll-in Shunt	Randomized Shunt
	HFrEF	HFrEF
n	49	101
Age	68.9±11.0	69.8 ± 11.1
Female, %	8.2%	16.8%
BMI, kg/m ²	31.1±5.8	29.1 ± 5.4
Therapies		
ICD, %	40.8%	88.10%
CRT, %	42.9%	48.50%
ARNI, %	59.2%	79.20%
SGLT2, %	22%	47.50%
Echo		
LVEF, %	28.2±6.7	30.0 ± 6.4
RVFAC, %	35.0±6.6	36.8 ± 6.9
TAPSE, mm	15.4±2.8	16.1 ± 3.4
Hemodynamics		
RAP, mmHg	11.6±4.6	8.9 ± 4.2
PAP mean, mmHg	31.6±8.7	25.6 ± 7.7
PCWP, mmHg	20.7±7.6	16.4 ± 6.6
PVR, Wood units	2.5±1.2	2.3 ± 1.3