

# Quantitative Analysis of Polyaromatic Hydrocarbons, including Benzo[a]Pyrene and Naphthalene, in Tobacco, Tobacco Products, Fiber-based Matrices, and Tobacco Derived Products using Gas Chromatography-Mass Spectrometry

|  |          |
|--|----------|
| <b>Purpose</b> .....   | <b>4</b> |
| <b>Applies to</b> .....  | <b>4</b> |
| <b>General information</b> .....                                   | <b>4</b> |
| Principle of the method .....                                      | 4        |
| Method scope, measurement range, and measurement uncertainty ..... | 6        |
| <i>Literature references</i> .....                                 | 8        |
| <i>Internal reference documents</i> (available upon request) ..... | 8        |
| <b>Risk assessment and safety instructions</b> .....               | <b>8</b> |
| Summarized risk assessment .....                                   | 8        |
| Substances hazardous to the environment:.....                      | 8        |
| Flammable substances:.....   | 8        |
| <b>Equipment</b> .....   | <b>9</b> |
| Apparatus and laboratory utensils .....                            | 9        |
| GC system: .....   | 9        |
| (b) (4) .....  | 9        |
| (b) (4) .....  | 10       |
| (b) (4) .....  | 10       |
| (b) (4) .....  | 12       |
| (b) (4) .....  | 12       |
| (b) (4) .....  | 13       |
| Other equipment: .....   | 14       |
| Chemicals, reagents and solvents.....                              | 14       |
| Check samples and reference materials .....                        | 15       |
| Preparation of standards .....                                     | 15       |

|  |           |
|--|-----------|
| General information .....                          | 15        |
| Preparations .....                                 | 15        |
| Comparison of old and new standards .....          | 20        |
| (b) (4) .....                                      | 21        |
| (b) (4) .....                                      | 21        |
| (b) (4) .....                                      | 21        |
| (b) (4) .....                                      | 21        |
| <b>Sample handling .....</b>                       | <b>21</b> |
| Sample storage and preparation .....               | 21        |
| Sample amount .....                                | 21        |
| <b>Analysis .....</b>                              | <b>21</b> |
| Calibration and verification of apparatus .....    | 21        |
| Sample stability .....                             | 23        |
| Analytical procedure .....                         | 23        |
| General information .....                          | 23        |
| Sample preparation .....                           | 23        |
| Solid phase extraction with Gilson SPE robot ..... | 24        |
| Special instructions .....                         | 24        |
| <b>Documentation .....</b>                         | <b>26</b> |
| <b>Data .....</b>                                  | <b>26</b> |
| Collection and storage of data .....               | 26        |
| How to find results and raw data .....             | 26        |
| Calculations .....                                 | 26        |
| Quality assurance .....                            | 27        |
| Control chart .....                                | 27        |
| Standard curve criteria .....                      | 28        |
| Duplicate and triplicate samples .....             | 28        |
| Confirmation .....                                 | 28        |
| Concentration above calibration curve .....        | 29        |
| Reporting of analysis results .....                | 29        |
| <b>Revision history .....</b>                      | <b>29</b> |
| <b>Person responsible .....</b>                    | <b>29</b> |
| <b>Validation report .....</b>                     | <b>30</b> |
| Supporting documentation for validation .....      | 30        |
| Specificity .....                                  | 30        |
| Repeatability .....                                | 34        |
| Precision within the laboratory .....              | 35        |
| Accuracy .....                                     | 37        |
| Limit of quantification (LOQ) .....                | 37        |
| Linearity .....                                    | 38        |

---

|   |    |
|---|----|
| Robustness.....                                     | 39 |
| Measurement range and measurement uncertainty ..... | 39 |
| Conclusion.....                                     | 40 |

***Purpose***

To determine the concentration of 18 polyaromatic hydrocarbons (PAHs) including benzo[a]pyrene (B[a]P) and naphthalene in tobacco, tobacco products, fiber-based matrices, and tobacco derived products with (b) (4) and head space (b) (4), and to determine the concentration for B[a]P with (b) (4).

***Applies to***

APRS

***General information*****Principle of the method**

(b) (4)

The method is a modified CORESTA method (b) (4). The modifications are:

(b) (4)

(b) (4)

Note: All reference documents and additional information stated “available upon request” are in Swedish. They are available upon request but need to be translated into English first.

**Method scope, measurement range, and measurement uncertainty**

(b) (4)

**Table 1.** (b) (4)

|         |  |  |  |  |
|---------|--|--|--|--|
| (b) (4) |  |  |  |  |
|---------|--|--|--|--|

**Table 2. Demarcations for GC-MS analysis of B[a]P.**

| <i>Analyte</i> | <i>Calibration range<br/>(ng/ml)</i> | <i>Measurement<br/>range (ng/g<br/>as is)</i> | <i>LOQ<br/>GC-MS<br/>(ng/g as is)</i> | <i>Measurement<br/>uncertainty<br/>GC-MS (%)<br/>Duplicate<br/>samples</i> |
|----------------|--------------------------------------|---|---------------------------------------|--|
|----------------|--------------------------------------|---|---------------------------------------|--|

(b) (4)

---

(b) (4)

*Literature references*

(b) (4)

(b) (4)

*Risk assessment and safety instructions*

(b) (4)

**Summarized risk assessment**

(b) (4)

Wear suitable protective clothing, protective gloves and goggles.

All work when preparing solutions for PAHs and internal standards are to be carried out in a fume cabinet, on a draw bench or places with spot extractor.

Protective gloves should be worn, along with protective goggles.

**Substances hazardous to the environment:**

(b) (4)

**Flammable substances:**


(b) (4)



---

## *Equipment*

### **Apparatus and laboratory utensils**



(b) (4)



Quality and Environmental  
Management System

Document Title

**Quantitative Analysis of PAHs**

Part of Process

Contract Analyses APS

Document Type

Method Description

Valid From

01/29/2019

Page

11(40)

Approved By

(b) (6)

Document Publisher

(b) (6)

*Table 3.* (b) (4)

(b) (4)

**GC parameters for headspace (naphthalene analysis)**

(b) (4)

**Table 4. Naphthalene: Quantifier and qualifier (SIM).**

| Name    | Abbreviation | Time<br>Segment<br>[min] | Quantifier | Qualifier |
|---------|--------------|--------------------------|------------|-----------|
| (b) (4) |              |                          |            |           |

(b) (4)

(b) (4)

(b) (4)

**Check samples and reference materials**

(b) (4)

**Preparation of standards**

(b) (4)

(b) (4)



(b) (4)

(b) (4)

(b) (4)

(b) (4)

#### Comparison of old and new standards

When preparing standards from a new primary standard, the difference in absolute peak area should not differ by more than 10%. Evaluate by analyzing calibration standard level 4 in triplicate and alternate the new and old standards when analyzing and performing the assessment based on the mean value. This is checked by setting the mean value of the old and new peak areas in Equation 1 below.

$$ABS|\bar{x}_{new}-\bar{x}_{old}| \leq 10\% \quad 1$$

Where;

$\bar{x}_{new}$  = average absolute peak area of new standard solution

$\bar{x}_{old}$  = average absolute peak area of old standard solution

## Preparation of other solutions

(b) (4)

## Sample handling

### Sample storage and preparation

(b) (4)

### Sample amount

The minimum sample amount is 0.01 g

## Analysis

### Calibration and verification of apparatus

Before a sample series can be run, prepare the instrument.

It is important that the instrument is checked before starting a sample sequence.

1. Before testing a sample series, the retention times must be checked using calibration standard 5. If the retention times have been changed resulting in the peaks falling outside their segment or close to the start and end, the data collection time segments must then be updated in the collection method.
2. Also check that area units for IS for calibration standard 5 is consistent with the area for IS at the last calibration. The area for IS must not have reduced by more than 50%. If the reduction is greater, investigate and correct as below:
  1. Replace liner
  2. Search for leaks with a leak detector
  3. Perform new autotune
  4. Cut about a 5 cm column on the injector side
  5. Contact the person responsible for the method

6. Replace column
  7. Clean/replace the ion source
- Inject calibration standard 5 again and evaluate if more measures need to be performed.

System performance is also checked when evaluating a batch. At start-up and at the end of each sample sequence, all calibration standards are also injected, which is evaluated as follows:

- (b) (4)
- Replace liner
- Perform new autotune
- Cut about a 5 cm column on the injector side
- Contact the person responsible for the method
- Replace column
- Clean/replace the ion source

A standard curve is generated in each sequence. Before samples are analyzed in a sample sequence (following equilibration of the system), the sample sequence is injected with all calibration standards, followed by a blank and two check samples that are evaluated as set out in DESC "Control Charts and Check Samples". Once all the samples have been analyzed, the check samples are injected again, and finally all calibration standards are also injected, (also applies when analyzing naphthalene).

Each sequence is built up according to:

(b) (4)

### **Sample stability**

The shelf life of prepared samples in vials with whole septum stored in the refrigerator is one week.

### **Analytical procedure**

#### General information

(b) (4)

#### Sample preparation

(b) (4)

(b) (4)



(b) (4)

## *Documentation*

### *Data*

#### **Collection and storage of data**

(b) (4)

#### **How to find results and raw data**

(b) (4)

#### **Calculations**

Quantification is performed by using internal standards, which means it is the area ratios of analyte and internal standard that are used as a response in the preparation of the calibration curve and the quantification of the samples.

The method includes a calibration table with linear calibration curves for each analyte. The calibration table has concentrations given in ng/ml. (b) (4)

(b) (4)

(b) (4)

## Quality assurance

### Control chart

For the internal quality control one check sample is injected, both before and after the unknown samples in the sequence, on every analysis occasion (b) (4)

. The analysis of the check sample after the unknown samples is used where troubleshooting/investigations are required.

(b) (4)

The reason that only one check sample is used is that the largest contribution to the drift is probably derived from the instrument and not from the sample preparation. It is therefore important that the criteria for the calibration residuals are met.

A remark is entered as a comment, for example, when new stock solutions are prepared, and when replacing the column.

### *Measures when verification is not approved:*

- Prepare new internal standard solution
- Replace septa
- Replace liner
- Cut the column at the inlet
- Recalibrate using existing calibration solutions
- Prepare new stock and/or calibration solutions and calibrate.

Standard curve criteria

Calibration of GC instrument with calibration standards is made for each sequence.

Concentrations are given in ng/ml. (b) (4)

If these criteria are not met, the measures below are taken in the following order:

- Higher calibration levels in the batch can and should be excluded, if there are no samples in the sequence that end up above the highest remaining calibration level. At least three calibration levels must be used.
- Re-analyze the sequence following a review/troubleshooting of the cause of the failure.  
Based on needs:
  1. Replace septa
  2. Replace liner
  3. Cut the column at the inlet
  4. Run an autotune on MS
- Re-prepare and re-analyze the sample
- Prepare new stock and/or calibration and internal standard solutions and calibrate.

Duplicate and triplicate samples

When analyzing with duplicate samples, the difference between the samples should not be greater than 3 standard deviations of what was measured when determining repeatability in the matrix in question. (b) (4)

When analyzing with triplicate samples, this difference between samples should not exceed 3.3 standard deviations.

Confirmation

(b) (4)

**Concentration above calibration curve**

Samples with concentrations above the calibration curve are mainly prepared using less sample weight. (b) (4)

(b) (4). Dilute enough to ensure that the response area of the analyte is estimated to fall within the upper half of the calibration range. The internal standard makes a correction for the dilution and there is therefore no need for an adjustment to the calculations.

**Reporting of analysis results**

(b) (4)

***Revision history***

| Version | Amendments |
|---------|------------|
| (b) (4) |            |

***Person responsible***

Director APS

***Validation report*****Supporting documentation for validation**

Calculations and all raw data used in MS Excel are available upon request.

The results are calculated as ng/g as is.

| Sample  | Description | Validation |
|---------|-------------|------------|
| (b) (4) |             |            |

(b) (4)  
(b) (4)

. The measurement range is applied to the different types of matrices in the validation such as tobacco flour, cigar and non-smoking products. An expanded relative measurement uncertainty has been calculated from the combined relative errors.

(b) (4)

**Specificity**

(b) (4)

(b) (4)

(b) (4)



(b) (4)

(b) (4)

**Repeatability**

(b) (4)


(b) (4)

*Table 11. Repeatability for GC-MS analysis.*

|                | Pooled RSD% | Max difference in % for<br>duplicate samples |
|----------------|-------------|--|
| <b>(b) (4)</b> |             |  |

**Precision within the laboratory**

Precision within the laboratory has been determined by three individuals who have analyzed three replicates of three tobacco matrices at six different timepoints where the extraction solution and standard curves have been varied. CRP1.1, CRP2.1 and CRP4.1 have been analyzed. In [Table 12](#) you can see the pooled RSD % for analysis with GC-MS/MS, after compilation, and corresponding data for the GC-MS analysis are found in [Table 13](#).

*Table 12. Pooled RSD % for GC-MS/MS analysis.**Table 13.* (b) (4)

|         |  | RSD % (pooled) |  |
|---------|--|----------------|--|
| (b) (4) |  |                |  |

**Accuracy**

(b) (4)

*Table 14. Accuracy data for GC-MS/MS analysis.*

(b) (4)

*Table 15. Accuracy data for GC-MS analysis.*

(b) (4)

**Limit of quantification (LOQ)**

LOQ is calculated according to  $(10 \times \text{the concentration found})/\text{S/N ratio}$  and the root mean square (RMS) algorithm was used in the (b) (4). Data from repeatability was used as supporting documentation, in the case of absence of substances in the tested matrices, LOQ has been calculated based on accuracy data (spiked levels). For most analytes, the LOQ/reporting limit has been set slightly higher than the estimated LOQ.

In Table 16, the limits of quantification for analysis with GC-MS/MS are listed and in Table 17 for analysis with GC-MS.

*Table 16. Limit of quantification for GC-MS/MS analysis.*

(b) (4)

*Table 17. Limit of quantification for GC-MS analysis.*

|  | LOQ GC-MS<br>(ng/g as is) |
|--|---------------------------|
|  |                           |

The LOD for the method is not calculated as it is not used for reporting.

### Linearity

(b) (4)

The pooled values for the relative standard uncertainty, when using a weighted (1/x) calibration curve, are included in the determination of the measurement uncertainty and are listed in the table below.

### Robustness

(b) (4)

### Measurement range and measurement uncertainty

(b) (4)

**Table 18. Measurement uncertainty when analyzed with GC-MS/MS in % of the detected concentration.**

|         | Single sample | Duplicate samples | Triplicate samples |
|---------|---------------|-------------------|--------------------|
| (b) (4) |               |                   |                    |

**Table 19. Measurement uncertainty when analyzed with GC-MS in % of the detected concentration.**

|         | Single sample | Duplicate samples | Triplicate samples |
|---------|---------------|-------------------|--------------------|
| (b) (4) |               |                   |                    |

**Conclusion**

The method is deemed fit-for-purpose.