

**Important Information for Augmentation Patients
About Mentor MemoryGel® Silicone Gel-Filled Breast
Implants
July 2009**

[GRAPHICS TO PROVIDE COVER]

Important Information for Augmentation Patients about Mentor MemoryGel® Silicone Gel-Filled Breast Implants July 2009

Table of Contents		Page No.
GLOSSARY		4
1.	WHAT SHOULD I CONSIDER IF I AM THINKING OF GETTING SILICONE GEL-FILLED BREAST IMPLANT AUGMENTATION?.....	11
1.1.	Breast Shape.....	11
1.2.	Silicone Gel-Filled Breast Implants.....	12
1.3.	Silicone Gel-Filled Breast Implants Eligibility.....	12
1.4.	Important Factors to Consider When Choosing Silicone Gel-Filled Implants (Precautions).....	13
2.	WHAT ARE POTENTIAL BREAST IMPLANT COMPLICATIONS?	14
2.1.	Rupture.....	15
2.2.	Capsular Contracture.....	18
2.3.	Additional Surgeries (Reoperations).....	18
2.4.	Implant Removal.....	19
2.5.	Unsatisfactory Results	19
2.6.	Pain.....	20
2.7.	Changes in Nipple and Breast Sensation.....	20
2.8.	Infection	20
2.9.	Hematoma/Seroma.....	20
2.10.	Breast Feeding.....	20
2.11.	Calcium Deposits in the Tissue Around the Implant.....	21
2.12.	Extrusion.....	21
2.13.	Necrosis	21
2.14.	Delayed Wound Healing.....	21
2.15.	Breast Tissue Atrophy/Chest Wall Deformity.....	21
2.16.	Lymphadenopathy	21
2.17.	Other Reported Conditions.....	22
3.	WHAT WERE THE RESULTS OF THE MENTOR CORE STUDY FOR AUGMENTATION AND REVISION-AUGMENTATION?.....	24
3.1.	Overview of Mentor Core Study.....	25
3.2.	3-Year Follow-up Rates in Augmentation Patients.....	25
3.3.	Benefits for Augmentation Patients.....	26
3.4.	3-Year Complication Rates in Augmentation Patients.....	27
3.5.	Main Reasons for Reoperation in Augmentation Patients	28
3.6.	Main Reasons for Implant Removal in Augmentation Patients.....	30
3.7.	Other Clinical Data Findings in Augmentation Patients.....	31
4.	WHAT ARE THE SURGERY CONSIDERATIONS FOR RECEIVING BREAST IMPLANTS?.....	32
4.1.	Surgical Considerations for Breast Augmentation	32
4.1.1.	Alternatives to Breast Augmentation with Silicone Breast Implants.....	32
4.1.2.	Choosing a Surgeon	32
4.1.3.	Implant Shape and Size.....	33
4.1.4.	Surface Texturing	33
4.1.5.	Implant Placement.....	33

4.1.6.	Incision Sites.....	34
4.1.7.	Additional Procedures at the Time of Breast Augmentation.....	35
4.1.8.	Palpability.....	35
4.1.9.	Surgical Setting and Anesthesia	35
4.1.10.	Postoperative Care.....	35
4.2.	Other Factors to Consider in Revision-Augmentation Surgery.....	36
5.	WHAT SORT OF FOLLOW-UP EXAMINATIONS SHOULD I HAVE?.....	36
5.1.	Breast Self-Examinations	36
5.2.	Screening for Silent Rupture.....	36
5.2.1.	Symptomatic Rupture.....	36
5.3.	Mammography	37
6.	WHAT ARE THE TYPES OF MEMORYGEL SILICONE GEL-FILLED BREAST IMPLANTS AVAILABLE FROM MENTOR?	37
7.	HOW SHOULD I REPORT PROBLEMS WITH MY IMPLANT?	38
8.	WHAT IS DEVICE TRACKING?	38
9.	WHAT ARE MENTOR’S PRODUCT REPLACEMENT POLICY AND LIMITED WARRANTIES?	39
10.	WHAT ARE OTHER SOURCES OF ADDITIONAL INFORMATION?.....	41
	ACKNOWLEDGMENT OF INFORMED DECISION	42
	INDEX	44
	REFERENCES	46

GLOSSARY

Areola	The pigmented or darker colored area of skin surrounding the nipple of the breast.
Asymmetry	Lack of proportion of shape, size, and/or position between the two breasts.
Autoimmune disease	A disease in which the body mounts an “attack” response to its own tissues or cell types. Normally, the body’s immune mechanism is able to distinguish clearly between what is a normal substance and what is foreign. In autoimmune diseases, this system becomes defective and mounts an attack against normal parts of the body, causing tissue injury. Certain diseases such as rheumatoid arthritis, lupus, and scleroderma are considered to be autoimmune diseases.
Axillary	Pertaining to the armpit area.
Biocompatible	The condition of being compatible with living tissues or systems without being toxic.
Biopsy	The removal and examination of tissues, cells, or fluid from the body.
Body Esteem Scale (BES)	A questionnaire that asks about a person’s body image. For females, the questionnaire asks about sexual attractiveness, weight concern, and physical condition.
Breast augmentation	A surgical procedure to increase breast size. For this document, it refers to placement of a breast implant. The first time a breast implant is placed to increase breast size, it is called primary augmentation. All subsequent times the implant is replaced, it is called revision-augmentation.
Breast implant	An internal artificial device or implant intended to replace the breast.
Breast mass	A lump or body in the breast.
Breast reconstruction	A surgical procedure to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality.
Calcification	Process of hardening by calcium salts.
Capsular contracture	<ul style="list-style-type: none">• A tightening of the tissue capsule surrounding an

implant, resulting in firmness or hardening of the breast and in squeezing of the implant if severe. Capsular contracture is classified by Grades. Grades III or IV are the most severe. Grade III often results in the need for additional surgery (reoperation) because of pain and possibly abnormal appearance. Grade IV usually results in the need for additional surgery (reoperation) because of pain and unacceptable appearance. Capsular contracture Grade II may also result in the need for additional surgery. Capsular contracture is a risk for implant rupture. Below is a description of each Grade.

- Grade I – Normally soft and natural appearance
- Grade II – A little firm, but breast looks normal
- Grade III – More firm than normal, and looks abnormal (change in shape)
- Grade IV – Hard, obvious distortion, and tenderness with pain

Capsule	Scar tissue that forms around the breast implant. Sometimes this capsule squeezes the implant, resulting in capsular contracture.
Capsulectomy	Surgical removal of the scar tissue capsule around the implant.
Capsulorrhaphy	Surgical stitching of a tear in the scar tissue capsule around the implant.
Capsulotomy (closed)	An attempt to break the scar tissue capsule around the implant by pressing or pushing on the outside of the breast. This method does not require surgery but is a known risk for rupture of the implant and is contraindicated.
Capsulotomy (open)	Surgical incision into the scar tissue capsule around the implant.
Congenital anomaly	An abnormal development in part of the body.
Connective tissue disease/disorder (CTD)	A disease, group of diseases, or conditions affecting connective tissue, such as muscles, ligaments, skin, etc. and/or the immune system. Connective tissue diseases (“CTDs”) that involve the immune system include autoimmune diseases such as rheumatoid arthritis, lupus, and scleroderma.
Contraindication	A use that is improper and should not be followed.

Failure to follow contraindications identified in the labeling could cause serious harm.

Contralateral	Opposite side.
Core Study	The primary clinical study of augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients that supported the approval of the premarket approval (PMA) application. Safety and effectiveness data are collected yearly through 10 years, with the follow-up from years 4 through 10 being performed as part of a postapproval Core Study.
Delayed wound healing	Delayed progress in the healing of an opened wound.
Displacement	Movement of the implant from the usual or proper place.
Epidemiological	Relating to the science of explaining the relationships of factors that determine disease frequency and distribution.
Extracapsular rupture	A type of rupture in which the silicone gel is outside of the scar tissue capsule surrounding the implant.
Extrusion	Skin breakdown with the pressing out of the implant through the surgical wound or skin.
Fibromyalgia	A disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue.
Fibrous tissues	Connective tissues composed mostly of fibers.
Granuloma	A lump or mass made of inflammatory cells surrounding a foreign substance due to longstanding inflammation.
Hematoma	A collection of blood within a space.
Hypertrophic scarring	An enlarged scar remaining after the healing of a wound.
Immune response	A bodily response to the presence of a foreign substance.
Infection	Invasion with microorganisms (for example, bacteria, viruses). An infection usually results in fever, swelling, redness, and/or pain.
Inflammation	The response of the body to infection or injury that is characterized by redness, swelling, warmth, pain, and/or loss of function.

Inframammary	Below the breast.
Inframammary fold	The crease at the base of the breast and the chest wall.
Inframammary incision	An incision made in the fold below the breast.
Inpatient surgery	A surgical procedure in which the patient is required to stay overnight in the hospital.
Intracapsular rupture	A type of rupture in which the silicone gel remains inside the scar tissue capsule surrounding the implant.
Lactation	The production and secretion of milk by the breast glands.
Low molecular weight silicones	Components of silicone of smaller molecular weight that may bleed (leak) out of silicone gel.
Lymphadenopathy	Enlargement of the lymph node(s).
Malposition	Implant malposition or displacement is when the implant is not in the correct spot in the breast. This could have been due to incorrect placement of the implant during the surgery or due to shifting of the implant position over time.
MRI	Magnetic resonance imaging. A radiographic examination that currently has the best ability to detect rupture of silicone gel-filled breast implants.
Mammary	Pertaining to the breast.
Mammography	<p>A type of x-ray examination of the breasts used for detection of cancer.</p> <p>A <u>screening mammogram</u> is an X-ray of the breast used to detect breast changes in women who have no signs or symptoms of breast cancer.</p> <p>A <u>diagnostic mammogram</u> is an X-ray of the breast that is used to check for breast cancer after a lump or other sign or symptom of breast cancer has been found. A diagnostic mammogram also may be used to evaluate changes found during a screening mammogram, or to view breast tissue when it is difficult to obtain a screening mammogram because of special circumstances, such as the presence of breast implants. A diagnostic mammogram takes longer than a screening mammogram because it involves more X-rays in order to</p>

obtain views of the breast from several angles.

Mammoplasty	Plastic surgery of the breast.
Mastopexy	Plastic surgery to move sagging breasts into a more elevated position.
Metastatic Disease	Spreading of cancer cells from the original site to other parts of the body.
Migration	Movement of silicone materials outside the breast implant.
Necrosis	Death of cells or tissues.
Outpatient surgery	A surgical procedure in which the patient is not required to stay in the hospital overnight.
Palpate	To feel with the hand.
Palpability	The ability to feel the implant.
Pectoralis	Major muscle of the chest.
Periareolar	Around the darkened or pigmented area surrounding the nipple of the breast.
Plastic surgery	Surgery intended for the improvement of appearance of the body.
Postoperatively	After surgery.
Primary breast augmentation	The first time a breast implant is placed for the purpose of breast augmentation.
Ptosis	Breast sagging that is usually the result of normal aging, pregnancy, or weight loss.
Reoperation	An additional surgery after your first breast implantation.
Revision-Augmentation	Refers to the correction or improvement of a primary augmentation. In the context of this document, it refers to surgical removal and replacement of breast implants that were placed originally for primary breast augmentation.
Rheumatological Disease/Disorder	A variety of diseases involving connective tissue structures of the body, especially the joints and fibrous tissue. These diseases are often associated with pain,

inflammation, stiffness, and/or limitation of motion of the affected parts. Can include autoimmune diseases. Fibromyalgia is a rheumatological disorder.

Rosenberg Self Esteem Scale	A questionnaire that measures self esteem. It consists of 10 statements related to overall feelings of self-worth or self-acceptance. The items are answered on a four-point scale ranging from strongly agree to strongly disagree.
Rupture	A tear or hole in the implant shell. Silicone implant ruptures may be silent or symptomatic. Ruptures can be intracapsular or extracapsular.
Saline	A solution that is made up of water and a small amount of salt.
Scar revision	A surgical procedure to improve the appearance of a scar.
Seroma	A build-up of the watery portion of the blood in a tissue location.
SF-36 Scale	A questionnaire intended to measure health-related quality of life. It includes questions that measure physical, mental, and social health. It assesses limitations in physical activities because of health problems; limitations in social activities because of physical or emotional problems; limitations in usual role activities because of physical health problems; bodily pain; general mental health; limitations in usual role activities because of emotional problems; vitality (energy and fatigue); and general health perceptions.
Silicone elastomer	A type of silicone that has elastic properties similar to rubber.
Silent rupture	A breast implant rupture without symptoms and which is not apparent except through appropriate imaging techniques such as MRI. Most silicone breast implant ruptures are silent. (see symptomatic rupture below)
Subglandular placement	Placement of a breast implant underneath and within the breast glands but on top of the chest muscle.
Submuscular placement	Placement of a breast implant wholly or partially underneath the chest muscle.
Surgical incision	A cut made to body tissue during surgery.

Symptom	Any perceptible change in the body or its functions that indicates disease or a phase of a disease.
Symptomatic	Any evidence or sign of disease or disorder reported by the patient.
Symptomatic rupture	A breast implant rupture that is associated with symptoms (such as lumps, persistent pain, swelling, hardening, or change in implant shape). Some silicone breast implant ruptures are symptomatic, but most are silent.
Tennessee Self Concept Scale (TSCS)	A questionnaire that evaluates how the patient sees herself and what she does, likes, and feels. The scale is intended to summarize an individual's feeling of self-worth, the degree to which the self-image is realistic, and whether or not that self-image is normal. It also measures the following aspects of how the patient feels about herself: moral-ethical, social, personal, physical, and family, identity, behavior, and self-satisfaction.

Important Information for Augmentation Patients about Mentor MemoryGel® Silicone Gel-Filled Implants

1. What Should I Consider if I am Thinking of Getting Silicone Gel-Filled Breast Implant Augmentation?

The purpose of this brochure is to help you in making an informed decision about having breast implants for augmentation (breast enlargement) or breast revision-augmentation (replacement) surgery. This brochure is not intended to replace consultation with your surgeon. This educational brochure is set up to provide you with information about risks and benefits of Mentor silicone gel-filled (MemoryGel®) breast implants. It is important that you read this entire brochure because you need to understand the risks and benefits and to have realistic expectations of the outcome of your surgery.


Please read this entire brochure carefully, and if you have any questions or there are things you do not understand, please discuss them with your surgeon before making any decisions. As part of your decision, both you and your surgeon will be required to sign the last page of this brochure to confirm your understanding of what you have read.

You should wait at least 1-2 weeks after reviewing and considering this information before deciding whether to have primary breast augmentation surgery. In the case of a revision-augmentation; however, your surgeon may find it medically necessary to perform surgery sooner.

1.1. Breast Shape

The breast consists of milk ducts and glands, surrounded by fatty tissue that provides its shape and feel. The chest muscle (pectoralis major muscle) is located beneath the breast. Factors such as pregnancy (when milk glands are temporarily enlarged), rapid weight loss, and the effects of gravity as you age, combine to stretch the skin, which may cause the breast to droop or sag.

It is important to realize that implants are used to make the breast larger. The implants alone may not adequately lift the breast, or correct the effects of pregnancy, weight loss, or skin stretching. Your surgeon may suggest additional procedures at the time of the breast augmentation, such as a procedure to help achieve improved breast lift (mastopexy).



1.2. Silicone Gel-Filled Breast Implants

A breast implant is a sac (implant shell) of silicone rubber (elastomer) filled with silicone gel, which is surgically implanted either under your breast tissue or under your chest muscle.



1.3. Silicone Gel-Filled Breast Implants Eligibility

Mentor MemoryGel Silicone Gel-Filled Breast Implants are indicated for females for the following uses (procedures):

- **Breast augmentation for women at least 22 years old.** Breast augmentation includes primary breast augmentation to increase the breast size, as well as revision surgery to correct or improve the result of a primary breast augmentation surgery.
- **Breast reconstruction.** Breast reconstruction includes primary reconstruction to replace breast tissue that has been removed due to cancer or trauma or that has failed to develop properly due to a severe breast abnormality. Breast reconstruction also includes revision surgery to correct or improve the result of a primary breast reconstruction surgery. (A separate patient brochure is available for and should be read for breast reconstruction.)

Contraindications

Breast implant surgery should not be performed in:

- Women with active infection anywhere in their body.
- Women with existing cancer or pre-cancer of their breast who have not received adequate treatment for those conditions.
- Women who are currently pregnant or nursing.

1.4. Important Factors to Consider When Choosing Silicone Gel-Filled Implants (Precautions)

- Breast implants are not lifetime devices, and breast implantation is likely not a one-time surgery. You will likely need additional unplanned surgeries on your breasts because of complications or unacceptable cosmetic outcomes. These additional surgeries can include implant removal with or without replacement, or they can include other surgical procedures. When you have your implants replaced (revision-augmentation), your risk of future complications increases compared to first time (primary) augmentation surgery, so you should also review the complication rates for revision-augmentation patients to see what future risk rates you may experience.
- Safety and effectiveness have not been established in patients with the following:
 - Autoimmune diseases (for example, lupus and scleroderma).
 - A weakened immune system (for example, currently taking drugs that weaken the body's natural resistance to disease).
 - Conditions that interfere with wound healing and blood clotting.
 - Reduced blood supply to breast tissue.
 - Radiation to the breast following implantation.
 - Clinical diagnosis of depression or other mental health disorders, including body dysmorphic disorder and eating disorders. Please discuss any history of mental health disorders with your surgeon prior to surgery. Patients with a diagnosis of depression, or other mental health disorders, should wait until resolution or stabilization of these conditions prior to undergoing breast implantation surgery.
- Many of the changes to your breast following implantation are irreversible (cannot be undone). If you later choose to have your implant(s) removed and not replaced, you may experience unacceptable dimpling, puckering, wrinkling, or other cosmetic changes of the breast, which can be permanent.
- Breast implants may affect your ability to breast feed, either by reducing or eliminating milk production.
- Rupture of a silicone gel-filled breast implant is most often silent. This means that neither you nor your surgeon will know that your implants have a rupture most of the time. In fact, the ability of a physical examination by a plastic surgeon who is familiar with breast implants to detect silicone breast implant rupture is 30%¹ compared to 89% for MRI.² You may need regular screening MRI examinations over your lifetime in order to determine if silent rupture is present. You should have your first MRI at 3 years after your initial implant surgery and then every 2 years thereafter. The cost of MRI screening may exceed the cost of your initial surgery over your lifetime. This cost, which may not be covered by your insurance, should be considered in making your decision.
- If implant rupture is noted on MRI, you should have the implant removed, with or without replacement.

- With breast implants, routine screening mammography for breast cancer will be more difficult. If you are of the proper age for mammography screening, you should continue to undergo routine mammography screening as recommended by your primary care physician. The implant may interfere with finding breast cancer during mammography. Because the breast and implant are squeezed during mammography, an implant may rupture during the procedure. More x-ray views are necessary for women with breast implants; therefore, you will receive more exposure to radiation. However, the benefit of having the mammogram to find cancer outweighs the risk of the additional x-rays. Be sure to inform the mammography technologist that you have implants.
- You should perform an examination of your breasts every month for cancer screening; however, this may be more difficult with implants. You should ask your surgeon to help you distinguish the implant from your breast tissue.
- You should perform an examination of your breasts for the presence of lumps, persistent pain, swelling, hardening, or change in implant shape, which may be signs of symptomatic rupture of the implant. These should be reported to your surgeon and possibly evaluated with an MRI to screen for rupture.
- After undergoing breast augmentation surgery (either primary or revision), your health insurance premiums may increase, your insurance coverage may be dropped, and/or future coverage may be denied. Treatment of complications may not be covered as well. You should discuss the complete extent of your insurance coverage with your insurance company before undergoing surgery.
- You should inform any other doctor who treats you of the presence of your implants to minimize the risk of damage to the implants.
- Mentor will continue its ongoing Core Study through 10 years to further evaluate the long-term safety and effectiveness of these products. In addition, Mentor has initiated a separate, 10-year postapproval study in the U.S. and Canada to address specific issues for which the Mentor Core Study was not designed to fully answer, as well as to provide a real-world assessment of some endpoints. The endpoints in the large postapproval study include long-term local complications, connective tissue disease (CTD), CTD signs and symptoms, neurological disease, neurological signs and symptoms, offspring issues, reproductive issues, lactation issues, cancer, suicide, mammography issues, and MRI compliance and results. Mentor will update its labeling as appropriate with the results of these two studies. You should also ask your surgeon if he/she has any available updated clinical information.

2. What are Potential Breast Implant Complications?

Undergoing any type of surgical procedure involves risks (some serious) such as the effects of anesthesia, infection, swelling, redness, bleeding, pain, and even death, which need to be balanced against the benefits of the breast augmentation surgery.

There are potential complications specific to breast implant surgery and breast implants, as described below. It should also be noted that the cited references include data from augmentation and/or reconstruction patients, as well as from a variety of manufacturers and implant models.

2.1. Rupture

Breast implants are not lifetime devices. Breast implants rupture when the shell develops a tear or hole. Rupture can occur at any time after implantation, but it is more likely to occur the longer the implant is implanted. The following things may cause your implant to rupture:

- damage by surgical instruments;
- stressing the implant during implantation and weakening it;
- folding or wrinkling of the implant shell;
- excessive force to the chest (for example, during closed capsulotomy), which is contraindicated);
- trauma;
- compression during mammographic imaging; and
- severe capsular contracture.

Breast implants may also simply wear out over time. Laboratory studies have identified some of the types of rupture for Mentor's product; however, it is not known whether these tests have identified all causes of rupture. These laboratory studies will continue postapproval.

Silicone gel-filled implant ruptures are most often silent. (MRI examination is currently the best method to screen for silent rupture.) This means that most of the time neither you nor your plastic surgeon will know if the implant has a tear or hole in the shell. The need to screen for rupture is why MRI is recommended at 3 years and then every 2 years thereafter. Sometimes there are symptoms associated with gel implant rupture. These symptoms include hard knots or lumps surrounding the implant or in the armpit, change or loss of size or shape of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast. If you are experiencing these symptoms, you should discuss the need for an MRI examination with your physician.

When MRI identifies a rupture, or if your surgeon determines that you have signs or symptoms of rupture, you should have the implant and any gel removed with or without replacement of the implant. It also may be necessary to remove the tissue capsule as well as the implant, which will involve additional surgery, with associated costs. If you have symptoms such as breast hardness, a change in breast shape or size, and/or breast pain, you should have an MRI to determine whether rupture is present.^{3,4}

There are also consequences of rupture. If rupture occurs, silicone gel may either remain within the scar tissue capsule surrounding the implant (intracapsular rupture), move outside the capsule (extracapsular rupture), or gel may move beyond the breast (migrated gel). There is also a possibility that rupture may progress from

intracapsular to extracapsular and beyond. There have also been health consequences reported in the literature. See below for details.

Rupture Information on Mentor MemoryGel Breast Implants

In Mentor's Core Study, rupture was assessed for patients who had scheduled MRIs to screen for silent rupture (i.e., part of the MRI cohort) and those who were not assessed for rupture by MRI (i.e., part of the non-MRI cohort).

- For primary augmentation patients in the MRI cohort, the rupture rate was 0.5% through 3 years. This means that through 3 years, 1 of every 200 augmentation patients had at least 1 ruptured breast implant.
- There was one primary augmentation patient in the Mentor Core Study with a suspected implant rupture detected via MRI, which has not been confirmed with examination of the implant following removal.
- For revision patients in the MRI cohort, the rupture rate was 7.7% after 3 years. This means that through 3 years, 1 of every 13 revision-augmentation patients had at least 1 ruptured breast implant.
- All of these implant ruptures were silent and were detected via MRI. There was one revision-augmentation patient who had removal of her breast implants after MRI and both implants were ruptured. The other implant ruptures have not yet been confirmed with removal and examination of the implant.
- The rupture rate beyond 3 years in Mentor's Core Study continues to be investigated. The Mentor Core Study will continue to monitor patients for rupture.

There were no ruptures reported in the non-MRI cohorts for either the primary augmentation or revision-augmentation patients of the Mentor Core Study through 3 years. Across all patients in the Mentor Core Study, of the 8 implants reported as ruptured, 4 showed intracapsular gel and 4 showed extracapsular gel on MRI (3 implants with extracapsular gel were in 2 revision-augmentation patients and 1 was in a primary reconstruction patient). For one of these implants with extracapsular gel, this was a confirmed case in which the device was explanted and the intracapsular gel rupture progressed into an extracapsular gel rupture as shown by MRIs at approximately 10 months and approximately 2 years. There were no cases of migrated gel.

Further rupture rate information on Mentor MemoryGel implants is provided from an unpublished European study known as the U.K. Sharpe and Collis Study. Silent rupture was assessed by a single MRI on 101 augmentation patients implanted with textured Mentor MemoryGel implants by one surgeon. The average age of the implants was approximately 9 years. Silent rupture was found in approximately 10% of these augmentation patients, which includes one patient for which the device was not explanted to confirm rupture. There were no cases of extracapsular rupture or migrated gel.

Additional information on rupture will be collected through Mentor's postapproval Core Study and large postapproval study.

Additional Information on Consequences of Rupture from Literature

Studies of Danish women evaluated with MRI involving a variety of manufacturers and implant models showed that about three-fourths of implant ruptures are intracapsular and the remaining one-fourth is extracapsular.⁵ Additional studies of Danish women indicate that over a 2-year period, about 10% of the implants with intracapsular rupture progressed to extracapsular rupture as detected by MRI.⁴ This means that for women with silicone gel rupture within the scar tissue capsule detected via MRI after 2 years, 1 in 10 of these women had progression of the gel outside the scar tissue capsule. Approximately half of the women whose ruptures had progressed from intracapsular to extracapsular reported that they experienced trauma to the affected breast during this time period or had undergone mammography. In the other half, no cause was given. In the women with extracapsular rupture, after 2 years, the amount of silicone seepage outside the scar tissue capsule increased for about 14% of these women. This means that for 100 women with silicone gel rupture outside the scar tissue capsule, the amount of gel outside the scar tissue capsule increased for 14 women 2 years later. This type of information pertains to a variety of silicone implants from a variety of manufacturers and implant models, and it is not specific to Mentor implants.

Below is a summary of information related to the health consequences of implant rupture, which have not been fully established. These reports were in women who had implants from a variety of manufacturers and implant models.

- Local breast complications reported in the published literature that were associated with rupture include breast hardness, a change in breast shape or size, and breast pain.⁴ These symptoms are not specific to rupture, as they also are experienced by women who have capsular contracture.
- There have been rare reports of gel movement to nearby tissues such as the chest wall, armpit, or upper abdominal wall, and to more distant locations down the arm or into the groin. This has led to nerve damage, granuloma formation (see glossary) and/or breakdown of tissues in direct contact with the gel in a few cases. There have been reports of silicone presence in the liver of patients with silicone breast implants. Movement of silicone gel material to lymph nodes in the armpit (axilla) also has been reported, even in women without evidence of rupture, leading to enlargement of the lymph nodes (lymphadenopathy).⁶
- Concerns have been raised over whether ruptured implants are associated with the development of connective tissue or rheumatic diseases and/or symptoms such as fatigue and fibromyalgia.^{7,8,9,10} A number of epidemiology studies have evaluated large populations of women with breast implants from a variety of manufacturers and implant models. These studies do not, taken together, support a significant association of breast implants with a typical, diagnosed rheumatic disease. Other than one small study,⁹ these studies do not distinguish whether the women had ruptured or intact implants.

2.2. Capsular Contracture

The scar tissue (capsule) that normally forms around the implant may tighten over time and compress the implant, making it feel firm and leading to what is called capsular contracture. Capsular contracture may be more common following infection, a build-up of the water portion of the blood (seroma), and a collection of blood (hematoma), and the chance of it happening may increase over time. Capsular contracture occurs more commonly in revision-augmentation patients than in primary augmentation patients. Because you may have your initial implants replaced, you should be aware that your risk of capsular contracture increases with revision-augmentation. Capsular contracture is a risk factor for implant rupture, and it is the most common reason for reoperation.

Symptoms of capsular contracture range from mild firmness and mild discomfort to severe pain, distorted shape of the implant, and ability to feel the implant (palpability). Capsular contracture is graded into 4 levels depending on its severity. Grades III or IV are considered severe and often additional surgery is needed to correct these grades:

- Grade I: the breast is normally soft and looks natural
- Grade II: the breast is a little firm but looks normal
- Grade III: the breast is firm and looks abnormal
- Grade IV: the breast is hard, painful, and looks abnormal

In Mentor's Core Study, for women receiving augmentation implants for the first time, the risk of severe capsular contracture was 8% through 3 years. This means that 8 out of every 100 women who received Mentor implants for primary breast augmentation had severe capsular contracture at least once during the first 3 years after receiving the implants.

For women receiving revision-augmentation implants, the risk of severe capsular contracture was 19% through 3 years. This means that 19 out of every 100 women who received Mentor implants for breast revision-augmentation had severe capsular contracture at least once during the first 3 years after receiving the implants.

Additional surgery may be needed in cases where pain and/or firmness are severe. This surgery ranges from removal of the implant capsule tissue, to removal and possible replacement of the implant itself. This surgery may result in loss of your breast tissue. Capsular contracture may happen again after these additional surgeries. Capsular contracture may increase the risk of rupture.³

2.3. Additional Surgeries (Reoperations)

You should assume that you will need to have additional surgeries (reoperations). In the Mentor Core Study, the reoperation rate was 15% for primary augmentation patients, which means that 15 out of every 100 women who received Mentor implants for primary augmentation had a reoperation during the first 3 years after receiving the implants. The reoperation rate was 28% for revision-augmentation patients, which means that 28 out of every 100 women who received Mentor implants for revision-augmentation had a reoperation during the first 3 years after receiving the implants.

Patients may decide to change the size or type of their implants, requiring additional surgery. Problems such as rupture, capsular contracture, asymmetry (lack of proportion of shape, size and/or position between the two breasts), hypertrophic scarring (irregular, raised scar), infection, and shifting can require additional surgery. Summary tables are provided in Section 3.5 that describe the reasons for performing additional surgeries in the Mentor Core Study. For women receiving primary augmentation implants, the three most common reasons for reoperation were severe capsular contracture, patient request for size/style change, and hematoma/seroma. For women receiving revision-augmentation implants, the three most common reasons for additional surgery were severe capsular contracture, patient request for style/size change, and biopsy.

2.4. Implant Removal

Because these are not lifetime devices, the longer you have your implants, the more likely it will be for you to have them removed for any reason, either because of dissatisfaction, an unacceptable cosmetic result, or a complication such as severe capsular contracture. Having your implants removed and replaced increases your chances of getting future complications.

For women receiving primary augmentation implants in Mentor's Core Study, 5% had their implants removed at least once through 3 years. Patient choice and severe capsular contracture were the most common reasons for implant removal. For women receiving revision-augmentation implants in Mentor's Core Study, 12% had their implants removed at least once through 3 years. The most common reasons were patient choice and severe capsular contracture.

Most women who have their implants removed have them replaced with new implants, but some women do not. If you choose not to replace your implants, you may have cosmetically unacceptable dimpling, puckering, wrinkling, and/or other potentially permanent cosmetic changes of the breast following removal of the implant. Even if you have your implants replaced, implant removal may result in loss of your breast tissue. Also, implant replacement increases your risks of future complications. For example, the risks of severe capsular contracture and reoperation increase for patients with implant replacement compared to first time placement. You should consider the possibility of having your implants replaced and its consequences when making your decision to have implants.

2.5. Unsatisfactory Results

Unsatisfactory results such as wrinkling, asymmetry, implant displacement (shifting), incorrect size, unanticipated shape, implant palpability, scar deformity, and/or hypertrophic scarring, may occur. Some of these results may cause discomfort. Pre-existing asymmetry may not be entirely correctable by implant surgery. Revision surgery may be recommended to maintain patient satisfaction, but carries additional considerations and risks. Selecting an experienced plastic surgeon may minimize, but not necessarily prevent, unsatisfactory results.

2.6. Pain

Pain of varying intensity and length of time may occur and persist following breast implant surgery. In addition, improper size, placement, surgical technique, or capsular contracture may result in pain. You should tell your surgeon about significant pain or if your pain persists.

2.7. Changes in Nipple and Breast Sensation

Feeling in the nipple and breast can increase or decrease after implant surgery. The range of changes varies from intense sensitivity to no feeling in the nipple or breast following surgery. While some of these changes can be temporary, they can also be permanent, and may affect your sexual response or your ability to nurse a baby. (See the paragraph on breast feeding below.)

2.8. Infection

Infection can occur with any surgery or implant. Most infections resulting from surgery appear within a few days to weeks after the operation. However, infection is possible at any time after surgery. In addition, breast and nipple piercing procedures may increase the possibility of infection. Infections in tissue with an implant present are harder to treat than infections in tissue without an implant. If an infection does not respond to antibiotics, the implant may have to be removed, and another implant may be placed after the infection is resolved (cleared up). As with many other surgical procedures, in rare instances, toxic shock syndrome has been noted in women after breast implant surgery, and it is a life-threatening condition. Symptoms include sudden fever, vomiting, diarrhea, fainting, dizziness, and/or sunburn-like rash. You should contact your doctor immediately for diagnosis and treatment if you have these symptoms.

2.9. Hematoma/Seroma

Hematoma is a collection of blood within the space around the implant, and a seroma is a build-up of fluid around the implant. Having a hematoma and/or seroma following surgery may result in infection and/or capsular contracture later on. Symptoms from a hematoma or seroma may include swelling, pain, and bruising. If a hematoma or seroma occurs, it will usually be soon after surgery. However, this can also occur at any time after injury to the breast. While the body absorbs small hematomas and seromas, some will require surgery, typically involving draining and potentially placing a surgical drain in the wound temporarily for proper healing. A small scar can result from surgical draining. Implant rupture also can occur from surgical draining if there is damage to the implant during the draining procedure.

2.10. Breast Feeding

Breast feeding difficulties have been reported following breast surgery, including breast reduction and augmentation. If your surgeon uses an incision around the colored portion surrounding the nipple (periareolar surgical approach), it may further increase the chance of breast feeding difficulties.

2.11. Calcium Deposits in the Tissue Around the Implant

Calcium deposits (calcifications) can form in the tissue capsule surrounding the implant. Symptoms may include and firmness. Deposits of calcium can be seen on mammograms and can be mistaken for possible cancer, resulting in additional surgery for biopsy and/or removal of the implant to distinguish calcium deposits from cancer. If additional surgery is necessary to examine and/or remove calcifications, this may cause damage to the implants. Calcium deposits also occur in women who undergo breast reduction procedures, in patients who have had hematoma formation, and even in the breasts of women who have not undergone any breast surgery. The occurrence of calcium deposits increases significantly with age.

2.12. Extrusion

Extrusion is when the breast implant comes through your skin. This may occur, for example, when your wound has not closed or when breast tissue covering your implants weakens. Radiation therapy has been reported to increase the likelihood of extrusion. Extrusion requires additional surgery and possible removal of the implant, which may result in additional scarring and/or loss of your breast tissue.

2.13. Necrosis

Necrosis is the death of cells or tissues. This may prevent or delay wound healing and require surgical correction, which may result in additional scarring and/or loss of your breast tissue. Implant removal may also be necessary. Factors associated with increased necrosis include infection, use of steroids, smoking, chemotherapy, radiation, and excessive heat or cold therapy.

2.14. Delayed Wound Healing

Some patients may experience a prolonged wound healing time. Delayed wound healing may increase the risk of infection, extrusion, and necrosis. Depending on the type of surgery or the incision, wound healing times may vary. Smoking may interfere with the healing process. You should contact your surgeon immediately if your wound does not heal within the period of time he/she has discussed with you.

2.15. Breast Tissue Atrophy/Chest Wall Deformity

The pressure of the breast implant may cause breast tissue thinning (with increased implant visibility and palpability) and chest wall deformity. This can occur while implants are still in place or following implant removal without replacement. Either of these conditions may result in additional surgeries and/or unacceptable dimpling/puckering of the breast.

2.16. Lymphadenopathy

Lymphadenopathy is a chronic enlargement of the lymph nodes. A lymph node is a round mass of tissue which makes cells as part of your immune system. The lymph nodes in the armpit (axilla) drain the breast area of fluid. Sometimes the enlarged lymph nodes are painful. If they become too large or painful, the lymph node(s) may

need to be surgically removed. Painful and/or enlarged lymph nodes should be reported to your doctor.

Literature reports associate lymphadenopathy with both intact and ruptured silicone breast implants. One study reported that armpit lymph nodes from women with both intact and ruptured silicone gel implants had abnormal tissue reactions, granulomas, and the presence of silicone.⁶ These reports were in women who had implants from a variety of manufacturers and implant models.

2.17. Other Reported Conditions

There have been reports in the literature of other conditions in women with silicone gel-filled breast implants. Many of these conditions have been studied to evaluate their potential association with breast implants. Although no cause and effect relationship has been established between breast implants and the conditions listed below, you should be aware of these reports. Furthermore, there is the possibility of risks, yet unknown, which in the future could be determined to be associated with breast implants.

- **Connective Tissue Disease (CTD)**

Connective tissue diseases include diseases such as lupus, scleroderma, and rheumatoid arthritis. Fibromyalgia, another CTD, is a disorder characterized by chronic pain in the muscles and soft tissues surrounding joints, with tenderness at specific sites in the body. It is often accompanied by fatigue. There have been a number of published epidemiological studies that have looked at whether having a breast implant is associated with having a typical or defined connective tissue disease. The study size needed to conclusively rule out a smaller risk of connective tissue disease among women with silicone gel-filled breast implants would need to be very large.^{3,8,9,10,11,12,13,14,15,16} The published studies taken together show that breast implants are not significantly associated with a risk of developing a typical or defined connective tissue disease.^{3,12,13,14} These studies do not distinguish between women with intact and ruptured implants. Only one study evaluated specific connective tissue disease diagnoses and symptoms in women with silent ruptured versus intact implants, but it was too small to rule out a small risk.⁹

- **CTD Signs and Symptoms**

Literature reports have also been made associating silicone breast implants with various rheumatological signs and symptoms such as fatigue, exhaustion, joint pain and swelling, muscle pain and cramping, tingling, numbness, weakness, and skin rashes. Scientific expert panels and literature reports have found no evidence of a consistent pattern of signs and symptoms in women with silicone breast implants.^{3,7,17,18,19} Having these rheumatological signs and symptoms does not necessarily mean you have a connective tissue disease; however, you should be aware that you may experience these signs and symptoms after undergoing breast implantation. If you notice an increase in these signs or symptoms, you should consider seeing a rheumatologist to determine whether these signs or symptoms are due to a connective tissue disorder or autoimmune disease.

- **Cancer**

Breast Cancer – Reports in the medical literature indicate that patients with breast implants are not at a greater risk than those without breast implants for developing breast cancer.^{20,21,22,23,24} Some reports have suggested that breast implants may interfere with or delay breast cancer detection by mammography and/or biopsy ; however, other reports in the published medical literature indicate that breast implants neither significantly delay breast cancer detection nor adversely affect cancer survival of women with breast implants.^{20,24,25,26,27}

Brain cancer – One recent study has reported an increased incidence of brain cancer in women with breast implants as compared to the general population.²⁸ The incidence of brain cancer, however, was not significantly increased in women with breast implants when compared to women who had other plastic surgeries. Another recently published review of four large studies in women with cosmetic implants concluded that the evidence does not support an association between brain cancer and breast implants.²⁹

Respiratory/lung cancer – One study has reported an increased incidence of respiratory/lung cancer in women with breast implants.²⁸ Other studies of women in Sweden and Denmark have found that women who get breast implants are more likely to be current smokers than women who get breast reduction surgery or other types of cosmetic surgery.^{30,31,32}

Cervical/vulvar cancer – One study has reported an increased incidence of cervical/vulvar cancer in women with breast implants.²⁸ The cause of this increase is unknown.

Other cancers – One study has reported an increased incidence of stomach cancer and leukemia in women with breast implants compared to the general population.²⁸ This increase was not significant when compared to women who had other types of plastic surgeries.

- **Neurological Disease, Signs, and Symptoms**

Some women with breast implants have complained of neurological symptoms (such as difficulties with vision, sensation, muscle strength, walking, balance, thinking or remembering things) or diseases (such as multiple sclerosis), which they believe are related to their implants. A scientific expert panel report found that the evidence for a neurological disease or syndrome caused by or associated with breast implants is insufficient or flawed.³

- **Suicide**

In several studies, a higher incidence of suicide was observed in women with breast implants.^{33,34,35,36} The reason for the observed increase is unknown, but it was found that women with breast implants had higher rates of hospital admission due to psychiatric causes prior to surgery, as compared with women who had breast reduction or in the general population of Danish women.³⁴

- **Effects on Children**

At this time, it is not known if a small amount of silicone may pass through from the breast implant silicone shell into breast milk during breastfeeding. Although there are no current established methods for accurately detecting silicone levels in breast milk,

a study measuring silicon (one component in silicone) levels did not indicate higher levels in breast milk from women with silicone gel-filled implants when compared to women without implants.³⁷

In addition, concerns have been raised regarding potential damaging effects on children born to mothers with implants. Two studies in humans have found that the risk of birth defects overall is not increased in children born after breast implant surgery.^{38,39} Although low birth weight was reported in a third study, other factors (for example, lower pre-pregnancy weight) may explain this finding.⁴⁰ This author recommended further research on infant health.

- **Potential Health Consequences of Gel Bleed**

Small quantities of low molecular weight (LMW) silicone compounds, as well as platinum (in zero oxidation state), have been found to diffuse (“bleed”) through an intact implant shell.^{3,41} The evidence is mixed as to whether there are any clinical consequences associated with gel bleed. For instance, studies on implants implanted for a long duration have suggested that such bleed may be a contributing factor in the development of capsular contracture³ and lymphadenopathy.⁶ However, evidence against gel bleed being a significant contributing factor to capsular contracture and other local complications, is provided by the fact that there are similar or lower complication rates for silicone gel-filled breast implants than for saline-filled breast implants. Saline-filled breast implants do not contain silicone gel and, therefore, gel bleed is not an issue for those products. Furthermore, toxicology testing has indicated that the silicone material used in the Mentor implants does not cause toxic reactions when large amounts are administered to test animals. It also should be noted that studies reported in the literature have demonstrated that the low concentration of platinum contained in breast implants is in the zero oxidation (most biocompatible) state.⁴² In addition, two separate studies sponsored by Mentor have demonstrated that the low concentration of platinum contained in its breast implants is in the zero oxidation (most biocompatible) state.

Mentor performed a laboratory test to analyze the silicones and platinum (used in the manufacturing process), which may bleed out of intact implants into the body. Over 99% of the LMW silicones and platinum stayed in the implant. The overall body of available evidence supports that the extremely low level of gel bleed is of no clinical consequence.

3. What Were the Results of the Mentor Core Study for Augmentation and Revision-Augmentation?

This section of this brochure summarizes the results of the Mentor Core Study conducted on Mentor’s silicone gel-filled breast implants for primary augmentation and revision-augmentation. The Mentor Core Study is the primary clinical study for this product. The results of the Mentor Core Study give you useful information on the experience of other women with Mentor silicone gel-filled implants. While the results cannot be used to predict your individual outcome, they can be used as a rough guide

of what you may expect. Your own complications and benefits depend on many individual factors.

As a note, supplemental safety information was also obtained from the Mentor Adjunct Study, the U.K. Sharpe/Collis Study, and the literature to help assess long-term rupture rate and the consequences of rupture for this product. The literature, which had the most available information on the consequences of rupture, was also used to assess other potential complications associated with silicone gel-filled breast implants. The key literature information is referenced throughout the Breast Implant Complications section above.

3.1. Overview of Mentor Core Study

The Mentor Core Study is a 10-year study to assess safety and effectiveness in augmentation, reconstruction, and revision (revision-augmentation and revision-reconstruction) patients. Patient follow-up is at 6 months, 12 months, 24 months, and annually through 10 years. Safety is assessed by complications, such as implant rupture, capsular contracture, and reoperation. Benefit (effectiveness) is assessed by patient satisfaction and measures of quality of life (QoL).

The Mentor Core Study consists of 1,007 patients, including 551 primary augmentation patients, 146 revision-augmentation patients, 251 primary reconstruction patients, and 59 revision-reconstruction patients. Of these patients, 202 primary augmentation patients, 57 revision-augmentation patients, 134 primary reconstruction patients, and 27 revision-reconstruction patients are in the MRI cohort, which means that they are assessed for silent rupture by MRI at years 1, 2, 4, 6, 8, and 10. The study is currently ongoing, with the results through 3 years reported in this brochure. Mentor will periodically update this brochure as more information becomes available. You should also ask your surgeon if he/she has any available updated clinical information.

Mentor's Core Study results indicate that the risk of at least one occurrence of any complication (including reoperation) at some point through 3 years after implant surgery is 37% for primary augmentation patients and 50% for revision-augmentation patients. The information below provides more details about the complications and benefits you may experience.

Described below are the benefits and complications reported in the Mentor Core Study for augmentation patients. The findings are described separately for primary augmentation and revision-augmentation patients.

3.2. 3-Year Follow-up Rates in Augmentation Patients

At the 3-year follow-up visit, data are reported for 88% of the eligible primary augmentation patients and 87% of the eligible revision-augmentation patients.

3.3. Benefits for Augmentation Patients

The Mentor Core Study measured a variety of outcomes that assessed the benefits of the implants. For augmentation, these outcomes included breast size change, satisfaction, and QoL measures. These outcomes were assessed before implantation and at 1, 2, and 3 years after surgery for those patients who still had their original implants and came back for follow-up visits.

Primary Augmentation Patients: For primary augmentation patients, 370 (67%) out of the original 551 patients were included in the analysis of cup size at 3 years. Of these 370 patients, 359 (97%) experienced at least one cup size increase. The average increase in circumferential chest size was 2.8 inches.

Mentor's satisfaction assessment was based on a single question of "Would the patient have this breast surgery again?" At 3 years, 456 (83%) of the 551 primary augmentation patients enrolled answered that question. Of these 456 patients, 445 (98%) stated to their surgeon that they would have the breast implant surgery again.

With regard to QoL measures at 3 years, an increase in self esteem was noted for 45% of patients after primary breast augmentation on the Rosenberg Self Esteem Scale. There was no change on the overall score of the Body Esteem Scale (BES), but the Sexual Attractiveness Subscale and the Chest Score of the BES increased. The SF-36 is a collection of scales assessing mental and physical health, and there was no change in the SF-36 after primary augmentation. The Tennessee Self Concept Scale (TSCS) is a survey completed by the patient that evaluates how the patient sees herself and what she does, likes, and feels. There was no change in the overall score for the TSCS.

Revision-Augmentation Patients: For revision-augmentation patients, 116 (79%) out of the original 146 patients were included in the analysis of circumferential chest size at 3 years. For these 116 patients, the average increase in circumferential chest size was 2.4 inches.

Mentor's patient satisfaction was based on a single question of "Would the patient have this breast surgery again?" At 3 years, 118 (81%) of the 146 revision-augmentation patients enrolled answered that question. Of these 118 patients, 111 (94%) stated to their surgeon that they would have the breast implant surgery again.

With regard to QoL measures at 3 years, no change in self esteem was noted following revision-augmentation surgery on the Rosenberg Self Esteem Scale or the BES. The SF-36 is a collection of scales assessing mental and physical health, and there were no changes in SF-36. The Tennessee Self Concept Scale (TSCS) is a survey completed by the patient that evaluates how the patient sees herself and what she does, likes, and feels. There was no change in the overall TSCS score.

3.4. 3-Year Complication Rates in Augmentation Patients

The 3-year complication rates are shown from the most common to the least common in Table 1 (primary augmentation) and Table 2 (revision-augmentation) below. The rates reflect the percentage of augmentation patients who experienced the listed complication at least once within the first 3 years after implantation. Some complications occurred more than once for some patients. The two most common complications experienced by primary augmentation patients within the first 3 years of implantation were reoperation (15.4%) and nipple sensation changes (10.4%).

**Table 1 -- 3-Year Complication Rates for MemoryGel Core Study Primary Augmentation Patients
N=551 Patients**

Key Complications	%
Reoperation	15.4
Capsular Contracture Grade III/IV	8.1
Implant Removal with Replacement with Study Device	2.8
Implant Removal without Replacement	2.3
Infection	1.5
Rupture (MRI Cohort) ¹	0.5
Other Complications occurring in ≥ 1% of patients²	
Nipple Complications ³	10.4
Scarring/Hypertrophic Scarring ³	6.7
Breast Mass ³	3.1
Hematoma ³	2.6
Ptosis (sagging) ³	2.3
Breast Sensation Changes ³	2.2
Breast Pain ³	1.7
Miscarriage ⁴	1.5
Trauma ⁵	1.3

1 - There was 1 patient with signs of rupture by MRI of one of her implants through the 3-year timepoint. This has not yet been confirmed with removal and visual inspection of the implant.

2 - The following complications were reported at a rate less than 1%: anaphylaxis, asymmetry, biopsy pending, bruising, deep vein thrombosis, granuloma, implant malposition/displacement, inflammation, lactation difficulties, new diagnosis of rheumatic disease (1 patient with Hashimoto's Thyroiditis, 1 patient with rheumatoid arthritis, and 1 patient with hypothyroidism), necrosis, placement damage (damage to breast implants during insertion, which were then removed while the patient was still on the operating table), position dissatisfaction, positive antinuclear antibodies negative for lupus, rash, suture reaction, seroma, and wrinkling.

3 - Mild occurrences were excluded.

4 - Preoperative miscarriage data were not collected.

5 - Lifted child and stroller; trauma sustained from motor vehicle accident; trauma to breast from fall; and first and second degree frostbite from ice bags placed on breasts the day after surgery to relieve operative pain.

The two most common complications experienced by patients within the first 3 years of revision-augmentation surgery were reoperation (28.0%) and capsular contracture Grade III/IV (18.9%). Notice that the rates for these two complications are higher than for primary augmentation. (For primary augmentation, reoperation was 15.4% and capsular contracture was 8.1%.)

**Table 2 -- 3-Year Complication Rates for MemoryGel Core Study Revision-Augmentation Patients
N=146 Patients**

Key Complications	%
Reoperation	28.0
Capsular Contracture Grade III/IV	18.9
Rupture (MRI Cohort) ¹	7.7
Implant Removal with Replacement with Study Device	6.5
Implant Removal without Replacement	5.9
Infection	1.4
Other Complications occurring in ≥ 1% of patients²	
Nipple Complications ³	10.5
Scarring/Hypertrophic Scarring ³	8.4
Breast Mass ³	6.6
Hematoma ³	2.8
Breast Sensation Changes ³	2.1
Seroma	2.1
Delayed Wound Healing ³	2.1
Wrinkling ³	2.1
Ptosis (sagging) ³	1.5
Breast Pain ³	1.5
Inflammation ³	1.4
Implant Malposition ³	1.4
Extrusion of Intact Implant	1.4

1 - Of the 4 patients who had signs of rupture on MRI, 1 patient had removal of her implants, which showed rupture of both of her implants. This occurred 2 years after she entered the Mentor Core Study as a revision-augmentation patient.

2 - The following complications were reported at a rate less than 1%: back and neck pain related to large implants, ectopic pregnancy, false positive for rupture on mammogram, granuloma, lactation difficulties, miscarriage, muscle spasm, new diagnosis of rheumatic disease (1 patient with rheumatoid arthritis), implant palpability/visibility, and trauma (blunt injury to left breast from being hit by fireworks).

3 - Mild occurrences were excluded.

3.5. Main Reasons for Reoperation in Augmentation Patients

There may be one or more reasons identified for having a reoperation (additional surgery after the primary or revision breast augmentation). Furthermore, there may be multiple surgical procedures (for example, implant removal with or without replacement, capsule procedures, incision and drainage, implant reposition, scar

revision, etc.) performed during a reoperation. In Mentor’s Core Study, there were 176 additional surgical procedures performed in 109 reoperations involving 83 primary augmentation patients.

Table 3 below provides the main reason for each reoperation in primary augmentation patients following initial implantation that were performed through 3 years. The most common reason for reoperation through 3 years in primary augmentation patients was because of capsular contracture (40 of 109 reoperations).

Table 3 -- Main Reasons for Reoperation in Primary Augmentation Patients through 3 Years

Reason for Reoperation	n
Capsular Contracture Grade II, III, IV	40
Patient Request for Style/Size Change	16
Hematoma/Seroma	12
Scarring/Hypertrophic Scarring	12
Biopsy	6
Asymmetry	5
Ptosis (sagging)	4
Infection	3
Delayed Wound Healing	2
Implant Malposition	2
Wrinkling	2
Breast Pain	1
Extrusion of Intact Implant	1
Necrosis	1
Suspected Rupture ¹	1
Tear in Capsule	1
Total	109

1 – The device was removed and found to be intact (not ruptured).

In Mentor’s Core Study, there were 105 additional surgical procedures performed in 58 reoperations involving 39 revision-augmentation patients. Table 4 below provides the main reason for each reoperation in revision-augmentation patients following initial implantation that were performed through 3 years. The most common reason for reoperation in revision-augmentation patients through 3 years was capsular contracture (23 of 58 reoperations).

Table 4 -- Main Reasons for Reoperation in Revision-Augmentation Patients through 3 Years

Reason for Reoperation	n
Capsular Contracture Grade II, III, IV	23
Patient Request for Style/Size Change	7
Biopsy	6
Hematoma/Seroma	5
Delayed Wound Healing	5
Scarring/Hypertrophic Scarring	3
Extrusion of Intact Implant	2
Implant Malposition	2
Asymmetry	1
Ptosis (sagging)	1
Infection	1
Wrinkling	1
Suspected Rupture ¹	1
Total	58

1 – The device was removed and found to be intact (not ruptured).

3.6. Main Reasons for Implant Removal in Augmentation Patients

The main reasons for implant removal among primary augmentation patients in the Mentor Core Study over the 3 years are shown in Table 5 below. There were 45 implants removed in 26 patients. Of these 45 implants, 24 were replaced. The most common reason for implant removal was patient request for style/size change (31 of the 45 implants removed).

Table 5 -- Main Reasons for Implant Removal in Primary Augmentation Patients through 3 Years

Reason for Removal	n
Patient Request for Style/Size Change	31
Capsular Contracture Grade III/IV	5
Breast Pain	2
Infection	2
Necrosis	2
Suspected Rupture ¹	1
Contralateral Explantation	1
Wrinkling	1
Total	45

1 – The device was removed and found to be intact (not ruptured).

The main reasons for implant removal among revision-augmentation patients in the Mentor Core Study over the 3 years are shown in Table 6 below. There were 30 implants removed in 18 patients. Of these 30 implants, 14 were replaced. The most common reason for implant removal was patient request (12 of the 30 implants removed).

Table 6 – Main Reasons for Implant Removal in Revision-Augmentation Patients through 3 Years

Reason for Removal	n
Patient Request for Style/Size Change	12
Capsular Contracture Grade III/IV	10
Patient Dissatisfied with Appearance	2
Asymmetry	1
Extrusion of Intact Implant	1
Scarring/Hypertrophic Scarring	1
Infection	1
Suspected Rupture ¹	1
Abnormal Mammogram	1
Total	30

1 – The device was removed and found to be intact (not ruptured).

3.7. Other Clinical Data Findings in Augmentation Patients

Below is a summary of clinical findings from Mentor’s Core Study with regard to connective tissue disease (CTD), CTD signs and symptoms, cancer, lactation complications, reproduction complications, and suicide. These issues, along with others, are being further evaluated as part of a Mentor large postapproval study involving patients followed through 10 years.

CTD Diagnoses

Three primary augmentation patients and one revision-augmentation patient in the Mentor Core Study were reported to have a new diagnosis of CTD according to a rheumatologist. These diagnoses were Hashimoto’s Thyroiditis at 2 years, two cases of rheumatoid arthritis at 2 and 3 years, and hypothyroidism at 2 years. It cannot be concluded that these CTD diagnoses were caused by the implants because there was no comparison group of similar women without implants.

CTD Signs and Symptoms

Data on over 100 self-reported signs and symptoms, including about 50 self-reported rheumatological symptoms, were collected. Compared to before having the implants, significant increases were found for fatigue, exhaustion, joint swelling, joint pain, numbness of hands, frequent muscle cramps, and the combined categories of fatigue, pain, and fibromyalgia -like symptoms in primary augmentation patients, and for joint pain in revision-augmentation patients. These increases were not found to be related to simply getting older over time. The Mentor Core Study was not designed to evaluate cause and effect associations because there is no comparison group of women without implants, and because other contributing factors, such as medications and lifestyle/exercise, were not studied. Therefore, it cannot be determined whether these increases were due to the implants or not, based on the Mentor Core Study. However, you should be aware that you may experience an increase in these symptoms after receiving breast implants.

Cancer

There were no primary augmentation patients with new diagnoses of breast cancer through 3 years in Mentor's Core Study. As previous breast cancer was an exclusion criteria for primary augmentation patients, there were no reports of breast cancer reoccurrence in this indication. There were no reports of new diagnoses or reoccurrence of breast cancer in revision-augmentation patients. There were no reports of other cancers, such as brain, respiratory, or cervical/vulvar.

Lactation Complications

Two (8%) of the 25 primary augmentation patients who attempted to breast feed following breast implantation in Mentor's Core Study through 3 years experienced difficulty with breast feeding. Of the 7 revision-augmentation patients who attempted to breast feed after receiving breast implants, 1 (14%) had difficulty breast feeding.

Reproduction Complications

Eight (1.5%) of the primary augmentation patients in Mentor's Core Study reported a miscarriage through 3 years. There were no reports of miscarriage in revision-augmentation patients.

Suicide

There were no reports of suicide in either the primary augmentation or revision-augmentation indications in Mentor's Core Study through 3 years.

4. What are the Surgery Considerations for Receiving Breast Implants?

This section provides a discussion of surgical considerations for breast augmentation .

4.1. Surgical Considerations for Breast Augmentation

4.1.1. Alternatives to Breast Augmentation with Silicone Breast Implants

For primary augmentation patients, alternatives may include:

- Accept your breasts as they are and have no surgery.
- Wear a padded bra or external prostheses.
- Have mastopexy surgery (breast lift) without an implant.
- Have surgery with saline implants.

For revision-augmentation patients, alternatives may include:

- No further revision
- Removal with or without replacement.

4.1.2. Choosing a Surgeon

When choosing a surgeon who is experienced with breast augmentation , you should know the answers to the following types of questions:

- How many breast augmentation implantation procedures does he/she perform per year?
- How many years has he/she performed breast augmentation procedures?

- Has he/she obtained training certification from Mentor to use its silicone gel-filled breast implants?
- Is he/she board certified, and if so, with which board?
- In which state(s) is he/she licensed to practice surgery? (Note that some states provide information on disciplinary action and malpractice claims/settlements to prospective patients, either by request or on the Internet.)
- What is the most common complication he/she encounters with breast augmentation?
- What is his/her reoperation rate with breast augmentation, and what is the most common type of reoperation he/she performs?
- Can he/she perform this surgery in a hospital, as well as in the surgeon's independent surgery center? (Note that hospitals require the demonstration of evidence of appropriate training in specific procedures before allowing surgeons to operate in their facilities.)

4.1.3. Implant Shape and Size

Depending on the desired shape you wish to achieve, you and your surgeon have implants with three different round profiles, or styles, from which to choose. Generally, the larger you want your cup size, the larger the breast implant the surgeon will consider (measured in cubic centimeters, or cc's), not in cup sizes, because this depends on the size and shape of the individual woman's chest.

Your surgeon will also evaluate your existing breast and skin tissue to determine if you have enough to cover the breast implant you are considering, or, in some cases, such as after pregnancy, too much extra skin. If you desire a breast implant size that is too large for your tissue, the surgeon may warn you that breast implant edges may be visible or palpable postoperatively. Also, excessively large breast implants may speed up the effects of gravity on the breast, and can result in droop or sag at an earlier age. A recent report indicates that larger sized implants (greater than 350cc) may be too large for many women, increasing the risk of developing complications such as implant extrusion, hematoma, infection, palpable implant folds, and visible skin wrinkling requiring surgical intervention to correct these complications.⁴³


4.1.4. Surface Texturing

Some studies suggest that surface texturing reduces the chance of severe capsular contracture,⁴⁴ while other studies do not.^{43,45} Mentor's Core Study did not show a difference in the likelihood of developing capsular contracture with textured implants compared to smooth-surfaced implants.

A textured implant may require a larger incision because the rougher textured surface makes it harder to place into the pocket without undue stress, which might damage the implant or decrease its durability.

4.1.5. Implant Placement

The breast implant can be placed either partially under the pectoralis major muscle (submuscular) or on top of



the muscle and under the breast glands (subglandular). You should discuss with your surgeon the advantages and disadvantages of the implant placement selected for you, as described in Table 7 below.

Table 7 – Comparison between Submuscular versus Subglandular Placement


Submuscular Placement	Subglandular Placement
Surgery may be longer	Surgery may be shorter
Recovery may be longer	Recovery may be shorter
May be more painful	May be less painful
Reoperation may be more difficult	May provide easier access for reoperation
Less visible and palpable implants	More visible and palpable implants
Less likelihood of capsular contracture ³	Greater likelihood of capsular contracture ^{43,45}
Easier imaging during mammography exam	More difficult imaging during mammography exam
May be preferable if you have thin or weakened breast tissue	May not be recommended if you have thin or weakened breast tissue.

4.1.6. Incision Sites

You should discuss with your surgeon the pros and cons for the incision site specifically recommended for you.

The incision size will be larger than for a saline breast augmentation. There are 3 common incision sites: under the arm (axillary), around the nipple (periareolar), or within the breast fold (inframammary).

- **Periareolar** - This incision is typically more concealed, but since it also involves cutting through the breast tissue, it is associated with a higher likelihood of breast feeding difficulties, as compared to the other incision sites.⁴⁶ Cutting through the tissue may increase the chance that there will be a change in breast and/or nipple sensation.
- **Inframammary** - This incision is generally less concealed than periareolar and associated with less breast feeding difficulties than the periareolar incision site. It is also the most commonly used incision site at the present time, and is felt to give the best access to and control of the breast implant pocket.
- **Axillary** - This incision is less concealed than periareolar and associated with less breast feeding difficulties than the periareolar incision site. If the incision is made under the arm, the surgeon may use a probe fitted with a miniature camera, along with minimally invasive (very small) instruments, to create a “pocket” for the breast implant. This approach is more difficult, and may increase the risk of damage to, and unexpected location of, the implant.



- **Umbilical** (belly button) - This incision site has not been studied in Mentor's Core Study and should not be used for a wide variety of reasons, including potential damage to the implant shell.

4.1.7. Additional Procedures at the Time of Breast Augmentation

Your surgeon will examine your breasts and help you make decisions to obtain the best result in your individual situation. In some cases, particularly after pregnancy or significant weight loss, implants alone may not address all of the issues, such as sagging or extra skin, affecting your breasts. This is particularly true when there is extra skin remaining from when the breasts were engorged with milk, or when you might have been carrying more weight.

In these situations, your surgeon may recommend a breast lift (mastopexy) to remove some of the extra skin, or to lift the breasts, at the time of implant placement. Mastopexy involves removing a strip of skin from under the breast or around the nipple to lift the nipple and breast location, and tighten the skin over the breast. Your surgeon will discuss the potential risks, and the location of the additional scars which might be required to lift your breasts or to remove the extra skin.

4.1.8. Palpability

Implants may be more palpable or noticeable if there is an insufficient amount of skin/tissue available to cover the implant and/or when the implant is placed subglandularly.

4.1.9. Surgical Setting and Anesthesia

Augmentation surgery is usually performed on an outpatient basis, in a specialized operating room which may be located in a hospital, a surgery center, or surgical suite in the surgeon's office. General anesthesia is commonly used, and local anesthesia with sedation is also an option. You should be sure to check with your surgeon and with the facility where the surgery will take place to become aware of the tests, presurgical examinations, and length of time you need to be without food or your routine medications prior to the surgical procedure.

4.1.10. Postoperative Care

You will probably feel somewhat tired and sore for several days following the operation, and your breasts may remain swollen and sensitive to physical contact for a month or longer. You may also experience a feeling of tightness in the breast area as your skin adjusts to your new breast size. The feeling in the breasts and nipple area also may be diminished during this time of swelling and immediate post surgery recovery. Other possible complications have been described above.

Postoperative care depends on each patient's situation, may involve the use of a special postoperative bra, compression bandage, or jog bra for extra support and positioning while you heal. Some surgeons may not want you to wear a bra at all for a period of time following the surgery. At your surgeon's recommendation, you will most likely be able to return to work within a few days, although for at least a couple of

weeks you should avoid any strenuous activities that could raise your pulse and blood pressure, or require strenuous use of your arms and chest. Your surgeon may also recommend breast massage exercises. Please be aware that smoking may interfere with the healing process and prolong your recovery.

Note: If you experience fever, do not feel well, or see noticeable swelling and/or redness or drainage in your implanted breast(s), you should contact your surgeon immediately.

4.2. Other Factors to Consider in Revision-Augmentation Surgery

Some revision surgeries require removal of an intact implant (for example, capsulotomy and pocket adjustments), while others do not require removal of the implant. Any device that has been removed during revision surgery should not be reimplanted. Mentor breast implants are “for single use only.”

5. What Sort of Follow-Up Examinations Should I Have?

5.1. Breast Self-Examinations

You should perform a breast self-examination monthly. This may be more difficult with an implant in place. In order to do this effectively, you should ask your surgeon to help you tell the difference between the implant and your breast tissue. Care should be taken not to squeeze the implant excessively. Any new lumps may be evaluated with a biopsy, as appropriate. If a biopsy is performed, care must be taken to avoid injuring the implant.

5.2. Screening for Silent Rupture

Because most ruptures of silicone breast implants are silent, in most cases, neither you nor your surgeon will be able to find evidence of rupture. Therefore, evaluation of your implants is recommended to screen for implant rupture. The best method of screening is currently MRI at a center with a breast coil, with a magnet of at least 1.5 Tesla. The MRI should be read by a radiologist who is familiar with looking for implant rupture.

It is recommended that your first MRI evaluation take place starting at 3 years after implant surgery and every 2 years, thereafter, even if you are experiencing no problems with your implant. If signs of rupture are seen on MRI, then you should have your implant removed, with or without replacement. More information on rupture is provided in Section 2 of this brochure. Your doctor should assist you in locating a radiology/screening center, as well as a radiologist who is familiar with the technique and equipment for proper MRI screening for silent rupture of your breast implant.

5.2.1. Symptomatic Rupture

Symptoms associated with rupture may include hard knots or lumps surrounding the implant or in the armpit, loss of size of the breast or implant, pain, tingling, swelling, numbness, burning, or hardening of the breast. If you notice any of these changes,

see your plastic surgeon so that he or she can examine the implants for rupture and determine whether you need to have an MRI examination to find out if your symptoms are due to rupture of the implant. If rupture has occurred, you should have your implant removed. More information on rupture is provided in Section 2 of this brochure.

You should monitor your breast implants for signs of symptomatic rupture when you check your breasts for lumps monthly. Examine your breast tissue by feeling for lumps. Then feel the breast implants. Move the implants around while looking in the mirror. Look for changes in shape, size, and feel of the implants. Know, and pay attention to, how the breast implants feel.

5.3. Mammography

The current recommendations for getting screening/preoperative mammograms are no different for women with breast implants than for those without implants. Mammography exams should be interpreted by radiologists experienced in the evaluation of women with breast implants. It is essential that you tell your mammography technologist that you have an implant before the procedure. You should request a diagnostic mammogram, rather than a screening mammogram, because more pictures are taken with diagnostic mammography. The technologist can use special techniques to reduce the possibility of rupture and to get the best possible views of the breast tissue. More information on mammography is provided in Section 1.4.


6. What are the Types of MemoryGel Silicone Gel-Filled Breast Implants Available from Mentor?

Mentor's silicone gel-filled breast implants, referred to as MemoryGel products, come in a variety of profiles and sizes. All currently available MemoryGel breast implants have either a textured shell or smooth surface shell (Table 8).

Table 8 -- Approved MemoryGel Implant Styles

Breast Implant Description	Size Range
Smooth, Round, Moderate Profile	100-800 cc
Textured Round, Moderate Profile	100-800 cc
Smooth, Round, Moderate Plus Profile	100-800 cc
Textured, Round, Moderate Plus Profile	100-800 cc
Smooth, Round, High Profile	125-800 cc
Textured, Round, High Profile	125-800 cc

The following diagrams illustrate the high, moderate plus, and moderate profiles.



Your surgeon will discuss the various styles and options of silicone breast implants that will help you achieve your desired result.

7. How Should I Report Problems with My Implant?

The Food and Drug Administration (FDA) requires that serious injuries (defined as those that need medical or surgical intervention to prevent permanent damage) be reported by hospitals if they are aware of the serious injuries. If you believe that you have experienced one or more serious problems related to your breast implants, you are encouraged to report the serious problem(s) through your health professional to the FDA. Although reporting by doctors or other health professionals is preferred, women may also report any serious problem directly through FDA's MedWatch voluntary reporting system. You can report by telephone to 1-800-FDA-1088; by FAX, use Form 3500 to 1-800-FDA-0178; electronically at <http://www.fda.gov/medwatch/index.html>; or by mail to MedWatch Food and Drug Administration, HF-2, 5600 Fishers Lane Rockville, MD 20857-9787. **Keep a copy of the MedWatch form completed by your doctor for your records.** The information reported to MedWatch is entered into databases to be used to follow safety trends (patterns) of a device and to determine whether further follow-up of any potential safety issues related to the device is needed.

8. What Is Device Tracking?

Silicone gel-filled breast implants are subject to Device Tracking by Federal regulation. This means that your physician will be required to report to Mentor the serial number of the device(s) you receive, the date of surgery, and information relating to the physician's practice. This information will be recorded on the Device Tracking Form supplied by Mentor with each silicone gel-filled breast implant.

Mentor strongly recommends that all patients receiving silicone gel-filled breast implants participate in Mentor's device tracking program. This will help ensure that Mentor has a record of each patient's contact information so that all patients, including you, can be contacted in the case of a recall or other problems with your implants that you should be made aware of. Please inform Mentor whenever your contact information changes.

9. **What Are Mentor's Product Replacement Policy and Limited Warranties?**

The following is a description of the assistance available from Mentor Lifetime Product Replacement Policy and the Mentor Advantage and Enhanced Advantage Limited Warranties.

The **Mentor Lifetime Product Replacement Policy** involves the free lifetime product replacement for its gel-filled and saline-filled breast implants, worldwide. When implant replacement is required and the Mentor Product Replacement Policy applies (see below), Mentor will provide the same or similar Mentor breast implant at no cost for the lifetime of the product. If a more expensive product is requested, Mentor will invoice the surgeon for the price difference.

The **Mentor Standard Advantage Limited Warranty** is free of charge to all patients who are implanted with Mentor gel-filled or saline-filled breast implants in the United States and Puerto Rico. When the limited warranty applies, Mentor provides the following:

- Financial assistance: For the first ten years following a breast implant procedure, Mentor will provide financial assistance up to \$1200 to help cover operating room, anesthesia, and surgical charges not covered by insurance. Financial assistance applies to covered events only (see below). Operating room and anesthesia charges will be given payment priority. In order to qualify for financial assistance, you will need to sign a Release Form.
- Free contralateral (opposite side) implant replacement upon surgeon request.
- Non-cancelable terms.

The **Mentor Enhanced Advantage Limited Warranty** is an optional limited warranty available for women who are implanted with Mentor gel-filled or saline-filled breast implants in the United States and Puerto Rico. To be eligible, the Mentor Enhanced Advantage Limited Warranty must be purchased for an enrollment fee of \$100 within 45 days from implantation. When the warranty applies, Mentor provides the following:

- Financial assistance: For the first ten years following a breast implant procedure, Mentor will provide financial assistance up to \$2400 to help cover operating room, anesthesia, and surgical charges not covered by insurance. Financial assistance applies to covered events only (see below). Operating room and anesthesia charges will be given payment priority. In order to qualify for financial assistance, you will need to sign a Release Form.
- Free contralateral implant replacement upon surgeon request.
- Non-cancelable terms.

With both the Mentor Standard Advantage and Mentor Enhanced Advantage Limited Warranties, it is important for you to also maintain your own records to ensure validation of your enrollment, as it is possible your surgeon may not retain your records for the entire duration of the limited warranty.

Products Covered

The Mentor Standard Advantage Limited Warranty coverage applies to all Mentor gel-filled and saline-filled breast implants that are implanted in the United States and Puerto Rico, provided they have been:

- Implanted in accordance with the Mentor package insert, current to the date of implantation, and other notifications or instructions published by Mentor; and
- Used by appropriately qualified, licensed surgeons, in accordance with accepted surgical procedures.

Events Covered

The Mentor Lifetime Product Replacement Policy, and the Standard Mentor Advantage and Enhanced Advantage Limited Warranties coverages apply to the following:

- Rupture due to localized stress, folding, manufacturing defect, patient trauma, or unknown cause.
- Other loss-of-shell integrity events, such as surgical damage may also be covered by these programs. Mentor reserves the right to determine if specific, additional events should be covered.

Events Not Covered

The Mentor Lifetime Product Replacement Policy and the Mentor Standard Advantage and Enhanced Advantage Limited Warranties coverages do not apply to the following:

- Removal of intact implants due to capsular contracture, or wrinkling.
- Loss of implant shell integrity resulting from reoperative procedures, open capsulotomy), or closed compression capsulotomy procedures.
- Removal of intact implants for size alteration.

Filing for Financial Assistance

- To file a Mentor Advantage claim for product replacement and/or financial assistance, the surgeon must contact the Mentor Product Evaluation Department at 1-866-250-5115 prompt #1 prior to replacement surgery.
- For financial assistance claims, a patient-specific Release form will be generated that you must sign and return.
- For either replacement or financial assistance claims, the surgeon must send the explanted, decontaminated Mentor breast implant(s) within six months of the date of explantation (implant removal) to:

Mentor Product Evaluation
3041 Skyway Circle North
Irving, Texas 75038-3540

- Upon receipt, review and approval of the completed claim, including receipt of the explanted product and your completion of a full general release, financial assistance will be issued.

This is a summary of the coverage of the Mentor Advantage and Enhanced Advantage Limited Warranties. It is an overview only and not a complete statement of the program. A copy of the complete Mentor Advantage and Enhanced Advantage Limited Warranties for saline-filled and silicone gel-filled breast implants may be obtained by writing or calling:

Consumer Affairs Department
Mentor Corporation
201 Mentor Drive
Santa Barbara, CA 93111
1-800-525-0245

A copy of the complete programs may also be obtained from your surgeon or by going to www.mentorcorp.com.

THESE ARE LIMITED WARRANTIES ONLY AND ARE SUBJECT TO THE TERMS AND CONDITIONS SET FORTH AND EXPLAINED IN THE APPLICABLE MENTOR LIMITED WARRANTIES. ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, BY OPERATION OF LAW OR OTHERWISE, INCLUDING BUT NOT LIMITED TO, IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS ARE EXCLUDED.

Mentor reserves the right to cancel, change, or modify the terms of the Mentor Advantage and Enhanced Advantage coverages. Any such cancellation, change, or modification will not affect the currently stated terms of the Mentor Advantage and Enhanced Advantage coverages for those already enrolled.

10. What are Other Sources of Additional Information?

Upon request, you will be provided with a copy of the package insert (Directions for Use). You can request a copy from your surgeon or from Mentor. The package insert has many undefined medical and technical terms because it contains information directed only to the surgeon.

For more detailed information on the preclinical and clinical studies conducted by Mentor, you are referred to the Summary of Safety and Effectiveness Data (SSED) for this product at <http://www.fda.gov/cdrh/breastimplants/>.

If you should decide to get breast implants, you will be given a device identification card with the style and serial number of your breast implant(s). This will be given to you right after your surgery. It is important that you keep a copy of this card because you may need to refer to that information at a later date.

For additional information or questions about Mentor breast implants, please call 1-800-MENTOR8.

Mentor Corporation

1-800-MENTOR8

www.mentorcorp.com

www.loveyourlook.com

Institute of Medicine Report on the Safety of Silicone Implants

www.nap.edu/catalog/9618.html

Food and Drug Administration

1-888-INFO-FDA or 240-276-3103

<http://www.fda.gov/cdrh/breastimplants/>

You can find important information in the FDA breast implant consumer handbook, which is available through the phone number or website provided above.

American Society of Plastic Surgeons

http://www.plasticsurgery.org/public_education/Silicone-Breast-Implant-Surgery.cfm



201 Mentor Drive
Santa Barbara
CA 93111 USA
(800) MENTOR-8
www.mentorcorp.com

ACKNOWLEDGMENT OF INFORMED DECISION

I understand that this patient brochure, “Important Information for Augmentation Patients About Mentor MemoryGel® Silicone Gel-filled Breast Implants,” is intended to provide the information regarding the risks and benefits of silicone gel-filled breast implants, both general and specific to Mentor’s MemoryGel products. I understand that silicone breast implant surgery involves risks and benefits, as described in this brochure. I also understand that the long-term (i.e., 10-year) safety and effectiveness of silicone gel-filled breast implants continue to be studied. I understand that reading and fully understanding this brochure is required, but that there also must be consultation with my surgeon.

By circling the correct response and signing below, I acknowledge:

- Y/N** I have had adequate time to read and fully understand this brochure;
- Y/N** I have had an opportunity to ask my surgeon any questions I may have about this brochure or any other issues related to breast implants or breast implant surgery;
- Y/N** I have considered the alternatives to silicone breast implants and have decided to proceed with silicone breast implant surgery;
- Y/N** I have had an adequate amount of time to review and understand this information, before scheduling my silicone breast implant surgery; and
- Y/N** I will retain this brochure, and I am aware that I may also ask my surgeon for a copy of this signed acknowledgment.

PATIENT (PRINT NAME)

SIGNATURE OF PATIENT*

DATED:

* A patient must be at least 22 years old for primary and revision breast augmentation with silicone breast implants.

By my signature below, I acknowledge that:

- My patient has had an adequate amount of time before making her final decision;
- All questions outlined above have been answered “Yes” by my patient;
- My patient has had an adequate amount of time before making her final decision; and
- Documentation of this Informed Decision will be retained in my patient’s permanent record.

SIGNATURE OF SURGEON

DATED:

Index

A

Areola, 4
Asymmetry, 4, 21, 22, 30, 32, 34, 35
Autoimmune disease, 4, 6, 9, 14, 25
Axillary, 4, 39

B

BES. *See* **Body Esteem Scale**
Biocompatible, 4, 27
Biopsy, 4, 21, 23, 25, 30, 32, 34, 41
Body Dysmorphic Disorder, 14
Body Esteem Scale, 4, 29
Breast augmentation, 4, 6, 9, 12, 13, 15, 16, 23, 27, 37, 39, 40
 primary augmentation, 4, 9, 12, 13, 14, 15, 17, 18, 19, 21, 27, 28, 29, 30, 31, 32, 34, 35, 36, 37
 revision-augmentation, 4, 6, 9, 12, 14, 15, 18, 19, 20, 21, 27, 28, 29, 30, 31, 32, 33, 34, 35, 36, 37, 41
Breast cancer, 8, 15, 25, 36
Breast mass, 4, 30, 31
Breast reconstruction, 4, 13

C

Calcification, 5, 23
Capsular contracture, 5, 16, 19, 20, 21, 22, 23, 26, 27, 30, 31, 32, 33, 34, 35, 38, 39, 45
Capsule, 5, 6, 7, 17, 18, 19, 20, 23, 32, 33
Capsulectomy, 5
Capsulorrhaphy, 5
Capsulotomy, 41, 45
Capsulotomy (closed), 5, 16
Capsulotomy (open), 5, 45
Congenital anomaly, 6
Connective tissue disease/disorder, 6, 16, 24, 25, 35
Contraindication, 6, 13
Contralateral, 6, 34, 44
Core Study, 6, 15, 17, 18, 20, 21, 27, 28, 30, 31, 32, 33, 34, 35, 36, 38, 39
CTD. *See* **Connective tissue disease/disorder**

D

Delayed wound healing, 6, 24, 31, 32, 34
Depression, 14
Displacement, 6, 8, 22, 30

E

Extrusion, 6, 23, 24, 31, 33, 34, 35, 38

F

Fatigue, 6, 10, 19, 24, 25, 36
Fibromyalgia, 6, 9, 24, 36
Fibrous tissues, 7

G

Granuloma, 7, 19, 24, 30, 32

H

Hematoma, 7, 19, 21, 22, 23, 30, 31, 32, 34, 38
Hypertrophic scarring, 7, 21, 22, 30, 31, 32, 34, 35

I

Immune response, 7
Infection, 7, 13, 16, 19, 21, 22, 23, 24, 30, 31, 32, 34, 35, 38
Inflammation, 7, 9, 30, 31
Inframammary fold, 7
Inframammary incision, 7
Inpatient surgery, 7

L

Lactation, 7, 16, 30, 32, 35, 36
Low Molecular Weight Silicones, 8
Lupus, 4, 6, 14, 24, 31
Lymphadenopathy, 8, 19, 24, 26

M

Malposition, 8, 30, 31, 33, 34
Mammogram, 8, 15, 23, 32, 35, 42
Mammography, 8, 15, 16, 18, 25, 39, 42
Mammoplasty, 8
Mastopexy, 8, 12, 37, 40
Mental Health Disorders, 14
Metastatic Disease, 9
Migration, 9, 17, 18
MRI, 8, 10, 14, 15, 16, 17, 18, 28, 30, 31, 41, 42

N

Necrosis, 9, 23, 24, 30, 33, 34

O

Outpatient surgery, 9, 40

P

Pain, 5, 6, 7, 9, 10, 11, 15, 16, 17, 19, 20, 22, 23, 24, 25,
30, 31, 32, 33, 34, 36, 39, 41
Palpability, 9, 19, 22, 24, 32, 40
Palpate, 9
Pectoralis, 9, 12, 38
Periareolar, 9, 23, 39
Plastic surgery, 8, 9
Postoperatively, 9, 38
Ptoisis, 9, 30, 31, 32, 34

R

Reoperation, 5, 9, 19, 20, 21, 28, 30, 31, 32, 33, 34, 37,
39
Rheumatological Disease/Disorder, 9
Rosenberg Self Esteem Scale, 10, 29
Rupture, 5, 8, 10, 11, 14, 15, 16, 17, 18, 19, 20, 21, 23,
27, 30, 31, 32, 33, 34, 35, 41, 42, 45
 extracapsular, 6, 10, 17, 18
 intracapsular, 7, 10, 17, 18
 primary augmentation, 17
 revision-augmentation, 17
 silent, 10, 11, 14, 16, 17, 18, 25, 28, 41

symptomatic, 10, 11, 15, 41, 42

S

Saline, 10, 27, 37, 39, 44, 45, 46
Satisfaction, 11, 22, 28, 29
Scar revision, 10, 32
Scleroderma, 4, 6, 14, 24
Screening, 8, 14, 15, 41, 42
Seroma, 10, 19, 21, 22, 23, 31, 32, 34
SF-36 Scale, 10, 29
Silent rupture, 14
Silicone elastomer, 10
Subglandular placement, 10, 38, 40
Submuscular placement, 11, 38
Surgical incision, 5, 11

T

Tennessee Self Concept Scale, 11, 29
Trauma, 4, 13, 16, 18, 30

W

Wrinkling, 14, 16, 21, 22, 31, 33, 34, 35, 38, 45

References

- ¹ Hölmich, L.R., et al. 2005a. The diagnosis of silicone breast implant rupture. Clinical findings compared to findings at MRI. *Ann Plast Surg* 54 (6): 583-9.
- ² Hölmich, L.R., et al. 2005b. The diagnosis of breast implant rupture: MRI findings compared to findings at explantation. 2005. *Eur J. Radiol.* 53: 213-25.
- ³ Bondurant, S., V.L. Ernster and R. Herdman, Eds. 2000. Safety of silicone breast implants. Committee on the Safety of Silicone Breast Implants, Division of Health Promotion and Disease Prevention, Institute of Medicine. Washington, D.C.: National Academy Press.
- ⁴ Hölmich, L.R., et al. 2004. Untreated silicone breast implant rupture. *Plast. Reconstr. Surg.* 114:204-214.
- ⁵ Holmich, L.R., et al. 2001. Prevalence of silicone breast implant rupture among Danish women. *Plast. Reconstr. Surg.* 108(4):848-858.
- ⁶ Katzin, W.E., et al. 2005. Pathology of lymph nodes from patients with breast implants: a histologic and spectroscopic evaluation. *Am J Surg Pathol.*29(4):506-11.
- ⁷ Berner, I., M., et al. 2002. Comparative examination of complaints of patients with breast-cancer with and without silicone implants. *Eur. J Obstet. Gynecol. Reprod. Biol.* 102:61-66.
- ⁸ Brown, S.L., et al. 2001. Silicone gel breast implant rupture, extracapsular silicone, and health status in a population of women. *J. Rheumatol.* 28:996-1003.
- ⁹ Hölmich, L.R., et al. 2003. Self-reported diseases and symptoms by rupture status among unselected Danish women with cosmetic silicone breast implants. *Plast. Reconstr. Surg.* 111:723-732.
- ¹⁰ Wolfe, F. and J. Anderson. 1999. Silicone filled breast implants and the risk of fibromyalgia and rheumatoid arthritis. *J. Rheumatol.* 26:2025-2028.
- ¹¹ Brinton, L.A., et al. 2004. Risk of connective tissue disorders among breast implant patients. *Am. J. Epidemiol.* 160(7):619-27.
- ¹² Janowsky, E.C., et al. 2000. Meta-Analyses of the Relation Between Silicone Breast Implants and the Risk of Connective-Tissue Diseases. *N. Engl. J. Med.* 342(11):781-90.
- ¹³ Lipworth, L.R.E., et al. 2004. Silicone breast implants and connective tissue disease: An updated review of the epidemiologic evidence. *Ann. Plast. Surg.* 52:598-601.
- ¹⁴ Tugwell, P., et al. 2001. Do silicone breast implants cause rheumatologic disorders? A systematic review for a court-appointed national science panel. *Arthritis Rheum.* (11):2477-84.
- ¹⁵ Weisman, M.H., et al. 1988. Connective-tissue disease following breast augmentation: A preliminary test of the human adjuvant tissue hypothesis. *Plast. Reconstr. Surg.* 82(4):626-30.
- ¹⁶ Williams, H.J., et al. 1997. Breast implants in patients with differentiated and undifferentiated connective tissue disease. *Arthritis and Rheumatism* 40(3):437-40.
- ¹⁷ Breiting, V.B., et al. 2004. Long-term health status of Danish women with silicone breast implants. *Plast. Reconstr. Surg.* 114:217-26.
- ¹⁸ Fryzek, J.P., et al. 2001. Self-reported symptoms among women after cosmetic breast implant and breast reduction surgery. *Plast. Reconstr. Surg.* 107:206-13.
- ¹⁹ Kjølner, K., et al. 2004. Self-reported musculoskeletal symptoms among Danish women with cosmetic breast implants. *Ann Plast Surg* 52(1):1-7.
- ²⁰ Brinton, L.A., et al. 2000. Breast cancer following augmentation mammoplasty (United States). *Cancer Causes Control.* 11(9):819-27. *J. Long Term Eff. Med. Implants.* 12(4):271-9.

-
- ²¹ Bryant, H., and Brasher, P. 1995. Breast implants and breast cancer--reanalysis of a linkage study. *N. Engl. J. Med.* 332(23):1535-9.
- ²² Deapen, D.M., et al. 1997. Are breast implants anticarcinogenic? A 14-year follow-up of the Los Angeles Study. *Plast. Reconstr. Surg.* 1997 99(5):1346-53.
- ²³ Herdman, R.C., et al. 2001. Silicone breast implants and cancer. *Cancer Invest.* 2001;19(8):821-32.
- ²⁴ Pukkala, E., et al. 2002. Incidence of breast and other cancers among Finnish women with cosmetic breast implants, 1970-1999. *J. Long Term Eff. Med. Implants* 12(4):271-9.
- ²⁵ Deapen, D., et al. 2000. Breast cancer stage at diagnosis and survival among patients with prior breast implants. *Plast Reconstr Surg.* 105(2):535-40.
- ²⁶ Jakubietz, M.G., et al. 2004. Breast augmentation: Cancer concerns and mammography – A literature review. *Plast. Reconstr. Surg.* 113:117e-122e.
- ²⁷ Miglioretti, D.L., et al. 2004. Effect of breast augmentation on the accuracy of mammography and cancer characteristics. *JAMA* 291(4):442-50.
- ²⁸ Brinton, L.A., et al. 2001b. Cancer risk at sites other than the breast following augmentation mammoplasty. *Ann. Epidemiol.* 11:248-56.
- ²⁹ McLaughlin, J.K. and L. Lipworth. 2004. Brain cancer and cosmetic breast implants: A review of the epidemiological evidence. *Ann. Plast. Surg.* 52(2):15-17.
- ³⁰ Cook, L.S. 1997. Characteristics of women with and without breast augmentation. *J. Amer. Med. Assoc.* 20:1612-7.
- ³¹ Fryzek, J.P., et al. 2000. Characteristics of women with cosmetic breast augmentation surgery compared with breast reduction surgery patients and women in the general population of Sweden. *Ann Plast Surg.* 45(4):349-56.
- ³² Kjølner K., et al. 2003. Characteristics of women with cosmetic breast implants compared with women with other types of cosmetic surgery and population-based controls in Denmark. *Ann Plast Surg.* 50(1):6-12.
- ³³ Brinton, L.A., et al. 2001a. Mortality among augmentation mammoplasty patients. *Epidemiol.* 12(3):321-6.
- ³⁴ Jacobsen, P.H., et al. 2004. Mortality and suicide among Danish women with cosmetic breast implants. *Arch. Int. med.* 164(22):2450-5.
- ³⁵ Koot, V., et al. 2003. Total and cost specific mortality among Swedish women with cosmetic breast implants: prospective study. *Br. J. Med.* 326(7388):527-528.
- ³⁶ Pukkala, E., et al. 2003. Causes of death among Finnish women with cosmetic breast implants, 1971-2001. *Ann. Plast. Surg.* 51(4):339-42.
- ³⁷ Lugowski, S.J., et al. 2000. Analysis of silicon in human tissues with special reference to silicone breast implants. *J. Trace Elem. Med. Biol.* 14(1):31-42.
- ³⁸ Kjølner, K., et al. 2002. Health outcomes in offspring of Danish mothers with cosmetic breast implants. *Ann. Plast. Surg.* 48:238-245.
- ³⁹ Signorello, L.B., et al. 2001. Offspring health risk after cosmetic breast implantation in Sweden. *Ann. Plast. Surg.* 46:279-286.
- ⁴⁰ Hemminki, E., et al. 2004. Births and perinatal health of infants among women who have had silicone breast implantation in Finland, 1967-2000. *Acta Obstet Gynecol Scand.* 83(12):1135-40.
- ⁴¹ Flassbeck, D.B., et al. 2003. Determination of siloxanes, silicon, and platinum in tissues of women with silicone gel-filled implants. 375(3):356-62 (for example, data from Patients B & C).
- ⁴² Stein, J., et al. 1999. In situ determination of the active catalyst in hydrosilylation reactions using highly reactive Pt(0) catalyst precursors. *J. Am. Chem. Soc.* 121(15):3693-3703. Chandra, G., et al. 1987. A convenient and novel route to bis(alkyne)platinum(0) and other platinum(0) complexes from Speier's hydrosilylation catalyst.

Organometallics. 6:191-2. Lappert, M.F. and Scott, F.P.A. 1995. The reaction pathway from Speier's to Karstedt's hydrosilylation catalyst. *J. Organomet. Chem.* 492(2):C11-C13. Lewis, L.N., et al. 1995. Mechanism of formation of platinum(0) complexes containing silicon-vinyl ligands. *Organometallics.* 14:2202-13.

⁴³ Henriksen, T.F., et al. 2005. Surgical intervention and capsular contracture after breast augmentation: a prospective study of risk factors. *Ann. Plast. Surg.* 54(4):343-51.

⁴⁴ For example: Seify, H., et al. 2005. Preliminary (3 years) experience with smooth wall silicone gel implants for primary breast augmentation. 54(3):231-5.

⁴⁵ For example: Kulmala, I., et al. 2004. Local complications after cosmetic breast implant surgery in Finland. *Ann. Plast. Surg.* 53(5):413-9.

⁴⁶ Hurst, N.M. 1996. Lactation after augmentation mammoplasty. *Obstet. Gynecol.* 87(1):30-4.