

Statistical Presentation  
to the Radiological Devices Panel Meeting  
March 5, 2001

Marina Kondratovich, Ph.D.  
Mathematical Statistician  
Division of Biostatistics, OSB/CDRH/FDA

# OUTLINE

- Primary and Secondary Hypotheses
- Results of MRMC-ROC Analysis  
(areas under ROC curves)
- Sensitivity and Specificity based 50% confidence rate
  - Averaged
  - Individual for each radiologist
- Sensitivity and Specificity based on location
  - Averaged
  - Individual for each radiologist
- Summary

# Hypotheses

**The Primary Hypotheses:**  
Radiologists using RapidScreen would increase the detection of *lung cancers 9-30 mm in size*.

## **The Secondary Hypotheses:**

- I. Radiologists using RapidScreen would increase their detection of *lung cancers that had previously been missed by two screening radiologists*.
- II. Radiologists using RapidScreen would increase their detection of *lung cancers 9-14 mm in size*.

## Reading conditions:

- Independent without CAD
- Sequential without CAD
- Sequential with CAD


The MultiReader MultiCase (MRMC) method was used for statistical analysis of these 3 reading conditions. The *area under ROC curve* ( $A_z$ ) was considered as a measure of diagnostic accuracy.

In this analysis, *only information about confidence rate* for each film was used. The information about locations of marks and recommended actions was not used in this type of analysis.

# Statistical Hypotheses


Efficacy of the RapidScreen = Superiority (in terms of Az)

For superiority, the following statistical hypotheses with clinical important difference  $\delta$  should be considered:



If, on the basis of the results from the study, the  $H_0$  is rejected, then  $H_a$  is true, i.e. we may conclude that the reading with CAD is superior to the reading without CAD.

The hypothesis  $H_0$  is rejected if the lower limit of the one-sided 95% CI (or the lower limit of the two-sided 90% CI) of the difference  $Az_{with\ CAD} - Az_{without\ CAD}$  is larger than  $\delta$ .



**Which reading condition is the baseline reading:**

Independent without CAD or  
Sequential without CAD?

Independent without CAD :

- Corresponds to current practice;
- Gives less biased point estimate of difference.

Sequential without CAD:

- Gives more biased point estimate of difference;
- Provides more sensitive probe of difference because this reading condition is more correlated with Sequential with CAD reading.

Conclusion: Both reading conditions should be used in comparison with the reading with CAD.

For example, Primary Hypothesis (all cancers):

$$Az \text{ for Sequential with CAD} = 0.865$$

$$Az \text{ for Independent w/o CAD} = 0.829$$

$$Az \text{ for Sequential w/o CAD} = 0.835$$

Difference between

Sequential with CAD and Independent w/o CAD is  
0.036 with 90% CI: (0.016; 0.058)

Difference between

Sequential with CAD and Sequential w/o CAD is  
0.030 with 90% CI: (0.019; 0.043)

It was demonstrated that

$$Az_{with \text{ CAD}} - Az_{without \text{ CAD}} > 1.9\%$$

## Results of the MRMC-ROC Study (areas under ROC curves)

	w/o CAD est. of Az	with CAD est. of Az	It was demonstrated $AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > \underline{\Omega}$
<b>all cancers (50 films)</b>	0.829	0.865	$AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > 3.6\%$
Priors (18 films)	0.723	0.744	$AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > 1.0\%$
	2.1%		
Currents	0.861	0.901	$AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > 1.5\%$
	4.0%		
Small 9-14 mm (38 films)	0.798	0.848	$AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > 2.6\%$
	5.0%		
Medium 15-19 mm	0.840	0.870	$AZ_{\text{with CAD}} - AZ_{\text{w/o CAD}} > 0.5\%$
	3.0%		
Large 20-27 mm	0.889	0.889	No improvement
	0.0%		

**Question:** whether these improvements in the areas under ROC curves (  $\Omega$  ) are of clinical value.

The area under ROC curve ( $A_z$ ) is sensitivity averaged over all specificities.

The meaning of  $\Omega$  is a such increase in averaged sensitivities which is clinically important.

# Sensitivity and Specificity Based on 50% Confidence Rate

Primary Hypothesis: All cancer cases

	Sensitivity	Specificity
<b>Independent without CAD</b>	0.703	0.800
<b>Sequential with CAD</b>	0.777	0.778

+0.074

-0.022

Comparing the two reading conditions

Without CAD

$Se_{without CAD}, Sp_{without CAD}$

With CAD

$Se_{with CAD}, Sp_{with CAD}$

If  $Se_{with CAD} > Se_{without CAD}$


$Sp_{with CAD} > Sp_{without CAD}$

then with CAD is clearly preferred.


If  $Se_{with CAD} > Se_{without CAD}$

$Sp_{with CAD} < Sp_{without CAD}$

?




## Comparing the Reading w/o CAD and Reading with CAD



Conclusion:  $NPV_{\text{with CAD}} > NPV_{\text{w/o CAD}}$   
 $PPV_{\text{with CAD}} = PPV_{\text{w/o CAD}}$

# Sensitivity and Specificity Based on 50% Confidence Rate



	<i>PPV</i> ↑	<i>PPV</i> ↓
<i>NPV</i> ↑	4	8
<i>NPV</i> ↓	2	1


## Sensitivity and Specificity based on location

Primary Hypothesis: All cancer cases

	<b>Sensitivity</b>	<b>Specificity</b>
<b>Independent without CAD</b>	0.662	0.759
<b>Sequential with CAD</b>	0.683	0.730

+0.021                      -0.029

# Sensitivity and Specificity Based on Location



	<i>PPV</i> ↑	<i>PPV</i> ↓
<i>NPV</i> ↑	4	3
<i>NPV</i> ↓	4	4