

1 patient and cancer characteristics in a
2 screening population. The relevant
3 characteristics are cancer size, breast
4 density, finding type, histologic type, and
5 palpability. And these are very important, as
6 Dr. Petrick and Dr. Gwise discussed earlier,
7 when designing studies. Whether you are doing
8 enriched or non-enriched studies, it's very
9 important to look at these characteristics in
10 the final population of patients.

11 With regard to cancer size, the
12 approximate distribution of cancer size on
13 screening mammograms is, approximately, as
14 follows: 35 percent less than or equal to 10
15 millimeters, 60 percent less than or equal to
16 15 millimeters and 75 percent less than or
17 equal to 20 millimeters. Larger cancers are
18 more readily identified and characterized on
19 mammography.

20 With regard to breast density,
21 approximately 10 percent of patients have
22 either almost entirely fatty breasts or

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 extremely dense breasts, and approximately 40
2 percent of patients have scattered
3 fibroglandular densities and 40 percent have
4 heterogeneously dense breasts. Greater breast
5 densities are associated with lower
6 sensitivity for breast cancer detection and a
7 higher incidence of interval development of
8 breast cancer following a negative mammogram.

9 With regard to finding type,
10 approximately, 30 to 40 percent will be
11 masses, 30 to 40 percent will be
12 microcalcifications, 10 to 20 percent will be
13 a combination of a mass and
14 microcalcifications and 10 to 20 percent will
15 be architectural distortion or focal
16 asymmetry.

17 With regard to histologic type,
18 approximately, 70 to 80 percent are invasive
19 cancers and 20 to 30 percent are Ductal
20 Carcinoma In Situ or DCIS.

21 By definition, patients who undergo
22 screening mammography are asymptomatic. The

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 retrospective studies have shown that
2 approximately two to five percent of patients
3 who undergo screening mammography actually
4 have symptoms that were unknown at the time of
5 the examination.

6 There are two types of mammography
7 devices currently on the market: Screen film
8 devices and digital devices. And there are
9 two types of digital devices: The DR and CR.

10 CR is a device that uses a photostimulable
11 phosphor.

12 There are two standard mammographic
13 projections that are obtained of each breast::
14 the craniocaudal (CC) view on the left side
15 and the mediolateral (MLO) oblique views on
16 the right side, and these are typically
17 displayed side-by-side right and left.

18 Breast cancer is detected on the
19 basis of four types of mammographic findings:

20 The characteristic morphology of a mass; the
21 shape and spatial configuration of
22 microcalcifications; distortion of breast

1 tissue architecture; and asymmetry between the
2 left and right breast.

3 Mammography is unique among imaging
4 tests as it must be performed and even
5 interpreted in accordance with the Mammography
6 Quality Standards Act or MQSA, but MQSA does
7 not apply to mammography CAD devices.

8 With regard to interpretation of
9 mammograms, the CC and MLO projections in each
10 breast are considered complementary and
11 necessary for interpretation. Mammography
12 studies are always interpreted by examination
13 of the CC views from each breast in the side-
14 by-side manner and likewise for the MLO views.

15 When a finding is identified on a
16 single view, whether it is the CC or the MLO,
17 the corresponding region on the complementary
18 view is examined in order to confirm the 3-
19 dimensionality of the finding. Comparison
20 should always be made to prior mammograms when
21 these are available.

22 With regard to reporting

Formatted: Line spacing: single

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 mammographic examinations, mammographic
2 characteristics and findings are typically
3 recorded according to the American College of
4 Radiology, that's the ACR, Breast Imaging
5 Reporting and Data System, that's the BI-RADS
6 system. BI-RADS is meant to standardize the
7 language and descriptions used in mammography
8 reports.

9 With regard to reporting location,
10 a finding should always be triangulated so
11 that it's three-dimensional location within
12 the breast is known. The ACR BI-RADS Atlas
13 recommends using a clock face for each breast
14 and divides the breast into anterior, middle
15 and posterior thirds on the CC and MLO views.

16 The BI-RADS system has final
17 assessment categories. These were developed
18 and standardized into seven categories that
19 correspond to the reporting requirements in
20 the MQSA.

21 Category 0 means that a patient
22 needs additional imaging evaluation and/or you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 need to obtain prior mammograms. Category 1
2 is definitely negative. Category 2 is
3 definitely benign findings. Category 3 is
4 when a finding is probably benign. An initial
5 short-term follow-up is suggested.

6 But it's very important to note
7 that according to the ACR BI-RADS Atlas, a
8 finding placed in Category 3 should have less
9 than a two percent chance of malignancy. The
10 BI-RADS Atlas also states that it is
11 inadvisable to render such an assessment that
12 is Category 3 when interpreting a screening
13 examination.

14 And, this point is important when
15 you are reading literature on reader
16 performance with mammography CAD devices. You
17 have to always consider how many patients are
18 put into this category, and how it is used in
19 the final analysis.

20 Category 4 is for suspicious
21 abnormalities when biopsy should be
22 considered. These are sometimes broken up

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 into 4A, B and C, based upon low, intermediate
2 and moderate concern. Lesions put into the 4C
3 Category should be or are expected to be a
4 malignancy.

5 Category 5 are highly suggestive of
6 malignancy, and these are lesions that should
7 have a probability greater than or equal to 95
8 percent of being cancer. Category 6 is for
9 lesions that have a known biopsy proof of
10 malignancy.

11 It's important to look at the
12 performance measures of mammographers. There
13 is great variability in the published
14 literature with regard to sensitivity that
15 ranges from 60 to 100 percent, as well as
16 specificity that ranges from 35 to 98 percent.

17 Using the data in the Breast Cancer
18 Surveillance Consortium, you see a sensitivity
19 of about 79 percent and a specificity of about
20 90 percent for that very large number of
21 screening mammograms.

22 Mammographic sensitivity is lowest

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in patients with dense breasts and for small
2 masses. Given that greater than 99
3 percent of patients who undergo screening
4 mammography do not have a cancer, the task of
5 the mammographer is essentially to find a
6 needle in a haystack. And a practicing
7 radiologist may, therefore, perform detection
8 and analysis very rapidly, that is both tasked
9 to perform almost simultaneously when reading
10 mammograms in a clinical setting.

11 Cancers that are visible on
12 mammograms may draw the radiologist's
13 attention, and they may be dismissed with or
14 without formal description in the radiology
15 report. But, if a finding draws the
16 radiologist's attention and is missed, this
17 does not constitute an error of detection, but
18 is more likely to be considered an error of
19 analysis.

20 Any device designed to reduce
21 radiologists' errors should obviously focus on
22 the types of cancers that radiologists tend to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 miss. But, it also might be beneficial to
2 detect cancers that are missed simply because
3 the radiologist is asleep at the wheel.
4 Either way, increased detection should always
5 be weighed against increased false positives.

6 Regarding false positive
7 mammograms, these greatly outnumber actual
8 breast cancers found. Approximately, 10
9 percent of patients who undergo screening
10 mammography will be recalled for diagnostic
11 mammography. At the same time, approximately,
12 0.4 percent of patients who undergo screening
13 mammography actually have a breast cancer.

14 Approximately 50 percent of all
15 women had at least one false positive
16 mammogram over 10 years of screening. False
17 positive mammograms can cause increased dose
18 exposure, biopsy, complications associated
19 with biopsy, and unnecessary anxiety.

20 What about false negative
21 mammograms? Approximately, 20 percent of
22 cancers are missed on screening mammography.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 It's important to look at what are these 20
2 percent of cancers. About 10 percent are
3 thought to actually be visible on the
4 mammogram and, approximately, 10 percent are
5 thought to be not visible.

6 Therefore, there is obviously, room
7 for improvement to capture approximately up to
8 an additional 10 percent, perhaps more, of
9 cancers that are otherwise visible, but go
10 undetected or misclassified on screening
11 mammography.

12 Of the 10 percent of missed cancers
13 that are visible on mammography and are
14 missed, approximately five percent are errors
15 of detection, and approximately five percent
16 are errors of analysis. And I would emphasize
17 these are approximations.

18 Compared with cancers that are not
19 missed, of the 10 percent of cancers that are
20 missed and are visible, most are masses,
21 architectural distortions or focal
22 asymmetries. They tend to be smaller in size,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 and they tend to be present in patients with
2 denser breasts.

3 How can we reduce radiologist
4 errors when interpreting mammograms? One
5 method that has been tried is double reading.

6 Double reading of screening mammograms, that
7 is reading by two radiologists, has been
8 advocated as a way to increase radiologist
9 detection of cancers.

10 And there are a lot of clinical
11 studies on this that have shown that double
12 reading of mammograms improves radiologist
13 detection by about five to 15 percent. But
14 typically, with an associated increase in
15 recall rate of about five to 10 percent,
16 unless consensus double reading is used with a
17 recall rate that is much lower.

18 Double reading of screening
19 mammograms can capture both errors of
20 detection and errors of analysis. And the
21 published literature does show that double
22 reading can capture a substantial portion of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the approximate 10 percent of cancers that are
2 visible but currently go undiagnosed.

3 So where did mammography CAD
4 devices come from? Well, CAD devices have
5 been developed as a potential replacement for
6 a double reader.

7 The intended use of current
8 commercially available mammography CAD devices
9 is to reduce errors of detection. That is,
10 current commercially available mammography CAD
11 devices attempt to capture the approximately
12 five percent of cancers that are visible but
13 are not detected.

14 The potential for improved
15 detection, as I noted earlier, should always
16 be weighed against the potential for false
17 positives interpretations by radiologists
18 using CAD devices. And just to give some
19 context, the screening cancer incidence is
20 approximately four per 1,000.

21 Current commercially available CAD
22 devices place at least, approximately, two

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 marks per patient, even assuming 100 percent
2 sensitivity for the CAD device and assuming
3 that marks are placed on both the CC and ML
4 view for each cancer. There will still be
5 about 249 false positive marks for every true
6 positive mark.

7 It's therefore important to measure
8 how easy or difficult it is for radiologists
9 to dismiss these false positive CAD marks.
10 It's also important to measure the effect of
11 false positive marks on potentially
12 distracting the radiologist from other
13 findings that may not be marked.

14 I just want to review mammography
15 CAD devices that have been approved by the
16 FDA. There are four mammography CAD systems
17 that have been approved through the PMA or
18 Pre-Market Approval application process. The
19 first approval order was issued in 1998. The
20 first-of-the-kind device was the subject of a
21 Radiological Devices Advisory Panel meeting
22 that was held on May 11, 1998, and the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 transcript and other information are all
2 available on-line.

3 All devices were first approved for
4 use with digitized versions of screened filmed
5 mammograms -- excuse me, of screen film
6 mammograms that were obtained for screening
7 purposes.

8 Further PMA supplements have been
9 approved over time to expand the use of CAD
10 devices to operate on digitized diagnostic
11 screen film mammograms, at least the CC and
12 MLO views, as well as mammograms obtained on
13 full field digital mammography or FFDM
14 devices. PMA supplements have also been
15 approved for modified software versions of CAD
16 algorithms.

17 The labeling of currently approved
18 mammography CAD devices have an indications
19 for use or IFU that is similar to the
20 following: They are intended to identify and
21 mark regions of interest on routine screening,
22 as well as the CC and MLO views of diagnostic

1 mammograms, in order to bring these regions of
2 interest to the attention of the radiologist
3 after the initial reading has been completed
4 and to assist the radiologist in minimizing
5 observational oversights by identifying areas
6 on the original mammogram that may warrant a
7 second review.

8 What about the data and information
9 that was provided for original approval of
10 mammography CAD devices? At the time of
11 original approval in 1998, there was limited
12 experience with use of mammography CAD devices
13 by radiologists in clinical practice.

14 The data that served as the basis
15 for approval for currently approved
16 mammography CAD devices typically included
17 four basic components: The first component
18 was standalone performance on missed cancers.

19 And missed breast cancers were identified by
20 obtaining prior mammograms from patients with
21 newly diagnosed cancer, that is, patients with
22 interval cancers; and determining if the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 cancers were actually visible in retrospect
2 and should have led to a clinical action.

3 Standalone performance on these
4 missed cancers was considered a surrogate for
5 the ability of the device to detect difficult
6 findings. And it was used to estimate the
7 maximum potential reduction of detection
8 errors if radiologists actually used the
9 device in practice.

10 The second component was standalone
11 performance on cancers detected at routine
12 screening mammography. And this was
13 considered a measure for the ability of the
14 device to detect more obvious and intermediate
15 level of difficulty findings.

16 The third component was standalone
17 performance on normal screening mammograms to
18 determine the rate of false positive CAD marks
19 on normal cases. The fourth component was
20 screening exams, and these databases may have
21 or may have not been enriched with some
22 cancers, and these were used to determine the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 potential increase in recall rate resulting
2 from use of the CAD devices.

3 What has been learned since the
4 original approval of these devices? Well,
5 there is a large body of literature out there
6 on standalone performance of mammography CAD
7 devices. There is also a large body of
8 literature on the subject of reader
9 performance testing of mammography CAD devices
10 where radiologist's performance is measured
11 both without and with use of the CAD device.

12 And these reader studies have
13 employed two different general designs, some
14 of which has been discussed previously:
15 Retrospective clinical performance testing and
16 prospective clinical performance testing.

17 The retrospective clinical studies
18 use radiologist interpretations that are not
19 part of actual clinical practice. The
20 prospective studies use radiologist
21 interpretations that are part of an actual
22 clinical practice. And the prospective

1 studies include two basic designs: The
2 sequential design and the historical control
3 design.

4 With the sequential design, the
5 images are presented to the radiologist
6 without CAD information. It requires
7 interpretation and then presents the same
8 images with CAD markings and allows the
9 radiologist to modify the assessment.

10 With the historical control design,
11 you compare radiologist performance over a
12 period of time without CAD devices to
13 radiologist performance over a period of time
14 after introduction of CAD devices. And it,
15 obviously, uses different patients.

16 Here are some of the key points
17 from the published literature which are
18 contained in more detail in the Panel briefing
19 document. Standalone testing has shown very
20 high sensitivity to mark calcifications but
21 much lower sensitivity to mark masses,
22 architectural distortions, or focal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 asymmetries. Standalone testing has also
2 shown a false positive mark rate of between
3 two and four marks per patient.

4 Reader performance testing has
5 shown conflicting results for detection of
6 invasive cancers. Reader performance testing
7 has also shown a trend toward CAD improving
8 radiologist's detection of calcifications,
9 especially DCIS. And reader performance
10 testing has shown an increase in recall rate
11 when using CAD devices in some studies, and
12 some of these are statistically significant
13 increases.

14 Now, I want to discuss some of the
15 clinical testing issues specific to
16 mammography CAD devices that will be the
17 subject of the questions that we have for the
18 Panel. And I tried to put on the top of the
19 slide the question numbers that correspond to
20 some of the comments and some of the questions
21 in the slides.

22 With regard to ground truth, this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 has been touched on previously, but this is
2 crucial for both standalone and reader
3 performance testing. Ground truth includes
4 whether or not the patient has a breast
5 finding, whether or not the patient has one or
6 more benign or malignant findings, the precise
7 location and extent of each finding on each
8 view, and the BI-RADS descriptors and final
9 assessment of each finding.

10 Ground truth for cancer is
11 determined by biopsy or surgery. Ground truth
12 for benign findings is determined by biopsy/
13 surgery or by one year follow-up mammogram.
14 And ground truth for normal is determined by a
15 one year follow-up mammogram.

16 Ground truth for the location and
17 extent, that is the lesion boundary of a
18 finding, is determined typically by a panel of
19 expert radiologists and can be annotated
20 either manually or digitally on an image.

21 What about standalone performance
22 testing? Standalone performance is highly

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 dependent on case selection, including the
2 following factors: the precise mammographic
3 characteristics in the case set, including
4 finding size; the pathologic type of lesions;
5 the number of masses versus
6 microcalcifications; and the distribution of
7 breast densities.

8 Standalone performance testing also
9 depends on the precise method of ground truth
10 determination for the location and extent of
11 disease; the precise scoring metric, for
12 example, per lesion, per CAD mark, per patient
13 versus per view; and the precise scoring
14 methodology.

15 For example, using overlap criteria
16 on the actual CAD mark or the actual region
17 identified by the CAD system itself; that is,
18 the algorithm segmentation. And the algorithm
19 segmentation may not be displayed to the user
20 typically.

21 Standalone performance can be done
22 using a larger database than used for reader

1 performance testing, and that is certainly a
2 potential advantage. And using a larger
3 database may allow meaningful stratified
4 analysis on clinical mammographic and
5 pathologic subgroups, which may be important
6 information, because it may influence user
7 confidence in the ability of the device to
8 detect findings in each of the important
9 subgroups.

10 You can also do -- use stratified
11 measures and this may include mammographic
12 finding types, such as masses,
13 microcalcifications, architectural distortions
14 and focal asymmetries, the different
15 pathologic types. You could do it by size.
16 You could do it by breast composition. In
17 particular, it might be important to look at
18 the number of small masses, that is, masses
19 less than 10 millimeters in size, in
20 particular in patients with dense breasts
21 since these are among the most difficult
22 findings to detect.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Both overall and stratified
2 standalone performance can be reported on a
3 per lesion, per view, per breast or per
4 patient basis. And one of the things that we
5 would like input from the Panel is given that
6 clinical actions following mammography are
7 finding-specific, what are the advantages and
8 disadvantages of each of the above reporting
9 measures?

10 Without standardized methodologies
11 for case selection, ground truth, scoring
12 metric, scoring methodology and reporting,
13 it's important to keep in mind that it may be
14 difficult, it may be invalid to compare
15 performance between devices from different
16 manufacturers or even different versions of a
17 device from the same manufacturer. And that's
18 another question that we do want some input
19 from the Panel.

20 It may also be difficult to account
21 for differences in detection location and
22 number of both true positives and false

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 positives in such comparisons.

2 Now, I would like to discuss reader
3 performance testing. While standalone
4 performance testing indicates how well the
5 device marks locations of interest, in the
6 absence of radiologist interaction, it does
7 not measure the safety or effectiveness of the
8 device for its intended uses and conditions of
9 use by a reader.

10 There are several types of reader
11 performance tests that are designed to
12 determine the impact of a CAD device on reader
13 performance. Some of these have been
14 discussed previously.

15 The prevalence of breast cancer in
16 a screening population, as I noted previously,
17 is, approximately, 0.4 percent. Following the
18 Least Burdensome approach, reader performance
19 testing may be accomplished using a so-called
20 enriched dataset when you enrich the dataset
21 with a significantly greater percentage of
22 patients with breast cancer. So, you have a

1 final population of patients with a prevalence
2 of cancer that is much higher than that in a
3 real screening population. And this is the
4 population of patients that you use for your
5 testing.

6 Ignoring prevalence, if the cancer
7 and non-cancer cases that are used for
8 enrichment, otherwise have clinical
9 mammographic and pathologic characteristics
10 that are typically seen in a screening
11 population, such testing simulates the so-
12 called field test, where a field test is
13 clinical assessment of a system in real
14 practice, you know, real life, real time.

15 However, the much higher prevalence
16 in enriched datasets can introduce bias, as
17 Dr. Gwise has previously discussed. And a
18 question that we're interested in looking at
19 is can reader performance testing using these
20 enriched datasets give an estimated measure
21 for device effectiveness as would be seen in
22 clinical practice?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 And this is also where, as I think
2 Dr. Gwise mentioned earlier, we're interested
3 in what effect size would be needed in order
4 to achieve this.

5 When enrichment is performed with
6 only difficult cases, the testing is usually
7 referred to as a stress test. In this
8 situation, the test dataset is enriched with
9 primarily or only difficult cases that
10 challenge the readers. Stress testing that is
11 done with only difficult cases will not
12 capture information about the effect of CAD on
13 cases that radiologists don't tend to miss.
14 And this may lead to an incomplete assessment.

15 Another thing we are interested in
16 getting input from the Panel is is stress
17 testing alone sufficient to measure safety and
18 effectiveness?

19 What about study endpoints for
20 reader performance testing? Both overall and
21 stratified reader performance, as I mentioned
22 previously, can be reported per lesion, per

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 view, per breast or per patient. And again,
2 given that the clinical action following
3 mammography are finding-specific, this is the
4 same question we had for standalone
5 performance, how critical is it to account for
6 reader location accuracy when you are looking
7 at reader performance testing?

8 What about stratified analysis?
9 Meaningful stratified analysis may include
10 breast density, finding size, finding type, or
11 histologic type. And again, this is something
12 we would like input from the Panel.

13 And now, I just want to touch on,
14 this is the very end of my talk, some other
15 issues for mammography CAD devices. Screen
16 film and Full Field Digital Mammography (FFDM)
17 devices have different spatial and contrast
18 resolution and different noise
19 characteristics.

20 Full Field Digital Mammography
21 systems also vary between one another in
22 spatial and contrast resolution because of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 differences of the types of solid state
2 detectors, the pixel sizes, and the quantum
3 and electronic noises. Different FFDM
4 manufacturers may use different technologies
5 and different image processing algorithms and
6 techniques.

7 Therefore, the standalone and
8 reader performance testing of CAD devices on
9 screen film images may differ from testing
10 results on FFDM images. This may also apply
11 to testing on different FFDM devices. And
12 again, a question that we would like input
13 from the Panel is is there a reason why
14 testing of CAD devices on a new or modified
15 image input, whether it be a new digitizer for
16 film or a new FFDM device, should that testing
17 be any different than the standalone and
18 reader performance testing already discussed?

19 Mammography CAD can be clinically
20 implemented as a second reader or as a
21 concurrent reader or perhaps other reader
22 paradigms. Second reading can only increase

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 reading time, while concurrent reading may
2 reduce reading time.

3 Another question we have for the
4 Panel is how are mammography devices currently
5 used in clinical practice? And then the last
6 bulleted point there, mammography CAD has been
7 reported with very high sensitivity for
8 calcifications. Is there a clinical role for
9 mammography CAD as a concurrent reader for
10 calcifications?

11 And I'll end the talk there.

12 CHAIRMAN GLASSMAN: Thank you very
13 much, Dr. Smith. Does the Panel have any
14 questions for Dr. Smith? Dr. Berry, yes?

15 DR. BERRY: So, Dr. Smith, it's
16 easy to see that CADs could improve
17 sensitivity. Are there any CADs that the
18 settings of the algorithms that improve
19 specificity as well?

20 DR. SMITH: As far as approved
21 indications for use, there is no approved
22 indication for use that would deal with that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 issue. But that certainly would be a possible
2 implementation.

3 DR. BERRY: So I don't understand.

4 Are devices that have been approved, do we
5 know what the sensitivity and specificity are
6 under specific circumstances?

7 DR. SMITH: I would have to ask you
8 if -- there are specific -- four different
9 specific devices. Each device has its own
10 performance testing and data that is provided
11 in a submission. The current CAD device
12 submissions have, essentially, not
13 exclusively, standalone performance testing
14 for those measures and testing with readers in
15 order to measure recall rates. But there is
16 no data available in current submissions for
17 reader performance testing that would give you
18 the measures that you are asking about.

19 DR. BERRY: Okay. I'm still
20 somewhat confused, but maybe I'll be
21 enlightened as the time goes along.

22 DR. SMITH: Well, I guess, if you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 are asking about the published literature,
2 that's different than what's in FDA
3 submissions.

4 DR. BERRY: So I'm trying to get a
5 feeling for is it sensitivity/specificity that
6 drives the clearance process? And are there
7 devices that show both, you know, improvements
8 in both? Is it necessary -- are there devices
9 that have improved on the basis of let's say a
10 statistically significant improvement in
11 sensitivity without loss of specificity?

12 DR. SMITH: All currently approved
13 devices, essentially, use the four components
14 that I described in the slide. It's strictly
15 standalone performance testing and reader
16 testing only to look at recall rates. So
17 there is no such data that has been used to
18 approve current devices.

19 DR. BERRY: Another question. So
20 just a follow-up. Is there any -- what is the
21 range of improvement in sensitivity or
22 specificity that one might see in these

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 devices?

2 DR. SMITH: Are you talking about
3 with a reader using the device?

4 DR. BERRY: Yes, yes.

5 DR. SMITH: There is a lot of
6 literature out there on that. And there has
7 been some conflicting results, depending on
8 the methodology used for the studies. There
9 have been basic -- two basic reader studies,
10 the sequential study, where readers read
11 without and then the same reader with CAD.
12 Those studies have not shown statistically
13 significant improvements, but the numbers of
14 cancers have typically been quite small.

15 The studies that have been
16 published have typically shown a trend toward
17 improved detection of cancers that manifest as
18 microcalcifications. But again, the numbers
19 of cancers are quite small, so they have not
20 been able to demonstrate statistically
21 significant results.

22 DR. BERRY: So there is currently

1 no benchmark then for sensitivity or
2 specificity for these products?

3 DR. SMITH: Well, that's one of the
4 reasons that we have the Panel here is to give
5 us some advice on this.

6 DR. BERRY: Thank you.

7 CHAIRMAN GLASSMAN: Dr. D'Orsi
8 next.

9 DR. D'ORSI: Hi. Is it correct to
10 assume that in the studies received by the
11 FDA, there is no statistically significant
12 difference in the data presented to the FDA
13 for AUC?

14 DR. SMITH: We have not gotten AUC
15 data with these submissions. Again, the data
16 that has been used in the currently approved
17 devices is standalone performance testing.

18 DR. D'ORSI: So there is no data on
19 that part?

20 DR. SMITH: That is correct.

21 CHAIRMAN GLASSMAN: The last
22 question.

1 DR. BOURLAND: I have a question
2 about the BCSC database. Who is -- has
3 accessibility to that database? For instance,
4 manufacturers. Are those images fully
5 characterized, so to speak, available in DICOM
6 form? And I'm sorry, the last part of the
7 question is what's the image quality assumed
8 for the images that, for instance, are inputs
9 to CAD? For instance, are those MQSA
10 compliant, so to speak, within the range or
11 outliers allowed?

12 DR. SMITH: This is a database
13 that, I think, is largely funded by the
14 National Cancer Institute. It currently uses
15 seven registries that link mammography and
16 pathologic findings. The database is
17 available -- the characteristics of the
18 patients and numerous other features of the
19 database are publicly available on-line. I do
20 not know about the availability of the actual
21 images.

22 They are all obtained from MQSA-

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 certified facilities. But I don't know if the
2 actual images themselves are available. I
3 presume most of the images that are in that
4 database are screen film. I don't know if
5 they have ever been digitized and available.

6 CHAIRMAN GLASSMAN: I want to move
7 on. We may have other questions, but we do
8 have a chance later on in our general
9 discussion to come back to these issues. So I
10 would like to now proceed with the first of
11 two Open Public Hearing sessions for today's
12 meeting. The second Open Public Hearing
13 session will follow the FDA presentation on
14 colon CADs this afternoon.

15 Ms. Wersto will now read a
16 statement prepared for Open Public Hearings.

17 EXEC. SEC. WERSTO: Both the Food
18 and Drug Administration, FDA, and the public
19 believe in a transparent process for
20 information gathering and decision making. To
21 ensure such transparency at the Open Public
22 Hearing session of the Advisory Committee

1 meeting, FDA believes that it is important to
2 understand the context of an individual's
3 presentation.

4 For this reason, FDA encourages
5 you, the Open Public Hearing speaker, at the
6 beginning of your written or oral statement to
7 advise the Committee of any financial
8 relationship that you may have with the
9 sponsor, their products, and if known, any of
10 their direct competitors.

11 For example, this financial
12 information may include a sponsor's payment of
13 your travel, lodging, or other expenses in
14 connection with your attendance at the
15 meeting.

16 Likewise, FDA encourages you at the
17 beginning of your statement to advise the
18 Committee if you do not have any financial
19 relationships. If you choose not to address
20 this issue of financial relationships at the
21 beginning of your statement, it will not
22 preclude you from speaking. Thank you.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 CHAIRMAN GLASSMAN: I would like to
2 remind the public observers, at this meeting,
3 that while this portion of the meeting is open
4 to public observation, public attendees may
5 not participate, except at the specific
6 request of the Chair.

7 I would like to ask, at this time,
8 that persons addressing the Panel come forward
9 to the microphone and speak clearly, so the
10 transcriptionist can make an accurate record
11 of what you say. Please provide an electronic
12 copy of your talk to the Executive Secretary
13 for use by the transcriptionist to help
14 provide this accurate record.

15 Prior to the meeting, we received
16 formal requests to speak during today's Open
17 Public Hearing session. Our first speaker is
18 Heang-Ping Chan, American Association of
19 Physicists in Medicine. And you will have
20 five minutes. It is my understanding that
21 after four minutes, you will get an orange
22 light there, so that you have a minute to sum

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 up. Thank you.

2 DR. CHAN: Good morning. I'm here
3 to present a statement on CAD on behalf of the
4 American Association of Physicists in Medicine
5 or AAPM. The statement was drafted by the
6 AAPM CAD Subcommittee. Here is a list of the
7 subcommittee members. I'm the Chair of the
8 subcommittee. Note that four of the members
9 listed in the footnote did not propose or
10 modify this statement.

11 Here is our statement. AAPM
12 recognizes that CAD, including Computer-Aided
13 Detection, Computer-Aided Diagnosis and more
14 broadly Computer-Assisted Image Analysis, will
15 be an indispensable part of diagnostic
16 medicine in the near future.

17 Although commercialization of CAD
18 systems is the bridge to clinical use, CAD
19 research in academia has been and will
20 continue to be a driving force for this
21 progress.

22 Although CAD systems have been

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 commercialized in several areas, the reported
2 performances are excellent in some areas, for
3 example, microcalcification detection on
4 mammograms, but modest in others. For
5 example, breast mass detection, lung nodule
6 detection, indicating that CAD is promising,
7 but the full potential of CAD has not yet been
8 realized. Thus, further research and
9 development of CAD should be strongly
10 encouraged.

11 Continued funding support for
12 research and development of CAD technologies
13 will be vital for improvement in current CAD
14 applications and development of CAD in new
15 areas.

16 FDA-approval of a CAD system should
17 be conditioned upon appropriate post-FDA-
18 approval prospective evaluations of the CAD
19 system in clinical practice. The majority of
20 prospective clinical trials today indicated
21 the promise of CAD in screening mammography,
22 despite some negative reports.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Standardization of study design and proper
2 evaluation technologies are clearly important
3 issues that need to be addressed.

4 Standardization of CAD evaluation
5 technologies includes, but is not limited to,
6 the following: clear definition of task,
7 patient population, reader training for
8 reading with CAD, without and with CAD reading
9 design, definition of truth, data analysis
10 methods, identification of biases and
11 variances, and endpoint for assessment of the
12 success or failure.

13 Journals should offer fair
14 opportunities for publishing rebuttal or
15 critical review of published studies. Quality
16 assurance procedures should be established for
17 CAD systems implemented in clinical use.

18 Radiologists should obtain training
19 on the proper interpretation of the
20 information provided by a specific CAD system
21 before using it clinically. Procedures should
22 be established to ensure that a CAD system is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 used as labeled.

2 Although an accurate CAD system can
3 potentially be used as a first reader, such
4 use should be first proven by properly
5 designed prospective studies to evaluate its
6 efficacy and approved by FDA. Vendors should
7 be prohibited from advocating off-label use of
8 a CAD system without proof and FDA-approval.

9 Thank you for your attention.

10 CHAIRMAN GLASSMAN: Thank you, Dr.
11 Chan. Next is Stephen Vastagh from NEMA.

12 MR. VASTAGH: Distinguished Chair,
13 Members of the Panel, Madam Secretary, FDA
14 staff, good morning. My name is Stephen
15 Vastagh. I am MITA staff liaison to the CAD
16 manufacturers. MITA would like to thank FDA
17 for holding this special Panel meeting on CAD.

18 I am employed by MITA, which is a
19 membership organization of the manufacturers.

20 Thus, I am and the organization is funded by
21 the dues of the manufacturers.

22 MITA has represented the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 manufacturers of medical imaging and therapy
2 systems for over 25 years. MITA has been on
3 the leading edge of technical issues. We
4 developed the early x-ray standards under the
5 leadership of a former director of the Office
6 of Device Evaluation.

7 Currently, we develop, maintain and
8 publish DICOM, the world's foremost standard
9 on imaging communication. Also, FDA experts
10 were along side MITA member experts in
11 international standardization of medical
12 devices.

13 Today, we are here on behalf of the
14 CAD group of MITA. In the interest of time, I
15 will not review this slide, which was our
16 organization and the CAD group in our
17 organization.

18 Our goal is to provide the MITA
19 members' views. Today, speakers from
20 Fujifilm, Hologic and tomorrow, from Phillips
21 Healthcare, GE Healthcare, and Medipattern
22 will present these views. Our presenters have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 significant experience in the CAD field.
2 These experts personally and their companies
3 institutionally are just as much compelled to
4 produce safe and effective CAD products as FDA
5 is compelled to review whether they are safe
6 and effective. The companies have a vested
7 interest that all their products be safe and
8 effective.

9 The scope of the CAD product is
10 defined by the claims the manufacturer makes.

11 These claims are interpreted into
12 manufacturers operating instructions to define
13 how the CAD product is to be used. MITA hopes
14 that one of the outcomes of this Panel meeting
15 is a common understanding of how CAD products
16 are to be used, and the relative roles of the
17 physicians and the CAD product.

18 We hope this will give confidence
19 to evaluate CAD products to the extent of
20 their claims and not require studies regarding
21 possible off-label uses. It is likely that
22 CAD technology is on the verge of increased

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 growth and diversification.

2 New technologies are developed to
3 accomplish similar functions; therefore, it is
4 important to evaluate the products on the
5 basis of their performance and not on the
6 basis of technology.

7 In view of the advance of CAD
8 technology, it is easy to reach a conclusion
9 that CAD should improve the performance of
10 radiologists. However, the manufacturers of
11 CAD products do not intend such function,
12 unless they explicitly claim so.

13 Given that physicians make
14 decisions, which will hopefully be made clear
15 during these two days, manufacturers should
16 not need to demonstrate improvement of
17 radiologist performance unless claimed.
18 Indeed, this is not a requirement for other
19 radiology devices.

20 Our presentations today and
21 tomorrow conclude with several
22 recommendations:

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 FDA to look at the state of the art
2 of CAD, consult with physicians and assess the
3 risk of CAD product submissions in light of
4 clinical experience.

5 Secondly, issue official guidance
6 concerning the separation between Class II and
7 Class III CAD devices and requirements. Use
8 the role and involvement of the physicians as
9 the guide to evaluate risk. Use a
10 collaborative guidance development process
11 that allows industry to provide input to the
12 process.

13 Recognize that CAD guidance
14 development generally takes two to three
15 years, as has been mentioned previously.
16 Also, recognize that industry is in great need
17 of a uniform and transparent interim solution
18 or an extraordinary expedited process for the
19 guidance development.

20 These recommended actions will help
21 introduce innovations and more advanced
22 technologies to patient care. With these

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 innovations in clinical practice, we can
2 expect new clinical studies to advance this
3 important science. These actions will also
4 help to differentiate CAD from products that
5 enhance and measure images, but are not CAD.

6 MITA looks forward to working with
7 the FDA on this important issue. Thank you
8 for your attention.

9 CHAIRMAN GLASSMAN: Thank you very
10 much. Next, Robert Nishikawa from the
11 University of Chicago.

12 DR. NISHIKAWA: Thank you very
13 much, and I thank you for the opportunity to
14 come and present some information here. This
15 slide just describes my activities with
16 Fujifilm Carestream and Hologic, but I'm here
17 representing myself and the views I'm going to
18 express here are those -- not necessarily
19 those of my colleagues, collaborators or
20 employer.

21 So evaluation of -- clinical
22 evaluation is difficult period. It's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 extremely difficult if cancer prevalence is
2 low. There is no perfect evaluation
3 methodology, some compromises need to be made
4 and this makes it difficult for deciding on
5 FDA-clearance and approval process.

6 Part of the problem, is there is a
7 view that the clinical data on CAD in
8 screening mammography is conflicting. What I
9 want to show here is that this is not true.
10 The clinical data actually presents a
11 consistent view that there is a benefit for
12 using CAD.

13 So there are 10 clinical studies.
14 I grouped them into two categories. One I'm
15 going to call longitudinal, which I think is
16 described as historical controls, and cross-
17 sectional, which is the sequential study.

18 If you look at the longitudinal
19 studies, you see the increase in the cancer
20 detection rates around one to two percent,
21 except for this one outlier here. The cross-
22 sectional studies present an increase whenever

1 cancer is detected at around 10 percent.

2 Oh, I see one more thing. And
3 also, if you look at the Gromet Study, that
4 actually measures sensitivity. It's about 8.2
5 percent, which is consistent with these
6 numbers. And the Fenton sensitivity increases
7 about 8.5 percent when you correct for a bias.

8 And it's also consistent with these studies.

9 So we developed in our lab a
10 simulation model to look at the effects of
11 these two different evaluation methods. We
12 start with, approximately, 50,000 women. We
13 screen them every year and all cancers grow at
14 the same rate.

15 So as one would expect, when
16 screening starts, there is a prevalence
17 effect, the cancer detection rate is high, but
18 it goes down to some steady state level. For
19 the same reason, we introduce CAD. There is a
20 prevalence effect, and then you hit the steady
21 state level.

22 But if you compare this curve to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 this curve, the difference is extremely small.

2 So if you're measuring changes in cancer
3 detection rate, it's very difficult to do. In
4 fact, if the interval cancer rate was 0, these
5 two lines would be collinear.

6 In a sequential reading method, you
7 are trying to measure this difference here and
8 it's actually quite doable.

9 Let's see, the previous example was
10 with minimizing all of the sources of
11 variability. You now add that cancers grow at
12 different rates, so different numbers of
13 cancers are present in the population each
14 screening year.

15 This does not take into account
16 variability. The radiologist will make this
17 graph even more noisy. So now, what we are
18 trying to do is measure some -- in
19 longitudinal studies, some time point over
20 here versus some time point over here. And
21 you can see depending on how the data is --
22 falls, you may get -- measure an increase, you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 may measure no change, or you may measure a
2 decrease in cancer detection rate.

3 But if you look at the difference
4 between the yellow line and the white line, it
5 is fairly consistent, even though the lines
6 are going up and down. So, that in the cross-
7 sectional methods or sequential method, you
8 are able to measure this difference, but in
9 the longitudinal or historical controls, you
10 can't measure a difference.

11 The reason here is fundamental.
12 See the -- the goal of CAD is not to find more
13 cancers. It's to find cancers earlier. So if
14 you measure CAD at detection rate, you are
15 making the assumption you are going to find
16 more cancers, which is not true.

17 And so using the longitudinal
18 method, I believe, is flawed. Measuring this
19 method, in fact, you're looking at how you can
20 decrease the false negative rate, which is
21 exactly what our CAD is trying to do. So this
22 is a -- measures the effectiveness of CAD more

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 directly than measuring the longitudinal --
2 than measuring a cancer detection rate and a
3 historical control method.

4 There is one publication that says
5 that, basically, CAD is bad to use, and I just
6 want to address this quickly. There was
7 actually six different results that I
8 summarized here from the study. I can't
9 discuss the reasons I have for all these
10 points, but I will point out, for example, the
11 sensitivity increase is not statistically
12 significant.

13 They quote a value of 4.5. This is
14 a biased estimate of the increase. I estimate
15 it's close to 8.5 percent. I don't know.
16 Sorry, 8.2 percent. I don't know if that's
17 statistically significant, but it's
18 substantially bigger than 4.5. I just
19 addressed the cancer detection rate to
20 increase. That's the wrong endpoint to use.

21 Many people have commented that the
22 high recall rate is because the readers were -

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 - the radiologists were inexperienced.

2 CHAIRMAN GLASSMAN: Excuse me, Dr.
3 Nishikawa, your time is up.

4 DR. NISHIKAWA: Okay. Okay.

5 CHAIRMAN GLASSMAN: Thank you.
6 Next is Dr. Carl Jaffe from the National
7 Cancer Institute.

8 DR. JAFFE: Thank you. I'm Dr.
9 Carl Jaffe, Branch Chief of Cancer Imaging
10 Program, Division of Cancer Treatment and
11 Diagnosis, NCI. I have no financial
12 disclosures to report, but I am speaking as an
13 individual, rather than for my colleagues or
14 my division.

15 The imaging test, I think we have
16 already had outlined, which is primarily the
17 issue of detection and classification. We
18 have already determined that expert observers
19 are good, but inconsistent. That effective
20 technology may reassure the best and
21 conceivably could lift the performance of
22 average performers. Man-machine systems, like

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 any graphic user interface, to require run-in.

2 The amount of time and the amount of
3 experience of that is still to be determined.

4 ROC analysis depicting the full range of
5 observer talents is very expensive and often
6 impractical.

7 We have already seen this graph
8 from Craig Dean's work on mammographers.
9 There is 110 mammographers nationwide
10 published in the 1990s. And it shows you that
11 a 10 percent false positive rate that the
12 sensitivity for observers ranges from 50 to
13 about 95 percent.

14 This tells you that people who
15 might volunteer for a reader study may already
16 be in the higher level. They may be a bias
17 sample set, because the actual clinical
18 practitioners are often falling well within a
19 range that is much below what might be
20 occurring in a confined reader environment,
21 that is in a testable environment.

22 The imaging test needed for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 improving critically important are for therapy
2 assessment. And quantitative monitoring of
3 change over time turns out to be as yet not
4 discussed in this particular meeting, but here
5 performance is mediocre to poor on the
6 observer's part. And hence, compromise work-
7 arounds like the RECIST technique that was
8 used, which are thought to be quantitative
9 are, in fact, semi-quantitative.

10 We can see that easily when you
11 look at a progress over time between two
12 images on the same patient and the observers
13 were unable to measure this consistently with
14 a range of over 114 percent variations on
15 that. We can also tell that ground truth is
16 almost unknowable in many of the situations,
17 particularly with lung, except in phantoms or
18 simulations.

19 Quantitative solutions are an
20 engineering exercise. They are not hard
21 science. So they must evolve over time and
22 must be proven by clinical use. One benchmark

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 is the comparative human performance, and that
2 is an effort that is being made by our
3 division under what is called the Lung Image
4 Database Consortium.

5 Validation is a clinical field
6 exercise. Clinical correlates and outcome
7 data are absolutely needed and, hence, the
8 need for big public databases that are
9 available for every developer of CAD systems
10 to test against.

11 You can see here, for instance, it
12 may be a little difficult to see the outline
13 of the observers, but a large dataset now has
14 become available in which the variation of
15 observers for the measuring of something that
16 is quite an obvious target become evident.
17 And we have seen a little bit of that earlier
18 in one of the presentations.

19 All of those, in a sense, are
20 truths of a form. What we are doing at NCI is
21 to try to help the -- and encourage the
22 development of effective CAD by developing

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 open public databases that contain the meta-
2 data and the clinical data that can be
3 correlated with it, so that new methods are
4 tested against something that at least the
5 performance becomes gradually evolutionarily
6 known.

7 For instance, in this variable
8 DICOM available dataset, the original first
9 collection was actually the virtual
10 colonoscopy set that was provided by the
11 Department of Defense Study of 2003.

12 So what are the unresolved
13 questions? It is how big must this referenced
14 image dataset be? What are the training and
15 testing components? And who possesses that
16 dataset? Who maintains it? Who makes it
17 available? And more importantly, who judges
18 performance and by what metric? Thank you.

19 CHAIRMAN GLASSMAN: Thank you very
20 much. Our next speaker is Steve Worrell from
21 Riverain Medical.

22 MR. WORRELL: Good afternoon. I do

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 not have slides today. I'll just be reading
2 from a prepared statement. My name is Steve
3 Worrell. I'm Vice President of Research at
4 Riverain Medical.

5 To date, Riverain Medical is the
6 only company to have demonstrated the safety
7 and effectiveness of a Computer-Aided
8 Detection device intended to identify and mark
9 regions of interest on frontal chest
10 radiographs. The device is intended for the
11 use as an aid only after the physician has
12 performed an initial interpretation of the x-
13 ray.

14 As a company who has committed
15 significant resources to maintain conformance
16 with the strictest Category 3 regulatory
17 requirements, we, obviously, have a strong
18 belief that Class III Regulation is necessary
19 to mitigate risk associated with today's chest
20 CAD devices.

21 PMA processes insure that CAD
22 devices are properly developed, tested,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 marketed, sold using approved indications of
2 use and clinically validated performance
3 claims. Chest CAD technology has evolved into
4 a complex and integral part of the practice of
5 medicine.

6 It's only through appropriate
7 controls that chest CAD has become and will
8 continue to be a trusted element of the
9 clinical work flow. Chest CAD is used to aid
10 physicians in discovering important findings
11 that would otherwise go undetected.

12 Riverain believes in the importance
13 of the PMA process and Class III Regulation of
14 chest CAD technology. Chest CAD is
15 appropriately classified as a Class III
16 medical device. It is part of a generic type
17 of device that is intended for a use which is
18 substantial, important in preventing
19 impairment to human health and its use can
20 prevent a risk to patients.

21 For example, false positives can
22 link to unnecessary CT exams. Currently CAD

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 devices are regulated as both Class II and
2 Class III. While this apparent inconsistency
3 in FDA Regulation can be explained, the
4 distinction between the regulatory
5 classifications is not clear.

6 In fact, the intended use and
7 design of Class II chest CAD devices are
8 encroaching in areas that deserve Class III
9 regulation. If the integrity in the PMA
10 process is to be preserved, a clear
11 distinction between Class II and Class III
12 chest CAD devices must be delineated by the
13 Agency.

14 As a Class III device, the Agency
15 has subjected Riverain chest CAD to
16 justifiably high degree of regulatory control.

17 As a consequence, Riverain maintains high
18 standards that affect all company activities
19 from the development, the manufacturing,
20 regulation of labeling, as well as promotion
21 and advertising.

22 Because of the importance of Class

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 III Regulation, Riverain advocates requiring a
2 reader study for initial PMA chest CAD
3 applications.

4 Further modifications could be
5 supported by standalone testing on independent
6 clinical test sets that demonstrate the
7 modified device is performing as good or
8 better than the approved device. This
9 approach is consistent with the Least
10 Burdensome provisions of the FDA Modernization
11 Act and provides reasonable assurance of
12 safety and effectiveness.

13 While we appreciate the need to
14 thoroughly evaluate devices prior to their
15 release, we urge the Agency not to regulate
16 CAD and prevent -- prevent superior products
17 from being released to the market.

18 In summary, Riverain believes the
19 following:

20 (1) The integrity of the PMA
21 process should be preserved.

22 (2) A reader study consistent with

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the intended use should be conducted for the
2 original PMA applications.

3 (3) Reasonable supporting evidence
4 for modifications to algorithms consistent
5 with the Least Burdensome provisions of the
6 FDA Modernization Act should be required.

7 Thank you for taking the time to
8 allow me to state my opinions.

9 CHAIRMAN GLASSMAN: Thank you very
10 much. Our next speaker is Akira Hasegawa of
11 Fujifilm Medical Systems USA.

12 DR. HASEGAWA: I'm Akira Hasegawa
13 from Fujifilm. The title of my talk is Risk
14 Assessment of CAD. As you know, there are
15 many different type of CAD. The different CAD
16 have different indication for use and IFU is
17 different. Risks are also different. Please,
18 note that standard -- the standard reading
19 procedure may be depending on organs and
20 modalities. Risk assessed by this comparison
21 may vary depending on organ and the
22 modalities.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 This slide shows the standard
2 reading procedure without CAD. A case comes
3 in, and the radiologist start searching
4 process -- searching suspicious area. Once he
5 identifies the suspicious area, he looks
6 closer and makes decisions. And it is
7 repeated until he goes through all area in the
8 case.

9 Here, we would like to insight that
10 radiologists are supposed to look at
11 everything, evaluate all suspicious areas, and
12 they should not overlook anything, even
13 without CAD.

14 Type 1. This is the standard
15 procedure without CAD. Again, even without
16 CAD, radiologists should not overlook
17 anything. However, radiologists are also
18 human being, and he may do some perceptual
19 oversight. To avoid such perceptual
20 oversight, after the completion of this
21 standard reading, radiologist can do optional
22 second read.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 And there is a CAD to support this
2 optional second read process. This CAD is
3 normally called Computer-Aided Detection, but
4 actually this CAD does not assist detection
5 process, it's here, but assist optional second
6 read process. Because this optional second
7 process after the standard reading procedure,
8 there is no effect to this standard reading
9 process.

10 Type 1 CAD is for optional second
11 read. It assist the users to do optional
12 second read. It helps users to reduce
13 potential oversight. The way it used as
14 recommended by the manufacturer, this type of
15 CAD does not affect users' initial reading.
16 And it does not provide any diagnostic
17 information to users.

18 Type 2. Again, this is the
19 standard reading procedure without CAD.
20 Again, radiologist should not overlook
21 anything without CAD. However, radio -- every
22 radiologist has -- sorry. Every radiologist's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 knowledge and experience are different. If
2 some radiologist does not have enough
3 knowledge and experience, he may make
4 cognitive oversight in decision making
5 process.

6 To avoid such a cognitive
7 oversight, there is a CAD to assist
8 radiologist to do decision making. Normally,
9 this is called as Computer-Aided Diagnosis,
10 because this CAD influences readers' decision
11 making process, which is a part of the
12 standard reading procedure.

13 Types of CAD is for interpretation.
14 It assist the users to make an
15 interpretation. It provides classification
16 information or diagnostic information to
17 users. It affect users' interpretation
18 process. The risk, it may affect users'
19 decision making ability negatively.

20 The Type 3. Again, this is the
21 standard reading procedure without CAD. There
22 is a CAD to assist radiologist to detect

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 suspicious areas. This type of CAD is CAD for
2 concurrent read, and it's normally called as
3 concurrent CAD. Obviously, this CAD affect
4 certain process in the standard reading
5 procedure.

6 The Type 3 CAD for concurrent read
7 it helps users to find the suspicious ROIs. It
8 affects users' initial searching process. It
9 does not provide diagnostic information to
10 users. Risk, it cause users satisfaction of
11 search. It may affect users searching
12 negatively.

13 CHAIRMAN GLASSMAN: Excuse me, Dr.
14 Hasegawa, but your time is up. Thank you.

15 DR. HASEGAWA: Oh, I thought we had
16 -- I have a 10 minutes?

17 CHAIRMAN GLASSMAN: You had five
18 minutes.

19 DR. HASEGAWA: Oh.

20 CHAIRMAN GLASSMAN: Does anyone on
21 the Panel have questions for any of the
22 speakers who have just presented? Dr. D'Orsi?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. D'ORSI: Dr. Hasegawa, are you
2 suggesting the use of CAD Type 1 and Type 2 on
3 screening exams?

4 DR. HASEGAWA: Actually, I not
5 really suggesting. It's -- my point is that
6 it's a depending IFU from manufacturers.

7 DR. D'ORSI: Thank you very much.
8 Could I do a question for Dr. Jaffe?

9 CHAIRMAN GLASSMAN: Go ahead.

10 DR. D'ORSI: Dr. Jaffe, the
11 database for the ACRIN Study, would that be
12 something that would be available to the FDA
13 for a sort of general testing set that could
14 be developed?

15 DR. JAFFE: ACRIN is a grantee of
16 ours. We have -- we are working out a complex
17 arrangement for the issue of their initiatives
18 for data-sharing and ones that we believe are
19 important for both parties, both NCI and the
20 field in general. We have done that through
21 the terms of award on the virtual colonoscopy
22 trial, which has just finished, and there will

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 be a substantial portion of that trial made
2 available.

3 The problem with the issue of the
4 DMIST trial is that it included devices as
5 well as CAD, as well as what could be CAD
6 components. When we get into the issue of
7 devices, it gets rather complex, because there
8 are multiple parties that have financial
9 interest in that. So that has been more
10 complex.

11 It is absolutely an area of great
12 interest to us, because we feel this field
13 really could develop much better if we had
14 some open databases that can act as
15 benchmarks. It puts more competition into the
16 system. So you don't actually have to file
17 with say the database and the operational
18 database, but it allows you to actually
19 develop your material, so you know how well
20 you are doing against other competing
21 procedures.

22 CHAIRMAN GLASSMAN: Thank you, Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Jaffe. Any other questions? If not, I want
2 to move on to our public speakers about
3 mammography CAD. First is Dr. Rachel Brem on
4 behalf of iCAD Medical. She is also a
5 grandparent with a new granddaughter, by the
6 way.

7 DR. BREM: Thank you very much.
8 I'm Rachel Brem. I have been involved in CAD
9 research for, approximately, a decade and have
10 been a Director with iCAD and am a Director,
11 but primarily I'm a practicing radiologist and
12 a Director of Breast Imaging at George
13 Washington University and personally interpret
14 over 10,000 mammograms a year.

15 What I would like to review today
16 is the scientific integrity, and the rigor of
17 the data presented for the original PMA
18 submission, as well as subsequent PMA
19 supplements that establish the safety and
20 efficacy of mammography CAD and that the
21 currently utilized paradigm of establishing
22 the safety and efficacy with the PMA is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 appropriate one as validated by the
2 literature.

3 The scientific data which we
4 presented with our original PMA submission
5 resulted in a 20 percent improvement in breast
6 cancer detection. And it was based on a large
7 scientific body of data, which not only showed
8 an improvement in breast cancer detection, but
9 that breast cancer with CAD was diagnosed 15
10 months earlier than without CAD.

11 The clinical utilization study of
12 nearly 4,000 patients demonstrated a
13 significantly -- a statistically insignificant
14 .5 percent increase in recall rate.

15 The data which we submitted for CAD
16 standalone performance utilized a large number
17 of breast cancers to speak to both the
18 stability of the data using ground truth as
19 well as the vast array of mammographic
20 presentations of breast cancer and 153 normals
21 to establish the false positive rate.

22 CAD -- A mammogram is a mammogram,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 regardless of how the image is obtained. In
2 the original PMA submission, which the base --
3 is the basis for establishing the safety and
4 efficacy of mammography CAD, requires
5 confirmatory studies that the standalone
6 performance of new image sources of
7 mammography demonstrate comparable
8 performance.

9 A mammogram is the same image,
10 regardless of whether it is a screen mammo --
11 film screen mammogram or digitally obtained.
12 And therefore, the confirmatory studies need
13 only show comparable performance with the
14 safety and efficacy based on the original PMA
15 submission.

16 The literature is full of studies
17 evaluating the improvement of breast cancer
18 detection, and whether it be sequential
19 studies that is on the same patient prior to
20 and after the implementation of CAD, which has
21 been reported in over 54,000 patients or
22 historical controls, that is the cancer

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 detection rate prior to and after the
2 implementation of CAD, which has been reported
3 in nearly 150,000 patients, have demonstrated
4 improvement in CAD detection from 4.5 to 19.5
5 percent. And regardless of the study design,
6 an approximate one percent increase in recall
7 rate.

8 And there is an obligatory need to
9 have an increase in recall rate in order to
10 identify more cancers. And as a practicing
11 physician, after an initial learning curve, I
12 can tell you that the false positive rate
13 becomes a minimal, if distractor at all. And
14 the increase in recall rate to patients, the
15 cost of the increase recall rate to patients
16 is not only acceptable, but based on my
17 patients, often welcome for the improved
18 detection of breast cancer.

19 CAD works, and the literature
20 supports that, whether it be in private
21 practice or in the academic setting, with an
22 experienced and novice mammographer, although

1 the impact is greater with a novice
2 mammographer, and in prospective and
3 retrospective trials.

4 Also, fatty and dense breast
5 mammography works, and this is important as
6 dense breasts not only hinder interpretation
7 of mammograms, but is a strong independent
8 risk factor for the development of breast
9 cancer.

10 And with regard to cancer size, CAD
11 works not only in large cancers, but in small
12 sub-centimeter cancers and even in cancers
13 that are 5 millimeters or less. With regard
14 to pathology, CAD works in the most common
15 invasive ductal carcinoma and the most
16 difficult to diagnose invasive lobular
17 carcinoma.

18 And make no mistake, ductal
19 carcinoma in situ is, in fact, breast cancer.

20 And although it is not yet invasive, if we
21 could diagnose all cancers in their in situ
22 phase, then we would be able to, essentially,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 cure breast cancer. Once a woman has an
2 invasive cancer, we can no longer assure her a
3 cure from breast cancer. And therefore, the
4 diagnosis of ductal carcinoma in situ is not
5 only desirable, but laudable.

6 So in summary, mammography works in
7 all clinical situations with various types of
8 readers, pathology size, and pathology and
9 size. And therefore, the currently utilized
10 paradigm of establishing the safety and
11 efficacy with the PMA submission and
12 subsequent confirmatory studies to establish
13 comparable performance of CAD is critical.

14 As a physician, as a woman, and as
15 a breast cancer survivor, the critical
16 technology of CAD is important for the optimal
17 diagnosis of breast cancer. Thank you.

18 CHAIRMAN GLASSMAN: Thank you. Our
19 next speaker is Mr. Julian Marshall from
20 Hologic/R2. Is Mr. Marshall here? Yes?

21 MR. MARSHALL: Good morning. In my
22 14 years at R2, now Hologic, I have had

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 countless opportunities to discuss CAD,
2 particularly mammography CAD, with
3 radiologists of all skill levels. The
4 conversation has become easier over time,
5 particularly in the U.S., where adoption has
6 been rapid, and an amounting body of evidence
7 has made it more and more apparent that
8 mammography CAD is clinically useful.

9 Over the last 10 years, hundreds of
10 papers have been written on mammography CAD.
11 They have varied widely in study design,
12 population, and controls. It seems that the
13 days of tiny academic studies is largely over,
14 replaced now by much larger studies with real
15 clinical basis authored by practicing
16 clinicians.

17 Yet, there is a surprising amount
18 of debate as to what truly will establish the
19 clinical efficacy of CAD. Early on, we were
20 told that CAD would not be established until
21 independent prospective studies were published
22 in peer reviewed journals. Subsequently, we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 became aware that it would be necessary for
2 those studies to demonstrate increase in
3 detection concordance with increase in workup
4 rate. And just in the last few years,
5 reduction in size, stage and age at detection.

6 By this time, a handful of
7 prospective clinical studies have been
8 published, most on the scale of 10 to 20,000
9 patients. All have demonstrated a cancer
10 detection rate increase, some concordant with
11 the increase in workup rate. Dr. Young's
12 paper is an outlier because CAD was used as an
13 adjunct to human double reading in a clinical
14 environment where any single radiologist can
15 trigger the recall.

16 In addition, please, note the dates
17 in red on this pie chart. Mammography CAD was
18 first approved in 1998, yet, it took three
19 years before the first independent prospective
20 study was published. Requiring such studies
21 prior to approval of a device will
22 dramatically delay putting these clinically

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 beneficial low-risk products into the hands of
2 clinicians.

3 When institutions implement CAD
4 into their clinical work flow, the incidence
5 of breast cancer doesn't suddenly increase.
6 Yet, the prospective studies all show an
7 increase in detection with use of CAD. It
8 wasn't until Dr. Cupples' published a paper in
9 2005, that it became clear why this was.

10 Dr. Cupples was first to point out
11 that the extra cancers detected were really
12 future years' cancers detected earlier. His
13 paper demonstrated a 164 percent increase in
14 detection rate of invasive cancers less than 1
15 centimeter and found those cancers in women
16 5.3 years younger.

17 He wrote that in multi-variable
18 analysis of invasive cancers early stage,
19 State 1, was strongly associated with
20 detection by CAD. Subsequent publications,
21 particularly by Dr. Nishikawa have discussed
22 this earlier detection effect in far more

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 detail.

2 Mammography CAD is an ancillary
3 source of information for the reading
4 radiologist, who remains the interpreter of
5 the images and the final arbiter. Should the
6 CAD device fail or results not be available,
7 the radiologist must still be able to finish
8 the reading in a timely manner.

9 As in all healthy debates,
10 mammography CAD papers and opinions exist that
11 indicate evidence to the contrary, such as the
12 heavily marketed paper by Fenton, et al. But
13 for each of those more critical studies, there
14 is also a lesser known body of published
15 commentary pointing out the significant
16 weaknesses in those studies.

17 There is clearly a mountain of
18 evidence indicating that mammography CAD has a
19 positive clinical impact, and we are pleased
20 by the recent publication of the paper by Dr.
21 Gromet, which included more than 118,000 cases
22 read with CAD and which showed statistically

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 significant performance improvement.

2 Fenton, et al, which is often used
3 as evidence that CAD does not work only
4 studied just over 31,000 cases with CAD, less
5 than 27 percent of the number used by Gromet.

6 Please, read the Gromet paper,
7 which the FDA has kindly provided to you.

8 In the last 10 years, the number of
9 women have confided that my doctor told me he
10 would not have found my breast cancer without
11 your system. I don't understand why a CAD
12 device that is already approved with the given
13 sensitivity/specificity cannot be upgraded to
14 a better algorithm with better performance
15 using the same assessment of standalone
16 performance. Thank you.

17 CHAIRMAN GLASSMAN: Thank you. I
18 was just told that I downgraded you to a
19 physician, Mr. Marshall. I apologize.

20 DR. MARSHALL: I can clarify that.

21 CHAIRMAN GLASSMAN: On the other
22 hand, our next speaker is a physician, Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Robyn Birdwell from the Society of Breast
2 Imaging.

3 DR. BIRDWELL: Thank you very much.

4 I am sponsored today by the Society of Breast
5 Imaging, but I come to you as a user of CAD
6 since the late 1990s, as a researcher in the
7 area and, in my opinion, CAD does make a
8 difference in reducing false negatives. This
9 is an opinion shared by other colleagues with
10 a modest increase in recall rate.

11 Radiologists overlook cancers. We
12 have heard many reasons for this. It happens
13 whether it is film screen or digital
14 mammograms. I turn to the DMIST study just to
15 show you that despite training, experience,
16 continuing medical education, review of
17 audits, we miss cancers.

18 It's a high volume, low prevalence
19 in the screening mammography world, four
20 cancers in every 1,000 examinations. And you
21 can see that having a tool to improve
22 sensitivity is a good idea.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 We miss cancers because they can't
2 be seen. We don't know what can help us for
3 that, at this time. Visible and detected by
4 the radiologist, but we have an error of
5 cognition. We see it, but we don't interpret
6 it correctly. This is not what CAD is
7 presently designed to help with. It's
8 visible, but simply overlooked by the
9 radiologist. It's an observational oversight
10 or an error in perception. This is where CAD
11 is helpful in our today's practice.

12 The performance and benefit has
13 been demonstrated to be equivalent, whether
14 the image source is film screen, digital
15 mammography or computed radiography. Again,
16 the present benefit, the so-called CADe, is in
17 making visible and overlooked cancers visible
18 to the radiologist by reminding us to look
19 again at this area.

20 How do we measure the performance
21 in the clinical setting? We must, as always,
22 focus on safety and efficacy. We assess

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 safety by recalls. We assess safety by biopsy
2 recommendations. Do we indeed recall more
3 patients because of CAD? Yes. If you want to
4 find more, you have to do more. Do we have
5 more biopsy recommendations? Yes. If you
6 want to find more, you need to do more.

7 Are we finding, however, additional
8 cancers that otherwise would have been
9 overlooked? Yes. We need to look at the
10 prospective sequentially read clinical
11 studies, you have heard much about most of
12 these today, looking at the conglomerate of
13 53,000 number of screens from studies done
14 both in the community as well as in the
15 academic world.

16 We see an overall increase of
17 cancers detected at 9.7 percent, an increase
18 in biopsies, and a moderate increase,
19 percentage increase in recalls. The numbers
20 here vary possibly because of the patient
21 populations, possibly because of the practice
22 numbers, maybe because of the type of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 radiologist reading these studies, but
2 overall, there is an increase.

3 This next slide, obviously, shows I
4 don't have CAD on my computer. This should
5 say historical studies in the first line. We
6 have heard about these as well. The Gur
7 Study, the Cupples, the very recent Gromet
8 Study, and then the much touted and difficult
9 to interpret Fenton Study.

10 But overall, we see again, at least
11 if you look at the low volume readers in the
12 Gur Study, an increase in cancer detection to
13 19.7 percent. The most recent study with
14 112,000 examination read pre- and,
15 approximately, the same number read after
16 initiation of CAD by Gromet, we see an
17 increase of 11 percent sensitivity. An
18 improvement in sensitivity of 11 percent.

19 And then the outlier of Fenton that
20 was a survey type study and bears its own
21 issues.

22 In conclusion, based on the USA

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 peer reviewed literature of prospective,
2 sequentially read or historically controlled
3 clinical studies, the preponderance of
4 evidence shows that currently FDA-approved CAD
5 systems and algorithm improvements are
6 efficacious, by increasing cancer detection,
7 and improving radiologist sensitivity, and
8 safe as looking at a modest increase in recall
9 rates concordant in most studies with the
10 increase in cancer detection.

11 The FDA review process for
12 mammography CAD it works as intended for women
13 and their caregivers. Thank you.

14 CHAIRMAN GLASSMAN: Thank you very
15 much. Next, Dr. Gillian Newstead for the
16 American College of Radiology.

17 DR. NEWSTEAD: Good morning. It's
18 a pleasure to be here, and I would like to
19 disclose that I do receive research as
20 indicated here, research support from Phillips
21 Medical, Bayer Health, and my spouse is a
22 stockholder in Hologic.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 I would like to move a little bit
2 from the mammography CAD into the broader
3 aspects of integration of the digital
4 enterprise for breast imaging, because it's
5 now very difficult for the healthcare provider
6 to assimilate clearly and effectively, all of
7 the different information presenting to us and
8 moving forward in breast cancer diagnosis and
9 screening. There is a vast increase in the
10 amount of information that is available to us
11 now, and radiologists are basically drowning
12 in large amounts of data.

13 Certainly, it is true with the
14 recommendations for MRI screening in high risk
15 women, in addition to mammography screening
16 for women with no known risk factors. And
17 moving into the diagnostic area of CAD as it
18 relates to breast cancer diagnosis, we do need
19 computer-assistance to view, manipulate and
20 analyze large complex imaging datasets, beyond
21 the mere interpretation and analysis of a two
22 view screening mammogram.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 In addition, if radiologists are
2 going to become very more involved in patient
3 care management, in terms of making real
4 diagnostic decisions that affect healthcare,
5 we need information regarding the patient's
6 medical records, clinical history, pathology
7 information, et cetera, which will allow us to
8 compose more detailed and management-type
9 decisions that would be important.

10 This is how I'm reading right now
11 at the University of Chicago. I have
12 multiple, in our reading room -- we have 27
13 monitors, all of them related to a different
14 kind of acquisition modality devices. We are
15 juggling mammograms, ultrasound, magnetic
16 resonance imaging, interventional studies,
17 reports, all of this partly from symptomatic
18 patients, partly from screen detected
19 mammography patients. We need to review
20 pathology to make management-type reports and
21 integrate the clinical findings.

22 Computers are going to be essential

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 for us in order to be able to do this. One
2 example I'm just going to show would be the
3 interpretation of breast MR images, where we
4 need computers to assist us in demonstrating
5 the kinetic and morphologic features of
6 lesions that we detect at screening.

7 The MRI Lexicon for breast MR
8 provides a very nice organization of reports
9 of findings that we will deliver in our
10 Lexicon reports for lesions that are detected.

11 It has been shown that radiologists are aided
12 by having computers depict the kinetic
13 characteristic lesions, rather than visually
14 estimating them during the reading process.

15 Computers can allow us to select
16 certain areas of thresholds and display the
17 enhancement characteristics of lesions that
18 are exhibiting kinetic enhancement beyond a
19 certain threshold, which is very helpful to us
20 when interpreted in context with the
21 morphology of the lesions.

22 Computers allow us to produce

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 volume metric assessments. This is important
2 in follow-up in the diagnostic setting,
3 display the enhancement curve of the lesion,
4 which you see here, and produce angiogenesis
5 maps, which depict again the dynamic
6 information appropriate to the particular
7 lesion detected.

8 So in general, computers are really
9 important in allowing us to progress from a
10 purely visual analysis of a 2-dimensional film
11 to a 3 and 4D analysis of advance technology
12 as we move forward, and how we incorporate
13 that will be important for this panel and
14 others to evaluate these devices and learn and
15 test them.

16 What does this mean for CAD
17 devices? For example, an MR, well, the
18 radiologist is using these devices as a tool
19 really to assist us in diagnosis. We are
20 using color mapping, angiogenesis mapping, and
21 many other display methods as a tool for
22 reaching an accurate diagnosis.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 This may help in depicting cancers,
2 that may be overlooked by the radiologist, and
3 morphologic assessment is going to be
4 important.

5 So in conclusion, I think computers
6 will increasingly become an active and
7 important play -- play an important role in
8 the diagnostic assessment and in the screening
9 assessment of lesions as technology proceeds
10 to give larger datasets some more advanced
11 information than we currently have. Thank you
12 very much.

13 CHAIRMAN GLASSMAN: Thank you very
14 much. Our next speaker is Dr. Robert
15 Nishikawa from the University of Chicago.

16 DR. NISHIKAWA: Thank you very
17 much. I already disclosed my financial
18 activities. So I showed this slide earlier.

19 I condensed the numbers down into this slide.
20 And the two, three important numbers are the
21 9.7 percent increase in sensitivity, 13
22 percent increase in recall rate and if you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 take the ratio of the pink to the green, you
2 get 1.3.

3 So let's suppose that's exactly
4 what the CAD does clinically, and we want to
5 be able to measure that in a reader study. So
6 those are the two points. And you can fit two
7 curves to those ROC curves, and you want to
8 compare the area under the curves.

9 Well, the area difference is only
10 .025. And if you try to design a reader study
11 to find that difference, it's a very huge
12 study. These are based on Obuchowski's work,
13 but she doesn't go down to that small of a
14 difference. So I don't know exactly what the
15 number of cases is, but I'm guessing even at
16 10 readers, it's probably a few thousand. I
17 don't think that's a practical observer study
18 to conduct.

19 So then, what is a reasonable
20 requirement for the FDA to approve or approve
21 CAD devices? Well, the goal of CAD is to
22 reduce the number of missed cancers by

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 radiologists. So then I think it's reasonable
2 to expect CAD to work something like double
3 reading with two independent readings. I'm
4 sorry. Independent double reading with two
5 radiologists, which in the literature shows,
6 approximately, 10 percent increase in
7 sensitivity and a comparable increase in
8 recall rate, basically, what you find in the
9 clinical studies to date.

10 So one may argue, that it's both
11 sensitivity and recall rate increase are just
12 shifting on the ROC curve. Well, that dotted
13 green curve is a single curve through those
14 two points. And in fact, if that was the true
15 curve, you would be shifting on the curve.
16 But that is not a valid ROC curve for people
17 who understand ROC analysis, those are the A
18 and B values and they are not consistent with
19 anything in the literature.

20 Also, if you look at the curve, it
21 says that radiologists, at the very bottom,
22 recall a lot of women before they find any

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 cancers, which is not found clinically.

2 Another argument is that reading
3 more aggressively, which is sliding on the --
4 it really moves you on the ROC curve, is
5 basically the radiologist changing their
6 threshold for recalling a patient. But what
7 CAD does, is help find cancers that were
8 missed. No matter how aggressively you read
9 the image, if you don't look at a certain part
10 of the image where the cancer is, you will
11 never find it. So a higher vigilance is what
12 CAD is giving you. And higher vigilance will
13 move you to a higher ROC curve.

14 I'm going to skip over that. So
15 what's a reasonable comparable increase in
16 recall rate? So here the yellow dot is a four
17 times increase in recall rate compared to the
18 increase in sensitivity. And you can fit an
19 ROC curve, that pink curve, which is
20 reasonable ROC curve for a radiologist, so you
21 can go up to a four times increase in recall
22 rate compared to sensitivity and still be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 considered operating on a different ROC curve.

2 So one reasonable endpoint might be
3 that the increase in recall rate to the
4 increase in sensitivity is less than 3, than
5 that's consistent with reading on two separate
6 ROC curves, which means you have an increase
7 in performance. You won't be able to measure
8 that in the study, because the differences are
9 very small.

10 So let me just comment on a couple
11 of other things that the Panel is considering.

12 Inclusion of benign cases in observer study.

13 I think this adds a complication to
14 interpretation of CAD. For CAD to be
15 effective, the radiologist needs to be
16 confident that the system can find cancers.

17 Since benign lesions that go to
18 biopsy appear malignant to radiologists, who
19 think that it's a cancer, then the computer
20 needs to point those out. But in scoring
21 those, since they are a benign lesion, scores
22 are false positive. So if the computer

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 doesn't point to it, the radiologist loses
2 confidence and doesn't use the computer
3 effectively. But if it does point to it, you
4 get penalized for pointing that out.

5 So this is a problem particularly
6 in an enriched test set. If you had the
7 normal prevalence, this wouldn't be a problem.

8 I also would like to comment on
9 evaluating CAD for different digital systems.

10 It's not reasonable to me to require PMA for
11 each -- for a CAD on each different digital
12 detector, since DMIST found no difference in
13 radiologists' performance between different
14 digital systems.

15 It's more reasonable to show CAD is
16 effective on one digital system and show
17 comparable performance on the other
18 manufacturers systems. Thank you.

19 CHAIRMAN GLASSMAN: Thank you. The
20 next speak is Dr. Mary Ellen Giger from the
21 University of Chicago. And she is a real
22 doctor. She is a Ph.D.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. GIGER: Thank you. I'm from
2 the University of Chicago. I have been
3 involved in CAD research for almost a quarter
4 of a century now. I'm representing myself. I
5 have research report from NIH and the Army,
6 and I am a shareholder and receive funding
7 from R2/Hologic.

8 Well, I'm concerned about the
9 timeliness and consistency of the translation
10 of CAD developments to clinical use. And, as
11 mentioned earlier, computers are increasingly
12 being incorporated into our lives. And the
13 progress of CAD depends on research funding,
14 careful clinical studies, and training with
15 the methods when introducing them into the
16 clinical arena.

17 Let's stop for a second and look at
18 the imaging chain. All these various
19 components ultimately lead to the radiologist
20 interpretation and decision. Many factors
21 affect the radiologist's performance levels,
22 including some such as shown here in the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 physical quality of the image.

2 Attempts to help in this area led
3 to new detector systems along the imaging
4 chain, as shown here where a screen film
5 system versus Full Field Digital Mammography
6 system was investigated.

7 Other areas along the imaging chain
8 affect the ability -- may help the ability of
9 the radiologist in terms of prior training and
10 knowledge, and help overcome some
11 interpretation conditions. These are
12 illustrated in orange here along the imaging
13 chain, which includes CAD, both in
14 quantitative image analysis and image display,
15 and how you ultimately output that CAD output
16 to the radiologist.

17 And how does the radiologist use
18 this? Well, the radiologist could read
19 without. They can perform double reading, or
20 they could read with CAD. So to get to the
21 point of my talk, and I was so worried about
22 time, I rushed through those probably too

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 quickly.

2 I'm just giving maybe a new way to
3 go in evaluating systems for the FDA. We know
4 that double reading improves detection
5 sensitivity with some increase in recall rate.

6 This occurs at the intended goal of double
7 reading is to improve sensitivity.

8 Currently, double reading by
9 radiologists is accepted, even though it has
10 never gone through any FDA-approval. From the
11 literature, the increase in sensitivity and
12 recall rate are similar for double reading and
13 single read with CAD.

14 With academia and commercial CAD
15 manufacturers, it might be useful to
16 demonstrate equivalency between single read
17 and CAD. Equivalency between single read with
18 CAD and radiologists double read, similar to
19 DMIST with the screen film mammography and
20 Full Field Digital Mammography. I am aware
21 that that study was looking for superiority,
22 however, for the FDA requirement, my

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 understanding is equivalency was accepted.

2 So the way to do this would be to
3 determine a performance standard for CAD based
4 on published data, require a CAD system to
5 meet that standard to be included in a
6 cooperative study, based on performance on
7 some test set. Perform the cooperative study,
8 similar to an ACRIN Study, but here the goal
9 will be equivalency and not superiority. And
10 allow only radiologists trained in CAD usage
11 to participate in this one big observer study.

12 After equivalency is demonstrated
13 via this multi-institutional study, as
14 improvements are made to a CAD system, it will
15 be necessary to demonstrate that the systems
16 performance meets or exceeds the earlier
17 specified performance standard.

18 One could use an independent
19 technology assessment institute, which would
20 be tasked with the performance assessment of
21 all new and improved CAD devices. The
22 institute would have a safe, sufficiently

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 large database with appropriate distributions
2 of cancer types to allow for random sampling
3 of cases for a CAD system that is undergoing
4 assessment.

5 Of course, one would have to worry
6 about the potential different types of false
7 positives from different systems. However,
8 for this large database, it should be large
9 enough so that one could randomly select a
10 subset of cases that match the distribution,
11 and this would help preserve the integrity of
12 the test set. Thank you.

13 CHAIRMAN GLASSMAN: Thank you very
14 much. Does anyone on the panel have any
15 questions for our last set of presenters? Dr.
16 D'Orsi?

17 DR. D'ORSI: Just one question for
18 Dr. Nishikawa and I think Dr. Brem. Do you
19 think that the use of CAD in secondarily
20 obtained digital images, i.e., from film, is
21 equal to the CAD operation from directly
22 attained information from direct digital? Do

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 you think that one requires a PMA? I agree
2 with you that the others do not, but what
3 about the difference between secondarily
4 derived digital data and primarily derived
5 digital data for CAD performance?

6 DR. NISHIKAWA: Are you talking
7 about digitizing the screen film images?

8 DR. D'ORSI: Correct.

9 DR. NISHIKAWA: I think that the
10 early overture is that they are probably
11 equivalent, but I am not 100 percent
12 convinced. But I don't think you necessarily
13 then have to go through a PMA to establish
14 efficacy.

15 DR. BREM: And I agree. I think
16 that the efficacy and safety was established
17 on digitized analog data and that, you know,
18 studies that show comparable performance is
19 what is needed for different image sources.

20 CHAIRMAN GLASSMAN: Dr. D'Orsi
21 again.

22 DR. D'ORSI: One question for Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Nishikawa. How did you arrive at the four
2 times increase recall to washout the increase
3 in sensitivity? Was this in order to obtain a
4 different ROC at a different level or was it
5 something else operative?

6 DR. NISHIKAWA: Right. So what I
7 tried to do was find -- so I have two points
8 in the ROC curve, and I'm moving one to the
9 right, increasing the recall rate, sensitivity
10 is the same. And I kept doing that until the
11 ROC curve became at least symmetric.

12 So if you look at the DMIST ROC
13 curves, they are all skewed to the left, which
14 is what you would expect. So in the earlier
15 green curve I showed you, was skewed to the
16 right, which is not a real curve. So once it
17 became symmetric, I said that's a possible ROC
18 curve for radiologists.

19 CHAIRMAN GLASSMAN: Question?

20 DR. BOURLAND: Yes. A question for
21 Dr. Giger. Are you suggesting an institute,
22 for instance, a radiological imaging center or

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701