



EXCOR[®] Pediatric VAD

Ventricular Assist Device

*with Stationary Driving Unit Ikus Rev. 2.1
for Pediatric Use*

Instructions for Use Rev. 3

Berlin Heart, Inc.
200 Valleywood, Suite A500
The Woodlands, TX 77380
USA

Tel. 281.882.3965
Fax 281.863.9701
Hotline: 866.249.0128

info@berlinheart.com
www.berlinheart.com

**HOTLINE**

866.249.0128

This instruction for use corresponds to the following product versions:

- *Ikus* software: from V 3.41 forward
- Laptop software: from V 3.41 V forward
- Laptop from CF30 forward

Manufacturer

**Berlin Heart GmbH
Wiesenweg 10
12247 Berlin
Germany**

Distributor

**Berlin Heart, Inc.
200 Valleywood, Suite A500
The Woodlands, TX 77380
USA**

Office Phone: 1.281.882.3965

Office Fax: 1.281.863.9701

Imprint

© 2011-05

Berlin Heart GmbH
Wiesenweg 10
12247 Berlin

All rights reserved.

This instruction for use (IFU) may not be reproduced, transmitted, transferred onto an electronic medium, stored in an electronic data processing system nor translated into another language, either in parts or as a whole, without the express prior written permission of *Berlin Heart GmbH*. Excluded are the sample copies contained in this instruction for use.

This instruction for use is intended for information purposes only. The contents of the instruction for use may be supplemented, modified or updated at any time without prior notice.

Any previous versions of the IFU become obsolete upon release of this version.

All translations of this instruction for use have been prepared and checked to the best of the translator's knowledge. The German version of this instruction for use shall prevail for all matters of interpretation within a legal context.

Licenses, trademarks and trade names

EXCOR® and  Berlin Heart® are trademarks of Berlin Heart GmbH. They are protected by the relevant legislation in Germany and other countries.

All other trademarks or brands mentioned in this instruction for use are acknowledged as being the property of the respective owners. These rights are specifically recognized and respected whenever such trademarks or brands are mentioned in this instruction for use.

Text and layout

Berlin Heart GmbH

Table of contents

1	Important safety information	13
1.1	Warnings	13
1.1.1	Storage and durability	13
1.1.2	Transport within the clinic	14
1.1.3	Device configurations.....	14
1.1.4	Procedural techniques - Ikus.....	15
1.1.5	Packaging and sterilization	17
1.1.6	Procedural techniques - pumps, cannulae, accessories	17
1.1.7	System.....	20
1.1.8	Procedures to minimize risk of thrombosis	21
1.1.9	Cleaning the components	21
1.1.10	Maintenance	21
1.1.11	Errors and corrective measures	22
1.1.12	Replacing the blood pump(s)	23
1.1.13	Driving blood pump(s) with the manual pump.....	23
1.1.14	Ambient conditions	24
1.1.15	Interaction with other procedures and therapies.....	25
1.2	Precautions	25
1.2.1	VAD placement technique.....	25
1.2.2	Ambient conditions	26
1.2.3	Caution while using the Ikus.....	26
1.2.4	Battery replacement and disposal.....	27
1.3	Obligations of the operator	27
2	General Information	29
2.1	Device description.....	29
2.2	Indications for use	29
2.3	Intended operation environment.....	29
2.4	Contraindications.....	30
2.5	Storage and durability	30
2.6	Battery replacement and disposal	30
2.7	Manufacturer's warranty.....	30
3	Description: blood pump, cannulae and accessories	33
3.1	EXCOR blood pumps	34
3.2	EXCOR cannulae	34
3.3	EXCOR accessories	35

4	Description: Ikus	37
4.1	Overview	37
4.2	Displays and operating elements	38
4.2.1	Connection panel.....	38
4.2.2	Display and operating panel.....	40
4.2.3	Power supply.....	40
4.3	Operating modes.....	41
4.3.1	Univentricular operation.....	41
4.3.2	Biventricular operation	41
4.4	Laptop computer with monitor program	42
4.5	Safety.....	42
4.5.1	Redundant design of pneumatic systems in univentricular operation.....	42
4.5.2	Redundant design of pneumatic systems in biventricular operation	42
4.5.3	Control computer with redundancy design	43
4.5.4	Battery operation.....	43
4.5.5	Manual pump.....	43
4.5.6	Password-protected user profiles (access passwords)	44
5	Instructions for use: Ikus	45
5.1	Start menu	46
5.1.1	Selecting an option in the start menu.....	46
5.1.2	Configuring user passwords.....	47
5.1.3	Saving data on USB stick.....	48
5.1.4	Changing date or time	48
5.2	Basic instructions for monitor program.....	49
5.2.1	Starting the monitor program	49
5.2.2	Shutting down the monitor program	49
5.2.3	Logging in and out of the monitor program.....	49
5.2.4	Logging in	49
5.2.5	Standard view – monitor program	50
5.2.6	Selecting monitor program options.....	51
5.2.7	Adjusting the parameter values	52
5.2.8	Browsing in the message window.....	53
5.3	Stopping the blood pump(s) and switching off the Ikus	53
5.3.1	Drive pause: stopping the Ikus temporarily.....	53
5.3.2	Pause left/ Pause right: stopping an individual blood pump	54
5.3.3	Drive OFF: switching the Ikus off.....	54
5.4	Switching between mains and battery operation	55
5.5	Changing over from univentricular to biventricular operation	58
5.5.1	Routine start-test when not in operation	58
5.6	Moving the Ikus.....	59
5.7	Transportation and packaging	60
5.7.1	Unloading the Ikus from the shipping crate.....	60
5.7.2	Loading the Ikus into the shipping crate	61

6	Commissioning the Ikus and setting parameters	65
6.1	Preparatory steps outside of the operating room	65
6.1.1	Connecting the tank unit.....	66
6.1.2	Switching on the Ikus	66
6.1.3	Starting the monitor program.....	67
6.1.4	Setting the test parameters	67
6.2	Intraoperative drive management	68
6.2.1	Disconnecting the tank unit from the Ikus.....	68
6.2.2	Selecting the operating mode (view Select operating mode)	68
6.2.3	Select the pump size	69
6.2.4	Select the cannula size	69
6.2.5	Display pump and cannula sizes	70
6.2.6	Setting the start-up parameters	71
6.2.7	Connecting the blood pump(s) to the Ikus.....	72
6.2.8	De-airing the blood pumps in single-step mode.....	72
6.2.9	Starting the blood pump (changing to standard view)	73
6.2.10	Checking the parameters when the pump is started and adjusting them	73
6.2.11	Switching from CPB support to VAD support.....	75
6.2.12	Possible complications.....	76
6.3	Postoperative drive management.....	77
6.3.1	After transfer to the ward	77
6.3.2	Follow-up treatment	77
7	Implantation: Preparations in the operating room	79
7.1	Preparing the components and materials required	79
7.2	Checking and adjusting the settings of the cable tie gun.....	79
7.3	Unpacking the sterile components.....	80
7.4	Moving the membrane to the end-of-diastole position.....	81
7.5	De-airing the blood pump.....	81
7.5.1	Inserting the de-airing needle	81
7.5.2	Rinsing and filling the blood pump	82
8	Implantation - surgical procedure	83
8.1	Cannula exit sites	84
8.2	Use of the cannula tunneling tip	85
8.3	Cannulae and connector set	85
8.4	Access	86
8.5	LV apex cannula	87
8.5.1	Anastomosis of inflow cannula with LV apex.....	87
8.5.2	Creating a transcutaneous tunnel for the LV apex cannula	88
8.6	Atrial cannula(e).....	88
8.6.1	Creating a transcutaneous tunnel for atrial cannula(e)	89
8.6.2	Anastomosis of atrial cannulae.....	89

8.7	Arterial cannula(e)	90
8.7.1	Creating a transcutaneous tunnel for arterial cannula	90
8.7.2	Anastomosis of the arterial cannula	91
8.8	Shortening the cannulae if necessary	92
8.9	Connecting the blood pumps to the cannulae	93
8.10	Removing the de-airing needle	94
8.11	Securing the connections.....	94
9	Implantation - anesthesia	97
10	Wound care and treatment	99
10.1	Removing the old dressings.....	100
10.2	Cleaning the blood pump.....	100
10.3	Cleaning of the wound	101
10.4	The new dressing.....	101
10.4.1	Preparing a new dressing	101
10.4.2	Applying a new dressing.....	102
10.5	Regular checks of blood pump(s) and cannulae.....	103
10.5.1	Visual inspection: pump filling and ejection	104
10.5.2	Visual inspection: deposits.....	104
10.5.3	Checks using the monitor program	105
10.5.4	Replacing the blood pump due to growth of the patient	105
11	Anticoagulation therapy	107
11.1	Before Implantation of the EXCOR	107
11.1.1	General considerations.....	107
11.1.2	Pre implantation	107
11.2	During Implantation - Cardiopulmonary Bypass	107
11.2.1	Cardiopulmonary Bypass (CPB)	107
11.2.2	Post CPB.....	107
11.3	Postoperative anticoagulation therapy	107
11.3.1	General Considerations	107
11.3.2	Starting anticoagulation therapy	107
11.3.3	Unfractionated heparin therapy (i.v.) Patient < 12 months	108
11.3.4	Unfractionated heparin therapy (i.v.) Patient ≥ 12 months.....	108
11.3.5	Thrombelastography (TEG®)	109
11.4	Low Molecular Weight Heparin	109
11.5	Oral Anticoagulation Therapy (only for patients ≥ 12 months of age who are taking a full oral diet)	109
11.6	Monitoring of Blood Count and Anticoagulation Status	109
11.7	Postoperative platelet inhibition therapy.....	110
11.7.1	Start of therapy.....	110
11.8	Adjunctive Medication	110

11.9	Anticoagulation Therapy	111
11.9.1	Therapeutic Heparin administration and adjustment	111
11.9.2	Oral Anticoagulation Therapy	112
12	Weaning and Explantation for BTR and BTT	113
12.1	Weaning Procedure	113
12.1.1	Introduction	113
12.1.2	Indication	113
12.1.3	Eligibility Criteria	113
12.1.4	Weaning Protocol	113
12.1.5	10 ml pump	114
12.1.6	25/ 30 ml pump	118
12.1.7	50/ 60 ml pump	122
12.1.8	Explantation Criteria	125
12.2	Explantation for BTR	126
12.2.1	Explantation with univentricular support	126
12.2.2	Explantation after biventricular support	126
12.3	Explantation for BTT	127
13	Error Messages and corrective measures	129
13.1	Pressure error / time error in system 1 (or in system 2 or 3)	130
13.2	Throttle valve error in system 1 (or system 2 or 3)	131
13.3	Please connect driving tube	131
13.3.1	Replacing a driving tube	131
13.4	Please check left / right pump and driving tube	132
13.5	Backup operation left/right	133
13.6	Error messages in emergency operating mode	133
13.6.1	UVAD, emergency operation!. Contact service immediately!	133
13.6.2	Emergency operation system 1 (or system 2 or 3). Contact Service now!	133
13.7	System 1 (or system 2) is defective!	134
13.8	System 3 (backup) is defective!	134
13.9	Alarm circuit fault: buzzer remains off (or on)	135
13.10	Backup computer faulty! Contact customer service.	135
13.11	Backup computer	135
13.11.1	... reports discrepancy in left/right pump output measurements!	135
13.11.2	... reports faulty measurements on the left (or right)	135
13.11.3	... reports an error in output measurement on the left (or right) pump	136
13.11.4	... reports faulty test	136
13.12	Measurement discrepancy in main computer (backup computer)	137
13.13	Parameter set update failure	137
13.14	Temperature sensors: <<8-digit binary code>>	137
13.15	Fault: <<16-digit binary code>> (<<type of fault>>)	138
13.16	Batteries discharged; battery operation not possible	138

13.17	Insufficient battery charge. Only limited battery operation	138
13.18	Error messages - Circuit breaker and internal battery fuse.....	139
13.19	Electronic malfunction. Contact customer service!	139
13.20	Acoustic alarm is not properly recognized	140
13.21	Error: no data/no reaction from the control computer.....	140
13.22	Left/right flow sensor defective. Notify Service!.....	141
13.23	Problem: <<Text>>	141
13.23.1	Self-test is not completed by passive computer!.....	141
13.24	Error messages during the start-up test.....	142
13.24.1	Battery test skipped (Battery problem!)	142
13.24.2	Additional messages during the start-up test.....	142
13.25	Discrepancy in pressure measurement: system 1 (or system 2 or 3)	143
13.26	Communication with Laptop failed	144
14	Troubleshooting and correcting faults	145
14.1	Replacing the blood pump(s)	146
14.1.1	Preparing a replacement blood pump	146
14.1.2	Replacing the right blood pump (RVAD/ BVAD)	147
14.1.3	Replacing the left blood pump (LVAD/ BVAD)	148
14.2	Restarting Ikus.....	149
14.3	Emergency pulse mode.....	150
14.3.1	Emergency pulse mode - switching the Ikus off.....	151
14.3.2	Ikus start-test following emergency pulse mode.....	151
14.4	Connecting the patient to a replacement Ikus	152
14.5	Driving blood pump(s) with the manual pump.....	154
14.6	Mains failure or breakdown of both control computers	155
14.7	Reading out the LOG files	156
14.8	Circuit breaker and battery fuse	157
15	Appendix	159
15.1	Overview: Product range and possible combinations	159
15.1.1	Blood pumps.....	159
15.1.2	Overview: Relationship: body weight – pump size.....	159
15.1.3	LV apex cannulae	160
15.1.4	Atrial cannulae	160
15.1.5	Arterial cannulae	160
15.1.6	Overview: Which cannulae should be used for which pump?	161
15.1.7	System accessories	161
15.1.8	Driving unit	161
15.1.9	Special components	161
15.1.10	Pump-cannula combinations	162
15.1.11	Blood pump combinations in biventricular mode	162
15.1.12	Relative systolic duration.....	162

15.2	Technical specifications	163
15.3	Symbols and tags	165
15.3.1	Ikus: identification plate	165
15.3.2	Ikus: connection panel	165
15.3.3	Symbol on numlock status LED.....	165
15.4	Sample copy: Ikus Incoming Checklist.....	166
15.5	Sample copy: EXCOR Implantation log.....	167
15.6	Sample copy: EXCOR pump log	170
15.6.1	Explanations on the pump log	170
15.6.2	Sample copy: EXCOR pump log	171
15.7	Sample copy: EXCOR parameter log.....	172
16	EMC tables	173
16.1	Essential Performance	173
16.2	Electromagnetic emissions.....	173
16.3	Electromagnetic immunity - part 1	174
16.4	Electromagnetic immunity - part 2.....	175
16.5	Recommended separation distances between portable and mobile RF communications equipment and the Ikus	176
	List of figures	177
	List of tables	179
	Index	181

Dear readers,

This Instructions for Use is intended for all medical personnel involved in caring for a patient who is being supported by an EXCOR® *Pediatric VAD* (referred to as EXCOR in this instruction for use).

The instruction for use provides information on the structure, principle of operation and application of the EXCOR in conjunction with the *Stationary Driving Unit Ikus* (referred to as *Ikus* in this instruction for use). To ensure patient safety and comfort, please read the instruction for use carefully.

Always make sure that only professional medical personnel who have been specifically trained in the use of the product are permitted to work with EXCOR.

Note: The recommendations in this manual are based on *Berlin Hearts* experience with the EXCOR. The decisions related to implantation, the components to be used, and patient care remain with the patients physicians.

The following pictograms and symbols are used in this instruction for use:



Indicates a hazardous situation which, if not avoided, **will** result in death or serious injury to the patient.



Indicates a hazardous situation which, if not avoided, **could** result in death or serious injury to the patient.



Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury to the patient and/ or damage to the device.



Notes are practices not related to personal injury. Possible damage to the device.



This symbol identifies measures and procedures which have proved useful and successful in conjunction with EXCOR and which we therefore recommend.



866.249.0128

This is the telephone number of the emergency hotline. The hotline desk is in operation 24 hours a day. This number is intended for use by medical personnel and should be used in cases of emergency only.



1. Individual steps of the instructions are numbered in sequential order.

Term definitions

Product life

How long the product can be used. No maintenance or repairs are performed after this date. With unsterile products, the product life starts on the day of (initial) shipment; with sterile products it starts on the day of implantation. All sterile products are for single use only.

Expiration date

How long the unused sterile product reliably maintains sterile. The device should not be used if the expiration date is past.

Maintenance interval

Interval in which the product needs to be serviced.

Warranty

In the case of justified claims, the company must choose to either repair or exchange the defective goods within a reasonable period of time. The buyer will be entitled to cancel or to reduce the order, in accordance with legal regulations, only when the defective goods should repeatedly fail or be deemed unacceptable. The buyer may not rectify the defects under any circumstances. The buyer is entitled to claims for defects only in accordance with item 5.5 of the General Business Terms and Conditions. (for additional information see section 5 of the Terms and Conditions)

Warranty period

All warranty claims expire after 12 months, calculated from the time that risk is transferred (for additional information see section 5 of the General Terms and Conditions).

The precise application on the individual components of the EXCOR will be described in section 1.3: Obligations of the operator, page 27.

Definition of the used font formates

Description	Meaning
bold, blue	software texts (messages and menus) except in headings and lists
<i>cursive</i>	proper names (except in headings and in registered trademarks)
"text"	quotation
<key>	key on the laptop keyboard
<<filler text>>	e.g. if texts in error messages are various
[dimension unit]	dimension units in tables; e.g. [mmHg]

1 Important safety information

1.1 Warnings



Before using EXCOR, read the Instructions for Use carefully.

Only qualified medical personnel trained specifically in the use of the system are permitted to work with EXCOR. Training courses can be arranged with *Berlin Heart, Inc.* Use by untrained personnel can pose a risk to the patient and the EXCOR.

Before starting the *Ikus*, make sure that a replacement *Ikus* is available in the hospital. If a replacement *Ikus* is not available, there is a risk that the patient cannot be cared for in the event of device malfunction.

The general rule is:

1 replacement *Ikus* if 1 or 2 systems are in use,

2 replacement *Ikus* if 3 or 4 systems are in use,

3 replacement *Ikus* if 5 or 6 systems are in use.

If more than 6 systems are in use the number of replacement *Ikus* has to be 1/2 of the active systems.

On the system EXCOR only use components of this system. Never use other components than those delivered by *Berlin Heart GmbH/ Berlin Heart Inc.* Otherwise the warranty is no longer valid.

The *Ikus* should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the *Ikus* should be observed to verify normal operation in the configuration in which it will be used.

The system *EXCOR Pediatric* and its components are permitted to be used only by prescription of the attending physician.

Unintended use can pose a risk to the patient and the EXCOR.

Do not use the EXCOR if there is any visible damage of the *Ikus* or any of its components.

If there is any malfunction of the *Ikus* while the driving unit is connected to the patient, the *Ikus* must immediately be replaced.

1.1.1 Storage and durability



The expiration date of each *EXCOR* product is found on the product labels located on both the outer and inner packaging. The pumps, cannulae and accessories must not be used after the expiration date and even not be re-sterilized. Otherwise there is a risk of patient infection.

An *EXCOR* blood pump may not be used on a patient for more than 1 year. After this it shall be replaced with new products.

1.1.2 Transport within the clinic

WARNING

To move the *Ikus* unit: push only, using the handle provided for this purpose. Avoid any sudden jerky motion. When passing over smaller obstructions, exercise extreme caution, pulling the *Ikus* unit backwards (i.e. handle first) across the obstruction if necessary.

To lift the *Ikus*: use only the lifting bars at the lower edge at each side of the unit to hold and lift it. Never attempt to lift the *Ikus* by its handle. The *Ikus* must always be lifted by at least 2 people, preferably 4.

Rolling the *Ikus* over sloping surfaces: ensure that the person pushing it is strong enough to push the *Ikus* in a controlled manner. The slope of the surface may not be steeper than 10° (exception: packing/unpacking of the *Ikus* into/ from the transport crate). Otherwise there is a risk of injury to the transporting persons or of damaging the *Ikus* !

If it is necessary to transport the patient within the clinic ensure that he is accompanied by a person trained to use the manual pump.




Fig. 1-1 Handle and lifting bar

1.1.3 Device configurations

WARNING

EXCOR was not designed to be used in combination with other systems, nor do any of the currently granted approvals allow for this. Use by untrained personnel poses a risk to the patient and to the EXCOR.

In univentricular operation (also for an RVAD): Always connect the driving tube of the blood pump to the red marked connector.

The units may only be operated with the disposable products and accessories specified in this document. Also see section 15.1.2: Overview: Relationship: body weight – pump size, page 159. Otherwise there is a risk of functional limitation and/or damage to the *Ikus*. Failure to observe this stipulation will invalidate all warranty agreements by *Berlin Heart Inc.*

The connection between the connector *External alarm* (Nurse call) and the internal alarm system of the clinic is not failsafe. The use of this feature does not release the user from supervising the *Ikus* and the displayed messages and alarms.

1.1.4 Procedural techniques - *Ikus*



Follow the numbered instructions exactly in their sequential order. Otherwise there is a risk of functional limitation or *Ikus* malfunction.

Protect the *Ikus* unit against dirt and contamination. Prevent foreign objects from falling or working their way into the connectors and ventilation slits. Keep all drive tube connectors covered at all times when not in use. Otherwise there is a risk of functional limitation and/or damage to the *Ikus*.

Place the *Ikus* driving unit on a firm and even surface.

Never place other objects on top of an *Ikus* driving unit.

When switching on the *Ikus*, always connect it to the mains power supply. This is the only way to ensure that the start test (see section 6.1.3: Starting the monitor program, page 67) is completed and that possible malfunctions are detected.

Only connect the mains cable of the *Ikus* driving unit to grounded mains power outlets. The supply voltage has to conform to the voltage requirements indicated on the *Ikus* identification plate. Only connect the mains cable of the *Ikus* driving unit to suitable power outlets complying with the electrical safety regulations of the country in which it is being used. Otherwise there is a risk of electrical shock or damage to the EXCOR system!

If there are several electrical devices in the vicinity of the patient, then these are to be connected to a central grounding point. The connector *Potential equalization* is to be used for this on the *Ikus*. Otherwise there is a risk of electric shock.

Switch on the *Ikus* 2 hours before use in order to sufficiently charge the batteries and to detect possible device errors during the start test. Important: During this time, always connect both tank units to the *Ikus* (see section 6.1.1: Connecting the tank unit, page 66)! Otherwise there is a risk that error messages are falsely generated.

After switching on the *Ikus*, pull the key out of the main switch (key switch) and store it in a safe place.

The *Ikus* power switch (toggle switch) must be switched on during the first startup and remain set to [I] position. Its position should never be changed. Otherwise there is a risk that the batteries will not be fully charged after battery operation.

The *Ikus* is designed for stationary operation and to be run on mains power (referred to as the mains in this instruction for use). Do not run it on battery operation unless this is absolutely necessary (e. g. when moving the patient within the clinic or during a mains failure).

Always take the battery capacity limits into consideration when using the system. When the system is working in battery operation, the patient must be accompanied by a person trained to use the manual pump. Thus the patient shall be guaranteed care in an emergency.

The air vents must not be covered or obstructed during operation. Otherwise the *Ikus* will not receive sufficient ventilation, which may lead to overheating. This may lead to malfunctions and a malfunction in device operation.

Only disconnect the *Ikus* from the mains if the charge level indicator shows that the batteries are completely charged (all yellow LEDs are lit).

Important: In order to prevent rapid and premature ageing of the batteries, the *Ikus* must always be run on the mains for at least 6 hours after using the battery operation. Only after this may the *Ikus* be run on battery operation again.

If the LEDs of the charge level indicator are blinking or the message **Batteries discharged - use power supply!** appears, immediately switch to mains operation! If the batteries are completely discharged (red LED is lit) and the drive continues to run on battery operation, there is the risk of a total malfunction of the *Ikus* and damage of the batteries.

Whenever the *Ikus* is running in battery operation, the patient must be accompanied by a person trained to use the manual pump. Thus the patient shall be guaranteed care in an emergency.

To prevent the batteries from aging quickly, every period of battery operation should be followed by at least 6 h of mains operation.

When the battery charge is low, the acoustic signal sounds every minute. The *Ikus* must be connected to the mains operation immediately.

If the batteries are discharged completely (red LED lights up), there is a danger of a total malfunction of the *Ikus* if battery operation is continued, and that the batteries will be damaged. If this happens, it can not be guaranteed that the *Ikus* will restart after connecting it to the mains.

Messages are only displayed when the monitor program is running. When the monitor program is shut down, the only indications that there is an error message are an acoustic signal and the fact that the indicator lamp on the *Ikus* handle lights up. There is no way of finding out what type of error message has been displayed.

When reading out the log files: Make sure that you always have the USB stick inserted and that there is sufficient capacity on the stick. Otherwise the LOG files might get lost as they are deleted from the hard disk as soon as they have been transferred onto a USB stick.

1.1.5 Packaging and sterilization



EXCOR blood pumps and cannulae are intended for single-use only. Otherwise there is a risk of infection.

The sterile components are sterilized using ETO and are packed in a double-layer sterile package. Check that the various layers of the sterile packaging are not damaged in any way before they are opened. Do not use the components if either of the sterile packages are damaged. The same applies to sterile components which have exceeded the expiration date as printed on the label. Otherwise there is a risk that the product is no longer sterile.

An aluminum-coated external packaging protects the *Carmeda*[®] *BioActive Surface (CBAS)* of the blood pump and its sterile packaging against fluctuations in relative humidity. Do not use blood pumps with damaged external packaging. Otherwise there is a risk that the CBAS coating may be compromised.

The following items are delivered in sterile condition: blood pumps, cannulae, driving tubes, de-airing set, de-airing hammer, tube connecting set., membrane set .

The external packaging and the outer surface of the outer sterile packaging are not sterile. These 2 packaging layers must be removed before the inner sterile packaging containing the product is handed over to the sterile field. Otherwise there is a risk that the sterile field will be contaminated.

1.1.6 Procedural techniques - pumps, cannulae, accessories



The preparation and use of blood pumps should only be performed by trained personnel. Surgical, nursing and perfusion personnel without experience in the use of EXCOR must complete the EXCOR Training Course which provides theoretical introduction and hands-on practical exercises in the operation of this system. The training program is organized and offered by *Berlin Heart, Inc.*

Only use sterile components which have been delivered in undamaged sterile condition (sterile packaging intact, expiration date not expired).

Only use blood pumps which have an undamaged aluminum-coated outer packaging.

The long-term storage conditions for all sterile products must be observed: temperature +15°C to 25°C, relative humidity: 35 % to 50 %. Store in a dry place! Otherwise there is a risk that the product is no longer sterile.

In order to prevent infection, use strict aseptic techniques during implantation and exercise extreme caution throughout the period of EXCOR cardiac support. Danger of infection!

The distal end of the cannulae can be trimmed. At least 5 cm (2 inches) of material without polyester velour covering should remain to allow visual inspection of the cannula/ titanium-connector junction. Otherwise there is a risk that possible deposits if formed, cannot be visualized.

Ensure proper placement of the cannulae, especially with respect to orientation of the LV apex cannula, to prevent suction of the myocardial wall.

Prior to initial operation of the blood pump(s) minimal initial start parameters have to be set on the laptop to ensure smooth transition from CPB to VAD support.
When connecting the blood pump(s) to the cannulae always observe the arrows on the inflow and outflow stubs. They show the blood flow direction. There is a risk of injury to the patient and severe pump malfunction if the titanium connectors on the end of the inflow and outflow stubs are not connected to the appropriate cannulae.

Do not touch or manipulate the blood pumps and cannulae with pointed or sharp-edged objects (surgical instruments, wire brushes, etc.). Otherwise there is a risk of blood pump and cannula leakage.

Creating a transcutaneous tunnel for the LV apex cannula: Always use cannula tunnelling tip, never use a sharp surgical instrument directly on the cannula.

If an *EXCOR connecting set* is required for implantation and the length of the tube part needs to be reduced, the tube part should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 1-1 Connector set: minimum length of connector tube



Follow exactly the instructions for using the de-airing set. Otherwise there is a risk of membrane damage.

Ensure that cannulae, blood pump(s) and driving tubes are not subject to external forces, like compression, traction or torsion forces, and are free of knots or sharp bends. Prevent the cannulae and connectors from being exposed to tensile forces. Otherwise there is a risk of obstruction of the air and blood flow.

When positioning the driving tubes follow hospital policies to mitigate the risk of adverse tubing and line incidents by routing the driving tubes in a clear pattern toward the feet and to the side.

Do not initiate cardiac support with the EXCOR blood pumps until the blood pumps have been completely de-aired. After connecting the cannulae, ensure removal of all air that is still in the atria or ventricle by performing single steps (**step left**, **step right**) with subsequent removal of the bubbles inside the pump via the de-airing needle. Otherwise there is a risk of embolism.

When removing the de-airing needle, never pull on the de-airing tube, but rather only on the de-airing needle.

Once the de-airing needle has been removed it cannot be re-inserted.

Rates < 60 bpm are intended to be used only for implantation and explantation. Never use the *Ikus* with a rate < 60 bpm without constant supervision.,

Under circumstances, the messages **Please check left pump and driving tube!** or **Please check right pump and driving tube!** are not generated with the 10 ml EXCOR blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

Secure each connection between blood pump and cannula with at least one cable tie as soon as the proper function of the EXCOR is established (see section 8.11: Securing the connections, page 94). Otherwise there is a risk of loose connections and inadequate blood supply to the patient..

At least every 4 hours, visually check that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, institute the appropriate corrective action.

Do not kink the drivelines. Otherwise there might not be sufficient pump output.

In no case should the cannulae either be kinked directly at the connector to the blood pump or at the transition area between velour and silicone.

Do not kink the cannulae needlessly. Otherwise there might not be sufficient pump output. Moreover, cannulae might be damaged.

Wound care and treatment: Before cleaning the wound (see 10.3: Cleaning of the wound, page 101), put on sterile disposable gloves, cap and mask.

Weaning: If the patient does not meet the eligibility criteria at any time during the weaning process: Resume pumping at rate prior to any weaning (initial rate, IR).

1.1.7 System



If a non-matching pump-cannula-combination (see section 15.1.10: Pump-cannula combinations, page 162) was chosen, use only the connector sets provided with the system in order to minimize the risk of clots at the junctions. Be aware of increased risk of thrombosis and hemolysis.

The cannula diameter may be adapted only once (either by using a staged cannula or a connector set.) Multiple staging could result in limited pump performance and compromised hemodynamics.

Do not install any additional software on the laptop. Otherwise there is a risk of damage to the original *Ikus* software. Risk of total malfunction of the *Ikus*!

Make sure that the <NumLk> key of the laptop is deactivated. The status LED on the laptop marked with a lock and/or a number (e.g. 1) should not be lit. Otherwise there is a risk of incorrect inputs.

Only use USB sticks included within the delivery of the EXCOR to store data. Do not use any other USB sticks with the laptop. Risk that a (wrong) USB stick is not recognized. If a (wrong) USB stick is recognized, then it may not be possible to save the data.

Never connect other USB devices to the USB port of the laptop than the delivered USB sticks. Otherwise there is a risk that the batteries of the *Ikus* will be discharged too fast.

Never connect wireless technology to the USB port of the laptop. Otherwise there is a risk of uncontrolled electromagnetic radiation which might interfere with other devices. The *Ikus* could also become more susceptible to emissions from other devices.

Prior to connecting and disconnecting the USB stick to the slot, the laptop must be switched off. Otherwise there is a risk that the USB stick will not be recognized. When removing the stick the stored data can be lost.

If the *Ikus* is operating in emergency pulse mode, immediately visually check whether the blood pump(s) is (are) filling and ejecting completely. If one pump is not filling and/or ejecting completely, the patient must be supported immediately using the manual pump (see section 14.5: Driving blood pump(s) with the manual pump, page 162). Otherwise there is a risk that the patient will not be supported sufficiently.

If the *Ikus* is operated by the backup system provide the patient immediately with a replacement *Ikus*.

Do not disconnect the *Ikus* from the mains power supply if the circuit breaker is triggered. Otherwise there is a risk that the driving unit immediately stops operating (see section 14.8: Circuit breaker and battery fuse, page 165).

Do not use water or fluids to cool the *Ikus*! Otherwise there is a risk of a short circuit or a malfunction of the device.

1.1.8 Procedures to minimize risk of thrombosis



Ensure complete filling/ejection of the pump.

When using staged cannulae or a connecting set, the pumping rate may not be greater than the respective value found in table 14-9, pg. 104, as the pump will not eject its full volume at higher rates.

At least every 4 hours, visually check of blood pump(s) for deposit formation.

1.1.9 Cleaning the components



Cleaning the pump and the drive line: Do not use any acetone or petroleum based products near the pump or drivelines. We recommend using only water or alcohol to clean the pump and the drive line.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the blood pump or drivelines as they may alter the surface of the product.

Cleaning the cannulae and transcutaneous exit site: Do not use any acetone or petroleum based products near the the cannulae and the transcutaneous exit site.

We recommend using chlorhexidine to clean the cannulae and transcutaneous exit site.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the cannulae and transcutaneous exit site as they may alter the surface of the product.

1.1.10 Maintenance



If the *Ikus* is not in operation, it requires maintenance every 6 months. If it is in operation, it requires maintenance after every 2000 hours of use. (In case of continuous operation, approx. 3 months)

If the ambient temperature is continuously above +30°C during operation, the maintenance interval or life of the batteries is reduced.

The *Ikus* shall only be serviced by *Berlin Heart GmbH/ Berlin Heart, Inc.* or those authorized by *Berlin Heart GmbH*. For this reason, this document does not contain any circuit or wiring diagrams.

Only replacement parts approved by the manufacturer may be used for repairs and servicing. Otherwise there is a risk of functional limitation or permanent damage of the *Ikus*.

1.1.11 Errors and corrective measures



Any time an error message has occurred, visually check that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles, then address the error message with the appropriate corrective action.

Check all information and messages in the message window of the monitor program at least every 4 hours. Take the necessary measures and (if required) notify the service department of *Berlin Heart, Inc.* The message window only shows a limited number of messages. Otherwise there is a risk that older information and messages can no longer be read and therefore corrective actions for older messages may no longer be possible.

If a message with the content **... Contact (customer) service (now)!** is displayed in the message field on the laptop, replace the *Ikus* immediately (see chapter 13: Error Messages and corrective measures, page 129).

If the emergency pulse mode is activated while the backup system is already active, the *Ikus* is no longer able to drive both pumps. In this case, the patient must immediately be supported using the manual pump (see section 14.5: Driving blood pump(s) with the manual pump, page 162). Otherwise there might not be sufficient pump output.

In order for a driving tube to be replaced, the pump must be stopped for a short time. If the left driving tube is being replaced in a driving unit providing biventricular support, the right pump must also be stopped while the driving tube is being replaced in order to avoid overloading of the pulmonary circulation (danger of pulmonary edema).

If the left pump is being replaced in a VAD providing biventricular support, the right pump must also be stopped while the pump is being replaced in order to avoid overloading the pulmonary circulation (danger of pulmonary edema).

Message **Left/right flow sensor defective Notify Service!**: Do not operate the *Ikus* without supervision! Otherwise an insufficient support of the patient might not be detected.

If the *Ikus* is operating in emergency pulse mode, the user must immediately visually check the blood pump(s) to determine whether the pump(s) are filling and ejecting completely. If one pump is not filling and/or ejecting completely the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152 and section 14.5: Driving blood pump(s) with the manual pump, page 154 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

If the emergency pulse mode is activated while the backup system is already active, the *Ikus* is no longer able to drive both pumps. In this case the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152 and section 14.5: Driving blood pump(s) with the manual pump, page 154 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

1.1.12 Replacing the blood pump(s)



When replacing a blood pump, follow the instruction given here. Otherwise the duration of the pump stop will be prolonged and the patient might suffer from inadequate support.

The blood pump may only be replaced under sterile conditions!

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs! These show the direction of the blood flow.

The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.

If the left pump is being replaced in a VAD providing biventricular support, the right pump must also be stopped while the pump is being replaced in order to avoid overloading the pulmonary circulation (danger of pulmonary edema).

1.1.13 Driving blood pump(s) with the manual pump



The use of the manual pump is only permitted for medical personnel trained in the use of it.

Pay attention to the colored markings on the driving tubes and on the connectors of the manual pump. Otherwise, there is a risk of lung edema.

Always keep manual pump attached to the *Ikus*. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

Call one or more persons to assist. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

The driving tubes and cannulae should be arranged in a bend-free position. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

When operating the manual pump with 1 hand, do not block the valves with your feet (see c in figure 14-3 and in fig. 14-4, page 144).

1.1.14 Ambient conditions



Protect the *Ikus* from exposure to moisture and wetness. Never store or operate the *Ikus* in a damp environment (e.g. bathroom, etc.). Otherwise there is a risk of functional limitation and/or *Ikus* malfunction.

In terms of electromagnetic compatibility (EMC) the *Ikus* is subject to special precautions! Avoid exposure to strong electromagnetic radiation (as generated by mobile/cell phones and cordless phones when switched on, electromagnetic security systems etc.), see chapter 16: EMC tables, page 173. Otherwise there is a risk of electromagnetic disturbances and fault-free functioning of the *Ikus* cannot be guaranteed.

When using a cell phone in the immediate environment of an *Ikus* in operation please make sure to keep a distance of at least 0.77 m. For further information please refer to chapter 16: EMC tables, page 173.

When using an RFID device in the immediate environment of an *Ikus* in operation please make sure to keep a distance of at least 1 m. For further information please refer to chapter 16: EMC tables, page 173.

Protect the *Ikus* against temperatures below +10°C and above +30°C; this includes extreme temperature changes and overheating (e.g. direct sunlight or from heaters). Otherwise there is a risk of functional limitation and/or *Ikus* malfunction.

If an ambient temperature of +30°C is continuously exceeded during operation, the lifetime of the batteries is reduced. Therefore, a person trained to use the manual pump should always be present in this case. This should ensure patient care in case of emergency.

Use the *Ikus* as far away as possible from environments containing flammable gases and use extreme caution. Otherwise there is a risk of explosion or gas ignition. The *Ikus* would be severely limited in function or malfunction altogether as a result of this damage.

Also see section 15.2: Technical specifications, page 163.

1.1.15 Interaction with other procedures and therapies



The following procedure is not possible:

- Magnetic resonance imaging
-

If EXCOR is used in interaction with other procedures and therapies, observe the movement of the membrane to determine whether the blood pump is filling and ejecting completely. If a pump is not filling and/ or ejecting completely, stop the interacting procedure or therapy and institute the appropriate corrective action.

In terms of electromagnetic compatibility (EMC) the *Ikus* is subject to special precautions! When exposing *Ikus* to the procedures and therapies listed below please observe EMC regulations given in chapter 16, page 173.

For the following procedures and therapies, the manufacturer does not expect any harmful interaction with the *Ikus* due to the general electromagnetic shielding of the device (see chapter 16, page 173). However, these procedures and therapies must only be applied after consultation with the treating physician.

- Radiotherapy
 - Nuclear diagnostics / nuclear therapy
 - Electro-stimulation therapy
 - Therapeutic ultrasonic treatment (e.g. lithotripsy)
 - External defibrillation
-

The following procedures and therapies have been tested in regard to their interaction with the *Ikus* and no harmful effects were found, however, these procedures and therapies must only be applied after consultation with the treating physician. Additionally the manufacturer does not guarantee that equivalent devices will not interfere.

- Diathermy
 - X-rays
 - Computed tomography
-

1.2 Precautions

1.2.1 VAD placement technique



Implantation - anesthesia: There should be an adequate supply of pre-matched stored blood, fresh frozen plasma and platelet concentrates available for immediate transfusion if required.

Implantation - anesthesia: Keep blood product transfusions to a minimum. Blood transfusions may lead to the development of antibodies, which are known to promote coagulation and inflammatory response.

The titanium connectors of the blood pumps have sharp edges designed to minimise the risk of clot formation at the junction. Be careful to avoid cutting yourself while connecting the pump and the cannulae.

1.2.2 Ambient conditions



The *Ikus* is intended solely for use in a hospital setting.

Before putting the *Ikus* into operation, check that the ambient conditions are suitable (see section 15.2: Technical specifications, page 163).

1.2.3 Caution while using the *Ikus*



When switching on the *Ikus* always make sure that the *Ikus* is switched on first and then the laptop and never vice versa! Otherwise there is a risk that during the start test error messages are falsely generated.

At least daily, the *EXCOR* cannulae should be inspected for signs of wear or damage. ADVICE: To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

At least every 4 hours, check visually that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, then take the appropriate corrective action.

Under certain circumstances, the message **left/right pump is not filling adequately** in some circumstances is not generated with the 10 ml *EXCOR* blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

After changing over to biventricular operation the device is operating in separate mode. All parameters are reset to the default parameters (see Tab. 14-3: Default standard parameters, page 152). The patient-customized parameters have to be adjusted again.

Do not switch the *Ikus* off unless the batteries are fully charged (i.e. all yellow charge indicator LEDs are on).

Replacing the blood pump due to growth of the patient: In children, plan to replace the pump(s) with a larger pump(s) in good time, to prevent the possibility of inadequate support due to an insufficient discharge rate.

Transport outside the clinic

CAUTION

The crate may only be transported as described in section 5.7: Transportation and packaging, page 60. Do not tilt or overturn the *Ikus* when it is packed inside the transport crate. Otherwise there is a risk that the *Ikus* is damaged or destroyed.

Always observe a resting period of 6 hours after each transportation before switching on the *Ikus*! The temperature of the *Ikus* should get adapted to the ambient temperature.

Keep all driving tube connectors covered at all times when not in use.

Also refer to section 5.7: Transportation and packaging, page 60.

1.2.4 Battery replacement and disposal

CAUTION

Only *Berlin Heart GmbH/ Berlin Heart, Inc.* service staff or persons authorized by the *Berlin Heart GmbH* service department may replace the batteries and dispose of them in accordance with the respective regulations.

1.3 Obligations of the operator

WARNING

Only qualified medical personnel trained specifically in the use of the system are permitted to work with EXCOR. Training courses can be arranged with *Berlin Heart, Inc.*

CAUTION

The operator (i.e. the hospital using the system) is responsible for instruction and care of the patient. The patient must be instructed on safety risks and cautionary measures (moisture, temperature, electromagnetic fields, etc.).

The operator is also responsible for adherence to the prescribed maintenance and service intervals (see section 1.2.4: Battery replacement and disposal, page 27).

A replacement *Ikus* and replacement equipment must always be available in the hospital.

If any of the components are damaged or if faults occur, inform *Berlin Heart GmbH/ Berlin Heart, Inc.* service department immediately. Do not use damaged components.

Only operate the *Ikus* with the components specified in this document. Never operate the *Ikus* with multiple-socket mains adapters or mains extension cables.

2 General Information

2.1 Device description

EXCOR is an extracorporeal, pneumatically driven ventricular assist device. It is designed to support the right and/or left ventricle when the native heart is unable to maintain normal blood flows and pressures even with help of drug therapy and intraaortic balloon counterpulsation. The device is designed for mid to long term mechanical support.

The EXCOR consists of 1 or 2 extracorporeal, pneumatically driven blood pumps and cannulae which connect the blood pump(s) to the atrium or ventricle and to the great arteries. The *Ikus* provides alternating air pressure to the blood pumps through driving tubes.

The blood pump is divided into an air chamber and a blood chamber by a multi-layer flexible polyurethane membrane. The alternating air pressure provided by the *Ikus* moves the membrane, thus filling and emptying the blood pump. Both the blood chamber and the polyurethane connectors are transparent to allow for visual detection of deposits and for monitoring the filling and emptying of the blood pump.

Valves (three-leaflet polyurethane valves) are located at the inlet and outlet positions of the blood pump connector stubs, thus ensuring the unidirectional blood flow.

Pulse rate, systolic drive pressure, diastolic suction pressure and the relative systolic duration can all be monitored and adjusted on the driving unit.

2.2 Indications for use

The EXCOR is intended to provide mechanical support as a bridge to cardiac transplantation for children. Pediatric patients with severe isolated left ventricular or biventricular dysfunction who are candidates for cardiac transplant and require circulatory support may be treated using the EXCOR.

For determination of the appropriate pump size, see section 15.1.2: Overview: Relationship: body weight – pump size, page 159.

2.3 Intended operation environment

Ikus is intended for use in a clinical setting. It can be used in any kind of hospital unit, e.g. OR, ICU, intermediate care unit or general care unit. It may be moved between clinical units using the built-in wheels, however in this case the patient must always be accompanied by a person trained in the use of the manual pump and emergency procedures. Thus, the patient shall be guaranteed care in case of an emergency.

Transporting the device during operation by any vehicles (e.g. ambulance, aircraft, etc.) is not allowed.

During movement of the device in operation within the clinic all electromagnetic compatibility precautions (EMC precautions) must be observed. See chapter 16, page 173. Otherwise there is a risk of electromagnetic disturbances and the fault-free operation of *Ikus* could not be guaranteed.

2.4 Contraindications

EXCOR patients with prosthetic aortic valves may have increased risk of thrombo-embolism due to the blood flow being shunted away from the native ventricular outflow. Patients whom cannot tolerate anticoagulation therapy should not be implanted.

2.5 Storage and durability



The expiration date of each EXCOR product is found on the product labels located on both the outer and inner packaging. The pumps, cannulae and accessories must not be used after the expiration date and even not be re-sterilized. Otherwise there is a risk of patient infection.

An EXCOR blood pump may not be used on a patient for more than 1 year. After this it shall be replaced with new products.

IMPORTANT: EXCOR must be stored at room temperature and be protected against extreme temperature fluctuations and moisture. Otherwise there is a risk of functional limitation and/or damage to the Ikus.

IMPORTANT: If the Ikus is not in use, run it once a month for 24 hours in order to ensure that all batteries are adequately charged. Refer to section 5.5.1: Routine start-test when not in operation, page 58. Otherwise there is a risk that the Ikus no longer functions correctly.

2.6 Battery replacement and disposal



Only *Berlin Heart GmbH/ Berlin Heart, Inc.* service staff or persons authorized by the *Berlin Heart GmbH* service department may replace the batteries and dispose of them in accordance with the respective regulations.

IMPORTANT: In the interest of achieving a definitive and comparable evaluation, the manufacturer assumes responsibility for both the disposal and analysis of explanted disposable articles.

2.7 Manufacturer's warranty

According to the General Terms and Conditions of the Berlin Heart GmbH the warranty is valid for 1 year.

All warranties apply only under the prescribed conditions of storage of the system, use in accordance with the instructions (intended use) and when the packaging is intact. This applies, in particular, to all sterile packaging and to the aluminum-coated outer packaging of the blood pump(s).

The warranty is no longer valid, if the Ikus has been opened or serviced by persons who are not members of the Berlin Heart GmbH/ Berlin Heart, Inc. service staff and/or who have not been authorized by the Berlin Heart GmbH service department to do so.

If an ambient temperature of +30°C is continuously exceeded during ongoing operation, the maintenance interval will be reduced.

Components	Product life (in sterile products starting from implantation)	Maintenance interval	Expiration date	Warranty
unsterile				
Ikus	max. 8 years	6 months or 2000 operating hours	x	1 year
Battery	exchange as needed	x	x	6 months
Manual pump	max. 6 years	yearly	x	1 year
sterile				
Driving tube	1 year (single-use product)	x	3 years	1 year
Blood pumps	1 year (single-use product)	x	3 years	1 year
cannulae	no limitation (single-use product)	x	3 years	1 year
Accessories	single-use product	x	3 years	first use

Tab. 2-1 Product life, maintenance interval, expiration date, warranty

3 Description: blood pump, cannulae and accessories

EXCOR is an extracorporeal electro-pneumatically driven ventricular assist device. It can be used for either univentricular or biventricular support. EXCOR is comprised of the following permanently active components:

- extracorporeal blood pump(s)
- inflow and outflow cannula(e)
- 1 driving tube for each blood pump
- *Ikus*




Fig. 3-1 EXCOR shown in situ as a biventricular assist device in pediatric application

Overview

The blood flows from the atrium or the ventricle through the inflow cannula into the blood chamber of the pump and then from this blood chamber through the outflow cannula into the aorta or into the pulmonary artery. A driving tube is used to connect the air chamber of the pump to the electro-pneumatic *Stationary Driving Unit Ikus*. *Ikus* generates the suction and driving pressures required to move the triple-layer membrane separating the blood chamber from the air chamber.

3.1 EXCOR blood pumps



Fig. 3-2 60 ml blood pump

EXCOR blood pumps have a transparent polyurethane (PU) housing which is divided into an air chamber and a blood chamber by a triple-layer membrane.

The blood chamber has an inflow and an outflow stub to which the inflow and outflow cannula, respectively, are connected. The pump stubs themselves are made of polyurethane, the end of each stub is fitted with a titanium connector to which the cannula will be connected. The valves located in the pump stubs keep the blood flowing in one direction. EXCOR blood pumps are available with three-leaflet valves made of polyurethane (10 - 60 ml stroke volume).

All surfaces of the pump coming into contact with the blood are coated with a *Carmeda® BioActive Surface (CBAS®)* coating. The transparent casing of the blood pump allows easy visual monitoring of the filling and emptying of the blood chamber.

The blood pump is equipped with a de-airing nipple which is used for de-airing the blood chamber when the pump is being commissioned.

The air chamber of the pump is equipped with a driving tube connector. This connector is used to connect the blood pump to the driving tube through which air is pumped from the *Ikus*. *Ikus* generates the suction and driving pressures required to move the blood pump's triple-layer membrane. A graphite powder layer is located between the membrane layers in order to minimize friction.

3.2 EXCOR cannulae

3 different types of cannulae are available for EXCOR in various sizes for each type:

- atrial cannulae (as inflow cannulae)
- LV apex cannulae (as inflow cannulae)
- arterial cannulae (as outflow cannulae)

The cannulae are made of tissue-friendly silicone. Polyester-velour suture rings enable convenient and safe anastomosis of the cannulae. The mid section of all cannulae is covered with polyester-velour in order to promote good ingrowth of the cannulae where they pass through the skin.

Some vascular cannulae have a shaping wire which allows the cannulae to be adapted to each individual patient's anatomic conditions.



Fig. 3-3 Cannula heads: 1) atrial cannula, 2) LV apex cannula, 3) arterial cannula

3.3 EXCOR accessories

The following *EXCOR* accessories are required in order to commission and operate *EXCOR*:

- 1 driving tube (PVC) for each blood pump
- 2 tank units
- 1 accessory (T00L-002) set which includes:
 - membrane set
 - de-airing set (2 x trocar, 2 x de-airing tube)
 - tube connecting set (cable ties, cable-tie gun)

There is enough material in 1 accessory set (T00L-002) to commission 2 *EXCOR* blood pumps.

4 Description: Ikus

4.1 Overview

The electro-pneumatic *Ikus* generates the suction and driving pressures required to drive the blood pump(s). The driving unit contains the pneumatic and electronic components as well as a laptop computer which serves as an interface to the operator.




Fig. 4-1 Ikus

Pneumatic systems

The *Ikus* has 3 pneumatic systems operating independently of each other. One pneumatic system is required for each blood pump, while the 3rd serves as an emergency backup. Each pneumatic system includes:

- a compressor
- pressure and suction limiters
- pressure and vacuum cylinders
- control electronic
- control valves

The compressor, the pressure and the suction limiters deliver constant pressure conditions in the pressure and vacuum cylinders. The control valves at the outlet of each cylinder allow optimal adjustment of the positive and negative (suction) pressure values.

Control computer

The *Ikus* system has 2 control computers operating independently of each other: the active (main) control computer and the backup control computer.

Laptop with monitor program

Messages and pressure graphs in the monitor program inform the user of the current status and working condition of the system. In addition, the laptop computer is used for commissioning the system and adjusting driving parameters. LOG files containing information on the system's operating status are recorded on the laptop's hard disk.

Manual pump

If there is no working *Ikus* available, the manual pump mounted on the *Ikus* can be used temporarily to drive the blood pump(s).

USB stick to store LOG files

On the last page of this instruction for use there are 2 USB sticks located in an envelope. The USB sticks are to be used for reading out and storing LOG files.

The 2nd USB stick serves as a backup, in the event that the first stick does not function properly. Do not under any circumstances use any other USB sticks on the *Ikus* laptop.

It is not permitted to use the USB sticks for purposes other than reading out and storing the LOG files on the *Ikus*. Never connect other USB devices (e.g. wireless technology) to the USB port of the laptop than the delivered USB sticks.

If both USB sticks are no longer available, contact *Berlin Heart, Inc.* immediately and request a replacement stick.

4.2 Displays and operating elements

4.2.1 Connection panel




Fig. 4-2 Connection panel



1 Driving tube connectors, left blood pump, red

The plug for the driving tube with the red marking is to be plugged to this connector.



2 Driving tube connectors, right blood pump, blue

The plug for the driving tube with the blue marking is to be plugged to this connector.



3 Button *Circuit breaker*

The *Ikus* is protected against overcurrent. The circuit breaker is resettable (see section 14.8: Circuit breaker and battery fuse, page 157).



4 Mains operation indicator

Illuminated when using power supply (default situation).



5 Power switch (toggle switch)

In this switch position (*[I]* position) the drive is operated through the mains. The batteries are charged at the same time.



In this switch position (*[0]* position), the mains power supply is interrupted. *Ikus* works in battery operation.



6 Mains power connector (with plug clip)

The *Ikus* is connected to the mains with the aid of the mains cable. The plug clip prevents an accidental loosening of the connection.

7 Main switch (key switch; 0 / I)

Drive switched off



Drive switched on



8 Connector *Potential equalization*

Should several electric devices be positioned in the vicinity of the patient, they are to be connected to a central grounding point. To this end, the connector *Potential equalization* is to be used on the *Ikus*.



9 Connector *External alarm (Nurse call)*

This connector serves as the link to the clinic's internal alarm system. A relay connects the internal alarm circuit with the connector *External alarm*. In the default configuration, the relay opens a pair of contacts whenever an alarm occurs, to set off an external alarm (via a 2-pin plug with a protective ground). Depending on the specifications of a specific hospital, the configuration of the connector *External alarm* can be modified so that the contact closes in case of an alarm. This setting modification is to be made by the *Berlin Heart, Inc.* service department.

Technical specifications: see section 15.2, page 163.

Main switch (key switch 0 / I) and power switch (toggle switch)



To prevent the batteries from being discharged by mistake, ensure that the power switch (toggle switch) is set to *[I]* position even if the *Ikus* is set to *[0]* position at the main switch (key switch).

While the *Ikus* is being operated, keep it connected to the stationary protective ground via the connector *Potential equalization*.

The main switch (key switch) is used to turn the driving unit on and off. The main switch (key switch) can only be operated with the key.

However, a connection to the mains is not established until the power switch (toggle switch) is set to *[I]* position as well. The power switch (toggle switch) is only meant to serve as an emergency measure to protect the driving unit in case of mains power disturbance. If the power switch (toggle switch) is set to *[0]* position, the *Ikus* is disconnected from the mains power supply and will continue to run automatically in battery operation. During this process, the batteries are discharged. When they are depleted, the *Ikus* will stop. For this reason, the power switch (toggle switch) must always be set to *[I]* position even if the *Ikus* is set to *[0]* position at the main switch (key switch).

4.2.2 Display and operating panel




Fig. 4-3 Display and operating panel

- | | | |
|---|--------------------------|--|
| 1 | Battery operation | The yellow LED lights up when the <i>Ikus</i> is running in battery operation |
| 2 | Alarm | The red LED lights up whenever there is an alarm message which has not been acknowledged by the monitor program |
| 3 | Mute alarm | Button to temporarily mute an alarm message. The message remains displayed on the laptop. |
| 4 | Operating hour counter | The number of all operating hours is displayed. When the <i>Ikus</i> was in operation for 2000 operating hours (maintenance interval) after the 1 st commissioning respectively after the last maintenance the display starts blinking and shows S -0.0h . The number of hours (i.e. S -10.0h) are displayed exceeding the maintenance interval. A maintenance sticker (on the right side under the laptop holder) contains the information when maintenance is due. |
| 5 | Battery charge indicator | 7 yellow LEDs indicate the battery charge state (see figure 4-3). The red LED lights up when the batteries are completely discharged! |



Fig. 4-4 Battery charge indicator with 7 yellow LEDs, red LED to the left of yellow LEDs

The LEDs of the running time display correspond to the maximum technically possible battery running time. However, the maximum permitted off-mains (battery) operating time is 30 minutes. After 30 minutes of battery operation, a warning message **Batteries discharged - use power supply!** is displayed. Immediately switch the *Ikus* over to mains operation. Failure to do so may result in failure of the device and damage to the battery system.

4.2.3 Power supply



The *Ikus* is designed for stationary operation and to be run on mains. Do not run it on battery operation unless this is absolutely necessary (e. g. when moving the patient within the clinic or during a mains failure).

To prevent the batteries from aging quickly, every period of battery operation should be followed by at least 6 h of mains operation.

Failure to do like mentioned above will severely reduce the capacity of the batteries and will greatly shorten the maximum off-mains (battery) operating time! This particularly applies to situations where the *Ikus* is operated at temperatures exceeding 30 °C. Always monitor the charging level display (see figure 4-3: Display and operating panel, page 40)!

Whenever the *Ikus* is running in battery operation, the patient must always be accompanied by a person trained to use the manual pump. Thus the patient shall be guaranteed care in an emergency.

Mains operation with integral battery charging function (mains operation)

Normally, the *Ikus* is powered by the mains. When the *Ikus* is in mains operation, the batteries which are completely or partially discharged are automatically recharged.

While the batteries are being charged, the battery charge indicator does not show the current charge correctly. The correct charge state is not indicated until the batteries have been fully charged. In this case, all yellow LEDs light up.

Battery operation

The *Ikus* has a rechargeable battery module (with 2 rechargeable batteries of 12V each) which can supply the system with power independently of the mains for a maximum of 30 minutes. Whenever the *Ikus* is being run in battery operation, an acoustic signal sounds at 10-minute intervals. When the battery charge is low, the acoustic signal sounds at 1-minute intervals.

4.3 Operating modes

4.3.1 Univentricular operation

- LVAD (Left Ventricular Assist Device)
- RVAD (Right Ventricular Assist Device)

In univentricular operation (also for an RVAD), the driving tube of the blood pump is always connected to the red marked connector.

4.3.2 Biventricular operation

Synchronous pulsing left/right

The systole cycles of the left and right blood pumps start simultaneously. Both pumps run at the same rate, the rate can only be adjusted via the left pump. The systolic pressure, diastolic pressure and the relative systolic duration can be set individually for each blood pump. Important: The rate of the right pump may not exceed that of the left pump.

Asynchronous pulsing left/right

The systole cycle of the right blood pump is started when the left pump switches to the diastole cycle. Both pumps run at the same rate, the rate can only be adjusted via the left pump. The systolic pressure, diastolic pressure and the relative systolic duration can be set individually for each blood pump. Important: The rate of the right pump may not exceed that of the left pump.

Separate mode left/right

Both blood pumps are operated completely independently of one another. All parameters can be adjusted as desired. Important: The rate of the right pump may not exceed that of the left pump.

4.4 Laptop computer with monitor program

NOTICE

If the laptop is switched off while the *Ikus* is in operation, the driving unit will continue to operate using the current parameter settings.

The laptop computer is integrated into the *Ikus* casing. The laptop, with its permanently installed monitor program, is used to start, set up and monitor the system.

The monitor program, when running, provides the user with a continuous stream of data on the system's condition, as well as information on events and function faults. At the same time, relevant data is stored on the hard disk for subsequent evaluation. A detailed description of the monitor program is given in chapter 5: Instructions for use: *Ikus*, page 45.



Fig. 4-5 Laptop




Fig. 4-6 Monitor program


4.5 Safety

4.5.1 Redundant design of pneumatic systems in univentricular operation

The blood pump is powered by 1 pneumatic system. The other 2 systems are redundant and serve as backups. If the active system fails, one backup system drives the pump. There is still an additional, redundant system available.

4.5.2 Redundant design of pneumatic systems in biventricular operation

Each blood pump is powered by 1 pneumatic systems. The 3rd system is redundant and serves as a backup. If one of the active systems fails, the 3rd (backup) system drives the corresponding pump. The performance of the *Ikus* remains the same. If 2 systems should fail, the remaining pneumatic system will drive both pumps. In this case, the *Ikus* is running in emergency operating mode.



Normal operation

Failure of 1 system

Failure of 2 systems

Fig. 4-7 Pneumatic system redundancy in biventricular operation

Normal operation	Failure of 1 system	Failure of 2 systems
System 1 left pump	System 1 failure	System 1 failure
System 2 right pump	System 2 right pump	System 2 failure
System 3 backup	System 3 left pump	System 3 left and right pump
	Message Backup operation left. Contact customer service. is displayed.	Message Emergency operation System 3. Contact service now! is displayed.

Emergency operating mode

If 2 pneumatic systems fail while in biventricular operation, the remaining system will drive both blood pumps. In this case, the *Ikus* will be operated in synchronous mode with a 250 mmHg systolic pressure, -100 mmHg diastolic pressure, 70 bpm and a relative systolic duration period of 40 %.

4.5.3 Control computer with redundancy design

The control computers are also designed to provide backup redundancy. The *Ikus* system has 2 independently operating control computers. The main computer controls the pneumatic systems and transmits important function-specific data to the laptop computer. The backup computer continuously compares its calculation results with those of the main computer. If it detects a difference, an error message is generated. For safety reasons, the control computers are located inside the *Ikus* casing and operate fully independently of the laptop.

Emergency pulsing

The final safety system of the drive's electronic system is the emergency pulse circuit board. Should both control computers fail, this hardware takes over control of the system. In order to preclude all possible sources of faults, the emergency pulse circuit board runs autonomously and can be influenced neither by the 2 control computers nor by the laptop.

In emergency pulse mode, the system will operate with the following settings:

Synchronous mode (biventricular)

	Systol. pressure [mmHg]	Diastol. pressure [mmHg]	Rate [bpm]	Relat. systole duration [%]
left	210	-40	70	40
right	150	-40	70	40

Tab. 4-1 Settings in emergency pulse mode

4.5.4 Battery operation

If the mains fails, the batteries will supply power to the system for at least 30 minutes. In this case, the battery operation indicator lights up. An acoustic alarm sounds at 10-minute intervals and a new message will be shown on the monitor program, showing the time the system has already been running in battery operation. The LEDs on the control panel of the handle show the charge state of the batteries.

4.5.5 Manual pump

If a working *Ikus* is no longer available, the blood pump(s) can also be temporarily operated (even in BVAD mode) by the manual pump supplied with the system (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

4.5.6 Password-protected user profiles (access passwords)

Only a user who has logged into the monitor program with the correct password will be allowed to change settings within the monitor program. When the system is delivered, a single default user profile is configured; up to 9 additional user profiles can be added. All system setting changes are logged specifically for the user who is carrying them out.

5 Instructions for use: Ikus

WARNING

Do not install any other software on this laptop.

Make sure that the <NumLk> key on the laptop is deactivated. The status LED on the laptop, identified by a lock symbol and/or a numeral (e. g. 1) should be off (see figures 5-1 to 5-3, page 45)

Only use the USB sticks provided with the device for saving data. Never connect any other USB stick to the laptop!

The laptop must be switched off before the USB stick is inserted in or removed from the USB port.

Check all the information and messages in the message window of the monitor program at least every 4 hours . If there is a message that indicates action must be taken, take the appropriate action and, if necessary, contact the *Berlin Heart, Inc.* emergency hotline.

When switching on the *Ikus*, always switch the *Ikus* on first and then the laptop and not vice versa.

NOTICE

For details on commissioning the system see chapter 6: Commissioning the *Ikus* and setting parameters, page 65.

The *Ikus* is controlled and monitored using the monitor program, which is permanently installed on the laptop. Only users who are properly logged in can use the monitor program to change system settings.



Fig. 5-1 LEDs below the mouse pad



Fig. 5-2 De-activate the *Caps Lock*.



Fig. 5-3 De-activate the *Numeric Lock*

5.1 Start menu

After the laptop is started, it will display the menu *Select language*.




Fig. 5-4 Select language

Afterwards the start menu will open up and offer the following options:



Fig. 5-5 Start menu

- | | |
|-------------------------------|---|
| 1. Start program | Select this option to start the monitor program. |
| 2. Entry codes | With this option password-protected user profiles (user no.: 1 digit; password: combination of max. 32 digits) for the monitor program can be edited. |
| 3. End | This option will shut down the operating system before switching off the laptop. The <i>Ikus</i> will continue to operate using the current parameter settings. |
| 4. Save data | This option allows the LOG files to be saved on a USB stick. This should only be done after having consulted the service department. |
| 5. Change date or time | This option enables the change of the system date and time on the laptop. |
| 6. Change language | Use this option to choose a different language for the monitor program user interface. |

5.1.1 Selecting an option in the start menu

After starting the monitor program (<1>, **1. Start program** in the start menu) the user has to log in with the user ID and password.

IMPORTANT: Before user profiles (<2>, **2. Entry codes**) can be entered, the initial user must first log in with user ID and password (as 4 space digit sequence) of the default user profile.

IMPORTANT: The default user profile is enclosed in the password envelope. This envelope can be found in the accessory box delivered with the *Ikus*.



Fig. 5-6 User ID



Fig. 5-7 Password

INSTRUCTION

1. Type in the number of the menu item (e. g. <1> for **1. Start program**). The option is activated immediately.

ADVICE

Keep the envelope containing the password in a safe place.

5.1.2 Configuring user passwords

A user profile (entry-code) is required to perform standard operating actions in the monitor program

IMPORTANT: The default user profile (user ID and password in password envelope as 4-digit code) is required in order to be able to manage user profiles.




Fig. 5-8 Entry codes

INSTRUCTION

1. If the monitor program is running: press <F10> to exit the monitor program. Confirm by pressing <X> or <1>.
2. From the start menu, select **2. Entry-codes** (<2>).
3. Enter the default user profile 4-digit code and confirm the input by pressing <Enter>. Now the password configuration program is started.
4. Use the arrow keys <↓>/<↑> to move the cursor to the desired user ID, then press <Enter> to confirm selection.
5. Enter the password (combination of max. 32 digits) and press <Enter> to confirm.
6. Repeat the entry and confirm with <Enter>.
7. If it is necessary, repeat steps 3 to 5 in order to configure additional user profiles.
8. To conclude and confirm all entries, press the <0> key. The passwords have now been assigned to the respective user IDs. The start menu is shown again.
9. In the start menu, select option **1. Start program** (<1>) to return to the monitor

program.

- Enter user ID and password, confirm by pressing <Enter>.

Other password configuration options

Key	Function
<ESC>	Discard all entries not yet confirmed with <0> and exit 2. Entry codes .
	Delete the selected user (press <Enter> to confirm)
<5>	Call up help text

Tab. 5-1 Other password configuration options

5.1.3 Saving data on USB stick

NOTICE

Consult the service department before proceeding with this option (refer to section 14.7: Reading out the LOG files, page 156).

5.1.4 Changing date or time



Fig. 5-9 Change date or time

INSTRUCTION

- If the monitor program is running: press <F10> to exit the monitor program and confirm by pressing <X> or <1>.
- From the start menu, select **5. Change date or time** (<5>).
- Move the cursor to the desired values using the <←> / <→> keys. Adjust the value with <↓>/<↑> (+/-1) and confirm with <Enter>.

Additional options changing date or time

Key	Function
<F6>	Switch between input field date and time fields
<F10>	Save settings and exit
<ESC>	Cancel and return to start menu

Tab. 5-2 Options changing date or time

5.2 Basic instructions for monitor program

5.2.1 Starting the monitor program

INSTRUCTION

1. The start menu is displayed on the screen. Select **1. Start program** (<1>).
2. Enter the user ID, then the password. Confirm the password with <Enter>.

5.2.2 Shutting down the monitor program

WARNING

Do not shut down the monitor program unless this is absolutely necessary (e. g. if new user profiles have to be set up). Restart the monitor program as soon as possible.

When the monitor program has been shut down:

- the system continues running with the currently set parameters
- no data on current events are recorded, the LOG files will be incomplete
- an acoustic signal and the indicator light on the handle control panel will alert the user of new messages. However there is no way to display the message that has been detected while the monitor program is shut down.

INSTRUCTION

1. Press <F10> and confirm by pressing <X> or <1>. The start menu is displayed on the laptop. The monitor program has now been shut down.

5.2.3 Logging in and out of the monitor program

IMPORTANT: This automatic forced logout does not occur after the first login during the very first commissioning of the *Ikus*.

IMPORTANT: The *Ikus* performs a self-test of the alarm circuit for each login procedure. If the self-test is completed successfully, the message **Acoustic alarm: OK** appears a few seconds later. **Important:** The self-test does not occur if, at the time of login, an alarm is pending that was muted with the button *Mute alarm* (see section 4.2.2: Display and operating panel, page 40).

IMPORTANT: If a problem arises during the self-test of the alarm circuit, the acoustic signal continues for test purposes and cannot be muted during this phase. Acknowledgment is only possible when a message appears (**Alarm circuit test failed - buzzer remains off!** or **Acoustic alarm is not properly recognized**).

To select options in the monitor program, to change parameters or scroll through the message window a user has to be logged into the monitor program. If no entries are made for several minutes, the monitor program will log the user out automatically. The *Ikus* continues to operate using the current parameter settings. The standard view is displayed.

5.2.4 Logging in



Fig. 5-10 User ID



Fig. 5-11 Password

INSTRUCTION

1. The monitor program is running. It now displays the field **User ID**. Enter user ID. The monitor program now displays the field **Password**.
2. Enter the password and confirm the password with <Enter>. The monitor program will now display the field **Log off**. The user now is logged into the monitor program. The alarm circuit test follows.

Logging out

INSTRUCTION

1. Select the menu item **Log off**, then press <Enter> to confirm. The monitor program now displays the field **User ID**. The user now is logged out of the monitor program. No further entries are possible until next login.

Cautionary measure

To ensure patient’s safety, the user always has to log out of the monitor program before leaving the vicinity of the system.

5.2.5 Standard view – monitor program

IMPORTANT: The displayed pressure graphs refer to the pneumatic pressures generated internally by the system and are given in millimeters of mercury [mmHg]. They do not display the patient’s blood pressure and cannot replace diagnostic measures by the medical personnel.

IMPORTANT: The monitor program also has 2 other views, which, however, are only displayed when the *Ikus* is being started or can be called up by the command **Drive pause**.

- view *Select operating mode*, see figure 6-3, page 68
- view *Pump size and single-step mode*, see figure 6-7, page 71




Fig. 5-12 Standard view of monitor program during start-test

- A Parameter table currently set parameters are displayed and can be changed here
- B Pressure graphs graphical display of the current pressure values in the respective pneumatic systems
- C Message window information on the system status (e. g. **test operation**) and messages (error messages and messages confirming error correction; see section 13: Error Messages and corrective measures, page 129).

Monitor program functions

Monitor and control system:

- If required: parameter adjustment
 - systolic pressure
 - diastolic pressure
 - rate
 - relative systolic duration (%Systole)
 - biventricular: operating mode (synchronous pulse; asynchronous pulse; separate mode left/right)
- View messages
- Acknowledge messages
- Pause drive
- Switch drive off

5.2.6 Selecting monitor program options

INSTRUCTION

1. Use the <←>/<→> keys to move the cursor to the desired field. A pop-up menu will be displayed. To move to the next line up or down, the <←> /<→> keys have to be pressed repeatedly.




Fig. 5-13 Selecting monitor program options

- with grey frame field is not activated
- black/ red frame field is activated
- colored background field is selected

Selecting pump size(s) and cannula sizes

INSTRUCTION

1. Use the <↓>/<↑> keys to move the cursor to the desired field. All options that can be selected for this field are displayed.
2. Pop-up menus in the parameter table: press <Enter> to confirm the selection.
3. Pop-up menus *Pump size* and *Cannula size*: after selecting the size, exit the field with <←>/<→>. The popup menu is closed and the selected option is displayed.




Fig. 5-14 Select pump size




Fig. 5-15 Select cannulae size

The fields highlighted in red in the pop-up menus *Select pump size* and *Select cannulae size* can be chosen. However, they are not the recommended option for each individual case (see section 15.1.10: Pump-cannula combinations, page 162 and section 15.1.11: Blood pump combinations in biventricular mode, page 162).

In individual cases, consideration is to be given as to whether a combination that is not recommended to be selected.

The final decision on the combination of blood pumps and cannulae is to be reached by the implanting surgeon, in consultation with *Berlin Heart, Inc., Clinical Affairs*.

5.2.7 Adjusting the parameter values

NOTICE

If the adjusted parameter was not confirmed with <Enter>, the parameter table will display the changed parameter until the automatic log off, but the *Ikus* will continue operating with the former value.

The changed parameter display can be canceled with <Esc>. The original value appears again in the parameter table.

To confirm the change, press <Enter>; the new value is displayed in the parameter table.

INSTRUCTION

1. Use the <<->/<->> keys to move the cursor to the desired field in the parameter table. The selected field is given a colored background.
2. Use the <↓>,<↑>/ <Bild-↓>,<Bild-↑> keys to adjust the value, then press <Enter> to confirm. The system will now operate using the new value.
3. Visually check that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles!

Parameter	Operation	Pressure [mmHg]		Rate	% Systole
		Systole	Diastole		
Left	L	200.0	-30.0	100.0	40.0
Right	R	170.0	-30.0	100.0	40.0

Fig. 5-16 Parameter table

Parameter	Range possible	<↓>/<↑> changes value by	<Bild-↓>/<Bild-↑> changes value by
Systolic pressure [mmHg]	60 to 350	2.5	25
Diastolic pressure [mmHg]	0 to -100	2.5	25
Rate [bpm]	30 to 150	1	10
Relative systolic duration [%]	20 to 70	1	10

Tab. 5-3 Parameter adjustment

5.2.8 Browsing in the message window

INSTRUCTION

1. Press the <7> key to move the cursor in the message window.
2. Scroll through the messages: either 1 message at a time with <↓>, <↑>; or 4 messages at a time with <Bild-↓>, <Bild-↑>.
3. Exit the message window by pressing either <Esc> or <Enter>. The field **Log off** is automatically activated and appears with a black background.

5.3 Stopping the blood pump(s) and switching off the Ikus




Fig. 5-17 Drive OFF (with stop options)

The pop-up menu *Drive OFF* offers the following options:

- **Drive pause**
- **Pause left**
- **Pause right** (only in biventricular operation)
- **Drive OFF**

5.3.1 Drive pause: stopping the Ikus temporarily

IMPORTANT: When restarting the system from a separate mode there can be an unintentional switch to the synchronous mode, although **L/R separate** continues to be displayed. The pressure graphs, however, distinctly indicate the synchronous mode.

INSTRUCTION

1. Select **Drive pause**, then press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The *Ikus* will stop. The view *Select operating mode* is displayed.




Fig. 5-18 Select operating mode

2. Select **Univentricular (UVAD)** or **Biventricular (BVAD)**, then confirm the selected operating mode with <Enter>. In biventricular mode the view *Pump size and single-step mode* appears. In univentricular mode, the seal test is performed first and then the view

Pump size and single-step mode appears.




Fig. 5-19 Pump size and single-step mode

- The cursor is located on the **OK** field. In order to restart the driving unit with the current settings, press the <Enter> key. Important: If the user leaves the **OK** field before pressing <Enter>, at least once **step left** or **step right**, respectively has to be performed again. Only now the **OK** field is highlighted again and it is possible to confirm **OK** by pressing <Enter> (see section 6.2.3: Select the pump size, page 69).

5.3.2 Pause left/ Pause right: stopping an individual blood pump

IMPORTANT: The option **Pause right** is only available in biventricular mode. The option is displayed in univentricular mode, but it is not activated. It is only possible to select **Pause left** (even for a RVAD).

➤ INSTRUCTION

- Select **Pause left** or **Pause right**, as required, then press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The selected pump will stop. The view *Pump size and single-step mode* is displayed.
- To restart the pump, at least once **step left** or **step right**, respectively has to be performed (see section 6.2.3: Select the pump size, page 69).

5.3.3 Drive OFF: switching the Ikus off



WARNING

The *Ikus* power switch (toggle switch) must always be in the [I] position even if the *Ikus* is set to [O] position at the main switch (key switch)!



CAUTION

The sequence of operations below must always be followed.

Do not switch the *Ikus* off unless the batteries are fully charged (i.e. all yellow charge indicator LEDs are on).

Keep all driving tube connectors covered at all times when not in use.

➤ INSTRUCTION

- Check the following: Are the batteries fully charged? (all yellow charge indicator LEDs should be lit). If not: Continue as instructed here up to and including step 4. Then switch off the laptop and leave the *Ikus* switched on until all yellow LEDs light up. Then switch the *Ikus* off using the main switch (key switch).
- In the monitor program, select the **Drive OFF** option and press <Enter> to confirm.
- Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The driving unit stops operation immediately and writes an operating LOG file.
- Wait until the LOG file has been completed. When the message **Switch drive off with main switch!** appears, press <F10> to shut down the monitor program. Confirm by pressing <X> or <1>.
- Select option **3. End** (<3>) in the start menu and switch off the laptop.
- Switch the *Ikus* off, when the batteries are fully charged. To do so, turn the main switch (key switch) to [O] position.

IMPORTANT: Always use the main switch (key switch) to switch off the *Ikus*.

IMPORTANT: If the *Ikus* is not in use, run it once a month for 24 hours in order to ensure that

all batteries are adequately charged. Refer to section 5.5.1: Routine start-test when not in operation, page 58. Otherwise there is a risk that the *Ikus* no longer functions correctly.

5.4 Switching between mains and battery operation



The *Ikus* is designed for stationary operation and to be run on mains. Do not run it on battery operation unless this is absolutely necessary (e. g. when moving the patient within the clinic or during a mains failure).

Whenever the *Ikus* is running in battery operation, the patient must always be accompanied by a person trained to use the manual pump.

Do not disconnect the *Ikus* from the mains unless the battery charge indicator shows that the batteries are full charged (all yellow LEDs light up). Important: In order to prevent premature discharge and aging of the batteries, the *Ikus* should be run in mains operation for at least 6 hours after every period of battery operation.

If the LEDs of the charge indicator flash on and off or if the message **Batteries discharged - use power supply!** is shown, switch to mains operation immediately!

If the batteries are discharged completely (red LED lights up), there is a danger of a total malfunction of the *Ikus* if battery operation is continued, and that the batteries will be damaged. If this happens, it can not be guaranteed that the *Ikus* will restart after connecting it to the mains.

When the battery charge is low, the acoustic signal sounds every minute. The *Ikus* must be connected to the mains operation immediately.



Always take along the mains cable when operating the *Ikus* in battery operation. In this way the system can be connected briefly to the mains again if it becomes necessary. If the patient has to be transported within the clinic choose a route as close to the power sockets as possible. Thus at any time the *Ikus* can be connected to the mains.



Fig. 5-20 Battery charge indicator with 7 yellow LEDs, red LED to the left of yellow LEDs

IMPORTANT: Whenever the *Ikus* is being run on battery power, an acoustic signal sounds at 10-minute intervals. These signals can be either acknowledged on the laptop or muted by pressing the button *Mute alarm* on the operating panel.

IMPORTANT: Battery charge indicator: 7 yellow LEDs indicate the battery charge state (see above). The red LED lights up when the batteries are completely discharged!

Switching over to battery operation**► INSTRUCTION**

1. Disconnect the *Ikus* from the mains. To do so, pull the plug from the socket. Never pull the cable from the driving unit.
2. Battery operation indicator lights up.
3. Observe all indicators and messages while the *Ikus* is running on battery power (see table 5-4, page 57).

Behavior of the *Ikus* with various battery charges

What?	When?	Meaning
The charge indicator shows the battery charge status	always This does not apply while the batteries are being charged (all yellow LEDs light up only <i>after</i> the batteries have been completely recharged)	charge status
Battery charge indicator lights up	during battery operation	battery operation
Message: battery operation , signal tone	every 10 minutes	OK ; can be either acknowledged on the laptop or muted by pressing the button <i>Mute alarm</i> on the operating panel
Message: Batteries discharged - use power supply!	after 30 minutes of operation, then at 10 minute intervals	The permissible runtime of 30 minutes on battery operation has been reached. Change to mains operation!
Message: Batteries discharged - use power supply!	Batteries have reached low charge state, 2-times per minute repeatedly	maximum possible runtime: only a few minutes. Switch to mains operation immediately!
yellow LED blinks (left)	Batteries are about to reach maximum possible runtime	Prewarning! Switch to mains operation immediately!
Message: Batteries discharged - use power supply! ...together with message (1-time): Fault: 0000 0011 1101 1000 (fault in power supply) red LED also lights up	Batteries are now no longer charged. (Batteries are completely dead) 2-times per minute repeatedly	Batteries are completely dead! Residual runtime: 0 minutes! Immediately switch to mains operation! Danger of total malfunction of ther <i>Ikus</i> !
Pump stands still - no further pump function; acoustic alarm is still audible	Batteries are complete emptied or there is a serious fault in the batteries	Immediately connect the <i>Ikus</i> to the mains. Until that supply the patient with the manual pump.
Acoustic alarm is audible; the displays are off and it is not possible to perform any action on the laptop	Batteries are defective or mains failure	Immediately connect the <i>Ikus</i> to the mains. Until that supply the patient with the manual pump.

Tab. 5-4 Displays and messages during battery operation mode

Switching over to mains operation

INSTRUCTION

1. Connect the *Ikus* to the mains.
2. The indicator *Mains operation* lights up. The indicator *Battery operation* in the handle goes out. The battery is now being charged.

5.5 Changing over from univentricular to biventricular operation

CAUTION

After changing over to biventricular operation the device is operating in separate mode. All parameters are reset to the default parameters (see Tab. 14-3: Default standard parameters, page 152). The patient-customized parameters have to be adjusted again.

INSTRUCTION

1. The pump which is to be connected must be already de-aired and connected to the cannulae.
2. In the monitor program, select the option **Drive pause** (see figure 5-17, page 53) and press <Enter> to confirm. Confirm the selection in the dialog window by pressing <X> or <1>. The active pump will stop. The view *Select operating mode* appears.
3. Open the driving tube connector identified by a blue marking.
4. When changing from RVAD to BVAD: disconnect the plug of the currently used pump from the connector marked red and connect it to the blue connector. To do so, take hold of the plug's release sleeve and pull the plug out of the socket.
5. Plug the driving tube of the de-aired new pump into the free driving tube connector. The sound of the plug snapping into place is clearly audible. Check that the plug is securely connected. To do so, grip the release sleeve above the grooved section and pull on it. Do not pull from the release sleeve, and never from the tube!
6. Select **Biventricular (BVAD)** and confirm with <Enter>. All following steps are the same as required for the start-up procedure. (see chapter 6: Commissioning the *Ikus* and setting parameters, page 65.)

5.5.1 Routine start-test when not in operation

IMPORTANT: If the *Ikus* is not in use, run it once a month for 24 hours in order to ensure that all batteries are adequately charged. Refer to section 5.5.1: Routine start-test when not in operation, page 58. Otherwise there is a risk that the *Ikus* no longer functions correctly.

INSTRUCTION

1. Procedure as described in section 6.1.1: Connecting the tank unit, page 66 and section 6.1.2: Switching on the *Ikus*, page 66.
2. Wait 24 hours.
3. Procedure as described in section 5.2.1: Starting the monitor program, page 49 and section 6.1.4: Setting the test parameters, page 67.
4. Observe and evaluate the curve representation.
5. Select the option **Drive OFF**, then confirm with <Enter>.
6. Confirm the selection in the dialog window by pressing <X> or <1>. The system immediately stops the operation and writes an operating LOG file.
7. Wait until the LOG file is complete (message: **Switch off drive with main key switch**).
8. End the monitor program and switch off the laptop.
9. All of the yellow battery LEDs are lighted: Switch off the *Ikus*. To do so, set the main switch (key switch) to [0] position.
Not all yellow LEDs are lit up: Leave the *Ikus* on the power supply until all of the yellow LEDs light up. Then set the main switch (key switch) to [0] position.
10. Disconnect both tank units from the *Ikus*.
11. Seal the driving tube connection sockets with seal plugs.

If the system detects an error in the start-up test**INSTRUCTION**

1. End the monitor program according to the instructions in the dialog.
2. Switch off the *Ikus* with the main switch (key switch).
Important: The power switch (toggle switch) remains on [I] position!
3. Important: Wait 10 minutes.
4. Restart the start-up procedure (see above).

If the system again detects an error in the start-up test or not all of the battery LEDs are illuminated after 24 hours, notify Berlin Heart:

HOTLINE

Notify Berlin Heart! 866.249.0128


5.6 Moving the Ikus

Fig. 5-21 Handle and lifting bar

Transport within the clinic**WARNING**

To move the *Ikus* unit: push only, using the handle provided for this purpose. Avoid any sudden jerky motion. When passing over smaller obstructions, exercise extreme caution, pulling the *Ikus* unit backwards (i.e. handle first) across the obstruction if necessary.

To lift the *Ikus*: use only the lifting bars at the lower edge at each side of the unit to hold and lift it. Never attempt to lift the *Ikus* by its handle. The *Ikus* must always be lifted by at least 2 people, preferably 4.

Rolling the *Ikus* over sloping surfaces: ensure that the person pushing it is strong enough to push the *Ikus* in a controlled manner. The slope of the surface may not be steeper than 10° (exception: packing/ unpacking of the *Ikus* into/ from the transport crate). Otherwise there is a risk of injury to the transporting persons or of damaging the *Ikus* !

If it is necessary to transport the patient within the clinic ensure that he is accompanied by a person trained to use the manual pump.

5.7 Transportation and packaging

CAUTION

Keep all driving tube connectors covered at all times when not in use.

To lift the *Ikus*: only use the lifting bars situated at the lower edge on each side of the unit to hold and lift it. Never attempt to lift the *Ikus* by the handle.

To transport the *Ikus*, only the supplied original transport crate must be used.

Make sure that the *Ikus* is standing firmly and securely in the crate and that the crate is closed and sealed properly.

The crate must always be transported as marked. Do not tip the *Ikus* driving unit stored in the transport crate or turn it upside down.

Always observe a resting period of 6 hours after each transportation before switching on the *Ikus*!

NOTICE

During storage and transportation, the power switch (toggle switch) should be in [I] position.

To transport the *Ikus*, only use the original transport crate supplied.

Ensure that the *Ikus* is placed firmly and securely inside the transport crate.

***Ikus* Accessories**

- 1 Instruction for use (supplied separately)
- 2 Tank units (supplied separately)
- 2 Keys
- 1 Power cable
- 1 Alarm connector
- 1 Envelope with password

Also see section 15.4: Sample copy: *Ikus* Incoming Checklist, page 166.

5.7.1 Unloading the *Ikus* from the shipping crate

NOTICE

Do not discard the shipping crate!

INSTRUCTION

1. Follow inspection steps per *Ikus Incoming Checklist*.
2. Place the crate with its back against a wall and/ or lock the brakes to prevent the crate from moving during unloading.
3. Open the door latches. Flip open the secure latch and turn counter clockwise to disengage latch from front door.
4. Swing front door open to the left.
5. Unlatch the ramp safety retainer.
6. Open the ramp while holding the black securing plate behind the ramp.
7. Open the front securing plate (is connected to the crate with a rope) and place the plate to the right of the crate.
8. Remove the white accessory box.
9. Disengage the *Ikus* front wheel brake (pulling the lever upwards).

10. Stand on the ramp and carefully pull out the *Ikus* by the handle bar.

For assistance please:



Notify Berlin Heart! 866.249.0128

5.7.2 Loading the *Ikus* into the shipping crate

During transport, the key remains in the main switch (key switch) in the [0] position.

Dimensions of the *Ikus* (W x H x D): 46 x 95 x 73 cm (approx. 18.5 x 37.5 x 29 inches) with laptop cover down

Weight of the *Ikus*: approx. 100.6 kg (approx. 219 lbs)

INSTRUCTION



Fig. 5-22 Empty shipping crate




Fig. 5-23 Prepare the laptop

1. Place the crate with its back against a wall and/ or lock the brakes to prevent the crate from moving during loading.
2. Open the door latches and swing front door open to the left. Unlatch the ramp safety retainer and open the ramp while holding the black door panel behind the ramp.
3. Make sure the laptop is closed and in parking position.



Fig. 5-24 Line up the *Ikus*



Fig. 5-25 Engage the brake

4. Line up the *Ikus* with the ramp. Push it upward on the handle bars; all the way into the crate.
5. Engage the *Ikus* front wheel brake (pushing the lever down).



Fig. 5-26 Place the accessory box



Fig. 5-27 Reattach the securing plate

6. Place the accessory box in front of the *Ikus* top of the bottom.
7. Reattach the securing plate:



Fig. 5-28 Lift up the ramp



Fig. 5-29 Fix the ramp

8. Lift up the ramp.
9. Fix the ramp with the safety retainer.




Fig. 5-30 Close the door




Fig. 5-31 Lock the crate

10. Close the front door and apply the latches.
11. Lock the crate with a cable binder.
12. The crate is now ready for storage/ shipping.



Fig. 5-32 Crate ready for storage/ shipping

6 Commissioning the *Ikus* and setting parameters

WARNING

Before starting the *Ikus*, make sure that a replacement *Ikus* is available in the hospital. If a replacement *Ikus* is not available, there is a risk that the patient cannot be cared for in the event of device malfunction.

The general rule is:

1 replacement *Ikus* if 1 or 2 systems are in use

2 replacement *Ikus* if 3 or 4 systems are in use

3 replacement *Ikus* if 5 or 6 systems are in use.

If more than 6 systems are in use the number of replacement *Ikus* has to be 1/2 of the active systems.

When switching on the system, make sure that the *Ikus* is first switched on and then the laptop and not vice versa!

CAUTION

Always observe a resting period of 6 hours after each transportation before switching on the *Ikus*! The temperature of the *Ikus* should get adapted to the ambient temperature.

Before putting the *Ikus* into operation, check that the ambient conditions are suitable (see section 15.2: Technical specifications, page 163).

Keep all driving tube connectors covered at all times when not in use.

NOTICE

This chapter describes the technical aspects of commissioning the system as well as perioperative and postoperative drive management. When commissioning the system, it is vital to observe the instructions given in chapter 7: Implantation: Preparations in the operating room, page 79.

The drive management procedures described here are intended as recommendations only. There is no substitute for careful patient observation and evaluation by the appropriate medical personnel.

6.1 Preparatory steps outside of the operating room

WARNING

Switch the *Ikus* on 2 hours before its intended use in order to adequately charge the batteries and so that a start-test can be performed to detect any possible faults in the device. During this period, always connect both tank units to the *Ikus* (see section 6.1.1: Connecting the tank unit, page 66)!

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the *Ikus* is set to [0] position at the main switch (key switch)! Otherwise the driving unit may stop completely because the batteries have become fully depleted.

The <NumLk> key on the laptop must be deactivated. The status LED on the laptop, identified by a lock symbol and/or a numeral (e. g. 1), should not light (see figures 5-1 to 5-3, page 45).




Fig. 6-1 Connection panel

6.1.1 Connecting the tank unit



Fig. 6-2 Tank unit

IMPORTANT: 2 separate tank units are used, each equipped with a driving tube.

IMPORTANT: The tank unit simulates a real operating situation during the warming-up period. If the test parameters have been set correctly (see section 6.1.4: Setting the test parameters, page 67), the tank units prevent the *Ikus* from generating false error messages.

INSTRUCTION

1. Remove the seal plugs from the driving tube connection socket.
2. Connect the plug of the 2 tank units to the driving tube connection sockets.

6.1.2 Switching on the *Ikus*

WARNING

The *Ikus* must always be connected to the power supply when it is switched on. This is the only way to ensure that the start-up test (see section 6.1.3: Starting the monitor program, page 67) is performed completely and possible malfunctions can be detected.

CAUTION

When switching on always take care to switch on the *Ikus* first and then the laptop. Never vice versa!

INSTRUCTION

1. Connect the *Ikus* to the mains. Secure the mains cable with the plug clip. Ensure that the mains power switch (toggle switch) is set to *[I]* position.
2. Turn the main switch (key switch) to the *[I]* position. The battery charge indicator will light up and the number of hours the driving unit has been operated will be displayed. The mains operation indicator lights up.

3. The menu *Select language* appears after the laptop is switched on.
4. Select the desired language by pressing the corresponding number key. It is not necessary to press <Enter> to confirm this selection.
5. The start menu is displayed on the laptop.

6.1.3 Starting the monitor program

INSTRUCTION

1. Select the **1. Start program** option (<1>) in the start menu.
2. Enter user ID and password, confirm with <Enter>. The *Ikus* will carry out a start-test.
3. Wait for the start-test to finish (it takes several minutes). Do not mute the acoustic signal. Exception: an error message is displayed. The messages in the message window will provide information on the status of the test. If the *Ikus* is found to be operating correctly, the view *Select operating mode* will be displayed next.

6.1.4 Setting the test parameters

NOTICE

When delivered, the *Ikus* has standard default parameters that it uses after each complete start-up process (see section 5.5.1: Routine start-test when not in operation, page 58 and section 14.3.2: *Ikus* start-test following emergency pulse mode, page 151). However, during the warm-up period with the tank unit, it is necessary to set the test parameters to prevent the *Ikus* from generating a false error message.

Systole [mmHg] left/right	Diastole [mmHg] left/ right	Rate [bpm]	Rel. systole duration [%] left/right	Operating type/mode
200	0	70	40	biventricular, synchronous mode

Tab. 6-1 Test parameters

INSTRUCTION

1. Even with planned univentricular operation: Navigate the cursor to **Biventricular (BVAD)**, confirm with <Enter>. The view *Pump size and single-step mode* appears.
2. Place the cursor with <↑> in the parameter table and navigate with <←>/<→> to the desired field in the parameter table. Adjust the value with <↓>, <↑> / <Bild-↓>, <Bild-↑>, then confirm with <Enter>.
3. Check that the batteries are sufficiently charged.
4. Disconnect the plug from the power supply and bring the *Ikus* immediately into operating room. Either acknowledge the message (**battery operation**) on the laptop or mute it by pressing the button *Mute alarm* on the operating panel.
5. In operating room: Reconnect the *Ikus* to the power supply.

If *Ikus* has detected a fault:

NOTICE

An overview of the messages that might occur during the start-test is given in section 13.24: Error messages during the start-up test, page 142.

The display **FAULT** in the parameter table indicates that the start-test has detected a fault. Such faults may be caused, for instance, by operating errors when the system was last shut down or during the start-up procedure.

INSTRUCTION

1. Shut down the monitor program as instructed in the dialog window.
2. Use the main switch (key switch) to switch off the *Ikus*. Important: The mains power switch (toggle switch) remains set to [I] position.
3. Wait for at least 10 minutes.
4. Commence the start procedure again (see section 6.1.3: Starting the monitor program, page 67).

If *Ikus* detects a fault again:

Do not use this *Ikus*.

HOTLINE

Notify Berlin Heart! 866.249.0128

6.2 Intraoperative drive management

6.2.1 Disconnecting the tank unit from the *Ikus*

INSTRUCTION

1. Select **Drive pause** then press <Enter> to confirm. Confirm the decision in the dialog window by pressing the <X> key or the <1> key. The *Ikus* will stop. The view *Select operating mode* is displayed.
2. Disconnect both tank units from the *Ikus*.
3. Seal the driving tube connection sockets with seal plugs.

6.2.2 Selecting the operating mode (view *Select operating mode*)




Fig. 6-3 Select operating mode

INSTRUCTION

1. In the operating room the *Ikus* should be reconnected to the mains immediately.
2. Select **Univentricular (UVAD)** or **Biventricular (BVAD)** with cursor, then confirm the selected operating mode with <Enter>.
3. In biventricular mode: the view *Pump size and single-step mode* appears (see figure 6-7, page 71).
In univentricular mode, a connector seal test is first performed: the *Ikus* checks whether the driving tube connector socket with the blue marking has been sealed.

Univentricular operation: connector seal test

NOTICE

The *Ikus* will repeat the connector seal test up to 3 times. If the test is still not successful, the system will switch itself off. Please contact the emergency hotline.

HOTLINE

Notify Berlin Heart! 866.249.0128

INSTRUCTION

1. The message **Please close right outlet** is shown. Check the driving tube connector identified by a blue marking. If the connector is open: use the seal plug to close it.
2. Move the cursor to **OK** and press <Enter> to confirm. The *Ikus* will test whether the connector is properly sealed. If it is, the view *Pump size and single-step mode* is displayed. If not, the *Ikus* will repeat the connector seal test.

6.2.3 Select the pump size

NOTICE

It is not possible to select a larger blood pump and/or a higher rate for the right pump than for the left one. The list of pump sizes available for the right pump is limited accordingly.

The sizes written in red ink in the pop-up menu *Select pump size* can be chosen. However, they are not the recommended option for each individual case (see figure 6-5, page 70).

In individual cases, consideration is to be given as to whether a combination that is not recommended is to be selected. The final decision on the combination of blood pumps and cannulae is to be reached by the implanting surgeon, in consultation with *Berlin Heart, Inc.*, Clinical Affairs.




Fig. 6-4 Select the pump size

INSTRUCTION

1. In the pop-up menu, select the desired pump size with <↓>, <↑> .
2. In biventricular mode: Using the cursor keys <←>, <→> , move to the field marked **right pump**, select the desired pump size with the cursor keys <↓>, <↑> .

6.2.4 Select the cannula size

NOTICE

The sizes written in red ink in the pop-up menu *Select the inflow and outflow cannula sizes* can be chosen. However, they are not the recommended option for each individual case (see section 15.1.10: Pump-cannula combinations, page 162 and section 15.1.11: Blood pump combinations in biventricular mode, page 162). In individual cases, consideration is to be given as to whether a deviant combination is to be selected.. The final decision on the combination of blood pumps and cannulae is to be reached by the implanting surgeon, in consultation with *Berlin Heart, Inc.*, Clinical Affairs.




Fig. 6-5 Cannula in for the left pump

INSTRUCTION

1. From the pop-up menu **cannula in** for the left pump, select the inflow cannula size for the left pump with <↓>,<↑>.
2. Move cursor with <←>/<→> to the pop-up menu **Canula out** for the left pump. With <↓>,<↑> select the outflow cannula size for the left pump.
3. In biventricular mode: Move the cursor with <←>/<→> to the pop-up menu **Canula in** for the right pump, select cannula size with <↓>,<↑>.
4. Move cursor with <←>/<→> to the pop-up menu **Canula out** for the right pump. With <↓>,<↑> select the outflow cannula size for the right pump.

6.2.5 Display pump and cannula sizes

It's possible to display the selected pump and cannula sizes on the laptop screen.

INSTRUCTION

1. Press <F10> to shut down the monitor program. Confirm the decision in the dialog window by pressing the <X> key or the <1> key.
2. The startmenu is displayed.
3. In the start menu, select the option **1. Start program** (<1>).
4. Enter user ID and password, confirm by pressing <Enter>.
5. In the message window the selected pump and cannula sizes are displayed.



Fig. 6-6 Display pump and cannula sizes

6.2.6 Setting the start-up parameters



To connect the blood pump(s), always set the start-up parameters!




Fig. 6-7 Pump size and single-step mode (with recommended start-up parameters)

INSTRUCTION

- Place the cursor with <↑> in the parameter table and navigate with <←>/<→> to the desired field in the parameter table. Adjust with <↓>, <↑> /<Bild-↓>, <Bild-↑> then confirm with <Enter>.
- Set the start-up parameters:

Systole [mmHg] left/ right	Diastole [mmHg] left/ right	Rate [bpm]	Rel. systol. duration [%] left/ right
100/ 80	-5/-5	30	40/40

Tab. 6-2 Start-up parameters

6.2.7 Connecting the blood pump(s) to the *Ikus*

WARNING

Do not kink either the driving tubes or the cannulae.

NOTICE

State of the blood pumps when they are initially connected: filled with sterile injectable saline, de-airing needle in place. To allow easier handling, the driving tubes are not connected until the inflow and outflow cannulae have been connected to the pump (see section 8.9: Connecting the blood pumps to the cannulae, page 93).

INSTRUCTION

1. Open the driving tube connector marked in red (univentricular) or both connectors (biventricular). To do so, pull the seal plugs out of the connector(s).
2. Connect the driving tube to the *Ikus*. To do so, push the plug of the driving tube into the connector. The sound of the plug snapping into place is clearly audible. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
3. In biventricular mode: observe the color of the markings.
4. In biventricular mode: repeat the procedure for the second pump.

Operating mode	<i>Ikus</i> connector
biventricular	LVAD: connector marked red RVAD: connector marked blue
univentricular	LVAD or RVAD: connector marked red

Tab. 6-3 Assignment: operating mode, blood pump, connector

6.2.8 De-airing the blood pumps in single-step mode

NOTICE

Each de-airing step (**Step left** / **Step right**) carries out half a pump cycle (systole or diastole), the 1th step being a diastole. Normally, several de-airing steps are required for each pump. In single-step mode, the pumps will operate using the pressures shown in the parameter table. It will not be possible to switch to the standard view unless at least 1 de-airing step has been completed for each connected pump.

INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. Move the cursor to the field marked **Step left**.
3. Lift the pump. The de-airing nipple is the highest point.
4. To trigger a single step, press the <Enter> key. If necessary, use the de-airing needle to vent the air from the pump (see section 7.5: De-airing the blood pump, page 81). After consulting the surgeon: If necessary, press <Enter> repeatedly to trigger further single steps until all air has been removed from the pump(s). If the the blood pump is not filling sufficiently, ensure there is sufficient preload and if necessary, increase the diastolic pressure.
5. In biventricular mode: Move the cursor to the field **Step right**. Repeat the procedure for the 2nd pump.




Fig. 6-8 Single-step mode

6.2.9 Starting the blood pump (changing to standard view)

WARNING

Do not start the pump(s) until all air has been removed.

Once the de-airing needle has been removed it cannot be re-inserted.

Only remove the de-airing needle after all air has been removed from the blood pump, the blood pump is running and the parameters have been adjusted (see section 6.2.10: Checking the parameters when the pump is started and adjusting them, page 73 and section 8.10: Removing the de-airing needle, page 94).

INSTRUCTION

1. Move cursor to the **OK** field and press <Enter> to confirm. The system now starts with the parameter values visible in the parameter table.

6.2.10 Checking the parameters when the pump is started and adjusting them

WARNING

In order to avoid air being sucked into the blood pump through the cannula anastomosis, adjust the parameters gradually. If air does enter the system, disconnect the driving tubes from the *Ikus* and de-air the system using the de-airing needle.

Continuously monitor all settings.

Once the de-airing needle has been removed it cannot be re-inserted.

NOTICE

If the pump is not filling adequately at this stage, increase the pre-load by adding volume from the CPB circuit. After adding volume, adjust the parameters on the laptop of the *Ikus* as described in the following table.

INSTRUCTION

1. Observe the left blood pump. Is the pump ejecting completely? If not: increase the left driving pressure if necessary.
2. Observe the right blood pump. Is the pump ejecting completely? If not: increase the right driving pressure if necessary.

Observe	Action / measure
Right pump Is the pump filling properly? (see below)	If not: check the filling pressure (central venous pressure; CVP) CVP too low: substitute volume CVP too high: increase suction pressure If no improvement occurs: check the position of the cannulae via echographic monitoring!
Left pump Is the pump ejecting properly?	If not: check mean arterial pressure (Guideline value: 70mmHg)

Observe	Action / measure
Compare left and right pump. Is left pump filling considerably worse than right pump?	If yes: increase suction pressure on left side If no improvement occurs: check the position of the cannulae via echographic monitoring!

Tab. 6-4 Pump filling criteria

Keep the following points in mind with regard to filling of the right pump:

The aim is to reduce the right ventricle’s load to a large extent but not completely. Signs that the RV load has been reduced completely are:

- filling of the pump depends largely on the respiratory cycle
- ventricle is empty/limp
- membrane stops abruptly during filling

Important: If the three above-mentioned phenomena are observed, do one of the following:

- reduce the diastolic pressure
- substitute volume

Adjusting parameters

►INSTRUCTION

1. Use the <←>/<→> keys to move the cursor to the desired field in the parameter table. The selected field is given a colored background.
2. Use the <↓>,<↑> or <Bild-↓>, <Bild-↑> keys to adjust the value, then press <Enter> to confirm the input.

Parameter	Range possible	<↓>/<↑> changes value by	<Bild-↓>/<Bild-↑> changes value by
Systolic pressure [mmHg]; driving pressure	60 to 350	2.5	25
Diastolic pressure [mmHg]; suction pressure	0 to -100	2.5	25
Rate [bpm]	30 to 150	1	10
Relative systolic duration [%]	20 to 70	1	10

Tab. 6-5 Parameter’s possible adjustments

In biventricular operation: adjusting the operating mode

To run the pumps in the asynchronous mode or separate mode instead of the synchronous mode the appropriate mode must be selected.

- asynchronous mode is recommended for patients who have a small thorax volume in comparison to the pump volume. In asynchronous mode, the intrathoracic blood volume remains unchanged.
- separate mode is useful, under some circumstances, for patients with intracardiac shunts.

►INSTRUCTION

1. Use the <←>/<→> keys to move the cursor to the field showing the current operating mode. A pop-up menu showing the available operating modes is opened (see table 6-3: Assignment: operating mode, blood pump, connector, page 72).
2. Select the desired operating mode with <↓>,<↑> and confirm with <Enter>. The system will now work in the selected mode.

Guideline values

The most important criteria when selecting drive parameters is that they ensure a good filling and emptying of the pump; the parameters must be set to achieve this goal.

NOTICE

The systolic driving pressure must be higher than the patient's physical systolic pressure. Important: If the systolic duration (% systole) is reduced or if very small cannulae are used, it may be necessary in some cases to select a higher value than recommended here.

The actual driving pressures achieved are influenced by the diameter of the cannulae used.

The following values are merely guideline values; they may not be appropriate in each individual case

Systolic pressure [mmHg], left/ right	Diastolic pressure [mmHg], left/ right	Rate [bpm]	Rel. systolic duration [%], left/ right
220/150	-40/-40	80	40/40

Tab. 6-6 Recommended guideline values for normal operation

ADVICE

Remove the de-airing needle after all air has been removed from the blood pump, the blood pump is running and the parameters have been adjusted (see section 6.2.10: Checking the parameters when the pump is started and adjusting them, page 73 and section 8.10: Removing the de-airing needle, page 94).

Important: Once the de-airing needle has been removed it cannot be re-inserted.

6.2.11 Switching from CPB support to VAD support

The aim here is to reduce the CPB flow and in doing so to shift the volume from the CPB to the patient (i.e. to the VAD).

WARNING

Secure the driving tubes and cannulae to the blood pump(s) as soon as the proper function of the EXCOR is established (see section 8.11: Securing the connections, page 94).

NOTICE

When using staged cannulae or a connecting set, do not set a pumping rate > 100 bpm, as the pump will not eject its full volume at higher rates. With these cannulae rates > 100 bpm are to be avoided.

INSTRUCTION

1. When the blood pump(s) starts to fill, reduce the CPB flow and gradually increase the EXCOR rate from an initial 30 bpm until CPB has been terminated and the required flow is achieved. Important: In doing so, make sure that the pump fills adequately, and if necessary regulate the driving pressure.
2. If necessary, adjust the systolic pressure, diastolic pressure and the systolic percent.

6.2.12 Possible complications

Decreased filling after stable filling conditions

If a good filling behavior was achieved at first (filling pressures LA/CVP < 10 mmHg and diastolic pulmonary artery pressure < 15 mmHg) with good drainage and nominal rate (normally 80 bpm), but the filling has over time, it usually will not help to increase the diastolic pressure.

Deterioration in the filling behavior despite stable inflow conditions may indicate hypovolemia or obstruction of the inflow cannula. The cause of deterioration in filling behavior must be identified and addressed.

NOTICE

Manipulations during implantation can severely influence the inflow temporarily – wait for the situation to stabilize before adjusting the values.

INSTRUCTION

1. Evaluate volume status and transfuse if necessary. Evaluate and if necessary correct the cannula position.

Pump filling deteriorates when thorax is closed

If atrial cannulation is used, a slight decrease in the filling may be observed in some cases when the thorax is closed. This may be caused by compression of the atria or a slight shift in the position of the cannulae.

INSTRUCTION

1. Evaluate volume status and transfuse if necessary. Important: Observe the effect volume replacement on the pump filling!
2. Increase suction pressure.

Distinct decrease in filling or generally poor inflow conditions on right side

INSTRUCTION

1. Make sure that there is no inflow obstruction.
2. If a suction pressure of less than -50 mmHg is necessary, increase the relative diastolic duration as an additional measure. At the same time, reduce the relative systolic duration. Important: Increase the driving pressure accordingly!

Incomplete ejection right/left

INSTRUCTION

1. Observe the arterial blood pressure, and at the same time observe the ejection movement of the pump membrane.
2. If complete emptying of the pump is no longer achieved, adjust the driving pressure accordingly. Important: Do not respond to extreme – temporary – increases in the arterial blood pressure (due to manipulation, catecholamine, etc.).

6.3 Postoperative drive management

NOTICE

The patient should receive the same treatment as is usual after any other major cardiac surgical procedure.

6.3.1 After transfer to the ward

If a good filling and stable ejection of the blood pump(s) is observed in the immediate post-operative period, it is normally not necessary to adjust the driving and suction pressures.

- Good filling means that the suction pressure is adequate.
- Stable ejection (at normal arterial blood pressure) means that the driving pressure is adequate.

WARNING

At least every 4 hours, visually check that the pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, appropriate measures are to be taken.

NOTICE

For further details on regular monitoring of pump(s) and cannulae, see section 10.5: Regular checks of blood pump(s) and cannulae, page 103.

6.3.2 Follow-up treatment

Guideline values and criteria for adjusting the parameter settings: see table 6-6: Recommended guideline values for normal operation, page 75.

It is only necessary to adjust the left driving pressure when

- the arterial blood pressure increases (e. g. after lifting sedation, when the patient wakes up)
- when the patient is mobilized (moving to an upright position, sitting, standing – in order to compensate for the additional hydrostatic pressure component).

7 Implantation: Preparations in the operating room

7.1 Preparing the components and materials required

NOTICE

Selection of blood pump(s): see section 15.1: Overview: Product range and possible combinations, page 159.

ADVICE

It is advantageous to provide a sterile table on which to place the prepared sterile components.

General (all sterile)

- 500 ml sterile injectable saline
- 2 small sterile basins
- 50 ml disposable syringe with luer lock connector
- suture (to secure the trocar to the de-airing nipple and the de-airing tube to the trocar)
- heavy scissors
- towel clamp, tube clamp
- other instruments and equipment as required for open-heart surgery

EXCOR components and accessories

- blood pump(s), each with a pump seal
- 1 driving tube for each blood pump
 - univentricular: driving tube, red (for LVAD and RVAD)
 - biventricular: 1 red driving tube and 1 blue driving tube
- inflow cannula(e) (atrial or LV apex cannula)
- outflow cannula(e)
- accessory set (T00L-002) for blood pumps with PU valves
 - membrane set
 - de-airing set (2 x trocar, 2 x de-airing tube)
 - tube connecting set (cable ties, cable-tie gun)

7.2 Checking and adjusting the settings of the cable tie gun

Before using the cable tie gun contained in the EXCOR *Tube connecting set* the accuracy of settings has to be checked and if necessary to be corrected.



Fig. 7-1 Cable tie gun

INSTRUCTION

1. Check if the following values are set:
 - coarse adjustment on STD (2)
 - fine adjustment on 5 (1)



2. In the case of deviations loosen the screw (4) and disassemble the locking cap (3).



3. Adjust the above-mentioned values with the adjusting wheels (6 and 5). Begin with adjusting wheel 6.



4. Assemble the locking cap (3) and secure it with the screw (4).



7.3 Unpacking the sterile components

WARNING

Only use sterile components which have been delivered in undamaged sterile condition (sterile packaging intact, expiration date not expired).

Only use blood pumps which have an undamaged aluminum-coated outer packaging.

INSTRUCTION

1. Pump: a non-sterile person opens the aluminum-coated package and removes the pump in its double sterile packaging.
2. The non-sterile person opens the outer sterile package.
3. A sterile person takes out the inner sterile package, opens it and places the components on the prepared sterile field.

7.4 Moving the membrane to the end-of-diastole position



Fig. 7-2 De-airing nipple and driving tube connector

INSTRUCTION

1. Pick up adapter tube, disposable syringe (membrane set) and the pump.
2. Connect the adapter tube to the disposable syringe.
3. Connect the free end of the adapter tube to the driving tube connector of the blood pump.
4. Remove all air from the air chamber of the pump. The blood pump membrane is now in the end-of-diastole position.
5. Seal the adapter tube with a tube clamp in order to keep the membrane in this position.

7.5 De-airing the blood pump

Prepare and place the following ready for use:

- blood pump(s) with pump seal(s)
- 1 de-airing set (trocar and a de-airing tube) for each blood pump
- 50 ml disposable syringe for each blood pump




Fig. 7-3 Pump with trocar in place (de-airing needle with inserted obturator)

7.5.1 Inserting the de-airing needle

WARNING

The membrane must be kept in the end-of-diastole position. Keep the clamped membrane set connected to the blood pump.

INSTRUCTION

1. Take hold of the trocar (de-airing needle with obturator) and remove the protective silicone cap.

2. Push the trocar as pictured above as far as it will go through the center of the blood pump's de-airing nipple. Never turn the trocar when inserting it, this increases the risk of removing a large piece of the silicone material in the de-airing nipple.
3. Remove the obturator.
4. Withdraw the de-airing needle by approx. 2 mm. Important: The tip of the cannula should still be visible in the blood chamber.
5. Use the suture to fix the de-airing needle to the de-airing nipple.
6. Remove the adapter tube from the pump.

7.5.2 Rinsing and filling the blood pump

ADVICE

Before commencing surgery, mark the points for the exit sites of the cannulae. The aim is to achieve a stable final position of the cannulae without exerting any tension on the skin. Caution: with biventricular support, 2 of the 4 cannulae will cross each other. This crossing point should be outside of the thorax as far as possible.

INSTRUCTION

1. Fill and empty the pump once or twice with sterile injectable saline.
2. Push the free end of the de-airing tube onto the trocar as far as it will go. Secure the de-airing tube to the trocar with a suture tie.
3. Fill the syringe with sterile injectable saline.
4. Connect the syringe to the stopcock end of the de-airing tube.
5. Slowly fill the pump with sterile injectable saline. Rock the pump back and forth to move any bubbles to the outflow stub.
6. Close the stopcock on the de-airing tube.
7. Tap the blood pump body gently in order to free all remaining bubbles. Remove all air from the pump through the outflow connector.
8. Use the seal caps to close the titanium cannula connectors.
9. Place the pump ready for connection with the connectors pointing up.

8 Implantation - surgical procedure

This chapter describes the product-specific measures to be observed when implanting an EXCOR blood pump. The setting of the *Ikus* parameters during and after implantation is described in chapter 6: Commissioning the *Ikus* and setting parameters, page 65. If the *Ikus* is brought into the OR, the *Ikus* should be prepared as described in section 6.1: Preparatory steps outside of the operating room, page 65.

Unless any specific instructions to the contrary are given, the same protocol as for any other major cardiothoracic surgical procedure should be followed. Implantation is accomplished using a CPB with bicaval cannulation. Implantation can be achieved with induced ventricular fibrillation or on a beating heart, hypothermia is usually not required.

WARNING

After implantation each cannulae and all connections must be inspected for it's solidity, safeness and tightness.

Do not start pump operation until the blood pump is completely free of air!

Do not touch or manipulate the blood pump with pointed or sharp-edged objects (e. g. surgical instruments)!

If a cannula is bent with flexible metal reinforcement to adjust it to the anatomical conditions: determine by visual inspection that the blood flow in the cannula is not restricted.

NOTICE

For the suture use an appropriate suture material. It should be a nonabsorbable monofilament, not traumatizing material.

ADVICE

For BVAD, carry out anastomosis of the cannulae in the following order:

- apical cannulation
 1. LV apex
 2. right atrium
 3. pulmonary artery
 4. aorta
 - atrial cannulation
 1. left atrium
 2. right atrium
 3. pulmonary artery
 4. aorta
-

8.1 Cannula exit sites




Fig. 8-1 Cannula position following implantation




Fig. 8-2 Suggested cannulae exit sites (Example: BVAD with LV apex cannulation)

8.2 Use of the cannula tunneling tip

The cannula tunneling tip is a sterile disposable product and is supplied with each cannula. Sizes available: see figure 8-3: Available sizes of cannula tunneling tips, page 85. Staged cannulae are supplied with 2 different tunneling tips.

INSTRUCTION

1. Push the cannula tunneling tip firmly into the distal end of the cannula.
2. Advance the forceps through the subcostal incision and the cannula tunnel into the mediastinum, so that the cannula tunneling tip can be gripped.
3. Use the forceps to firmly grip the flat end piece, pull it through the cannula tunnel and the skin incision and position it.
4. Carefully remove the tunneling tip from the cannula by bending it back and forth.

Refer to the respective cannula type as described in sections 8.3 to 8.6 of the instruction for use to determine the sequence of cannulae anastomosis and tunneling.




Fig. 8-3 Available sizes of cannula tunneling tips




Fig. 8-4 Use of cannula tunneling tip

8.3 Cannulae and connector set

To avoid damages of cannulae careful attention should be paid to the following safety precautions.

WARNING

During implantation the *Cannula Tunneling Tip* (provided with each cannula) should be used during implantation of the *EXCOR* system.

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
- After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.

If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

If the connector set is being used: Secure each connection between blood pump and cannula with at least one cable tie.

If an *EXCOR connecting set* is required for implantation and the length of the tube part needs to be reduced, the tube part should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 8-1 Connector set: minimum length of connector tube

⚠ WARNING

If replacement of an *EXCOR* blood pump is required, the following procedures should be observed:

- The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.
 - If a connecting set needs to be cut for a pump replacement, ensure that there will be sufficient length of the tube part remaining to meet the minimum length recommendations.
-

Do not kink the drivelines. Otherwise there might not be sufficient pump output

Do not kink the cannulae needlessly. Otherwise there might not be sufficient pump output. Moreover, cannulae might be damaged.

At least daily, the *EXCOR* cannulae should be inspected for signs of wear or damage. **ADVICE:** To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

In no case should the cannulae either be kinked directly at the connector to the blood pump or at the transition area between velour and silicone.

8.4 Access

➤ INSTRUCTION

1. Median sternotomy. Make sure that there is absolutely no bleeding.
2. Insert standard cardiopulmonary bypass cannulae (bicaval cannulation).
3. Initiate extracorporeal circulation.
4. Place a vent in the left atrium, if necessary.

8.5 LV apex cannula

Refer to section 8.2: Use of the cannula tunneling tip, page 85.

8.5.1 Anastomosis of inflow cannula with LV apex

WARNING

During anastomosis of the LV apex cannula, make sure that the cannula head is facing in the right direction: the long side of the head should be parallel to the lateral wall. This prevents the ventricular lateral wall from being sucked into the tip of the cannula. After the cannula head has been placed, its position can be checked by means of the flow direction arrow on the cannula body (except LV apex cannulae C10A-030, C14A-040, C18A-020). The arrow is aligned with the long side of the cannula head (see figure 8-6: Ideal position of the LV apex cannula, page 87).

INSTRUCTION

1. If indicated, initiate ventricular fibrillation as needed.
2. Apical excision of the LV: The ideal implant position of the LV cannula is slightly off-center of the LV apex toward the lateral wall. The distance from LAD/ septum to the center of the excised muscle core is about 2 cm for children.
3. We recommend to excise a circular apical core with a diameter slightly smaller than the size of the cannula head.
4. Start with muscle core incision on the side away from the septum/ LAD (see b in figure 8-6, page 87) to avoid septal injury.
5. Check left ventricle for thrombi and excise the excess trabeculae.




Fig. 8-5 Anastomosis of LV apex cannula




Fig. 8-6 Ideal position of the LV apex cannula



a Long side of LV apex cannula head

Fig. 8-7 Head of LV apex cannula

8.5.2 Creating a transcutaneous tunnel for the LV apex cannula

WARNING

Always use the cannula tunneling tip provided (see section 8.2: Use of the cannula tunneling tip, page 85) to advance the cannula through the prepared transcutaneous tunnel. Never use a sharp surgical instrument directly on the cannula.

Make sure that the blood pump and cannulae come to rest in a stable position without tension.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
- After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.
- If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

The skin incision must be slightly smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible, the cannula exit sites should be on different planes (see fig. 8-2, page 84).

INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Incise the pericardium widely in a lateral direction. Prepare the cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Tunnel the LV apex cannula through the transcutaneous passage by using a pair of forceps to firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision..
Important: Do not rotate the cannula while pulling it through the tunnel. At the end of this procedure, the apex of the heart should be in its native position without torsion.
4. Terminate ventricular fibrillation if necessary.

8.6 Atrial cannula(e)

Refer to section 8.2: Use of the cannula tunneling tip, page 85.

ADVICE

For atrial cannulae supplied with a forming wire, the transcutaneous tunnel should be created and the cannula advanced through the tunnel and skin incision prior to the anastomosis.
For all other atrial cannulae, the sequence is arbitrary.

8.6.1 Creating a transcutaneous tunnel for atrial cannula(e)



If possible, always use the cannula tunneling tip provided (see section 8.2: Use of the cannula tunneling tip, page 85) to advance the cannula through the prepared transcutaneous tunnel.

If it is necessary to apply a clamp directly to the cannula in order to pull the cannula through the skin, the following procedures should be observed:

- Position the clamp at the distal end of the cannula
- After the cannula has been pulled through the skin, cut off and discard the part of the cannula where the clamp was applied.
- If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.

Care must be taken to ensure that the cannulae come to rest in a stable position free of tension.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision. Important: Do not rotate the cannula while pulling it through the tunnel.

The incision must be slightly smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible the cannula exit incisions should be on different planes.

INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Prepare the cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Using a pair of dressing forceps, tunnel the cannula through the transcutaneous tunnel. Important: Do not rotate the cannula while pulling it through the tunnel.

8.6.2 Anastomosis of atrial cannulae

Right atrium



Create the anastomosis laterally, directly above the tricuspid valve.


a) closed technique

Fig. 8-8 Cannulation of right atrium




Fig. 8-9 Suture technique, right atrium

INSTRUCTION

1. Make a running (purse-string) suture with monofilament, secured with pledgets at 4 positions.
2. Place 4 single U-sutures secured with pledgets on each side of the purse string suture.
3. Make a sufficiently long incision inside of the suture circle and extend it as required.
4. Push the cannula down on the sutures, at the same time slightly reduce the venous inflow to the CPB while inflating the lung in order to prevent negative pressure in the left atrium.
5. Remove all air from the cannula and use a tube clamp to clamp the cannula below the anastomosis.

b) open technique with bicaval cannulation

With bicaval cannulation, the right atrial cannula can be inserted in an open technique.

Left atrium

The procedure for anastomosis of the left atrium corresponds to the procedure applied to the right atrium.

ADVICE

Place anastomosis at the junction of the right upper pulmonary vein and the left atrium. The atrial wall is the recommended implantation location. The pulmonary vein should be left intact.

8.7 Arterial cannula(e)

Refer to section 8.2: Use of the cannula tunneling tip, page 85.

ADVICE

For cannulae supplied with a forming wire, the transcutaneous tunnel should be created and the cannula advanced through the tunnel and skin incision prior to the anastomosis.

8.7.1 Creating a transcutaneous tunnel for arterial cannula**WARNING**

Care must be taken to ensure that the blood pump and cannulae come to rest in a stable position.

Do not touch or manipulate the silicone cannulae with pointed or sharp-edged objects (e. g. surgical instruments).

Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision. Important: Do not rotate the cannula while pulling it through the tunnel.

The incision must be smaller than the cannula diameter (to ensure good ingrowth) but large enough to prevent skin necrosis.

Plan the cannula exit sites appropriately. Leave an adequate bridge of skin and subcutaneous tissue between the cannula exit incisions to prevent breakdown and necrosis of the skin and tissue. If possible the cannula exit incisions should be on different planes (see fig. 8-2, page 84).

► INSTRUCTION

1. Prepare the transcutaneous tunnel. Ensure that the incision is large enough.
2. Prepare cannula tunnel by blunt dissection. Important: Do not tunnel transperitoneally.
3. Using a pair of forceps, firmly grip the flat end piece of the tunneling tip and pull it through the cannula tunnel and the skin incision. Important: Do not rotate the cannula while pulling it through the tunnel.

8.7.2 Anastomosis of the arterial cannula

Aorta

► INSTRUCTION

1. Tangentially clamp the ascending aorta and make a longitudinal opening of a length which is suitable for the cannula diameter. If necessary, offset the incision laterally to the right by up to 45°.
2. Anastomose the cannula using ten teflon-backed double-reinforced individual monofilament (e. g. 4-0 EB) U-sutures. (If simpler conditions are encountered, a running suture can be made instead.)
3. Remove all air from the cannula and use a tube clamp to clamp the cannula below the anastomotic site. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.




Fig. 8-10 Anastomosis of the aortic cannula

Pulmonary artery

► INSTRUCTION

1. Make a longitudinal incision of a size suitable for the cannula diameter in the pulmonary artery.
2. Anastomose the cannula using 10 teflon-backed, double-reinforced individual monofilament (e. g. 4-0 EB) U-sutures. (If simpler conditions are encountered, a

- running suture can be made instead.)
- Remove all air from the cannula and use a tube clamp to close it below the anastomosis. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge




Fig. 8-11 Cannulation of the pulmonary artery

8.8 Shortening the cannulae if necessary

WARNING

If an *EXCOR* connecting set is required for implantation and the length of the tube part needs to be reduced, the tube part should be cut but only to achieve the following minimum lengths:

Part Number	Diameter Reduction	Minimum Length
A12-016	16 to 12 mm	90 mm
A09-012	12 to 9 mm	75 mm
A06-009	9 to 6 mm	60 mm

Tab. 8-2 Connector set: minimum length of connector tube

INSTRUCTION

- Cut the cannulae to the required length. Make the cut perpendicular to the cannula axis and ensure that the cut is straight.
- Make sure that the lengths of the 2 cannulae leading to the same pump match. It must be possible to connect the cannulae to the pump without having to exert any tension.

8.9 Connecting the blood pumps to the cannulae

⚠ WARNING

Ensure that cannulae, blood pump(s) and driving tubes are not subject to external forces and are free of kinks or sharp bends.

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs. These show the direction of the blood flow.

Type of support	Anastomosis of inflow cannula to	Points upwards...
LVAD	apex	blood chamber
RVAD	atrium	air chamber
LVAD	atrium	air chamber


Tab. 8-1 Anastomosis and direction of the blood chambers

NOTICE

Finally, the driving tube is connected to the *Ikus*. The *Ikus* is started and the parameters are gradually adjusted (see section 6.2.10: Checking the parameters when the pump is started and adjusting them, page 73).

➤ INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. Release the tube clamps, flush the cannulae and then use tube clamps to clamp the cannulae below the exit sites. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
3. First connect the inflow cannula to the pump, then connect the outflow cannula. When doing so, add sterile injectable saline with a bulb syringe in order to connect the pump air free. Be careful to avoid damaging the gloves and the inner cannula (lumen) and pump surfaces.
4. Release the tube clamps, de-air the pump(s) and the cannulae.
5. Connect the driving tube to the blood pump. Biventricular: use the red driving tube for the left blood pump and the blue driving tube for the right blood pump. Univentricular: always use the red driving tube.



- a right pump (air-chamber pointing upwards)
- b inflow cannula from right atrium
- c outflow cannula to pulmonary artery
- d left pump (blood-chamber pointing upwards)
- e inflow cannula from LV apex
- f outflow cannula to ascending aorta

Fig. 8-12 Final position of the blood pumps, for example: BVAD with LV apex cannulation

8.10 Removing the de-airing needle

WARNING

When removing the de-airing needle, never pull on the de-airing tube, but on the de-airing needle itself.

Before removing the de-airing needle, be sure that the de-airing tube is secured to the de-airing needle. Important: Once the de-airing needle has been removed it cannot be re-inserted.

NOTICE

Do not remove the de-airing needle until all air is removed, the blood pump is running, all parameters have been adjusted and the chest has been closed. (see section 6.2.10: Checking the parameters when the pump is started and adjusting them, page 73).

INSTRUCTION

1. Cut the suture material between the de-airing needle and the de-airing nipple (see image 1 in figure 8-13, page 94). Important: Leave the ligature around the de-airing nipple (see image 2 in figure 8-13, page 94).
2. Pull the de-airing needle out of the de-airing nipple.




Fig. 8-13 Removing the de-airing needle

After the patient has been weaned from the CPB and the proper function of the EXCOR is established, the connections of the driving tubes and cannulae to the blood pump(s) have to be secured.

8.11 Securing the connections

WARNING

All connections have to be secured by at least 1 cable tie. 2 cable ties may be used. Exception: connection between drive line and drive line connector of the blood pump: 1 cable tie only!


INSTRUCTION

1. Pick up the Tube connecting set.
2. Secure the following connections:
 - inflow cannula on the connector
 - outflow cannula on the connector
 - drive line on the drive line connector (1 cable tie only!)

3. The 1. cable tie must be positioned exactly on the groove profile of the connector (1).
Important: the heads of the cable ties have to be directed away from the patient's body.
4. Fasten the cable ties by the cable tie gun. Important: pay attention to 7.2, page 79.



5. A 2nd cable tie can be used optionally. If a 2nd cable tie shall be used (2) it has to be positioned above the 1st cable tie. **IMPORTANT:** the heads of the cable tie straps should both be staggered and directed away from the patient's body.



6. If an EXCOR Connecting set is required for implantation after that secure also those connections with cable ties. Proceed thereby as described in the instruction steps 3 to 5.

9 Implantation - anesthesia

The following risk factors should be closely monitored for anesthetic and hemodynamic management:

- right heart function during LVAD implantation
- coagulopathy
- renal insufficiency
- abnormal reactions to inotrope administration
- pulmonary hypertension

CAUTION

There should be an adequate supply of pre-matched stored blood, fresh frozen plasma and platelet concentrates available for immediate transfusion if required.

Keep blood product transfusions to a minimum. Blood transfusions may lead to the development of antibodies, which are known to promote coagulation and inflammatory response.

ADVICE

Medication for right ventricular afterload reduction should be available for use in the operating room (nitric oxide NO, phosphodiesterase inhibitor, prostaglandin, etc)

Auto-transfusion equipment (e. g. Cellsaver) should be available for use in the operating room.

For patients with an LVAD, start ventilation with nitric oxide or administer the appropriate medication to treat pulmonary hypertension and reduce afterload for right ventricle 15 minutes before weaning from the CPB. This can help to prevent or lower the risk of right ventricular failure.

Monitoring procedure

Intraoperative monitoring should include the same monitoring procedures applied during major cardiothoracic surgery:

- central venous line
- Swan-Ganz catheter (if appropriate)
- arterial line
- ECG
- pulseoximetry
- central temperature monitor
- urine catheter

Additional recommended monitoring procedures

- cardiac output calculation (if appropriate)
- intraoperative transesophageal echocardiogram (inflow cannula position, heart valve function, intracardial shunts, volume status)
- right heart function in case of LVAD

Any other monitoring processes can be used (e. g. neurological monitoring) at the anesthesiologist's discretion.

10 Wound care and treatment

Cannula exit sites should be treated like open wounds. The patient's wounds should always be attended to by a small group of nurses in the inpatient area.

The only way to ensure there is a minimum risk of infection is to provide good wound care.

WARNING

Before cleaning the wound (see 10.3: Cleaning of the wound, page 101), put on sterile disposable gloves, cap and mask.

Cleaning the pump and the drive line: Do not use any acetone or petroleum based products near the pump or drivelines. We recommend using only water or alcohol to clean the pump and the drive line.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the blood pump or the drive line as they may alter the surface of the product.

Cleaning the cannulae and transcutaneous exit site: Do not use any acetone or petroleum based products near the cannulae and the transcutaneous exit site.

We recommend using chlorhexidine to clean the cannulae and transcutaneous exit site.

IMPORTANT: Do not use any corrosive or colored solutions or organic solvents to clean the cannulae and the transcutaneous exit site as they may alter the surface of the product.

NOTICE

Do not stick bandages to the cannulae. Over time, remnants of adhesive contaminate the cannulae and increase the risk of infection.

Do not use any adhesive on the velour coating of the cannula as it is difficult to remove and may adversely manipulate the cannula.

Do not use organic solvents near the EXCOR Pediatric such as petroleum ether or turpentine oil, as they could damage the cannulae and the pumps. The plastic parts must not get in contact with chlorinated hydrocarbon (e.g. chloroform), thinners (e.g. acetone, naphtha, toluol, xylene, heptane) or similar compounds.

Do not mark or write on the plastic parts.

Material required (with biventricular access):

- Sterile dressing tray
- Disinfectant i.e. 2% chlorhexidine solution
- Clean gloves
- Mask
- Sterile gloves and towel
- *Metalline*[®] drain compress
- 2X2 gauze, 4X4 gauze
- Adhesive dressing (i.e. *Mepore*[®])
- Adhesive remover
- Non sting barrier film sticks
- Abdominal pads
- Tape
- Tubular bandage (i.e. *Burnnet*)



Fig. 10-1 Materials for dressing change

How often to change the dressing

If the wound is dry and not infected:

- POD 1- once a day
- POD 11-28 every second day, if the wound is dry and not infected
- POD >28 twice a week, if the wound is dry and not infected

If the wound shows signs of infection: clean wound and change dressing twice a day

10.1 Removing the old dressings**► INSTRUCTION**

1. Unpack all the material required to dress the wound and place this within reach on a sterile sheet.
2. Put on disposable gloves, remove old dressings.
3. Take off the disposable gloves, put on the sterile gloves.
4. Remove old dressing using no-touch technique.
5. Examine the places where the cannulae pass through the skin and if changes are apparent take appropriate measures if necessary.
6. Use adhesive remover to remove any adhesive dressing.
Important: adhesive remover (depending on contents) might damage cannula and the pump, use only on skin.

10.2 Cleaning the blood pump

Fig. 10-2 Cleaning the blood pump

► INSTRUCTION

1. Cleanse the exposed cannula and the pump head with disinfectant (i.e. 2% chlorhexidine)

- solution) then place on sterile towel.
2. Observe cannulae and cannulae exit sites.
3. Remove gloves.

10.3 Cleaning of the wound

INSTRUCTION

1. Hand hygiene, prepare sterile dressing tray, put on sterile gloves. If assistance is necessary notify Berlin Heart.
2. 4X4 gauze soaked in 2% chlorhexidine cleanse each cannula exit site in a circular motion outward to a radius of approximately 10 cm.
3. Using a new soaked 4X4 repeat 2 more times beginning at the exit site and clean in larger circles each time.



Fig. 10-3 Cleanse each cannula exit site

4. Wrap 4X4 gauze soaked in 2% chlorhexidine around cannula and gently cleanse with back/forth motion.
5. Repeat with each cannula exit site.
6. Cleanse entire cannula (upper and bottom side).
7. 4X4 gauze soaked in 2% chlorhexidine solution.
8. Starting at the exit site moving down cannula approximately 1 cm from exit site.
9. Repeat for each cannula exit site.
10. Allow chlorhexidine to dry completely.



Fig. 10-4 Cleanse with back/forth motion



Fig. 10-5 Cleanse entire cannula

10.4 The new dressing

10.4.1 Preparing a new dressing

INSTRUCTION

1. Apply non sting barrier film to skin around cannulae. Non sting barrier prevents skin maceration around cannula exit sites.



Fig. 10-6 Non sting barrier film

10.4.2 Applying a new dressing

► INSTRUCTION

1. Wrap a Metalline drain compress around each cannula (from right to left, slit always facing upwards (see figure 10-7).
2. Attach the Metalline drain compresses above the cannulae using sterile bandages. First secure the outer compresses, then the inner compresses (see figure 10-8).



Fig. 10-7 Metalline drain compress



Fig. 10-8 Secure with a sterile bandage

3. Pass a gauze compress folded lengthwise beneath the 2 left cannulae. The open end of the folded compress should point in the direction of the wound. Pull the cannulae into place by tugging the compress slightly (see figure 10-9).
4. Fold the left end of the compress upwards, diagonally to the right and secure with a sterile bandage (see figure 10-9).
5. Fold the right end of the compress upwards, diagonally to the left and secure with a sterile bandage (see figure 10-11).



Fig. 10-9 Gauze compress under the cannulae



Fig. 10-10 Fold the left end of compress and secure



Fig. 10-11 Fold the right end of compress and secure

6. Repeat this procedure for the 2 right cannulae. In this way, the 4 cannulae are padded so that they do not press on the skin or wound (see figure 10-12).

7. Cover the entire wound broadly with gauze compresses (see fig. 10-13).



Fig. 10-12 Cannulae are padded



Fig. 10-13 Cover with sterile gauze compresses

8. Secure the upper part of the dressing with a sterile bandage (see figure 10-14).
 9. Finally, seal the dressing at the left and right side, below the cannulae and between the individual cannulae with strips of adhesive bandage (e. g. Leukoplast), see figure 10-15.



Fig. 10-14 Secure with a sterile bandage



Fig. 10-15 Seal with strips of adhesive bandage

10. Place tubular bandage (i.e. *Burnnet*) around patient (see figure 10-16).
 11. Tie in front to secure dressing.



Fig. 10-16 Tubular bandage

10.5 Regular checks of blood pump(s) and cannulae

Frequency of inspection: every 4 hours



Everyone involved in caring for an EXCOR patient must be trained to carry out a visual check, to evaluate the filling behavior of the blood pump(s) and to detect deposits.



At least daily, the EXCOR cannulae should be inspected for signs of wear or damage. ADVICE: To avoid needless kinking of the cannulae use a mirror for inspection of the bottom side of the blood pump.

At least every 4 hours, check visually that the blood pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, then take the appropriate corrective action.

Under certain circumstances, the message **left/right pump is not filling adequately** in some circumstances is not generated with the 10 ml EXCOR blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

10.5.1 Visual inspection: pump filling and ejection

The filling and ejection behavior of a blood pump is optimal when the membrane surface is completely smooth at the end-of-systole and end-of-diastole positions. Check visually that the pump(s) is (are) filling and ejecting completely over a period of several pump cycles. If a pump is not filling and/ or ejecting completely, take the appropriate corrective action.

Cautionary measures

For all blood pumps: check the position and condition of the driving tube and the cannulae (inflow deterioration due to kinks in cannulae/driving tubes is rather rare).

For all blood pumps: check the membrane movement.

Medical examination of patient

Check CVP, mean arterial pressure and adjust therapy if necessary.

Check the volume status:

- amount of bleeding
- increased urine output (use of diuretics?)
- tamponade
- Important: Increasing the suction pressure will not bring about any distinct improvement if there is not sufficient volume available.

LVAD: observe the functions of the right ventricle.

Adjusting the parameter values

Only adjust the parameters if the measures listed above have no effect or in case of:

- Mobilization of patient: adjust the systolic pressure, both left and right. When pressures have increased, do not reduce these again, even when the patient is lying down.
- Signs of low cardiac output: the membrane is moving properly while at the same time a decrease in urine output, lactate increase and dyspnea (shortage of breath) can be observed. In this case, increase the rate and adjust other settings as required.

► INSTRUCTION

1. Use the <<->/<->> keys to move the cursor to the desired field in the parameter table. The selected field is given a colored background.
2. Use the <↓>,<↑> and <Bild↓>, <Bild↑> keys to adjust the value, then press <Enter> to confirm the input. The system will now operate using the new settings.

Cautionary measure

Confirm each changed parameter value by pressing <Enter>. The system does not take over the new, changed value until it has been confirmed with <Enter>.


10.5.2 Visual inspection: deposits

ADVICE

Enter all the changes to the parameter values into the parameter log.
(see section 15.7: Sample copy: EXCOR parameter log, page 172).

Check the blood pump(s) for visible deposits (fibrin, clots) every 4 hours. If deposits develop, check the pump(s) every hour.

Checking the pump areas which come in contact with blood



- 1 transition inflow cannula – inflow connector
- 2 inflow stub in front of inflow valve (only for pumps with PU valves)
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 (only for pumps with PU valves) outflow stub behind outflow valve
- 11 transition outflow connector – outflow cannula

Fig. 10-17 Diagram of EXCOR blood pump (top view of blood chamber)

ADVICE

During the visual check, first clean the blood pump then illuminate the blood chamber with a flashlight. This makes it easier to detect deposits. Enter all of the findings into the blood pump log. (see section 15.6: Sample copy: EXCOR pump log, page 170).

Cautionary measures

Initial signs of deposits: check anticoagulation therapy and adjust therapy if necessary.
Floating deposits inside the pump: replace the pump!

10.5.3 Checks using the monitor program

Record all drive parameters and adjust if necessary.
Objective: the blood pump(s) must fill and eject completely in each pumping cycle, the diastolic pressure should be as low as possible.

ADVICE

Record the parameter values once a day.

To record the parameters use the sample copy in section 15.7: Sample copy: EXCOR parameter log, page 172.

10.5.4 Replacing the blood pump due to growth of the patient

CAUTION

In children, plan to replace the pump(s) with a larger pump(s) in good time, to prevent the possibility of inadequate support due to an insufficient discharge rate.

The pump selected at the time of transplantation may not be adequate for the entire period of cardiac support. Growth and/or weight gain can result in the patient not receiving adequate support. Use the chart in section 15.1.3: Overview: relationship of body weight – pump size, page 122, to plan, in good time, which pump(s) the patient may need to change over to. This chart is for guideline purposes only and is not binding for each individual case. This decision must be taken by the surgeon in consultation with Berlin Heart GmbH.



Notify Berlin Heart! 866.249.0128

The blood pump(s) must be replaced as described in section 14.1: Replacing the blood pump(s), page 146.

11 Anticoagulation therapy

11.1 Before Implantation of the EXCOR

11.1.1 General considerations

Patients with an EXCOR system must be maintained on anticoagulation therapy.

Anti-Xa levels should be specific to the drug being used, either unfractionated heparin or enoxaparin.

The TEG® may be useful in managing unfractionated heparin and antiplatelet therapy. Please contact *Berlin Heart*, Clinical Affairs for further information.

11.1.2 Pre implantation

The following laboratory tests should be considered prior to implantation.

- Platelet Function Studies, INR, PTT, fibrinogen, antithrombin III, and platelet count to establish a baseline. Assessment for thrombophilia by measuring Protein C, S, Factor V Leiden, Prothrombin 20210 defect, as well as Heparin Induced Thrombocytopenia (HIT) is recommended.

11.2 During Implantation - Cardiopulmonary Bypass

11.2.1 Cardiopulmonary Bypass (CPB)

Use unfractionated heparin as per institutional protocol for cardiopulmonary bypass.

11.2.2 Post CPB

Completely reverse heparin with protamine sulphate as per institutional protocol.

The goal post-CPB is to achieve normal (institution specific) coagulation parameters (INR, PTT, fibrinogen, platelet count).

In the early post-operative period, the possibility of surgical bleeding, GI bleeding, internal bleeding in the retro-peritoneum or other bleeding diathesis is possible and must be monitored.

If the patient is bleeding despite normal coagulation parameters consider:

- Von Willebrand's
- Surgical bleeding

11.3 Postoperative anticoagulation therapy

11.3.1 General Considerations

Primary tests used to evaluate anticoagulation in the patient include antifactor Xa levels and/or PTT.

11.3.2 Starting anticoagulation therapy

During the first 24 hours following implantation, no anticoagulants should be administered.

Approximately 24 - 48 hours after implantation, commence unfractionated heparin therapy (i.v.) if the following criteria are met:

- Platelet count >20,000/ μ l
- Normal Platelet Function Studies
- Minimal bleeding in infants and young children.

11.3.3 Unfractionated heparin therapy (i.v.) Patient < 12 months

- Initial dose 15 IU/kg/hour.
- Do not use a bolus
- After 6 hours if the patient does not have increased bleeding, increase the heparin infusion to 28 IU/kg/hour (therapeutic dose).

6 hours after increasing the heparin to the therapeutic dose, obtain a PTT and an antifactor Xa level.

If the anti factor Xa level is desired range (0.35-0.5 U/ml) and the PTT is in the therapeutic range (institution dependent), then either the PTT or anti factor Xa level may be used to follow the heparin therapy.

If the anti factor Xa level is <0.35 U/ml or >0.5 U/ml, increase or decrease the heparin infusion, respectively until the anti factor Xa level is the therapeutic range (see Tab. 11-1: Unfractionated Heparin adjusted to maintain an anti factor Xa level of 0.35 to 0.50 U/ml., page 111).

Anti factor Xa levels should be obtained daily. Important: hyperbilirubinemia may result in falsely low anti factor Xa levels. If anti Xa levels do not correlate with the PTT in this setting, consider using the PTT to monitor heparin therapy.

Antithrombin should be >70%. If the antithrombin is <70%, treat according to institutional protocol.

11.3.4 Unfractionated heparin therapy (i.v.) Patient ≥ 12 months

Initial dose 10 IU/kg/hour.

Do not use a bolus.

After 6 hours if the patient does not have increased bleeding, increase the heparin infusion to 20 IU/kg/hour (therapeutic dose).

6 hours after increasing the heparin to the therapeutic dose, obtain a PTT and an anti factor Xa level.

If the anti factor Xa level is desired range (0.35-0.5 U/ml) and the PTT is in the therapeutic range (institution dependent), then either the PTT or anti factor Xa level may be used to follow the heparin therapy.

If the anti factor Xa level is < 0.35 U/ml or > 0.5 U/ml, increase or decrease the heparin infusion, respectively until the anti factor Xa level is the therapeutic range (see Tab. 11-1: Unfractionated Heparin adjusted to maintain an anti factor Xa level of 0.35 to 0.50 U/ml., page 111).

Anti factor Xa levels should be obtained daily. Important: hyperbilirubinemia may result in falsely low anti factor Xa levels. If anti Xa levels do not correlate with PTT in this setting, consider using the PTT to monitor heparin therapy.

Antithrombin should be >70%. If the antithrombin is <70%, treat according to institutional protocol.

NOTICE

If during standard unfractionated heparin therapy:

1. Platelet count is < 40,000/μl revert to the Stage I heparin dose for continuous infusion (see Tab. 11-1, page 111)
2. Platelets <20,000/ul discontinue heparin and consider evaluation for heparin induced thrombocytopenia (HIT).

If the anti factor Xa or PTT is too low or too high during heparin therapy, never use a bolus of heparin or protamine. Instead, increase or decrease the heparin dose, IU/hour, as required (see Tab. 11-1, page 111).

11.3.5 Thrombelastography (TEG®)

TEG® analysis may be useful in managing the anticoagulation and anti-platelet therapy. Please contact *Berlin Heart Inc.*, Clinical Affairs for further information.

11.4 Low Molecular Weight Heparin

At 48 hours following surgery if all bleeding has stopped, the creatinine is within normal limits, and the patient is hemodynamically stable, switching from unfractionated heparin to low molecular weight heparin (LMWH) is recommended.

- Patient < 3 months start administration of Enoxaparin at 1.5 mg/kg subcutaneously every 12 hours.
- Patient > 3 months start administration of Enoxaparin at 1 mg/kg subcutaneously every 12 hours.
- Stop heparin infusion and administer LMWH (subcutaneously) simultaneously.
- Obtain the first anti factor Xa level at 4 hours after the 2nd LMWH dose is administered. See Tab. 11-2, page 111 for monitoring and dosing.
- Anti factor Xa therapeutic range: 0.6 to 1.0 U/ml.
- Anti factor Xa should be monitored along with platelet count, and creatinine
- When using LMWH, monitor Anti factor Xa daily. Once the Anti Factor Xa level is in the therapeutic range at a stable dose, monitor twice a week for 2 weeks, and then weekly.

11.5 Oral Anticoagulation Therapy (only for patients ≥ 12 months of age who are taking a full oral diet)



This section only applies to patients ≥ 12 months. Oral anticoagulation in children < 12 months of age is not recommended due to difficulties with monitoring the warfarin effect.

When the patient's condition has been fully stabilized (e.g. hemodynamically stable, no evidence of bleeding, etc), switch to oral anti-coagulation therapy with a vitamin K antagonist (target INR: 2.7 to 3.5), with an initial loading dose of 0.2 mg/kg/day. Do not exceed maximum loading dose of 5mg/day. The INR must be checked daily in the first 4 weeks, twice a week for the next 4 weeks (if INR is stable), and once a week thereafter (see Tab. 11-3: Warfarin loading dose to maintain an INR of 2.7 - 3.5 (6Monagle, P, et al.), page 112 and Tab. 11-4: Warfarin Maintenance Dosing for Day 5 and longer to maintain INR 2.7-3.5, page 112).

Until the target INR is achieved, simultaneous administration of warfarin and heparin is necessary (approximately 4 days). Once the target INR is achieved, heparin therapy can be discontinued. If the INR decreases to = 2.7, administer LMW heparin immediately and then q12h until an INR of > 2.7 is achieved. (Table 5, Appendix 2) If INR is 2.0- 2.7 use an enoxaparin dose of 0.5 mg/kg targeting an anti factor Xa level of 0.3-0.5, if INR is <2.0 use an enoxaparin dose of 1 mg/kg targeting an anti factor Xa level of 0.5 - 1.0.

When unable to achieve a stable INR with warfarin, LMWH should be used instead. Discontinue the warfarin and administer LMWH as per previously discussed age related dosing (see Tab. 11-2, page 111).

11.6 Monitoring of Blood Count and Anticoagulation Status

Monitoring the anticoagulation status as well as infection risk, and renal and hepatic function is important and should be monitored with the following frequency:

- Daily while on UFH, twice a week while on enoxaparin/coumadin for 4 weeks then once week : Fibrinogen, D-dimer, aPTT, PT/INR, Platelet Count, TEG®, Antithrombin, WBC, HgB, HCT, BUN/SCr, AST/ALT, bilirubin T/D, prealbumin, CRP.
- While on UFH obtain anti factor Xa level daily.

- While on enoxaparin obtain anti factor Xa daily until in therapeutic range and on a stable dose, then twice a week for two weeks and then weekly.

If infection is suspected, appropriate measures must be taken immediately (antibiotic therapy, adjustment of the anticoagulation and platelet inhibition therapy) and increased monitoring of the coagulation system. In addition, in the setting of hemodynamic instability, organ dysfunction, and inadequate anticoagulation daily monitoring should be performed until any of these issues are resolved

11.7 Postoperative platelet inhibition therapy

As individual patient responses vary to the anti-platelet agents, the optimum dosage for each patient will be that which minimizes both the risk of thromboembolic complications when the dose is too low and the risk of hemorrhagic complications when the dose is too high. Acetylsalicylic acid (ASA) and dipyridamole are the anti-platelet agents recommended.

11.7.1 Start of therapy

Dipyridamole

At 48 hours after surgery, start dipyridamole, 4mg/kg/day p.o. divided into 4 doses (1 mg/kg Q6) (maximum dose 15mg/kg/day). If the following are present:

- All bleeding has stopped, AND
- The patient is hemodynamically stable AND,
- Platelet studies do not show significantly decreased function,
- Platelet count is > 40,000/ μ l,

Acetylsalicylic Acid

At 4 days post implantation, following the removal of all drainage tubes, start acetylsalicylic acid (ASA) 1mg/kg/day p.o., divided into 2 doses (0.5 mg/kg Q 12), if the following are present:

- Platelet studies show platelet inhibition in the presence of AA < 70%

The ASA dose should split and be administered two times daily (0.5 mg/kg Q 12) due to the short half life and the high turnover of the platelets (approximately 10 % new platelets per day).

11.8 Adjunctive Medication

The inflammation parameters (Tissue factor pathway inhibitor, prothrombin fragment 1-2, fibrinogen, Factor VIII) for patients on ventricular assist device support are often elevated above normal. Accordingly, the physician may choose to administer the following medications at his/her discretion to facilitate the overall anticoagulation/anti-platelet management of the patient:

- Omega-3 fatty acids (e.g. DHA/EPA), have been shown to have an anti-inflammatory effect and also decrease premature activation of platelet membrane. Omega-3-fatty acids are composed of long chain polyunsaturated long chain carbons. Only alpha-linolenic acid (ALA) of the omega-3 family is truly essential.

Antioxidants (Vitamin C and E) also have been shown to have an anti-inflammatory effect, and may be considered.

11.9 Anticoagulation Therapy

11.9.1 Therapeutic Heparin administration and adjustment

NOTICE

This table assumes the site therapeutic PTT is 60 to 85 seconds (⁶Monagle, P, et al.). Each site should use their hospital calculated therapeutic range.

Stage	Description	Anti factor Xa [u/ml]/PTT	Infusion	Hold heparin	Rate Change [%]	Repeat PTT
I	Initial Dose (first 6 hours)					
	Infant < 12 mo		15 IU/kg			
	Child ≥12mo		10 IU/kg			
II	Therapeutic Dose					
	Infant < 12 mo		28 IU/kg/h			after 6h
	Child ≥12mo		20 IU/kg/h			after 6h
III	Adjustment					
		<0.1/<50	0	0	+15%	4h
		0.1-0.34/50-60	0	0	+10%	6h
		0.35-0.50/60-85	0	0	0	next day
		0.51-0.70/86-95			-10%	6h
		0.71-0.89/96-120		30 min.	-10%	4h
		= 0.90/ >120		60 min.	-15%	4h

Tab. 11-1 Unfractionated Heparin adjusted to maintain an anti factor Xa level of 0.35 to 0.50 U/ml.

Anti Factor Xa level U/ml?	Hold Next Dose?	Dose Change?	Repeat Anti Factor Xa?
< 0.35	no	increase dose by 25%	4h after next dose
0.36 - 0.45	no	increase dose by 15%	4h after next dose
0.46 - 0.59	no	increase dose by 10%	4h after next dose
0.6 - 1.0	no	no	4h after next dose
1.1 - 1.25	no	decrease dose by 20%	4h after next dose
1.26 - 1.5	no	decrease dose by 30%	4h after next dose
1.6 - 2.0	yes for 3h	decrease dose by 40%	Before next dose then 4h after next dose
> 2.0	yes, until anti factor Xa level is <0.5 U/ml	decrease dose by 50%	Before next dose is administered, if >0.5U/ml (therapeutic level), do not give next enoxaparin dose & repeat anti Xa level in 12h. When level <0.5 U/ml, administer 50% original dose.

Tab. 11-2 Enoxaparin, low molecular weight heparin dosing (⁶Monagle, P, et al.)

11.9.2 Oral Anticoagulation Therapy

Stage	INR	Action
Day 1	1.0 - 1.8	0.2 mg/kg orally
Day 2-4	1.1 - 1.3	repeat day 1 loading dose
	1.4 - 1.9	50% of day 1 loading dose
	2.0 - 3.0	50% of day 1 loading dose
	3.1 - 3.5	25 % of day 1 loading dose
	> 3.5	hold dosing until INR is < 3.5

Tab. 11-3 Warfarin loading dose to maintain an INR of 2.7 - 3.5 (⁶Monagle, P, et al.)

Stage	INR	Action
Maintenance : = Day 5 and long term	1.1 - 1.9	increase dose by 40 -50%
	2 - 2.4	increase dose by 10%
	2.7 - 3.5	no change
	3.6 - 4.0	administer next dose at 50% then restart at 20% less maintenance dose
	4.1- 5.0	hold one dose then 20% less maintenance dose

Tab. 11-4 Warfarin Maintenance Dosing for Day 5 and longer to maintain INR 2.7-3.5

INR 2.7 to 3.5	use only warfarin p.o.
INR = 2.7	use warfarin plus enoxaparin as outlined in section 5 until INR = 2.7

Tab. 11-5 Drugs and Dose for specific INR range

12 Weaning and Explantation for BTR and BTT

12.1 Weaning Procedure

12.1.1 Introduction

This document summarizes the clinical guideline for weaning and explantation of the EXCOR. The decision to wean the EXCOR should be made cautiously after careful review of all available clinical and laboratory data. This document should be considered a guideline only. As always treatment must be individualized to each patient based on his/her unique clinical circumstances.

It is important to recognize that prolonged pump stoppage and operation of the device at lower beat rates is not recommended due to the risks of blood stagnation and thrombus formation. This risk increases with the smaller blood pumps (e.g. 10, 25 and 30 ml devices) where the luminal sizes and flow rates are the lowest. Therefore, a size-based guideline has been developed to test the adequacy of the native circulation without a prolonged pump stoppage using a combination of gradual weaning, brief pump stoppages, careful anticoagulation monitoring, invasive hemodynamic testing, and a brief afterload challenge. It is not recommended that weaning proceed unless all parameters especially those pertaining to anticoagulation have been fully optimized. This protocol reflects the most recent understanding of the safest possible weaning strategy based on the collective US and European experience to date. Consultation with *Berlin Heart, Inc.* prior to weaning and explantation is strongly recommended.

12.1.2 Indication

Weaning may be considered in children supported with the EXCOR judged to have sufficient evidence of myocardial recovery to provide adequate systemic perfusion independent of VAD support.

12.1.3 Eligibility Criteria



Continuous reassessment of eligibility criteria is critical to reducing the risks associated with weaning of VAD support. At all times each of the weaning criteria must be satisfied in order to proceed with the weaning protocol.

Special attention must be taken to ensure the patient's anticoagulation status remains within the targeted range.

Weaning of the EXCOR may be considered in subjects who meet the following eligibility criteria:

- LVEDD within normal limits (<98th percentile, or Z-score of +2)
- EF = 45% (i.e. no less than mild dysfunction)
- Lactate <3 mmol/L
- No clinical evidence of thromboembolism or bleeding
- Anticoagulation markers within target parameters

12.1.4 Weaning Protocol



Rates < 60 bpm are intended to be used only for implantation and explantation. Never use the *Ikus* with a rate < 60 bpm without constant supervision.

If the patient does not meet the eligibility criteria at any time during the weaning process: Resume pumping at rate prior to any weaning (initial rate, IR).

The weaning protocol can be divided into 5 steps and generally takes one week to complete.

- Day 0 (and throughout the weaning process). Confirmation of eligibility criteria for weaning.
- Day 0. Acute weaning challenge
- Day 1-4. Graduated weaning challenge with non-invasive assessment (echo).
- Day 5. Pump stoppage with invasive hemodynamic assessment with afterload challenge.
- Day 6. Pump stoppage with invasive hemodynamic assessment in OR (full anticoagulation).

This size-based weaning protocol accounts for physiologic differences in heart rate and stroke volume observed in children of varying ages.

12.1.5 10 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	50 bpm
total weaning interval	Difference between initial rate and explantation rate: TWI = IR – WR	TWI	Please enter: IR ____ bpm – WR 50 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table 12-2.

Tab. 12-1 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 50\ bpm + 0.75 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 50\ bpm + 0.50 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 50\ bpm + 0.25 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 12-2 Reduced rate day 1 to day 3

10 ml pump Weaning Sequence

10 ml pump Weaning Sequence	
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 3-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After a total time of 10 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After a total time of 20 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

10 ml pump Weaning Sequence	
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. rattle, clapping) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 3-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. rattle, clapping) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 3-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (50 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

10 ml pump Weaning Sequence		
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. After 5 minutes, reduce the pump rate step-wise from WR (50 bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 4. After a total time of 10 min at 30 bpm, stop the pump for 3 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 3 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 3 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 6-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

12.1.6 25/ 30 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	40 bpm
total weaning interval	Difference between initial rate and explanation rate: $TWI = IR - WR$	TWI	Please enter: IR ____ bpm – WR 40 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table 12-4.

Tab. 12-3 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 40\ bpm + 0.75 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 40\ bpm + 0.50 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 40\ bpm + 0.25 \times TWI\ (___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 12-4 Reduced rate day 1 to day 3

25/ 30 ml pump Weaning Sequence

25/ 30 ml pump Weaning Sequence	
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 5-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After a total time of 10 min at 30 bpm, stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 5-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. After a total time of 20 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 4. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

25/ 30 ml pump Weaning Sequence	
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (___ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. patty cake) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (___ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (___ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. patty cake) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 15-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (40 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

25/ 30 ml pump Weaning Sequence		
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. After 5 minutes, reduce the pump rate step-wise from WR (50 bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. 4. After a total time of 30 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 15 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 5 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 20-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR. 	<p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p> <p>_____</p>
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p>	<p><input type="checkbox"/> NO - STOP</p> <p><input type="checkbox"/> YES - Proceed</p> <p>MD: _____</p>

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

12.1.7 50/ 60 ml pump

The individual weaning progress is based upon the following parameters:

Parameter	Explanation	Abbr.	Value
initial rate	rate prior to any weaning	IR	Please enter: IR = _____ bpm
weaning rate	lowest rate achieved during weaning process, depends on pump size	WR	30 bpm
total weaning interval	Difference between initial rate and explanation rate: $TWI = IR - WR$	TWI	Please enter: IR ____ bpm – WR 30 bpm = TWI ____ bpm
reduced rate	rate resumed at the end of day 1 to 3	RR ₁ to RR ₃	Please refer to table 12-6.

Tab. 12-5 Important parameters for weaning progress

Reduced rate (RR _x)	Calculation
RR ₁	Please enter: $RR_1 = WR\ 30\ bpm + 0.75 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₂	Please enter: $RR_2 = WR\ 30\ bpm + 0.50 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$
RR ₃	Please enter: $RR_3 = WR\ 30\ bpm + 0.25 \times TWI\ (\ ___ \text{ bpm}) = ____ \text{ bpm}$

Tab. 12-6 Reduced rate day 1 to day 3

50/ 60 ml pump Weaning Sequence

50/ 60 ml pump Weaning Sequence	
Day 0	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer unfractionated heparin (UFH) 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from IR (_____ bpm) to 30 bpm in increments of 5 bpm q5 min. After 5 minutes at 30 bpm, reassess LV size and function. After an additional 5 minutes (i.e. total time = 10 min at 30 bpm), stop the pump for 5 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 5-minute pump stop, reconnect pump to Ikus and resume pump speed at IR(____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD:_____</p>
Day 1	<p>After confirmation of eligibility criteria, the following steps should be performed sequentially under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise by from the IR (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 10 min at 30 bpm, stop the pump for 10 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 10-minute pump stop, reconnect pump to Ikus and resume pumping at rate RR₁ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD:_____</p>
Day 2	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. After 5 minutes, reduce the pump rate step-wise from RR₁ (_____ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. After a total time of 15 min at 30 bpm, stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. After 15-minute pump stop, reconnect pump to Ikus and resume pumping at RR₂ (____ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD:_____</p>

50/ 60 ml pump Weaning Sequence	
Day 3	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce the pump rate step-wise from RR₂ (___ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. ambulate) as clinically appropriate, where possible 4. After a total time of 30 min at 30 bpm, stop the pump for 20 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After 20-minute pump stop, reconnect pump to Ikus and resume pumping at RR₃ (___ bpm).
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>
Day 4	<p>After confirmation of eligibility criteria, the following steps should be performed under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 2. After 5 minutes, reduce pump rate step-wise from RR₃ (___ bpm) to 30 bpm in increments of 5 bpm q 5 min. After 5 minutes at 30 bpm, reassess LV size and function. 3. Initiate exercise with gentle age-appropriate play tasks (e.g. ambulate) as clinically appropriate, where possible. 4. After a total time of 30 min at 30 bpm, stop the pump for 30 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while Ikus is disconnected. 5. After a 30-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to Ikus and resume pumping at WR (30 bpm). If the patient does not meet all criteria, reconnect Ikus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

50/ 60 ml pump Weaning Sequence	
Day 5	<p>After confirmation of eligibility criteria, the following steps should be performed in the cath lab under echo guidance ¹:</p> <ol style="list-style-type: none"> 1. Obtain standard access for RHC (if possible with out sedation). 2. Administer UFH 75 units/kg x _____ kg = _____ mg IV x 1 [max 5000 units]. 3. Assess LV size and function to obtain data for comparison. 4. Stop the pump for 15 min and reassess LV size and function. During pump stoppage, use the manual pump to pump twice q30 seconds to minimize the risk of thrombus formation due to blood stagnation while lkus is disconnected. 5. After 15 minutes, initiate norepinephrine infusion at 0.01 mcg/kg/min IV gtt titrated to MAP 20% above baseline x 5 min. While doing so, proceed pumping manually twice q30 seconds. 6. If LV size and function acceptable, proceed pumping manually twice q30 seconds for 15 min. While doing so, reassess LV size & function, and record RAP, PAP, PCWP and MVS. 7. After 30-minute pump stop: If the patient meets all eligibility criteria, reconnect pump to lkus and resume pumping at 50 bpm until the actual surgical procedure of explantation takes place. If the patient does not meet all criteria, reconnect lkus and resume pumping at IR.
	<p>Did the patient satisfy all 5 eligibility criteria throughout the period of weaning and pump stoppage?</p> <p><input type="checkbox"/> NO - STOP <input type="checkbox"/> YES - Proceed MD: _____</p>

¹ TTE unless echo windows insufficient. The last weaning increment may be less than 5 bpm if the wean interval is not a multiple of 5.

12.1.8 Explantation Criteria

NOTICE

ASA and dipyridamole should be discontinued 24-hours prior to device explantation; coumadin/Enoxaparin should be transitioned back to unfractionated heparin (titrated to therapeutic levels).

Milrinone 0.75 µg/kg/min should be started 12 hours prior explantation. ACE inhibitor, β-Blocker and Spirinolactone should be not stopped.

In the operating room, explantation should be considered if the following criteria are met with the pump stopped for 20 minutes (after anticoagulation has been established in the target range for cardiopulmonary bypass):

- LVEDD less than 98th percentile (Z-score less than +2)
- EF ≥45% (i.e. no more than mild ventricular dysfunction)
- Normotensive on only Milrinone (no other inotropes)
- Lactate <3 mmol/L
- LVEDP < 12 mm Hg
- Resting CI of > 2.8 L/min/m²

Surgery should be performed without Cardiopulmonary Bypass. Control all bleeding immediately during and post implantation.

12.2 Explantation for BTR

12.2.1 Explantation with univentricular support

The procedure is analogous to that used after BTT (see 12.3: Explantation for BTT, page 127). Sew over all anastomosis areas where cannulae were placed.

12.2.2 Explantation after biventricular support

Stopping the right pump

INSTRUCTION

1. Select **Pause right** (see figure 12-1), then press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The right pump will stop. The view *Pump size and single-step mode* is shown (see figure 6-7, page 71). The cursor is located on the **OK** field.

Parameter	Operation	Pressure [mmHg]		Rate	% Systole
	Normal	Systole	Diastole		
Left	L	200.0		0.0	40.0
Right	R	170.0		0.0	40.0

Buttons at the bottom: Alarm off, L/R separate, **Pause right** (highlighted), Log off.

Fig. 12-1 Pause right

2. Unplug the driving tube of the right pump from the connector on the *Ikus*. Use the seal plug to seal the connector.
3. To confirm the **OK** selection, press <Enter>. The *Ikus* continues running. The screen shows the standard view.

Switching the *Ikus* off

WARNING

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [O] position!. Otherwise there is a risk that the drive may fail in future due to the *Ikus* batteries being totally discharged.

Always follow the above sequence of operations. Always use the key switch to switch off the *Ikus*.

Do not switch the *Ikus* off unless the batteries are fully charged. Leave the *Ikus* switched on until all yellow LEDs light up, then switch off the *Ikus* with main switch (key switch).

Keep all driving tube connectors covered at all times when not in use.

INSTRUCTION

1. Put the patient on cardiopulmonary bypass (CPB).
2. Disconnect the driving tubes and connect both tank units to the *Ikus*.
3. Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.
4. Next in the monitor program, select the option **Drive OFF** (see figure 12-2, page 127) and press <Enter> to confirm.
5. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops operation immediately and writes an operating log.
6. Disconnect the driving tube(s) from the connector(s). To do so, take hold of the plug's release sleeve and pull the plug out of the connector.

7. Use the seal plugs to seal the driving tube connector sockets.
8. Wait until the log has been completed. When the message **Switch drive off with main switch!** appears, press <F10> to shut down the monitor program. Confirm by pressing the <X> key or the <1> key.
9. Select **3. End** (<3>, see figure 12-3, page 128) in the start menu and switch off the laptop.
10. Switch the *Ikus* off, provided that the batteries are fully charged. To do so, turn the key switch to [0] position.

12.3 Explantation for BTT

NOTICE When planning and timing the transplantation, be aware that massive adhesions may exist in the transplant recipient.

Preparing the donor organ

ADVICE Leave adequate lengths of the aorta and the pulmonary artery attached to the donor organ in order to be able to continue using those parts of the original vessels used for anastomosis of the VAD cannulae.

Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.

Switching the *Ikus* off

WARNING The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [0] position! Otherwise there is a risk that the drive may fail due to the *Ikus* batteries being totally discharged.

CAUTION Always follow the above sequence of operations. Always use the key switch to switch off the *Ikus*.

Do not switch the *Ikus* off unless the batteries are fully charged. To do this leave the *Ikus* switched on until all yellow LEDs light up, then switch the *Ikus* off using the key switch.

INSTRUCTION

1. Put the patient on cardiopulmonary bypass.
2. Disconnect the driving tubes and connect both tank units to the *Ikus*.
3. Leave the *Ikus* running with the tank units until the patient is stable on CPB and the blood pumps have been explanted.
4. Next in the monitor program, select the **Drive OFF** option and press <Enter> to confirm (see figure 12-2).




Fig. 12-2 Drive OFF

5. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops operation immediately and writes an operating log.
6. Disconnect the driving tube(s) from the connector(s). To do so, take hold of the release

- sleeve and pull this out of the connector.
7. Use the seal plugs to seal the driving tube connectors.
 8. Wait until the log has been completed. When the message **Switch drive off with main switch!** appears, press <F10> to shut down the monitor program. Confirm by pressing the <X> key or the <1> key.
 9. Select **3. End** (<3>, see figure 12-3) in the start menu and switch off the laptop.
 10. Switch the *Ikus* off, provided that the batteries are fully charged. To do so, turn the key switch to [0] position.

```
          Berlin Heart <R>
          Ikus2000 <R> Rev. 2.1
          Build 2009.06
          Copyright <C> 1997-2009

1. Start Program
2. Entry codes
3. End
4. Save data
5. Change date or time
6. Change language

Input:
```

Fig. 12-3 Start menu

Removing the VAD cannulae

► INSTRUCTION

1. Clamp off the cannulae.
2. Disconnect the pump from the cannulae.
3. Remove the cannulae. Sew over the anastomosis areas of the atrium.

The remaining procedure is the same as for any primary orthotopic heart transplantation.

13 Error Messages and corrective measures

This chapter describes all *Ikus* error messages and explains what measures should be taken if an error does occur.

Whenever a message is displayed, always follow the exact instructions provided in this instruction for use.

Keep calm!

It is necessary to be very observant as long as the cause of a message has not been corrected.

HOTLINE

Notify Berlin Heart! 866.249.0128

NOTICE

Some errors immediately re-trigger an alarm as long as they are still active even after being acknowledged on the laptop. Also in this case first mute the alarm on the handle in order to prevent permanent retriggering of the acoustic alarm. Before acknowledging the alarm wait for the error message and then take appropriate action.

When an error message occurs, the following happens:

- An acoustic signal (continuous tone) is emitted.
- The indicator light on the control panel of the handle lights up.
- A red border is shown around the field **Alarm off** in the monitor program display.
- In the message window of the monitor program, a text is displayed informing about the time of occurrence, the type of fault and the corrective measures which must be taken. Important: Always observe the instructions! In addition, some of the more complex error messages contain an 8- or 16-digit binary code which enables the service department to identify the exact cause of the fault.

What to do when an error message is shown

INSTRUCTION

1. Check the status of the patient.
2. Observe the filling and ejection behavior of the blood pump (visual check) over several pump cycles!
3. Carry out the appropriate measures and acknowledge the message.




Fig. 13-1 Example of an error message in the message window

Acknowledging an error message

WARNING

Messages are only displayed when the monitor program is running. When the monitor program is shut down, the only indications that there is an error message are an acoustic signal and the fact that the indicator lamp on the *Ikus* handle lights up. There is no way of finding out what type of error message has been displayed.

If the monitor program has been shut down, no data on current events are recorded and the LOG files will be incomplete. Therefore do not shut down the monitor program unless this is absolutely necessary (e. g. if it is necessary to set up new user profiles)! Restart the monitor program as soon as possible after this!

The monitor program shows some of the error messages with an additional reference to system 1, 2 or 3. The system descriptions left and right refer to the internal arrangement of the pneumatic systems and not to the left or right pump.

This refers to the following connections specifically:

- system 1 (left): the system connected to the red connector
- system 2 (right): the system connected to the blue connector
- system 3 (backup): backup system

All messages, together with their time of occurrence, are recorded in the log file. In the case of some messages, the cause of the fault may correct itself automatically after a short time. In this case, a corresponding message is shown (e. g. **Please check left driving tube and pump** or **Left: driving tube/ pump OK**).

INSTRUCTION

1. Press the button on the control panel of the handle to mute the error message. This switches off the acoustic signal temporarily. The **Alarm off** field lights up red in the monitor program.
2. Read the displayed message(s) carefully. Observe the instructions provided in the message text.
3. Take corrective measures immediately!
4. Acknowledge the message in the monitor program. To do so, move the cursor to the **Alarm off** field and then press <Enter> to confirm. Otherwise the acoustic alarm will sound again after 10 minutes at latest. Exceptional cases: see **Please check left pump and driving tube!** or **Please check right pump and driving tube!** in section 13.4, page 132 and **Temperature sensors: <<8-digit binary code>>** in section 13.14, page 137.

13.1 Pressure error / time error in system 1 (or in system 2 or 3)

Pressure fault in system 1 (left), please wait ... or
Pressure fault in system 2 (right), please wait ... or
Pressure fault in system 3 (backup), please wait ... or

Time error system 1 (left), please wait ... or
Time error system 2 (right), please wait ... or
Time error system 3 (backup), please wait ...

The air volume pumped per cycle by the *Ikus* or the pressure required to pump it has changed. The *Ikus* checks whether there is an internal system fault (e. g. a compressor has broken down) or whether there is an external error (e. g. driving tube not connected or driving tube leak).



WARNING

The blood pump of the corresponding system is stopped for the duration of this test (for approx. 10 seconds).

If the *Ikus* has detected a pneumatic system fault...

the *Ikus* switches to backup operation and the message **Backup operation left/right. Contact customer service.** is shown (see section 13.5: Backup operation left/right, page 133), together with a request to inform the service department.



HOTLINE

Notify Berlin Heart! 866.249.0128

If the *Ikus* has detected an external fault ...

it displays the message **Please connect driving tube!** (see section 13.3: Please connect driving tube, page 131). Provide the patient immediately with a replacement *Ikus*.

If the main computer and the backup computer arrive at different results in the test phase ... the backup computer generates the message **Backup computer reports faulty test** (see section 13.11.4: ... reports faulty test, page 136).

Follow the instructions exactly as described in each section.

13.2 Throttle valve error in system 1 (or system 2 or 3)

[Throttle fault in system 1 \(left\), please wait ...](#) or
[Throttle fault in system 2 \(right\), please wait ...](#) or
[Throttle fault in system 3 \(backup\), please wait ...](#)



HOTLINE

Notify Berlin Heart! 866.249.0128

In an active system, the *Ikus* has detected a throttle valve error (throttle: internal module for pneumatics test) and switches to backup operation (see section 13.5: Backup operation left/right, page 133). Provide the patient immediately with a replacement *Ikus*. A throttle is an internal assembly for the pneumatic test.

13.3 Please connect driving tube

[Please connect driving tube!](#)

INSTRUCTION

1. Inspect the driving tube and the connectors.
2. If a plug is not seated correctly: re-insert it correctly. To do so, grip the release sleeve and pull it out of the connector, then plug it back in. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
3. If the driving tube is defective: replace it.

13.3.1 Replacing a driving tube

The driving tubes (red/blue) have a maximum life of 1 year (see section 15.2: Technical specifications, page 163) and must be replaced after 1 year.



WARNING

In order for a driving tube to be replaced, the pump must be stopped for a short time. If the left driving tube is being replaced in a driving unit providing biventricular support, the right pump must also be stopped while the driving tube is being replaced in order to avoid overloading of the pulmonary circulation (danger of pulmonary edema).

Material

- 1 driving tube, red or blue
- 1 tube connecting set (cable tie, cable-tie gun), from accessory set

INSTRUCTION

1. If required log into the monitor program by entering user ID and password and confirm with <Enter>.
2. Select the option **Drive pause** and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The *Ikus* will stop.
3. Carefully cut the cable tie on the defective driving tube.
4. As soon as the pump has stopped, remove the defective driving tube from the pump.
5. Connect the new driving tube to the blood pump. To do so, carefully push the end of the driving tube onto the driving tube connector.

6. Remove the defective driving tube from the *Ikus* connector socket. To do so, take hold of the release sleeve and pull this out of the connector.
7. Connect the new driving tube to the connector, which is now free. The sound of the plug snapping into place is clearly audible. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
8. The view *Select operating mode* appears on the monitor. Select **Univentricular (UVAD)** or **Biventricular (BVAD)**, then confirm the selected operating mode with <Enter>. In biventricular mode the view *Pump size and single-step mode* appears. In univentricular mode, the connector seal test is performed first and then the view *Pump size and single-step mode* appears.
9. To confirm the **OK** selection, press <Enter>.
10. The system starts up again using the defined parameters.
11. Check whether the pump is filling correctly and, if necessary, adjust the *Ikus* parameters.
12. Secure the pump end of the driving tube with a cable tie strap. Important: Only the cable ties and cable tie guns provided should be used. See 8.11, page 94.

13.4 Please check left / right pump and driving tube

Please check left pump and driving tube! or
Please check right pump and driving tube!

The *Ikus* has detected an excessively deviating flow.



WARNING

Rates < 60 bpm are intended to be used only for implantation and explantation. Never use the *Ikus* with a rate < 60 bpm without constant supervision.

Under certain circumstances, the messages **Check left pump and driving tube** or **Please check right pump and driving tube!** are not generated with the 10 ml *EXCOR* blood pump due to the low volume of air which is moved in the pump. Therefore in pumps of this size, pay special attention to the movement of the membrane and ensure that each pump fills and empties completely.

IMPORTANT: First mute the message on the drive unit, do not acknowledge yet in the monitor program. Otherwise, alarms will occur by mistake. When muted, the alarm will be audible again in an 1-minute interval.

INSTRUCTION

1. Mute the message on the drive unit.
2. Inspect the driving tube and the cannulae: are kinks blocking the flow? Correct the positions to ensure an unimpaired flow.
3. Inspect the driving tube and the plugs. If a plug is not positioned correctly: re-insert it correctly. If the driving tube is defective: replace it (see section 13.3.1: Replacing a driving tube, page 131).
4. If necessary, adjust the parameters. Wait until the message **Left: Pump output measurement activated** or **Right: Pump output measurement activated** appears.
5. If necessary, correct the position of the cannulae.
6. Assess the hemodynamic status of the patient (volumes, MAP, PAP, CVP, etc.)

If the message appears again:

INSTRUCTION

1. Mute the message on the drive.
2. Confirm the parameter values. For that the cursor must be in the parameter table.
3. Assess the patient's hemodynamic status (volume, MAP, PAP, CVP...)

4. Monitor the membrane movement and make sure that the pump(s) are filling and ejecting completely. If this is quite correct, confirm the message in the monitor program.
5. If the message appears again, then the parameters can be modified slightly to normalize the flow.
6. Assess the hemodynamic status of the patient (volumes, MAP, PAP, CVP, etc.)
7. Monitor the membrane movement and make sure that the pump(s) are filling and ejecting completely. If this is functioning properly, then acknowledge the message in the monitor program.
8. If the message appears again, check if it is necessary to adjust the cannulae.

If the message appears again:



Notify Berlin Heart! 866.249.0128

13.5 Backup operation left/right

[Backup operation \(right\). Contact customer service](#) or
[Backup operation \(left\). Contact customer service](#)

The left (or right) pneumatic system has failed. The respective pump is now being powered by the backup system.

Provide the patient immediately with a replacement *Ikus*.



Notify Berlin Heart! 866.249.0128

13.6 Error messages in emergency operating mode

13.6.1 UVAD, emergency operation!. Contact service immediately!

[UVAD, emergency operation! Contact service immediately!](#)

Appears in univentricular mode



No more system is available as a redundancy. If the only remaining intact pneumatic system fails, there is a risk that the *Ikus* will stop running altogether.

The blood pump is being driven by the last intact pneumatic system.

The pneumatic system continues running with the currently set parameters. The parameters can still be adjusted if necessary. However, it is not possible to switch over to biventricular operation.

Provide the patient immediately with a replacement *Ikus*.

React and take appropriate measures immediately!



Notify Berlin Heart! 866.249.0128

13.6.2 Emergency operation system 1 (or system 2 or 3). Contact Service now!

[Emergency operation System 1. Contact service now!](#) or
[Emergency operation System 2. Contact service now!](#) or
[Emergency operation System 3. Contact service now!](#)

Appears in **biventricular mode**

Both blood pumps are being driven by the last intact pneumatic system.



There is no longer a redundant backup system. If the only remaining intact system fails, there is a risk that the *Ikus* will stop running altogether.

The pneumatic system will now operate with fixed parameters (synchronous mode, systolic pressure 250 mmHg, diastolic pressure -100 mmHg, 70 bpm, relative systolic duration 40 %). It is not possible to change these settings.

Possible causes:

- 2 of the 3 pneumatic systems have developed faults.
- The *Ikus* driving unit has been running on battery power for too long (Error message: **Emergency operation due to empty batteries: Risk of total failure!**). It was no longer possible to establish reliable conditions. In order to ensure that this is not due to a control circuit defect, the *Ikus* has switched over to the backup system and the backup control computer. The message **Backup computer started! Contact customer service.** is displayed.

Provide the patient immediately with a replacement *Ikus*.

React and take appropriate measures immediately!



Notify Berlin Heart! 866.249.0128

13.7 System 1 (or system 2) is defective!

System 1 (left) is defective! or
System 2 (right) is defective!

The respective pneumatic system has developed a fault. The backup system (system 3) is activated. If this error occurs during backup operation, *Ikus* will be running in emergency operating mode. In this case, *Ikus* has no more redundancy.

Immediately



Provide the patient with a replacement *Ikus*.



Notify Berlin Heart! 866.249.0128

13.8 System 3 (backup) is defective!

System 3 (backup) is defective!

The backup pneumatic system is detected as a fault during backup operation.

The *Ikus* runs in emergency operating mode. Provide the patient immediately with a replacement *Ikus*.

Immediately



1. Provide the patient with a replacement *Ikus*.



Notify Berlin Heart! 866.249.0128

13.9 Alarm circuit fault: buzzer remains off (or on)

Alarm circuit test failed - buzzer remains off! or
Alarm circuit test failed - buzzer remains on!



If this is the case, the *Ikus* will not generate an acoustic signal in an alarm situation, or it will generate a wrong signal. Observe the messages displayed in the message window carefully and look out for the visual alarm signal in the display and operating panel. Do not operate the *Ikus* without supervision!

A fault in the alarm circuit is discovered during the self-test for the alarm circuit or when an alarm situation occurs. Depending on the type of fault, the message appears:

Alarm circuit test failed - buzzer remains off or
Alarm circuit test failed - buzzer remains on or
Acoustic alarm is not properly recognized

Provide the patient immediately with a replacement *Ikus*.



Notify Berlin Heart! 866.249.0128

13.10 Backup computer faulty! Contact customer service.

Backup computer faulty! Contact customer service.



Depending on when the fault occurs, the error message **Processor down. Inform service!** may not appear due to unfavorable resynchronization time response. In this case the drive will provide a visual and acoustic signal, however, no error message will be displayed.

One of the processors has failed. If the main computer fails, the message **Backup computer started! Contact customer service.** also appears. Provide the patient immediately with a replacement *Ikus*.

Immediately...



Notify Berlin Heart! 866.249.0128

13.11 Backup computer

13.11.1 ... reports discrepancy in left/right pump output measurements!

Backup computer: discrepancy in left pump output measurement or
Backup computer: discrepancy in right pump output measurement

The main computer has detected a flow error. This message only appears in conjunction with the message in section 13.4: Please check left / right pump and driving tube, page 132. Immediately take the measures described there.



Notify Berlin Heart! 866.249.0128

13.11.2 ... reports faulty measurements on the left (or right)

Backup computer: reports faulty measurement left! or
Backup computer: reports faulty measurement right!

The main and backup computers show different results.

If this message appears only once,
no further measures are necessary.

If the message appears again:

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of a malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

HOTLINE

Notify Berlin Heart! 866.249.0128

13.11.3 ... reports an error in output measurement on the left (or right) pump

Backup computer reports error: left pump output measurement or
Backup computer reports error: right pump output measurement

The backup computer has detected a flow error, but the main computer has not.

INSTRUCTION

1. Check the filling and emptying of the blood pump(s).
2. Check the plausibility of the measured values.
3. Reset all of the parameter values. In addition, log out of the monitor program and then log back into it. Navigate the cursor with <<-> / <-> to the desired field, adjust the value with <↓>, <↑>/ <Bild ↓>, <Bild-↑>, then confirm with <Enter>. The system works with the new value.

If the message appears again:

HOTLINE

Notify Berlin Heart! 866.249.0128

13.11.4 ... reports faulty test

Backup computer reports faulty test

The main computer cannot end the test phase to check the system after the error message

Pressure fault in system 1 (left), please wait ... or
Pressure fault in system 2 (right), please wait ... or
Pressure fault in system 3 (backup), please wait ... or

Time error system 1 (left), please wait ... or
Time error system 2 (right), please wait ... or
Time error system 3 (backup), please wait ... (see section 13.1, page 130).

IMPORTANT: The blood pump being driven by the tested system stops for approx. 10 seconds during the test phase.

INSTRUCTION

1. Check the status of the patient.
2. If the *Ikus* switches to backup operating mode or this message repeatedly appears, inform Service immediately.
3. Provide the patient immediately with a replacement *Ikus*.



HOTLINE

Notify Berlin Heart! 866.249.0128

13.12 Measurement discrepancy in main computer (backup computer)

[Measurement discrepancy in main computer!](#) or
[Measurement discrepancy in backup computer!](#)



WARNING

Between the error message [Measurement discrepancy in main computer!](#) and the 2nd (real) error message there can be a delay of several seconds. At any case, wait for both error messages.

Only the active processor has detected an error. This message only appears in conjunction with an additional, actual error message. Immediately take all of the necessary measures for this second message.

If the message appears one time at the backup computer,
it is a pure information message.

If the message appears several times
contact the service. Provide the patient immediately with a replacement *Ikus*.



HOTLINE

Notify Berlin Heart! 866.249.0128

13.13 Parameter set update failure

[Parameter set update failure](#)

The *Ikus* cannot store the changed parameter values in the internal system memory. The driving unit is operating with the changed values, but if a reset were to be performed, the old values would be valid again.

INSTRUCTION

1. Check the parameter values regularly. After a reset, if the *Ikus* continues to work with the old parameter values: Re-adjust the parameter values.

If the message appears again:
Provide the patient immediately with a replacement *Ikus*.



HOTLINE

Notify Berlin Heart! 866.249.0128

13.14 Temperature sensors: <<8-digit binary code>>

[Temperature sensors: <<8-digit binary code>>](#)



WARNING

Do not use water or other liquids to cool the *Ikus*! Otherwise there is a risk of short circuit and/or device malfunction.

IMPORTANT: An alarm was triggered by one of the sensors and an user muted it (mute interval: 10 min). During this mute interval now another sensor generates another alarm which is also muted by the user. This means that the remaining mute time of the first triggered alarm will be extended for another full mute interval.

INSTRUCTION

1. Determine whether internal or external influences have caused the driving unit to over-heat. Is it exposed to direct heat from external sources? Is the ambient temperature too high?

Fault: <<16-digit binary code>> (<<type of fault>>)

2. If possible, remedy the situation (move the *Ikus* away from the heater, etc.). Provide adequate ventilation. Acknowledge the message. The *Ikus* takes a few minutes to cool down. Acknowledge the message repeatedly if necessary.

Usually, overheating of the *Ikus* is due to external factors such as direct thermal radiation (e. g. direct sunlight or from heaters). Overheating may also be caused by an internal fault, but this rarely occurs.

If external influences can be excluded as factors causing the message...



Notify Berlin Heart! 866.249.0128

IMPORTANT: When passing on the error code to the service department directly (by telephone or fax): remember to state all 8 digits!

If necessary the service department will request to read out the LOG files and send a copy to Berlin Heart, Inc. (see section 14.7: Reading out the LOG files, page 143).

13.15 Fault: <<16-digit binary code>> (<<type of fault>>)

Fault: <<16-digit binary code>> or
<<type of fault>>



Notify Berlin Heart! 866.249.0128

IMPORTANT: When passing on the error code to the service department directly (by telephone or fax): remember to state all 8 digits!

If necessary the service department will request to read out the LOG files and send a copy to Berlin Heart, Inc. (see section 14.7: Reading out the LOG files, page 143).

13.16 Batteries discharged; battery operation not possible

ABORT: Batteries discharged! No Battery operation possible!

The batteries have become discharged during mains operation by a battery malfunction.



Continue to operate the *Ikus* on mains! There is a risk of total malfunction if battery operation is used!

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).



Notify Berlin Heart! 866.249.0128

If this message appears during the start test

it is a pure information message. By all means continue to operate the *Ikus* on mains! After a sufficient charging time in mains operation the message **Battery charge OK**. appears.

13.17 Insufficient battery charge. Only limited battery operation

Insufficient battery charge. Only limited battery operation.

Keep *Ikus* connected to mains! Danger of total shutdown after brief battery operation.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. If the message **Battery charge OK.** does not appear within 6 hours provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

If the batteries reach a sufficient state of charge the message **Battery charge OK.** appears.

13.18 Error messages - Circuit breaker and internal battery fuse

DANGER:Battery fuse test failed! - Check circuit breaker!

No battery mode possible! Contact service.

or

DANGER:Internal battery fuse test failed!

No battery mode possible! Contact service.

The following applies to both error messages:



Ikus must be kept connected to the power supply! Danger of total shutdown in battery mode!

INSTRUCTION

1. Try to reset the circuit breaker on the connection panel (press button *Circuit breaker*).
2. If the button *Circuit breaker* can be pressed, the message **Battery fuse test OK.** will appear within 10 minutes. After this the error has been eliminated and battery operation is possible again.
3. If the button *Circuit breaker* cannot be pressed or triggers again, the internal battery fuse has failed. In this case the replacement *Ikus* must be connected immediately (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152).

Also see section 14.8: Circuit breaker and battery fuse, page 157.

13.19 Electronic malfunction. Contact customer service!

Electronic malfunction. Contact customer service!

This message appears if the internal power supply of the electronic equipment is faulty.



Danger of total malfunction of the *Ikus* in the event of an additional error! Take immediate action!

INSTRUCTION

1. Provide patient with a replacement *Ikus* and switch off the malfunctioning *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of malfunction of the

defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).



Notify Berlin Heart! 866.249.0128

13.20 Acoustic alarm is not properly recognized

Acoustic alarm is not properly recognized

If the message **Acoustic alarm: OK** appears within 8 seconds after the error message, the *Ikus* is working perfectly. No further measures are necessary.

If the message **Acoustic alarm: OK** does not appear within 8 seconds, the alarm circuit is defective. A possible error might not have been detected.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Important: Do not operate the *Ikus* without supervision!
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of a malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).



Notify Berlin Heart! 866.249.0128

13.21 Error: no data/no reaction from the control computer

Error: no data from Master. or
Error: no reaction from Master.

In the event of the simultaneous malfunction of both control computers or the malfunction of the power supply, the *Ikus* cannot generate a specific error message in the message window.

INSTRUCTION

1. Check the status of the patient.
2. Evaluate the malfunction scenario:
 - no alarm, but the *Ikus* continues to function: no communication between the control computers and the laptop
 - optical and acoustic alarm, and the *Ikus* is working in emergency mode: both control computers are defective
 - The *Ikus* is halted; acoustic alarm only: Defective power supply
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152).
If no replacement *Ikus* is available and the *Ikus* is running: Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
If no replacement *Ikus* is available and the defective *Ikus* is halted: Support the patient with the manual pump. (see section 14.5: Driving blood pump(s) with the manual pump, page 154).



Notify Berlin Heart! 866.249.0128

13.22 Left/right flow sensor defective. Notify Service!

Left flow sensor fault. Contact customer service. or
Right flow sensor fault. Contact customer service.

The corresponding flow sensor is defective. Although the *Ikus* continues to run, the excessively low flow would not be detected - possibly due to a kink in a cannula or the driving tube.

WARNING

Do not operate the *Ikus* without supervision! Otherwise an insufficient support of the patient might not be detected.

NOTICE

In biventricular mode there is a flow alarm for one of the blood pumps. The user stops the pump operation of the other blood pump and restarts it. This affects that during the restart the flow alarm will be deleted, whether its reason continues to exist or not.

IMPORTANT: If the defective *Ikus* is not replaced, this message appears repeatedly at 10-minute intervals.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of a malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

HOTLINE

Notify Berlin Heart! 866.249.0128

13.23 Problem: <<Text>>

Problem: <<Text>>

If this message appears during the start-up test:

see section 13.24: Error messages during the start-up test, page 142.

If this message appears outside of the start-up test,

a serious problem exists.

The *Ikus* has detected a serious problem.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. If possible: provide the patient with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of a malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

HOTLINE

Notify Berlin Heart! 866.249.0128

13.23.1 Self-test is not completed by passive computer!

Self-test is not completed by passive computer!

The passive processor was unable to end the self-test of the alarm circuit.

If, within 8 seconds, the message Alarm circuit test OK appears, the *Ikus* is working perfectly. No further measures are necessary.

If the message Alarm circuit test OK does not appear within 8 seconds, the alarm circuit is defective. A possible error might not have been detected.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and emptying of the blood pump(s), and monitor the function of the *Ikus*. Do not operate the *Ikus* without supervision!
3. Provide the patient immediately with a replacement *Ikus* and switch off the malfunctioning *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152). If no replacement *Ikus* is available: Support the patient, if necessary (in the event of a malfunction of the defective *Ikus*), by means of the manual pump! (see section 14.5: Driving blood pump(s) with the manual pump, page 154).

HOTLINE

Notify Berlin Heart! 866.249.0128

13.24 Error messages during the start-up test

In addition to the error messages listed below there are additional messages possible that are described in chapter 12. Always take those measures corresponding to the indicated messages.

The error messages listed here can only occur during the start test. At any case wait for the end of the start test. Afterwards the corresponding measures are to be initiated.

13.24.1 Battery test skipped (Battery problem!)

Battery test skipped (Battery problem!)

The charge level of the batteries is too low to permit battery operation.

INSTRUCTION

1. Operate the *Ikus* on mains. Battery operation is only possible if all of the yellow LEDs are illuminated.

13.24.2 Additional messages during the start-up test

If the *Ikus* detects an error during the start-up test, one of the following messages appears in the message window depending on the nature of the error:

- Problem: batteries have very different charges.**
- Problem: battery controller: batteries are discharged**
- Problem: Charge unit fault. Contact customer service.**
- Problem: Laptop fault. Contact customer service.**
- Problem: Power relay. Contact customer service.**
- Problem: Mains sensor fault. Contact customer service.**
- Problem: Mains voltage. Check power and switch.**
- Problem: Power pack fault. Contact customer service.**
- Problem: WR2 fault / relay board**
- Problem: Relay 1 faulty. Contact customer service**
- Problem: Relay 2 faulty. Contact customer service**
- Problem: WR1 not switching. Contact customer service.**
- Problem: 0000 0011 0011 0111(fault in power supply),** see section : **IMPORTANT:** When passing on the error code to the service department directly (by telephone or fax): remember to state all 8 digits!, page 138.

Restart the *Ikus* after the start-up test if one of the above messages appears**INSTRUCTION**

1. According to display: If necessary: enter <7> to read any messages; go back by pressing <Enter>. Important: enter <1> or <x> to exit the monitor program. A 2. window appears.
2. 2. window: According to the display, select the option **End**. The monitor program writes a log file and is terminated. Wait until the message **Switch off drive with main switch!** appears.
3. Switch off the main switch (key switch) and then switch it back on again.

If the *Ikus* now ends the start-up test without a message:

Start up the *Ikus*.

If one of the above messages appears again:

Do not start up the *Ikus*. Provide the patient immediately with a replacement *Ikus*.

HOTLINE

Notify Berlin Heart! 866.249.0128

13.25 Discrepancy in pressure measurement: system 1 (or system 2 or 3)

Discrepancy in pressure measurement: system 1 (left) or
Discrepancy in pressure measurement: system 2 (right) or
Discrepancy in pressure measurement: system 3 (backup)

This message is a pure information message.

INSTRUCTION

1. Check the status of the patient.
2. Check the filling and ejection behavior of the blood pump(s). Do not operate the *Ikus* without supervision!
3. Compare nominal values with the values actually generated. In the event of discrepancies from high set values, if necessary reduce rate accordingly and/or vary relative systolic duration, lower systolic driving pressure (e.g. from 210 mmHg to 200 mmHg), and/or raise the diastolic driving pressure parameters (e.g. from -40 mmHg to -45 mmHg) if appropriate.
4. Scroll down message window to check if message has previously appeared.
5. Continue to observe message window to see if message reappears.

The information message **Discrepancy in pressure measurement: system 1 (2 or 3)** appears when the compressed air being produced at a constant rate by the compressor in the present pump cycle is unable to meet the compression levels required to achieve the defined parameters.

The message can have the following causes:

- Extreme parameter values have been set (extremely high values, extremely low values)
- Restricted pneumatic performance (maintenance required)
- System malfunction

If the message appears several times within 24 hrs:**NOTICE**

Damage to the *Ikus* is technically impossible when this message appears!

Avoid a further increase in the rate or, alternatively, the driving pressure for parameters set at high values. Otherwise the *Ikus* might perform a test phase, which means that the blood pump of the corresponding system is briefly stopped for the duration of this test

The information message **Discrepancy in pressure measurement: system 1 (2 or 3)** is not accompanied by an acoustic signal.



Notify Berlin Heart! 866.249.0128

13.26 Communication with Laptop failed

Communication with Laptop failed!

This message informs the user that the communication between the control computers and the monitor program has been temporarily interrupted. This occurs if the monitor program is exited, for example to administrate user IDs and passwords. Since this message only appears when the communication is re-established, this message does not require user intervention.

14 Troubleshooting and correcting faults



Notify Berlin Heart! 866.249.0128

Problem	Cause of problem / action to be taken
Visible <i>Ikus</i> faults	Notify 866.249.0128
Visible blood pump faults	Replace the pump, see section 14.1: Replacing the blood pump(s), page 146.
Deposits in the pump	Initial deposits: check anticoagulation status and adjust therapy if necessary. If floating deposits are detected (may cause thromboembolic complication): replace the pump, see section 14.1: Replacing the blood pump(s), page 146..
Pump is filling or ejecting blood incorrectly	Assess the condition of the patient and the hemodynamic status. If necessary, adjust the system parameters. Check whether it may be necessary to manipulate the cannulae, see section 13.4: Please check left / right pump and driving tube, page 132.
<i>Ikus</i> : the graph display stops moving, parameters cannot be adjusted	Possible causes <ul style="list-style-type: none"> • faulty communications between control computer and laptop • batteries not supplying enough current • the electronics (main and backup control computers) have failed What to do? Switch the laptop off and then back on again, wait for the start-up procedure to be completed, then start the monitor program. Important: The <i>Ikus</i> continues running with the set parameters.
The graphs remain frozen	The <i>Ikus</i> is operating in emergency pulse mode, see 14.3: Emergency pulse mode, page 150. Notify 866.249.0128 Restart the <i>Ikus</i> after consulting the service department staff. See section 14.2: Restarting <i>Ikus</i> , page 149.
Acoustic and visual alarm from the <i>Ikus</i> , message Error: no data from Master or Error: no reaction from Master	Possible causes <ul style="list-style-type: none"> • simultaneous malfunction of both control computers • power supply malfunction What to do? Assess the condition of the patient and the hemodynamic values. Notify 866.249.0128 immediately.

Problem	Cause of problem / action to be taken
Pump stands still - no further pump function; acoustic alarm is still audible	<p>Possible causes</p> <ul style="list-style-type: none"> complete failure of the <i>Ikus</i> batteries are complete emptied or serious fault in the batteries <p>What to do?</p> <p>Immediately connect the <i>Ikus</i> to the mains. Until that supply the patient with the manual pump. See also section 5.4: Switching between mains and battery operation, page 55 and 14.5: Driving blood pump(s) with the manual pump, page 154.</p> <p>Notify Berlin Heart! 866.249.0128</p>

Tab. 14-1 Possible problems

14.1 Replacing the blood pump(s)

WARNING

When replacing a blood pump, follow the instruction given here. Otherwise the duration of the pump stop will be prolonged and the patient might suffer from inadequate support.

The blood pump may only be replaced under sterile conditions!

When connecting the blood pump(s), pay attention to the direction of the arrows on the inflow and outflow stubs! These show the direction of the blood flow.

The cable tie covering the *EXCOR* cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula.

BVAD: If the left pump is being replaced, the right pump must also be stopped while the pump is being replaced. Otherwise there is the risk of pulmonary edema.

NOTICE

If the replacement pump has a larger volume than the one being replaced,

- the use of a connector set must be considered
- the corresponding parameter in the view *Pump size and single-step mode* must be updated

IMPORTANT: When 2 blood pumps need to be replaced, replace the right blood pump in the first place, subsequently replace the left blood pump.

IMPORTANT: Sedate the patient if necessary and administer a bolus of Heparin according to the the anticoagulation protocol.

14.1.1 Preparing a replacement blood pump

Material

- 1 replacement blood pump of appropriate type and size
- 1 driving tube, red or blue
- 1 accessory set (for blood pumps with PU valves) with tube connecting set;

IMPORTANT: Only the cable ties and cable tie guns provided should be used.

► INSTRUCTION

1. Bring membrane to the end-of-diastole position, position de-airing needle, rinse and fill pump with sterile injectable saline (see section 14.5: Driving blood pump(s) with the manual pump, page 154).
2. Connect the driving tube to the respective driving tube connector of the pump.
3. Place the pump, ready for connection, with the titanium connectors pointing upwards.

14.1.2 Replacing the right blood pump (RVAD/ BVAD)

Material

- 1 prepared replacement blood pump (see section 14.1.1: Preparing a replacement blood pump, page 146)
- 1 tube connecting set (cable tie, cable-tie gun), included in the accessory set. Only the cable ties and cable tie guns provided should be used.

Stopping the right blood pump and detaching the blood pump from Ikus

► INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. The cable tie covering the EXCOR cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula. Check cannulae immediately to make sure they are not damaged.
3. If necessary log into the monitor program by entering user ID and password, confirming the password with <Enter>.
4. BVAD: Reduce rate of left blood pump to 30 bpm. Use <←>/ <→> to navigate cursor to the respective field of the parameter table, then use <↓> to adapt value. Confirm with <Enter>.
5. In the monitor program, select the option **Pause left** respectively **Pause right** and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The right blood pump will stop.
RVAD: **Pause left**
BVAD: **Pause right**
The view *Pump size and single-step* mode is displayed.
6. As soon as the right pump has stopped, clamp off the cannulae beneath the right pump to be replaced and slide the cannulae off the pump. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
7. Check cannulae for visible deposits. If necessary, remove these deposits carefully.
8. Remove the driving tube of the pump to be replaced from the connector. To do so, take hold of the release sleeve and pull this out of the connector.

Connect new right blood pump to the Ikus

► INSTRUCTION

1. Fill the free ends of the cannulae with sterile saline solution. Make sure that all air has been removed. Connect the prepared replacement pump to the cannulae.
2. Plug the new driving tube into the freed connector. The plug snaps into place clearly audible.
3. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
4. Release the tube clamps from the cannulae.

Starting the Ikus

► INSTRUCTION

1. Move the cursor to the field **step left** (RVAD) respectively **step right** (BVAD).
2. RVAD: Confirm **Step left** with <Enter> to trigger a single step.
BVAD: Confirm **Step right** with <Enter> to trigger a single step.
3. If any air bubbles are visible remove them via the de-airing needle. When all air has been completely removed from the left pump: remove the de-airing needle.
4. Move cursor to the **OK** field and press <Enter> to confirm. The driving unit starts up again using the defined parameters.
5. Check whether the pump is filling correctly and, if necessary, adjust the parameters.
6. Secure all connections with cable ties. See 8.11, page 94.

14.1.3 Replacing the left blood pump (LVAD/ BVAD)

▲ WARNING

BVAD: If the left pump is being replaced, the right pump must also be stopped while the pump is being replaced. Otherwise there is the risk of pulmonary edema.

Material

- 1 prepared replacement blood pump (see section 14.1.1: Preparing a replacement blood pump, page 146)
- 1 tube connecting set (cable tie, cable-tie gun), included in the accessory set. Only the cable ties and cable tie guns provided should be used.

Stopping the right left blood pump and detaching the blood pump from Ikus

► INSTRUCTION

1. Bring the patient into the Trendelenburg position.
2. The cable tie covering the EXCOR cannula on the stub of the blood pump should be removed carefully. Important: never use a sharp instrument, for example, a scalpel or scissors, to remove the cable tie. This may cause damage to the cannula. Check cannulae immediately to make sure they are not damaged.
3. If necessary log into the monitor program by entering user ID and password, confirming the password with <Enter>.
4. In the monitor program, select the option **Pause left** respectively **Drive pause** and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The left blood pump respectively both blood pumps will stop.
LVAD: **Pause left**. The view *Pump size and single-step* mode is displayed.
BVAD: **Drive pause**. The view *Select operating mode* is displayed.
5. As soon as the right pump(s) has/ have stopped, clamp off the cannulae beneath the left pump to be replaced and slide the cannulae off the pump. If it is necessary to clamp any other part of the cannula that is not covered with velour, cover the part of the cannula that will be clamped with a gauze sponge.
6. Check cannulae for visible deposits. If necessary, remove these deposits carefully.
7. Remove the driving tube of the left pump to be replaced from the connector. To do so, take hold of the release sleeve and pull this out of the connector.

Connect new left blood pump to the Ikus

► INSTRUCTION

1. Fill the free ends of the cannulae with sterile saline solution. Make sure that all air has been removed. Connect the prepared replacement pump to the cannulae.
2. Plug the new driving tube into the freed connector. The plug snaps into place clearly audible.
3. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!

4. Release the tube clamps from the cannulae.

Starting the Ikus

► INSTRUCTION

1. Move the cursor to the field **step left**.
2. Confirm **Step left** with <Enter> to trigger a single step.
3. If any air bubbles are visible remove them via the de-airing needle. When all air has been completely removed from the left pump: remove the de-airing needle.
4. Move cursor to the **OK** field and press <Enter> to confirm. The driving unit starts up again using the defined parameters.
5. Check whether the pump is filling correctly and, if necessary, adjust the parameters.
6. Secure all connections with cable ties.

14.2 Restarting Ikus

▲ WARNING

Do not switch off and restart the *Ikus* unless the service consultant requests to do so (e. g. in emergency operating mode).

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [O] position!. Otherwise there is a risk that the drive may fail in future due to the *Ikus* batteries being totally discharged.

If the graphs in the monitor program are frozen (not moving) and the parameters cannot be adjusted even after the laptop and the monitor program have been restarted, then the *Ikus* is operating in emergency pulse mode. In this case, do not proceed as described here, but as described in section 14.3: Emergency pulse mode, page 150 instead.

► INSTRUCTION

1. Support patient with a replacement *Ikus* (see section 14.4: Connecting the patient to a replacement Ikus, page 152) or manual pump (see section 14.5: Driving blood pump(s) with the manual pump, page 154).
2. Use the seal plugs to seal the driving tube connectors.
3. Switch off the *Ikus* driving unit. In the monitor program, select the **Drive OFF** option and press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> or the <1> key. Important: If it is not possible to select the **Drive OFF** option, the *Ikus* is running in emergency pulse mode. In this case, proceed as explained in section 14.3: Emergency pulse mode, page 150.
4. Wait until the log has been completed. If the message **Switch off drive with main switch!** appears, turn the key switch to the [O] position.
5. Switch off the laptop.
6. Switch *Ikus* on again immediately. To do so, turn the key switch to the [I] position.
7. Switch the laptop on. Select the **1. Start program** option (<1>). Enter user ID and password, confirming the password with <Enter>.
8. Check all parameters and re-adjust them if necessary.
9. For univentricular operation, unplug the connector marked in red, for biventricular operation, unplug both connectors. To do so, pull the seal plug(s) out of the respective connector(s).
10. Disconnect the driving tube(s) from the replacement *Ikus* or the manual pump and connect it to the *Ikus*. Important: Observe the colored markings. The sound of the plug snapping into place is clearly audible.

11. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
12. Move cursor to the **OK** field and press <Enter> to confirm. The standard view is shown. The system will operate with the current parameter settings.
13. Switch off the replacement *Ikus* driving unit (see section 5.3.3: Drive OFF: switching the *Ikus* off, page 54).

14.3 Emergency pulse mode



If the *Ikus* is operating in emergency pulse mode, the user must immediately visually check the blood pump(s) to determine whether the pump(s) are filling and ejecting completely. If one pump is not filling and/or ejecting completely the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152 and section 14.5: Driving blood pump(s) with the manual pump, page 154 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

If the emergency pulse mode is activated while the backup system is already active, the *Ikus* is no longer able to drive both pumps. In this case the patient must be supported immediately with the replacement *Ikus*. Use the manual pump while securing the replacement *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152 and section 14.5: Driving blood pump(s) with the manual pump, page 154 resp.). Otherwise there is the risk that the patient will not be supported sufficiently.

IMPORTANT: In emergency pulse mode a controlled shut down is not possible.

IMPORTANT: In emergency pulse mode it is not possible to acknowledge the acoustic alarm.

If the graphs in the monitor program are not moving and the parameters cannot be adjusted even after the laptop and the monitor program have been restarted, then the *Ikus* is operating in emergency pulse mode. Both control computers have failed and are no longer communicating with each other. The emergency pulse board has taken over control of the left and the right pneumatic system.

Replace the *Ikus* with a backup *Ikus* if possible. If no replacement *Ikus* is available, the *Ikus* will continue to support the patient in emergency pulse mode until a replacement *Ikus* is ready.

In emergency pulse mode the system operates with the following settings:

Synchronous mode (biventricular)

	Systol. pressure [mmHg]	Diastol. pressure [mmHg]	Rate [bpm]	relat. systole duration [%]
left	210	-40	70	40
right	150	-40	70	40

Tab. 14-2 Settings in emergency pulse mode

Always immediately ...



Notify Berlin Heart! 866.249.0128

14.3.1 Emergency pulse mode - switching the *Ikus* off



Only proceed as described below if a replacement *Ikus* is available to assume the support of the patient as described here.

The *Ikus* power switch (toggle switch) should always be in the [I] position, even if the main switch (key switch) is in the [O] position!. Otherwise there is a risk that the drive may fail in future due to the *Ikus* batteries being totally discharged.

INSTRUCTION

1. Support the patient by a replacement *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152).
2. Use the seal plugs to seal the driving tube connectors on the *Ikus*.
3. Switch off the *Ikus*. To do so, turn the key switch to [O] position.
4. Switch off the laptop.

14.3.2 *Ikus* start-test following emergency pulse mode



Do not reconnect the original *Ikus* to the patient until the *Berlin Heart, Inc.* service department has evaluated the LOG files or has serviced the driving unit.

Wait for at least 5 minutes after switching the *Ikus* off while in emergency pulse mode. Otherwise, only the service staff will be able to restart it.

Never connect other USB devices (e.g. wireless technology) to the USB port of the laptop than the delivered USB sticks.

INSTRUCTION

1. After 5 minutes: Insert the USB stick into the USB port and switch on the *Ikus*. To do so, turn the key switch to the [I] position. The battery charge indicator will light up and the number of hours the driving unit has been operated to date will be displayed. The mains operation indicator lights up.
2. Switch the laptop on. The menu *Select language* appears.
3. Select the desired language by pressing the corresponding number key. It is not necessary to press <Enter> to confirm this selection. The start menu is displayed on the laptop.
4. Select the **1. Start program** option (<1>). Enter user ID and password, confirming the password with <Enter>. The *Ikus* will carry out a start-test.
5. Wait for the start-test phase to finish (this takes a few minutes). Do not mute the acoustic signal. The messages in the message window inform the user of the current status of the test. If the system is found to be operating correctly, the view *Select operating mode* will be displayed next.
6. Select **Drive OFF**, then press <Enter> to confirm.
7. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops operation immediately and writes an operating log.
8. Wait until the log has been completed (message: **Switch off drive with main switch!**). Important: Do not switch the *Ikus* off yet.
9. Shut down the monitor program. Press <F10> and confirm by pressing the <X> key or the

- <1> key. The start menu is displayed on the laptop.
- 10. Insert the USB stick in the laptop before the laptop is switched on. In the start menu, select **4. Save data** (<4>). The LOG files are copied onto the USB stick. After that the start menu appears again.
- 11. Switch off laptop. Take the USB stick out of the port (but never when the laptop is switched on).
- 12. Send the LOG files by e-mail to service@berlinheart.com.

14.4 Connecting the patient to a replacement *Ikus*

CAUTION

The *Ikus* must always be connected to the power supply when it is switched on. This is the only way to ensure that the start-up test is performed completely and possible malfunctions can be detected.

This is necessary if ...

- maintenance is required
- the *Ikus* is defective

Switching the replacement *Ikus* on

INSTRUCTION

1. Prepare the replacement *Ikus* and connect it to the mains. Secure the mains cable with the plug clip. Ensure that the power switch (toggle switch) is set to **[I]** position.
2. Use the seal plugs to seal both driving tube connectors.
3. Switch on the replacement *Ikus*. To do so, turn the key switch to the **[I]** position. The battery charge indicator will light up and the number of hours the driving unit has been operated to date will be displayed. The mains operation indicator lights up.
4. Switch the laptop on. The menu *Select language* appears.

Setting the parameter values of the replacement *Ikus*

NOTICE

If no parameter values are entered into the replacement *Ikus*, the replacement *Ikus* starts up with the following default parameter values (standard parameters):

Systole [mmHg] left/ right	Diastole [mmHg] left/ right	Rate [bpm]	Rel. diast. duration [%] left/ right	Operation mode
210/ 130	-40/ -20	80	40/40	biventricular, synchronous and separate

Tab. 14-3 Default standard parameters

INSTRUCTION

1. Select the desired language by pressing the corresponding number key. It is not necessary to press <Enter> to confirm this selection. The start menu is displayed on the laptop.
2. Select the **1. Start program** option (<1>). Enter user ID and password, confirming the password with <Enter>. The *Ikus* will carry out a start-test.
3. Wait for the start-test phase to finish (this takes a few minutes). Do not mute the acoustic signal. The messages in the message window inform the user of the current status of the test. If the driving unit is found to be operating correctly, the view *Select operating mode* will be displayed next.
4. Select **Univentricular (UVAD)** or **Biventricular (BVAD)**, then confirm the selected operating mode with <Enter>.

5. In biventricular mode, the view *Pump size and single-step mode* is shown. In univentricular mode, a connector seal test is first performed (taking approx. 10 seconds). The *Ikus* checks whether the driving tube connector with the blue marking has been sealed. Then the view *Pump size and single-step mode* is displayed.
6. Transfer all the parameters from the original *Ikus* to the replacement *Ikus*.
7. Move cursor to the field **step left** and press <Enter> to confirm.
Biventricular: Move cursor to the field **step right** and press <Enter> to confirm.
8. Move the cursor to the **OK** field. Important: Do not confirm **OK** yet.

The default settings for systole, diastole and relative systole duration depend on whether the pump is registered as the left or right pump in the monitor program.

Univentricular: Default setting as for left pump.

Connecting the blood pump(s) to the replacement *Ikus*

► INSTRUCTION

1. For univentricular operation, unplug the connector marked in red, for biventricular operation, unplug both connectors. To do so, remove the seal plugs from the respective connector(s).
2. If possible, log into the monitor program of the original *Ikus*. Select **Drive OFF**, then press <Enter> to confirm. Respond to the prompt in the dialog window by pressing the <X> key or the <1> key. The system stops immediately. If it is not possible to select the **Drive OFF** option, proceed with triggering single steps (**Step left/ Step right**, see instructions in section 14.1: Replacing the blood pump(s), page 146) without stopping the *Ikus*.
3. As soon as the original *Ikus* has stopped, remove the driving tube(s) from it. To do so, take hold of the release sleeve and pull this out of the connector(s).
4. Connect the driving tube(s) to the replacement *Ikus*.
Important: Observe the colored markings. The sound the plug snapping into place is clearly audible. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. Do not pull from the release sleeve, and never from the tube!
5. To confirm the **OK** selection, press <Enter>. The replacement *Ikus* will start up using the defined parameter settings.
6. Check whether the pump is filling correctly and, if necessary, adjust the parameters.
7. Switch off the original *Ikus* (see section 5.3.3: Drive OFF: switching the *Ikus* off, page 54).

14.5 Driving blood pump(s) with the manual pump




Fig. 14-1 Patient on *Ikus*



Fig. 14-2 Patient on manual pump

This is necessary if ...

- the power supply to the *Ikus* cannot be ensured
- the *Ikus* has to be restarted (e.g. emergency operating mode) and there is no replacement *Ikus* available

WARNING

The use of the manual pump is only permitted for medical personnel trained in the use of it.

Pay attention to the colored markings on the driving tubes and on the connectors of the manual pump. Otherwise, there is a risk of lung edema.

Always keep manual pump attached to the *Ikus*. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

Call one or more persons to assist. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

The driving tubes and cannulae should be arranged in a bend-free position. Otherwise in an emergency situation the adequate support of the patient is not guaranteed.

When operating the manual pump with 1 hand, do not block the valves with your feet (see c in figure 14-3 and in fig. 14-4, page 144).

NOTICE

Seal the connector(s) on the *Ikus* immediately after removing the driving tube(s) in order to avoid contaminants from entering the system.

IMPORTANT: In biventricular mode: the blood pumps are driven asynchronously by the manual pump.

INSTRUCTION

1. The patient is lying down.
2. Disconnect the driving tube(s) from the *Ikus*. To do so, take hold of the release sleeve and pull this out of the connector.
3. Connect the driving tube(s) to the manual pump.
Important: Observe the colored markings.
4. Check that the plug is securely connected. To do so, grip the plug body above the release sleeve and pull on it. (see f in figure 14-4, page 155:) Do not pull from the release sleeve, and never from the tube!
5. Pump steadily and rhythmically at roughly 60 to 80 strokes per minute. Important: Move the piston so far that the membrane reaches its final position. The piston need not necessarily be moved to its end position.
6. Perform a visual check of the blood pump to verify that the membrane is moving and that blood is being pumped.




Fig. 14-3 Manual pump

Fig. 14-4 Plug on the driving tube




Fig. 14-5 Examples to operate the manual pump

The manual pump can be operated with both hands or with one hand (placing the pump between the feet). Alternating between two-handed or one-handed pumping, as well as using the left or right hand, is allowed. When doing so, care of the patient must remain ensured.

14.6 Mains failure or breakdown of both control computers

INSTRUCTION

1. Assess the condition of the patient.
2. Check the filling and ejection behavior of the blood pump(s).
3. If possible: ensure support of patient with a replacement *Ikus* and switch off the faulty *Ikus* (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152).

If both control computers fail at the same time or if there is a mains failure, the *Ikus* cannot generate any specific error messages in the message window. An acoustic alarm sounds and the indicator lamp on the handle lights up. Depending on the type of fault, the message **Error: no data from Master**. or **Error: no reaction from Master** appears in the message window.

Ikus is running in emergency pulse mode (see section 14.3: Emergency pulse mode, page 150).



Notify Berlin Heart! 866.249.0128

14.7 Reading out the LOG files

This is necessary if it is not possible to clearly identify function faults even after consultation with the service department.



If the user exits the monitor program, it is not possible to identify any incoming messages. For this reason, always start the monitor program again immediately after saving the data.

Use 1 of the 2 USB sticks provided with the device to save the data. Do not connect any other USB device to the laptop (e.g. wireless technology).

Make sure that you always have the USB stick inserted and that there is sufficient capacity on the stick. Otherwise the LOG files might get lost as they are deleted from the hard disk as soon as they have been transferred onto a USB stick.

INSTRUCTION

1. Press <F10> to exit the monitor program and confirm the intention in the dialog window by pressing the <X> key or the <1> key. The start menu is displayed. The *Ikus* will continue to operate using the current parameter settings.
2. Switch off the laptop. Insert the USB stick into the port (never do this while the laptop is switched on:). Switch on the laptop again.
3. Select the **4. Save data** option in the start menu. The LOG files are saved onto the USB stick. After completion the start menu appears.
4. Switch off the laptop and remove the USB stick (never do this while the laptop is still switched on!)
5. Switch the laptop on again. To return to the monitor program select the **1. Start program** option in the start menu.
6. Enter user ID and password, confirming the password with <Enter>.
7. Send the LOG files by e-mail to service@berlinheart.com.



Fig. 14-6 Laptop CF30 with inserted USB stick



Fig. 14-7 USB stick with extended plug (operating position)

Fig. 14-8 USB- stick with retracted plug (position to transport and store the stick)

14.8 Circuit breaker and battery fuse

The *Ikus* battery pack is protected against excessive current.

There is an internal battery fuse that protects against excessive load current. From outside it is unamenable and after triggering it has to be replaced by the service team.

A second fuse (the circuit breaker) protects against excessive charging current. The circuit breaker is located on the connection panel and it is resettable after one-time triggering (button *Circuit breaker* flips out) by the operator.

Activated circuit breaker or internal battery fuse in mains operation

WARNING

Never disconnect the *Ikus* from the mains when the circuit breaker or the internal battery fuse are activated. This will cause the driving unit to stop immediately.

INSTRUCTION

1. Ensure that the *Ikus* is connected to the mains.
2. Check if the circuit breaker was triggered. Push the button *Circuit breaker* back in place to produce power supply again. Important: Only press the button briefly. Never keep the button pressed for a longer period, because otherwise the retriggering of the circuit breaker would not be detected.
3. If the circuit breaker is reactivated, do not press the button in again. Never keep the button pressed for a longer period of time. If possible, support the patient with a replacement *Ikus*. (see section 14.4: Connecting the patient to a replacement *Ikus*, page 152).

HOTLINE

Notify Berlin Heart! 866.249.0128

Activated circuit breaker or internal battery fuse in battery operation mode

WARNING

The *Ikus* stops immediately. The blood pumps are no longer being driven.

Whenever the *Ikus* is running in battery operation, the patient must always be accompanied by a person trained to use the manual pump.

NOTICE

If the circuit breaker or internal battery fuse is activated in battery operation, the *Ikus* generates an acoustic alarm.

INSTRUCTION

1. Check if the resettable circuit breaker was triggered.

- If yes: immediately push it back in place. Start the *Ikus* again just in case that it does not happen automatically.
2. If the circuit breaker is triggered again, immediately ensure the support of the patient with the manual pump.



Notify Berlin Heart! 866.249.0128

15 Appendix

15.1 Overview: Product range and possible combinations

15.1.1 Blood pumps

Designation	Article number	Inflow / outflow [mm]
Blood pump PU valves 10 ml in/out ø 6 mm	P10P-001	6/6
Blood pump PU valves 25 ml in/out ø 9 mm	P25P-001x01	9/9
Blood pump PU valves 30 ml in/out ø 9 mm	P30P-001x01	9/9
Blood pump PU valves 50 ml in/out ø 12 mm	P50P-001	12/12
Blood pump PU valves 60 ml in/out ø 12 mm	P60P-001	12/12

Tab. 15-1 Blood Pump

15.1.2 Overview: Relationship: body weight – pump size




Fig. 15-1 Relationship: body weight - pump size

15.1.3 LV apex cannulae

Designation	Article number	Length of head [mm]	Overall length [mm]	Lumen diameter [mm]
Apex cannula for infants	C14A-040	14	220	5
Apex cannula for small children	C18A-020	18	250	6
Apex cannula for children, staged	C22A-004	28	270	12, 9; head 9
Apex cannula	C27A-001	38	265	12

Tab. 15-2 LV apex cannulae

15.1.4 Atrial cannulae

Designation	Article number	Length of head [mm]	Length of corpus [mm]	Lumen diameter [mm]
Atrial cannula for infants	C15V-040	15	200	5
Atrial cannula for small children	C19V-020	19	250	6
Atrial cannula for children, staged (with mandrin)	C22V-004	22	280	9, 12; head 9
Atrial cannula for children, staged (with mandrin)	C25V-004	25	280	9, 12; head 9
Atrial cannula (with mandrin)	C22V-002	22	330	12
Atrial cannula (with mandrin)	C26V-002	26	330	12

Tab. 15-3 Atrial cannulae

15.1.5 Arterial cannulae

Designation	Article number	Head angle [°]	Length of corpus [mm]	Lumen diameter [mm]
Arterial cannula for infants	C80G-040	80	200	5
Arterial cannula for small children	C80G-021	80	250	6
Arterial cannula for children, staged	C60G-004	60	280	9, 12; head 9
Arterial cannula for children, staged	C85G-004	85	280	9, 12; head 9
Arterial cannula	C60G-002	60	330	12
Arterial cannula	C85G-002	85	330	12

Tab. 15-4 Arterial cannulae

15.1.6 Overview: Which cannulae should be used for which pump?

Pump: connector \varnothing [mm]	Which pump?	Cannula: lumen \varnothing [mm] where cannula joins pump	Which inflow cannula?	Which outflow cannula? (arterial cannula)
6	P10P-001	5/ 6	C15V-040 (AT) C19V-020 (AT) C14A-040 (AP) C18A-020 (AP)	C80G-040 C80G-021
9	P25P-001x01 P30P-001x01	9	C22V-004 (AT;SC) C25V-004 (AT;SC) C22A-004 (AP;SC)	C60G-004 (SC) C85G-004 (SC)
12	P50P-001 P60P-001	12	C22V-004 (AT;SO) C25V-004 (AT;SO) C22V-002 (AT) C26V-002 (AT) C22A-004 (AP;SO) C27A-001 (AP)	C60G-004 (SO) C85G-004 (SO) C60G-002 C85G-002
Explanation:	AT atrial cannula AP apex cannula SO staged (stepped diameter) cannula, original diameter		SC staged (stepped diameter) cannula, diameter after cutting to size	

Tab. 15-5 Which cannulae for which pump?

15.1.7 System accessories

Designation	Article number
Accessory set for blood pumps with PU valves (membrane set, de-airing set and tube connecting set)	T00L-002
Driving tube, red; length: 200 cm	L20H-002
Driving tube, blue; length: 200 cm	L20H-003
Tank unit	1600422

Tab. 15-6 System Accessories

15.1.8 Driving unit

Designation	Article number
EXCOR® Stationary Driving Unit <i>Ikus</i> (115V/ 60Hz) - SW 3.41	D03I-111

Tab. 15-7 Driving unit

15.1.9 Special components

Designation	Article number
Connector set for cannulae \varnothing 6 mm to \varnothing 9 mm	A06-009
Connecting set for cannulae \varnothing 9 mm to \varnothing 12 mm	A09-012
Cannula tunnelling tip	attached to cannula

Tab. 15-8 Special components

15.1.10 Pump-cannula combinations

Cannulation		Blood pumps				
∅ inflow cannula	∅ outflow cannula	10 ml	25 ml	30 ml	50 ml	60 ml
5 mm	5 mm	130 bpm				
6 mm	5 mm	130 bpm				
6 mm	6 mm	130 bpm	80 bpm	65 bpm		
9 mm	6 mm		100 bpm	90 bpm		
9 mm	9 mm		130 bpm	130 bpm	130 bpm	105 bpm
12 mm	9 mm				130 bpm	105 bpm
12 mm	12 mm				130 bpm	125 bpm

Tab. 15-9 Pump-cannula combinations



Pump-cannula combinations in which not every parameter combination is recommended (pump rate, % systole, systolic and diastolic pressure,) can lead to incomplete filling and emptying of the blood pump.

Rate value (bpm)

The value indicated is the upper threshold for pump rates. Values that are below the upper threshold are within the acceptable range. Values that are higher than the upper threshold are in a questionable range. The threshold values have been determined (in vitro) taking a mean arterial blood pressure of 120 mmHg as a basis.

Rate value (bpm)

Red marked values displayed on the laptop: These parameter combination (pump rate, % systole, systolic and diastolic pressure) for these pump-cannula combination can lead to incomplete filling and emptying of the blood pump(s). Observe the filling behavior of the blood pump(s)!

in biventricular mode

The lower value of both pump rates (corresponding to the pump sizes used) must also be considered. The higher of the 2 pump rates should be disregarded.

15.1.11 Blood pump combinations in biventricular mode

The following combinations are recommended:

- left pump 10 ml - right pump 10 ml (10 ml/ 10 ml)
- left pump 30 ml - right pump 25 ml (30 ml/ 25 ml)
- left pump 60 ml - right pump 50 ml (60 ml/ 50 ml)

Check whether a blood pump combination that is not recommended has been selected for the patient. The final decision on the combination of blood pumps and cannulae is to be reached by the implanting surgeon, in consultation with *Berlin Heart, Inc* Clinical Affairs.

15.1.12 Relative systolic duration

The relative systolic duration is adjustable in the range of 20% and 70%. The upper and lower threshold (20-30% and 60-70%) are marked in red on the laptop. For these values it cannot be guaranteed that the activated pressure parameters are achievable for each single case.

15.2 Technical specifications

Product	Electro-pneumatic extracorporeal ventricular assist device <i>EXCOR® Pediatric VAD with Stationary Driving Unit Ikus</i>
Manufactured by:	Berlin Heart GmbH Wiesenweg 10 12247 Berlin Germany
Classification	Class 3
Overall system (except sterile products)	
Ambient temperature in operation	+10 °C ... +30 °C with restrictions of the battery performance up to +35°C
Ambient temperature, transportation and storage	-10 °C ... +50 °C; 6 h warming-up period before commissioning after transport
Max. permitted ambient magnetic field strength	10 A/m
Relative humidity of environment	45 % to 75 %
Ambient atmospheric pressure	max. 2000 m (6562 ft) above MSL (mean sea level)
Pump	
Dimensions	Refer to product data sheets
Material	Casing and membranes: polyurethane Driving tube adapter: polyoxymethylene Connectors: titanium
Coating of blood contact surfaces	<i>Carmeda® BioActive Surface (CBAS®)</i>
Max. period of use	max. 1 year
Cannulae	
Dimensions	Refer to product data sheets
Material	Silicone, partially reinforced with plastic webbing, partially encased with suture-suitable polyester velour; some equipped with flexible metal reinforcement: wire 2 mm, circular steel Rd 1.4301; apex cannula with a titanium alloy shell
For all sterile products	
Long-term storage conditions	Temperature: +15°C to 25°C Relative humidity: 35 % to 50 % Store in a dry place!

<i>Ikus</i>	
Dimensions (W x H x D)	46 x 95 x 73 cm with laptop cover down (approx. 18.5 x 37.5 x 29 inches) 46 x 120 x 73 cm with laptop cover open (approx. 18.5 x 47.5 x 29 inches)
Driving tube (red/blue)	max. 1 year period of use
Weight	100.6 kg (approx. 219 lb)
Input voltage	115 VAC
Frequency	60 Hz
Power drawn	575 VA
Mains fuse	5 A
Mains cable	10 A, hospital grade
Connector <i>External alarm</i>	electical data: max 1 A, 24 V insulation specifications: 2.5 mm/ 4 mm clearance and creepage distance between alarm contact and 24 V extra low voltage inside the device (coil-sided) insulation test voltage: 500 V
Protection class	IPX1 (protection against touching live parts not tested, tested safety from vertically dripping water)
Pump rate	30 to 150 bpm
Systolic pressure:	60 to 350 mmHg
Diastolic pressure	-100 to 0 mmHg
Pressure display accuracy	± 10 %
Relative systolic duration	20 % to 70%
Off-mains operating time	max. 30 minutes
Battery charging time	6 h
Maintenance interval	2000 operating hours (at the latest after 6 months). In the event of permanently higher ambient temperatures than recommended, the maintenance intervals can shorten drastically.
Product life <i>Ikus</i>	max. 8 years

Tab. 15-10 Technical specifications

15.3 Symbols and tags

15.3.1 *Ikus*: identification plate



Fig. 15-2 Applied part type B, according to DIN EN 60601-1: 1990 + A1:1993 + A2:1995

IPX1

IPX1 is a tag on the identification plate.

EC529 and the most recent IEC60529 is an international standard that describes a system to classify the degrees of protection provided by enclosures of electronic equipment. IPX1 indicates tested safety from vertically dripping water. For a medical product IPX1 indicates drip-proof, a higher than ordinary level of protection from drips, leaks, and spills.

15.3.2 *Ikus*: connection panel




Fig. 15-3 *Ikus* main switch (key switch) in [0] position; the device is switched off



Fig. 15-4 *Ikus* main switch (key switch) in [I] position; the device is switched on




Fig. 15-5 Connector *Potential equalization*, symbol according to DIN EN 60601-1: 1990 + A1:1993 + A2:1995

15.3.3 Symbol on numlock status LED



Fig. 15-6 Example of symbol for numlock status LED on laptop

15.5 Sample copy: EXCOR Implantation log



Berlin Heart®

Implantation Record Form

EXCOR® VAD



This form applies **only** to USA and Canada



Please fill out the form (3 pages), and fax it to Berlin Heart, Inc. *immediately* after implantation (fax: 866.540.5026).

After replacing a blood pump, please fill out the "Pump Replacement" section (page 1), list the supplies used on page 2/3, and fax (3 pages) to Berlin Heart Inc. (fax: 866.540.5026)

Hospital	City/Country
-----------------	---------------------

Patient data (for Berlin Heart registry)

Patient's initials	Sex m <input type="checkbox"/> / f <input type="checkbox"/>	Age	Body size [cm]	Weight [kg]
Patient-No. (BH Site No. followed by the patient No. ie: 004-103)	IABP pre-op n <input type="checkbox"/> y <input type="checkbox"/>		ECMO pre-op n <input type="checkbox"/> y <input type="checkbox"/> since _____ (Date)	
	On transplantation list n <input type="checkbox"/> y <input type="checkbox"/> since _____ (Date)			
Ischemic CMP <input type="checkbox"/>	Idiopathic CMP <input type="checkbox"/>	Acute Myocarditis <input type="checkbox"/>	Postcardiotomy <input type="checkbox"/>	
Acute Myocardial Infarction <input type="checkbox"/>	Congenital <input type="checkbox"/>		Other <input type="checkbox"/>	
PAP mean [mmHg]	CVP [mmHg]	MAP [mmHg]	LVEF % FS %	
CI [l/min/m ²]	NYHA	LVEDP [mmHg]	LVEDD [mm]	
Creatinine [mg/dl]	Total Bilirubin [mg/dl]	Platelet count [/μl]	Leukocytes [/μl]	

Implantation

Date	Surgeon					
Type BVAD <input type="checkbox"/> LVAD <input type="checkbox"/> RVAD <input type="checkbox"/>	Access medial <input type="checkbox"/> lateral <input type="checkbox"/>		Left-sided cannulation atrial <input type="checkbox"/> apical <input type="checkbox"/>			
LVAD	Pump type: PU valve <input type="checkbox"/> Tilting-disk valve <input type="checkbox"/>					
	Pump size: 10 ml <input type="checkbox"/> 25 ml <input type="checkbox"/> 30 ml <input type="checkbox"/> 50 ml <input type="checkbox"/> 60 ml <input type="checkbox"/> 80 ml <input type="checkbox"/>					
RVAD	Pump type: PU valve <input type="checkbox"/> Tilting-disk valve <input type="checkbox"/>					
	Pump size: 10 ml <input type="checkbox"/> 25 ml <input type="checkbox"/> 30 ml <input type="checkbox"/> 50 ml <input type="checkbox"/> 60 ml <input type="checkbox"/> 80 ml <input type="checkbox"/>					

Pump replacement

Left pump <input type="checkbox"/>	Reason for replacement _____
Date	Location of deposit inflow <input type="checkbox"/> outflow <input type="checkbox"/> pump chamber <input type="checkbox"/>
Right pump <input type="checkbox"/>	Reason for replacement _____
Date	Location of deposit inflow <input type="checkbox"/> outflow <input type="checkbox"/> pump chamber <input type="checkbox"/>



Berlin Heart®

Implantation Record Form

EXCOR® VAD

This form applies **only** to USA and Canada

Please record the lot numbers of the used EXCOR® components used and the components to be kept as back-up and fax to Berlin Heart, Inc. *immediately* after implantation (fax: 866.540.5026).

Hospital/City		Date of Implantation	
Patient ID (BH Site No. followed by the patient No. ie: 004-103)			
Ikus-No.		Ikus hours of operation	
Replacement Ikus		Ikus hours of operation	
Ikus-No.		Ikus hours of operation	
Item	Lot-No.		Article No.
	used	b/u	
EXCOR Blood Pumps with PU valves			
10 ml in/out Ø 6 mm			P10P-001
25 ml in/out Ø 9 mm			P25P-001x01
30 ml in/out Ø 9 mm			P30P-001x01
50 ml in/out Ø 12 mm			P50P-001
60 ml in/out Ø 12 mm			P60P-001
80 ml in/out Ø 12 mm			P80P-001***
EXCOR Blood Pumps with Tilting-disk valves			
50 ml in/out Ø 12 mm			P50M-001***
60 ml in/out Ø 12 mm			P60M-001***
80 ml in/out Ø 12 mm			P80M-001***
80 ml out/in Ø 12 mm (in/out exchanged)			P80M-005***
80 ml in/out Ø 16 mm			P80M-003***
80 ml out/in Ø 16 mm (in/out exchanged)			P80M-004***
EXCOR Apex Cannulas			
Ø 5 mm, L 22 cm (Apex cannula for infants)			C14A-040
Ø 6 mm, L 25 cm (Apex cannula for small children)			C18A-020
Ø 12/9 mm, L 27 cm (Apex pediatric cannula, staged)			C22A-004
Ø 12 mm, L 26,5 cm (Apex cannula, one-piece)			C27A-001
Ø 16 mm, L 33 cm (Apex cannula)			C41A-050***
EXCOR Atrial Cannulas			
Ø 5 mm, L 20 cm, head 15 mm (Atrial cannula for infants)			C15V-040
Ø 6 mm, L 25 cm, head 19 mm (Atrial cannula for small children)			C19V-020
Ø 12/9 mm, L 28 cm, head 22 mm (Atrial pediatric cannula, staged)			C22V-004
Ø 12/9 mm, L 28 cm, head 25 mm (Atrial pediatric cannula, staged)			C25V-004
Ø 12 mm, L 33 cm, head 22 mm (Atrial cannula)			C22V-002
Ø 12 mm, L 33 cm, head 26 mm (Atrial cannula)			C26V-002
Ø 12 mm, L 33 cm, head 30 mm (Atrial cannula)			C30V-002

*** Not available for general use in the US and Canada

+++ Not available for general use in the US

Berlin Heart, Inc.
200 Valleywood, Suite B400
The Woodlands, TX 77380
www.berlinheart.com



Berlin Heart®

Implantation Record Form EXCOR® VAD



This form applies **only** to USA and Canada



Item	Lot-No.		Article No.
	used	b/u	
EXCOR Arterial Cannulas			
Ø 5 mm, L 20 cm (Arterial cannula for infants)			C80G-040
Ø 6 mm, L 25 cm (Arterial cannula for small children)			C80G-021
Ø 12/9 mm, L 26 cm (Graft-adapter pediatric cannula, staged)			C00P-004+++
Ø 12/9 mm, L 28 cm, 85° (Arterial pediatric cannula, staged)			C85G-004
Ø 12/9 mm, L 28 cm, 60° (Arterial pediatric cannula, staged)			C60G-004
Ø 12 mm, L 33 cm, 60° (Arterial cannula)			C60G-002
Ø 12 mm, L 33 cm, 85° (Arterial cannula)			C85G-002
Ø 12 mm, L 26 cm (Graft-adapter cannula)			C00P-001+++
Ø 16/12 mm, L 36 cm, 85° (Arterial cannula, staged)			C85G-050+++
Ø 16 mm, L 26 cm (Graft-adapter cannula)			C00P-050+++
Connecting Set for Cannulas			
Ø 6/9 mm			A06-009
Ø 9/12 mm			A09-012
Ø 12/16 mm			A12-016***
Accessories			
Accessory set Tilting-disk valves			T00L-001***
Accessory set PU-valves			T00L-002
Driving tube, red Ø 6/8 mm, L 2 m			L20H-002
Driving tube, blue Ø 6/8 mm, L 2 m			L20H-003
Tank unit			1600422


*** Not available for general use in the US or Canada

+++ Not available for general use in the US

Date	Signature
-------------	------------------

15.6 Sample copy: EXCOR pump log

15.6.1 Explanations on the pump log



- 1 transition inflow cannula - inflow connector
- 2 only on pumps with PU valves: inflow stub in front of inflow valve
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 only on pumps with PU valves: outflow stub behind outflow valve
- 11 transition outflow connector - outflow cannula


Fig. 15-7 Numbering of the checkpoints



To briefly describe the findings, use the following letter codes:

- p - small punctual deposit
- P - large punctual deposit
- a - small area of deposit
- A - large area of deposit
- f - small strand
- F - large strand
- t - small thrombus
- T - large thrombus
- ~above the respective letter indicates floating deposits

Example: Plotting of the deposits



- 1 small laminar deposit
- 2 small suture on the inflow valve
- 3 small specks
- 4 larger laminar deposit


Fig. 15-8 Plotting of the deposits

Example: Notation with letter code

Datum date	Zeit time	Name Sign.	Linke Pumpe/ Left pump		ml: 50ml No.: 0815										
					1	2	3	4	5	6	7	8	9	10	11
01.01.	8:00	z.B.			a		F			p		A			

Fig. 15-9 Notation for letter code

15.6.2 Sample copy: EXCOR pump log



- 1 transition inflow cannula - inflow connector
- 2 only on pumps with PU valves: inflow stub in front of inflow valve
- 3 inflow valve
- 4 inflow stub behind inflow valve
- 5 area between inflow and outflow stubs
- 6 remaining area of blood chamber
- 7 transition blood chamber - membrane (directly above the reinforcement ring)
- 8 outflow stub in front of outflow valve
- 9 outflow valve
- 10 only on pumps with PU valves: outflow stub behind outflow valve
- 11 transition outflow connector - outflow cannula

Fig. 15-10 EXCOR blood pump with checkpoint numbers

Patient:

Lot No.: Pump left

Lot No.: Pump right

Datum date			Zeit time			Name Sign.			Linke Pumpe/ Left pump				Rechte Pumpe/ Right pump																							
									ml:				ml:																							
									No.:				No.:																							
									1	2	3	4	5	6	7	8	9	10	11	1	2	3	4	5	6	8	7	9	10	11						

16 EMC tables

16.1 Essential Performance

The following essential performance was verified in the electromagnetic immunity tests: The driving unit must drive the EXCOR blood pumps with the set parameters. This was controlled by monitoring the following acceptance criteria:

- The pump rate may not deviate by more than 10% or 50 ms (the higher value is valid).
- The relative systolic duration may not deviate by more than 10% or 50 ms (the higher value is valid).
- The driving pressure in engaged condition at the end of the systole or diastole may not deviate by more than 10% or 20 mmHg (the higher value is valid).
- No failure alarms may occur.
- No undesired change may occur between mains and battery operation.
- In battery operation the battery alarms must occur correctly.
- No switch to emergency pulse mode may occur.

16.2 Electromagnetic emissions

The *Ikus* is intended for use in the electromagnetic environment specified below. The customer or the user of the *Ikus* should ensure that *Ikus* is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The <i>Ikus</i> uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The <i>Ikus</i> is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Tab. 16-1 Emissions characteristics

16.3 Electromagnetic immunity - part 1

The *Ikus* is intended for use in the electromagnetic environment specified below. The customer or the user of the *Ikus* should ensure that it is used in such an environment

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec	< 5 % U_T (>95 % dip in U_T) for 0.5 cycle 40 % U_T (60 % dip in U_T) for 5 cycles 70 % U_T (30 % dip in U_T) for 25 cycles < 5 % U_T (>95 % dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the <i>Ikus</i> requires continued operation during power mains interruptions, it is recommended that the <i>Ikus</i> be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	100 A/m	The <i>Ikus</i> can be used up to 1 m from power cables carrying up to 100 A.
NOTE U_T is the a.c. mains voltage prior to application of the test level.			

Tab. 16-2 Electromagnetic immunity - part 1

16.4 Electromagnetic immunity - part 2

The *Ikus* is intended for use in the electromagnetic environment specified below. The customer or the user of the *Ikus* should ensure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz outside ISM bands ^a	10 V	Portable and mobile RF communications equipment should be used no closer to any part of the <i>Ikus</i> , including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 0.35 \sqrt{P}$
	10 V _{rms} 150 kHz to 80 MHz in ISM bands ^a	10 V	$d = 1.2 \sqrt{P}$
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	30 V/m 80 MHz to 6 GHz	$d = 0.4 \sqrt{P}$ 80 MHz to 800 MHz $d = 0.77 \sqrt{P}$ 800 MHz to 6 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^c should be less than the compliance level in each frequency range. ^d Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Tab. 16-3 Electromagnetic immunity - part 2

- a The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.
- b The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in these frequency ranges.
- c Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the *Ikus* is used exceeds the applicable RF compliance level above, the *Ikus*

should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the *Ikus*.

d Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 10V/m.

16.5 Recommended separation distances between portable and mobile RF communications equipment and the *Ikus*

The *Ikus* is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the *Ikus* can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the *Ikus* as recommended below, according to the output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m			
	150 kHz to 80 MHz outside ISM bands $d = 0.35 \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 0.4 \sqrt{P}$	800 MHz to 6 GHz $d = 0.77 \sqrt{P}$
0.01	0.04	0.12	0.04	0.08
0.1	0.11	0.38	0.13	0.24
1	0.4	1.2	0.4	0.8
10	1.1	3.8	1.3	2.4
100	3.5	12	4	7.7

Tab. 16-4 Separation distance depending on frequency of transmitter

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

- NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
- NOTE 2 The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.
- NOTE 3 An additional factor of 10/3 has been incorporated into the formulae used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.
- NOTE 4 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

List of figures

Fig. 1-1	Handle and lifting bar.....	14
Fig. 3-1	EXCOR shown in situ as a biventricular assist device in pediatric application	33
Fig. 3-2	60 ml blood pump	34
Fig. 3-3	Cannula heads: 1) atrial cannula, 2) LV apex cannula, 3) arterial cannula	34
Fig. 4-1	Ikus.....	37
Fig. 4-2	Connection panel.....	38
Fig. 4-3	Display and operating panel	40
Fig. 4-4	Battery charge indicator with 7 yellow LEDs, red LED to the left of yellow LEDs	40
Fig. 4-5	Laptop	42
Fig. 4-6	Monitor program.....	42
Fig. 4-7	Pneumatic system redundancy in biventricular operation	42
Fig. 5-1	LEDs below the mouse pad	45
Fig. 5-2	De-activate the Caps Lock.....	45
Fig. 5-3	De-activate the Numeric Lock.....	45
Fig. 5-4	Select language.....	46
Fig. 5-5	Start menu	46
Fig. 5-6	User ID.....	47
Fig. 5-7	Password.....	47
Fig. 5-8	Entry codes	47
Fig. 5-9	Change date or time	48
Fig. 5-10	User ID.....	49
Fig. 5-11	Password.....	49
Fig. 5-12	Standard view of monitor program during start-test	50
Fig. 5-13	Selecting monitor program options	51
Fig. 5-14	Select pump size.....	51
Fig. 5-15	Select cannulae size	52
Fig. 5-16	Parameter table.....	52
Fig. 5-17	Drive OFF (with stop options)	53
Fig. 5-18	Select operating mode.....	53
Fig. 5-19	Pump size and single-step mode.....	54
Fig. 5-20	Battery charge indicator with 7 yellow LEDs, red LED to the left of yellow LEDs	55
Fig. 5-21	Handle and lifting bar.....	59
Fig. 5-22	Empty shipping crate.....	61
Fig. 5-23	Prepare the laptop	61
Fig. 5-24	Line up the Ikus	61
Fig. 5-25	Engage the brake	61
Fig. 5-26	Place the accessory box.....	62
Fig. 5-27	Reattach the securing plate.....	62
Fig. 5-28	Lift up the ramp	62
Fig. 5-29	Fix the ramp	62
Fig. 5-30	Close the door.....	62
Fig. 5-31	Lock the crate.....	62
Fig. 5-32	Crate ready for storage/ shipping	63
Fig. 6-1	Connection panel.....	66
Fig. 6-2	Tank unit.....	66
Fig. 6-3	Select operating mode.....	68
Fig. 6-4	Select the pump size	69
Fig. 6-5	Cannula in for the left pump	70
Fig. 6-6	Display pump and cannula sizes	71
Fig. 6-7	Pump size and single-step mode (with recommended start-up parameters)	71
Fig. 6-8	Single-step mode.....	73
Fig. 7-1	Cable tie gun.....	79
Fig. 7-2	De-airing nipple and driving tube connector.....	81
Fig. 7-3	Pump with trocar in place (de-airing needle with inserted obturator)	81
Fig. 8-1	Cannula position following implantation	84

Fig. 8-2	Suggested cannulae exit sites (Example: BVAD with LV apex cannulation)	84
Fig. 8-3	Available sizes of cannula tunneling tips	85
Fig. 8-4	Use of cannula tunneling tip	85
Fig. 8-5	Anastomosis of LV apex cannula	87
Fig. 8-6	Ideal position of the LV apex cannula	87
Fig. 8-7	Head of LV apex cannula	87
Fig. 8-8	Cannulation of right atrium	90
Fig. 8-9	Suture technique, right atrium	90
Fig. 8-10	Anastomosis of the aortic cannula	91
Fig. 8-11	Cannulation of the pulmonary artery	92
Fig. 8-12	Final position of the blood pumps, for example: BVAD with LV apex cannulation	93
Fig. 8-13	Removing the de-airing needle	94
Fig. 10-2	Cleaning the blood pump	100
Fig. 10-1	Materials for dressing change	100
Fig. 10-3	Cleanse each cannula exit site	101
Fig. 10-4	Cleanse with back/forth motion	101
Fig. 10-5	Cleanse entire cannula	101
Fig. 10-6	Non sting barrier film	102
Fig. 10-7	Metalline drain compress	102
Fig. 10-8	Secure with a sterile bandage	102
Fig. 10-9	Gauze compress under the cannulae	102
Fig. 10-10	Fold the left end of compress and secure	102
Fig. 10-11	Fold the right end of compress and secure	102
Fig. 10-12	Cannulae are padded	103
Fig. 10-13	Cover with sterile gauze compresses	103
Fig. 10-14	Secure with a sterile bandage	103
Fig. 10-15	Seal with strips of adhesive bandage	103
Fig. 10-16	Tubular bandage	103
Fig. 10-17	Diagram of EXCOR blood pump (top view of blood chamber)	105
Fig. 12-1	Pause right	126
Fig. 12-2	Drive OFF	127
Fig. 12-3	Start menu	128
Fig. 13-1	Example of an error message in the message window	129
Fig. 14-1	Patient on Ikus	154
Fig. 14-2	Patient on manual pump	154
Fig. 14-3	Manual pump	155
Fig. 14-4	Plug on the driving tube	155
Fig. 14-5	Examples to operate the manual pump	155
Fig. 14-6	Laptop CF30 with inserted USB stick	157
Fig. 14-7	USB stick with extended plug (operating position)	157
Fig. 14-8	USB stick with retracted plug (position to transport and store the stick)	157
Fig. 15-1	Relationship: body weight - pump size	159
Fig. 15-2	Applied part type B, according to DIN EN 60601-1: 1990 + A1:1993 + A2:1995	165
Fig. 15-3	Ikus main switch (key switch) in [0] position; the device is switched off	165
Fig. 15-4	Ikus main switch (key switch) in [I] position; the device is switched on	165
Fig. 15-5	Connector Potential equalization, symbol according to DIN EN 60601-1: 1990 + A1:1993 + A2:1995	165
Fig. 15-6	Example of symbol for numlock status LED on laptop	165
Fig. 15-7	Numbering of the checkpoints	170
Fig. 15-8	Plotting of the deposits	170
Fig. 15-9	Notation for letter code	170
Fig. 15-10	EXCOR blood pump with checkpoint numbers	171

List of tables

Tab. 1-1	Connector set: minimum length of connector tube.....	18
Tab. 2-1	Product life, maintenance interval, expiration date, warranty	31
Tab. 4-1	Settings in emergency pulse mode.....	43
Tab. 5-1	Other password configuration options	48
Tab. 5-2	Options changing date or time	48
Tab. 5-3	Parameter adjustment.....	53
Tab. 5-4	Displays and messages during battery operation mode.....	57
Tab. 6-1	Test parameters	67
Tab. 6-2	Start-up parameters.....	72
Tab. 6-3	Assignment: operating mode, blood pump, connector.....	72
Tab. 6-4	Pump filling criteria	74
Tab. 6-5	Parameter's possible adjustments	74
Tab. 6-6	Recommended guideline values for normal operation.....	75
Tab. 8-1	Connector set: minimum length of connector tube.....	86
Tab. 8-2	Connector set: minimum length of connector tube.....	92
Tab. 8-1	Anastomosis and direction of the blood chambers	93
Tab. 11-1	Unfractionated Heparin adjusted to maintain an anti factor Xa level of 0.35 to 0.50 U/ml.	111
Tab. 11-2	Enoxaparin, low molecular weight heparin dosing (⁶ Monagle, P, et al.)	111
Tab. 11-3	Warfarin loading dose to maintain an INR of 2.7 - 3.5 (⁶ Monagle, P, et al.)	112
Tab. 11-4	Warfarin Maintenance Dosing for Day 5 and longer to maintain INR 2.7-3.5.....	112
Tab. 11-5	Drugs and Dose for specific INR range.....	112
Tab. 12-1	Important parameters for weaning progress	114
Tab. 12-2	Reduced rate day 1 to day 3	114
Tab. 12-3	Important parameters for weaning progress	118
Tab. 12-4	Reduced rate day 1 to day 3	118
Tab. 12-5	Important parameters for weaning progress	122
Tab. 12-6	Reduced rate day 1 to day 3	122
Tab. 14-1	Possible problems	146
Tab. 14-2	Settings in emergency pulse mode.....	150
Tab. 14-3	Default standard parameters	152
Tab. 15-1	Blood Pump	159
Tab. 15-2	LV apex cannulae.....	160
Tab. 15-3	Atrial cannulae	160
Tab. 15-4	Arterial cannulae.....	160
Tab. 15-5	Which cannulae for which pump?.....	161
Tab. 15-6	System Accessories	161
Tab. 15-7	Driving unit.....	161
Tab. 15-8	Special components	161
Tab. 15-9	Pump-cannula combinations.....	162
Tab. 15-10	Technical specifications	164
Tab. 16-1	Emissions characteristics.....	173
Tab. 16-2	Electromagnetic immunity - part 1	174
Tab. 16-3	Electromagnetic immunity - part 2	175
Tab. 16-4	Separation distance depending on frequency of transmitter.....	176

Index

A

access passwords	44
accessories	
description - blood pump, cannulae, ~	33
EXCOR components and ~	79
EXCOR Pediatric ~	35
system ~	161
accessory set	
for blood pumps with PU valves	79, 161
for commissioning and operation	35
adverse events	27
alarm	
external ~ connector	164
visual signal display	40
ambient temperature	21, 24
in operation	163
anticoagulation	
status	105, 145
therapy	107
asynchronous pulse	
in biventricular operating mode	51
seperate pulsing left/right	41

B

battery	
~ charge indicator	40
~ charge indicator with 7 yellow LEDs	40
~ charging time	164
~ operation indicator	40
~ power operation	41
activated circuit breaker	157
circuit breaker	157
replacement and disposal	27
switching over to ~ operating mode	56
toggle between mains and ~ operation	55
battery operation not possible	138
biventricular	
~ operation, adjusting the operating mode	74
assist device, shown in place	33
dysfunction	29
emergency operating mode in ~ mode	133
emergency pulse mode	43
explantation after ~ support	126
operating mode	41
overview, assignment	72
test parameter	67
biventricular operation	
pneumatic system redundancy in ~	42
pneumatic system redundancy in ~, figure	42
blood pump	
biventricular operation	41
combinations in biventricular mode	162
connecting the ~ to the cannulae	93

connecting to the replacement Ikus	153
consists of	29
de-airing	81
EXCOR Pediatric ~	34
explantation	113
operation period	13, 30
preparing a replacement ~	146
product range	159
regular checks of ~	103
replacing	146
rinse and fill the ~	82
stopping	53
stopping an individual ~	54
univentricular operation	41
BTT, bridge to transplant	127
BVAD see also biventricular	
BVAD, biventricular	
changing from RVAD to ~	58
manual pump in ~ mode	43
selecting operating mode	53

C

cannula tunneling tip	85
available sizes of ~	85
use of ~	85
cannulae	
apex ~	87, 160
arterial ~	160
atrial ~	88, 160
connecting the blood pumps to the ~	93
description	33
EXCOR Pediatric ~	34
outflow ~	90
regular checks of the ~	103
removing the VAD ~	127
shortening	92
technical specifications	163
types of ~	34
Which ~ should be used for which pump?	161
cardiac transplantation	29
Carmeda® BioActive Surface	17, 34, 163
cautionary measures	
adjusting the parameter values	104
EMC	24, 25
filling behavior of the blood pump(s)	104
monitor program	50
visual inspection - deposits	105
checks	
~using the monitor program	105
regular ~ of blood pumps and cannulae	103
Cleaning	21
components	
and accessories	79
permanently active ~	33
preparing the ~ and materials required	79
connection panel	37, 38, 66, 165
connector seal test	68

contraindications	30	what to do?	129
control computer	37	error messages	
breakdown of both ~s	155	acoustic alarm is not properly recognized	140
emergency operating mode	150	alarm circuit fault - acoustic alarm cannot be ac-	
redundant design	43	tivated/deactivated	135
counter showing number of operating hours	40	backup computer...	
D		... is reporting abnormal measurements on the left	
de-airing		(or right)	135
~ hammer	35, 79	... is reporting the error - flow measurement on the	
~ needle, figure	81	left (or right)	136
~ nipple, figure	81	...is reporting a checking error	136
~ set	35, 79, 81, 161	...is reporting flow measurements deviation on the	
~ step	72	left (or right)	135
inserting the ~ needle	81	backup operation left/right	133
removing the ~ needle	94	batteries discharged	138
deposits		check the left (or right) driving tube and pump ..	132
~ in the pump	145	computer faulty - contact customer service	135
floating	145, 170	connecting the drive tube	131
large thrombus	170	emergency operating mode	133
large, laminar	170	biventricular mode	133
large, speck-like	170	error - no data/no reaction from the control computer	
plotting of the ~	170	140
small thrombus	170	left/ right flow sensor defective. Contact Service	
small, laminar	170	department!	141
small, speck-like	170	pressure error/ time error in system 1 (or in system 2	
visual inspection	104	or 3)	130
dressing		problem - <<text>>	141
applying a new ~	102	replacing a drive tube	131
removing the old ~	100	system 1 (or system 2) is defective	134
driving pressure		system 3 is defective	134
possible range	74	temperature sensors - <<8-digit binary code>> ...	137
systolic	29	throttle valve error in system 1 (or system 2 or 3)	131
driving tube	33	131	
as an active component	33	updating of parameter block failed	137
connector	34, 39	error messages during the start-up test	
connector, figure	34, 81	additional messages during the start-up test	142
function	33	battery operation test skipped due to battery problem	
in univentricular operation	41	142
on the tank unit	66	Essential Performance	173
replacing	131	events, adverse	27
E		examination, medical ~ of patient	104
echocardiogram, intraoperative trans-esophageal	97	EXCOR Pediatric	79
electro-stimulation therapy	25	parameter log, sample copy	172
EMC		pump log, sample copy	171
special precautions required	24, 25	EXCOR Pediatric HDE Study	146
tables	173	explantation	113
emergency operating mode	43, 150	~ after biventricular support	126
error messages	133	~ after univentricular support	126
Ikus start-test following ~	151	external alarm	
restarting Ikus in ~	149	~ connector	164
equipotential bonding connector	66	~ connector, figure	66
error message		H	
~s and corrective measures	129	hours of use	21
acknowledging	129		

I

identification plate	15, 165
Ikus	
ambient conditions	26
commissioning the ~ and setting parameters	65
description, Stationary Driving Unit Ikus	37
designation	11
device description	29
error messages and corrective measures	129
instructions for use	45
markings and symbols	165
moving the ~	59
product life time	164
recommended separation distances	176
software version	1
technical specifications	163
transport	27
~ation and packaging	60
troubleshooting and correcting faults	145
implantation	
anesthesia	97
preparations in the operating room	79
surgical procedure	83
infection	99
risk of ~	99
to prevent ~	18
inspection	
filling behavior of blood pump	104
frequency of ~	103
of pump areas which come in contact with blood	105
visual ~ - deposits	104
instruction and care of the patient	27
interaction with other procedures and therapies	26

K

key switch (main switch)	
and power switch (toggle switch)	39
figure	66
restart the Ikus	143
switching on the Ikus	66
switching the Ikus off	54, 126, 127

L

laptop	
~ computer with monitor program	37, 42
software	1
start menu	46
use of ~	45
version	1
LOG files, reading out the ~	156
LVAD	
~ anastomosis	93
~ anesthesia	97
~ assignment	72
term	41

M

magnetic resonance imaging	25
main switch (key switch)	165
~ and power switch	39
figure	66, 165
mains power operation	
activated circuit breaker in ~	157
mains power cable	55
switching over to ~	58
toggle between ~ and battery operation	55
maintenance	27, 152
~ interval	164
manual pump	
driving a blood pump with the ~	154
figure	155
purpose	38, 43
membrane set	35, 79, 161
message window	67, 129, 151
browsing in the ~	53
messages, see error messages	
mobilization	104
monitor program	45
basic instructions	49
figure	42
functions	51
laptop computer with ~	37, 42
logging in and out of the ~	49
shutting down the ~	49
standard view	50
start menu	46
starting	49
monitoring	
~ procedure	97

N

nuclear	
~ diagnostics	25
~ therapy	25
numlock	
deactivated at the start	45, 65
Symbol on ~ status LED	165

O

obturator	81
operating hours, see hours of use	
operating manager, obligations of ~	27
operating mode	72
biventricular operation	41
emergency	43
selecting the ~	68
synchronous pulse	51
univentricular operation	41

P		S	
packaging		safety risks	27
and sterilization	17	sample copy	
and transportation	60	Blood pump replacement	172
sterile ~	81	EXCOR Implantation log	167
parameter		EXCOR Pediatric parameter log	172
adjusting the ~ values	52, 104	EXCOR Pediatric pump log	171
adjustment	53	Ikus Incoming Checklist	166
adjustment, possible range	74	selecting pump size(s) and cannula sizes	51
checking the ~s when the pump is started and adjusting them	73	self-test of the alarm circuit	49, 142
commissioning the Ikus and setting ~s	65	separate mode	
EXCOR Pediatric parameter log	172	biventricular operating mode	51
set the test ~s	67	left/right	41
start-up ~s, recommended	71	single-step mode	
table	50, 52	figure	73
test ~s	67	removing air from the blood pump(s) in ~	72
password		view - pump size and ~	71
~-protected user profiles	44	standard default parameters	67, 152
envelope	47	start-up parameters, see parameter	
options in the start menu	46	start-up test	
platelet concentrates	25, 97	error messages during ~	59, 142
product range and possible combinations	159	following emergency operating mode	151
power supply	40	routine ~ when not in operation	58
battery operation	41	sterilization and packaging	17
failure	145, 154	storage	
mains ~	41	and durability	27
mains power failure or breakdown of both control computers	155	Ikus, ambient temperature	163
power switch, see toggle switch		suction pressure, diastolic	29
precautions, general	25	possible range	74
pulsing		synchronous pulse	51, 67
asynchronous, left/right	41	in emergency operating mode	43
independent, left/ right	41	synchronous pulsing left/right	41
synchronous, left/ right	41		
pump log		T	
explanations on the EXCOR Pediatric ~	170	tank unit	65
sample copy EXCOR Pediatric ~	171	connecting	66
pump size		disconnect the ~ from the Ikus	68
pop-up menu for selcting the ~	69	technical specifications	162
relationship between bodyweight and ~	159	test parameters	
selecting the ~	69	setting the ~	67
view - pump size and single-step mode	71	table	67
pump, see blood pump	29	time error, see error messages	
		toggle switch (power switch)	66
		during transportation	60
		main switch and power switch	39
		restart the Ikus	149
		switch off the Ikus	54, 126, 127
		switch on the Ikus	66
		transplantation	
		orthotopic heart ~	128
		transport	27
		~ation and packaging	60
		outside the clinic	27
		proceeding after ~	27, 65
		technical specifications	163
		within the clinic	59
R			
radiotherapy	25		
replacing blood pump	146		
restart Ikus (in emergency operation mode)	149		
routine start-test when not in operation	58		
RVAD			
~ anastomosis	93		
~ assignment	72		
changing from ~ to BVAD	58		
stopping an individual blood pump	54		
term	41		

treatmen and wound care	99
Trendelenburg position	72, 93, 147, 148
trocar	
as part of the de-airing set	35, 79
figure	81
to de-air the blood pump	81

U

ultrasonic treatment	25
univentricular	
operation	41
overview, assignment	72
which driving tube?	79
user profile	44
UVAD see also univentricular	
UVAD, univentricular	
change over from ~ to biventricular operation	58
connector seal test	68
default settings	153
explantation after ~ support	126
redundant design of pneumatic systems	42

V

valves	
control ~	37
three-leaflet polyurethane ~	29
visual inspection	
deposits	104
filling and ejecting of the blood pump	22
filling behavior of blood pump	104
regular ~ of blood pumps and cannulae	103

W

weaning	128
wound care	99
and treatment	99

X

X-rays	25
--------------	----

