

## SUMMARY BASIS OF APPROVAL

Reference No.: 95-1588

Proper Name: Human Immunodeficiency  
Virus Type 1 [HIV-1]

Applicant: Cambridge Biotech Corp.    Trade Name: HIV-1 Western Blot  
1500 East Gude Drive  
Rockville, MD 20850-5307

### I. INDICATIONS FOR USE

The HIV-1 Western Blot Kit is an *in vitro* qualitative assay for the detection and identification of antibodies to Human Immunodeficiency Virus Type 1 (HIV-1) contained in human serum, plasma, or urine. It is intended for use as a more specific test with specimens found to be repeatedly reactive using a screening procedure, such as an Enzyme Immunosorbent Assay (EIA).

The urine application of this HIV-1 Western Blot Kit is not intended for initial screening or testing of blood donors. The product includes the following restrictions:

- The HIV-1 Western blot urine procedure is intended for use with urine samples only and must not be used with other bodily fluids.
- The assay must be performed in strict accordance with these instructions to obtain accurate, reproducible results.
- Although a Positive result may indicate infection with HIV-1, a diagnosis of Acquired Immunodeficiency Syndrome (AIDS) can be made only if an individual meets the case definition of AIDS established by the Centers for Disease Control.
- Individuals may present incomplete banding patterns due to the natural history of AIDS or other immunodeficiency states, e.g.:
  - Infants born to HIV-1 infected mothers, but who are uninfected, may display incomplete patterns as passively acquired maternal antibodies begin to disappear.
- Do not use this kit as the sole basis of diagnosis of HIV-1 infection.
- A Negative result does not exclude the possibility of HIV-1 infection.
- The Cambridge HIV-1 Western Blot Kit is a biological product which, although highly consistent, does display variation from lot to lot. Examples of these variations include bands which have a slightly wavy or slanted appearance, small artifacts within the banding area, and a light smearing pattern across a

set of strips. These are considered normal assay variations which infrequently affect assay interpretation. However, if they do interfere with the assay interpretation, the assay should be invalidated and repeated.

- The HIV-1 Urine Western Blot Controls are required when testing urine specimens with the Cambridge Biotech HIV-1 Western Blot Kit. Each Western Blot kit intended for use for testing urine is shipped with a matched set of urine controls.

## II. BRIEF DESCRIPTION OF TEST

The HIV-1 Western Blot Kit is an enzyme-linked immunosorbent blot technique ("Western Blot") intended to be used to detect antibodies to Human Immunodeficiency Virus Type 1 (HIV-1) in human serum, plasma, or urine. Its primary purpose is for use as an additional, more specific test on human serum, plasma or urine specimens found to be repeatedly reactive using a licensed HIV-1 antibody detection test.

The Cambridge Biotech HIV-1 Western blot kit is manufactured by Cambridge Biotech Corporation from HIV-1 propagated in an H9/HTLV-III<sub>B</sub> T-Lymphocyte cell line. Partially purified HIV-1 virus is inactivated by treatment with psoralen and ultraviolet light, and detergent disruption. Specific HIV-1 proteins are fractionated according to molecular weight by electrophoresis of the partially purified, inactivated HIV-1 virus preparation on a polyacrylamide slab gel in the presence of sodium dodecylsulfate (SDS). The separated HIV-1 proteins are electrotransferred from the gel to a nitrocellulose membrane which is then washed, blocked (to minimize nonspecific immunoglobulin binding), and packaged.

Individual nitrocellulose strips are incubated with specimens and controls. During incubation, if HIV-1 antibodies are present in the specimen, they will bind to the viral antigens present on the nitrocellulose strips. The strips are washed again to remove unbound material. Visualization of the human immunoglobulins specifically bound to the HIV-1 proteins is accomplished *in situ* using a series of reactions with goat anti-human IgG conjugated with biotin, avidin conjugated with horseradish peroxidase (HRP), and the HRP substrate, 4-chloro-1-naphthol. If antibodies to any of the major HIV-1 antigens are present in the specimen in sufficient concentration, bands corresponding to the position of one or more of the following HIV-1 proteins (p) or glycoproteins (gp) will be seen on the nitrocellulose strip: p17, p24, p31, gp41, p51, p55, p66, gp120, gp160 (the number refers to the apparent molecular weight in kilodaltons).

The procedural steps of the HIV-1 Western blot [originally designed for use with only serum and plasma] have been slightly modified to allow use of the test kit with urine. The modified steps include sample dilution and conjugate incubation times only. Each HIV-1 Western Blot Kit intended for use with urine is shipped

with a set of HIV-1 Urine Controls. The urine controls replace the serum controls and are used with the other test components whenever urine specimens are tested.

### Components of the Test Kit

1. NITROCELLULOSE STRIPS - Each NITROCELLULOSE STRIP contains separated, bound antigenic proteins from partially purified, inactivated HIV-1, in sufficient quantity to detect human antibodies. Bovine protein is present as a blocking agent. Strips are consecutively numbered.
2. NON-REACTIVE SERUM CONTROL - Normal human serum, non-reactive for antibodies to HIV-1 and non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide and 0.005% thimerosal as preservatives.
3. WEAKLY REACTIVE SERUM CONTROL - Inactivated human serum containing a low titer of antibodies to HIV-1 antigens and non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide and 0.005% thimerosal as preservatives.
4. STRONGLY REACTIVE SERUM CONTROL- Inactivated human serum containing a high titer of antibodies to HIV-1 antigens and non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide and 0.005% thimerosal as preservatives.
5. WASH BUFFER - Supplied as a 20x concentrate. When diluted, this reagent contains 0.02 M tris, 0.1 M NaCl, 0.3% Tween 20, and 0.005% thimerosal as a preservative, at pH 7.4.
6. BLOTTING BUFFER - Supplied as a 10x concentrate. When diluted, this reagent contains 0.02 M tris, 0.1M NaCl, heat-inactivated normal goat serum, and 0.1% thimerosal as a preservative, at pH 7.4.
7. CONJUGATE 1 - Biotinylated Goat Anti-human IgG (heavy and light chain) antibodies. Contains 0.002% thimerosal as a preservative.
8. CONJUGATE 2 - Avidin conjugated horseradish peroxidase. Contains 0.01% thimerosal as a preservative.
9. SUBSTRATE A - 7.8 mM solution of 4-chloro-1-naphthol in an alcohol solution.
10. SUBSTRATE B - Aqueous hydrogen peroxide solution (0.02%) in citrate buffer.
11. BLOTTING POWDER - nonfat dry milk (sterile).

Packaged separately, but shipped with each Western blot kit intended for use with urine:

12. **NEGATIVE URINE CONTROL** - Human urine negative for antibodies to HIV-1 antigens. A serum sample from the source was non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide as a preservative.
13. **LOW POSITIVE URINE CONTROL** - Inactivated human urine positive for antibodies to HIV-1 antigens. A serum sample from the source was non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide as a preservative.
14. **HIGH POSITIVE URINE CONTROL** - Inactivated human urine containing a high titer of antibodies to HIV-1 antigens. A serum sample from the source was non-reactive for hepatitis B surface antigen and antibodies to HCV. Contains 0.1% sodium azide as a preservative.

### III. MANUFACTURING AND CONTROLS

#### A. Manufacturing and Controls

The HIV-1 Western Blot is manufactured at Cambridge Biotech Corporation facilities and the Urine Controls are manufactured at Calypte Biomedical facilities described below.

Raw materials intended for use in the product are subjected to appropriate quality control evaluation before they are accepted for use in manufacturing.

All kit components are tested for sterility or monitored for low bioburden content consistent with the procedures provided in the manufacturer's PLA. Final performance testing of the Western blot kit measures potency, reproducibility, sensitivity, and specificity.

#### 1. Manufacturing of HIV-1 viral lysate.

HIV-1 viral antigens are purified from HIV-1 virus following propagation in a Human T-cell line (designated H9/HTLV-III<sub>B</sub>) under defined conditions. The viral material and cellular debris are separated from intact cells by cross-flow filtration and concentrated by filtration through a molecular weight cut-off system. The resulting HIV-1 virus concentrate is inactivated by the addition of psoralen and exposure to ultraviolet light. Additional processing via centrifugation, sucrose gradient purification, and detergent lysis yields a final HIV-1 viral lysate which is tested for purity, identity and strength.

2. Manufacturing of nitrocellulose strips, finished components and final kit assembly.

Production of the nitrocellulose strips involves fractionating specific HIV-1 proteins according to molecular weight by electrophoresis on a polyacrylamide slab gel in the presence of sodium dodecylsulfate (SDS). The separated HIV-1 proteins are electrotransferred from the gel to a nitrocellulose membrane which is then washed, blocked (to minimize non-specific immunoglobulin binding), and packaged. Finished nitrocellulose strips must meet in-house specifications for consistency, sensitivity, and specificity using an in-house panel.

The HIV-1 positive and negative serum controls have been shown to be non-reactive for HBsAg, HIV-1, HIV-2, and antibody to HCV and are evaluated against a reference control for functionality and tested for sterility or bioburden. The positive control is heat-inactivated and tested for infectivity.

The Conjugate 1 is prepared by a manufacturer of goat anti-human biotinylated conjugate. Final titration, formulation, bottling, and labeling are performed at Cambridge Biotech. The Conjugate is tested for functionality and sterility or bioburden.

The Conjugate 2 is prepared by a manufacturer of avidin conjugated horseradish peroxidase. Final titration, formulating, bottling, and labeling are performed at Cambridge Biotech. The Conjugate is tested for functionality and sterility or bioburden.

Eight kit components, wash buffer, blotting buffer, substrate A, substrate B, blotting powder, negative urine control, low positive urine control, and high positive urine control are formulated and packaged by a vendor. Each component is subjected to testing at Cambridge Biotech Corporation prior to acceptance of the component.

3. Calypte Biomedical is the sole source of urine controls.

## B. Stability

Stability testing was performed on individual kit components at the recommended storage temperature (2-8°C) and at elevated temperatures (25°C, 37°C, or in some cases 45°C). At regular time points, the components were substituted into a Master Kit Lot to assess their performance over time. The proposed stability dating is based on real-time data of at least three lots of components stored at the recommended temperature (2-8°C). Using this method, expiration dating of 12 months has been assigned to all components except Conjugate 1 and 2, the Negative Urine Control, and the High Positive Urine Control, which have been

assigned 9 months, and the Low Positive Urine Control which has been assigned 6 months.

The urine controls are packaged separately from, but attached to, the rest of the Western blot kit. The master serum kit and the urine control set have the same expiration date, which is based on the shortest dated component in both parts of the kit. The master serum kit and the urine control set are matched, and this kit type for use in testing urine specimens is assigned a new lot number, and lot released by FDA, so that the master serum kit reagent lots and the urine reagent lots may not be interchanged.

### C. Methods of Validation

All components of the HIV-1 Western blot and urine control kit are routinely monitored by the following methods:

1. NITROCELLULOSE STRIPS - Strip lots are tested by Western blot for the presence of required proteins. Each lot of lysate used for strip preparation is tested for purity and potency by Coomassie gel, Western blot, reverse phase chromatography, and peptide mapping.
2. NON-REACTIVE SERUM CONTROL - New lots of negative control are functionally compared to the approved reference. Testing is also performed for hepatitis B surface antigen, antibodies to HCV, pH and sterility. After vialing, testing is repeated for functionality, including a sensitivity and specificity panel, fill volume, physical appearance, and bioburden.
3. WEAKLY REACTIVE SERUM CONTROL - New lots of weakly reactive control are functionally compared to the approved reference. Testing is also performed for hepatitis B surface antigen, antibodies to HCV, pH and sterility. After vialing, testing is repeated for functionality, including a sensitivity and specificity panel, fill volume, physical appearance, and bioburden.
4. STRONGLY REACTIVE SERUM CONTROL - New lots of strongly reactive control are functionally compared to the approved reference. Testing is also performed for hepatitis B surface antigen, antibodies to HCV, pH and sterility. After vialing, testing is repeated for functionality, including a sensitivity and specificity panel, fill volume, physical appearance, and bioburden.
5. WASH BUFFER - Testing for functionality, pH, conductivity, bioburden, and fill volume is performed on every lot.
6. BLOTTING BUFFER - Testing for functionality, pH, conductivity, bioburden, and fill volume is performed on every lot.
7. CONJUGATE 1 - After receipt of the conjugate concentrate from the supplier, coarse and fine titration is performed. After dilution, the material is

tested for sterility, functionality, and pH. Vialled conjugate concentrate is evaluated for functionality, fill volume and bioburden.

8. CONJUGATE 2 - After receipt of the conjugate concentrate from the supplier, coarse and fine titration is performed. After dilution, the material is tested for sterility, functionality, and pH. Vialled conjugate concentrate is evaluated for functionality, fill volume and bioburden.
9. SUBSTRATE A - Testing for functionality, pH, conductivity, bioburden, and fill volume is performed on every lot.
10. SUBSTRATE B - Testing for functionality, pH, conductivity, bioburden, and fill volume is performed on every lot.
11. BLOTTING POWDER - Testing for functionality, bioburden, and fill weight is performed on every lot.
12. NEGATIVE URINE CONTROL - New lots of negative control from a source whose serum is non-reactive for hepatitis B surface antigen and antibodies to HCV are tested for bioburden and by Western blot (no bands present).
13. LOW POSITIVE URINE CONTROL - New lots of low positive control from a source whose serum is non-reactive for hepatitis B surface antigen and antibodies to HCV are tested for bioburden and by Western blot (weak gp160 band and other bands may be present).
14. HIGH POSITIVE URINE CONTROL - New lots of high positive control from a source whose serum is non-reactive for hepatitis B surface antigen and antibodies to HCV are tested for bioburden and by Western blot (gp160, gp120, gp41 and p24 bands and other bands may be present).
15. KITS - Each lot is evaluated by an internal sensitivity and specificity panel, and the FDA panel. Each kit lot designated for use with urine is also tested for lot release using a urine release panel. A representative sample of each kit is submitted to CBER for lot release testing and review of protocols summarizing pertinent product testing prior to distribution. — master lots of the HIV-1 Western blot kits were submitted for evaluation in support of licensing.

D. Labeling

The labeling, including immediate container and package labels and the package insert (directions for use), have been reviewed for compliance with 21 CFR 610.60; 610.61, 610.62, and 809.10, and were found to be satisfactory. The package insert states "The Cambridge Biotech HIV-1 Western Blot Kit is an *in vitro* qualitative assay for the detection and identification of antibodies to Human Immunodeficiency Virus Type 1 (HIV-1). This more specific assay is used as a supplemental test with specimens (serum, plasma, or urine) that tested repeatedly reactive using a

screening procedure (Enzyme Immunosorbent Assay ). Urine specimens are to have been tested with the Calypte/Seradyn Sentinel HIV-1 Urine EIA Kit." The product trade name, Cambridge Biotech HIV-1 Western Blot Kit, is not known to conflict with another biologic or device trade name.

E. Establishment Inspection

An establishment inspection was conducted by the FDA district office and CBER (October 6 - 17 and November 7, 1997) at the Cambridge Biotech Corporation facility in Rockville, MD (1500 E. Gude Dr., \_\_\_\_\_). Office of Compliance, CBER, determined on April 30, 1998 that all facility issues were resolved satisfactorily.

F. Environmental Impact Analysis Report

A detailed EIAR was filed by Cambridge Biotech Corporation for the two facilities responsible for the manufacture of the HIV-1 Western blot kit. There was a finding that no significant environmental impact resulted from manufacture of this product. A summary of procedures taken by the manufacturer to protect the environment is stated as follows.

1. Manufacturing of HIV-1 viral lysate.
  - a. Liquid waste from propagation of viruses is decontaminated by autoclaving or chemical treatment and released into the sewer.
  - b. Solid waste is decontaminated by autoclaving, boxed and labeled as non-infectious waste, and removed for incineration.
  - c. The BSL3 facility is decontaminated with \_\_\_\_\_ between campaigns, or at least annually, for HEPA filter changing and certification.
  - d. \_\_\_\_\_ is used to neutralize residual \_\_\_\_\_ gas.
  - e. Air compressors for the self-contained breathing apparatus have oil and water separators. Waste is labeled as hazardous and removed by a certified waste hauler.
2. Manufacturing of nitrocellulose strips, all finished components (except for urine controls) and kit.
  - a. All liquid chemical waste from nitrocellulose strip manufacture is removed as hazardous chemical waste according to U.S. EPA and Maryland requirements for \_\_\_\_\_

- b. Methanol vapor is exhausted through a \_\_\_\_\_ filter and out the laboratory exhaust stacks. Methanol exhaust is routinely monitored and exhaust is below the state of Maryland Air Quality Regulations units for total volatile organics.
- c. The human defibrinated plasma used to manufacture negative controls is certified as non-reactive for HBsAg, antibodies to HCV and antibodies to HIV-1.
- d. The HIV-1 positive human plasma used to manufacture the positive controls is heat inactivated and tested to confirm lack of infectivity. The human plasma is also confirmed to be antibody positive by Western blot assay and ELISA. Waste is decontaminated by autoclaving, labeled as non-infectious, and removed for incineration.
- e. All wastes from kit components are handled as hazardous chemical wastes in accordance with U.S. EPA, and Maryland EPA. Waste is removed and incinerated.
- f. Employee training meets the requirements of the OSHA Hazard Communication Standard (29 CFR 1910.1200), and the Bloodborne Pathogen Standard (29 CFR 1910.1030).

#### IV. BIOLOGICAL PRINCIPLES OF THE TEST

In the Cambridge Biotech HIV-1 Urine Western Blot the HIV-1 viral antigens are produced by propagation of HIV-1 virus in an H9/HTLV-III<sub>B</sub> T-Lymphocyte cell line. The purified virus is inactivated by treatment with psoralen and ultraviolet light, and detergent disruption. Specific HIV-1 proteins are fractionated according to molecular weight by electrophoresis of the purified, inactivated HIV-1 virus preparation on a polyacrylamide slab gel in the presence of sodium dodecylsulfate (SDS).

The separated HIV-1 proteins are electrotransferred from the gel to the nitrocellulose membrane which is then washed, blocked (to minimize nonspecific immunoglobulin binding), cut into individual strips, and packaged.

Human urine specimens or urine controls are mixed with a diluted blotting buffer and incubated with individual nitrocellulose strips. During the specimen incubation period, if HIV-1 antibodies are present in the specimen or control, they will bind to the viral antigens present on the nitrocellulose strips. The strips are washed to remove unbound material. Visualization of the human immunoglobulins specifically bound to HIV-1 proteins is accomplished *in situ* using a series of reactions with goat anti-human IgG conjugated with biotin, avidin conjugated with horseradish peroxidase (HRP), and the HRP substrate, 4-chloro-1-naphthol. If

antibodies to any of the major HIV-1 antigens are present in the specimen in sufficient concentration, bands corresponding to the position of one or more of the following HIV-1 proteins (p) or glycoproteins (gp) will be seen on the nitrocellulose strip: p17, p24, p31, gp41, p51, p55, p66, gp120, gp160 (the number refers to the apparent molecular weight in kilodaltons).

## V. PRECLINICAL DATA

### A. Preclinical studies

#### 1. Failure Mode Analysis

The procedure for urine testing using the Cambridge HIV-1 Western Blot kit was developed by determining the optimum dilution for the sample incubation step (1:2) and the optimum incubation times for Conjugate 1 and 2 (120 minutes each). The time range for sample incubation previously determined for serum was not altered and was validated in the urine procedure. The wash steps and substrate incubation time also remain as described for serum and were validated for the urine procedure. The composition and preparation of all kit reagents for use in the urine procedure are as described for serum.

#### 2. Reproducibility

To evaluate the reproducibility of the Western blot kits when testing urine specimens, a study was performed by testing a panel of strongly reactive, weakly reactive and negative specimens on 3 different HIV-1 Western blot kit lots. All specimens tested performed as expected.

#### 3. Stability of HIV-1 antibodies in urine specimens

Antibody to HIV-1 in urine specimens were shown to be stable for 12 months when stored at 2-8°C in the presence of Stabilur™. When the specimens are stored at 2-8°C in the absence of preservative, the specimens were shown to be stable for 4 months.

Freezing of urine specimens intended for testing with the Western blot procedure for urine is not recommended.

#### 4. Effects of urine conditions

Potential effects from a variety of urine conditions were evaluated. In such cases as urine pH, creatinine, glucose and osmolality the Western blot results were unchanged.

## 5. Analytical Sensitivity [Urine Specimens]

Sensitivity of the Cambridge HIV-1 Western blot was evaluated on urine specimens obtained from HIV-1 seropositive subjects. Of the originally tested urine samples, 10 had matched serum samples. The specimens were obtained from patients with AIDS and patients who were asymptomatic and symptomatic for infection with HIV-1. The specimens were tested at dilutions of 1:2, 1:10, 1:50, 1:500 and 1:1000 until the major viral bands were no longer visible on the Western blot. (See Section VI Clinical Data.)

## PERFORMANCE CHARACTERISTICS

### PERFORMANCE STUDIES

Three studies were conducted to evaluate the performance of the Cambridge Biotech HIV-1 Western Blot Kit using a urine assay procedure. The performance was evaluated by comparing the results of urine specimens to the results of paired serum specimens tested with a licensed HIV-1 Western Blot Kit.

One study (Study 1) evaluated 696 archived urine specimens. The specimens were from low risk (N=200), high risk (N=37) and HIV-1 positive (N=377) populations. The HIV-1 positive populations included patients symptomatic (N=55) and asymptomatic (N=87) for HIV-1 infection, AIDS patients (N=115) and HIV-1 positive subjects from foreign sites (N=120) whose clinical status was unknown. Other specimens (N=82) were obtained and evaluated from subjects with medical conditions unrelated to HIV-1 infection that might result in antibodies cross-reactive with HIV-1 proteins.

The Cambridge Biotech HIV-1 Urine Western Blot results compared to serum Western Blot results are presented in Table E. Two additional studies (Study 2, Study 3) evaluated 1,240 prospectively collected urine specimens. Study 2 evaluated specimens from subjects whose HIV-1 clinical status was unclassified (N=197), subjects who were HIV-1 negative but at high risk of HIV-1 infection (N=51) and subjects with non-HIV related medical conditions (N=1). Study 3 evaluated low risk (N=315), high risk (N=303) and HIV-1 positive (N=175) populations, including AIDS patients. The HIV-1 positive populations included patients symptomatic (N=38) and asymptomatic (N=36) for HIV-1 infection and AIDS patients (N=101). Other specimens (N=198) were also obtained from subjects with unrelated medical conditions that might result in assay interference.

In the three studies combined, 1,936 paired urine and serum specimens collected from multiple geographical locations within the United States and from foreign sites were evaluated at four testing laboratories throughout the United States. The status of the subject was based upon the paired serum result or documented clinical status of the subject.

Two additional special studies were conducted using specimens from Study 3 to assess the performance of the Cambridge Biotech HIV-1 Western Blot kit. One evaluation involved Western blot testing of urine specimens paired to serum EIA non-reactive Western blot indeterminate specimens (N=109). The second evaluation involved urine Western blot testing of urine EIA repeatedly reactive specimens paired to serum EIA non-reactive specimens from uninfected individuals (N=114).

## SENSITIVITY STUDIES

### Sensitivity In HIV-1 Seropositive Individuals

The sensitivity of the Cambridge Biotech HIV-1 Western Blot Kit using urine was evaluated by comparing the urine results to the results obtained from testing paired serum specimens collected from individuals who were HIV-1 seropositive and from individuals clinically diagnosed as AIDS patients. The results of this study is shown in Table E.

In the combined studies, the Cambridge Biotech HIV-1 Western Blot Kit using urine obtained from the 215 patients clinically diagnosed with AIDS identified 213 of 215 (99.1%) patients as Western Blot positive. There were two (2) false negative urine Western Blot results in the AIDS population. The Cambridge Biotech HIV-1 Western Blot Kit using urine obtained from HIV-1 positive symptomatic, asymptomatic and unclassified groups correctly identified 533 of 533 (100%) patients as Western Blot positive. In the combined population of AIDS patients and other HIV-1 positive patients tested in this study, the Cambridge Biotech HIV-1 Western Blot Kit using urine correctly identified 746 of 748 (99.7%) patients as positive.

**Table E**

**Urine and Paired Serum Western Blot Results for Confirmed HIV-1 Seropositive Individuals and AIDS Patients (N=748)**

Risk Group	N	Serum Western Blot Results <sup>f</sup>			Urine Western Blot Results		
		Pos <sup>i</sup>	Neg	Ind	Pos	Neg	Ind
AIDS <sup>a</sup>	215	215	0	0	213	2	0
Sympt <sup>b</sup>	93	93	0	0	93	0	0
Asympt <sup>c</sup>	123	122	0	1 <sup>h</sup>	123	0	0
Un-classified <sup>d, e</sup>	317	296 <sup>g</sup>	0	1 <sup>h</sup>	317	0	0
<b>Total</b>	<b>748</b>	<b>726</b>	<b>0</b>	<b>2</b>	<b>746</b>	<b>2</b>	<b>0</b>

<sup>a</sup> One hundred and fifteen (115) specimens from Study 1, 100 specimens from Study 3.

<sup>b</sup> Fifty-five (55) specimens from Study 1, 38 from Study 3.

<sup>c</sup> Eighty-seven (87) specimens from Study 1, 36 from Study 3.

<sup>d</sup> One hundred twenty (120) specimens from Study 1, 197 from Study 2.

<sup>e</sup> The clinical status of these HIV-1 positive subjects was unknown.

<sup>f</sup> A licensed serum HIV-1 Western Blot Kit was used when testing serum specimens.

<sup>g</sup> Twenty (20) of the 316 specimens were from Uganda and were not confirmed by Western blot. The specimens were confirmed by a second manufacturer's EIA and by agglutination.

<sup>h</sup> The specimen did not meet the required band intensity criterion for a positive on serum Western blot and therefore was discordant with urine Western blot. However, the patient was known positive by previous clinical diagnosis.

<sup>i</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

## Sensitivity in High Risk Populations

The sensitivity was also determined in 391 individuals at high risk of HIV-1 infection but of unknown HIV status. Of the 391 high risk subjects tested, 327 were substance abusers and 64 were prostitutes, bisexuals, homosexuals, and other individuals with acknowledged risk factors. The results obtained are provided in Table F.

The results of testing urine specimens from these high risk populations showed that seventeen (17) of seventeen (17) (100%) urine Western blot specimens were correctly identified as positive when compared to the paired serum Western blot results. Of the twenty (20) urine Western blot positives, three (3) urine specimens were paired to serum EIA non-reactive specimens (urine false positives). While the significance of urine positivity in the absence of serum reactivity is not known, the results for these three samples must be classified as false positive in the absence of follow-up testing or clinical information to resolve the infection status of these individuals (see Table F, footnote f).

Sixty-nine (69) specimens were urine EIA repeatedly reactive and urine Western blot negative and were paired to serum EIA non-reactive specimens.

One (1) specimen was urine EIA repeatedly reactive and urine Western blot indeterminate and was paired to a serum EIA non-reactive specimen.

Three hundred (300) urine specimens were EIA non-reactive and Western blot negative and were paired to serum EIA non-reactive specimens. The one urine EIA non-reactive specimen that was Western blot indeterminate was paired to a serum EIA non-reactive specimen.

In this study the sensitivity of the urine Western blot was 100% (17 of 17) for seropositive individuals. In these high risk populations, the specificity of the urine Western blot for EIA repeatedly reactive urine specimens was 94.5% (69 of 73) for seronegative individuals.

Table F

Comparison of Cambridge Biotech HIV-1 Western Blot Results  
Using Urine and Paired Serum for High Risk Populations  
(N=391)

Risk Group	Serum				Urine			
	EIA	N	Western Blot <sup>b,d</sup>		EIA	N	Western Blot	
			Pos	N			Pos <sup>e</sup>	N
High	RR	17	17	17	RR	90	20 <sup>f</sup>	20 <sup>f</sup>
			0	0			1 <sup>a</sup>	1 <sup>a</sup>
			0	0			69	69
	NR	374	0	0	NR	301	0	0
			123	123			1 <sup>c</sup>	1 <sup>c</sup>
			230	230			300 <sup>g</sup>	300 <sup>g</sup>
			20 <sup>h</sup>	20 <sup>h</sup>			0	0
			1 <sup>i</sup>	1 <sup>i</sup>			0	0

<sup>a</sup> Non viral bands

<sup>b</sup> A licensed HIV-1 Western Blot Kit was used when testing serum specimens.

<sup>c</sup> p24 only

<sup>d</sup> Serum EIA non-reactive specimens from Study 1 (N=20 drug abusers) were not tested by Western blot.

<sup>e</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

<sup>f</sup> Three (3) urine specimens were Western blot positive and paired to serum EIA non-reactive specimens (urine false positives).

<sup>g</sup> Two hundred twenty-nine (229) were correctly classified as negative when compared to the paired serum Western blot results.

<sup>h</sup> NT - Not Tested

<sup>i</sup> UNR - Unreadable

## Frequency of Virus Specific Bands in High Risk Populations

The frequency of virus specific bands and interpretation of results for urine specimens tested with the Cambridge Biotech HIV-1 Western Blot Kit from these high risk populations are presented in Table G.

The results show that in these combined high risk populations of 391 specimens, only one of the 301 (0.3%) urine EIA non-reactive specimens demonstrated any viral bands. The band present was p24, resulting in an indeterminate interpretation. Twenty (20) of 90 (22.2%) urine EIA repeatedly reactive specimens were identified as positive on the basis of the presence of the gp160 band. Seventeen (17) of the 20 (85%) urine Western blot positive specimens were paired to serum Western Blot confirmed HIV-1 positive specimens. One (1) urine specimen was paired to a serum Western Blot negative specimen and two (2) were paired to serum Western Blot indeterminate specimens. One (1) of the urine EIA repeatedly reactive specimens was indeterminate due to the presence of non-viral bands on the blot. Sixty-nine (69) urine EIA repeatedly reactive specimens were correctly classified as Western blot negative based on the paired serum results. The serum EIA results were non-reactive.

**Table G**

**Frequency of Virus Specific Bands and Interpretation of Results of Urine Specimens from a High Risk Population Tested by the Cambridge Biotech HIV-1 Western Blot Kit (N=391<sup>a</sup>)**

Urine EIA	Urine WB		Frequency of Virus Specific Bands <sup>b</sup>								
		N	p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
NR <sup>c</sup> N=301	Pos <sup>g</sup>	0	0	0	0	0	0	0	0	0	0
	Ind	1	0	1	0	0	0	0	0	0	0
RR <sup>d</sup> N=90	Pos <sup>g</sup>	20	1	8	8	16	8	1	16	17	20 <sup>f</sup>
	Ind	1 <sup>e</sup>	0	0	0	0	0	0	0	0	0

<sup>a</sup> Thirty seven (37) specimens from Study 1, 51 specimens from Study 2, 303 specimens from Study 3.

<sup>b</sup> Band patterns for negative specimens do not appear in this table. By definition, negative specimens show no reactivity.

<sup>c</sup> NR indicates non-reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>d</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>e</sup> Specimen with non-viral bands.

<sup>f</sup> Three (3) of the 20 urine specimens had only a gp160 or gp120 and gp160 band present. All 3 were from Study 1.

<sup>g</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

## Frequency of Virus Specific Bands in AIDS Patients

Specimens from 215 AIDS patients were tested. Table H presents the frequency of viral specific bands observed and interpretation of results for these AIDS patients.

The results show that 213 of 215 (99.1%) specimens from AIDS patients were positive on the basis of the presence of a gp160 band when tested with the Cambridge Biotech HIV-1 Western Blot Kit.

**Table H**

### Frequency of Virus Specific Bands and Interpretation of Results of Urine Specimens from AIDS Patients (N=215<sup>a</sup>)

Urine EIA	Urine WB	N	Frequency of Virus Specific Bands <sup>b</sup>								
			p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
RR <sup>c</sup> N=215	Pos <sup>e</sup>	213	50	104	130	199	114	44	166	210	213 <sup>d</sup>
	Neg	2	0	0	0	0	0	0	0	0	0
	Ind	0	0	0	0	0	0	0	0	0	0

<sup>a</sup> One hundred and fifteen (115) specimens from Study 1, 100 specimens from Study 3.

<sup>b</sup> Band patterns for negative specimens do not appear in this table. By definition, negative specimens show no reactivity.

<sup>c</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>d</sup> Fifteen (15) of the 213 urine specimens had only a gp160 or gp120 and gp160 band present (11 from Study 1, 4 from Study 3).

<sup>e</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

## Frequency of Virus Specific Bands in Non-AIDS HIV-1 Positive Populations

The frequency of virus specific bands in non-AIDS HIV-1 populations (N=533) was also determined. The frequency of virus specific bands in HIV-1 symptomatic, asymptomatic and unclassified HIV-1 positive patients (120 foreign specimens from Study 1, 197 HIV-1 positive patients from Study 2) is given in Table I.

All 533 specimens were serum EIA repeatedly reactive and 531 serum Western blot positive, 2 were serum Western blot indeterminate. The results in Table I demonstrate the presence of a gp160 band in 533 of 533 of the HIV-1 positive urine specimens tested with the Cambridge Biotech HIV-1 Western Blot Kit, classifying all of the urine specimens as Western blot positive. There were no urine specimens that were negative or indeterminate in this population.

Table I

**Frequency of Virus Specific Bands and Interpretation of Results of  
Urine Specimens from Non-AIDS HIV-1 Positive Populations  
(N=533<sup>a</sup>)**

Urine EIA	Urine WB		Frequency of Virus Specific Bands <sup>b</sup>								
			N	p17	p24	p31	gp41	p51	p55	p66	gp120
RR <sup>d</sup>	Pos <sup>f</sup>	533	136	301	319	462	317	118	433	512	533 <sup>e</sup>
N=533	Ind	0	0	0	0	0	0	0	0	0	0

<sup>a</sup> Two hundred and sixty two (262) specimens from Study 1, 197 from Study 2, 74 from Study 3.

<sup>b</sup> Band patterns for negative specimens do not appear in this table. By definition, negative specimens show no reactivity.

<sup>d</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>e</sup> Fifty (50) of the 533 urine specimens had only a gp160 or gp120 and gp160 band (37 from Study 1, 8 from Study 2, 5 from Study 3).

<sup>f</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

### Performance Using a Low Intensity Criterion

The performance of Western blot testing when using an intensity criterion for the gp160 band was evaluated by testing a subset of urine samples from known HIV-1 positive patients (Table J). One hundred and seventy-two (172) urine samples were retested. The results were interpreted by comparing the intensity of the gp160 band of the specimen with that of the LOW POSITIVE URINE CONTROL. The testing was performed, and results were interpreted at two different sites. A description of the known HIV-1 positive urine specimens that were retested and the results obtained at both sites are presented in Tables J and K.

The results show that the 2 sites read the same specimen as indeterminate. The sites differed on their interpretation for a second specimen (Site 1 read the sample as positive, Site 2 read the specimen as indeterminate). Therefore, the frequency of indeterminate results in a sample population of 172 HIV-1 EIA RR urine specimens ranged from 0.6% (1 of 172 specimens read at Site 1) to 1.2% (2 of 172 specimens read at Site 2).

The sensitivity of the Western blot using the intensity criterion associated with the LOW POSITIVE URINE CONTROL ranged from 98.3% (169/172) to 98.8% (170/172).

**Table J**  
**Description of HIV-1 Positive Specimens Retested**  
**On The Cambridge Biotech HIV-1 Western Blot Kit**

Study (N)	Study Group	Clinical Condition ((Samples Collected (N)))	N Tested
1 (100)	HIV-1 Seropositive Women	AIDS (26) HIV-1 Symptomatic (38) HIV-1 Asymptomatic (36)	25 <sup>a</sup> 37 <sup>b</sup> 35 <sup>c</sup>
2 (75)	AIDS Patients	AIDS (75)	75
Total (175)		175	172

- <sup>a</sup> One of the 26 AIDS patient specimens had insufficient volume to test by Western blot.
- <sup>b</sup> One of the 38 HIV-1 positive symptomatic patient specimens had insufficient volume to test by Western blot.
- <sup>c</sup> One of the 36 HIV-1 positive asymptomatic patient specimens had insufficient volume to test by Western blot.

**Table K**  
**Western Blot Results for Known HIV-1**  
**Positive Patients Using a gp160 Band Intensity Criterion.**

N	No Intensity Criterion			Using gp160 Band Intensity Criterion Western Blot Results					
	Prior <sup>a</sup> Western Blot Result			Site 1			Site 2 <sup>g</sup>		
Tested	Pos	Ind	Neg	Pos	Ind	Neg	Pos	Ind	Neg
172 <sup>b</sup>	171	N/A	1 <sup>c</sup>	170	1 <sup>d</sup>	1 <sup>e</sup>	169	2 <sup>f</sup>	1 <sup>e</sup>

- <sup>a</sup> "Prior" indicates data from Table E: analysis without an intensity criterion.
- <sup>b</sup> Three of the original 175 urine specimens had insufficient volume to test by Western blot.
- <sup>c</sup> Specimen was borderline EIA repeatedly reactive (S/CO range 0.972 to 1.916), negative on Western blot in Table E.
- <sup>d</sup> Specimen is +/- gp160/+gp120 only. The EIA testing associated with Western blot testing using Intensity Criterion was 1.284 S/CO. The specimen was previously urine Western blot false negative (Table E).
- <sup>e</sup> Specimen was urine EIA non-reactive and false negative (S/CO=0.825).
- <sup>f</sup> One Indeterminate was the specimen identified in footnote d.
- <sup>g</sup> The blots were interpreted approximately 7 days after processing.

## SPECIFICITY STUDIES

### Specificity in Low Risk Groups

The specificity of the Cambridge Biotech HIV-1 Western Blot Kit was assessed by testing specimens from 515 EIA seronegative subjects at low risk for HIV-1 infection. The subjects were insurance applicants. Insurance applicants are presumed to be at low risk for HIV-1 infection. The results obtained from testing paired urine and serum specimens from low risk uninfected individuals by Western blot are provided in Table L.

The results show that in this low risk population, the specificity of the Cambridge Biotech HIV-1 Urine Western Blot was 100% (515 of 515 urine specimens were Western blot negative). There were no (0%) urine Western blot indeterminate or positive specimens.

**Table L**

**Comparison of Cambridge Biotech HIV-1 Western Blot Results Using Urine and Paired Serum Specimens from Low Risk Populations (N=515)**

Risk Group	N	Serum			Urine			
		Western Blot Results <sup>c</sup>			EIA RR	Western Blot Results		
		Pos	Neg	Ind		Pos <sup>d</sup>	Neg	Ind
Low <sup>a</sup>	200 <sup>b</sup> 315	NT 0	NT 284	NT 31	0 1	0 0	200 315	0 0

<sup>a</sup> Two hundred (200) specimens from Study 1, 315 specimens from Study 3.

<sup>b</sup> The 200 archived serum specimens (Study 1) were EIA NR and were not tested by Western blot.

<sup>c</sup> A licensed HIV-1 Western Blot Kit was used when testing serum specimens.

<sup>d</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

**Frequency of Virus Specific Bands in Low Risk Groups**

The frequency of virus specific bands and interpretation of results for urine specimens tested with the Cambridge Biotech HIV-1 Western Blot Kit for these low risk groups are presented in Table M.

The results show that in these populations of 515 low risk subjects, none of the urine specimens were positive or indeterminate when tested with the Cambridge Biotech HIV-1 Western Blot Kit. All 515 serum specimens were EIA non reactive; 200 serum specimens were not tested by Western blot; of 315 serum specimens that were tested by Western blot, 284 were negative, 31 were indeterminate, and none were positive. All 515 urine specimens were identified as Western blot negative, including the one urine specimen that was EIA initially reactive. This demonstrates the high specificity of the Cambridge Biotech HIV-1 Western Blot Kit for urine in low risk populations.

Table M

**Frequency of Virus Specific Bands and Interpretation of Results of  
Urine Specimens from Low Risk Populations  
Tested by the Cambridge Biotech HIV-1 Western Blot Kit  
(N=515<sup>a</sup>)**

Urine EIA	Urine WB	N	Frequency of Virus Specific Bands <sup>b</sup>								
			p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
NR <sup>c</sup> N=514	Pos	0	0	0	0	0	0	0	0	0	0
	Ind	0	0	0	0	0	0	0	0	0	0
RR <sup>d</sup> N=1 <sup>e</sup>	Pos	0	0	0	0	0	0	0	0	0	0
	Ind	0	0	0	0	0	0	0	0	0	0

<sup>a</sup> Two hundred (200) specimens from Study 1, 315 specimens from Study 3.

<sup>b</sup> Band patterns for negative specimens do not appear in this table. By definition, negative specimens show no reactivity.

<sup>c</sup> NR indicates non-reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>d</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>e</sup> The one urine specimen initially reactive in the EIA was not repeat tested.

**Frequency of Virus Specific Bands in Other Groups**

Two additional special studies of paired urine and serum specimens were collected for evaluation by Western blot. The first evaluation involved Western blot testing of urine specimens paired to serum EIA non-reactive Western blot indeterminate specimens (N=109). The purpose of the study was to demonstrate the specificity of the urine Western blot for those samples that are repeatedly reactive on the urine EIA. The results are shown in Table N.

In this evaluation, 4 of 109 (3.7%) urine specimens paired to serum Western blot indeterminate specimens were urine Western blot indeterminate. None were positive. One hundred and five (105) of 109 (96.3%) were negative. These results show the specificity of the Cambridge Biotech HIV-1 Western Blot Kit for samples from uninfected individuals who are repeatedly reactive (false positive) on the urine EIA.

Table N

**Frequency of Virus Specific Bands and Interpretation of Results of  
Urine Specimens Paired to Serum Western Blot Indeterminate Specimens from  
Uninfected Individuals  
Tested by the Cambridge Biotech HIV-1 Western Blot Kit  
(N=109<sup>a</sup>)**

Urine EIA	Urine WB		Frequency of Virus Specific Bands								
		N <sup>d</sup>	p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
NR <sup>b</sup> N=66	Pos	0	0	0	0	0	0	0	0	0	0
	Ind	1	0	0	0	0	0	0	0	1	0
RR <sup>c</sup> N=43	Pos <sup>e</sup>	0	0	0	0	0	0	0	0	0	0
	Ind	3	0	2	0	0	0	0	1	0	0

<sup>a</sup> All 109 specimens were from Study 3.

<sup>b</sup> NR indicated non-reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>c</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

<sup>d</sup> Includes specimens with non-viral bands.

<sup>e</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

The second evaluation involved urine Western blot testing of urine EIA repeatedly reactive specimens paired to serum EIA non-reactive specimens from uninfected individuals (N=114). The purpose of this study was to assess the utility of the urine Western blot for individuals who are not infected but whose urine specimens are repeatedly reactive using the HIV-1 EIA. The frequency of virus specific bands in each group of urine specimens tested is demonstrated in the Table O.

Of the 114 urine EIA repeatedly reactive specimens paired to serum EIA non-reactive specimens, 109 of 114 (95.6%) were Western blot negative, 5 of 114 (4.4%) were Western blot indeterminate and none were positive. This demonstrates the ability of the Cambridge Biotech HIV-1 Western Blot Kit to resolve urine EIA repeatedly reactive specimens from uninfected individuals as negative or indeterminate.

Table O

**Frequency of Virus Specific Bands and Interpretation of Urine EIA False Positive Urine Specimens Tested by the Cambridge Biotech HIV-1 Western Blot Kit (N=114<sup>a</sup>)**

Urine EIA	Urine WB	N	Frequency of Virus Specific Bands								
			p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
RR <sup>b</sup> N=114	Pos	0	0	0	0	0	0	0	0	0	0
	Ind	5	0	5	0	0	0	0	0	0	0

<sup>a</sup> All 114 specimens were from Study 3.

<sup>b</sup> RR indicates repeatedly reactive in the Calypte Biomedical HIV-1 Urine EIA.

**Specificity in Subjects with Medical and Other Conditions**

The specificity of the Cambridge Biotech HIV-1 Western Blot Kit was evaluated using urine from individuals with medical conditions unrelated to HIV-1 and from individuals with potential interfering substances in their urine. The results of testing urine specimens from these individuals with the Cambridge Biotech HIV-1 Urine Western Blot Kit are provided in Table P.

The results obtained from testing urine specimens collected from patients with diseases and potentially interfering substances showed that 2 urine specimens were urine EIA repeatedly reactive and Western blot positive and were paired to serum EIA repeatedly reactive and Western blot positive specimens.

For subjects with non-HIV-1 medical conditions, the specificity of the urine Western blot for EIA repeatedly reactive urine specimens in this study was 90.1% (100/111).

Table P

**Comparison of Cambridge Biotech HIV-1 Western Blot Results Using Urine and Paired Serum in Populations with Disease and Potentially Interfering Substances (N=281)**

Group	N	Urine EIA Results		Urine Western Blot Results		
		RR	NR	Pos <sup>j</sup>	Neg	Ind
Autoimmune <sup>a</sup>	25	7	18	0	25	0
Kidney/Liver <sup>b</sup>	59	32	27	0	55	4
STD <sup>c</sup>	37	5	32	0	37	0
Urine Cond. <sup>d</sup>	47	22	25	1	45	1
Pregnant <sup>e</sup>	63	25	38	1	59	3
Neoplasms <sup>f</sup>	35	17	18	0	35	0
Multiple Transfusions <sup>g</sup>	13	5	8	0	10	3
Multiparous <sup>h</sup>	2	0	2	0	2	0
<b>Total</b>	<b>281</b>	<b>113</b>	<b>168</b>	<b>2<sup>m</sup></b>	<b>268</b>	<b>11<sup>k</sup></b>

<sup>a</sup> Twenty (20) specimens from Study 1, 5 specimens from Study 3.

<sup>b</sup> Twenty (20) specimens from Study 1, 39 specimens from Study 3.

<sup>c</sup> STD = Sexually transmitted disease. Twenty two (22) specimens from Study 1, 15 specimens from Study 3.

<sup>d</sup> Twenty (20) specimens from Study 1, 27 specimens from Study 3.

<sup>e</sup> Sixty three (63) specimens from Study 3.

<sup>f</sup> Thirty five (35) specimens from Study 3.

<sup>g</sup> One (1) specimen from Study 2, 12 specimens from Study 3.

<sup>h</sup> Two (2) specimens from Study 3.

<sup>j</sup> No intensity criterion associated with the LOW POSITIVE URINE CONTROL was used during this analysis.

<sup>k</sup> These specimens were all urine EIA repeatedly reactive.

<sup>m</sup> These two specimens were true serum positive as confirmed.

## ANALYTICAL SENSITIVITY

### Dilution of Paired Urine and Serum Specimens

Ten (10) paired urine and serum specimens from HIV-1 positive individuals were tested in a dilution study using the Cambridge Biotech HIV-1 Western Blot Kit. The urine specimens were tested at dilutions of 1:2, 1:10, 1:50, 1:500 and 1:1,000. The serum specimens were tested at dilutions of 1:101, 1:1,000, 1:5,000, 1:25,000 and 1:50,000. The urine specimens were tested according to the urine procedure. The serum specimens were tested according to the manufacturer's package insert instructions. The major viral bands (gp160, gp120, gp41 and p24) were observed for presence or absence on Western blot at each dilution. The results are reported as the last dilution at which a band is visible on Western blot without comparison to an intensity criterion.

A side by side comparison of the urine and serum dilutions are presented in Table Q. The comparison shows that the gp160 band can be observed at a higher dilution than the gp120, gp41 or p24 bands for both urine and serum specimens tested in dilution series.

The difference in the analytical sensitivity between the licensed Serum Western Blot Kit and the Urine Western Blot Kit for antibodies to different proteins ranged from a ratio of greater than 50 to greater than 25,000.

**Table Q**  
**Comparison of Dilutions of HIV-1 Positive Paired Urine & Serum Specimens Tested on Cambridge Biotech Corporation's HIV-1 Western Blot Kit**

**Last Dilution At Which Major Viral Band is Present**

**Major Viral Bands<sup>^</sup>**

No.	Specimen No.	Clinical Classification*	gp160		gp120		gp41		p24	
			Urine	Serum	Urine	Serum	Urine	Serum	Urine	Serum
1	SP00000604	AIDS	1:500	>1:50,000	1:10	>1:50,000	1:10	>1:50,000	-	-
2	SP00000606	AIDS	1:50	>1:50,000	1:2	>1:50,000	1:2	>1:50,000	-	1:1,000
3	SP00001709	HIV-1 Symptomatic	1:10	>1:50,000	-	>1:50,000	-	>1:50,000	-	1:10,000
4	SP00001752	HIV-1 Symptomatic	1:50	>1:50,000	1:2	>1:50,000	1:2	>1:50,000	1:10	>1:50,000
5	SP00001783	HIV-1 Asymptomatic	1:50	>1:50,000	1:2	>1:50,000	-	1:10,000	-	1:101
6	CL1	F-TL	1:50	>1:50,000	-	>1:50,000	-	>1:50,000	-	>1:50,000
7	CL35	F-TL	1:10	>1:50,000	-	1:25,000	-	1:10,000	-	1:10,000
8	CL83	F-TL	1:1,000	>1:50,000	1:50	>1:50,000	1:10	>1:50,000	>1:1,000	>1:50,000
9	175	F-TZ	1:500	>1:50,000	1:50	>1:50,000	1:10	>1:50,000	1:1,000	>1:50,000
10	3.030751	F-Abidj	1:10	>1:50,000	-	>1:50,000	-	1:25,000	1:10	>1:50,000

\* F-TL = Foreign specimen from Thailand

\* F-TZ = Foreign specimen from Tanzania

• F-Abidj = Foreign specimen from Abidjan

<sup>^</sup> (-) indicates these bands were not present for the original serum specimen diluted 1:101 according to the serum Western Blot procedure and for the original urine specimen diluted 1:2 according to the Western Blot procedure for urine.

## REPRODUCIBILITY

The reproducibility of the Cambridge Biotech HIV-1 Western Blot Kit was evaluated by testing a panel of urine specimens at three (3) geographically separate sites. A panel of 12 specimens of defined viral reactivity was provided to the three sites for evaluation.

The panel consisted of specimens strongly reactive and weakly reactive for antibodies to HIV-1 and specimens non-reactive for antibodies to HIV-1. Each panel member was tested in duplicate on three different lots of the Cambridge Biotech HIV-1 Western Blot Kit. The testing was performed by at least two different operators at each of 3 sites over multiple days.

In addition to the 12 specimens, a high positive urine control, low positive urine control and a negative urine control were provided for testing with each kit.

The results of this analysis demonstrate the reproducibility of the Cambridge Biotech HIV-1 Western blot for urine specimens with HIV-1 antibody activity to the gp160 viral gene product at the limit of visual detection. The combined results of testing are provided in Table R.

Table R

**Number of Replicates of Western Blots With Reactive Bands  
(% Reactive Replicates)**

Spec/ Res <sup>a</sup>	# of Reps <sup>h</sup>	p17	p24	p31	gp41	p51	p55	p66	gp120	gp160
HPC <sup>b</sup>	24	23 (96)	24 (100)	24 (100)	24 (100)	15 (63)	13 (54)	24 (100)	24 (100)	24 (100)
LPC <sup>c</sup>	24	3 (13)	24 (100)	23 (96)	17 (71)	6 (25)	0 (0)	21 (88)	24 (100)	24 (100)
NC <sup>d</sup>	24	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)
1 Neg <sup>e</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)
2 Neg <sup>e</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)
3 S. Pos <sup>f</sup>	48	16 (33)	2 (4)	48 (100)	20 (42)	28 (58)	18 (38)	48 (100)	48 (100)	48 (100)
4 S. Pos <sup>f</sup>	40	40 (100)	40 (100)	40 (100)	40 (100)	34 (85)	22 (55)	40 (100)	40 (100)	40 (100)
5 W. Pos <sup>g</sup>	40	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	14 (35)	0 (0)	40 (100)	40 (100)
6 W. Pos <sup>g</sup>	40	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	40 (100)	40 (100)
7 W. Pos <sup>g</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	28 (58)	48 (100)	48 (100)
8 W. Pos <sup>g</sup>	48	0 (0)	46 (96)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	48 (100)	48 (100)
9 W. Pos <sup>g</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	2 (4)	48 (100)	48 (100)
10 W. Pos <sup>g</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	48 (100)	48 (100)
11 W. Pos <sup>g</sup>	48	0 (0)	24 (50)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	48 (100)	48 (100)
12 W. Pos <sup>g</sup>	48	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	0 (0)	48 (100)	48 (100)

<sup>a</sup> Spec/Res=Specimen ID#/Expected Result

<sup>b</sup> HPC=High Positive Urine Control

<sup>c</sup> LPC=Low Positive Urine Control

<sup>d</sup> NC=Negative Urine Control

<sup>e</sup> Neg=Negative Specimen

<sup>f</sup> S Pos=Strong Positive Specimen

<sup>g</sup> W Pos=Weak Positive Specimen

<sup>h</sup> Repls=Replicates

VII. PACKAGE INSERT

The package insert for the Cambridge Biotech HIV-1 Western Blot Kit is attached.

SUMMARY BASIS OF APPROVAL  
Reference #95-1588


Signatures:

  
\_\_\_\_\_  
Kori Francis, Co-Chairperson  
Licensing Review Committee

4-30-98  
Date

  
\_\_\_\_\_  
Indira Hewlett, Ph.D., Co-Chairperson  
Licensing Review Committee


4/27/98  
Date

  
\_\_\_\_\_  
Edward Tabor, M.D., Director  
Division of Transfusion Transmitted Diseases

4/17/98  
Date

  
\_\_\_\_\_  
Michael Calabro, Ph.D., Regulatory Coordinator  
Division of Blood Applications

5/7/98  
Date

  
\_\_\_\_\_  
Mary L. Gustafson, Director  
Division of Blood Applications

5-26-98  
Date