



**U.S. FOOD & DRUG
ADMINISTRATION**

Bacteriological Analytical Manual
Chapter 19a: Detection of Cyclospora
and Cryptosporidium from Fresh
Produce: Isolation and Identification
by Polymerase Chain Reaction (PCR)
and Microscopic Analysis
October 2004 Edition

Table of Contents

Authors	3
1. Materials and Equipment	3
2. Reagents	3
3. Optional Reagents, Supplies, and Equipment	4
4. Wash Procedure for Fresh Produce	4
5. Isolation of Parasitic Contaminants from Juices, Cider and Milk	6
6. Isolation of Parasitic Contaminants from Large Volumes of Water	6
7. Slide Preparation and Microscopic Analysis – Cyclospora cayetanensis	6
8. Slide Preparation and Microscopic Analysis & Cryptosporidium spp.	7
9. PCR Analysis	8
References	15

Note

Updated methods are available:

- (2017) Chapter 19b: [Detection of Cyclospora cayetanensis in Fresh Produce using real-time PCR](#)
- (2020) Chapter 19c: [Dead-end Ultrafiltration for the Detection of Cyclospora cayetanensis from Agricultural Water](#)

Authors

Palmer A. Orlandi (ret), Christian Frazar, Laurenda Carter, and Dan-My T. Chu
 Contact: [Sonia Almeria](#) or [Mauricio Durigan](#)

1. Materials and Equipment

- a. BagPage[®]+ filter bags (400 ml) and Bag Clips (*Interscience*, St Nom, France)
- b. Envirochek[™] sampling capsule, 1 µm nominal (Pall Gelman Laboratory)
- c. Rocker platform
- d. Rotary shaker
- e. Rotating wheel
- f. 150 ml analytical filter units (Nalgene, Cat No. 130-4045)
- g. 25 mm disposable filter funnels (Whatman Biosciences)
- h. Dynal L10 tubes Prod. No. 740-03)
- i. Dynal MPC[®]-1 (Prod. No 120.01)
- j. Dynal MPC[®]-S (Prod. No 120.20)
- k. Vacuum manifold (Vac-Man[®] Laboratory Vacuum Manifold 20-sample capacity; Promega)
- l. FTA Filters, Classic Card of 960circle pre-printed card formats (Whatman Biosciences)
- m. Photo laminating sheets (Scotch[®])
- n. Conical 250 ml centrifuge tubes
- o. Sorvall RT7 Plus refrigerated centrifuge or equivalent (to centrifuge 250 ml conical centrifuge tubes)
- p. Epi-illuminated Fluorescence Microscope-equipped with the following: UV 1A filter block (excitation Filter, EX 365/10; Dichroic mirror, DM 440; Barrier Filter, BA-400; or equivalent. Optics for differential interference contrast (DIC). For *Cryptosporidium* spp.: appropriate filter for viewing fluorescein isothiocyanate (FITC) conjugated oocysts. (*Table 1*)
- q. Glass microscope slides and cover slips
- r. Blotting paper
- s. Heating block for incubation at 56°C
- t. Single punch (6 mm diameter) hole puncher
- u. Thin-walled 0.65 ml PCR tube (PGC)
- v. ART pipette tips (Molecular BioProducts)
- w. PTC-200 DNA Engine (MJ Research) or comparable thermal cycler.
- x. Horizontal gel electrophoresis apparatus and power supply.
- y. Polaroid camera or digital imaging system to capture ethidium bromide-stained gels
- z. Polaroid Type 667 film
- aa. UV Transilluminator

2. Reagents

- a. Water
 - i. Deionized water (for washing produce [dH₂O])
 - ii. Sterile deionized water (for PCR procedure)
- b. Envirochek[™] elution buffer (0.01M Tris, pH 7.4 ; 0.001 M EDTA; 1% SDS)

- c. Silicone vacuum grease
- d. Albumin, bovine (BSA) (Sigma, A-7030)
- e. Celite (Sigma, C-8656)
- f. Polyvinyl polypyrrolidone (PVPP) (Sigma, P-6755)
- g. NET buffer: 0.1 M Tris, pH 8.0, 0.15M NaCl, 0.001M EDTA
- h. NET-BSA buffer-NET buffer containing 1% (w/v) BSA
- i. 20% celite in NET-BSA buffer (w/v)
- j. 10% PVPP (w/v) in dH₂O
- k. Dynabeads anti-Cryptosporidium Kit (DynaL Biotech Inc, USA)
- l. Hydrofluor Combo Detection Kit for *Cryptosporidium* and *Giardia* (Strategic Diagnostics, Inc.)
- m. 0.1 N HCl
- n. 0.1 N NaOH
- o. Immersion oil
- p. Clear fingernail polish, slide compound, paraffin wax or equivalent
- q. FTA Purification Buffer (Whatman Biosciences)
- r. FTA filter wash buffer: 0.01 M Tris, Ph 8.0; 0.1 mM EDTA
- s. DNA Primers-See Table 1 in PCR Section
- t. HotStartTaq™ Master Mix Kit (Qiagen)
- u. 0.5x Tris-acetate-EDTA buffer (0.5x TAE)
- v. Molecular biology-grade agarose (BioRad)
- w. Ethidium bromide
- x. 6x gel loading solution
- y. 100 bp and 25 bp DNA ladders (Invitrogen)
- z. *Vsp I* restriction endonuclease (Promega)
- aa. *Dra II* restriction endonuclease (Hoffman-La Roche)
- bb. NuSieve® 3:1 agarose (Biowhitaker Molecular Applications)

3. Optional Reagents, Supplies, and Equipment

- a. LightCycler®-FastStart DNA Master SYBR Green I kit (Roche Diagnostics)
- b. LightCycler® Capillaries (Roche Diagnostics)
- c. LightCycler® System (Roche Diagnostics)

4. Wash Procedure for Fresh Produce

This procedure can be used to analyze for potential contaminants on fresh leafy produce (lettuce, herbs, etc) and berries (e.g. raspberries) and may be applicable to other fresh produce.

- a. Place produce to be analyzed (10-25 g of fresh leafy produce or 50g of fresh berries) in a BagPage®+ filter bag, add 100 ml dH₂O and seal with Bag Clip.
- b. Place on rocker platform and gently agitate for 30 minutes at room temperature, inverting the bag after 15 minutes.
- c. Decant supernatant from the BagPage®+ filter bag into clean 50 ml conical centrifuge tubes and centrifuge for 20 minutes at 2,000xg.

Isolation of *Cyclospora* from Fresh Produce Washes

- a. Aspirate supernatants (without disturbing debris pellets) to a volume not to exceed 45 ml when combined.
- b. Suspend pellets in remaining supernatants and combine.
- c. Add 2.5 ml 20% (w/v) celite in NET-BSA suspension (ensure that celite is thoroughly suspended prior to its addition to sample washes). Mix samples on a rotating wheel at room temperature for 15 minutes.
- d. Add 1.0 ml of 10% PVPP suspension (w/v). Mix samples on a rotating wheel at room temperature for 15 minutes.
- e. Prepare a 150 ml analytical filter unit by removing the membrane filter (0.45 µm or 0.2 µm) but retaining the grade 4 filter backing. Attach to a vacuum source.
- f. Pre-wet the analytical filter with a small volume (~10ml) of NET buffer.
- g. Decant celite/PVPP-containing sample wash into a prepared 150 ml analytical filter unit and vacuum filter. Ensure that the liquid passes through the filter septum to remove particulates and the celite particles from the suspension as it is filtered. The adsorbant celite should prevent filter clogging.
- h. Rinse the container with 10 ml NET to recover as much of the residual celite-containing sample and decant into the filter unit. Then rinse the celite and particulate material trap onto the filter with an additional 10 ml volume of NET.
- i. Prior to FTA filtration in **step 4i**, save ~10% of the filtered sample for microscopic examination.
- j. Prepare filter funnel unit(s) containing FTA filter disk and attach to vacuum manifold.
- k. Under vacuum, pre-wet the FTA filter assembly.
- l. Slowly decant filtrate from **step 4i** into filter funnel unit while under vacuum until entire sample has passed through FTA filter.
- m. In succession while filter funnel unit is still attached to vacuum manifold, rinse filter twice with 10 ml of FTA purification buffer and twice with 10 ml of FTA filter wash buffer.
- n. Remove filter funnel unit from vacuum manifold, disassemble unit and dry FTA filter disk on 56°C heating block.

Isolation of *Cryptosporidium* spp. from Fresh Produce Washes

- a. Aspirate supernatants (without disturbing debris pellets) to a volume not to exceed 10 ml when combined.
- b. Suspend pellets in remaining supernatants and combine. **NOTE:** pellet volume should not exceed 0.5 ml packed volume as it will interfere with subsequent steps that employ immunomagnetic adsorption of *Cryptosporidium* oocysts.
- c. Follow directions accompanying Dynabead anti-*Cryptosporidium* Kit for immunomagnetic separation (IMS) of *Cryptosporidium* oocysts from fresh produce washes using recommend tubes and magnetic capturing devices.
- d. Elute captured *Cryptosporidium* oocysts in a 0.1 ml volume of 0.1N HCl for 5 minutes at room temperature.
- e. Neutralize acid eluates with 0.01 ml 1N NaOH. Save ~10% of the sample for microscopic examination.
- f. Dilute with 10 ml NET buffer.
- g. Proceed as in **Section 4A, steps j-n.**

5. Isolation of Parasitic Contaminants from Juices, Cider and Milk

This procedure can be used to analyze for contaminants in liquid samples such as orange juice, apple juice, apple cider, milk and milk products.

Isolation of *Cyclospora* from Juices, Cider and Milk

- a. Adjust a 25 ml aliquot of liquid sample to pH 8.0.
- b. Add an equal volume of NET buffer and mix well.
- c. Proceed as in **Section 4A, steps c-n**

Isolation of *Cryptosporidium* spp. from Juices, Cider and Milk

A 10 ml volume of cider, juice or milk product is directly sampled by IMS using the Dynabead anti-*Cryptosporidium* Kit along with the recommend tubes and magnetic capturing devices as in **Section 4B, steps c-f**.

6. Isolation of Parasitic Contaminants from Large Volumes of Water

This procedure is designed to isolate contaminants from a designated water source (stream, river, reservoir, standing water, runoff, etc).

- a. Place a single Envirochek™ sampling capsule in line with the water source to be sampled. Ensure that the flow rate does not exceed the standards established by the manufacturer.
- b. Collect a water sample from 10 L of flow through.
- c. Remove contaminants captured by the filter using 125 ml of elution buffer. Filters should be treated with elution buffer, sealed and agitated on a rotating wheel (moderate speed) for at least 5-10 minutes following the manufacturer's recommendations. Decant filter rinse into a 250 ml conical centrifuge tube.
- d. Repeat **step 6c** and combine filter rinses.
- e. Centrifuge the sample at 1,500-2,000 x g for 20 minutes allowing the centrifuge to coast to a stop (do not use the brake)
- f. Proceed as described in **Sections 4A and 4B** for the isolation of *Cyclospora* and *Cryptosporidium* spp. respectively.

7. Slide Preparation and Microscopic Analysis – *Cyclospora cayetanensis*

Cyclospora oocysts emit a cobalt blue autofluorescence with the UV-1A emission filter or blue-green with broader emission spectra filters under ultraviolet illumination. Prepare slides in duplicate and examine slides under ultraviolet illumination as described below.

Laboratories should use a microscope reticle capable of measuring organism in the 8-10 μm range to verify oocyst size when organisms are recovered. Compare morphological characteristics of presumptive oocysts to those in a known standard.

Slide Preparation

Centrifuge the volume set aside for microscopic analysis (**Section 4A, step i**), at $1,500 \times g$ for 10-15 min at 4°C . Aspirate supernatant to within 0.5 ml of debris pellet. Uniformly suspend pellet material by gentle, repeated pipetting.

- a. Apply silicone vacuum grease to edge of cover slip
- b. Place 10 μl of suspended debris to a clean glass microscope slide and prepare a wet mount using pre-greased cover slip

Microscopy

- a. Examine slide under UV light at $400\times$ magnification. *Cyclospora* oocysts autofluoresce cobalt blue. **Examine slide at multiple planes under UV. Occasionally, debris in slide preparations make it difficult to only view the slide on one plane.**
- b. Switch from epifluorescence microscopy to bright field microscopy or differential interference contrast microscopy. Determine oocyst size of any presumptive oocysts at $1000\times$ magnification. Compare to standards. Confirm internal structures of presumptive *Cyclospora* oocysts as compared to a standard.
- c. Seal cover slips to the glass slides of presumptive positives with fingernail polish, slide compound or paraffin wax.
- d. Document presumed positive samples with photographs taken at multiple planes.

8. Slide Preparation and Microscopic Analysis & *Cryptosporidium* spp.

Microscopic examination of produce washings and other liquid samples for the presence of *Cryptosporidium* spp. oocysts will be conducted using commercially available immunomagnetic bead separation (IMS) kits and immunofluorescence labeling kits.

- A. An aliquot ($\sim 10\mu\text{l}$) of the IMS-derived sample **Section 4B, step e**, is FITC-labeled using the Hydrofluor Combo Detection Kit for *Cryptosporidium* and *Giardia* (Ensys Inc., Research Triangle Park, NC) per the manufacturer's instructions and examined in conventional DIC and epifluorescence mode.

Table 1: Parameters for Epifluorescence Microscopy†

INCIDENT LIGHT
Light Source-Mercury Vapor
200W, 100W, or 50W

Excitation Filter	Dichroic Filter	Barrier Filter	Red Suppression Filter
KP500	TK510	K510 or K530	BG38
FITC	TK510	K530	BG38
D. Tungsten-Halogen 50 and 100 W			
KP500	TK510	K510 or K530	BG38
FITC	TK510	K530	BG38

†Taken from protocol provided with HYDROLFUOR-Combo Detection Kit for Giardia cysts and Cryptosporidium oocysts (Ensys Inc., Research Triangle Park, NC)

9. PCR Analysis

The molecular detection of *Cyclospora* spp. and *Cryptosporidium* spp. is independently accomplished using nested PCR protocols. The differential identification of *Cyclospora cayetanensis* from other closely related non-human pathogenic parasites (*i.e.* *Eimeria* spp.) employs a nested multiplex PCR assay. This assay can be accomplished using either a conventional thermal cycler with heated lid or a real-time PCR platform using the Roche LightCycler®.

The detection of *Cryptosporidium* spp. also involves nested PCR amplification. However, differentiation and speciation of *Cryptosporidium* spp. requires further analysis by a restriction fragment length polymorphism (RFLP) assay. Please note the following: *C. parvum* genotype I has been renamed *C. hominis*; *C. parvum* genotype II (bovine strain) is now referred to as *C. parvum*.

A. DNA Primers

Table 2: DNA Primer Sequences for Cyclospora-specific PCR Amplification†

Primer Designation	Primer Specificity	Primer Sequence (5'-3')	Amplicon Size (bp)	Designated Application
F1E (forward)	<i>Cyclospora</i> and <i>Eimeria</i> spp.	TACCCAATGAAAACAGTTT	636	Primary Amplification
R2B (reverse)	<i>Cyclospora</i> and <i>Eimeria</i> spp.	CAGGAGAAGCCAAGGTAGG	636	Primary Amplification
CC719 (forward)	<i>C. cayetanensis</i>	GTAGCCTCCGCGCTTCG	298	Nested Amplification

PLDC661 (forward)	<i>C. cercopithecii</i> , <i>C. colobi</i> , <i>C. papionis</i>	CTGTCGTGGTCATCGTCCGC	361	Nested Amplification
ESSP841 (forward)	<i>Eimeria spp.</i>	GTTCTATTTTGTGGTTTCTAGGACCA	174	Nested Amplification
CRP999 (reverse)	<i>Cyclospora</i> and <i>Eimeria spp.</i>	CGTCTCAAACCCCCTACTGTCTG		Nested Amplification

†All primer sequences were derived from the published sequences for the 18S rRNA genes of the respective organisms.

Table 3: DNA Primer Sequences for Cryptosporidium Genus-specific PCR Amplification†

Primer Designation	Primer Specificity	Primer Sequence (5'-3')	Amplicon Size (bp)	Designated Application
ExCry1 (forward)	<i>Cryptosporidium</i> spp.	GCCAGTAGTCATATGCTTGTCTC	844	Primary Amplification
ExCry2 (reverse)	<i>Cryptosporidium</i> spp.	ACTGTAAATAGAAATGCCCCC	844	Primary Amplification
NesCry3 (forward)	<i>Cryptosporidium</i> spp.	GCGAAAAAACTCGACTTTATGGAAGGG	590-593	Nested Amplification
NesCry4 (reverse)	<i>Cryptosporidium</i> spp.	GGAGTATTCAAGGCATATGCCTGC	590-593	Nested Amplification

†All primer sequences were derived from the published sequences for the 18S rRNA genes of the respective organisms.

B. General Sample Preparations for Primary PCR Amplifications

- Punch marked triplicate areas (6 mm diameter) from dried FTA filter disk using a single punch hole-puncher. **Decontamination of the hole-puncher is not necessary as cross contamination between samples from the hole-puncher has been found to be negligible. However, the individual researcher may wish to swab the punch with ethanol between sample disks if it is deemed appropriate.**
- Insert filter punches snugly into bottom of 0.65 ml thin-walled PCR tubes.
- Dispense 50 µl HotStartTaq™ Master Mix into each PCR tube.
- Prepare reagent master mix (see *Table 4*) with the appropriate forward and reverse DNA primers (see *Tables 2 and 3*) and dispense into each PCR tube.
- All PCR analyses must include positive and negative controls (see *Table 5*).
- Mix tubes with gentle tapping.
- Follow the appropriate thermal cycling protocol (*Table 6 and 7*) for primary PCR amplification.

Table 4: General PCR Conditions for Primary PCR Amplification

Component		Volume (µl)*	Final Concentration
FTA Filter Disk (DNA Template)			
HotStartTaq™ Master Mix		50.0	†
Reagent Master Mix	MgCl ₂ , 25 mM	2.0	2.0 ‡
	Forward Primer, 10 µM	2.0	0.2 µM
	Reverse Primer, 10 µM	2.0	0.2 µM
	Sterile deionized water	44.00	

*100 µl total volume

†Final concentrations for components in the HotStartTaq™ Master Mix are as follows: 200 µM of each dNTP, 1.5 mM MgCl₂ and 2.5 U HotStartTaq™ DNA Polymerase

‡Final MgCl₂ concentration is that contributed by both the HotStartTaq™ Master Mix and 25 mM MgCl₂stock

Table 5: Controls for PCR amplifications

Control Type	Condition/Organism	
Negative Control-1	Reagent blank-no filter	
Negative Control-2	Reagent blank + unspotted, washed filter	
†Positive Controls:	<i>Cyclospora</i> Analysis:	<i>C. cayetanensis</i>
†Positive Controls:	<i>Cyclospora</i> Analysis:	‡† <i>Cyclospora</i> spp. (NHP)
†Positive Controls:	<i>Cyclospora</i> Analysis:	* <i>Eimeria</i> spp.
†Positive Controls:	<i>Cryptosporidium</i> Analysis:	† <i>C. hominis</i> (formerly know as <i>C. parvum</i> genotype 1)
†Positive Controls:	<i>Cryptosporidium</i> Analysis:	<i>C. parvum</i> , (formerly know as <i>C. parvum</i> genotype II (Bovine))
†Positive Controls:	<i>Cryptosporidium</i> Analysis:	† <i>C. baileyi</i>
†Positive Controls:	<i>Cryptosporidium</i> Analysis:	† <i>C. serpentis</i>

†Whenever possible, positive control FTA filters should be spotted with at least 10³ organisms

‡Non-human primate-derived *Cyclospora* spp.

†Not routinely available.

*Most available *Eimeria* spp. are suitable.

Table 6: PCR Thermal Cycling Parameters for Cyclospora and Eimeria spp.

	Step	Number of Cycles	Temperature and Time
Primary PCR	Initial Activation	1	95°C; 15 min
	Amplification	35	Denaturation: 94°C; 30 sec
	Amplification	35	Annealing: 53°C; 30 sec
	Amplification	35	Extension: 72°C; 90 sec
	Final Extension	1	72°C; 10 min
†Nested Multiplex PCR	Initial Activation	1	95°C; 15 min
	Amplification	25	Denaturation: 94°C; 15 sec
	Amplification	25	Annealing: 66°C; 15 sec

†This is a stringent, 2-step amplification program (simultaneous annealing and extension at 66°C). Likewise, it does not require a final extension step.

Table 7: PCR Thermal Cycling Parameters for Cryptosporidium spp.

	Step	Number of Cycles	Temperature and Time
Primary PCR	Initial Activation	1	
	Amplification	40	Denaturation: 94°C; 45 sec
	Amplification	40	Annealing: 53°C; 75 sec
	Amplification	40	Extension: 72°C; 45 sec
	Final Extension	1	72°C; 7 min
Nested PCR	Initial Activation	1	95°C; 15 min
	Amplification	35	Denaturation: 94°C; 25 sec
	Amplification	35	Annealing: 65°C; 25 sec
	Amplification	35	Extension: 72°C; 25 sec
	Final Extension	1	72°C; 7 min

Conventional Nested Multiplex PCR Amplification for the Differential Identification of Cyclospora and Eimeria spp.

1. Dispense 25 µl of HotStartTaq™ Master Mix into each PCR tube.
2. Prepare reagent master mix (*Table 8*) and dispense into all tubes
3. Complete reaction samples with the addition of the desired volume (1-3 µl) of primary amplicon solution.
4. Be sure to include all positive and negative controls as in the primary amplification reactions.
5. Mix tubes with gentle tapping.
6. Follow the appropriate thermal cycling protocol listed in *Table 6*.

Table 8: Assay Conditions for the Conventional Nested Multiplex PCR Amplification of Cyclospora and Eimeria spp.

Component	Volume (µl)*	Final Concentration
HotStartTaq™ Master Mix	25.0	†
MgCl ₂ , 25 mM	1.0	2.0 mM‡
CC719 (forward primer), 10 µM	1.0	0.2 µM
PDCL661 (forward primer), 10 µM	1.0	0.2 µM
ESSP841 (forward primer), 10 µM	1.0	0.2 µM
CRP999 (reverse primer), 10 µM	1.0	0.2 µM
Sterile deionized water	19.00	
DNA Template (primary amplicon)	1.0	

*50 µl total volume

†Final concentrations for components in the HotStartTaq™ Master Mix are as follows: 200 µM of each dNTP, 1.5 mM MgCl₂ and 2.5 U HotStarTaq™ DNA Polymerase

‡Final MgCl₂ concentration is that contributed by both the HotStartTaq™ Master Mix and 25 mM MgCl₂ stock.

OPTIONAL-Real Time Multiplex PCR Amplification for the Differential Identification of Cyclospora and Eimeria spp. using the Roche LightCycler® System

- a. Place the requisite number of glass capillaries in a pre-chilled cooling block (with accompanying centrifuge adapters)
- b. Prepare a reagent master mix (*Table 9*) and dispense into individual 0.65 ml PCR tubes

- c. Complete reaction samples with the addition of the 1 µl of primary amplicon solution.
- d. Be sure to include all positive and negative controls as in the primary amplification reactions.
- e. Mix tubes with gentle tapping.
- f. Dispense reaction mixture into glass capillaries, cap, and centrifuge briefly (5-10 sec at 3000 rpm) in a bench-top microcentrifuge.
- g. Transfer glass capillaries to LightCycler® carousel.
- h. Follow the appropriate thermal cycling protocol for (*Table 10*).
- i. Real time confirmation of pathogen can be made by melting curve analysis (*see Table 11*)
- j. For final confirmation, samples can be recovered from each glass capillary.
- k. Uncap, invert capillary into a 0.65 ml PCR tube containing 5 µl of gel loading solution, and briefly centrifuge (5-10 sec at 3000 rpm) in a bench-top microcentrifuge.
- l. Following instruction for agarose gel electrophoresis (**Part 9, Section F, steps b-f**)

Table 9: Assay Conditions for the Real-time Nested Multiplex PCR Amplification of Cyclospora and Eimeria spp.

Component	Volume (µl)*	Final Concentration
LightCycler®-FastStart DNA Master SYBR Green	2.0	-
MgCl ₂ , 25 mM	1.6	3.0 mM
CC719 (forward primer), 10 µM	1.0	0.5 µM
PDCL661 (forward primer), 10 µM	1.0	0.5 µM
ESSP841 (forward primer), 10 µM	1.0	0.5 µM
CRP999 (reverse primer), 10 µM	1.0	0.5 µM
Sterile deionized water	11.40	-
DNA Template (primary amplicon)	1.0	-

*20 µl total volume

Table 10: LightCycler® Thermal Cycling Parameters for Real-time Multiplex PCR Amplification of Cyclospora and Eimeria spp.

Step	Number of Cycles	Temperature and Time	Comments
Hot Start	1	95°C; 10 min	
Amplification	30	Denaturation: 95°C; 15 sec	
	30	Annealing: 66°C; 15 sec	Single Fluorescence Acquisition
Melting Curve Analysis	1	95°C; 15 sec	
	1	65°C; 15 sec	
	1	98°C; 0.1°C/sec	Continuous Fluorescence Acquisition

Table 11: Expected Results for Real-time Multiplex PCR Amplification of Cyclospora and Eimeria spp. by Melting Curve Analysis

Primer Designation	Primer Specificity	Primer Sequence (5'-3')	Amplicon Size (bp)	Amplicon T _m (°C)
CC719	<i>C. cayetanensis</i>	GTAGCCTTCCGCGCTTCG	298	85°C
PLDC661	<i>C. cercopitheci</i> , <i>C. colobi</i> , <i>C. papionis</i>	CTGTTCGTGGTCATCGTCCGC	361	91°C
ESSP841	<i>Eimeria</i> spp.	GTTCTATTTTGTGGTTTCTAGGACCA	174	81°C

Nested PCR Amplification for the Differential Identification of *Cryptosporidium* spp.

- Dispense 25 µl of HotStartTaq™ Master Mix into each PCR tube.
- Prepare reagent master mix (Table 12) and dispense into all tubes
- Complete reaction samples with the addition of the desired volume (1-3 µl) of primary amplicon solution.
- Be sure to include all positive and negative controls as in the primary amplification reactions.
- Mix tubes with gentle tapping.
- Follow the appropriate thermal cycling protocol listed in Table 7.

Table 12: Assay Conditions for Nested Amplification of *Cryptosporidium* spp.

Component	Volume (µl)*	Final Concentration
HotStartTaq™ Master Mix	25.0	†
MgCl ₂ , 25 mM	1.0	2.0 ‡
NesCry3 (forward primer), 10 µM	1.0	0.2 µM
NesCry4 (reverse primer), 10 µM	1.0	0.2 µM
Sterile deionized water	21.0	
DNA Template (primary amplicon)	1.0	

*50 µl total volume

†Final concentrations for components in the HotStartTaq™ Master Mix are as follows: 200 µM of each dNTP, 1.5 mM MgCl₂ and 2.5 U HotStarTaq™ DNA Polymerase

‡ Final MgCl₂ concentration is that contributed by both the HotStartTaq™ Master Mix and 25 mM MgCl₂ stock

Agarose Gel Electrophoresis

- Mix 10 µl of nested amplification product with 2-3 µl of gel loading solution.
- Load samples into wells of a 1.5% agarose gel prepared with 0.5 x TAE containing 0.2 µg/ml ethidium bromide. Include at least one lane containing 100 bp DNA ladder to approximate the size of any amplicon present.
- Run the gel at 125 volts (constant voltage) for at least 30 min.
- PCR products on the agarose gel can be visualized by using a UV transilluminator. Photograph the gel to have a permanent record of the results using a Polaroid Type 667 film (or a digital system, if you decide to include that in the material & methods).
- The primary amplicon from primer pair F1E/R2B for *Cyclospora* PCR may not be visible; therefore, only product from the nested reaction should be electrophoresed.
- Predicted sizes of PCR amplicons from *Cyclospora* spp., *Eimeria* spp. and *Cryptosporidium* spp. are listed in Tables 1 and 2

Restriction Fragment Length Polymorphism (RFLP) Analysis of *Cryptosporidium* spp. Nested PCR Amplicons to Determine the Presence of *C. parvum* Oocysts and Distinguish *C. hominis* from *C. parvum*

- a. A 590 bp (genotype I) or a 593 bp product (genotype II) following nested PCR is a presumptive positive for the presence of *Cryptosporidium* spp.
- b. The restriction patterns resulting from the digestion of the nested amplicon with restriction endonucleases *Vspl* and *Drall* will distinguish between *C. parvum* and *C. hominis* (*Vspl* digestion) and *C. parvum* from *C. baileyi* and *C. serpentis* (*Drall* digestion).
- c. Combine 15 µl of the *Cryptosporidium* nested PCR amplicon with one unit *Vspl*, 2.0 µl of 10x enzyme buffer, and 0.2 µl BSA solution (include with enzyme). Adjust final volume to 20 µl with sterile dH₂O.
- d. For digestion with *Drall*, combine 15 µl of the *Cryptosporidium* nested PCR amplicon with one unit *Drall*, 2.0 µl of 10x enzyme buffer, and adjust final volume to 20 µl.
- e. Positive controls for *C. parvum* and other species of *Cryptosporidium* must be digested in the same manner and alongside test samples.
- f. Incubate digestion samples for at least 2 hr at 37°C.
- g. Analyze samples by gel electrophoresis using a 3% NuSieve gel prepared with 0.5x TAE and 0.2% ethidium bromide.
- h. Mix 10-15 µl of nested amplification product with 2-3 µl of gel loading solution and load into wells of gel. Include at least one lane containing 25 bp DNA ladder to estimate the size of restriction fragments present.
- i. Run the gel at 125 volts (constant voltage) for at least 45 min.
- j. RFLP patterns can be visualized on the agarose gel by using a UV transilluminator. Photograph the gel to have a permanent record of the results using a Polaroid Type 667 film (or a digital system).
- k. Predicted banding patterns confirmatory for the presence of *Cryptosporidium* spp. are listed in *Table 13*.

Table 13: Evaluation of Nested PCR Amplification and RFLP Analysis to Differentiate *Cryptosporidium* spp.

Organism	PCR Amplicon		RFLP Digestion Products (bp)	
	Primary (bp)	Nested (bp)	<i>Vspl</i>	<i>Drall</i>
<i>C. hominis</i> (formerly <i>C. parvum</i> , genotype I)	844	593	503 and 90	-
<i>C. parvum</i> (formerly <i>C. parvum</i> , genotype II)	840	590	-	-
<i>C. baileyi</i>	831	579	-	295 and 284 [†]
<i>C. serpentis</i>	836	583	-	298 and 284 [†]
<i>C. muiris</i>	-	-	-	-
<i>C. wrairii</i>	-	-	-	-

[†]Indistinguishable within an agarose gel

References

- Centers for Disease Control. 1997. Outbreak of cyclosporiasis--Northern Virginia-Washington D.C.-Baltimore, Maryland, Metropolitan area, 1997. *Morbid. Mortal. Weekly Rep.* **46**:689-691.
- Centers for Disease Control. 1998. Update: outbreak of cyclosporiasis--Ontario Canada, May 1998. *Morbid. Mortal. Weekly Rep.* **47**:806-809.
- Eberhard, M.L., A.J. da Silva, B.G Lilley, and N.J. Pieniazek. 1999. Morphological and molecular characterization of new *Cyclospora* species from Ethiopian monkeys: *C. cercopithecii* sp.n., *C. colobi* sp.n., and *C. papionis* sp.n. *Emerg. Infect. Dis.* **5**:651-658.
- Herwaldt, B.L. 2000. *Cyclospora cayetanensis*: a review, focusing on the outbreaks of cyclosporiasis in the 1990's. *Clin. Inf. Dis.* **31**:1040-1057.
- Herwaldt, B.L., M.-L. Ackers, and the *Cyclospora* Working Group. 1997. An outbreak in cyclosporiasis associated with imported raspberries. *N. Engl. J. Med.* **336**:1548-1556.
- Herwaldt, B.L., M.J. Beach and the *Cyclospora* Working Group. 1997. The return of *Cyclospora* in 1997: another outbreak of cyclosporiasis in North America associated with imported raspberries. *Ann. Intern. Med.* **130**:210-220.
- Ho, A.Y., A.S. Lopez, M.G. Eberhart, R. Levenson, B.S. Finkel, A.J. da Silva, *et al.* 2002. Outbreak of cyclosporiasis associated with imported raspberries, Philadelphia, Pennsylvania, 2000. *Emerg. Inf. Dis.* **8**:783-788.
- Jinneman, K.C., J.H. Wetherington, W.E. Hill, A.M. Adams, J.M. Johnson, B.J. Tenge, N-L. Dang, R.L. Manger, and M.M. Wekell. 1998. Template preparation for PCR and RFLP of amplification products for the detection and identification of *Cyclospora* sp. and *Eimeria* spp. oocysts directly from raspberries. *J. Food Prot.*, **61**:1497-1503.
- Jinneman, K.C., J.H. Wetherington, W.E. Hill, C.J. Omiescinski, A.M. Adams, J.M. Johnson, B.J. Tenge, N-L. Dang, R.L. Manger, and M.M. Wekell. 1999. An oligonucleotide-ligation assay for the differentiation between *Cyclospora* and *Eimeria* spp. polymerase chain reaction amplification products. *J. Food Prot.*, **62**:682-685.
- Koumans, E.H.A. D.J Katz, J.M. Malecki, S. Kumar, S.P. Wahlquist, M.J. Arrowood, *et al.* 1998. An outbreak of cyclosporiasis in Florida in 1995: a harbinger of multistate outbreaks in 1996 and 1997. *Am J. Trop. Med. Hyg.* **59**:235-242.
- Lopez, A.S., D.R. Dodson, M.J. Arrowood, P.A. Orlandi, A.J. da Silva, J.W. Bier, *et al.* 2001. Outbreak of cyclosporiasis associated with basil in Missouri in 1999. *Clin. Infec. Dis.* **32** :1010-1017.
- Lopez, F.A., J. Manglicmot, T.M. Schimidt, C.Yeh, H.V. Smith, and D.A. Relman. 1999. Molecular characterization of *Cyclospora*-like organisms from baboons. *J. Infect. Dis.* **179**:670-676.

- Orlandi, P.A. and K.A. Lampel. 2000. Extraction-free, filter-based template preparation for the rapid and sensitive PCR detection of pathogenic parasitic protozoa. *J. Clin. Microbiol.* **38**:2271-2277.
- Orlandi, P.A., D.-M. T. Chu, J.W. Bier, and G.J. Jackson. 2002. Parasites and the food supply. *Food Technology* **56**:72-81.
- Orlandi, P.A., L. Carter, A.M. Brinker, A.J. da Silva, K.A. Lampel, and S.R. Monday. 2003. Targeting single nucleotide polymorphisms in the 18S rRNA gene to differentiate *Cyclospora* spp. and *Eimeria* spp. by multiplex PCR. *Appl. Environ Microbiol.* Submitted for publication.
- Ortega, Y., C.R. Roxas, R.H. Gilman, N.J. Miller, L. Cabrera, C. Taquiri, and C.R. Sterling. 1997. Isolation of *Cryptosporidium parvum* and *Cyclospora cayetanensis* from vegetables collected in markets of an endemic region in Peru. *Am. J. Trop Med. Hyg.* **57**:683-686.
- Ortega, Y., C.R. Sterling, R.H. Gilman, V.A. Cama, and F. Diaz. 1993. *Cyclospora* species-a new protozoan pathogen of humans. *N. Engl. J. Med.* **328**:1308-1312.
- Ortega, Y., R.H. Gilman, and C.R. Sterling. 1994. A new coccidian parasite (Apicomplexa:Eimeriidae) from humans. *J. Parasitol.* **80**:625-629.
- Pieniazek, N.J. and B.L. Herwaldt. 1997. Reevaluating the molecular taxonomy: is human-associated *Cyclospora* a mammalian *Eimeria* species? *Emerg. Inf. Dis.* **3**:381-383.
- Relman, D.A., T.M. Schmidt, A. Gajadhar, M. Sogin, J. Cross, K. Yoder, O. Sethabutr, and P. Echeverria. 1996. Molecular phylogenetic analysis of *Cyclospora*, the human intestinal pathogen, suggests that it is closely related to *Eimeria* species. *J. Infect. Dis.* **173**:440-445.
- Sherchand, J.B., J.H. Cross, M. Jimba, S. Sherchand, M.P. Shrestha. 1999. Study of *Cyclospora cayetanensis* in health care facilities, sewage water and green leafy vegetables in Nepal. *Southeast Asian J. Trop. Med. Public Health* **30**:58-63.
- Slifko, T.R., H.V. Smith, and J.B. Rose. 2000. Emerging parasite zoonoses associated with water and food. *Int. J. Parasitol.* **30**:1379-1393.
- Sturbaum, G.D., Y.R. Ortega, R.H. Gilman, C.R. Sterling, L. Cabrera, and D. Klein. 1998. Detection of *Cyclospora cayetanensis* in wastewater. *Appl. Environ. Microbiol.* **64**:2284-2286.
- Sturbaum, G.D., C. Reed, P.J. Hoover, B.H. Jost, M.M. Marshall, and C. R. Sterling. 2001. Species-specific, nested PCR-restriction fragment length polymorphism detection of single *Cryptosporidium parvum* oocysts. *Appl. Environ. Microbiol.* **67**:2665-2668.
- Wurtz, R. 1994. *Cyclospora*: a newly identified intestinal pathogen of humans. *Clin. Inf. Dis.* **18**:620-623.

- Yoder, K.E., O. Sethabutr, and D.A. Relman. 1996. PCR-based detection of the intestinal pathogen *Cyclospora*, pp. 169-176. In D.H. Persing (ed.), PCR protocols for emerging infectious diseases, a supplement to diagnostic molecular microbiology: principles and applications. ASM Press, Washington, DC.