

OCT 31 2000

## 510(K) SUMMARY OF SAFETY AND EFFECTIVENESS

FOR

**BOSTON RGP Lenses Wet Shipped In Boston Advance Comfort Formula  
Conditioning Solution And Stored For Up To 30 Days****1. SUBMITTER INFORMATION:**

Polymer Technology  
Global Vision Care  
1400 N. Goodman Street  
Rochester, New York 14603-0450

**2. CONTACT PERSON:**

	Debra Ketchum
	Manager, Regulatory Affairs
Address:	1400 North Goodman Street
	P.O. Box 30450
	Rochester, New York 14603-0450
Telephone No.:	(716) 338-8638
Fax No.:	(716) 338-0702
E-mail Address:	dketchum@bausch.com

**3. DEVICE IDENTIFICATION:**

Classification Name:	Rigid Gas Permeable (hydrophobic) Contact Lens Material
Proprietary Name:	BOSTON Contact Lens Materials
Common Name:	fluoro silicone acrylate rigid gas permeable contact lens material

**4. PREDICATE DEVICE:**

Paragon HDS (paflucocon B) and Menicon Z (tisilfocon B) Rigid Gas Permeable Contact Lenses wet shipped and up to 30 days storage have been selected as the Predicate Devices.

**5. DESCRIPTION OF THE DEVICE:**

Boston RGP Contact Lenses are composed of silicone acrylate or fluoro silicone acrylate copolymers wet shipped in Boston Advance Comfort Formula Conditioning Solution and stored for up to 30 days.

**6. INDICATIONS FOR USE:**

BOSTON XO (hexafocon A), BOSTON EO (enfluocon B), BOSTON ES (enfluocon A), BOSTON 7 (satafocon A) and BOSTON RXD (itabisfluorofoccon A) RGP Contact Lenses are indicated for daily wear for the correction of refractive ametropia (myopia, hyperopia, astigmatism and presbyopia) in aphakic and not-aphakic persons with nondiseased eyes. The lens may be disinfected using a chemical disinfecting system only.

BOSTON IV (itafocon B), BOSTON II (itafocon A) RGP Contact Lenses are indicated for daily wear for the correction of visual acuity for non-aphakic persons with myopia, hyperopia, or keratoconus and for the correction of corneal astigmatism up to 4.00 diopters. The lens is disinfected using a chemical (not heat) disinfection system recommended in the labeling.

**7. DESCRIPTION OF SAFETY AND SUBSTANTIAL EQUIVALENCE:**

The applicant performed stability, compatibility, and microbiology testing on *Boston RGP Contact Lenses wet shipped in Boston Advance Comfort Formula Conditioning Solution and stored for up to 30 days*. This testing in conjunction with toxicology testing on the shipping case plastics supports the claim of substantial equivalence to Menicon Z and Paragon HDS lenses wet shipped and stored for up to 30 days.

**Toxicology on the Shipping Lens Case:**

**Acute Ocular Irritation:**

Acute Ocular Irritation test was performed on the shipping lens case and produced no ocular irritation.

**Systemic Injection**

The material from the shipping lens case meets the requirements of the Systemic Injection Test and is considered non-toxic.

**Stability/Compatibility:**

Silicone acrylate (SA) and fluorosilicone acrylate (FSA) tinted rigid gas permeable contact lenses were subjected to a thirty-day soak in Boston Advance Comfort Formula Conditioning Solution according to the lens compatibility protocol. The average changes for each parameter (diameter, base curve and power), relative to the initial measurements were determined and compared to the DRAFT ISO/CD 8321-2: Optics

and Optical Instruments – Contact Lenses- Part 1: Specifications for rigid corneal and scleral contact lenses.

After soaking in the contact lens carrying cases at room temperature for thirty days the silicone acrylate and fluoro silicone acrylate rigid gas permeable contact lenses were determined to be physically compatible with Boston Advance Comfort Formula Conditioning Solution.

#### Microbiology

A bioburden study was completed. A set of test lenses was cleaned with Boston Laboratory Lens Cleaner and subjected to bioburden testing. In this test, two sets of lenses were tested to validate the storage in Boston Advance Comfort Formula after 30 days. One set was stored dry (control) and the other set was stored in Boston Advance Comfort Formula. This established the "cleanliness" of the test samples prior to entering the stability study. Testing showed that the colony forming units (CFU) per lens was less than 10. The acceptance criteria is 100 CFU per lens.

#### **8. SUBSTANTIAL EQUIVALENCE**

*Boston XO RGP Contact Lenses wet shipped in Boston Advance Comfort Formula Conditioning Solution and stored for up to 30 days is substantially equivalent to the currently marketed Paragon HDS (paflucocon B), approved in P870024, (S042), and Menicon Z (tisilfocon B) Rigid Gas Permeable Contact Lens 510(k) Premarket Notification No. K972443.*



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

OCT 31 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Debra Ketchum  
Manager, Regulatory Affairs  
Polymer Technology  
1400 N. Goodman St.  
P.O. Box 30450  
Rochester, NY 14603-0450

Re: K002025

Trade Name: Boston XO (hexafocon A), Boston EO (enfluocon B), Boston ES  
(enfluocon A), Boston 7 (satafocon A), Boston RXD  
(itabisfluorofoccon A), Boston IV (itafocon B), Boston II  
(itafocon A) Boston RGP Contact Lenses (Wet Shipped and up to 30 day  
storage in Boston Advance Comfort Formula Conditioning Solution)

Regulatory Class: II  
Product Code: 86 HQD  
Dated: October 4, 2000  
Received: October 5, 2000

Dear Ms. Ketchum:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Ms. Debra Ketchum

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-6413. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in cursive script that reads "A. Ralph Rosenthal".

A. Ralph Rosenthal, M.D.  
Director  
Division of Ophthalmic and Ear,  
Nose and Throat Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Polymer Technology  
1400 North Goodman Street  
P.O. Box 30450  
Rochester, NY 14603-0450

**Indications for Use Statement**

510(k) Number (if known): K002025

Device Name: BOSTON XO (hexafocon A) RGP Contact Lens  
BOSTON EO (enfluocon B) RGP Contact Lens  
BOSTON ES (enfluocon A) RGP Contact Lens  
BOSTON 7 (satafocon A) RGP Contact Lens  
BOSTON RXD (itabisfluorofoccon A) RGP Contact Lens  
BOSTON IV (itafocon B) RGP Contact Lens  
BOSTON II (itafocon A) RGP Contact Lens

*Indications for Use:*

BOSTON XO (hexafocon A), BOSTON EO (enfluocon B), BOSTON ES (enfluocon A), BOSTON 7 (satafocon A) and BOSTON RXD (itabisfluorofoccon A) RGP Contact Lenses are indicated for daily wear for the correction of refractive ametropia (myopia, hyperopia, astigmatism and presbyopia) in aphakic and not-aphakic persons with nondiseased eyes. The lens may be disinfected using a chemical disinfecting system only.

BOSTON IV (itafocon B), BOSTON II (itafocon A) RGP Contact Lenses are indicated for daily wear for the correction of visual acuity for non-aphakic persons with myopia, hyperopia, or keratoconus and for the correction of corneal astigmatism up to 4.00 diopters. The lens is disinfected using a chemical (not heat) disinfection system recommended in the labeling.

(PLEASE DO NOT WRITE BELOW THIS LINE-CONTINUE ON ANOTHER PAGE IF NEEDED)  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use

X

OR

Over-The-Counter-Use \_\_\_\_\_

Myra Smith  
(Division Sign-Off)

Division of Ophthalmic Devices

510(k) Number K002025