

OCT 11 2000

## 510(k) Summary

*Page 1 of 3*  
*K002345*

### Surgi-Vision Urethral Coil

Common/Classification Name: Accessory to  
Magnetic Resonance Diagnostic Device, 21 CFR 892.1000

Surgi-Vision, Inc.  
9250 Rumsey Road, Suite 100  
Columbia, MD 21045

Contact: Nancy E. Taylor, Prepared: July 28, 2000

#### A. LEGALLY MARKETED PREDICATE DEVICES

The **Surgi-Vision Urethral Coil** is substantially equivalent to the Medrad MRInnervu Endorectal MRI Coils, which were cleared for marketing on June 26, 1995, in premarket notifications K952232 and K952235.

#### B. DEVICE DESCRIPTION

The **Surgi-Vision Urethral Coil** is a specialty coil for use in MRI imaging of the anatomical regions surrounding the urethra. The signals picked up by the coil are conducted through a small coaxial cable to a connection with the standard surface coil connector for GE MRI systems. The coil and cable are completely sealed inside the insulating layer.

#### C. INTENDED USE

The **Surgi-Vision Urethral Coil** is recommended for high-resolution Magnet Resonance Imaging of the female urethra. The single use, disposable Urethral Coil was designed for insertion into the urethra of the patient during MRI scans in order to obtain improved image quality in the anatomical regions surrounding the urethra. The unique coil shape facilitates orientation of the coil to the anatomy. This product is to be used with a 1.5T MRI machine.

#### **D. SUBSTANTIAL EQUIVALENCE SUMMARY**

K002345

The **Surgi-Vision Urethral Coil** has similar but not identical indications for use as the legally marketed predicate device. However, the intended use, providing intercavitary reception of MRI signals for the purpose of creating an MRI image of nearby anatomy, is the same.

The **Surgi-Vision Urethral Coil** has the same technological characteristics as the predicate device. Both proposed and predicate devices have an electronic matching circuit, a connecting coaxial cable, and an intercavitary probe with a radiofrequency receiving coil. However, there are differences in construction and design that make it necessary to provide performance data to assure substantial equivalence. Such performance data are available and do demonstrate substantial equivalence.

#### **E. TECHNOLOGICAL CHARACTERISTICS**

See Section D, above.

#### **F. TESTING**

Surgi-Vision carried out testing and/or analysis of the **Surgi-Vision Urethral Coil** that addressed the following issues:

- (1) Possibility of excess RF heating;
- (2) Possibility of increased susceptibility of patients to peripheral nerve stimulation;
- (3) Imaging performance; and
- (4) Mechanical testing.

The results of the heating experiments demonstrate that there is no excess heating when SVUC is positioned in a phantom that is representative of clinical conditions. The change in temperature observed during use of the SVUC is not significantly different than that observed without the coil. The calculations done to determine current leakage by the MRI pulsed gradient field demonstrate that there is no possibility of increased susceptibility of patients to nerve stimulation. Finally, the imaging performance demonstrates enhanced resolution of

K002345

urethral regions. The results of the testing demonstrate that there are no safety problems for imaging of a patient using the **Surgi-Vision Urethral Coil** if the instructions for use are followed.

The imaging performance was evaluated in a porcine model as shown in Exhibit VI. The porcine model was selected for its body mass and similarity of cutaneous and subcutaneous tissue structure. Comparison images were done using a surface coil. The SVUC images were superior to the surface coil images taken with the same imaging parameters.

Mechanical test demonstrated tensile results that exceeded the one pound minimum. These results are described in Exhibit XI.

#### **G. CONCLUSIONS**

This pre-market submission has demonstrated Substantial Equivalence as defined and understood in the Federal Food Drug and Cosmetic Act and various guidance documents issued by the Center for Devices and Radiological Health.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

OCT 11 2000

Nancy E. Taylor  
CEO/President  
Surgi-Vision, Inc.  
9250 Rumsey Rd. Suite 100  
Columbia, MD 21045

Re: K002345  
Surgi-Vision Urethral Coil  
Dated: August 1, 2000  
Received: August 2, 2000  
Regulatory class: II  
21 CFR 892.1000/Procode: 90 MOS

Dear Ms. Taylor:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment ~~date of the Medical Device Amendments~~, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination ~~assumes compliance with the Current Good Manufacturing Practice requirements~~, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4639. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

Daniel G. Schultz, M.D.  
Captain, USPHS  
Acting Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure (s)

## STATEMENT OF INDICATIONS FOR USE

510(K) Number (if known): K002345

Device name: Surgi-Vision Urethral Coil

### Indications for Use:

The Surgi-Vision Urethral Coil is recommended for high-resolution Magnetic Resonance Imaging of the female urethra. The single use, disposable coil has been designed to be inserted in the urethra of the patient during MRI scans in order to obtain improved image quality in the anatomical regions surrounding the urethra. The unique coil shape facilitates orientation of the coil to the anatomy. This product is to be used with a 1.5T MRI machine.

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

\_\_\_\_\_  
Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use Δ  
(Per 21 CFR 801.109)

OR

Over-the-Counter Use \_\_\_\_\_

*David A. Legman*  
(Division Sign-Off)  
Division of Reproductive, Abdominal, ENT,  
and Radiological Devices  
510(k) Number K002345