

MAR 10 2000

K993805

510(K) SUBMISSION FOR NITRI-CARE NITRILE POWDER-FREE STERILE MEDICAL EXAMINATION GLOVE
SUBMISSION DATE: 1999-11-08

SUMMARY OF SAFETY AND EFFECTIVENESS

A. INFORMATION

1. SUBMITTER'S

Name: BEST MANUFACTURING COMPANY

Address: 579 Edison Street
Menlo, GA 30731 USA

Telephone Number: 706 862 2302

Contact Person: David C. Young

Date Summary Prepared: 1999-11-08

2. NAME OF DEVICE

Trade or Proprietary Name: NITRI-CARE Nitrile Powder-Free Sterile
Medical Examination Glove

Common or Usual Name Sterile Nitrile Powder-Free Patient
Examination Glove

Classification Name: Patient Examination Glove

3. PREDICATE DEVICE

N-DEX Nitrile Powder-Free Medical
Examination Glove, K992170

IDENTIFICATION NAME,
NUMBER

4. DESCRIPTION OF DEVICE

a. How the device functions:

Nitrile rubber films form an excellent barrier to body fluids and
bloodborne pathogens.

b. Scientific concepts that form the basis for the device:

The nitrile rubber is water tight under normal conditions of use. It's
tensile properties cause it to conform to the hand, allowing fine movement
necessary for treatment. The absence of natural rubber latex in the
product yields no latex protein allergens.

c. Physical and performance characteristics such as design, materials, and physical properties:

Nitrile rubber is known to create a superior barrier to bloodborne
pathogens and body fluids.

5. STATEMENT OF INTENDED USE, INCLUDING DESCRIPTION OF THE DISEASE OR CONDITIONS THAT THE DEVICE WILL ADDRESS

This is a disposable device, intended for medical purposes, that is worn on the examiner's hand to prevent contamination between the patient and examiner. Powder-free examination gloves are suitable in situations where powder is not desirable. Sterile gloves are suitable where a sterile examination glove is required.

6. EXPLANATION OF SIMILARITIES OR DIFFERENCES TO PREDICATE DEVICE

- The proposed device is identical to the predicate examination glove K992170 device, except for the following:
The proposed device is labelled "Sterile".

B. IF SE DECISION BASED ON PERFORMANCE DATA:

1. DISCUSSION OF NON-CLINICAL TESTS

<u>Specification</u>	<u>Proposed</u> NITRI-CARE Nitrile Powder-Free Sterile Medical Examination Glove	<u>Predicate</u> N-DEX Nitrile Powder-Free Medical Examination Glove
Performance Standards	ASTM	ASTM
Watertightness	ASTM	ASTM

2. DISCUSSION OF CLINICAL TESTS

<u>Specification</u> <u>Safety</u>	<u>Proposed</u>	<u>Predicate</u>
Rabbit Irritation	Passes	Passes
Guinea Pig Sensitization	Passes	Passes
Modified Draze Test (Human Study)	Passes	Passes

DESCRIPTION OF SUBJECTS

For the Modified Draze Test, 200 human subjects were used. The criteria for inclusion in the study was as discussed in the study, pages 6 and 7, paragraphs 3.11 and 3.12 (see Section N: Human and Animal Testing).


DISCUSSION OF SAFETY OR EFFECTIVENESS DATA OBTAINED
(with specific reference to adverse effects and complications)

See Section N: Human and Animal Testing, page 4 "SUMMARY", of the
Modified Draze Test.

3. CONCLUSIONS DRAWN FROM NONCLINICAL AND CLINICAL TESTS
THAT DEMONSTRATE SAFETY AND EFFECTIVENESS, AND
Performance => PREDICATE PRODUCT

The NITRI-CARE Nitrile Powder-Free Sterile Medical Examination
Glove has been carefully compared to a legally marketed device in the
510(k). The data summaries indicate that the proposed product meets or
exceeds accepted scores for the predicate product in both nonclinical tests
and satisfies the requirements for a safe and effective "powder-free"
medical glove.

Pursuant to 21 C.F.R. 807.87 (j), I David C. Young, Director, Regulatory
Affairs and Quality Assurance, certify that to the best of my knowledge
and belief and based upon the data and information submitted to me in the
course of my responsibilities as the Director, Regulatory Affairs and
Quality Assurance, for the Best Manufacturing Company, and in reliance
thereupon, the data and information submitted in this premarket
notification are truthful and accurate and that no facts material to a review
of the substantial equivalence of this device have been knowingly omitted
from this submission.



David C. Young, Director, Regulatory Affairs & Quality Assurance

1999-11-08

DATE



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR 10 2000

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. David C. Young
Director, Regulatory Affairs
& Quality Assurance
Best Manufacturing Company
579 Edison Street
Menlo, Georgia 30731-0008

Re: K993805
Trade Name: NITRI-CARE Nitrile Powder-Free Sterile
Medical Examination Glove, Blue, Cherry-
Flavored
Regulatory Class: I
Product Code: LZA
Dated: February 15, 2000
Received: February 22, 2000

Dear Mr. Young:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any


Page 2 -Mr. Young

the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4690. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

K993805

510(K) SUBMISSION FOR NITRI-CARE NITRILE POWDER-FREE STERILE MEDICAL EXAMINATION GLOVE
SUBMISSION DATE: 1999-11-08

INDICATIONS FOR USE

Applicant: **Best Manufacturing Company**

510(k) Number (if known) *

Device Name: **NITRI-CARE Nitrile Powder-Free Sterile Medical Examination**
Glove, Blue, Cherry-flavored

The **NITRI-CARE Nitrile Powder-Free Sterile Medical Examination Glove** is a disposable device intended for medical purposes; that is worn on the examiner's hand, to prevent contamination between the patient and examiner (21 CFR 880.6250).

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH Office of Device Evaluation (ODE)


Prescription Use _____
Per 21 CFR 801.109

OR

Over-The-Counter X

* For a new submission, do NOT fill in the 510(k) number blank.

PAGE 3 OF 11



(Division Sign-Off)

Division of Dental, Infection Control,
and General Hospital Devices

510(k) Number K993805