

## Summary of Safety and Effectiveness for the Irrigating Cannula

*Submitted by*  
CA-Guard, LTD.  
610 E. Olympia Ave  
Punta Gorda, FL 33950  
Phone: (941) 833-1102

**Contact Person:** Al Weisenborn  
**Device Trade Name:** Irrigating Cannula  
**Common Name:** Cannula  
**Classification Name:** Endoscopes and Accessories, per 21 CFR §876.1500

### Identification of a Legally Marketed Predicate Device

The Irrigating Cannula is substantially equivalent to the cannula manufactured and marketed by the United States Surgical Corporation as part of the Versaport® Trocar System pursuant to 510(k) K954108.

### General Description

The Irrigating Cannula is a non-toxic, sterile, single use, disposable surgical access port. The device is intended to be inserted percutaneously into the peritoneal cavity to provide access for endoscopic instruments. The cannula is intended to replace the cannula that is supplied with the United States Surgical Versaport® 11mm Trocar System.

The proximal end of the cannula is equipped with a female luer connector. The luer communicates with 6 longitudinal grooves on the surface of the inner member of the cannula. The inner member is covered with polyolefin shrink tubing. An array of twenty-four (24) holes deliver infusate to the compromised tissue surfaces of the trocar wound.

### Intended Use

The CA-Guard Irrigating Cannula is intended to provide an access port to body cavities for endoscopes and endoscopic accessories when intraoperative irrigation or infusion of trocar wound surfaces is desired. The cannula has been tested for compatibility with lidocaine and marcaine. Saline and water may also be used as irrigants. The CA-Guard Irrigating Cannula is inherently a needleless system and therefore, minimizes the possibility of "needle sticks" for healthcare workers, in accordance with the Needlestick Safety and Prevention Act of 2000.

The CA Guard Irrigating Cannula is designed to be mechanically compatible with all components of the United States Surgical Versaport® 11mm Trocar System using the 5 mm PLUS SEAL.

### **Summary of Technological Characteristics**

The Irrigating Cannula was compared to the predicate device using 14 points of comparison and found to be equivalent.

### **Summary of Performance Data**

Performance of the device was characterized and compared to that of the predicate utilizing 5 tests. Additionally, the cannula was tested for compatibility with the indicated infusates. The tests demonstrated the Irrigating Cannula is substantially equivalent to the cannula manufactured and marketed by the United States Surgical Corporation as part of the Versaport® Trocar System pursuant to 510(k) K954108. The materials of construction have been carefully selected for their long history of biocompatibility.

Since the Irrigating Cannula embodies technological characteristics essentially identical to those of the predicate device, we believe the device is safe and effective and that it performs as well as or better than the predicate device. The device has been designed and developed utilizing design control methods in compliance with the QSR. The Irrigating Cannula will be manufactured per specifications and good manufacturing practices that ensure the device is safe and effective for its intended use.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

MAY 24 2002

CA-GUARD, LTD.  
c/o Mr. Al Weisenborn  
Official Correspondent  
19526 East Lake Drive  
MIAMI FL 33015

Re: K020614  
Trade/Device Name: Irrigating Cannula  
Regulation Number: 21 CFR 876.1500  
Regulation Name: Endoscope and accessories  
Regulatory Class: II  
Product Code: 78 KOG  
Dated: February 21, 2002  
Received: February 25, 2002

Dear Mr. Weisenborn:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure

**Indications for Use**Page 1 of 1510(k) Number (if known): K020614Device Name: Irrigating Cannula

Indications for Use:

The CA-Guard Irrigating Cannula is intended to provide an access port to body cavities for endoscopes and endoscopic accessories when intraoperative irrigation or infusion of trocar wound surfaces is desired.

*David A. Segura*  
(Division Sign-Off)  
Division of Reproductive, ~~Abdominal~~,  
and Radiological Devices  
510(k) Number K020614

PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE IF NEEDED

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ☒  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use ☐

(Optional Format 1-2-96)