

UNITED STATES OF AMERICA
 FOOD AND DRUG ADMINISTRATION
 CENTER FOR DEVICES AND RADIOLOGICAL HEALTH
 MEDICAL DEVICES ADVISORY COMMITTEE

OBSTETRICS AND GYNECOLOGY DEVICES PANEL

68TH MEETING

THURSDAY,
 JUNE 3, 2004

The Panel met at 8:30 a.m. in the Whetstone Room of the Gaithersburg Holiday Inn, Two Montgomery Village Avenue, Gaithersburg, Maryland, Dr. Kenneth L. Noller, Chair, presiding.

PRESENT:

KENNETH L. NOLLER, M.D., Panel Chair
 SUSAN M. ASCHER, M.D., Temporary Voting Member
 ANDREW I. BRILL, M.D., Temporary Voting Member
 CAROL L. BROWN, M.D., Member
 LAWRENCE A. CRUM, Ph.D., Temporary Voting Member
 RALPH B. D'AGOSTINO, Ph.D., Temporary Voting Member
 MICHAEL P. DIAMOND, M.D., Non-Voting Member
 EVELYN R. HAYES, Ph.D., Member
 PAUL J.A. HILLARD, M.D., Member
 GRACE M. JANIK, M.D., Temporary Voting Member
 KLEIA R. LUCKNER, J.D., M.S.N., Consumer Representative
 HUGH MILLER, M.D., Member
 ANNE C. ROBERTS, M.D., Temporary Voting Member
 THADDEUS V. SAMULSKI, Ph.D., Temporary Voting Member
 STEPHEN B. SOLOMON, M.D., Non-Voting Member
 JONATHAN W. WEEKS, M.D., Member
 BRANDFORD J. WOOD, M.D., Temporary Voting Member
 JOYCE WHANG, Ph.D., Panel Executive Secretary
 MICHAEL T. BROWN

This transcript has not
 been edited and FDA
 makes no representation
 regarding its accuracy

FDA REPRESENTATIVES:

NANCY BROGDON, Director, Div. Of Reproductive,
Abdominal and Radiological Devices
COLIN POLLARD, Chief, Obstetrics and Gynecology
Devices Branch
JULIA A. CORRADO, M.D., Medical Officer
KATHRYN S. DAWS-KOPP
NOEL DEL MUNDO, M.D.
BRUCE A. HERMAN
LOREN A. ZAREMBA, Ph.D.

SPONSOR REPRESENTATIVES:

ROB NEWMAN, M.S., R.A.C.
BOBBIE S. GOSTOUT, M.D.
GINA K. HESLEY, M.D.
ELIZABETH A. STEWART, M.D.
CLARE M.C. TEMPANY, M.D.
KOBI VORTMAN, Ph.D.

I-N-D-E-X

| | |
|--|-----|
| Introductory Remarks | 4 |
| Open Public Hearing | 14 |
| Presentation by Sponsor | 14 |
| Presentation by FDA | 88 |
| Panel Discussion | 163 |
| Open Public Hearing | 264 |
| Panel Deliberations and Vote | 292 |
| Adjourn | 363 |

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

P R O C E E D I N G S

(8:25:23 a.m.)

DR. NOLLER: Everyone take their seats, please. We have a very full day so I want to get started exactly on time. My name is Ken Noller, and I'd like to call the meeting to order. This is the Meeting of Obstetrics and Gynecology Devices Panel. I request that everyone in attendance please sign in. If you have not done so, please go out and sign in at the front desk now.

I also note for the record that the voting members present constitute a quorum as required by 21 CFR Part 14. I'm going to ask the panel members to introduce themselves. Let's start at this end, please.

MS. MOONEY: Mary Lou Mooney. I'm the Vice President of Clinical Regulatory and Quality for SenoRx, and I'm the Industry Rep to the panel.

MS. LUCKNER: Kleia Luckner, Hospital Administrator, Toledo, Ohio, and I am the Consumer Rep.

DR. D'AGOSTINO: Ralph D'Agostino from

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Boston University, Biostatistician.

2 DR. BRILL: Andrew Bill. I am a Professor
3 OB-GYN, University of Illinois.

4 DR. HILLARD: Paula Hillard, Professor of
5 OB-GYN and Pediatrics, University of Cincinnati.

6 DR. DIAMOND: Michael Diamond, Professor
7 OB-GYN, Wayne State University, Detroit Michigan.

8 DR. ROBERTS: Anne Roberts, Professor of
9 Radiology, University of California - San Diego.

10 DR. NOLLER: I'm Ken Noller, Professor and
11 Chair of Tufts University OB-GYN.

12 DR. WHANG: I'm Joyce Whang. And I'm an
13 FDA Reviewer and the Executive Secretary for this
14 panel.

15 DR. BAILEY: I'm Mike Bailey. I'm also a
16 Reviewer in the OB-GYN Devices group, and I'm an
17 Assistant Executive Secretary.

18 DR. BROWN: Hi. Carol Brown, I'm a Panel
19 Member. I am an Assistant Professor at Cornell Weill
20 Medical College, OB-GYN and a GYN Oncologist at
21 Memorial Sloan-Kettering Cancer Center.

22 DR. CRUM: I'm Larry Crum from the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 University of Washington. I'm Director of the Center
2 for Industrial and Medical Ultrasound at the
3 University of Washington.

4 DR. JANIK: Grace Janik, Clinical
5 Professor at the Medical College of Wisconsin,
6 Reproductive Endocrinologist.

7 DR. SAMULSKI: Thad Samulski, Duke
8 University Medical Physics.

9 DR. HAYES: Evelyn Hayes, Professor of
10 Nursing, University of Delaware.

11 DR. ASCHER: Susan Ascher, Radiologist,
12 Georgetown University Hospital.

13 DR. WOOD: Bradford Wood, Interventional
14 Radiologist, National Institutes of Health.

15 MS. BROGDON: I'm Nancy Brogdon. I'm not
16 a member of the panel. I'm the Director of FDA's
17 Division of Reproductive, Abdominal, and Radiological
18 Devices.

19 DR. SOLOMON: Steve Solomon from
20 Department of Radiology, Johns Hopkins.

21 DR. NOLLER: Thank you. For the press,
22 the FDA press contact is Colin Pollard who is sitting

1 here in the front row. I don't expect that we'll have
2 any super controversial outbursts today, but we would
3 like everyone to please be courteous, turn off your
4 cell phones, and if you have anything to say, wait
5 until you're recognized and then come to the table.
6 For the audience and the panel members I will
7 recognize people before they speak. Our Executive
8 Secretary has some things to read into the Minutes.

9 DR. WHANG: There will be OB-GYN Devices
10 Panel on July 26th and 27th, so the remaining panel
11 meeting date for this year is October 25th to 26th.

12 We are pleased to introduce a new voting
13 member to this panel, Dr. Paula Hillard of the
14 Department of Obstetrics and Gynecology and the
15 Department of Pediatrics at the University of
16 Cincinnati, College of Medicine.

17 Today we will have eight temporary voting
18 members, Drs. Ascher, Brill, Crum, D'Agostino, Janik,
19 Roberts, Samulski and Wood. And I will now read into
20 the record the appointments to temporary voting status
21 signed by Daniel Schultz, M.D., the Acting Director
22 for the Center of Devices and Radiological Health.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 "Pursuant to the authority granted under
2 the Medical Devices Advisory Committee Charter dated
3 October 27th, 1990, and amended August 18th, 1999, I
4 appoint the following individuals as voting members of
5 the Obstetrics and Gynecology Devices Panel for this
6 meeting on June 3rd, 2004; Susan M. Ascher, M.D.,
7 Andrew I. Brill, M.D., Lawrence A. Crum, Ph.D., Ralph
8 B. D'Agostino, Ph.D., Grace M. Janik, M.D., Kenneth E.
9 Najarean, M.D., Anne C. Roberts, M.D., Thaddeus V.
10 Samulski, Ph.D., Bradford J. Wood, M.D.

11 For the record, these people are special
12 government employees and are consultants to this
13 panel. They have undergone the customary conflict of
14 interest review, and they have reviewed the material
15 to be considered at this meeting."

16 I will now read the conflict of interest
17 statement for this meeting. "The following
18 announcement addresses conflict of interest issues
19 associated with this meeting, and is made a part of
20 the record to preclude even the appearance of an
21 impropriety. To determine if any conflict existed,
22 the Agency reviewed the submitted agenda, and all

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 financial interests reported by the committee
2 participants. The Conflict of Interest statutes
3 prohibit special government employees from
4 participating in matters that could affect their or
5 their employer's financial interests. However, the
6 Agency has determined that participation of certain
7 members and consultants, the need for whose services
8 outweighs the potential conflict of interest involved
9 is in the best interest of the government. Therefore,
10 full waivers have been granted for Dr. Susan Ascher
11 and Anne Roberts, and limited waivers have been
12 granted for Drs. Michael Diamond and Steven Solomon
13 for their interest in firms that could potentially be
14 affected by the panel's recommendations.

15 Dr. Ascher's waiver involves a contract to
16 her employer funded for less than \$100,000 per year
17 with a competing firm. Dr. Roberts' waiver involves
18 a stockholding in a competing firm in which the value
19 is between \$15,001 and \$25,000. Dr. Diamond's
20 limited waiver involves a contract to his institution
21 for the sponsor study in which he had no involvement
22 in data generation or analysis, and for which total

1 funding to the institution was less than \$100,000.
2 Dr. Solomon's limited waiver involves a contract to
3 his institution for the sponsor study in which he had
4 no involvement in data generation or analysis, and for
5 which funding to the institution is unknown.

6 The waivers of Dr. Ascher and Dr. Roberts
7 allow them to participate fully in today's
8 deliberations. The limited waivers for Dr. Diamond
9 and Dr. Solomon allow them to participate in the panel
10 discussions, but exclude them from voting.

11 Copies of these waivers may be obtained
12 from the Agency's Freedom of Information Office, Room
13 12A-15 of the Parklawn Building. We would like to
14 note for the record that the Agency took into
15 consideration other matters regarding Drs. Diamond
16 and Solomon. They reported current interests with
17 firms at issue, but in matters that are not related to
18 today's agenda. The Agency has determined, therefore,
19 that these individuals may participate fully in the
20 panel's deliberations.

21 In the event that the discussions involve
22 any other products or firms not already on the agenda

1 for which an FDA participant has a financial interest,
2 the participant should excuse him or herself from such
3 involvement and the exclusion will be noted for the
4 record.

5 With respect to other participants, we ask
6 in the interest of fairness that all persons making
7 statements or presentations disclose any current or
8 previous financial involvement with any firm whose
9 products they may wish to comment on.

10 Transcripts for today's meeting are
11 available from Neal R. Gross and Company of
12 Washington, D.C., at (202) 234-4433, and videos are
13 available from FDA Live at (301) 984-0001, or FDA
14 Advisory Committee.com at (800) 627-8171.

15 Any presenters to the panel who have not
16 already done so should provide FDA with a hard copy of
17 your remarks, including overheads. Michelle Byrnes
18 will collect these from you at the podium.

19 DR. NOLLER: Thank you. First, I'd like
20 to introduce Colin Pollard, Chief of Obstetrics and
21 Gynecology Devices Branch of the Food and Drug
22 Administration.

1 DR. POLLARD: Thank you, Dr. Noller, and
2 I just have a few brief comments to kick off our panel
3 meeting today. I want to welcome all the panel
4 member, and thank you very much for coming from near
5 and far.

6 Today you'll be looking at a PMA for a
7 high-intensity focused ultrasound system, really a new
8 surgical modality that uses conventional MR Imaging
9 for pre-op treatment planning and MR thermal mapping,
10 really a new feature of MR technology for interactive
11 treatment feedback. And treatment of uterine fibroids
12 is the very first indication that's coming before this
13 center in a PMA. The technology, obviously, looks
14 capable of many other clinical applications and the
15 center is currently working on a plan to optimize our
16 regulatory review approach.

17 As you'll hear later in our presentation,
18 we put together something of a designer review team
19 for this PMA drawing from all parts of our center,
20 especially from the technical side, and as we look
21 around the table here I see several familiar faces,
22 but lot of new faces. And really, we've put together

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 something of a designer panel, as well, and so we're
2 very much looking for your input.

3 This PMA, the center granted expedited
4 review to based on unique features and advantages.
5 The FDA review is still ongoing, but we consider it
6 quite appropriate at this stage to hear the panel
7 input even as we continue to work our way through many
8 review issues. And finally, that we feel this
9 technology pushes the traditional envelope of clinical
10 management, and if the panel gets to that point, we'll
11 be definitely looking for input regarding training and
12 labeling, and credentialing and that sort of thing.

13 So with those initial comments, Dr.
14 Noller, I turn it back to you. Thank you.

15 DR. NOLLER: Thank you. Let me just ask
16 that Drs. Miller and Weeks introduce themselves,
17 please.

18 DR. MILLER: Hugh Miller from Arizona.

19 DR. NOLLER: And what do you do, Dr.
20 Miller?

21 DR. MILLER: I'm a Maternal-Fetal Medicine
22 Specialist.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 MR. WEEKS: Jonathan Weeks from
2 Louisville, Kentucky, Maternal-Fetal Medicine, Norton
3 Health Care.

4 DR. NOLLER: Thank you both. First, let
5 me ask before we open the public hearing, is there
6 anyone present who will be speaking, request speaking
7 at this session of the public hearing? All right. So
8 if there's no one at this time, I will not read the
9 conflict statement then. We'll move right ahead to
10 the presentation by the sponsor.

11 I'd like to introduce Rob Newman from
12 InSightec. The sponsor has been granted one hour and
13 15 minutes for their presentations. I ask the panel
14 members to hold all questions until the end of the
15 presentation.

16 MR. NEWMAN: Good morning, Chairman
17 Noller, and thank you very much, ladies and gentlemen
18 of the panel and the audience. I'm Rob Newman. I'm
19 from InSightec in Dallas, Texas. My trip here has
20 been paid for by my company. I am a member of the
21 sponsor.

22 I'd like to introduce other members of our

1 team today. Dr. Elizabeth Stewart is an Associate
2 Professor of Gynecology from Harvard Brigham & Women's
3 Hospital. Dr. Clare Tempany is a Professor of
4 Radiology at Brigham & Women's. They are Co-PIs of
5 the study. Dr. Stewart is the Lead PI for the study.

6 Also, Kobi Vortman, the President of
7 InSightec is here with us today. Karin Coyne, Senior
8 Research Scientist from MEDTAP International, who has
9 helped us with the quality of life work, and some of
10 the biostatistics. Kathy McDermott from MedTrials,
11 our CRO. We also have a guest here, Dr. Bobbie
12 Gostout, who is Assistant Professor from the Mayo
13 Clinic in Rochester, and Dr. Gina Hesley, an
14 Instructor of Radiology, who are Co-PIs from the Mayo
15 Clinical site.

16 This is an outline of our discussion
17 today. I'll give a brief introduction. Dr. Stewart
18 will talk about, in general, an overview of uterine
19 fibroids and their application here. I'll give a
20 brief overview of the device description. If you had
21 a chance to review the video in the package, I think
22 that will cover some of it, and there was quite a bit

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 of material in the panel package, so I won't go over
2 that in any great detail.

3 Dr. Clare Tempany will talk about a review
4 of MR anatomy, and what's commonly seen on MR that may
5 be a little bit more than what some of you see in a
6 regular clinical practice. She'll also discuss the
7 treatment development process. Dr. Stewart will talk
8 about clinical design trial results. I'll cover some
9 elements of training, in addition to what was in the
10 panel package, and then Dr. Stewart will summarize.

11 The indications for use for this device is
12 its for use in pre or peri-menopausal women with
13 symptomatic fibroids. The fibroids to be treated must
14 be visible on non-contrast MRI and should enhance on
15 contrast MR.

16 Outside the U.S., the system has received
17 CE Mark in Europe in 2002. Its commercially available
18 in Europe, Israel and Japan. In the U.S., the only
19 applications are investigational. We have treated
20 approximately 600 women worldwide for uterine
21 fibroids. And I'd like to introduce Dr. Stewart, who
22 will introduce the topic.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. STEWART: Mr. Chairman, panel members
2 and guests, my travel expenses were paid by InSightec
3 today. As Mr. Newman said, I serve as a Clinical
4 Trial Investigator for the company, and am a
5 Consultant for the company, but abide by the Harvard
6 Medical School ethical guidelines that limit
7 consulting when an investigator is involved in
8 clinical research.

9 I want to start today by talking about the
10 important problem of uterine fibroids. As everyone in
11 this room probably knows, this is a very important
12 clinical problem for women. That are very common
13 tumors and the prevalence rates vary from anywhere
14 from about 20 percent of women to being affected, to
15 more recent estimates looking at high-risk populations
16 by ultrasound where the prevalence of clinically
17 detectible fibroids appears to be in the range of 75
18 percent.

19 Most of the discussion regarding uterine
20 fibroids centers around cost and the costs are
21 substantial for a healthcare system. It's estimated
22 that the cost for hysterectomy alone in the U.S. along

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 per year is in excess of \$2 billion. This is really
2 the tip of the iceberg because it doesn't even start
3 to take into account other surgical options, non-
4 surgical options, medical options and various
5 alternative treatments that women seek to try to
6 control their symptoms.

7 I think it's important to note also that
8 there is information regarding productivity in women
9 with menorrhagia and so this is probably an under-
10 estimate for the kind of women that we're seeing in
11 our study who have clinically significant fibroids.
12 The estimation from 2000 was that lost productivity
13 due to menorrhagia or excessive menstrual flow is in
14 the range of \$1,600 per woman per year.

15 I think it's important to realize that
16 fibroids are a common source of morbidity for women.
17 They cause a lot of symptoms that tend to cluster in
18 several different areas. Menorrhagia or excessive
19 menstrual flow is an extremely important problem due
20 to fibroids. And for women, this really limits their
21 ability to carry out their work or their interactions
22 with their families. There are many women that spend

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 up to two weeks every month with significant menstrual
2 bleeding, and there are many women who have such
3 significant menstrual bleeding that they cannot attend
4 to any other activity for an hour or more without
5 having to stop to deal with changes in sanitary
6 protection.

7 Pain and discomfort are significant
8 symptoms related to uterine fibroids. Many clinically
9 significant fibroids are in the range of a three,
10 four, five month pregnant uterus, and this gives women
11 significant symptoms in terms of urinary frequency,
12 urgency bladder discomfort, pelvic discomfort.

13 These symptoms have been shown to
14 significantly impair health-related quality of life,
15 and in several studies there have been demonstrations
16 that women with uterine fibroids have significantly
17 lower health-related quality of life than population
18 norms. Uterine fibroids have also been linked to time
19 away from work and other activities that are important
20 to the economic system. And the Rand Corporation
21 estimated that medical therapy may fail to control the
22 symptoms in approximately two-thirds of women, so we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 do need better therapies for uterine fibroids.

2 There are treatment options for uterine
3 fibroids, but I think if you look at the range of
4 options available for uterine fibroids in contrast to
5 the woman who has a normal uterus and her options
6 available for endometrial ablation, the contrast is
7 clear. This panel has approved many devices for
8 endometrial ablation, and many of those are restricted
9 to women who have a structurally normal uterus.

10 Hysterectomy is a good solution for
11 uterine fibroids. It is very effective in solving the
12 symptoms, but it does have a significant morbidity
13 associated with it, and a significant time away from
14 work and family. For many women today to have the six
15 week recovery for a major surgery is something that
16 they cannot incorporate into their work and their
17 family.

18 Myomectomy is an option for women who have
19 a desire to retain their uterus but want resolution of
20 their fibroid symptoms. Clearly, there are some women
21 that are amenable to minimally invasive Myomectomies
22 if the fibroid is in the right position at the serosal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 or the mucosal surface. However, again this modality
2 is a surgical modality, and can have significant
3 recovery associated with it.

4 Uterine Artery Embolization has been an
5 important option that has been added in the past
6 decade for women with uterine fibroids. It has
7 significantly decreased recovery time and fewer
8 complications than hysterectomy. However, this
9 modality is associated with pain and fever post-
10 operatively, and there's increasing attention to the
11 fact that there is an age-related impairment of
12 ovarian function and this may be particularly an issue
13 for certain groups of women.

14 Thermally ablative therapies have been
15 tried previously for uterine fibroids. Many people
16 have had experience with either myolysis or
17 cryomyolysis, and there's a small experience with RF-
18 ablation. These techniques have not really made it
19 into the general gynecologist armamentarium, probably
20 because of a lack of thermal monitoring.

21 With these prior therapies, there was no
22 gauging of temperature, and so you couldn't tell had

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 you established a sufficient temperature to destroy
2 the tissue. If not, you probably decreased your
3 efficacy. Or if you exceeded the temperature goal,
4 potentially you were injuring normal tissue and
5 causing problems with adhesions or other follow-up.

6 Again, we do have drug therapies, but they
7 tend to fall into two broad categories. Drugs such as
8 oral contraceptives and progestins are widely used to
9 control fibroid symptoms, but they tend to not be
10 efficacious in the long-term.

11 On the other hand, GnRH agonists are very
12 effective, but their side effects are significant, and
13 their cost is significant, and so these drugs really
14 haven't been great long-term choices for women with
15 uterine fibroids.

16 We see that there's a spectrum of options
17 available for uterine fibroids that for women with
18 severe disease or who require a definitive solution,
19 hysterectomy is still a choice. But many women are
20 sitting down here with expected management and dealing
21 with significant levels of symptomatology because they
22 fear the surgical invasiveness of the other options,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 or because they cannot, again, take the time and the
2 recovery that's necessary to undergo a very invasive
3 option.

4 We think MRI guided focused ultrasound
5 surgery will be a very important option to offer
6 women. It will give them the symptom relief that they
7 require with significantly less invasiveness than many
8 of the other options.

9 There are several unique things that are
10 important to know about MRI guided focused ultrasound.
11 It is a non-invasive, rather than a minimally invasive
12 surgery. There is no surgical incision. There is no
13 probe that goes into the fibroid. It is able to be
14 accomplished as an out-patient procedure. Again, it
15 serves the uterus and is uterine sparing.

16 The other important issue is that it is a
17 fibroid-specific therapy. Unlike something like
18 uterine artery embolization that targets the entire
19 uterus, the fibroid is specifically targeted so that
20 there is no impact on the myometrium or the
21 endometrium.

22 Again, the real time feedback on

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 temperature gives you precise thermal ablation, and
2 this is very important both for the optimization of
3 safety and efficacy. Again, you can know that your
4 temperature is getting to a therapeutic level and
5 causing tissue destruction, and yet remaining in a
6 safe range. And we have found that this procedure
7 does not preclude or complicate future treatment
8 options.

9 I will return the presentation to Mr.
10 Newman, who will talk a little bit more about the
11 device.

12 MR. NEWMAN: Thank you, Dr. Stewart. I'll
13 just briefly review some of the key points of the
14 device itself. As I said, much of this material is in
15 the panel package, so I won't belabor the issues.

16 MR guided focused ultrasound is really a
17 combination of two things, the idea of focused
18 ultrasound as a source of thermal energy, and MR to
19 plan and control the treatment in progress. There's
20 two main components; one is the patient table and the
21 electronics that's attached to the MR system. In the
22 top of the patient table is the transducer, and

1 there's a water bath here. The patient lies on top of
2 that. The energy is transmitted through the abdominal
3 wall and focuses on a point inside the body.

4 Out next to the operator console of the MR
5 is the control console for the focused ultrasound.
6 Here's the regular MR console, and they sit side-by-
7 side so that you can see your work on both systems
8 during the treatment. Once the treatment begins, all
9 of the control of the treatment and all the
10 observation of the patient images is done from the
11 ExAblate workstation. Next slide, please.

12 Just a brief history of focused
13 ultrasound. Although this may be one of the first
14 times that many of you have heard about it, focused
15 ultrasound is a technology that's been around for a
16 long time. There are publications as early as the
17 1930s. We didn't invent this. We're just kind of the
18 latest people to carry on a long line of research in
19 this. The Fry Brothers in the 40s and 50s did a lot
20 of work on this looking at focused ultrasound in the
21 brain and other places in the body.

22 Lele carried on work looking at several

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 tumors. There's also some carried on with some work
2 using focused ultrasound for acoustic hemostasis and
3 other applications. In 1993, Hynynen, Cline and
4 others wrote the first paper on the combination of MR
5 with focused ultrasound using MR for the thermal
6 imaging. And then 1995 to present is ExAblate, the
7 development of the device we're discussing today.

8 The transducer, a little of the physics of
9 the transducer. The transducer lies here. The energy
10 passes through the skin and intervening tissue to
11 focus at a point. The energy is highly focused. It's
12 kind of like taking a magnifying glass and focusing
13 the sun's energy, so you can put your hand above the
14 magnifying glass, below the magnifying glass, and it
15 isn't until you get right at the point that the energy
16 is highly concentrated.

17 The density in the far field, the energy
18 is attenuated and absorbed along the beam path so
19 while it's highly concentrated here, it falls off with
20 distance in the far field.

21 The focused ultrasound energy propagates
22 through tissue and skin. It's blocked by air, such as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 in bowel or the rectum behind the focal point, and
2 it's absorbed by bone, such as the pubic bone or the
3 sacrum in the far field.

4 This is just brief picture showing the
5 patient lies on top of it, so here's the transducer
6 underneath. Here's the overall beam path for this
7 entire volume, and in the blue is the focal path for
8 a single sonication.

9 The one thing that's different about this
10 as a source of thermal energy is we ablate one small
11 piece at a time, as opposed to a cryoprobe where you
12 create a large two, or three, or four centimeter
13 lesion in one go, or RF ablation. We build up the
14 treatment from a series of these individual
15 sonications that are approximately 25 by 25 by 10
16 millimeters, so you're ablating about a half a cubic
17 centimeter at a time. A single sonication takes about
18 20 seconds, and the target is to raise the tissue in-
19 between 65 and 85 degrees Centigrade. If you raise
20 tissue above 57 degrees Centigrade for one second,
21 it's ablative.

22 There's a rapid fall-off. It's a highly

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 focused transducer, so there's a rapid fall-off with
2 distance, so just a very few millimeters away from the
3 focal spot you're back at normal body temperature.
4 There's also a combining effect here because tissue's
5 sensitivity to temperature is very time-dependent, so
6 when you look at one pixel, it's the time temperature,
7 it's the product of time and temperature that dictates
8 whether you've had ablation or reversal heating of
9 that point. So when do a treatment, the physician
10 draws a region of treatment around the area to be
11 treated, and then the system tiles it, if you will,
12 with a series of these jellybeans, so that basically
13 you draw a region of treatment and the system figures
14 out how many jellybeans are in the jar. So how many
15 will it take to completely cover this volume, so you
16 can do a single layer, you can do multiple layers, and
17 here's what it looks like in the horizontal plane, or
18 looking down on it from above.

19 The treatment is controlled by MR
20 thermometry. You're doing MR continuously throughout
21 the treatment, so it isn't like you do a single
22 planning image or a stereotactic plan before. You're

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 using MR continuous throughout the energy delivery, so
2 approximately every three seconds you're acquiring an
3 MR image. The accuracy in vivo in fibroid tissue is
4 about 3 degrees Centigrade. And what we're doing is
5 we're measuring change in temperature with MR. We
6 can't measure absolute temperature, but we measure
7 change relative to body core temperature.

8 We take the information. We use this to
9 tell us -- we can see the focal spot. We can tell
10 where the energy is being delivered in three
11 dimensions, and we can quantify the temperature to get
12 this time/temperature information from each pixel.

13 At the end of each sonication, using this
14 time/temperature information, we can calculate the
15 volume of tissue that exceeded the dose, and we can
16 use this information to plan the next sonication.

17 This is just a quick picture showing
18 during a single sonication we're acquiring an image
19 every three seconds here over a 15 second sonication
20 in this case, so you can see at 1.7 seconds, you can
21 see the spot start to show up on the MR image. At the
22 end we take this information, calculate the volume of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 tissue that was ablated, and we can draw one of these
2 time versus temperature histograms here or maps, so we
3 can where this little red cursor - it's hard to see in
4 this slide - but there's a little cursor here, and you
5 can see the time/temperature history out to 98
6 seconds.

7 If we would move that cursor somewhere
8 here in the background away from the focal point,
9 you'll just see some bouncing around, plus or minus a
10 few degrees of normal body temperature.

11 We've done extensive thermal modeling in
12 both 2D and 3D looking at a simulation of energy along
13 the beam path to quantify the absorption of the
14 energy, and to look at the tissue characteristics.
15 We've done this to explore boundary conditions, to
16 look at what type of dosimetry would be appropriate
17 for maximum effectiveness, and to minimize thermal
18 damage outside the treatment volume. And to really
19 simulate things that we can't really do in vivo to
20 look at kind of worst case scenarios of energy
21 delivery and absorption that we wouldn't be able to do
22 in vivo.

1 Pre-clinical evaluation was extensive
2 where we looked at transducer designs, looking at
3 transducer power, verification of ability to control
4 the focal spot. There's a lot of work done in
5 cavitation. Some of you may be familiar with focused
6 ultrasound in other applications, such as lithotripsy.
7 In that application, you're trying to generate
8 cavitation. The whole point is to generate a very
9 high energy shockwave to shatter a stone, such as a
10 kidney stone. In our application, we only want
11 thermal effects, and we want no cavitational effects,
12 so there's a lot of design in the system and the use,
13 limitations on the use to make sure that we limit our
14 effects to thermal effects.

15 There was testing of the biocompatibility
16 and a lot of animal testing in both -- for both our
17 system and in the literature. There are several dozen
18 publications over the last 15 years on the use of
19 thermal imaging in MR.

20 Next I'd like to introduce Dr. Clare
21 Tempany, who will give us an overview of MR anatomy
22 for treatment planning.

1 DR. NOLLER: If I could interrupt for just
2 a second. Our support personnel, could we have the
3 temperature turned down a little bit, whoever is doing
4 that.

5 DR. TEMPANY: Thank you. Good morning,
6 Mr. Chairman, panel members, and guests. My name is
7 Clare Tempany. I, too, am a Clinical Trial
8 Investigator at the Brigham. My trip and
9 accommodations have been paid for by the company. I
10 work as a consultant like Dr. Stewart for the
11 company, and work within the Harvard Medical School
12 Guidelines for Conflict of Interest and Ethics in
13 Research.

14 What I'd like to do for you today is two
15 things. I'd like to introduce you to the MR imaging
16 anatomy, display of anatomy and pathology that's used
17 in this trial. It's used routinely in clinical
18 imaging today, and then walk you through a typical
19 clinical treatment.

20 Female pelvic MRI has become a very
21 powerful diagnostic tool, and it's been available now
22 to us in radiology for over 15 years. It exclusively

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 displays the female pelvic anatomy as you see in these
2 what are called T2-weighted images for you on this
3 slide. On the left you see a sagittal view, and on
4 your right is a coronal view. And on the left, you
5 can see the anatomy of the uterus and cervix displayed
6 with the substructure of the zonal architecture of the
7 uterus displayed with the layers delineated for you.
8 And on the right you see the same thing with the ovary
9 on either side.

10 Many of you are more familiar perhaps with
11 pelvic ultrasound, and these are images of patients
12 with fibroids where you can see an enlarged uterus
13 here in the center, and then you see a slightly
14 different appearing uterus in the right side here.
15 The texture and tissue characterization of ultrasound
16 is somewhat limited to either solid or cystic, where
17 we can see the differences here with the cystic
18 component on the right.

19 A little bit of MR anatomy and how we
20 visualize these fibroids before we determine whether
21 they're eligible for treatment or not, selected images
22 here. Now we're going to walk through several planes

1 just to show you the display of the anatomy, and on
2 the left you can see an axial view with the patient
3 lying prone. The blue line represents the sagittal
4 image on your right, and all of the relevant
5 structures will be labeled, obviously, but you can see
6 a very typical uterine leiomyoma sitting here in the
7 center. It's classically a typical one that has a
8 very low signal intensity or it's black, and it has a
9 very sharp border. This is what we call a cookie-
10 cutter sharp border which delineates and
11 differentiates this from say adenomyosis, which will
12 not have such a sharp border.

13 The coronal plane here you can now see
14 nicely posteriorly as delineated up here on the blue
15 line, but way in deep at the back of the pelvis here
16 and the woman is standing in front of us, you can see
17 the sacral nerves coming down here posteriorly,
18 coming down along the lateral aspect of the pelvis to
19 exit through the sacrosciatic notch. And as we come
20 forward, you can see more anteriorly now. We're
21 coming into the uterus. We see the large fibroid. We
22 can see its relationship to the bladder. It's very

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 easy to understand some of the symptomatology this
2 patient has experienced when you see images like this
3 with a large uterine fibroid pressing on the bladder.

4 Now axial planes in a typical treatment
5 position now with the patient is lying again prone,
6 and you can see the uterine fibroid sitting here. We
7 see the anterior skin there, and you can see the
8 direction of the beam as you will see in a minute.
9 And there's the fibroid. These are the anterior
10 rectus muscles here anteriorally, and posteriorally we
11 see the fat, the bowel, and the sacral nerves.

12 We've learned a lot about uterine fibroid
13 or leiomyoma imaging over the years with MRI, and done
14 many pathological correlation studies, and have
15 determined that there are many types of fibroids, as
16 you've known in the clinical world for many years.
17 And these can be seen and characterized well in MRI.
18 And to just summarize some of them here for you where
19 you see about five different varieties described.

20 The top two are probably the typical ones
21 that we would treat in this trial, or have treated in
22 this trial, and these consist of the classic leiomyoma

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 which is a fiber muscular stroma. It's of low signal
2 in every imaging sequence we have. In other words,
3 it's black, it's easy to see.

4 A different type is what we call
5 hypercellular, where it's also known as the white
6 fibroid. It appears at high-signal intensity on T2-
7 weighted images. Those are the ones we've treated.

8 The other group will represent ones that
9 we wouldn't treat, which are non-enhancing leiomyomas
10 basically, once that have already undergone
11 spontaneous degeneration or necrosis in vivo, and
12 obviously, of varying patterns also.

13 Just to show you some more examples of the
14 range of the types of appearances of fibroids, here's
15 a woman who's had very significant fibroid burden.
16 Everything with an F on it is clearly a fibroid here.
17 This is a coronal T2-weighted image, as if she's
18 standing in front of us, her urinary bladder in white,
19 and you see how this may appear like a five-month
20 gravida uterus.

21 On the right side we see a different
22 patient with multiple fibroids and an unusual

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 appearing one here posteriorally that's already
2 undergone degeneration. This is a large cystic
3 degenerated leiomyoma, and we know it's degenerated
4 because we have post contrast images here in the
5 middle that's after the injection of intravenous
6 Gadolinium, and it shows no evidence of enhancement or
7 it stays black with no perfusion; thus, indicating
8 that it's necrotic. So let's just move into a
9 treatment process now.

10 Much of that imaging will occur prior to
11 the patient's being determined as eligible for the
12 trial, and we have identified, selected the patient
13 and identified the target treatment, and this is what
14 now happens on the day. So starting the night before,
15 the patient will receive written guidelines about the
16 therapy and what to expect during the treatment. She
17 will review that. She will have prepared the
18 abdominal wall, removing abdominal hair from the
19 umbilicus down to below the pubic bone. This is
20 important because we want the skin to be as smooth as
21 possible, and to not interfere with any coupling or
22 cause decoupling of the ultrasound beam as it will

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 transmit through the skin.

2 She remains NPO from midnight because we
3 use intravenous conscious sedation, and clearly don't
4 want to have any problems with food. So we have the
5 patient then come in the next morning. I meet with
6 the patient. We review the treatment guidelines with
7 her. We review what sort of sensations or experiences
8 she may feel during the treatment. We develop the
9 communication ritual to tell her when we're going to
10 do a sonication, she tells me what she feels, and we
11 sort of discuss all of that communication issue before
12 we go in the room at all.

13 I also then consent her for administration
14 of intravenous conscious sedation per our hospital
15 guidelines. Once that is done, the IV line is sited,
16 the Foley Catheter is placed. We use a Foley Catheter
17 clearly to control the bladder during the procedure to
18 make sure that the bladder stays empty. As you know,
19 when the bladder fills, the uterus moves, and treating
20 a moving target is clearly difficult, so we use the
21 Foley Catheter to control that.

22 At the same time in parallel, the room

1 check is going on. There's a phantom checking of the
2 system occurring, and after all of that is done, the
3 patient then comes into the room, is positioned on the
4 table in the coil, her vital signs monitoring devices
5 are placed in position, obviously her pulse ox, blood
6 pressure cuff, et cetera. The nurse will remain in
7 the room with her at all times, and both she and the
8 nurse will have a small little sonic button in their
9 hands which will allow them to terminate an individual
10 sonication should the patient experience an unusual
11 severe pain. She has full control of the therapy
12 itself at the time.

13 So here's just some pictures. You can see
14 this is the MRI magnet, this is the table, patient
15 sitting getting ready to go into position. She then
16 turns over and lies prone, positioning the pelvis over
17 the transducer. The transducer, as you've seen
18 already, is in the table surrounded by degassed water
19 so she lowers the skin down onto the water bath
20 basically with a gel pad also in side it, and she
21 makes direct contact with the skin into the water.

22 A little bit about our conscious sedation

1 and monitoring of the patient during the procedure.
2 We use standard intravenous conscious sedation
3 medications at our site. We use Versed and Fentanyl.
4 These are administered to provide a combination of
5 both analgesia and sedation. It's clearly important
6 that the patient's anxiety and any claustrophobia that
7 she may be experiencing in the magnet be aided by the
8 administration of these medications. Patients,
9 obviously, may experience positional pain lying on
10 their stomach in the magnet for the duration of the
11 procedure. And again, the analgesic effect is useful
12 for that. And we obviously want to try to reduce any
13 pain from sonication so we use the Fentanyl.

14 Typical doses that have been used in the
15 procedures, and these are the total doses, range from
16 as little as 25 mics of Fentanyl to 250 mics, to .25
17 to 5 of Versed. These are both intravenously. We
18 also give patients an oral non-steroidal anti-
19 inflammatory at the very beginning of the procedure.
20 Usually, typically 75 milligrams of Voltaren has been
21 used.

22 The medication then is given as required.

1 Before we start any treatment, we will give a very
2 small incremental initial dose of Versed and Fentanyl,
3 and then depending on how the patient is feeling,
4 responding during the therapy, we will give further
5 doses during the procedure, so that's why the ranges
6 are quite wide here. Some patients require very
7 little, some patients require a little bit more.

8 So let's just start now with treatment
9 planning. The patient is positioned on the table, and
10 you can see the transducer. And this is a good
11 positioning on your left here, as opposed to the one
12 on the right where the transducer is too high. And
13 you can see this is a very large field of view image
14 here. The uterus is really too low, and we would have
15 to angle too steeply to treat that, so we clearly can
16 readjust the transducer and the patient at this stage
17 before we start going any further.

18 We will then take three planes of pelvic
19 MRI images, as you've just seen, an axial, sagittal,
20 and coronal to again define our target, to allow me to
21 draw the contour to target volume superimposed on the
22 fibroid at that point, and those images are coming up

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 in a minute. We'll show you how we do that.

2 Just to remind you that in the trials, we
3 have protocol treatment guidelines. Single treatments
4 initially were limited to 120 minutes. The maximal
5 thermal dose per fibroid was limited to less than 100
6 CCs, and you could see treatment for all fibroids, if
7 more than one was treated, was a total of 150 Ccs.
8 We have a maximum of four fibroids that could be
9 treated in any one setting.

10 The protocol treatment guidelines
11 delineate a little bit further in detail here for you,
12 and the schema on the right really explains it all
13 nicely. The large black circle is a fibroid. The
14 smaller one on the inside the region of treatment, the
15 ROT, is the circle that I would draw as the sub-volume
16 in the fibroid. We have to work with the guidelines,
17 obviously, remaining within 15 millimeters of the
18 outer serosal lining, and 15 millimeters of the
19 endometrial lining. And so this clearly restricted
20 somewhat the volume of the fibroid that we could
21 actually treat during the initial safety and efficacy
22 evaluations.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Here are some pictures of the same thing.
2 You can see the treatment plan. Now the sonication
3 grid has been overlaid and you can see these
4 jellybeans, as Mr. Newman has already referred to, and
5 you can see them overlaid here on the images. Before
6 we do anything now, the next thing to do is to walk
7 through each of these sonications and determine is the
8 beam path going to be safe, and will it remain within
9 the guidelines. So we work through this system here
10 where we see the beam path on each and every one of
11 these. And there are some images now, just to show
12 you how that's done. You can see the passage of the
13 beam going through in green here, and the focal point
14 is delineated there on the sagittal view, the axial,
15 and the coronal.

16 What can be in the way? Well, things that
17 certainly can be in the way that we can identify
18 relatively easily are things like scars that would be
19 in the skin from prior surgeries, clearly things that
20 are in the skin such as scar can cause defocusing of
21 the ultrasound beam as it's passing through, cause
22 local heating of the skin, and something that we try

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 to avoid at all costs. And so it's fairly simple to
2 do this, we simply identify the scar ahead of time,
3 and then using the roll and tilt mechanism of the
4 transducer, we can angle around that area.

5 The same thing with bowel loops. Bowel
6 loops are relatively easy to see here. You can see
7 them on this T2-weighted image. They're the dark
8 structures up at the top, and you can see again, we
9 can angle either up and around the bowel loop, or if
10 necessary, just simply not treat that area, clearly
11 not go anywhere close to the bowel loop. The
12 sonication can simply be deleted.

13 Same thing here if we're looking at the
14 distal field. We can evaluate the location of the
15 sciatic nerves, and we can determine whether or not
16 the beam is going to pass through, and angle and roll,
17 and tilt again to avoid it.

18 Okay. So now we're ready to go. The
19 first thing we do is the geometric accuracy, and so
20 this is when a low powered sonication will be
21 delivered, and the very first set of images will come
22 up as that's being delivered, and this is a cropped

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 down view here of the face map, and you can see that
2 we're determining the accuracy; first of all, the
3 visibility of the sonication. Can I see it at all?
4 And you can hardly see it because it's covered by the
5 red cross, but underneath that little red cross is a
6 white dot, and that's the first sonication that's
7 being delivered. And the initial assessment is good,
8 now I see it. Is it in the right place? And so here
9 you can that it's off by about 5 millimeters, so we
10 will readjust all of the anatomical and adjust the
11 geometric alignment so that the green overlies the
12 red, and that they are absolutely concurrent.

13 The next step then is to move into a
14 therapeutic sonication dose, so we increase the power
15 up to typically 100, 140 watts, and we start the
16 actual procedure with therapeutic doses being
17 delivered. We compare this as it's going. We modify
18 the treatment parameters as necessary. As you'll see
19 in a minute, we're constantly looking at the feedback
20 mechanisms of the thermal imaging, to determine if
21 we've achieved a therapeutic dose or not.

22 Throughout the procedure I'm in constant

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 communication with the patient. This is very
2 important because there's a one-on-one communication
3 between myself, the patient, and the nurse in the
4 room, but I will talk to her, tell her we're about to
5 start a sonication at the beginning of one, and then
6 at the end of that ask her if she's experienced any
7 sensations, and if she has any concerns.

8 These are some examples now of the typical
9 dose profile. On the left, you see a sagittal view or
10 a long axis of a sonication, so you see the jellybean
11 shape. This is a short axis view where you see it on
12 end. And then these are three incidents of things
13 that could happen, so if you look at the bottom left,
14 we have a sonication that's achieving a thermal dose
15 that's probably too hot. The temperature, you
16 probably can't read that, I'm sure I can't either -
17 it's 100 degrees is what that one has reached, so
18 clearly, that's a little too hot. So what we do in
19 that situation is to back down on the power before we
20 go any further, so we wouldn't continue to treat
21 without changing parameters once we've seen that.

22 Similarly, in the opposite direction, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 next one demonstrates a sonication that achieved a 50
2 degree temperature, and that's too cold, so the first
3 step then would be to increase the power to bring it
4 up to a therapeutic dose, which we like to see between
5 60 and 80 degrees.

6 Calcifications occur in fibroids, as you
7 know, very frequently. They can be small punctate
8 little pieces of calcium within a fibroid, or you can
9 have a very more densely, heavily calcified one. The
10 latter patient doesn't usually get into the trial
11 because we can identify that in imaging, and a dense
12 rim of calcification precludes treatment using this
13 treatment modality. But small punctate calcifications
14 are impossible to see ahead of time, and this is what
15 may happen.

16 The ultrasound beam will be reflected off
17 the sonication, will simply not achieve any
18 therapeutic dose, so we simply move onto the next
19 location, delete that sonication, so to speak, and
20 don't treat that specific area. So an overview of the
21 treatment cycle is seen here for you for an individual
22 sonication. Before anything happens, the MR scanning

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 starts. Then the sonication is delivered, then a
2 tissue cooling period occurs, and at all times the
3 images are being acquired, and this total is about 2
4 minutes. So this is an extremely interactive
5 treatment.

6 As you've seen already, the entire beam
7 path is checked prior to delivery of sonication,
8 irregularities of skin, bowel, and beam path are
9 evaluated. We have multiple tools available to avoid
10 critical structures, things that we would not want to
11 have the beam pass through, and we use each and every
12 image to modify the next sonication, so it's a very
13 iterative process, so we're learning from the last
14 sonication what to do for the next sonication. And we
15 do this with the MR imaging that's continuously
16 occurring during the procedure.

17 So there are some safety issues,
18 obviously, where motion is a problem, if patients were
19 to move during the procedure, as I've already said,
20 heat treating a moving target is not good. So we
21 obviously have prevented that now by the Foley
22 Catheter placement with the bladder being controlled.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We also obviously coach the patient that she should
2 not move her pelvis. She's lying prone, and if any of
3 you have ever had an MRI scan, you know that we strap
4 patients in on the table, that there's a coil around
5 it, so it's fairly restrictive. There's not a lot of
6 room to maneuver and to move around, so these are
7 things obviously to our advantage.

8 We also use the restraint strap which is
9 strapped around the outer pelvis to hold the patient
10 on the table. The sedation somewhat helps also, but
11 clearly she can still move if she really so desires.
12 We monitor this with both sets of images. The real
13 time images being acquired during the sonication are
14 very easy to see motion, because it's like watching a
15 movie. You're sort of seeing a cine loop, so to
16 speak, so you can see changes if she was to move her
17 skin or her spine.

18 We also place fiducials at the beginning
19 of the imaging sequences, and these are the little red
20 marks you see here. And those are monitored
21 carefully, as well, to ensure that they don't change
22 in position over the procedure.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 The outcome assessment while the patient
2 is still on the table, essentially as we're going
3 we're developing this blue map in the center here
4 which represents all of the therapeutic sonications,
5 that have been delivered, and have been therapeutic;
6 in other words, reached the goal temperature delivery.
7 And so the blue is the area that we will expect to see
8 the necrosis. And at the end of the procedure, we
9 confirm this by injecting Gadolinium and evaluating
10 the necrotic tissue. And as you can see, nicely maps.
11 The blue area here is now seen as the black area here,
12 which is the non-enhancing or necrotic tissue.

13 Again, some other images from the end of
14 a treatment. Typical treatments look like this. They
15 can range from relatively small sub-volumes to
16 slightly larger volume here, with the areas of
17 necrosis seen in the area of treatment. And I thank
18 you for your attention, and pass the podium back to
19 Dr. Stewart, who will continue with the clinical trial
20 design.

21 DR. STEWART: Thank you. In moving on to
22 clinical trials in fibroids, that can be quite a

1 daunting task. As the Duke Evidence-Based Practice
2 report has shown on, despite the fact of a wealth of
3 clinical experience with uterine fibroids, this isn't
4 a lot of good evidence on which to base therapy.

5 We were fortunate in going into our
6 feasibility study having information from an in vitro
7 model using a rodent model, and using ultrasound
8 guided high-intensity focused ultrasound that showed
9 treatment with this energy modality was feasible for
10 uterine fibroids. And we wanted to get several
11 important things out of our feasibility study.

12 First of all, we wanted to make sure that
13 this was a safe treatment for women. We also wanted
14 to confirm our targeting accuracy. As Clare has
15 discussed, the feedback we get from the MR is
16 important, and we are depending on the non-enhancing
17 volume representing the tissue that we have
18 successfully ablated, so we did want to get pathologic
19 confirmation of this ablation. And this is actually
20 something that hasn't been done with previous
21 therapies, such as myolysis, cryomyolysis, or even
22 uterine artery embolization. And we wanted to take

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 this information to help refine our pivotal trial
2 design.

3 The study design was that it was an open
4 trial for women who were scheduled for hysterectomy.
5 They were to undergo MRI guided focused ultrasound
6 three to thirty days in advance of their
7 hysterectomies. In your panel packet, it appears that
8 there are two distinct studies that our center and St.
9 Mary's in London has described in one area, and then
10 the other three sites are described in another area.
11 However, because women were reluctant to go through
12 treatment and hysterectomy, recruitment in the
13 original cohort suffered, and so as time went on,
14 these other sites began recruiting patients, as well.
15 And then, in fact, the Israeli National Health Service
16 made hysterectomy optional for that group of patients.
17 They felt that it was unethical to require women to
18 undergo this therapy and then not have the option of
19 opting out of definitive therapy, so our trial design
20 changed somewhat midstream, but we followed all of
21 these patients, and reported them together.

22 We were able to confirm our pathological

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 information, and this is a diagram from our
2 manuscript. This is the treated fibroid, and this is
3 the pre-MR imaging that shows Gadolinium going
4 throughout the fibroid indicating good perfusion.
5 This is the post-treatment Gadolinium MRI where you
6 see a large area of non-enhancement. And then this is
7 the hysterectomy specimen. You can see on gross
8 examination that there is a clear lesion that
9 corresponds to the targeted area. And on microscopic
10 exam, there appear to be coagulative necrosis
11 corresponding to this area.

12 We were also able to confirm that there is
13 a relationship between the targeted volume, the non-
14 enhancing volume, and the pathologically correlated
15 area of tissue destruction. In this particular
16 fibroid from our St. Mary's site, you see the thermal
17 dose volume in A, the B is a little bit bigger, the
18 non-enhanced volume, and the pathologic area confirmed
19 more closely to this non-perfused volume.

20 We did find that the non-perfused volume
21 in general over-estimated the amount of tissue
22 destruction, but we found that in all cases the area

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of targeting was confined to the treated fibroid.
2 There was one case where microscopic evidence of
3 sonication was seen at the serosal border; however, in
4 retrospect, it appears that that was incorrectly
5 targeted. This was one of the cases where the bladder
6 filled and the target moved, and is a reason why we
7 adopted a Foley Catheter with our pivotal trial
8 treatment.

9 So we were able to confirm pathologically
10 that the tissue that we thought we destroyed was
11 destroyed. We also were very pleased with our results
12 in terms of patient treatment. All but one patient
13 were able to be treated as an out-patient. There was
14 a single hospitalization overnight for control of
15 nausea. There was no post embolization syndrome.
16 There, in fact, was very little pain in women
17 undergoing this protocol. And most of the patients
18 that we saw were not even taking over-the-counter
19 medications at the time we saw them within 72-hours of
20 their treatment.

21 The one safety issue that we did see in
22 this initial protocol was there, there was a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 significant incidence of infection seen post-
2 hysterectomy. They did not occur between the focused
3 ultrasound and the hysterectomy, but following the
4 hysterectomy. And we stopped the trial at the time we
5 saw the first three infections. We reviewed our
6 procedures. At that time, we did change our protocol
7 to institute prophylactic antibiotics. And once those
8 were instituted, we didn't see further significant
9 infections. And our pivotal protocol did not have
10 prophylactic antibiotic use.

11 We also used the information in the trial
12 to mitigate the adverse events we saw. We found early
13 on that paying attention to the skin in various forms
14 was important. Initially, patients were not shaving
15 and there were small skin burns at the area where
16 there may have been loss of coupling of the ultrasound
17 to the skin. We also, again, found the importance of
18 mapping the scars, and incorporating those into
19 treatment planning. Because scar tissue is very
20 similar to fibroid tissue, some of the energy would
21 stop at that point and patients would be
22 uncomfortable, and so we used a lot of the information

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 from this feasibility study to define the optimal
2 treatment protocol to embark on our pivotal study.

3 The major issue when embarking on the
4 pivotal study was the selection of a control group,
5 and there are always issues with picking the perfect
6 control group. And it's especially important, I
7 think, to put this in the context of the times. At
8 the time that the selection was going on, it was
9 December, 2001. Although uterine artery embolization
10 today might appear to be the best alternative, as a
11 control group, this was not really possible at that
12 time. There were no embolic agents that had received
13 FDA approval at that time. And with extensive
14 negotiations with the FDA and the investigators, we
15 looked at the other alternatives. And we felt that
16 looking at a surgical option would really give us
17 important safety information. It was important to
18 have a contemporaneously recruited control group, and
19 not to depend in historical controls.

20 Again, abdominal myomectomy in many ways
21 appears to be an important option. The issue for this
22 group of patients was that many of them may not be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 symptomatic. They would be pursuing treatment to
2 attain fertility. They would also tend to be younger
3 than the symptomatic patients that we were seeing.
4 And with our group, we specifically wanted to recruit
5 women with a threshold level of symptomatology.
6 Therefore, we decided that although no control group
7 was perfect, that abdominal hysterectomy would be the
8 best alternative.

9 With our knowledge of difficulties in
10 recruitment and our pivotal study, and also
11 information we were gaining from the experience with
12 uterine artery embolization trials, at this time many
13 groups were trying to perform randomized trials
14 between conventional surgical therapies and uterine
15 artery embolization. And no one succeeded in having
16 sufficient enrollment, so in that group of patients
17 there were generally case series or parallel controls.
18 And again, this is the study design that we settled
19 on.

20 The hysterectomy group and the focused
21 ultrasound group were enrolled in parallel. They met
22 the same inclusion and exclusion criteria, and both

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 received the same six month follow-up. We also chose
2 to separate the sites TAH and focused ultrasound so
3 that you did not have investigator bias channeling
4 good prognosis patients into focused ultrasound, and
5 bad prognosis patients into hysterectomy, so the sites
6 were all separated. And with our power calculations
7 we found that a 3-2 ratio would give us the desired
8 number of patient's in each arm.

9 The inclusion criteria included women who
10 were not pursuing future pregnancy. We felt it was
11 not ethical to treatment women who desired future
12 fertility until we had information regarding the
13 efficaciousness of this treatment. They were all pre-
14 menopausal or peri-menopausal women. They did have
15 both clinical exam and MRI consistent with fibroids.
16 The fibroids needed to be visible on contrast MR, and
17 feasible for treatment.

18 We also chose to have a minimum symptom
19 severity score, so that they had to score over 40
20 points on a scale of 100 to be included in this
21 protocol. The exclusion criteria were fairly obvious.
22 Women who could not undergo MR were not included.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Women with excessive uterine sizes in excess of 24
2 weeks, or women that were too heavy to fit in to the
3 MRI equipment were excluded. We also excluded anyone
4 with an undiagnosed pelvic mass, or other worrisome
5 pelvic pathology.

6 The primary hypothesis for our pivotal
7 study was that we would see at least a 10 point
8 improvement in the uterine fibroid symptom and quality
9 of life symptom severity score. This is the only
10 validated quality of life score specific for uterine
11 fibroids. And we felt that in our treated group, we
12 would have at least 50 percent of our patients
13 achieving this goal.

14 We realized that the treatment modality
15 would likely not be as effective as hysterectomy given
16 the limitations, but we felt that this was an
17 important landmark in demonstrating the efficacy.

18 We also evaluated several important
19 secondary hypotheses. We wanted to look at the
20 significant clinical complications in both arms to
21 compare safety. We wanted to look at the trajectory
22 of recovery, and also the costs involved.

1 For those of you not familiar with the
2 uterine fibroid symptom and quality of life measure,
3 this again is a disease-specific validated measure.
4 It was developed specifically for uterine fibroids and
5 it has two different parts. The symptom severity
6 score, which you'll see in this presentation referred
7 to as the SSS, has eight questions that relate to the
8 fibroid specific symptoms, pain, bleeding and bulk.

9 There is also a component to the health-
10 related quality of life which has six different sub-
11 scales as is common with all quality of life
12 questionnaires. And this questionnaire was developed
13 from an ethnically diverse set of focus groups to
14 really get input of fibroid patients, and what they
15 felt their significant symptoms were.

16 Also during the validation process, this
17 was correlated with the SF-36, which is really the
18 standard measurement of quality of life, as well as a
19 menorrhagia questionnaire indicating its comportance
20 with symptoms of menstrual blood loss.

21 For those of you not familiar with the
22 questionnaire, you'll see that the symptom severity

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 score addresses issues such as heavy bleeding during
2 your menstrual period, passing blood clots. It also
3 looks at bulk related symptoms, feelings of tightness
4 and pressure, frequency of urination or nocturia or
5 feeling fatigued. And patients are asked to rate
6 their symptoms on a five point Likert scale from not
7 at all to a very great deal.

8 This is data from the initial validation
9 of this questionnaire that you'll see the two parts
10 are divided here to the symptom severity score, and
11 these are the sub-scales of the health related quality
12 of life. One of the first things you'll notice is
13 that there's an inverse relationship between them.
14 For symptom severity score, the women in blue who are
15 women with uterine fibroids, have a higher score, so
16 higher scores mean higher symptoms. Whereas, with the
17 health-related quality of life, the normal women tend
18 to have higher scores and impaired related quality of
19 life is reflected in a lower score.

20 It's also interesting to note the absolute
21 levels of the symptom severity score. In this study,
22 looking at women with symptomatic fibroids, the mean

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 score was 44; whereas, the mean score for normal women
2 was 23 or about a 20 point difference between the two
3 groups.

4 This clinical difference was the primary
5 reason we selected our 10 point difference between
6 treatment success and treatment failure; that if 20
7 points represents the difference between women with
8 fibroids and normal, getting 50 percent relief of
9 symptoms appears to be an appropriate clinical end-
10 point.

11 There were also standard methodologic
12 reasons to choose this. That 10 points is very
13 similar to the standard deviation in the population.
14 The standard error of the mean and gives a moderate
15 effect size, as well.

16 We did not depend only on one outcome. We
17 also looked at additional efficacy measures. We used
18 the SF-36 which gives standard health-related quality
19 of life. We looked at several measures of disability
20 days, some assessment of an overall treatment effect,
21 and also patient's treatment satisfaction.

22 This is a schematic drawing of the pivotal

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 study design that at the screening visit we perform
2 the MR prior to the treatment, as well as the symptom
3 screening. A hematocrit ruled out serious anemia, and
4 during the treatment visit we again got information
5 regarding symptomatology.

6 We took seriously that this was a new
7 technology, and that there wasn't a lot of experience
8 with follow-up, so we have everyone come back for a
9 physical exam within a week so that we would not miss
10 important issues that arose, so patients came back and
11 did have a hematocrit and a physical exam at that
12 time.

13 The one month and the three month follow-
14 up were generally by phone, but then there was a full
15 visit at six months with a physical exam and MR exam,
16 and again complete testing.

17 The pivotal study design was originally
18 designed to have outcomes at six months. However,
19 later we have extended follow-up so that we're now
20 seeing patients who are continuing on at 12 months, 24
21 months and 36 months. And again, getting information
22 on quality of life, as well as MR exams at that time.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We wanted to try to capture significant
2 clinical complications, and what we did at this time
3 was we went to the literature. The paper by Dicker,
4 et al, arose out of the collaborative study of
5 sterilization. And they felt it was important at that
6 time to try to define characteristics that could be
7 used to compare treatment.

8 We used their criteria, but tried to
9 update it both for the change and length of stay that
10 has occurred since the 1970s, and also some of the
11 differences that we would potentially see with this
12 new therapy included as additions or things like
13 discharged going to a rehabilitation facility,
14 discharged with either a catheter or a drain, or also
15 various interventional treatments that may not qualify
16 under their definition of surgical procedures.

17 While this would seem to favor picking up
18 complications from hysterectomy, I think it's
19 important to remember that if there had been
20 inappropriate targeting and significant injury of
21 adjacent structures, these complications would have
22 been seen and picked up if the treatment had

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 significant side effects in that way.

2 So moving onto the results of the trial,
3 as we talked about earlier for the pivotal trial,
4 there were separate sites for hysterectomy and MRI
5 guided focused ultrasound. There were three U.S.
6 sites and several through Europe and Israel. There
7 were also hysterectomy groups and about half of the
8 enrollment for both arms came from the U.S., and half
9 from out of the U.S.

10 There was fairly equal distribution of
11 patients through the sites. There wasn't a primary
12 site that contributed all of the patients. And we
13 looked at the demographics between the patients
14 undergoing focused ultrasound, and the patients
15 undergoing hysterectomy. We knew that since this was
16 not a randomized trial, there were likely to be some
17 differences. We did find them similar in age, and
18 fairly typical for women with fibroids. There was a
19 statistically different finding in body mass index
20 with the women undergoing hysterectomy being somewhat
21 heavier. And both groups of women had significantly
22 elevated symptom severity scores.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 As you'll recall in the validation study,
2 the women with fibroids typically had scores in the
3 40s, and both of our groups this mean score was over
4 60. And again, there was a difference between these
5 two groups with women who had elected definitive
6 therapy for hysterectomy having a somewhat higher
7 score.

8 There were more black women in the
9 hysterectomy group. Again, probably a relationship of
10 site selection, but all women in both groups were pre-
11 menopausal by and large.

12 There were some differences in co-
13 morbidities. The women undergoing hysterectomy were
14 more likely to have diabetes and hypertension, and
15 the women undergoing focused ultrasound were more
16 likely to have thyroid disease. As you'll see later
17 on, we looked at these differences between the focused
18 ultrasound group and the hysterectomy group to see if
19 these differences affect the treatment outcome.

20 We did perform an intention to treat
21 analysis, so that every patient who received focused
22 ultrasound is included, and so our denominator in the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 slide you'll see is 109 patients. There were three
2 withdrawals from the study less than six months, and
3 11 patients were non-evaluable. We did, however, do
4 calculations for both evaluable patients and intention
5 to treat patients, and they were similar.

6 The characteristics of the fibroid
7 patients were consistent with women who had
8 symptomatic fibroids. The average uterine volume was
9 approximately 600 Ccs, but there were clearly a number
10 of women who had uteruses in the range of 1,000 cubic
11 centimeters or more. The average total fibroid load,
12 meaning calculating the volume of the fibroids without
13 the myometrium was in the range of 300 to 400 cubic
14 centimeters. And patients had an average of two to
15 three fibroids, but as many as 12. And although one
16 to four fibroids could be treated during this
17 protocol, in a average most women got one treated.

18 We excluded from treatment fibroids that
19 were amenable to either hysteroscopic or laparoscopic
20 myomectomy, so although these say submucosal and
21 subserosal, they were probably more accurately
22 classified as partially submucosal or partially

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 subserosal with a large intramural component. And we
2 also looked at differences in location when we
3 assessed treatment outcome.

4 . So when looking at the treatment
5 parameters again with the intention to treat patients,
6 the baseline fibroid volume was about 300 Ccs. The
7 non-perfused volume at the end of treatment was in the
8 range of 68 cubic centimeters, so we had approximately
9 24 percent of the fibroid that had been treated during
10 this protocol. That at six months, there had been a
11 decrease in size from about 330 to 295. This
12 percentage of shrinkage is similar to the non-perfused
13 volume. Again, it's not a large absolute number, but
14 it is proportional to the amount targeted for
15 treatment.

16 Looking at our primary efficacy and the
17 symptom severity score, again we hypothesized that at
18 least 50 percent of our patients would have a 10 point
19 improvement. We were substantially in excess of that.
20 Over 70 percent of our patients reached this targeted
21 improvement, and this was statistically significant.

22 We also found that in fact the symptom

1 severity score at entry was in the range of 60. By
2 three months there was already clear evidence of a
3 treatment effect with a mean treatment level going
4 down to 41, and then some continued improvement
5 between three months and six months. And you can also
6 see here, this is the criteria we set for entry, so
7 many women at the three or the six month time point
8 would not have had symptoms sufficient to qualify for
9 enrollment if they had come at that point in time to
10 seek treatment.

11 This is the distribution of changes in
12 symptom severity score, so again this line indicates
13 the threshold for success, or 10 points or more.
14 These are the patients who had no improvement, or one
15 to ten points of improvement, so everyone from here
16 over is a treatment success.

17 The mean patient improvement, however, was
18 about two and a half times what we had predicted and
19 the mean treatment improvement was approximately 24
20 points. There were, however, some patients who
21 improved as much as 60 points in symptom severity.

22 When we turn to look at the health-related

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 quality of life subscales, these parallel the changes
2 that we saw in symptom severity score. Because of the
3 inverse relationship these lines go up rather than
4 down, so again you see a significant change or marked
5 change between baseline and three months, and then
6 some improvement from three to six months.

7 We use the SF-36 to be able to compare
8 more accurately the patients between the focused
9 ultrasound and the hysterectomy arm. What we see
10 again in the focused ultrasound group is the same
11 pattern of improvement that already at one month
12 you're seeing improvement in some scales, continued
13 improvement at three months, and stabilization from
14 three to six months. In contrast, the women who
15 underwent hysterectomy had marked impairment in some
16 of their functioning at one month, and it took them
17 three months to six months to get back to where they
18 were and, in fact, to note improvement following the
19 treatment.

20 The significant difference is in terms of
21 disability between the two groups. I think it's
22 important to note not only the differences between the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 groups, but the absolute level for the focused
2 ultrasound patients. When looking at the days -- this
3 is follow-up at one month following treatment. There
4 were only 1.4 days of missed work on average for the
5 women in the focused ultrasound group. Whereas, women
6 undergoing hysterectomy clearly has much more short-
7 term disability with 18 days. And parallel the days
8 that women with focused ultrasound were kept from
9 their normal activities averaged about three days.
10 And they again spent only about a day and a half in
11 bed, so these numbers demonstrate the significant
12 improvement and short-term recovery seen with this
13 treatment.

14 We also looked at resource utilization
15 through six months. Because of our different sites in
16 different countries we didn't bring this down to
17 dollars, but looked at encounters with the healthcare
18 system. This takes into account not only all of the
19 scheduled study visits for the MRI guided focused
20 ultrasound patients, but for those patients that
21 elected additional therapy, or went on to additional
22 procedures. All of those resource utilizations are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 captured.

2 We found that there was a significantly
3 different length of stay. Only 1 percent of our
4 focused ultrasound patients stayed more than five
5 hours post treatment. They also had substantially
6 fewer provider encounters, fewer additional
7 procedures, and fewer diagnostic tests.

8 We looked at a logistic regression model
9 to see if our baseline differences affected outcomes,
10 so in the model we included not only the things that
11 differed between our groups, such as race and BMI, but
12 also looked at other variables of interest, such as
13 age, country of treatment, fibroid location, percent
14 non-perfused volume. And the only predictor of success
15 was baseline symptom severity score. In other words,
16 the most highly symptomatic patients were the patients
17 that improved the most.

18 We also looked at patient satisfaction and
19 asked patients were they satisfied with their
20 treatment, was it effective in eliminating their
21 symptoms, and would you recommend this to a friend?
22 And again, over 70 percent of women answered

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 affirmatively to these three questions.

2 We did continue to follow patients beyond
3 the pivotal study, and attempted to bring patients in
4 for follow-up between six and twelve months. Again,
5 we start with our intent to treat population of 109.
6 We found that 91 patients continued on, 9 patients
7 declined to be included in the follow-up. They had
8 enrolled for a six month trial, and elected not to
9 come back, and 9 were withdrawn, which left us with 82
10 evaluable patients at 12 months.

11 We found in following this group that 23
12 patients had gone on to alternative therapies, and
13 four patients had elected and were offered additional
14 focused ultrasound treatments. Both of these groups
15 of patients are then included as treatment failures in
16 our 12-month analysis.

17 So the original study was, indeed,
18 designed for six month follow-up and we did contact as
19 many patients as we could to come back. Because of
20 the date that we started to do this, there was some
21 lag, so although it's reported as 12-month follow-up,
22 the actual mean follow-up was approximately 14 months.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 The success rates do not look as promising at this
2 point. If we look at our intent to treat group,
3 there's only approximately a 38 percent success rate.
4 And again, the patients who declined to come back for
5 us to follow-up or chose alternative treatments are
6 included here. And if you look at our evaluable
7 patients, it is slightly higher at 51.

8 I think what's notable is that there were
9 a substantial number of women who were still improved
10 with the mean treatment being targeted at
11 approximately 20 percent of their fibroid load. The
12 other thing that is interesting about the results at
13 12 months is that we still could see significant
14 decrements in the treatment parameters as measured by
15 the symptom severity score, so that at baseline again,
16 we're coming in at about 61 points, and going down to
17 points in the mid to high 30s at six months and twelve
18 months.

19 Part of the issue with the twelve month
20 data may be that fibroid symptoms returned. This is
21 clearly a common problem in the literature, and is
22 well described for myomectomy. Again, it appears to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 be an issue that may be applicable to uterine artery
2 embolization, as well. But again, many of the studies
3 with uterine artery embolization also have relatively
4 short-term follow-up. And I think that our original
5 treatment parameters were aimed at the maximization of
6 safety and, therefore, may not have optimally targeted
7 the amount of fibroids to get sustained treatment.

8 However, there were still significant
9 patient satisfaction with treatment success at 12
10 months. Again you see in blue the six month data, and
11 the twelve month data in yellow, so the patients were
12 still very happy with the treatment option that they
13 had pursued.

14 Turning our attention to safety, I think
15 it's important, first of all, to note what we did not
16 see; that many devices that are approved have
17 significant complications. In many case series, there
18 have been patient deaths or urgent unintended
19 procedures. There were none of those in this
20 treatment. There were no bowel injuries. There were
21 no hospitalizations for pain control or post
22 embolization syndrome. So compared to some concerns

1 that we had at the beginning, we were very happy that
2 there were not severe safety issues that we
3 encountered.

4 Looking at a strict definition of adverse
5 events, we found that 19 percent of patients in the
6 focused ultrasound group had no adverse events
7 compared with 1 percent in the total abdominal
8 hysterectomy. We chose for this protocol because of
9 its novel technology to strictly define adverse events
10 more similar to what you would see in a drug study
11 than a typical device study. We knew that this was a
12 device that didn't have clear predicates, and we
13 wanted to make sure that we were not missing adverse
14 events.

15 However, we found that when we looked at
16 device or procedure-related serious adverse events, we
17 still did very well with only 2 percent of MRI guided
18 focused ultrasound patients having serious adverse
19 events compared to 13 percent in our contemporaneously
20 enrolled group.

21 We found that the body systems in which
22 adverse events were found to be similar in most cases.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 On average, women undergoing focused ultrasound had
2 about two adverse events versus four for the total
3 abdominal hysterectomy.

4 We also wanted to define what we thought
5 prospectively would be device or treatment-related
6 adverse events. We that non-significant events might
7 include fever or pain in the treatment area, swelling
8 or firmness in the treatment area, or minor skin
9 burns. However, we felt that either skin burns that
10 caused ulceration or any kind of nerve damage should
11 be termed significant anticipated events, and that we
12 were especially looking for these events as treatment
13 unfolded.

14 Again, as we saw in our feasibility study,
15 there was a substantial decrease in the amount of pain
16 patients had both during this procedure and post
17 procedure, compared to some alternative therapies.
18 Interoperatively, the patients reported on average
19 mild discomfort and mild to moderate pain. And then
20 at post procedure, their levels of both pain and
21 discomfort were significantly closer to no pain at all
22 than to mild.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 There were patients who had some severe
2 pain during the procedure. And as Dr. Tempany
3 discussed, we do have the ability to redose pain
4 medications during the procedure. Only one patient,
5 however, related her pain as severe post procedure.
6 And again, we found that narcotic use following the
7 procedure was very rare, and that even over-the-
8 counter medication use was rare in the days following
9 treatment.

10 We wanted to look at adverse events again
11 to see if our baseline differences and the co-morbid
12 conditions or the demographics affected these
13 outcomes. Clearly, there were some co-morbid
14 conditions in the hysterectomy group that may have
15 made them more likely to experience complications. We
16 found, however, that in controlling for this, the odds
17 ratios still showed that there was significantly
18 increased risk of dermatologic, gastrointestinal CNS
19 and pain adverse events in the group undergoing
20 hysterectomy compared to the group undergoing focused
21 ultrasound.

22 Again, we wanted to look at the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 significant clinical complications to make sure that
2 we captured significant events, and again use the
3 literature to prospectively define this. We found
4 that the patients undergoing hysterectomy were more
5 likely to have a significant clinical complication
6 with about 46 percent of the group in the hysterectomy
7 group having an adverse event, as opposed to 12
8 percent in the focused ultrasound group.

9 One of the interesting comparisons is
10 looking at fever and antibiotic use, given that the
11 patients in the focused ultrasound group did not
12 receive prophylactic antibiotics; whereas, the
13 patients undergoing hysterectomy traditionally did.
14 And still, the incidents of fever and antibiotic use
15 started after the prophylactic antibiotics for
16 presumed infection were lower. The transfusion rate
17 was also low. There were no unintended surgical
18 procedures, no discharges with appliances. There were
19 several rehospitalizations, but none requiring
20 interventional treatment, and no death or life-
21 threatening events.

22 We also found that there were differences

1 in -- there were significant differences in clinical
2 complications. And our most serious adverse events
3 included device-related adverse events; that we found
4 that there were several instances of leg pain, which
5 again we had identified as an anticipated event.
6 There were also some skin burns, although most were
7 first and second degree burns that resolved easily.

8 Our most important device-related, and our
9 only device-related SAE involved a patient who had a
10 treatment where there was injury to the sacral nerves.
11 I think this is the case that pointed out to us the
12 importance of having patients talk to us about their
13 pain and discomfort. And that this patient did not
14 receive pain medication at her request, and was noted
15 at post treatment time to have weakness and nerve
16 conduction studies confirmed injury.

17 However, by 12 months she has resumed a
18 high level of physical activity, and in fact has run
19 a marathon since her treatment. She had significant
20 symptom improvement, and continues to be a part of our
21 study.

22 We also put in a number of steps to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 mitigate the risk of nerve damage. As Dr. Tempany
2 talked, there are several things that can be done in
3 the treatment planning and the use of feedback from
4 the patient. Since we instituted these measures,
5 there's been no significant nerve injury, and the
6 incidence of nerve injury has been minimal. So we
7 also did a number of simulations that allowed us to
8 look at this issue. And so again, we had one event
9 since learning from this important case.

10 Looking at the serious adverse events,
11 again we classified everyone that was hospitalized as
12 having an adverse event, an SAE, even if it was felt
13 not to be device or procedure-related; again, sacral
14 nerve injury, nausea. And then there were four women
15 that went on to additional therapy, which we felt was
16 really progression of disease and not device-related.

17 There is one complication on commercial
18 treatment that resulted in a patient death. A single
19 patient in one of our outside the U.S. sites had a
20 pulmonary embolism following commercial treatment.
21 This was investigated by the local M&M committee, and
22 it was felt that her death was not related to the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 procedure. And in fact, in retrospect it turns out
2 she had several important thrombotic risk factors that
3 had not been identified.

4 We do have a continued access protocol,
5 and have continued to treat patients. It's very
6 similar to our pivotal study with mild changes in the
7 treatment parameters that allow slightly increased
8 treatment time and treatment volume. And we've been
9 enrolling patients in this protocol since April of '03
10 with 89 patients treated to-date. The adverse events
11 in this group have been significantly less than in our
12 pivotal study, and indicate that our mitigation steps
13 have been successful. And we don't have enough of
14 these patients to six months to comment on efficacy,
15 but the three month efficacy appears similar to the
16 pivotal study.

17 So in summary, we only had one device-
18 related SAE, and a low incidence of adverse events.
19 We confirmed that this treatment can be safely
20 performed as an out-patient, and have learned from our
21 experience to design a safer study protocol.

22 We also found that we met our primary

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 efficacy point with a significant margin. We had a
2 much lower symptom severity score than we had
3 predicted, and all of the measures of improvement tend
4 to move together to show patient improvement.

5 I'll turn the program over to Rob to talk
6 about the training.

7 MR. NEWMAN: I'd just like to speak
8 briefly to amplify on the information that's in the
9 panel packet about the training program. We believe
10 that this is truly a non-invasive surgical
11 alternative. This is a scalpel of sorts, a non-
12 invasive one, but it is a scalpel. The physician
13 controls the delivery of therapy, and the system
14 provides the ability for real-time interactive control
15 of that looking at the results from the treatment
16 itself.

17 The system works only a 1.5T MRI system.
18 We believe that this is necessary. It's the current
19 state-of-the-art for pelvic imaging for assessment of
20 anatomy and pathology. And it also gives us the image
21 quality that we need for accurate temperature
22 measurements. These symptoms have a high level of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 service, and are in wide use throughout the medical
2 community.

3 This system will only be used under the
4 direct supervision of trained physicians. This is not
5 something that would be used by anybody else. We
6 believe that the gynecology and radiologic expertise
7 is required, and the nursing requirements for these
8 kinds of treatments are similar to what is currently
9 being used in hospitals for regular interventional
10 radiological interventional control, so there's
11 nothing unique about that part of the treatment.

12 The training for all installations will
13 include the entire team, doctors, the MR technologists
14 and nursing. It's divided into two phases. One is
15 the system operation, the technology side of it. The
16 other part will be the clinical issues, which will be
17 covered by preceptorships at clinical sites involving
18 topics of patient selection, treatment planning,
19 anesthesia, adverse event management and those kinds
20 of things.

21 First, treatments will be supervised. And
22 on our system, every sonication on our system is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 recorded and kept in a file, so we have a log of every
2 treatment we've ever done. This allow us both to
3 review prior treatments if you had an adverse event,
4 or if you've had something interesting. It also
5 builds us a continuously growing teaching file that we
6 can use for future sites.

7 Just a brief overview, the kinds of things
8 would be what you expect we would cover in the
9 classroom part of it on system components, and the
10 physiology, device, protocol development, and we would
11 follow this up with training after the procedures have
12 begun at a specific site.

13 InSightec has a continued commitment to
14 studying MR guided focused ultrasound. We think that
15 this is -- there's an ongoing process here, a lot we
16 can learn. As we've described before, we have the
17 continued access protocol is in progress. We've
18 treated 89 of 250 patients, and we intend to complete
19 that 250 patients and collect three-year follow-up
20 data on them to look at -- to gather more data on
21 safety and efficacy of the system. And will provide
22 us a lot of information on improvements in treatment

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 planning, and ways to make it more effective. And we
2 also have additional studies ongoing outside the
3 United States, and will include the analysis of that
4 in our development of future features.

5 DR. STEWART: So in summary, I think we've
6 demonstrated to you that the device that we're
7 presenting has a low risk of serious adverse events.
8 We were very careful to try to capture all events that
9 occurred, and to report as completely as we could to
10 make sure that this novel technology did not have any
11 unintended side effects that we were missing.

12 One of the important issues with this
13 technology is that it is fibroid-specific. And I
14 think that that has benefits beyond what we've
15 demonstrated today. The risk of complications is
16 significantly lower than hysterectomy. And I think if
17 we had chosen other control groups, we would have
18 probably been able to demonstrate significant
19 differences with other treatment modalities.

20 We've had a very low incidence of device-
21 related events. And because this technique employs
22 conscious sedation rather than anesthesia, there is

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 also a decreased risk of anesthesia-related events.

2 We have seen a clinically significant
3 improvement in these patients. Patients are very
4 vocal about voicing their improvement with this
5 treatment, and we have been able to capture that by a
6 number of different modalities. We designed our study
7 and well-exceeded both our primary and our secondary
8 end-points. And to be able to gain this kind of
9 improvement without surgical incision, without major
10 disability I think is a major step forward. The fact
11 that these procedures can be performed as out-patients
12 is important, as is the fact that it preserves the
13 uterus.

14 Many women, I think, with fibroids tend to
15 live with their symptoms rather than go through some
16 of the treatment options. Some women have significant
17 disability that they put up with day in and day out
18 because of their concerns regarding invasive
19 therapies. And I think MRI guided focused ultrasound
20 surgery gives us an important new choice, and an
21 important choice to help reduce the symptoms of
22 uterine fibroids for women. Thank you.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. NOLLER: Thank you. We very much
2 appreciate the sponsor staying within their time
3 limit.

4 Our next presentation will be by the FDA.
5 By the clock we are using up here, it is now 10:13.
6 We will take a break until 10:30 by this clock, for 17
7 minutes, and then the FDA will make their
8 presentation.

9 (Whereupon, the proceedings in the above-
10 entitled matter went off the record at 10:07:53 a.m.
11 and went back on the record at 10:26:03 a.m.)

12 DR. NOLLER: Okay. We'll reconvene now,
13 please. And again, I'll ask the panel to hold its
14 questions until after the FDA presentation. At that
15 time we will I think have about 30 minutes to
16 formulate and ask some questions. I'd like to
17 introduce Kathryn Daws-Kopp, who will lead us through
18 the FDA presentation.

19 MS. DAWS-KOPP: Good morning, ladies and
20 gentlemen, distinguished panel members and guests.
21 I'm Kathy Daws-Kopp, the Lead Reviewer for FDA on this
22 PMA. My presentation will give a brief overview --

1 DR. NOLLER: Excuse me. Turn to the sound
2 up. We can't hear her.

3 MS. DAWS-KOPP: Okay. Good morning. I'm
4 Kathy Daws-Kopp, the Lead Reviewer for FDA on this
5 PMA. My presentation will give a brief overview of
6 FDA's review process on this PMA to orient you for the
7 remainder of the FDA presentations.

8 You may notice as we go through our
9 presentation that you'll be hearing some of the same
10 things that the company said. Our intention is to
11 focus on the issues we felt were important in our
12 review of the file.

13 I'll start off by describing the history
14 of regulatory interactions with the company, and I'll
15 describe components of the device from a regulatory
16 perspective. I'll provide a list of the PMA review
17 team, and briefly discuss what we did in reviewing the
18 PMAs, and I'll follow that with a list of some major
19 issues that are still ongoing with this review, some
20 of which are part of the panel discussion questions.
21 I'll close with an agenda of the remaining FDA topics
22 and presenters.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 This is a brief overview of the history of
2 FDA review on this device. The sponsor first came to
3 FDA with a feasibility study in 2000. That file was
4 reviewed by another branch in FDA, the General Surgery
5 Devices Branch, who consulted with us on the file.

6 In late 2001, our branch took over review
7 and the sponsor came to us to discuss a pivotal study.
8 The study was given conditional approval in March of
9 2002, followed by full approval in May. We worked
10 with the company on the protocol, and the study
11 includes as you've heard both U.S. and foreign sites.

12 In 2003, when they had completed
13 enrollment of the pivotal trial, the sponsor requested
14 permission to conduct a continued access study which
15 allows the company to continue to enroll patients
16 while they're working on preparation of a PMA, and
17 while the PMA review is ongoing.

18 For a number of reasons, the proposed
19 protocol for the continued access study differ
20 somewhat from the pivotal study. The continued access
21 study was given conditional approval in June of 2003,
22 and full approval in August. We received the PMA

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 submission on January 27th, 2004, and I'd just like to
2 note that that file received expedited review status.

3 The ExAblate system is made up of the
4 following basic components; patient table, operator
5 workstation, software, equipment cabinet. The patient
6 table is a standard MR table that has been modified to
7 house the ultrasound transducer and associated
8 equipment, and was already described by the company.

9 It should be noted that the MR system is
10 a commercially available GE device, the Signa 1.5T MRI
11 system is not commercially approved for thermography
12 at the site. Software in the ExAblate device uses MR
13 information from the GE device for mapping and
14 targeting, as well as these new thermography
15 functions.

16 This is the indication for use the company
17 has already presented, but we'd like to go over this
18 again. ExAblate is intended for use in pre and peri-
19 menopausal women with symptomatic uterine fibroids.
20 Patients must have a uterine size of less than 24
21 weeks, and be family complete. The fibroid or
22 fibroids to be treated must be visible on non-contrast

1 MR and should enhance on contrast MR imaging.

2 This is a list of the review team. As you
3 can see, a number of people have been involved in the
4 review of this PMA application in the areas of
5 clinical, statistical, epidemiology, MRI, ultrasound
6 software, bioresearch monitoring, patient labeling,
7 human factors, and manufacture.

8 This slide lists the things that we look
9 at during our review. For software and hardware we
10 look at safety and effectiveness. Examples of safety
11 issues for software and hardware include electric
12 shock, EMI shielding, and unintended burns. Examples
13 of effectiveness are adequate targeting and thermal
14 dose delivery.

15 We specifically look at requirements in
16 testing. We check to see that the device is designed
17 to do what the sponsor or manufacturer says it will
18 do. And we look to see that they do tests that check
19 to see that it works the way it's supposed to.

20 For bioresearch monitoring, we look at
21 study execution, including recordkeeping and informed
22 consent administration, as examples. For

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 manufacturing, we look at compliance with design
2 controls both included in inspection. Bioresearch
3 monitoring inspects clinical sites, as well as any
4 records related to the conduct of the trial at the
5 sponsor's facility. Manufacturing connects an
6 inspection at the manufacturing facilities.

7 Bioresearch monitoring inspection is
8 common for clinical trials, but is not required. A
9 pre-approval manufacturing inspection is required.
10 Drs. Corrado and Del Mundo will address clinical and
11 statistical reviews during their presentations.

12 This is a list of our current major
13 ongoing issues. This is not a comprehensive list of
14 all issues. We are still discussing the thermal
15 accuracy of the system with the company. Dr. Loren
16 Zaremba will discuss this further in his presentation.
17 We're still discussing adverse events that occurred,
18 and appropriate medications to employ in response to
19 these events. This will be discussed further by Dr.
20 Noel Del Mundo. We will also discuss how the
21 treatment in control groups differed, which Dr.
22 Corrado will be discussing in her talk.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 A pre-approval inspection is required, as
2 I mentioned. FDA is working to get this inspection
3 completed in a timely manner. Review of the labeling
4 for a device is an integral part of the scientific
5 review; however, we do not complete our review of
6 labeling until we have finished the rest of our review
7 of the file. These last two items, inspection and
8 labeling will not be discussed further by other
9 presenters today.

10 The rest of FDA's presentation will
11 proceed as follows. Dr. Corrado will provide a
12 summary of the clinical study and results. Dr.
13 Zaremba will discuss the MR thermal mapping review.
14 Bruce Herman will discuss the ultrasound-related
15 review concerns, and Dr. Del Mundo will close FDA's
16 presentation with a safety analysis discussion that
17 will cover what we have considered most significant
18 adverse events. Thank you for your time and
19 attention, and I will now turn the floor over to Dr.
20 Corrado.

21 DR. CORRADO: Thanks a lot, Kathy. Good
22 morning, everybody. I'm Julia Corrado, and I'm a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 member of the review team.

2 You have all already heard about the
3 clinical trials of ExAblate from Dr. Stewart and Dr.
4 Tempany, and I am going to be covering some of the
5 same material, but I'm going to try to give an FDA
6 perspective on that material. And I will try very
7 hard to avoid unnecessary redundancy.

8 I'm going to be starting with a brief
9 description of the feasibility study. I will then
10 describe in more detail the pivotal clinical study,
11 and the aspects of that study as you see here. And
12 finally, I will give a very, very brief synopsis of
13 the continued access study.

14 I'd just like to say who the -- normally
15 we don't spend much time talking about the feasibility
16 study at panel meetings, but this one was especially
17 important because it signaled to us a couple of
18 aspects of this treatment that we really wanted to
19 scrutinize closely when it came to the pivotal study.

20 This feasibility study was prospective.
21 It was non-randomized. It was conducted at two
22 centers, and I'll just digress for a second. Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Stewart described five centers. There was an IDE
2 pivotal study that was conducted under FDA approval,
3 and that was conducted at a center in the U.S. and one
4 in Britain. And I'm speaking just about that
5 feasibility study in my next couple of slides.

6 It was a pre-hysterectomy study. The
7 women who volunteered were scheduled for hysterectomy,
8 but they agreed to undergo the ExAblate procedure
9 approximately a month prior to hysterectomy. And we
10 approved the study for 15 subjects and 13 subjects
11 received treatment.

12 The objectives were already described by
13 Dr. Stewart. There were, in general, two types of
14 tissue effect that are noted from ExAblate. I won't
15 speak about them further, but there is a thermal
16 coagulative necrosis and then there is an ischemic
17 necrosis. The difference is that the thermal
18 coagulative necrosis is caused by direct heating, and
19 the ischemic necrosis results from lack of blood flow
20 to surrounding tissue following heating.

21 In the summary of the feasibility study,
22 the pathologist from Brigham & Women's described the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 tissue effect as follows; that the volume of necrosis
2 was sometimes larger than the treated area. That's a
3 very important point that I'm going to be emphasizing.
4 The treatment effect consists of bland and highly
5 uniform coagulative-type necrosis with relatively
6 sharp outline, scattered interstitial hemorrhage, and
7 variable amounts of acute inflammation consisting
8 mostly of neutrophils.

9 The next point also should be noted, and
10 that is that the only abnormality noted in the
11 myometrium outside of the fibroid, this was beyond the
12 fibroid capsule, was microscopic coagulative necrosis
13 extending one to two millimeters beyond the fibroid.
14 This is the only case where we saw this effect, that
15 there was a treatment effective beyond the fibroid
16 capsule. But nevertheless, we thought it was
17 important, as I'll describe further.

18 The purpose of the next slide is to
19 illustrate something I just hinted at, and that is
20 that the volume of effected tissue is different from
21 the thermal dose volume; that is, the volume that was
22 actually targeted. And there are two volumes that we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 can talk about from the feasibility study. One is the
2 non-profuse volume immediately following treatment.
3 This is on, I believe, T1-weighted images with
4 Gadolinium enhanced MRI. But also from this
5 population, most of these women underwent hysterectomy
6 so we also have volumes from hysterectomy specimens.
7 And what I'd like you to notice here is that there is
8 a consistent -- the non-profuse volume and the volume
9 from histology are consistently greater than the
10 thermal dose volume, which led us to feel that we
11 wanted to be cautious in how the pivotal clinical
12 study was conducted because we did not want to get
13 injuries resulting from tissue necrosis beyond the
14 targeted area.

15 As always, as we would expect during any
16 kind of a clinical study of an investigational device,
17 problems were encountered during treatment. For
18 example, several patients received what was described
19 as sub-optimal treatment due to excessive fat layers
20 within the beam path. And in one case, the portion of
21 the fibroid that the clinician wanted to treat was too
22 close to intestine, and that limited treatment in that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 case. In three cases, patients did not receive
2 treatment due to tissue aberration and scar in the
3 beam path that caused the patient to experience pain.

4 FDA, of course, always looks closely at
5 adverse events and clinical trials, and we saw the
6 following. But before I go into these adverse events,
7 let me just note that despite that enhanced volume
8 effect that I have described, we did not see any
9 evidence of thermal injury to tissue adjacent to the
10 uterine serosa, and this is one of the types of
11 adverse events that we always watch very closely in
12 devices that treat uterine pathology, so we did not
13 see any such adverse events.

14 What we did see was bleeding post
15 ExAblate, two first degree skin burns, a couple of
16 cases of nausea and vomiting, and some post-
17 hysterectomy adverse events that we would not be able
18 to argue were related to the treatment. They were
19 probably related to the hysterectomy.

20 As Dr. Stewart mentioned, there is also
21 feasibility data from outside of the United States.
22 And interestingly, in this study although 56 patients

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 received the ExAblate treatment, only four of those
2 patients elected to undergo hysterectomy, and that was
3 as of 14-month follow-up. So there is relatively less
4 hysterectomy data from this feasibility study
5 population.

6 The next couple of slides I'm not going to
7 spend much time on, but I would just like to say that
8 they demonstrate a trend, at least, towards non-
9 profuse volume being greater than the thermal dose
10 volume, although it was not uniform as it was in the
11 smaller feasibility study conducted at Brigham &
12 Women's and at St. Mary's in London.

13 In the feasibility studies that were
14 conducted in Israel, again this was not conducted
15 under FDA IDE regulation. However, there was one
16 adverse event that in hindsight we probably under-
17 appreciated at the time, and that was a case of
18 sciatica post treatment. This patient had symptoms as
19 of three weeks following her treatment, which at that
20 time were described as improving, and at that time she
21 was referred to a neurologist. I'm going to at least
22 allude to this adverse event later in my discussion.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com