

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

UNITED STATES OF AMERICA
FOOD AND DRUG ADMINISTRATION

+ + + + +

CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

+ + + + +

ORTHOPEDIC AND REHABILITATION DEVICES PANEL

+ + + + +

THURSDAY, SEPTEMBER 8, 2005

+ + + + +

The meeting was held in Salons A, B, and C
of the Hilton Washington, D.C. North, 620 Perry
Parkway, Gaithersburg, Maryland, at 8:34 a.m., Dr.
Sanjiv H. Naidu, Acting Panel Chairperson, presiding.

PRESENT:

SANJIV H. NAIDU, M.D. Ph.D., ACTING PANEL CHAIR

BRENT A. BLUMENSTEIN, Ph.D., DEPUTIZED VOTING MEMBER

CHOLL W. KIM, M.D., Ph.D., VOTING MEMBER

JAY D. MABREY, M.D., DEPUTIZED VOTING MEMBER

MICHAEL B. MAYOR, M.D., DEPUTIZED VOTING MEMBER

HARRY B. SKINNER, M.D., Ph.D., DEPUTIZED VOTING MEMBER

PAMELA ADAMS, INDUSTRY REPRESENTATIVE

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

PRESENT (Continued):

CONNIE F. WHITTINGTON, MSN, RN, CONSUMER
REPRESENTATIVE

JANET L. SCUDIERO, M.S., EXECUTIVE SECRETARY

MARK MELKERSON, M.S., FDA

NEAL R. GROSS
COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

C O N T E N T S

	<u>PAGE</u>
Appointment to Temporary Voting Status	5
Conflict of Interest Statement	6
Introductions	11
Post Market Study Design, Susan Gardner, Ph.D. ...	13
Division Update, Glenn Stiegman	19
Open Public Comment:	
Susan Krasny, Ph.D., R.A.C.	26
William Maloney, III, M.D.	32
<u>Birmingham Hip Resurfacing System, P040033</u>	
Sponsor Presentation, Marcos Velez-Duran	46, 66
Derek McMinn, FRCS	47
Tim Band	61
George DeMuth	77
Neal Defibaugh	86
FDA Presentation, John Goode, M.S.	117
Chang S. Lao, Ph.D.	151
Panel Deliberations and FDA Questions	168
Open Public Comment:	
Craig Thomas, M.D.	290

P R O C E E D I N G S

(8:34 a.m.)

MS. SCUDIERO: Good morning. We're ready to being this meeting of the Orthopedic and Rehabilitation Devices Panel.

I am Jan Scudiero, the Executive Secretary of this panel and a reviewer in the Division of General Restorative and Neurological Devices.

First, the usual housekeeping matters. If you haven't already done so, please sign the attendance sheets that are on the tables by the door.

Information on today's agenda and for panel meeting minutes and transcripts as well is there.

The next tentatively scheduled meeting for this panel on November 3rd and 4th is canceled because there's no agenda item ready for panel review.

Upcoming panel meetings are announced on an Advisory Panel Web site, the Federal Register and in the telephone information line. Please monitor the Web site for future meeting announcements.

Finally, as a courtesy to others in the

1 room, please, turn off or put on silence your cell
2 phones during the meeting.

3 Thank you.

4 Dr. John Kirkpatrick is unable to be with
5 us today.

6 I will now read into the record three
7 agency statements prepared for this meeting. They are
8 the appointment of temporary panel chair statement,
9 the appointment of temporary voting member statement,
10 and the conflict of interest statement.

11 First, I appoint Sanjiv H. Naidu, M.D.,
12 Ph.D., a voting member of the Orthopedic and
13 Rehabilitation Devices Panel as Acting Panel Chair for
14 the September 8th and 9th, 2005 meeting of the panel,
15 and this was signed by Dr. Daniel G. Schultz,
16 Director, Center for Devices and Radiological Health,
17 on September 7th.

18 The appointment to temporary voting
19 status: pursuant to the authority granted under the
20 Medical Devices Advisory Committee Charter dated
21 October 27th, 1990, and amended April 20th, 1995, I
22 appoint the following as voting members of the

1 Orthopedic and Rehabilitation Devices Panel for the
2 duration of this meeting on September 8th, 2005:
3 Brent A. Blumenstein, Ph.D., Jay D. Mabrey, M.D.,
4 Michael B. Mayor, M.D., and Harry B. Skinner, M.D.
5 Ph.D.

6 For the record, these people are special
7 government employees and are consultants to this panel
8 or another panel under the Medical Devices Advisory
9 Committee. They have undergone the customary conflict
10 of interest review and have reviewed the material to
11 be considered at this meeting.

12 The conflict of interest statement: the
13 Food and Drug Administration is convening today's
14 meeting of the Orthopedic and Rehabilitation Device
15 Panel of the Medical Devices Advisory Committee, under
16 the authority of the Federal Advisory Committee Act,
17 FACA, of 1972.

18 The panel meetings provide transparency
19 into the agency's deliberative processes. With the
20 exception of the industry rep. all members of the
21 committee are special government employees, or SGEs,
22 or regular federal employees from other agencies and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 are subject to federal conflict of interest laws and
2 regulations.

3 Consequently, in the interest of
4 transparency and in the spirit of disclosure, the
5 following information on the status of this panel's
6 compliance with federal ethics and conflict of
7 interest laws as covered by, but not limited to those
8 found at Part 18, U.S. Code 208, and Part 21, U.S.
9 Code, Section 355(n)(4). This is being provided to
10 the participants in today's meeting and to the public.

11 FDA has determined that the members of
12 this panel are in compliance with the federal ethics
13 and conflict of interest laws, including, but not
14 limited to, Part 18 U.S. Code Section 208, and Part
15 21, U.S. Code Section 355(n)(4). Under Part 18 U.S.
16 Code Section 208, applicable to all government
17 agencies and Part 21 U.S. Code Section 355(n)(4),
18 applicable to FDA, Congress has authorized FDA to
19 grant waivers to special government employees who have
20 limited financial conflict when it is determined that
21 the agencies need for the particular individual's
22 services outweighs his or her potential financial

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

17071 224-4433

www.nealrgross.com

1 conflict of interest.

2 Members who are special government
3 employees at today's meeting, including special
4 government employees appointed as temporary voting
5 members, have been screened for potential financial
6 conflicts of interest of their own, as well as those
7 imputed to them, including those of their employer,
8 spouse, or minor child, related to the discussions of
9 today's meeting. These interests may include
10 investments, consulting, expert witness testimony,
11 contracts, grants, CRADAs, teaching or speaking and
12 writing, patents and royalties, and primary
13 employment.

14 Today's agenda involves the review of a
15 pre-market approval application, PMA, for a hip
16 resurfacing system intended to relieve hip pain and
17 improve hip function in patients who have adequate
18 bone stock and are at risk of requiring more than one
19 hip joint replacement over their lifetimes.

20 This is a particular matters meeting
21 during which specific matters related to the PMA will
22 be discussed.

1 In accordance with Part 18, U.S. Code
2 Section 208(b)(3), a full waiver has been granted to
3 Dr. Michael Mayor. Dr. Mayor's waiver involves a
4 patent licensed to assist a company of a competitor.
5 He receives less than \$15,001 in royalties for the
6 patent.

7 Copies of each acknowledgement and consent
8 to disclosure statement signed by each participant at
9 today's meeting who received a conflict of interest
10 waiver along with the statement will be available for
11 review at the registration table during the meeting
12 and will be includes as part of official meeting
13 transcript.

14 A copy of the written conflict of interest
15 statement may be obtained by submitting a written
16 request to the agency's Freedom of Information Office,
17 Room 12A-30 of the Parklawn building.

18 Lastly, Ms. Pamela Adams is the industry
19 rep. acting on behalf of all related industry and is
20 employed by Etex Corporation, Incorporated. In the
21 event that the discussions involves any other products
22 or firms not already on the agenda for which an FDA

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 224-4433

www.nealrsgross.com

1 participant may have a financial interest, all meeting
2 participants are reminded that they're required by
3 Part 18, U.S. Code 208, to exclude themselves from
4 such deliberations and announce their exclusion for
5 the record.

6 Finally, in the interest of public
7 transparency with respect to all other participants,
8 we ask that they publicly disclose prior to making any
9 statements any current or previous financial
10 involvement with any firm whose product they may wish
11 to comment upon.

12 Thank you.

13 I would now like to turn the meeting over
14 to Dr. Naidu, our Acting Chairman for the day.

15 PANEL CHAIRPERSON NAIDU: Good morning. My
16 name is Sanjiv Naidu. I am the Acting Chairperson for
17 the Orthopedic and Rehabilitation Devices Panel today.
18 I'm Professor of Orthopedic Surgery at Penn State
19 College of Medicine, and I'm also a materials
20 scientist.

21 At this meeting the panel will be making a
22 recommendation to the FDA on the approvability of the

1 pre-market application, P040033, for the Smith &
2 Nephew Birmingham resurfacing hip system, the BHR
3 system. It is a hip joint, metal on metal, semi-
4 constrained, hybrid prosthesis, cemented femoral
5 component, uncemented acetabular component. It is
6 intended to relieve hip pain and improve hip function
7 in patients who have adequate bone stock and are at
8 risk of requiring more than one hip joint replacement
9 over their lifetimes.

10 Before we begin, I would like to ask our
11 distinguished panel members who are generously giving
12 their time to help FDA in the matter being discussed
13 today and the other FDA staff seated at this table to
14 introduce themselves. Please state your name, your
15 area of expertise, your position, and affiliation.

16 Mr. Melkerson, if you could start off.

17 MR. MELKERSON: I'm Mark Melkerson. I'm
18 the Acting Director for the Division of General
19 Restorative and Neurological Devices.

20 DR. MAYOR: Michael Mayor, orthopedic
21 surgeon, Professor of Orthopedics at the Dartmouth-
22 Hitchcock Medical Center and the Dartmouth Medical

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 School. I'm co-director of the Thayer Engineering
2 School, Dartmouth Biomedical Engineering Laboratories.

3 DR. BLUMENSTEIN: Brent Blumenstein, a
4 biostatistician working out of Seattle.

5 DR. MABREY: Jay Mabrey, orthopedic
6 surgeon, Medical Director of the Orthopedic Motion
7 and Sports Performance Laboratory at Baylor University
8 Medical Center and also Chief of Orthopedics at Baylor
9 University Medical Center in Dallas.

10 DR. KIM: I'm Choll Kim. I'm the
11 Assistant Professor of Orthopedic Surgery at the
12 University of California, San Diego. I'm the Director
13 of the Spine Research Lab there.

14 DR. SKINNER: My name is Harry Skinner.
15 I'm Professor and Chair of Orthopedic Surgery at the
16 University of California, Irvine.

17 MS. WHITTINGTON: My name is Connie
18 Whittington. I'm an orthopedic clinical nurse
19 specialist at Piedmont Hospital where I serve as the
20 Coordinator for Orthopedic Research.

21 MS. ADAMS: I'm Pamela Adams. I'm Chief
22 Operating Officer at Etex Corporation.

1 PANEL CHAIRPERSON NAIDU: Thank you.

2 I would like to note that for the record
3 that the voting members present constitute a quorum as
4 required by 21 CFR Part 14.

5 There will be two brief presentations
6 before the main agenda topic. First is Dr. Susan
7 Gardner who will speak on post market study design.

8 DR. GARDNER: Okay. Good morning. I'm
9 going to spend just a few minutes telling you about an
10 important programmatic change that has taken place in
11 the center.

12 The essence of the change is a move of the
13 condition of approval studies program from the Office
14 of Device Evaluation to the Office of Surveillance and
15 Biometrics.

16 Briefly, the Office of Surveillance and
17 Biometrics is involved both in premarket and post
18 market activities. We're involved in the premarket
19 review because the statisticians and the
20 epidemiologists of the center are in the Office of
21 Surveillance and Biometrics, and we also have a major
22 role in the post market in that all the adverse events

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

1 and the post market monitoring tools, the medical
2 device reporting system, and the MedSun program are in
3 OSB.

4 We're also responsible for analyzing the
5 data to come in on these post market monitoring tools
6 and characterizing the risk and working with the rest
7 of the center to identify post market problems and
8 take action on those. We're responsible for
9 coordinating the center response to health care
10 professionals in risk communication, and we're
11 responsible for interpretation of the medical device
12 reporting regulation.

13 The legislative 21 CFR 814.82 says that
14 post approval requirements can include continuing
15 evaluation and periodic reporting on the safety,
16 effectiveness, and reliability of the device for its
17 intended use, and this is the basis of the condition
18 of approval studies program.

19 The impetus for the change came from a
20 study, an internal evaluation that we did in CDRH of
21 our CoA Program. We look at this in about 2000, 2001,
22 and we went back and we looked at all of the PMAs that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE. N.W.
WASHINGTON D.C. 20005-1701

(202) 234-4433

www.nealrgross.com

1 had been approved from 1998 through the year 2000.
2 There were 127 PMAs, and 45 of those had condition of
3 approval orders.

4 And what we found unfortunately is that we
5 were really unable to find some of these studies, and
6 we realized this was because we had no standardized
7 tracking procedures for tracking the results of these
8 studies. As one might expect there had been a
9 turnover in lead reviewers as people naturally moved
10 through the organization or changed jobs, and also
11 that in ODE with their focus on premarket, it was
12 really difficult for them under their current
13 resources to continue to follow these studies and give
14 them the attention that they need.

15 So we came up with a plan to change the
16 program. The point of the change or the goal of the
17 change was to make sure that we could obtain this post
18 market information at this critical period when the
19 device enters the market and we could continue to
20 assess the safety and effectiveness as the device
21 moves from the clinical trial into the real world use,
22 and this obviously would allow us to better

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 characterize the risk-benefit profile and add to our
2 ability to make sound scientific decisions as we
3 continue to monitor these devices.

4 We actually did a pilot that went for
5 about two years. So when we officially changed the
6 program on January 1st of 2005, we were fairly well
7 prepared, and in addition to that a number of studies
8 already were being monitored in OSB. We have
9 developed and instituted an automatic tracking system
10 that is up and working, and not only do we have
11 studies being tracked from January 2005. We've gone
12 back to pick up studies that have been approved
13 earlier. So we're also following those.

14 The point of the tracking system obviously
15 is to acknowledge to industry when studies are
16 received and to follow up if we don't receive the
17 information that we're supposed to have.

18 Another fundamental change is we have
19 added an epidemiologist to the PMA review team when we
20 think that we're going to have or it looks likely that
21 we'll have a condition of approval study if the device
22 is approved. The epidemiologist is tasked with the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 development of the post market monitoring plan during
2 the premarket review working with the rest of the PMA
3 team. The epidemiologist has the lead for developing,
4 and let me emphasize we're really looking to develop
5 well formulated post market questions and make sure
6 that these studies are important if we're going to ask
7 industry to do them.

8 The epi person will have the lead in the
9 design of the study protocol and in the evaluation of
10 study progress and results after approval, and so when
11 the results of these studies come into OSB, we will
12 then look at the results and then turn to our
13 colleagues in the premarket arena and go back to the
14 PMA team, review what the progress has been and see if
15 there's anything else that we need to do.

16 So with these changes, why do we think we
17 will do better? Well, first of all, we think that
18 with a real emphasis on working with industry to make
19 sure that the questions that we're asking in the
20 premarket arena are really important fundamental
21 questions, and we have a good study protocol design.
22 Everybody will be motivated, first of all, to get the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

1001 234-4433

www.nealrgross.com

1 studies done, and secondly, to evaluate them and make
2 sure that the results are important.

3 For everybody acknowledgement of the work
4 that you do and feedback on the studies motivate
5 people to do it. So, again, industry won't feel that
6 their work is falling into a black hole, and we will
7 be able to interact with them to make sure that the
8 results are on target.

9 We will be posting the status of the
10 studies on CDRH Web site, and also if studies are not
11 done, we do have the ability under Section 522 to
12 mandate a post market study, but again, we are hoping
13 if people work together to have a good study and a
14 good protocol that we won't have to go there.

15 What does this mean to the Advisory Panel?

16 Well, first of all, let me emphasize it certainly
17 doesn't change the standards for making sure that the
18 product is reasonably safe and effective before it's
19 approved.

20 However, during the approval process we
21 will sometimes lay out deliberately and sometimes, of
22 course, it will come up naturally in discussion

1 questions that deal with the post market piece of the
2 approval process, and your input, your suggestions of
3 possible approaches for the post market pieces or
4 questions that you're concerned about will be really
5 important to us and certainly will be taken under
6 consideration if the device is approved and we decide
7 to do a condition of approval study.

8 We are also committed to coming back to
9 you; either FDA or industry, to come back to you and
10 give you the results of condition of approval studies
11 after the product is approved.

12 Any questions?

13 (No response.)

14 DR. GARDNER: Thanks.

15 PANEL CHAIRPERSON NAIDU: Thank you, Dr.
16 Gardner.

17 Now, Mr. Glenn Stiegman will give us a
18 division update on the activities since 2004 panel
19 meeting.

20 MR. STIEGMAN: Hi. My name is Glenn
21 Stiegman. This is the panel update for the Orthopedic
22 and Rehabilitation Panel of September 8th and 9th.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 Action on items from previous Advisory
2 Committee meetings. From the June 2nd, 2004 meeting,
3 P040006 from the DePuy Spine, the Charite artificial
4 disc. This device was approved October 26th of 2004
5 for spinal arthroplasty in patients with degenerative
6 disk disease at one level from L4 to S1.

7 A post approval study to further document
8 the incidence of complications, such as migration and
9 subsidence is ongoing.

10 Other significant approvals that have
11 occurred since the last panel update, two ceramic on
12 ceramic hips, one on December 17th, 2004, the Smith &
13 Nephew Reflections ceramic acetabular system, which is
14 indicated for patient requiring primary total hip
15 arthroplasty due to noninflammatory arthritis, such as
16 osteoarthritis, avascular necrosis and traumatic
17 arthritis.

18 On May 4, 2005, the DePuy Orthopedics, the
19 Dura Option ceramic hip system for patients with
20 noninflammatory degenerative joint disease, such as
21 osteoarthritis, avascular necrosis, congenital hip
22 dysplasia, and post traumatic arthritis.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 734-4433

www.nealrgross.com

1 On December 3rd, 2004, P010029, Savient
2 Pharmaceuticals, which has recently been bought out by
3 Ferring Pharmaceuticals, the Nuflexxa one percent
4 sodium hyaluronate for treatment of pain and
5 osteoarthritis of the knee in patients who have failed
6 to respond adequately to conservative nonpharmacologic
7 therapy and simple analgesics. This particular device
8 is made by bacterial fermentation instead of being
9 extracted from chicken products like other hyaluronate
10 products.

11 An HDE, the humanitarian device exemption,
12 H030009 from Synthes USA, the vertical expandable
13 prosthetic titanium rib, which is indicated for the
14 treatment of thoracic insufficiency syndrome. It's
15 going to immature patients. TIS is defined as
16 inability of the thorax to support normal respiration
17 or lung growth.

18 Other significant 510(k) clearances,
19 February 18th, 2005, Blackstone Medical, we cleared a
20 laminoplasty fixation system for holding bone graft in
21 place during laminoplasty procedures.

22 From the Zimmer Trabecular Technology, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 trabecular metal osteonecrosis interventional implant
2 which is indicated to treat patients with
3 osteonecrosis of the femoral head.

4 Guidances that have been published since
5 the last panel update, the clinical data presentations
6 for orthopedic device applications and the clinical
7 trial considerations, vertebral augmentation devices
8 to treat spinal insufficiency fractures.

9 The Division of General, Restorative, and
10 Neurological Devices staffing changes. The new or
11 permanent orthopedic staff in ORDB, Ronald Jean, Dr.
12 Kristin Mills, John Holden, and Dr. Khan Li have all
13 joined our staff since the last panel update on a
14 permanent or new basis.

15 Also, no longer with the Orthopedic
16 Devices Branch or FDA, Barbara Zimmerman is now the
17 Deputy Director of the Division of Cardiovascular
18 Devices. Dr. Celia Witten is now the Director of the
19 Office of Cellular Tissue and Gene Therapies at CBER.
20 Dr. Martin Yahiro is moving on to industry on
21 September 17th. Genevieve Hill is returning to school
22 and Dr. Michael Schlosser is returning to private

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 practice.

2 Current division staff changes, Mark
3 Melkerson is the Acting Division Director on a four-
4 month detail. Ted Stevens and Barb Buch are both
5 Acting Deputy Division Directors for DGRND. Myself,
6 Glenn Stiegman, is the new Branch Chief of Orthopedic
7 Devices Branch, and Aric Kaiser is currently on a
8 four-month detail on the Office of Surveillance and
9 Biometrics.

10 And lastly, in conclusion, we'd like to
11 take a moment of silence to remember one of our
12 bright, young, highly esteemed orthopedic device
13 reviewers, Jonathan Lim who passed away last week.
14 More than just a colleague, Jon was a good friend to
15 many and will be sorely missed. Jon's spirit will
16 always live in the pleasant memories that he created
17 in our hearts.

18 (Pause in proceedings.)

19 MR. STIEGMAN: Thank you.

20 PANEL CHAIRPERSON NAIDU: Thank you, Mr.
21 Stiegman.

22 We will now proceed with the open public

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 hearing portion of the meeting. Prior to the meeting,
2 two people asked to speak in the open public hearing.

3 They will speak in order of their request to speak.

4 We ask that you speak clearly into the
5 microphone as the transcriptionist is dependent on
6 this means of providing an accurate record of this
7 meeting.

8 Please state your name and the nature of
9 any financial interest you may have in this or any
10 other medical device company. Ms. Scudiero will now
11 read the open public hearing statement.

12 MS. SCUDIERO: Both the FDA and the public
13 believe in a transparent process for information
14 gathering and decision making. To insure such
15 transparency at the open public hearing session of the
16 Advisory Committee meeting, FDA believes it is
17 important to understand the context of any
18 individual's presentation.

19 For this reason, FDA encourages you, the
20 open public hearing or industry speaker, at the
21 beginning of your written or oral statement to advise
22 the committee of any financial relationship you may

1 have with the sponsor, its products and, if known, its
2 direct competitors.

3 For example, this information may include
4 the sponsor's payment of your travel, lodging or other
5 expenses in connection with your attendance at the
6 meeting.

7 Likewise, FDA encourages you at the
8 beginning of your statement to advise the committee if
9 you do not have any such financial relationships. If
10 you choose not to address the issue of financial
11 relationships at the beginning of your statement, it
12 will not preclude you from speaking.

13 I would like to note for the record that
14 prior to the meeting three letters with general
15 comments on the agenda item were received. Copies of
16 these letters were given to the panel and the PMA
17 sponsor this morning. These are now part of the
18 record for this meeting.

19 Dr. Naidu.

20 PANEL CHAIRPERSON NAIDU: Thank you.

21 The first open public hearing presented is
22 Dr. Susan Krasny, Vice President of the Orthopedic

1 Surgical Manufacturers Association.

2 Dr. Krasny, you have ten minutes.

3 DR. KRASNY: Thank you.

4 I am the Senior Director of Regulatory
5 Affairs and Clinical Affairs with Stryker Spine, and I
6 have no financial interest in the sponsor today.

7 Good morning. My name is Susan Krasny,
8 and I speak here today representing the Orthopedic
9 Surgical Manufacturers Association, OSMA, of which I
10 am Vice President.

11 OSMA, a trade association of over 30
12 member companies, welcomes this opportunity to provide
13 general comments at today's Orthopedic Advisory Panel
14 meeting. OSMA's comments should not be taken as an
15 endorsement of the products being discussed today. We
16 ask instead that our comments be considered during
17 today's panel deliberations. These comments represent
18 the careful compilation of member companies' views.

19 OSMA was formed over 40 years ago and has
20 worked cooperatively with the FDA, the American
21 Academy of Orthopedic Surgeons, the American Society
22 for Testing Materials and other professional medical

1 societies and standard development bodies. This
2 collaboration has helped to insure that orthopedic
3 medical products are safe, of uniform high quality,
4 and supplied in quantities sufficient to meet national
5 needs.

6 Association membership currently includes
7 over 30 companies who produce over 85 percent of all
8 orthopedic implants intended for clinical use in the
9 United States.

10 OSMA has a strong and vested interest in
11 insuring the ongoing availability of safe and
12 effective medical devices. The deliberations of the
13 panel today and the panel's recommendation to the FDA
14 will have a direct bearing on the availability of new
15 products.

16 We make these comments to remind the panel
17 of the regulatory burden that must be met today. We
18 urge the panel to focus its deliberations on the
19 product safety and effectiveness based on the data
20 provided.

21 The FDA is responsible for protecting the
22 American public from drugs, devices, food and

1 cosmetics that are either adulterated or unsafe or
2 ineffective. However, FDA has another role, that to
3 foster innovation.

4 The Orthopedic Devices Branch is fortunate
5 to have available a staff of qualified reviewers,
6 including Board certified orthopedic surgeon to
7 evaluate the types of applications brought before this
8 panel. The role of this panel is also very important
9 to the analysis of the data in the manufacturer's
10 application and to determine the availability of new
11 and innovative products in the U.S. marketplace.

12 Those of you on the panel have been
13 selected based on your expertise and training. You
14 also bring the view of practicing clinicians who treat
15 patients with commercially available products.

16 OSMA is aware that you have received
17 training from the FDA on the law and regulation, and
18 we do not intend to repeat that information today. We
19 do, however, want to emphasize two points that may
20 have a bearing on today's deliberations:

21 One, reasonable assurance of safety and
22 effectiveness, and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

1/20/21 2:44:33

www.nealrgross.com

1 Two, valid scientific evidence.

2 There is a reasonable assurance that a
3 device is safe when it can be determined that the
4 probable benefits outweigh the probable risks. Some
5 important caveats associated with this oversimplified
6 statement include valid scientific evidence and proper
7 labeling, and that safety data may be generated in the
8 laboratory, in animals and in humans.

9 There is reasonable assurance that a
10 device is effective when it provides a clinically
11 significant result. Again, labeling and valid
12 scientific evidence play important roles in this
13 determination.

14 The regulation and the law clearly state
15 that the standard to be met is reasonable assurance of
16 safety and effectiveness. "Reasonable" is defined as
17 moderate, fair, and inexpensive.

18 The regulation states that well controlled
19 investigations shall be the principal means to
20 generate the data used in the effectiveness
21 determination.

22 The following principles are cited in the

1 regulation as being recognized by the scientific
2 community as essential in a well controlled
3 investigation.

4 One, a study protocol;

5 Two, method of selecting patients;

6 Method of observation and recording
7 results;

8 And, four, comparison of results with
9 control.

10 The panel has an important job today. You
11 must listen to the data presented by the sponsor,
12 evaluate the FDA presentations, and make a
13 recommendation about the approvability of the
14 sponsor's application.

15 We speak for many applicants when we ask
16 for your careful consideration. Please keep in mind
17 that the standard is reasonable assurance, balancing
18 the benefits with the risks. The regulatory standard
19 is not proof beyond a shadow of a doubt.

20 When considering making recommendations
21 for further studies, remember that the FDA takes these
22 recommendations seriously, often as a consensus of the

1 panel as a whole and they may delay the introduction
2 of a useful product or result in burdensome and
3 expensive additional data collection.

4 Therefore, you play an important role in
5 reducing the burden of bringing new products that you
6 and your colleagues use in treating patients to the
7 market. Please be thoughtful in weighing the
8 evidence. Remember that the standard is a reasonable
9 assurance of safety and effectiveness and that there
10 is a legally broad range of valid scientific evidence
11 to support the determination.

12 OSMA thanks the FDA and the panel for the
13 opportunity to speak today. Our association trusts
14 that its comments are taken in the spirit offered, to
15 help the FDA decide whether to make a new product
16 available for use in the U.S. marketplace.

17 OSMA are present in the audience and are
18 available to answer to answer questions any time
19 during the deliberations today.

20 Thank you.

21 PANEL CHAIRPERSON NAIDU: Thank you, Dr.
22 Krasny.

1 The next presenter is Dr. William Maloney,
2 III.

3 Dr. Maloney.

4 DR. MALONEY: Thank you, and good morning.

5 My name is Bill Maloney, and I'm the
6 professor and Chairman of Orthopedics at Stanford
7 University School of medicine.

8 And while we're getting up my
9 presentation, I'll do my conflict disclosure. I
10 design hip implants for Zimmer, for which I receive
11 royalties. I am not involved in any surface
12 replacement design. I design knee implants for Wright
13 Medical for which I receive royalties. I do not own
14 stock in any orthopedic implant companies. My travel
15 expenses have been reimbursed by Wright Medical, and
16 currently my time is not being compensated to come
17 here.

18 So I'm here to make some comments on this
19 submission, and for those of you on the panel who do
20 not know me, in addition to being professor and
21 Chairman of Stanford University, I'm currently Chief
22 of Joint Replacement Service at Stanford.

1 Prior to this I was the Chief of the Joint
2 Replacement Service at Washington University School of
3 Medicine. I'm a member of the North American Hip
4 Society and International Hip Society. I led the
5 American Joint Replacement Registry effort in trying
6 to get that off the ground in this country and chaired
7 that committee for two years, and I'm currently a
8 member of the Quality Improvement Program for Medicare
9 working with Dr. David Hunt.

10 I want to say quite clearly what I'm not
11 here to do. I'm not here to indict Smith Nephew.
12 It's a well established orthopedic company, a great
13 reputation. I'm not here to indict Total Resurfacing
14 Arthroplasty or the specific implant, and I'm
15 certainly not here to indict Derek McMinn, who is a
16 colleague of mine and certainly a well known
17 arthroplastic surgeon.

18 What I am here to do is to make some
19 comments on the study methodology which I think
20 appears to me as a problem, to talk a little bit about
21 conflict of interest. Obviously that's an important
22 thing in this country and one as a chairman of an

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

1/20/21 234.4433

www.nealrgross.com

1 academic department at a major medical school I deal
2 with on a regular basis, and make a comment on
3 compelling medical need.

4 When you look at studies, there are a
5 variety of studies in the orthopedic literature. Most
6 of them are retrospective. There's prospective data
7 collection with retrospective data review, like
8 registry as the Mayo Clinic Registry would be a good
9 example of that.

10 There's a prospective study. In a
11 prospective study you generate a hypothesis. You
12 design the study to test the hypothesis. You get IRB
13 approval. You recruit patients, and you consent the
14 patients. You collect the data, and you analyze the
15 data.

16 Then you have an IDE study, which is a
17 prospective study that goes beyond those requirements.

18 There's site monitoring. There's concurring data
19 verification against source documentation. There's
20 adverse event documentation and reporting, and there's
21 FDA inspections, akin to an IRS audit.

22 When you look at the norms for an IDE

1 study protocol, these have been well established. The
2 protocol is predetermined. The safety and efficacy
3 endpoints are well defined. The control group is
4 homogeneous with the treatment group. There's well
5 established inclusion and exclusion criteria. There's
6 predefined follow-up intervals and direct patient
7 evaluation. The conclusions are based on a
8 predetermined, valid statistical plan, and there's
9 substantial patient site monitoring which is
10 mandatory.

11 These trials are usually multi-center
12 under a common protocol. There's contemporaneous
13 accountability for all adverse events and complete
14 documentation of all protocol deviations, regardless
15 of severity.

16 The data that we're currently or you're
17 currently reviewing in this current PMA submission is
18 a combination of retrospective data and data collected
19 prospectively and reviewed retrospectively. In that
20 submission it states there were not predefined follow-
21 up time windows, standardized clinical evaluations,
22 adverse event report forms or standardized

1 radiographic evaluations.

2 So clearly, that doesn't meet the standard
3 that we're used to in the evaluation of Class 3
4 devices.

5 When you look at critical study elements,
6 the survivorship data is from a single surgeon. It's
7 from Derek McMinn's data, and the justification of
8 that is as follows. The Oswestry Outcomes Center
9 records information regarding complications, deaths or
10 revisions, but the data is not verifiable. You cannot
11 verify the primary source data. Therefore, the only
12 verifiable data is Dr. McMinn's data.

13 They further go on to state, "Comparison
14 of data from the Oswestry Outcomes Center database to
15 the data from Mr. McMinn's clinical records provides
16 additional assurance that an accurate, up-to-date
17 survivorship data is known."

18 But the problem here is "additional." The
19 only data we have here on survivorship that can be
20 verified is Dr. McMinn's data, and it's not done
21 concurrently and not by an independent study monitor.

22 The single surgeon series, as we all know,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 has significant disadvantages. Derek is an expert
2 surgeon. He's been doing this operation since 1989.
3 He is the pioneer in this operation and deserves a lot
4 of credit, but the question is: were his results be
5 optimal to the surgeon community at the U.S. at large?

6 Clearly, I don't think that's the case.
7 This does not represent U.S. surgeon data or U.S.
8 patients, and this is a new operation in the United
9 States. We have no experience with it, and surgical
10 training is going to be critical.

11 When we look at the specific cohorts in
12 this study, there are three. There's an X-ray cohort.

13 There's the Oswestry cohort, and there's the McMinn
14 cohort.

15 The Oswestry Outcome Center, which it's
16 unclear how that's funded, is a self-administered
17 patient questionnaire. There's no direct patient
18 contact, which is, of course, required in ID studies.

19 They report an outcome measure which we're not
20 familiar with, an OSHIP score which appears to be a
21 combination of a Harris HIP score and a Merle
22 D'Aubigne score, and presumably there's patient

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(707) 234-4433

www.nealrgross.com

1 consent in this study, although it's not well
2 documented.

3 In the McMinn cohort, it appears to have
4 no functional outcome data. It is specifically just
5 survivor data, as the outcome data was not available
6 for review of this submission, and it doesn't appear
7 that those patients in the McMinn cohort were actually
8 consented for a research study.

9 They certainly were consented for surgery.
10 The PMA submission clearly states that they were
11 given options as it relates to hip arthroplasty, but
12 that's fundamentally different than being consented
13 for a research study.

14 What about the radiographic analysis?
15 This, again, appears to be primarily a retrospective
16 radiographic analysis. Clearly at the beginning of
17 this PMA station they stated that there's no
18 radiographic protocol.

19 If you look at the numbers in this cohort,
20 there's 124 hips. Two hips were excluded because of
21 patient death in one patient. Four hips were revised.
22 So they left with 118. Ten of those we then lost to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 follow-up. So you're down to 108, and then 19 hips
2 had no immediately post-op X-ray. So there's 89 hips
3 which were potentially valuable.

4 It goes on to state there were immediate
5 post op films on 89 of 108 procedures with five
6 radiographs, but the sponsor stated that these films
7 were of low quality, portable films and unusable for
8 the purposes of precise postoperative comparisons.

9 Therefore, baseline films for the purpose
10 of comparisons were made in each of the 108 cases in
11 the postoperative time period, usually within three
12 months, but eight of the 108 procedures had baseline
13 evaluations performed at the time points ranging from
14 110 to 860 days.

15 So in some cases the baseline X-rays were
16 quite a bit of time after the index procedure. In a
17 prospective IDE study, these would all be protocol
18 violations and could be potentially excluded. So you
19 could actually end up with a cohort in the X-ray study
20 of zero if you had a strict analysis of the data.

21 The next topic I want to briefly discuss,
22 and it has already been a significant topic here this

1 morning is one of conflict of interest. In this
2 country, we are under great scrutiny as it relates to
3 our conflict, and as a department chairman, again, I
4 deal with this on a regular basis.

5 Conflict of interest exists when an
6 investigator and/or his or her family has the
7 potential for financial gain. Financial conflicts in
8 human subject research at least at our institution are
9 of special concern, and I'm sure they are across this
10 country.

11 How, Derek is the product designer. He's
12 the pioneer here, and he's one of the co-founders of
13 Midland Medical Technology, who we do not have
14 disclosure in the PMA submission as to what the exact
15 conflict is here. This company was purchased by Smith
16 and Nephew for 67 million pounds, and that's a fair
17 sum of money in anybody hands, and it will pay another
18 33 million pounds for the Food and Drug Administration
19 or if the Food and Drug Administration allows this
20 procedure into the United States.

21 The company's history, this has been an
22 extremely successful company, and I wish I was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 involved with it. It was incorporated in 1996, and in
2 1997 they issued about 100,000 shares, 60,000 of which
3 were immediately transferred to Sky Fall, Limited and
4 30,000 to Maldeney, Limited. Who owns those we have
5 no idea, and who has gotten the 66 million pounds
6 needs to be addressed because we have to do that in
7 this country, and the playing field needs to be level.

8 Now, Derek has got a five-year contract
9 based on the Daily Mail with his new employers. He
10 told the Daily Mail, "I will stay on until I retire
11 and drop. I am completely addicted to this product."

12 Now, if that's the case, if he's an
13 employer of the sponsor, and again, we need to have
14 clarification there, it appears that the data in this
15 PMA is entirely from the clinical practice of a Smith
16 & Nephew employee.

17 In this country and for the FDA for the
18 IDE studies, you rely primarily on the individual's
19 IRB to oversee issues of patient protection. This
20 includes the study consent form and disclosure of
21 financial conflict of interest.

22 At my institution I would not be able to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 do an IDE study that involved a Zimmer product. I
2 would not be able to do a clinical study that involved
3 a Zimmer hip product. It's simply not mitigatable.

4 The academic community in this country has
5 weighed in on this issue. The NIH has weighed in on
6 this issue, and the federal government has weighed in
7 on this issue. Significant financial conflicts cannot
8 be mitigated. It doesn't mean that the data is bad,
9 and it doesn't mean that the investigator is
10 dishonest. It just means that the appearance of
11 significant financial conflicts currently are not felt
12 to be mitigatable.

13 The last point relates to compelling
14 medical need. Now, some people have said, well, these
15 products are great, we need to have them on the
16 market, but clearly with resurfacing arthroplasty
17 that's not the case. Conventional hip replacement is
18 a good operation. It has an extremely long track
19 record, and currently in this country there are three
20 IDEs being performed and two PMAs that are pending
21 before the FDA.

22 Corin and DePuy and Wright Medical all

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrross.com

1 have IDE studies in this country. So clearly, we're
2 going to have a little bit higher quality data to
3 review in the relatively near future.

4 So how do I summarize my evaluation and
5 those of us who are involved in the American Joint
6 Replacement Registry?

7 This is a PMA device. It's the first of
8 its kind in the United States. It's supported by
9 data, and this is a quote from the PMA submission,
10 from essentially one surgeon, Derek McMinn, who is a
11 very qualified surgeon and very expert surgeon.

12 That source of data is also the product
13 designer and the founder of the company and the things
14 they implant. IDE conventions are not followed. The
15 data set does not reflect the high bar set by the FDA
16 for approval of Class III orthopedic devices, and I
17 remind you that the last time this type of data was
18 accepted at least that I could find was back when the
19 Mittelmeier hip was approved. This was a ceramic-on-
20 ceramic total hip replacement that had excellent
21 clinical results in Europe, principally by Dr.
22 Mittelmeier. There was no U.S. IDE. It was a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

7071 234-4433

www.nealrgross.com

1 clinical disaster in this country, and it doomed
2 ceramic-to-ceramic hip replacements in this country
3 for years.

4 This would be precedent setting in terms
5 of how all Class III devices are evaluated for
6 approval in this country, and I'll tell you I was
7 involved, and I was in the orthopedic implant side.
8 If this goes forward, this would be the last IDE done
9 in this country. There's no reason to do an IDE if
10 retrospective data from another country is acceptable.

11 And maybe that's the way we should go, but
12 clearly, it's a significant change in our current
13 standards. Maybe the bar is too high, and maybe it is
14 too burdensome, and it probably is, and I probably
15 need to have some changes made, but this is a big jump
16 from what we're currently used to, and it potentially
17 is a significant patient safety issue.

18 Thank you for your attention.

19 PANEL CHAIRPERSON NAIDU: Thank you, Dr.
20 Maloney.

21 Is there anyone else in the room who would
22 like to address the panel at this point? If so,

1 please raise your hand and come forward and state your
2 name, affiliation and whether you have any involvement
3 in a medical device firm.

4 (No response.)

5 PANEL CHAIRPERSON NAIDU: I don't see
6 anyone. Please note that there will be a second open
7 public session in the afternoon. If anyone else would
8 like to address the panel about today's agenda topic,
9 you may speak in the afternoon.

10 We will now proceed to the sponsor
11 presentation for the Birmingham hip resurfacing
12 device. We will then have a short break and proceed
13 with the FDA presentation.

14 After lunch the panel will deliberate on
15 the approvability of the PMA. Before the panel votes
16 on the approvability of the PMA, there will be a
17 second open public hearing and FDA and sponsor
18 summations.

19 I would like to remind public observers at
20 this meeting that while this meeting is open for
21 public observation, public attendees may not
22 participate except at the specific request of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 224-4433

www.nealrsgross.com

1 panel.

2 We will begin with the sponsored
3 presentation. The first Smith & Nephew presenter is
4 Mr. Marcos Velez-Duran, Vice President of Clinical and
5 Regulatory Affairs and Quality Assurance. He will
6 introduce the other Smith and Nephew presenters.

7 Mr. Duran.

8 MR. VELEZ-DURAN: Yes. Good morning. My
9 name is Marcos Velez-Duran, and I am the VP of
10 Regulatory and Clinical Affairs and Quality for Smith
11 & Nephew. I'm here today to present to you a summary
12 of the data presented in the pre-market approval for
13 the Birmingham hip resurfacing product.

14 I will take this opportunity to thank the
15 panel and FDA for the time and effort that they have
16 put in into the review of this PMA. In particular, I
17 would like to thank FDA for embracing the spirit of
18 the regulation. As demonstrated by the then
19 interactively review in this PMA and concerning the
20 data that will be presented today as meeting the
21 requirement of valid scientific evidence as outlined
22 under the law.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 The presenters for this morning as me,
2 Marcos Velez; Mr. Derek McMinn will present design
3 history; Mr. Tim Band will present on preclinical
4 information. I will come back and present the summary
5 of the clinical data. George DeMuth will present on
6 the statistical issues associate with this PMA and
7 data analysis, and Neal Defibaugh will review the
8 labeling, plus approval studies and plan for training.

9 In addition, we have a number of people
10 supporting our efforts today that will get up and
11 answer questions as necessary, that is, Dr. Cecil
12 Rorabeck, Professor James Richard from the Oswestry
13 Center, Mr. Joseph Daniels, Professor Anthony
14 Unsworth, Dr. Roger Rogerson, Marie Marlow, Dr. Marc
15 Thomas, and Sally Maher.

16 Now I would like to introduce Mr. McMinn
17 to give the next presentation.

18 DR. MCMINN: Thank you, Marcos.

19 Good morning, ladies and gentlemen. I'm
20 Derek McMinn. I'm an orthopedic surgeon from
21 Birmingham, England, and I'm the co-inventor of the
22 Birmingham hip resurfacing. I do have a financial

1 interest in the Birmingham hip resurfacing. I'm a
2 consultant for Smith & Nephew, and I'm a non-Executive
3 Director of Smith & Nephew.

4 For the nonclinicians, I just want to do a
5 very quick review of the indications for hip
6 arthroplasty. On the left there you see a normal hip,
7 and on the right you can see an osteoarthritic hip
8 with loss of some articular cartilage.

9 Next.

10 When you look at an arthritic femoral
11 head, there is the retained cartilage, and there is
12 the variable on the femoral head. So the difference
13 between a pain free hip and a very painful hip with a
14 disabled patient is the loss of a few millimeters of
15 articular cartilage.

16 Now, for that lesion, we take a really
17 aggressive approach to it. We resect the femoral head
18 and neck, insert a cup in the stem to the shaft of the
19 femur called a total hip replacement. Why do we do
20 that? Because we've all been trained to do it and
21 because for the vast majority of patients it works
22 extremely well.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrimes.com

1 Next.

2 Here's a patient of mine, a total hip
3 replacement freshly done on the left. A few years
4 later unfortunately the patient has developed
5 loosening of the femoral component with substitutes.
6 It had to be revised.

7 Next.

8 This is another patient of mine with two
9 hip replacements and both sides have failed for
10 different reasons. On this side, a typical total hip
11 replacement. Where is the polyethylene, and that is
12 caused a linear pattern of osteolysis and loosening of
13 the socket.

14 On the other side we've got a hybrid total
15 hip replacement and an uncemented cup, and again,
16 there's been wear of the polyethylene, and this
17 patient has nasty pelvic osteolysis, and on the femur
18 side it's removed most of Zone 7 in the femur.

19 The problem that we all face now is
20 largely that set-up with pelvic osteolysis thanks to
21 wear of polyethylene, and the problem is made worse
22 the younger and more active the patient is.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrross.com

1 Next.

2 So to summarize the problems with total
3 hip replacement, there's an incidence of dislocation.

4 In our country it's three to four percent. There's a
5 leg lengthening issue with total hip replacement,
6 which is a common cause for patient unhappiness.
7 There's acetabular and femoral loosening which is
8 reducing with time, but particularly acetabular and
9 femoral osteolysis is a big problem. There's the
10 inevitable stress shielding of inserting a stem into
11 the shaft of the femur and a revision surgery, as we
12 all know, is difficult and expensive and in young
13 patients may be recurrent.

14 So why can't we just replace the worn out
15 acetabular and femoral articular cartilage? The
16 answer is we can, and that's a typical model for a hip
17 resurfacing where we put a thin shell over the ball
18 and a thin shell into the socket to replace the worn
19 out articular cartilage.

20 Is that new? Not at all. This was Sir
21 John Charnley's first attempt at hip arthroplasty, and
22 he used Teflon shells both on the socket and the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 femoral head, and here's some removed two years with
2 complete wear through of the socket and the head.

3 Then came the 1970s generation of his
4 resurfacings using the materials available at that
5 time, polyethylene on the socket, metal shell on the
6 femoral head.

7 Here's a case of ours, loose cup, a common
8 complication of these cemented resurfacing models.

9 Here's another one I revised, loose
10 socket, femoral component solid. When we sliced that
11 component there are holes in the femoral head. Why
12 are they there?

13 Here you can see an entry hole at the
14 head-neck junction.

15 Next.

16 Here you can see a granuloma on histology
17 and on polarized light microscopy this head is stuffed
18 full of polyethylene debris. So the polyethylene
19 debris problem has loosened the cup, and it's starting
20 to erode the femoral head.

21 Is hip resurfacing new in the United
22 States? Not at all. There's some extremely good

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

1/20/21 234 4433

www.nealrgross.com

1 innovators in the United States. The first two models
2 of hip resurfacing are cemented, and the problems are
3 largely as I outlined already in the previous models.

4 The second two models are uncemented, and
5 they did, indeed, address the issue satisfactorily of
6 component loosening. But of course, as you can see
7 from those sample radiographs, they did not address
8 the problem of tremendous osteolysis because of
9 polyethylene debris generated as an inevitable result
10 of the large head articulating on polyethylene.

11 So to summarize the 1970s and '80s
12 resurfacing, there's a small instance of femoral neck
13 fracture and collapsed femoral heads. That was
14 manageable. It wouldn't have finished resurfacing in
15 the '70s and '80s. The thing that finished it was the
16 large head articulating on polyethylene gave massive
17 debris of polyethylene, and that loosened cemented
18 cups and cause osteolysis with cementless components.

19 So the conclusion of that was that hip
20 resurfacing was not a viable operation with varying
21 materials available at the time. Our experience in
22 Birmingham with the '70s and '80s resurfacings

1 unfortunately was no better than anyone else's, and
2 this is a survivorship curve with the cases we did,
3 and you can see that we've got a 50 percent revision
4 rate at six to seven years.

5 As a trainee I put some of these in, and
6 then as an orthopedic surgeon in the late '80s, I got
7 to revise most of these cases. So I was seeing
8 patients turning up in my clinics with resurfacings
9 that had failed because of polyethylene debris
10 generated from large femoral heads.

11 Also in my clinics, I was following up my
12 predecessor's cases of metal-metal total hip
13 replacements, and the first total hip replacement of a
14 metal-metal nature in our hospital was done in 1966.
15 So I was seeing large headed implants, metal-metal,
16 that were surviving well over 20 years.

17 To summarize the metal-metal total hip
18 replacement experience in use since 1960, osteolysis is
19 rare and severe osteolysis is exceedingly rare. Peter
20 Walker showed a few decades ago that equatorial
21 bearing must be avoided. Polar bearings work well,
22 and large heads work well.

1 There you can see the difference between
2 an equatorial bearing with the head jamming with
3 increasing load and the successful polar bearing where
4 the ball is slightly smaller than the inside of the
5 cup.

6 Now, it eventually dawned on me that what
7 we needed to do to resurrect hip resurfacing as a
8 viable operation was to take the concept of
9 resurfacing and put with it the large headed metal-
10 metal that had been clinically proven since 1960. And
11 interestingly, that thought didn't just occur to me.
12 It occurred also to the late Professor Heintz Wagner
13 from Germany, and without knowing what each other was
14 up to, we ended up both putting in our first metal-
15 metal hip resurfacings in February 1991.

16 So to summarize our first six years'
17 experience with metal-metal, we started with a
18 loosening problem. We then had to change the
19 component design, but cementless fixation of the
20 socket was best. Cemented fixation of the femoral
21 head was best.

22 The operation eventually was extremely

1 reliable. Hip function was good, and it satisfied all
2 of those goals that we had talked about in the '70s,
3 namely that where failure occurred revision to a total
4 hip replacement was easy because the femoral canal
5 was not violated.

6 We had three cups available in those early
7 systems. There are now 23 cups available in the
8 Birmingham hip resurfacing system, and we designed
9 that and this first use was 1997. The implant has a
10 porous surface on the cup with a hydroxyapatite
11 coating.

12 And there you can see a single layer of
13 beads integrally cast with the substrate metal, and
14 that means that you don't have to heat the implant by
15 centering to glue on the beads. That has two effects.

16 The beads are highly unlikely to come off, and
17 second, because it's not subjected to the heat of
18 centering it does not destroy the carbide
19 microstructure of the metal, which you can clearly see
20 and which we believe is very important for
21 satisfactory metal-metal long-term bearing
22 performance.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 224-4433

WWW.NEALRGROSS.COM

1 This is a case that I had to remove for
2 hematogenous infection a number of years after
3 implantation, and you can see excellent bone ingrowth
4 and on growth onto that socket.

5 This is another case, and you can see good
6 ingrowth into the porous surface. So this device
7 appears to be acting as intended.

8 On the femoral side, we've gone for a
9 cemented component because that's what I've done for
10 the last 13 and a half years, and that's, therefore,
11 what we did with the Birmingham resurfacing.

12 On the femur, we've got microinterlock of
13 Cancellous bone into the peripheral femoral head bone
14 and a cementless stem. So we have a specific design
15 to line-to-line contact of implant on bone, and as the
16 implant is inserted, the high pressure drives low
17 viscosity cement into the peripheral femoral head
18 giving good microinterlock in succession.

19 On the right I'm not sure if you can see
20 it, but that's blown up in an attempt to show you that
21 that's a tapered stem put into a parallel hole, and
22 from this point down the stem has a gap between it and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the bone, and I did that deliberately to try and avoid
2 distal loading and proximal stress healing of the
3 implant.

4 This is the first Birmingham hip
5 resurfacing that I carried out, that was carried out
6 ever on the 30th of July 1997 in a 38 year old man. I
7 can't quite read it from there. That is his postop X-
8 ray, one year and two-year X-ray.

9 Next.

10 Five-year X-ray. Now, what can we tell
11 from this and other X-rays? Well, he's got no
12 heterotopic ossification. He's got no radiolucent
13 lines. That's not too surprising. Bach's
14 radiographic study published from Melvin in the
15 Journal of Bone and Joint Surgery showed that they had
16 no radiolucent lines.

17 There's no gross migration of this
18 component, but what about minor migration? Because we
19 know from our Swedish colleagues that early minor
20 migration can herald late loosening. Can we tell
21 about minor migration from a plain X-ray with an
22 accuracy of plus-minus three millimeters? No.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Can we tell about stress shielding? No,
2 because you need a 30 percent change in bone density
3 to be able to see it on a plain x-ray. So we need
4 special studies. So we've done the RSA with our
5 Swedish colleagues from the Carl Linsek Institute in
6 Stockholm, and there's also been an RSA study
7 published in the Journal of Bone and Joint Surgery
8 also from Oxford, England.

9 And on the femoral side, there is no
10 measurable migration out to two years. On the
11 acetabular side, there's .2 millimeters of migration
12 of the socket within the first two months, and then
13 there's no further migration out to two years.

14 If we're looking at density, we need to do
15 a longitudinal Dexter study. This has been done by
16 Kishita (phonetic) and again published in the Journal
17 of Bone and Joint Surgery, and his control group was a
18 standard proximate porous coated cementless stem, and
19 on the cementless stem the bone density in Zone 7 goes
20 down. That's what we all expect.

21 Interestingly, with the Birmingham
22 resurfacing in Zone 7, the bone density goes up. So

1 not only are we retaining bone in the femur, but when
2 we load that femoral head normally in the absence of
3 pain, the density in the upper femur goes up.

4 So hip resurfacing is not new. Metal-
5 metal total hip replacement is not new. The new thing
6 is combining metal-metal with resurfacing.

7 Now, all we need to know is: is the metal
8 line exposure to patients from a metal-metal
9 Birmingham hip resurfacing any different to a
10 contemporary metal-metal total hip replacement?

11 We've done a number of studies. We've
12 included some of these for your interest. This is
13 whole blood cobalt analyzed by high resolution ICP
14 mass spec, which is the most accurate means we have of
15 measuring these.

16 In our one-year Birmingham, with 50 and
17 54 heads, the levels and mean and standard deviations
18 are shown. The one year metasuls are shown, and they
19 were 28 millimeter heads and the five-year Birmingham
20 are show, and there's no difference between the large
21 headed Birmingham metal-metal, and the 28 millimeter
22 metasul metal-metal.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 For chromium, again, there's no difference
2 between the large headed Birminghams and the 28
3 millimeter metasuls.

4 Now, we were interested in the total
5 amount of metal lines produced. So we've got 24-hour
6 output of cobalt on a lot of patients, and here is a
7 longitudinal study of the Birminghams at two years,
8 and the amount of metal line produced per day in the
9 Birmingham with a 50 and 54 head is no different to
10 the metal line produced in the 28 millimeter metasuls.

11 At five year the Birminghams metal line
12 production is no different to the 28 millimeter
13 metasul production.

14 So to summarize, hip resurfacing is not
15 new. Metal-metal total hip replacement is not new.
16 With the new part of combining metal-metal with
17 resurfacing, I've shown you that fixation of the femur
18 and the acetabular component of the Birmingham hip
19 resurfacing is good as demonstrated by two published
20 RSA studies in the Journal of Bone and Joint Surgery.

21 The loading of the proximal femur is favorable as
22 judged by DEXA and, again, published in the Journal of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrross.com

1 Bone and Joint Surgery.

2 I've just shown you that the metal line
3 exposure to patients is no different with a metal-
4 metal Birmingham hip resurfacing compared to a 28
5 millimeter headed metasul bearing total hip
6 replacement.

7 So the conclusion is that the Birmingham
8 resurfacing device is based upon the lessons learned
9 from previous resurfacing designs and historic metal-
10 metal bearings.

11 Thank you.

12 I'd now like to pass over to Tim Band, who
13 was a young, fresh faced metallurgist when I first met
14 him, but I've grown him then by totally unreasonable
15 demands in the last ten years.

16 MR. BAND: Thank you, Derek.

17 Good morning. My name is Tim Band. I'm
18 an employee of Smith & Nephew. I'm going to present
19 on the device description and preclinical testing.

20 Before I do that, with the Chairman's
21 permission, I'd like to present the panel with a
22 physical example of the component for your physical

1 consideration.

2 PANEL CHAIRPERSON NAIDU: That would be
3 great. Thank you.

4 MR. BAND: The Birmingham hip resurfacing
5 device, as we've heard, is a device comprising of two
6 components, the femoral head and acetabular cup. Both
7 components are produced by high carbon cast cobalt
8 chromium material. The femoral head is available in
9 six sizes of four millimeter increments. These are
10 described by the external diameter of the femoral
11 component and are in sizes between 38 millimeters and
12 58 millimeters.

13 The femoral head achieves fixation and
14 stability through the use of bone cement, and on the
15 inside of the femoral head, there are six recesses to
16 assist in the stability.

17 There are 12 sizes of acetabular cup which
18 means there are two sizes of cup per femoral head
19 component. Again, they're described by their external
20 diameter, and they start from 42 millimeters in
21 diameter through to 66 millimeters.

22 So, for example, for a 38 millimeter

1 femoral head component, you could use either a 44 or a
2 46 millimeter diameter acetabular cup in the event of
3 any mismatch between the femoral device and the
4 recently prepared acetabular.

5 The acetabular cup is a porous coated, HA
6 coated, cementless device, and they're also a
7 dysplasia option and bridging cup with screws to
8 provide primary stability.

9 So to summarize the technical
10 characteristics of this system, we can say that the
11 material is a high carbon, as cast cobalt chromium
12 material. The specification for the components allows
13 the femoral head to be smaller than the acetabular cup
14 to provide a polar bearing, as described by Mr. McMinn
15 in his lessons learned presentation.

16 As I showed you, the range of devices I
17 described are appropriate to meet the anatomical sizes
18 of presenting patients. The acetabular cup is an
19 uncemented HA porous coated device, and the femoral
20 component achieves stability with the six recesses
21 inside the femoral head and the use of bone cement.

22 So both components are produced from

1 cobalt chrome molybdenum alloy conforming to ASTM F-
2 75, an ISO 5832 specifications. And the
3 biocompatibility of this material has been proven by
4 its benign clinical use since as early as the 1930s.

5 The first specification for this material
6 grade and composition was actually described in the
7 early 1960s, but interestingly the material grade has
8 remained fairly unchanged to the present day.

9 So on the slide here you can see an
10 example of a first generation metal-metal bearing
11 which had survived for over 30 years in vivo. Its
12 failure mode was typical of this type of device, which
13 was for tissue growth around the smooth cobalt chrome
14 surface.

15 We know, of course, today that this is not
16 a very good fixation surface, but the slide on the
17 right is a micrograph. This is a section of the
18 component which has then been polished and chemically
19 etched to reveal the metallurgical phases in the
20 microscope here.

21 And as we can see, it's a biphasic
22 material, which means there are two microstructural

1 phases in the material. The light background is the
2 matrix, which is essentially cobalt, chrome, and
3 molybdenum maximums, and the small, dark particulate
4 component is a precipitate carbide which is made of
5 chromium and molybdenum carbide, which forms during
6 the solidification of this alloy during the investment
7 casting process.

8 The Birmingham hip resurfacing which is
9 shown below has its microstructure shown on the right-
10 hand side, and as you can see, its specification was
11 based upon the extensive forensic study of those first
12 generation metal-metal bearings. The microstructures
13 are comparable.

14 So to talk about our preclinical studies
15 which were submitted to the FDA and summaries have
16 been made in your panel packs, all of the studies were
17 conducted in accordance with the FDA guidance
18 documents and were provided for the components, the
19 beaded surface and the HA coating. The component
20 testing included wear testing, friction testing,
21 femoral standard fatigue testing and kinematics were
22 assessed by simulating range of motion.

1 Metallographic microstructure metrology examinations
2 were also carried out to characterize the device.

3 The preclinical studies carried out on the
4 beaded surface included static shear, shear fatigue
5 and static tensile strength testing, while on the
6 substrate yield ultimate tensile elongation and
7 abrasion testing were carried out.

8 Finally, the preclinical studies on the
9 hydroxylapatite coating included environmental
10 stability, coating thickness, static shear and tensile
11 strength, and analysis of the chemical and
12 crystallographic characteristics of the coating.

13 So, in summary, the device which was
14 designed from the lessons learned that Mr. McMinn
15 presented in his earlier presentation have all been
16 evaluated on the preclinical data, and the results
17 confirm that the device showed performance intended in
18 vivo.

19 Thank you.

20 I'd now like to hand over to my colleague,
21 Mr. Marcos Velez-Duran.

22 MR. VELEZ-DURAN: Marcos Velez-Duran with

1 Smith & Nephew, and I'm to present the summary of the
2 clinical data and clinical studies.

3 This product was first available in Europe
4 in the U.K. specifically in July of 1997. Since then
5 it has been currently marked in 23 countries. Some of
6 those countries include Canada, Australia, Japan, and
7 all over Europe.

8 There have been more than 33,000 implants
9 implanted worldwide at the time of the PMA submission.

10 The evidence of safety and effectiveness
11 presented in this PMA is based on a consecutive series
12 of 2,385 cases, surgeries that occur from July of 1997
13 to 2004. The safety data that's presented in the PMA,
14 and you have available in your packet, is a review of
15 all 2,385 cases in this study. The effectiveness data
16 is based on a total of the first consecutive series of
17 1,626.

18 The effectiveness data is based on an
19 independent review and follow-up by the Oswestry
20 Center Registry, and those are surgeries from July
21 1997 to March of 2002.

22 The radiographic study that you also have

NEAL R. GROSS
COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON D.C. 20005-3701

1 in your packet is based on a total of 124 cases.
2 Those are the first consecutive 124 cases conducted by
3 Mr. McMinn.

4 The study design has been described as
5 prospective, and it's a consecutive series of like I
6 mentioned, 2,385 procedures on the BHR, and no further
7 design changes are in the series. So 100 percent of
8 the products in this series are of the same design.

9 It's a single center by surgeries
10 conducted by Mr. Derek McMinn and their Birmingham
11 hospital using a proven surgical technique.

12 There is five years' follow-up on the
13 series, and we include safety, survivorship, pain or
14 function assessment, and also patient satisfaction.

15 The strength of this data is that one is
16 consecutive clinical series, and for the safety
17 assessment there was a 100 percent audit of all cases
18 at patient records and at the Oswestry Center.

19 On the effectiveness side, the prospective
20 registry was independent of Mr. McMinn, was performed
21 by the Oswestry Outcomes Center and include, like I
22 mentioned previously, 1,626 cases, and in addition to

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the patient self-assessment of pain and function,
2 there is also a patient satisfaction question.

3 Also, as part of effectiveness, there's an
4 independent radiographic evaluation at five years
5 using a prospective protocol on 124 cases.

6 There have been some questions about the
7 comparability of the U.K. and the U.S. in terms of
8 patient populations and practice of medicine, and we
9 will surely talk about this later in this meeting.
10 Our observation is that there are similar target
11 populations in both countries, and joint surgeries are
12 performed in the same matter in both countries and
13 similar in hospital procedures.

14 There are three specific cohorts that are
15 mentioned in the PMA. There is the X-ray cohort,
16 Oswestry cohort, and McMinn cohort. I would like to
17 make the clarification that they're all related to the
18 very same patient, the total of 2,385 cases. It just
19 happened that the X-ray cohort and the Oswestry cohort
20 were patients that were followed by the Oswestry
21 Center registry, but all three cohorts were included
22 in the safety data.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 The effectiveness measurement, the primary
2 effectiveness measure for this PMA, incident by
3 survivorship and secondary effectiveness measures
4 include the Oswestry modified Harris Hip Score, which
5 we will refer to through the presentation as OSHIP
6 patient satisfaction.

7 In addition, there is the five-year
8 radiographic assessment.

9 Safety measures include the primary safety
10 measure is the number of revisions and percentage of
11 the population study, as well as all other adverse
12 events.

13 The study population of these cases
14 include mostly men with osteoarthritis.

15 In terms of accountability, we have an
16 excellent patient follow-up at five years of 90.8
17 percent, and at five-year survivorship of 98.4
18 percent. The survivorship is also available and is
19 consistent with publish reports in British Journal of
20 Bone and Joint. You can see series from Australia,
21 other parts of Europe and also U.K. The survivorship
22 estimate is consistent to the one we observe in this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 PMA.

2 In addition, the Oswestry Center followed
3 a total of 5,000 cases, including Mr. McMinn. So 140
4 surgeons collected 3,374 cases in over five years, a
5 maximum of five years to survivorship is 96.3 percent.

6 So it is consistent with what we see with Mr. McMinn.

7 So the success, other clinicians are able to repeat
8 the same success.

9 In terms of clinician rates as compared to
10 a comparative group presented in the PMA for the BHR
11 on all cohorts or the specific group that's X-ray,
12 Oswestry cohort, it's the primary population group for
13 effectiveness. You can see that the percent of
14 revision for BHR compares very well with published
15 articles on existing total hip replacements.

16 The Oswestry modified Harris Hip Score was
17 developed by the Oswestry Outcomes Center. It is
18 validated and compares well to the Harris Hip Score.
19 It's a patient self-assessment of pain and function,
20 and it took questions and scored similar to the Harris
21 Hip Score, with a difference that flexion and
22 extension questions are different from the Harris Hip

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

1/20/2024 4433

www.nealrgross.com

1 Score in that the OSHIP score asks relevant questions
2 to patients and activity of daily living.

3 There's also patient self-assessment of
4 satisfaction that's very simple and we will present it
5 later in this presentation.

6 To compare the Harris Hip Score to the
7 Oswestry hip score, I put together the two scoring
8 systems, and as you can notice, the only difference is
9 in the section of shoes and socks and deformity and
10 range of motion. Those have been captured in the
11 Oswestry score in the area of movement, and it equals
12 13 points. So at the end of the day both scoring
13 systems are up to 100 points.

14 The result of the OSHIP scores is as
15 follows. There is a baseline. There is an outreach
16 of 59.8 points and at five years, 95 points, a
17 significant improvement over time.

18 And if you were to define success based on
19 the scoring system as patients that have greater than
20 or equal to 80 points at five years, 93.2 percent of
21 those patients are considered successful. Even if we
22 were to define success on the total score as greater

1 than or equal to 90 points, at five years our success
2 rate would be 85.5 percent.

3 The second, in the patient satisfaction,
4 these are the questions that were asked of the
5 patients. You can see they're very straightforward
6 questions, no room for misinterpretation of the
7 question.

8 And at five years, patient satisfaction
9 was 99.5 percent of extremely pleased or pleased with
10 the operations.

11 On the radiographic data, there were
12 predefined failure and success criteria. The
13 predefined failure criteria is presented here as
14 presence of the incomplete or complete radiolucencies
15 or a radiolucency in all zones and migration of
16 components greater than two millimeters or a change in
17 acetabular orientation of greater than or equal to
18 five degrees.

19 The five-year radiographic success result
20 is 97.2 percent. There were three out of the 108 that
21 were radiographic failures, representative of 2.8
22 percent, and we want to note that no radiographic

1 failure in the series studied in revision.

2 Those radiographic findings are consistent
3 with what Mr. McMinn presented previously on RSA
4 studies, which are the best ways of assessing
5 migration and other radiographic observations.

6 The primary safety outcome with revision
7 and of the total, 2,385 cases there were 27 revisions
8 and probably as previously explained, that represented
9 1.1 percent of the total population. The majority of
10 the revisions were on the femoral neck fracture, and
11 they occurred as we can see by the average time to
12 revision, early, followed by infection and collapsible
13 head, which happen between two and three years on an
14 average.

15 We wanted to spend a little time
16 explaining the method of data collection for our first
17 event. We collected all first events in patient
18 charts and the OSHIP questionnaires. There was no
19 differentiation at the time of the collection between
20 clinical observation and actual event. The collection
21 of the adverse event was conducted by separate
22 consulting group. Mr. McMinn was not included in that

1 assessment or review of the chart.

2 Some hips have more than one adverse event
3 over time, and adverse event documented multiple times
4 were treated as separate events. It was a very
5 comprehensive review of everything that was in the
6 patient chart.

7 In terms of result and device related
8 events for the overall cohort, the AVN, actually all
9 of the device related events were at one percent or
10 less. For AVN, there's one percent for which were at
11 grade the one year post op, and 31 out of the 35 total
12 observations of AVN were interoperative observations.

13 On the femoral hip collapse, there were 15
14 total reported events. Eight were evacuated. Seven
15 were interoperative observations, again, femoral neck
16 fractures less than one percent, component migration
17 less than one percent.

18 In summary, we have a very large series of
19 hips, 2,305, followed for a long term, a long-term
20 follow-up for five years. Excellent patient follow-up
21 of 90.8 percent at five years.

22 There's an independent assessment of pain

1 and function by a separate research center that was
2 presented and not related to Mr. McMinn. There's an
3 assessment of patient satisfaction, independent
4 critique or evaluation, confirmation of the PMA there
5 resolved by other BHR series, and literature result
6 for a total hip replacement.

7 In the effectiveness side, the
8 survivorship was 98.4 at five years. The Oswestry
9 score was an average of 95 at five years; patient
10 satisfaction, 99.5, or extremely pleased with the
11 operation at five years; radiographic evaluation, 99.2
12 success at five years. In terms of safety, revisions
13 are 27, which represent only 1.13 percent of the total
14 population and very low instance of the adverse event,
15 like I mentioned, all of them less than or equal to
16 one percent.

17 In conclusion, we have talked before in
18 the presentations before me that the BHR device design
19 was based upon the lessons learned from previous
20 resurfacing design and historical metal-on-metal
21 bearings; that the clinical data confirms that the
22 device should perform as is intended in vivo, and now

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 the clinical data offer a reasonable assurance of
2 safety and effectiveness.

3 Thank you.

4 George.

5 MR. DeMUTH: Hello. I'm George DeMuth.
6 I'm a consultant to Smith & Nephew, and I do not have
7 an interest, financial interest, in the product or the
8 company.

9 So I'm going to provide some statistical
10 commentary and go to the next slide.

11 I had a few background comments, and then
12 I want to touch on the accountability and efficacy in
13 terms of surviving OSHIP and have some conclusions.
14 Go ahead.

15 I want to just start with some
16 interpretational issues and we'll come back to them
17 later, but the OSHIP has actually been collected
18 prospectively, but we're analyzing it down the way
19 now, and that constitutes some independent follow-up,
20 but we do have the single site issue and trying to
21 figure out what to compare against, which is part of
22 what I think the discussion will be today.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrsgross.com

1 The FDA raised questions about pooling of
2 the data. My first reaction was it's a single site,
3 single physician. People are being handled the same
4 way. So I'd like to be inclusive, but it's also of
5 interest to see if we can see if some of the subsets
6 of patients are behaving differently. They're better
7 or worse than expected.

8 And I want to take this as a brief note to
9 thank the FDA because they provide some nice
10 commentary and, I think, added to the analysis.

11 So the data I'm going to present will
12 actually be for unilateral hips in the common tables.

13 I don't have an n here. It was actually more than
14 1,100, and it will be the X-ray and Oswestry cohorts.

15 Go to the next slide.

16 And the reason I want to do the X-ray and
17 Oswestry cohorts is because there we have the data
18 follow-up coming off the OSHIP data and we can use
19 that for censoring. So in this case we know we've
20 gotten contact from the patient and have a much better
21 feel for whether they, you know, are being aware or we
22 know what they're doing, and for the McMinn cohort,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 224-4433

www.nealrgross.com

1 those patients we would have to use a cutoff date, and
2 that would be conservative.

3 So we'll focus on the group that we know
4 the best about.

5 Go ahead.

6 All right. So this is actually the OSHIP
7 thing. I want to start and talk about here on the
8 bottom of the slide because in terms of survival we
9 want to know what their last visit is, and they may
10 have some intermediate missing OSHIP data, but we want
11 to know at the end how they're doing.

12 So we had a very nice, 91 percent follow-
13 up at five years across the cohorts, and very good,
14 88, greater than 85 percent three and four. So I
15 think in terms of survival in that kind of follow-up
16 we're very happy about the data.

17 Now, there's some missing data at the one
18 and two year data on the OSHIP insert, some missing
19 baseline that doesn't show up here as well, and so in
20 terms of evaluating the OSHIP, we have to have that in
21 mind.

22 Now, sort of back up to the top, there

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 224-4433

www.nealrross.com

1 wasn't a flag in the Oswestry Outcomes Center that
2 said patients were discontinued. So what we did is
3 just to try to get some idea. That was to look at
4 whether patients had missed their last two expected
5 visits, and that doesn't even mean that they won't
6 come back in, but in this case, it was 84 out of the
7 1,626 hips. That's just slightly more than five
8 percent.

9 So I think in terms of the survival
10 endpoint, that's good as well. So let's go ahead.

11 I don't want to spend a lot of time here,
12 but this just a report. Really the rate here for this
13 cohort, the unilateral X-ray and Oswestry cohorts:
14 98.3 percent survival, 1.7 percent revision rate. The
15 upper bound is 3.1 percent of revision rate. So it
16 looks solid as we expect with a lot of patients.

17 So let's go.

18 So without trying to do some of these
19 covariate analyses, I think there will be some
20 discussion later about P values. I'm just presenting
21 both here. None of these are significant for the
22 cohort and the gender and the age at five years.

1 There's, you know, a very small survival difference.

2 I wanted to pick up and just look at the
3 baseline data, people that had baseline and didn't.
4 It's not different, but the people that didn't have
5 baseline available were just slightly worse, 1.7
6 percent -- 1.3 percent. Sorry.

7 Let's go ahead.

8 We looked at the diagnosis. This is the
9 reason for resurfacing. We saw a significant or
10 marginally significant effect for AVN. Nothing else
11 was different. All of these P values are pair-wise
12 comparisons of the osteoarthritis group. So I don't
13 think that's an unexpected result.

14 Let's continue.

15 This just shows the all hips, unilateral
16 versus bilateral. Bilateral is still a little bit
17 better, but again, good across the board.

18 Go ahead.

19 So the X-ray and Oswestry cohort, I think,
20 is a good focus here. It's the best follow-up
21 information we have, and it looks very good in terms
22 of the data, the file we have and the revision rates,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 less than two percent at five years, 3.1 percent upper
2 bound.

3 There's a limited number of revisions.
4 Probably doesn't have a lot of power to see subgroup
5 effects, and there was one significant effect in terms
6 of AVN.

7 So go ahead.

8 This is the OSHIP data, and this is
9 observed data, and I don't want to spend long here
10 because you've seen it, but a big increase from 60 to
11 mid-90s. That's basically on average 35 to 36 points
12 as an increase, and it's just very large. So let's go
13 to the next slide.

14 I think this probably to me is a little
15 bit better because you see the majority of patients
16 less than 70, and then everybody is piled up to be
17 greater than 90 or greater than 80. So it looks like
18 you have a very strong response in the OSHIP.

19 And so this summary is just a big
20 response, a good response for five years.

21 We did model the OSHIP two-plus year
22 results. That's basically figuring if they had a two-

1 year point, you took that one or you took the next one
2 after it. Almost everybody has that, and there were
3 some significant effects there. I think we just have
4 a lot of power because now we're talking about, you
5 know, 1,100 hips and effects on the order of three
6 points, you know, would be significant, three to four
7 points.

8 And in terms of a 35 point increase, I'm
9 not sure that we have to worry about that a lot.

10 So if you touch on this missing data issue
11 back here, a couple of cuts of it and a couple of ways
12 to look at it. One, it summarized patients that had
13 missing baseline, looked at their post treatment, and
14 you just compared that. You wrote sort of
15 descriptively to the other patients. Those
16 differences are less than three points, around two
17 points, indicates that those patients may have had a
18 slightly lower mean.

19 Didn't repeat a measured mile, but there's
20 very little difference there, but that's probably
21 dominated by the observed data.

22 And then the other point is that if you

1 look at the two-plus data, it looks like, again, that
2 has very high value. So I don't think you can rule
3 out an association here, but it looks small in respect
4 to the big improvement we're seeing on the average.

5 Let's go ahead.

6 This is just patients and just to point
7 out and we saw it before; all most everybody is
8 pleased, and this coincides with high survival, low
9 revision rates and the very big improvement in OSHIP.

10 So let's go ahead.

11 This is going to touch back. I hope to
12 capture here a compendium of the issues that will be
13 raised later on about the study design, a single
14 investigator, no control. There are no a priori
15 sample size. I think some of these I answered, and
16 then we have a very large sample size. I think the
17 OSHIP was very sensitive to reasonably small
18 differences relative to the treatment size.

19 I think we had very good follow-up in
20 that. We got very tight confidence bounds around the
21 survival. I think the survival is a very objective
22 endpoint here, and so I think that's favorable.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We have external radiographic review. It
2 helps, I think that the efficacy data is being
3 collected in terms of OSHIP externally and an
4 independent safety review. So I think it's important
5 to put this study sort of in the context of the other
6 studies that we had, the published literature where we
7 see similar survival rates and nice or, you know, high
8 survival rates and high OSHIP satisfaction look
9 internally consistent.

10 So second page.

11 There aren't too many comments about this.

12 There is a question about the study design for the
13 validation work, and I don't have too many comments.
14 Just say there were some usually high correlations.
15 There was a difference in the mean OSHIP and HHS. At
16 least it's in the direction that the OSHIP is lower,
17 which would make it seem more conservative.

18 I think, too, we have to go back and look
19 at OSHIP in terms of patient satisfaction, that
20 consistency as well.

21 All right. Go ahead.

22 So in summary, back to sort of what I

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 224-4433

www.nealrsgross.com

1 thought were data points here, there's a large number
2 of procedures. We have, I think, very good follow-up.

3 Survival rate in the X-ray and Oswestry cohort for
4 the hips of 83. -- 98.3 percent, sorry. Large OSHIP
5 improvements. Patient satisfaction look very good,
6 and to the extent we did a sensitivity analysis, our
7 evaluations of the missing baseline data, missing data
8 don't lead me to believe it should interfere with a
9 major interpretation of a very large increase from
10 baseline.

11 This brings me back to the slide we've
12 been building on. I think we've started with the
13 lessons learned, and there's a history about how the
14 device is what it is and how we should have
15 expectations of good survival for this device. In the
16 preclinical data, overview of the clinical data and
17 now some statistical analyses showing, you know, very
18 good effects and a lot of data with good follow-up.

19 So I'll turn it back to Neal.

20 Thanks very much.

21 MR. DEFIBAUGH: Hello. My name is Neal
22 Defibaugh, and I'm an employee of Smith & Nephew.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 244-4433

www.nealrsgross.com

1 I would like to spend just a few minutes
2 talking about some approval considerations, namely our
3 proposed labeling, training and post approval study.

4 Next slide, please.

5 We believe the data that we have shown you
6 today demonstrates that the BHR device is intended or
7 can be used for patients who are at risk of requiring
8 more than one hip joint replacement over their
9 lifetimes, and while it's impossible to predict
10 exactly who may require a future revisions, some
11 factors that are known to increase the risk of
12 revision or include patients of a young age and
13 anybody less than 55 years and/or high activity level.

14 Next slide, please.

15 Therefore, specific indications for use
16 that we believe the data supports include
17 noninflammatory degenerative joint disease such as
18 osteoarthritis, AVN, dysplasia, DDH, as well as
19 inflammatory DJD, such as rheumatoid arthritis.

20 Next slide please.

21 Obviously there are contraindications to
22 the use of the BHR device. General complications you

1 would see for any type of hip replacement include
2 things like infection, patients who are skeletally
3 immature, and patients with conditions that would
4 compromise the implant stability or the post operative
5 recovery period.

6 Review of clinical data also indicates
7 that very clearly patients with inadequate bone stock,
8 such as severe osteopenia, patients with avascular
9 necrosis, greater than 50 percent involvement of the
10 femoral head, and patients with cysts, multiple cysts
11 in the femoral head greater than one centimeter.
12 These types of conditions are not appropriate for use
13 of the BHR system.

14 Additionally, out of an abundance of
15 caution, given the metal-metal nature of the device,
16 we believe that women of child bearing age and
17 patients with renal failure should be contraindicated.

18 Just briefly, I'll touch on training.
19 It's our intention to send a group of what we call
20 core surgeons to view live surgery in Birmingham,
21 U.K., as well as receive additional lecture and
22 workshop interactions. These surgeons would then be

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the ones who would come back to the U.S. and
2 ultimately after approval train other surgeons who
3 have an interest in the BHR system.

4 This training would be held in North
5 America, and there would also be additional
6 information available for surgeon resources, such as
7 surgical techniques and web site information with
8 advice and technique information.

9 FDA requested us to provide a post
10 approval study protocol, and although we do not
11 believe necessarily that a post approval study is
12 needed due to the long-term follow-up and large
13 patient population we already have, we certainly
14 submitted a protocol that was based upon a template of
15 a previously approved device, Class III hip device
16 that we had provided the FDA earlier. So that's the
17 template for this, generally a prospective,
18 nonrandomized survivorship study with clinical
19 radiographic evaluations through five years.

20 Next slide.

21 In years six through ten, we will do
22 postcard follow-up. Any explants returned to the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 sponsor would be analyzed as appropriate, and we would
2 make appropriate progress reports to FDA obviously.

3 So to round out the theme we've been
4 building on, the device design of BHR is based upon
5 lessons learned from previous resurfacing designs and
6 metal-metal bearings. Extensive preclinical testing
7 confirms the device should perform as intended in
8 vivo. The clinical data confirms there's a reasonable
9 assurance that the safety and effectiveness of the
10 device.

11 The statistical analysis demonstrates
12 robustness of the efficacy results, and we believe,
13 finally, that the clinical results support the
14 proposed labeling.

15 I would like to thank the panel and FDA
16 for all of their efforts on behalf of this PMA.

17 Thank you.

18 PANEL CHAIRPERSON NAIDU: Thank you.

19 I'd like to thank the sponsor and the
20 representatives for their presentations.

21 I'd like to address the panel now.
22 Remember that you'll have time for questions this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 afternoon, but nevertheless, if you have extensive
2 questions for the sponsor at this point, please ask it
3 so that they'll have some time to prepare for the
4 afternoon session.

5 Does anybody on the panel have any
6 questions at this time for the sponsor? Dr. Mabrey.

7 DR. MABREY: With regards to the surgical
8 training of the core surgeons and then the subsequent
9 training of surgeons within the United States, if you
10 could provide us some information on the length of
11 that training, the number of cases that they would
12 see, whether or not observation of the trained
13 surgeons within the United States would be performed,
14 and how they would be evaluated as to their
15 proficiency in implanting the device; whether or not
16 that would be a requirement for using the device or
17 simple attendance at the training course.

18 PANEL CHAIRPERSON NAIDU: Thank you, Dr.
19 Mabrey.

20 Anybody from sponsor willing to tackle
21 that question?

22 MR. VELEZ-DURAN: I just wanted to make

1 sure. So you want us to address now or later, but
2 we're ready to address it now.

3 PANEL CHAIRPERSON NAIDU: Yeah, we have a
4 little time. So we can go ahead and we've got 15
5 minutes before we break.

6 MR. VELEZ-DURAN: To discuss the issues of
7 training we have brought a couple of people that are
8 working very closely on the plan, and I would defer
9 that question to one of my colleagues, and also we
10 have some clinicians that have also participated in
11 that discussion.

12 So Marc.

13 DR. THOMAS: My name is Dr. Marc Thomas,
14 and I'm an employee with Smith & Nephew.

15 The training of the core group of surgeons
16 will initially be a two-day, on-the-ground course
17 associated with didactic issue and hands-on component
18 which will be in the form of sawbones and also viewing
19 live surgery.

20 After that course that they will attend,
21 we had to time it within an acceptable period of when
22 the people or the surgeons will get the device in

1 their hands. Before that occurs, we will have another
2 meeting where, once again, there will be follow-up
3 lectures to solidify the information for that core
4 group of surgeons.

5 There will be another opportunity for a
6 hands on component, and there will be another
7 opportunity to view our surgery to see any tips or
8 tricks and have an opportunity for a last question and
9 answer session with one of the designing surgeons from
10 the U.K.

11 At that stage if the FDA allows, they will
12 be allowed access to the device and work on their own
13 proficiency through performing the surgery on their
14 own patients for a period of time until they feel that
15 they have the self-proficiency to stand up and teach
16 their colleagues.

17 DR. MABREY: Will there be an opportunity
18 for cadaveric dissection as well or will the entire
19 procedure be conducted on sawbones?

20 And just for the audience's clarification,
21 sawbones are basically plastic representations of the
22 bones without soft tissue attachments versus

1 cadaveric, which includes all of the soft tissue and
2 muscle attachments.

3 DR. THOMAS: There will be an overture in
4 the United States for use of the instrumentation on
5 cadaveric specimens.

6 DR. MABREY: Thank you.

7 MR. VELEZ-DURAN: I would also like to
8 take the opportunity to introduce Dr. Rogerson. He is
9 a U.S. clinician from Wisconsin, and he has reviewed
10 also our plan for training and we'd like his comments
11 on that.

12 DR. ROGERSON: Thank you.

13 I'd like to preface this that I am an
14 orthopedic surgeon in Madison, Wisconsin, specializing
15 in the treatment of shoulder, hip, and knee from
16 arthroscopy to joint replacement. I've been very
17 active in the Arthroscopy Association, presently on
18 the board of directors, and vice chair of the
19 Orthopedic Learning Center where we do cadaveric
20 courses and training for arthroscopic surgeons
21 learning shoulder, ship, and knee arthroscopy.

22 I have no financial interest in this

1 whatsoever, although my air fare and hotel room are
2 being reimbursed by Smith & Nephew, and my main
3 interest in this proceeding is in trying to help
4 facilitate the introduction of this technology in the
5 United States. Some of that is a selfish, personal
6 reason. I have some arthritis brewing in my hip, and
7 I have been very interested in looking at alternative
8 to conventional total hip arthroplasty, as a lot of my
9 patients. And I have personally send about 30
10 patients over to Europe to have Dr. Kunda Schmidt
11 perform the Birmingham replacement, and I have been
12 basically amazed at the results of these patients
13 coming back and their activity level.

14 I've also visited Dr. Schmalzeried,
15 Michael Mont and staff here in the United States. So
16 I have no particular allegiance to Smith & Nephew, but
17 I do feel that the technology that is being presented
18 to you is extremely important, and I hope that it at
19 some point will make its way through the FDA.

20 In terms of the training program, which I
21 think is imperative because when you look at this
22 replacement; it is a different animal than a

1 conventional total hip replacement. You actually do
2 not have the exposure that you'd have by cutting off
3 the head and neck to evaluate the acetabulum.

4 The soft tissue dissection is critical and
5 the learning curve for this procedure is going to be
6 much greater for surgeons in any setting for this
7 procedure.

8 Smith & Nephew, as you've seen, has had
9 extensive experience with this procedure worldwide.
10 They've developed a well structured, educational
11 experience that has previously been employed in other
12 countries, and what it really consists of is graduated
13 release of this prosthesis into the market so that it
14 basically guards against a full release, a full launch
15 with poor training and ultimately poor results.

16 The core group of surgeons that would go
17 to Birmingham for the two-day intensive training would
18 be subjected to educational didactic lectures, which
19 would go through the history of the development, the
20 tribology, the metallurgy, the surgical technique, the
21 indications and contraindications for the procedure,
22 and particularly the surgical technique, and that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 would be augmented by surgical observation of Dr.
2 McMinn or Dr. Treacy, and also work with sawbones in
3 Europe would have a hands-on laboratory in that
4 setting.

5 That would hopefully be about two months
6 before the introduction of the device on the market,
7 and then right before the introduction of the device
8 on the market, there would be another intensive
9 training session for that core group of surgeons, and
10 they would go to a learning facility, again have
11 didactic lectures, but more importantly, more exposure
12 to the surgical procedure both in terms of DVD/CD
13 education. Sometimes you can actually see more. My
14 experience at the learning center is that although the
15 cadaveric dissection is critical, sometimes a well
16 developed and particularly well executed DVD of the
17 procedure where you can zoom in and show the critical
18 parts of that procedure are equally as valuable as the
19 cadaveric dissection.

20 So during that, right before the release
21 those surgeons would be, again, educated as to the
22 indications, contraindications, but more importantly

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 have a very fresh experience in terms of the operative
2 exposure and operative technique.

3 They would then go back to their local
4 setting and start performing these operations on
5 patients. The desire of the company is to have field
6 representatives who are very experienced in this
7 technique, go out to this core group of physicians and
8 observe and help in the first ten cases and to make
9 sure that they are getting through their learning
10 curve, and when that learning curve has -- and that
11 may vary depending on the core group of surgeons --
12 but when that learning curve has been attained, then
13 those core group of surgeons would become the teachers
14 of the next group that would occur at the full launch.

15 So that when it comes time for the full
16 launch for the rest of the United States, there would
17 be a number of regional facilities where you have this
18 core group of surgeons that have excellent experience.

19 The participants would then go to a lab setting in
20 the United States where they would go through a very
21 similar procedure that the core group of doctors went
22 through in England, do that in the United States and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 then have the ability to have monitoring come into
2 their OR with the field representatives, but also have
3 the opportunity to visit the core group of surgeons,
4 observe surgery and get the same type of course that
5 the core surgeons got in Europe.

6 So I think that from my standpoint in
7 terms of learning how to do a procedure that Smith &
8 Nephew has thought this out. There seems to be a well
9 structured program that has, I think, been effective
10 worldwide up to this point, and I think that it will
11 be even more defined and meticulous in the United
12 States because of the wider launch that will occur
13 here.

14 DR. MABREY: How many surgeries do you
15 anticipate being part of this initial learning curve
16 for your 15 core surgeons?

17 DR. ROGERSON: I think that the level of
18 expertise in terms of comfort level in the learning
19 curve, I would think that you would need to have
20 probably 30 surgeries under your belt before you would
21 start to get through the learning curve.

22 That could vary depending on the surgeon

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 and what his experience has been. For example, I've
2 up to this point been doing a reasonable number of
3 metal on metal, big femoral head stem prostheses. So
4 I've got very good experience with the acetabular
5 components of a joint replacement when it's a metal on
6 metal.

7 So the learning curve of the acetabular
8 might be easier for one surgeon versus another, but I
9 think that it would be hard to put a number of cases
10 on it, but I think the comfort level would be
11 documented by the company and by the field
12 representatives when somebody has gotten through their
13 learning curve.

14 PANEL CHAIRPERSON NAIDU: Thank you. Dr.
15 Blumenstein.

16 DR. BLUMENSTEIN: How many deaths were in
17 these cohorts?

18 MR. VELEZ-DURAN: I'm sorry. Could you
19 repeat the question?

20 DR. BLUMENSTEIN: How many patients died?

21 MR. VELEZ-DURAN: If you'll just wait a
22 minute, we have the data. We're looking for it.