

1 collision.

2 DR. ROBINSON: Okay.

3 MR. YANG: The air bag situation we model
4 as impact testing, not as fatigue testing.

5 DR. ROBINSON: I have one more for Dr.
6 Cunningham.

7 CHAIRMAN WHALEN: Please.

8 DR. ROBINSON: Dr. Cunningham, towards the
9 end of your presentation you had a group of patients
10 that were explanted and then reimplanted. I think
11 there was 60 of them if I remember correctly. I may
12 be a little bit off on that. I may have missed it.

13 Is it too early to comment on those 60
14 that have been reimplanted?

15 DR. CUNNINGHAM: You're referring to the
16 revision group?

17 DR. ROBINSON: I believe so. It was
18 towards the end of your presentation, yeah.

19 DR. CUNNINGHAM: We can discuss the
20 revision group.

21 MR. PURKAIT: Before we start the revision
22 group, one question to answer Ms. Dubler about that,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 no, we have not had that incorporated in the patient
2 level. We would be doing that, the information about
3 Betadine use.

4 And to answer your, Doctor -- I can see
5 your --

6 DR. ROBINSON: Robinson.

7 MR. PURKAIT: -- Robinson, about the
8 safety margin, the question that you have on fatigue,
9 I'd just like to conclude that by saying the model
10 that you have shown here is fatigue testing only to
11 show the safety margin in terms of the load factor
12 that it can take over a time period.

13 The question that you have is data impact.
14 We have tested against a model called 35 miles per
15 hour collision. If somebody had the amount of impact
16 energy, it would be on the chest, whether it would
17 withstand the breast or not, and we found that our
18 product, it takes about three times more than impact
19 energy to cause rupture.

20 So I just want to clarify that.

21 DR. ROBINSON: That's what I was looking
22 for. Thanks.

1 MR. PURKAIT: Thank you.

2 To answer the question on the other areas,
3 I'd like to call Dr. Gene Poggio to show some of the
4 information, and then I'll have Dr. Cunningham explain
5 the clinical data on that.

6 DR. POGGIO: This actually connects with
7 what I mentioned at the beginning when I said there
8 was one exception to when we discontinued the patient
9 to explanation. We did for all of the analyses except
10 for these analyses where we actually used that as the
11 data, if you will, as a baseline for the next set of
12 patients.

13 So I'll run through this focusing on the
14 saline perspective part, but revision patients are
15 basically defined by the FDA as patients that are
16 replacing their original implant regardless of whether
17 your original implant was for augmentation or
18 reconstruction.

19 And in the saline prospective study -- and
20 I must apologize here. The 196 is actually the number
21 of devices. It's 124 patients, and the 215 in the
22 large, simple trial is, indeed, the number of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 patients.

2 Next slide.

3 This slide showed the complication rates
4 in the LST for revision patients with 95 percent
5 confidence intervals, in the "large simple trial," I
6 should say, and this is on a per patient basis at 12
7 months.

8 And the next slide.

9 And now we're looking on a per device
10 basis at 36 months, and we did actually so that we
11 have the estimated rates for the major complications
12 here with 95 percent confidence intervals, and we
13 tested whether there was a significant difference
14 between -- and I'm sorry. This is for prior, where
15 the previous implant was for augmentation. The next
16 slide is the reconstruction.

17 We compared whether there was a
18 statistically significant difference between these
19 estimated rates and the rates per the original
20 augmentation, and there was no significant difference
21 with the exception of explantation, which was somewhat
22 higher.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Next.

2 And then for reconstruction, again, the
3 same five major complications with estimated rates,
4 complication rates, and 95 percent confidence
5 intervals, and again tested for significant difference
6 between this, and in this case the original
7 reconstruction, and we have two significant
8 differences, one up and one down.

9 And then we also have information on
10 effectiveness, and I think I'd rather turn that over
11 to Dr. --

12 MR. PURKAIT: -- show the complication
13 rate and then Dr. Anderson will show the
14 effectiveness.

15 MR. POGGIO: Okay.

16 DR. CUNNINGHAM: It was of interest to me
17 to try to determine or theorize why these rates of
18 clinically significant changes in these patient
19 cohorts over their primary indication for
20 implantation, and to come up with a clinical story
21 that explains it.

22 The explantation group, which is higher

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 after the revision, the causes for explantation, the
2 second explantation, are capsular contracture,
3 deflation, and infection, and I think what we're
4 seeing here is a group of patients who had an
5 intervention to try to solve a problem. They had a
6 capsular contracture. They had some other problem and
7 had an implant placed in an attempt to solve that
8 problem.

9 And I think this increased explantation
10 rate shows that the problem was attempted to be
11 solved, but was not, in effect, solved. So a patient
12 who in a clinical situation, as we were discussing
13 earlier, might have the signs and symptoms of a
14 cellulitis, you try to deal with it with intravenous
15 or oral antibiotics. It does not resolve. You offer
16 the patient the choice.

17 The choices are: we remove your device,
18 let you have no device for a period of time while you
19 heal, and then replace the device. That's one option
20 that we offer patients.

21 Another option that we offer patients is
22 we can go in, we can take out the infected device, we

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 can clean the space, we can put you on oral
2 antibiotics or IV antibiotics, and we can try and save
3 you a third operation or another operation and put the
4 implant back in at this time.

5 And I think what we're seeing here is that
6 there are times where we try to do that and were not
7 successful.

8 With respect to the reconstruction
9 patients, I think the decreased number of capsular
10 contracture Baker Grade III or IV that we see after
11 reimplantation indicates the opposite story, but this
12 is a group of patients where we were able to
13 successfully treat a problem, namely, that of capsular
14 contracture, by operating on the patient, taking out
15 the scar capsule contracture, dividing through the
16 scar capsule, whatever, and that this decreased number
17 indicates we've been successful.

18 With respect to deflation, it's a more
19 difficult question to answer for me clinically as to
20 why this group has a slightly greater risk of
21 deflation than when they originally had their device
22 placed, and I feel that I have to come back to the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 sense that it's a more adverse environment.

2 Perhaps if the surgery was done for a
3 capsular contracture and the contracture continued to
4 exist for some reason, we know that that might be
5 associated with a greater incidence of rupture.

6 Now, the bigger question is: how do
7 patients who have already undergone a frustrating
8 experience? They've had high hopes of success. The
9 operation has been done, and it hasn't worked out.
10 They've had to have another implant placed. So a very
11 significant question is: how do they take this? How
12 do they respond to it? What's their satisfaction
13 level?

14 And I would like to ask Dr. Anderson to
15 spend a second or two talking about that.

16 DR. ANDERSON: We had used the breast
17 evaluation questionnaire to assess patient
18 satisfaction in the saline perspective study. So we
19 decided to look at patient satisfaction in this
20 revision group of augmentation patients on the three
21 aspects, size, shape, and firmness.

22 And as you can see, despite the fact that

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 they experienced revision procedures, they still
2 reported very high levels of satisfaction at the 36
3 months follow-up with respect to size, shape, and
4 firmness.

5 CHAIRMAN WHALEN: Dr. Bandeen-Roche.

6 DR. BANDEEN-ROCHE: I believe --

7 CHAIRMAN WHALEN: Excuse me one second.

8 On that same issue?

9 PARTICIPANT: No.

10 CHAIRMAN WHALEN: We're going to someone
11 else and then I'll get to you.

12 DR. BANDEEN-ROCHE: I believe my question
13 series is for Dr. Anderson: concerns about the
14 quality of life data, and certainly include the lack
15 of a control group.

16 And so first of all, just correct me if
17 I'm wrong, but my understanding is you really don't
18 have any true quality of life data for the
19 augmentation patients. It's body appearance and self-
20 esteem data rather than quality of life, is it not?

21 DR. ANDERSON: I suppose you could
22 characterize it as body image.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. BANDEEN-ROCHE: Okay. In terms of the
2 slight increase of self-esteem, I'll just voice my
3 concerns, and I just want you to respond. If women
4 went into the surgery at a particularly low point of
5 their self-esteem, which is certainly reasonable,
6 slight increase would be consistent with regression to
7 the mean, wouldn't it?

8 DR. ANDERSON: I'm not a statistician,
9 but, yes, I understand the concept you're talking
10 about, and that's probably true.

11 Are you referring to the Tennessee self-
12 concept scale?

13 DR. BANDEEN-ROCHE: Yes, I am.

14 DR. ANDERSON: Okay. With respect to that
15 scale, we've acknowledged that that scale is maybe not
16 the best assessment to have been utilized, and that's
17 one of the reasons I didn't present it in my
18 presentation, even though we did achieve clinical
19 significance. We didn't know if the results were --
20 I mean statistical significance -- we didn't know if
21 the result were clinically meaningful.

22 DR. BANDEEN-ROCHE: Okay. Thank you.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 My next question goes to the
2 reconstruction patients. So you noted increases in
3 quality of life on the FLIC. I'm just going to state
4 something, and just correct me.

5 My impression is that those increases are
6 not distinguishable from increases that could have
7 happened anyway just because they've recovered from
8 their cancer.

9 DR. ANDERSON: It's very possibly true
10 that they would have over a period of three years
11 adjusted to their cancer and shown an improvement. In
12 my clinical experience, however, I do see that there
13 is a tremendous amount of satisfaction with
14 restoration of the breast in these patients.

15 DR. BANDEEN-ROCHE: Thank you.

16 CHAIRMAN WHALEN: Ms Brinkman.

17 MS. BRINKMAN: Yeah. Along that same
18 vein, I'm interested in your Beck depression inventory
19 because in reconstruction patients, you say, you know,
20 that their scores have decreased, but, I mean, is that
21 a decrease, and how do you know the difference whether
22 it's a decrease due to the fact that they've finished

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 cancer therapy treatment or whether it's actually the
2 fact that they've had an implant?

3 How do you separate that?

4 DR. ANDERSON: Well, you know, it was
5 relatively nice to see that they weren't a real
6 depressed group of patients to start with, which is
7 consistent with one of my studies, which looked at
8 psychological adjustment in breast reconstruction
9 patients, and I suppose that it is theoretically
10 possible that levels of depression would have
11 decreased over time in these patients.

12 Again, I relate to my clinical experience.
13 These patients are overwhelmingly satisfied and
14 pleased to wake up from surgery with a breast mound.

15 CHAIRMAN WHALEN: Dr. Burkhardt.

16 DR. BURKHARDT: I wasn't going to open the
17 door to Betadine, Dr. Cunningham, but it's already
18 been opened, and I think I have to walk through it.

19 My recollection is that Mentor initially
20 sent out a flyer to users of breast implants saying
21 that Betadine was a problem with the integrity of the
22 implant, and that was from a study in which Betadine

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 had been placed within the device, and the problem at
2 that time was a valve failure.

3 And then there was some furor about it
4 within the plastic surgery community, and as I recall,
5 Mentor did another study or perhaps a parallel study
6 with the implant immersed in Betadine solution and
7 found no problems with that in vitro.

8 Now, is my memory of that correct?

9 MR. PURKAIT: Some of them, correct, and
10 some of them -- if I may have your indulgence, I'd
11 like to kind of go and kind of give this.

12 When we looked into the Betadine, this was
13 brought to our attention by many different surgeons.
14 They are the one who called us and said, "Look. Maybe
15 you should take a look at it because some of the
16 implants are showing failure because of some reason we
17 do not know."

18 Well, when we started looking into their
19 information and the data, we realized that there was
20 a large amount of Betadine was used with our implants
21 in all conditions, whether it has been soaked or put
22 inside the cavity or put inside the implant.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 We initiated our internal investigation
2 and studies to understand was there any relationship
3 between the Betadine and the implant failure. We have
4 done several experiments three different times, one
5 with the solution inside, one with the soaking, and
6 also to match up the acidity of this particular one to
7 make sure that the acidity doesn't have anything to do
8 with it.

9 So the series of experiments have
10 indicated that, one, either the implant does fail even
11 with the contact, and when we learned that --

12 DR. BURKHARDT: I'm sorry. I didn't
13 understand when you said that. That the implant does
14 fail?

15 MR. PURKAIT: Failed, yes. The implant
16 failed. When I said "failed," it means the pads tends
17 to come out of the shell or the shell itself, the
18 surface of the shell looks like getting weaker, and
19 you can easily probably break through that. Those are
20 the kind of observations we have seen.

21 And when we saw that, we realized that any
22 way we want to do it -- in other words, if we can go

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 back and probably do some more experiments to tease
2 out exactly what conditions and what time, we kind of
3 felt that this is our responsibility to contact the
4 agency with that information, and we did so, and at
5 that point in time, the agency and ourselves decided
6 that we should put that immediately with the patient
7 information, with information that with even the
8 slightest contact will provide or will probably fail
9 or show the loss of integrity of the implant in the
10 future.

11 So that's where it is.

12 DR. BURKHARDT: For Dr. Cunningham, were
13 the failures that you observed with the Betadine
14 irrigation, were they valve failures or were they the
15 usual fold flaw failure, or do you know?

16 MR. PURKAIT: We have the failure, the
17 deflation, as we have shown you before. We only saw
18 two fold flaw failure. We did see some failure
19 because of the tear of the shell. That's the largest
20 number.

21 Now, it's very difficult sometimes to
22 exactly identify the tear was already there or it was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the loss of the shell thickness for some reason. So
2 we couldn't really identify that very well.

3 DR. BURKHARDT: Thank you very much.

4 CHAIRMAN WHALEN: Dr. Morykwas.

5 DR. MORYKWAS: I just wanted to ask a
6 follow-up question of Dr. Cunningham on the infection
7 that I brought up before, and just to simplify things,
8 I'll say an early infection is one that is apparent
9 before the incision is totally healed, and a late
10 infection is one after the incision is healed.

11 Do you have the percentage of early versus
12 late, and was it at all correlated to the surface type
13 of the implant, the smooth versus textured?

14 And I guess as a follow-up, did you use
15 any of the partially textured implants?

16 MR. PURKAIT: The last question first.
17 Partially textured implant was not used to understand
18 that phenomenon, but as far as the breakdown of those,
19 we'll talk to Dr. Gene Poggio to see if he can tease
20 out the information for you.

21 DR. MORYKWAS: Okay.

22 DR. POGGIO: I can answer part of your

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 question pretty readily. Looking overall, across
2 augmentation and reconstruction, and I have it
3 separately for those if you'd like, in the first year
4 -- this is using the Kaplan-Meier estimates and
5 looking at changes from year to year. So how much
6 happened in the first year and then how much
7 increased; did it increase between the end of the
8 first year and the end of the second year?

9 So infection overall, 2.8 percent in the
10 first year, 0.5 percent in the second year, and 0.16
11 percent in the third year. So it's almost all in the
12 first year.

13 DR. CUNNINGHAM: But in terms of teasing
14 out, I mean you're asking for a time frame of one or
15 two weeks, and the first interval follow-up were data
16 reported as four to six weeks, and these would
17 presumably be detected earlier than that, but they
18 would be, you know reported as they accumulated.

19 I think in clinical practice they occur
20 most frequently within the first two weeks, and it's
21 a little bit different than a wound infection without
22 a device because most plastic surgeons, as Dr.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Burkhardt has teased out, try to use antibiotics,
2 Betadine, whatever they can, to make this complication
3 go as close to zero as possible.

4 So there are times where the -- plus
5 giving prophylactic antibiotics -- so there are times
6 where the time course is shifted or delayed out into
7 the future than in a wound where there is, you know,
8 no prophylactic antibiotic, not as much irrigation,
9 but clinically my impression is that they almost all
10 are most apparent within the first two to three weeks,
11 and it's very rare that you see a late complication
12 associated with, say, some dental procedure or some
13 other seeding.

14 Here's the time to occurrence. Four to
15 six weeks is 53 percent, which is the majority. Six
16 months is 24 percent. Twelve months is 17 percent for
17 reconstruction patients, and then in terms of the
18 infection by device, the textured device, both the
19 SILTEX and the SPECTRUM, are more likely -- it looks
20 like about 85 to 90 percent -- are more likely to have
21 an infection.

22 CHAIRMAN WHALEN: Dr. Change.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. CHANG: I'm going to quote Dr.
2 Burkhardt again because it is my concern that on the
3 face of it, although in the ten year retrospective
4 study patient satisfaction rate was high, 93 percent,
5 could you help me out? And I presume you're looking
6 at infection, significant capsular contracture and
7 deflation rate. Could you help me out and in your own
8 words, how would you explain to a consumer, to a
9 patient that despite one in four complication rate at
10 ten years, that this is, indeed, a safe product?

11 Overall complication rate is 27 percent.

12 DR. CUNNINGHAM: I can speak to, you know,
13 what I see clinically, and perhaps Dr. Anderson can
14 speak to that as well.

15 First of all, no plastic surgeon wants a
16 surprised or unhappy patient, particularly when we're
17 doing elective surgery, aesthetic surgery. So I think
18 one of the ways to explain the fact that despite a one
19 in four risk for complication or reoperation patients
20 are generally satisfied goes to the degree to which
21 they are informed.

22 If I as a plastic surgeon whitewash the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 possible complications and have the kind of risk of
2 complications that we've demonstrated today, I'm going
3 to have a lot of very unhappy, surprised patients
4 saying, "How could this happen to me?"

5 Whereas if I go as far as I can to stress
6 what the risks are, make clear what things could go
7 wrong, and make clear that the patient understands
8 that and we're not pretending that it's not going to
9 happen to them, we're saying it could happen to you;
10 I think that sets an expectation set that makes
11 anything that does occur more acceptable to a patient
12 and not something that's going to diminish their
13 overall result.

14 But I think when we ask the patients to
15 rate themselves on a strongly dissatisfied,
16 dissatisfied, satisfied, very satisfied -- I can't
17 remember what the fifth one was -- the vast majority
18 of them also when we asked them would you do it again,
19 the 90 percent range was, yes, they would.

20 So I think part of it is they're well
21 informed, and so they tolerate the complications.

22 CHAIRMAN WHALEN: Ms. Dubler.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MS. DUBLER: I'd like to pursue that just
2 a little bit because I have an epistemological
3 problem. In the context of research where IRBs have
4 to weigh the risk and benefit, there's an articulate,
5 although not the majority of scholars who argue that
6 the informed consent process can cure the defects in
7 the risk-benefit ratio; that it doesn't really matter
8 what the risk-benefit ratio is. If you can tell the
9 patient and the patient can make an informed choice,
10 it solves your problem.

11 But in the context of the FDA's finding
12 that something is safe and effective, I'm not sure
13 that an informed consent response solves the
14 underlying problem. So we have a 27.6 overall
15 complication rate.

16 I don't know. I think women want it. I
17 think they're satisfied. Your data seems to show
18 that, but I don't know how we can find it safe. Now,
19 maybe that's a question for the company. Maybe it's
20 a question for the FDA, but that's a problem for me.

21 DR. CUNNINGHAM: Well, I think --

22 (Applause.)

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. CUNNINGHAM: You know, I think we
2 wouldn't want to try to make the assessment for the
3 patient. We are constantly asking them what their
4 assessment is. That's the basis of the data that Dr.
5 Anderson presented. It's the basis of the data that
6 I presented.

7 I think there are two different risk-
8 benefit sets. There's obviously a different set for
9 the augmentation patient versus the one for the
10 reconstruction patient.

11 I think another thing that we haven't
12 really talked about extensively here is that we're not
13 just counting pieces of chalk marking on a blackboard.
14 The complications that we talk about, some of them are
15 significant and impact definitely on a patient's life,
16 and some of them are things that they can control and
17 determine themselves.

18 For instance, the patient who wants to
19 change their size is recorded as a reoperation, but
20 it's not an obligatory reoperation. It's something
21 that they choose to do to improve their result.

22 So it's not exactly the same model as some

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 of the others that you might be looking at because the
2 patient is in so much control of the process relative
3 to, you know, you need a heart valve, you need a
4 pacemaker, this you need. It's a different situation.

5 A lot of these operations are things that
6 patients choose to do.

7 CHAIRMAN WHALEN: Dr. Burkhardt.

8 DR. BURKHARDT: Yes. Dr. Cunningham, I
9 forget exactly what the figures were. Seventy-two of
10 88 implants that were removed were replaced at the
11 same time, and the problem was apparently size. So
12 rather than being a concern about the integrity of the
13 device or anything, it does raise some question about
14 how the size is picked in the first place.

15 How did you do that in your study?

16 DR. CUNNINGHAM: I think we left the size
17 determination to the surgeon's individual practice.
18 There's no way that the company could help the surgeon
19 decide what size the patient would need for the
20 patient to be happy, and I think certainly I have seen
21 in my own clinical practice where a patient might come
22 in with one set of expectations before they have any

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 surgery, and then as time goes on, their expectations
2 and their desires change after the operation.

3 So there are times when patients come in
4 and demonstrate or discuss a certain look or
5 appearance, and then after the surgery they say, you
6 know, "I would like to enhance that further," and
7 that's part of what we see when the implants are
8 changed, particularly for a larger size.

9 DR. BURKHARDT: I guess my point would be
10 that as I understand it, then that doesn't reflect
11 deficiency of any sort on the implant, only in the
12 decision making process as to the size originally.

13 DR. CUNNINGHAM: I think what you're
14 reflecting is a communication issue or a change in
15 communication or change in desire on the part of the
16 patient, not an implant related problem.

17 DR. BURKHARDT: Thank you.

18 MS. DUBLER: Could I follow up with that,
19 please?

20 CHAIRMAN WHALEN: Yes.

21 MS. DUBLER: I just want to follow up this
22 discussion because I think it's very interesting, and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that is the Betadine discussion and perhaps the
2 reoperation to put in a larger prosthesis reflects on
3 the quality of the practice of the surgeons more than
4 it does on the device itself.

5 So, again, it's hard for me to separate
6 out how much of this reoperation and overall
7 complication rate is related to the device and how
8 much of it is related to surgical patterns of
9 practice, and how we allocate or understand that I
10 think makes a difference in this context.

11 DR. CUNNINGHAM: I take that as more of a
12 statement than a question.

13 MS. DUBLER: Okay.

14 MS. BRINKMAN: I just have one question.
15 You're talking about your ten year study, your follow-
16 up. You're using the SEER data for that?

17 MR. PURKAIT: That's what the cancer
18 patients primary.

19 MS. BRINKMAN: And that sample size is a
20 little over 200? That's about it, of saline?

21 MR. PURKAIT: No.

22 DR. CUNNINGHAM: No, the ten year data

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 that I presented is not the SEER data. The SEER data
2 was another part of the way of looking at the
3 reconstructed patients and what their failure rate
4 was, their deflation rate.

5 MS. BRINKMAN: But in the SEER data, the
6 saline implant population number was about a couple
7 hundred for the reconstruction?

8 DR. CUNNINGHAM: That data set, I think,
9 was complicated by the fact that it was a
10 retrospective study. A lot of charts were looked at,
11 and there were definite incidences, and perhaps Dr.
12 Poggio or one of the others can pull that out, where
13 devices were recorded as implants, but clearly as you
14 look back over the medical record, they were not ever
15 meant to be permanent devices. They were soft tissue
16 expanders, not implants.

17 So I think the large number of those that
18 were soft tissue expanders and not implants kind of
19 clouds that whole data set for us and it makes us very
20 hard for us to interpret what is the actual failure
21 rate in the SEER data for reconstruction patients.

22 CHAIRMAN WHALEN: Are there other members

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 of the panel who have questions for the sponsor?

2 PARTICIPANT: One quick one.

3 CHAIRMAN WHALEN: Well, I'm asking only
4 because I've been asked to move the rest of the
5 questions to a later point in time, but if indeed,
6 we're done, we're done. But if there are others, we
7 are going to move on now to the FDA presentation, and
8 there will be time later for more sponsor questions,
9 which I would project, looking at the schedule, should
10 be some time before 2:00 a.m.

11 (Laughter.)

12 CHAIRMAN WHALEN: So we will move on to
13 the FDA presentation, and I thank the sponsor for
14 their presentation.

15 For all my fellow panel members, as the
16 FDA is coming up for their presentation, please be
17 aware it is past 5:00 p.m. We now are on overtime,
18 which means that Jim Dillard will thank us twice on
19 Friday instead of once for our work.

20 PARTICIPANT: -- other questions later on?

21 CHAIRMAN WHALEN: Not at the table, sir,
22 because the FDA will be coming up, but please if you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 could stay in this vicinity, that would be good
2 because there will be questions, and also you have a
3 summation period later on.

4 (Pause in proceedings.)

5 DR. BERKOWITZ: I'll present the FDA
6 presentation of the Mentor saline filled and SPECTRUM
7 saline filled breast implants.

8 I'm David Berkowitz, the lead reviewer,
9 and I will give an overview of the status of the
10 preclinical testing, and then I'll finish with a one
11 slide summary of the medical device reports for the
12 Mentor prostheses.

13 And then we'll hear from Sahar Dawisha,
14 who is the clinical reviewer and will review the
15 clinical results.

16 And then we'll hear from Phyllis
17 Silverman, who is the statistical reviewer.

18 To describe the device first, the saline
19 filled device is available in six styles. The styles
20 are determined by two things. One is the shape, like
21 the round, the profile or the contour, and the other
22 is the nature of the surface. The surface is either

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 smooth or textured, and SILTEX is the Mentor name for
2 textured. So SILTEX implies textured, and the
3 surfaces are either smooth, SILTEX or SILTEX PT, which
4 is a partially textured device.

5 The saline filled device has a diaphragm,
6 an anterior diaphragm valve, and of course, it's
7 filled with physiological saline, and both devices,
8 obviously the shells are made from silicone
9 elastomers.

10 The SPECTRUM device differs from the
11 saline filled device in that it can be postoperatively
12 adjusted. The volume can be postoperatively adjusted,
13 and when the desired volume is reached, the little
14 valve for postoperative filling can be removed under
15 local anesthesia.

16 The SPECTRUM device has a posterior kink
17 plug valve, and like the saline filled, obviously the
18 filler is also saline.

19 The indications for use are augmentation,
20 reconstruction, asymmetry, ptosis, aplasia, hypoplasia
21 of the breast, replacement, and combined breast and
22 chest wall deformities.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Mentor has done extensive chemical
2 analyses on the device. They've looked at the levels
3 of the volatiles, extractables, and metals, and these
4 are important because, one, they characterize the
5 materials present in the device and, secondly, they
6 determine what is there and how quickly those things
7 can diffuse out of the device to cause either local or
8 systemic toxicity.

9 These are some of the toxicology testing.
10 The pharmacokinetics testing came from the literature,
11 but also relied upon the chemical determinations sine
12 by knowing how much is present and how much could leak
13 out, we know what the dangers are.

14 So, in fact, it turns out in terms of
15 systemic toxicity, even if all of the low molecular
16 weight components present in the device leaked out
17 immediately, it would still be a wide margin of safety
18 between the levels, say, the blood levels obtained and
19 the toxic levels.

20 The middle group of things that were
21 determined are all the, I think, quite commonly done
22 biocompatibility things that are tested on most

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 devices, and I just want to mention at the end that
2 the immunotoxicity is important for these devices, and
3 it was quite extensive.

4 That is, in addition to doing simply the
5 hematology and counting cells, they enumerated, you
6 know, the B and the T cells and the T cell subsets.
7 They enumerated the killer cells, for example, and
8 they also estimated the effects of implanted shell
9 material on killer cell activity and on things like
10 the mixed leukocyte response.

11 So various aspects of the immunology were
12 broadly tested.

13 The remaining, these toxicological things
14 have to do with mutagenesis and in bacterial testing
15 and in mammalian cells, and finally culminating in a
16 two year rat carcinogenicity study which demonstrated
17 no carcinogenicity.

18 The company also did a reproductive
19 toxicology and teratology study, and that was also
20 negative, which would have been expected.

21 The mechanical testing, on the other hand,
22 is not complete. Mentor has done some mechanical

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 testing on each of these topics, and we have also
2 received some recent data which have not yet been
3 reviewed. But with all the testing we now have, we
4 think that the data are not yet sufficient to make a
5 recommendation for all the implant styles proposed in
6 the PMA.

7 So we are still going back and forth with
8 Mentor on the mechanical testing.

9 Finally now, to change topics, I'd like to
10 put up a slide summarizing the medical device
11 reporting for Mentor, and this slide summarizes the
12 medical device reports that FDA has received for the
13 Mentor saline filled breast implants during the last
14 three year period.

15 The first column, that is the Maude
16 reports -- oops, I missed the slide -- the first
17 column, the Maude reports, the Maude system received
18 reports directly from patients, health care providers,
19 practitioners, and from manufacturers.

20 The second column lists the five most
21 frequently reported adverse events that are reported
22 in summary form by Mentor on a quarterly basis. So

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 these are the sums of the reports for the last -- for
2 over the three years shown above.

3 So that's all I'll say about this, and I
4 think now I'll ask Dr. Dawisha to come and begin the
5 clinical report.

6 DR. DAWISHA: Good afternoon. I am Sahar
7 Dawisha, a Medical Officer in the Division of General
8 and Restorative Devices, and I will be presenting
9 FDA's clinical perspective of the information provided
10 in the Mentor Corporation's saline filled breast
11 implant PMA.

12 The clinical studies reported in the PMA
13 are summarized on this slide and consist of a
14 retrospective assessment of implant removal from the
15 SEER data base, a one year large, simple trial, or
16 LST; the saline prospective study, or SPS; and the
17 Mentor retrospective study.

18 The SEER and LST were conducted in
19 response to suggestions from FDA in 1994 on the type
20 of information needed for PMA approval submission.

21 The SPS is a prospective clinical study which was
22 approved by FDA in 1995 after all augmentation and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 some reconstruction patients had already been
2 enrolled.

3 Because the saline prospective study
4 contains the main safety and effectiveness
5 information, I will focus on this study summarizing
6 the SEER and LST only briefly. I will not be
7 discussing the Mentor retrospective study because the
8 patient population in this study is highly selected
9 and because data ascertainment bias severely limits
10 the conclusions drawn from this study.

11 The sponsor funded a retrospective
12 analysis of implant removal in a breast cancer
13 population cohort from the surveillance epidemiology
14 and end results or SEER program of the National Cancer
15 Institute because they were having difficulty
16 enrolling reconstruction patients in their clinical
17 studies. Women with a diagnosis of breast cancer in
18 the years of 1983, '85, '87, and '89 with any type of
19 breast implant, including silicone gel filled, saline
20 filled, and tissue expanders, were asked to respond to
21 a questionnaire regarding implant removal.

22 The results of this study are shown here.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The nonresponse rate was 20 percent overall. Of the
2 1,375 total implants removed -- I'm sorry. Of the
3 1,375 total implants, 32 percent were removed.

4 Of the 252 saline implants, 43 percent
5 were removed.

6 There was information provided based on
7 the reason for removal, and excluding the 28 saline
8 implants removed as part of planned reconstruction --
9 these are the tissue expanders that Dr. Cunningham was
10 referring to -- the reasons for saline implant removal
11 are shown.

12 Capsular contracture constituted the
13 single most common reason for implant removal, 35
14 percent of implant removal.

15 The large, simple trial was designed as a
16 prospective study of a large number of patients
17 followed only for the safety endpoints of capsular
18 contracture, infection, rupture, deflation, and
19 explantation for a total of one year. The sample size
20 of 3,000, and 5,000 patients was proposed by the
21 sponsor to estimate complication rates with a
22 precision of one to two percent.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The LST results at one year are shown here
2 on a by patient basis. The analysis method used here
3 is the Kaplan-Meier survival analysis, which shows a
4 risk of a first occurrence of a complication, along
5 with a 95 percent confidence interval, which is shown
6 in parentheses.

7 Not that for the total group, which is in
8 this column, the confidence interval are all within
9 one to two percent as proposed by the sponsor.
10 However, the intervals for the reconstruction and
11 revision patients are much larger, and in some cases
12 there was insufficient information to estimate the
13 proportion.

14 There were a total of 2,373 patients
15 enroll with the majority as augmentation, and the
16 follow-up rate at one year was approximately 47
17 percent.

18 Of the four complications studied here,
19 capsular contracture, Baker Grade III or IV was
20 generally the complication encountered with the
21 greatest overall frequency.

22 Furthermore, you can see that for the two

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 complications in which there was sufficient
2 information, which would be explant and capsular
3 contracture Grade III or IV, the revision patients
4 generally have rates between those of augmentation and
5 reconstruction.

6 With the exception of this study and the
7 implants in the SPS in which there was replacement and
8 follow-up information, the sponsor has not collected
9 safety and effectiveness information on revision
10 patients.

11 You'll be asking the panel questions to
12 discuss the revision indication.

13 Before I discuss the SPS in detail, I
14 would like to show you the implant style studied in
15 Mentor Corporation's clinical studies, as well as
16 those not studied for which the sponsor is seeking
17 approval.

18 Note that the sponsor is no longer
19 manufacturing implants with an oval shape or with a
20 leaf valve. The implants with a contour profile shape
21 have a greater contouring than those that are
22 contoured, and the major difference here is the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 partially textured devices. These implants appear to
2 be unique from all others in that they are textured
3 only on the posterior surface. The texturing is more
4 pronounced. The posterior textured layer is an
5 additional vulcanized layer, and this layer is made
6 from a different, softer silicone.

7 The sponsor has been asked to clarify how
8 this new texturing differs from the type of texturing
9 in their clinical studies and to explain whether and
10 how the clinical performance can be inferred from this
11 new texturing method.

12 The saline prospective study was initiated
13 in 1993 and approved in 1995 after augmentation and
14 summary construction patient enrollment. The study is
15 a prospective, open label, multi-center study with
16 three years of total follow-up for patients seeking
17 primary augmentation and primary reconstruction.

18 Safety was based on local complications,
19 and effectiveness was based on breast dimension
20 changes, patient satisfaction, and quality of life
21 measures.

22 The sponsor collected lactation and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

1 reproduction history at baseline and connective tissue
2 disease, or CTD, symptomatology and breast conditions
3 at baseline and at follow-up.

4 A sample size of 1,200 to 1,500 patients
5 with at least 20 percent of these as reconstruction
6 was proposed to estimate the 95 percent confidence
7 interval precision for complications.

8 The patient disposition at three years is
9 shown here. Of the 1,265 augmentation patients
10 enrolled, approximately 76 percent completed their
11 three year visit. For reconstruction, the completion
12 rate was 66 percent at three years.

13 Of the patients who were withdrawn, the
14 majority for augmentation were lost to follow-up. For
15 the majority of reconstruction patients who were
16 withdrawn, the majority were explanted. The 15
17 patient deaths reported in this study were not implant
18 related. For the 49 augmentation and 75
19 reconstruction patients who underwent explantation,
20 subsequent complications are not included in the
21 Kaplan-Meier complication rates to follow.

22 The three year cumulative Kaplan-Meier

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 rates of first occurred and 95 percent confidence
2 intervals for selected complications are shown here on
3 a per patient basis for the patients in the saline
4 prospective study. You can see that the largest 95
5 confidence intervals are plus or minus three points
6 for augmentation and plus or minus five points for
7 reconstruction.

8 Note that the capsular contracture shown
9 here includes both Baker Grade III or IV and Baker
10 Grade unknown or unreported, and that the category of
11 any complication here includes reoperation.

12 The cumulative risk of a first occurrence
13 of a complication is 43 and 73 percent, respectively,
14 for augmentation and reconstruction. Although only
15 the three year rate is shown here, the cumulative rate
16 of first occurrence of any complication increases over
17 time and has not leveled off by three years of follow-
18 up.

19 The cumulative risk of at least one
20 reoperation for any reason over the three year period
21 is 13 percent for augmentation and 40 percent for
22 reconstruction, and these rates, as well, are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 increasing over time.

2 In general, the complication rates for
3 reconstruction are three times greater than for
4 augmentation.

5 The cumulative rates of first occurrence
6 of leakage/deflation, implant removal, breast pain,
7 wrinkling, and nipple changes, which includes both
8 loss of nipple sensation and intense nipple
9 sensitivity, are shown here as well.

10 The most common types of reoperation
11 procedures performed through three years is shown here
12 based on the number of procedures. Percentages do not
13 sum to 100 because I have omitted infrequently
14 performed procedures from this table.

15 There were 358 and 353 reoperation
16 procedures performed in the augmentation and
17 reconstruction patients, respectively, through three
18 years. For the category of removal with replacement,
19 I combined the following categories reported in the
20 PMA: implant size exchange, secondary augmentation,
21 replacement, and revision.

22 Scar/wound revision includes skin

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 adjustment and wound management. Capsule related
2 procedures include capsulotomy and capsulectomy.

3 The most commonly performed procedure for
4 augmentation patients was removal with replacement, 32
5 percent of the procedures performed in these patients.
6 And for reconstruction, it was a capsule procedure, 28
7 percent of the procedures.

8 Not shown here are the nine implant
9 removals without replacement and augmentation and 40
10 in reconstruction patients.

11 The information shown on this slide was
12 provided subsequent to the PMA submission at the
13 agency's request, and it shows the reasons for implant
14 removal through three years on a by implant basis. If
15 an implant was reported to have been removed for
16 multiple reasons, the hierarchy for categorization
17 into this table is shown in the footnote below the
18 table.

19 Cosmesis includes asymmetry, ptosis,
20 wrinkling, and scarring. Of the 136 augmentation and
21 116 reconstruction implants removed over the three
22 years of follow-up, other than a patient request for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 a size or shape change, the single most common reason
2 for augmentation implant removal was due to
3 leakage/deflation.

4 Twenty-three percent of all removed
5 augmentation implants were due to leakage/deflation.
6 Infection and capsular contracture constituted the
7 most common reasons for reconstruction implant
8 removal, each at 26 percent of all reconstruction
9 implant removal.

10 For both augmentation and reconstruction,
11 if you were to take the complications and sum those,
12 you would see that the majority of implants were
13 removed due to a complication rather than due to
14 patient request for a size or shape change.

15 In an effort to characterize the
16 complication rate in revision patients, the sponsor
17 was asked to provide the cumulative Kaplan-Meier first
18 occurrence complication rate on a by implant basis for
19 those implants which were removed and replaced during
20 the study and for which there was follow-up
21 information.

22 This table summarizes this information

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 through three years of follow-up. Because the sample
2 size and follow-up is lower than for the primary
3 implantation in the SPS, the confidence interval for
4 this table are much wider than reported for the
5 primary implantation group.

6 Note that the capsular contracture here,
7 as well as in the other table I showed you, includes
8 both Baker Grade III/IV and Baker grade unknown or
9 unreported, and the any complication category here
10 includes reoperation.

11 The risk of a first occurrence of any
12 complication for this group is similar for these
13 implants compared to primary implantation. However,
14 for the major complications of reoperation, implant
15 removal, capsular contracture, and leakage deflation,
16 the rates are higher than for primary implantation and
17 lower than for primary reconstruction, which is
18 similar to the revision complication rate reported in
19 the LST.

20 For most other complication the rates are
21 similar or lower than for primary implantation.
22 You'll be asked to address a revision indication in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the panel questions.

2 The sponsor performed Cox regression
3 analysis for several patient at implant variables
4 which they showed you and for the complications of
5 infection, capsular contracture, deflation,
6 reoperation, and removal, and selected associations
7 are shown on this slide.

8 There were no associations with infection.
9 Implants with leaf valves had two times higher
10 capsular contracture than those with diaphragm valves.
11 Recall that the sponsor is no longer manufacturing
12 implants with leaf valves.

13 Surgical pocket irrigation with Betadine
14 was associated with a three and a half times greater
15 risk of deflation than without, and implants with
16 SPECTRUM valves were associated with a twofold higher
17 risk of both implant removal and reoperation than
18 those without.

19 With respect to other safety issues, the
20 sponsor collected breast cancer and connective tissue
21 disease information at baseline and at follow-up. Of
22 note, there were two augmentation patients who

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 developed breast cancer during the course of the
2 study, one patient two years and the other patient
3 five months after implantation. At any follow-up
4 visit complaints of fatigue, generalized achiness
5 and/or joint pain were reported infrequently in
6 patients without such a report at baseline.

7 There were six confirmed and 31
8 unconfirmed cases of connective tissue diseases
9 reported over the course of the study. The six
10 confirmed cases are shown here with the indications
11 shown as well.

12 There were two patients with
13 osteoarthritis and one with an undetermined arthritis
14 and one with ankylosis spondylitis in the
15 reconstruction group.

16 In the augmentation patients, there was
17 one patient with systemic lupus erythematosus and
18 rheumatoid arthritis reported during the course of the
19 study.

20 Without a control group of sufficient
21 numbers of similar types of patients followed for the
22 same duration, conclusion regarding the association of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the CTDs shown here with these implants cannot be
2 made.

3 With respect to effectiveness, the
4 augmentation patients and delayed reconstruction
5 patients experienced an increase in breast size.
6 Augmentation patients were generally satisfied, and
7 they experienced small, but statistically significant
8 improvements in one of the ten subscales of the
9 multidimensional body self-relations questionnaire, or
10 MBSRQ, and small but statistically significant
11 improvements in the Tennessee self-concepts scale.

12 Reconstruction patients experienced
13 statistical improvements in the functional living
14 index of cancer scale, or FLIC, and immediate
15 mastectomy patients experienced improvements in the
16 Beck depression inventory.

17 There were no statistical improvements in
18 the MBSRQ or in the Tennessee self-concept scale for
19 reconstruction patients.

20 Recall that the SPS was initiated in 1993
21 prior to FDA approval in 1995. Shortly after FDA
22 approval, the sponsor was informed that continued

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 follow-up beyond three years would be advisable. The
2 sponsor contacted patients who were still
3 participating in the SPS at that time to solicit
4 continued follow-up in the form of a yearly postcard
5 assessing for deflation.

6 Of the 1,045 augmentation patients in the
7 SPS at that time, 519 or 50 percent agreed to the four
8 to ten year follow-up. Of these 519, 362 patients
9 returned postcards, a 70 percent response rate, and in
10 these 362, there were 36 deflations reported or a rate
11 of ten percent.

12 Of the 375 reconstruction patients in the
13 SPS at this time, 186, or 45 percent, agreed to the
14 four to ten year follow-up. Of the 186, 144 returned
15 their postcards, an 86 percent response rate, and
16 deflation was reported in 17, or 12 percent, of these
17 patients.

18 You'll be asked in the panel questions to
19 address the duration and type of follow-up information
20 needed to fully characterize the long term safety of
21 these implants.

22 In summary, the cumulative risk of a first

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 occurrence -- of a first complication is 43 percent
2 for augmentation and 73 percent for reconstruction; is
3 increasing with time; and has not leveled off by three
4 years.

5 Cumulative complication rates of
6 reoperation and removal have not leveled off as well
7 at three years.

8 Although cumulative local complication
9 rates are increasing, the types of local complications
10 are well characterized, and the rates are precisely
11 defined. In augmentation patients, most reoperations
12 are implant removal. For both augmentation and
13 reconstruction, most implants are removed due to a
14 complication rather than due to a patient request for
15 a size or shape change.

16 Breast size benefits were realized for
17 augmentations and quality of life changes were
18 evident, but small. For reconstruction patients,
19 quality of life measures generally improved.

20 You will be asked to discuss these safety
21 and effectiveness issues in the panel questions to
22 follow.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 Thank you.

2 I'd like to now introduce Ms. Phyllis
3 Silverman, who will be discussing the statistical
4 perspective.

5 MS. SILVERMAN: Good afternoon, or perhaps
6 I should say good evening. I'm Phyllis Silverman, the
7 statistical reviewer for the Mentor PMA.

8 The statistical sections of this PMA are
9 well written, comprehensive, and address nearly all of
10 the requests put forth in the draft guidance for
11 breast implants. The sponsor's PMA contains safety
12 and effectiveness data from five studies.

13 Since the saline prospective study is the
14 only one that utilizes the device in question,
15 includes all of the endpoints of interest, and
16 fulfills the recommended follow-up, I consider it to
17 be the primary study, with the others lending various
18 degrees of support.

19 Because of the approximate 50 percent loss
20 to follow-up with the large, simple trial, the ability
21 to draw meaningful conclusions from this trial is
22 limited. Therefore, my comments will focus on the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 saline prospective study.

2 Because there was no control group and
3 therefore no claims of superiority or equivalence, the
4 safety and effectiveness results for this device must
5 be evaluated by way of descriptive statistics.
6 Complication rates, implant survival curves, and
7 effectiveness parameters must be evaluated from a
8 clinical perspective.

9 As a statistician, my role is not to judge
10 the acceptability of these rates, but to evaluate the
11 validity of the data presentation as well as point out
12 any weaknesses in the study design and analysis. I
13 will start with some comments on sample size.

14 Because there were no null and alternative
15 hypotheses for the primary endpoints, hence making
16 statistical power a non-issue, the adequacy of the
17 sample size was determined by the desired precision
18 around the estimates of complication and reoperation
19 rates. The larger the sample size, the smaller the
20 width of the 95 percent confidence intervals which are
21 used to represent the precision.

22 We wanted to insure that the width of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 confidence intervals would be no more than about plus
2 or minus four and a half percent when rates were high,
3 for example, a 40 or 50 percent complication rate, and
4 only about one to two percent when rates were low.

5 The sponsor's enrollment of 1,265
6 augmentation patients and 425 reconstruction patients
7 resulted in a three year accountability sufficient to
8 meet this precision. Therefore, I feel the sample
9 size was adequate.

10 This brings me to the Kaplan-Meier curve.
11 The sponsor used Kaplan-Meier curves to estimate the
12 occurrence of complications and adverse events. This
13 technique allows women who were not followed for the
14 entire three years to contribute information to the,
15 quote, survival curve for the time that they were in
16 the study. They either experienced the event in
17 question or they are, quote, censored at their last
18 follow-up, which means they are dropped from the
19 denominator at that point.

20 I feel that this is the best technique
21 that the sponsor could have used for this type of
22 data. There are, however, three weaknesses with this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 methodology as applied to these particular data. The
2 information collected reflects prevalence and not
3 incidence, and thus, new cases of infection or
4 contracture, for example, could not be distinguished
5 from continuing cases.

6 Therefore, survival curves are based on
7 the time to the first occurrence of each complication
8 and multiple occurrences could not be analyzed.

9 Secondly, because a patient explanted or
10 revised will be censored from the table and not be in
11 the pool to experience other complications, there is
12 the issue of competing risks which can add an
13 uncalculated bias to these rates.

14 And, thirdly, with the exception of
15 deflation, explant, and reoperation, the exact time of
16 onset of a complication could not be known, but would
17 generally have been noted at the next scheduled
18 follow-up.

19 This interval censoring as it is called
20 can add an additional unknown bias to the data.
21 Therefore, the curves are not as exact as if one were
22 measuring an endpoint like mortality in days.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 A summary of the Kaplan-Meier rates for
2 the primary safety endpoints is reiterated on the next
3 slide. It can be seen from the table that the
4 complication rates in the reconstruction group are
5 about three times what they are in the augmentation
6 group. Because of considerable difference between the
7 augmentation and reconstruction cohorts, these rates
8 must be considered separately by indication and
9 evaluated from a risk-benefit perspective.

10 Is a three year explant rate of almost 27
11 percent or a re-op. rate of 40 percent acceptable for
12 re-com. patients? As you can see from the slide, the
13 95 percent confidence intervals were plus or minus two
14 percent or less for the augmentation patients and plus
15 or minus three to five percent for the reconstruction
16 patients. This is consistent with the guidelines.

17 Now, I would like to discuss some possible
18 biases with the data. There are several sources of
19 possible bias with this data. With three year follow-
20 up missing for approximately 25 percent of the cohort,
21 there could be a non-respondent bias in that women who
22 were having problems were more likely to return for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 follow-up than those who were happy with their
2 implants. This would result in an overestimation of
3 complication rate, or there could be the opposite
4 scenario. Patients with complications were not
5 returning for follow-up because they were seeking help
6 elsewhere. This would result in an under estimation
7 of complication rates.

8 A key assumption for the Kaplan-Meier
9 analysis is that the censoring distribution is
10 independent of the survival distribution. What this
11 means in English is that whether or not a patient
12 returns for follow-up should be unrelated to their
13 level of satisfaction with their implants.

14 Since we do not know to what degree this
15 is true or the reasons for patients not returning, we
16 cannot ascertain this bias. We can only acknowledge
17 that there probably is some, and the complication
18 rates must be evaluated with this in mind.

19 Because many of the complications are
20 self-reported, there is likely also to be some recall
21 bias, especially with the reporting of connective
22 tissue disease and the rheumatology screening. This

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 bias could go either way. That is, it could inflate
2 or decrease the rates, depending on how a patient's
3 memory compares to reality.

4 This is a weakness of a study design where
5 follow-ups are infrequent and medical conditions are
6 not always confirmed by a physician.

7 The last bias I wish to discuss is
8 investigator or site bias. With a study design of 153
9 sites, it is virtually impossible to justify pooling
10 on a statistical basis, and the sponsor did not
11 attempt it.

12 Although there is always the possibility
13 for difference in follow-up or results among sites, I
14 feel that any site or investigator bias would probably
15 be minimal, especially compared to some of the other
16 variables that emerged as related to outcome.

17 For example, surgical approach, valve
18 type, and implant shape are significantly associated
19 with contracture, and valve type is also associated
20 with explant and reoperation.

21 Incision size and use of Betadine
22 irrigation was significantly associated with

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 deflation. The sponsor presented an extensive
2 analysis of co-variables, such as valve type, implant
3 shape, laterality, incision size and surface type by
4 use of Cox regression.

5 Surface type was not significantly
6 associated with contracture, the very thing it was
7 meant to reduce.

8 Of the three biases discussed, my belief
9 is that the nonrespondent bias is of most concern, and
10 that the other two are probably minimal. This brings
11 me to effectiveness.

12 The sponsor presented a very thorough
13 analysis of effectiveness by way of descriptive
14 statistics resulting from numerous surveys
15 administered and objective breast measurements.

16 In addition, before and after comparisons
17 of some effectiveness endpoints showed statistically
18 significant changes. However, I question the
19 interpretation of the phrase "statistically
20 significant increase in breast size." It does not
21 appear to mean anything from a clinical perspective.

22 The data must be looked at in the broader

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 sense. The data indicate to me that even though there
2 were some dissatisfied patients, breast implants are
3 overall effective from both a physical and emotional
4 perspective.

5 Conclusion. In summary, I found the data
6 analysis presented in this PMA to be comprehensive.
7 The sponsor's analysis was consistent with the
8 methodologies laid out in the guidance. The
9 complication rates must not be taken as exact, but
10 rather as estimates subject to the biases discussed
11 earlier.

12 I would like to close just by presenting
13 a few more statistics. Because there could be
14 multiple complications per patient, and even
15 correlations between adverse events, for example,
16 contracture and pain, I would like to leave you with
17 the complication free rates at one, two, and three
18 years. These rates are not subject to the problem of
19 competing risk and would be of particular interest to
20 a prospective patient in making an informed decision.

21 Although some complications are more
22 serious than others, the data show the complications

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 are, indeed, frequent. Approximately 57 percent of
2 augmentation patients are complication free at three
3 years, as opposed to only 27 percent of reconstruction
4 patients. If breast implants are deemed acceptable
5 for market, women must be presented with these figures
6 so they can make an informed decision from their own
7 personal risk-benefit perspective.

8 Thank you for your attention.

9 CHAIRMAN WHALEN: Thank you, Ms. Silverman
10 and the entire team.

11 Well, our next order of business would be
12 to have the FDA entertain questions from the panel.
13 We have a sort of unscheduled break that we must take
14 because apparently part of this room is not reserved
15 for this block of time. So we're all going to get a
16 little bit closer, too.

17 So if we would please take a 15 minute
18 break while they resize this room and hopefully that
19 will be sufficient time for them to do what they have
20 to do.

21 (Whereupon, at 5:53 p.m., a recess was
22 taken, to reconvene at 6:15 p.m., in the same place.)

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 E-V-E-N-I-N-G S-E-S-S-I-O-N

2 (6:15 p.m.)

3 CHAIRMAN WHALEN: I would like to ask that
4 at the real conclusion of the FDA presentation that
5 Dr. Berkowitz review for us the FDA questions that are
6 going to be posed to the panel as the next step.

7 DR. BERKOWITZ: Question one, while the
8 sponsor provided no long term clinical data on their
9 implant, fatigue testing and fold flaw testing
10 provides some information in the long term rupture
11 leakage of the implants. Please comment on the
12 sponsor's methodology and results for each of these
13 tests.

14 CHAIRMAN WHALEN: Dr. Berkowitz, if you
15 could just read all of the questions, we're not going
16 to go into the deliberation upon each one just yet.

17 DR. BERKOWITZ: All right. Question two,
18 given the data for augmentation patients in the SPS
19 and other data provided by the sponsor, is there
20 reasonable insurance as defined in 21 CFR 860.7 that
21 the product is both safe and effective for
22 augmentation patients?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Question three says given the data for
2 reconstruction patients in the SPS and the other
3 information provided by the sponsor, is there
4 reasonable assurance as defined in 21 CFR 860.7 that
5 the product is both safe and effective for
6 reconstruction patients.

7 Question four, with the exception of the
8 LST one year follow-up and the implants in SPS in
9 which there was continued follow-up after
10 explantation, the sponsor has not collected safety and
11 effectiveness information for the cohort of revision
12 patients. Yet the sponsor proposes revision as an
13 indication for use. Given that this cohort typically
14 represents at least 30 percent of the patients
15 presenting for breast implantation, please discuss
16 whether there is sufficient safety and effectiveness
17 data to include revision as an indication and whether
18 the sponsor should evaluate the safety and
19 effectiveness for revision patients as a condition of
20 approval.

21 Please also comment on the information
22 that would be useful to collect in a post approval

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 study.

2 Question five, given that the sponsor's
3 data show increasing cumulative rates per year for the
4 majority of complications for both reconstruction and
5 augmentation patients, please comment on whether there
6 is adequate follow-up data to characterize the
7 frequency and types of long term adverse events.

8 Please address the following pertaining to
9 long term adverse events: (a) the minimum duration of
10 follow-up; (b) the type of visit, i.e., active or
11 passive, and (c) which types of complications should
12 be assessed.

13 Question six, the sponsor's SPS study was
14 not designed to provide information in the following
15 long term issues of pertinence to women with implants:
16 one, the interference on the ability of screening
17 mammography to detect tumors in breasts with implants;
18 two, the interference with lactation; and, three, the
19 effects on offspring from women with implants.

20 Please discuss whether the sponsor should
21 evaluate these issues as a condition of approval. If
22 so, please discuss the appropriate methods for

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 addressing these issues.

2 And the last question is: given the
3 heterogeneity of surgical practices and post operative
4 management of mammary implantation, please comment on
5 the important issues which should be included in
6 physician training.

7 CHAIRMAN WHALEN: Thank you.

8 We now, therefore, proceed to the panel
9 discussion and review of FDA's seven questions, and we
10 will start that off by having three of the panel
11 members in specific areas as lead reviewers make
12 comments in their areas of expertise. Those three
13 will be Dr. Li in mechanical testing, Dr. Burkhardt on
14 the clinical study, and Dr. Blumenstein on statistical
15 considerations.

16 First, for mechanical testing, Dr. Li.

17 DR. LI: Thank you.

18 Let me first say as an overall comment
19 that it appears that mechanical failure of this device
20 in the form of leakage or rupture is one of the
21 primary reasons for revision and reoperation, and this
22 is purely a mechanical failure in my view of either

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 some combination of material, design and environment,
2 and it's unfortunate that that's actually the one area
3 that was the most incomplete in your PMA.

4 The FDA has done a nice job summarizing
5 their comments regarding your testing in the
6 deficiency letter I believe you have received, and in
7 general I agree with virtually all of your comments,
8 but let me highlight a couple of them, I think, that
9 are more important perhaps than the others.

10 One is, I think, you need to provide data
11 for all models that you intend to market and not just
12 the ones that you have selected. This goes for all
13 the sets of testing.

14 The other is that most of the data, as I
15 understand it, was done with the material of
16 construction that you refer to as PTC, RTC, silicone.
17 Yet the final products are made with Sytech (phonetic)
18 silicone, and although you make some comparisons of
19 basic material property similarities, I believe this
20 argument is insufficient to merely make a material
21 swap in those raw material properties alone.

22 For instance, one of your own data points

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 suggests that, in fact, Sytech and PTC are not
2 equivalent in terms of if you look at one of your
3 numerous tables that showed the elongation and break
4 strength of Sytech versus PTC on dry heat. The Sytech
5 silicone actually has improved properties of
6 elongation and break strength. However, if you do
7 that same comparison and the components are gamma
8 sterilized, the order is reversed, and in fact, the
9 Sytech is less strong and has less elongation.

10 So certainly I think your claim that the
11 materials are equivalent is not supported.

12 That raises an odd issue. It appears
13 through the literature that you have one particular
14 Model 1600 which apparently from my reading may or may
15 not be gamma sterilized, which is a little confusing
16 to me. I'm not sure when you choose to gamma
17 sterilize it and when you choose not to gamma
18 sterilize it, and this might be important because
19 basically all of the gamma sterilized material
20 properties are significantly less than those that were
21 heat sterilized. **

22 So my questions would be: why do you do

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 this? How do you decide when they're gamma sterilized
2 and when they're not? And more importantly, does the
3 physician at the end of the line know when he buys a
4 1600 component if it's gamma sterilized or not and
5 that there's, in fact, a material property difference?

6 In a more general sense, the thing that I
7 was most taken back is these implants have been around
8 for decades and the leakage and deflation has always
9 been a key complaint or indication for revision, but
10 none of the testing actually directly addressed this
11 particular mode of failure other than your fatigue
12 test, which as you describe it, is a catastrophic
13 failure under extremely high loads and high cycle,
14 which is some indication of leakage, but certainly not
15 any mirror of what happens clinically apparently in
16 these leakage phenomena.

17 So on the fatigue testing, there is a high
18 variation in your results. If I can just quote a
19 couple of your numbers in general, for your Model 2600
20 of the 175 milliliters volume and an 80 pound load,
21 your cycles to failure varied from 3,000-something to
22 over 32,000-something, a factor of ten from the best

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 to the worst.

2 If you take that same model and do the 325
3 milliliter sample at 75 pounds, again you get about a
4 factor of ten difference, from 16,000 for the worst to
5 126,000 for the best.

6 And lastly, if you take the Model 2400 at
7 the same volume and load, it actually fails at 850
8 cycles rather than the thousands mentioned previously.
9 So although your average graph looks very good, if you
10 factor in the actual cycle fatigue for each individual
11 product, there's at least a factor of ten from best to
12 worst for every component, and then the gamma
13 sterilized version, the 1600 which would be expected
14 to have the worst values, is not done.

15 You do do a lifetime survivorship in a
16 couple of different ways. I don't dispute the
17 methodology. However, everybody should be reminded
18 that that safety factor is for that particular test.
19 So if the end use was, in fact, that kind of cyclic
20 high speed loading in your test rig, then that safety
21 factor would, in fact, be appropriate.

22 But I think it's undoubtedly true that

1 that particular mode of failure is not what happened
2 clinically. So I think it's not supported at all that
3 a safety factor is, in fact, carried over to the
4 clinic.

5 A technical detail, I think, that carries
6 through all of this is the variation of properties and
7 final results as a function of the percent of fill.
8 I couldn't see actually on your reports that I saw how
9 you exactly filled each one. Were they filled to the
10 same volume? Were they filled to the same pressure?
11 I wasn't that exactly sure how that was done.

12 And also in subsequent device tests, there
13 might be cases where the worst case scenario is a
14 device that's under inflated, and in another test
15 scenario, the worst case scenario might be where it's
16 slightly over inflated, and I see no addressing of the
17 issue of inflation percentage at all, and this
18 actually might be one of the surgically related
19 phenomena that Dr. Cunningham alluded to.

20 A bothersome thing on the load thing is
21 this. It seems to be a belief by yourself and Dr.
22 Cunningham that how the surgeon puts it in, in fact,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 makes a difference on the outcome. In fact, there's
2 even a general agreement that somehow that provides
3 higher stresses. Maybe there's folds; maybe there
4 isn't; maybe it's puncture; maybe it isn't. But
5 absolutely none of those particular factors are
6 addressed in any of your testing components, so
7 basically remains anecdotal even after 30-plus years
8 of use.

9 The other tests are -- I guess I don't
10 know what to do with them. You do a static and
11 dynamic rupture test, which is either dropping of
12 weights on something or just squeezing it until it
13 breaks, and those are interesting kind of device
14 tests, but I'm not actually sure how clinically
15 relevant either one of those particular tests are
16 unless you're going to tell me car accidents are
17 actually one of the reasons that some of those devices
18 fail.

19 The abrasion test is even more peculiar.
20 I'm not actually sure what the clinical consequence of
21 where is. Are you projecting that the clinical
22 consequence of where is that the device thins and,

1 therefore, is more likely to rupture or are you
2 worried about where and the fact that it creates some
3 kind of particulate debris that goes on to cause some
4 kind of systemic effect?

5 But in either case neither one of those
6 particular issues is addressed, and you also use a
7 Tabor abrasion test for which panel they stick a flat
8 piece of membrane on a device and then rub against it
9 a very roughened surface. In the crudest sense,
10 sometimes it's a piece of sandpaper, and again, it's
11 a relatively crude test. They're only testing a
12 portion of the device, and again, I don't really know
13 what to do with the device, nor am I sure that
14 particular method of where it is, again, clinically
15 relevant.

16 The tear test let's me get into the area
17 of retrieval analysis. Depending on what numbers we
18 see, the number of devices that were deflated varies
19 somewhere from maybe three percent to, you know, some
20 relatively much higher number at the end. So even if
21 it's a five percent deflation rate and you've
22 implanted 500,000 of these, there's 25,000 retrieved

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 deflated devices somewhere, and I only see a report
2 that you provide that looks at ten of these devices
3 that were retrieved for either deflation and/or wear,
4 which is a little confusing because you didn't tell me
5 which of the ten were retrieved for deflation and
6 which were wear, and I didn't know if the ones that
7 you called wear were also deflated.

8 So be that as it may, you did look at
9 where -- tried to assess where the flaws were that
10 caused the leakage, and the short answer is they
11 appear to be everywhere. They may or may not be
12 around folds. Some of your retrieved implants had
13 folds, permanent folds. Others did not.

14 In most cases the cracks that you
15 associated with leakage were nowhere near the folds.
16 There were a couple that were on the folds, but I
17 don't know if this is just a statistical chance that
18 that's where the crack and the fold happened to meet,
19 but certainly you can have folds without leaks, and
20 certainly you could have cracks in areas without
21 folds.

22 So your tear testing was a cruder sense.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 The cracks that you identified in retrieval run the
2 size of hundreds of microns, very small, sometimes
3 even microscopic sizes, but your tear test is
4 relatively gross when we take a big piece of material
5 and you just try to pull it apart.

6 So I'm not quite sure of the relevance of
7 this more macroscopic tear test to the crack
8 initiation or propagation that you've identified as
9 failure modes in the retrieved implants.

10 So, again, you've done a lot of tear
11 testing, but I actually have no idea how to relate
12 that to the clinical situation, and you also tore it
13 only in one direction, and multi-directional tears and
14 then assessments of you might have survived the tear
15 test, but I didn't see you look at the samples as
16 closely as you looked at the retrieved devices to see
17 if, in fact, you created creases or cracks that may
18 not have failed as a tear, but may have caused
19 pinholes or whatever bit enough to cause a leakage.

20 I'm almost done. Bear with me.

21 The next to the last item is this issue of
22 fold flaw. I think it's certainly a reasonable

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 hypothesis that somehow that these permanent folds
2 that end up in this device are somehow related to the
3 failure, but it's an interesting thing.

4 If you take a brand new implant and try to
5 fold it and you just fold it in half and let it go,
6 the fold doesn't stay there. It goes away, but just
7 because it's a piece of nice, resilient rubber, but in
8 these retrieved devices, that fold line is often rigid
9 and hard, which indicates either a chemical and/or
10 structural change in material along that fold. That's
11 why that fold is permanent as opposed to if you take
12 a brand new implant, fold it up, do some kind of
13 fatigue test. The chemistry and the structure of the
14 two folds, I believe, are completely different.

15 So I'm not quite sure that the fold flaw
16 tests, as difficult as they are to run, again, have a
17 clinical relevance.

18 And then an item that I didn't put much
19 weight to until today's discussion is what for lack of
20 a better term let me call reverse diffusion. It seems
21 as if somehow the inside of these bags get infected
22 somehow and microorganisms find their way inside, and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 perhaps a mechanism for that is there's got to be a
2 way for the bacteria to get into the device.

3 One of the earlier speakers even suggested
4 that the valve, in fact, might be two way in the sense
5 that you can let liquid in as well as liquid out, but,
6 again, I see none of your testing that addresses that.

7 So in summary, I would say that the
8 testing you have done actually has been rather
9 extensive as far as the number of samples and effort
10 that you've put in, but unfortunately hasn't really
11 helped me answer the question at all, is that will
12 this implant leak; how often will it leak; where will
13 it leak; and there's this big mystery in my mind of
14 the mechanism of why the implants fail more in
15 reconstruction than they do in the original
16 augmentation.

17 There's got to be a biomechanical reason.
18 This is a mechanical failure, and you ought to be able
19 to define how that happens and develop a test to
20 address that issue.

21 So in the end of it all, I believe that
22 the FDA is correct in saying that the tests are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 incomplete and they should be completed, and
2 additionally I think the reoperation rate due to
3 inflation -- I find if I carry it over to the devices
4 I normally work with are alarmingly high.

5 To have a device fail in the two to three
6 year period mechanically is extremely surprising, and
7 it's also amazing to me that it seems to be tolerated
8 as just something that you just have to live with in
9 these implants, and I don't really see how the design
10 and the materials change or the testing really
11 addresses that issue.

12 Let me stop there.

13 CHAIRMAN WHALEN: Thank you, Dr. Li.

14 I should point out to all the panel
15 members, lead reviewers and the rest of the panel,
16 that this is the appropriate juncture when, if there
17 are any further questions for the sponsor or any
18 specific questions to the FDA presenters that these
19 questions be raised.

20 That being said, there were probably too
21 numerous to count questions, but some of them were
22 rhetorical. Some of them were comments, and some of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 them were questions that perhaps an answer is needed,
2 and so with that preface specifically for you, Dr. Li,
3 among those questions you raise or any others that you
4 have, would you like to direct any question
5 specifically to sponsor or FDA at this juncture?

6 DR. LI: Let's see. Well, I guess I would
7 like to -- well, I'm not exactly sure because the
8 question is kind of broad ended. I'm struggling with
9 asking kind of a non -- all of my specific questions
10 might be kind of trivial, and the big question I'm not
11 sure we can get into.

12 CHAIRMAN WHALEN: Well, there also will be
13 a closing summation, about ten minutes for each, the
14 sponsor and the FDA, that we'll get to eventually
15 where anything that has been raised during this
16 discussion can be addressed, although no new data will
17 be raised.

18 DR. LI: Let me ask one general question
19 then.

20 MR. PURKAIT: (Inaudible.)

21 CHAIRMAN WHALEN: Well, the point I was
22 making is if I were to say to you right now could you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 answer each of his questions, we'd all have to get
2 cots because it might take the night.

3 DR. LI: Well, let me ask you one general
4 question then. Do you believe that with all the
5 testing that you have provided that you can a priori
6 determine what the leakage will be? Because you
7 clearly have some idea of things that you think are
8 important, for instance, surgical placement, just to
9 bring up another issue, where the size of the incision
10 that you at least feel anecdotally are related to the
11 performance of the device.

12 Yet I was frustrated by that because you
13 don't have a hard number, biomechanical data that
14 says, you know, when you make the incision this small,
15 the force goes up 30 percent and the stress goes up
16 this high and this leads to this and this leads to
17 that.

18 I really kind of -- I don't see that
19 particular sequential kind of argument that we
20 normally apply to a device failure and every other
21 medical device that have been in being applied to this
22 particular device.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 So what makes you -- what gives you the
2 confidence or the belief that if you change, for
3 instance, as an example, from PTC to Sytech or from
4 something to partially textured, which you have not
5 clinically proven; what leads you factually to
6 believe, other than your personal belief, but in terms
7 of data that would you believe that the leakage rate
8 is going to be the same, smaller or bigger than your
9 previous device?

10 MR. PURKAIT: I think -- let me go try to
11 see. You have about -- I don't know -- 15 different
12 areas that you have questioned, which is quite
13 interesting because you get me going for the next two
14 hours I probably can do that, but let me see if I
15 could summarize and try to answer those, the simple
16 question first and then go to the complex one.

17 You asked a question about the differences
18 between the sterilization of dry heat versus gamma.
19 Right now we all do dry sterilization. So one issue
20 about gamma is out right now.

21 Now, the question was: why did you do the
22 gamma nd dry heat and how would one --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. LI: If it's out, you don't have to
2 answer the question.

3 MR. PURKAIT: Okay. Good. So that was
4 simple.

5 DR. LI: It's in the application though,
6 as I understand it, right?

7 MR. PURKAIT: The second question you had
8 about the PTC active versus --

9 DR. LI: Just to clarify that, I raised
10 the issue of that particular one because it was
11 highlighted and takes up many pages in your PMA
12 application. So you're now withdrawing that
13 particular --

14 MR. PURKAIT: No withdrawing. You see any
15 manufacturer operation always have an optional
16 sterilization procedure because you cannot rely upon
17 one particular type of sterilization. We do qualify
18 both dry heat and gamma sterilization, and as we have
19 established the process, validation of the dry heat,
20 we have converted all of the sterilization to dry
21 heat.

22 In case that there be a need in the future

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 we probably will do gamma sterilization, but this is
2 not a withdrawal of the gamma sterilization.

3 DR. LI: Well if you're going to do it at
4 all, then I think you need to answer the questions.

5 MR. PURKAIT: The gamma sterilization and
6 the dry heat sterilization, we have compared the data,
7 and that has been submitted in the PMA, and if you
8 look carefully in the PMA, you will see that the gamma
9 sterilization does reduce to some extent the
10 mechanical properties in comparison to dry heat.

11 However, the range of the properties, what
12 we see by the gamma, is far superior to the expected
13 results that we believe we are going by, such as the
14 ASTM standards. For example, if you have a 350
15 percent elongation, that's what we kind of maintain.
16 Our product shows consistently over 350 percent
17 elongation.

18 You might question, well, what does this
19 350 percent elongation mean, you know. Remember in a
20 body when you put this thing in a cavity, we have a
21 liquid elastomer which goes up to 700 percent
22 elasticity. Now, here we are testing for 350 percent.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 In a body probably under all loads we know that we can
2 measure, this probably wouldn't extend more than 20
3 percent to 30 percent.

4 So, you know, to look into the proper
5 perspective, the elastomer elasticity in this case is
6 much superior or much more higher than is reported in
7 the body.

8 DR. LI: Well, let me stop you here. This
9 is my general problem, I guess, with the testing.
10 First of all, the ASTM methods for those of you who
11 aren't into ASTM methods are proposed standard ways of
12 doing tests, and there's not an ASTM method yet that
13 I've read that doesn't have a disclaimer in there that
14 if you meet these standards, it has nothing to do with
15 projected clinical performance.

16 So if you followed the ASTM standard,
17 you're basically telling me you're following a
18 standard test, but as the ASTM itself says, it's not
19 performance related.

20 And this is the general problem or concern
21 I've got with all of your testing.

22 MR. PURKAIT: I understand that.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 DR. LI: It performs to some standard, but
2 I can't make the connection to the clinical case where
3 maybe 19 percent are failing by leakage. Your tests
4 suggest that.

5 MR. PURKAIT: I agree with you. The ASTM
6 is not our Bible; that we follow ASTM, therefore,
7 everything is good. The ASTM is a standard that's
8 accepted across the country, across all product lines,
9 all devices that exist. So we do follow their
10 standard.

11 At the same time, we supplement many other
12 tests to show that not only do we meet ASTM. We also
13 have other tests to show that we go beyond that. So
14 ASTM is not the only study that we do to say that our
15 product --

16 DR. LI: I understand that, but at the end
17 of the day, you still have a 19-plus percent, up to a
18 19 or more percent leakage rate.

19 MR. PURKAIT: Now you talk to my heart.
20 If you look at the in vitro versus in vivo situation,
21 unfortunately we are at a loss to exactly simulate
22 what happened in the body process in vitro.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 For example, you have seen in our data
2 today that if somebody used Betadine, if somebody used
3 different incision size, if somebody used a bilateral,
4 if somebody used a different valve type, there are
5 some clinical indications that will occur. A clinical
6 problem will occur that it cannot replicate every time
7 in vitro.

8 I'm not saying that we're not going to try
9 for that. However, at this point in time we took into
10 consideration the best we can, and we continue to
11 study that all the time.

12 DR. LI: I'm not saying that you didn't do
13 your best. What I'm saying is that there's a
14 disconnect for me between the data you generated and
15 the prediction of in vivo performance. I mean I'm not
16 disputing the hard work that you put into it or you
17 sincerity in doing --

18 MR. PURKAIT: Well, I'm not going to argue
19 on that. I'm trying to make the point that some of
20 the test conditions, what we have used, does have some
21 real meaning behind it.

22 That includes one of the areas you also

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 have addressed the fact that our fold flaw or the
2 leakage things -- let me address the leakage one.

3 To understand the deflation or leakage in
4 a body, we considered that there are three ways,
5 primarily three ways that it can fail: rupture of the
6 shell, the valve failure, or maybe fold flaw, or maybe
7 other reasons in the clinical.

8 The rupture in the shell, we try five
9 testings, such as fatigue, the static rupture, the
10 static impact, and so forth. For the valve
11 competence, we have three different tests for valves.
12 We have valve burst test. We have valve special test.
13 We have valve -- the flow properties test.

14 For the fold flaw, we believe that the
15 fold flaw test is very unpredictable because you
16 cannot predict where and when, how the fold will be
17 formed. That has a lot to do with how it has been
18 implanted and how these devices have been put there.

19 And the other question I think previously
20 asked about, the special inferences in different
21 locations, initially all of these implants are folded
22 and put inside the cavity and then has been placed and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 then been inflated.

2 So the pressure generated there, whether
3 you put it in submuscular or subglandular, is really
4 determined by what location you're putting and,
5 regardless whatever the pressure is, we always test
6 for the worst case condition. We always test for the
7 extreme conditions. So, therefore, we believe that
8 even if it is within the range, it will maintain the
9 properties.

10 DR. LI: Maybe one last -- I'll try to
11 make it a last response.

12 If I were to take your data at completely
13 face value, I think I would walk away with the
14 impression that this device is near bulletproof. I
15 mean you have to go in your testimony to rather
16 extreme conditions to get rupture, fatigue. Nothing
17 happens in the fold flaw.

18 Now, if I look at all of your tests, it
19 actually looks extremely good, except for the --

20 MR. PURKAIT: It is good.

21 DR. LI: But you get 19 percent failure
22 rate, and over 40 percent of them are revised.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 MR. PURKAIT: But over 40 percent of
2 those --

3 DR. LI: Which is an enormously high
4 number.

5 MR. PURKAIT: If you break it down into
6 cosmetic versus non-cosmetic --

7 DR. LI: Okay. So 20 percent. Take half.

8 MR. PURKAIT: Okay.

9 DR. LI: It's still a high number in three
10 years. So that's the disconnect that I'm going after,
11 right?

12 I mean your data looks excellent, right?
13 I mean if I just looked at your data in and of itself,
14 I would say from a materials and design standpoint it
15 looks excellent. Okay? In the absence of any
16 clinical data, you know, I'd probably have a
17 completely different view, right? But the problem is
18 I do have a clinical --

19 CHAIRMAN WHALEN: Just to focus upon this
20 for a moment --

21 DR. LI: Yeah.

22 CHAIRMAN WHALEN: -- in the interest of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the flow of things, I don't think we're talking about
2 mechanical testing per se anymore. We are talking
3 about a highly clinically significant issue.

4 DR. LI: Well, I think they have to be
5 linked to be meaningful.

6 CHAIRMAN WHALEN: Indeed, and we're going
7 to proceed to other clinical issues. I guess
8 refocusing, is there something specific that you would
9 like to inquire about in terms of other mechanical
10 testing that could have or should have been done?

11 DR. LI: Well, I guess, for instance, why
12 haven't you looked at the effect of percent fill on
13 the results as an example?

14 MR. PURKAIT: How do you mean?

15 DR. LI: In other words, doing a fold --
16 pick a test. Pick a dynamic test, fatigue, fold,
17 whatever.

18 MR. PURKAIT: Let's say you mention
19 about --

20 DR. LI: And then do an under fill -- do
21 a 20 percent under filled, 20 percent over filled, and
22 then --

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 MR. PURKAIT: Well, we don't suggest
2 anybody to under fill. In our label copy, we clearly
3 say, please, please, please don't.

4 DR. LI: So are you going to say that
5 never happens?

6 MR. PURKAIT: I don't know, but we say
7 that's what is supposed not to happen. We can't
8 control this. We test in the nominal volume. That's
9 the way they come in there.

10 You mentioned something. Just to clarify,
11 the 2600 model, 2400 model, they were failed about
12 3,332 cycles, 16,118, those failures if you're looking
13 back in the data was intentionally done to understand
14 at what pressure and at what load that we can make a
15 failure so that we can make a model. Those does not
16 show a premature failure of those materials or those
17 devices.

18 DR. LI: That wasn't clear in the
19 application then in that case.

20 DR. CUNNINGHAM: If I could address your
21 under fill issue from a clinical point of view, there
22 is a body of information within plastic surgery

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 medical literature which I think would have almost all
2 plastic surgeons feel very strongly that they should
3 not under fill these devices because of the
4 predilection for possible folds, fold flaws.

5 So I think most plastic surgeons know that
6 there is a small range of fill which these devices are
7 meant to perform within, and in fact, the way you
8 determine what size implant you want to use is fairly
9 precise because there are sizers which are connected
10 with a tube and have the same footprint and dimensions
11 as the implant which you place in, fill with saline
12 until you reach the look or appearance that you think
13 is appropriate for that patient and conforms with the
14 discussions that you've had with that patient, and
15 that gives you the amount of saline and allows you to
16 choose which device to use.

17 So that surgeons are able to choose within
18 that narrow fill range exactly which device is
19 appropriate for the patient. So I think in clinical
20 practice the real world is in this case perhaps more
21 precise than the kind of range of testing that you
22 might think we would have to do.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 DR. LI: Well, in their own literature, in
2 their PMA document, they said that -- I don't remember
3 the exact phrase -- but that 20 percent over fill
4 would be something that they would consider as an
5 upper end of what might happen clinically.

6 So at 20 percent over fill, why weren't
7 things tested at 20 percent over fill if that's what
8 they stated as an over fill potential level?

9 MR. PURKAIT: I'm not sure I recall that,
10 but maybe you are referring to the SPECTRUM product
11 where you can go for adjustment purposes as we allow
12 in our particular product.

13 DR. LI: Right.

14 MR. PURKAIT: Yeah, that is okay for the
15 SPECTRUM, not for the regular fixed volume one.

16 DR. LI: Well, the other medical -- again,
17 I hate to harp upon other medical devices, but
18 typically in these there's a zone where you want the
19 surgeon to be or the physician to be in the
20 implantation of this device, but for reasons either by
21 skill or by necessity, the person's anatomy or
22 something beyond the physician's control, they can't

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 always hit that exact target, and they might have to
2 do something, make a decision to go out of that
3 extreme.

4 It happens. Right? It happens. You know
5 you don't want it to happen. Most of the time it
6 doesn't happen, but it happens. Right?

7 And you're faced with a case here where
8 you have a very high number of failures and an
9 alarmingly little analysis of those retrieved devices,
10 right? And with the absence of that information of
11 how the device actually deflates, I don't really know
12 how you can discount any possible mechanism.

13 CHAIRMAN WHALEN: Any other questions, Dr.
14 Li?

15 DR. LI: I think I'm done.

16 DR. ALLEN: Will I have an opportunity to
17 respond?

18 CHAIRMAN WHALEN: Actually at the
19 summation period if you wish, yes, but no.

20 And I'll get to your question in a moment.
21 We're going to do the lead reviewers first and then
22 we're going to go to general questions from the panel,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 but we're going to -- okay. Dr. Burkhardt is here for
2 the clinical study.

3 DR. BURKHARDT: That was all just on
4 question one?

5 CHAIRMAN WHALEN: Actually we've not done
6 question one yet.

7 (Laughter.)

8 DR. BURKHARDT: How do you want me to do
9 this? Do you want me just to address question one?

10 CHAIRMAN WHALEN: No, no, no. This is
11 just a general clinical study review, not with
12 reference to any of the particular seven questions.

13 DR. BURKHARDT: Oh.

14 CHAIRMAN WHALEN: An overview.

15 DR. BURKHARDT: Okay. My thoughts on the
16 clinical trial are, first, that -- my thoughts on this
17 whole thing, I think is what you're driving at.

18 CHAIRMAN WHALEN: Yes, sir.

19 DR. BURKHARDT: The question of systemic
20 illness and second generation problems, the reports
21 are experiential. We don't have any scientific data
22 on that, and what we have doesn't support it, and I

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 think they have to be disregarded.

2 What we have to worry about is local
3 complications, and I'm not enough of an engineer to
4 understand what happened, why this incidence was as
5 high as it was, but my understanding is that if you
6 eliminate needle punctures, which are there, and valve
7 failures, that the major problem here is fold flaw
8 failure, and my understanding of that -- you correct
9 me if I'm incorrect -- is that that occurs because of
10 abrasion, internal abrasion at the end of a fold.

11 In other words, it's not material fatigue.
12 Am I correct about that? Because that's what I've
13 been told.

14 MR. PURKAIT: To some extent fold flaw,
15 that what we know of, we can speak for, could occur,
16 could fail in at least three or four different ways.
17 One of the mechanisms may be abrasion. The other
18 would be, as I mentioned, a localized stress
19 concentration. It's a creep factor.

20 If you have a two fold (phonetic) and some
21 weight has been there for a long time, the material
22 crimps, and that might give a pinhole, and it will

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 fail.

2 And the other probably would be a multiple
3 fold that causes abrasion with different surfaces, and
4 that might be a problem there.

5 DR. BURKHARDT: Thank you.

6 But the underlying problem is that they're
7 going to fold. If you take an oval, three dimensional
8 or round three dimensional thing with an oval cross-
9 section and you stand it on its end, which is what we
10 see in these patients when they stand up, for
11 instance, they're going to get folds in them. The
12 material is not perfectly elastic, and there's no way
13 that I know of that you can get around that.

14 DR. LI: But, Dr. Burkhardt, perhaps to
15 clear it up, there's some pictures that were difficult
16 to see in black and white, but the FDA gave me the
17 color versions. These are ten retrieved devices that
18 Mentor supplied a photograph. They did a very nice
19 job on these particular set of ten.

20 But if you look through these photographs
21 of ten, the cracks which may have been the leak are
22 delineated in black, and on Figure 6 there, you

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 actually see running across the horizon a white line
2 that's a fold flaw, and you notice the two cracks that
3 probably caused the leak are nowhere near it.

4 And if you look at the ten all the way
5 through, more often than not the cracks and the pin
6 holes that they identified were nowhere near the fold.

7 Now, in two cases, I think they were, but
8 in the other eight they were not. So my point is this
9 generation of these small cracks is not mimicked in
10 any of their testing that I've seen.

11 DR. BURKHARDT: Well, these are implants
12 that are removed, right? And --

13 DR. LI: For leakage.

14 DR. BURKHARDT: Yeah, for leakage, and
15 what you're saying is that there's no permanent fold
16 there that you can see now.

17 DR. LI: Well, on Figure 6 there's kind of
18 a ghostly white line that runs across the horizon. I
19 think it's Figure 6, the page I handed you. Now, some
20 of them you'll see like a white line that goes across
21 a gray background. That white line is the permanent
22 fold, and the black lines that they've delineated are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 where the cracks and pin holes are, and you'll see
2 that oftentimes, most times that white line -- and
3 sometimes it runs in all different directions -- is
4 nowhere near the black lines that they identify as the
5 source of the crack or the leak.

6 So in other words, this fold flaw thing
7 still, after 30 years, may or may not be the reason
8 these things leak based on the data that they've
9 supplied.

10 DR. BURKHARDT: Well, I see some big black
11 lines on here, but these were not cracks in the
12 implant. Am I correct?

13 DR. LI: Yes. Yes, the short ones are
14 all --

15 DR. BURKHARDT: Well, now let's find out
16 about that before we decide that.

17 DR. CUNNINGHAM: I'm not sure exactly what
18 picture --

19 DR. BURKHARDT: You've got big, long black
20 lines in these pictures. Were those the cracks?

21 DR. LI: The long ones are wear. The
22 short ones are cracks. If you see a really long black

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 line -- and correct me if I'm wrong -- but if I know
2 your nomenclature, there is some implants that have
3 very long black lines in them. Those are wear lines.

4 Other ones are very short black lines.
5 The short black lines are cracks and holes.

6 DR. BURKHARDT: Okay, and so I guess I
7 don't understand the point.

8 DR. LI: Well, the point is it is not so
9 simple as you would like to make it, that if there's
10 a fold it's bad, and if there's no fold, it's good,
11 right? Because some of those have fold --

12 DR. BURKHARDT: Why not?

13 DR. LI: Because you've got the data right
14 there. If folds and leaks were directly associated,
15 the black line, the small black marks to indicate
16 holes should be right on all of those ghostly white
17 lines that are folds, and they are not.

18 DR. BURKHARDT: I think I'm not
19 sophisticated enough to follow that line of reasoning.
20 If the black lines are wear lines, which I would
21 interpret as being the convex lines of the wear --

22 DR. LI: So, for instance, we're looking

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 at Figure 6 in your report. This line here, this off
2 white, that's a fold line as I understand it, and
3 these two lines are where the cracks are.

4 If the fold line was a source of cracks,
5 these crack lines should be right on that line, and
6 they are not.

7 DR. BURKHARDT: The black lines are
8 described as wear lines, not cracks.

9 DR. LI: "Location of cracks in relation
10 to the wear lines." Right. There's cracks and wear
11 lines, right.

12 DR. BURKHARDT: There's a red line for a
13 crack and black lines for the wear lines.

14 DR. LI: The problem is on this one
15 unless --

16 CHAIRMAN WHALEN: In the microphone,
17 please.

18 DR. LI: I guess the problem on this one
19 is unless Mentor can say that the wear lines are not
20 leak lines, I assume that when you had a wear line it
21 may or may not be the source of a leak; is that
22 correct?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 So some of those wear lines may leak and
2 some of those wear lines may not; is that correct?

3 MR. PURKAIT: That's correct.

4 DR. LI: Okay.

5 MR. PURKAIT: But I'm not sure which
6 picture and what you are talking there because it's
7 hard for me to really --

8 DR. LI: I understand. I'm just trying to
9 fill in for Dr. Burkhardt how I was looking at those
10 photos.

11 DR. BURKHARDT: Okay. Where were we?
12 Nobody under fills these implants with any knowledge.
13 You can't really legislate physician behavior, but --

14 DR. LI: But you could design for it.

15 DR. BURKHARDT: Well, maybe you can design
16 for it, and again, I don't know. Maybe these could be
17 improved, and that's your jurisdiction more than mine.

18 The leakage rate that is generally
19 reported from fold flaw failure is in the five to ten
20 percent range, and I don't remember the slides that
21 you showed that brought it up to 19 or why it got to
22 be so high.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Questions were raised about the texture to
2 implant. I think it would be worthwhile for everybody
3 to know why these implants were textured and what gave
4 rise to the origin of the textured implant because it
5 didn't have anything to do with easy insertion through
6 a small incision.

7 A number of years ago a polyurethane
8 covered implant came out from another company that was
9 called the MIM (phonetic). It was widely accepted in
10 plastic surgery and short term had a very small
11 incidence of hardness or capsular contracture. Long
12 term there's some question about what happened. Most
13 of us think they all got hard.

14 At that time the companies that were
15 producing the silicone implants, the silicone shell
16 implants could not use a polyurethane covered implant
17 because it was patented. So they had to try to do
18 something to compete economically with the success of
19 the polyurethane covered implant, and the response was
20 to texturize the surface of the silicone implant.

21 Most of us who are in the field thought
22 that that wouldn't accomplish anything at all, but

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com