

FOOD AND DRUG ADMINISTRATION

+ + + + +

MEDICAL DEVICES ADVISORY COMMITTEE

+ + + + +

EAR, NOSE AND THROAT DEVICES PANEL

+ + + + +

MEETING

+ + + + +

Friday, July 21, 2000

+ + + + +

This transcript has not
been edited and FDA
makes no representation
regarding its accuracy

The panel met 9:15 a.m. in the Goshen Room
of the Gaithersburg Holiday Inn, Two Montgomery
Village Avenue, Gaithersburg, Maryland, Dr. Carl A.
Patow, Chairman, presiding.

PRESENT:

CARL A. PATOW, M.D., Chairman

ALEXA CANADY, M.D., Temporary Voting Member

WILLIAM H. DUFFELL, JR., Ph.D., Industry Rep.

HOWARD FRANCIS, M.D., Temporary Voting Member

A. JULIANNA GULYA, M.D., Temporary Voting Member

LINDA J. HOOD, Ph.D., Temporary Voting Member

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

PRESENT (Continued):

ANJUM KAHN, M.D., Voting Member

PAUL R. KILENY, Ph.D., Voting Member

ROSS J. ROESER, Ph.D., Temporary Voting Member

CLOUGH SHELTON, M.D., Voting Member

GAYLE E. WOODSON, M.D., Voting Member

SARA M. THORNTON, Executive Secretary

FDA REPRESENTATIVES:

NANCY C. BROGDON

TERI M. CYGNAROWICZ

I. SIDNEY JAFFEE, M.D.

JAMES K. KANE, Ph.D., CCC-A

VASANT MALSHET, Ph.D.

ALFRED MONTGOMERY, M.D.

MORRIS WAXLER, Ph.D.

SPONSOR REPRESENTATIVES:

PATTI L. ARNDT, CCC-A

DERALD E. BRACKMANN, M.D. FACS

KIARA A. EBINGER, M.S.

WILLIAM E. HITSSELBERGER, M.D.

MARTYN L. HYDE, Ph.D.

STEVE OTTO, M.A.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

SPONSOR REPRESENTATIVES (Continued):

ROBERT SHANNON, Ph.D.

STEVE STALLER, Ph.D.

CHRIS VAN DEN HONERT, Ph.D.

RONALD E. WEST, M.B.A.

PUBLIC SPEAKERS:

HENRY J. ILECKI, Ph.D.

DONNA McLAUGHLIN

GAIL UMPHREY

C-O-N-T-E-N-T-S

PAGE

Introductory Remarks, Sara Thornton	5
Conflict of Interest Statement	11
Public Comment:	
Gail Umphrey	14
Donna McLaughlin	20
Henry J. Ilecki, Ph.D.	27
Division Update, Nancy Brogdon	34
Branch Update, Morris Waxler, Ph.D.	38
Presentation of PMA P000015	39
Sponsor's Presentation:	
Ron West, M.B.A.	39
Patti L. Arndt, CCC-A	43, 88
Martyn L. Hyde, Ph.D.	51
Derald E. Brackmann, M.D., FACS	57
William E. Hitselberger, M.D.	67
Kiara A. Ebinger, M.S.	70
FDA Presentation:	
Morris Waxler, Ph.D.	129
Teri Cygnarowicz	129
I. Sidney Jaffee, M.D.	131
James K. Kane, Ph.D.	134
Additional Comments by the Sponsor, Patti L. Arndt	144
Committee Deliberations	146
Closing Comments by the Sponsor, Ron West	209

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

P-R-O-C-E-E-D-I-N-G-S

(9:15 a.m.)

CHAIRMAN PATOW: I'd like to welcome you all and up front say how much I appreciate all of the hard work that the panelists have put in in reviewing the materials for today and also industry and FDA for their efforts in preparation for today.

We are very much going to try to stay on schedule, and so I will be mentioning before each speaker the time allotted for those talks so that we can keep on schedule. On the other hand, we want to make sure that we get all of the appropriate data so that we can make some good decisions.

At this time I'd like to introduce Sara Thornton, Executive Secretary, for introductory remarks.

MS. THORNTON: Good morning and welcome to the second day of the Ear, Nose, and Throat Devices Panel meeting.

Before we proceed with today's agenda, I just have a few announcements. I'd like to remind everyone that there's a sign-in sheet for attendance

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 record in the registration area just outside the
2 meeting room. All handouts for today's meeting are
3 available at the registration table.

4 If you have any messages for the panel
5 members and FDA participants, information or special
6 needs, you should direct that through Ms. Ann Marie
7 Williams or Ms. Carol Coy, who are available in the
8 registration area.

9 If you should need an assistive listening
10 device, please see Ms. Williams or Ms. Coy.

11 The phone number for calls to the meeting
12 area is (301) 948-8900.

13 Lunch for the panel will be in the
14 farthest area of the restaurant. We have tables
15 reserved back there during lunchtime for the panel and
16 for the FDA.

17 In consideration of the panel, the sponsor
18 and the agency, we ask that those of you with cell
19 phones and pagers either turn them off or put them on
20 vibration mode while in this room.

21 Lastly, will ^{**}all meeting participants
22 please speak into the microphone and give your name

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 clearly so that the transcriber will have an accurate
2 recording of your comments? The panel, I believe,
3 need do that just a few times, and then they will be
4 clear on who is responsible for what comment.

5 Now, at this time before I ask the panel
6 to introduce themselves, I'd like to extend a special
7 welcome and introduce to the public the panel and the
8 FDA staff, two panel consultant members who are with
9 us for the first time today.

10 Dr. Howard Francis. Dr. Francis is an
11 Assistant Professor with the Division of Neurotology
12 and Skull Base Surgery, Department of Otolaryngology-
13 Head and Neck Surgery at the Johns Hopkins University
14 School of Medicine in Baltimore.

15 And Dr. Linda Hood is a Professor at the
16 Kresge Hearing Research Laboratory of the South
17 Department of Otorhinolaryngology, Louisiana State
18 University Health Science Center in New Orleans,
19 Louisiana.

20 I'd like you to know that the panel
21 consumer representative "has been prevented from
22 attending this meeting due to illness.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 I would like to extend a special welcome
2 as I introduce to you Dr. Alexa Canady. Dr. Canady is
3 Chief of Neurosurgery and Peter Schotanus Endowed
4 Professor in pediatric neurosurgery at the Children's
5 Hospital of Michigan and Vice Chairman of Neurosurgery
6 at Wayne State University, Detroit. Dr. Canady has
7 been a consultant on the Neurological Devices Panel of
8 FDA's Medical Devices Advisory Committee since 1994
9 and currently serves as the panel chair.

10 Dr. Canady, we realize that your
11 preparation and attendance was above and beyond the
12 call of duty, and we are grateful for your willingness
13 to participate with us today.

14 I'll now ask the others at the panel table
15 to introduce themselves, starting with Dr. Duffell.

16 DR. DUFFELL: I'm Bill Duffell. I'm the
17 industry rep. I am Vice President of Clinical
18 Research and Regulatory Affairs for Cyberonics, Inc.
19 in Houston, Texas, and beginning August 1, I'll be
20 with Gambro BCT in a similar capacity in Lakewood,
21 Colorado.

22 DR. GULYA: I am Julie Gulya, Clinical

1 Professor of otolaryngology, head and neck surgery at
2 the George Washington University, and I'm also
3 Director of the clinical trials program at the NIDCD.

4 DR. SHELTON: Clough Shelton. I'm a
5 professor at University of Utah and an otologist.

6 DR. KAHN: I'm Angie Kahn who's in private
7 practice in otolaryngology in Silver Spring, Maryland.
8 I'm affiliated with George Washington University and
9 Uniformed Services Health Sciences as Associate
10 Clinical Professor.

11 CHAIRMAN PATOW: I'm Carl Patow. I'm the
12 Executive Director for the Health Partners Institute
13 for Medical Education. It's a large, nonprofit
14 educational institution associated with a managed care
15 organization in Minneapolis.

16 I'm also on the clinical faculty of the
17 University of Minnesota.

18 DR. KILENY: I'm Paul Kileny, Professor of
19 Otolaryngology at the University of Michigan Medical
20 School and Director of Otology at the University of
21 Michigan Health System. **

22 DR. WOODSON: I'm Gayle Woodson. I'm

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Professor of Otolaryngology at the University of
2 Tennessee, Memphis, and in about a month I'll be
3 moving to Gainesville where I'll be Professor at the
4 University of Florida, Gainesville.

5 DR. ROESER: I'm Ross Roeser. I'm a
6 Professor at the University of Texas, Dallas, in the
7 Program and Communications Sciences and Disorders.
8 I'm also the Director of the Callier Center for
9 Communication Disorders, which is a component of the
10 University of Texas, Dallas, a large center in Dallas
11 specializing in communications disorders, and I'm a
12 Clinical Professor at the University of Texas
13 Southwestern Medical Center in the Department of
14 Otorhinolaryngology-Head and Neck Surgery.

15 MS. BROGDON: I'm Nancy Brogdon. I'm not
16 a member of the panel. I'm FDA's liaison to the
17 panel. I'm the Acting Director of the Division of
18 Ophthalmic and ENT Devices.

19 MS. THORNTON: Thank you very much, panel.

20 On behalf of the FDA, I wish to extend our
21 sincere appreciation to the ^{**}panel for the time they've
22 taken from their busy schedules to prepare for and

•
NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 participate in this meeting today.

2 Thank you, Dr. Patow.

3 I'd like to now read the conflict of
4 interest statement for this meeting.

5 The following announcement addresses
6 conflict of interest issues associated with this
7 meeting and is made a part of the record to preclude
8 even the appearance of impropriety. To determine if
9 any conflict existed, the agency reviewed and
10 submitted an agenda, and all financial interests
11 reported by the committee participants.

12 The conflict of interest statutes prohibit
13 special government employees from participating in
14 matters that could affect their or their employer's
15 financial interests. However, the agency has
16 determined that participation of certain members and
17 consultants, the need for whose services outweigh the
18 potential conflict of interest involved, is in the
19 best interest of the government.

20 We would like to note for the record that
21 the agency took into consideration certain matters
22 regarding Drs. Paul Kileny and Clough Shelton. These

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 panelists reported past interests in firms at issue,
2 but in matters that are not related to today's agenda.

3 Therefore, the agency has determined that
4 they may participate fully in today's deliberations.

5 In the event that the discussions involve
6 any other products or firms not already on the agenda
7 for which an FDA participant has a financial interest,
8 the participant should excuse him or herself from such
9 involvement, and the exclusions will be noted for the
10 record.

11 With respect to all other participants, we
12 ask in the interest of fairness that all persons
13 making statements or presentations disclose any
14 current or previous financial involvement with any
15 firm whose products they may wish to comment upon.

16 I'd like to read the appointment to
17 temporary voting status for today's meeting.

18 Pursuant to the authority granted under
19 the Medical Devices Advisory Committee charter, dated
20 October 27, 1990, and as amended August 18th, 1999, I
21 appoint the following individuals as voting members of
22 the Ear, Nose, and Throat Devices Panel for this

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 meeting on July 21st, 2000: Dr. Howard Francis; Dr.
2 Julianna Gulya; Dr. Linda Hood; Dr. Ross Roeser; Dr.
3 Alexa Canady.

4 For the record, these individuals are
5 special government employees and consultants to this
6 panel or other panels under the Medical Devices
7 Advisory Committee. They have undergone the customary
8 conflict of interest review and have reviewed the
9 materials to be considered at this meeting.

10 Signing for Dr. Feigal, Linda Kahan. Dr.
11 Feigal is Director of the Center for Devices and
12 Radiological Health. This is dated 7/11/2000.

13 Thank you, Dr. Patow.

14 CHAIRMAN PATOW: Thank you.

15 At this point, I'd like to read a charge
16 to the panel regarding confidentiality.

17 I'd like to remind the panel that we're
18 not to discuss any PMAs under consideration with
19 anyone else, including FDA staff and other panel
20 members. For our own protection, we must be very
21 cautious about the perception of bias and conflict of
22 interest that can arise at a public meeting attended

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 by members of industry who may be in market
2 competition with teach other.

3 To that end I would caution you against
4 having extended conversations with individuals who are
5 not on the panel, conversations that might be
6 misinterpreted by others as demonstrating favoritism
7 or bias.

8 At this point, I'd like to go directly
9 then to our open public hearing session. We have
10 three individuals who have asked to speak in the open
11 public hearing session, and I'd ask that each of them
12 limit their comments to ten minutes or less.

13 The first speaker is Gail Umphrey, and I'd
14 ask that each speaker announce, please, their
15 affiliation and also who has paid or supported them to
16 this conference.

17 Thank you.

18 Is Gail Umphrey here today? Yeah.

19 MS. UMPHREY: Okay. My name is Gail
20 Umphrey. This is my husband, Varn.

21 I received the implant in 1994, and I'm
22 very honored to be here and to be participating in the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 usage of the ABI.

2 I just want to say -- give you some facts
3 on what it was like before I had the implant. I had
4 almost full hearing until about the age of 23. I was
5 totally deaf for two years before the implant; slowly
6 lost my hearing from 1976 to 1992.

7 Without hearing for the two years, I tried
8 raising my four children, which was very difficult
9 being totally deaf; very difficult communicating, a
10 lot of frustration. The kids had to -- I did mainly
11 lip reading, all lip reading. I've never had sign
12 language training. So it would be a lot of
13 frustration. If I couldn't understand, they would
14 have to write a note. A lot of frustration.

15 So having this problem, a lot of people,
16 including my family, friends, would leave me out of
17 conversation because it was so difficult. That was
18 sad for me. Like during even holidays I would find
19 myself grabbing a magazine and going in another room
20 and reading a magazine while my family was enjoying
21 the holidays. **

22 I know they didn't intend to do it, but it

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 was difficult. You know, very frustrated to
2 understand on both parts.

3 I felt as though my family, they were
4 taking care of me instead of I was taking care of my
5 family. And during those two years, you know, my
6 children just grew up more or less like without a mom.
7 I couldn't do my part as being a full mom and wife.

8 I just went on day to day surviving and
9 really not living. I would find myself if I was in a
10 supermarket and I would see a friend, an acquaintance
11 that I wouldn't see on a day to day basis and didn't
12 know that I had all of these problems, I would avoid
13 them. I would go the opposite way before they would
14 notice me.

15 Now with the implant I talk to everyone,
16 and I'm not avoiding a situation. Without the
17 implant, during the time I was deaf, I would never ask
18 for help from, let's say, somebody in a store or what
19 have you. Now I do. I have no problem with that.

20 I have absolutely wonderful hearing with
21 this ABI. I use the phone." I'm able to call home and
22 check on the kids, and that was not there for two

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 years, of course.

2 Well, you can say that. Now, you know,
3 we're able to go out with friends and continue like I
4 did before I lost my hearing.

5 Can you help me? I have all of these
6 notes sitting here.

7 Okay. I also was not aware of my
8 surroundings. I think just going out and walking to
9 the mailbox several times, and I mean my kids have
10 even done it during the time of no hearing, the car
11 would drive up the driveway, and of course I'd be
12 walking up the driveway and not even be aware that
13 there was a car behind me. Now there is no problem.

14 I hear birds. I hear cars and everything.

15 I think to sum up all of this, the two
16 years without hearing, it was really the worst time of
17 my life for me and also for my family. I was just
18 surviving. That's all. In fact, at the end of the
19 day I couldn't wait until it was bedtime so that the
20 day was over with.

21 And I look back now and I wonder how I
22 even got through the two years. It sounds scary to me

1 even to think of no sounds, absolutely no sound for
2 two years.

3 With the NF-2, you know, NF-2 is not the
4 greatest. It just seems to have what I have found to
5 be able to live with it. I could deal with it. And
6 my husband, he's by my side all through this.

7 Being able to hear with the ABI has given
8 me my life back, and this is really what I would like
9 to see for others to be able to get that sound, get
10 their hearing back, to be able to live life to its
11 fullest because there is so much that you miss without
12 hearing.

13 I felt normal again. I felt extremely
14 happy. There isn't anything that I do not do now.
15 Without the ABI, I would not go to wedding functions,
16 different entertainment functions. My family would
17 leave, let's say, to go to a movie. I would stay home
18 because, of course, deafness; you can't hear a movie.

19 I go to the movie now, and I can
20 understand that movie. So I'm part of my family.

21 They would do many things. They weren't
22 intentionally leaving me out, but I would stay home in

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 my home. I felt safe in my home.

2 The ABI to me is a miracle.

3 CHAIRMAN PATOW: Are there any questions
4 by the panel members?

5 MS. UMPHREY: The only problem with the
6 ABI is the microphone. It's really not compatible
7 with a microphone.

8 CHAIRMAN PATOW: I understand.

9 Just for the record, if I could ask you
10 are you affiliated with any corporation?

11 MS. UMPHREY: Oh, yes. Cochlear
12 Corporation asked me to be here only because -- and I
13 also asked them, you know, but this had happened, and
14 we knew it was going to happen for quite some time.
15 I wanted them to. I wanted to be here.

16 CHAIRMAN PATOW: And did they pay --

17 MS. UMPHREY: I wanted you to hear what I
18 had to say.

19 CHAIRMAN PATOW: I appreciate that.

20 MS. UMPHREY: Yes, they did.

21 CHAIRMAN PATOW: They paid your way here.

22 MS. UMPHREY: They paid my expenses, my

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 flight.

2 CHAIRMAN PATOW: Thank you very much.

3 MS. UMPHREY: Okay. Thank you.

4 CHAIRMAN PATOW: Thank you.

5 Our second speaker this morning is Donna

6 McLaughlin.

7 MS. McLAUGHLIN: Do they know who I am or
8 should I --

9 CHAIRMAN PATOW: If you could.

10 MS. McLAUGHLIN: Just checking.

11 Good morning.

12 CHAIRMAN PATOW: Good morning.

13 MS. McLAUGHLIN: And I bring you greetings
14 from the great State of South Carolina. It's such an
15 honor for me to be among such a fine group of people
16 this morning, and I'm glad to be here.

17 I'm here to share with you the story of a
18 miracle that occurred in my life last year. In 1988
19 I was diagnosed with neurofibromatosis 2, a genetic
20 disorder of the nervous system. It is estimated that
21 this conditions occurs in one in every 40,000 births
22 and is found on genes, chromosome number 22.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 In 1988, my diagnosis included bilateral
2 acoustic neuromas. The nine tumors on both of my
3 auditory nerves of either side of my brain stem, as
4 well as nine other tumors in random spots between the
5 meninges into my brain.

6 Since 1988 I've had three tumors removed,
7 including two which were the acoustic neuromas, the
8 last acoustic neuroma being removed 16 months ago.
9 When I awoke from the long surgery, I was profoundly
10 deaf.

11 In January prior to the surgery, Dr. Gary
12 Jackson of Nashville, Tennessee informed my husband
13 and I that I was a candidate for an auditory brain
14 stem implant. The implant took place on March the
15 24th, 1999, the same day that the last tumor was
16 removed.

17 I spent 60 days following the surgery in
18 total deafness. It was such a scary and frightening
19 time for my family and myself. In the surgery my
20 facial nerves suffered some trauma. My taste buds
21 suffered some shock. I've also experienced problems
22 with depression, my equilibrium, tear ducts, saliva

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 glands, and swallowing.

2 I have learned to function, but it has
3 been a battle every day. To some folks, they would
4 have said that my battle would have been too much. If
5 you can imagine, even chocolates taste bad to me.

6 (Laughter.)

7 MS. McLAUGHLIN: I've lost 65 pounds since
8 March of '99. The weight loss within itself is a
9 miracle.

10 (Laughter.)

11 MS. McLAUGHLIN: But not the miracle that
12 I've come to talk to you about today. The real
13 miracle came for me on June the 1st, 1999 when my ABI
14 was activated. You see, I can now somewhat hear.
15 Some folks call it the marvels of technology and some
16 call it modern medicine. You can call it what you
17 like. I simply call it my miracle.

18 My faith tells me that miracles come from
19 God, and I believe with all my heart that I received
20 a miracle.

21 Some of you may also ask what do I hear.
22 Now that's a hard thing for me to answer, that I hear

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 funny things. Let me share to you.

2 I can hear the bell on the car when the
3 door is left open and the keys are in it. I can hear
4 the sound of a Kleenex being crumpled up. I can hear
5 the wind blowing in the trees. I can hear the blinker
6 on the car when I've left the turn signal on too long.
7 I can hear the frogs croaking in our pond. I can hear
8 the crickets chirping on a still night, and I can hear
9 my dogs barking.

10 These are just a few of the things I hear,
11 and they don't sound exactly like they once sounded
12 like, but I'm hearing something, and my brain is
13 learning exactly what I'm hearing when I hear it.

14 I'm still in the process of learning. The
15 most humorous experiences happened with my ABI, and I
16 think it does us all good to look at our life
17 experiences and get a little laugh every now and then.
18 So I'm going to share with you this little experience
19 that I had.

20 After I got hooked up for the first time
21 with my ABI, my hubby and I ventured down the street
22 in Nashville, Tennessee for lunch. I went in to eat

1 lunch, and I went into the ladies room to wash my
2 hands, and I washed my hands, and I tore off a paper
3 towel, and the paper two went "zip."

4 Well, I stood there, and I thought, "Hey,
5 I don't remember when I've ever heard a paper towel
6 being torn off," or if I ever had heard a paper towel,
7 what the sound was like. So I tore off another one.

8 (Laughter.)

9 MS. McLAUGHLIN: I sat there, and I kind
10 of got tickled at myself, and I thought, "What are you
11 doing?"

12 And in walks this lady, and she looked at
13 me like, "What are you doing, lady?" And I could
14 vision myself standing in a huge pile of paper towels
15 just standing there listening to the zip.

16 So I began to get teary eyed and think
17 about how blessed this little device was really going
18 to make my life in the years to come, and I didn't
19 tear off anymore paper towels. I just went back to
20 eat my lunch.

21 I got ahead of myself. Excuse me. As the
22 lunch has gone by, I've had numerous tune-ups, as my

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 children affectionately call them, at the Wilkerson
2 Center in Nashville, and this is where I go to the
3 audiologist, Susan Amberg, there, and she works with
4 me patiently for hours at a time, and makes
5 adjustments to the tones and volumes of my device.

6 At my first reading without the aid of my
7 ABI, I could comprehend the lip reading only 55
8 percent of the time. Now, with the aid of my device
9 and my lip reading skills, I am able to comprehend 99
10 percent of the time.

11 Today I look at life totally different.
12 My friends tell me my self-esteem has improved
13 immensely. I look for the blessings in my life, and
14 I'm here to tell you that my ABI has truly been a
15 blessing.

16 Since my surgery I have been blessed with
17 many wonderful things. Let me share with you a few of
18 them. I was blessed with the ability to hear the
19 crowds' applause when my son, Sam Roland, won his
20 first Tennessee walking horse national celebration
21 down in Shelbyville, Tennessee back in September.

22 I was able to attend my eldest daughter's

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 wedding in November and know when she and my new son-
2 in-law said, "I do."

3 I'm able to hear every afternoon my middle
4 child Celia's dog bark and let me know that she's home
5 from work.

6 I have been blessed with the ability to
7 sit in church on some Sunday mornings and look at my
8 husband and say, "I believe they're singing 'The Old
9 Rugged Cross,' right?"

10 I'm really blessed. My hope is that God
11 will grant me the resolve to use my talents and my
12 gift of hearing to the benefit of others. Until my
13 activation I knew no one involved in the National Ear
14 Foundation, the Cochlear Corporation, or any of the
15 panels -- any of the members of this panel with the
16 Food and Drug Administration. I have no knowledge of
17 any one of you knowing either me or my family.

18 However, in some wonderful way I think the
19 Lord has brought us all together. He has allowed me
20 to hear again, and you are a natural part of my
21 miracle, and I thank you from the bottom of my heart.

22 CHAIRMAN PATOW: I want to thank you for

1 your comments this morning. What I need to know from
2 you is two things. Are you affiliated with any
3 company? And has anyone paid your way to come here
4 today?

5 MS. McLAUGHLIN: Thank you.

6 MR. WEST: Can I help you? He needs to
7 know if we paid your way here today, Cochlear
8 Corporation.

9 MS. McLAUGHLIN: Not that I'm aware of.
10 You invited me.

11 (Laughter.)

12 MS. McLAUGHLIN: I don't know that you
13 did.

14 MR. WEST: Well, I think we did take care
15 of her air fare.

16 CHAIRMAN PATOW: Thank you very much.

17 And thank you for your comments.

18 MS. McLAUGHLIN: Thank you very much.

19 CHAIRMAN PATOW: Our third speaker this
20 morning is Henry Ilecki, Ph.D.

21 DR. ILECKI: My name is Henry Ilecki. I
22 am the Director for Audiology Practice in Industry and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Private Practice at the American Speech-Language-
2 Hearing Association, ASHA.

3 ASHA has supported or is supporting my
4 appearance here this morning, and I have no other
5 company affiliations.

6 Members of the Ear, Nose, and Throat
7 Devices Panel, good morning, and thank you for the
8 opportunity to offer general commentary on the use of
9 auditory brain stem implants and the management of
10 patients with profound hearing loss secondary to
11 bilateral surgical lesions of the acoustic nerve.

12 The American Speech-Language-Hearing
13 Association is a professional and scientific
14 organization that represents over 98,000 audiologists,
15 speech-language pathologists, and hearing and speech
16 scientists. The association encourages the
17 development, evaluation, and implementation of
18 procedures, programs, and technologies holding promise
19 in the areas of identification, evaluation, and
20 treatment of individuals with hearing loss and related
21 disorders.

22 The work of the Food and Drug

1 Administration in this regard has been highly
2 constructive and beneficial to the public health, and
3 there is every anticipation that results of this
4 record will be extended into the newly expanding realm
5 of auditory brain stem implants.

6 Following the diagnosis of acoustic
7 neuromas, patient management decisions are made that
8 are based on a variety of factors. These include the
9 size of the tumors, the overall health of the patient,
10 and the feasibility of hearing preservation.

11 When a decision is made to proceed with a
12 surgical approach that will result in complete
13 deafness, the only means to provide the patient with
14 hearing postoperatively is by means of the auditory
15 brain stem implant, or ABI.

16 While clinical experience with this
17 technology is relatively limited compared, for
18 example, to the research findings of experiences of
19 Cochlear implantees, preliminary indications suggest
20 that ABI recipients generally receive benefits of
21 sound detection and discrimination that are similar to
22 those afforded by the first generation of Cochlear

1 implants.

2 ASHA applauds this technology and
3 recognizes the hope it provides to persons having to
4 undergo the hearing debilitating effects of some of
5 the surgical remedies to excise acoustic neuromas. In
6 this promising and emerging technology, ASHA
7 recommends the following areas of investigation be
8 included during the FDA's deliberative process.

9 In reviewing desired outcomes, the
10 question should be posed as to what extent ABIs alter
11 recipients' perception of disability and their sense
12 of quality of life. What differences, if any, are
13 there between recipients of ABIs and individuals
14 choosing not to have the procedure or different
15 procedures? How do ABI recipients fare compared to
16 Cochlear implant recipients?

17 The second broad area of study recommended
18 by ASHA concerns the vast array of existing assistive
19 technologies in the marketplace available to the
20 myriad users of conventional hearing aid devices.

21 Assistive listening devices range in size
22 and cost from the simple strap mounted, battery

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 operated telephone amplifier, to FM or infrared
2 devices used to enhance the signal to noise ratio of
3 stage material broadcast, to compatibly equipped
4 hearing aid users in the audiences of concert and
5 lecture halls.

6 The extent to which ABI utilizes existing
7 assistive technology in special listening
8 circumstances needs to be documented and, where
9 lacking, development encouraged.

10 A particular concern to this clinical
11 population, that is, persons who are bilaterally
12 deafened, is their post implantation access to
13 assistive technology of the alerting variety. It is
14 essential that implantees be able to benefit from such
15 potential life saving devices as smoke, fire, and
16 carbon monoxide detectors and security alarms, as well
17 as convenience announcement devices interfacing with
18 doorbells or knockers, sleep alarms, et cetera.

19 ASHA urges the panel to consider a general
20 recommendation recognizing the essential role
21 performed by the audiologist as a critical hearing
22 care professional in candidacy consideration and the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 rehabilitative process.

2 Decisions relevant to the implantation
3 process should relate to auditory status and auditory
4 processing information derived through comprehensive
5 pre, peri, and post implantation audiological
6 evaluation performed by an audiologist.

7 Consider that even when the surgery spares
8 the cochlear nerve anatomically, it is frequently
9 compromised physiologically. Thus, an intact nerve
10 following tumor excision may be unresponsive.
11 Audiological intraoperative monitoring incorporating
12 the auditory brain stem response should, therefore, be
13 an essential aid to the implantation decision.

14 Certainly in the area of Cochlear
15 implants, but as well in all device based forms of
16 intervention, the critical component to successful
17 patient outcome has been shown to be dependent upon
18 regular intensive and quality post surgical device
19 orientation, counseling, and rehabilitation by the
20 audiologist. Thus, crucial areas of audiologists'
21 participation in an ABI practice protocol would
22 include determination of candidacy, preoperative

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 counseling, and post operative mapping, and otologic
2 rehabilitation.

3 In addition, speech-language pathologists
4 may also be involved in speech and voice assessment
5 and treatment of ABI patients in need of such
6 services.

7 In concluding, the review and regulation
8 of medical devices for safety and efficacy, such as
9 hearing aids and implantable devices, is a critically
10 important function of the Food and Drug
11 Administration. It is the view of the American
12 Speech-Language-Hearing Association that in evaluating
13 emerging technologies and applications, the Food and
14 Drug Administration recognizes and promotes the value
15 of an audiological component to insure the eventual
16 clinical acceptance, utility, and successful outcomes
17 with auditor brain stem implants.

18 Thank you for the opportunity of
19 addressing this panel.

20 CHAIRMAN PATOW: Thank you.

21 I want to thank each of our speakers this
22 morning for their valuable comments, and I'm certain

1 that the information they've shared will be very
2 helpful to the panel.

3 At this point, we'll proceed to our open
4 committee discussion session. Nancy Brogdon will at
5 this time give the Division update.

6 MS. BROGDON: Good morning. I have
7 several announcements to make. Most of these I said
8 yesterday, but we have new people here on the panel
9 and in the audience today.

10 First, I'd like to let you know that Dr.
11 Ralph Rosenthal, who is our Division Director in the
12 Division of Ophthalmic and ENT Devices, is working
13 temporarily in our Center Director's office on Health
14 Care Financing Administration issues. He's expected
15 back in a few months, but I'm the Deputy Director, and
16 that's why I'm sitting here today.

17 Secondly, I'd like to announce that our
18 Office Director has been selected -- this Office
19 Director is over the six reviewing divisions who
20 review the whole spectrum of medical devices, and this
21 new Director is Dr. Bernard Statland.

22 He was here yesterday visiting, but was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 unable to be here today.

2 Dr. Statland received his M.D. degree and
3 his Ph.D. in biochemistry from the University of
4 Minnesota. He did residencies at the University of
5 Copenhagen and the University of Minnesota Hospitals.

6 He also served in the Public Health
7 Service in New Orleans and at the NIH Clinical Center.

8 Dr. Statland is a clinical pathologist,
9 and he's held a number of positions, including Medical
10 Director and CEO of laboratories at the North Shore
11 Long Island Jewish Health System, and he has run his
12 own consulting firm in Minneapolis. He has many
13 publications in several areas of interest, and we're
14 happy to have him on board.

15 I'd also like to announce that Mr. Harry
16 Sauberman, who is the Chief of the ENT Branch, is
17 working currently in the Office of Device Evaluation
18 on special projects. Among those is partnering with
19 governments of other countries.

20 In the meantime, Dr. Morris Waxler is the
21 Acting Chief of the ENT Branch. Dr. Waxler's
22 experience as a Branch Chief and as a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 neuropsychologist made him well suited for selection
2 to be Acting Director of this branch.

3 I'd like to introduce a new reviewer in
4 the ENT Branch, Dr. James Kane. Dr. Kane is an
5 audiologist. He has a B.S. in speech and hearing
6 science from California University of Pennsylvania.
7 He has his Master's and Ph.D. in audiology from the
8 University of Pittsburgh, and he did a post doc. at
9 the University of Pittsburgh Medical School.

10 He's practiced for 22 years both privately
11 and in the Veterans Administration, and he's held
12 various supervisory positions. Jim is a Fellow of the
13 American Academy of Audiology, and we welcome him to
14 the branch.

15 One last item. As I mentioned yesterday,
16 we have three voting members from this panel who will
17 have completed their four-year terms in October and
18 probably are not likely to attend another meeting as
19 a voting member, and presented certificates to Drs.
20 Woodson and Duffell yesterday.

21 Today I'd like to read a letter from
22 Commissioner Jane Henney, Commissioner of Food and

1 Drugs, to Dr. Shelton.

2 "Dear Dr. Shelton:

3 "I would like to express my deepest
4 appreciation for your efforts and guidance during your
5 term as a member of the Ear, Nose, and Throat Devices
6 Panel of the Medical Devices Advisory Committee. The
7 success of this committee's work reinforces our
8 conviction that responsible regulation of consumer
9 products depends greatly on the participation and
10 advice of the nongovernmental health community.

11 "In recognition of your distinguished
12 service to the Food and Drug Administration, I am
13 pleased to present you with a certificate."

14 Signed Dr. Jane Henney.

15 That completes my announcements.

16 CHAIRMAN PATOW: Thank you.

17 I also would like to express my
18 appreciation to each of the members who have completed
19 their term for the hard work that they've put in in
20 reviewing the PMAs and other documents, and in
21 participating in these panels.

22 At this time we'll have the branch update.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Morris Waxler, Ph.D., will present.

2 DR. WAXLER: Good morning again. I'd like
3 to introduce our branch. Our branch is Karen Baker,
4 who's an R.N., scientific reviewer; Teri Cygnarowicz,
5 who is an audiologist; Dr. Sid Jaffee, Sidney Jaffee,
6 who's our Medical Officer; Dr. James Kane, our
7 audiologist, whom you've already met; Dr. Vasant
8 Malshet, who's our toxicologist; and Dr. Alfred
9 Montgomery, who is not here.

10 (Laughter.)

11 DR. WAXLER: He was here just a moment.
12 He slipped out on me.

13 In addition, we have reviews often from
14 other folks at CDHR, including Dr. Brian Beard and Dr.
15 Bill Regnault and Dr. Victor Krauthammer. I won't
16 give the list of all of those who have helped us, but
17 in this particular case, I would like to acknowledge
18 Dr. Rhonda Ballum (phonetic) from the Office of Orphan
19 Product Development. I think I got that right, and
20 for her special role in this application.

21 The Cochlear Corporation received a grant
22 from the Orphan Products Division to study this area,

1 and we appreciate your efforts.

2 That's it for me.

3 CHAIRMAN PATOW: Thank you, Dr. Waxler.

4 We'll begin now our consideration of PMA
5 P000015, the nucleus auditory brain stem implant,
6 presented by Cochlear Corporation.

7 I understand that there are seven
8 presenters scheduled. We will have an hour for their
9 presentation, and we will plan to end, therefore, a
10 couple of minutes after 11 o'clock.

11 MR. WEST: My name is Ron West, and I'm
12 President of Cochlear Corporation, the sponsor of PMA
13 000015, and this morning I'm going to give you just
14 some context of the history of the development of this
15 auditory brain stem device and introduce our speakers.

16 We do have seven speakers. Mine will be
17 brief, and we should be able to finish on time.

18 The history of the auditory brain stem
19 implant project really goes back over 20 years, and
20 the first auditory brain stem implant was performed by
21 Dr. William House and Dr. William Hitselberger in
22 1979. It was a single channel device with a ball type

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 electrode and was successful, and that patient is
2 still benefitting from that implantation today.

3 Some five or six years later, the House
4 Institute applied for an IDE with the ENT Branch, and
5 that was granted, and they implanted a series of 25
6 patients, again, with a single channel device that
7 utilized a platinum plate electrode.

8 The early cases were percutaneous
9 connection, and it was later changed to a
10 transcutaneous connection.

11 Next slide, please.

12 Then in the early '90s, there was a
13 collaboration that began between Cochlear Corporation
14 and our design people in Sydney, Australia, the
15 Huntington Medical Research Institute, and the House
16 Ear Institute to develop a multi-channel device. An
17 initial series of three patients were implanted at the
18 House Ear Institute during that period of '92 to '93,
19 and then we filed for an IDE to expand that program in
20 1994.

21 That multi-channel device was based on the
22 nucleus 22 stimulator design and utilized an eight

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 platinum disk array. That program led over the
2 intervening seven to eight years to a premarket
3 approval application that was filed in March, and that
4 application is based on the nucleus 24 auditory brain
5 stem implant and involves some 90 subjects studied for
6 safety and 60 subjects studied for effectiveness.

7 That's a very brief overview of what's
8 been a long journey to this day, and I'd like to take
9 a moment just to thank the ENT Branch for their
10 efforts and review of this PMA application, as well as
11 the panel.

12 I know what it takes to write one of these
13 applications, and I think it's equally or maybe more
14 onerous to review them. They're massive in nature,
15 and I applaud your efforts.

16 I would also like to acknowledge the
17 Orphan Device Office's support of this program. Dr.
18 Waxler mentioned that, but I think it's important that
19 everyone understand that we might not be here today
20 without their support.

21 There were ^{**}actually two grants. The
22 initial one was to the House Ear Institute and

1 involved, I think, \$300,000 over three years, and then
2 there was a second grant to Cochlear Corporation of a
3 similar amount. So without that investment, I don't
4 know that this PMA would have been possible.

5 So thank you, Rhonda. We really
6 appreciate that support.

7 Now, if I may -- next slide, please -- I'd
8 like to introduce our speakers.

9 Ms. Patti Arndt, who's our Manager of
10 Clinical Studies, will give an overview of the PMA.

11 And then Professor Martyn Hyde, who is our
12 statistical consultant and is Professor of
13 Biostatistics at the University of Toronto, will cover
14 the statistical design.

15 And then Dr. Brackmann and Dr. William
16 Hitselberger, surgical consultants to Cochlear from
17 the House Ear Institute, will cover device safety in
18 surgery.

19 And then Ms. Kiara Ebinger, a senior
20 clinical study specialist for Cochlear will cover the
21 device effectiveness data."

22 And then we'll round out the session with

1 Ms. Arndt again covering the proposed labeling, post
2 market surveillance recommendations, and conclusion.

3 So without further ado, I'd like to ask
4 Patti to come up and give the overview.

5 MS. ARNDT: Thank you very much.

6 Just wait a second for the slide here.

7 What I'd like to do is just to describe
8 the device that you'll be discussing and deliberating
9 this morning and explain to you some of the bases for
10 our recent PMA submission.

11 The device, as you know, is the nucleus 24
12 auditory brain stem implant. The N24 ABI is intended
13 for use in individuals who are 12 years of age or
14 older who have been diagnosed with neurofibromatosis
15 Type 2. That's hard to say this early in the morning.

16 The device is implanted during the
17 recipient's first or second side tumor removal
18 surgery.

19 Next slide, please.

20 Just to give you a reminder, NF2 has a
21 number of auditory manifestations. Certainly the
22 primary symptom is the presence of bilateral

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 vestibular schwannomas. These tumors generally grow
2 progressively over time, and individuals lose their
3 hearing at the same time. So NF2 recipients
4 demonstrate a progressive hearing loss as their tumors
5 grow.

6 By the time that their second tumor is
7 removed surgically, they are totally bilaterally
8 profoundly deaf. Because the connection between the
9 cochlea and the first brain stem nucleus, first
10 auditory brain stem nucleus is disrupted, these
11 patients aren't viable candidates for either hearing
12 aids or cochlear implants as treatments. Really their
13 only possible treatment, their only access to sound is
14 through electrical stimulation of the cochlear nucleus
15 via the ABI.

16 Next slide.

17 So the device that we're asking you to
18 consider this morning is composed of a number of
19 pieces. It is a system for patient use. The primary
20 piece of this system is the implanted portion, the
21 nucleus 24 auditory brain stem implant, which consists
22 of a receiver stimulator and a 21 electrode brain stem

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 array.

2 The external equipment worn by the patient
3 includes the SPrint body worn speech processor and
4 headset. The SPrint is programmed to implement the
5 spectral peak or SPEAK encoding strategy. Part of the
6 headset is an adhesive retainer disk which allows
7 patients to wear their headsets following the removal
8 of the internal magnet from the device, which allows
9 them ultimately to have MRIs, which are an important
10 diagnostic tool in this population.

11 There are a family of speech processor
12 accessories which allow these patients access to
13 assistive listening devices, as was mentioned by the
14 ASHA representative earlier.

15 And then lastly there are device
16 programming systems, both a portable system, a desktop
17 system, and some software that allow audiologists to
18 program the devices.

19 As you can see from the stars on the
20 slide, all but two of these components have been
21 previously approved by FDA^{**} or cleared for commercial
22 use as part of the nucleus 24 Cochlear implant system.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Next slide, please.

2 Just to remind you, the nucleus 24
3 receiver stimulator component is identical to the
4 receiver stimulator used in the nucleus 24 Cochlear
5 implant, which was released by FDA in June of 1998.

6 This particular receiver stimulator allows
7 for some flexibility in programming options for these
8 patients. The electronics have a capacity for high
9 rate stimulation. The implant is capable of
10 implementing multiple speech processing strategies.
11 It can be programmed to implement a number of
12 stimulation modes, three monopolar modes, a variable
13 bipolar mode, and a common ground mode.

14 As I mentioned before, the magnet that is
15 placed kind of in the center of the receiver coil can
16 be removed for MRI, and then lastly, this particular
17 technology allows for some comprehensive telemetry,
18 which allows the audiologist to verify the
19 functionality of the device and includes a tool called
20 neural response telemetry.

21 Next slide, please.

22 This is just a picture of the device.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 You've probably all seen it in a lot of detail. I'll
2 just point out to you the magnet in the middle of the
3 receiver coil on the right-hand side of the slide, and
4 that generally is removed for these patients prior to
5 implantation, which, again, allows them to undergo
6 MRIs.

7 Next slide.

8 This is a little more detail on the brain
9 stem electrode array. It features 21 platinum disk
10 electrodes, each seven millimeters in diameter. The
11 electrodes are carried or organized in a three-by-
12 seven matrix and are carried on a silicone pad, which
13 is eight and a half by three millimeters.

14 And then lastly, the electrode is backed
15 by a T-shaped piece of mesh, and what the mesh does is
16 it allows connected tissue to grow through the mesh
17 following implantation, which ultimately fixes the
18 electrode pad against the cochlear nucleus.

19 Next slide.

20 So what is the basis of our request today?

21 As you know,** we have submitted U.S.
22 clinical trial data collected under IDE G930077 using

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

www.nealrgross.com

1 the nucleus 22 auditory brain stem system. This is a
2 system that uses the same receiver stimulator
3 component as the commercially available nucleus 22
4 Cochlear implant, which is our previous generation
5 Cochlear implant system, and it features an eight
6 channel brain stem electrode.

7 We have also submitted safety and
8 effectiveness data that was obtained as part of a
9 European clinical trial on 27 subjects, and this
10 receiver-stimulator component, again, is the nucleus
11 22 component, but this time it's coupled with a 20 or
12 21 channel brain stem electrode.

13 As Ron mentioned, the numbers in the U.S.
14 trial are rather large. We've studied 90 subjects for
15 safety and 60 for effectiveness.

16 Next slide.

17 In addition, we've submitted to FDA quite
18 extensive laboratory testing, which verifies the
19 function of the nucleus 24 ABI.

20 And then lastly we have clinically
21 validated all of the components of the N24 implant.
22 So the receiver-stimulator component, which has been

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 surgically implanted in over 10,000 Cochlear
2 implantations worldwide, and then finally the 21
3 electrode array was evaluated as I just mentioned in
4 Europe and has now received the European CE Mark.

5 Next slide.

6 All of that is summarized in this
7 difficult to see slide, but the point of it is that
8 the nucleus 24 ABI, which is described in the lower
9 bottom right-hand box, has evolved from two different
10 Cochlear implant systems, a previous generation
11 system, our current system, and two different versions
12 of an ABI system, both of which have been based on the
13 nucleus 22 technology, one with an eight electrode
14 array, and the other with a 20 or 21 electrode array.

15 Next slide.

16 So why are we in such an unusual
17 situation?

18 I think we need to kind of recall the
19 context that NF2 puts us in when we start studying
20 ABIs. NF2 is a disease that occurs at a very low
21 incidence and low prevalence. There just aren't a lot
22 of people around to receive ABIs and to study it. In

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 this day and age our Cochlear implant and ABI
2 technology is rapidly changing. We really are making
3 very exciting improvements all the time.

4 So in that context, it's very hard to
5 recruit enough investigational subjects to continually
6 evaluate new technology.

7 Why do we think it's important that you
8 approve the nucleus 24?

9 First of all, it's our desire to provide
10 these patients with the current generation technology.
11 This technology has features that are very well suited
12 to ABI recipients and needed by them. Perhaps even
13 more than our Cochlear implant recipients, ABI
14 patients can benefit from the capacity of the N24
15 technology, its increased flexibility for programming,
16 and ultimately its upgradability.

17 And then lastly, we strongly believe that
18 we have shown the safety and effectiveness of the
19 nucleus 24 receiver, the 21 electrode array, and
20 certainly the SPEAK encoding strategy through the
21 other mechanisms that I've just mentioned.

22 Thank you.

1 DR. HYDE: Good morning. I'm Martyn Hyde.
2 I'm someone who has no financial interest in the
3 company or in the outcome of today's proceedings, but
4 I am being reimbursed for expenses incurred in
5 attending this meeting.

6 I'd like to give you a brief, very brief
7 outline of the experimental design and statistical
8 methods.

9 Next slide, please.

10 I'm going to touch upon these points.

11 Sorry. Can you go back one.

12 The study design itself; the disposition
13 of the study cohort; the analytical approach; some
14 issues of power and sample size; and a few simple
15 conclusions.

16 Next slide.

17 The study design was a familiar one to
18 you, I believe, the single subject, repeated measures
19 design replicated in 60 subjects. The key features of
20 this design, of course, are that old treatment
21 conditions can apply on each subject. So we can make
22 within subject comparisons.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Most importantly about that, we can make
2 within subject significance tests. We can do
3 statistical significance tests within individuals
4 using the binomial distribution, and then we achieve
5 generalizability of those single subject results using
6 group aggregate statistics.

7 Next slide.

8 The disposition of the cohort is
9 summarized here. Ninety-two subjects within F2
10 received the implant, and two of those subjects
11 deceased for reasons unrelated to the ABI. That
12 leaves 90 subjects who yielded the safety data. There
13 was a subgroup of 13 subjects who had no auditory
14 precept on initial activation, and 17 subjects, other
15 subjects, yielded no data.

16 There was a variety of reasons for that,
17 the most common reason being that the subjects were
18 too ill to attend the follow-up sessions.

19 This leaves an effectiveness cohort of 60
20 subjects providing effectiveness data, and these were
21 obtained at eight sites, but I should remark that
22 there was an overwhelming preponderance of subjects,

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 a substantial number of subjects at one site, more
2 than half of the subjects.

3 Next slide, please.

4 The general analytical approach is that
5 the subjects who lacked auditory precepts constitute
6 a distinct and separate sub-cohort. They are
7 qualitatively distinct from the balance of the
8 effectiveness cohort. So they are separated and
9 described.

10 The second point, the missing data, they
11 appear in my opinion to be unrelated to probable
12 outcome, and so I believe that there is no necessity
13 for adjusting the outcome data for bias.

14 I believe, therefore, that the
15 effectiveness cohort of 60 should be considered to be
16 representative of the target population, and in that
17 cohort we used completely standard descriptive or
18 inferential statistical methods to detect and quantify
19 the treatment effects.

20 Next slide, please.

21 We targeted outcome data at six months,
22 but some subjects did not have six month data

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 available at the time of analysis. They offered three
2 month data. Therefore, we believe that because it's
3 possible that they might not have achieved the
4 asymptotic performance levels at six months, this
5 outcome data may be slightly conservative.

6 The actual outcome measures were various
7 and may be divided into four main areas: the sound
8 and speech recognition test; the test of lip reading
9 enhancement; and two secondary measures, which are
10 subjective performance ratings and subjective benefit
11 ratings.

12 Next slide, please.

13 The specific statistical methods that we
14 used, we used single subject significance tests of
15 improvement from chart score for the sound recognition
16 tests. These tests are, of course, one-tailed because
17 you can't get worse than chance score, and they were
18 based on the binomial model.

19 The second thing, we used single subject
20 significance tests if performance changed with device
21 activation and without device activation. These tests
22 are conservatively two-tailed, and of course, based on

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 the binomial distribution as well.

2 We summarized the cohort data by using
3 proportions of subjects who had statistically
4 significant improvements or achieved specific
5 performance levels. I just want to point out that
6 that is, in fact, an extremely conservative approach.
7 It is unusual to have statistical significance in an
8 individual, and simply to count the proportion of
9 people who achieved that is a conservative
10 representation of the outcome.

11 And then we used other statistics that are
12 always required by FDA in these kinds of designs,
13 typical ones such as just group means, medians, and
14 standard deviations.

15 Next slide, please.

16 A few comments on power. Power, of
17 course, as you all know, is the probability of
18 detecting a treatment effect when it genuinely exists,
19 and we usually need to prove that we can achieve a
20 power level of .8, .9, with an alpha of less than .05.

21 Group size is usually computed to achieve
22 the target power for a minimum practically significant

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 treatment effect size.

2 In this single subject design, the
3 significance tests are possible in individuals. So
4 the single subject power is dictated mainly by the
5 number of items in the test list. That's the basic
6 feature of the binomial model, and the power depends
7 purely on the number of items in the test list.

8 But we also require high group power in
9 order to confirm the consistency of effects across
10 subjects and to imply that the findings are
11 generalizable.

12 Next slide, please.

13 What we found was highly significant
14 improvements for several of the primary measures and
15 in many individual subjects. For several of the
16 primary measures the actual group power that we found
17 was extremely high, such that the likelihood of the
18 observed results being obtained by chance if there
19 were no treatment effects is close to zero, virtually
20 negligible, massive statistical significance at the
21 group level.

22 So I believe that the power of the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 experiment was clearly sufficient at both the single
2 subject and the group levels. The cohort sample size
3 and effect consistency were also, I believe,
4 sufficient to yield reasonable generalizability of the
5 findings.

6 So my simple conclusions, in my opinion,
7 the study sample is representative and is sufficient.
8 the statistical approach and methods used were valid
9 and were appropriate, and the proposed performance
10 claims are conservative and are fully supported by the
11 data.

12 Thank you.

13 DR. BRACKMANN: Hello. I'm Derald
14 Brackmann. I'm an otologist, neurotologist, and
15 practice at the House Ear Clinic in Los Angeles. I am
16 on the Medical Advisory Board of Cochlear Corporation.
17 They have paid my expenses to this meeting. I have no
18 financial interest in the company other than that.

19 Safety of a device, of course, is of
20 utmost importance. By our oath we are committed to do
21 no harm, and I know that is the FDA's primary concern
22 as well, that devices must prove to be safe as well as

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 effective.

2 But safety is of paramount importance.
3 I'd like to address that.

4 Next slide, please, Jennifer.

5 The quick review -- I know this is
6 redundant for most of you. You've already heard that
7 the external receiver stimulator is identical to the
8 currently used Cochlear implant devices.

9 Cochlear implants, of course, are not
10 appropriate for patients with NF2 because in most
11 cases in removal of their tumor the auditory never,
12 the pathway to the higher centers, is destroyed. So
13 you may think of this as just like a Cochlear implant,
14 except that the interface with the auditory pathway is
15 at the next substation or next way station, the
16 cochlear nucleus.

17 So it's identical in many ways except that
18 the electrode stimulates the next higher pathway in
19 the auditory chain.

20 Next, please.

21 Fortunately, nature was kind to patients
22 and to us in that the cochlear nucleus lies on the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 roof of a natural pathway from the brain stem. It
2 lies in the roof of the lateral recess of the fourth
3 ventricle. So the electrode can be placed onto the
4 surface of the cochlear nucleus without the necessity
5 of putting anything into the brain itself.

6 So there's a natural pathway which will
7 hold the electrode in position, and the electrode is
8 designed to fit exactly into that natural pathway so
9 that upon insertion, it will be in surface contact
10 with the cochlear nucleus and thus stimulate the
11 auditory pathway.

12 Next, please.

13 So this is obviously not a real tumor, but
14 a simulation of a tumor that shows the lateral recess
15 of the fourth ventricle. The flocculus lies
16 posteriorly, and you locate the eighth nerve, follow
17 the eighth nerve into the lateral recess and then
18 insert the electrode over the surface of the cochlear
19 nucleus in the lateral recess of the fourth ventricle.

20 And as I said, the electrode is designed
21 to fit in that it has a carrier of dacron which has
22 some surface adhesion, which will hold the electrode

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 in place.

2 Next please.

3 The total number of patients that have
4 been implanted is 92. Two succumbed to the disease
5 prior to stimulation, and so that we have safety data
6 on the remaining 90 subjects.

7 I do want to spend a minute telling about
8 those two patients because those, of course, that's
9 the ultimate complication, and I do want to address
10 that briefly.

11 One patient died of serratia meningitis.
12 This is a Gram-negative meningitis. He had done very
13 well in surgery. In fact, he was discharged from the
14 hospital, seemed to be doing well. He developed a
15 headache, was readmitted, and initially had negative
16 spinal fluid cultures.

17 In infection was being harbored in the
18 ventricle, in the lateral ventricle. He then had a
19 sudden turn for the worst and positive cultures for
20 serratia. Serratia meningitis has an overall
21 mortality of about 80 percent. It's extremely
22 virulent, extremely difficult to treat. We had a

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 limited autopsy. There was no evidence of infection
2 or loculation of any pus or infection around the
3 implant. We do not believe that his death or his
4 meningitis and death was in any way related to the
5 auditory brain stem implant.

6 The other death was a patient who had an
7 early massive tumor. He had an early stroke, a brain
8 stem stroke, never recovered from that. He never
9 gained consciousness so he could be stimulated.
10 Again, there was no evidence of any relationship to
11 the auditory brain stem implant.

12 So we believe that these two deaths were
13 related to disease and not to the auditory brain stem
14 implant.

15 Next slide, please.

16 Complications for the remaining 90
17 patients we divided into three groups. As you see
18 here, minor complications were those that could be
19 resolved by either changing the external device. Some
20 of them just resolved spontaneously. Some required
21 reprogramming, things of that nature. All of the
22 minor complications were resolved with these measures.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Major complications were those that
2 required hospitalization, revision or removal of the
3 implant.

4 And finally, a complication -- it could be
5 called a complication, but as I'll show you in a
6 minute, there were a number of patients who did not
7 receive auditory precepts from the device, and they
8 did not require removal or any treatment, but they did
9 not benefit from the device because they did not
10 stimulate.

11 Next slide, please.

12 First I'll address the non-stimulation
13 patients. There's a little confusion in the data
14 because there are actually 14 who did not stimulate at
15 first hook-up, and there's 13 that you see at the top
16 of this graph, and then there's non-stimulation with
17 the first ABI, who subsequently performed well with
18 the second ABI. So if you take those two groups,
19 there were 14 who did not stimulate at initial hook-
20 up.

21 Two patients had good stimulation for a
22 period of time, in one case a couple of years, and

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 then had deterioration of the device. In these cases,
2 in one of them, CT identified that the electrode had
3 migrated, and the other the cause was unknown, but
4 there was non-stimulation initially on 14 of the
5 patients.

6 Next, please.

7 Two major complications. In one case the
8 skin broke down over the receiver stimulator. A flap
9 was placed over this, but that flap also broke down,
10 and the patient had poor skin turgor. It was deemed
11 that it was not reasonable to continue to try to cover
12 it, and the patient elected and the physicians
13 concurred that it was best to remove the device.

14 The electrode was cut in the mastoid
15 cavity. The electrode pad was not removed from the
16 brain stem. There have been no further complications.
17 The skin healed, and there have been no further
18 complications related to that.

19 Another patient developed a skin flap
20 infection. A decision was made to revise the flap.
21 During the revision, the electrode was moved, and
22 intraoperatively electrically evoked response.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Audiometry responses could not be obtained, and it was
2 again elected to remove the device, and the skin
3 healed and there have been no further sequelae, but
4 the patient, of course, has no benefit from the
5 device.

6 All of the other complications that are
7 listed there are minor and temporary. At the time of
8 initial stimulation, one patient had fluid beneath the
9 flap which subsequently resolved. One patient had
10 some dizziness, and one electrode was found to produce
11 that. That electrode was programmed out of the
12 device, and that was resolved.

13 The headache resolved spontaneously, and
14 deprogramming two electrodes resulted in relief of the
15 light headedness/dizziness with the device use.

16 So all of these complications, except the
17 two that required explanation, have been resolved
18 either expectantly or with device reprogramming.

19 Next, please.

20 Some considered device complications. You
21 know, the magnets are removed, and nobody has yet
22 mentioned, but this is compatible with MRI. We have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 done hundreds of MRIs on these patients. So there's
2 no internal ferromagnetic material.

3 The way the device is held in place is
4 that the patients glue a metallic disk over the
5 implant, and two patients have had some irritation
6 from the glue. One has solved that by using adhesive
7 on either side of the disk. The other just by
8 changing it every night, not leaving it on for three
9 or four days the way some patients do. So that
10 problem has been resolved.

11 Two patients had some sound quality
12 changes, one after a surgery for the other side, and
13 that was reprogrammed, and the problem resolved.

14 One patient heard some clicking and
15 popping, and we offered to reprogram it, and he said,
16 no, it went away, and there was one patient who had
17 pain on electrode, one electrode or two electrodes
18 actually, and they were programmed out and that
19 problem solved.

20 So next slide, please.

21 So when you address, again, going back,
22 addressing the initial question that I posed, is it

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 safe. We believe it is. The extrusion rate is very
2 similar to Cochlear implants. That's a known
3 complication of Cochlear implants. The flaps need to
4 be designed well, but no problems, no harm was done
5 except that the patient doesn't hear after explanation
6 of the device.

7 ~~Next, please.~~ Next, please.

8 All but two of the medical surgical
9 complications were minor. They were all resolved with
10 reprogramming or with doing something with the device
11 of a minor nature. The non-stimulators are a problem
12 that is apparent.

13 Next, please.

14 The European data supports our experience.
15 In Europe, there have been on ABI related neurological
16 problems reported. They have not had major
17 complications attributed to the ABI, and so their
18 experience has been similar to ours.

19 Next, please.

20 So to summarize, I think these things have
21 already been said. Seventeen, point, eight percent of
22 the subjects do not receive auditory precepts. In the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 majority of those, it's distortion of the cochlear
2 nucleus at the time of surgery. These patients in our
3 series in the United States, two-thirds of the
4 patients were done at the time of second tumor
5 removal. These patients typically hold onto their
6 hearing and tumors as long as possible.

7 These are huge tumors by the time the
8 second tumor is removed, and there's distortion in the
9 area of the brain stem and cochlear nucleus which
10 makes device placement difficult in some cases.

11 The other complications are minor and have
12 been resolved with reprogramming. There have been no
13 life threatening complications or risks that we have
14 been able to associate with the auditory brain stem
15 implant. There have been no device failures over an
16 eight-year period.

17 So I conclude that the device is safe and
18 that no patient has been harmed by its placement.

19 Thank you.

20 DR. HITSELBERGER: I'm Dr. Hitselberger,
21 and I don't have any business connections with
22 Cochlear or anybody else.

1 (Laughter.)

2 DR. HITSELBERGER: They paid for my way to
3 come here.

4 I'm going to talk about the neurologic
5 evaluation. Most of this stuff, I think, you all
6 know. One in 40,000, I just computed that. That
7 means that there are 6,250 patients in the United
8 States with NF2. Whether they all come to have an
9 implant I don't know.

10 Bilateral eighth nerve tumors, and the
11 symptoms from these tumors vary depending on the
12 location of the tumor. As you can see, some are along
13 the spine. Brain tumors: loss of balance, headaches,
14 and also you can get meningiomas and even mixed
15 tumors, with the average life expectancy, 40 years.

16 Next one.

17 And what we wanted to do was just evaluate
18 the patients neurologically to see that there wasn't
19 any change in the status that could be attributed to
20 the ABI and not to the disease. It's not easy to do,
21 but I think in general we were able to accomplish
22 that.

1 I felt there were eight parameters that
2 should be evaluated to really ascertain whether or not
3 the implant itself was responsible for neurologic
4 deterioration.

5 Next slide.

6 And these were obviously the cerebellar
7 function. It's right near the cerebellum, sensory
8 function, descending and ascending tracts from the
9 cerebrum; extraocular motion; symmetry of palate;
10 ninth and tenth nerves; swallowing;
11 sternocleidomastoid; all of these are nine, ten, and
12 11, which are near there, and just as a general check,
13 see if there's any disparity in the reflexes and the
14 intraocular pressure, i.e., is there any increased
15 intracranial pressure from whatever cause.

16 In the 80 cases that were evaluated, we
17 found no ABI that we could feel for sure were related
18 to the ABI.

19 Next slide.

20 Eighty of the 90. There were no changes
21 vis-a-vis the parameters I just went over for this.

22 Do you have that slide on the -- yeah, let

1 me just show you.

2 Let me just show you a little bit about --
3 go ahead. Oh, press the button? This is what the
4 tumors look like, and as you can see, after you take
5 the tumor out, there's a tremendous amount of
6 distortion of the brain, and this is why we
7 probably -- there may be intrinsic deficit in the
8 nucleus itself, which it's amazing that it works as
9 well as it does.

10 Next slide.

11 Now, this is just a post-op to show you
12 what the ABI looks like on a scan done post-op. We
13 can ascertain exactly where it is. We look at all of
14 these patients to make sure that it's where we said it
15 is.

16 Next slide.

17 And that's another one just to show the
18 same thing.

19 That's really all I have to say. Thank
20 you very much.

21 MS. EBINGER: Good morning. I'm Kiara
22 Ebinger. I'm Senior Clinical Study Specialist with

1 Cochlear Corporation.

2 The effectiveness of the nucleus auditory
3 brain stem implant system in adults and young adults
4 with NF2 is supported by results from 60 ABI
5 recipients in the United States and 27 in Europe.
6 Each of these 87 individuals either have suffered or
7 will suffer complete bilateral deafness as a result of
8 acoustic tumor removal.

9 Participants in this investigation were
10 individuals 12 years of age or older who were
11 diagnosed with NF2. There were no preoperative
12 audiological criteria due to the life threatening
13 nature of CPA tumors and the elimination of residual
14 hearing by the tumor removal surgery.

15 Subjects were implanted with the ABI
16 during the first or second side tumor removal surgery,
17 and previous gamma knife therapy in the vicinity of
18 the cochlear nucleus was a contraindication to
19 participation in this trial.

20 Next slide, please.

21 During this investigation 92 subjects were
22 implanted with the ABI. Unfortunately two of those

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 subjects died from causes unrelated to the ABI prior
2 to device activation, leaving a safety sample of 90
3 subjects.

4 Of those 90, 30 were not included in the
5 effectiveness sample for the following reasons.
6 Thirteen individuals did not receive auditory precepts
7 at initial device activation. Ten subjects were
8 unable to attend the three or six month evaluation
9 usually due to health concerns or additional surgical
10 procedures. Six subjects had been recently implanted
11 and had not yet reached the three month evaluation,
12 and one subject was explanted prior to his three month
13 evaluation.

14 Sixty stimulable subjects with a minimum
15 of three months of device experience and whose data
16 were reported to us by February 9th, 1999, were
17 submitted in support of device effectiveness.

18 Next slide.

19 A few biographic characteristics of the
20 sample are shown on this slide. On average,
21 participants were implanted at about 33 years of age,
22 and one-third of them were implanted at the time of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 first side tumor removal surgery, and about two-thirds
2 of the group was female.

3 Next slide, please.

4 The study design and investigational
5 protocol are reviewed briefly here. As I mentioned,
6 there were no preoperative audiological criteria for
7 the trial. Comprehensive neurological and
8 audiological evaluations were conducted at initial
9 device activation, three, six, nine, and 12 months
10 after activation, and annually thereafter.

11 An extensive battery of recorded dependent
12 measures was used to assess a variety of auditory
13 skills, and patient questionnaires were also used to
14 assess patient satisfaction and perceived benefit from
15 the device.

16 Major findings from the trial are
17 summarized in the next few slides. Due to the fact
18 that patients with acoustic tumors have usable hearing
19 preoperatively, device effectiveness couldn't be
20 evaluated by comparing pre to post-op performance.

21 Instead, device effectiveness was
22 evaluated using the binomial statistic to compare each

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 subject's performance to chance performance on each of
2 a variety of speech perception tests.

3 Measures of lip reading enhancement also
4 were used to assess device effectiveness. To quantify
5 the improvement in speech perception when auditory
6 cues provided by the ABI were used in conjunction with
7 visual cues obtained from lip reading, the binomial
8 statistic was used to compare performance in those two
9 conditions.

10 The six month test interval was chosen as
11 the submission interval in support of device
12 effectiveness because it represented a considerable
13 amount of device use. As described previously,
14 patients with NF2 often have many health problems and
15 sometimes were unable to attend follow-up
16 appointments. So if six month data were not available
17 for a given subject, three month data were
18 substituted. This was felt to be a conservative
19 treatment since it represented less experience with
20 the device as opposed to more experience with the
21 device.

22 Although all 60 subjects were evaluated on

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 every measure -- sorry -- at three or six month
2 intervals, a few were not evaluated on every test. So
3 the number of patients tested varied slightly from 60.

4 Next slide.

5 Auditory performance was assessed using a
6 variety of recorded tests administered without visual
7 cues. Chance performance levels for each test are
8 indicated by the dotted lines on this slide.

9 The sound effects recognition test, or
10 SERT, is a closed set measure of an individual's
11 ability to identify common environmental sounds. This
12 test reflects an individual's connection to and
13 awareness of the auditory world around him. Chance
14 performance for the SERT is 25 percent.

15 As shown here, the group mean was well
16 above chance, at 54 percent correct, with individual
17 scores ranging from 13 to 83 percent. Eighty-two
18 percent of the subjects scored significantly above
19 chance, indicating that the ABI facilitates an
20 important connection between ABI recipients and the
21 world around them.

22 The MTS, or monosyllable trochee spondee

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 test is a closed set measure of word identification
2 that also can be scored to reflect correct
3 identification of the stress pattern of a test item.
4 Chance performance for stress pattern scoring is 33
5 percent. Mean performance for this group of ABI
6 recipients was well above chance at 75 percent
7 correct, with scores ranging from 21 to 100 percent.
8 Eighty-eight percent of the subjects scored above
9 chance on this measure.

10 The investigational protocol also included
11 two closed set measures of word identification, the
12 MTS-Word and NU-CHIPS tests. On the MTS-Word test, 80
13 percent of recipients scored significantly above the
14 8.3 percent chance level, with a group mean of 35
15 percent correct. MTS-Word scores ranged from eight to
16 88 percent correct for individuals.

17 For the NU-CHIPS words, chance performance
18 is 25 percent correct, and the mean NU-CHIPS score for
19 this group was 43 percent, with 67 percent of the
20 sample scoring significantly above chance. NU-CHIPS
21 scores ranged from four to 78 percent for each
22 individual.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Performance across these two measures of
2 word recognition are consistent and demonstrate that
3 ABI recipients are able to possess word identification
4 skills.

5 CID sentences were administered to assess
6 subject's open sentence recognition using sound alone.
7 This is obviously an extremely challenging task for
8 ABI recipients, and in general subjects didn't perform
9 as well on this measure, with a mean score of four
10 percent correct.

11 Next slide.

12 Individual results for CID sentences,
13 sound alone, are shown here. As you can see, many
14 subjects did score zero on this test. However, scores
15 ranged as high as 58 percent correct. Nine subjects
16 scored significantly above chance, and two subjects
17 scored over 50 percent correct.

18 Overall, these results are encouraging and
19 are consistent with those previously reported for
20 single channel Cochlear implant recipients.

21 Next slide.

22 Results of lip reading enhancement are

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 displayed here. Normal hearing individuals commonly
2 integrate visual cues with auditory cues in order to
3 understand a message, and obviously hearing impaired
4 individuals do the same, relying even more heavily on
5 visual cues.

6 For this reason, measures of the benefit
7 received from assimilating auditory cues with visual
8 cues are a very accurate indicator of everyday
9 functioning for a hearing impaired individual.

10 Three speech perception tests were
11 administered in multiple conditions, sound alone,
12 vision alone, and sound plus vision. The three
13 measures used to evaluate lip reading were the Iowa
14 medial vowel and consonant tests and the CUNY
15 sentences test, all of which are shown here.

16 For all three of these measures, sound
17 alone performance is represented by the gold bars,
18 vision alone by the bluish bars, and sound plus vision
19 by the lavender bars.

20 Lip reading enhancement is defined as the
21 improvement of scores when performance using the ABI
22 sound with lip reading or sound plus vision was

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealgross.com

1 compared to performance using lip reading cues alone
2 or vision alone.

3 On this slide that's represented by the
4 improvements seen between the blue bars and the
5 lavender bars.

6 Overall scores for the Iowa medial vowels
7 were higher than those for the consonants, probably
8 because vowels are easier to hear and to lip read.
9 Using sound alone, mean score for the vowels test was
10 25 percent correct. Using vision alone, the mean
11 improved to 65 percent, and when sound was used in
12 conjunction with vision, the mean score further
13 improved to 71 percent, with individual scores ranging
14 from 19 to 96 percent correct.

15 Fifty-five out of the 57 subjects scored
16 significantly above chance when using the ABI with lip
17 reading.

18 Because vowels are relatively easy to lip
19 read and many of these subjects lip read very well,
20 it's not surprising that just 11 percent of the
21 subjects demonstrated statistically significant
22 improvements with lip reading since the vision alone

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 scores were very good.

2 For the medial consonant test, the middle
3 three bars on this chart, the mean score in the sound
4 alone condition was 19 percent. Mean performance
5 improved to 38 percent using vision alone and to 52
6 percent when visual and auditory cues were combined.
7 Sound plus vision scores ranged from 30 to 85 percent
8 correct. One-third of the group showed statistically
9 significant enhancements in lip reading when using the
10 ABI for this test.

11 Fifty-eight subjects were tested with the
12 CUNY sentences using sound alone, the mean CUNY
13 sentence score was four percent and ranged from zero
14 to 57 percent correct. Using lip reading alone, the
15 mean increased to 31 percent, and when sound and lip
16 reading were used together, the performance mean
17 further improved to 54 percent, ranging in individuals
18 from seven to 96 percent correct.

19 Next slide.

20 Individual results for the CUNY sentences
21 are shown here. Vision alone scores are represented
22 by the lavender portion of each bar, and sound plus

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 vision is represented by the maroon portion of each
2 bar.

3 You can see just by the amount of maroon
4 on this slide that the ABI does provide a good deal of
5 benefit to these patients. Specifically, using the
6 ABI in conjunction with lip reading, 85 percent of
7 subjects improved significantly over their lip reading
8 alone scores.

9 In addition, when using both sound and
10 vision, every subject scored significantly above
11 chance. These results are very exciting since the
12 understanding of sentences when using both sound and
13 vision is a fundamental component of everyday
14 communication.

15 Next slide.

16 This slide summarizes results for CUNY
17 sentence scores over several years of device use.
18 Again, vision alone scores are represented by lavender
19 and sound plus vision by maroon. As you can see,
20 although subject numbers do decrease over time,
21 performance does remain very stable.

22 Notice that the vision alone scores or the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 lavender portion stays essentially the same, while the
2 auditory component of performance appears to improve
3 slightly over time. Without the ABI, this group of
4 ABI patients' communicative abilities would be
5 reflected just by the lavender portion of those bars.

6 Next slide.

7 As additional support for the efficacy of
8 the ABI system, European clinical trial results were
9 included in our PMA application and are summarized
10 very briefly here. Seventeen clinical trial patients
11 were implanted with the nucleus 22 21-electrode
12 device. Data are included for another ten pilot
13 subjects who are implanted with devices that differed
14 only very slightly from the clinical trial device.

15 Demographic characteristics of this group
16 were very similar to what we saw in the U.S. trial,
17 with the exception that there were more males in the
18 European group.

19 Statistical grouping of this data was
20 precluded by the fact that testing was conducted in
21 six different languages and using different methods
22 and materials.

1 Next slide.

2 Overall, the results of the European trial
3 were very consistent with the findings of the U.S.
4 trial. As shown here, the ABI provided most
5 recipients with the ability to recognize environmental
6 sounds, identify words in a closed set, achieve open
7 set sentence understanding in a few cases, and enhance
8 lip reading with auditory cues.

9 Next, please.

10 I think we can all grasp these performance
11 results from an intellectual perspective, but most of
12 us really can't conceive how these functional gains
13 really impact the quality of life for individuals with
14 NF2. So in an effort to better understand these
15 results from the patient's perspective, the
16 investigational protocol also included patient
17 questionnaires.

18 Next slide.

19 Performance questionnaire was administered
20 to patients at each evaluation beginning at three
21 months. This questionnaire assessed issues of daily
22 device use and perceived helpfulness in a variety of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 listening situations.

2 The next couple of slides summarize some
3 of the more salient questionnaire items. Data are
4 shown for 51 of the 60 subjects who completed the
5 questionnaire at either their three or six month eval.

6 Next slide.

7 When asked about daily device use,
8 respondents reported using the speech processor
9 between zero and 17 hours per day, with a mean for all
10 subjects of just over seven hours a day. Eighty
11 percent of the subjects who reported little or no
12 daily use were first side subjects with usable hearing
13 in the other ear.

14 Subjects who were implanted at the time of
15 second side tumor removal reported using the device
16 for an average of about ten hours per day.

17 Next slide.

18 Subjects were asked to rate the
19 helpfulness of the ABI in a variety of specific
20 listening situations using a six point rating scale,
21 where a rating of one indicated no help at all, and a
22 rating of six indicated very helpful.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 Overall subjects rated the ABI most
2 helpful in situations where auditory input could be
3 supplemented with visual cues. These situations,
4 including lip reading in a one-on-one situation,
5 listening to a familiar voice with lip reading in
6 quiet, identifying environmental sounds, and general
7 social interaction.

8 Next slide.

9 Subjects reported that the ABI was least
10 helpful for strictly auditory tasks where lip reading
11 cues were not available. These situations included
12 listening to the radio, listening to music, and
13 listening to either an unfamiliar or a familiar voice
14 and noise without the aid of lip reading.

15 In summary, results of the performance
16 questionnaire demonstrate that recipients rated the
17 ABI as being helpful in many social situations,
18 particularly when the information from the ABI could
19 be supplemented with visual cues.

20 Next slide.

21 The performance questionnaire reflected a
22 relatively short duration of experience with the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 device as patients completed it at either the three or
2 six month evaluation. In order to better assess long
3 term device use and patient outcomes, a final
4 questionnaire was administered. Responses to this
5 questionnaire reflected up to six years of device use.

6 Forty-four subjects returned the final
7 questionnaire prior to closure of the databases. In
8 general, despite fairly modest performance outcomes,
9 such as small number of recipients with open set
10 speech perception, patients reported enhanced quality
11 of life and improved communicative function when using
12 the ABI.

13 Next slide.

14 Seventy-five percent of respondents
15 reported wearing their devices daily. The remaining
16 25 percent were not regular users for a variety of
17 reasons, including implantation during first side
18 tumor removal with usable hearing in the other ear,
19 other acute health concerns or surgeries, and in a few
20 cases, lack of perceived benefit.

21 Eighty percent of respondents reported
22 that they received benefit from their ABI. When asked

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 about the decision to receive the implant, 84 percent
2 felt that the decision to get the implant was the
3 right one, and 73 percent reported that they would
4 recommend an ABI to others.

5 These results indicate that the majority
6 of respondents are satisfied with the ABI and with the
7 benefits they receive.

8 Next slide.

9 In conclusion, performance results
10 demonstrate that the ABI allow recipients to recognize
11 environmental sounds, identify stress patterns in
12 words, enhance their lip reading abilities, and in a
13 few cases, achieve open set, auditory only sentence
14 understanding.

15 Perhaps more importantly, results of
16 patient questionnaires confirm that these auditory
17 benefits can lead to significant and meaningful
18 improvement in quality of life for these individuals.

19 The results of the clinical trial as
20 reported to FDA in detail and summarized here support
21 the effective application of the nucleus 24 ABI in
22 individuals with NF2.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 these days.

2 It's important, we believe, to clearly
3 communicate the relatively high non-stimulation rate
4 as well as other possible adverse effects.

5 The claims that we proposed at our speech
6 perception claims in this section are based on the
7 individual results for the 60 effectiveness subjects.
8 We're also requesting permission to propose and use
9 four claims, I believe, that are based on
10 questionnaire data coming from the written responses
11 of 44 subjects.

12 Next slide.

13 The areas where we're asking for claims
14 are just Kiara has described. We believe that the
15 benefits of the ABI pertain to the identification of
16 environmental sounds, lip reading enhancement, open
17 set sentence recognition, as well as self-reported
18 benefit satisfaction and information regarding device
19 use.

20 Next slide.

21 The identification of environmental sound
22 claims that we've proposed are here. The first claim

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 is an aggregation of the single subject results
2 basically counting up the number of subjects who
3 demonstrated statistically significant performance
4 improvements compared to chance using the binomial
5 model as Martyn explained.

6 So 82 percent of the subjects, 49 of 60,
7 scored significantly above chance, which was 43
8 percent, on a recorded closed set test of
9 environmental sound identification. So most of the
10 sample is experiencing benefit, statistically
11 significant benefit in this area.

12 The second claim attempts to attach to
13 that some magnitude of the effect. So using the ABI,
14 subjects recognized 54 percent of common environmental
15 sounds on average and 65 percent of the sample, 39 of
16 60 subjects, recognized 50 percent or more of the
17 sounds. So not only do we want to tell patients that
18 they're very likely to improve; we want to give them
19 some kind of an index as to how much improvement to
20 expect.

21 Next slide. **

22 This is the same sort of format for the

1 lip reading enhancement claims. Eighty-five percent
2 of the tested subjects, 49 of 58, demonstrated
3 statistically significant improvements in open set
4 sentence understanding when using the ABI in
5 conjunction with lip reading.

6 Secondly, the average sentence recognition
7 score improved from 31 percent for lip reading alone
8 to 54 percent when subjects combined auditory
9 information from the ABI with lip reading.

10 Next slide.

11 This final claim concerns open set
12 sentence recognition. It is quite modest. Using
13 sound alone, just 12 percent of study participants
14 scored greater than ten percent on a difficult open
15 set test of sentence understanding.

16 We believe that a claim like this is
17 important, again, in order to give ABI patients some
18 context, something to which they can compare the
19 results that they're likely to hear about Cochlear
20 implant patients.

21 Next slide.

22 These are our questionnaire based claims.

1 Firstly, 61 percent of subjects, nine of 31, who
2 received the device following removal of a second side
3 tumor reported using the ABI on a daily basis for ten
4 or more hours.

5 Eighty percent, 35 of 40 of the
6 respondents reported receiving benefit from the ABI,
7 and 84 percent, 37 of 44, indicated that the decision
8 to receive the ABI was the right one, and lastly, 73
9 percent, 32 of 44, of the respondents would recommend
10 an ABI to others.

11 Next slide.

12 The last section of the insert that I'd
13 like to discuss pertains to training requirements. As
14 Dr. Brackmann and Dr. Hitselberger described earlier,
15 this isn't an easy surgical procedure. We believe
16 that physicians should be quite experienced in both
17 tumor removal procedures and Cochlear implant surgery.

18 Obviously they should be thoroughly
19 familiar with the anatomy of the fourth ventricle and
20 the ABI surgical procedure itself.

21 We strongly believe that the ABI project
22 is a team project, and that the implanting surgeon

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 work with an experienced team of professionals. These
2 professionals should include a neurosurgeon, a
3 neurologist, an audiologist, and lastly, an
4 electrophysiologist.

5 Next slide.

6 The implanting physician and ABI team will
7 be asked to attend a manufacturer sponsored training
8 program, and then secondly, to kind of further the
9 training program that they attend, we will be asking
10 that the newly implanting ABI surgeon invite a
11 designated consultant, someone identified by us to be
12 present and support his or first ABI surgery.

13 Next slide.

14 Lastly, with respect to post market
15 surveillance, PMS programs, as you know, are intended
16 to address long term questions regarding the safety
17 and effectiveness of a device. We don't believe that
18 PMS is indicated for the ABI system at this time.

19 The safety issues, we believe, are very
20 well characterized in the current sample. As you
21 know, there is a relatively high rate of non-
22 stimulations. However, the other complications are

•
NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 generally minor and are characteristic of Cochlear
2 implant surgery.

3 The longitudinal data that Kiara
4 summarized that is reported in quite more detail in
5 our submission to FDA clearly supports the stability
6 of the effectiveness outcomes over time. I'd just
7 like to remind you that this IDE has been open since
8 1994. Some of these patients have been implanted six,
9 seven years. They've contributed a lot to this
10 investigation and, to be honest, they're kind of tired
11 of coming.

12 The sample includes over 200 person-years
13 of experience, which again we believe will represent
14 the longitudinal effects of this device.

15 That concludes our presentation. We'll be
16 very happy to entertain your questions.

17 CHAIRMAN PATOW: Thank you very much.

18 I'd like to thank the presenters for their
19 very clear and organized presentations.

20 We now have some time for panel questions.

21 Dr. Hood? *

22 DR. HOOD: I had just a couple of

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 questions relative to the insert and the numbers. I'm
2 unclear about the environmental sounds.

3 CHAIRMAN PATOW: That would be a great
4 idea. If we could have the presenters come up to the
5 table where the microphones are, that would save some
6 time.

7 DR. HOOD: From the data presented, I
8 understand that the environmental sounds has a chance
9 level of 25 percent, and in the labeling it talks
10 about subjects scoring significantly above chance of
11 43 percent. I'm wondering if it's a different measure
12 that has a different chance level that's being used
13 for this claim.

14 MS. ARNDT: I'm so sorry, Linda. I didn't
15 quite hear the first part of the question.

16 DR. HOOD: Okay.

17 MS. ARNDT: Could you tell me one more
18 time?

19 DR. HOOD: It has to do with the
20 environmental sounds and the therapeutic claim that
21 talks about a chance score of 43 percent. In the data
22 that we were presented relative to the SERT, the

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 environmental sounds has a chance level of 25 percent.

2 MS. ARNDT: Right.

3 DR. HOOD: I'm wondering if there are
4 different measures contributing to this.

5 MS. ARNDT: No, I think you've identified
6 an error. Chance performance on the SERT is 25
7 percent.

8 DR. HOOD: Okay.

9 MS. ARNDT: So the results actually looked
10 better than we claimed.

11 DR. HOOD: Okay, and just one other minor
12 point along that line. On the questionnaire results,
13 you mentioned 61 percent of the subjects, but nine of
14 31, and I wondered if that's just another calculation?

15 MS. ARNDT: No, there we're very
16 specifically talking about device use in patients who
17 received -- I'm sorry, Linda. Tell me one more time.
18 I should have brought my presentation up.

19 DR. HOOD: Okay.

20 MS. ARNDT: You're referring to the?

21 DR. HOOD: The questionnaire results with
22 an n of 44 subjects.

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 MS. ARNDT: Yes.

2 DR. HOOD: And you're talking about the
3 fact that 61 percent of the subjects receiving the
4 device following removal of a second side tumor report
5 using it on a daily basis for ten or more hours, and
6 the number that you have describing the subjects that
7 contributed to that I think is just another dropped
8 number.

9 MS. ARNDT: It actually isn't. This
10 denominator is an n of 44. We had 44 recipients who
11 sent us back their questionnaires. Of these 44, 31
12 received the device following the removal or at the
13 time of the removal of their second neuroma. The
14 remaining patients received the device at the removal
15 of their first side tumor, which means they still had
16 normal or near normal hearing in the contralateral
17 ear.

18 So we were trying to describe device use
19 in the population of patients who replied on this
20 questionnaire for whom you would expect them to use
21 the device. They are profoundly deaf.

22 Does that make sense?

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS

1323 RHODE ISLAND AVE., N.W.

WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 We've addressed the other issue in another
2 section of the package insert called the clinical
3 considerations, and there we have described to
4 potential recipients that if they do receive an ABI at
5 the time that their first tumor is removed, because
6 their hearing is likely to be normal, that they're not
7 likely to use it very much until the time that the
8 second tumor is removed.

9 DR. HOOD: Okay.

10 DR. BRACKMANN: But it is 31 of 44, which
11 is 66, not nine of 31.

12 MS. ARNDT: Oh, I see. Well, if we have
13 an error, we'll correct that.

14 DR. HOOD: Okay.

15 MS. ARNDT: Long road to a little house.

16 DR. HOOD: Could I ask one other question?

17 CHAIRMAN PATOW: Yes.

18 DR. HOOD: This is relative to the magnet,
19 and am I clear that magnets are not used with this,
20 that the internal device -- that none of the patients
21 have the magnet on the internal device?

22 DR. BRACKMANN: That is correct. You have

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

1 the option. It's a removable magnet.

2 We do have a subset of patients who have
3 only bilateral acoustic tumors. There are families in
4 which the phenotype only produces bilateral acoustic
5 tumors so that you would have the option perhaps when
6 both tumors are removed of leaving the magnet in
7 place, but we have not done that.

8 We recommend removal of the magnet in all
9 cases. None of the device is, therefore,
10 ferromagnetic, and so it is MRI compatible with
11 removal of the magnet.

12 DR. HOOD: Okay, and I also just wanted to
13 be clear that the 60 patients that were entered into
14 the effectiveness data all are using the adhesive
15 method of connection rather than the magnetic.

16 MS. ARNDT: That's right.

17 DR. HOOD: Thank you.

18 CHAIRMAN PATOW: Other questions? Dr.
19 Woodson.

20 DR. WOODSON: Yes. This is Dr. Woodson.

21 Dr. Brackmann, you mentioned that in most
22 of the cases of the non-stimulation it was probably

NEAL R. GROSS

COURT REPORTERS AND TRANSCRIBERS
1323 RHODE ISLAND AVE., N.W.
WASHINGTON, D.C. 20005-3701

(202) 234-4433

www.nealrgross.com

1 due to cochlear nucleus distortion, which makes me
2 wonder if there were other factors, and I'm wondering
3 could you tell at the time of surgery that they were
4 not going to stimulate? Is there any factor that you
5 could identify maybe on the pre-op scans or any way of
6 identifying prior to putting an implant in that it's
7 going to be futile for a particular patient?

8 DR. BRACKMANN: Of the 14 non-stems, there
9 were seven in which an electrically about the EABR
10 could not be elicited intraoperatively. So that would
11 have predicted that they would have been non-
12 stimulators.

13 On the other hand, we have two or three
14 that had no intraoperative EABRs who were very good
15 performers. So we have felt insecure in not offering
16 it to patients despite the fact that the EABR is not
17 obtainable.

18 In the other seven where EABRs were
19 obtained intraoperatively, we have to believe that it
20 was device migration. In a couple of those cases,
21 we've identified a very "patchless, large, lateral
22 recess of the fourth ventricle which did not hold the