

at

AT

1

DEPARTMENT OF HEALTH AND HUMAN SERVICES
FOOD AND DRUG ADMINISTRATION
CENTER FOR DEVICES AND RADIOLOGICAL HEALTH

* * * * *

OPHTHALMIC DEVICES PANEL
101ST MEETING

* * * * *

Friday, July 20, 2001

9:35 a.m.

Main Conference Room
9200 Corporate Boulevard
Office of Device Evaluation
Rockville, Maryland

MILLER REPORTING COMPANY, INC.
735 8th Street, S.E.
Washington, D.C. 20003-2802
(202) 546-6666

PARTICIPANTS

Joel Sugar, M.D., Chairperson
Sara M. Thornton, Executive Secretary

VOTING MEMBERS

Michael Grimmett, M.D.
Janice M. Jurkus, O.D.
Alice Y Matoba, M.D.
Jose S. Pulido, M.D.
Jayne S. Weiss, M.D.

CONSULTANTS DEPUTIZED TO VOTE

Karen Bandeen-Roche, Ph.D.
Timothy B. Edrington, O.D.
Timothy T. McMahon, O.D.
Barry A. Weissman, O.D., Ph.D.
Karla S. Zadnik, O.D., Ph.D.

INDUSTRY REPRESENTATIVE

Marcia S. Yaross, Ph.D.

FDA

A. Ralph Rosenthal, M.D.

C O N T E N T S

Call to Order: Joel Sugar, M.D.	4
Introductory Remarks: Sara M. Thornton	4
Conflict of Interest Statement: Sara M. Thornton	7
Appointment to Temporary Voting Status: Sara M. Thornton	9
Open Public Hearing	
Professor Brien A. Holden, Ph.D.	10
Associate Professor Deborah F. Sweeney, Ph.D.	21
Dr. D. James Kerr	29
Division Update: A. Ralph Rosenthal, M.D.	36
Branch Updates	
Everette T. Beers, Ph.D.	37
Donna R. Lochner	41
James F. Saviola, O.D.	42
PMA 0010019	
Sponsor Presentation	
Alicia M. Plesnarski, RAC	48
John McNally, O.D., F.A.A.O.	54
Scott Robirds, O.D., F.A.A.O.	82
FDA Presentation	
James F. Saviola, O.D.	131
Myra K. Smith	131
Bernard P. Lepri, O.D., M.S., M.Ed.	131
Committee Deliberations	
Janice M. Jurkus, O.D., M.B.A.	142
Alice Y. Matoba, M.D.	148
Open Public Hearing	
Professor Brien A. Holden, Ph.D.	217
Closing Comments	220
Voting Options Read	225
Panel Recommendation	226

P R O C E E D I N G S

Call to Order

DR. SUGAR: I would like to call this meeting of the Ophthalmic Devices Panel to order. We will have introductory remarks from Sara Thornton.

Introductory Remarks

MS. THORNTON: Good morning and welcome to the 101st Meeting of the Ophthalmic Devices Panel. Before we proceed with today's agenda, I have a few short announcements to make. Bear with me. I would like to remind everyone to sign in on the attendance sheets in the registration area just outside the room here. All the handouts for today's meetings are available at the registration table.

Messages for the panel members, the FDA participants, information or special needs should be directed through Ms. Annemarie Williams or Mr. Demarc Thompson who are available in the registration area. If you need the phone number for someone to reach you out there, it is 301 443-8011.

In consideration of the panel, the sponsor and the agency, we ask that those of you with cell

1 phones and pagers either turn them off or put them
2 on vibration mode while in this room.

3 We ask that all panel meeting participants
4 speak into the microphone and give your name
5 clearly so that the transcriber will have an
6 accurate recording of your comments.

7 The next Ophthalmic Devices Panel Meeting
8 will be on Friday September 21, 2001. All
9 available information for that meeting will be on
10 the FDA Advisory Committee website within the next
11 few weeks. Should the September meeting be held
12 here, we will be pleased to be able to invite you
13 back to enjoy new carpeting and thorough painting
14 that have taken place in your absence.

15 Now, at this time, I would like to extend
16 a special welcome and introduce to the public, the
17 panel and the FDA staff four new panel consultants
18 who are with us today for the first time.

19 Dr. Timothy Edrington is a Professor of
20 Optometry and Chief of the Cornea and Contact Lens
21 Service at the Southern California College of
22 Optometry in Fullerton, California.

23 Dr. Timothy McMahon is as Associate
24 Professor of Optometry in the Department of
25 Ophthomology and Visual Sciences at the University

1 of Illinois at Chicago.

2 Dr. Barry Weissman is Professor of
3 Ophthalmology and Chief of the Contact Lens Service
4 of the Jules Stein Eye Institute and Department of
5 Ophthalmology at the UCLA School of Medicine.

6 Dr. Karla Zadnik is an Associate Professor
7 of Optometry and Physiological Optics at Ohio State
8 University College of Optometry and a Glenn A. Frye
9 Endowed Professor since 1999.

10 We greet you as special government
11 employees and welcome you to the panel table today.

12 To continue, will the remaining panel
13 members please introduce themselves beginning with
14 --

15 DR. SUGAR: Ralph, we can start with you.

16 DR. ROSENTHAL: Ralph Rosenthal. I am the
17 Division Director.

18 DR. WEISS: Jayne Weiss, panel member.

19 DR. GRIMMETT: Michael Grimmitt, Bascom
20 Palmer Eye Institute, Miami, Florida.

21 DR. MATOBA: Alice Matoba, Baylor College
22 of Medicine, Houston, Texas.

23 DR. JURKUS: Jan Jurkus, Illinois College
24 of Optometry in Chicago.

25 DR. SUGAR: Joel Sugar, University of

1 Illinois at Chicago.

2 DR. PULIDO: Jose Pulido, University of
3 Illinois, Chicago.

4 DR. BANDEEN-ROCHE: Karen Bandeen-Roche,
5 Johns Hopkins University, Baltimore.

6 DR. YAROSS: Marcia Yaross, Allergan,
7 Irvine, California and industry representative to
8 the panel.

9 MS. THORNTON: Thank you, panel. I would
10 like to note for the record and with regret that
11 Ms. Lynn Morris, our panel consumer representative
12 cannot be with us today. Earlier this week, she
13 fell and broke her ankle and is doing her best to
14 rest comfortably at home. We wish her well and
15 look forward to having her with us at the next
16 meeting.

17 I am your executive secretary, Sara
18 Thornton.

19 **Conflict of Interest Statement**

20 MS. THORNTON: I will now read the
21 conflict of interest statement for the record. The
22 following announcement addresses conflict of
23 interest issues associated with this meeting and is
24 made part of the record to preclude even the
25 appearance of an impropriety.

1 To determine if any conflict existed, the
2 agency reviewed the submitted agenda for this
3 meeting and all financial interests reported by the
4 committee participants. The conflict of interest
5 statutes prohibit special government employees from
6 participating in matters that could affect their or
7 their employer's financial interest.

8 However, the agency has determined that
9 participation of certain members and consultants,
10 the need for those services outweighs the potential
11 conflict of interest involved, is in the best
12 interest of the government.

13 Therefore, a waiver has been granted for
14 Dr. Karla Zadnik for her financial interest in a
15 firm at issue that could potentially be affected by
16 the panel's recommendations. The waiver allows
17 this individual to participate fully in today's
18 deliberation. A copy of this waiver may be
19 obtained from the agency's Freedom of Information
20 Office, Room 12A-15 of the Parklawn Building.

21 We would like to note for the record that
22 the agency took into consideration other matters
23 regarding Drs. Karen Bandeen-Roche, Timothy
24 Edrington, Timothy McMahon, Barry Weissman and
25 Karla Zadnik. These individuals reported past or

1 current interest in firms at issue but in matters
2 that are not related to today's agenda.

3 The agency has determined, therefore, that
4 they may participate fully in all panel
5 deliberations. In the event that the discussions
6 involve any other products or firms not already on
7 the agenda for which an FDA participant has a
8 financial interest, the participant should excuse
9 him or herself from such involvement and the
10 exclusion will be noted for the record.

11 With respect to all other participants, we
12 ask, in the interest of fairness, that all persons
13 making statements or presentations disclose any
14 current or previous financial involvement with any
15 firm whose products they may wish to comment upon.

16 **Appointment to Temporary Voting Status**

17 MS. THORNTON: I would like to read now
18 the appointment to temporary voting status.
19 Pursuant to the authority granted under the Medical
20 Devices Advisory Committee charter dated October
21 27, 1990 and as amended August 18, 1999, I appoint
22 the following individuals as voting members of the
23 Ophthalmic Devices Panel for this meeting on July
24 20, 2001; Dr. Karen Bandeen-Roche, Dr. Timothy
25 Edrington, Dr. Timothy McMahon, Dr. Barry Weissman,

1 Dr. Karla Zadnik.

2 For the record, these individuals are
3 special government employees and consultants to
4 this panel or other panels under the Medical
5 Devices Advisory Committee. They have undergone
6 the customary conflict of interest review and have
7 reviewed the material to be considered at this
8 meeting. Signed Dr. David W. Feigal, Jr.,
9 Director, Center for the Devices and Radiological
10 Health, June 28, 2001.

11 Thank you, Dr. Sugar.

12 DR. SUGAR: Thank you, Sally.

13 We will now move to the Open Public
14 Hearing. We have thirty minutes so I presume each
15 of the three presenters will limit themselves to
16 ten minutes and will start out their presentation
17 with a statement of any financial interest and who
18 is sponsoring their appearance here.

19 Dr. Holden?

20 Open Public Hearing

21 DR. HOLDEN: Thank you, Mr. Chairman and
22 thank you for the opportunity of saying a few
23 words. I have documented in the notes that were
24 distributed our commercial linkages. They include
25 royalty arrangements and intellectual property

1 development with a variety of corporations
2 including Ciba Vision.

3 The government of Australia makes it
4 mandatory for us, when we develop intellectual
5 property, to receive royalties and those royalties
6 are distributed according to the contract with the
7 government. We, indeed, have collaborative
8 projects which are intellectual property and
9 royalty generating with the organizations listed on
10 the slides.

11 What brought me to Washington by way of
12 money was money from my own organization. I am not
13 sponsored to speak here nor do I have any shares in
14 any company other than Tulstra, Australia. I have
15 never bought shares in the ophthalmic industry so I
16 don't gain to benefit in that way. My organization
17 certainly does gain to benefit from both consulting
18 and collaborative money.

19 [Slide.]

20 The main reason I am here is because this
21 is an extremely historic occasion for the
22 consideration of extended wear and the genius of
23 Otto Wichterle back in the '60's both predicted and
24 worked on daily exposables and extended wear in his
25 time.

1 [Slide.]

2 I thought I would show you what thirty
3 years of extended-wear research has done to me, Mr.
4 Chairman.

5 [Slide.]

6 The problem we have had is microbial
7 keratitis. It is the only serious adverse event as
8 defined by Stulka and others in the literature as
9 it can lead to significant loss of vision.

10 [Slide.]

11 It was another genius, Montague Ruben from
12 Moorfields Eye Hospital, that blew the whistle on
13 extended wear back in the early '70's in particular
14 in soft-lens extended wear for aphakic patients.

15 [Slide.]

16 There are a variety of studies that I have
17 listed in my handout. The classic one in 1989 was
18 Poggio, Glynn and Schine and colleagues where
19 ulcerative keratitis in extended wear was at the
20 rate of 21 per 10,000 people or 1 in 500 people,
21 five times greater than with daily wear.

22 [Slide.]

23 Many may not know that ten years later,
24 another landmark paper was published by Cheng et
25 al. in Lancet showing 1 in 500 or 20 per 10,000 was

1 still the norm for extended wear of current
2 hydrogen lenses at that time, although the
3 situation for daily wear seemed to have been
4 improved as there was an 8.3 times less risk with
5 current daily-wear lenses.

6 [Slide.]

7 In our studies over the last ten years or
8 so with low-Dk soft extended wear that have
9 involved about 2,278 eye years and about 1,000
10 wearers, we find a much higher rate of microbial
11 keratitis in those wearers at 2.5 times, Poggio, in
12 fact, in these prospective case-control studies.

13 [Slide.]

14 In fact, for every 191 patient years of
15 extended wear of current hydrogels, we find a case
16 of microbial keratitis. In fact, the survival
17 analysis shows us that, over time, the number of
18 people being affected is quite substantial. The
19 last point there is a patient of ours who we have
20 been following for five years, one of 37 patients,
21 in fact, who has just recently had microbial
22 keratitis.

23 [Slide.]

24 However, despite the issue of microbial
25 keratitis, the loss of best-corrected visual acuity

1 has relatively been an underplayed variable in our
2 understanding of these issues. In fact, if you
3 look at Cheng and Nilsson's paper, the loss of two
4 lines of best-corrected visual acuity is at a rate
5 of 1 in 40,000 contact-lens wearers.

6 With extended wear, it averages about 1 in
7 12,000 contact-lens wearers.

8 [Slide.]

9 This is in comparison, if you like, to
10 LASIK where the loss of two lines best-corrected
11 visual acuity is 1 in 32, some 300 to 1,000 times
12 higher than it is with contact lenses,

13 [Slide.]

14 Of course, that is not unexpected as
15 contact lenses have only really one really sight-
16 threatening side effect whereas LASIK, in fact, has
17 a number.

18 [Slide.]

19 If we look at the intraoperative
20 complications, the postoperative complications and,
21 indeed --

22 [Slide.]

23 -- the summary of the two that is recorded
24 in the handouts that I have prepared, an average of
25 3.2 percent is the literature figure these days for

1 significant complications following LASIK leading
2 to a 313 per 10,000 complication rate.

3 [Slide.]

4 Not coincidentally is the loss of two or
5 more lines of best-corrected visual acuity with
6 LASIK is also recorded in the literature averaging
7 around 3.1 percent or 1 in 32 people.

8 [Slide.]

9 Turning back to contact lenses, we knew
10 from 15 to 20 years ago and the research by many
11 people from the United States, Europe and
12 Australia, that the major problems that we have had
13 with infections have been somewhat related to the
14 sickness of the epithelium continuing with current
15 materials on extended wear.

16 The closed-eye environment is virtually
17 anoxic with current lenses for extended wear
18 leading to a thin, attenuated poorly metabolizing
19 epithelium. The adherence of pathogenic bacteria
20 is increased and if patient is in the circumstance
21 where they introduce massive numbers of bacteria,
22 infection can, indeed, result.

23 [Slide.]

24 So the hypothesis starting out some 15
25 years ago or maybe even earlier than that, 20 years

1 ago, was that hypoxia would provide a healthy
2 epithelium and better able to resist for the eye
3 infection.

4 [Slide.]

5 In 1994, George Mertz and I published what
6 we thought was necessary to avoid hypoxia with
7 contact lenses.

8 [Slide.]

9 Since both the Bausch & Lomb and Ciba
10 Vision lenses have been released for experimental
11 and clinical use around the world, in fact, this
12 data from Fonn shows the overnight swelling
13 response with high Dk soft is very low compared
14 with the lenses that are on market at the present
15 time.

16 [Slide.]

17 Perhaps more importantly, the ongoing
18 clinical indicators, particularly microcysts, show
19 that, compared with low Dk soft lenses, high Dk
20 soft lenses have virtually no microcyst response.

21 [Slide.]

22 A colleague of mine, Eric Papas, has shown
23 that, as you increase oxygen transmissibility to
24 the levels we see today, limbal redness actually
25 disappears.

1 [Slide.]

2 In fact, although the lighting is poor
3 here, we would see that vascularization of the
4 peripheral cornea, when patients are refitted with
5 high Dk soft, those vessels unfill, if I can use
6 that term.

7 [Slide.]

8 Of course, Dwight Cavanagh and colleagues
9 have documented with human epithelial cells the
10 decrease in adherence of Pseudomonas with the wear
11 of higher oxygen-permeability contact lenses.

12 [Slide.]

13 So what is our situation with regard to
14 the risk of microbial keratitis? We have been
15 looking at about 1,000 eye years of patients with
16 microbial keratitis being our number-one
17 requirement for these studies. As yet, we have
18 found no cases of microbial keratitis over these
19 1,000 eye years.

20 [Slide.]

21 When we look at the survival analysis of
22 the two, we are only at the p equals 0.09 stage for
23 significance of difference, but there is obviously
24 a difference in trend. That is promising, but it
25 is not conclusive.

1 [Slide.]

2 When we pool the data from B&L and Ciba
3 Vision premarketing and research studies --

4 [Slide.]

5 -- we get an eye-wearer figure of around
6 3,000 eye years.

7 [Slide.]

8 When we look at the figure from Cheng et
9 al., 48 microbial keratitises out of 24,000 Dutch
10 contact-lens as opposed to 0 out of 3,000. That
11 also, indeed, looks promising.

12 [Slide.]

13 In the marketplace, there has been an
14 influx of contact-lens wear of high Dk extended
15 wear. As yet, there is one report that we have
16 received and we are monitoring these things as
17 closely as we can, in the last week, in fact, of
18 microbial keratitis in the 55,000 wearers in the
19 United Kingdom.

20 [Slide.]

21 In Australia, high Dk soft has been on
22 marketplace for 24 months.

23 [Slide.]

24 In the first year, it captured 5 percent
25 of contact-lens wearers and it is actually doubling

1 every 12 months with two-thirds of the patients on
2 30 nights extended wear and one-third on daily
3 wear.

4 [Slide.]

5 Currently, 13 percent of all new patients
6 and 30 percent of refits are wearing high Dk soft
7 lenses, according to the data recently published by
8 Wood and Morgan.

9 [Slide.]

10 So, indeed, the penetration rate in
11 Australia of the contact-lens market is around 13
12 percent.

13 [Slide.]

14 There have been four events of microbial
15 keratitis in Melbourne reported recently seen at
16 the Victorian Eye and Ear Hospital. All were 16 to
17 22-year-old males. Maybe swimming was a factor.
18 Two occurred with each lens type on the market.
19 Three of them were culture positive but none were
20 Pseudomonas. Two of them resolved to 20/25 and two
21 had no effect on vision.

22 [Slide.]

23 If we look at that rate, we are talking
24 about 1 in 16,000 wearer years being an indicated
25 figure for microbial keratitis, MK, in Australia.

1 [Slide.]

2 Globally, there are about 250,000 wearers
3 of high Dk soft with about 175,000 patient years.
4 There are 9 MK case reports that we received. Four
5 have led to one line loss of acuity, three no
6 effect and two we don't have the data, one in
7 Italy, one in France, one in the U.S., four in
8 Melbourne, one in the U.K. and one in Norway.

9 At that rate, 9 in 175,000 wearer years
10 looks fairly promising compared with the previous
11 experience.

12 [Slide.]

13 If we take the worst case for Victoria,
14 20,000 high Dk soft-lens wearers in Victoria, four
15 that we know about and, perhaps, four that we
16 don't, we are still looking at a factor of some six
17 times less microbial keratitis per wearer year than
18 in low Dk soft lenses.

19 [Slide.]

20 So, globally, that is very promising.
21 What is even more promising is that there are yet
22 to be reported any cases of loss of visual acuity
23 of two lines or more of best-corrected visual
24 acuity in the 175,000 wearer years that have so far
25 existed around the planet.

1 [Slide.]

2 So where do we go from here? Microbial
3 keratitis is the only contact-lens serious adverse
4 event that is likely to occur with high Dk soft.
5 High Dk soft looks very promising but we need
6 continued postmarket surveillance targeted at the
7 annualized incidence of microbial keratitis
8 especially recording visual outcome. Such studies
9 need to collect that data.

10 [Slide.]

11 In addition, the world needs a gold
12 standard, properly controlled, scientifically valid
13 benchmark study of the prevalence and relative risk
14 of microbial keratitis and with colleagues around
15 the world, we are undertaking such studies.

16 Thank you very much for your attention

17 DR. SUGAR: Thank you, Dr. Holden.

18 Dr. Deborah Sweeney will now give the next
19 presentation.

20 DR. SWEENEY: Good morning.

21 [Slide.]

22 Thank you for this opportunity. I have no
23 commercial interest in any ophthalmic industry and
24 Professor Holden has already outlined the
25 commercial linkages of the CRCERT and CCLRU which I

1 am employed by and CRCERT and CCLRU have provided
2 the funding for my attendance here today.

3 What I hope to do briefly is talk to you
4 about what we feel is the value or the worth of the
5 development of these new high-Dk silicone materials
6 and what that means to both our patients and us as
7 practitioners.

8 [Slide.]

9 In surveys conducted at the CCLRU of
10 nearly 1500 patients that have either been wearing
11 contact lenses or are interested in contact-lens
12 wear, when questioned about their preferred mode of
13 wear, we can see overwhelmingly that patients are
14 interested in being able to wear their lenses where
15 they can sleep in a modality either in extended
16 wear or a continuous-wear basis.

17 [Slide.]

18 Other surveys have recorded what we know
19 as practitioners to be the case that contact-lens
20 wearers are very interested in being able to see in
21 the morning on awakening. 79 percent of the
22 patients in this survey had considered refractive
23 surgery but had not elected to have the procedure
24 and 65 percent report that their contact-lens care
25 and routine interferes with their lifestyle.

1 [Slide.]

2 A group of educators earlier this year
3 looked at the area of the patient wanting to feel
4 normal, to be without any correctional vision
5 problem and so rated their impression of the
6 average patient's desire for achieving continuous
7 vision with freedom of spectacles and over
8 70 percent rated this desire as very high.

9 [Slide.]

10 So, as our patients want a variety of
11 factors or needs to be met from their vision
12 correction, the primary two are comfort,
13 particularly with contact lenses -- they want to be
14 unaware of these lenses -- they want a no-fuss and
15 no-bother modality of vision correction.

16 Together with the LV Prasad Institute in
17 India and the CCLRU, we have conducted a number of
18 prospective clinical studies on both neophytes and
19 experienced in a range of modalities from spectacle
20 wear, daily wear, daily disposable, conventional
21 extended wear and of continuous wear.

22 [Slide.]

23 As part of these studies, as well as
24 collecting the clinical data, we also survey our
25 patients of their attitudes and administer

1 questionnaires regarding their attitudes and
2 satisfaction with both continuous wear, their
3 previous lens experience and their attitudes to
4 LASIK.

5 [Slide.]

6 In this group of 80 patients that have
7 experienced continuous wear for an average 12
8 months, when we ask these patients what they liked
9 most about being able to wear lenses on a 30-night
10 wear schedule, overwhelmingly, the main reason for
11 liking this modality is the convenience that it
12 offers as well as their ability to see in the
13 morning and comfort.

14 [Slide.]

15 This issue of convenience and what it
16 offers to our patients, when we look at the
17 different modalities here, is quite obvious. Here,
18 in daily wear, 30 percent, roughly, of the
19 patients, convenience is rated as the most likable
20 thing of their schedule and that rises extremely
21 high to when we get to 30-night continuous wear and
22 we see a rating of over 85 percent.

23 [Slide.]

24 When asked about their overall
25 satisfaction with the modality, 80 percent of our

1 patients rate their satisfaction when asked on a 1
2 to 100 scale where 100 indicates excellent
3 satisfaction as over 85 percent.

4 [Slide.]

5 They also, when asked to rate various
6 aspects of both convenience, safety, vision comfort
7 and how their eyes appear, the appearance or lack
8 of redness that is discernable with continuous
9 wear, they all give very high satisfaction ratings
10 for the performance of these lenses with this
11 modality.

12 [Slide.]

13 We have also asked a group of patients
14 that were previous daily wearers and have since
15 moved to continuous wear to look back at their
16 previous daily-wear experience and compare overall
17 satisfaction, convenience, vision, comfort, comfort
18 at end of day and just how clean their lenses feel.
19 For all these attributes, the patients rate their
20 overall satisfaction or their experience in
21 continuous wear as being significantly better than
22 their daily-wear experience.

23 [Slide.]

24 When we asked our patients what the
25 disadvantages, if any, were of 30 nights and

1 wearing lenses on this schedule, 50 percent of our
2 patients reported they saw no disadvantage in
3 wearing lenses in this way. We still have the
4 remaining problem, 13 percent rated dryness and
5 discomfort.

6 [Slide.]

7 Having experienced extended wear or
8 continuous wear for an average 27 months, the
9 majority of the patients, now 92 percent, want to
10 be able to sleep in their lenses either for
11 continuous-wear purposes, and that is over, now, 70
12 percent of patients or at least on an extended-wear
13 basis.

14 [Slide.]

15 In the studies that we conducted LVP and
16 CCLRU, our patients are on a 30-night schedule.
17 However, they are encouraged to remove their lenses
18 for an overnight break or temporarily for a clean,
19 rinse and reinsertion as needed. We also allow
20 them to use unit-dose saline for morning and night
21 if they wish.

22 When we look at the patient's success or
23 their ability to be able to achieve 30 nights wear
24 by looking at the number of nights of consecutive
25 wear which they achieve, we see in 82 percent of

1 all visits, patients are able to wear their lenses
2 consecutively for 28 to 30 nights and a further 12
3 percent are able to wear them for 21 to 28 nights
4 without needing any removal.

5 [Slide.]

6 This data here is from the 12-month visit
7 where we look at the percentage of patients who do
8 not remove their lenses at all for an overnight
9 schedule or an overnight removal outside their
10 schedule. At this visit, 68 percent do not take
11 their lenses out for an overnight break.

12 14 percent are taking them out once and a minor
13 percentage are taking them out more than once for
14 an overnight break during their 30-night schedule.

15 [Slide.]

16 As well as monitoring the number of
17 overnight removals across time and, as you can see
18 here, this does not change across the 30 months
19 that we have monitored these patients and it
20 averages that 71 percent of our patients are able
21 to achieve 30 nights of continuous wear with no
22 overall break.

23 We also rate, or collect information about
24 how many times they temporarily remove their lenses
25 for a quick rub and rinse as well as overnight

1 removal. Again, across this 30-month period, there
2 is no change and it averages 53 percent that are
3 not needing to take their lenses out at all during
4 that period.

5 [Slide.]

6 Collectively, the CCLRU clinicians have
7 been involved in development of extended wear and
8 extended-wear research now for over 25 years not
9 only with the conventional Dk materials, also with
10 the high-Dk rigid materials as well as silicon
11 elastomer and silicon hydrogels.

12 Despite this vast experience, we still, as
13 clinicians, feel uncomfortable about using low-Dk
14 extended wear even in the clinical trials that we
15 conduct, and the reasons, primarily, are because of
16 the problems with hypoxia, safety and infection and
17 the concerns of ocular redness. It is for these
18 reasons that we value the development of these new
19 high-Dk materials.

20 [Slide.]

21 Our patients are very enthusiastic about
22 both the convenience and freedom from spectacles
23 that continuous wear offers. When we surveyed over
24 200 of our patients and asked, "Have you ever
25 considered refractive surgery to permanently

1 correct your vision correction?" we found that 69
2 percent have.

3 [Slide.]

4 However, after they have worn continuous
5 wear, 30 nights continuous wear, we ask again those
6 143 patients that had considered refractive surgery
7 what they now prefer as a means of permanent vision
8 correction. Now, only 39 percent of those original
9 143 patients are considering refractive surgery and
10 the others are happy to remain with continuous wear
11 as their vision correction option.

12 [Slide.]

13 So, in summary, I would just like to say
14 that we believe the value and the worth of this
15 development of continuous wear and high-Dk
16 materials for our patients is that it offers a very
17 convenient modality for permanent vision correction
18 and for practitioners, the decreased hypoxic effect
19 is, of course, of great interest to us.

20 Thank you.

21 DR. SUGAR: Thank you, Dr. Sweeney.

22 The next presentation will be by Dr. James
23 Kerr.

24 DR. KERR: Good morning. I am in private
25 optometric practice in Saskatoon, Saskatchewan,

1 Canada. I do not now work for Ciba Vision nor have
2 I ever. They paid my way here, but I have no other
3 financial interest in Ciba and they have had no
4 input into my remarks.

5 They are based on my clinical experience
6 with this product and that clinical experience
7 began when I was involved in the Canadian clinical
8 trials of the Focus Night and Day lens beginning in
9 March of 1999. I fitted twelve patients according
10 to that protocol which involved using a
11 competitor's contact lens in one eye and the Focus
12 Night and Day Lens in the other.

13 It became immediately obvious to me that
14 this new product was superior to anything we had
15 used before, so much so that at the conclusion of
16 this study, all of the patients continued to wear
17 the Focus Night and Day lens on a 30 day-and-night
18 continuous-wear schedule.

19 The lens was then approved for 30-day
20 continuous wear in Canada in June of 1999. Since
21 that time, our office has ordered over 900 six-
22 month supplies of this lens. This represents
23 approximately 500 different patients. The majority
24 of these patients wear the lens on a 30-day
25 continuous-wear cycle.

1 Ciba then expanded the parameters of the
2 lens in April of this year. Up until that time,
3 the lens had one base curve and limited powers and
4 that limited the fittings that we could do with
5 this lens. When they expanded the parameters, we
6 had a much wider range of fitting, much wider range
7 of powers.

8 Since that time, the lens has simply taken
9 off in our practice. We have five doctors
10 prescribing it. All five are involved now.
11 Something like 75 percent of our 14-day disposable
12 lens wearers choose to switch to this lens when
13 they are advised of the features and benefits of
14 the lens.

15 Many, if not most, are skeptical partly
16 because we have always discouraged extended wear
17 and partly because most patients, in spite of our
18 opposition to extended wear, they have either
19 intentionally or otherwise slept with their lenses,
20 convention lenses, and they find them to either
21 stick or fog up or both. After learning that this
22 new product does not do this, most patients are
23 interested and, after a trial period of one month,
24 they are usually very enthusiastic.

25 Our experience to date with this lens has

1 been as follows: the lenses are now very
2 comfortable when they fit properly and, with only
3 two base curves, there are limitations to the
4 fitting. But we expect good comfort with the lens.
5 We have seen no corneal edema at all. We have seen
6 less surface deposition than most other lenses.

7 We have seen no neovascularization. Most
8 patients find the lenses do not dry as much as
9 other lenses. I have an asterisk here. I live in
10 Saskatchewan which is something like the Sahara
11 desert so most new products, most contact-lens
12 products, if they are dry at all, we have big
13 trouble with them. With this one, we have had no
14 difficulty there at all.

15 I haven't seen, to this date, any lens
16 stick to the cornea. Most patients wake in the
17 morning and either blinking or installation of a
18 wetting drop renders the lenses immediately
19 comfortable. We expect less limbal injection and
20 whiter eyes than any daily-wear lens presently
21 available. I think you have seen reference to
22 that. This is the third time in a row. These eyes
23 are whiter than any product we have had before.

24 We have seen no giant papillary
25 conjunctivitis to this point and we have seen no

1 cases of microbial keratitis or other significant
2 infections. I, personally, have seen two patients
3 develop a contact-lens-induced acute red eye. Both
4 of these patients responded well to topical
5 antibiotic steroid drops and were able to resume
6 continuous wear of the contact lenses without any
7 loss of vision or recurrence to this date.

8 Other complications have been minor but
9 include lens coating, dryness, lens awareness and
10 mucin balls. Our patients' acceptance of this lens
11 has been a surprise to us. Because of previous
12 product failures, there is a natural skepticism and
13 resistance to the concept of continuous wear. But,
14 as more and more patients are successful, the
15 demand is truly amazing.

16 The benefits are obvious as it affords to
17 patient who is handicapped by refractive error to
18 live their lives in a less complicated way with far
19 less risk of adverse events than any other form of
20 correction.

21 This correction is adjustable so that as
22 their eyes change, so can the correction. It is
23 safe. It is reversible and the complications that
24 do arise can be easily managed.

25 It is not without risk. I think it would

1 be unreasonable to expect that we will not see
2 ocular infections, perhaps some serious ones. But
3 it is risk management that we must consider.

4 We know these lenses pass more oxygen
5 along the cornea to maintain its natural resistance
6 to disease. We know that certain ocular pathogens
7 do not adhere to the corneal epithelium and they do
8 with conventional hydrogels. We know that
9 decreasing chemical damage to the epithelium caused
10 by current multipurpose solutions will increase the
11 cornea's resistance to disease.

12 We know that compliance with current
13 disposable protocols and cleaning protocols
14 increases the likelihood of infection with
15 conventional hydrogels. It is intuitive that such
16 significant improvements will decrease the risk of
17 ocular complications of contact-lens wear compared
18 to current systems. This certainly seems to be
19 born out in my clinical experience.

20 The real risk management, however, is in
21 comparing 30-day continuous-wear Focus Night and
22 Day lens to refractive surgery. I practice in a
23 city that is well advanced in refractive surgery.
24 We have three laser centers in a city of 200,000
25 people. Two have been operating for six years or

1 more and there is a very high public use of and
2 demand for this form of continuous vision.

3 In my own practice, I have over 600
4 patients who have had refractive surgery and the
5 results have been truly outstanding. But the
6 statistics do hold out and, in my practice, I now
7 have 18 to 20 patients who have had complications
8 resulting in permanently reduced best corrected
9 vision.

10 When this happens, it is, indeed,
11 permanent and irreversible. Since the
12 incorporation of the Focus Night and Day lens into
13 my practice, I have gone from sending ten patients
14 a month for refractive surgery to sending two
15 patients a month. None of the patients who have
16 chosen the contact-lens path have lost any vision
17 and, indeed, we have been able to adjust their
18 correction to provide optimum vision.

19 I feel that offering this alternative form
20 of continuous wear has, therefore, resulted in me
21 preventing vision loss in fifteen to twenty
22 patients who may have otherwise have opted for
23 refractive surgery while still providing them with
24 continuous vision.

25 The Focus Night and Day lens is the first

1 real improvement we have had in contact-lens
2 technology in over a decade. It has become an
3 important tool in our practice and I expect its use
4 to continue to grow to the point where we use this
5 sort of product in every contact-lens application.

6 I also believe it gets us much closer to
7 the point that when we remove a lens, we throw it
8 away, whether it is a single-day, 30-day or,
9 perhaps, some day, a year or more. Focus Night and
10 Day lenses truly revolutionized our practice. We
11 think it is a shame that this product is not yet
12 available for citizens of this country.

13 Thank you.

14 DR. SUGAR: Thank you.

15 If any panelists have questions for the
16 previous three presenters, we have a minute or two
17 to allow that.

18 Seeing none, we will move on to the open
19 committee discussion and the Division Update by Dr.
20 Rosenthal.

21 **Division Update**

22 DR. ROSENTHAL: I just have one issue to
23 announce to you, Mr. Chairman and Panel, and that
24 is that Nancy Brogden, the Deputy Director of this
25 Division, has been promoted to Director of the

1 Division of Diagnostic, Radiologic, Abdominal
2 Devices. That includes the rest of the body except
3 the eyes.

4 DR. SUGAR: And ears.

5 DR. ROSENTHAL: And ears, and nose and
6 throat. David Whipple has been appointed the
7 Deputy Director of our Division.

8 DR. SUGAR: Thank you.

9 We will now have the Branch updates. Dr.
10 Beers.

11 **Branch Updates**

12 DR. BEERS: Good morning. I am Everette
13 Beers. I am Acting Chief of Diagnostic and
14 Surgical Devices Branch. There have been no
15 personnel changes in the Branch since we updated
16 last November. We have approved a PMA and cleared
17 some 510(k)s.

18 I do want you to be aware that all
19 approvals and clearances, additional information
20 can be obtained on the FDA websit, fda.gov/cdrh.
21 You have to figure it out yourself from there.
22 Sometimes we can get there and sometimes we can't.

23 For PMA approvals, P930016, Supplement 12,
24 VISX LASIK Hyperopic Astigmatism for up to 5.00
25 diopter sphere and up to 3.00 diopter cylinder was

1 approved April 27, this past April 27.. That's all
2 we had on PMA approvals.

3 For 510(k)s, I did want you to be aware of
4 some of the 510(k)s. This panel does not see
5 510(k)s but those are our less risky devices but we
6 frequently have very forward-looking and cutting-
7 edge-technology types of devices in these areas
8 that we call 510(k)s.

9 The first one is K01199 cleared in June
10 2001. That is the Bausch & Lomb Proview Eye
11 Pressure Monitor, formerly the Fresco Phosphene
12 Tonometer. It was cleared for over-the-counter
13 home use. The tonometer is utilized on the closed
14 eyelid and requires a subjective response of the
15 perception of the phosphene which is placed on the
16 eye, on the closed eye. When you see a phosphene,
17 a little spring tells you what your eye pressure
18 reading is.

19 Another one I wanted you to be aware, that
20 we have put up on the CDRH website a Keratome LASIK
21 Guidance. Previously, keratomes have not been
22 allowed to state that they were for use with LASIK.
23 Now, we are saying that it is allowed to say that
24 the keratome can be used for LASIK. Keratomes are
25 class I devices. Lasers are class III devices. It

1 gets all very confusing, but we have put a guidance
2 up there that has changed previously when class I
3 devices were not allowed to advertise or be
4 indicated for class III indication of LASIK.

5 Finally, in the 510(k) area, wave-front
6 analysis autorefractometers, or aberrometers, are
7 exempt with limitations. Exemption means that you
8 do not have to submit a 510(k) for these prior to
9 marketing.

10 The product code for these devices is NCF.
11 These exempt 510(k) devices do not have to submit
12 premarket notification to FDA but, in accordance
13 with Section 513(i)(1)(E) of the Food, Drug and
14 Cosmetic Act, these exempt aberrometers must carry
15 the warning in their labeling that the safety and
16 effectiveness of using the data from this device,
17 whatever it is, have not been established for
18 determining treatments involving higher-order
19 aberrations of the eye such as coma and spherical
20 aberrations.

21 You can also see 510(k) K000637 for the
22 limitations on this device. If you are not
23 familiar with Wave Front autorefractometers, I
24 think I mentioned in the handout for the open
25 session a little bit more about these devices. In

1 general, they use a laser beam reflected from the
2 retina to determine distortions through the entire
3 visual system of the eye.

4 These aberrations include sphere, cylinder
5 and axis and the higher-order aberrations such as
6 coma and spherical aberrations. So they are used
7 like regular refractometer to get your sphere,
8 cylinder and axis. You can also get some other
9 readings, but they are not allowed to use those to
10 do those higher-order aberrations for refractive
11 treatments.

12 Finally, I think most of you are aware of
13 our LASIK websit, www.fda.gov/cdrh/lasik. We have
14 had tens of thousands of hits on that website and
15 we found that it has been very useful to consumers
16 as well as many practitioners.

17 Finally, I know that some of you, as
18 practitioners and also being on the panel,
19 occasionally receive questions from consumers
20 regarding something that is up with the FDA. You
21 really should forward those consumer questions to
22 the Office of Health and Industry Programs, the
23 Division of Small Manufacturers Assistance, DSMA.

24 You can have them call 800 638-2041 or
25 they can send in an e-mail request to

1 dsma@cdrh.fda.gov. That information should also be
2 in the open session handout that is in your
3 package.

4 Are there any questions?

5 DR. SUGAR: Go ahead, Marcia.

6 DR. YAROSS: Not a question, but I would
7 really like to commend the agency on the guidance
8 document on the LASIK indication for keratomes. I
9 think that that was really the clearest application
10 of least burdensome and I believe it is much
11 appreciated.

12 DR. BEERS: Thank you.

13 DR. SUGAR: Thank you, Dr. Beers.

14 Now Donna Lochner will speak for the
15 Intraocular and Corneal Implants Branch.

16 MS. LOCHNER: I would like to announce the
17 PMA approvals since the last panel meeting. First,
18 Staar Surgical P000026 AquaFlow Collagen Glaucoma
19 Drainage Device, Model CGDD-20, was approved on
20 July 12. This PMA was reviewed by the panel in
21 November, 2000.

22 The next two PMAs were not reviewed by the
23 panel because we felt that there were no new issues
24 of safety and effectiveness presented. The first
25 one, Pharmacia P990080 for the CeeOn Edge Foldable

1 UV-Absorbing PC IOL Model 911A was approved on
2 April 5.

3 Last, Anika Therapeutics, P00046, which
4 was a licensing PMA in which Bausch & Lomb provided
5 reference rights to P810025 which is Amvisc sodium
6 hyaluronate was approved April 18. This means that
7 Anika has approval to distribute and manufacturer
8 the Amvisc sodium hyaluronate under their label.

9 At this time, Anika did not request
10 distribution under the Anika label. Instead, they
11 received approval for Staar Surgical Company to
12 distribute the product as Staarvisc II sodium
13 hyaluronate.

14 That concludes my updates.

15 DR. SUGAR: Thank you.

16 Next, Jim Saviola will talk on the
17 Vitreoretinal and Extraocular Devices Branch.

18 DR. SAVIOLA: Thank you, Dr. Sugar. Good
19 morning, everybody. There are a few clearances and
20 PMA approvals that I wanted to inform you about
21 this morning. I had neglected in my prepared
22 remarks to mention a website that we were involved
23 in developing recently. I thank Dr. Beers for
24 jarring my memory on that.

25 About two months ago, there was a website

1 posted on the Center web that addresses questions
2 and answers regarding purchasing contact lenses
3 from the Internet or from other sources other than
4 from an eye-care practitioner. That is something
5 that we developed in response to inquiries we were
6 getting regarding prescription dispensing and
7 things of that nature. So if people are interested
8 in that, I would refer you to our website for that.

9 In the class II area for 510(k)
10 clearances, the first area I would like to discuss
11 is orthokeratology lens clearances. On February
12 28, 2001, we cleared the Paragon Fluroperm 151 for
13 daily-wear orthokeratology. That K number was
14 010109. The labeling for that product includes
15 reference to a previous Paragon study involving the
16 Fluroperm 60 material, so that is where that data
17 came from for that new clearance.

18 Polymer Technology received a clearance
19 for the Polymer Tech Boston EO lens and also for
20 the Polymer Tech Boston Equalens II, both on
21 February 16. Those K numbers are K003932 and 3933.
22 Those two had included references to the Contex
23 AirPerm clinical study that was conducted by
24 Contex.

25 Polymer had also received a clearance for

1 their XO material back in August of 2000. So, with
2 all these new clearances, there are now a total of
3 six orthokeratology lenses cleared. As you see,
4 some of them are based on original clinical studies
5 such as the Context AirPerm and also the Paragon
6 Fluroperm 60. Others are using data as a reference
7 within the context of determining equivalency and
8 being able to do that in class II which is
9 something you can't really do in class III.

10 In the lens-care product area, I told you
11 last meeting about Opti-Free EXPRESS Multi-Purpose
12 Disinfecting Solution manufactured by Alcon and how
13 they received a clearance for the No-Rub care
14 directions. Their first clearance was in July,
15 2000. Those are for lenses replaced for 30 days or
16 less followed by a second clearance in October of
17 2000 to remove the 30-day limitation to include
18 soft lenses prescribed on any replacement schedule.

19 We now have two more care products that
20 have received a "no-rub" clearance for lenses
21 replaced 30 days or less. K003252 cleared on
22 February 21 for Allergan Complete Multipurpose
23 solution and K003345 cleared March 26 for Ciba
24 Vision's AO Sept One-Care peroxide solution which
25 has a surfactant added to the peroxide.

1 With all of these clearances, there is
2 still wording in the labeling to advise users that
3 additional products or procedures such as rubbing
4 their lenses may be recommended by the eye-care
5 practitioner.

6 In the class II area, I had one PMA to
7 inform you of, and that is Vistakon (lenefilcon a)
8 soft hydrophilic contact lens which was approved on
9 February 16, 2001. That is now indicated for daily
10 wear and for extended wear up to seven days. As
11 Everette mentioned, some of these products are not
12 reviewed by the panel and seven-day extended lenses
13 are one of those that, in class III, we do not
14 refer for panel recommendation and review.

15 I neglected to report down the PMA number
16 for that, so if anybody is curious, I can ask Stan
17 Rakowsky or any of the other representatives of
18 Vistakon in the audience this morning.

19 Thank you very much.

20 DR. SUGAR: Thank you. Are there any
21 questions from the panel of the Division Chiefs?
22 Branch Chiefs; sorry. Sorry for the promotion.

23 Dr. Pulido?

24 DR. PULIDO: Dr. Saviola, have any of the
25 orthokeratology lenses come before panel?

1 DR. SAVIOLA: Not for the daily-wear
2 indication. For extended wear, overnight wear, we
3 plan to take those, at least the first one, for
4 panel.

5 DR. SUGAR: Dr. Yaross.

6 DR. YAROSS: I did have one question for
7 Ms. Lochner. Can you provide an update on the
8 status of the reevaluation of the age indications
9 for IOLs?

10 MS. LOCHNER: For those that may not know
11 about this issue, FDA has been doing basically --
12 we have been doing some research, basically a
13 metaanalysis of the literature as well as working
14 with Dr. Apple's group on postmortem globes with
15 the Academy on their outcomes base to compile data
16 to present a case for lowering the age indication
17 for IOLs to adults instead of age 60 and over.

18 We have actually completed substantially
19 the body of the work, the actual analysis, and this
20 is being prepared in the hopes of publication.
21 Right now, this publication, this draft
22 publication, is being reviewed by the authors
23 within FDA, the Academy and Dr. Apple's group to
24 follow.

25 Then we plan to submit this article for

1 publication. It is our hope that once this article
2 is published, sponsors could use this as a
3 reference of valid safety and effectiveness data to
4 support a lowered age indication.

5 So, in summary, the article has been
6 drafted. It is in the "being-reviewed" stage prior
7 to publication.

8 DR. YAROSS: Thank you.

9 DR. SUGAR: Jim?

10 DR. SAVIOLA: There was another answer I
11 forgot to give, too, Dr. Pulido. The very first
12 ortho-K lens we had for Contex we did refer out for
13 panel homework assignment to one of the panel
14 members.

15 DR. BEERS: Regarding the LASIK indication
16 for keratomes, I should mention that the keratomes
17 must meet certain requirements before they are
18 allowed to use that LASIK indication. So you
19 should look at the website at that guidance to
20 determine whether or not their keratome meets those
21 requirements.

22 DR. SUGAR: Thank you, Chiefs.

23 PMA P010019

24 We are now going to move on to the
25 discussion of the PMA at hand today, PMA P010019.

1 This will begin with a one-hour presentation by the
2 sponsor. I would like to remind everybody to state
3 their name before speaking so that the scribes can
4 have this in the transcript.

5 MS. PLESNARSKI: Good morning.

6 [Slide.]

7 My name is Alicia Plesnarski. I am a
8 regulatory specialist in the Global Regulatory
9 Affairs Group of Ciba Vision Corporation. I have
10 been with the company for about ten years and, for
11 most of that time, I have been the regulatory
12 project leader on the project team for the PMA
13 device.

14 Today, I am very proud and excited to be
15 here as part of this team. We are here to present
16 and discuss Ciba Vision's PMA P010019 for See3
17 (lotrafilcon A) Soft Contact Lenses with an
18 indication for up to 30-night extended wear.

19 [Slide.]

20 My presentation will be brief. I will
21 talk a little bit about our company and the PMA
22 device and then introduce the rest of the team.

23 First of all, Ciba Vision is a eye-care
24 unit of Novartis. We began as a small start-up
25 company in 1980 and have grown into a global

1 corporation involved in research and development,
2 manufacturing and marketing of ophthalmic products.
3 Our corporate headquarters are located in the
4 suburbs of Atlanta, Georgia and today our company
5 makes and markets contact lenses, lens-care
6 products and intraocular lenses and we maintain
7 FDA-registered manufacturing facilities on three
8 continents.

9 [Slide.]

10 Regarding the PMA device under
11 consideration today, the See3 lotrafilcon A soft
12 contact lens for up to 30-night extended wear is
13 classified as a class III medical device. The lens
14 material, lotrafilcon A, is a 24 percent water, 76
15 percent fluorosilicon-containing hydrogel which is
16 surface treated.

17 As a low water, nonionic polymer, this
18 lens material falls into FDA group 1 and, while the
19 lens has many physical and optical characteristics
20 that are similar to other soft contact lenses, one
21 extraordinary feature of this lens material is that
22 it has oxygen permeability of 140 delivering an
23 oxygen transmissibility of about 175 for a -3.00
24 diopter lens with a center thickness of 80 microns.

25 [Slide.]

1 To date, the lenses haven't been marketed
2 in the U.S. Ciba Vision has obtained FDA 510(k)
3 daily-wear marketing clearance in May of 1997.
4 Outside the U.S., the lenses are marketed under the
5 trade name Focus Night and Day in product packaging
6 that bears the CE mark.

7 In early 1999, a global market
8 introduction began and the product was launched in
9 many countries in the European Union, in Canada and
10 in Australia. Today, the lens has over 250,000
11 wearers in over 40 countries.

12 [Slide.]

13 Now, while soft contact lenses have been
14 on the market for over thirty years, much of the
15 exciting advancement in contact-lens material
16 properties has occurred more recently. The early
17 '90's marked the beginnings of a strong commitment
18 and targeted initiatives by industry to develop 30-
19 night continuous wear as a safe and effective
20 vision-correction option.

21 In terms of research and development of
22 next-generation contact-lens materials, the
23 progress in this area can be followed in the dozens
24 of scientific articles published in professional
25 journals regarding high-Dk lenses and extended

1 wear.

2 [Slide.]

3 Regarding our PMA and development of the
4 SEE3 soft contact lens, this project was initiated
5 in the early '90's and our goal was to develop and
6 market a noninvasive, safe and effective and
7 convenient 30-night extended-wear soft contact
8 lens. We are talking now of a period of about ten
9 years and, over the course of lens development,
10 there have also been some significant developments
11 in the regulatory area for medical devices.

12 With passage of the Medical Devices
13 Directives and CE marketing requirements in the
14 European Union and revision to the FDA GMP
15 Regulation to include design controls, the SEE3
16 lotrafilcon A lens became one of our first projects
17 to proceed under a formalized design-control system
18 compliant to both the FDA quality-system regulation
19 and ISO 9001 Quality Systems requirements for
20 design controls.

21 Before we move ahead, I wanted to mention
22 some commonly used terms you will be hearing this
23 morning. SEE3, lotrafilcon A and Focus Night and
24 Day are the project name, lens material name and
25 trade name for the contact lenses. The phrases

1 extended wear and continuous wear may be used
2 interchangeably and we mean no differences in these
3 phrases.

4 [Slide.]

5 In just a moment, we are going to be
6 moving on to the clinical findings, but I did want
7 to mention the lens has undergone a comprehensive
8 series of nonclinical testing to support product
9 safety. Some of those tests are listed on the
10 slide, but the actual list of testing exceeds those
11 recommended by the FDA in 1989 and 1994, Contact
12 Lens Guidance Documents and includes additional
13 physical-chemical testing, biocompatibility studies
14 as well as analysis of worn lenses.

15 Wherever possible, the methods conformed
16 to the applicable ISO or ANSI standards for
17 contact-lens testing.

18 [Slide.]

19 The important findings from all
20 nonclinical testing are that the lenses are
21 nontoxic and biocompatible. They are stable and
22 compatible with lens-care solutions. They have
23 material properties which are consistent with or
24 better than other soft contact lenses and these
25 properties remain unchanged after lens wear.

1 The results of all nonclinical tests
2 support the safety of lotrafilcon A lenses for
3 their intended use.

4 [Slide.]

5 At this point, I would like to introduce
6 the rest of our team. Presenting today, and up
7 next, will be Dr. John McNally who will provide
8 information on the clinical study design and
9 results. After John, Dr. Scott Robirds will talk
10 about product labeling and our proposed postmarket
11 study protocol.

12 [Slide.]

13 Also with us today and available to help
14 with questions and other information are Dr. Curtis
15 McKenney from our Research Clinic who has been on
16 the SEE3 project since its beginnings and Dr. Gary
17 Cutter, a biostatistician who worked with us on a
18 consultant basis regarding study design and
19 statistical analysis.

20 In addition, the president of our lens
21 business, Stuart Heap, is also here with us today.
22 On behalf of Ciba Vision, Stuart authorized payment
23 of our travel expenses to Washington and we are
24 hopeful he is going to do the same for our return
25 tickets back to Atlanta this afternoon. Stuart

1 will have some closing remarks later today.

2 That concludes my presentation and I thank
3 you for your time and attention. Up next, Dr. John
4 McNally.

5 DR. McNROBALLY: Thank you, Alicia.

6 [Slide.]

7 Good morning. My name is John McNally.
8 Today, I have the pleasure of presenting the
9 culmination of over a decade of extended-wear
10 research carried out by many hands from around the
11 world. I started my own interest in extended-wear
12 research some twenty-five years ago in the
13 laboratories of Dr. Mandell at U.C. Berkeley School
14 of Optometry.

15 I have since been with Ciba Vision for
16 twenty years, continuous and extended years I might
17 say, serving in various clinical, regulatory and
18 research management positions. I am currently the
19 head of continuous-wear research programs.

20 [Slide.]

21 This morning, I will briefly touch on some
22 of the background information regarding the
23 product. Then I will provide an overview of the
24 results of the clinical trial and provide some
25 comments in response to questions we have received

1 from the reviewers thus far.

2 [Slide.]

3 Here, I would like to reemphasize three of
4 the distinguishing properties of the lens material
5 that may be of importance for our discussions this
6 morning. Of course, the high oxygen permeability,
7 the low water content, the nonionic nature of the
8 material and the modulus which, for the panel's
9 reference, is higher than many soft lenses on the
10 market but is not unlike a number of contact lenses
11 that have been on the market for many years.

12 [Slide.]

13 The oxygen permeability of lotrafilcon A
14 is due to the siloxane content of the material.
15 Unlike hemabased hydrogels which require increases
16 in the water content to increase the oxygen
17 permeability, as shown in the curve on the bottom
18 of this illustration, it is obvious to see that the
19 lotrafilcon A polymer, shown here in the upper left
20 quadrant of the graph, is a departure from that
21 principle and a clear breakthrough in terms of
22 oxygen permeability.

23 [Slide.]

24 Critical to the unique nature and
25 performance of this polymer as well was the

1 discovery of the requirement for a continuous
2 hydrogel phase allowing the movement of ions
3 through the lens which is then responsible for or
4 related to the lens movement and the maintenance of
5 ion mobility. We have included a paper describing
6 this work in your panel packet.

7 [Slide.]

8 In the early phases of clinical
9 development, we studied several of the important
10 performance outcomes required for successful
11 extended wear, namely overnight corneal swelling,
12 bacterial colonization, lens-surface cleanliness.
13 In the panel packet, we have included summaries of
14 this work or published articles, when available.

15 I will briefly review the results of these
16 three.

17 [Slide.]

18 In a study of overnight corneal swelling
19 published by Fonn and coworkers, the SEE3 lens
20 produced a mean corneal swelling of 2.7 percent
21 overnight compared to 8.7 percent for the Acuvue
22 control clearly demonstrating one advantage of the
23 increased oxygen transmissibility.

24 [Slide.]

25 In a study of bacterial colonization of

1 the lenses during wear published by Keay and
2 coworkers, SEE3 lenses were aseptically removed
3 from the eye after 30 nights of continuous wear and
4 compared to Acuvue lenses sampled after six nights
5 of continuous wear. There were no significant
6 differences in the number of sterile samples, as
7 shown here, nor in the amount of types of bacteria
8 found, thereby showing no increased bacterial
9 colonization over the 30-day period.

10 [Slide.]

11 In a clinical study conducted at Ciba
12 Vision, lenses were retrieved for analysis of
13 protein buildup. SEE3 lenses were retrieved after
14 30 nights of continuous wear and Acuvue lenses
15 after six nights of continuous wear.

16 In this and similar studies, the SEE3 lens
17 made of the nonionic lotrafilcon polymer shows
18 remarkably less protein buildup than the control
19 lens, in this case that of lotrafilcon which is an
20 ionic polymer.

21 [Slide.]

22 After these early studies and prior to the
23 launch of the product in 1999, we completed an
24 international safety and effectiveness trial. The
25 rates of adverse events in that trial are presented

1 here. No statistical difference was found between
2 the SEE3 lens after 30 nights of continuous wear
3 and the control lens, Acuvue, at six nights of
4 continuous wear.

5 These rates for adverse events are similar
6 to the rates found in the U.S. trial that I will
7 discuss in just a few minutes.

8 [Slide.]

9 As you have heard, we launched the product
10 internationally in 1999 and currently there are
11 approximately 2.5 million lenses in the
12 marketplace. From that, we estimate that we have
13 approximately 250,000 wearers representing a
14 cumulative experience of approximately 100,000
15 patient years.

16 These numbers are updated from those
17 included in your packet and these are the current
18 numbers and represent our best knowledge. We have
19 had five cases of potential infectious keratitis
20 reported to us. I use the word "potential" because
21 of differing definitions by practitioners around
22 the world. But, nonetheless, these were severe
23 adverse events.

24 Based upon this information, our best
25 estimate for infectious keratitis is 5 in 100,000

1 patient years although we realize that this will
2 not hold up to epidemiological scrutiny.

3 Earlier this year, we added a second base
4 curve based upon feedback from the marketplace as
5 well as our findings in our clinical trials. We
6 also added plus lenses and high minus lenses at the
7 same time.

8 [Slide.]

9 Overall, feedback from the international
10 marketplace has found the Focus Night and Day
11 product to offer a desirable alternative for those
12 seeking the convenience of around-the-clock vision
13 correction.

14 We have also had numerous anecdotal
15 reports of less dryness and less redness from
16 wearers. The lens offers flexibility both in terms
17 of wearing regimen as well as the ability to easily
18 adjust refractive correction as required and it has
19 been particularly well received in the higher
20 refractive powers.

21 [Slide.]

22 So now to the results of the safety and
23 effectiveness study in the United States. After
24 briefly reviewing the study design, I will present
25 the key results and the conclusions and some

1 elements of clarification required to address the
2 reviewer's questions that we have received.

3 The evidence presented in the PMA packet
4 is in support of the indications being sought, in
5 particular the wearing schedule indication of up to
6 30 nights of continuous wear and the reduction in
7 contact-lens dryness symptoms.

8 [Slide.]

9 The objective of the study was stated as
10 follows: to determine whether the SEE3 lens when
11 worn for up to one month extended wear and replaced
12 on a monthly basis performed as well as the Acuvue
13 control lens when worn for up to one week extended
14 wear and replaced on a weekly basis.

15 [Slide.]

16 This was one of the largest prospective
17 contact-lens studies conducted to date in support
18 of safety and effectiveness. It was a one-year
19 open-label randomized controlled clinical trial
20 involving 59 investigative sites. As I mentioned,
21 there were differences in both the wearing schedule
22 and replacement frequency with SEE3 being worn for
23 up to a month extended wear and replaced monthly
24 and the control lens weekly extended wear and
25 replaced weekly.

1 Additionally, during the study, the SEE3
2 lens was available in a single base curve whereas
3 the control lens was available in multiple base
4 curves.

5 [Slide.]

6 The primary safety endpoint was
7 infiltrates of grade 3 or greater or any
8 infiltrates with overlying fluorescein staining.
9 This is a conservative endpoint as contact-lens
10 infiltrates are not usually infectious in nature
11 and rarely lead to reduction in visual acuity.
12 However, this endpoint may serve as a threshold
13 surrogate for an infectious ulcerative keratitis
14 or, as it is commonly referred to in the contact-
15 lens industry, microbial keratitis.

16 Microbial keratitis is a rare corneal
17 complication and is therefore prohibitive to study
18 in a premarket trial and is better suited to
19 postmarket evaluation such as we will propose
20 later.

21 [Slide.]

22 The primary effectiveness endpoints were
23 the visual acuity and the wearing time achieved
24 with the contact lenses.

25 [Slide.]

1 The sample size for the study was based
2 upon a noninferiority statistical design. This
3 type of design allowed us to test at the alpha 0.05
4 level whether the SEE3 lens was worse than the
5 control by a specified amount referred to as the
6 equivalence margin.

7 For the safety endpoint just discussed,
8 the equivalence margin was set at 5 percent and the
9 estimated endpoint rates were set at 8.6 percent
10 for reasons discussed in the clinical protocol and
11 report that you have received.

12 A noninferiority study design has the
13 advantage that we specifically set out to prove
14 that you are not different by a certain amount
15 unlike the statistical outcome from many studies
16 where equivalence is claimed because a difference
17 wasn't detected.

18 The null hypothesis, then, is that the
19 rates are different by 5 percent or more and
20 noninferiority would be demonstrated by
21 statistically rejecting this null hypothesis.
22 Although this study design preceded the draft FDA
23 extended-wear guidance for extended-wear lenses, it
24 closely aligns with the statistical principles in
25 that guidance and the examples provided as well.

1 [Slide.]

2 Under these assumptions just discussed,
3 the sample size at 80 percent power would be 389
4 subjects per group. In order to have greater than
5 80 percent power, we increased our sample size to
6 500. This provided us a robust study design that
7 would provide 87 percent power at the estimated
8 rate of 8.6 percent as shown on this graph.

9 It also would provide adequate power
10 across a wide range of potential outcomes as shown
11 here. We were also satisfied at the gut-feel
12 clinician level with the maximum SEE3 rate we might
13 observe in the trial, shown in the bottom row, and
14 still reject the null hypothesis.

15 [Slide.]

16 Our final enrollment target was set at
17 700. We included a 15 percent allowance for the
18 possible inability to fit all subjects with a
19 single base-curve parameter in the SEE3 product and
20 made a further 20 percent adjustment for dropouts
21 that may occur over a year's period of observation.

22 [Slide.]

23 So, to the results. Today, I will discuss
24 the enrollment and the accountability, the
25 discontinuations, the primary safety endpoint and

1 adverse events and the effectiveness endpoints of
2 visual acuity and wearing time and the results
3 regarding contact-lens dryness.

4 [Slide.]

5 697 SEE3 and 698 control subjects were
6 enrolled. 39 SEE3 and 17 control subjects were
7 unable to be dispensed. The difference between the
8 two here is due to 20 SEE3 subjects that did not
9 achieve an acceptable fit at this time. In total,
10 658 SEE3 eyes and 681 control eyes were dispensed
11 equating to 1,316 and 1,362 eyes, respectively.

12 [Slide.]

13 The demographics were representative of
14 the contact-lens-seeking population and the two
15 groups were nearly identical. Subjects were
16 actually dispensed in the power ranges you see
17 listed here, +6.00 to -6.00 for SEE3 and +4.50 to -
18 6.50 for Acuvue. Approximately 95 percent of the
19 subjects in each group were myopic, as is typical
20 of the current contact-lens-wearing population.

21 [Slide.]

22 The groups were also well matched in terms
23 of previous contact-lens-wear experience as shown
24 in this chart. It is important to note that
25 approximately 60 percent of the subjects were new

1 to extended wear and, thus, were not established
2 veterans.

3 [Slide.]

4 A larger proportion of the SEE3 group, 175
5 versus 102, were discontinued from the study, many
6 for reasons we had foreseen, as I will discuss
7 next. Accountability was excellent in the study
8 with complete data available on 96 percent of the
9 subjects' dispensed lenses.

10 [Slide.]

11 The four biggest differences and reasons
12 for discontinuations are seen here; discomfort,
13 lens fit, biomicroscopy and acuity with contact
14 lenses. Let me discuss each one of these
15 separately.

16 [Slide.]

17 The difference in discontinuations for
18 discomfort was largely driven by the difference in
19 the first week and overall in the first month.
20 After the first month, the rates were similar, as
21 you can see graphically depicted on the bottom of
22 the slide.

23 [Slide.]

24 The same is true for discontinuations for
25 lens fit, as you can see in the chart, and again in

1 the graph at the bottom. All but one of the SEE3
2 discontinuations for lens fit were due to
3 unacceptably flat or loose-fitting lenses.

4 [Slide.]

5 Discomfort and fit discontinuations are
6 very likely related. With the SEE3 lens, when the
7 lens is too flat relative to the cornea, the lens
8 edge will lift or even buckle, as you see here.
9 This, of course, is an extreme case and this
10 subject would likely have been discontinued for
11 lens fit.

12 However, when the edge lift is more subtle
13 or sporadic, then it may not be observed by the
14 investigator during biomicroscopy but is evident to
15 the wearer by lid sensation or discomfort at the
16 area of edge lift. We have addressed this both in
17 the fitting guide, as Dr. Robirds will explain, and
18 in a subsequent development of a second base curve.

19 [Slide.]

20 In response to reviewer questions, we
21 examined various factors concerning the discomfort
22 discontinuations. We found no correlation to
23 corneal curvature, refractive power or lens fit
24 with only a very slight trend towards steeper
25 corneas and towards higher myopia. This lack of

1 correlation can perhaps be explained by the fact
2 that, on the steepest corneas, the lens fit was
3 obviously flat and these subjects were discontinued
4 for lens fit.

5 This would leave, then, only a trend for
6 these other factors as they relate to discomfort
7 discontinuations.

8 [Slide.]

9 Reviewers were also interested in the
10 investigators' decisions to discontinue subjects
11 for biomicroscopy and especially asked about the
12 severity of the findings. Four SEE3 subjects were
13 discontinued at the first event, one subject for a
14 peripheral ulcer or CLPU in the second week of
15 continuous wear with the infiltrate and staining
16 grades as shown here.

17 Two of these subjects were discontinued
18 for infiltrative keratitis, where I have IK listed,
19 because the event occurred in the first month of
20 wear and the investigator recommended against
21 continuing. One subject with a previous history of
22 Thygeson's was discontinued shortly into the study
23 due to a reoccurrence felt by the investigator to
24 be unrelated to the product.

25 Three subjects were discontinued because

1 the event listed in the table in your report was
2 their second event. One was a peripheral ulcer and
3 two with infiltrative keratitis. The other
4 discontinuations for biomicroscopy were four with
5 papillary conjunctivitis, three of the four from
6 one investigative site, and all subjects with
7 previous history of papillary conjunctivitis, and
8 five other subjects for early microacyst rebound or
9 dimple veiling.

10 [Slide.]

11 Two control subjects discontinued for
12 events because they occurred in the second week,
13 one for infiltrative keratitis and one for herpes
14 keratitis. This latter was listed in the report
15 table as intraepithelial keratitis. One subject
16 was discontinued at the second occurrence of
17 episcleritis and one additional control subject,
18 with a peripheral ulcer, was discontinued for
19 "other ulcer" and more appropriately should have
20 been included in this listing.

21 [Slide.]

22 Several questions were raised regarding
23 discontinuations for contact-lens acuity. As
24 mentioned in the report, we encountered an issue in
25 our packaging design causing a small percentage of

1 lenses to adhere to the package, thus distorting
2 the optics as you see here in this photo. Although
3 we adjusted the packaging design, we did not
4 replace the clinical inventory.

5 All of the discontinuations for acuity
6 were in the first three months of the study. After
7 that, if a wearer experienced substandard vision
8 when they put in a new lens, they simply replaced
9 it with another one from the pack.

10 A question was also asked about engraving
11 and deposits and I will address that at this time
12 since you can see that the lens is engraved in this
13 photograph. The surface treatment of the lens is
14 applied after the engraving and, therefore, the
15 engraving presents no problems for tear-film
16 deposits.

17 [Slide.]

18 So, as I have explained, we found that the
19 majority of the differences in discontinuations
20 occurred in the first month and many for the fit
21 reasons we had foreseen. We have provided guidance
22 for this in the labeling as we will discuss.

23 [Slide.]

24 Now the results regarding the primary
25 safety endpoint.

1 [Slide.]

2 3.1 percent of the control group and 5
3 percent of the SEE3 group experienced an endpoint
4 infiltrate. You remember this is infiltrates grade
5 3 or greater or infiltrates with any overlying
6 staining. These unadjusted rates were not
7 statistical different.

8 [Slide.]

9 From these rates, we performed a survival
10 or life-table analysis that would account for all
11 subjects' time in the study and allow us to better
12 estimate annualized rates for the safety endpoint.
13 This life-table graph is in the report and shows
14 the survivors or, as we say, those not voted off
15 the island for experiencing an endpoint infiltrate.

16 [Slide.]

17 From that analysis, we obtained the
18 estimated annualized rates of 3.3 percent for the
19 control and 6.1 percent for the SEE3 lens. As
20 pointed out in the report, this is a conservative
21 estimate for the control rate since two control
22 peripheral ulcers were not included in this
23 statistical analysis.

24 One ulcer occurred at six months but was
25 treated by a non-study ophthalmologist over the

1 holiday season. The ulcer was later confirmed by
2 the investigator by the presence of a corneal scar.
3 However, since no data regarding the infiltrate,
4 itself, was provided, we did not include this
5 subject in this calculation.

6 A second ulcer was seen in another subject
7 in the control group at the 12-months visit.
8 However, since the visit occurred at 378 days,
9 including it in the life-table analysis would
10 greatly overestimate the control rate since so few
11 subjects were still in the study at that time, at
12 378 days.

13 Still, based on the noninferiority test I
14 outlined earlier in this presentation, we calculate
15 the p-value to be 0.0465 allowing us to reject the
16 null hypothesis and demonstrate noninferiority.

17 [Slide.]

18 In response to the reviewers' questions,
19 we examined various factors concerning the
20 incidence of endpoint infiltrates. We found no
21 correlation to refractive power, corneal curvature
22 or lens fit.

23 [Slide.]

24 In the clinical report, we characterized
25 the endpoint and analyzed various risk factors. I

1 note here only a few observations. The infiltrates
2 were mostly paracentral and limbal and few were
3 central. Subjects with a history of a previous
4 event were at a higher risk for having an event and
5 there was a trend for a higher rate in smokers
6 although it was not statistically significant in
7 our study.

8 [Slide.]

9 Our findings, then, in the primary safety
10 endpoint were as follows: SEE3 was found to be
11 noninferior to the control by the equivalence
12 margin defined in advance and no subjects lost best
13 corrected acuity with any endpoint infiltrate. We
14 will provide guidance from our findings in the
15 labeling.

16 [Slide.]

17 Now I will briefly cover overall adverse-
18 event rates and discuss other eyes that required
19 treatment during the course of the study. The
20 primary safety endpoints just discussed were all
21 considered adverse events and thus are also
22 included in the overall rates that follow.

23 [Slide.]

24 In line with the draft guidance for
25 extended-wear lenses, we classified adverse events

1 as serious, significant or nonsignificant using
2 examples provided in that guidance. Roughly, these
3 categories can be thought of as follows: serious
4 adverse events are potentially sight-threatening
5 events. For contact lenses, this would be optical
6 axis or infectious ulcers. Based on the guidance,
7 we also included any events in this category that
8 had the presence of any anterior chamber reaction.

9 Significant events are not directly sight
10 threatening but are usually treated to preclude
11 potential escalation or other sequelae.

12 Nonsignificant events are those that are typically
13 managed through temporary removals of the lens or
14 other palliative procedures.

15 [Slide.]

16 It is probably most meaningful to look at
17 the cumulative rates. The rates are listed on this
18 chart with the cumulative rate, shown on the
19 bottom, being 9.4 percent for SEE3 and 8.3 percent
20 for the control. If you remember, these rates were
21 very similar to those that I presented earlier for
22 the international safety and effectiveness trial.

23 Neither this cumulative rate nor any of
24 the rates shown here were statistically different.
25 Further details regarding events were included in

1 the report in the panel packet.

2 [Slide.]

3 In summary, regarding adverse events, we
4 found no difference in incidence between the two
5 groups, no cases of microbial keratitis and no loss
6 of best-corrected visual acuity. The rates of
7 these events will be provided in the labeling.

8 [Slide.]

9 A statistical difference was found in the
10 proportion of eyes requiring management for
11 contact-lens-induced papillary conjunctivitis or
12 CLPC on this chart. For the SEE3 subjects, 1 of
13 the 59 investigative sites reported 7 of the 30
14 papillary conjunctivitis subjects. All of these
15 seven had had a previous history or CLPC.

16 [Slide.]

17 We found no correlation with surface
18 deposits, lens fit, corneal curvature or refractive
19 power. However, the location in the early onset in
20 the SEE3 group suggest a possible mechanical origin
21 for these cases. Subjects with a previous history
22 of papillary conjunctivitis were at higher risk in
23 this trial and, in the labeling, we addressed the
24 potential increased risk of CLPC.

25 [Slide.]

1 Now I will discuss the primary
2 effectiveness endpoints.

3 [Slide.]

4 The visual acuity results can be
5 summarized as shown here. I have already mentioned
6 the best-corrected visual-acuity results. The
7 visual acuity with contact lenses worn remained
8 within two lines of baseline for 98 percent of the
9 evaluations over the course of the study. 83
10 percent of these evaluations were 20/20 or better
11 and, although not shown here, 99 percent were 20/30
12 or better.

13 In approximately 90 percent of the
14 evaluations, subjects rated the vision 8 or higher
15 on a 10-point scale.

16 [Slide.]

17 We evaluated wearing time in several ways.
18 First, we collected the prescribed wearing time.
19 This was the wear schedule assigned by the
20 investigator based on the case history and clinical
21 findings for each subject throughout the course of
22 the study and this was recorded at each visit.

23 Next, with the assistance of a diary, we
24 also collected data from the subjects about the
25 time of each removal and the reasons for that

1 removal. In this section, I will also address a
2 question from the reviewers regarding the
3 relationship of wearing time to adverse events.

4 [Slide.]

5 Here I think it is best to summarize the
6 results first since some clarification is needed
7 based upon review comments. Regarding the
8 prescribed wearing time, no subjects were
9 permanently prescribed less than a full indication
10 in either group. However, prescribed wearing times
11 were temporarily reduced in order to manage signs
12 or symptoms.

13 [Slide.]

14 Two tables about prescribed wearing time
15 were included in the text and I will briefly
16 explain the data from the one-month visit table.
17 This table shows 91.2 percent of the SEE3 subjects
18 and 93.9 percent of the control subjects were
19 assigned the full indication at this one-month
20 visit. The remaining subjects were temporarily
21 assigned shorter wearing schedules to manage
22 whatever was happening at that particular time.

23 Over the course of the study, the
24 prescribed wearing time at scheduled visits was at
25 full indication for more than 90 percent of the

1 subjects in the study at the time of the visit.

2 [Slide.]

3 From the subjects, we learned that lenses
4 were removed overnight for a variety of reasons
5 including, of course, the scheduled removal but
6 also for symptoms or problems the subject was
7 experiencing or as needed for sickness or other
8 demands in their life. The average wearing time was
9 based upon the period between overnight removals.

10 [Slide.]

11 Although not highlighted in the text,
12 Table 13A and 13B, the trend analysis profile,
13 recorded the average wearing time over the course
14 of the study. You can see that after the first
15 month, the average wearing time for SEE3 was 26 to
16 27 days which consisted primarily of many subjects
17 at 30 nights and a smaller group temporarily at
18 shorter times.

19 [Slide.]

20 This chart presented in the report is a
21 compilation of all reported wearing intervals from
22 all subjects, and this includes months 1 through
23 12. For the completed subjects, 67.2 percent
24 represents the percentage of time the 22 to 31 was
25 recorded as the maximum wearing time in a month.

1 Over the course of a year, a single subject may be
2 counted in several different categories on this
3 chart depending on how they were wearing the lenses
4 in that month.

5 Remember that all subjects wore the lens
6 for full indication except for temporary periods.
7 Basically, 67.2 represents the overall patient
8 months on the completed subjects recorded at 22 to
9 31 nights and 88.1 percent represents the number of
10 patient months at continuous intervals greater than
11 seven months.

12 As you would expect, in the discontinued
13 group, the wearing times were not as long since
14 they were having difficulties with the lenses and
15 ultimately discontinued. There may be alternative
16 ways of representing the wearing time data in the
17 labeling that the panel may prefer.

18 [Slide.]

19 We were asked to evaluate the relationship
20 between wearing time and adverse events for the
21 SEE3 lenses. This study design does not allow us
22 to do that in a direct dose-response fashion since
23 we did not have groups assigned at various wearing
24 schedules. Remember that over 90 percent were
25 described at the full indication over the course of

1 the study.

2 However, to address the root concern of
3 the question, we looked at our data in several
4 different ways. We calculated, for each subject,
5 their own individual average of reported
6 consecutive nights slept in the lens and looked for
7 a relationship between this and adverse events. We
8 found no increased risk with the increased average
9 consecutive nights slept in the lens.

10 We also looked at the consecutive months
11 of 30-night wear prior to the event; that is,
12 whether the subject had worn the lenses for one or
13 two or ten consecutive months of the 30-night
14 regimen. Here we also found an increased risk for
15 the increasing months at the full indication.

16 Finally, we looked at the consecutive days
17 of wear in a given lens at the time of the event;
18 that is, whether the lens had been worn for five or
19 ten or 20 consecutive nights at the time of the
20 event and we also found no increased risk with
21 increasing nights of continuous wear of a given
22 lens.

23 [Slide.]

24 We also summarized the reasons for the
25 shorter wearing time reported by the subjects. In

1 addition to the temporary reductions prescribed by
2 the investigators, the reasons for unscheduled
3 overnight removals were summarized in table 25 and
4 26 and temporary daytime removals in tables 27 and
5 28.

6 For both groups, the main reasons for
7 unscheduled overnight removals were eyes needed
8 rest, irritation or allergy. For daytime removals,
9 for both the test and control group, the main
10 reasons were to clean or for irritation or for
11 dryness. The multiple other reasons are listed in
12 the tables in your report.

13 [Slide.]

14 In summary on wearing time, the prescribed
15 and reported wearing schedules were predominantly
16 the full indication. Symptoms, problems or
17 lifestyle needs led to a temporary reduction in
18 wearing time and, as analyzed in the SEE3 group, we
19 were unable to show evidence of increased adverse
20 events with increased wearing time.

21 [Slide.]

22 As the final part of my section, I will
23 present the results supporting the finding of less
24 contact-lens-induced dryness. This finding was
25 supported by data gathered in the case history, in

1 the subjective questionnaires as well as the other
2 subject diaries.

3 [Slide.]

4 Text table 19 in the report showed the
5 symptoms reported to the investigator at each
6 visit. While dryness remains the most often
7 reported symptom in both groups, we found a
8 statistically relevant decrease in the reported
9 symptoms of dryness with SEE3 in both the completed
10 and discontinued subjects which you see highlighted
11 here.

12 [Slide.]

13 In the subject questionnaire, we also had
14 statistically fewer reports of dryness upon
15 awakening with the SEE3 lens. I have graphically
16 represented the data from all patients presented in
17 table 17A through D of the report and you can
18 clearly see the shift towards less problems with
19 dryness.

20 [Slide.]

21 Finally, and probably most important of
22 the three, in the completed subjects, we found
23 fewer unscheduled overnight removals for reasons of
24 dryness in the SEE3 group compared to the control
25 group as now highlighted here in the text table 25.

1 This difference was statistically significant with
2 p equals 0.02.

3 [Slide.]

4 Based upon the consistency of these
5 findings from three different sources, the SEE3
6 lenses demonstrated reduced dryness symptoms
7 compared to the control. This is consistent with
8 our international experience as well.

9 [Slide.]

10 So, as my concluding slide, we feel that
11 the scientific evidence presented in this PMA
12 application provides reasonable assurance that
13 safety and effectiveness have been demonstrated for
14 the requested indications.

15 [Slide.]

16 I will now turn the presentation to Dr.
17 Scott Robirds who will discuss the proposed
18 labeling and the postmarket protocol.

19 Thank you.

20 DR. ROBIRDS: Thanks, John.

21 [Slide.]

22 Good morning. I head up the Global
23 Clinical Affairs Group at Ciba Vision. I have had
24 the pleasure of working in clinical research and
25 regulatory affairs at Ciba for the past fifteen

1 years. Prior to working at Ciba, I was an
2 associate in a contact-lens practice for three
3 years. That is as far back as I care to go.

4 I will be presenting a proposed product
5 labeling and a summary of the postapproval study
6 for the SEE3 PMA.

7 [Slide.]

8 I will start with the product labeling.
9 As many of you know, the FDA has provided the
10 contact-lens industry with guidance documents that
11 provide helpful direction related to product
12 labeling. The guidance for package inserts,
13 practitioner fitting guide and patient instruction
14 booklet are very comprehensive and the majority of
15 the proposed labeling for the SEE3 product is
16 consistent with these guidance documents.

17 So I am just going to focus on the
18 elements of the proposed labeling that are unique
19 to our product and, in some cases, are a departure
20 from the published FDA guidance for contact-lens
21 labeling.

22 As you can see here, the product name is
23 Focus Night and Day and the product description
24 portion of the labeling, a summary of lens
25 properties, is presented. This is the same

1 information presented earlier by Ms. Plesnarski.

2 [Slide.]

3 The product description also includes the
4 proposed approval range of lens parameters and the
5 parameters that will be initially available. You
6 can see the available parameters are a diameter of
7 13.8, base curves of 8.4 and 8.6, and powers will
8 range from +6.00 to -10.00 in either quarter or
9 half steps dependent upon the power selected.

10 [Slide.]

11 All of this information will be present in
12 the package insert. The practitioner fitting guide
13 will present only the available lens parameters
14 seen here in gold.

15 [Slide.]

16 There are four indication statements that
17 were submitted with the proposed labeling. I will
18 just go through these. The first one is fairly
19 straightforward and deals with vision correction
20 and states that, "Focus Night and Day soft contact
21 lenses are indicated for the optical correction of
22 refracted ametropia in phakic or aphakic persons
23 with nondiseased eyes with up to approximately 1.50
24 diopters of astigmatism."

25 [Slide.]

1 The second relates to wearing time. Here
2 we stated that, "Focus Night and Day may be worn
3 for daily or extended wear for up to 30 nights as
4 recommended by the eye-care professional."

5 [Slide.]

6 The third indication relates to
7 replacement intervals and lens-care systems. Here
8 we say that, "Lenses should be replaced every month
9 and when removed between replacements must be
10 cleaned and disinfected with a chemical
11 disinfecting system before reinsertion."

12 [Slide.]

13 The fourth indication states that, "Focus
14 Night and Day lenses may reduce dryness symptoms
15 that are present with regular hydrogel soft
16 lenses." This claim is driven by findings in our
17 FDA study as presented earlier by Dr. McNally.

18 In the reviewers' comments, there was a
19 concern that this claim may be interpreted as being
20 applicable for use with patients having aqueous
21 tear deficiency or other pathological dry-eye
22 conditions. Our intent was not to claim an
23 indication of dry-eye relief in patients with
24 pathological dry eyes but to claim a reduction in
25 dryness symptoms secondary to routine contact-lens

1 wear. However, we acknowledge that we should
2 clarify this in the language of the claim.

3 Another review commented that since our
4 testing was confined to only one type of control
5 lens, the data could not support such a broad claim
6 for all regular hydrogen soft lenses. This is a
7 valid comment and the claim should be modified to
8 account for this.

9 [Slide.]

10 In the warning section, a portion of the
11 standard language about ulcerative keratitis has
12 been deleted, as you can see here. Our rationale
13 for deleting this section is that we have not found
14 this wording to be fully applicable to Focus Night
15 and Day as we have just heard in Dr. McNally's
16 presentation.

17 We would propose adding a statement such
18 as, "The incidence of microbial keratitis with
19 extended-wear lenses is approximately 20 per
20 10,000," or, alternatively, "Not all individuals
21 can wear lenses for up to 30 nights continuously.
22 Individual wearing times should be determined in
23 consultation with your eye-care practitioner."

24 We believe these types of statements would
25 address one of the reviewer's requests that a

1 statement be added that not all patients can
2 tolerate continuous wear.

3 [Slide.]

4 In the precaution section, we have added
5 three statements that relate to suitability as a
6 contact-lens candidate. The added text is seen
7 here in yellow. These precautions are added to the
8 standard contact-lens labeling as a result of our
9 FDA trial and relate to patients with a history of
10 acute inflammatory reactions, giant papillary
11 conjunctivitis or ocular allergies.

12 Subjects with histories of these
13 conditions were at a higher risk for repeated
14 occurrence of the condition compared to subjects
15 without such histories.

16 [Slide.]

17 Now, in the adverse-event section, we have
18 added a chart that calls out annual rates for
19 selected events as seen during the FDA trial.
20 Placing results of clinical trials in product
21 labeling is routine for pharmaceutical agents in
22 many medical devices. However, no other contact-
23 lens labeling contains this type of information.
24 We have proposed a listing of corneal inflammatory
25 events that occurred in the trial presented in

1 order of most frequent to less frequent.

2 Additional event types and/or rates could
3 be added to this section as proposed by the panel
4 reviewers.

5 [Slide.]

6 In the wearing-schedule section, we have
7 added a chart that identifies the average achieved
8 wearing time for those who completed the one-year
9 FDA trial. As you have heard, virtually all
10 subjects in the Focus Night and Day group were
11 prescribed 30 nights extended wear throughout the
12 study duration but because of symptoms or simply
13 lifestyle requirements, mid-month removals did
14 occur and it was our goal to present information in
15 the labeling that reflected the wearing experience
16 of the subjects in the trial.

17 But, as mentioned earlier, there may be
18 alternative ways of presenting information
19 regarding wearing time that the panel recommends
20 for this section.

21 [Slide.]

22 Also in the wearing-schedule section, we
23 emphasized the importance of close monitoring
24 during the first month of 30-night extended wear.
25 We have added to the standard labeling that

1 patients should be monitored closely during the
2 first month of 30-night continuous wear. If
3 problems occur during this first month, the patient
4 may not be suitable for the full 30-night wearing
5 schedule.

6 This addition is made because many of the
7 problems were noted in the first month of the FDA
8 trial as described earlier. This statement is
9 another alternative to address the reviewers'
10 request that a statement be added that not all
11 patients can tolerate continuous wear.

12 [Slide.]

13 In the lens-fit assessment section of the
14 practitioner fitting guide, we had added a
15 statement about lens-edge standoff and a separate
16 statement about reduced comfort as often being the
17 only signal of a loose-fitting lens. These
18 statements are included in the section that also
19 describes the characteristics of a well-fitting
20 lens or tight-fitting lens and communicates that
21 the lens should demonstrate a satisfactory push-up
22 test and have 0.1 to 0.5 millimeters of movement
23 with the blink.

24 These statements should improve early
25 fitting performance without encouraging fitting

1 practices that result in excessively tight-fitting
2 lenses.

3 [Slide.]

4 So, in summary, product labeling is an
5 essential tool used to distribute key safety and
6 effectiveness information to practitioners and
7 patients. The indication of up to 30 nights
8 extended wear is an important change from the
9 current six-night extended-wear products and,
10 therefore, warrants modifications for current
11 extended-wear labeling.

12 As you have heard, we have taken key
13 information from our FDA study and have used that
14 data to modify the parts of the labeling that you
15 see here.

16 [Slide.]

17 At this point, I would like to talk
18 briefly about the postapproval evaluation that we
19 have submitted.

20 [Slide.]

21 We believe that the preapproval clinical
22 trials have given reasonable assurance that Focus
23 Night and Day is safe and effective as indicated
24 for up to 30 nights extended wear. The high oxygen
25 permeability and biocompatible nature of these

1 lenses are the primary reason for this clinical
2 success. However, certain important low-incidence
3 events such as microbial keratitis require a large
4 trial to determine the event rate.

5 For example, our study of over 650
6 subjects for a year had no cases of microbial
7 keratitis allowing us to conclude that the rate of
8 MK is no greater than approximately 45 in 10,000.
9 But a postapproval evaluation will allow us to
10 increase our confidence that the actual rate is
11 much lower than this.

12 So, to address this important issue, we
13 are working with the agency to design a
14 postapproval evaluation. The questions we have
15 chosen to ask during this evaluation are, number
16 one, is the annualized microbial keratitis rate
17 greater than 20 per 10,000 in Focus Night and Day
18 wearers and, number two, is there vision loss in
19 any case of microbial keratitis that is equal to
20 two or more lines of Snellen acuity.

21 [Slide.]

22 With respect to study design and
23 rationale, both case-control and prospective study
24 designs were considered as alternative approaches
25 for this postapproval evaluation. The case-control

1 study design was ruled out based on our current
2 international marketing experience which predicts
3 that it would be difficult to get sufficient
4 numbers of microbial keratitis cases in a timely
5 fashion.

6 The small number of cases of which we have
7 been informed suggest that a case-control study may
8 require a long case-collection phase in order to
9 get sufficient case numbers.

10 During the February '98 Ophthalmic Panel
11 Meeting, Dr. Schein recommended that these studies
12 should be observational with simple, important
13 outcomes recorded. He stated there should be no
14 doubt that something significant happened.

15 Also at the November, 2000 Ophthalmic
16 Panel Meeting, Dr. Bullimore commented that, in
17 addition to simply observing the rates of important
18 events such as MK, that the loss of best visual
19 acuity should be noted as with refractive-surgery
20 studies.

21 Finally, it was also mentioned at that
22 November, 2000 panel meeting that a dedicated
23 effort should be made to standardize the definition
24 of microbial keratitis for the purpose of
25 consistently counting endpoint events.

1 So these were our parameters in designing
2 our postmarket evaluation.

3 [Slide.]

4 This is a summary of our proposed design.
5 You can see, we have selected a single-group
6 observational study design consisting of 2000 Focus
7 Night and Day wearers all prescribed 30 nights
8 extended wear.

9 Between 100 and 200 clinical sites will
10 participate and the observational period will be
11 for one year. Our endpoints will be microbial
12 keratitis and the loss of best visual acuity of two
13 lines or greater after resolution of MK or other
14 inflammation.

15 [Slide.]

16 An important element of this evaluation
17 would be the use of an independent board selected
18 from ophthalmologists, optometrists,
19 epidemiologists, et cetera, to define the endpoints
20 prior to the evaluation and to review cases during
21 the observation period.

22 Information will be collected from the
23 wearer and the practitioner at three points;
24 baseline, 6 and 12 months. Simple accountability
25 data and a questionnaire will be completed if there

1 has been no event. With an infiltrative event
2 occurring during the observational period, detailed
3 information about the course and outcome of the
4 information will be recorded.

5 The independent board will review these
6 reports and determine if the event was a microbial-
7 keratitis endpoint. So we feel this study design
8 would provide important information about rare
9 inflammatory events in a timely fashion giving us
10 an early warning signal of an unexpectedly high
11 incidence of microbial keratitis and would also
12 allow us to gather more information for future
13 improvements aimed at reducing even minor
14 inflammation.

15 [Slide.]

16 So, at this point, I would like to make a
17 few closing comments. We believe international
18 premarket FDA studies and international market
19 experience have given reasonable assurance that
20 Focus Night and Day is safe and effective for the
21 proposed indications.

22 The proposed labeling adds significant new
23 information about product performance which should
24 enlighten practitioners and patients about the
25 benefits and risks of this product. We have

1 pointed out a few areas in the labeling that could
2 be further enhanced. We feel we have made a pretty
3 good start at the modifications required for this
4 new indication.

5 [Slide.]

6 The proposed postapproval evaluation will
7 add additional information about significant sight-
8 threatening events in a relatively short period
9 after market launch. This evaluation, plus
10 information from the U.S. Medical Device Reporting
11 System and global postmarket vigilance, will allow
12 timely and sufficient monitoring of product
13 performance.

14 This concludes the sponsor's presentation.
15 Thank you for your attention.

16 DR. SUGAR: Thank you. We have about
17 forty minutes until the proposed lunchtime. We
18 have scheduled fifteen minutes for questions for
19 the sponsor. Sometimes it takes longer, sometimes
20 it takes shorter, than that. Then the FDA
21 presentation.

22 I don't see many people squirming so my
23 proposal is that we work through until lunchtime
24 and get as far along in the program as we can,
25 unless someone has strong objections, even weak

1 objections.

2 Then what I would like to do is ask the
3 panel for questions of the sponsor. Alice?

4 DR. MATOBA: I am Alice Matoba. I have a
5 question on the slide from Dr. McNally's
6 presentation on study design, sample size and
7 power. I know you expected an incidence of 8.6
8 percent for endpoint infiltrates in the Acuvue
9 group and you selected an n of 500. So that gave
10 you a power base of 7 percent.

11 But this table shows the power increasing
12 as your expected incidence increases. I think that
13 is incorrect. In fact, you found only an incidence
14 of 3.1 percent in your Acuvue group and that would
15 tend to decrease the power of your study. So my
16 question is, have you reassessed the numbers to see
17 whether, indeed, you did have adequate power to
18 detect a 5 percent difference.

19 DR. McNALLY: Our goal was to determine if
20 we had a difference from 5 percent. Indeed, we did
21 detect a difference from 5 percent

22 DR. MATOBA: It was not significant.

23 DR. McNALLY: It was significant because
24 we were trying to be significantly different than 5
25 percent in a noninferiority design. With a p of

1 0.046, which is -- we set the alpha at 0.05, so
2 0.046 being less than 0.05 says that we are
3 statistically different than 5 percent.

4 DR. MATOBA: Okay. Let me just rephrase
5 it, then. You were looking at a phenomenon that
6 you expected to have an incidence of 8.3 percent.
7 So you selected an n of 500. But if, actually, the
8 phenomenon had an incidence of only 3 percent,
9 would you not expect to need a greater n to detect
10 a significant difference between your product and
11 the product you are comparing it to?

12 DR. McNALLY: This might have to be one
13 that we refer to the statisticians but these
14 calculations here were based on the assumptions to
15 prove a 5 percent difference. The power does
16 increase as the rates go lower for this
17 noninferiority test.

18 DR. MATOBA: But my second question, can
19 you explain --

20 DR. McNALLY: Can we give this maybe the
21 statistician because they can maybe explain it a
22 little better than a clinician.

23 DR. CUTTER: I am Gary Cutter. I am a
24 consultant to Ciba Vision. You are correct about
25 power for a test of no difference between two

1 groups. As the event rate goes down, you would
2 need more sample size if you were doing a null
3 hypothesis of equality between the two groups.
4 This is a noninferiority hypothesis so the
5 difference of 5 percent was fixed and, therefore,
6 as the event rate is smaller, you actually increase
7 power because you have a bigger fixed difference
8 relative to a smaller standard error.

9 You are not trying to get a difference
10 between two groups which now has a smaller standard
11 error. You have a fixed difference with a smaller
12 standard error. Your power actually goes up. This
13 study was designed with the 8.6 event rate from
14 prevalence data and we then were being conservative
15 if, in fact, the rate would be lower in an
16 incidence study where you follow the patients.

17 I think those calculations are correct.
18 It does have to do with the uniqueness, maybe, of
19 the noninferiority study. But it was what I think
20 we were attempting to do.

21 DR. MATOBA: Thank you.

22 DR. SUGAR: Dr. Bandeen-Roche, do you want
23 to comment?

24 DR. BANDEEN-ROCHE: Yes. This is Karen
25 Bandeen-Roche. I just wanted to comment that I

1 agree with the explanation just given. I also want
2 to follow up with a related question which was
3 that, certainly, the Acuvue rate was much smaller
4 than what was projected. You just said you used a
5 prevalence rate, but do you have any other
6 comments? Were you surprised by the Acuvue rate
7 being 3 percent rather than the projected 8
8 percent? How do you account for that?

9 DR. McNALLY: I can make a clinical
10 comment on it. When we started the study, or when
11 we designed the study which I guess would have been
12 '98 or so, we looked for studies in the literature
13 to say what is the expected rate with extended
14 wear. There weren't that many studies. We could
15 only really find one published study by Levy which
16 said a 12 percent rate. Then we were looking for
17 staining over the top which we explained in the
18 report, which gave us the 8.6.

19 We found no other studies other than our
20 own studies which we didn't have extensive year-
21 long studies so there was very little data at that
22 time. I think there has been a lot more data
23 generated in the literature since then.

24 But two panels or three panels ago when we
25 discussed extended wear, one thing that did show

1 out in the study is we thought that the cumulative
2 rate of serious and significant and nonsignificant
3 types of adverse events was about the 10 percent
4 level. I don't know if you recall this. We came
5 out with these 8.something and 9.something rates.

6 So the overall rates, I think, were pretty
7 much in line with what was more in the literature
8 at the time, I think, which was about a 10 percent
9 rate of adverse events. The infiltrate rate,
10 especially as we defined it, we really didn't have
11 a lot go on. We took our best estimate and tried
12 to design a study that would give us that
13 flexibility from 12 percent to 2 percent to still
14 have enough power to perform the noninferiority
15 test.

16 DR. SUGAR: Dr. Pulido?

17 DR. PULIDO: Jose Pulido. I guess you
18 chose the 5 percent because you didn't want to have
19 a 50 percent increase over what was already out
20 there. If the one was 10 percent, you didn't want
21 greater than 15 percent; right?

22 DR. McNALLY: That's right. If we took
23 8.6 and you said a 5 percent, that is not even a
24 two times increase.

25 DR. PULIDO: Correct.