

JUN 30 2000

K000860

Medisys PLC
Futura Safety Syringe 510(k) Notification

510(k) Summary

Futura Safety Syringe

Medisys PLC

Prepared March 15, 2000

Product Name: Futura Safety Syringe

Manufacturer: A Commercial Manufacturing site has not been selected.

Generic Name: Piston syringe with needle

Classification Name: Piston Syringe with needle

Contact Person: Sheila W. Pickering Ph.D.
2081 Longden Circle
Los Altos, California 94024
Telephone/Fax 650 969 6114

A. Legally Marketed Predicate Device

The Futura device is substantially equivalent to the following predicate devices with regard to device features, specifications, and intended use.

| Sponsor | Predicate Device |
|--------------------------------|---------------------|
| New Medical Technology, Ltd. | NMT Safety Syringe |
| Retractable Technologies, Inc. | VanishPoint Syringe |

B. Device Description

The Futura Safety Syringe is a 3cc piston type hypodermic syringe with an automated needle retraction system. It is a sterile, non-toxic, non-pyrogenic, retractable syringe designed to provide a safe and reliable method for intramuscular and subcutaneous injection of drugs and/or fluids while helping to provide protection from accidental needlestick.

Futura has an integral hypodermic needle available in 5/8", 1", and 1-1/2" lengths with needle gauges of 21G, 23G, 25G, 26G, 27G and 28G. Indications for use, dimensions, materials and operation are substantially equivalent to the predicate devices, i.e. the NMT Safety Syringe™ and Retractable Technologies' VanishPoint™ syringes.

000036

C. Intended Use

The device is intended to be used for subcutaneous and intramuscular uses.

The Futura Safety Syringe is designed to provide a safe and reliable method of injecting medication and fluids into patients and also helps to protect the user from potential needlesticks. The Futura Safety Syringe functions as a conventional hypodermic syringe except for its ability to retract the contaminated needle safely inside the syringe immediately after the completion of the patient injection. Complete delivery of the syringe contents activates the retracting mechanism.

D. Substantial Equivalence

The following tables show the basis for substantial equivalence

| Product Name (K number) | Predicate Devices | | | Submission Device | SE |
|----------------------------|--|---|--|--|-----|
| | New Medical Technology NMT Safety Syringe K982431 | Retractable Technologies VanishPoint Syringe K946219 | | Medisys Futura Safety Syringe | |
| Intended Use | <p>The device is intended to be used for subcutaneous and intramuscular uses.</p> <p>The function of the NMT Safety Syringe is to provide a safe and reliable method of injecting medication into a patient that also protects the user from potential needlesticks.</p> <p>The Futura Safety Syringe functions as a conventional hypodermic syringe except for its ability to retract the contaminated needle inside the syringe immediately after the completion of the patient injection. Complete delivery of the syringe contents activates the retracting mechanism.</p> | <p>To retract and contain the contaminated needle after injection and render the syringe inoperable for further injections.</p> | | <p>The device is intended to be used for subcutaneous and intramuscular uses.</p> <p>The function of the Futura Safety Syringe is to provide a safe and reliable method of injecting medication into a patient that also helps to protect the user from potential needlesticks.</p> <p>The Futura Safety Syringe functions as a conventional hypodermic syringe except for its ability to retract the contaminated needle safely inside the syringe immediately after the completion of the patient injection. Complete delivery of the syringe contents activates the retracting mechanism.</p> | Yes |
| Users/Site of Use | Hospitals, clinics, laboratories | Hospitals, clinics, laboratories | | Hospitals, clinics, laboratories | Yes |

Medisys PLC
Futura Safety Syringe 510(k) Notification

| Product Name (K number) | Predicate Devices | | | SE |
|---------------------------------|---|--|--|-----|
| | New Medical Technology NMT Safety Syringe K982431 | Retractable Technologies POP-N-LOCK Syringe K946219 | Submission Device Medisys Futura Safety Syringe | |
| Syringe type | Not provided in labeling | Plunger, antistick with hypodermic needle | Plunger, antistick with hypodermic needle | Yes |
| Volume | 3cc | 3cc -40M | 3cc | Yes |
| Length of barrel and hub | Not provided in labeling | 3.47 inches | 4.19 inches | Yes |
| Diameter | Not provided in labeling | 0.40 inches | 0.42 inches | Yes |
| Materials | Not provided in labeling | Polypropylene rubber, stainless steel, epoxy, lubrican | Polypropylene rubber, stainless steel, epoxy, lubricant | Yes |
| Available needle gauge sizes | 20 and 25G 1 and 1½ in | 25G, ½ in 22G, 1½ in 23G, 1 in. 22G, 1in. 21G, 1½ in | 28, 5/8 in 27, ½ in 26, 1 and 1½ in 25, 1 and 1½ in 23, 1 and 1½ in 21, 1 and 1½ in | Yes |
| Needle cover color | Not provided in labeling | Colorless | Colorless | Yes |

000039

E. Performance Data

Design verification testing was conducted to confirm that the device met all functional specifications and applicable standards. Also, a simulated use evaluation with healthcare professionals demonstrated that the Futura device was substantially equivalent in performance to the predicate devices.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JUN 3 0 2000

Medisys PLC
C/O Sheila W. Pickering, Ph.D.
Regulatory Affairs Consultant
for Medisys PLC
2081 Longden Circle
Los Altos, California 94024

Re: K000860
Trade Name: Futura Safety Syringe
Regulatory Class: II
Product Code: MEG
Dated: May 23, 2000
Received: May 30, 2000

Dear Dr. Pickering:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Good Manufacturing Practice for Medical Devices: General (GMP) regulation (21 CFR Part 820) and that, through periodic GMP inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2 - Dr. Pickering

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4692. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski
Director
Division of Dental, Infection Control,
and General Hospital Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

FDA Submission Cover Sheet

510(k) Number (if known): Not applicable *K 000 860*

Device Name: Futura Safety Syringe

Indications For Use:

The device is intended to be used for subcutaneous and intramuscular use.

The function of the Futura Safety Syringe is designed to provide a safe and reliable method of injecting medication and fluids into patients and also helps to protect the user from potential needlesticks.

The Futura Safety Syringe functions as a conventional hypodermic syringe except for its ability to retract the contaminated needle safely inside the syringe immediately after the completion of the patient injection. Complete delivery of the syringe contents activates the retracting mechanism.

B. Golden

(Division Sign-Off)

Division of Dental, Infection Control,
and General Hospital Devices

510(k) Number *K 000 860*

(PLEASE DO NOT WRITE BELOW THIS LINE - CONTINUE ON ANOTHER PAGE
IF NEEDED)

Concurrence Of CDRH, Office Of Device Evaluation (ODE)

Prescription Use _____

OR

Over-The-Counter Use

(Per 21CFR 801)