

JUL 27 2000

## SUMMARY OF SAFETY AND EFFECTIVENESS

**Sponsor:** Biomet Inc.  
Airport Industrial Park  
P.O. Box 587  
Warsaw, IN 46581-0587

**Contact Person:** Tracy J. Bickel  
(219) 372-1761

**Device(s):** New Bio-Moore Endo Head, Taper Adapter **k002106**

**Classification:** Hip joint femoral (hemi-hip) metallic cemented or uncemented prosthesis  
(888.3360)

**Device Description:** The New Bio-Moore Endo Heads are Co-Cr-Mo shells without the factory assembled, titanium inserts. They are designed to mate with Biomet taper type I trunnions through the use of large Ti-6Al-4V cylindrical taper adapters (the object of this 510k). The surgeon assembles the chosen taper adapter and Endo Head based on trial components at the time of surgery. The taper adapter is placed on the trunnion of the hip femoral and impacted followed by the impaction of the head onto the taper adapter. As with the original design, the New Bio-Moore Endo Heads are intended to articulate with a patient's natural acetabulum.

### Intended Use:

The New Bio-Moore Endo Heads are intended for use in:

- 1) non-inflammatory degenerative joint disease, including osteoarthritis and avascular necrosis
- 2) rheumatoid arthritis
- 3) correction of functional deformities
- 4) revision procedures where other devices and treatments have failed
- 5) treatment on non-union, femoral neck and trochanteric fractures of the proximal femur with head involvement

The Bio-Moore Endo Heads are intended for impaction either as cemented or press-fit femoral component and are for single use implantation.

**Potential Risks:** The potential risks associated with this device are the same as with any joint replacement device. These include, but are not limited to:

Reaction to bone cement	Blood vessel damage	Bone fracture
Deformity of the joint	Soft tissue imbalance	Infection
Cardiovascular disease	Delayed wound healing	Hematoma
Fracture of the cement	Metal sensitivity	Dislocation
Implant loosening/migration	Fracture of the components	Excessive wear
Tissue growth failure	Nerve damage	

000054



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

JUL 27 2000

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Ms. Tracy J. Bickel  
Regulatory Specialist  
Biomet, Inc.  
P.O. Box 587  
Warsaw, Indiana 46581-0587

Re: K002106  
Trade Name: Bio-Moore Endo Head, Taper Adapter  
Regulatory Class: II  
Product Code: JDI and JDG  
Dated: July 5, 2000  
Received: July 12, 2000

Dear Ms. Bickel:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general control provisions of the Act. The general control provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Ms. Tracy J. Bickel

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

*Donna R. Witten*

*SR* Celia M. Witten, Ph.D., M.D.  
Director  
Division of General, Restorative and  
Neurological Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

510(k) Number (if known): K002106  
Device Name: New Bio-Moore Endo Head, Taper Adapter  
Indications for Use:

The New Bio-Moore Endo Heads are intended for use in:

- 1) non-inflammatory degenerative joint disease, including osteoarthritis and avascular necrosis
- 2) rheumatoid arthritis
- 3) correction of functional deformities
- 4) revision procedures where other devices and treatments have failed
- 5) treatment of non-union, femoral neck and trochanteric fractures of the proximal femur with head involvement

The Bio-Moore Endo Heads are intended for impaction on either a cemented or press-fit femoral hip component and are for single use implantation.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use X  
(Per 21 CFR 801.109)

OR

Over-The-Counter Use \_\_\_\_\_  
(Optional Format 1-2-96)

Danner R. Kechner  
(Division Sign-Off)  
Division of General Restorative Devices  
510(k) Number K002106