

**ADMINISTRATIVE INFORMATION**

Manufacturer Name:	MacroPore, Inc. 6740 Top Gun Street San Diego, CA 92121
Official Contact:	Kenneth K. Kleinhenz Director of Regulatory Affairs Telephone (858) 458-0900 Fax (858) 458-0994

**DEVICE NAME**

Classification Name:	Plate, Cranioplasty, Non-Alterable
Trade/Proprietary Name:	MacroPoreNS CraniLoc

**ESTABLISHMENT REGISTRATION NUMBER**  
2031733

**DEVICE CLASSIFICATION AND PRODUCT CODE**

As shown in 21CFR 882.5330 Non-Alterable Cranioplasty Plates are intended for use in craniotomy procedures and are classified as Class II. Non-Alterable Cranioplasty Plates have been assigned Product Code GXN.

**INTENDED USE**

MacroPoreNS CraniLoc is indicated for the fixation of bone flaps after a craniotomy.

**DEVICE DESCRIPTION****Design Characteristics**

MacroPoreNS CraniLoc is a resorbable, macroporous implant manufactured from polylactic acid (PLA). The MacroPoreNS CraniLoc is composed of a threaded male and tapped female disk that screw into each other. The MacroPoreNS CraniLoc utilizes two mating disks that incrementally engage each other via a threaded shaft, resulting in a compression force that stabilizes adjacent bone fragments without the need for external fixators (screws, tacks, etc.). The MacroPoreNS CraniLoc is designed to fit into an drilled hole or osteotomy cut and exert vertical pressure to stabilize flap movement in both the vertical and horizontal directions. Typically, 3 or more MacroPoreNS CraniLocs will be used to secure a bone flap. Once engaged to the desired compression force, the shaft of the MacroPoreNS CraniLoc can be cut to make the shaft flush with the anatomy. Various manual instruments (drivers, excision devices, etc.) are used in conjunction with the MacroPoreNS CraniLoc to assist in the installation process.

The MacroPoreNS CraniLoc comes in two basic designs, one for use in burr holes and another style for use in the osteotomy cut. The design features between the burr hole version and the osteotomy line version are virtually identical except that the burr hole models have a round stem while the osteotomy line models have a semi-round stem with two flattened sides. The MacroPoreNS CraniLoc is provided in various sizes to fit burr holes ranging from 8mm to 19mm in diameter. The MacroPoreNS CraniLoc is also designed to fit into the osteotomy line where 1-2mm of clearance is available. The pore size of the upper and lower disks ranges from 500 microns to 3000 microns in diameter, with pores distributed uniformly throughout. The thickness of the MacroPoreNS CraniLoc's disks ranges from 0.50 mm to 2.0 mm according to the anatomical region to be treated. The outer diameter of the MacroPoreNS CraniLoc shaft ranges from 1mm to 10mm. The shaft of the MacroPoreNS CraniLoc's male component may be provided as either solid or cannulated, depending on strength requirements. The larger shaft diameters may also be porous with holes to allow for in growth of bone and tissues.

### Material Composition

The MacroPoreNS CraniLoc is fabricated from polylactic acid (PLA).

### In Vitro Testing

Mechanical testing of the MacroPoreNS CraniLoc demonstrates that the device is substantially equivalent to the predicate. Test results indicate that the mechanical properties of the MacroPoreNS CraniLoc are substantially equivalent to the mechanical properties of the predicate device under indication for use conditions.

### EQUIVALENCE TO MARKETING PRODUCT

The MacroPoreNS CraniLoc shares indications and design principles with the following predicate devices which have been determined by FDA to be substantially equivalent to pre-amendment devices: Aesculap CranioFIX Titanium Clamp System and Walter Lorenz RapidFlap Cranial Clamp, Class II medical device that were cleared for marketing in the United States under K972332 and K991029 respectively.

The MacroPoreNS CraniLoc is composed of the same raw material and processed in the same manner as previously FDA cleared medical devices from MacroPore

- K983360 MacroPore Protective Sheet
- K972913 MacroPore (Protego System) plates, screws, and Protective Sheet
- K000992 MacroPoreDX Distractor
- K994158 MacroPoreOS Protective Sheet

**Indications For Use**

The MacroPoreNS CraniLoc and the predicate device share substantially equivalent indications for use as they both are indicated for the fixation of bone flaps after a craniotomy. Specifically, the MacroPoreNS CraniLoc is indicated for the fixation of bone flaps after a craniotomy.

**Design and Materials**

The MacroPoreNS and the predicate device share substantially equivalent design features as both devices are composed of two disks that are attached through a threaded stem. Both devices secure bone in an identical manner as both the MacroPoreNS CraniLoc and the predicate device press the bone cover and the vault of the cranium between two disks. Both the MacroPoreNS CraniLoc and the predicate device utilize two mating disks that incrementally engage each other via a threaded shaft, resulting in a compression force that stabilizes adjacent bone fragments without the need for external fixtures (screws, tacks, etc.). The MacroPoreNS CraniLoc and the predicate also share design characteristics of fitting into a drilled burr hole or fitting into the osteotomy cut. Both the MacroPoreNS CraniLoc and the predicate share principles of operation as they both exert vertical pressure to stabilize flap movement in both the vertical and horizontal directions. Both the MacroPoreNS CraniLoc and the predicate utilize various manual instruments (drivers, excision devices, etc.) to assist in the installation process.

## SUMMARY : TABLE OF SUBSTANTIAL EQUIVALENCE

The MacroPoreNS CraniLoc is substantially equivalent to the Aesculap CranioFIX and the Walter Lorenz RapidFlap devices in the following respects:

	Subject Device			Predicate Device			Related Device		
	MacroPoreNS MacroPoreNS CraniLoc	Aesculap CranioFIX (K972332)	Walter Lorenz RapidFlap (K991029)	MacroPore Protective Sheet (Protego System) (K972913)	Codman Ethisorb Dura Patch (K991413)	Neuroregen Neurotube (K983007)			
<b>Intended Use</b>	MacroPoreNS CraniLoc is indicated for the fixation of bone flaps after a craniotomy.	Aesculap's CranioFIX Titanium Clamp System is intended for use in refixation of cranial bone flaps after craniotomy.	RapidFlap is indicated for the re-attachment of the bone flap after a craniotomy.	MacroPore Protective Sheet is intended for use in trauma and reconstructive procedures in the midface and craniofacial skeleton: 1. Comminuted fractures of the naso-ethmoidal and infraorbital areas 2. Comminuted fractures of the frontal sinus wall 3. Trauma of the midface or craniofacial skeleton 4. Reconstructive procedures of the midface or craniofacial skeleton. The system is not intended for use in the mandible and/or for full load bearing procedures.	The Codman Ethisorb Dura Patch is an absorbable, synthetic implant for bridging defects of the dura mater.	The Neurotube is intended for single use in patients with an injury to a peripheral nerve, in which the nerve gap is more than or equal to 8mm but less than or equal to 3cm. The nerve gap may be created primarily at the time of injury or created secondarily at the time of exploration of failed primary repair.			
<b>Design</b>	Two disks that are attached through a threaded stem	Two disks that are attached through a threaded stem	Two disks that are attached through a threaded stem	Plates, screws, protective sheet mesh, and tacks of various shapes and sizes.	Sheet	2.3mm x 4cm long Tube			
<b>Material</b>	Poly (L-lactide-co-D,L-lactide) 70:30, amorphous	Titanium	Titanium	Poly (L-lactide-co-D,L-lactide) 70:30, amorphous	Polyglactin 910 (90% polyglycolide, 10% L-lactide) and polydioxanone.	Polyglycolic acid (PGA)			
<b>Product Code</b>	GXN	GXN	GXN	HRS and HWC	GXQ	JXI			



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAY - 2 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

Mr. Kenneth K. Kleinhenz  
Director of Regulatory Affairs  
MacroPore, Inc.  
6740 Top Gun Street  
San Diego, California 92121

Re: K002334  
Trade/Device Name: MacroPoreNS CraniLoc  
Regulation Number: 882.5250 and 882.5360  
Regulatory Class: II  
Product Code: GXR and HBW  
Dated: January 31, 2001  
Received: February 1, 2001

Dear Mr. Kleinhenz:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the current Good Manufacturing Practice requirement, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic (QS) inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

Page 2 - Mr. Kenneth K. Kleinhenz

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or at (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsmamain.html>".

Sincerely yours,

A handwritten signature in black ink, appearing to read "Neil R. P. Ogden", followed by a small flourish.

Celia M. Witten, Ph.D., M.D. *for*

Director

Division of General, Restorative and

Neurological Devices

Office of Device Evaluation

Center for Devices and

Radiological Health

Enclosure

Device Name: MacroPoreNS CraniLoc

K002334

**Indications for Use:**

MacroPoreNS CraniLoc is indicated for the fixation of bone flaps after a craniotomy.

NRO for crani  
(Division Sign-Off)  
Division of General, Restorative  
and Neurological Devices

510(k) Number K002334

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NECESSARY)

---

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use

~~X~~

OR

Over-The-Counter Use

\_\_\_\_\_