

AUG 29 2001

K003535

## **Solo-SAFE™ Releasable 510(k) Summary**

### **Submitter of the 510(k)**

Spectrum Biotech Inc.  
3990 Sunnycrest Drive  
North Vancouver, British Columbia  
V7R 3C9  
Tel: 604.987.4884  
Fax: 604.985.5103

### **Regulatory Consultant:**

Ian McDougall, B.Sc.  
821 East 17<sup>th</sup> Street  
North Vancouver, British Columbia  
V7L 2X2  
Tel: 604.831.9311  
Fax: 604.987.0268  
E-mail: imcdougall@medscape.com

Pursuant to Section 510(k) of the Federal Food, Drug and Cosmetic Act, Safety and Effectiveness information is enclosed for the following device:

1. **Trade/Proprietary Name:**  
Solo-Safe™ Safety Syringe  
Size – 3cc
2. **Predicate Device:** BD Safety-Lok K920321 and K924072
3. **Common/Usual Name:**  
Safety Syringe, Hypodermic Syringe, Piston Syringe with Safety Mechanism
4. **Classification Name:**  
Piston Syringe with integral needle
5. **Classification:**  
Class II per 21 CFR 880.5860
6. **Product Code:**  
MEG– Anti-stick Syringe
7. **Classification Panel:**  
Panel 80 (General Hospital Use)

8. **Predicate Device Information:**

A claim of substantial equivalence is made to: BD Safe-Lok K920321 and K924072

**Device Description**

The Solo-SAFE™ Syringe is a piston syringe with integral needle combined with a sharps injury prevention mechanism. The Solo-SAFE™ Syringe is a sterile, non-toxic, non-pyrogenic, retractable syringe designed to provide a safe and reliable method for intramuscular injection of drugs and/or fluids while helping to provide a protection from accidental needlestick. The injury prevention mechanism allows the user to actively pull the contaminated needle back into the barrel of the syringe after use. The user accomplishes this by pulling back on the plunger which has connected to the Solo-SAFE™ needle after being fully depressed. By pulling back on the plunger the user is able to retract the Solo-SAFE™ needle within the barrel of the syringe.

**Substantial Equivalence**

The following tables show the basis for substantial equivalence:

	<b>Predicate Device (K920321)</b>	<b>Submission Device</b>	<b>SE</b>
	<b>BD Safety Lok 3cc Syringe</b>	<b>Solo- SAFE 3cc Syringe</b>	
<b>Indications for Use</b>	Indicated for IM injection	Indicated for IM injection	YES
<b>Safety Feature</b>	Active safety feature, manually activated by user	Active safety feature, manually activated by user	YES
<b>Sites of Use/Users</b>	Hospitals, clinics, laboratories	Hospitals, clinics, laboratories	YES
<b>Syringe Type</b>	Plunger, antistick with hypodermic needle	Plunger, antistick with hypodermic needle	YES
<b>Volume</b>	3cc	3cc	YES
<b>Available Needle Ga</b>	25G x 5/8 inches 23G x 1 inches 22G x 1 ½ inches 22G x 1 inch 21G x 1 ½ inches	21G x 1 ½ inches	YES
<b>Materials</b>	Polypropylene, rubber, stainless steel, adhesive, lubricant	Polypropylene, rubber, stainless steel, lubricant, adhesive	YES

**Statement of Intended Use**

The Solo-SAFE Syringe is indicated for use in the administration of an intramuscular (IM) injection.

The Spectrum 3cc Solo-SAFE Syringe primary intended use is for injection of medications and fluids into the patient. The Solo-SAFE has an active safety feature that is safe and effective for intramuscular (IM) injections only. Its secondary intended use is the safety feature of the device which helps prevent sharps injuries when using the device for its primary intended use. Solo-SAFE is a sterile, single use, disposable syringe.

**Performance Data**

Design verification testing was conducted to confirm that the device met functional and specifications. Also a simulated use study with healthcare professionals demonstrated that the Solo-SAFE syringe was substantially equivalent in performance to its predicate device.



AUG 29 2001

Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

C/O Mr. Ian McDougall  
Regulatory Consultant  
Spectrum Biotech, Incorporated  
821 East 17<sup>th</sup> Street  
North Vancouver, British Columbia  
V7L 2X2

Re: K003555  
Trade/Device Name: Solo-Safe Safety Syringe, Size-3CC  
Regulation Number: 880.5860  
Regulatory Class: II  
Product Code: MEG  
Dated: August 11, 2001  
Received: August 16, 2001

Dear Mr. McDougall:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements

concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4613. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,



Timothy A. Ulatowski  
Director  
Division of Dental, Infection Control  
and General Hospital Devices  
Office of Device Evaluation  
Center for Devices and  
Radiological Health

Enclosure

K003555

**Indications for Use**

*The Solo-SAFE Syringe is indicated for use in the administration of an intramuscular (IM) injection.*

*The Spectrum 3cc Solo-SAFE Syringe primary intended use is for injection of medications and fluids into the patient. The Solo-SAFE has an active safety feature that is safe and effective for intramuscular (IM) injections only. Its secondary intended use is the safety feature of the device which helps prevent sharps injuries when using the device for its primary intended use. Solo-SAFE is a sterile, single use, disposable syringe.*

*Patricia Curran*

(Division Sign-Off)

Division of Dental, Infection Control,  
and General Hospital Devices

510(k) Number K003555