

MAR 29 2001

K004051

II. 510(k) Summary

Submitted by: Eagle TWE, Inc.
2090 Tucker Industrial Road, Suite A-5
Tucker, GA 30084

Contact Person: Keith F. Heintz
Phone (864) 375-1817
Fax (864) 375-0505
E-mail Origrethnk@aol.com

Date Prepared: December 21, 2000

Proprietary Name: Sheer Thin Aromas Yellow Latex Condom

Common Name: Condom

Classification Name: Condom (Rubber) Contraceptive

Predicate Device: Condom with Flavor
Line One Laboratories
510(k) Document Control Number K970633

Description of the Device: This condom is made of a natural rubber latex sheath, which completely covers the penis with a closely fitting membrane. This condom is a straight walled, nipple end, yellow colored condom that uses a dusting powder and is lubricated with silicone oil and a banana flavor.

Intended Use of the Device: This latex condom has the same intended use as the predicate condom. The condom is used for contraception and for prophylactic purposes to help prevent pregnancy and the transmission of sexually transmitted diseases.

Technological Characteristics: The latex condom has the same technological characteristics as the predicate condom identified above. The design of the condom is in accordance with the ASTM Latex Condom Standard D3492 and the condom is made out of natural latex rubber.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

MAR 29 2001

Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

Mr. Keith F. Heintz
Vice President
Eagle TWE, Inc.
2090 Tucker Industrial Road
TUCKER GA 30084

Re: K004051
Sheer Thin Aromas Yellow Male Latex Condom
Dated: December 21, 2000
Received: December 29, 2000
Regulatory Class: II
21 CFR §884.5300/Procode: 85 HIS

Dear Mr. Heintz:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.


Please be advised that, as of March 25, 1998, labeling for latex condoms (21 CFR §884.5300 and §884.5310) must comply with Use Labeling for Latex Condoms: Expiration Dating, 21 CFR §801.435. Therefore, an expiration date, supported by test data developed under the conditions specified in §801.435(d), must be displayed prominently and legibly on condom labeling. For condoms with spermicidal lubricant, the effective shelf life of the spermicide must be compared with the shelf life of the condom and labeled with the earlier of the two expiration dates. Although supporting data is not to be provided in your 510(k) submission, §801.435(j) requires that you maintain this data and that it be available for inspection by FDA.

Page 2 - Mr. Keith Heintz

Furthermore, §801.435(e) requires that if your real-time test data fails to confirm the shelf life estimated by the methods in §801.435(d), then you must relabel all products to reflect the actual shelf life. Condoms are not to be labeled with an expiration date that gives a shelf life more than five years.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at (301) 594-4616. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR §807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597, or at its Internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in dark ink, appearing to read "Daniel G. Schultz", with a stylized flourish at the end.

Daniel G. Schultz, M.D.
Acting Director, Division of Reproductive,
Abdominal, and Radiological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Indications for Use Statement

510(k) Number:

K004051

Device name:

Sheer Thin Aromas Yellow Latex Condom

Indications for Use:

The Sheer Thin Aromas condom is used for contraception and for prophylactic purposes to help prevent pregnancy and the transmission of sexually transmitted diseases.

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER PAGE IF NEEDED)

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use _____
(Per 21 CFR 801.109)

OR Over-The-Counter-Use ✓ _____
(Optional Format 1-2-96)

David G. Segerson
(Division Sign-Off)
Division of Reproductive, Abdominal, ENT,
and Radiological Devices
510(k) Number K004051