

JAN 25 2002

K011074

Section II - Summary of Safety and Effectiveness

(1) Contact Information

Vincent Cutarelli
Sr. Vice President, Regulatory Affairs and Quality Assurance
Telephone: (949) 595-5424
FAX: (949) 595-4766
email: vcutarelli@endocare.com

(2) Company Information

Endocare, Inc.
7 Studebaker
Irvine, CA 92618
Telephone: (949) 595-4770
FAX: (949) 595-4766

(3) Device Name

Cryocare® Surgical System

(4) Device Description

The Cryocare™ Surgical System consists of a control unit that operates one to eight single-use, disposable CryoProbes. The system utilizes inert argon gas as the cooling agent and is currently available in 1, 4 and 8-CryoProbe configurations. The control unit is software-controlled and operates off standard 110/230 VAC wall power. A 486 IBM-compatible microprocessor serves as the host computer and a screen displays the status of the system. System control is accomplished either directly through keys on the console itself (e.g., 1-probe system) or through a remote control keypad (e.g., 4 and 8-probe system). The CryoProbes operate on the Joule-Thompson principle and the freezing capacity is limited only to the distal tip of the probe. The CryoProbes incorporate a thermocouple to measure temperatures at the probe tip. The thermocouple is mounted inside each CryoProbe tip and its signal is used to monitor and control some operations of the system. The control unit can also control one to eight independent TempProbes™ to monitor temperatures in surrounding tissues. The temperature probes are standard T-type needle thermocouples.

RECEIVED
DEC 19 8 07 AM '01
FROM: JRNH/00E/ENDOCARE

(5) **Indications for Use**

The Cryocare® Surgical System is intended for use in open, minimally invasive or endoscopic surgical procedures in the areas in general surgery, urology, gynecology, oncology, neurology, dermatology, ENT, proctology, pulmonary surgery and thoracic surgery. The system is designed to freeze/ablate tissue by the application of extreme cold temperatures including prostate and kidney tissue, liver metastases, tumors, skin lesions, and warts. In addition, the system is intended for use in the following indications:

General Surgery

- Destruction of warts or lesions
- Palliation of tumors of the oral cavity, rectum and skin
- Ablation of leukoplakia of the mouth, angiomas, sebaceous hyperplasia, basal cell tumors of the eyelid or canthus area, ulcerated basal cell tumors, dermatofibromas, small hemangiomas, mucocele cysts, multiple warts, plantar warts, hemorrhoids, anal fissures, perianal condylomata, pilonidal cysts, actinic and seborrheic keratoses, cavernous hemangiomas, recurrent cancerous lesions

Urology

- Ablation of prostate tissue in cases of prostate cancer and benign prostatic hyperplasia

Gynecology

- Ablation of malignant neoplasia or benign dysplasia of the female genitalia

Oncology

- Ablation of cancerous or malignant tissue
- Ablation of benign tumors
- Palliative intervention

Neurology

- Freezing of nerve tissue in pain management/cryoanalgesia

Dermatology

- Ablation or freezing of skin cancers and other cutaneous disorders

Proctology

- Ablation of benign or malignant growths of the anus or rectum
- Ablation of hemorrhoids

Thoracic Surgery

- Ablation of arrhythmic cardiac tissue
- Ablation of cancerous lesions

(6) Predicate or Legally Marketed Device(s)

Cryocare® Surgical System
CryoGen Cryosurgical System 3
Wallach Painblocker WA5000 Cryosurgical System

(7) Substantial Equivalence

The Cryocare® Surgical System is substantially equivalent to the Cryocare® Surgical System that was determined to be substantially equivalent on April 10, 1998 (reference K980110), the CryoGen Cryosurgical System 3 that was determined to be substantially equivalent on December 26, 2000 (reference K003050) and the Wallach Painblocker WA5000 Cryosurgical System that was determined to be substantially equivalent on April 30, 1986 (reference K854334).

(8) Technological Characteristics

The technological characteristics of the Cryocare® Surgical System are the same as those of the predicate devices. These devices are substantially equivalent in terms of design, materials, principle of operation, product specifications and sterilization.



Food and Drug Administration
9200 Corporate Boulevard
Rockville MD 20850

JAN 25 2002

Endocare, Inc.
Vincent Cutarelli
Sr. Vice-President, Regulatory Affairs and Quality Assurance
7 Studebaker
Irvine, California 92618

Re: K011074
Trade Name: Cryocare® Surgical System
Regulation Number: 878.4350
Regulation Name: Cryosurgical Unit and Accessories
Regulatory Class: II
Product Code: GEH
Dated: November 26, 2001
Received: December 18, 2001

Dear Mr. Cutarelli:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.


Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (Sections 531-542 of the Act); 21 CFR 1000-1050.

Page 2 – Mr. Vincent Cutarelli

This letter will allow you to begin marketing your device as described in your Section 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 21 CFR Part 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4659. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>

Sincerely yours,

for Mark N. Melnick

Celia M. Witten, Ph.D., M.D.
Director
Division of General, Restorative
and Neurological Devices
Office of Device Evaluation
Center for Devices and
Radiological Health

Enclosure

Indications For Use

510(k) Number: K011074

Device Name: Cryocare™ Surgical System

Indications for Use: The Cryocare® Surgical System is intended for use in open, minimally invasive or endoscopic surgical procedures in the areas in general surgery, urology, gynecology, oncology, neurology, dermatology, ENT, proctology, pulmonary surgery and thoracic surgery. The system is designed to freeze/ablate tissue by the application of extreme cold temperatures including prostate and kidney tissue, liver metastases, tumors, skin lesions, and warts. In addition, the system is intended for use in the following indications:


General Surgery

- Destruction of warts or lesions
- Palliation of tumors of the oral cavity, rectum and skin
- Ablation of leukoplakia of the mouth, angiomas, sebaceous hyperplasia, basal cell tumors of the eyelid or canthus area, ulcerated basal cell tumors, dermatofibromas, small hemangiomas, mucocele cysts, multiple warts, plantar warts, hemorrhoids, anal fissures, perianal condylomata, pilonidal cysts, actinic and seborrheic keratoses, cavernous hemangiomas, recurrent cancerous lesions

Urology

- Ablation of prostate tissue in cases of prostate cancer and benign prostatic hyperplasia

Concurrence of CDRH, Office of Device Evaluation (ODE):


(Division Sign-off)
Division of General Restorative
and Neurological Devices
510(k) Number K011074

Prescription Use: X
(Per 21 CFR 801.109)

K011074

Gynecology

- Ablation of malignant neoplasia or benign dysplasia of the female genitalia

Oncology

- Ablation of cancerous or malignant tissue
- Ablation of benign tumors
- Palliative intervention

Neurology

- Freezing of nerve tissue in pain management/cryoanalgesia

Dermatology

- Ablation or freezing of skin cancers and other cutaneous disorders

Proctology

- Ablation of benign or malignant growths of the anus or rectum
- Ablation of hemorrhoids

Thoracic Surgery

- Ablation of arrhythmic cardiac tissue
- Ablation of cancerous lesions

Concurrence of CDRH, Office of Device Evaluation (ODE):

for Mark N. Melkerson
(Division Sign-Off)
Division of General, Restorative
and Neurological Devices
510C Number K011074

Prescription Use: X
(Per 21 CFR 801.109)