

AUG 15 2002



755 8<sup>th</sup> Court, Suite #4  
PO Box 650790  
Vero Beach, FL 32965-0790  
Tel: 561-569-5955  
Fax: 561-569-4430  
e-mail: [sales@zanderivf.com](mailto:sales@zanderivf.com)

**510(k) Summary****20 April 2001**

**Submitted by:** Zander Medical Supplies, Inc.  
755 8<sup>th</sup> Court, Suite 4  
P.O. Box 650790  
Vero Beach, FL 32965-0790  
(561) 569-5955  
(561) 569-4430 Fax  
Contact: Friedel MW Zander, President/CEO  
E-mail: [fzander@zanderivf.com](mailto:fzander@zanderivf.com)

**Manufactured by:** Biotronics  
Unit 1, Morton Farm  
Morton Eye, Leominster  
HR6 ODP U.K.  
+44 1568 612-402  
+44 1568 616-088 Fax  
Contact: Steve Goode, Managing Director  
E-mail: [sales@biotronics.net](mailto:sales@biotronics.net)

**Prepared by:** Friedel MW Zander  
Zander Medical Supplies, Inc.

**Trade Name:** Biotronics DB1 Embryo Freeze  
**Common Name:** DB1

**Class:** II                      **CFR:** 884.6120                      **Procode:** 85 MQG

**Predicate Device:** There is currently no Predicate Device for a freezer of this sort.

**Device Description:**

The DB1 EmbryoFreeze has been developed to meet the requirements of specialists working in both human and veterinary embryo transfer. This portable freezer can be used in most locations with its own internal battery and is easy to program and use, with a throughput of up to 100 straws per cycle. The DB1 is a "solid state" design without moving components. There are no solenoid valves, pumps or motors. The straw plate is precisely controlled and the temperature is balanced against the cold liquid nitrogen (LN<sub>2</sub>) chamber.

## **Intended Use:**

The DB1 EmbryoFreeze is used for freezing biological materials such as embryos, oocytes, etc., in cryostraws for cryopreservation using Liquid Nitrogen (LN2).

## **Technical Details:**

### **Power Supply**

Internal Battery(s):	12V-DC 7Ah. Sealed Lead Acid. 1 or 2 installed.
External DC supply:	12V-DC 10A fuse. Maximum 14V-DC.
Mains Power:	85V-AC to 264V-AC 47 to 440 Hz or 120V-DC to 370V-DC 2.5A.
Run time on internal battery:	Approximately 3 hours on one battery / 6 hours on two batteries.

### **Fuses**

Mains:	5A (5A spare)
Heater:	5A
12V-DC Input:	10A (internal)
Internal Battery:	10A (internal)
Internal Power Supply:	4A (internal)

### **LN2 Supply**

LN2 Capacity:	Approximately 2.5 litres.
Static Holding Time:	Up to 2 hours.
LN2 Consumption:	3 hour use from 2.5 litres (after initial boil off).

### **Case**

Size:	L=405mm x W=270mm x H=295mm (including feet).
Weight:	9.5 Kgs with one battery, additional battery add 2.65 Kgs.
Material:	Polyurethane high density foam.

### **Controls**

Display:	240 x 120 pixels, 30 x 16 characters.
Back Light:	Fluorescent EL tube.
Buttons:	Push button with green LED.
Alpha Wheel:	Optical rotary encoder.
On/Off Button:	Push on and push off.

### **Connections**

Mains:	Filtered and fused (plus spare fuse).
12V-DC:	Internally fused, centre positive.
Printer:	Standard 25 way parallel printer interface (HP PCL Compatible).
Serial:	Standard 9 pin RS232 interface. Use null modem lead.
Service:	12 way of which 3 are COM, N/O & N/C alarm relay outputs.

### **General Details:**

Cooling & Heating Rates:	Minimum rate: 0.1°C Maximum rate: 5.0°C
Temperature Range:	Minimum temperature: -80°C Maximum temperature: 150°C
Timed Hold Range:	Minimum timed hold: 1 minute Maximum timed hold: 99 minutes Setting resolution: 1-minute steps. Timing resolution: 0.01 seconds.
Control Resolution:	0.01°C. Data logging and display 0.1°C
Aim Temperature Resolution:	1°C
Maximum Program Lines:	8
Number of Stored Programs:	8
Password Characters:	4
Data Log Frequency:	Log recorded every 10 seconds when running a program.

### **Label Examples:**

Importer & Distributor:



Manufacturer:

A sample of the Manufacturer's labeling can be found in the User Manual, page 4. It is screen-printed onto the case, with the exception of the CE Label, the serial number on the back and mains warning label on the base. These are self-adhesive standard labels.

### **General Details:**

Cooling & Heating Rates:	Minimum rate: 0.1°C Maximum rate: 5.0°C
Temperature Range:	Minimum temperature: -80°C Maximum temperature: 150°C
Timed Hold Range:	Minimum timed hold: 1 minute Maximum timed hold: 99 minutes Setting resolution: 1-minute steps. Timing resolution: 0.01 seconds.
Control Resolution:	0.01°C. Data logging and display 0.1°C
Aim Temperature Resolution:	1°C
Maximum Program Lines:	8
Number of Stored Programs:	8
Password Characters:	4
Data Log Frequency:	Log recorded every 10 seconds when running a program.

### **Label Examples:**

Importer & Distributor:



Manufacturer:

A sample of the Manufacturer's labeling can be found in the User Manual, page 4. It is screen-printed onto the case, with the exception of the CE Label, the serial number on the back and mains warning label on the base. These are self-adhesive standard labels.



Food and Drug Administration  
9200 Corporate Boulevard  
Rockville MD 20850

AUG 15 2002

Mr. Friedel M. Zander  
President /CEO  
Zander Medical Supplies, Inc.  
755 8<sup>th</sup> Court, Suite 4  
P.O. Box 650790  
VERO BEACH FL 32965-0790

Re: K011327  
Trade/Device Name: DB1 Embryo Freeze  
Regulation Number: 21 CFR 884.6120  
Regulation Name: Assisted reproduction accessories  
Regulatory Class: II  
Product Code: 85 MQG  
Dated: July 30, 2002  
Received: Aug 1, 2002

Dear Mr. Zander:

We have reviewed your Section 510(k) premarket notification of intent to market the device referenced above and have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act) that do not require approval of a premarket approval application (PMA). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.

If your device is classified (see above) into either class II (Special Controls) or class III (PMA), it may be subject to additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 898. In addition, FDA may publish further announcements concerning your device in the Federal Register.

Please be advised that FDA's issuance of a substantial equivalence determination does not mean that FDA has made a determination that your device complies with other requirements of the Act or any Federal statutes and regulations administered by other Federal agencies. You must comply with all the Act's requirements, including, but not limited to: registration and listing (21 CFR Part 807); labeling (21 CFR Part 801); good manufacturing practice requirements as set forth in the quality systems (QS) regulation (21 CFR Part 820); and if applicable, the electronic product radiation control provisions (sections 531-542 of the Act); 21 CFR 1000-1050.


This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801), please contact the Office of Compliance at one of the following numbers, based on the regulation number at the top of this letter:

8xx.1xxx	(301) 594-4591
876.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4616
884.2xxx, 3xxx, 4xxx, 5xxx, 6xxx	(301) 594-4616
892.2xxx, 3xxx, 4xxx, 5xxx	(301) 594-4654
Other	(301) 594-4692

Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21 CFR Part 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers, International and Consumer Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its Internet address <http://www.fda.gov/cdrh/dsma/dsmamain.html>.

Sincerely yours,



Nancy C. Brogdon  
Director, Division of Reproductive,  
Abdominal, and Radiological Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



755 8<sup>th</sup> Court, Suite #4  
PO Box 650790  
Vero Beach, FL 32965-0790  
Tel: 561-569-5955  
Fax: 561-569-4430  
WATS-Line: 1-800-820-3029  
e-mail: [sales@zanderIVF.com](mailto:sales@zanderIVF.com)

510(k) NUMBER (IF KNOWN) : K011327

DEVICE NAME : DB1 EmbryoFreeze

INDICATIONS FOR USE :

The DB1 EmbryoFreeze is used for freezing biological materials such as embryos, oocytes, etc., in cryostraws for cryopreservation using Liquid Nitrogen (LN2).

(Division Sign-Off)

Division of Reproductive, Abdominal,  
and Radiological Devices

510(k) Number K011327

(PLEASE DO NOT WRITE BELOW THIS LINE – CONTINUE ON ANOTHER  
PAGE IF NEEDED).

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-the-Counter-Use \_\_\_\_\_  
(Optional Format 1-2-96)