

JUL 17 2001

K011499



## SUN BIOMEDICAL LABORATORIES, INC.

604 VPR CENTER, 1001 LOWER LANDING ROAD, BLACKWOOD, NJ 08012  
Tel. 856-401-1080 Fax. 856-401-1090

### 510(k) CONTENT SUMMARY

1. Name of Manufacturer:

Sun Biomedical Laboratories, Inc.  
604 VPR Center, 1001 Lower Landing Rd.  
Blackwood, NJ 08012

2. Trade Name: Visualine® Amphetamine DipStrip Assay

3. Common Name:

An in-vitro immunoassay test by visual color comparison for the detection of Amphetamine in human urine samples. **This test is to be used for professional use only.**

4. Regulation # and Classification:

Reg. #862-~~3170~~<sup>3160</sup>, Class II Device

5. Test Description:

The Visualine® Amphetamine DipStrip test is based on the principle of antigen-antibody complexation and is used for the analysis of Amphetamine in urine samples. The assay utilizes a competitive immunochromatographic technique involving a sample of test urine delivered in a sample well on the device that holds the porous membrane. When the drug is present in the urine test sample, the drug or drug metabolite competes for the limited antibody sites on the colored microspheres. When an adequate amount of drug is present, it will fill the limited antibody binding sites. This will prevent attachment of the colored microspheres to the probe site on the membrane. Therefore, a positive urine sample will inhibit the formation of precipitin at the probe site.

A reference or control line with a secondary antibody reaction is added to the membrane strip to indicate that the sample is properly wicking on the membrane. This control line should always be present. A negative urine sample will produce two colored lines and a positive urine sample will show only one, the control line.



DEPARTMENT OF HEALTH & HUMAN SERVICES

Public Health Service

Food and Drug Administration  
2098 Gaither Road  
Rockville MD 20850

JUL 17 2001

Ming Sun Ph.D.  
President  
Sun Biomedical Laboratories, Inc.  
604 VPR Center  
1001 Lower Landing Road  
Blackwood, NJ 08012

Re: 510(K) Number: K011499  
Trade/Device Name: Visualine® Amphetamine Dip Strip Test Kit  
Regulation Number: 862.3170  
Regulatory Class: II  
Product Code: DKZ  
Dated: May 11, 2001  
Received: May 15, 2001

Dear Dr. Sun:

We have reviewed your Section 510(k) notification of intent to market the device referenced above and we have determined the device is substantially equivalent (for the indications for use stated in the enclosure) to legally marketed predicate devices marketed in interstate commerce prior to May 28, 1976, the enactment date of the Medical Device Amendments, or to devices that have been reclassified in accordance with the provisions of the Federal Food, Drug, and Cosmetic Act (Act). You may, therefore, market the device, subject to the general controls provisions of the Act. The general controls provisions of the Act include requirements for annual registration, listing of devices, good manufacturing practice, labeling, and prohibitions against misbranding and adulteration.


If your device is classified (see above) into either class II (Special Controls) or class III (Premarket Approval), it may be subject to such additional controls. Existing major regulations affecting your device can be found in the Code of Federal Regulations, Title 21, Parts 800 to 895. A substantially equivalent determination assumes compliance with the Current Good Manufacturing Practice requirements, as set forth in the Quality System Regulation (QS) for Medical Devices: General regulation (21 CFR Part 820) and that, through periodic QS inspections, the Food and Drug Administration (FDA) will verify such assumptions. Failure to comply with the GMP regulation may result in regulatory action. In addition, FDA may publish further announcements concerning your device in the Federal Register. Please note: this response to your premarket notification submission does not affect any obligation you might have under sections 531 through 542 of the Act for devices under the Electronic Product Radiation Control provisions, or other Federal laws or regulations.

Page 2

This letter will allow you to begin marketing your device as described in your 510(k) premarket notification. The FDA finding of substantial equivalence of your device to a legally marketed predicate device results in a classification for your device and thus, permits your device to proceed to the market.

If you desire specific advice for your device on our labeling regulation (21 CFR Part 801 and additionally 809.10 for in vitro diagnostic devices), please contact the Office of Compliance at (301) 594-4588. Additionally, for questions on the promotion and advertising of your device, please contact the Office of Compliance at (301) 594-4639. Also, please note the regulation entitled, "Misbranding by reference to premarket notification" (21CFR 807.97). Other general information on your responsibilities under the Act may be obtained from the Division of Small Manufacturers Assistance at its toll-free number (800) 638-2041 or (301) 443-6597 or at its internet address "<http://www.fda.gov/cdrh/dsma/dsmamain.html>".

Sincerely yours,

A handwritten signature in dark ink, appearing to read "Steven Gutman". The signature is fluid and cursive, with the first name "Steven" and last name "Gutman" clearly distinguishable.

Steven I. Gutman, M.D., M.B.A.  
Director  
Division of Clinical Laboratory Devices  
Office of Device Evaluation  
Center for Devices and Radiological Health

Enclosure



# SUN BIOMEDICAL LABORATORIES, INC.


604 VPR CENTER, 1001 LOWER LANDING ROAD, BLACKWOOD, NJ 08012  
Tel. 856-401-1080 Fax. 856-401-1090

Page 1 of 1

510 (k) Number: K011499

DEVICE NAME: Visualine® Amphetamine Dip Strip Test

INDICATIONS FOR USE: The Visualine® Amphetamine Dip Strip Assay is used for qualitative testing for the presence of Amphetamines in urine samples at 1000 ng/ml. This test provides only a preliminary screening result; a more specific alternative method should be used to confirm the test result. This test is intended for use by medical professionals.

  
(Division Sign-Off)  
Division of Clinical Laboratory Devices  
510(k) Number K011499

Concurrence of CDRH, Office of Device Evaluation (ODE)

Prescription Use ✓  
(Per 21 CFR 801.109)

OR

Over-The-Counter-Use \_\_\_\_\_  
(Optional Format 1-2-96)