

K954570

510(k) SUMMARY

APR 30 1996

SUBMITTER: Ortho Diagnostic Systems Inc
1001 U.S. Highway 202
P.O. Box 350
Raritan, NJ 08869-0606

CONTACT: Joanne Harris
Tel.: (908) 218-8404
FAX.: (908) 218-8168

DEVICE NAME:

ORTHO™ ImmunoCount Flow Cytometry System

ORTHO ImmunoCount Flow Cytometry System
includes:

ORTHO CytoronAbsolute™ Laser Flow Cytometer
(optional AutoBioSampler)

ORTHO ImmunoCount II Software

ORTHO TRIO Monoclonal Antibodies (Murine)

ORTHO TRIO Control
(FITC/PE/CyP)
Monoclonal Antibodies (Murine)

ORTHO TRIO CD4/CD8/CD3
(OKT™4A/OKT8/OKT3)
(FITC/PE/CyP)
Monoclonal Antibodies (Murine)

ORTHO TRIO CD16/CD19/CD3
(NK/OKB™19A/OKT3)
(FITC/PE/CyP)
Monoclonal Antibodies (Murine)

Ortho-Count Calibration Kit
(currently under review, Ref. No. K-935720)

PREDICATE:

FACScan™ Flow Cytometer
Simulset Software
Simultest™ reagents K93292

Simultest Control IgG₁/IgG_{2a}
(IgG₁ FITC/IgG_{2a}PE)

Simultest LeucoGATE™
(CD45/CD14[Anti-HLe 1/Leu™M3])

Simultest CD3/CD4 (Leu-4/3a)

Simultest CD3/CD8 (Leu-4/Leu-2a)

Simultest CD3/CD16+CD56
(Leu-4/11c+19)

Simultest T and B Cell Test
CD3/CD19 (Leu4/Leu12)

Hematology Analyzers

Coulter T 540 K896873

Coulter JT2 K874383

AGROS K895127

DATE: September 29, 1995

DEVICE DESCRIPTION:

ORTHO ImmunoCount Flow Cytometry System consists of ORTHO CytoronAbsolute Laser Flow Cytometer, ImmunoCount II Software, ORTHO TRIO Monoclonal Antibodies for lymphocyte immunophenotyping, and Ortho-Count Calibration Kit for calibration and verification of system calibration.


INTENDED USE:

ORTHO ImmunoCount Flow Cytometry System is intended to be used for lymphocyte immunophenotyping. Results may be reported as either percent positive cells, or as absolute counts of lymphocytes and lymphocyte subsets.

TECHNOLOGICAL CHARACTERISTICS:

ORTHO ImmunoCount Flow Cytometry System is substantially equivalent to the combination of flow cytometer, lymphocyte immunophenotyping reagents, and hematology analyzer when used for determination of percentages and absolute counts of lymphocytes and lymphocyte subsets.

The traditional method for determination of absolute counts of lymphocyte subsets combines lymphocyte immunophenotyping reagents and a flow cytometer for determination of lymphocyte subset percentages, with a hematology analyzer for determination of total lymphocyte count, to obtain absolute counts of lymphocyte subsets.



PERFORMANCE DATA:

ORTHO ImmunoCount Flow Cytometry System was compared to the traditional method for determining absolute counts of lymphocytes and lymphocyte subsets in a clinical study conducted at four laboratories. The instruments and reagents used in the clinical evaluation of the ORTHO ImmunoCount Flow Cytometry System are:

FACScan™ Flow Cytometer with Simulset Software
- Becton Dickinson Immunocytometry Systems

Simulset™ Monoclonal Antibody Immunophenotyping Reagents
-Becton Dickinson Immunocytometry Systems

Simulset Control IgG₁/IgG_{2a} (IgG₁ FITC/IgG_{2a} PE)
Simulset LeucoGATE™ (CD45/CD14 [Anti-HLe-1/Leu™-M3])
Simulset CD3/CD4 (Leu-4/3a)
Simulset CD3/CD8 (Leu-4/Leu-2a)
Simulset CD3/CD16+CD56 (Leu-4/11c+19)
Simulset T and B Cell Test - CD3/CD19 (Leu4/Leu12)

Coulter® JT2 Automated Hematology Analyzer
-Coulter Electronics

Coulter T Series Automated Hematology Analyzer
-Coulter Electronics

Argos Hematology Analyzer
-ABX France

Blood specimens were collected from 200 apparently healthy donors. Normal donors were defined as any person who was not clinically diagnosed with an immunodeficiency disease or hematological malignancy. The donors were male (44.5%) and female (55.5%) with an age range of 21 years to 62 years. Three geographically distinct clinical sites in the United States were used to determine absolute lymphocyte counts and four geographically distinct clinical sites in the United States were used to determine the percent positive stained cells. A total of 200 normal whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and analyzed in the ORTHO ImmunoCount Flow Cytometry System for percent positive stained cell determination. A total of 151 normal whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and analyzed in the ORTHO ImmunoCount Flow Cytometry System for absolute lymphocyte count determination.

TABLE A

IMMUNOCOUNT SYSTEM LYMPHOCYTE SUBSET PERCENT POSITIVE					
NORMAL DONORS N = 200					
ImmunoCount	Mean%	Range%	Predicate	Mean%	Range%
CD3 ⁺ CD4 ⁺	45.5	20 - 68	CD3 ⁺ CD4 ⁺	45.2	20 - 65
CD3 ⁺ (4/8/3)	74.5	36 - 101	CD3 ⁺ (3/4)	73.2	29 - 91
CD3 ⁺ CD8 ⁺	24.5	9 - 48	CD3 ⁺ CD8 ⁺	26.9	9 - 55
CD3 ⁺ (4/8/3)	74.5	36 - 101	CD3 ⁺ (3/8)	73.5	31 - 90
CD16 ⁺ CD3 ⁻	11.6	1 - 58	CD(16&56) ⁺ CD3 ⁻	13.1	2 - 65
CD3 ⁺ (16/19/3)	75.1	39 - 89	CD3 ⁺ (3/16&56)	73.6	31 - 91
CD19 ⁺	13.3	3 - 51	CD19 ⁺	13.3	3 - 28
CD3 ⁺ (16/19/3)	75.1	39 - 89	CD3 ⁺ (3/19)	73.3	30 - 89
Mean CD3 ⁺	74.8	38 - 95	Mean CD3 ⁺	73.4	30 - 90
IMMUNOCOUNT SYSTEM LYMPHOCYTE SUBSET ABSOLUTE COUNTS					
NORMAL DONORS N = 151					
ImmunoCount	Mean (cells/μl)	Range (cells/μl)	Predicate	Mean (cells/μl)	Range (cells/μl)
CD3 ⁺ CD4 ⁺	855	257 - 1785	CD3 ⁺ CD4 ⁺	900	250 - 1940
CD3 ⁺ (4/8/3)	1417	477 - 3372	CD3 ⁺ (3/4)	1472	480 - 3530
CD3 ⁺ CD8 ⁺	481	112 - 1588	CD3 ⁺ CD8 ⁺	551	100 - 1750
CD3 ⁺ (4/8/3)	1417	477 - 3372	CD3 ⁺ (3/8)	1477	470 - 3500
CD16 ⁺ CD3 ⁻	221	19 - 1536	CD(16&56) ⁺ CD3 ⁻	266	60 - 2280
CD3 ⁺ (16/19/3)	1428	512 - 3423	CD3 ⁺ (3/16&56)	1476	470 - 3530
CD19 ⁺	255	35 - 826	CD19 ⁺	270	50 - 840
CD3 ⁺ (16/19/3)	1428	512 - 3423	CD3 ⁺ (3/19)	1469	480 - 3540
Mean CD3 ⁺	1423	494 - 3398	Mean CD3 ⁺	1473	475 - 3525
ISUM	1905	667 - 4241	lymphocyte	2010	700 - 4159

Blood specimens were collected from 119 HIV-antibody positive donors. HIV-antibody positive donors were defined as any individual positive for the HIV antibody inclusive of AIDS or AIDS Related Complex (ARC) according to the current criteria for these diseases, Morbidity Mortality Weekly Report, Vol. 41/No. 31, May 8, 1992). The donors were male (88.1%) and female (11.9%) with an age range of 3 years to 75 years. Three geographically distinct clinical sites in the United States were used to determine absolute lymphocyte counts and four geographically distinct clinical sites in the United States were used to determine the percent positive stained cells. A total of 119 HIV-antibody positive whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and analyzed in the ORTHO ImmunoCount Flow Cytometry System for percent positive stained cell determination. A total of 89 HIV-antibody positive whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and analyzed in the ORTHO ImmunoCount Flow Cytometry System for absolute lymphocyte count determination.

TABLE B

IMMUNOCOUNT SYSTEM LYMPHOCYTE SUBSET PERCENT POSITIVE					
HIV ANTIBODY POSITIVE DONORS N = 119					
ImmunoCount	Mean%	Range%	Predicate	Mean%	Range%
CD3 ⁺ CD4 ⁺	17.7	0 - 55	CD3 ⁺ CD4 ⁺	17.6	0 - 56
CD3 ⁺ (4/8/3)	80.1	18 - 115	CD3 ⁺ (3/4)	78.1	24 - 96
CD3 ⁺ CD8 ⁺	57.3	7 - 89	CD3 ⁺ CD8 ⁺	58.6	16 - 91
CD3 ⁺ (4/8/3)	80.1	18 - 115	CD3 ⁺ (3/8)	78.3	21 - 96
CD16 ⁺ CD3 ⁻	8.8	0 - 45	CD(16&56) ⁺ CD3 ⁻	10.6	1 - 50
CD3 ⁺ (16/19/3)	80.9	23 - 95	CD3 ⁺ (3/16&56)	78.4	14 - 97
CD19 ⁺	10.4	2 - 42	CD19 ⁺	10.8	2 - 40
CD3 ⁺ (16/19/3)	80.9	23 - 95	CD3 ⁺ (3/19)	77.4	22 - 95
Mean CD3 ⁺	80.5	21 - 103	Mean CD3 ⁺	78.1	20 - 96
IMMUNOCOUNT SYSTEM LYMPHOCYTE SUBSET ABSOLUTE COUNTS					
HIV ANTIBODY POSITIVE DONORS N = 89					
ImmunoCount	Mean (cells/μl)	Range (cells/μl)	Predicate	Mean (cells/μl)	Range (cells/μl)
CD3 ⁺ CD4 ⁺	260	0 - 1029	CD3 ⁺ CD4 ⁺	313	0 - 1400
CD3 ⁺ (4/8/3)	1078	42 - 2986	CD3 ⁺ (3/4)	1262	50 - 3650
CD3 ⁺ CD8 ⁺	758	39 - 1961	CD3 ⁺ CD8 ⁺	931	40 - 2810
CD3 ⁺ (4/8/3)	1078	42 - 2986	CD3 ⁺ (3/8)	1267	40 - 3600
CD16 ⁺ CD3 ⁻	93	5 - 572	CD(16&56) ⁺ CD3 ⁻	144	20 - 720
CD3 ⁺ (16/19/3)	1079	34 - 2927	CD3 ⁺ (3/16&56)	1267	30 - 3650
CD19 ⁺	131	13 - 554	CD19 ⁺	170	20 - 870
CD3 ⁺ (16/19/3)	1079	34 - 2927	CD3 ⁺ (3/19)	1248	40 - 3650
Mean CD3 ⁺	1078	38 - 2957	Mean CD3 ⁺	1261	40 - 3638
ISUM	1303	118 - 3691	lymphocyte	1568	200 - 4600

Linear regression analyses of absolute lymphocyte counts and percent positive cells obtained using the ORTHO ImmunoCount Flow Cytometry System and the traditional (predicate) method are presented below.

TABLE C

IMMUNOCOUNT SYSTEM VS. PREDICATE METHOD LINEAR REGRESSION					
NORMAL N = 200 AND HIV ANTIBODY POSITIVE N = 119					
ImmunoCount	Predicate	Graph	R	Slope	Intercept
CD4 ⁺ (4/8/3)	CD4 ⁺ (3/4)	3a	1.00	1.00	0
CD3 ⁺ (4/8/3)	CD3 ⁺ (3/4)	3b	0.83	0.87	12
CD8 ⁺ (4/8/3)	CD8 ⁺ (3/8)	3c	0.98	1.01	-2
CD3 ⁺ (4/8/3)	CD3 ⁺ (3/8)	3d	0.85	0.90	9
CD16 ⁺ CD3 ⁺ (16/19/3)	CD16 ⁺ CD3 ⁺ (3/16+65)	3e	0.91	0.90	0
CD3 ⁺ (16/19/3)	CD3 ⁺ (3/16+56)	3f	0.93	0.88	11
CD19 ⁺ (16/19/3)	CD19 ⁺ (3/19)	3g	0.88	0.98	0
CD3 ⁺ (16/19/3)	CD3 ⁺ (3/19)	3h	0.91	0.91	9
Mean CD3 ⁺	Mean CD3 ⁺	3i	0.91	0.91	9

TABLE D

IMMUNOCOUNT SYSTEM VS. PREDICATE METHOD LINEAR REGRESSION					
NORMAL N = 151 AND HIV ANTIBODY POSITIVE N = 89					
ImmunoCount	Predicate	Graph	R	Slope	Intercept
CD4 ⁺ (4/8/3)	CD4 ⁺ (3/4)	3A	0.94	0.89	25
CD3 ⁺ (4/8/3)	CD3 ⁺ (3/4)	3B	0.89	0.84	119
CD8 ⁺ (4/8/3)	CD8 ⁺ (3/8)	3C	0.92	0.78	43
CD3 ⁺ (4/8/3)	CD3 ⁺ (3/8)	3D	0.89	0.84	112
CD16 ⁺ CD3 ⁺ (16/19/3)	CD16 ⁺ CD3 ⁺ (3/16+65)	3E	0.94	0.76	7
CD3 ⁺ (16/19/3)	CD3 ⁺ (3/16+56)	3F	0.88	0.84	123
CD19 ⁺ (16/19/3)	CD19 ⁺ (3/19)	3G	0.89	0.86	10
CD3 ⁺ (16/19/3)	CD3 ⁺ (3/19)	3H	0.89	0.85	124
Mean CD3 ⁺	Mean CD3 ⁺	3I	0.89	0.84	118
ISUM	Lymphocyte	3J	0.88	0.86	87
All lymphocyte subsets	All lymphocyte subsets	3K	0.96	0.89	3

WITHIN LABORATORY REPRODUCIBILITY

Five normal whole blood specimens were used in a within-laboratory reproducibility study at three independent laboratories. The normal whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and run in the ORTHO ImmunoCount Flow Cytometry System. The samples were run diluted (0.5X WBC), concentrated (4X WBC) and undiluted to simulate low, high and normal samples. The leukocyte-depleted blood was prepared using a Sepracell™ unit (Baxter Fenwall). Buffy coat preparations were prepared from the donor samples, counted and adjusted to a 4X concentration. Serial dilutions were made from the 4X concentrate of each sample's buffy coat and diluting this concentrate with leukocyte depleted whole blood from the same sample. All samples were run in replicates of 10. The coefficients of variations (CV) for normal, concentrated and diluted samples were compared and are contained in Table E. The coefficient of variation was generated from variability between the ten replicates from each of the three laboratories. The low CV values demonstrate strong within-laboratory reproducibility for each lymphocyte subset. The CV values for normal and low counts for CD16+ and CD19+ values were higher due to the low counts for these populations. The higher CVs for CD16+ and CD19+ values do not significantly affect the ImmunoSum (ISUM) CV.

TABLE E			
ORTHO IMMUNOCOUNT SYSTEM			
WITHIN-LABORATORY REPRODUCIBILITY			
COEFFICIENT OF VARIATION OF ABSOLUTE LYMPHOCYTE COUNTS			
N(LOW) = 139; N(NORMAL) = 142; N(HIGH) = 146			
	LOW	NORMAL	HIGH
ORTHO ImmunoCount System	CV	CV	CV
CD3 ⁺ (4/8/3)	5.2	4.9	4.3
CD3 ⁺ CD4 ⁺	6.2	5.2	4.5
CD3 ⁺ CD8 ⁺	6.4	6.4	4.5
CD16 ⁺ CD3 ⁺	16.3	10.1	6.0
CD3 ⁺ (16/19/3)	5.0	4.9	4.3
CD19 ⁺	12.0	10.0	5.7
ISUM Lymph Count	5.0	4.9	4.2

BETWEEN LABORATORY REPRODUCIBILITY

Five normal whole blood specimens were compared between three independent laboratories. The whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and analyzed in the ORTHO ImmunoCount Flow Cytometry System. Low, normal and high count preparations were made according to the method described in the Within-Laboratory Reproducibility section. All samples were run in replicates of ten. The data collected from all laboratories are shown in Table F. The coefficient of variation was generated by comparing the three laboratories absolute lymphocyte subset means from the 10 replicates. The low CV values demonstrate acceptable between-laboratory reproducibility at all three laboratories. The CV values for normal and low counts for CD16+ and CD19+ values were higher due to the low counts for these populations. The higher CVs for CD16+ and CD19+ values do not significantly affect the ImmunoSum (ISUM) CV.

TABLE F ORTHO IMMUNOCOUNT SYSTEM BETWEEN-LABORATORY REPRODUCIBILITY COEFFICIENT OF VARIATION OF ABSOLUTE LYMPHOCYTE COUNTS				
LOW COUNT PREPARATION : N = 139				
ORTHO ImmunoCount System	Lab 1 Mean Cells/ μ L	Lab 2 Mean Cells/ μ L	Lab 3 Mean Cells/ μ L	Between-site CV
CD3 ⁺ (4/8/3)	750	740	736	0.9
CD3 ⁺ CD4 ⁺	481	480	473	0.9
CD3 ⁺ CD8 ⁺	231	221	224	2.2
CD16 ⁺ CD3 ⁺	146	122	160	13.1
CD3 ⁺ (16/19/3)	754	748	740	0.9
CD19 ⁺	116	98	121	10.3
ISUM Lymph count	1011	963	1009	2.6
NORMAL COUNT PREPARATION: N = 142				
ORTHO ImmunoCount System	Lab 1 Mean Cells/ μ L	Lab 2 Mean Cells/ μ L	Lab 3 Mean Cells/ μ L	Between-site CV
CD3 ⁺ (4/8/3)	1494	1468	1445	1.6
CD3 ⁺ CD4 ⁺	958	953	933	1.3
CD3 ⁺ CD8 ⁺	456	434	437	2.6
CD16 ⁺ CD3 ⁺	363	290	309	11.5
CD3 ⁺ (16/19/3)	1502	1515	1489	0.8
CD19 ⁺	221	212	250	8.5
ISUM Lymph count	2067	2002	2028	1.6
HIGH COUNT PREPARATION: N = 146				
ORTHO ImmunoCount System	Lab 1 Mean Cells/ μ L	Lab 2 Mean Cells/ μ L	Lab 3 Mean Cells/ μ L	Between-site CV
CD3 ⁺ (4/8/3)	5961	5718	5663	2.7
CD3 ⁺ CD4 ⁺	3786	3687	3667	1.7
CD3 ⁺ CD8 ⁺	1835	1693	1683	4.8
CD16 ⁺ CD3 ⁺	1548	1334	1347	8.4
CD3 ⁺ (16/19/3)	5990	5929	5737	2.2
CD19 ⁺	962	844	999	8.3
ISUM Lymph count	8417	8028	7959	3.0

SPECIMEN AGE CRITERIA

A total of 64 whole blood specimens, normal (34 donors) and HIV-antibody positive (30 donors), were tested at three external clinical sites over an extended time period. The samples were collected in EDTA anticoagulant tubes. Samples were prepared and tested within 24 hours, 48 hours and 72 hours after collection. The absolute counts for each ORTHO TRIO Monoclonal Antibody, run in the ORTHO ImmunoCount Flow Cytometry System, are compared in below. The data demonstrates that samples processed within 24 to 72 hours of collection produce comparable results. CI refers to confidence interval of the mean slope.

TABLE G

ORTHOMUNOCOUNT SYSTEM ABSOLUTE LYMPHOCYTE COUNTS OVER A 72-HOUR TIME PERIOD NORMAL DONORS N=34; HIV-ANTIBODY POSITIVE DONORS N=30					
ORTHOMunocount System	Mean Count 24 Hours (Cells/ μ L)	Mean Count 48 Hours (Cells/ μ L)	Mean Count 72 Hours (Cells/ μ L)	Mean Slope	CI
CD3 ⁺ (4/8/3)	1367	1388	1362	-1.5	2.3
CD3 ⁺ CD4 ⁺	615	631	619	-0.4	1.0
CD3 ⁺ CD8 ⁺	664	665	653	-1.2	1.3
CD16 ⁺ CD3 ⁺	188	199	249	0.8	1.4
CD3 ⁺ (16/19/3)	1360	1356	1338	-3.6	2.4
CD19 ⁺	207	206	192	-0.0	0.7
ISUM Lymph Count	1755	1761	1779	-2.9	2.6

LINEARITY

Four normal whole blood specimens were used in a linearity study at Ortho Diagnostic Systems. The whole blood specimens were stained with ORTHO TRIO Monoclonal Antibodies and run in the ORTHO ImmunoCount Flow Cytometry System. The samples were run diluted (0.5X, 0.25X, 0.125X, 0.065X), concentrated (5X, 4X, 2X) and undiluted (1X) to yield low, high, and normal cell counts. The leukocyte depleted blood was prepared using a SepracellTM unit (Baxter Fenwall.) Buffy coat preparations were prepared from the donor sample, counted and diluted to a 5X concentration. Serial dilutions were made from the 5X concentrate of each sample's buffy coat and diluting this concentrate with leukocyte depleted whole blood from the same unit. Table H shows the correlation coefficient of the linear regression analysis of absolute counts for each marker and ISUM.

TABLE H

ORTHO IMMUNOCOUNT SYSTEM LINEAR RANGE OF ABSOLUTE COUNTS NORMAL DONORS N=4		
ORTHO ImmunoCount System	R	Range Cells/μL
CD3 ⁺ (4/8/3)	1.00	75 - 9829
CD3 ⁺ CD4 ⁺	1.00	0 - 6300
CD3 ⁺ CD8 ⁺	1.00	0 - 2952
CD16 ⁺ CD3 ⁻	0.99	16 - 1796
CD3 ⁺ (16/19/3)	1.00	59 - 10390
CD19 ⁺	1.00	19 - 2507
ISUM Lymph Count	1.00	95 - 14694

CONCLUSION

- The performance of the ImmunoCount Flow Cytometry System was equivalent to the predicate method for the determination of absolute lymphocyte counts and lymphocyte subset percentages and absolute counts CD3+, CD3+CD4+, CD3+CD8+, CD16+CD3-, and CD19+.
- ORTHO ImmunoCount Flow Cytometry System demonstrated excellent within and between laboratory reproducibility for low, normal, and high absolute counts for total lymphocytes and lymphocyte subsets.
- ORTHO ImmunoCount Flow Cytometry System, when used to assay samples at extended time points after collection (48 through 96 hours after collection), showed that samples assayed up to 72 hours after collection produced equivalent absolute counts for total lymphocytes and for each of the lymphocyte subsets CD3+, CD3+CD4+, CD3+CD8+, CD16+CD3-, and CD19+.
- ORTHO ImmunoCount Flow Cytometry System demonstrated excellent linearity when compared against dilution over a wide range of absolute counts for each lymphocyte subset and ISUM.