

K961968

FEB - 4 1997

May 17, 1996

SUMMARY OF SAFETY AND EFFECTIVENESS

SUBMITTED BY:

Virginia C. Weinknecht
Becton Dickinson Microbiology Systems
P O Box 243
Cockeysville, MD 21030-0243

NAME OF DEVICE:

Trade Name:	BBLCRYSTAL™ Gram Positive ID System
Common Name/Description:	Miniaturized Microorganism ID System
Classification Name:	Microbiology - Discs, Strips, and Reagents, Microorganism Differentiation

PREDICATE DEVICES:

BioMerieux Vitek, Inc., API Staph (K813614)
BioMerieux Vitek, Inc., API 20 Strep (K813610)
BioMerieux Vitek, Inc., API Coryne (K910304)

DEVICE DESCRIPTION:

INTENDED USE: The BBLCRYSTAL™ Gram Positive (GP) Identification (ID) System is a miniaturized identification method employing modified conventional, fluorogenic, and chromogenic substrates. It is intended for the identification of frequently isolated aerobic gram-positive bacteria from clinical specimens.

INDICATIONS FOR USE: Use of the BBLCRYSTAL™ Gram Positive Identification System is indicated when the aerobic gram-positive organisms described in the attached table have been isolated in pure culture from clinical specimens in a clinical laboratory, and identification of the microorganisms is desired.

TABLE B-I: Taxa List for BBLCRYSTAL™ Gram Positive Identification System

<i>Actinomyces pyogenes</i>	<i>Lactococcus garvieae</i>	<i>Staphylococcus kloosii</i>
<i>Aerococcus species</i> (includes <i>A. urinae</i> and <i>A. viridans</i>)	<i>Lactococcus lactis</i> ssp <i>cremoris</i>	<i>Staphylococcus lentus</i>
<i>Aerococcus urinae</i>	<i>Lactococcus lactis</i> ssp <i>hordniae</i>	<i>Staphylococcus lugdunensis</i>
<i>Aerococcus viridans</i>	<i>Lactococcus lactis</i> ssp <i>lactis</i>	<i>Staphylococcus pasteurii</i> ¹
<i>Alloisococcus otitidis</i> ¹	<i>Lactococcus raffinolactis</i>	<i>Staphylococcus saccharolyticus</i>
<i>Arcanobacterium haemolyticum</i> ¹	<i>Lactococcus species</i> (includes <i>L. lactis</i> ssp <i>cremoris</i> , <i>L. lactis</i> ssp <i>hordniae</i> , <i>L. lactis</i> ssp <i>lactis</i> and <i>L. raffinolactis</i>)	<i>Staphylococcus saprophyticus</i>
<i>Bacillus brevis</i>	<i>Leuconostoc citreum</i>	<i>Staphylococcus schleiferi</i> (includes <i>S. schleiferi</i> ssp <i>caguli</i> and <i>S. schleiferi</i> ssp <i>schleiferi</i>)
<i>Bacillus cereus</i>	<i>Leuconostoc lactis</i>	<i>Staphylococcus sciuri</i>
<i>Bacillus circulans</i>	<i>Leuconostoc mesenteroides</i> ssp <i>mesenteroides</i>	<i>Staphylococcus simulans</i>
<i>Bacillus coagulans</i>	<i>Leuconostoc pseudomesenteroides</i>	<i>Staphylococcus vitulus</i>
<i>Bacillus licheniformis</i>	<i>Leuconostoc species</i> (includes <i>L. citreum</i> , <i>L. lactis</i> , <i>L. mesenteroides</i> ssp <i>mesenteroides</i> and <i>L. pseudomesenteroides</i>)	<i>Staphylococcus warneri</i>
<i>Bacillus megaterium</i>	<i>Listeria grayi</i> ¹	<i>Staphylococcus xylosum</i>
<i>Bacillus pumilus</i>	<i>Listeria ivanovii</i> ssp <i>ivanovii</i>	<i>Stomatococcus mucilaginosus</i>
<i>Bacillus species</i> (includes <i>B. brevis</i> , <i>B. circulans</i> , <i>B. coagulans</i> , <i>B. licheniformis</i> , <i>B. megaterium</i> , <i>B. pumilus</i> and <i>B. sphaericus</i> , <i>P. alvei</i> , <i>P. macerans</i>)	<i>Listeria monocytogenes</i>	<i>Streptococcus acidominimus</i>
<i>Bacillus sphaericus</i>	<i>Listeria murrayi</i>	<i>Streptococcus agalactiae</i>
<i>Bacillus subtilis</i>	<i>Micrococcus kristinae</i>	<i>Streptococcus anginosus</i>
<i>Corynebacterium aquaticum</i>	<i>Micrococcus luteus</i>	<i>Streptococcus bovis</i> (includes <i>S. bovis</i> I and <i>S. bovis</i> II)
<i>Corynebacterium bovis</i>	<i>Micrococcus lylae</i>	<i>Streptococcus constellatus</i>
<i>Corynebacterium diphtheriae</i> (includes <i>C. diphtheriae</i> ssp <i>gravis</i> , <i>C. diphtheriae</i> ssp <i>mitis</i> and <i>C. diphtheriae</i> ssp <i>intermedius</i>)	<i>Micrococcus roseus</i>	<i>Streptococcus cricetus</i> ¹
<i>Corynebacterium genitalium</i>	<i>Micrococcus sedentarius</i>	<i>Streptococcus crista</i>
<i>Corynebacterium jeikeium</i>	<i>Micrococcus species</i> (includes <i>M. kristinae</i> , <i>M. luteus</i> , <i>M. lylae</i> , <i>M. roseus</i> , <i>M. sedentarius</i>)	<i>Streptococcus dysgalactiae</i>
<i>Corynebacterium kutscheri</i>	<i>Oerskovia species</i> (includes <i>O. turbata</i> , and <i>O. xanthineolytica</i>)	<i>Streptococcus equi</i> (includes <i>S. equi</i> ssp <i>equi</i> and <i>S. equi</i> ssp <i>zooepidemicus</i>)
<i>Corynebacterium propinquum</i>	<i>Paenibacillus alvei</i>	<i>Streptococcus equi</i> ssp <i>equi</i>
<i>Corynebacterium pseudodiphtheriticum</i>	<i>Paenibacillus macerans</i>	<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>
<i>Corynebacterium pseudogenitalium</i>	<i>Pedilococcus demnosus</i>	<i>Streptococcus equinus</i>
<i>Corynebacterium pseudotuberculosis</i>	<i>Pedilococcus parvulus</i>	<i>Streptococcus gordonii</i>
<i>Corynebacterium renale</i> group	<i>Pedilococcus pentosaceus</i>	<i>Streptococcus intermedius</i>
<i>Corynebacterium species</i> (includes <i>C. aquaticum</i> , <i>C. bovis</i> , <i>C. kutscheri</i> , <i>C. propinquum</i> , <i>C. pseudodiphtheriticum</i> , <i>C. pseudotuberculosis</i> , <i>C. renale</i> group, <i>C. striatum</i> and <i>C. ulcerans</i>)	<i>Pedilococcus species</i> (includes <i>P. demnosus</i> , <i>P. parvulus</i> , and <i>P. pentosaceus</i>)	<i>Streptococcus milleri</i> group (includes <i>S. anginosus</i> , <i>S. constellatus</i> and <i>S. intermedius</i>)
<i>Corynebacterium striatum</i>	<i>Rhodococcus equi</i>	<i>Streptococcus mitis</i>
<i>Corynebacterium ulcerans</i>	<i>Rothia dentocariosa</i> ¹	<i>Streptococcus mitis</i> group (includes <i>S. mitis</i> and <i>S. oralis</i>)
<i>Enterococcus avium</i>	<i>Staphylococcus aureus</i>	<i>Streptococcus mutans</i>
<i>Enterococcus casseliflavus/gallinarum</i>	<i>Staphylococcus auricularis</i>	<i>Streptococcus mutans</i> group (includes <i>S. cricetus</i> , <i>S. mutans</i> and <i>S. sobrinus</i>)
<i>Enterococcus durans</i>	<i>Staphylococcus capitis</i> (includes <i>S. capitis</i> ssp <i>capitis</i> and <i>S. capitis</i> ssp <i>ureolyticus</i>)	<i>Streptococcus oralis</i>
<i>Enterococcus faecalis</i>	<i>Staphylococcus caprae</i>	<i>Streptococcus parasanguis</i>
<i>Enterococcus faecium</i>	<i>Staphylococcus carnosus</i>	<i>Streptococcus pneumoniae</i>
<i>Enterococcus hirae</i>	<i>Staphylococcus cohnii</i> (includes <i>S. cohnii</i> ssp <i>cohnii</i> and <i>S. cohnii</i> ssp <i>ureolyticus</i>)	<i>Streptococcus porcinus</i>
<i>Enterococcus raffinosus</i>	<i>Staphylococcus cohnii</i> ssp <i>cohnii</i>	<i>Streptococcus pyogenes</i>
<i>Enterococcus solitarius</i>	<i>Staphylococcus cohnii</i> ssp <i>ureolyticus</i>	<i>Streptococcus salivarius</i>
<i>Erysipelothrix rhusiopathiae</i>	<i>Staphylococcus epidermidis</i>	<i>Streptococcus salivarius</i> group (includes <i>S. salivarius</i> and <i>S. vestibularis</i>)
<i>Gardnerella vaginalis</i>	<i>Staphylococcus equorum</i>	<i>Streptococcus sanguis</i>
<i>Gemella haemolysans</i>	<i>Staphylococcus felis</i>	<i>Streptococcus sanguis</i> group (includes <i>S. crista</i> , <i>S. gordonii</i> , <i>S. parasanguis</i> and <i>S. sanguis</i>)
<i>Gemella morbillorum</i>	<i>Staphylococcus gallinarum</i>	<i>Streptococcus sobrinus</i>
<i>Gemella species</i> (includes <i>G. haemolysans</i> , and <i>G. morbillorum</i>)	<i>Staphylococcus haemolyticus</i>	<i>Streptococcus uberis</i>
<i>Globicatella sanguis</i>	<i>Staphylococcus hominis</i>	<i>Streptococcus vestibularis</i>
<i>Helicobacterium kunzli</i>	<i>Staphylococcus intermedius</i>	<i>Turicella ostioides</i> ¹

¹ These taxa have <10 unique CRYSTAL profiles in the current database.

PRODUCT DESCRIPTION:

The main component of the BBLCRYSTAL™ GP ID System is the BBLCRYSTAL GP panel assembly, consisting of the CRYSTAL base and lid. The BBLCRYSTAL lid consists of 29 dehydrated biochemical/chromogenic/fluorogenic substrates and one fluorogenic negative control, on the ends of plastic prongs. The CRYSTAL base consists of 30 matching wells; its design allows inoculation of all 30 wells in a single step by pouring the suspension of pure culture into the target area and tilting the base. The test inoculum rehydrates the dried substrates and initiates test reactions.

The pure culture suspension is prepared by picking several small colonies of the same morphology from media such as Trypticase® Soy Agar with 5% Sheep Blood or Columbia Agar with 5% Sheep Blood, or alternatively selective media such as Phenylethyl Alcohol Agar with 5% Sheep Blood or Columbia CNA Agar with 5% Sheep Blood. A standardized suspension of this culture is prepared in the BBLCRYSTAL™ ANR, GP, RGP, N/H ID Inoculum Fluid provided. The suspension is added to the target area of the panel base, which the user then rocks back and forth to inoculate all the wells contained in the base.

After the base/lid assembly has been incubated for 4 hours at 35-37°C, the assembly is placed on the BBLCRYSTAL Panel Viewer and the color reactions are visually compared to the BBLCRYSTAL GP Color Chart provided. Each reaction is scored as a positive (+) or negative (-) and recorded on the BBLCRYSTAL GP Report Form.

After all results are read, a 10-digit numerical profile is calculated by assigning a value of 4, 2, or 1 to each positive reaction. (Negative reactions are scored as "0".) The values for each column are then added together to obtain the 10 digit Profile Number.

The BBLCRYSTAL ID System Electronic Codebook is loaded into the user's PC and the appropriate database is selected. Then the Profile Number and results of any off-line tests are entered, and the Codebook gives one of the following three results:

- (a) a definitive ID;
- (b) a tie between two or more species; or
- (c) no ID possible with data submitted.

In the case of a definitive ID or a tie between two or more potential ID's, the user can access the statistics for that ID as well as background information for the species identified.

In the case where no ID is possible, the Codebook suggests that the user perform a purity check of the test isolate. If culture purity has been confirmed, then it is likely that (i) the test isolate is producing atypical Crystal reactions (which may also be caused by procedural errors), (ii) the test species is not part of the intended taxa, or (iii) the system is unable to identify the test isolate with the required level of confidence. Once user error has been ruled out, the Codebook suggests that additional testing must be done to establish an identification.

PERFORMANCE DATA:

Clinical Correlation:

In a study conducted at four clinical laboratories, the performance of the BBLCRYSTAL Gram Positive ID System was evaluated against a combination of the API Staph, API 20 Strep, and API Coryne Identification Systems, and conventional methodologies. Fresh, routine isolates arriving in the clinical laboratory, as well as previously identified isolates of the clinical trial sites' choice were utilized in the study. A total of 735 gram positive aerobic isolates were tested; 90% (668/735) of these isolates were correctly identified (including supplemental testing) using the BBLCRYSTAL™ GP ID System; 7.6% (56/735) were incorrectly identified; and 1.5% (11/735) yielded a "No Identification" result.

Reproducibility:

At the same four clinical laboratories, reproducibility of the BBLCRYSTAL™ GP ID System was established by testing ten (10) Quality Control strains in triplicate on three days. Evaluations were performed of individual and overall reproducibility of substrate reactions, of QC organism reactions, and of inter- and intra-laboratory reproducibility.

Overall reproducibility was calculated as 96.7%. Reproducibility of individual substrate reactions ranged from 79.2 to 100%; individual QC organism reaction reproducibility ranged from 91.4 to 99.8%; and individual site reproducibility from 95.5 to 97.8%.