

DEC 17 1996

IN 50115 P1072  
**RÜSCH.**  
INTERNATIONAL  
**Group Regulatory Affairs**  
A Subsidiary of Teleflex Incorporated (USA)

Tall Pines Park  
Jaffrey, NH 03452  
(603) 532-7706  
FAX (603) 532-8211 or 6108

**510(k) SUMMARY**  
[As required by 21 CFR 807.92]

1. **Submitter and Contact Person**

Neil R. Armstrong  
Group Regulatory Affairs Manager  
Rüsch International  
Tall Pines Park  
Jaffrey NH 03452  
Tel: (603) 532 7706  
Fax: (603) 532 8211

2. **Device Name**

Trade Name

Rüsch Simplastisc Foley Catheterization Set

Common Name

Urological Catheter

Classification Name

Catheter, Urological 78 KOD  
21 CFR 876.5130; Class II medical device

3. **Comparison Devices**

Rüsch Simplastisc Foley Catheter - Preamendment  
Urotec/Franklin Soft Simplastisc Catheters - K851684  
Inmed Foley Catheterization Tray - K832363

4. **Description of Device**

The Rüsch Simplastisc Foley Catheterization Set contains a pack of lubricating jelly and a Rüsch Simplastisc Foley Catheter.

The Rüsch Simplastisc Foley Catheters used in the kit will be available:

- in 2 way and 3 way formats
- in various sizes from 12 - 26 F.G.
- in two hardnesses
- with various sizes of balloon
- with various tip and eye configurations.

Common features of all the catheters are: a clear polyvinyl chloride tube, a radiopaque stripe of BaSO<sub>4</sub> loaded PVC is fully encapsulated in the tube wall, a latex balloon attached by adhesive bonding, a funnel connected to the main lumen and a luer activated valve for filling and emptying the balloon.

5. **Intended Use**

The Rüsç Simplastic Foley Catheterization Set contains:  
 one Rüsç Simplastic Foley Catheter  
 one pack of lubricating jelly.

The Rüsç Simplastic Foley Catheter is intended to be used to pass fluids to and from the urinary tract. A pack of lubricating jelly is intended to assist insertion of the catheter through the urethra.

6. **Summary of Technological Characteristics**

The Rüsç Simplastic Foley Catheterization Set is manufactured from the same materials and by the same processes (including sterilization) as the Rüsç Simplastic Foley Catheter (preamendment) and Urotec/Franklin Soft Simplastic Catheters (K851684), and contains the same "gel" pack as the Inmed Foley Catheterization Tray (K832363).

Sterile catheters have been biocompatibility tested in accordance with ISO 10993 and the FDA "Blue Book Memo" #G95/1.

7. **Summary of Performance Data**

Laboratory bench testing has been completed to section #XI of the FDA "Draft Guidance for the Content of Pre-market Notifications for Conventional and Antimicrobial Foley Catheters":

- drainage lumen flow rate
- balloon resistance to rupture
- pullout resistance of inflated balloon
- maintenance of balloon inflation over extended time
- manufacturing tolerances for tip, balloon and shaft
- deflation after period of inflation
- biocompatibility testing for patient contacting materials