

K964116

DEC 24 1996

SUMMARY OF SAFETY AND EFFECTIVENESS**Quinidine Method for Bayer Technicon Immuno 1[®] System**

Listed below is a comparison of the performance between the Immuno 1 Quinidine method and a similar device that was granted clearance of substantial equivalence (Syva EMIT[®] 2000 Quinidine Assay, Behring Diagnostics Inc.). The information used in the Summary of Safety and Effectiveness was extracted from the Immuno 1 Quinidine method sheet and the Syva EMIT[®] 2000 Quinidine Assay insert.

INTENDED USED

This *in vitro* method is intended to quantitatively measure Quinidine, an anti-arrhythmic drug, in human serum or plasma (heparin) using Syva EMIT[®] 2000 Assay on a *Technicon Immuno-1[®]* system. Measurements of quinidine are used in the diagnosis and treatment of quinidine overdose and in monitoring serum levels of quinidine to ensure appropriate therapy.

METHOD	Immuno 1 Quinidine		Syva EMIT [®] 2000 (predicate Device)	
Part No.	Reagents	T01-3678-51	Reagents	4Q019UL
	Calibrators	T03-3680-01	Calibrators	4Q109UL
Minimum Detectable Conc.	0.06 µg/mL		0.25 µg/mL	
Precision (Between-Run)	1.2 µg/mL	5.5%	1.5 µg/mL	6.5%
	3.5 µg/mL	4.8%	3.4 µg/mL	5.2%
	6.0 µg/mL	5.3%	5.7 µg/mL	5.6%
Correlation	$y = 1.02x - 0.16$			
	where			
	y = Immuno 1 Quinidine method			
	x = Syva EMIT [®] 2000 Quinidine Assay*			
	n = 25			
	r = 0.991			
	S _{yx} = 0.148 µg/mL			

*This assay was performed on Roche COBAS FARA II[®] Instrument.

Gabriel J. Muraca, Jr.

Gabriel J. Muraca, Jr.
Manager Regulatory Affairs
Bayer Corporation
511 Benedict Avenue
Tarrytown, New York 10591-5097

10/10/96
Date

* Please Note: Additional # patients requested, DATA provided in 11/26/96 memo. No significant change in correlation. R013.

12/9/96