

# BECKMAN

## Summary of Safety & Effectiveness

### IMAGE™ Immunochemistry System Kappa (KAP) and Lambda (LAM) Light Chain Reagents

1.0 **Submitted By:**

K904260

Annette Hellie  
Regulatory Specialist, Product Submissions  
Beckman Instruments, Inc.  
200 S. Kraemer Blvd., W-337  
Brea, California 92822-8000  
Telephone: (714) 993-8767  
FAX: (714) 961-4457

FEB 10 1997

2.0 **Date Submitted:**

24 October 1996

3.0 **Device Name(s):**

3.1 **Proprietary Names**

IMAGE™ Immunochemistry System Kappa Light Chain (KAP) Reagent  
IMAGE™ Immunochemistry System Lambda Light Chain (LAM) Reagent

3.2 **Classification Name**

Immunoglobulin (light chain specific) immunological test system (21 CFR § 866.5550)

4.0 **Predicate Device(s):**

IMAGE System Reagent	Predicate	Manufacturer	Docket Number
IMAGE System Kappa (KAP) Reagent	Beckman Kappa Light Chain (KAP) Reagent	Beckman Instruments, Inc.	K884276A K902484
IMAGE System Lambda (LAM) Reagent	Beckman Lambda Light Chain (LAM) Reagent	Beckman Instruments, Inc.	K884597A K902484

Beckman Instruments, Inc.

**5.0 Description:**

The IMMAGE Immunochemistry System KAP and LAM Reagents in conjunction with Beckman Calibrator 1, are intended for use in the quantitative determination of kappa and lambda light chain concentrations respectively in human serum and urine samples on Beckman's IMMAGE Immunochemistry System.

**6.0 Intended Use:**

The IMMAGE Immunochemistry System Kappa (KAP) Reagent, when used in conjunction with Beckman IMMAGE™ Immunochemistry Systems and Beckman Calibrator 1, is intended for the quantitative determination of kappa light chains (free and bound) by rate nephelometry.

The IMMAGE Immunochemistry System Lambda (LAM) Reagent, when used in conjunction with Beckman IMMAGE™ Immunochemistry Systems and Beckman Calibrator 1, is intended for the quantitative determination of lambda light chains (free and bound) by rate nephelometry.

**7.0 Comparison to Predicate(s):**

The following table shows similarities and differences between the predicates identified in Section 4.0 of this summary.

Reagent	Aspect/Characteristic	Comments
SIMILARITIES		
IMMAGE System KAP and LAM Reagents	Initial Analytic Range	Same as Beckman Kappa and Lambda Light Chain reagents
	Nephelometric methodology	
	Antibody source (goat)	
DIFFERENCES		
IMMAGE System KAP and LAM Reagents	Antigen excess testing solution	IMMAGE KAP & LAM have antigen excess testing solution included in the reagent cartridge, while the Beckman Kappa and Lambda require off-line preparation of the solution.
	Antibody concentration	IMMAGE KAP and LAM have a higher antibody concentration than the Beckman Kappa and Lambda reagents

## 8.0 Summary of Performance Data:

The data in the Premarket Notification on safety and effectiveness supports a finding of substantial equivalence to chemistry test systems already in commercial distribution. Equivalence is demonstrated through method comparison, stability, and imprecision experiments that relate results obtained from the Beckman Reagents to the IMAGE System Reagents.

**Method Comparison Study Results**  
**IMAGE Kappa (KAP) and Lambda (LAM) Reagents**

Analyte	Sample Type	Slope	Intercept	r	n	Reference Method
IMAGE KAP Reagent	serum	1.037	-9.16	0.987	220	Beckman KAP Reagent
	urine	0.970	1.36	0.980	103	
IMAGE LAM Reagent	serum	1.020	-23.5	0.991	236	Beckman LAM Reagent
	urine	1.009	-0.08	0.954	40	

**Stability Study Results**

Reagent	Product Claim
IMAGE KAP & LAM	24 month shelf-life 14 day open container stability 14 day calibration stability

**Estimated Within-Run Imprecision**

Sample	Mean (µg/mL)	S.D. (µg/mL)	%C.V.	N
<b>Kappa Light Chain (KAP) Serum</b>				
Level 1	540	15.6	2.9	80
Level 2	1113	24.1	2.2	80
Level 3	2377	61.2	2.6	80

Sample	Mean (µg/mL)	S.D. (µg/mL)	%C.V.	N
<b>Kappa Light Chain (KAP) Urine</b>				
Level 1	2.33	0.059	2.5	30
Level 2	9.34	0.236	2.5	30
Level 3	19.2	0.25	1.3	30

Sample	Mean (µg/mL)	S.D. (µg/mL)	%C.V.	N
<b>Lambda Light Chain (LAM) Serum</b>				
Level 1	247	5.1	2.1	80
Level 2	718	14.0	1.9	80
Level 3	1233	37.2	3.0	80

Sample	Mean (µg/mL)	S.D. (µg/mL)	%C.V.	N
<b>Lambda Light Chain (LAM) Urine</b>				
Level 1	8.69	0.204	2.4	30
Level 2	17.1	0.20	1.2	30
Level 3	30.0	0.89	3.0	30

This summary of safety and effectiveness is being submitted in accordance with the requirements of the Safe Medical Device Act of 1990 and the implementing regulation 21 CFR 807.92.